Engineering Conferences International ECI Digital Archives

Single-use Technologies II: Bridging Polymer Science to Biotechnology Applications

Proceedings

5-7-2017

XDR-500 MO—Single-use fermentor for microbial processes

Jozsef Vasi GE Healthcare, Sweden

Andreas Castan GE Healthcare, Sweden, andreas.castan@ge.com

Sandeep Kristiansson GE Healthcare, Sweden

Follow this and additional works at: http://dc.engconfintl.org/biopoly_ii Part of the <u>Materials Science and Engineering Commons</u>

Recommended Citation

Jozsef Vasi, Andreas Castan, and Sandeep Kristiansson, "XDR-500 MO—Single-use fermentor for microbial processes" in "Single-use Technologies II: Bridging Polymer Science to Biotechnology Applications", kta Mahajan (Genentech, Inc., USA) Gary Lye (University College London, UK) Regine Eibl-Schindler (Zurich University of Applied Science, Switzerland) Eds, ECI Symposium Series, (2017). http://dc.engconfintl.org/biopoly_ii/13

This Abstract and Presentation is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Single-use Technologies II: Bridging Polymer Science to Biotechnology Applications by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.



Xcellerex[™] XDR-500 MO single-use fermentor for microbial processes

József Vasi, Andreas Castan and Sandeep Kristiansson GE Healthcare Bio-Sciences AB, Björkgatan 30, 75184 Uppsala, Sweden

Abstract and introduction

Single-use bioreactors are well-established in mammalian cell-based biopharmaceutical productions. However, adoption of single-use technologies for microbial/yeast-based production has been slow due to the unique requirements of fermentation processes. The typical challenges include supporting high-cell density cultivations, high oxygen transfer capacity, and efficient heat removal.

Here, we present the single-use XDR-500 MO stirred-tank fermentor system specifically designed for microbial/yeast applications. The performance of the fermentor was evaluated using an *E. coli* fed-batch process producing a domain antibody (dAb) as a model system. The process was originally developed for a lab-scale conventional stainless steel (SS) system. This work describes the strategy applied to transfer and scale-up this process to the single-use bioreactor (SUB) without the need for extensive re-optimization due to technology change. Moreover, process economy of dAb production based on single-use and conventional technology was compared *in silico*.

Cell densities and dAb expression levels similar to the lab-scale conventional bioreactor were achieved in XDR-500 MO, confirming feasibility of single-use technology even in microbial processes.

Materials and methods

Fed-batch process for *E. coli* RV308 clone expressing a dAb, originally developed for lab-scale conventional SS bioreactor, was transferred and scaled up to XDR-500 MO SUB. Mineral salt medium complemented with yeast extract in combination with 60% (w/v) glycerol substrate feed was used for cultivations. Process parameters critical for dAb expression were kept constant between technologies, while others were adapted to single-use technology (Table 1). Dissolved oxygen level was cascade controlled by stirrer speed, airflow, and oxygen enrichment of airflow in the single-use system. Substrate feeding strategy was adapted to the SUB and three different rates were evaluated. Medium and feed solutions were prepared in Xcellerex XDUO mixer bags and sterile filtered into the to single-use bioreactor. All cultures have been induced at an OD_{600nm} of 80.






Table 1. Process transfer and scale-up strategy

Process parameters				
Kept constant	Adjusted			
pH and DO set points	DO control strategy			
Medium and feed composition	Medium and feed preparation p			
Induction conditions	Substrate feed profile			
Temperature set point				

GE, the GE Monogram, and Xcellrex are trademarks of General Electric Company. © 2017 General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information. GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden. For local office contact information, visit gelifesciences.com/contact 29254432 AA 03/2017

procedures

Experimental results

Improved process understanding during transition to SUB indicated substrate overfeeding in the SS bioreactor. Evaluation of three reduced feed regimes showed that even applying as low as 60% of the original feed rate to XDR-500 MO was sufficient to maintain cell growth, as shown for XDR-Run 3 in Figure 1, and product formation was found to be at comparable levels to the original process. The reduced substrate feed also led to decreased oxygen consumption rate.

Levels of expressed dAb achieved in XDR-500 MO after process transfer and scale-up were reproducible and comparable with the SS bioreactor (Fig 2) even without process optimization. Furthermore, it was shown that dAb concentration increased linearly after induction (Fig 3), indicating a straightforward strategy to improve product titer.




Fig 1. Cell growth in XDR-500 MO at studied feed rates. Cell growth even at lowest feed rate was comparable with stainless steel bioreactor.


Process economy

Comparison of process economy between single-use and conventional SS production facilities for dAb manufacturing revealed up to twice as high facility utilization when producing in a single-use system (Table 2).


According to our model, at 4–15 annual batches (e.g., for toxicology or early clinical phase studies), it is more costeffective to manufacture in a single-use facility compare with a conventional facility (Fig 5).

Table 2. Comparison of process economy for dAb production between XDR-500 MO SUB and SS bioreactor (scenarios for single and multiple product facilities were evaluated)

Parameter	Scenario 1 – Single product facility		Scenario 2 – Multi product facility	
Type of Bioreactor	SS	SUB	SS	SUB
Maximum number of batches	100	150	67	135
Relative cost comparison	1×	1.47×	1×	1.40×
Relative facility utilization	1×	1.5×	1×	2×



densities (Fig 4).



and stainless steel reactors. Product concentrations were determined after 5 h induction time.

Fig 2. Comparison dAb expression levels obtained in single-use XDR-500 MO





Fig 5. Total production cost to number of batches produced annually in single-use and stainless steel facilities.



Critical process parameters like dissolved oxygen level, pH, and temperature were controlled efficiently within their targeted range under the entire cultivation. Temperature was maintained at $37^{\circ}C \pm 0.1^{\circ}C$ even when metabolic heat generation peaked at 38 000 BTU/h. Instead of increasing headspace pressure, DO was maintained with oxygen enrichment, enabling cell growth to high cell



Conclusions

- XDR-500 MO allowed robust control of challenging process parameters (DO, pH, temperature).
- The model process successfully transferred and scaled up from SS bioreactor to SUB without the need for extensive modifications.
- Cell growth and dAb production level were comparable across technologies/scales even with the implemented changes.
- Higher production flexibility and facility utilization can be achieved with single-use systems.
- Single-use systems are viable alternatives to conventional SS bioreactors for microbial processes.