

5-6-2018

Novel transposase tools for cell-line engineering

Ferenc Boldog

ATUM (formerly DNA2.0), USA, fboldog@atum.bio

Jeremy Minshull

ATUM (formerly DNA2.0), USA

Mark Welch

ATUM (formerly DNA2.0), USA

Kate Caves

ATUM (formerly DNA2.0), USA

Follow this and additional works at: <http://dc.engconfintl.org/ccexvi>



Part of the [Engineering Commons](#)

Recommended Citation

Ferenc Boldog, Jeremy Minshull, Mark Welch, and Kate Caves, "Novel transposase tools for cell-line engineering" in "Cell Culture Engineering XVI", A. Robinson, PhD, Tulane University R. Venkat, PhD, MedImmune E. Schaefer, ScD, J&J Janssen Eds, ECI Symposium Series, (2018). <http://dc.engconfintl.org/ccexvi/53>

This Abstract and Presentation is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Cell Culture Engineering XVI by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

Novel Transposases for Cell-Line Development & Engineering

Ferenc Boldog*, Sowmya Rajendran, Maggie Lee, Molly Hunter, Candice Talsalt, Lynn Webster, Vivi Truong, Elizabeth Hart, Nicolay Kulikov, Thomas Purcell, Mark Welch, Miles Scotcher**, Claes Gustafsson*** and Jeremy Minshull
 ATUM (formerly DNA2.0), Newark, CA USA. *fboldog@atum.bio, **miles@atum.bio, ***claes@atum.bio



Abstract

The generation of stable cell lines is critical for the commercial production of protein therapeutics. Current methodologies to introduce engineered recombinant genes into production strains relies on homologous recombination, a method limited by poor integration rates (<<1%), concatemer formation and common transgene rearrangements. To address this limitation, ATUM has identified and engineered two orthologous transposases, designated Leap-In 1 (from frog) and Leap-In 2 (from insect), that stably integrate synthetic transposons into the host genome with close to 100% integration frequency and an absolute minimum of concatemers or other genetic instabilities. Each synthetic transposon can encode up to 4 unique open reading frames (ORFs) flanked by inverted terminal repeats recognized by the cognate transposase. Each ORF is under discrete expression control thereby enabling tuning of expression ratios by adjustment of each ORF independently. Leap-In transposases speed cell recovery times, enabling the generation of high expressing stable pools (up to 5g/L) in as little as 2 weeks post transfection.

Leap-In transposase benefits

1	2	3
Application-optimized expression constructs Maximize expression levels using ATUM's proprietary gene sequence and vector optimization tools. Adjust and tune expression levels of your gene(s) by testing combinations of vector elements. Express multichain proteins (such as bispecific antibodies) or multiple genes within a pathway .	Shortened stable pool recovery times Clonal productivity distribution is characteristically higher and more uniform in stable pools. Fewer clones need to be screened to enable identification of highly productive clonal cell lines. Initiate process development earlier with more representative pool derived reagents.	Stable integration combined with structural integrity Transposases integrate the entire, intact transposon, thereby maintaining its structural integrity , and expression balance . Clonal cell lines created using Leap-In transposases exhibit genetic stability and stable productivity over 60 population doublings. Clonal cells lines created using Leap-In transposases exhibit genetic stability and stable productivity over 60 population doublings. Control integration copy number by manipulating dose of transposase and synthetic transposon. Foot-print free excision of transposon from host genome.

Leap-In transposase process

Legend:
 ▼ Transposon Recognition Site
 ▭ Inverted Terminal Repeat
 ▭ Gene of interest
 ● Transposase protein

- Transposon and transposase (as mRNA) co-transfected into target cell
- Broad host cell range
- ATUM uses 2 in-house cell lines, both with full traceability
 - HD-BIOP3 GS null (Horizon Discovery)
 - DG44 DHFR null (Prof. Lawrence Chasin)
- Single copies of entire transposon are independently integrated by a "cut and paste" mechanism into multiple sites across the host genome
- 48 hours post transfection, no mRNA or transposase are detectable


Process flow:
 Stable pools → Stable pool ranking → Clonal cell line generation → Stability testing → Research Cell Bank

1 Application - optimized expression constructs

- Up to 4 open reading frames under independent expression control.
- No size limit for synthetic transposon.

- Multiple, parallel selection stringencies.
- High productivities in stable pools up to 5g/L and clonal up to 6.5g/L.

Adjustment of expression ratios for 2-, 3- and 4-ORF vectors




2 Shortened timelines

Transposases speed cell pool recovery post-transfection.

Stable pool	Variable Parameters		
	GS promoter	Transposase	MSX (µM)
1	mPGK	YES (wild type)	50
2	mPGK	NO	50
3	Attenuated-1	YES (wild type)	0
4	Attenuated-1	NO	0
5	negative control	NO	0


Upper graph - Under typical selection conditions (strong GS promoter, 50 µM MSX), cell pools created using transposases recover viability >7 days faster than cell pools created by random integration.
 Lower graph - Under drug-free selection conditions (attenuated GS promoter), cell pools created using transposases recover viability in ~3 weeks, whereas cell pools created by random integration fail to recover.

Screening and selection of clones expressing transposon-based antibody construct




3 Stable integration combined with structural integrity

Transposase	Transposon based integration	Random integration	TG/TG fusion
+	24/25	1/25	2
+	23/24	1/24	1
-	No site identified	No site identified	23




An antibody expression construct was transfected into HD-BIOP3 GS null CHO1 cells (Horizon Discovery) in the presence (+ transposase) or absence (- transposase) of Leap-In transposase. Targeted Locus Amplification was used to analyze the structural integrity of expression cassette integrations.

Clonal cell lines generated using Leap-In transposases exhibit genetic stability and stable productivity over 60 population doublings.




Nine CHO GS KO clones, derived from the transposase-mediated integration of an antibody expression cassette, were grown in 24-deep well plate for 60 generations (population doublings (PD)) in the presence (-Gln) and absence (+Gln) of selection. Productivity and integrated copy number of clones after 60 PD were compared to T0. No significant change in productivity or copy number under the with glutamine (+Gln) and minus glutamine (-Gln) growth conditions were observed compared to T0.

1:1 correlation between copy number and specific productivity in transposase mediated DG44 stable pools.




Stable pools derived from the transposase-mediated integration of an antibody expression cassette in DG44 cells, were grown under fed-batch conditions. Productivity and copy number were measured. The correlation between copy number and specific productivity indicates that each integrated copy is intact and functional.

Foot-print free excision of transposon from host genome.




Studies are currently ongoing to demonstrate that integration of single copies of transposons is possible.

Excision efficiency of wild type and engineered hyperactive Leap-In@1 transposases



Hyperactive transposase variants designed using ProteinGPS® demonstrate foot-print free excision. Variants were ranked by excision efficiency measured as per figure on left. Different concentrations of hyperactive transposase mRNA demonstrated a dose response. Hyperactive transposase variants tested show an approximately 80-fold higher activity compared to wild type in this assay.

Typical timeline of cell line development at ATUM



- Stable pool-derived protein is available as early as 2 weeks post transfection.
- Clonal cell line without stability information, available as early as 12 weeks post-transfection.
- RCBs (Research Cell Banks) released as early as 20 weeks post-transfection.