

9-25-2018

# Continuous Process Performance Enhancements for 50-500 L S.U.B.s

Ben Madsen

*Thermo Fisher Scientific, USA*, [ben.madsen@thermofisher.com](mailto:ben.madsen@thermofisher.com)

Jordan Cobia

*Thermo Fisher Scientific, USA*

Nephi Jones

*Thermo Fisher Scientific, USA*

Kevin Mullen

*Thermo Fisher Scientific, USA*

Follow this and additional works at: [http://dc.engconfintl.org/sut\\_iii](http://dc.engconfintl.org/sut_iii)

 Part of the [Engineering Commons](#)

---

## Recommended Citation

Ben Madsen, Jordan Cobia, Nephi Jones, and Kevin Mullen, "Continuous Process Performance Enhancements for 50-500 L S.U.B.s" in "Single-Use Technologies III: Scientific and Technological Advancements", Weibing Ding, Amgen Martina Micheletti, University College London Robert Repetto, Pfizer Eds, ECI Symposium Series, (2018). [http://dc.engconfintl.org/sut\\_iii/46](http://dc.engconfintl.org/sut_iii/46)

This Abstract and Presentation is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Single-Use Technologies III: Scientific and Technological Advancements by an authorized administrator of ECI Digital Archives. For more information, please contact [franco@bepress.com](mailto:franco@bepress.com).



# Continuous Process Performance Enhancements for 50-500 L S.U.B.s


September 25, 2018

Ben Madsen, Ph.D.  
Senior Process Development Engineer, Advanced Technologies  
Thermo Fisher Scientific


- Perfusion Overview and Workflow
  - Mass transfer limitations
  - Mixing limitations
- S.U.B. Enhancements
- Results
  - Mass transfer performance
  - Cell culture case study
  - Cell culture mass transfer models



# Perfusion vs. Fed-batch



- Feed addition from D3-end
- Volume change in reactor
- High titers achievable
- Large production vessels ( $\geq 2000\text{L}$ )



- Continuous filtration for  $\geq 21$  days
- Cell retention usually via filtration (0.2-0.65 $\mu\text{m}$ )
  - Concentrated Fed-batch uses 30-50 kDa
- Used for:
  - Unstable products
  - Intensified seed train
  - Cell banking

# Perfusion Advantages/Challenges

- Advantages


- Smaller reactors
- Healthier culture with less waste = predictable product quality
- Higher cell densities
- Consistency from steady state production
- Less frequent bioreactor turnaround
- Clean harvest stream for DSP
- Continual harvest ensures some success despite late-stage contamination
- Cost savings

- Challenges

- What is steady state?
- Process control and integration
- More complex setup (Regulatory)
- High validation costs
- Filter fouling leads to culture shutdown
- Higher gas usage = foaming risk
- Sensor stability



# Perfusion Work Flow



- Thermo Fisher Scientific Offerings:

- Cell line development
- Media optimization
- Media hydration
- Sterile media hold
- **Bioreactor production vessels**
- Product pooling
- Waste containment
- Chromatography resins
- Bulk storage and final fill




# Overcoming Legacy S.U.B. System Limitations

- Custom 500 L S.U.B. configuration
  - Standard impeller (25.1 cm diameter)
  - Standard laser DHS (980 pores, 0.368 mm Ø)
  - 2 ea frit (1 standard, 1 additional)
- Standard Dynamic  $k_L a$  testing
  - 1 g/L pluronic
  - 3.5 g/L HEPES
  - DHS = **15 slpm**
  - Agitation = **150 RPM** ( $65 \text{ W/m}^3$ )

		Additional Frit (slpm)		
		10	20	40
Standard Frit (slpm)	10	36.5	37.6	38.1
	20	37.6	40.5	42.6
	40	39.6	40.2	41.9


- Limited benefit from running frits at extreme flow rates even at max RPM
  - Combined air sparge of 95 slpm (0.19 VVM) = ~42/hr
  - Second option would be to increase DHS sparging rate



# S.U.B. Mixing Enhancements

- Implement upsized impeller (1000 L impeller in 500 L S.U.B., etc.)

S.U.B. Size		50	100	250	500
Max Motor RPM		200	200	150	150
Standard Impeller	Diameter (cm)	11.11	14.6	20	25.1
	Speed at 20 W/m <sup>3</sup> (RPM)	182.5	145.8	117.1	101.1
	Max PIV at max RPM (W/m <sup>3</sup> )	26.3	51.6	42.0	65.4
Upsized Impeller	Diameter (cm)	14.6	20	25.1	32.1
	Speed at 20 W/m <sup>3</sup> (RPM)	115.7	86.3	80.2	67.1
	Max PIV at max RPM (W/m <sup>3</sup> )	<b>103.2</b>	248.9	130.8	223.7




All systems meet objective of 100 W/m<sup>3</sup>

Note: Safely within the qualified design space of the 1000L/2000L S.U.M. (rated up to 350 RPM).




# S.U.B. BPC Perfusion Enhancements

1. Enhanced DHS for each S.U.B. size
  - Increased DHS efficiency (smaller pore/bubble)
  - Boost capacity (0.2 VVM @ <15m/s GEV)
  - Frit optional in all S.U.B. designs
2. Upsized impeller in each S.U.B. size (n+1)
3. SU foam sensor for automated antifoam dispensing
4. Port for cell bleed – maintain cell density in perfusion cultures
5. 1" ID ATF/TFF connections available



# Characterization Results with Enhanced Sparging and Mixing




# 50L Perfusion S.U.B. Case Study Setup

	Components/Conditions	Setting
Cells/Media	Cell Line	CHO-DP12 (Supplied by Repligen)
	Seed Density	$0.4 \times 10^6$ cells/mL
	Target Cell Density	Various up to maximum
	Media Composition	CD OptiCHO media + 100 ng/mL LR3 + 4 mM Glutamax
	Base	N/A
	Antifoam	Antifoam C (10,000 ppm simethicone stock)
	Glucose	450 g/L stock
Bioreactor	Bioreactor W.V.	40 L
	Temp	37°C
	DO Setpoint	30% ( $O_2$ /air through DHS)
	pH	$7.0 \pm 0.2$ ( $CO_2$ )
	Air	Headspace sweep at 2 slpm
	Agitation	107-184 RPM 20-100 $W/m^3$
	Antifoam dosing	Sensor-based antifoam addition
	Glucose Feed	Maintain 1-3 g/L with constant feed
Cell Bleed	Sensor-based cell bleed (ABER)	
ATF	ATF System	ATF6SU
	ATF Filter	F6:RF02PES 0.2 micron PES HF
	Filter SA	$2.53 \text{ m}^2$
	Perfusion rate	Day 0-3: None
		Day 3-26: variable depending on cell density/nutrient demands, 1-4 VVD (40-160 L/day)
	Flux	0.67-2.67 LMH
	ATF Rate	19.2 LPM
	Shear Rate	$2264 \text{ s}^{-1}$



# Perfusion Process Automation

- Standard
  - DO (TruFluor DO)
  - pH
  - RTD
- Foam sensor linked to antifoam pump
- Cell density sensor linked to cell bleed pump
- S.U.B. load cells linked to media pump (Heaviside equation)




$$\text{Feed} = 2 \cdot \text{Target} \cdot \left( \frac{1}{1 + \exp^{-2 \cdot k \cdot x}} \right)$$



# Cell density balance (probe-less)

$$V_R \cdot \frac{dX_R}{dt} = V_R \cdot \mu \cdot X_R - Q_B \cdot X_R$$
$$\mu = \frac{\ln(X_2) - \ln(X_1)}{t_2 - t_1} + \frac{Q_B}{V_R}$$
$$Q_B = \left( \mu - \frac{\ln(X_2/X_1)}{t_2 - t_1} \right) \cdot V_R$$




## Nomenclature


## Equation

- Change in cell density = effective growth – cell bleed
- Rearrange, integrate, solve for  $Q_B$  to obtain bleed rate
- Best to average growth rates over multiple days

- $V_R$  = reactor working volume
- $X_R$  = cell density
- $dX_R/dt$  = cell density change over time
- $\mu$  = cell growth rate
- $Q_B$  = reactor bleed rate



# Glucose Mass Balance Equations



Cell Growth:

$$V_R \frac{dX_R}{dt} = V_R \cdot \mu \cdot X_R - Q_B \cdot X_R$$

Overall Glucose Consumption:

$$\frac{dG_R}{dt} = Q_M \cdot G_M + Q_G \cdot G_G - Q_P \cdot G_P - Q_B \cdot G_R - X_R \cdot V_R \cdot q_G$$

If

$$\frac{dG_R}{dt} = 0 \text{ and } \frac{dX_R}{dt} = 0$$


Then

$$Q_G = \frac{[D \cdot G_P + \mu \cdot G_R + X_R \cdot q_G - (D + \mu) \cdot G_M] \cdot V_R}{(G_G - G_M)}$$


Nomenclature

$Q$  = flow  
 $G$  = [glucose]  
 $\mu$  = growth  
 $q_G$  = gluc consumption  
 $V$  = volume  
 $X$  = [cells]  
 $D$  = perfusion rate

# Cell Culture Results – 50 L S.U.B. Cells




# Cell Culture Results – 50 L S.U.B. Culture Parameters






# Cell Culture Results – 50 L S.U.B. with Antifoam Control Loop



# Cell Culture Results – 50 L S.U.B. Gassing, Agitation



# Mass Balance Models – Cells = f(gassing, mixing)

- Can we build an equation to model sustainable cell concentration based on mixing and gassing?

- $OUR = q_{O_2} \cdot C_X$

- $OTR = k_L a \cdot (C_{O_2}^* - C_{O_2})$


$$C_{O_2} = pp_{O_2,air} \cdot SP$$

$$C_{O_2}^* = \left( \frac{F_{O_2} \cdot pp_{O_2} + F_{air} \cdot pp_{O_2,air}}{F_{TOT}} \right)$$

$$k_L a = \alpha \cdot PIV^\beta \cdot vvm^\gamma$$


- If  $OUR = OTR$

- $\therefore q_{O_2} \cdot C_X = \alpha \cdot PIV^\beta \cdot vvm^\gamma \left( \frac{F_{O_2} \cdot pp_{O_2} + F_{air} \cdot pp_{O_2,air}}{F_{TOT}} - pp_{O_2,air} \cdot SP \right)$



# Mass Balance Models – $\text{CO}_2 = f(\text{acids, base})$

- $\text{pH} = \text{baseline pH} - \text{acid buildup} + \text{base buildup}$
- Acid/base accumulations
  - $\text{CO}_2$  concentration – controllable ( $k_L a_{\text{CO}_2}$ )
  - Lactate
  - Ammonium
  - Ammonia
- $k_L a_{\text{CO}_2} = \alpha \cdot PIV^\beta \cdot vvm^\gamma$
- $\therefore \text{CO}_2 = \frac{\text{pH}}{k_L a_{\text{CO}_2} \cdot (\text{acid buildup}) \cdot (\text{base buildup})}$



# Conclusions

- S.U.B. enhancements lead to 3-4 fold increases in oxygen mass transfer from legacy vessels
- Scalability of S.U.B.s
- High degree of automation vital to perfusion success
- Achieved  $>260 \times 10^6$  cells/mL in cell culture
  - Room for more oxygen mass transfer at all scales
- Highly automated process controlling DO, pH,  $dCO_2$
- Predictive mass transfer models based on culture data
- Further work to be done optimizing media and exchange rates