

# Plutonium Isotopic Analysis System for Plutonium Samples Enriched in $^{238}\text{Pu}$ in EP 60/61 and Fuel-Clad Containers

## Volume 3 Part 2: Software Listings

Wayne D. Ruhter

July 1, 1991

Prepared under the auspices of WestinghouseHanford Company in partial  
fulfillment of Memorandum Purchase Order MJL-XMV-012823.

LAWRENCE LIVERMORE NATIONAL LABORATORY  
University of California • Livermore, California • 94551



Available to DOE and DOE contractors from the Office of Scientific and Technical Information  
P.O. Box 62, Oak Ridge, TN 37831 Prices available from (615) 576-8601, FTS 626-8401

Available to the public from the National Technical Information Service, U.S. Department of Commerce  
5285 Port Royal Rd., Springfield, VA 22161 •

MASTER

MAR 16 1992

DISTRIBUTION OF THIS DOCUMENT IS UNLIMITED

## **DISCLAIMER**

**This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency Thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.**

## **DISCLAIMER**

**Portions of this document may be illegible in electronic image products. Images are produced from the best available original document.**

#### DISCLAIMER

Work performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under contract number W-7405-ENG-48.

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial products, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

## PREFACE

This software manual is addressed to the Westinghouse Hanford's Radioisotope Power Systems Facility personnel (programmers and supervisors) who maintain the software on the Pu-238 isotopic analysis system. The document is divided into two parts. Part 1 describes the computer codes that control the system, analyze the spectral data, and determine the relative plutonium abundances. Part 2 contains the software listings of the analysis codes.

In addition to this manual, two other manuals complete the description of the Pu-238 isotopic analysis system. The Pu-238 Isotopic Analysis Users Manual (UCRL-CR-107467, Vol. 1) explains how to operate and calibrate the system. The Pu-238 Isotopic Analysis Hardware Manual (UCRL-CR-107467, Vol. 2) describes the system components and provides service and maintenance information. In that volume are detailed instructions on required settings for the nuclear instrument modules.

All commercial items mentioned in this manual are assumed to be functioning correctly for the purposes of system operation. Users are referred to individual equipment manufacturers' manuals for details of operation, troubleshooting, and maintenance of this commercial equipment.

-----  
NOTE: If any of the Pu-238 Isotopic Analysis codes are subsequently changed, users are asked to print out the changed listing and replace the listing now in this manual. In that way, this document will be kept current.  
-----

CONTENTS

PART 2

Chapter 11

11. SUBROUTINE LISTINGS . . . . . 11- 1

Chapter 12

12. PULIB LISTINGS . . . . . 12- 1

Chapter 13

13. CHANGE PROGRAM LISTINGS . . . . . 13- 1

C\*\*\*\*\*

C ACQR38 (ACQR90 subroutine for PU238 instrument)

C This routine is a special adaptation of the ACQR90 routine for the Pu238 isotopic analysis instrument. Routine is instrument specific to Canberra Series 90 data collection. Acquire data on the Series 90 Set up to send 'AS' (suspend acquisition) and 'AR' (resume acquisition) to VAX

C SSJ 20-JUL-79

C Rev. 19-jun-84 lc

C Rev. 25-sep-84 cms

C Rev. 26-OCT-84 MJC

C Rev. 31-OCT-84 MJC

Split ACQR90 and CMPLTN subroutines Deleted LOGICAL \* 1 TKEY & added it CONTRL.CMN

C Rev. 21-DEC-84 WDR

C Rev. 03-JAN-85 MPK

Modified for the solid isotopic instrument Restored DATTIM option argument Added VAX comm for suspend and resume

C Rev. 15-JAN-85 MPK

C Rev. 21-JAN-85 MPK

C Rev. 04-FEB-85 MPK

Disabled T key during MC runs Check VAXUP in comm calls Add sleep time in ADC/terminal polling loop Put END= and ERR= returns in terminal read after suspend, to prevent operator from bombing the program

C Rev. 08-MAR-85 MPK

C Rev. 27-JAN-86 WDR

Include CVAX common Test for rotation before turning on ADCs to make sure it is rotation is enabled

C Rev. 02-APR-86 WDR

C Rev. 05-JUL-87 RDP

C Rev. 15-JUL-87 RDP

C Rev. 18-AUG-87 RDP

Modified for Micro VAX II Modified for uVAX II w/ ADAC Drivers Using Single Key Input for VAX Included calls to control Load/Unload, Rotation, and SET lights.

C REV. 04-SEP-87 RDP

Include a resume fcn, if assay is suspended due to rotation error.

C Rev. 19-Jan-88 WDR

C Rev. 30-Aug-89 WDR

C Rev. 02-May-90 WDR

C Rev. 08-May-91 WDR

Sleep call changed to real parameter Modified for Pu238 instrument Removed COMM calls for Pu238 instrument Added test for sample type in elevator motion check

C\*\*\*\*\*

C SUBROUTINE ACQR90

C-----

C\*\*\*\*\*

C Real Time Interrupt Variables for Keyboard Interrupt

C COMMON/RTI/ INKEY,INCHAR  
C INTEGER INKEY  
C CHARACTER INCHAR\*1

C LOGICAL\*1 EFLG, RFLG, ROTERR, ELERR  
C CHARACTER ACHAR\*1  
C INTEGER\*2 ON6,OFF6  
C INTEGER\*4 ISTAT,IMASK

C The following common block files are generic for all instruments

```

C
C   INCLUDE 'CASSAY.CMN'
C   INCLUDE 'CDATE.CMN'
C   INCLUDE 'CONTRL.CMN'
C   INCLUDE 'CLUNS.CMN'
C   INCLUDE 'CVAX.CMN'
C   INCLUDE 'ADAC.CMN'
C   INCLUDE 'CTIM.CMN'

```

```

C
C   The following common block files are instrument specific for
C   gamma-ray instruments

```

```

C       INCLUDE 'CONT90.CMN'
C       INCLUDE 'CSPECT.CMN'

```

```

C   Disable LOAD/UNLOAD sample capability; enable rotation

```

```

C       CALL ROTMON(RFLG,EFLG)
C       IMASK = 512
C       ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
C       CALL ERROR_CHECK(ISTAT,IOSB)

```

```

C   We are generating live data; set the live flag
C   Reset the termination flag

```

```

C       LIVE = .TRUE.
C       TKEY = .FALSE.

```

```

C   Clear the time
C   Erase the existing spectrum unless this is a non-first cycle
C   of a cumulative run

```

```

12   CONTINUE
C       CALL CTIM(IADC,,IERR)
C       CALL CHK90('CTIM ',IERR)
C       CALL CDAT(IADC,,IERR)
C       CALL CHK90('CDAT ',IERR)

```

```

C   Preset (clock) time
C   Start collecting

```

```

13   CONTINUE
C       CALL PSET(IADC,,1,1,TIMEPR,0.0,0,0)

```

```

C   Check to see that sample is rotating

```

```

C       RFLG = .TRUE.
C       EFLG = .TRUE.
C       CALL ROTMON(RFLG,EFLG)
C       TYPE *, 'EFLG IN ACQR38', EFLG, MATTYP(1)
C       IF (RFLG) GOTO 14

```

```

1003  TYPE 1003
C       FORMAT('/', ' !! Sample not rotating. Enable rotation !!')
C       TYPE 1004, 'Press ENTER to retry (type A and a <CR> to ABORT) ->'
C       READ (5,1006,END=13,ERR=13) ACHAR
C       IF ((ACHAR.EQ.'A').OR.(ACHAR.EQ.'a')) THEN
C           TYPE *, ' Abort Process -- Problem with Rotation Motor'
C           ABORT = .TRUE.
C           GOTO 400

```



```

ENDIF
1004 FORMAT( #, ' ', A)
1006 FORMAT( B0A1)
14 CONTINUE
IF (.NOT. EFLG .AND. MATTYP(1).EQ. 'E') THEN
TYPE 1113
1113 FORMAT( //, ' !! Sample elevator not working. Enable elevator !!')
TYPE 1004, 'Press ENTER to retry (type A and a <CR> to ABORT) ->'
READ (5, 1006, END=13, ERR=13) ACHAR
IF ((ACHAR.EQ. 'A').OR.(ACHAR.EQ. 'a')) THEN
TYPE *, ' Abort Process -- Problem with Sample Elevator'
ABORT = .TRUE.
GOTO 400
ENDIF
ENDIF
16 CONTINUE
CALL CMPLTN
C
C Now turn on ADCs
C
15 CONTINUE
CALL COLL(IADC, , IERR)
CALL CHK90('COLL ', IERR)
C
C Check the status of the ADC
C If it is idle, acquisition has completed
C
50 I = INQADC(IADC, ISTAT, MEMSIZ, IGAIN, IRANGE, IOFF, INPUTS)
CALL CHK90('INQADC', I)
IF (I .EQ. 0) GO TO 100
CALL SLEEP(0.5)
C
C Check sample rotation and elevator; if stopped, terminate measurement
C
ROTERR = .FALSE.
ELERR = .FALSE.
CALL ROTMON(RFLG, EFLG)
IF (RFLG .AND. EFLG) GO TO 51 !No problem
CALL DATTIM(1) !Rotation stopped; suspend meas.
IF (.NOT. RFLG) THEN
TYPE 1005, NOWTIM
1005 FORMAT(' SAMPLE NOT ROTATING. MEASUREMENT SUSPENDED AT ', AB//)
ACHAR = 'S'
ROTERR = .TRUE. !Suspended by Rotation Error
GOTO 55
ELSE
C
TYPE *, EFLG, MATTYP(1)
IF (MATTYP(1) .EQ. 'E') THEN
TYPE 1016, NOWTIM
1016 FORMAT(' SAMPLE ELEVATOR STOPPED. MEASUREMENT SUSPENDED AT ', AB//)
ACHAR = 'S'
ELERR = .TRUE. !Suspended by elevator problem
GOTO 55
ENDIF
endif
C
Acquisition not complete
C Check for quit (Q) or suspend (S) or terminate (T) key on console
C
51 CONTINUE

```

```

IF (INKEY.NE.O) THEN
    ACHAR = INCHAR      ! SET KEY VALUE
    INKEY = 0          ! RESET INKEY FLAG
ELSE
    GOTO 50
ENDIF

```

```

L
C Terminal input seen - check whether action required
C

```

```

52 CONTINUE
IF (ACHAR .EQ. 'Q') GOTO 55
IF (ACHAR .EQ. 'S') GOTO 55
IF (ACHAR .NE. 'T') GOTO 50

```

```

C
C If we are here, T key was seen. Allow termination unless MC
C run is in progress.
C

```

```

IF (.NOT. (MP.OR.MB)) GO TO 55
TYPE 1002

```

```

1002 FORMAT(/20X,' T option ignored during meas. control;'
1          'use Q to quit.//')
GOTO 50

```

```

C
C Quit (Q) or suspend (S) or terminate (T) key recognized
C Turn on the keyboard and stop the analyzer
C

```

```

55 CONTINUE
CALL ADCOFF(IADC,,IERR)
CALL CHK90('ADCOFF',IERR)
IF (ACHAR .EQ. 'Q') GOTO 70
IF (ACHAR .EQ. 'T') GOTO 60
IF (.NOT. RFLG .OR. (.NOT. EFLG .AND.
+ MATTYP(1).EQ.'E')) GOTO 56

```

```

C
C Suspended by S key
C

```

```

TYPE 1008,NOWTIM
1008 FORMAT(/20X,' Acquisition suspended by S key at ',AB//)
TYPE *

```

```

56 TYPE 1058

```

```

1058 FORMAT('!To resume, press RETURN (press A to abort) ->')

```

```

57 IF (INKEY.NE.O) THEN          !Wait for Key Hit
    INKEY = 0
    IF ((INCHAR.EQ.'A').OR.(INCHAR.EQ.'a')) THEN
        IF (.NOT. RFLG) THEN
            TYPE*, ' Abort Process -- Problem with Rotation Motor'
        ELSE
            TYPE*, ' Abort Process -- Problem with Sample Elevator'
        ENDIF
        ABORT = .TRUE.
        GOTO 75          !exit acquisition and clean up
    ELSE
        GOTO 58          !resume
    ENDIF

```

```

ELSE
    CALL ROTMON(RFLG,EFLG)
    IF (RFLG) THEN
        CALL SLEEP(2)          !Wait a sec to make sure
        CALL ROTMON(RFLG,EFLG)
    ENDIF

```

```

        IF (RFLG.AND.ROTERR) THEN      !If rotation has restarted
            ROTERR=.FALSE.
            GOTO 58
        ELSE
            IF ((EPLG.and. ELERR) .and. RFLG) THEN
                ELERR = .FALSE.
                GOTO 58
            ENDIF
            ENDIF
            GOTO 57
        ENDIF
58      CONTINUE                        !Resume
1001   FORMAT(A1)
        CALL DATTIM(1)
C      IF (VAXUP) CALL COMM('AR')      !*Communications - assay resumed
        TYPE 1007,NOWTIM
1007   FORMAT(//20X,' Acquisition resumed at ',A8//)
        GO TO 15
C
C   Terminated by T key
C   Stop cycling
C
        60   CALL DATTIM(1)
            TYPE 1019,NOWTIM
1019   FORMAT(//20X,' Acquisition terminated by T key at ',A8//)
            TKEY = .TRUE.
            STPCYC = .TRUE.
            GO TO 200
C
C   Aborted
C
        )   ABORT = .TRUE.
            TYPE 1000,NOWTIM
1000   FORMAT(//20X,' Acquisition aborted by Q key at ',A8//)
75     CONTINUE
C   Turn off rotation, enable LOAD/UNLOAD
C   Enable LOAD/UNLOAD
400    IMASK = 512
        CALL = BCLRW(IMASK,H_DATA,IOSB,,)
        CALL ERROR_CHECK(ISTAT,IOSB)
        GO TO 500
C
C
100    CONTINUE
C
C   Acquisition stopped
C   Retrieve the live (0) and clock (1) times
C
- 200  CALL FROMCA(IADC,2,0,SPCTRM,IERR,)
        CALL CHK90('FROMCA',IERR)
        TIMELV = SPCTRM(1)
        TIMECL = SPCTRM(2)
C
C
500    RETURN
        END

```

```
C*****
C
C ADATE
C
C     Set the date and time for assay
L     16-Jan-85 MPK
C
C*****
C
C     SUBROUTINE ADATE
C
C     INCLUDE 'CDATE.CMN'
C
C     CALL DATTIM(1)
C     ASDATE=NOWDAT
C     ASTIME=NOWTIM
C     RETURN
C     END
```



```

C*****
C
C AEND
C
C      Do VAX communication at end of assay.
C      Message type 'AE' is sent unless the assay was terminated
L      early, in which case 'AT' is sent. 'AQ' is not sent,
C      because it is generated elsewhere.
C
C      Enable sample UNLOAD and disable sample rotation.
C
C      16-Jan-84 MPK
C      25-Aug-87 WDR
C      11-Nov-87 WDR  Added test for autocytle
C
C*****
C
C      SUBROUTINE AEND
C
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CVAX.CMN'
C      INCLUDE 'CASSAY.CMN'
C
C      INCLUDE 'ADAC.CMN'
C
C      Test for autocytle and number of cycles.  If not done, don't do
C      do the following
C
C          IF (.NOT.AUTOCY)GO TO 5
C          IF (ICYCLE.LT.NUMCYC)GO TO 7
C
C      Turn off rotation, enable LOAD/UNLOAD
5      IMASK = 4096
C      Send HIGH to reset
C      ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
C      CALL ERROR_CHECK(ISTAT,IOSB)
C      Return to LOW
C      ISTAT = BCLRW(IMASK,H_DATA,IOSB,,)
C      CALL ERROR_CHECK(ISTAT,IOSB)
C      IMASK = 2048
C      Send HIGH to reset
C      ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
C      CALL ERROR_CHECK(ISTAT,IOSB)
C      Return to LOW
C      ISTAT = BCLRW(IMASK,H_DATA,IOSB,,)
C      CALL ERROR_CHECK(ISTAT,IOSB)
C Turn Off Rotation
C      IMASK = 1
C      ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
C      CALL ERROR_CHECK(ISTAT,IOSB)
C Enable LOAD/UNLOAD
C      IMASK = 512
C      CALL = BCLRW(IMASK,H_DATA,IOSB,,)
C      CALL ERROR_CHECK(ISTAT,IOSB)
7      RETURN
C      END

```

C\*\*\*\*\*

C  
C ASOPT

C Responds to the assay option selected by the operator  
C from the menu. Calls to perform the assay, print the  
L results, and write data to files are performed here.  
C The data is sent to the VAX from this routine also.

C  
C CMS 10-MAY-84  
C REV.11-SEP-84 CMS  
C REV 11-NOV-84 WDR Removed WRDAT since already done for LLNL  
C instruments. Livermore systems always write  
C data to disk.  
C Rev 29-Nov-84 MPK Replace list-directed with  
C formatted I/O  
C Rev 04-Jan-85 MPK Remove VAX communication that is  
C handled in ASSAY - change rest  
C to meet specs  
C Rev. 18-JAN-85 MPK Add calls to ADATE and AEVAX  
C Rev. 21-JAN-85 MPK Check VAXUP in comm calls  
C Rev. 12-FEB-85 MPK Redo logical flags  
C Rev. 20-FEB-85 MPK Added assay type 'A','S' for log  
C Rev. 29-APR-85 MPK Send 'MS' message to VAX  
C Rev. 02-APR-86 WDR Modified for Micro VAX II  
C Rev. 13-Jul-87 RDP Modified for uVax II

C\*\*\*\*\*

C  
C SUBROUTINE ASOPT

C The following common block files are generic for all instruments

C  
C INCLUDE 'CDIAG.CMN'  
C INCLUDE 'CONTRL.CMN'  
C INCLUDE 'CLUNS.CMN'  
C INCLUDE 'CBACKG.CMN'  
C INCLUDE 'CVAX.CMN'  
C INCLUDE 'CASSAY.CMN'  
C INCLUDE 'CFILES.CMN'  
C INCLUDE 'CLOG.CMN'  
C INCLUDE 'CMSG.CMN'

C  
C The following common block files are instrument specific for  
C gamma and x-ray instruments

C  
C INCLUDE 'CSPECT.CMN'

C  
C ASFLG = .TRUE.  
C RTYP = ' Assay '

C  
C Check that measurement control has been performed  
C within the required time limits and without  
C error

C  
C CALL MCCHK  
C IF(ABORT) GO TO 500

C  
C Get set up parameters from the dialog

```

C  Filename returned in common if writing to disk
C
      CALL ADATE
      CALL DIALOG
      IF (ABORT) GOTO 500

L  Announce assay to VAX and to operator
C
      TYPE 1000
1000  FORMAT (/' Assay beginning')
C      IF (VAXUP) CALL COMM('SA')      !***communications - start assay
      100  CALL ASSAY
          CALL AEND
          IF(ABORT)GO TO 500

C
C  Print results of assay - if diagnostic error flag is set
C  write result with error message
C
      CALL PRINT1      !instrument specific

C
C  Check if VAX is up - send assay results to VAX
C  If VAX is down and writing raw data - write a temporary file
C
      IF (.NOT.VAXUP) GOTO 400
C  300  CALL COMM('RA')      !***communications-send results
      300  CONTINUE
C      IF(VAXUP.AND.RAWDTA)CALL COMM('RD')!***communications-send raw data
          ICHAR = 0

C
C  Write assay results and ancillary information to the assay
C  log file 'ASLOGF.INS'

400  LVFLAG = VAXUP
      ASTYPE(1) = 'A'
      ASTYPE(2) = 'S'
      CALL WRITLG

500  RETURN
      END

```



```
C
C*****
```

```
C ASSAY38 (ASSAY routine for Pu238 isotopic instrument)
```

```
C Routine performs assay runs for assay,background,
C measurement control, calibration, from disk options.
C It then calls the analysis routines and the data
C diagnostic routines.
```

```
C CMS 26-APR-84
```

```
C Rev. 05-OCT-84 CMS ADDED ABORT=.TRUE. IF ERROR IN OPENING FILE
```

```
C Rev. 22-OCT-84 MJC Changed 'BACKGR.DAT' to BKFIL in opening DKLUN1
```

```
C Rev. 01-NOV-84 wdr adapted to LIVERMORE needs for solid isotopic
C instrument
```

```
C Rev. 17-DEC-84 WDR modified for solid isotopic instrument
```

```
C Rev. 03-JAN-85 MPK Restored DATTIM option argument
```

```
C Revised COMM calls to follow specs
```

```
C Rev. 18-JAN-85 MPK Removed calls to send 'AE' and 'AT'
C messages, which will be handled elsewhere.
C Sending 'AQ' is still done here.
C Removed dead code at 501.
```

```
C Rev. 13-FEB-85 MPK Redo logical flags
```

```
C Rev. 20-FEB-85 MPK Remove Ref to CS90.CMN
```

```
C Continue line that was being cut off
```

```
C Rev. 02-APR-86 WDR Modified for Micro VAX II
```

```
C Rev. 30-Aug-89 WDR Modified for Pu238 instrument
```

```
C*****
```

```
C SUBROUTINE ASSAY
C LOGICAL * 1 IANS
```

```
C The following common block files are generic for all instruments
```

- INCLUDE 'CONTRL.CMN'
- INCLUDE 'CDATE.CMN'
- INCLUDE 'CLUNS.CMN'
- INCLUDE 'CDIAG.CMN'
- INCLUDE 'CASSAY.CMN'
- INCLUDE 'CTIM.CMN'
- INCLUDE 'CFILES.CMN'
- INCLUDE 'CLOG.CMN'
- INCLUDE 'CVAX.CMN'
- INCLUDE 'CMCONS.CMN'
- INCLUDE 'CBACKG.CMN'

```
C The following common block files are instrument specific for
C gamma ray instruments
```

- INCLUDE 'CSPECT.CMN'
- INCLUDE 'CONT90.CMN'
- INCLUDE 'CONTRS.CMN'

```
C Set spectrum pointer to guarantee fresh data
```

```
C IFIRST = -500 !instrument specific
```

```
C Check if assay is from disk - open file for analysis
```

```

C
IF(AC) GO TO 200          !current data in MCA
IF(.NOT. FROMDK) GO TO 130
OPEN(UNIT=DKLUN,NAME = FILNAM,ACCESS = 'DIRECT',TYPE = 'OLD',
1   RECORDSIZE = 128,ERR = 504)
GO TO 200
504 CALL ERRMSG(1,FILNAM)
   ABORT = .TRUE.
   GO TO 500

C
C Read data from disk file or read current data
C to retrieve information needed for analysis
C
200 CALL DATA90(0,4,SUM)    !DATA90 is instrument specific
    TIMELV = SPCTRM(1)      !Livetime
    TIMECL = SPCTRM(2)      !Real time
    IF(AC ) GO TO 180
    GO TO 184                !Assay data from disk

C
C Do the actual data collection
C
130 CALL DATTIM(1)
    CALL ACQR90              !collect data - instrument specific
    CALL DATTIM(1)
    IF(.NOT. ABORT)GO TO 180
C   IF (VAXUP) CALL COMM('AQ')    !*communications-assay quit
    GO TO 500
180 CALL WRDAT          ! write data to disk file - write ADC #1
    CLOSE(UNIT=DKLUN)
184 CALL P238GRP        !instrument specific routine
    CLOSE(UNIT=DKLUN)
500 RETURN
    END

```

```
C*****  
C  
C ASYP238  
C  
C main program for the SRP_FU238 program  
C  
C MPK 21-Nov-84  
C*****  
C  
C PROGRAM ASYP238  
C CALL EXEC  
C END
```

C  
C\*\*\*\*\*

C AUTCYC

C Responds to the operator request for assay in autocyte  
L mode where assays are collected repeatedly for the number  
C of times input in the dialog mode. This routine also  
C responds to autocyte from disk

C CMS 10-MAY-84

- C Rev. 25-SEP-84 CMS
- C Rev. 16-NOV-84 WDR modified standard deviation calculation
- C Rev. 04-JAN-85 MPK changed VAX COMM to conform to specs
- C Rev. 18-Jan-85 MPK Add calls to ADATE and AEVAX
- C Add logic to handle 'T' option
- C Rev. 22-Jan-85 MPK Check VAXUP in COMM calls
- C correct statistics
- C Rev. 12-Feb-85 MPK Redo logical flags
- C Rev. 20-Feb-85 MPK Declared XSSQ, XMEAN, XSUM as REAL\*8
- C Added assay type 'A', 'S' for log
- C Rev. 02-May-90 WDR Removed AEVAX and COMM calls

C\*\*\*\*\*

```
C  
C SUBROUTINE AUTCYC  
C REAL*8 XSSQ, XMEAN, XSUM  
C LOGICAL * 1 SUMMRY
```

C The following common block files are generic for all instruments

```
C  
C INCLUDE 'CONTRL.CMN'  
C INCLUDE 'CVAX.CMN'  
C INCLUDE 'CFILES.CMN'  
C INCLUDE 'CDIAG.CMN'  
C INCLUDE 'CBACKG.CMN'  
C INCLUDE 'CLUNS.CMN'  
C INCLUDE 'CASSAY.CMN'  
C INCLUDE 'CLOG.CMN'  
C INCLUDE 'CQSUMS.CMN'
```

C The following common block files are instrument specific for  
C gamma ray instruments

```
C  
C INCLUDE 'CSPECT.CMN'  
C  
C RTYP = 'Autocycl'  
C AUTOCY= .TRUE.  
C ASFLG = .TRUE.  
C  
C XMEAN = 0.0  
C STDEV = 0.0  
C VAR = 0.0  
C XSQ = 0.0  
C ICYCLE = 0
```

C Check that measurement control has been performed within  
C the required time and without error

```

C
    IF(FROMDK) GO TO 100
    CALL MCCHK
    IF(ABORT) GO TO 500
100  CALL DIALOG
    IF (ABORT) GO TO 500
    CALL NEWNAM(0)
C
C Collect assays for the number of times returned from dialog
C
200  DO 10 I = 1, NUMCYC
    CALL ADATE
    IF (VAXUP) CALL COMM('SA')          !*communications - start assay
    CALL ASSAY
    IF(ABORT)GO TO 500
    ICYCLE = I
C
C If autocycle from disk - don't send results to VAX
C
    IF(FROMDK) GO TO 30
C
C Write assay log entry
C
35  LVFLAG = VAXUP
    ASTYPE(1) = 'A'
    ASTYPE(2) = 'S'
    CALL WRITLG
C
C Print results
C
30  CALL PRINT1                          !instrument specific routine
    IF (.NOT.TKEY) GO TO 32
    ICYCLE = ICYCLE - 1                  !Reset ICYCLE
    GO TO 15
C
C Accumulate statistics
C
32  XMEAN = XMEAN + ASANS
    XVAL(I)= ASANS                      !save answer for standard deviation calc
C
C Create next filename
C
    IF(I .NE. NUMCYC) CALL NEWNAM(1)
10  CONTINUE
C
C Calculate average and standard deviation from assay results
C
15  IF (ICYCLE.LE.1)GO TO 500
    XMEAN = XMEAN/ICYCLE
    DO 50 J = 1,ICYCLE
    XSQ = XSQ + (XVAL(J) - XMEAN)**2
50  CONTINUE
    STDEV = SQRT(XSQ)/(ICYCLE-1)
C
    NUMCYC = ICYCLE
    SUMMRY = .TRUE.
    CALL HEADER(SUMMRY)
    WRITE(LOUT, 1000)ICYCLE,XMEAN, STDEV
1000 FORMAT(' Average from ',I3,' runs = ',F10.4,
1      ' standard deviation = ',F10.4)

```

500

CALL AEND  
RETURN  
END

C C C C C C C C C

```
SUBROUTINE CALEFF(EFCOF,EPK,EMVLOG,DEFLOG)
```

```
DIMENSION EFCOF(2)
```

```
IF(EPK .GT. EFCOF(12))GO TO 20 !Test vs. Xover energy.
```

```
DEFLOG = EFCOF(7)
```

```
DO 10, I = 8,11
```

```
DEFLOG = DEFLOG + EFCOF(I) * EMVLOG**(I-7)
```

```
GO TO 50
```

```
DEFLOG = EFCOF(1)
```

```
DO 40, I = 2,6
```

```
DEFLOG = DEFLOG + EFCOF(I) * EMVLOG**(I-1)
```

```
RETURN
```

```
END
```

```
C SUBROUTINE CALGPM (EPK,VALUE,PCERR)
C Last edit 6-Mar-85 (J.B.N.)
C Last edit 27-Aug-85 (R.G.) To self- absorption correction
C Edit 12-Sep-88 (r.g.) SSMU
C ~~~~~
```

```
INCLUDE 'GRPLGCM.PMS'
```

```
INTEGER*2 MTRL
DATA MTRL/'GE'/
```

```
ELOG = ALOG (.001 * EPK)
SMU = 0.0
GAMAS = VALUE
IF (NF1 .EQ. 0)GO TO 10
CALL CALMU (NF1,SMU,EPK,ELOG,MTRLZ,CMPOS,NSMBL,EDGE,COFMU)
```

```
10 SSMU = SMU
C IF (RUNFLG(4) .EQ. 0)THEN
C IF (((DCNST(1).EQ.0) .AND. (DCNST(12).GT.0)) .OR.
1 (GEOM - DCNST(15) .LE. 0))THEN
VALUE = 1.0 !Skip efficiency corrections, but allow
IF (SMU .GT. .001) GAMAS = GAMAS *SMU / (1. -EXP(-SMU))
SSMU = 2.*SMU !To increase the error assoc. with value
GO TO 20 !absorber corrections.
ENDIF
```

```
C [DONT change DEPTH in common-- GPM will use 1/2 depth.]
CALL GPM (EPK,ELOG,GEOM,SURFC,DEPTH,SMU,VALUE,GAMAS,DCNST)
IF (EPK .GT. 500.)GO TO 20
```

```
C CORRECT FOR GE "DEAD LAYER" ATTENUATION.
```

```
CALL CALMU (1,XMU,EPK,ELOG,MTRL,SMU,NSMBL,EDGE,COFMU)
GAMAS = GAMAS * EXP (XMU)
20 IF (NF2 .EQ. 0)GO TO 30 !No absorber present.
ATN = 0.0
```

```
CALL CALMU (NF2,AMU,EPK,ELOG,NABS,ABSRB,NSMBL,EDGE,COFMU)
```

```
C ADJUST ABSORBER ATTENUATION OF PHOTONS TO ACCOUNT FOR GAMMAS
C TRAVERSING THE ABSORBER AT SKEWED ANGLES.
```

```
C WRITE (*,*)'AMU',EPK,AMU
C AMU = AMU * VALUE
C IF (AMU .LT. 6.9)ATN = EXP (AMU)
C GAMAS = GAMAS * ATN
```

```
C CALCULATE ERROR IN ATTENUATION USING AN ERROR OF 2%
C IN ABSORBER THICKNESSES
```

```
C ADRA = 2.*AMU
C ADRA = 1.3
C PCERR = SQRT(PCERR**2 + SSMU**2 + ADRA**2)
```

```
30 VALUE = GAMAS
RETURN
END
```



SUBROUTINE CALIB

CC

C GRPANEL sub to calculate GAIN, ZERO, & standard reference points  
corrected for linearity. Called once for each grouping.

By J.B.Niday Last edit 25-Sept-87 (JBN)  
CC

INCLUDE 'GRPLGCM.PMS'  
  
INCLUDE 'CONTRL.CMN'  
INCLUDE 'CLUNS.CMN'  
INCLUDE 'CDIAG.CMN'

DIMENSION XC(2),TREFCH(2)  
DATA LP/-1/  
DATA LP/6/

IF(NUSHPS .NE. 0) THEN !Parameter specification has precedence.  
NOC = 1  
EXP1 = ATAIL  
GO TO 10  
ENDIF

IF(SHAPC(1) .NE. 0)THEN !Shape constants on file.  
NOC = 0  
TSHPS = SHAPC(5)  
ELSE  
NOC = 2  
ENDIF

10 BTAIL = TSHPS \* GAIN  
EXP2=BTAIL  
GAN = GAIN  
ZRO = ZERO

IF (LONGPR)WRITE(LOUT,1005)KTGRP+1,EST,EEND  
1005 FORMAT(///' Group',I3,' includes',F7.1,' to',F7.1,' keV')

CC

C NGAIN = 1 is a (byte) flag indicating fixed gain; otherwise, two  
C internal reference peaks are to be used to establish energy scaling.  
C If the first reference peak's energy = 0, this is interpreted to request  
C calculation of the energy of channel zero from ZERO & GAIN

IF(NGAIN .GT. 0)GO TO 150  
  
IF((STD(1).EQ.OLSTD1) .AND. (STD(2) .EQ. OLSTD2))GO TO 150  
  
OLSTD1 = STD(1)  
OLSTD2 = STD(2)

1010 IF (LONGPR)WRITE(LOUT,1010)STD(1),STD(2)  
+ FORMAT(' Recalibrate Gain using energy standards at',F9.3,  
' and',F9.3,' keV')  
GSQ = GAIN \* GAIN

DO 100 J=1,2

```
TREFCH(J) = 0.0
XC(J) = 0.
IF(STD(J) .EQ. 0)GO TO 100
```

```
IF(NOC .EQ. 0) THEN
    TFWHM = SQRT(SHAPC(1) + SHAPC(2) * STD(J))
ELSE
    TFWHM = FWHM
ENDIF
```

```
IF(TFWHM .EQ. 0) THEN
    T = 0.3
    IF(GAIN .LT. 0.15) T = 0.03 !For low-noise detectors.
    TFWHM = SQRT(T +.002*STD(J)) !Estimate FWHM.
ENDIF
```

```
ALFA=-2.7726*GSQ/TFWHM**2
HW = SQRT(.69315/-ALFA)
XLIM = 20. * HW
REF=(STD(J)- ZERO)/GAIN
```

```
C Look for calibration peaks using the approx. GAIN & ZERO
DO 90 K = 1,6
IF(XLIM .GT. (MXDP-10))XLIM = MXDP-10
NCH = XLIM + 0.5
ST1 = REF - XLIM/2.0
IST1 = ST1 + 0.5
ST1 = IST1
1014 IF (LONGPR)WRITE(LOUT,1014)NCH,IST1,STD(J)
    + FORMAT(' CALIB searching',I4,' chans from',I5,' for ',
    + F7.1,' keV peak')
```

```
15 CALL RDCHNS(Y,IST1,NCH,INPDEV)
NS = REF - ST1 - 2.0 * HW + 0.5
IF(NS .LT. 1)NS = 1
NE = REF - ST1 + 2.0 * HW + 0.5
IF(NE .GT. NCH)GO TO 90
1015 IF(K .GT. 1)TYPE 1015,NS+IST1-1,NE+IST1-1
    FORMAT(' Raise limits to chans:',I5,' and',I5)
```

```
C Locate highest & lowest points of STD peaks.
17 CALL MAXVAL (NS,NE,JPK,YMAX,Y)
IF((JPK .GT. NS) .AND. (JPK .LT. NE))GOTO 25
HW = 1.5 * HW
GO TO 15
25 CALL MINVAL (1,JPK,MS,BG1,Y)
CALL MINVAL (JPK,NCH,ME,BG2,Y)
YTST = YMAX * 0.9
IF((JPK .LE. MS) .OR. (JPK .GE. ME))GO TO 90
IF((YTST .GT. BG1) .AND. (YTST .GT. BG2))GO TO 35
GO TO 80
```

```
35 CALL AVE (MS,ME,SUMY,AV,Y)
YS = 0.0
DBG = BG2 - BG1
DO 40 I = MS,ME
YS = YS + Y(I)
Y(I) = Y(I) - BG1 - DBG * YS/SUMY
IF(Y(I) .LE. 0)Y(I) = 0.1
40 CONTINUE
```

```

YTST = Y(JPK)*0.5
M = 0
DO 45, I = JPK+1,ME
M = M+1
IF(Y(JPK+M) .LE. YTST)GO TO 50
CONTINUE
GO TO 80

50 IF((JPK-M) .LT. MS)GO TO 80

55 IF(NUSHPS .NE. 0)GOTO 60
SG = SHAPC(1)/GSQ + SHAPC(2) *(ST1/GAIN +ZERO/GSQ)+.462
ALFA=-2.7726/SG
EXP1=EXP(SHAPC(3) +STD(J) *SHAPC(4))

C DETERMINE PEAK POSITIONS OF REFERENCE PEAKS

60 CALL GFIT(JPK,M,ALFA,EXP1,EXP2,PKPO,PHT,ER,Y)
IF(ER .EQ. -2.)GO TO 78

TREFCH(J) = PKPO + ST1 - 1.
CALL LINCOR(CLIN,TREFCH(J),XC(J))

GO TO 100

78 TYPE 1078,ER
1078 FORMAT(' Error in GFIT; flag =',F4.0)

80 TYPE 1080,(JPK+IST1-1),NCH,IST1
1080 FORMAT(' Reject peak at chan',I5,' in',I4,' chan grp',
+ ' beginning at',I5)
XLIM = XLIM * 1.5 !Widen search limits & retry.
90 CONTINUE

TYPE 1090,KTGRP+1
IF(IDFLG(3) .EQ. 0)WRITE(LP,1090),KTGRP+1
1090 FORMAT(' CALIBRATION FAILURE -- ON Group #',I3,
+ '; will use the preexisting GAIN')
GO TO 150

100 CONTINUE

TMP1 = STD(1)
IF(TMP1 .EQ. 0) TMP1 = ZERO !Use E of chan 0.
GAIN = (STD(2) - TMP1) / (TREFCH(2) - TREFCH(1))
GAN = (STD(2) - TMP1) / (XC(2) - XC(1))
IF(STD(1) .NE. 0) THEN
ZERO = STD(1) - GAIN * TREFCH(1)
ZRO = STD(1) - GAN * XC(1)
ZROC = TREFCH(1) - STD(1)/GAIN
ENDIF
IF ((ABS(ZROC) .GT. 5.) .AND. (ABS(GAIN-0.25).GT. .0025))THEN
WRITE (LOUT,1140)
1140 FORMAT(' ***WARNING: Zero-Channel intercept or gain incorrect.
+ Please center 152.7-keV peak at 612 and the 766.4-keV peak at
+ channel 3063. Check stabilizers')
DGFLB(2) = .true.
ABORT = .TRUE.
END IF

```

```
IF(INTFLG .NE. 0)INTFLG = INTFLG + 1      !Flag the new GAIN
```

```
150 ST = (EST - ZERO) / GAIN
END = (EEND - ZERO) / GAIN
IST = ST + 0.5
IEND = END + 0.5
C Correct for system nonlinearity by correcting the reference
C energy value
REFCH = (IST + IEND) / 2.
CALL LINCOR(CLIN,REFCH,PP)
REFEN = ZRO + GAN * PP
IF (LONGPR)WRITE(LOUT,1155) IST,EST,EEND,ST,REFEN,ZERO,GAIN
1155 FORMAT(' St-Ch Energy End energy start Refen Zero',
+ ' Gain'/I6,6F10.4)
C TYPE 9999, IST,EST,EEND,ST,REFEN,ZERO,GAIN
C9999 FORMAT (I5,6F10.4)
```

```
RETURN
```

```
END
```

```
C - - - - -
C CORRECT FOR ANY EFFECT OF SYSTEM NONLINEARITY ON ENERGY
C MEASUREMENT (CLIN(I) VALUES =0 FOR LINEAR RESPONSE).
```

```
SUBROUTINE LINCOR(CLIN,CH,X)
```

```
DIMENSION CLIN(2)
```

```
PP = CH
PKP = PP
DO 10, I=1,4
PKP = PKP * CH
10 PP = PP + CLIN(I) * PKP
X = PP
RETURN
END
```

C\*\*\*\*\*

C CHK90

C

C Routine is instrument specific for gamma ray instrument

C using Series 90

C Check error return on Series 90 call

C

C CALLING ARGUMENTS:

C ROUTIN = Routine from which the error was returned

C IERR = Error number

C

C SSJ 01-NOV-79

C REV. 14-May-84 lc

C\*\*\*\*\*

C

C SUBROUTINE CHK90 (ROUTIN,IERR)

C

C REAL\*8 ROUTIN

C

C IF(IERR .GE. 0) GO TO 100

C TYPE 1000,IERR,ROUTIN

1000 C FORMAT(//20X,'\*\*\*\* Error number ',I3,' returned from ',A6//)

C

100 C RETURN

C END

```

C*****
C CMPLTL      (CMPLTN subroutine for Livermore use)
C
C      Estimate completion time for the assay,
C      given the current time & the preset time
C
C      SSJ 25-OCT-82
C      Rev. 20-Jun-84 lc
C      REV. 26-OCT-84 MJC Split from ACQR90
C                          Deleted argument to CALL DATTIM
C      REV. 13-DEC-84 WDR Modified to determine start time of
C                          spectrum and pass results using CTIM.CMN
C      REV. 03-JAN-85 MPK Restored DATTIM option argument
C      Rev. 02-APR-86 WDR Modified for Micro VAX II
C*****

```

```

C
C      SUBROUTINE CMPLTN

```

```

C      INTEGER*2 HOURS,MINS
C      LOGICAL*1 ENDTIM(5)
C      CHARACTER*8 WHEN

```

```

C
C      The following common block files are generic for all instruments

```

```

C      INCLUDE 'CDATE.CMN'
C      INCLUDE 'CMONTH.CMN'
C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CTIM.CMN'

```

```

C      Get current date & time
C      CALL DATTIM(1)

```

```

C      Get start time in Julian days and seconds since midnite

```

```

C      CALL JULDAT(NOWDAT,NOWTIM,SDAT,STIM)

```

```

C      Get seconds since midnight

```

```

C      SEC = SECNDS(0.0)

```

```

C      Add preset time in seconds

```

```

C      Take modulo seconds in a day

```

```

C      SEC = SEC + TIMEPR
C      DAY = 24. * 3600.
C      NDAYS = SEC / DAY
C      SEC = AMOD(SEC,DAY)

```

```

C      Compute hours & minutes

```

```

C      HOURS = SEC / 3600.
C      MINS = AMOD(SEC,3600.) / 60.

```

```

C      Encode time & zero-fill spaces

```

```

C      ENCODE(5,1000,ENDTIM) HOURS,MINS

```

```

C      FORMAT(I2,':',I2)

```

```

C      IF(ENDTIM(1) .EQ. ' ') ENDTIM(1) = '0'

```

```

C      IF(ENDTIM(4) .EQ. ' ') ENDTIM(4) = '0'

```

C  
C Determine when - today, tomorrow, or how many days  
C

IF(NDAYS .EQ. 0) WHEN = ' Today'

IF(NDAYS .EQ. 1) WHEN = ' Tomorrow'

IF(NDAYS .GT. 1) WHEN = ' days'

IF(NDAYS .GT. 1) ENCODE(2,1004,WHEN) NDAYS

1004 FORMAT(I2)

C  
TYPE 1005, TIMEPR, NOWTIM, NOWDAT, WHEN, ENDTIM

1005 FORMAT(/' Starting a run of ',F7.0,' s at ',A,' ('A,')',

X 3X,'Ends ',A,' at ',5A1/)

C  
RETURN

END

```

C
C*****
C
C CURDTA
C
C     Responds to the operator request for assay from
C     current data in MCA
C
C     CMS 10-MAY-84
C     REV.11-SEP-84 CMS
C     Rev. 08-Jan-85  MPK     remove reset of AC flag
C                             (now done by INTFLG)
C     Rev. 02-Apr-86  WDR     Modified for Micro VAX II
C
C*****
C
C     SUBROUTINE CURDTA
C
C     The following common block files are generic for all instruments
C
C     INCLUDE 'CONTRL.CMN'
C     INCLUDE 'CASSAY.CMN'
C     INCLUDE 'CFILES.CMN'
C     INCLUDE 'CLUNS.CMN'
C     INCLUDE 'CVAX.CMN'
C     INCLUDE 'CBACKG.CMN'
C
C     The following common block files are instrument specific for
C     gamma ray and x-ray multichannel instruments
C
C     INCLUDE 'CSPECT.CMN'
C
C     RTYP = 'From MCA'
C     AC = .TRUE.
C     CALL DIALOG
C     IF(ABORT)GO TO 500
C     CALL ASSAY
C     IF(ABORT) GO TO 500
C     CALL PRINT1
C     RETURN
C     END
C
500

```

!instrument specific routine



```

C*****
C DATA90
C      Return Region of Interest data from the Series 90
C
C      Calling arguments:
C      NFIRST - First channel to retrieve
C      NLAST  - Last channel to retrieve
C      SUM    - Returns the sum of the channels
C
C
C      SSJ 03-JUL-79
C      REV. 10-May-84 LC
C      REV. 26-OCT-84 MJC  CHANGED LAST TO MEMSIZ FROM 4095
C      REV. 02-APR-86 WDR  Modified for Micro VAX II
C      REV. 05-JUL-87 RDP  Modified for uVAX II w/ ADAC Drivers
C*****
C
C      SUBROUTINE DATA90 (NFIRST,NLAST,SUM)
C
C
C      The following common block files are generic for all instruments
C
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CLUNS.CMN'
C
C      The following common block files are instrument specific for
C      (gamma ray) instruments
C
C      INCLUDE 'CSPECT.CMN'
C      INCLUDE 'CONT90.CMN'
C
C      Check the validity of the input channel numbers
C
C      IF(NFIRST .LT. 0 .OR. NLAST .GE. 4095) GO TO 500
C      NUMCHN = NLAST - NFIRST + 1
C      IF(NUMCHN .GT. 128) GO TO 500
C
C      Check if it is necessary to get new data or whether the data
C      exists in /cspect/ already
C
C      IF(NFIRST .GE. IFIRST .AND. NLAST .LE. IFIRST+255) GO TO 50
C
C      IF(FROMDK) GO TO 20
C
C      Get the data from the Series 90
C      Get 256 channels
C
C      LAST = NFIRST+127
C      IF(LAST .GT. MEMSIZ) LAST = MEMSIZ
C      TYPE 1200,NFIRST,LAST
C1200  FORMAT(' Debug - FROMCA',215)
C      CALL FROMCA(IADC, LAST, NFIRST, SPCTRM, IERR, )
C      CALL CHK90('FROMCA', IERR)
C      IFIRST = NFIRST
C      GO TO 50
C

```

```

C Get data directly from a disk file
C Check which records the first and last channels are in
C
20 IRECF = NFIRST/128 + 1
   IRECL = NLAST/128 + 1
   IF(IRECF .LT. IRECL) GO TO 25
C
C First and last channels are in the same record
C Read the whole record
C
   READ(DKLUN'IRECF) SPCTRM
   IFIRST = (IRECF-1) * 128
   GO TO 50
C
C Get half the data from the first record
C and half from the second
C
25 READ(DKLUN'IRECF) (DUM,I=1,128),(SPCTRM(I),I=1,128)
   READ(DKLUN'IRECL) (SPCTRM(I),I=129,256)
   IFIRST = (IRECF-1) * 256 + 128
C
C Sum the channels asked for
C
50 J = NFIRST - IFIRST + 1
   SUM = 0.0
   DO 60 I=1,NUMCHN
     SUM = SUM + SPCTRM(J)
60 J = J + 1
   GO TO 600
C
C Error in input parameters
500 TYPE 1000,IADC,NFIRST,NLAST
800 FORMAT(//20X'**** DATA90 INPUT ERROR - ADC#',I4/
X          '          NFIRST = ',I5,'          NLAST = ',I5,' ****'/)
C
600 CONTINUE
D WRITE(LOUT,1600) NFIRST,NLAST,SUM
D1600 FORMAT(2I6,F10.0)
RETURN
END

```

```

C*****
C DATTIM
C
C Retrieve, verify or change the system date & time
C
C SSJ 25-JAN-79
C REV. 13-SEP-83 LC
C Rev 02-Jan-85 MPK Include option to sync to date and
C time received from VAX and stored in
C the CVAX common. Adapt to the NSR/FMF
C generic program. Add automatic operator
C input of date and time if date and time
C are found to be unreasonable.
C
C Rev 14-Jul-87 RDP Modified for uVAX II
C Rev 18-Sep-87 RDP Remove Operator from initializing date
C and time.(IOPT=0)
C
C Options:
C
C 0 - Get initial date and time from operator input
C and give to monitor.
C
C 1 - Update date and time from monitor.
C
C 2 - Update and let operator verify or change.
C
C 3 - Sync to VAX date and time (note that these
C must have been set up by a call to VAX
C communication immediately before calling
C DATTIM).
C
C*****

```

SUBROUTINE DATTIM(OPTION)

INTEGER OPTION

The following common block files are generic for all instruments

```

INCLUDE 'CDATE.CMN'
INCLUDE 'CMONTH.CMN'
INCLUDE 'CVAX.CMN'

```

```

CHARACTER*9 INPSTR
CHARACTER*23 DATE_TIME

```

LOGICAL\*1 VFD

INTEGER\*4 STATUS, SYSTIM(2), SYS#BINTIM

INTEGER\*2 SYS#SETIME

EXTERNAL SYS#BINTIM, SYS#SETIME

```

DATA AMONTH /'JAN','FEB','MAR','APR','MAY','JUN',
1           'JUL','AUG','SEP','OCT','NOV','DEC'/
DATA ALOWDT /0./           !Julian date for 01-Jan-1900
DATA AHIDAT /36525./       !Julian date for 31-Dec-99

```

Set internal option variable. If it is zero, force operator input. If it is 3, sync to date and time from VAX.

IOPT=OPTION

```
IF (IOPT.LE.0) GOTO 100
IF (IOPT.GE.3) GOTO 300
VFD=.FALSE.
```

```
C
C Update alphanumeric date and time.
```

```
10 CALL DATE (NOWDAT)
CALL TIME (NOWTIM)
```

```
C
C If the date and time have been operator verified, return now.
```

```
IF (VFD) GOTO 500
```

```
C
C For option 1, return if date and time from VAX are reasonable;
C otherwise, make the operator input them by setting the internal
C option to 0.
```

```
IF (IOPT.NE.1) GOTO 40
CALL JULDAT (NOWDAT,NOWTIM,ADAT,ATIM)
VFD=.TRUE.
```

```
1 IF ((ADAT.GE.ALWDT).AND.(ADAT.LE.AHIDAT).AND.
(ATIM.GE.0.0).AND.(ATIM.LE.86399.0)) GOTO 500
```

```
IOPT=0
GOTO 100
```

```
40 CONTINUE
```

```
C
C Print out date and time for operator verification.
```

```
1000 TYPE 1000,NOWDAT,NOWTIM
FORMAT(/1X,A,2X,A)
VFD=.TRUE.
```

```
C Allow operator to change date.
```

```
100 WRITE(*,1003)
1003 FORMAT(' /'#$Enter Date in the form 01-FEB-85 ->')
IF (IOPT.NE.0) WRITE(*,1006)
```

```
1006 FORMAT('$ (or RETURN for no change)')
READ(5,1014,END=100,ERR=100) INPSTR
```

```
C! READ(5,1014,END=100,ERR=100) NCHAR,INPSTR
```

```
C! 1014 FORMAT(Q,A)
```

```
1014 FORMAT(A)
```

```
IF ((IOPT.NE.0).AND.(INPSTR.EQ.' ')) GOTO 200
```

```
IF (INPSTR(3:3).NE.'-') GOTO 100
```

```
C! IF (NCHAR.NE.0) GOTO 100
```

```
IF (IOPT.EQ.0) GOTO 110
```

```
GOTO 200
```

```
C
C We appear to have a date from the operator. Call JULDAT to
C parse it, and check its components for reasonableness.
```

```
110 VFD=.FALSE.
CALL JULDAT(INPSTR,NOWTIM,ADAT,ATIM)
IF (IYEAR.GE.84 .AND. IYEAR.LE.99) GOTO 120
TYPE 1024
```

```
1024 FORMAT (/20X,'**** INCORRECT YEAR NUMBER ****'/)
GOTO 100
```

```
120 IF (IMONTH.GE.1 .AND. IMONTH.LE.12) GOTO 130
TYPE 1025
```

```
1025 FORMAT (/20X,'**** INCORRECT MONTH NAME ****'/)
```

```

GOTO 100
130 IF (IDAY.GE.1 .AND. IDAY.LE.31) GOTO 170
TYPE 1020
1020 FORMAT (//20X,'**** INCORRECT DAY NUMBER ****/')
GOTO 100
170 DATE_TIME(1:8)=INPSTR(1:8)
J = INDEX(DATE_TIME,'-')+1
IF (J.GT.3) GOTO 175
TYPE 8000
8000 FORMAT(//20X,' **** INCORRECT FORMAT 04-APR-88 ****')
GOTO 100
175 J = INDEX(DATE_TIME(J+1:),'-')+1+J
DATE_TIME(J:J)='1'
DATE_TIME(J+1:J+1)='9'
DATE_TIME(J+2:J+2)=INPSTR(J:J)
DATE_TIME(J+3:J+3)=INPSTR(J+1:J+1)
C
C Allow operator to change time.
C
200 WRITE(*,1004)
1004 FORMAT(' '//#$Enter Time in the form 01:47:59 ->')
IF (IOPT.NE.0) WRITE(*,1006)
READ(5,1014,END=200,ERR=200) INPSTR
IF ((IOPT.NE.0).AND.(INPSTR.EQ.' ')) GOTO 500
IF (INPSTR(3:3).NE.':') GOTO 200
IF (IOPT.EQ.0) GOTO 210
GOTO 300
C
C We appear to have a time from the operator. Call JULDAT to
C parse it, and check its components for reasonableness.
C
10 VFD=.FALSE.
CALL JULDAT (NOWDAT,INPSTR,ADAT,ATIM)
IF (IHOUR.GE.0 .AND. IHOUR.LE.23) GO TO 220
TYPE 1034
1034 FORMAT (//20X,'**** INCORRECT HOUR NUMBER ****/')
GO TO 200
220 IF (IMIN.GE.0 .AND. IMIN.LE.59) GO TO 230
TYPE 1036
1036 FORMAT (//20X,'**** INCORRECT MINUTE NUMBER ****/')
GO TO 200
230 IF (ISEC.GE.0 .AND. ISEC.LE.59) GO TO 270
TYPE 1038
1038 FORMAT (//20X,'**** INCORRECT SECOND NUMBER ****/')
GO TO 200
270 CONTINUE
C
C Load DATE_TIME character string with time and pad with blanks
DATE_TIME(12:12)=' '
DATE_TIME(13:20)=INPSTR(1:8)
DATE_TIME(21:23)=' '
C
TYPE 9000,DATE_TIME
9000 FORMAT(1X,23A)
STATUS = SYS#BINTIM(DATE_TIME,SYSTIM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
STATUS = SYS#SETIME(SYSTIM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
C
C Repeat update and verification.
C

```

```
IOPT=2
GOTO 10
```

```
C
C Sync to the VAX time.
C
DO CALL JULDAT(VAXDAT,VAXTIM,ADAT,ATIM)
C Get VAX time from comm
C! CALL COMM('TS')
C Load DATE_TIME character string with time and date (modified date
C to meet Date String requirements.
DATE_TIME(1:9) = VAXDAT
DATE_TIME(10:11)= DATE_TIME(8:9)
DATE_TIME(8:9) = '19'
DATE_TIME(12:12)=' '
DATE_TIME(13:20)= VAXTIM
DATE_TIME(21:23)=' '
C TYPE 9000,DATE_TIME
C Update Time
STATUS = SYS#BINTIM(DATE_TIME,SYSTEM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
STATUS = SYS#SETIME(SYSTEM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
500 RETURN
C
C
END
```

```

C*****
C DIALGS      (DIALOG subroutine for solid isotopic system)
C
C      Assay and Background dialogue with operator for the solid
C isotopic instrument
C
C      ***This routine needs to be modified for the communications
C      package ***
C
C      SSJ 17-AUG-79
C      REV. 31-May-84 lc
C      REV. 05-Sept-84 TES 10:28
C      REV. 01-NOV-84 WDR CREATE FILENAME FOR CURRENT DATA
C      REV. 11-OCT-84 MJC: CHANGED ACCEPT TO READ WITH ERR TO BRANCH
C                          BACK TO 'ENTER NUMBER OF CYCLES/CALIBRATIONS'
C      REV. 29-NOV-84 MPK Replace list-directed with formatted I/O
C      REV. 19-DEC-84 WDR Modified for solid isotopic instrument
C      REV. 14-FEB-85 MPK      Rewrite to do "real" VAX comm
C                          Replace all ACCEPTS with READS having
C                          END and ERR exits
C                          Use YESNO to ask yes/no questions
C      REV. 20-FEB-85 MPK      Check ABORT after calling GETSMP
C      Rev. 23-OCT-85 WDR      Modified communications with the VAX
C                          installed 'LO' message to VAX
C      Rev. 21-NOV-85 WDR      Get Cadmium absorber thickness when
C                          calibrating
C      Rev. 02-Apr-86 WDR      Modified for Micro VAX II
C      Rev. 13-Jul-87 RDP      Modified for uVax II
C                          Convert BYTE Strings to CHARACTER strings
C      Rev. 13-Aug-87 RDP      Preliminary Hardware and Acquisition
C                          check added before performing ASSAY.
C      Rev. 02-May-90 WDR      Removed COMM call for Pu238 instrument

```

```

C*****

```

```

C      SUBROUTINE DIALOG

```

```

C      LOGICAL*1 YESNO,MTL
C      CHARACTER*15 CURDTA
C      REAL*4 MONTH
C      INTEGER SAMLEN

```

```

C      The following common block files are generic for all instruments

```

```

C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CFILES.CMN'
C      INCLUDE 'CONFIL.CMN'
C      INCLUDE 'CVAX.CMN'
C      INCLUDE 'CLOG.CMN'
C      INCLUDE 'CCAL.CMN'
C      INCLUDE 'CBACKG.CMN'
C      INCLUDE 'CMCONS.CMN'
C      INCLUDE 'CMC.CMN'
C      INCLUDE 'CDATE.CMN'
C      INCLUDE 'CMONTH.CMN'
C      INCLUDE 'P238COM.CMN'

```

```

C      The following common block files are instrument specific for

```

```

C gamma ray instruments
C
C      INCLUDE 'CSPECT.CMN'
C
C      DATA 'SAMPLN' / 79 /
C      DATA CURDTA / 'CURDTA.' : ' /
L
C Get operator and sample ID
C
C      CALL GETOPS(MTL)
C      IF(ABORT)RETURN
C
C Get filename if from disk
C
C      8      IF( .NOT.(AC .AND. WRITNG .OR. FROMDK))GO TO 85
C      9      TYPE 1059
C 1059     FORMAT(' '//Enter the full filename [ FILENM.EXT ] -> ')
C      READ (5,1001,END=9,ERR=9)FILNAM
C      IF(FILNAM(1:1) .EQ. ' ') GO TO 450
C
C Get number of cycles for autocycling
C
C      85     NUMCYC = 1
C           ICYCLE = 1
C           IF(.NOT. AUTOCY )GO TO 245
C      86     TYPE 1020
C 1020     FORMAT('/'Enter number of cycles -> ')
C 1013     FORMAT(I)
C           READ(5,1013,END=86,ERR=86)NUMCYC
C           IF(NUMCYC .LT. 0) GO TO 400
C           IF(NUMCYC .EQ. 0) GO TO 86
C           IF(FROMDK)GO TO 500
C
C Skip rest of code if no data is to be taken
C
C      245    IF( FROMDK)GO TO 500
C           IF(WRITNG .OR. AUD)GO TO 295
C
C Create a dummy filename for current data
C This is required for Livermore solution analysis code
C Get start date for data
C
C 247     FILNAM = CURDTA
C           CALL NEWNAM(0)
C           GOTO 350
C
C Set filename to sample id if writing or VAX is not up
C
C 295     IF(.NOT. WRITNG .AND. VAXUP)GO TO 350
C           IF (AC) GO TO 500
C           K = 1
C           DO 300 I = 1, SAMLEN
C             FILNAM(K:K) = SAMPID(I:I)
C             K = K + 1
C 300     CONTINUE
C           CALL NEWNAM(0)
L Final instructions for inserting the sample
C
C 350     IF(AC .OR. FROMDK) GO TO 500

```



C Perform Preliminary Check on Acquisition System

CALL PRECHK

IF (ABORT) GO TO 400

IF (BACKG) THEN

TYPE 1057

057       FORMAT(/'#Remove any sample, then press RETURN' ->')

ACCEPT 1001

1001       FORMAT(80A)

ENDIF

CALL ADATE

GO TO 500

C

C Abort the measurement because of negative input

C

400       TYPE 1019

1019       FORMAT(/' Measurement aborted by operator input')

C

450       ABORT = .TRUE.

C

500       RETURN

END



C\*\*\*\*\*

C ERRMSG

C

C Read appropriate error message from the  
C error message file and output to user terminal  
C Calling arguments:

C INUM = Message calling argument  
C ARG = ASCII argument passed to be printed (15A1)

C File format:

C INUM (Message call number) = I4  
C IARG (ASCII argument passed) = A1  
C MESSAG (ASCII error message) = 72A1

C

C 21-Apr-84 LC

C REV. 12-OCT-84 MJC: CLOSED DKLUN1 IF ARG = ' '  
C Rev. 19-Apr-85 MPK Remove redundant initialization  
C of DKLUN and DKLUN1  
C Rev. 29-Apr-85 MPK Put MESSAG and ICHAR in CMSG common  
C Rev. 02-Apr-86 WDR Modified for Micro VAX II  
C Rev. 13-Jul-87 RDP Modified for uVax II  
C Convert BYTE strings to CHARACTER strings

C\*\*\*\*\*

C

C

C SUBROUTINE ERRMSG(INUM,ARG)

C

C

C CHARACTER ARG\*15,IARG\*1

C

C The following common block files are generic for all instruments

C INCLUDE 'CFILES.CMN'  
C INCLUDE 'CLUNS.CMN'  
C INCLUDE 'CMSG.CMN'

C

C

C Open error message file  
C Locate proper error message

C

10 OPEN (UNIT=DKLUN1,NAME=ERRFIL,TYPE='OLD',ERR=25)  
10 READ(DKLUN1,1000,ERR=20,END=20)NUM,IARG,ICHAR,  
X (MESSAG(I),I=1,ICHAR)  
1000 FORMAT(I4,A1,Q,72A1)  
IF (NUM .EQ. -999)GO TO 20  
11 IF (NUM .EQ. INUM) GO TO 15  
GO TO 10

C

C Print error message

C

15 TYPE 1015, NUM, (MESSAG(I),I=1,ICHAR)  
1015 FORMAT(/5X,'\*\*\*\* ERROR # ',I3/5X,72A1,' \*\*\*\*')  
IF(IARG .NE. ' ')GO TO 17  
CLOSE(UNIT=DKLUN1)  
GO TO 30  
17 TYPE 1016, ARG  
1016 FORMAT(7X,A15)  
CLOSE(UNIT=DKLUN1)  
GO TO 30

C

C Error number not found in file

C

20 TYPE 1020, INUM, ERRFIL

1020 FORMAT(/5X, ' Error # ', I4, ' not located in ', A15)

CLOSE(UNIT=DKLUNI)

GO TO 30

C

C Error opening file

C

25 TYPE 1025, ERRFIL

1025 FORMAT(/5X, '\*\*\*\* UNABLE TO OPEN ', A15, ' \*\*\*\*')

C

C

30 RETURN

END



INTEGER\*4 ISTATUS,VKID,PBID

C System Routines

INTEGER\*4 SMG#CREATE\_PASTEBOARD  
INTEGER\*4 SMG#CREATE\_VIRTUAL\_KEYBOARD  
INTEGER\*4 SMG#ENABLE\_UNSOLICITED\_INPUT  
INTEGER\*4 SMG#DISABLE\_UNSOLICITED\_INPUT  
INTEGER\*4 SMG#DELETE\_VIRTUAL\_KEYBOARD  
INTEGER\*4 SMG#DELETE\_PASTEBOARD

EXTERNAL GET\_KEY

INCLUDE '(\$SYSSRVNAM)'

-----  
LOGICAL \* 1 NOPT,SUPV,AUTSND  
CHARACTER OPT(25)\*4, OPTIN\*4,OPTR\*3  
CHARACTER SAVFIL\*15, ACHAR\*1  
INTEGER OPTLEN,MENPRM,ON6,OFF6  
EQUIVALENCE (OPTR(1:1),OPTIN(2:2))

The following common block files are generic for all instruments

INCLUDE 'CFILES.CMN'  
INCLUDE 'CLUNS.CMN'  
INCLUDE 'CDIAG.CMN'  
INCLUDE 'CVAX.CMN'  
INCLUDE 'CLOG.CMN'  
INCLUDE 'CPASS.CMN'  
INCLUDE 'CMEAS.CMN'  
INCLUDE 'CMONTH.CMN'  
INCLUDE 'CDATE.CMN'  
INCLUDE 'CONTRL.CMN'  
INCLUDE 'CMC.CMN'  
INCLUDE 'CBACKG.CMN'  
INCLUDE 'CONFIL.CMN'  
INCLUDE 'CASSAY.CMN'  
INCLUDE 'CMCONS.CMN'  
INCLUDE 'CCAL.CMN'  
INCLUDE 'p238com.cmn'

The following common block files are instrument specific for  
gamma-ray instruments

INCLUDE 'CSPECT.CMN'  
INCLUDE 'CONT90.CMN'

DATA MAXASY,LOGREC /100,22/  
DATA JSW /"44/, ON6 /"100/, OFF6 /"177677/  
DATA MCA /1/, IADC /1/, MEMSIZ /4096/ !CHANGE FOR INSTRUMENT  
DATA MCMAX,MCREC /100,22/  
DATA DKLUN /1/, DKLUN1 /2/, LOUT /6/  
DATA MENPRM /25/, OPTLEN /4/  
DATA ASLOG //ASYLOG.WHS//  
DATA CONFIL //PARMTR.WHS//  
DATA BKFIL //BACKGR.WHS//  
DATA MCLOG //MCLOGF.WHS//  
DATA NAMSFL //CONSTA.WHS//  
DATA ERRFIL //ERRMSG.INS//  
-----

```

C ** Single Key Interrupt for VAX input
C Create PasteBoard
      ISTATUS = SMG$CREATE_PASTEBOARD(PBID,'SYS#INPUT')
      IF (.NOT.ISTATUS) CALL LIB$SIGNAL(%VAL(ISTATUS))
C
C Create Virtual Keyboard
      ISTATUS = SMG$CREATE_VIRTUAL_KEYBOARD(VKID,'SYS#INPUT')
      IF (.NOT.ISTATUS) CALL LIB$SIGNAL(%VAL(ISTATUS))
C
C Enable Virtual Input
      ISTATUS = SMG$ENABLE_UNSOLICITED_INPUT(PBID,GET_KEY,VKID)
      IF (.NOT.ISTATUS) CALL LIB$SIGNAL(%VAL(ISTATUS))
C
C-----
C Announce startup
C
      TYPE *,' Starting WHC_ASY238 program Version 1.02'
C
C Initialize logical flags and hardware for the instrument
C If problem with the hardware - abort flag is set
C
      CALL START
      IF (ABORT) GOTO 500
C
C Reset default settings for the Assay options
C Argument 0 resets default flags and parameters
C Argument 1 permits operator to select flags and parameters
C
      CALL RESET(0)
C
C set values for 8232 stabilizer
C
      CALL STABSET
C
C Get time sync from VAX and set date and time
C
C! Get sync time from VAX and set time and date.
C! If VAX is down test if time and date are reasonable.
C! if not, get time and date from user.
      IOPT = 3
      CALL COMM('TS')
      IF (.NOT.VAXUP) GO TO 15
      CALL DATTIM(IOPT)
C
C Initialize menu option list
C
      15 DO 20 I = 1, MENPRM
          OPT(I) = ' '
      20 CONTINUE
C
C Close line printer if in use
C
      30 IF (LOUT.NE.6) GOTO 100
      CLOSE(UNIT=LOUT,ERR=100)
C
C Open a non-existent file on the system disk to clear
C the directory, in case floppies are switched by the
C operator; applies to noswap
C
      100 OPEN(UNIT = DKLUN, NAME = 'NOFILE.DAT', TYPE = 'OLD',

```

```

1      ERR = 105)
      CLOSE(UNIT=DKLUN)
105   CONTINUE

C
C Check status of hardware - determines if hardware is okay
  to run
L
C
C 110   ABORT = .FALSE.
      CALL STATUS(0)           !instrument specific routine
      IF (ABORT) GOTO 500

C
C Re-initialize logical flags
C
      CALL INTFLG

C
C Get system date/time
C
120   CALL DATTIM(1)

C
C Prompt operator for menu option
C Menu options are:
C 202  A   - Assay
C 206  MB  - Measurement control - Bias
C 208  MP  - Measurement control - Precision
C 210  H   - Help - display menu
C 212  AC  - Assay Current data in MCA
C 214  AD  - Assay data from Disk
C 216  AU  - AUTOCYCLE
C 218  AUD - AUTOCYCLE from Disk
C 222  C   - Calibration
C 224  CD  - Calibration from Disk
C 226  D   - Default
C 228  LA  - List Assay log
C 230  LM  - List Measurement control log
C 232  ST  - Status
C 234  R   - Read data from disk
C 236  W   - Write data to disk
C 238  OU  - change OUTPUT listing device
C 240  HS  - Help Supervisor mode
C 500  X   - eXit

C
      TYPE 1003, NOWDAT, NOWTIM
1003  FORMAT(' /,1X,A9,2X,A8,3X,'Enter OPTION ( H to see Menu) -> ',#)

C
C Check LN monitors while waiting for operator input.
C Read operator input - if error in input, output a
C message and prompt again
C
C
125   CALL LNMN38
      IF (ABORT) GOTO 120

C
C-----
C
C Test if Key has been press...
      IF (INKEY.NE.0) THEN
          ACHAR = INCHAR           ! Get Character
          INKEY = 0               ! Reset INKEY Flag
          TYPE 2001,ACHAR        ! Echo Character

```



```

2001     FORMAT('+',A,#)
        ELSE
          CALL SLEEP(0.5)
          GOTO 125           ! Return to loop
        ENDIF

```

```

-----
1004     READ (5,1004,END=120,ERR=120) OPTR
        FORMAT(A3)
        OPTIN(1:1) = ACHAR
        IF (OPTIN(1:1).EQ.' ') GOTO 120

```

C  
C Menu searches a table of options available and matches the  
C input to a menu option in the table - if there is not a  
C match, a flag is set and the operator is prompted again.  
C When a match is found menu returns a number of the routine  
C to be called.

```

        CALL MENU(OPTIN,MNUM,NOPT)
        IF (ABORT) GOTO 130
        IF (NOPT) GOTO 120
        GOTO 140
130     CALL ERRMSG(102)
        ABORT = .FALSE.
        GOTO 110
140     IF (MNUM .GT. 5) GO TO 211
        GOTO (202,204,206,208,210) MNUM

```

C  
C The following is a list of the subroutines that can  
C be executed from the menu option list - these menu options  
C do not require a password to access them

```

202     CALL ASOPT
        GO TO 260
204     CONTINUE           !Background option removed
206     CALL MCBIAS
        GO TO 260
208     CALL MCPREC
        GO TO 260
210     SUPV=.FALSE.
        CALL HELP(SUPV)
        GO TO 260

```

C  
C These menu options require a password

```

211     CALL PASS(OK)
        IF(.NOT. OK) GO TO 110
        MPNUM = MNUM - 5
        GO TO (212,214,216,218,220,222,224,226,228,230,
1         232,234,236,238,240,500) MPNUM
212     CALL CURDTA
        GO TO 260
214     AD = .TRUE.
        GO TO 2201

```

```

C
218  AUD = .TRUE.
      FROMDK = .TRUE.
C
216  CALL AUTCYC
      GO TO 260
L
220  CONTINUE          !Background from disk removed
      GO TO 260
2201 CALL FRDSK
      GO TO 260
C
224  CALDSK = .TRUE.
      FROMDK = .TRUE.
C
222  CONTINUE          !CALIBRATION OPTION REMOVED
C 222  CALL CALIBS
      GO TO 260
C
226  CALL RESET(1)
      GO TO 260
C
228  AUTSND=.FALSE.
      CALL LSLOG(AUTSND)
      GO TO 260
C
230  AUTSND=.FALSE.
      CALL LMCLOG(AUTSND)
      GO TO 260
C
232  CALL STATUS(1)          !instrument specific routine
      GO TO 260
C
C Prompt for filename for reading and writing data to/from disk
C
234  SAVFIL=FILNAM
231  TYPE 2005
2005  FORMAT('/$Enter full filename [FILENM.EXT] -> ')
      READ (5,2006,END=231,ERR=231) FILNAM
2006  FORMAT(A15)
      IF(FILNAM(1:1) .EQ. ' ') GO TO 260
      IF (OPTIN(1:1).EQ.'W') GOTO 236
      CALL PUT90(IADC,0)          !instrument specific routine
      GO TO 237
236  CALL GET90(IADC)          !instrument specific routine
237  FILNAM=SAVFIL
      IADC = IADCSV
      GO TO 260
C
238  TYPE 2007
2007  FORMAT(' /$Enter 6 for line printer, 5 for terminal -> ')
      READ (5,2008,END=238,ERR=238) LOUT
2008  FORMAT(I)
      IF(LOUT .EQ. 0) LOUT = 5
      IF(LOUT .NE. 6 .AND. LOUT .NE. 5) GO TO 238
      GO TO 110
L
240  SUPV = .TRUE.
      CALL HELP(SUPV)

```

```

      GO TO 110
C
260  IF(ABORT) GO TO 30
C
C If the Vax is up, send all files that have not yet been sent
      IF(.NOT. VAXUP) GO TO 30
      TYPE 2010
2010  FORMAT(/' Files being sent to VAX ')
      AUTSND =.TRUE.
      CALL LSLOG(AUTSND)
C
C
      GO TO 110
500  CONTINUE
C
C Exit from program
C
C      CALL COMM('DS')
C-----
C *** Remove Single Key Input
C Disable Unsolicited Virtual Keyboard Input
      ISTATUS = SMG#DISABLE_UN SOLICITED_INPUT(PBID)
      IF (.NOT.ISTATUS) CALL LIB#SIGNAL(%VAL(ISTATUS))
C
C Delete Virtual Keyboard
      ISTATUS = SMG#DELETE_VIRTUAL_KEYBOARD(VKID)
      IF (.NOT.ISTATUS) CALL LIB#SIGNAL(%VAL(ISTATUS))
C
C Delete Pasteboard
      ISTATUS = SMG#DELETE_PASTEBOARD(PBID)
      IF (.NOT.ISTATUS) CALL LIB#SIGNAL(%VAL(ISTATUS))
C-----
      RETURN
      END

```

SUBROUTINE FIT

C  
 CCC  
 C FIT modified to use MATINV and linear peak heights.  
 C Last edit (R.G.) Statement 464 correction  
 C Last edit for VAX 13-OCT-87  
 L Modified for Pu238 instrument WDR 1-Sep-89  
 C  
 CCC  
 C

```

INCLUDE 'GRPLGCM.PMS'
INCLUDE 'GRPRFL.PMS'

```

C The following common is required for the instrument

```

INCLUDE 'CONTRL.CMN'

```

```

BYTE NFLG1,NFLG2,NFLG3,NFLG4,NFLG5
EQUIVALENCE (NFLG1,NFLG(1)),(NFLG2,NFLG(2))
EQUIVALENCE (NFLG3,NFLG(3)),(NFLG4,NFLG(4)),(NFLG5,NFLG(5))
EQUIVALENCE (KK,NWRDS(2)) !A fudge to transfer KK to OUT.

```

```

DIMENSION AX(MXFP)
EQUIVALENCE (AX,VM)
DOUBLE PRECISION APHA(MXFP,MXFP),XSUM(MXFP)

```

```

IF(MAXITR .EQ. 0)MAXITR = 10 !Max iterations may be specified
LFLG2 = 0
LFLG4 = 0
FIXFLG2 = 0
FIXFLG4 = 0
ITFLG=0 !INIT ITERATION COMPLETION FLAG
ASLP = SHAPC(2)/GAIN

```

C C C C C C C C C C  
 C Five NFLG's flag the freedom of ALPHA, short-term tail size &  
 C slope, and long-term tail size & slope, respectively.  
 C Reiteration comes back to here.

```

1 KK = NFLG1 + NFLG2 + NFLG3 + NFLG4 + NFLG5

SG = -2.7726 / ALPHA
CX = SQRT(-ALPHA)
XPTS = NDPTS

```

C CALCULATE PEAK TAILING INTENSITIES AT THE GROUP BOUNDARIES.  
 C THESE VALUES USED LATER TO ADJUST NET COUNTS OF PEAK  
 C GROUPING.

```

BG(1) = 0.
BG(2) = 0.
TLS = EXP4/GAIN
IF ((EN(1) .GT. 4500.) .AND. (TLS .LT. .1) .AND. (TLS .GT. 0.))
+ GOTO 56 !Assume Alpha Spectrum where long term tail = BG

```

```

DO 55 J = 1,NP
AMP1 = EXP1
AMP2 = EXP3
XI = PKPOS(J) + XTRAPL !Correct for remote bkgds.

```

```

DO 50 L = 1,2
CC = 0.0
IF(GAMA(J) .EQ. 0.)GO TO 10
WL = GAMA(J) * CX
CC = CX**XI
EXPNT = ALPHA * XI * XI
EXPNT1 = 0.
IF(EXPNT .GT. -14.)EXPNT1 = EXP(EXPNT)
CALL BWF2(WL,CC,EXPNT,EXPNT1)
10  BG(L)=BG(L)+PKHT(J)*(CC+AMP1*EXP(-EXP2*XI)
1   * (1.- EXP(.4*ALPHA*XI*XI)) + AMP2*EXP(-EXP4*XI))
AMP1 = 0.
AMP2 = 0.
XI = XPTS - (PKPOS(J) + XTRAPH)
50  CONTINUE
55  CONTINUE

56  YP = 0.
DBG = BG(2) - BG(1)
CALL AVE(1,NDPTS,YN SUM,AV, YNET)

C   ZERO COEFFICIENT ARRAY

DO 60 J=1,MXFP
XSUM (J) = 0.
DO 60 I=1,MXFP
60  APHA(I,J)=0.
DO 61 J=1,NP
61  SUM(J) = 0.0
CHISQ = 0.

FORM MATRIX OF LINEAR EQUATIONS

DO 200 I = 1,NDPTS
DO 65 JK = 1,MXFP
65  AX(JK) = 0.
XI = I
FX = 0.0
DALFA = 0.0
DBETA = 0.0
DGAMA = 0.0
DDELTA = 0.0
DCHI = 0.0
YP = YP + YNET(I)
DELY(I) = YNET(I) + BG(1) + YP * DBG / YNSUM
RM(I) = DELY(I)
TYNET = YNET(I)
IF (TYNET .LE. 0) TYNET = TYNET + 1.
EY = ABS(TYNET) + AVEBG + 1.0
IF (EY .LT. 0.) EY = 1.0
EY = SQRT(EY) + 0.001*EY           !Channel width uncertain

C   Prepare for profile plotting only if requested.
IF(MPLGS(4) .EQ. 0)GO TO 75
KPHF = 0           !Clear ctr. for free peak profiles.
SVYNET(I) = DELY(I)
ERYNET(I) = EY
DO 70 JK = 1,MXNP           !Clear profile arrays.
70  KKB(JK) = 0
PKPRFL(JK) = 0.

```

```
LBF = 0          !Buffer #
XEN(I) = EST + (I-1)*GAIN
```

```
C DETERMINE WEIGHTING FACTOR
```

```
75 WT = 1.0 / EY
   KPK = KK
```

```
C CALCULATE COEFFICIENTS FOR EACH OF "NP" PEAKS
```

```
DO 170 J = 1, NP
DLTAX = XI - PKPOS(J)
DX = DLTAX
DLX = DLTAX
DLXSQ = DLTAX * DLTAX
FX2 = 0.0
FX3 = 0.0
CB1 = 0.
CB2 = 0.
EXPNT = 0.0
PHT = PKHT(J)
ALFA = -2.7726 / (ASLP * PKPOS(J) + SG)
EXPNT = ALFA * DLXSQ
EXPNT3 = 0.
IF (EXPNT .GT. -14.) EXPNT1 = EXP(EXPNT)
```

```
C CALCULATE GAUSSIAN COMPONENT OF PEAK
```

```
IF (DLTAX .LT. 0) EXPNT3 = EXP (0.4 * EXPNT)
FX1 = EXPNT1
IF (GAMA(J) .LE. 0) GO TO 80
```

```
C COMPUTE MODIFIED SHAPE FOR X-RAY PEAK
```

```
CX = SQRT(-ALFA)
WL = GAMA(J) * CX
CC = CX * DLTAX
CALL BWF2(WL, CC, EXPNT, EXPNT1)
FX1 = CC
```

```
80 IF (NFLG1 .GT. 0) THEN !Is peak width parameter free?
   DALFA = DALFA + PHT * DLXSQ * EXPNT1 !Yes
```

```
ENDIF
```

```
IF (DLTAX .GE. 0) GO TO 130
```

```
EXPNT = EXP2 * DLTAX
```

```
IF ((EXPNT + 10.) .LT. 0) GO TO 110
```

```
EXPNT2 = EXP (EXPNT)
```

```
CB = PHT * EXPNT2 * (1.0 - EXPNT3)
```

```
CB1 = EXP1 * EXP2 * CB
```

```
FX2 = EXP1 * CB
```

```
IF (NFLG2 .GT. 0) THEN !Is short term tailing amplitude free?
```

```
   DBETA = DBETA + CB !Yes
```

```
ENDIF
```

```
IF (NFLG3 .LE. 0) GO TO 110
```

```
C SHORT TERM TAILING SLOPE FREE
```

```
DGAMA = DGAMA + DLTAX * CB
```

```
110 IF (EXP3 .LE. 0) GO TO 120
```

```
EXPNT = EXP4 * DLTAX
```

```
IF ((EXPNT + 10.) .LT. 0) GO TO 120
```

```
CB = PHT * EXP(EXPNT) * (1. - EXPNT3)
```

```
CB2 = EXP3 * EXP4 * CB
```

```
FX3 = EXP3 * CB
```

```

IF(NFLG4 .GT. 0) THEN !Is long term tailing amplitude free?
  DDELT = DDELT + CB !Yes
ENDIF

```

```

IF(NFLG5 .LE. 0)GO TO 120

```

```

L LONG TERM TAILING SLOPE FREE

```

```

DCHI = DCHI + DLTAX * CB
120 IF (NFLG3+NFLG5 .EQ. 0) THEN
  IF (NFLG2 .EQ. 0) FX = FX + FX2
  IF (NFLG4 .EQ. 0) FX = FX + FX3
ELSE
  FX = FX + FX2 + FX3
ENDIF

```

```

130 IF(NYFLG(J) .GT. 0)THEN !Is peak size freed?
  KPK = KPK + 1 !Yes, count free param.
  IF(MFLGS(4) .NE. 0)KPHF = KPHF + 1 !Count the buffer
  L=KPK
  VALU=FX1*WT
  GO TO 140
ENDIF

```

```

C Peak size is not free. [NYFLG is not positive.]

```

```

IF (HIGHT(J) .GE. 0.) THEN
  FX = FX + FX1 *PHT !Peak size is fixed absolutely.
  IF(MFLGS(4) .EQ. 0) GO TO 160
  MFLGS(4) = 2 !Flag: a profile stored at MXNP
  LBF = MXNP !Assign 1 buf. to all fixed peaks
  PKPRFL(MXNP) = PKPRFL(MXNP) + FX1*PHT + FX2 + FX3 !Store
  IF(ITFLG .EQ. 0) GO TO 160
  L = MXNP
  GO TO 150
ENDIF

```

```

C Peak size is fixed relative to another peak.

```

```

LPK = -NYFLG(J) !Find index of primary peak.
IF (HIGHT(LPK) .GT. 0.) THEN
  FX = FX-HIGHT(J)*PKHT(LPK)*FX1 !Primary is fixed also.
  GOTO 160 !Skips plotting!
ENDIF

```

```

VALU=-HIGHT(J)*FX1 *WT
L = KK

```

```

IF(LPK-1 .GT. 0) THEN
  DO 135 M=1,LPK-1 !Skip stg. space already used up.
  IF(NYFLG(M) .GT. 0)L= L+1
  IF(NXFLG(M) .GT. 0)L=L+1
135 CONTINUE

```

```

ENDIF
L= L+1 !Locn. indexing profile of primary .

```

```

140 AX(L) =VALU + AX(L)
IF(MFLGS(4) .EQ. 0)GO TO 160 !No plots

```

```

Store peak profile of each free peak.
PKPRFL(L) = PKPRFL(L) + FX1 * PHT + FX2 + FX3
LBF = KPHF

```

IF(ITFLG .EQ. 0)GO TO 160

Iterations are finished, clean up plotting.

IF(NYFLG(J) .LT. 0) THEN !Is peak size freed?

IF(KKB(L) .EQ. 0)GO TO 160 !No. Look for index of relative

LBF = KKB(L) !Substitute this profile.

GO TO 150

ENDIF

DO 148, K = J,NP

IF(-NYFLG(K) .EQ. J)GO TO 153

CONTINUE

PRFILS(I,LBF) = PKPRFL(L)

GO TO 160

KKB(L) = LBF

C C C C C

PEAK POSITION FREE?

IF (NXFLG(J) .GT. 0)THEN

KPK = KPK + 1 !Yes. Count the free parameter

AX(KPK) = (-2.\*ALFA \* DX \*EXPN1 \*PHT -CB1-CB2) \* WT

ENDIF

SUM(J) = SUM(J) + FX1 \* PHT + FX2

RM(I) = RM(I) - FX1 \* PHT - FX2 - FX3

CONTINUE

L = 0

IF (NFLG1 .GT. 0) THEN

L = L + 1

AX(L) = DALFA \* WT

ENDIF

IF (NFLG2 .GT. 0) THEN

L = L + 1

AX(L) = DBETA \* WT

ENDIF

IF (NFLG3 .GT. 0) THEN

L = L + 1

AX(L) = DGAMA \* EXP1 \* WT

ENDIF

IF(NFLG4 .GT. 0)THEN

L = L + 1

AX(L) = DDELTA \* WT

ENDIF

IF(NFLG5 .GT. 0) THEN

L = L + 1

AX(L) = DCHI \* EXP3 \* WT

ENDIF

DELY(I) = (DELY(I)-FX) \* WT

DO 190 KM = 1,KPK

XSUM(KM) = XSUM(KM) + AX(KM)\*DELY(I)

DO 190 KL = 1,KPK

APHA(KM,KL) = APHA(KM,KL) + AX(KM)\*AX(KL)

RM(I) = RM(I) \* WT

CHISQ = CHISQ + RM(I)\*\*2

CONTINUE



C Matrix is now loaded.

```
IF(ITER .EQ. 1)GO TO 210
RCHISQ = CHISQ/(NDPTS-KPK)
IF (LONGPR)WRITE (LOUT,1208)RCHISQ
)8 FORMAT(' REDUCED CHI-SQUARE VALUE IS ',E15.6)
```

210 IF(ITFLG .EQ. 1) GO TO 500  
CALL MATINV(APHA,MXFP,KPK,DET)

```
IF(DET .EQ. 0.) THEN
TYPE *,'MATRIX IS SINGULAR'
ABORT = .TRUE.
RETURN
ENDIF
```

```
POS = NDPTS/2 !For FWHM at mid-energy
IF (NP .EQ. 1) POS = PKPOS(1) !At peak position
FWHM = GAIN * SQRT(ASLP*POS + SG - .462)
EX2 = EXP2/GAIN
EX4 = EXP4/GAIN
```

9040 IF (LONGPR)WRITE (LOUT,9040)  
FORMAT(/,' ITERATION FWHM EXP1 EXP2 EXP3 EXP4')

9030 IF (LONGPR)WRITE (LOUT,9030)ITER,FWHM,EXP1,EX2,EXP3,EX4

9030 FORMAT(I3,6X,F9.5,4F8.4)

9010 IF (LONGPR)WRITE (LOUT,9010)  
FORMAT(/,' CHANNEL KEV PEAK HEIGHT')

```
DO 280 J=1,NP
```

```
PKP=PKPOS(J)+OFFSET
```

```
PKE =(PKP - REFCH)*GAIN + REFEN
```

```
IF (LONGPR)WRITE (LOUT,9020)PKP,PKE,PKHT(J)
```

```
CONTINUE
```

9020 FORMAT(2F9.2,F14.3)

```
DO 310 L = 1,KPK
```

```
DEL(L) = 0.
```

D IF(APHA(L,L) .LT. 0)TYPE 1310,L

D1310 FORMAT(/' \*\*\* Warning - Negative matrix element at L=',I3/)

```
DM(L) = SQRT(ABS(APHA(L,L)))
```

```
DO 310 J=1,KPK
```

310 DEL(L) = DEL(L) + XSUM(J)\*APHA(L,J)

C

C INCREMENT THE TRIAL VALUES AND CHECK FOR COMPLETION OF ITERATIVE  
C PROCESS.

```
DO 320 J = 1,5
```

```
IF(NFLG(J) .LE. 0) ERR(J+2) = 0.
```

320 CONTINUE

```
L = 0
```

```
LFLG = 0
```

```
IF (NFLG(1) .LE. 0)GO TO 410
```

C INCREMENT PEAK WIDTH PARAMETER

```
L = L + 1
```

```
ERR(3) = DM(L)
```

```
ADEL = ABS(DEL(L))
```

```

ALPHA = ALPHA + DEL(L) * (1. + ADEL/(ALPHA - ADEL))
IF ((ADEL + .01 *ALPHA) .GT. 0) LFLG = 1
410 IF (NFLG(2) .LE. 0)GO TO 425
L = L + 1
ERR(4) = DM(L)
L
INCREMENT SHORT TAILING AMPLITUDE

IF (NFLG3+NFLG5 .EQ. 0) GOTO 419
IF ((EXP1 + DEL(L)) .LE. 0) DEL(L) = -.5 * EXP1
ADEL = ABS(DEL(L))
EXP1 = EXP1 + DEL(L) * (1. - ADEL/(EXP1 + ADEL))
IF ((ADEL - DM(L) - .01 * EXP1) .GT. 0) LFLG = 1
GOTO 425

419 IF(DEL(L) .GT. 0.)GO TO 420
IF (LONGPR)TYPE 1419
1419 FORMAT(' SHORT TERM TAILING AMPLITUDE NEGATIVE--SET TO 0.001')
DEL(L) = 0.001
LFLG2 = LFLG2 + 1
IF(LFLG2 .GT. 4) THEN
    NFLG2 = 0
    LFLG = 1
    FIXFLG2 = 1
ENDIF
420 IF (ABS(EXP1 - DEL(L)) - .5 * DM(L) .GT. 0) LFLG = 1
EXP1 = DEL(L)

425 IF (NFLG(3) .LE. 0)GO TO 435
C
INCREMENT SHORT TERM TAILING SLOPE

L = L + 1
ERR(5) = DM(L)
IF ((EXP2 + DEL(L)) .LE. 0)DEL(L) =- .5 * EXP2
ADEL = ABS(DEL(L))
EXP2 = EXP2 + DEL(L) * (1. - ADEL / (EXP2 + ADEL))
IF ((ADEL - DM(L) - .02* EXP2) .GT. 0) LFLG = 1
435 IF(NFLG(4) .LE. 0)GO TO 440
C
INCREMENT LONG TERM TAILING AMPLITUDE

L = L + 1
ERR(6) = DM(L)
IF (NFLG3+NFLG5 .EQ. 0) GDTO 437
IF ((EXP3 + DEL(L)) .LE. 0) DEL(L) = -.5 * EXP3
ADEL = ABS(DEL(L))
EXP3 = EXP3 + DEL(L) * (1. - ADEL/(EXP3 + ADEL))
IF ((ADEL - DM(L) - .01 * EXP3) .GT. 0) LFLG = 1
GOTO 440

437 IF(DEL(L) .GT. 0.)GO TO 438
IF (LONGPR)TYPE 1437
1437 FORMAT(' LONG TERM TAILING AMPLITUDE NEGATIVE--SET TO .0001')
DEL(L) = 0.0001
LFLG4 = LFLG4 + 1
IF (LFLG4 .GT. 4) THEN
    NFLG4 = 0
    LFLG = 1
    FIXFLG4 = 0

```

```

ENDIF
438 IF (ABS(EXP3 - DEL(L)) - 0.5 * DM(L) .GT. 0) LFLG = 1
    EXP3 = DEL(L)
440 IF(NFLG(5) .LE. 0)GO TO 445
    INCREMENT LONG TERM TAILING SLOPE

    L = L + 1
    ERR(7) = DM(L)
    IF((EXP4 + DEL(L)) .LE. 0)DEL(L) = -0.5 * EXP4
    ADEL = ABS(DEL(L))
    EXP4 = EXP4 + DEL(L) * (1. -ADEL / (EXP4 + ADEL))
    IF((ADEL - DM(L) - 0.02 * EXP4) .GT. 0) LFLG = 1
445 K = L
    DO 460 J = 1,NP
    IF (NYFLG(J) .LE. 0)GO TO 450 !If fixed height

C GET NEW PEAK HEIGHTS

    K = K + 1
    PKHT(J) = DEL(K)
    IF (DEL(K) .LT. 0.0) PKHT(J) = -1.0 !Damper
450 IF (NXFLG(J) .LE. 0)GO TO 460 !If fixed posn.

C INCREMENT PEAK POSITIONS
    K = K + 1
    ADEL = ABS(DEL(K))
    PKPOS(J) = PKPOS(J) + DEL(K) * (1. - ADEL/(1.+ADEL))
    IF((ADEL - DM(K) - .01) .GT. 0)LFLG = 1
460 CONTINUE

REESTABLISH PEAK ENERGY AND INTENSITY RELATIONSHIPS
C SPECIFIED IN INPUT.

    DO 470, J = 1,NP
    IF(NYFLG(J) .GE. 0)GO TO 465
    IF(HIGHT(J) .GT. 0)GO TO 464
    IF(HIGHT(J) .EQ. 0)THEN
        PKHT(J) = 0 !Eliminated peak.
        GO TO 465
    ENDIF
    L = -NYFLG(J)
    PKHT(J) = -HIGHT(J) * PKHT(L)
    GO TO 465
464 PKHT(J) = PKHT(J) *HIGHT(J) / SUM(J)

465 IF(NXFLG(J) .GE. 0)GO TO 470
    L = -NXFLG(J)
    PKPOS(J) = PKPOS(L) -(EN(L) - EN(J))/GAIN
470 CONTINUE
    IF(LFLG .EQ. 0)ITFLG = 1
    IF (ITER .GE. MAXITR) ITFLG = 1
    IF(ITFLG .EQ. 0)ITER = ITER + 1

C ITERATE CALCULATION AGAIN

    GO TO 1

500 IF(ITER .GE. MAXITR.AND. LONGPR)TYPE 1500
1500 FORMAT(' MAXIMUM ITERATIONS REACHED.')
```

```

C      Check for missing peaks. If found, set flags & retry fitting.
      L = 0
      DO 530, J=1,NP
      IF(PKHT(J) .GE. 0) GO TO 530
      L = L + 1
      HIGHT(J) = 0.
      PKHT(J) = 0.
      NYFLG(J) = -J
      NXFLG(J) = 0
      SUM(J) = 0.
530    CONTINUE
      IF(L .EQ. 0)      RETURN

      ITFLG = 0
      ITER = 1
      GO TO 1          !Try again.

      END

```

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

```

```

C
C      VERSION 29-OCT-82

      SUBROUTINE BWF2(WL,X,EXPNT,EXPN1)

      DIMENSION CA(5),CB(6)
      DATA CA/.729281,.0344557,.127639,-.0325239,.00642257/
      DATA CB/.933301,-4.45603,1.16191,-.141124,-.0352464,.00744387/

      C1 = 1. - .5642*WL
      C2 = .15915 * WL
      C3 = .25 * WL**2
      C4 = .6366 * WL
      RT = C1 + C3
      X = ABS(X)

      IF (X .GE. 2.2) GOTO 10
      XX = X * X
      AX = 1.
      SUM = 0.0
      DO 4 I = 1,5
      AX = AX * XX
-4    SUM = SUM + CA(I) * AX
      X = SUM * EXP (-XX)
      GO TO 20
-10   XLOG = ALOG (X)
      SUM = CB(1)
      DO 12 I = 2,6
      SUM = SUM + CB(I) * XLOG ** (I-1)
12    X = EXP (SUM)
20    AA = (C1 - C2*EXPNT + C3*(1. + 2.*EXPNT))*EXPN1 + C4*X
      X = AA / RT

      RETURN
      END

```

C\*\*\*\*\*

C  
C FRDSK

C Responds to operator selection of assay from disk and  
C background from disk

C  
C CMS 10-MAY-84  
C REV.11-SEP-84 CMS  
C Rev 29-Nov-84 MPK Replace list-directed with  
C formatted I/O  
C Rev. 12-FEB-85 MPK Redo logical flags  
C Rev. 04-Apr-86 WDR Modified for Micro VAX II  
C Rev. 13-Jul-87 RDP Modified for uVAX II  
C Convert BYTE strings to CHARACTER strings

C\*\*\*\*\*

C  
C SUBROUTINE FRDSK  
C CHARACTER SAVFIL\*15

C The following common block files are generic for all instruments

C  
C INCLUDE 'CLUNS.CMN'  
C INCLUDE 'CONTRL.CMN'  
C INCLUDE 'CBACKG.CMN'  
C INCLUDE 'CFILES.CMN'

C  
C FROMDK = .TRUE.  
C RTYP = 'Fromdisk'  
C IF(BD) TYPE 1000  
C 1000 FORMAT (' Background from disk')  
C IF(AD) TYPE 1001  
C 1001 FORMAT (' Assay from disk')

C Get parameter set ups from dialog

C  
C CALL DIALOG  
C IF(ABORT)GO TO 500

C Save original filename

C  
C SAVFIL = FILNAM

C Perform assay

C  
C CALL ASSAY  
C IF (ABORT) GOTO 500

C Restore filename

C  
C FILNAM = SAVFIL

C Print results

C  
C CALL PRINT1 !instrument specific routine  
C 500 RETURN

END

1944-1945  
1946-1947  
1948-1949  
1950-1951  
1952-1953  
1954-1955  
1956-1957  
1958-1959  
1960-1961  
1962-1963  
1964-1965  
1966-1967  
1968-1969  
1970-1971  
1972-1973  
1974-1975  
1976-1977  
1978-1979  
1980-1981  
1982-1983  
1984-1985  
1986-1987  
1988-1989  
1990-1991  
1992-1993  
1994-1995  
1996-1997  
1998-1999  
2000-2001  
2002-2003  
2004-2005  
2006-2007  
2008-2009  
2010-2011  
2012-2013  
2014-2015  
2016-2017  
2018-2019  
2020-2021  
2022-2023  
2024-2025

-----  
SUBROUTINE GETIME(NDEV,LIVTIM,RLTIM,DEDTIM)

CC

L A sub for picking live,real, & start times out of first  
C two channels of spectrum.

CC

REAL\*4 LIVTIM,RLTIM

DIMENSION BUF(2)

CALL RDCHNS (BUF,1,2,NDEV)

LIVTIM = BUF(1)

RLTIM = BUF(2)

DEDTIM = 100. \*(RLTIM - LIVTIM)

RETURN

END

SUBROUTINE GET\_KEY(PBID,VKID)

C\*\*\*\*\*

C Routine gets input from the virtual keyboard buffer.

C It sets a char detected flag (INKEY = 1) and the character

C detected (INCHAR). The character is echoed to the display.

INCLUDE '(\$TRMDEF)'

COMMON/RTI/ INKEY,INCHAR  
 INTEGER INKEY  
 CHARACTER INCHAR\*1,INPUT\*4  
 INTEGER\*4 STATUS,VKID,PBID,LENGTH,TIMEOUT

INTEGER\*4 SMG\$READ\_STRING

C Read Key Input

LENGTH = 1  
 TIMEOUT = 3

STATUS = SMG\$READ\_STRING(VKID,  
 INPUT,  
 LENGTH,  
 TRM\$M\_TM\_NOECHO+TRM\$M\_TM\_CVTLOW,  
 TIMEOUT)

IF (.NOT.STATUS) THEN  
 INCHAR= '  
 INKEY = 0  
 CALL LIB\$SIGNAL(%VAL(STATUS))

ELSE  
 INCHAR= INPUT(1:1)  
 INKEY = 1

ENDIF

END



```

C*****
C
C  GETOPS
C      Get operator and sample ID
C      Output argument
C
C      MTL (LOGICAL*1)      .TRUE. if material type and
C                          isotopics are needed
C
C      31-May-84 lc
C      Rev 14-Feb-85 MPK      Put in "real" VAX communications
C                          Supply end and err exits for
C                          terminal input
C      Rev 08-Oct-85 MPK      Adapted to neutron coinc counters
C                          GETOP and GETSMP combined
C      Rev 21-Oct-85 MPK      Pad "canned" sample ID properly
C      Rev 23-OCT-85 WDR      Revised for solution assay system
C      Rev 04-Apr-86 WDR      Modified for Micro VAX II
C*****
C
C      SUBROUTINE GETOPS(MTL)
C
C      The following common block files are generic for all instruments
C
C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CONTRL.CMN'
C
C      LOGICAL*1 YESNO,MTL
C      CHARACTER BSAMP*4,FSAMP*4
C      DATA BSAMP /'BKG '/
C      DATA FSAMP /'FOIL'/
C
C      Get operator ID; set abort flag if it is null
C
C      10      TYPE 1000
C 1000      FORMAT (' '//Enter operator ID (up to 8 chars)',
C      1          ' or RETURN to escape -> ')
C      READ (5,1002,END=10,ERR=10) OPERID
C 1002      FORMAT(A80)
C      ABORT = OPERID.EQ.' '
C      IF (ABORT) RETURN
C
C      See whether sample ID, material type, and isotopics are needed
C
C      IF (FROMDK) RETURN
C      IF (AC) RETURN
C      IF (BACKG) GO TO 100
C      IF (.NOT.MP) GO TO 120
C      IF (MP .OR. MB)GO TO 120
C
C      Supply "canned" sample ID for background or meas control
C
C 100      MTL=.FALSE.
C          SAMPID=FSAMP
C          IF (BACKG) SAMPID=BSAMP
C          SAMPID=' '
C          RETURN
C

```

C Get sample ID; set abort flag if it is null

C

120 MTL=.TRUE.

20 TYPE 1001

1001 FORMAT (' /' Read sample barcode or enter sample ID  
1 /' (RETURN to escape) -> ' )  
READ (5,1002,END=20,ERR=20) SAMPID  
ABORT = SAMPID.EQ.' '  
RETURN  
END

```

C*****
C GET90
C
C Routine specific to gamma ray instruments using
C Canberra Series 90
C Get a spectrum from the Series 90 and write to disk;
C Read chunks of 256 channels at a time
C
C SSJ 28-JUN-79
C Rev. 13-Jun-84 lc
C REV. 12-Dec-85 WDR Put start time and date in channels
C 3 and 4
C Rev. 19-Feb-85 MPK Change ACCEPT to READ with ERR and
C END exits
C Rev. 20-MAR-85 WDR Modified to do FROMCA rather than READ80
C Rev. 04-APR-86 WDR Modified for Micro VAX II
C Rev. 05-JUL-87 RDP Modified for uVAX II w/ ADAC Drivers
C*****
C
C SUBROUTINE GET90(NADC)
C
C LOGICAL*1 ANS,SMALL
C
C The following common block files are generic for all instruments
C
C INCLUDE 'CONTRL.CMN'
C INCLUDE 'CLUNS.CMN'
C INCLUDE 'CFILES.CMN'
C INCLUDE 'CASSAY.CMN'
C INCLUDE 'CTIM.CMN'
C INCLUDE 'CBACKG.CMN'
C INCLUDE 'CDATE.CMN'
C
C The following common block files are instrument specific for
C gamma-ray and x-ray multichannel instruments
C
C INCLUDE 'CSPECT.CMN'
C INCLUDE 'CONT90.CMN'
C
C Initialize for writing
C
C NSPEC = MEMSIZ
5 NRECS = NSPEC/128
C KTIMES = 0
C SMALL = .FALSE.
C DUM = -999.0
-D WRITE(5,5568)NSPEC,NRECS,NADC
D5568 FORMAT(' NSPEC = ',15,'NRECS = ',15,'NADC = ',15)
C
C Open the output file on disk first
C If any problems occur (such as the program bombing
C due to no disk), we do not want the Series 90
C to be opened up yet
C
C Allocate channel, free it when through
C
15 CONTINUE
OPEN(UNIT=DKLUN1,NAME=FILNAM,ACCESS='DIRECT',TYPE='NEW',

```

```

X          RECORDSIZE=128,MAXREC=NRECS+1,ERR=20)
  IREC = 1
  GO TO 25
C
C Problem in opening the file
20      CALL ERRMSG(1,FILNAM)
C
C Free channel IOCHAN
C
C      CALL IFREEC(IOCHAN)
C
C      KTIMES = KTIMES + 1
C      IF(KTIMES .LE. 1) GO TO 15
C      TYPE 1005
1005     FORMAT(/'#Enter A to abort writing spectra to disk',
X         '      -> ')
C      READ (5,1007,END=15,ERR=15)ANS
1007     FORMAT(A1)
C      IF (ANS .EQ. 'A') WRITNG = .FALSE.
C      IF (ANS .EQ. 'A') GOTO 500
C      GO TO 15
C
C Read one 256-channel block
C Check for errors in the read
C
25      ICHNF = 0
C      ICHNL = 127
30      CONTINUE
C      CALL FROMCA(NADC,ICHNL,ICHNF,SPCTRM,IERR,)
C      CALL CHK90('FROMCA',IERR)
C      NUMCHN = IERR
C      IERR = IERR-128
C      IF(IERR .EQ. 0 .AND. IREC .EQ.1) GO TO 32
C      IF(NUMCHN .NE. 128) TYPE 1002,NUMCHN
1002     FORMAT(/20X,'****', 16,' CHANNELS READ - EXPECTED 128 ****')
C      IF(IREC .GT. 1) GO TO 35
C      GO TO 500
C
C
C First record
C Check to see if file is big enough
C Change number of records for actual spectrum size
C
32      TYPE 1003,NSPEC,NSTART
1003     FORMAT(/15,' Channels starting at channel # ',15/)
C      IF(NSPEC/128 .GT. NRECS) GO TO 33
C      NRECS = NSPEC/128
C      GO TO 35
C
C This file is not large enough
C Go on to close the file and the series 90
C and start over
C
33      SMALL = .TRUE.
C      GO TO 100
C
C In the case of a cumulative run, insert the cumulative
C live and clock times in place of SPCTRM(1) & (2)
C
35      IF((.NOT. CUMLTV) .OR. (IREC .NE. 1)) GO TO 50

```

```
SPCTRM(1) = TIMELV
SPCTRM(2) = TIMECL
```

```
C
C      Insert the start time and date in SPCTRM(3) and (4)
```

```
C      Make sure it is the first record
```

```
50      IF (IREC .NE. 1) GO TO 55
        SPCTRM(3) = STIM
        SPCTRM(4) = SDAT
        TYPE *, 'SDAT IN GET90', SDAT
```

```
C      Write this 128 channel block to disk
```

```
55      WRITE(DKLUN1'IREC)SPCTRM
        IREC = IREC + 1
        ICHNF = ICHNL + 1
        ICHNL = ICHNL + 128
        IF(IREC .LE. NRECS) GO TO 30
```

```
C      Close the ADC and the disk file
C      Write current background data to disk
```

```
100     IF(SMALL) GO TO 5
        WRITE(DKLUN1'IREC)DUM,BK,SGBK,TIMEBK,BKDAT,BKTIM,
X         BKOP,BKSFIL,SAMPID,OPERID,SYSID,ASDATE,ASTIME,
X         MATTYP(1)
```

```
101     CLOSE(UNIT=DKLUN1)
        CALL IFREEC(IOCHAN)
```

```
        NBLKS = IREC
        TYPE 1006,NBLKS,FILNAM
```

```
1006     FORMAT(/' Spectrum retrieval completed --',
X         I6,' blocks written to ',15A)
```

```
500     RETURN
        END
```

```
C*****
C HEADR38      (HEADER subroutine for PU-238 isotopic system)
C
C      Print assay for PU-238 isotopic instrument
C
C      Calling arguments:
C          SUMMRY = Logical flag for writing summary for cycles
C
C      07-Jun-84 lc
C      04-Sept-84      TES  16:38
C      11-SEP-84      CMS
C      06-NOV-84      WDR      Modified for solid isotopic instrument
C      27-APR-90      WDR      Modified for Pu-238 isotopic instrument
C
C*****
```

```
C
C      SUBROUTINE HEADER(SUMMRY)
C
C      LOGICAL*1  SAMPLE,SUMMRY
C
C      The following common block files are generic for all instruments
C
C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CDATE.CMN'
C      INCLUDE 'CFILES.CMN'
C      INCLUDE 'CLUNS.CMN'
C      INCLUDE 'CONFIL.CMN'
```

```
C
C      Print the header showing date, time, operator, etc.
C
C      WRITE(LOUT,1999)
C 1999  FORMAT(/1X,80(1H*))
C
C      WRITE(LOUT,1065)RTYP
C 1065  FORMAT(1X,'PU-238 isotopic gamma assay - ',T37,'Run Type: ',
C      X      AB)
C
C      IF(.NOT. SUMMRY .AND. (AUTOCY .OR. MP)) WRITE(LOUT,1007)ICYCLE
C 1007  FORMAT(/' Cycle: ',I5)
C      IF(SUMMRY) WRITE(LOUT,1077) NUMCYC
C 1077  FORMAT(/19X,' Summary for ',I4,' Cycles')
C
C      SAMPLE = LIVE .AND. .NOT. FROMDK
C      IF(SAMPLE) WRITE(LOUT,1001) SAMPID,ASDATE,ASTIME
C      IF(.NOT. SAMPLE) WRITE(LOUT,1006) FILNAM,ASDATE,ASTIME
C 1001  FORMAT(/' Sample ID:      ',A12,T35,'Measurement Date: ',A9,2X,
C      X      AB)
C 1006  FORMAT(/' Input file:    ',A15,T35,'Measurement Date: ',A9,2X,
C      X      AB)
C
C      WRITE(LOUT,1005) OPERID,NOWDAT,NOWTIM
C 1005  FORMAT(' Operator ID:    ',A8,T35,'Current Date:      ',A9,2X,AB)
C
C      IF(WRITNG) WRITE(LOUT,1002) FILNAM
C 1002  FORMAT(' Data file:      ',A15)
C
```

```
WRITE(LDUT,1066) CONFIL, FILDAT, FILTIM
1066  X   FORMAT(' Constants File: ',A15,T35,' Constants Date: ',A9,2X,
           AB)
C
C
WRITE(LDUT,1003) TIMELV,TIMECD,REMARK
003  X   FORMAT(' Live time (sec):',F8.0,T35,' Clock time (sec):',F8.0,
           //1X,B0A1)
C
C
WRITE(LDUT,1999)
C
500  RETURN
      END
```

```

C
C*****
C
C  HELP38      (HELP subroutine for Pu-238 isotopic instrument)
C
C  Outputs to the terminal the list of menu options for the
C  operator.  If the supv flag is true, supervisor options
C  are listed.
C
C  Calling arguments:
C  SUPV = Logical flag passed from calling program
C         If SUPV = .TRUE. then operator can select from
C         supervisor options
C         If SUPV = .FALSE. then operator is not allowed
C         to select from supervisor options
C
C  CMS 26-APR-84
C  REV. 07-Sep-84 CMS
C  REV. 26-JUN-85 WDR      Modified for solids isotopic system
C                          Removed B and BD options
C  REV. 04-APR-86 WDR      Modified for Micro VAX II
C*****
C
C  SUBROUTINE HELP(SUPV)
C  IMPLICIT INTEGER (A - Z)
C  CHARACTER*1 REPLY
C
C  The following common block files are generic for all instruments
C
C  INCLUDE 'CLUNS.CMN'
C
C  Output the list of options - does not include supervisor options
C
C 1000 WRITE(5,1000)
C      FORMAT(/
C  X      ' A - Assay'//
C  X      ' MB - Measurement control - bias'//
C  X      ' MP - Measurement control - precision'//
C  X      '      Options during data acquisition:'//
C  X      ' Q - Quit assay - no results '//
C  X      ' S - Suspend assay temporarily'//
C  X      ' T - Terminate assay with results'//)
C
C  If supervisor mode output the rest of the options
C
C  IF(SUPV)THEN
C  TYPE 900,'Press RETURN to continue ->'
900  FORMAT(1X,A,#)
C  READ (5,920)REPLY
920  FORMAT(A1)
C  WRITE(5,1001)
1001  FORMAT(/ ' SUPERVISOR OPTIONS '//
C  X      '      '//
C  X      ' AC - Current data in MCA '//
C  X      ' AD - From disk'//
C  X      ' AU - Autocycle'//
C  X      ' AUD- Autocycle from disk'//
C  X      ' D - Default'//
C  X      ' LA - List assay log'//

```



```
X      / LM - List measurement control log'/  
X      / OU - Change output listing device'/  
X      / ST - Status'/  
X      / R  - Read data from disk'/  
X      / W  - Write data to disk'/  
X      / X  - Exit from program'/)
```

```
ENDIF
```

```
RETURN  
END
```

```
C  
C
```

```

C*****
C INIT90
C
C Routine specific to initializing Series 90
C
C SSJ 13-AUG-79
C
C REV 16-FEB-84 LC
C REV 22-OCT-84 MJC          CHANGED 'PUT MCA IN REMOTE . . . '
C                            TO 'PUT SERIES 90 ON LINE . . . '
C                            IN LABEL 1020
C Rev 19-Feb-85 WDR          Change ACCEPT to READ with END and
C                            ERR exit
C Rev 04-Apr-86 WDR          Modified for Micro VAX II
C Rev 05-Jul-87 RDP          Modified for uVAX II w/ ADAC Drivers
C

```

```

C*****
C
C SUBROUTINE INIT90
C
C LOGICAL*1 ANS
C INTEGER IARY(12)
C
C The following common block files are instrument specific for
C gamma ray instruments
C
C INCLUDE 'CONT90.CMN'
C
C DATA MCA/1/
C
C Try initializing the series 90 a maximum of 3 times
C
C ITIMES = 0
3 ITIMES = ITIMES + 1
  CALL INIT(MCA,,IERR)
  IF(IERR .EQ. 0) GO TO 10
  IF(ITIMES .GE. 3) GO TO 8
C
C Some problem
C Check if MCA is not in remote (ierr = -1)
C
  IF(IERR .NE. -1) GO TO 3
  TYPE 1020
- 1020 FORMAT(' //Put Series 90 on line - hit //return// when ready ',
  X          '-> ')
  READ (5,1001,END=3,ERR=3)
- 1001 FORMAT(A1)
  GO TO 3
C
C Print error msg and terminate
C
8 CALL CHK90('INIT ',IERR)
  TYPE 1030
1030 FORMAT('//20X, '**** SERIES 90 INIT UNSUCCESSFUL 3 TIMES'//
  X          10X, 'PROGRAM EXITING'//
  X          10X, 'CHECK YOUR "AN" HANDLER ****')
  RETURN
C

```

C Detach any tasks

C

10 DO 11 I=1,3

11 CONTINUE

C

20 RETURN

END

## SUBROUTINE INITLP

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C
C      GRPANL sub to begin line printer output of parameters & file names.
C      Called from READIN (before any peak groupings)
L
C      By J.B.Niday
C      Last edit for VAX 14-Mar-88 (J.B.N.)
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

      INCLUDE 'GRPLGCM.PMS'

      INCLUDE 'CFILES.CMN'
      INCLUDE 'CLUNS.CMN'
C
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

      BYTE IBUF(80)
      EQUIVALENCE (WM,IBUF),(WM(30),LTH)

      BYTE DETFNM(MXFNM),SHPFNM(MXFNM)
      EQUIVALENCE (DETFNM,RM),(SHPFNM,RM(MXFNM+1))

      DIMENSION XC(2)
      DATA LP/-1/
      DATA LP/6/
      - - - - -
      DEDTIM = 100. * (RLTIME-TIMEL)/RLTIME !Temp arg. only.

      Tabulate some of initialization stuff on printer
      CALL VERS(1)
      WRITE(LOUT,1010)(DTE(J),J=1,9),(TBUF(I),I=1,8),FILNAM,
+ (ICFIL(I),I=1,12)
1010 FORMAT(/30X,'Analysis on ',9A1,' at ',8A1,' of '//
+ ' Spectral data file: ',A,'Parameter control file: ',12A1)

      DO 5, I = 1,LTH
      IF (IBUF(I) .NE. '/') GO TO 5
      II = I
      GO TO 12
5      CONTINUE
      GO TO 15
12      WRITE(LP,1012)(IBUF(I),I=II,LTH)
1012 FORMAT(4X,' Using run-time control options: ',30A1)

15      WRITE(LP,1015)(EXPTID,DTZERO,IZYR)
      WRITE(LP,1020)(SMPID(I),I=1,12),(DSNAME(I),I=1,6),DELTAT

      IF((DETFNM(1) .NE.0).AND.(DETFNM(1) .NE. ' ')) THEN
          WRITE(LP,1030)(DETFNM(I),I=1,MXFNM)
      ELSEIF ((DFILNM(1) .NE.0).AND.(DFILNM(1) .NE. ' ')) THEN
          WRITE(LP,1035)(DFILNM(I),I=1,MXFNM)
      ENDIF

      IF (SHAPC(1) .EQ. 0) GO TO 20
      IF((SHPFNM(1) .NE.0).AND.(SHPFNM(1) .NE. ' ')) THEN
          WRITE(LP,1040)(SHPFNM(I),I=1,MXFNM)
      ELSEIF ((SFILNM(1) .NE.0).AND.(SFILNM(1) .NE. ' ')) THEN
          WRITE(LP,1041)(SFILNM(I),I=1,MXFNM)

```

ENDIF

```
20      WRITE(LP,1622)(BKGDFL(I),I=1,MXFNM-1)
        WRITE(LP,1625)TIMEL,DEDTIM
        WRITE(LP,1630)DSTTIM,ICTYR,DSTPTM,1STPYR
15      FORMAT(/' Expt: ',I2A1,9X,'Zero time: ',F12.6,' (',I4,')')
1020    FORMAT(' Sample: ',I2A1,6X,' Counted on: ',6A1,6X,'Decay',
+        ' (days to midpt.): ',1PE11.5)

1030    FORMAT(' Efficiency data updated from file: ',30A1)
1035    FORMAT(' Efficiency data originally obtained from file: ',30A1)

1040    FORMAT(' Shape parameters updated from file: ',30A1)
1041    FORMAT(' Shape parameters originally obtd. from file: ',30A1)

1622    FORMAT(' Background peak file used: ',29A1)
1625    FORMAT(' Live time',F9.2,' minutes with',F7.3,'% dead time')
1630    FORMAT(' Starting day',F9.4,' (',I4,') and ending on',F9.4,
+        ' (',I4,')')
        WRITE(LP,1635)MCHBUF,XTXBUF
1635    FORMAT(' Input text:',5X,32A1/1X,70A1)

        IF(RUNFLG(3) .GT. 0) WRITE (LP,1638)
1638    FORMAT(/' NOTE: Control file used was for a DIFFERENT detector',
+        '-system!!!!!!'//)
        IF(RUNFLG(4) .EQ. 0) WRITE (LP,1640)
1640    FORMAT(/' Photon calculations were not requested'//)

C      SUMMARIZE DATA FOR CALCULATION OF PHOTON INTENSITIES

        WRITE(LP,1643)DCNST(13),DCNST(14),(DSNAME(I),I=1,6)
        WRITE(LP,1645)GEOM,GEOM+DCNST(15),SURFC,DEPTH,RHO,SMPWT,WT100,
+        ANSFCT

1643    FORMAT(/' Detector parameters: Height =',F7.3,' cm.; Radius =',
+        F7.3,' cm. for ',6A1)
1645    FORMAT(' Sample midplane to detector window =',F7.3,' cm.;',
+        ' total dist. used =',F7.3,' cm.'/' Sample surface area =',F6.2,
+        ' ; Depth =',F6.3,' cm.; Density =',F7.3,' g/cm3'/' Weight =',
+        F11.6,'g; 100% weight =',F12.6,'g; Normalization factor =',
+        1PE12.4)

660    IF(NF1 .EQ. 0)GO TO 680
1660    FORMAT (/ ' Sample composition: ' / ' Element Amount(%)')
        WRITE(LP,1660)
        DO 675, K=1,NF1
        WRITE(LP,1675)MTLZ(K),CMPOS(K)
675    CONTINUE
1675    FORMAT(3X,A2,F12.4)

680    WRITE(LP,1680)
1680    FORMAT (/ ' Absorbers used: ',#)
        IF(NF2 .EQ. 0) THEN
            WRITE(LP,1681)
            GO TO 700
        ENDIF
31     FORMAT(' NONE'//)
        WRITE(LP,1682)
1682    FORMAT(/' Element & amount (g/cm2)')
        DO 685, K=1,NF2
```

```
685     WRITE(LP,1675)NABS(K),ABSRB(K)
700     IF(KTGRP.EQ.0)WRITE(LP,1001),GAIN,ZERO
1001    FORMAT(' Initialize at approx. GAIN of',F9.6,
+ ' keV/channel.; & chan: ZERO =',F9.4,' keV//')
      RETURN
      END
```

C\*\*\*\*\*

C  
C INTFLG -

C Initialize logical flags to be used in the assay  
C procedure. This is called when the system is first  
C initialized and after each option has been exercised.

C  
C CMS 7-MAY-84

C Rev 2-Oct-84 CMS Added CCAL.CMN Include statements

C REV 11-OCT-84 MJP Added common block CCAL

C Rev 08-Feb-85 MPK Initialize all flags that need  
C initialization

C Rev 13-Jun-85 MPK Change for NCC instruments

C Rev 19-Jul-85 MPK Make DGFLG an array; add CONDIR  
C and BAD flags

C Rev 04-Apr-86 WDR Modified for Micro VAX II

C\*\*\*\*\*

C  
C SUBROUTINE INTFLG

C The following common block files are generic for all instruments

C  
C INCLUDE 'CONTRL.CMN'  
C INCLUDE 'CBACKG.CMN'  
C INCLUDE 'CCAL.CMN'  
C INCLUDE 'CVAX.CMN'  
C INCLUDE 'CDIAG.CMN'  
C INCLUDE 'CMC.CMN'  
C INCLUDE 'CASSAY.CMN'

C initialize logical flags

C  
C ABORT = .FALSE.  
C BACKG = .FALSE.  
C FROMDK = .FALSE.  
C AUTO CY = .FALSE.  
C OLDDATA = .FALSE.  
C CUMLTV = .FALSE.  
C STPCYC = .FALSE.  
C ASFLG = .FALSE.  
C AC = .FALSE.  
C AUD = .FALSE.  
C AD = .FALSE.  
C MB = .FALSE.  
C MP = .FALSE.  
C TKEY = .FALSE.  
C CHANG = .FALSE.  
C CHANG1 = .FALSE.

C  
C BD = .FALSE.

C  
C CAL = .FALSE.  
C CALDSK = .FALSE.

C  
C VAXUP = .FALSE.  
C NOGO = .FALSE.

C

```
DO 100 I=1,10  
100 DGFLG(I)=.FALSE.  
.CONDIT = .FALSE.  
BAD = .FALSE.  
MCERR = .FALSE.
```

C

```
RETURN  
END
```



```
C*****
C JULDAL      (Livermore version of JULDAT)
C
C      Using ASCII strings for date, time, compute Julian day
C      From Sue Johnson's Julday
```

```
L      Calling arguments:
C      IDAT - Input ascii string for date 30-apr-84
C      ITIM - Input ascii string for time 13:54:17
C      XDAT - Output Julian day
C      XTIM - Output time of day in seconds
```

```
C      12-APR-83 LC
C      REV. 02-MAY-84 LC      for generic package
C      REV. 31-AUG-84 CMS
C      REV. 13-DEC-84 WDR      made day 0 1/1/1900
C                               modified days calculation
C                               Note: Day 0 1/1/1900
```

```
C      Rev. 11-Oct-84 mjc:      IDAT is accepted as 8 characters (1-May-84)
C                               or 9 characters (10-May-84)
C      Rev. 30-Oct-84 mpk      Fix error in 11-Oct correction -
C                               tighten code
C      Rev. 31-Oct-84 mpk      Tighten code some more, and handle
C                               time format more flexibly. This is
C                               done by copying date and time into
C                               an auxiliary array, and editing the
C                               array to convert all special characters
C                               into commas. Then the DECODE can use
C                               the comma-delimited integer field
C                               feature.
```

```
C      Also, correct Julian date algorithm
C      to use integer arithmetic.
C      In the correct algorithm, day 0 comes
C      out to 30-OCT-79, not 28-OCT as
C      stated previously.
```

```
C      Rev. 01-Nov-84 mpk      Leave spaces in date/time string
C                               rather than editing to commas.
C      Rev. 02-Nov-84 mpk      Latest installment in space saga:
C                               Ignore leading spaces and replace
C                               trailing spaces with commas; this
C                               makes it work right for time inputs
C                               like 11:23 (seconds omitted).
```

```
C      Rev. 10-DEC-85 WDR      Modified to calculate number of days
C                               from 1/1/1900 rather than as described
C                               above.
```

```
C      Rev. 3-Apr-86 WDR      Modified for Micro VAX II
```

```
C*****
```

```
C      SUBROUTINE JULDAT(IDAT,ITIM,XDAT,XTIM)
```

```
C      CHARACTER XMONTH*3
C      CHARACTER*1 IDAT(9),ITIM(8), DATTIM(23)
C      LOGICAL*1 TERMF8
C      DIMENSION NDAYS(12)
```

```
C      DATA NDAYS /31,28,31,30,31,30,31,31,30,31,30,31/
```

```
C      The following common block files are generic for all instruments
```

```

C      INCLUDE 'CMONTH.CMN'
C
C      Decode ASCII strings
C
      DO 5 I=1,9
5      DATTIM(I)=IDAT(I)
      DO 6 I=1,8
6      DATTIM(I+9)=ITIM(I)
D      TYPE 1998,DATTIM
C
C      Replace special chars ('-' or ':') and trailing spaces with commas
C
      TERMFG = .TRUE.
      DO 7 I=1,17
      IF (DATTIM(I).NE.' ') GO TO 71
      IF (TERMFG) GO TO 7
      GO TO 72
71     IF (DATTIM(I).GE.'0' .AND. DATTIM(I).LE.'9') GO TO 73
      IF (DATTIM(I).GE.'A' .AND. DATTIM(I).LE.'Z') GO TO 73
72     DATTIM(I)=','
      TERMFG = .TRUE.
      GO TO 7
73     TERMFG = .FALSE.
7     CONTINUE
D      TYPE 1998,DATTIM
D1998  FORMAT (' EDITED DATE/TIME: ',17A1)
      IDAY = 0
      XMONTH = ' '
      IYEAR = 0
      IHOOR = 0
      IMIN = 0
      ISEC = 0
      DECODE (9,1000,DATTIM,ERR=75) IDAY,XMONTH,IYEAR
1000  FORMAT (I3,A3,1X,I4)
75    DECODE (8,1001,DATTIM(10),ERR = 76) IHOOR,IMIN,ISEC
1001  FORMAT (2I5,I2)
76    CONTINUE
C
C      Compare month to common
C
      IMONTH = 0
      DO 9 I=1,12
9      IF(XMONTH.EQ.AMONTH(I))IMONTH = I
D      TYPE 1999,IMONTH,IDAY,IYEAR,IHOOR,IMIN,ISEC
D1999  FORMAT(' DATE:',3I3,' TIME:',3I3)
C
C      Convert time to seconds
C
      XTIM = ISEC + 60. * (IMIN + 60. * IHOOR)
C
C      Convert date to days
C
      IF (IYEAR .GT. 1900) IYEAR = IYEAR - 1900
      NMONTH = IMONTH - 1
10     XDAT = AINT (365.25 * IYEAR)
C      1 + AINT (30.6001 * NMONTH)
      2 + IDAY

```

```
C      3      + XTIM/86400.      !uncomment to include fractional day
C
C      Add the number of days up to this month
C
C      IF (NMONTH.EQ.0)GO TO 16
C      DO 15 J=1,NMONTH
13      XDAT = XDAT + NDAYS(J)
16      IF(NMONTH .LT. 3)GO TO 500
      YR = IYEAR
C
C      Include an extra day for leap year
C
      IF(AMOD(YR,4.).EQ. 0.)XDAT = XDAT + 1.
C
500      RETURN
      END
```

```

C*****
C
C  LGMENU
C      Lists menu options for the list assay log and list
C      measurement control log
C      Options are:
C          Return to main menu
C          Entries not sent to VAX
C          Entries between two dates
C          All entries
C          'N' entries
C
C      Returns number of option selected
C
C      Calling arguments:
C      Input:  LGNAME - name of log being listed, 23A1
C      Output: NANS - menu option number, I
C
C      CMS 14-MAY-84
C      REV 10-SEP-84 mjc
C      Rev 24-Jan-85 MPK      Add menu option 0 - "return to main menu"
C      Rev 25-Jan-85 MPK      Add code to handle variable-length name
C      Rev 04-Feb-85 MPK      Reset control/O before menu output
C      Rev 04-Apr-86 WDR      Modified for Micro VAX II
C      Rev 04-Apr-90 WDR      Modified for Pu-238 instrument; 1st
C                             option searches for sample ID
C
C*****

```

```

C
C      SUBROUTINE LGMENU(NANS,LGNAME)
C
C      LOGICAL * 1  LGNAME(23)
C
C      The following common block files are generic for all instruments
C
C      INCLUDE 'CLUNS.CMN'
C
C      DO 100 L=1,22
C      IF (LGNAME(L+1).EQ.0) GO TO 110
100  CONTINUE
C      L=23
C      110  TYPE 1000, (LGNAME(I),I=1,L)
1000  FORMAT('/' Menu of Options for Listing ',23A1)
C      TYPE 1001
1001  FORMAT('/' 0.  Return to main menu'
C      X      // 1.  List entries with specified sample ID'
C      X      // 2.  List entries between two dates'
C      X      // 3.  List all entries'
C      X      // 4.  List n entries')
10    TYPE 1002
1002  FORMAT('/' '$ Enter option number -> ')
C      READ(5,1003,END=10,ERR=10)NANS
1003  FORMAT(I)
C      IF(NANS .LT. 0 .OR. NANS .GT. 4) GO TO 10
C      RETURN
C      END

```

```
C
C*****
C
C LHEAD
C Outputs header information for the assay log list option
C
C CMS 14-MAY-84
C Rev 23-Aug-84 mjc
C Rev 11-Oct-84 cms Changed title output from 'Data'
C to 'Assay'
C Rev 04-Apr-86 WDR Modified for Micro VAX II
C Rev 05-Aug-87 WDR Modified format statements for
C character output
C*****
```

```
C
C SUBROUTINE LHEAD
```

```
C The following common block files are generic for all instruments
```

```
C INCLUDE 'CLUNS.CMN'
C INCLUDE 'CLOG.CMN'
C INCLUDE 'CDATE.CMN'
```

```
C CALL DATTIM(1)
```

```
WRITE(LOUT,1099)NOWDAT,NOWTIM
```

```
1099 FORMAT(/'1',10X,'Assay log as of ',A9,2X,A8///
X 1X,'Entry Date Time Sample Operator File ',
X ' Cycles Results 240 Err'/22X,'Id',50X,
X 'Flg',1X,'Flg'//)
```

```
RETURN
```

```
END
```

C\*\*\*\*\*

C  
C LNMN38

C This subroutine checks the signals from the Liquid Nitrogen  
C Liquid Level Monitors (Ortec 729A) in the Pu238 isotopic  
U instrument through the ADAC 1616 CCI interface board. If  
C bit #5 in the data register at address ZWA1: goes low,  
C this indicates a low liquid nitrogen level in detector.

C REV. 06-AUG-87 RDP Changes made to ADAC interface.

C\*\*\*\*\*

C SUBROUTINE LNMN38

C INTEGER\*4 VALUE, ISTAT, ISTATE  
C INTEGER\*2 IBUFF(2)

C The following common block is generic to all instruments

C INCLUDE 'CONTRL.CMN'

C Common Containing ADAC Channel Information

C INCLUDE 'ADAC.CMN'

C ABORT = .FALSE.

C First, clear data register

C VALUE = 96  
C ISTAT = BTSTW(VALUE, ISTATE, C\_DATA, IOSB, , , )  
C ESTATUS = ERROR\_CHECK(ISTAT, IOSB)  
C IF (.NOT. ESTATUS) TYPE\*, 'PROBLEM WITH ADAC (#)'

C Now test bit #5

C VALUE = 32  
C ISTAT = BTSTW(VALUE, ISTATE, C\_DATA, IOSB, , , )  
C WRITE (\*, \*) 'ISTATE', ISTATE  
C ESTATUS = ERROR\_CHECK(ISTAT, IOSB)  
C IF (.NOT. ESTATUS) TYPE\*, 'PROBLEM WITH ADAC (#5)'  
C IF (ISTATE.EQ.1) GOTO 500

C Bit #5 is low, send warning message

C TYPE 900  
900 FORMAT(' ')  
C TYPE 999  
999 FORMAT('/' '!!!WARNING!!!')  
C TYPE 1000  
1000 FORMAT(' !!!Liquid Nitrogen level low in detector dewar!!!')  
C TYPE 1500  
1500 FORMAT(' Press RETURN to continue ')  
C ACCEPT 1600  
1600 FORMAT(80A1)  
C ABORT = .TRUE.  
500 RETURN  
C END

C\*\*\*\*\*

C  
C LMCLOG

C List the measurement control log entries

C Options for listing:

- C All entries
- C 'N' entries
- C Entries not yet sent to VAX
- C Entries within a pair of dates
- C Entries are listed from current date

C Calling argument:

- C AUTO If .FALSE., the routine runs in interactive mode and prompts for listing options.
- C If .TRUE., there is no operator interaction or listing; all entries not yet sent to the VAX are sent.

C CMS 14-MAY-84

- C REV. 7-SEP-84 mjc
- C REV. 18-OCT-84 MJC closed LOUT if it equaled 6 (line printer)
- C Rev. 23-Jan-85 MPK Remove date syntax checking, and use JULDAT to check dates for correctness
- C Update communications to use pseudo-message types 'LB' (log bias) and 'LP' (log precision)
- C Rev. 24-Jan-85 MPK Allow multiple listings on a single entry
- C Add end and error exits to all terminal input
- C Rev. 25-Jan-85 MPK Change 'carr ret' in prompt to 'RETURN'
- C Allow escaping from date mode
- C Rev. 28-Jan-85 MPK Correct and clean up indexing and listing logic
- C Rev. 30-Jan-85 MPK Add automatic mode for sending entries to VAX
- C Add additional error checks for log file access
- C Check VAXUP when sending to VAX
- C Update MCVFLG in file when an entry is sent to the VAX successfully
- C Rev. 15-Feb-85 MPK Use YESNO to ask yes/no questions
- C Rev. 10-Jul-85 MPK Adapt to new log format
- C In option for listing entries not sent to VAX, do not list background entries
- C Rev. 04-Apr-86 WDR Modified Micro VAX II
- C Rev. 02-May-90 WDR Removed COMM calls for Pu238 instrument

C\*\*\*\*\*

C SUBROUTINE LMCLOG(AUTO)

C CHARACTER\*1 INPDTE(9)  
 C INTEGER\*2 NANS,ENTNUM  
 C INTEGER\*2 LPTR,FPTR,IPTR,INTPTR(100)  
 C LOGICAL\*1 VAXFLG,DATFLG,STRDFL  
 C LOGICAL\*1 AUTO,YESNO

C The following common block files are generic for all instruments

```
INCLUDE 'CLUNS.CMN'  
INCLUDE 'CLOG.CMN'  
INCLUDE 'CDATE.CMN'  
INCLUDE 'CFILES.CMN'  
INCLUDE 'CMC.CMN'  
INCLUDE 'CMONTH.CMN'  
INCLUDE 'CONTRL.CMN'  
INCLUDE 'CASSAY.CMN'  
INCLUDE 'CVAX.CMN'  
INCLUDE 'CDIAG.CMN'  
INCLUDE 'CMEAS.CMN'
```

Open log file MCLOGF for reading

```
IERR=1  
OPEN (UNIT=DKLUN,NAME=MCLOG,TYPE='OLD',ACCESS='DIRECT',  
1 RECORDSIZE=MCREC, ERR=90)
```

Read first record to get start and end pointers

```
150 IERR=2  
READ (DKLUN'1,ERR=90) LPTR,FPTR  
IERR=5  
IF (LPTR.LT.2 .OR. LPTR.GT.MCMAX) GO TO 90  
IF (FPTR.LT.2 .OR. FPTR.GT.MCMAX) GO TO 90  
IF (AUTO) GO TO 350
```

Message to operator to select a listing option

```
CALL LGMENU(NANS,'Measurement control log')  
IF (NANS.EQ.0) GO TO 500  
VAXFLG=.FALSE.  
DATFLG=.FALSE.
```

Option 201 goes to entries not yet sent to VAX  
Option 202 goes to entries between two dates  
Option 203 goes to all entries  
Option 204 goes to N entries

```
GO TO (201,202,203,204)NANS
```

List entries not yet sent to VAX

```
201 VAXFLG=.TRUE.  
GO TO 203
```

List entries between two dates

```
202 DATFLG=.TRUE.
```

Set flag for start date and get start date

```
800 STRDFL=.TRUE.  
TYPE 1003
```

1003 FORMAT('/ /' \$ Enter start date (most recent) [09-AUG-84] -> ( )

```
READ (5,1004,END=800,ERR=800) NCHAR,INPDTE
```

1004 FORMAT(Q,9A)

```
IF (NCHAR.EQ.0) GO TO 150  
CALL JULDAT(INPDTE,NOWTIM,SDAT,XTIM)
```



GO TO 820

C

C Reset flag for start date and get stop date

C

```
810 STRDFL=:FALSE:
    TYPE 1005
1005 FORMAT(' '//* Enter stop date (oldest) [02-AUG-84] -> ')
    READ (5,1004,END=810,ERR=810) NCHAR,INPDTE
    IF (NCHAR.EQ.0) GO TO 150
    CALL JULDAT(INPDTE,NOWTIM,EDAT,XTIM)
```

C

C Check the start or stop date for validity,

C using the numbers JULDAT saves in IYEAR, IMONTH, and IDAY

C

```
820 IF (IYEAR.GE.84 .AND. IYEAR.LE.99) GO TO 830
    TYPE 1020
1020 FORMAT (//20X,'**** INCORRECT YEAR NUMBER ****')
    GO TO 850
830 IF (IMONTH.GE.1 .AND. IMONTH.LE.12) GO TO 840
    TYPE 1021
1021 FORMAT (//20X,'**** INCORRECT MONTH NAME ****')
    GO TO 850
840 IF (IDAY.GE.1 .AND. IDAY.LE.31) GO TO 860
    TYPE 1022
1022 FORMAT (//20X,'**** INCORRECT DAY NUMBER ****')
850 IF (STRDFL) GO TO 800
    GO TO 810
860 IF (STRDFL) GO TO 810
    IF (SDAT.GE.EDAT) GO TO 203
    TYPE 1023
1023 FORMAT (//20X,'**** START DATE IS EARLIER THAN STOP DATE ****')
    GO TO 800
```

C

C List N entries

C

```
204 TYPE 1006,MCMAX-1
1006 FORMAT(' '* How many entries do you want listed ? (1, ->
    1 I3,' maximum) -> ')
    READ (5,1000,END=205,ERR=205) N
1000 FORMAT (I4)
    IF (N .EQ. 0) GO TO 150
    IF (N .GE. 1 .AND. N .LT. MCMAX) GO TO 220
```

C

C Error message for improper input

C

```
205 TYPE 1009
1009 FORMAT(' '//
    1 ' **** ANSWER WAS NON-NUMERIC, NEGATIVE, OR TOO LARGE ****')
    GO TO 204
```

C

C List all entries

C

```
203 N = MCMAX - 1
```

C

C Output header for listing

C

```
220 CALL MCHEAD
```

C

C Main loop for listing entries

C

```
LGNUM = 1
IPTR=LPTR
```

```
DO 300 I = 1, N
IERR=2
READ (DKLUN/IPTR,ERR=90)
1   ASDATE, ASTIME, SAMPID, OPERID,   MTYP, MCRUNS,
2   VALUE,   SIGMA,   VAR, ARCHV1, ARCHV2, DGFLG,
3   BIASOK, PRECOK, BACKOK, MCERR, MCVAX
IF (.NOT.DATFLG) GO TO 230
CALL JULDAT(ASDATE,ASTIME,XDAT,XTIM)
IF (XDAT.GT.SDAT) GO TO 240
IF (XDAT.LT.EDAT) GO TO 400
230 IF ( VAXFLG .AND. (MCVAX .OR. MTYP(1).NE.'M') ) GO TO 240
WRITE (LOUT,1002) LGNUM,
1   ASDATE,ASTIME,
2   OPERID,MTYP,VALUE,SIGMA,VAR,
3   MCRUNS,MCVAX,MCERR
1002 FORMAT (I4,1X,A6,1X,A5,1X,A8,1X,2A,1PE13.5,
1   1PE11.3,E14.6,I3,1X,L,1X,L)
INTPTR(LGNUM) = IPTR
LGNUM = LGNUM + 1
240 IF (IPTR.EQ.FPTR) GO TO 400
IPTR=IPTR-1
IF (IPTR.LT.2) IPTR=MCMAX
300 CONTINUE
```

If the line printer is open, close it in order to finish printing

```
400 IF(LOUT .EQ. 6)CLOSE(UNIT = LOUT)
```

Prompt for sending entries to VAX

```
IF (.NOT.YESNO('Send entry to VAX')) GO TO 150
IF (YESNO('All files not yet sent')) GO TO 350
```

```
380 TYPE 1014
```

```
1014 FORMAT(' '//# Type number of entry (RETURN to quit).=>' ')
READ (5,1000,END=380,ERR=380) ENTNUM
IF (ENTNUM.EQ.0) GO TO 150
```

Check if number entered is in the range of entries listed

```
IF (ENTNUM .GE. 1 .AND. ENTNUM .LT. LGNUM) GO TO 390
```

Error message for number less than one or greater than maximum entries listed

```
TYPE 1015,ENTNUM
```

```
1015 FORMAT('/' **** ENTRY NUMBER',I4,' WAS NOT LISTED ****'/)
GO TO 380
```

Read the entry from the file and send it to the VAX

```
390 IERR=2
```

```
READ (DKLUN/INTPTR(ENTNUM),ERR=90)
1   ASDATE, ASTIME, SAMPID, OPERID,   MTYP, MCRUNS,
2   VALUE,   SIGMA,   VAR, ARCHV1, ARCHV2, DGFLG,
3   BIASOK, PRECOK, BACKOK, MCERR, MCVAX
IF (MTYP(2).NE.'B') GO TO 391
```

```
CALL COMM('RB')           !*communications - log bias
```

```

GO TO 393
391 IF (MTYP(2).NE.'P') GO TO 392
C CALL COMM('RP') !*communications - log precision
393 IF (.NOT.VAXUP) GO TO 150
IF (MCVAX) GO TO 380
MCVAX=.TRUE.
IERR=3
WRITE (DKLUN'INTPTR(ENTNUM),ERR=90)
1 ASDATE, ASTIME, SAMPID, OPERID, MTYP, MCRUNS,
2 VALUE, SIGMA, VAR, ARCHV1, ARCHV2, DGFLG,
3 BIASOK, PRECOK, BACKOK, MCERR, MCVAX
GO TO 380
392 TYPE 1030,MTYP
1030 FORMAT (//20X,'**** ENTRY TYPE ',2A1,
1 ' CANNOT BE SENT TO VAX ****'/)
GO TO 380

C
C All entries not sent to VAX yet will be sent here
C
350 IPTR=LPTR
C
355 IERR=2
READ (DKLUN'IPTR,ERR=90)
1 ASDATE, ASTIME, SAMPID, OPERID, MTYP, MCRUNS,
2 VALUE, SIGMA, VAR, ARCHV1, ARCHV2, DGFLG,
3 BIASOK, PRECOK, BACKOK, MCERR, MCVAX
IF (MCVAX) GO TO 361
IF (MTYP(2).NE.'B') GO TO 363
C CALL COMM('RB') !*communications - log bias
GO TO 364
363 IF (MTYP(2).NE.'P') GO TO 361
CALL COMM('RP') !*communications - log precision
364 IF (VAXUP) GO TO 362
IF (AUTO) GO TO 500
GO TO 150
362 MCVAX=.TRUE.
IERR=3
WRITE (DKLUN'IPTR,ERR=90)
1 ASDATE, ASTIME, SAMPID, OPERID, MTYP, MCRUNS,
2 VALUE, SIGMA, VAR, ARCHV1, ARCHV2, DGFLG,
3 BIASOK, PRECOK, BACKOK, MCERR, MCVAX
361 IF (IPTR.EQ.FPTR) GO TO 365
IPTR=IPTR-1
IF (IPTR.LT.2) IPTR=MCMAX
GO TO 355
365 IF (AUTO) GO TO 500

C
C Check if more are to be sent
C
IF (YESNO('Send more')) GO TO 380
GO TO 150

C
C ** Log file handling errors **
C
90 CALL ERRMSG(IERR,MCLOG)
C
Clean up and exit
L
500 CLOSE (UNIT=DKLUN,ERR=501)
501 RETURN

```

END

## SUBROUTINE LODGRP

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C
C Sub to load a specified peak grouping and remove its background:
C Optionally, fix any specified energy & intensity relationships
L between peaks. In the case of single peaks specified for simple
C integration only, sum around the centroid and return.
C By J.B.Niday Edited by R. Gunnink 1-June-84
C Last edit 13-Sep-85 (R.G.) PKHT in DO 200 loop.
C Last edit 17-Jun-87 (J.B.N.)
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

```

```
INCLUDE 'GRPLGCM.PMS'
```

```
INCLUDE 'CONTRL.CMN'
```

```
INCLUDE 'CLUNS.CMN'
```

```

C Check for size of group (EDIGRP cannot do)
OFFSET = IST - 1.0
NDPTS = IEND-IST+1
CHANS = NDPTS
IF(NDPTS .GT. MXDP) THEN
    TYPE 1002,KTGRP,IST,IEND,MXDP
    EST = 0.
    GO TO 400
ENDIF
1002 FORMAT (' Group',I3,', (channels',I5,' to',I5,') exceeds limit',
+ ' of',I5,' data points--REJECT IT!')

```

```
C *** Run-time modification of shape params that change with energy
```

```

GSQ = GAIN * GAIN
IF(NUSHPS .NE. 0)GO TO 13
IF (SHAPC(1) .EQ. 0)GO TO 10 !Flag for SHAPC file
AVEN = 0.5*(IST+IEND) * GAIN + ZERO
SG = (SHAPC(1) + SHAPC(2)*EST)/GSQ +.462
FWHM = SQRT(SG) * GAIN
ATAIL = EXP(SHAPC(3) + AVEN * SHAPC(4))
BTAIL = SHAPC(5) * GAIN
EXP3 = SHAPC(7)
EXP4 = SHAPC(8) * GAIN
GO TO 16

```

```

10 IF (INTFLG .NE. 0) GO TO 13
EST = 0. !Set to abort the group!
TYPE 1010,KTGRP
1010 FORMAT(' *** Cannot fit peaks without some shape',
+ ' information!!!') Abort group #',I4)
IF(IDOFLG(3) .EQ. 0) WRITE(LOUT,1010)KTGRP
GO TO 400

```

```

13 IF(SHAPC(1) .EQ. 0) SHAPC(2) = .0017
C Set a default for slope of peak width parameter.
POS = NDPTS/2. !For FWHM at mid-energy of grouping.
IF (NP .EQ. 1) POS = EN(1)/GAIN - OFFSET !At peak position
SG = FWHM**2/GSQ - SHAPC(2)*POS + .462
BTAIL = TSHP5 * GAIN
EXP4 = TSHP8 * GAIN

```

```

16 ALPHA = -2.7726/SG
EXP1 = ATAIL
EXP2 = BTAIL

DO 19, I=1, NP
Convert intrinsic widths from eV to channel values.
GAMA(I) = GAMA(I) * .001/GAIN
19 CONTINUE

C Determine peak grouping location and read in data, after first
C determining the backgrounds to be subtracted.
C Background areas were specified for channels adjacent to the
C group data, somewhat separated ("remote"), or both. In the latter
C case, the slope for extrapolation is to be determined.

XTRAPL = 0      !Store extapolation dist. for remote cases.
XTRAPH = 0
NUMR = 0
NUMS = BKLO / GAIN + 0.5      !Note: These values may be zero.
NUME = BKHI / GAIN + 0.5

IF((NUMS .LE. 0) .AND. (BKLO .NE. 0)) NUMS = 1
C Modified to protect against high gain cases (J.B.N. 24-Jan-85)
IF((NUME .LE. 0) .AND. (BKHI .NE. 0)) NUME = 1

C Convert bkgrd slopes to fractions/channel.
SLPLO = 0.01 * SLP(1) * GAIN      !Note: SLP may be zero
SLPHI = 0.01 * SLP(2) * GAIN

C Detn. average background level in front of peak grouping.
C Correct sloping backgd to channel nearest data points

IF(RBKEN(1) .NE. 0) GO TO 30
ISTLO= IST-NUMS      !No remote background used on low E side.

CALL RDCHNS(Y,ISTLO,NUMS,INPDEV)

CALL AVE (1,NUMS,SM,AVBG1,Y)

TLS = EXP4/GAIN
IF (EN(1) .LT. 4500.) GO TO 36
IF ((TLS .LT. 0.1).AND.(TLS .GT. 0.)) AVBG1 = 0.01      !Assume Alpha Spectr

GO TO 36

C A remote low E bkgrd has been specified.
30 ISTLO = ((RBKEN(1) - ZERO) / GAIN) + 0.5
NUMR = RBKWID(1) /GAIN + 0.5
NPTS = IST - ISTLO

CALL RDCHNS (Y,ISTLO,NPTS,INPDEV)

CALL AVE(1,NUMR,SM,AVBGLO,Y)

IF(BKLO .NE. 0) THEN      !Both bkgds given; calc. slope
CALL AVE(NPTS-NUMS,NPTS,SM,AVBG1,Y)
ADIST = NPTS - NUMR/2. - NUMS/2.
SLPLO = (AVBG1-AVBGLO) / (AVBGLO * ADIST)
SLP(1) = SLPLO / (.01*GAIN)      !Preserve SLP(1) for printout in
ELSE

```

```

XTRAPL = NPTS - NUMR/2. !Remote bkgd only!
BGLO = AVBGLO * (1. + SLPLO * XTRAPL)
SLPLO = SLPLO * AVBGLO !Finally, convert to cts/chan
NUMS = NUMR !To fix BPTS & carry forward for OUT
GO TO 40
ENDIF

36 BGLO = AVBG1 * (1. + SLPLO * NUMS/2.)
SLPLO = SLPLO * AVBG1 !Finally, convert to cts/chan

C Now detn avg. bkgd. level on high E side of group
40 NUMR = 0
IF(RBKEN(2) .NE. 0) GO TO 43

CALL RDCHNS(Y,IEND+1,NUME,INPDEV)

CALL AVE (1,NUME,SM,AVBG2,Y)
GO TO 48

C A remote bkgd specified on high E side.
43 ISTHI = ((RBKEN(2) - ZERO) / GAIN) + 0.5
NUMR = RBKWID(2) / GAIN + 0.5
NDXHR = ISTHI - IEND
NPTS = NUMR + NDXHR - 1

CALL RDCHNS(Y,IEND+1,NPTS,INPDEV)

CALL AVE(NDXHR,NPTS,SM,AVBGHI,Y)

IF(BKHI .NE. 0) THEN
CALL AVE(1,NUME,SM,AVBG2,Y) !Both bgds used
ADIST = NPTS - NUMR/2. - NUME/2.
SLPHI = (AVBGHI - AVBG2) / (AVBGHI*ADIST)
SLP(2) = SLPHI / (.01 * GAIN)
ELSE

XTRAPH = NPTS - NUMR/2. !Remote bkgd only!
BGHI = AVBGHI * (1. - SLPHI * XTRAPH)
SLPHI = SLPHI * AVBGHI
NUME = NUMR
GO TO 50
ENDIF

48 BGHI = AVBG2 * (1. - SLPHI * NUME/2.)
SLPHI = SLPHI * AVBG2

C Remove background continuum from peak grouping; return net cts.

50 CALL RDCHNS(Y,IST,NDPTS,INPDEV)
52 CALL BKGRD(1,NDPTS,BGLO,BGHI,SLPLO,SLPHI,NPTS,SMY,Y,YNET)
C TYPE *, 'BGLO', BGLO, BGHI, SLPLO, SLPHI
BPTS = NUMS + NUME
AVEBG = (BGLO + BGHI) / 2.
IF(SMY .LT. 0)SMY = 0.
ER= SQRT(SMY + AVEBG*CHANS*(1. + CHANS/BPTS))
ERR(1) = ER !Used only if SMY equal to 0.0
ER = ER / (SMY+1.) * 100.
IF(ER .GT. 100.) SMY = 0.0 !Less than 50% confidence.
IF(SMY .NE. 0)ERR(1) = ER

```

PKPO = (EN(1) - REFEN)/GAIN + REFCH

IF((SMY.EQ.0.) .OR. (ERR(1).GE.50)) THEN

IF(INTFLG .NE. 0)GO TO 92 !Scheduled for integration.

WRITE (LOUT,1068)KTGRP,ERR(1),SMY

IF(IDFLG(3) .EQ. 0)WRITE (LOUT,1068)KTGRP,ERR(1),SMY

GO TO 92 !Even though area was to be fitted.

ENDIF

1068 FORMAT(/' !! For group #',I3,' Error=',F9.3,' SumY=',F10.3/

+ ' Will skip fitting & just integrate group!')

IF(INTFLG .EQ. 0) GO TO 100 !Peak scheduled to be fitted.

C C C C C C C C C

Simply locate peak for integration.

CALL MAXVAL(1,NDPTS,JPK,YMAX,YNET)

IF((JPK .EQ. 1) .OR. (JPK .EQ. NDPTS))GO TO 80

Y2 = 0.5 \* YMAX

M = 0

DO 70, J = JPK,NDPTS

M = M + 1

IF(Y(JPK+M) .LE. Y2)GO TO 75

70 CONTINUE

GO TO 80

75 IF((JPK-M) .GE. 1)GO TO 85

80 PKPO = JPK

TYPE 1080,KTGRP

1080 FORMAT(' For group',I3,' peak position will be inexact')

GO TO 90

CALL GFIT(JPK,M,ALPHA,ATAIL,BTAIL,PKPO,PHT,ER,YNET)

FWHM = GAIN\* SQRT(-2.7726/ALPHA)

IF(ER .EQ. -2.)GO TO 80

C Assume slow variation within group of any non-linearity effect.

90 PKPO = PKPO + OFFSET

92 INTFLG = 1 !Turns on flag, in case peak scheduled for fitting.

IF (LONGPR)WRITE (LOUT,1090)IST,IEND,PKPO,SMY,ERR(1)

1090 FORMAT(' By simple integration of region:',I5,' to',I5,/9X,

+ 'Peak chan =',F8.2,'; Net counts =',F12.0,' +- ',F7.2,'%')

C Transfer results of integration now to simplify storage in OUT

PKPOS(1) = PKPO

SUM(1) = SMY

ENJ(1) = REFEN + (PKPO-REFCH) \* GAIN

GO TO 400

C C C C C C C C C C

100 DO 200 J = 1,NP

PKPOS(J) = (EN(J) - REFEN)/ GAIN + REFCH - OFFSET

NPK = PKPOS(J) + 0.5

PKHT(J) = YNET(NPK)

IF(HIGHT(J).LE.0.) GOTO 102

PKHT(J) = HIGHT(J) / SQRT(-3.14/ALPHA)

WL = GAMA(J) \* SQRT(-ALPHA)

IF (WL .GT. 0.) PKHT(J) = PKHT(J)\*(1.-.5642\*WL+.25\*WL\*\*2)

GO TO 200

C Adjust relative intensities for local variation of shapes



```

102 K = -NYFLG(J)
S1 = SHAPC(1)
IF(S1 .EQ. 0.)S1 = SG * GSQ
HIGHT(J)=HIGHT(J)*SQRT((S1+SHAPC(2)*EN(K))/(S1+SHAPC(2)*EN(J)))
RATIO = 0.0
IF ((GAMA(K) .EQ. 0.) .AND. (GAMA(J) .EQ. 0.)) GOTO 200
IF ((GAMA(K) .GT. 0.) .AND. (GAMA(J) .GT. 0.)) GOTO 200
C X ray and Gamma Ray are being related. Make height adjustment.
WL = (GAMA(J)+GAMA(K))*SQRT(-ALPHA) !One of GAMA's = 0.0
RATIO = 1.-.5642*WL +.25*WL**2 !X ray/Gamma peak height ratio.
IF (GAMA(K) .GT. 0.0) RATIO = 1. / RATIO
HIGHT(J) = HIGHT(J) * RATIO
200 CONTINUE
FPKHT=0.

```

```

C FIX THE INTENSITY RELATIONSHIPS THAT HAVE BEEN SPECIFIED
C IN THE INPUT. FIRST READJUST PEAK HEIGHT
C ESTIMATES FOR INTERFERENCES FROM NEIGHBORING PEAKS.

```

```

DO 250 J=1,NP
IF(NYFLG(J) .GT. 0)GO TO 250
IF(HIGHT(J) .GE. 0)GO TO 250
K = -NYFLG(J)
DX = PKPOS(J) - PKPOS(K)
FCT = -HIGHT(J) * EXP(ALPHA * DX**2)
IF (HIGHT(K) .LE. 0.) PKHT(K) = PKHT(K) / (1. + FCT)
PKHT(J) = -HIGHT(J) * PKHT(K)
250 CONTINUE

```

```

DO 300 J = 1,NP
IF(NYFLG(J) .LT. 0)GO TO 300
L = J - 1
CONST = 0.
IF(L .LE. 0)GO TO 290
C Is neighbor related?
270 IF(NYFLG(L) .GE. 0)GO TO 280
IF(-NYFLG(L) .EQ. J)GO TO 290
280 DX = PKPOS(J) - PKPOS(L)
EXPN = EXP(ALPHA * DX**2)
PKHT(J) = PKHT(J) - PKHT(L) * EXPN *(1.- CONST * EXPN)
290 IF(L .GT. J)GO TO 300
CONST = 0.5
L = J + 1
IF(L .LE. NP)GO TO 270
300 CONTINUE

```

```

C FIX THE ENERGY RELATIONSHIPS THAT HAVE BEEN
C SPECIFIED IN THE INPUT.

```

```

DO 390 J =1,NP
IF(NXFLG(J) .GE. 0)GO TO 390
L = -NXFLG(J)
PKPOS(J) = PKPOS(L) -(EN(L) - EN(J))/GAIN
390 CONTINUE
ITER = 1

```

```

400 RETURN
END

```

C\*\*\*\*\*

C  
C LSL068 (LSLOG subroutine for Pu-238 isotopic assay system.  
C Includes listing of isotopic abundances.)

C List the assay log entries  
C Options for listing:  
C All entries  
C 'N' entries  
C Entries not yet sent to VAX  
C Entries within a pair of dates  
C Entries are listed from current date

C Calling argument:  
C AUTO If .FALSE., the routine runs in interactive mode  
C and prompts for listing options.  
C If .TRUE., there is no operator interaction or listing;  
C all entries not yet sent to the VAX are sent.

C CMS 14-MAY-84

C Rev 27-SEP-84 mjc  
C Rev 25-sep-84 cms  
C Rev 18-OCT-84 MJC closed LOUT if it equaled 6 (line printer)  
C Rev 25-Jan-85 MPK Remove date syntax checking, and  
C use JULDAT to check dates for  
C correctness

C Update communications to use pseudo-  
C message type 'LA' (log assay)  
C Allow multiple listings on a single  
C entry  
C Add end and error exits to all  
C terminal input

C Allow escaping back to menu from  
C most input prompts

C Rev 28-Jan-85 MPK Correct and clean up indexing and  
C listing logic

C Rev 29-Jan-85 MPK Add automatic mode for sending  
C entries to VAX  
C Add additional error checks for  
C log file access  
C Check VAXUP when sending to VAX  
C Update LGFLAG in file when an  
C entry is sent to the VAX successfully

C Rev 15-Feb-85 MPK Use YESNO to ask yes/no questions  
C Rev 20-Feb-85 CMS Added error bar and assay type variables  
C to log entries

C Rev 27-Feb-85 MPK Check that log entry has proper type  
C before sending to VAX

C Rev 20-Mar-85 MPK Correct interchange of variables  
C in write statement

C Rev 08-Jul-85 MPK Adapt to new log format

C Rev 10-Jul-85 MPK In option for listing entries not  
C sent to VAX, list only "AS" entry type

C Rev 19-Sep-85 WDR Added isotopic abundances to the output

C Rev 04-Apr-86 WDR Modified for Micro VAX II

C Rev 13-Jul-87 RDP Modified for uVAX II

C Convert BYTE strings to CHARACTER strings

C Rev 20-Mar-90 WDR Modified for Pu-238 system; changed  
C LVFLAG to C240FG

```
C
C *****
C
C SUBROUTINE LSLOG(AUTO)
```

```
CHARACTER*1 INPDTE(9)
CHARACTER*8 SMPID
INTEGER*2 NANS,ENTNUM
INTEGER*2 LPTR,FPTR,IPTR,INTPTR(100)
LOGICAL*1 VAXFLG,DATFLG,STRDFL,NAMFLG
LOGICAL*1 AUTO,YESNO,CYCBUF(3)
```

```
C
C
C The following common block files are generic for all instruments
```

```
C
C
C INCLUDE 'CLUNS.CMN'
C INCLUDE 'CLOG.CMN'
C INCLUDE 'CDATE.CMN'
C INCLUDE 'CFILES.CMN'
C INCLUDE 'CASSAY.CMN'
C INCLUDE 'CMONTH.CMN'
C INCLUDE 'CONTRL.CMN'
C INCLUDE 'CVAX.CMN'
C INCLUDE 'CDIAG.CMN'
C INCLUDE 'CMEAS.CMN'
C INCLUDE 'CABUND.CMN'
```

```
C
C Open log file ASYLOG for reading
```

```
C
C IERR=1
C OPEN (UNIT=DKLUN,NAME=ASLOG,TYPE='OLD',ACCESS='DIRECT',
C 1 RECORDSIZE=LOGREC,ERR=90)
```

```
C
C Read first record to get start and end pointers
```

```
C
C 150 IERR=2
C READ (DKLUN'1,ERR=90) LPTR,FPTR
C IERR=5
C IF (LPTR.LT.2 .OR. LPTR.GT.MAXASY) GO TO 90
C IF (FPTR.LT.2 .OR. FPTR.GT.MAXASY) GO TO 90
C IF (AUTO) GO TO 350
```

```
C
C Message to operator to select a listing option
```

```
C
C CALL LGMENU(NANS,'Assay log')
C IF (NANS.EQ.0) GO TO 500
C VAXFLG=.FALSE.
C DATFLG=.FALSE.
C NAMFLG=.FALSE.
```

```
C
C Option 201 goes to entries not yet sent to VAX
```

```
C Option 202 goes to entries between two dates
```

```
C Option 203 goes to all entries
```

```
C Option 204 goes to N entries
```

```
C
C GO TO (201,202,203,204)NANS
```

```
C
C List entries WITH SPECIFIC SAMPLE ID
```

```

201     TYPE 1101
1101    FORMAT(' '//& Enter sample ID to search for ->')
        READ (5,1102,END=800,ERR=800)SMPID
1102    FORMAT(A)
        IF (SMPID.EQ.' ') GO TO 150
        NAMFLG = .TRUE.
        GO TO 203

C
C List entries between two dates
C
202     DATFLG=.TRUE.

C
C Set flag for start date and get start date
C
800     STRDFL=.TRUE.
        TYPE 1003
1003    FORMAT(' '//& Enter start date (most recent) [09-AUG-84] -> ')
        READ (5,1004,END=800,ERR=800) NCHAR,INPDTE
1004    FORMAT(Q,9A)
        IF (NCHAR.EQ.0) GO TO 150
        CALL JULDAT(INPDTE,NOWTIM,SDAT,XTIM)
        GO TO 820

C
C Reset flag for start date and get stop date
C
810     STRDFL=.FALSE.
        TYPE 1005
1005    FORMAT(' '//& Enter stop date (oldest) [02-AUG-84] -> ')
        READ (5,1004,END=810,ERR=810) NCHAR,INPDTE
        IF (NCHAR.EQ.0) GO TO 150
        CALL JULDAT(INPDTE,NOWTIM,EDAT,XTIM)

C Check the start or stop date for validity,
C using the numbers JULDAT saves in IYEAR, IMONTH, and IDAY
C
820     IF (IYEAR.GE.84 .AND. IYEAR.LE.99) GO TO 830
        TYPE 1020
1020    FORMAT (//20X,'**** INCORRECT YEAR NUMBER ****')
        GO TO 850
830     IF (IMONTH.GE.1 .AND. IMONTH.LE.12) GO TO 840
        TYPE 1021
1021    FORMAT (//20X,'**** INCORRECT MONTH NAME ****')
        GO TO 850
840     IF (IDAY.GE.1 .AND. IDAY.LE.31) GO TO 860
        TYPE 1022
1022    FORMAT (//20X,'**** INCORRECT DAY NUMBER ****')
850     IF (STRDFL) GO TO 800
        GO TO 810
860     IF (STRDFL) GO TO 810
        IF (SDAT.GE.EDAT) GO TO 203
        TYPE 1023
1023    FORMAT (//20X,'**** START DATE IS EARLIER THAN STOP DATE ****')
        GO TO 800

C
C List N entries
C
204     TYPE 1006,MAXASY-1
1006    FORMAT(//& How many entries do you want listed ? (1,
1         13,' maximum) -> ')
        READ (5,1000,END=205,ERR=205) N

```

```
1000  FORMAT (I4)
      IF (N .EQ. 0) GO TO 150
      IF (N .GE. 1 .AND. N .LT. MAXASY) GO TO 220
```

```
C
C Error message for improper input
```

```
205  TYPE 1009
1009  FORMAT(' ' /
1     ' **** ANSWER WAS NON-NUMERIC, NEGATIVE, OR TOO LARGE ****' )
      GO TO 204
```

```
C
C List all entries
```

```
203  N = MAXASY - 1
```

```
C
C Output header for listing
```

```
220  CALL LHEAD
```

```
C
C Main loop for listing entries
```

```
      LGNUM = 1
      IPTR=LPTR
```

```
      DO 300 I = 1, N
      IERR=2
```

```
      READ (DKLUN,IPTR,ERR=90)
```

```
1     ASDATE, ASTIME, SAMPID, OPERID, MATTYP, ASTYPE,
2     LGFILE, ICYCLE, NUMCYC, ASANS, ERRBAR, DGFLG,
3     BIASOK, PRECOK, BACKOK, BAD, CONdit, C240FG,
4     (WTPCT(J), PCT(J),J=1,6)
```

```
      IF (.NOT.DATFLG) GO TO 228
      CALL JULDAT(ASDATE,ASTIME,XDAT,XTIM)
      IF (XDAT.GT.SDAT) GO TO 240
      IF (XDAT.LT.EDAT) GO TO 400
```

```
228  IF (.NOT.NAMFLG)GOTO 230
      IF (SMPID.EQ.SAMPID)GOTO 232
      GOTO 240
```

```
230  IF ( VAXFLG .AND. (ASTYPE(1).NE.'A') ) GO TO 240
      ENCODE (3,1001,CYCBUF) NUMCYC
```

```
1001  FORMAT (I3)
```

```
231  IF (CYCBUF(1).NE.' ') GO TO 232
```

```
      CYCBUF(1)=CYCBUF(2)
      CYCBUF(2)=CYCBUF(3)
      CYCBUF(3)=' '
      GO TO 231
```

```
232  WRITE (LOUT,1002) LGNUM,
```

```
1     ASDATE,ASTIME,
2     SAMPID,OPERID,LGFILE,
3     ICYCLE,CYCBUF,ASANS,C240FG,BAD,(237+IJ,WTPCT(IJ),IJ=1,2),
4     WTPCT(5),wtpct(3),WTPCT(6),wtpct(4),ASANS
1002  FORMAT (I4,'.',A9,1X,A5,1X,A8,1X,A8,1X,A13,14,'/',,3A,
1     '1PE10.3,4X,L,2X,L,/,4X,2(' %',I3,'=',OPF7.4),
2     ' %240=',F7.4,' %241=',F7.4,' %242=',F7.4,' %Am=',F7.4,/,
3     4X,'Specific power =',1PE10.3)
```

```
      INTIPTR(LGNUM) = IPTR
```

```
      LGNUM = LGNUM + 1
```

```
240  IF (IPTR.EQ.FPTR) GO TO 400
```

```
      IPTR=IPTR-1
```

```
      IF (IPTR.LT.2) IPTR=MAXASY
```

300 CONTINUE

C If the line printer is open, close it in order to finish printing

400 IF (LOUT.EQ.6) CLOSE (UNIT=LOUT)

C Prompt for sending entries to VAX

GOTO 150

IF (.NOT.YESNO('Send entry to VAX')) GO TO 150

IF (YESNO('All files not yet sent')) GO TO 350

380 TYPE 1014

1014 FORMAT(' '// \$ Type number of entry (RETURN to quit) -> ')

READ (5,1000,END=380,ERR=380) ENTNUM

IF (ENTNUM.EQ.0) GO TO 150

C Check if number entered is in the range of entries listed

IF (ENTNUM .GE. 1 .AND. ENTNUM .LT. LGNUM) GO TO 390

C Error message for number less than one or  
C greater than maximum entries listed

TYPE 1015,ENTNUM

1015 FORMAT('/' '\*\*\*\* ENTRY NUMBER',I4,' WAS NOT LISTED \*\*\*\*'/)

GO TO 380

C Read the entry from the file and send it to the VAX

390 IERR=2

READ (DKLUN'INTPTR(ENTNUM),ERR=90)

1 ASDATE, ASTIME, SAMPID, OPERID, MATTYP, ASTYPE,

2 LGFILE, ICYCLE, NUMCYC, ASANS, ERRBAR, DGFLG,

3 BIASOK, PRECOK, BACKOK, BAD, CONDIR, C240FG,

4 (WTPCT(I), PCT(I),I=1,6)

IF (ASTYPE(1).NE.'A') GO TO 392

CALL COMM('RA') !\*communications - log assay

IF (.NOT.VAXUP) GO TO 150

IF (LVFLAG) GO TO 380

LVFLAG=.TRUE.

IERR=3

WRITE (DKLUN'INTPTR(ENTNUM),ERR=90)

1 ASDATE, ASTIME, SAMPID, OPERID, MATTYP, ASTYPE,

2 LGFILE, ICYCLE, NUMCYC, ASANS, ERRBAR, DGFLG,

3 BIASOK, PRECOK, BACKOK, BAD, CONDIR, C240FG,

4 (WTPCT(I), PCT(I),I=1,6)

GO TO 380

392 TYPE 1030,ASTYPE

1030 FORMAT ('//20X,'\*\*\*\* ENTRY TYPE ',2A1,

1 ' CANNOT BE SENT TO VAX \*\*\*\*'/)

GO TO 380

C All entries not sent to VAX yet will be sent here

350 IPTR=LPTR

355 IERR=2

READ (DKLUN'IPTR,ERR=90)

1 ASDATE, ASTIME, SAMPID, OPERID, MATTYP, ASTYPE,

2 LGFILE, ICYCLE, NUMCYC, ASANS, ERRBAR, DGFLG,

```

3      BIASOK, PRECOK, BACKOK,   BAD, CONDIR, C240FG,
4      (WTPCT(I),                PCT(I),I=1,6)
C      IF (LVFLAG) GO TO 361
C      IF (ASTYPE(1).NE.'A') GO TO 361
C      CALL COMM('RA')           !*communications - log assay
C      IF (VAXUP) GO TO 362
C      IF (AUTO) GO TO 500
C      GO TO 150
C362   LVFLAG=.TRUE.
362   CONTINUE
      IERR=3
      WRITE (DKLUN,IPTR,ERR=90)
1      ASDATE, ASTIME, SAMPID, OPERID, MATTYP, ASTYPE,
2      LGFILE, ICYCLE, NUMCYC, ASANS, ERRBAR, DGFLG,
3      BIASOK, PRECOK, BACKOK,   BAD, CONDIR, C240FG,
4      (WTPCT(I),                PCT(I),I=1,6)
361   IF (IPTR.EQ.FPTR) GO TO 365
      IPTR=IPTR-1
      IF (IPTR.LT.2) IPTR=MAXASY
      GO TO 355
365   IF (AUTO) GO TO 500
C
C   Check if more are to be sent
C
      IF (YESNO('Send more')) GO TO 380
      GO TO 150
C
C   ** Log file handling errors **
C
90    CALL ERRMSG(IERR,ASLOG)
C
C   Clean up and exit
C
500   CLOSE (UNIT=DKLUN,ERR=501)
501   RETURN
      END

```

```

C*****
C  MCBIAS
C      Measurement control check for instrument bias
C          Generic package
C
C      05-Jun-84 1c
C
C      Rev. 07-Jun-84 1c
C      REV. 05-Sep-84 TES
C      REV. 25-Sep-84 CMS
C      Rev. 08-Jan-85 MPK      Move call to VAX comm and add
C                              call to DATTIM so that VAX
C                              gets start time
C
C      Rev. 18-JAN-85 MPK      Add calls to ADATE and AEVAX
C                              remove reference to VAXCOM
C
C      Rev. 22-Jan-85 MPK      Check VAXUP in COMM calls
C
C      Rev. 12-Feb-85 MPK      Redo logical flags
C
C      Rev. 15-Feb-85 MPK      Use YESNO to ask yes/no questions
C
C      Rev. 20-Feb-85 MPK      Changed BIACT=2.58 to BIACT=3.00
C
C      Rev. 07-Mar-85 MPK      Included CVAX common
C
C      Rev. 04-Apr-86 WDR      Modified for Micro VAX II
C
C      Rev. 02-May-90 WDR      Removed COMM and AEVAX for Pu238
C                              instrument
C

```

```

C*****

```

```

C      SUBROUTINE MCBIAS

```

```

C      LOGICAL*1 SUMMRY,YESNO
C      REAL*8 WORD

```

```

C      The following common block files are generic for all instruments

```

```

C      INCLUDE 'CMCONS.CMN'
C      INCLUDE 'CMEAS.CMN'
C      INCLUDE 'CLUNS.CMN'
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CMC.CMN'
C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CDATE.CMN'
C      INCLUDE 'CVAX.CMN'

```

```

C      DATA WORD // BIAS //
C      MB = .TRUE.

```

```

C      Set constants for bias warning limit and action limit
C      HISTSD comes from the parameter file

```

```

C      BIWARN = 1.96 * HISTSD
C      BIACT  = 3.00 * HISTSD
C      RTYP  = 'MC Bias '

```

```

C      Initialize repeat count
C      Assay

```

```

C      SUMMRY = .FALSE.
C      ITIMES = 0
C      ICYCLE = 1
C      CALL DIALOG

```



IF(ABORT) GO TO 500

C

C Announce measurement control bias to VAX

C

CALL ADATE

IF (VAXUP)CALL COMM ('SB') !\*communications - start bias run

CALL NEWNAM(0)

5 CALL ASSAY

IF(ABORT) GO TO 500

VALUE = ASANS ! Instrument specific variable here.

SIGMA = ERRBAR

C

C Initialize bias check pass flag

C Get difference between measured value & reference standard value

C

BIASOK = .TRUE.

DIFF = VALUE - BIREF

DIFF = ABS(DIFF)

C

C Difference less than warning limit

C

IF(DIFF .LT. BIWARN) GO TO 250

C

C Difference equals or exceeds warning limit

C but is less than action limit

C

IF(DIFF .GE. BIACT) GO TO 200

CALL ERRMSG(77,' BIAS WARNING ')

GO TO 250

C

C Difference equals or exceeds action limit

200 BIASOK = .FALSE.

CALL ERRMSG(78,' ')

ITIMES = ITIMES + 1

C

250 CALL PRINT1

WRITE(LOUT,1200)VALUE,BIREF,DIFF

1200 FORMAT(/' Current value = ',F10.4/' Reference value = ',F10.4/

X ' Difference = ',F10.4)

IF(BIASOK .OR. ITIMES .GT. 1) GO TO 300

IF (YESNO('Do you want to repeat bias check'))GO TO 5

C

C Send to VAX

C Write to measurement control log

C IDCAL is to be defined for specific instrument

C

300 CALL AEND

C300 CONTINUE

VAR = BIREF

MCRUNS = 1

IF(.NOT. BIASOK) MCERR = .TRUE.

MTYP(1) = 'M'

MTYP(2) = 'B'

C

IF (VAXUP) CALL COMM('RB') !\*communications - bias results

MCVAX = VAXUP

CALL WTMCLG

C

500 RETURN

END



```

C   have been checked.
C   Also initialize flags for run expired, and run out of limit.
C   A run that is missing from the log will appear to be an expired run.
C
C   LBKG = .TRUE. !Set to TRUE for solid isotopic instrument
C   LPREC = .FALSE.
C   LBIAS = .FALSE.
C   LBKTIM = .FALSE. !Set to FALSE for solid isotopic instrument
C   LPRTIM = .TRUE.
C   LBITIM = .TRUE.
C   LBKERR = .FALSE.
C   LPRERR = .FALSE.
C   LBIERR = .FALSE.

```

```

C
C   Open measurement control log file and read first record to obtain
C   number of records in file.
C

```

```

C   OPEN (UNIT=DKLUN,NAME=MCLOG,ACCESS='DIRECT',TYPE='OLD',
C   1      RECORDSIZE=MCREC,ERR=105)
C   GO TO 110
105  CALL ERRMSG(1,MCLOG)          !** open error
C   GO TO 650
110  READ (DKLUN'1,ERR=115) LPTR,FPTR
C   GO TO 120
115  CALL ERRMSG(2,MCLOG)          !** read error
C   GO TO 650

```

```

C
C   Check that number of records in file is within range
C

```

```

120  IF (LPTR.LT.2 .OR. LPTR.GT.MCMAX) GO TO 125
C   IF (FPTR.GE.2 .AND. FPTR.LE.MCMAX) GO TO 200
125  CALL ERRMSG(5,MCLOG)
C   GO TO 650

```

```

C
C   Convert the current date and time to "julian" days and seconds
C

```

```

200  CALL JULDAT(NOWDAT, NOWTIM, XNDAT, XNTIM)

```

```

C
C   Read through the MC log looking for records of runs that have not
C   been checked yet. Stop when all three run types (bias, precision,
C   and background) have been checked, or when the entire log has been
C   read.
C

```

```

C   IPTR = LPTR
205  READ (DKLUN'IPTR,ERR=115)
C   1      ASDATE, ASTIME, SAMPID, OPERID, MTYP, MCRUNS,
C   2      VALUE, SIGMA, VAR, ARCHV1, ARCHV2, DGFLG,
C   3      BIASOK, PRECOK, BACKOK, MCERR, MCVAX

```

```

C   IF (LBIAS) GO TO 210          !bias already checked
C   IF (MTYP(2).NE.'B') GO TO 210 !not a bias record

```

```

C
C   Check a bias record. Save the error flag, and copy the time limit
C   for a bias run. Also set the flag to remember that a bias record
C   was checked.
C

```

```

C   LBIAS=.TRUE.
C   LBIERR=MDEF LG
C   TIMLIM=BITIME
C   GO TO 230

```

```

C
210  IF (LPREC) GO TO 220                !precision already checked
      IF (MTYP(2).NE.'P') GO TO 220    !not a precision record
C
C  Check a precision record.
      LPREC=.TRUE.
      LPRERR=MCEFLG
      TIMLIM=PRTIME
      GO TO 230
C
220  IF (LBKG) GO TO 240                !background already checked
      IF (MTYP(2).NE.'G') GO TO 240    !not a background record
C
C  Check a background record
      LBKG=.TRUE.
      LBKERR=MCEFLG
      TIMLIM=BKTIME
C
C  Common code to do measurement control checks.
C  Find the interval (in days) between the current time and the
C  time of the run in the log. Compare this to the allowed time
C  interval, and set a flag for the appropriate run type to
C  indicate whether the run is current or expired.
C
230  CALL JULDAT(ASDATE, ASTIME, XDAT, XTIM)
      DIF = (XNDAT+XNTIM/86400.) - (XDAT+XTIM/86400.)
      LTIME=DIF.GT.TIMLIM
      IF (MTYP(2).EQ.'B') LBITIM=LTIME
      IF (MTYP(2).EQ.'P') LPRTIM=LTIME
      IF (MTYP(2).EQ.'G') LBKTIM=LTIME
C
C  If all three MC types have been checked, stop looking
C
240  IF (LBKG.AND.LPREC.AND.LBIAS) GO TO 300
      IF (IPTR.EQ.FPTR) GO TO 300
      IPTR = IPTR -1
      IF (IPTR.LT.2) IPTR=MCMAX
      GO TO 205
C
C  Check flags for runs expired or out-of-limit, and generate
C  appropriate messages
C
300  IF (LBITIM) TYPE 1000
- 1000  FORMAT (' **** MEASUREMENT CONTROL BIAS RUN',
1        ' SHOULD BE MADE ****')
      IF (LPRTIM) TYPE 1001
- 1001  FORMAT (' **** MEASUREMENT CONTROL PRECISION RUN',
1        ' SHOULD BE MADE ****')
      IF (LBKTIM) TYPE 1002
1002  FORMAT (' **** MEASUREMENT CONTROL BACKGROUND RUN',
1        ' SHOULD BE MADE ****')
      IF (LBIERR) TYPE 1010
1010  FORMAT (' **** LAST MEASUREMENT CONTROL BIAS RUN',
1        ' WAS OUT OF LIMIT ****')
      IF (LPRERR) TYPE 1011
- 1011  FORMAT (' **** LAST MEASUREMENT CONTROL PRECISION RUN',
1        ' WAS OUT OF LIMIT ****')
      IF (LBKERR) TYPE 1012

```

```
1012  FORMAT (' **** LAST MEASUREMENT CONTROL BACKGROUND RUN',  
1      ' WAS OUT OF LIMIT ****')
```

C

C Set global flags to indicate whether meas control is OK

C

```
      BIASOK=.NOT.(LBITIM.OR.LBIERR)
```

```
      PRECDK=.NOT.(LPRTIM.OR.LPRERR)
```

```
      BACKOK=.NOT.(LBKTIM.OR.LBKERR)
```

C

C If a run out of limit was found, get the operator's permission

C before continuing

C

```
      IF ( .NOT.(LBIERR.OR.LPRERR.OR.LBKERR) ) GO TO 650
```

```
      IF (YESNO('Do you want to continue')) GO TO 650
```

C

C Abort exit

C

```
640  ABORT = .TRUE.
```

C

C Non-abort exit

C

```
650  CLOSE (UNIT=DKLUN)
```

```
      RETURN
```

```
      END
```

```
C
C*****
C
C MCHEAD
C Outputs header information for the measurement control
C log list option
C
C CMS 14-MAY-84
C
C REV. 28-Aug-84 MJC
C REV. 04-Apr-86 WDR Modified for Micro VAX II
C*****
```

```
C SUBROUTINE MCHEAD
```

```
C The following common block files are generic for all instruments
```

```
C INCLUDE 'CLUNS.CMN'
C INCLUDE 'CLOG.CMN'
C INCLUDE 'CDATE.CMN'
C INCLUDE 'CMC.CMN'
```

```
C CALL DATTIM(1)
```

```
C WRITE(LOUT,1099)NOWDAT,NOWTIM
```

```
1099 FORMAT('1',10X,'Measurement control log as of ',A,
X 2X,A//1X,'Ent Date Time Operator Ty ',
X ' Result Fractional Standard No VX Er '
X /44X,'Error Value Run Fg Fg'/)
```

```
C RETURN
```

```
C END
```

C\*\*\*\*\*

C MCPREC

C Measurement control check for precision

C 06-Jun-84 lc

C Rev. 07-Jun-84 lc

C REV. 06-Sep-84 TES

C REV. 25-Sep-84 CMS

C Rev. 08-Jan-85 MPK

Move call to VAX comm and add  
call to DATTIM so that VAX  
gets start time

C Rev. 18-Jan-85 MPK

Add calls to ADATE and AEVAX  
Remove reference to to VAXCOM

C Rev. 22-Jan-85 MPK

Check VAXUP in COMM calls

C Rev. 14-FEB-85 WDR

Limits for 5 and 15 runs were reversed

C Rev. 19-Feb-85 MPK

Correct statistics

C Rev. 20-Feb-85 MPK

Use YESNO to ask yes/no question

C Rev. 21-Feb-85 MPK

Declared XVAL,XSIG,VALUE,XMEAN as REAL\*8

C Rev. 07-Mar-85 MPK

Declare YESNO as Logical\*1

C Rev. 18-Apr-85 MPK

Include CVAX common

C Rev. 18-Apr-85 MPK

Remove REAL\*8 declaration of VALUE  
(causing misalignment of CRESLT common)

C Rev. 04-Apr-86 WDR

Modified for Micro VAX II

C Rev. 02-May-90 WDR

Removed COMM and AEVAX for Pu238  
instrument

C Rev. 03-May-90 WDR

Add call to AEND

C Rev. 03-May-90 WDR

C\*\*\*\*\*

C SUBROUTINE MCPREC

C LOGICAL \* 1 SUMMRY,YESNO

C REAL\*8 XVAL(25), XSIG(25), XMEAN

C REAL \* 8 WORD

C REAL \* 4 PRWARN(2), PRACT(2)

C The following common block files are generic for all instruments

C INCLUDE 'CMCONS.CMN'

C INCLUDE 'CMEAS.CMN'

C INCLUDE 'CLUNS.CMN'

C INCLUDE 'CONTRL.CMN'

C INCLUDE 'CASSAY.CMN'

C INCLUDE 'CDATE.CMN'

C INCLUDE 'CMC.CMN'

C INCLUDE 'CVAX.CMN'

C DATA WORD /'PRECISIO'/

C MP = .TRUE.

C RTYP = 'MC PRECN'

C Initialize - set up values for the precision warning limits  
C and action limits - these values are dependent on the NRUNS  
C selected in the parameter file. First make sure NRUNS is  
C equal to 5 or 15.

C IF (NRUNS.NE. 5 .AND. NRUNS.NE. 15)NRUNS = 5

C IF(NRUNS .EQ. 15) GO TO 150

```

C
C Set limits for 5 runs
C
      PRWARN(1) = 0.12
      PRWARN(2) = 2.79
      PRACT(1) = 0.05
      PRACT(2) = 3.72
      GO TO 155
C
C NRUNS equals 15
C
150  PRWARN(1) = 0.40
      PRWARN(2) = 1.87
      PRACT(1) = 0.29
      PRACT(2) = 2.24
C
155  ITIMES = 0
5    XMEAN = 0.0
      SIGSQ = 0.0
      SUMMRY = .FALSE.
C
C Dialog with operator
C
      CALL DIALOG
      IF(ABORT) GO TO 500
      NUMCYC = NRUNS           !For printout in Header summary
C
C
C Assay, print single assay results
C Accumulate results
C
      ICYCLE = 0
      DO 10 I=1,NRUNS
      ICYCLE = ICYCLE + 1
      CALL ADATE
      CALL ASSAY
      IF(ABORT) GO TO 500
      CALL NEWNAM(0)
      VALUE = ASANS
      SIGMA = ERRBAR
      CALL PRINT1
      XVAL(I) = VALUE
      XSIG(I) = SIGMA
      XMEAN = XMEAN + VALUE
      ABSIG = SIGMA * VALUE
      SIGSQ = SIGSQ + ABSIG*ABSIG
10  CONTINUE
C
C Average results
C Square mean sigma
C
      XMEAN = XMEAN/NRUNS
      SGMEAN = SIGSQ/NRUNS
C
C Accumulate variance
C
      SUM2 = 0.0
      DO 15 I=1,NRUNS
      VAR = (XVAL(I)-XMEAN) * (XVAL(I) -XMEAN)
      SUM2 = SUM2 + VAR

```



15 CONTINUE

```
C
C Compute square of standard deviation
C Compute chi square (standard deviation **2 / mean sigma **2)
C
      STDEV2 = SUM2/(NRUNS-1)
      CHISQ = STDEV2/SGMEAN
C
C Set up variables to write to MC log file.
C
      VALUE = CHISQ
      SIGMA = STDEV2
      VAR = SGMEAN
      MCRUNS = NRUNS
C
C Check for chi square limits
C
      PRECOK = .TRUE.
      IF(CHISQ .GT. PRWARN(1) .AND. CHISQ .LT. PRWARN(2))
X          GO TO 250
C
C Warning limits exceeded
C
      IF(CHISQ .LT. PRACT(1) .OR. CHISQ .GT. PRACT(2))
X          GO TO 200
      CALL ERRMSG(75,' ')
      GO TO 250
C
C Action limits exceeded
C
700      PRECOK = .FALSE.
      CALL ERRMSG(76,' ')
C
C Print results
C
250      SUMMRY = .TRUE.
      CALL HEADER(SUMMRY)
      WRITE(LOUT,1500)
1500      FORMAT(//3X,'Cycle',T15,'Result',T30,'Sigma',T50,'Abs. Sigma'/)
      DO 50 I=1,NRUNS
          WRITE(LOUT,1501)I,XVAL(I),XSIG(I),XSIG(I)*XVAL(I)
1501      FORMAT(5X,I3,T13,F10.4,T28,F8.4,T48,F8.4)
          50 CONTINUE
      STDEV = SQRT(STDEV2)
      WRITE(LOUT,1502)XMEAN,SGMEAN,STDEV,CHISQ
1502      FORMAT(// Average result = ',T25,F10.4/' Mean Sq sigma = ',T25,
X          F10.4/' Standard deviation = ',T25,F10.4/' Chisquare = ',
X          T25,F10.4)
      ITIMES = ITIMES + 1
      IF(PRECOK .OR. ITIMES .GT. 1) GO TO 300
      IF (YESNO('Do you want to repeat precision check?')) GO TO 5
C
C Send results to VAX
C Write to measurement control log
C
300      CALL AEND
      IF(.NOT. PRECOK) MCERR = .TRUE.
      MTYP(1) = 'M'
      MTYP(2) = 'P'
C
      IF(VAXUP) CALL COMM ('RP') !communications - precision results
```

MCVAX = VAXUP

CALL WTMCLG

500 RETURN

END

C\*\*\*\*\*

C MENU

C

C Searches the table 'MENU.TXT' to locate a match  
C with the option entered by the operator. If no  
C match is found a flag is set and the routine  
L returns to the Executive.

C

C Calling arguments and format:

C OPTIN (Option entered by operator) A4

C MNUM (Number of entry in table) I2

C NOPT (Logical flag set if no match) L

C

C

C CMS 25-APR-84

C

C\*\*\*\*\*

C

C

SUBROUTINE MENU(OPTIN,MNUM,NOPT)

LOGICAL \* 1 NOPT

CHARACTER OPT(25)\*4

CHARACTER OPTIN\*4,OPTTBL\*4

INTEGER OPTNUM,MENPRM,OPTLEN

C

C The following common block files are generic for all instruments

C

INCLUDE 'CFILES.CMN'

INCLUDE 'CLUNS.CMN'

INCLUDE 'CONTRL.CMN'

C

DATA MENPRM /25/, OPTLEN /4/

NOPT = .FALSE.

C

C Open the file 'MENU.TXT' for reading

C

OPEN(UNIT=DKLUN,NAME='MENU.TXT',TYPE='OLD',ERR=100).

GOTO 20

100 CALL ERRMSG(1,'MENU')

ABORT = .TRUE.

GOTO 70

C

C Read the menu file and store data in the internal menu

C

20 READ(DKLUN,10,END=15,ERR=200) OPTNUM,OPTTBL

10 FORMAT(I3,A4)

GOTO 505

200 CALL ERRMSG(2,'MENU')

GOTO 70

505 IF (OPTNUM .GT. MENPRM .OR. OPTNUM .LT. 1) GOTO 20

OPT(OPTNUM) = OPTTBL

GOTO 20

15 CLOSE(UNIT=DKLUN)

C

C Look for the option match in the table

C

DO 60 I = 1, MENPRM

IF (OPTIN.NE.OPT(I)) GOTO 60

MNUM = I

GOTO 70

60 CONTINUE

C  
C  
C

Unrecognized option - flag set and return

TYPE \*, '\*\*\*\* OPTION NOT RECOGNIZED - TYPE H TO SEE MENU \*\*\*\*'  
NOPT = .TRUE.

70 RETURN  
END

SUBROUTINE MLR (MM, M, N, A, Y, X, B, R, W, V)

INTEGER\*2 MM, M, N

REAL A(MM,N), Y(M), X(N), B(MM,N), R(M), W(N), V(N)

INTEGER MLRI

COMMON /MLRCC/ MLRI

-----  
MLR (MINIMIZE LINEAR RESIDUALS) COMPUTES A SOLUTION  
TO THE (POSSIBLY) OVERDETERMINED LINEAR SYSTEM

$$A X = Y$$

IN THE SENSE THAT THE LENGTH OF THE RESIDUAL VECTOR

$$R(X) = Y - A X$$

IS MINIMIZED, USING THE METHOD OF HOUSEHOLDER REFLECTIONS  
WITH ITERATIVE IMPROVEMENT.

-----  
CALLING SEQUENCE:

COMMON /MLRCC/ MLRI

MLRI = <INTEGER VALUE>

CALL MLR (MM, M, N, A, Y, X, B, R, W, V)

WHERE:

MLRI TELLS WHICH "ENTRY" TO USE:

MLRI = 1: MLR

MLRI = 2: MLRA USE MLR AGAIN WITH SAME MATRIX A, NEW  
NEW RIGHT-HAND SIDE Y.

MLRI = 3: MLRE COMPUTE DIAGONAL OF INVERSE OF NORMAL  
MATRIX, A'A.

MLRI = 4: MLRS COMPUTE FULL INVERSE OF NORMAL MATRIX.

MLRI IS SET TO 1 IN A DATA STATEMENT IN MLR .

MM IS THE NUMBER OF ROWS SPECIFIED IN THE DIMENSION STATEMENT  
FOR A.

M IS THE NUMBER OF EQUATIONS (ROWS).

N IS THE NUMBER OF UNKNOWNNS (COLUMNS).

A IS THE TWO-DIMENSIONAL ARRAY OF COEFFICIENTS.

Y IS THE ONE-DIMENSIONAL ARRAY OF CONSTANTS (M WORDS).

X IS A ONE-DIMENSIONAL ARRAY FOR THE UNKNOWNNS (N WORDS).

B IS A TEMPORARY STORAGE ARRAY OF MM\*N WORDS.

IN CASE MLRI=4, B IS REPLACED BY THE FULL INVERSE OF THE  
NORMAL MATRIX. B IS ASSUMED TO BE DIMENSIONED EXACTLY AS A.

R IS A ONE-DIMENSIONAL ARRAY FOR THE RESIDUALS (M WORDS).

IN CASE MLRI=3, R IS REPLACED BY THE DIAGONAL ELEMENTS  
OF THE INVERSE OF THE NORMAL MATRIX (N WORDS).

W IS A TEMPORARY STORAGE ARRAY OF N WORDS.

V IS A TEMPORARY STORAGE ARRAY OF N WORDS.

INPUT ARGUMENTS: MLRI, MM, M, N, A, Y.

WORKING ARGUMENTS: B, R, W, V.

OUTPUT ARGUMENTS: X, R.

STANDARD FORTRAN FUNCTIONS REFERENCED: DBLE, SIGN, SQRT

OTHER ROUTINES CALLED: MLRERR

---

CAUTION:

1. A CALL WITH MLRI.GT.1 ASSUMES THAT:
  - A. MLR HAS PREVIOUSLY BEEN CALLED WITH MLRI=1, AND
  - B. THE CONTENTS OF ARRAYS A, B, W, V HAVE NOT BEEN  
DISTURBED SINCE THE PREVIOUS CALL.  
(THERE IS NO WAY TO CHECK FOR THIS.)
2. MLR MAY NOT BE CALLED AGAIN WITH MLRI.GT.1  
IF THE MLRI=4 OPTION IS USED.

---

STANDARD FORTRAN VERSION WRITTEN BY F. N. FRITSCH (LLL).  
DATE LAST CHANGED: 24 MAY 1978.

---

NOTE.. THE USER MAY MONITOR THE ERROR FLAG JFLAG BY INCLUDING  
THE FOLLOWING COMMON BLOCK DECLARATION IN THE CALLING  
PROGRAM. (JFLAG.EQ.0 FOR NORMAL RETURN. )  
(JFLAG.NE.0 FOR ERROR RETURN. )

```
INTEGER JFLAG  
COMMON /MLRERF/ JFLAG
```

.....DECLARE LOCAL VARIABLES.

```
INTEGER I, IN, IP, J, JM, K  
REAL RE, REOLD, TEMP, YE  
DOUBLE PRECISION RD, XD
```

```
DATA MLRI /1/
```

.....CHECK FOR INVALID INPUTS.

```
JFLAG = -1  
IF (MLRI.LT.1 .OR. MLRI.GT.4) GO TO 510  
JFLAG = 0  
IF (M.GT.MM) JFLAG = 55  
IF (N.GT.M) JFLAG = 56  
IF (N.LT.1) JFLAG = 57  
IF (JFLAG.NE.0) GO TO 510
```

```
GO TO (1, 100, 300, 400), MLRI
```

```

C-----
C
C...CASE MLRI = 1 ( MLR:
C ( NORMAL USAGE -- COMPLETE SOLUTION.
1 CONTINUE
C
C INITIALIZE B TO A .
DO 20 I=1,M
  DO 10 J=1,N
    B(I,J) = A(I,J)
10 CONTINUE
20 CONTINUE
C
C.....TRIANGULARIZATION PHASE.
C
DO 70 I=1,N
  TEMP = 0.
  DO 30 K=I,M
    TEMP = TEMP + B(K,I)**2
30 CONTINUE
C TEST FOR SINGULARITY.
IF (TEMP .EQ. 0.) GO TO 500
TEMP = SQRT(TEMP)
W(I) = SIGN (TEMP, B(I,I))
B(I,I) = B(I,I) + W(I)
V(I) = - B(I,I) * W(I)
IF (I .EQ. N) GO TO 80
IP = I+1
DO 60 J=IP,N
  TEMP = 0.
  DO 40 K=I,M
    TEMP = TEMP + B(K,I)*B(K,J)
40 CONTINUE
  TEMP = TEMP/V(I)
  DO 50 K=I,M
    B(K,J) = B(K,J) + TEMP*B(K,I)
50 CONTINUE
60 CONTINUE
70 CONTINUE
C
80 CONTINUE
C-----
C
C...CASE MLRI = 2 ( MLRA:
C ( RE-USE MLR WITH SAME A , DIFFERENT Y .
100 CONTINUE
C
C INITIALIZE X TO ZERO.
DO 110 I=1,N
  X(I) = 0.
110 CONTINUE
C INITIALIZE R TO Y AND COMPUTE LENGTH OF Y .
YE = 0.
DO 120 I=1,M
  R(I) = Y(I)
  YE = YE + Y(I)**2
120 CONTINUE
C ALREADY DONE IF INPUT VECTOR IS ZERO.
IF (YE .EQ. 0.) RETURN

```

REOLD = YE

C  
C.....FORWARD SOLUTION PHASE.

C  
130 CONTINUE  
COMPUTE VECTOR Z .  
DO 160 I=1,N  
TEMP = 0.  
DO 140 K=I,M  
TEMP = TEMP + B(K,I)\*R(K)  
140 CONTINUE  
TEMP = TEMP/V(I)  
DO 150 K=I,M  
R(K) = R(K) + TEMP\*B(K,I)  
150 CONTINUE  
160 CONTINUE

C  
C.....BACK-SUBSTITUTION PHASE.

C  
COMPUTE CORRECTION, USING R FOR TEMPORARY STORAGE.  
DO 190 I = 1, N  
IN RUNS FROM N DOWN TO 1.  
IN = N - I + 1  
R(IN) = R(IN) / W(IN)  
IF (IN .EQ. 1) GO TO 190  
DO 180 K = 2, IN  
J RUNS FROM IN-1 DOWN TO 1.  
J = IN - K + 1  
R(J) = R(J) + R(IN)\*B(J,IN)  
180 CONTINUE  
190 CONTINUE

78-05-24

C  
C ADD CORRECTION TO X : (DOUBLE PRECISION)

DO 200 I=1,N  
XD = X(I)  
RD = R(I)  
XD = XD - RD  
X(I) = XD  
200 CONTINUE

C  
C.....ITERATIVE IMPROVEMENT PHASE.

C  
DOUBLE PRECISION RESIDUAL CALCULATION.  
DO 230 I=1,M  
RD = Y(I)  
DO 220 J=1,N  
RD = RD - DBLE(A(I,J))\*DBLE(X(J))  
220 CONTINUE  
R(I) = RD  
230 CONTINUE  
C COMPUTE LENGTH OF NEW RESIDUAL VECTOR.  
RE = 0.  
DO 240 I=1,M  
RE = RE + R(I)\*\*2  
240 CONTINUE  
C DONE WITH ITERATIVE IMPROVEMENT IF LENGTH OF R NOT DECREASED.  
IF (RE .GE. REOLD) GO TO 270  
REOLD = RE  
GO TO 130



```
270 CONTINUE
C TEST FOR NUMERICAL SINGULARITY.
  IF (RE .LT. YE) RETURN
  JFLAG = 58
  GO TO 510
```

```
-----
C...CASE MLRI = 3 ( MLRE:
C ( COMPUTE DIAGONAL OF (A'A)-INVERSE = (T'T)-INVERSE
```

```
300 CONTINUE
C
  DO 350 I = 1, N
C.....COMPUTE I-TH ROW OF T-INVERSE AND STORE IN R(I), R(I+1),
C      ..., R(N).
    R(I) = -1. / W(I)
    IF (N .LE. I) GO TO 330
    IP = I+1
    DO 320 J = IP, N
      TEMP = 0.
      JM = J-1
      DO 310 K = I, JM
        TEMP = TEMP + R(K)*B(K,J)
310      CONTINUE
      R(J) = TEMP / W(J)
320      CONTINUE
C.....COMPUTE I-TH DIAGONAL ELEMENT OF (T-INVERSE) X (T-INVERSE)'
C      AND STORE IN R(I).
330      CONTINUE
      TEMP = 0.
      DO 340 J = 1, N
        TEMP = TEMP + R(J)**2
340      CONTINUE
      R(I) = TEMP
350      CONTINUE
C
    RETURN
```

```
-----
C...CASE MLRI = 4 ( MLRS:
C ( COMPUTE FULL INVERSE OF A'A = T'T .
```

```
400 CONTINUE
C
C.....FORM (T-INVERSE) IN B .
  DO 430 I = 1, N
    B(I,I) = -1. / W(I)
    IF (N .LE. I) GO TO 440
    IP = I+1
    DO 420 J = IP, N
      TEMP = 0.
      JM = J-1
      DO 410 K = I, JM
        TEMP = TEMP + B(I,K)*B(K,J)
410      CONTINUE
      B(I,J) = TEMP / W(J)
420      CONTINUE
430      CONTINUE
C.....COMPUTE (T-INVERSE) X (T-INVERSE)' .
```

```

440 CONTINUE
      DO 470 I = 1, N
        DO 460 J = I, N
          TEMP = 0.
          DO 450 K = J, N
            TEMP = TEMP + B(I,K)*B(J,K)
          450 CONTINUE
          B(I,J) = TEMP
          B(J,I) = TEMP
        460 CONTINUE
      470 CONTINUE

```

```

      RETURN

```

---

```

C...ERROR RETURNS.

```

```

500 CONTINUE
      JFLAG = 59

```

```

510 CONTINUE
      FOLLOWING CALL ASSUMES THAT MLRERR IGNORES RE AND YE
      UNLESS JFLAG = 58 .

```

```

      CALL MLRERR (RE, YE)

```

```

      RETURN

```

---

```

      END
      SUBROUTINE MLRERR (RE, YE)
      REAL RE, YE

```

```

      INTEGER JFLAG
      COMMON /MLRERF/ JFLAG

```

```

      INTEGER OUB
      COMMON /MLRUNIT/ OUB

```

---

```

      M L R   E R R O R   R O U T I N E

```

```

      THIS ROUTINE IS CALLED IF, AND ONLY IF, AN ERROR IS ENCOUNTERED
      IN MLR .

```

---

```

      CALLING SEQUENCE:

```

```

          COMMON /MLRERF/ JFLAG
          JFLAG = <NONZERO INTEGER VALUE>
          CALL MLRERR (RE, YE)

```

```

      WHERE:

```

```

      JFLAG INDICATES WHICH ERROR HAS OCCURRED. ALLOWABLE VALUES ARE:
      JFLAG = -1: ILLEGAL VALUE OF MLRI .

```

```
C          JFLAG = 55:  M > MM .
C          JFLAG = 56:  N > M .
C          JFLAG = 57:  N < 1 .
C          JFLAG = 58:  SINGULAR SYSTEM DETECTED IN MLRA .
C          [LENGTH(R) .GE. LENGTH(Y)].
C          JFLAG = 59:  SINGULAR MATRIX DETECTED DURING TRIANGULAR-
L          IZATION PHASE.
C          [AN ELEMENT OF W IS ZERO.]
```

```
C          RE, YE ARE THE SQUARED LENGTHS OF R AND Y , RESPECTIVELY,
C          AND ARE IGNORED UNLESS JFLAG=58 .
```

```
C          OUB DETERMINES THE DISPOSITION OF ERROR MESSAGES:
C          OUB < 0:  NO PRINTED MESSAGES.
C          OUB = 0:  TELETYPE OUTPUT (DEFAULT).
C          OUB > 0:  OUTPUT TO LOGICAL UNIT OUB .
```

```
-----
C          WRITTEN BY F. N. FRITSCH (LLL).
C          DATE LAST CHANGED:  10 OCTOBER 1975.
C          -----
```

```
C          INTEGER MLRI
C          COMMON /MLRCC/ MLRI
C          C.....DECLARE LOCAL VARIABLES.
C          INTEGER LOUT, JFL
C          DATA OUB /0/
C          C.....SET UP LOUT FOR 11/70.
C          LOUT = OUB
C          TELETYPE OUTPUT IF OUB=0 .
C          IF (OUB .EQ. 0) LOUT = 5
C          IF (LOUT .GT. 0) WRITE (LOUT, 2000) JFLAG
C          C.....HANDLE SPECIAL CASE JFLAG = -1 .
C          IF (JFLAG .NE. -1) GO TO 1000
C          IF (LOUT .GT. 0) WRITE (LOUT, 2001) MLRI
C          RETURN
C          C.....CHECK FOR ILLEGAL VALUE OF JFLAG .
C          1000 CONTINUE
C          IF (JFLAG.LT.55 .OR. JFLAG.GT.59) GO TO 1060
C          IF (LOUT.LE.0) RETURN
C          JFL = JFLAG - 54
C          GO TO (1010, 1020, 1030, 1040, 1050), JFL
C          C.....CASE JFL = 1 (JFLAG = 55).
C          1010 CONTINUE
C          WRITE (LOUT, 2005)
C          WRITE (LOUT, 2010)
C          RETURN
L          C.....CASE JFL = 2 (JFLAG = 56).
C          1020 CONTINUE
```

```
WRITE (LOUT, 2005)
WRITE (LOUT, 2020)
RETURN
```

```
C
C.....CASE JFL = 3 (JFLAG = 57).
030 CONTINUE
WRITE (LOUT, 2005)
WRITE (LOUT, 2030)
RETURN
```

```
C
C.....CASE JFL = 4 (JFLAG = 58).
1040 CONTINUE
WRITE (LOUT, 2040) RE, YE
RETURN
```

```
C
C.....CASE JFL = 5 (JFLAG = 59).
1050 CONTINUE
WRITE (LOUT, 2050)
RETURN
```

```
C.....ENDCASE (JFL)
```

```
C
C.....ILLEGAL JFLAG VALUE. (LOGICALLY, THIS CASE CANNOT OCCUR.)
1060 CONTINUE
IF (LOUT .LT. 0) LOUT = 59
WRITE (LOUT, 2060)
CALL EXIT
```

```
C
C-----
```

```
C
C.....FORMATS.
```

```
000 FORMAT (/24H *** MLR ERROR NUMBER ,I4,6H ***)
2001 FORMAT (6X,23HILLEGAL VALUE OF MLRI: ,I5)
2005 FORMAT (6X,38HSPECIFICATION ERROR IN SUBROUTINE MLR:)
2010 FORMAT (8X,36H M GREATER THAN MM - STORAGE OVERLAP )
2020 FORMAT (8X,48H N GREATER THAN M - MORE UNKNOWN THAN EQUATIONS )
2030 FORMAT (8X,14H N LESS THAN 1 )
2040 FORMAT (6X,27HMATRIX IS SINGULAR IN MLR: /
* 8X,49H SUM OF SQUARES OF RESIDUALS R(I) GREATER THAN OR,
* 42H EQUAL TO SUM OF SQUARES OF VARIABLES Y(I) /
* 8X,18H SUM OF R(I)**2 = ,E15.8,
* 5X,18H SUM OF Y(I)**2 = ,E15.8)
2050 FORMAT (6X,27HMATRIX IS SINGULAR IN MLR: /
* 8X,24H AN ELEMENT OF W IS ZERO )
2060 FORMAT (6X,38HILLEGAL ERROR NUMBER -- JOB TERMINATED)
```

```
C
END
```

```

C*****
C
C NEWNMS      (NEWNAM subroutine for PU-238 isotopic instrument)
C
C      Create a new filename for an autocycle run
C      NEWNMS is a special adaptation of NEWNAM for the solid
C      isotopic instrument.
C
C      Calling arguments:
C      OPTION = 0  create the name with .INS as the extension
C
C      If autocycle mode:
C
C      OPTION = 0  create the name with .IO1 where 'I' is
C                  the shortened extension for the instrument
C      OPTION = 1  update the name to .IO2, .IO3, etc.
C
C      If user inputs the extension, assume that the extension
C      is correct for the first file and do not change it
C
C
C      SSJ 19-DEC-79
C
C      REV. 29-Jun-82  SSJ
C      REV. 07-Sep-84  CMS
C      REV. 21-NOV-84  WDR      Modified for solid isotopic instrument
C      REV. 14-FEB-85  WDR      Check for null file name, and
C                               don't touch it if it is.
C
C      REV. 04-APR-86  WDR      Modified for Micro VAX II
C      REV. 13-JUL-87  RDP      Modified for uVAX II
C                               Convert BYTE strings to CHARACTER strings
C*****
C
C      SUBROUTINE NEWNAM(OPTION)
C
C      CHARACTER ACHAR*1
C      INTEGER*2 OPTION
C
C      The following common block files are generic for all instruments
C
C      INCLUDE 'CFILES.CMN'
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CCAL.CMN'
C
C      If the file name is null, bug out now
C
C      IF (FILNAM.EQ.' ') RETURN
C
C      Check whether option is 0 (create initial extension)
C      or 1 (increment the extension)
C
C      IF(OPTION .EQ. 1) GO TO 200
C
C      Rules:
C      The first blank or dot ends the filename proper
C      - only six characters may be in the filename.
C      If a colon appears in letter 3 or 4, six more letters
C      may appear in the name.

```

```

C      All special characters will be converted to a '*'
C
C      Initialize letter pointer K
C      Initialize count of letters in filename, J, not counting device
C
      K = 0
      J = 0
C
C      Look for blank, dot, or colon
C      Count up to 6 characters after the colon
C      Colon must be in letters 3 or 4
C
10      K = K + 1
      ACHAR = FILNAM(K:K)
      IF(ACHAR .EQ. '.') GO TO 80
      IF(J .EQ. 8 .OR. ACHAR .EQ. ':') GO TO 100
      IF(ACHAR .NE. ':') GO TO 20
C
C      Found:
C      Make sure it is the 3rd or 4th character and not repeated
C
      IF(J .NE. 0 .AND. (K .EQ. 3 .OR. K .EQ. 4)) GO TO 15
      GO TO 20
15      J = 0
      GO TO 10
C
C      Check for alpha-numeric characters
C      Replace anything bad with *
C
20      IF((ACHAR .GE. 'A' .AND. ACHAR .LE. 'Z') .OR.
X          (ACHAR .GE. '0' .AND. ACHAR .LE. '9') )
X          GO TO 25
      FILNAM(K:K) = '*'
25      J = J + 1
      GO TO 10
C
C      User typed the extension
C
80      GO TO 500
C
C      Insert dot
C      Insert extension 'INS' if not autocycle
C      If autocycle, insert extension 'IOO'
C
100     IF(K .GT. 11) STOP 'Filename problem'
      FILNAM(K:K) = '.'
      FILNAM(K+1:K+1) = 'W'
      IF (AUTOCY) GO TO 300
      IF (CAL) GO TO 300
      FILNAM(K+2:K+2) = 'H'
      FILNAM(K+3:K+3) = 'B'
      GO TO 310
300     FILNAM(K+2:K+2) = 'O'
      FILNAM(K+3:K+3) = 'I'
C
      . Insert blanks at end
L
310     K = K + 4
150     FILNAM(K:K) = ' '

```

```
K = K+1
IF(K .LE. 15) GO TO 150
GO TO 500
```

```
C
C Update the previous filename
C First search for the dot
```

```
L
200 K = 0
210 K = K + 1
IF(K .GT. 15) STOP 'No dot in filename'
IF (FILNAM(K:K) .EQ. '.') GO TO 220
GOTO 210
```

```
C
C Update the last digit
C If it goes beyond 9 (:), update the previous digit
```

```
C !Reset character after dot
C 220 FILNAM(K+1:K+1)= CHAR(ICHAR(FILNAM(K+1:K+1)) - 1)
220 K = K + 3 !It was bumped for ADC #2 data
225 IF (FILNAM(K:K) .EQ. '.') GOTO 500
FILNAM(K:K) = CHAR(ICHAR(FILNAM(K:K)) + 1)
IF (FILNAM(K:K) .NE. ':') GOTO 500
FILNAM(K:K) = '0'
K = K - 1
GO TO 225
```

```
C
500 RETURN
END
```

SUBROUTINE OUT

C  
CC  
C Sub to accumulate results of peak fitting, calculate photon  
C emission rate of source, tabulate data & results, and store them;  
if requested.  
L Last edit for VAX 13-Oct-87  
C Modified for Pu238 instrument WDR 1-Sep-89  
C  
CC

INCLUDE 'GRPLGCM.PMS'  
INCLUDE 'GRPGMS.PMS'  
C The following common is required for the instrument

INCLUDE 'CONTRL.CMN'  
INCLUDE 'CLUNS.CMN'

C  
DIMENSION YSUM(MXNP),RES(MXDP)  
EQUIVALENCE (KK,NWRDS(2))  
EQUIVALENCE (YSUM,SUM),(RES,DELY)

DIMENSION QFT(100),SDOF(MXNP),TCTS(MXNP),ENDM(MXNP)  
DIMENSION QFTL(MXPKS)  
EQUIVALENCE (WM(100),ENDM)  
BYTE BKSMBL,EFMBL

CHARACTER\*8 TBKGD  
CHARACTER\*22 TSLPS

CHARACTER\*6 TWID  
CHARACTER\*37 TINT  
CHARACTER\*13 TNGY  
CHARACTER\*66 TOUT  
CHARACTER\*9 ALPJ,STAMP,STSLP,TLAMP,TLSLP

EQUIVALENCE (TWID,TOUT(1:6)),(TOUT(10:22),TNGY)  
EQUIVALENCE (TOUT(29:),TINT)

C  
C DATA LP/6/  
C DATA LP/-1/ !For VAX  
C DATA EFMBL/' '/

IF(RUNFLG(4) .EQ. 0) EFMBL = '\*' !Assuming det. efficiency = 1

C  
C IOFLG(2) is set when last grouping done (preparing summary)  
C IOFLG(3) is set to inhibit line printer output  
C INTFLG is set for single peaks that were simply integrated  
C

IF(IOFLG(2) .NE. 0)GO TO 300

IOFLG(3) = 1 ! NO LINE PRINTER  
IF(IOFLG(3) .NE. 0)GO TO 30

ALPJ = 'FIXED'  
STAMP = 'FIXED'  
STSLP = 'FIXED'  
TLAMP = 'FIXED'  
TLSLP = 'FIXED'



```

IF(FIXFLG2 .NE. 0)STAMP = 'FIXED *'
IF(FIXFLG4 .NE. 0)TLAMP = 'FIXED *'
IF(NFLG(1) .NE. 0)ALPJ = 'FREE'
IF(NFLG(2) .NE. 0)STAMP = 'FREE'
IF(NFLG(3) .NE. 0)STSLP = 'FREE'
IF(NFLG(4) .NE. 0)TLAMP = 'FREE'
IF(NFLG(5) .NE. 0)TSLP = 'FREE'

```

```
WRITE(LOUT,1010),KTGRP,EST,EEND,IST,IEND
```

```
TSLPS = '(as prespecified)'
```

```

IF (BKLO .EQ. 0) THEN
  TBKGD = ' REMOTE '
  WRITE(LOUT,1011)TBKGD,RBKWID(1),NUMS,'low'
  WRITE(LOUT,1013)RBKEN(1)
  WRITE(LOUT,1012)SLP(1),TSLPS
  GO TO 15
ENDIF

```

```

  TBKGD = 'Adjacent'
  IF (RBKEN(1) .EQ. 0) THEN
    WRITE(LOUT,1011)TBKGD,BKLO,NUMS,'low'
    WRITE (LOUT,1012)SLP(1),TSLPS
  ELSE
    TSLPS = 'determined internally'
    WRITE(LOUT,1011)TBKGD,BKLO,NUMS,'low'
    WRITE (LOUT,1012)SLP(1),TSLPS
    WRITE(LOUT,1014)RBKWID(1),RBKEN(1)
  ENDIF

```

```

TSLPS = '(as prespecified)
IF (BKHI .EQ. 0) THEN
  TBKGD = ' REMOTE '
  WRITE(LOUT,1011)TBKGD,RBKWID(2),NUME,'high'
  WRITE(LOUT,1013)RBKEN(2)
  WRITE(LOUT,1012)SLP(2),TSLPS
  GO TO 20
ENDIF

```

```

  TBKGD = 'Adjacent'
  IF (RBKEN(2) .EQ. 0) THEN
    WRITE(LOUT,1011)TBKGD,BKHI,NUME,'high'
    WRITE (LOUT,1012)SLP(2),TSLPS
  ELSE
    TSLPS = 'determined internally'
    WRITE(LOUT,1011)TBKGD,BKHI,NUME,'high'
    WRITE (LOUT,1012)SLP(2),TSLPS
    WRITE(LOUT,1014)RBKWID(2),RBKEN(2)
  ENDIF

```

```

1010  FORMAT(// ---- Analyzed peak grouping No.',I3,' from',F9.3,
+     / ' to',F9.3,' keV'// by reading channels',I5,' to',I5,' and',
+     / ' subtracting background determined over')
1011  FORMAT (3X,A,' area of',F8.3,' keV(',I3,' chans) on the ',A,
+     / ' E side,'$)
1012  FORMAT (4X,'extrapolated with a slope of',F9.4,'%keV; ',A)
1013  FORMAT('+ beginning at',F8.3,' keV')
1014  FORMAT(5X,'using also a REMOTE area of',F8.3,' keV starting',

```

```

+ ' at',F8.3,'keV')
20 IF(NGAIN .NE. 0)GO TO 25
WRITE(LOUT,1020)STD(1),STD(2)
GO TO 30
200 FORMAT(/' Reference peak energies used for gain:',F8.3,
1 ' and',F8.3)
25 WRITE(LOUT,1025)
1025 FORMAT (/ ' GAIN IS FIXED')

IF(INTFLG .GT. 1)WRITE(LOUT,1030),GAIN,ZERO
1030 FORMAT(' GAIN was redetermined as',F8.4,' and ZERO as',F8.3)
30 IF (INTFLG .EQ. 0)GO TO 40
TQFTL = 0.
CHAN = PKPOS(1)
EPK = ENJ(1)
CTS = SUM(1)
PCERR = ERR(1)
C Estimate FWHM for summed peaks, if needed.
IF(FWHM .LE. 0) FWHM = SQRT (.3 + .002*EPK)
NP = 1
IF(IOFLG(3) .EQ. 0) WRITE(LOUT,1031)KTGRP
1031 FORMAT(/' - Using simple integration for group #',I3)

IF(RUNFLG(4) .NE. 0) THEN
IF (LONGPR) WRITE (LOUT,1035)
IF(IOFLG(3).EQ.0) WRITE(LOUT,1035)
ELSE
IF (LONGPR)WRITE (LOUT,1036)
IF(IOFLG(3).EQ.0) WRITE(LOUT,1036)
ENDIF

1035 FORMAT (/ ' Peak Channel Energy Counts ',
+ ' Gammas/min Pcterr')
1036 FORMAT (/ ' Peak Channel Energy Counts ',
+ ' Counts/min * Pcterr')

GO TO 155

40 IF(IOFLG(3) .NE. 0) GO TO 43
WRITE (LOUT,1042) ALPJ,STAMP,STSLP,TLAMP,TLSLP
IF(NUSHPS .NE. 0) THEN
WRITE(LOUT,1040)
ELSE
WRITE (LOUT, 1041)
ENDIF
IF ((FIXFLG2 + FIXFLG4) .NE. 0) WRITE(LOUT,1043)
1040 FORMAT (' These fitting parameters initialized as specified ',
+ ' in the control file.')
1041 FORMAT (' Initial values of fitting parameters taken from',
+ ' shape parameter file.')
1042 FORMAT(' Parameter type: ',9X,'Amplitude Slope'/
+ ' Alpha',23X,A/' Short term tail',13X,A,A/
+ ' Long term tail',14X,A,A/)
1043 FORMAT(' Starred parameters turned negative 5 times & were'
+ ' fixed at a small positive value')

43 SG = -2.7726/ ALPHA
POS = NDPTS/2 !For FWHM at mid-energy.

```

```
IF (NJ .EQ. 1) POS = PKPOS(1) !At peak position
SG = SG + ASLP*POS - .462
ALFA = -2.7726/SG
FWHM = SQRT(SG)
TSHP5 = EXP2/GAIN
TSHP8 = EXP4/GAIN
```

```
44 DO 44 I=1,100
QFT(I)=0.0
QFT(50)=1.
DX=0.
L=50
```

```
DO 46 I=51,100
DX=DX+1.
L=L-1
QFT(I)=SQRT(EXP(ALFA*DX*DX))
QFT(L)=QFT(I)
IF(QFT(I)-.01)48,46,46
46 CONTINUE
```

```
48 ASUM = 0.0
DO 50 J=1,NP
SDOF(J)=0.0
SUMR(J) = 0.0
```

```
C CALCULATE PORTION OF PEAK AREAS OUTSIDE BOUNDARIES
```

```
IF (EXP2 .NE. 0.0) YSUM(J)=YSUM(J)+PKHT(J)*EXP1/EXP2*
+ EXP(-EXP2*PKPOS(J))
```

```
FRACT1 = 0.
FRACT2 = 0.
IF(GAMA(J) .LE. 0) GOTO 49
```

```
C GAMA(J)=GAMA(J)*.001/GAIN
FRAC1 = ATAN(2. * PKPOS(J)/GAMA(J))
FRACT1 = 1. - .6366 * FRAC1
DELT = NDPTS - PKPOS(J)
FRAC2 = ATAN(2. * DELT/GAMA(J))
FRACT2 = 1. - .6366 * FRAC2
YSUM(J) = YSUM(J) * (1. + FRACT1 + FRACT2)
```

```
49 ASUM = ASUM + YSUM(J)
50 CONTINUE
```

```
TSUM=0.0
```

```
DO 60 I = 1,NDPTS
YT = YNET(I) + AVEBG
IF (YT .LE. 0.0) YT = 1.0
RES(I) = RM(I) * SQRT (YT)
RR = RM(I)**2
TSUM = TSUM + RR
```

```
DO 55 J=1,NP
NPK=PKPOS(J)+0.5
K=I-NPK+50
IF(K.GT.100)GOTO 55
IF(K .LE. 0)GOTO 55
SDOF(J)=SDOF(J)+QFT(K)
SUMR(J)=SUMR(J)+QFT(K)*RR
55 CONTINUE
```

```

60     CONTINUE

      L = KK
      DO 70 J=1,NP
      VM(J) = 0.      !Will use to re-sort diagonal matrix elem.
      ENDM(J) = 0.
70     TCTS(J) = 0.
      DO 80 J = 1,NP
      IF(NYFLG(J) .LE. 0)GO TO 73
      L = L + 1
      VM(J) = DM(L)
      I = J
      GO TO 75
73     I = -NYFLG(J)
75     TCTS(I) = TCTS(I) + YSUM(J)
      IF(NXFLG(J) .LE. 0)GO TO 80
      L = L + 1
      ENDM(J) = DM(L)
80     CONTINUE

      DOF = NDPTS - L
      QFIT = TSUM / DOF
      RMSD = 1. + 100.*(DOF * (QFIT - 1.0) / ASUM)
      FWHM = FWHM * GAIN

      IF (LONGPR)WRITE (LOUT,1090)KTGRP,ITER,QFIT,FWHM,ERR(3),GAIN,ZERO,
1     EXP1,ERR(4),TSHPS,ERR(5),RMSD,EXP3,ERR(6),TSHPS,ERR(7)
1090  FORMAT(/12X,'Results from group',I3,' (' ,I3,' iterations)'/
+ ' QFIT',F8.3,'; FWHM',F8.5,'+/- ',F6.4,'; KEV/CHANNEL',
+ F9.5,'; ZERO',F9.3,/
+ 20X,' EXP1',F9.5,'+/- ',F7.4,' EXP2',F9.5,'+/- ',F7.4,' NQFIT',
+ F6.3,8X,' EXP3',F9.5,'+/- ',F7.4,' EXP4',F9.5,'+/- ',F7.4)

      IF(IDFLG(3) .NE. 0)GO TO 150
      WRITE(LOUT,1095)
1095  FORMAT('OPEAK ENERGY WIDTH PEAK-ENERGY PEAK-',
+ 'INTENSITY',18X,'AREA'/' NO. keV eV',
+ 7X,2('SPECS',14X),8X,'RATIO'/)
      DO 110 I=1,NP
      TOUT = '
      TENJ = ENJ(I)
      IF(TENJ .LT. 0)TENJ = -TENJ !Remove X-RAY flag.
      TGAMA = GAMA(I)*GAIN/.001
      IF (GAMA(I) .NE. 0) WRITE (TWID,1100).GAMA(I)*GAIN/.001

      IF(NXFLG(I).LT.0) THEN
          WRITE (TNGY,1101), 'FIXED VS # ',-NXFLG(I)
      ELSEIF (NXFLG(I) .EQ. 0) THEN
          TNGY = ' FIXED'
      ELSE
          TNGY = ' FREE'
      ENDIF

      IF (NYFLG(I) .GT. 0) THEN
          TINT = ' FREE'
          GO TO 105
      ELSE IF (HIGHT(I) .GE. 0) THEN

```

```

WRITE (TINT,1103)' Fixed at ',HIGHT(I)
ELSE
WRITE (TINT,1104) -NYFLG(I),YSUM(I)/YSUM(-NYFLG(I))
ENDIF
105 WRITE (LOUT,1105) I,TENJ,TOUT
110 CONTINUE

1100 FORMAT(F6.1)
1101 FORMAT(A,I2)
1103 FORMAT(A,F10.2,' Counts')
1104 FORMAT('Fixed vs # ',I2,13X,F9.5)
1105 FORMAT(I4,F10.3,A)

WRITE(LOUT,1090)KTGRP,ITER,QFIT,FWHM,ERR(3),GAIN,ZERO,EXP1,
1 ERR(4),TSHP5,ERR(5),RMSD,EXP3,ERR(6),TSHP8,ERR(7)

150 IF(RUNFLG(4) .NE. 0) THEN
IF (LONGPR)WRITE (LOUT,1150)
IF(IOFLG(3) .EQ. 0)WRITE(LOUT,1150)
ELSE
IF (LONGPR)WRITE (LOUT,1151)
IF(IOFLG(3) .EQ. 0)WRITE(LOUT,1151)
ENDIF

1150 FORMAT (/ ' Peak Channel Energy Error Counts ',
+ ' Gammas/min Pcterr Qfit')
1151 FORMAT (/ ' Peak Channel Energy Error Counts ',
+ ' Counts/min * Pcterr Qfit')

155 DO 230 J = 1,NP
IF(INTFLG .NE. 0)GO TO 175
SG = -2.7726/ALPHA + ASLP * PKPOS(J)
ALFA = -2.7726/SG
PCERR = 0.0
L = J
IF(NYFLG(J) .LT. 0)L = -NYFLG(J)
IF(PKHT(L) .LE. 0)GO TO 160

TQFTL = SUMR(L)/(SDOF(L)-1.)
IF(TQFTL .LT. 1.) TQFTL = 1.0
PCERR=SQRT(TCTS(L)/YSUM(J))*(100.*VM(L)*SQRT(TQFTL))/PKHT(L)

160 CTS = YSUM(J)
CHAN = PKPOS(J) + OFFSET
ERRNGY = 0.
IF(NXFLG(J) .GT. 0)ERRNGY = ENDM(J) * GAIN

175 EPK = REFEN + (CHAN - REFCH) * GAIN

C SUBTRACT KNOWN BACKGROUND PEAKS FOR THIS DETECTOR

BKSMBL = ' '
DO 180, JJ = 1,NBKPKS
IF(BKPKEN(JJ) .EQ. 0)GO TO 190
IF(ABS(EPK-BKPKEN(JJ)) .GT. FWHM)GOTO 180 !Window = FWHM keV
BKSMBL = 'b'

```

```

CTS1 = CTS
CTS = CTS - BKPKCT(JJ)*TIMEL
IF(CTS .LE. 0.)GOTO 179
PCERR = PCERR * CTS1 / CTS
IF (PCERR .LT. 100.) GO TO 180
CTS = 0.0
PCERR = 0.01 * PCERR * CTS1
180 CONTINUE

C THIS PORTION OF CODING CALCULATES ABSOLUTE PHOTON INTENSITIES

190 GAMAS =0.
VALUE=CTS
IF (CTS .EQ. 0.0) VALUE = 2.0 * PCERR

C CALGPM is always called for absorber corrections as well as efficiency
C If efficiency corrns. not made, counts/min will be corrected.
CALL CALGPM(EPK,VALUE,PCERR)

IF(CTS .EQ. 0.) PCERR = 200.

GAMAS=VALUE
GAMAS = (GAMAS / TIMEL)
C Note that TIMEL is LIVTIM in minutes

IF(INTFLG .NE. 0)GO TO 206
IF(ERRNGY .EQ. 0.)GO TO 201
IF (LONGPR)WRITE (LOUT,1200)J,CHAN,EPK,ERRNGY,CTS,BKSMBL,GAMAS,PCERR,TQF
IF(IOFLG(3).EQ.0)WRITE(LOUT,1200)J,CHAN,EPK,ERRNGY,CTS,BKSMBL,
+ GAMAS,PCERR,TQFTL
JO FORMAT (I5,F9.3,F10.3,F8.3,F12.0,A1,1PE13.4,OPF9.3,F8.3)
GO TO 210
201 IF (LONGPR)WRITE (LOUT,1201)J,CHAN,EPK,CTS,BKSMBL,GAMAS,PCERR,TQFTL
1201 FORMAT (I5,F9.3,F10.3,8X,F12.0,A1,1PE13.4,OPF9.3,F8.3)
IF(IOFLG(3).EQ.0)WRITE(LOUT,1201)J,CHAN,EPK,CTS,BKSMBL,GAMAS,
+ PCERR,TQFTL
GO TO 210

206 IF (LONGPR)WRITE (LOUT,1206)J,CHAN,EPK,CTS,BKSMBL,GAMAS,PCERR
1206 FORMAT (I5,F9.3,F10.3,F12.0,A1,1PE13.4,OPF9.3,F8.3)
IF(IOFLG(3).EQ.0)WRITE(LOUT,1206)J,CHAN,EPK,CTS,BKSMBL,GAMAS,PCERR

C For identification purposes it will be important to have
C peak arrays ordered by energy.
210 NG = NGAMS+1
IF(NGAMS .EQ. 0)GO TO 225
IF(EPK .GE. ENGY(NGAMS))GO TO 225
DO 215, JK=NGAMS,1,-1
IF(EPK .GT. ENGY(JK))GO TO 216
215 CONTINUE
JK = 0
C Shuffle arrays to make a slot at JK+1
216 DO 220, I=NGAMS,JK+1,-1
ENGY(I+1) = ENGY(I)
PKS(I+1) = PKS(I)
PKE(I+1) = PKE(I)
CHNL(I+1) = CHNL(I)
CNTS(I+1) = CNTS(I)
GAMS(I+1) = GAMS(I)

```

```

PCTER(I+1) = PCTER(I)
QFTL(I+1) = QFTL(I)
BKFLGS(I+1) = BKFLGS(I)
220 CONTINUE
NG = JK+1

225 CHNL(NG) = CHAN
ENGY(NG) = EPK
CNTS(NG) = CTS
GAMS(NG) = GAMAS
PCTER(NG) = PCERR
QFTL(NG) = TQFTL
BKFLGS(NG) = BKSMBL
NGAMS = NGAMS+1
DELT = FWHM
PKE(NG) = EPK + DELT !Upper bound of energy window
IF (NG .EQ. 1) GOTO 228 !Don't apply to first peak
DEN = EPK - ENGY(NG-1) !Distance from previous peak
IF (DEN .GT. 2.*FWHM) GOTO 228 !Distance O.K.
DELT = .5*DEN !Use narrower energy window
PKE(NG-1) = ENGY(NG-1) + DELT !Re-adjust upper bound of previous peak
228 PKS(NG) = EPK - DELT !Lower bound
230 CONTINUE

IF(RUNFLG(4) .EQ. 0) THEN
    IF (LONGPR)WRITE (LOUT,1232) !Explanatory footnote on photon
    IF(IOFLG(3) .EQ. 0)WRITE(LOUT,1232)
1232 FORMAT('/ * Detector efficiency was assumed to be unity. '//)
ENDIF

232 IF(IOFLG(4) .NE. 0)GO TO 300
IF(INTFLG .NE. 0)GO TO 300

IF (LONGPR)WRITE (LOUT,1230)
IF(IOFLG(3) .EQ. 0)WRITE(LOUT,1230)
1230 FORMAT('/ CHANNEL ENERGY INDEX YNET
+ 'RECONSTRUCT RESIDUALS STD DEV')
YP = 0.
DBG = BG(2)-BG(1)
CALL AVE(1,NDPTS,TAREA,AV,YNET)
DO 240 I=1,NDPTS
    IICH= I + IST - 1
    EENG = (FLOAT(IICH) - REFCH)* GAIN + REFEN
    YP = YP + YNET(I)
    YR = YNET(I) + BG(1) + DBG * YP/TAREA
    IF (LONGPR)WRITE (LOUT,1240)IICH,EENG,I,YNET(I),YR,RES(I),RM(I)
    IF(IOFLG(3) .NE. 0)GO TO 240
    WRITE(LOUT,1240) IICH,EENG,I,YNET(I),YR,RES(I),RM(I)
1240 FORMAT(I5,F12.3,I6,2X,3F13.2,F10.2)
240 CONTINUE

300 IF(IOFLG(2) .EQ. 0)GO TO 400 !Not finished yet.
C Prepare a summary after last grouping.
IF(KTGRP .LE. 1)GO TO 320
IF (LONGPR)WRITE (LOUT,1300)KTGRP
IF(IOFLG(3) .EQ. 0)WRITE(LOUT,1300),KTGRP
70 FORMAT('/ SUMMARY OF PEAK RESULTS ON',I3,' GROUPS')

IF(RUNFLG(4) .NE. 0) THEN
    IF (LONGPR)WRITE (LOUT,1305)

```

```

        IF(IOFLG(3) .EQ. 0) WRITE(LOUT,1305)
    ELSE
        IF (LONGPR)WRITE (LOUT,1306)
        IF(IOFLG(3) .EQ. 0) WRITE(LOUT,1306)
    ENDIF

1305  FORMAT (/ ' Peak Channel      Energy      Counts      ',
+ ' Gammas/min  Pcterr   Qfit')
1306  FORMAT (/ ' Peak Channel      Energy      Counts      ',
+ ' Counts/min * Pcterr   Qfit')

    DO 310, I= 1,NGAMS
    IF (LONGPR)WRITE (LOUT,1206)I,CHNL(I),ENGY(I),CNTS(I),BKFLGS(I)
+ ,GAMS(I),PCTER(I),QFTL(I)
    IF(IOFLG(3) .NE. 0)GO TO 310
    WRITE(LOUT,1206),I,CHNL(I),ENGY(I),CNTS(I),BKFLGS(I),GAMS(I),PCTER(I),
+ QFTL(I)
310  CONTINUE

    IF(RUNFLG(4) .EQ. 0)THEN
        IF (LONGPR)WRITE (LOUT,1232)
        IF(IOFLG(3) .EQ. 0)WRITE(LOUT,1232)
    ENDIF

C      Convert to fractional error, after adding 0.5% for LDMTRX
C      (Even if evaluation not requested at this time)
320  DO 325 J=1,NGAMS
325  FR CER(J) = .01 * SQRT(.25 + PCTER(J)**2)

400  RETURN
    END

```



SUBROUTINE P238AB

P238AB

Version 1.1  
Author: Wayne Ruhter  
Date: 14-Oct-88  
Date: 01-Feb-90 Modified efficiency values  
Date: 09-Mar-90 Fixed bug in Pa-233 correction to 376-keV peak  
Date: 14-Mar-90 Added Specific heat calculations  
Date: 16-Mar-90 Added decay corrections for precision and bias measurements

\*\*\*\*\*

Abstract:

This program reads in gammas/min generated by the GRPANL program that are corrected for Cd absorption and absorption by the EP60/61 stainless steel container. This routine estimates the thickness (g/cm<sup>2</sup>) of the Pu sample, and small changes in the shape of the detector efficiency curve from that of a nominal 20% coax detector. The calculations are made using the measured intensities of the following 12 peaks: 152, 208, 238, 375, 376, 413, 583, 727, 743, 766, 786, 860 keV. The equations involve 8 unknowns: the dpm's of 228Th, 238Pu, 239Pu, 241Pu, and 241Am if present, the Pu sample thickness, the change in the slope of the efficiency curve, and the change in the curvature of the efficiency curve. The equations are non-linear in form, requiring an iterative approach to their solution. The following boundary constraints are also applied:  
.01 < Pu (g/cm<sup>2</sup>) < 40.

The following common blocks are specific to the P238AB program

INCLUDE 'GRPLGCM.PMS'  
INCLUDE 'BRPGMS.PMS'  
INCLUDE 'P238CCM.CMN'

The following common blocks are generic to all instruments

INCLUDE 'CONTRL.CMN'  
INCLUDE 'CLUNS.CMN'  
INCLUDE 'CABUND.CMN'  
INCLUDE 'CASSAY.CMN'  
INCLUDE 'DTIM.CMN'

CHARACTER RSPFILE\*20

DIMENSION ENRG(12),BR(5,12),PUMU(12),ANSR(6),XSUM(6),INDX(5),  
+ TEMP(12),TEMP2(12),DAT(5,12),XML(12),TIMEP(6),WT(12),  
+ ER(6),ER1(5),BRMS(6),GAMZ(6),PCT1(6),PCT2(6),SPHT(6),SPROW(5),  
+ ACOEF(325,16)



```

C           !Efficiency curvature will be turned on later.
NCURV = 0
C           !Initially assume 4 isotopic components
NM = 5
NUM= 5
NP = 12
C           !Fraction of sample considered to be Pu.
PUGPSC=3.
PUFRAC = 1.
FRACX = 1.-PUFRAC
C
DO 2 J = 1,12
DO 2 K = 1,5
C           !Transfer branching coef. into arrays.
BR(K,J) = DAT(K,J)
CONTINUE
C
C Read in gamma intensity information
C
OPEN (UNIT=INDEV,NAME=STGFIL,TYPE='OLD',ERR=5)
GO TO 6
CLOSE(INDEV)
WRITE (*,9999)STGFIL
WRITE (NDEV1,1003)STGFIL
9999 FORMAT(1X,15A1)
1003 FORMAT(1X,A,' File not found')
GO TO 4
ABORT = .TRUE.
RETURN
6
CALL RDGAMS
C
A376 = 0.
A311 = 0.
KK = 0
DO 10 J = 1,NGAMS
C IF (ABS(ENGY(J)-376.6).LE.0.5)THEN
C     A376 = GAMS(J)
C
C END IF
C IF (ABS(ENGY(J)-311.9).LE.0.5)THEN
C     A311 = GAMS(J)
C
C END IF
DO 8 K = 1,NP
DIFE = ENGY(J) - ENRG(K)
8 IF (ABS(DIFE).LE.0.6)GOTO 9
GOTO 10
9
KK = KK + 1
CNTS(KK) = GAMS(J)+1.
IF (KK.EQ.7)PCTER(J)=1.
WT(KK) = 1./(PCTER(J)*CNTS(KK))
IF (KK .NE.5)THEN
IF (PCTER(J).GE.2.)ABORT=.TRUE.
ENDIF
10
CONTINUE
C
C Check peak areas, if any are zero abort analysis
C
IF (ABORT)THEN
TYPE *,('Insufficient data for analysis. Analysis Aborted.')
GOTO 500

```

```

END IF
      !Get initial estimate of 228Th, 238Pu, 239Pu, and 241Pu
C      ATNUM = 237.
      DO 16 J = 1,NM
      K = INDX(J)
      ELOG = ALOG(ENRG(K))
      PUMAS = PUGPSC * PUMU(K)
      CORR = 1./PUMAS*(1.-EXP(-PUMAS))*EXP(SHAPC(12))
C      !Total correction factor
      + *(ELOG-5.3375)+SHAPC(13)*(ELOG**2-28.4889))
C      !Initial answers in dpm.
      GPM = CNTS(K)
      ATNUM = ATNUM+1.
      IF (J.EQ.3)ATNUM =241.
      IF (J.EQ.4)ATNUM = 241.
      IF (J.EQ.5)ATNUM= 228.
      CNVRT= 1.6603E-21*ATNUM*T1HLF(J)/.693147
      ANSR(J) = GPM/ BR(J,K) /CORR
      ANSR(J) = ANSR(J) * CNVRT
C      WRITE (*,*)'ANSR',J,ANSR(J),GPM,BR(J,K),CORR,CNVRT
16     CONTINUE
      CNTS375 = CNTS(4)
C
      IF (LONGPR)WRITE (LOUT,9102)
9102   FORMAT ('      PU      DELPU      ESLP      ',
      + 'DELSLP      ECURV      DELCURV')
      IF (LONGPR)
      + WRITE (LOUT,9104) PUGPSC,DELPU,ESLP,DELSLP,
      + ECURV,DELCRV
C
      ELOG = ALOG(311.9)
      PUMAS = PUGPSC * .476
      CORR = 1./PUMAS*(1.-EXP(-PUMAS))*EXP(SHAPC(12))
C      !Total correction factor
      + *(ELOG-5.3375)+SHAPC(13)*(ELOG**2-28.4889))
      CNTS311 = A311/CORR/.38
      ELOG = ALOG(ENRG(5))
      PUMAS = PUGPSC * PUMU(4)
      CORR = 1./PUMAS*(1.-EXP(-PUMAS))*EXP(SHAPC(12))
C      !Total correction factor
      + *(ELOG-5.3375)+SHAPC(13)*(ELOG**2-28.4889))
      CNTS376=CNTS311*CORR*.00655
      CNTS(4) = CNTS375-CNTS376
      DO 15 J=1,16
      DO 15 K=1,325
15     ACDEF(K,J)=0.
      DO 30 J = 1,NP
      PUMAS = PUGPSC *(PUMU(J)*PUFRAC +XMU(J)*FRACX)
      PUEXP = EXP(-PUMAS)
C      !Pu absorption
      PUCCF = (1.-PUEXP)/PUMAS
      ELOG = ALOG(ENRG(J))
C      !det. eff
      EFF1 = EXP(SHAPC(12))*(ELOG-5.3375)+SHAPC(13)*
1     (ELOG**2-28.4889))
C      IF (ENRG(J).EQ.152.48)EFF1=EFF1/1.05
      EN100 = ENRG(J) - 208.
C      !Change from nominal curve.
      EFF2 = (1. +ESLP*EN100 +ECURV*EN100**2)
C      WRITE (6,*)'EFF',ENRG(J),EFF1,EFF2,PUCCF

```

```

TEMP(J) = PUCOF *EFF1 *EFF2
TEMP2(J) = EFF1 *EFF2
C                                     !Composite correction
CORR = TEMP(J) *WT(J)
C
ASUM = 0.0
ATNUM = 237.
DO 20 K = 1,NM
ATNUM = ATNUM+1.
IF (K.EQ.3)ATNUM=241.
IF (K.EQ.4)ATNUM=241.
IF (K.EQ.5)ATNUM= 228.
CNVRT= 1.6603E-21*ATNUM*T1HLF(K)/.693147
C                                     !Calculate coefficient
ACDEF(J,K) = BR(K,J) *CORR/CNVRT
C                                     !Total dpms.
ASUM = ASUM + BR(K,J) *ANSR(K)/CNVRT
20 CONTINUE
C
CORR = CORR * ASUM
NCMP = NUM
C                                     !Pu not to be evaluated
C
1 IF (NPU .EQ. 0) GOTO 23
NCMP = NCMP + 1
PUG = PUGPSC
C                                     !Pu coefficient
ACDEF(J,NCMP)=CORR/PUCOF/PUG*(PUEXP-(1.-PUEXP)/PUMAS)
C                                     !Eff. slope not determined
23 IF (NSLP .EQ. 0) GOTO 30
NCMP = NCMP + 1
C                                     !Eff. slope coefficient
ACDEF(J,NCMP) = CORR/EFF2 *EN100
C                                     !Eff. curvature not to be determined
IF (NCURV .NE. 1) GOTO 30
NCMP = NCMP + 1
C                                     !Eff. curvature coefficient...
ACDEF(J,NCMP) = ACDEF(J,NCMP-1) *EN100
30 CONTINUE
C
C
DO 40 J = 1,NCMP
XSUM(J) = 0.0
DO 40 K = 1,NCMP
40 ALPHA(J,K) = 0.0
C
DO 42 J = 1,NP
DO 41 K = 1,NCMP
C
XSUM(K) = XSUM(K) +ACDEF(J,K) * CNTS(J) * WT(J)
C
DO 41 L = 1,NCMP
41 ALPHA(K,L) = ALPHA(K,L) +ACDEF(J,K) *ACDEF(J,L)
42 CONTINUE
C
C                                     !Invert matrix
CALL MATINV (ALPHA,B,NCMP,DET)
C
L IF (DET .EQ. 0.)THEN
WRITE (*,*)'MATRIX IS SINGULAR IN "P235AB" ROUTINE'

```

```
ABORT = .TRUE.  
END IF
```

```
DO 44 K = 1,NCMP  
ANSR(K) = 0.0  
DM(K) = APHA(K,K)  
DO 43 L = 1,NCMP
```

```
!Calculate answers  
ANSR(K) = ANSR(K) +XSUM(L) *APHA(K,L)  
CONTINUE
```

```
DO 48 J = 1,NP  
RM(J) = CNTS(J) * WT(J)  
DO 46 K = 1,NCMP
```

```
!Residual  
RM(J) = RM(J) - ACOEF(J,K) * ANSR(K)
```

```
CONTINUE
```

```
Determine goodness-of-fit
```

```
SUMRM = 0.  
DOF = NP - NCMP  
DO 49 J = 1,NP  
SUMRM = SUMRM + RM(J) * RM(J)  
QFIT = SUMRM/(DOF-1.)  
ICNVRG = 1  
NCMP = NUM
```

```
DELPU = 0.0
```

```
!Pu was not determined
```

```
IF (NPU .EQ. 0) GOTO 54  
NCMP = NCMP + 1  
DELPU = ANSR(NCMP)  
ADEL = ABS(DELPU)  
DELP = DELPU *(1.-ADEL/(PUGPSC+ADEL+1.))  
IF ((ABS((DELPU+OLDEL)/DELPU)) .LT. .1) DELP = 0.5*DELP
```

```
!Increment Pu answer.
```

```
PUGPSC = PUGPSC + DELP  
OLDEL = DELPU
```

```
!Check upper limit
```

```
IF (PUGPSC .LT. 100.) GOTO 52
```

```
!Maximum value
```

```
PUGPSC = 100.  
GOTO 54
```

```
!Check lower limit
```

```
IF (PUGPSC .GT. .002) GOTO 53
```

```
PUGPSC = .002  
GOTO 54
```

```
!Not converged
```

```
IF (ABS(DELPU) .GT. .1) ICNVRG = 0
```

```
!Eff. slope not determined
```

```
IF (NSLP .EQ. 0) GOTO 55
```

```
NCMP = NCMP + 1  
DELSLP = ANSR(NCMP)
```

```
!Absolute value of ESPL change
```

```
ADEL = ABS(DELSLP)
```

```
!Increment slope
```

```
ESLP=ESLP+DELSLP*(1.-ADEL/(ABS(ESLP)+.001+ADEL))
```

```

C
55 DELCRV = 0.0
C                                     !Eff. curvature not determined
C                                     IF (NCURV .NE. 1) GOTO 56
C                                     NCMP = NCMP + 1
C                                     DELCRV = ANSR(NCMP)
L                                     !Absolute value of ECURV change
C                                     ADEL = ABS(DELCRV)
C                                     !Increment
C                                     ECURV = ECURV + DELCRV*(1.-ADEL/(ABS(ECURV)+.00001 +ADEL))
C                                     IF (ADEL .GT. .000002) ICNVRG = 0
56 CONTINUE
C
C                                     IF (LONGPR)WRITE (LOUT,9031) QFIT
C                                     IF (LONGPR) WRITE (LOUT,9104) PUGPSC,DELPV,
+                                     ESLP,DELSLP, ECURV,DELCRV
9104 FORMAT (2F8.3,4F10.6)
C                                     IF (LONGPR)WRITE (LOUT,9130)
9130 FORMAT (/ ' Energy Gammas Wt. Fctr. Residual Tot. Eff. ',
+ ' Det. Eff. ')
C                                     IF (LONGPR)WRITE (LOUT,9090), (ENRG(J),CNTS(J),WT(J),RM(J),TEMP(J),
+ TEMP2(J), J=1,NP)
9090 FORMAT (OPF9.3,F11.0,F9.6,F8.1,P2E12.3)
C                                     !Has converged
C                                     IF ((NSLP .EQ. 1) .AND. (ICNVRG .EQ. 1)) GOTO 90
C                                     IF (ICNVRG.EQ.1)GOTO 97
C                                     !O.K. to calc. ESLP
C                                     IF (ABS(DELPV).LT. 0.1)NSLP =1
C
C                                     ITER = ITER + 1
C                                     Inspect 583-keV residual, if still large, reduce weighting
C                                     IF (ITER.GE.2 .AND. RM(7).GT. 5.0)WT(7) = WT(7)/1000.
C
C                                     IF (ITER .LT. 20) GOTO 18
C                                     WRITE (LOUT,9110)
9110 FORMAT ( ' WARNING: DID NOT CONVERGE IN "P238AB" ROUTINE ')
C
90 IF (NCURV .NE. 0) GOTO 92
C GOTO 92
C NPU = 0
C NCURV = 1
C GOTO 18
92 IF (NPU .EQ. 0) GOTO 96
C                                     !Test Pu boundary conditions
C NPU = 0
C GOTO 18
C
96 IF (NCURV .EQ. 3) GOTO 97
C NCURV = NCURV + 1
C GOTO 18
C
97 CONTINUE
C
C ATNUM = 237.
C SUM1 = 0.
C SUM2 = 0.

```

C The results of the least-squares solution of the set of equations used are  
C in units of relative disintegration rates of each of the isotopic

C components. These results are now decay corrected (if requested), and  
C reinterpreted for reporting purposes.

```
C
      DELT = 0.
      IF (MP .DR. MB)DELT = 1440.*(SDAT-32852.)
DO 70 J = 1,4
      IF (J.EQ.4)GOTO 72
      ATNUM=ATNUM+1
      IF (J.EQ.3)ATNUM=ATNUM+1
C
C                                     !Atomic number
72      ER(J) = 100.*SQRT(DM(J)) /ANSR(J)
C
C                                     !Statistical error
C      GRMS(J) = 1.6603E-21*ANSR(J)*ATNUM*T1HLF(J)/.693147
      GRMS(J) = ANSR(J)
C
C                                     !Milligrams of J
      DCYNST = .693147 / T1HLF(J)
C
C                                     !Decay constant
      GRMZ(J) = GRMS(J) * EXP (DCYNST * DELT)
C
C                                     !Grams at zero time
C                                     !Only add up Pu
      IF (J.GT.3)GOTO 70
      SUM1 = SUM1 + GRMS(J)
C
C                                     !Sum of grams
      SUM2 = SUM2 + GRMZ(J)
C
C                                     !Sum at zero time
```

```
70      CONTINUE
C
C Calculate Pu-240/Pu-239 ratio, if PU240<0
```

```
C
      C240 = 0
      P240 = 0.
      P240T = 0.
      IF (PU240.LT.0)THEN
          C240FB = .TRUE.
          C240 = 1
          GRMS40 = .07411*GRMS(2) + 2.9804*GRMS(3)
          SUM1 = SUM1 + GRMS40
          DCYNST = .693147 / T1HLF(6)
C
C                                     !Decay constant
          GRMZ40 = GRMS40 * EXP (DCYNST * DELT)
          SUM2 = SUM2 + GRMZ40
          ER240 = .1
      ELSE
```

```
C
C Correct totals for missing Pu-240 abundance
```

```
C
          C240FB=.FALSE.
          IF (PU240.EQ.0)THEN
              P240T = 2.0
              ELSE
              P240T = PU240
              ENDIF
          P240T = P240T *EXP(.693147 *DELT/ T1HLF(6))
          ER240 = 0.01
          P240 = P240T
```

```
      ENDIF
      G239 = GRMS(2)
```

```
C
      POW = 0.0
      BCPDW = 0.0
```



```

C
C ... CORRECT TOTALS FOR PU-242 ABUNDANCE
      IF (PU242.EQ.0)THEN
          P242T = 0.1
      ELSE
          P242T = PU242
      ENDIF
      SUMN = 1. - 0.01 * (P240T +P242T)
      SUM1 = SUM1 / SUMN
      SUM2 = SUM2 / SUMN
      ER242 = 0.1
      P242 = P242T
C
DO 78 J = 1,4
      !Percent abundance of Jth isotope
      WTPCT(J) = GRMS(J) / SUM1 * 100.
      !Percent abundance at zero time
      PCT2(J) = GRMZ(J) / SUM2 * 100.
      !Specific power due to Jth isotope
      SPPOW(J)= .01* SPHT(J)*PCT2(J)
      POW = POW + SPPOW(J)
C
      !Abundance relative to 239Pu
      GRMS(J) = GRMS(J) / 6239
      GRMZ(J) = GRMZ(J) / 6239
      IF (MP .OR. MB)THEN
          GRMS(J) = GRMZ(J)
          WTPCT(J) = PCT2(J)
      ENDIF
78 CONTINUE
      IF (C240.EQ.1)P240 = GRMS40/SUM1 * 100.

      SP240 = .01 * P240 * SPHT(6)
      POW = POW + SP240

      SP242 = .01 * P242 * .1146
      POW = POW + SP242
C
C The following set of calculations propogate the statistical uncertainties
C of the original results to the quantities being reported. The treatment
C acknowledges the covariant dependancy of the normalized isotopic
C abundances following the procedures reported in J. of INMM, vol. 12, No.2,
C (1983)
C
      SBPOW = 0.0
      XPCT = (.01*P240 *ER240 *SUM2)**2 + (.01*P242 *ER242 * SUM2)**2
DO 82 J = 1,4
      !Count and zero time fractions of J component
      FRAC1 = .01 * WTPCT(J)
      FRAC2 = .01 * PCT2(J)
      !Start with Pu240 component. Treat as independent variable.
      ER1(J) = (FRAC1 *P240 *ER240/WTPCT(2))**2 +
1          (FRAC1 * P242 * ER242/WTPCT(2))**2
      PCT(J) = (FRAC2 *P240 *ER240/PCT2(2))**2 +
1          (FRAC1 * P242 * ER242/PCT2(2))**2
      IF (J .LT. 4) XPCT = XPCT + ALPHA(J,J)
      NK = 3
C
      !If component is Am, set new limit.
      IF (J .EQ. 4) NK = 4
C
      !Do K-row by L-column integration

```

```

DO 81 K = 1,NK
C                                     !Assume it is a diagonal element (ie. K=J)
FK1 = 1.-FRAC1
FK2 = 1.-FRAC2
IF (K .NE. J .OR. J .EQ. 4) THEN
C                                     !K not= J or component is Am
    FK1 = -FRAC1
    FK2 = -FRAC2
ENDIF
IF (K .EQ. 4) THEN
C                                     !Component is Am.
    FK1 = 1.0
    FK2 = 1.0
ENDIF

DO 80 L = 1,NK
C                                     !Assume it is a diagonal element (ie. L=J)
FL1 = 1.-FRAC1
FL2 = 1.-FRAC2
IF (L .NE. J .OR. J .EQ. 4) THEN
C                                     !L not= J or component is Am
    FL1 = -FRAC1
    FL2 = -FRAC2
ENDIF
IF (L .EQ. 4) THEN
C                                     !Component is Am.
    FL1 = 1.0
    FL2 = 1.0
ENDIF

C                                     !Integrate sigma squared.
ER1(J) = ER1(J) + FK1 *FL1 *APHA(K,L)
PCT(J) = PCT(J) + FK2 *FL2 *APHA(K,L)
C                                     !Also integrate the specific power error.
IF (J .EQ. 4) SGPOW = SGPOW + (SPPOW(K)-POW)**2 *APHA(K,L)
80 CONTINUE
81 CONTINUE
C                                     !Calculate the uncertainties
ER1(J) = 100.*SQRT(QFIT*ER1(J))/(SUM1*FRAC1)
PCT(J) = 100.*SQRT(QFIT*PCT(J))/(SUM2*FRAC2)
82 CONTINUE
C                                     !Uncertainty in specific power
SGPOW = SGPOW + (SP240 * ER240)**2
SGPOW = 100.*SQRT(SGPOW) / (SUM2*POW)

C
C
C ... Open ASCII results file; first create file name from SAMPID

DO 85 II = 1,20
IF (SAMPID(II:II).EQ. ' ')GOTO 87
85 CONTINUE
87 IF (II.LE.0)THEN
    TYPE *,'Error in SAMPID; Results file ID = DUMMY.RST'
    RSFFILE = 'DUMMY.RST'
ELSE
    RSFFILE = SAMPID
    RSFFILE(II:II+4)='.RST'
ENDIF
OPEN (UNIT=DKLUN1,FILE=RSFFILE,STATUS='NEW',ACCESS='SEQUENTIAL')
WRITE (LOUT,9015)

```

```

9015  FORMAT(1H1,80(1H*))//)
      WRITE (LOUT,9017)
      WRITE (DKLUN1,9017)
9017  FORMAT(' PU-238 gamma-ray ')
      WRITE (LOUT,9020)SAMPID
      WRITE (DKLUN1,9020)SAMPID
9020  FORMAT(' Assay results for sample: ',A,/)
      WRITE (LOUT,9022)
9022  FORMAT(1X,80(1H*))//)
C
      WRITE (LOUT,9031) QFIT,FWHM
      WRITE (DKLUN1,9031) QFIT,FWHM
9031  FORMAT (' QFIT =',F7.2,15X,'FWHM at 766.4 keV =',F5.1,' keV',/)
      WRITE (LOUT,9032)
      WRITE (DKLUN1,9032)
9032  FORMAT (15X,'ISOTOPIC ANALYSIS AT')
      IF (MB .OR. MP)THEN
      WRITE (LOUT,9042)
      WRITE (DKLUN1,9042)
9042  FORMAT(19X,'ZERO TIME (1/1/1990) ',
+ /9X,' RATIOS          WT.PCT.    PCTERR    Milliwatts/(gm to
+tal Pu)')
      ELSE
      WRITE (LOUT,9033)
      WRITE (DKLUN1,9033)
9033  FORMAT(19X,'COUNT TIME ',
+ /9X,' RATIOS          WT.PCT.    PCTERR    Milliwatts/(gm to
+tal Pu)')
      ENDIF
      K = 0
      DO 105 J=1,3
      K = K +1
      IF (J.EQ.3)THEN
      K = K+1
      IF (C240.EQ.1)THEN
      WRITE (LOUT,9037)P240,SP240
      WRITE (DKLUN1,9037)P240,SP240
      ELSE
      IF (PU240.GT.0)THEN
      WRITE (LOUT,9036)P240,SP240
      WRITE (DKLUN1,9036)P240,SP240
      ELSE
      WRITE (LOUT,9038)P240,SP240
      WRITE (DKLUN1,9038)P240,SP240
      ENDIF
      ENDIF
      ENDIF
      WRITE (LOUT,9034) (237+K,BRMS(J),WTPCT(J),ER1(J),SPPQW(J))
      WRITE (DKLUN1,9034) (237+K,BRMS(J),WTPCT(J),ER1(J),SPPQW(J))
105  CONTINUE
9034  FORMAT (' PU',I3,' =',F10.7,10X,F11.5,F7.2,7X,F8.3)
      IF (PU242 .GT. 0.)THEN
      WRITE (LOUT,9039)P242,SP242
      WRITE (DKLUN1,9039)P242,SP242
      ELSE
      WRITE (LOUT,9041)P242,SP242
      WRITE (DKLUN1,9041)P242,SP242
      ENDIF
      WRITE (LOUT,9035)BRMS(4),WTPCT(4),PCT(4),SPPQW(4)

```

```
WRITE (DKLUN1,9035)GRMS(4),WTPCT(4),PCT(4),SPPOW(4)
9035 FORMAT (' Am241 =',F10.7,10X,F11.5,F7.2,7X,F8.3)
9036 FORMAT (' Pu240 Operator input value ',F11.5,14X,F8.3)
9037 FORMAT (' Pu240 Calculated value ',4X,F11.5,14X,F8.3)
9038 FORMAT (' Pu240 Default value ',4X,F11.5,14X,F8.3)
39 FORMAT (' Pu242 Operator input value ',F11.5,16X,1PE10.3)
9041 FORMAT (' Pu242 Default value ',4X,F11.5,16X,1PE10.3)
WRITE (LOUT,9040)POW,SGPOW
WRITE (DKLUN1,9040)POW,SGPOW
9040 FORMAT(45X,' Total =',F8.3,' +/-',F4.2,'%')
CLOSE(UNIT=DKLUN1)
```

```
C
C Pass results to the system. Isotopic values passed through CABUND.
ASANS = POW
ERRBAR = SGPOW/100.
WTPCT(5) = P240
WTPCT(6) = P242
500 RETURN
END
```

SUBROUTINE P238GRP

Last edit for VAX 18-Jan-85
Edited for Micro VAX II 30-Dec-86
Edited for Pu238 SRP project 30-Aug-89

Originated by R.Gunnink & W.D.Ruhter
Modified by J.B.Niday to use overlays, COMMON blocks pre-prepared
by EDIGRP, for easy loading, consistent COMMON definitions in all
sections and programs, and optional automatic operation.
Further modifications included reading in the peak groupings separately
from the other parameters, in order to allow multiple groups & also
saving the photon results for later nuclide identification in GEVAL.

A version of GRPANL modified for SRP Pu238 analysis system.

INCLUDE 'GRPLGCM.PMS'
INCLUDE 'GRPGMS.PMS'
INCLUDE 'GRPEVC.PMS'

INCLUDE 'CONTRL.CMN'
INCLUDE 'CFILES.CMN'
INCLUDE 'CASSAY.CMN'
INCLUDE 'CDIAG.CMN'
INCLUDE 'CLUNS.CMN'

DATA NROW/0/,NCOL/0/,NCMP/0/

INDEV = 1
INPDEV = 4
NDEV = 3

Read file of pre-edited control info & open the spectral file.

CALL READ38

NGAMS = 0 !Init counter.

IOFLG(5) = 0

IOFLG(2) = 0

CALL RDGRP !Read in a block of group specifications

IF(EST .EQ. 0.)IOFLG(2)=1

IF(IOFLG(2) .NE. 0)GO TO 200 !A summary is needed now.

CALL CALIB

KTGRP = KTGRP + 1

CALL LOGGRP

Look at peak position and FWHM results for the 766.4-keV
peak for measurement control purposes. Report problems.

IF (KTGRP, EQ.9)THEN !766.4-KEV PEAK

TYPE \*, 'FWHM IN P238SRP',FWHM

FW766 = FWHM !Keep this value for later

IF (ABS(PKPOS(1) - 3063.) .GT. 3.)THEN

WRITE (LOUT,1100)

FORMAT(' 766.4-keV peak outside limit. Check stabilizer:')

```

        DGFLG(4) = .TRUE.
        ENDIF
        IF (FWHM .GT. 2.4 )THEN
1150      WRITE (LOUT,1150)
+      FORMAT(' !WARNING! 766.4-keV energy resolution exceeds
        2.4 keV. Check count rate or detector.')
        DGFLG(3) = .TRUE.
        ENDIF
    ENDIF
    IF(EST .EQ. 0.)GO TO 50          !Abort the illegal group!
    IF(INTFLG .NE. 0)GO TO 200

C
    CALL FIT
    IF (ABORT)RETURN
    IF (KTGRP. EQ. 2)THEN
C      SHAPC(1) = FWHM**2 - 208. * SHAPC(2)
C      WRITE (*,*)'FWHM',FWHM,SHAPC(1)
    ENDIF

C
    IF(MFLGS(4) .EQ. 0)GO TO 200
C      Write profiles into ASCII files for later plotting.
C      IF(IOFLG(5) .EQ. 2)CALL PLOTFL

200    CALL OUT

        IF(IOFLG(2) .EQ. 0)GO TO 50
        CLOSE(INPDEV)

        IF(IOFLG(5) .EQ. 0)GO TO 220
        IF(IOFLG(5) .EQ. 2)TYPE 1205
1205    FORMAT(/' Run program HIST to obtain plots from *.TMP files')

C      For peak fitting only, RUNFLG(4)=0, RUNFLG(5)=0, old KALPFG=0
C      For photon calculation after fitting, RUNFLG(4)=1, old KALPFG =1
C      For full evaluation at once RUNFLG(4) & RUNFLG(5)=1, old KALPFG=3
C      For evaluation based on counts only RUNFLG(4)=0,RUNFLG(5)=1

C
220    CALL SVPKS

C      Continue with nuclide evaluation
    PRPK(1) = DCNST(17)
    PRPK(2) = DCNST(18)
    NROW = 0
    NCOL = 0
    NCMP = 0
    NPRGMS = 0
    MGMS = 0
    NSPEC = 0

C      Perform Pu238 analysis on stored results file
C      Pass FWHM at 766 keV to this subroutine for reporting
    FWHM = FW766
    CALL P238AB
    RETURN
    END

```

```

C *****
C PASS
C
C     Check the password in order to allow operator to
C     enter supervisory mode
C
C     CALLING ARGUMENTS:
C         OK = Logical flag indicating validity of operator password
C
C     SSJ 01-Jul-82
C     REV. 08-May-84 lc
C     REV. 19-APR-85 MPK      Remove logical*1 OK
C                             (Redundant with /CPASS/
C     REV. 04-APR-86 WDR      Modified for Micro VAX II
C     REV. 06-AUG-87 RDP      Modified to accept password w/o
C                             echo to screen. All character are
C                             converted to UPPERCASE.
C *****

```

```

C     SUBROUTINE PASS(OK)

```

```

C     CHARACTER          WORD*10
C     INTEGER JSW,BIT12,NOT12,CR,LF
C     REAL*8 PWORD

```

```

C     INTEGER*4          STATUS,VKID,LENGTH
C     INTEGER*4          SMG#CREATE_VIRTUAL_KEYBOARD
C     INTEGER*4          SMG#READ_STRING
C     INTEGER*4          SMG#DELETE_VIRTUAL_KEYBOARD

```

```

C     INCLUDE '($TRMDEF)'

```

```

C     The following common block files are generic for all instruments

```

```

C     INCLUDE 'CPASS.CMN'

```

```

C     EQUIVALENCE (PWORD,WORD)

```

```

C     DATA CR,LF /"15","12/"

```

```

C     OK = .FALSE.

```

```

C     Prompt operator for password

```

```

C     TYPE 1000

```

```

1000  FORMAT(/' You have selected a supervisory option',
X     ' which requires a password.'/)

```

```

C     Store up to 8 characters

```

```

C     STATUS = SMG#CREATE_VIRTUAL_KEYBOARD(VKID,'SYS#INPUT')
C     IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))

```

```

C     STATUS = SMG#READ_STRING(VKID,
C     WORD,
C     'Enter password ->',
C     10,
C     TRM#M_TM_NOECHO+TRM#M_TM_CVTLOW)
C     IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))

```

```
20      TYPE 1004
1004    FORMAT(1X)
```

```
C
      STATUS = SMG#DELETE_VIRTUAL_KEYBOARD(VKID)
      IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
```

```
C  Check if operator wanted to escape
```

```
C      IF(WORD(1:1).EQ.' ') GOTO 500
```

```
C  Compare with all of the passwords
```

```
C      DO 30 I=1,3
          IF(PWORD.EQ.PASSWD(I)) GOTO 40
```

```
30      CONTINUE
```

```
C  Improper password
```

```
C  Inform operator
```

```
C      TYPE 1002
```

```
1002    FORMAT(//10X,'**** IMPROPER PASSWORD - NO ENTRY',
X        ' TO SUPERVISORY MODE ****'/)
      GO TO 500
```

```
C  Password is OK - Enter supervisory mode
```

```
C 40      OK = .TRUE.
```

```
C 500      RETURN
          END
```



## SUBROUTINE PLOTFL

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

```

```

C      A module to prepare ASCII files for external plot routines.

```

```

C      Last edit for VAX 1-Oct-85

```

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

```

```

      INCLUDE 'GRPLGCM.PMS'

```

```

      INCLUDE 'GRPRFL.PMS'

```

```

      DATA PRFNAM/'P','R','F','L','O','O','.',',','T','M','P',19*' ',0/

```

```

      DATA NETFNM/'Y','N','E','T','O','O','.',',','T','M','P',19*' ',0/

```

```

      DATA REMFNM/'R','M','O','O','.',',','T','M','P',21*' ',0/

```

```

      ENCODE(2,1020,PRFNAM(5))KTGRP

```

```

      IF(PRFNAM(5) .EQ. ' ')PRFNAM(5) = 'O'

```

```

      ENCODE(2,1020,NETFNM(5))KTGRP

```

```

      IF(NETFNM(5) .EQ. ' ')NETFNM(5) = 'O'

```

```

      ENCODE(2,1020,REMFNM(3))KTGRP

```

```

      IF(REMFNM(3) .EQ. ' ')REMFNM(3) = 'O'

```

```

1020  FORMAT(I2)

```

```

      K= 0

```

```

      IF(MFLGS(4) .NE. 2) GO TO 20

```

```

      Transfer fixed peak buffer for plots

```

```

      K = 1

```

```

      DO 10 I = 1,NDPTS

```

```

      PRFILS(I,KPHF+1) = PRFILS(I,MXNP)

```

```

20  OPEN(NDEV,FILE=PRFNAM,STATUS='NEW')

```

```

      DO 50 L = 1,KPHF+K

```

```

      WRITE(NDEV,1025)L,EXPTID,SAMPID,EST,EEND

```

```

1025  FORMAT(' Profile',I3,' for ',2(1X,12A1),'; E range',F9.3,

```

```

      + ' to',F9.3' keV')

```

```

      WRITE(NDEV,1030)NDPTS

```

```

1030  FORMAT(I4)

```

```

      WRITE(NDEV,1040)(XEN(I),PRFILS(I,L),I=1,NDPTS)

```

```

1040  FORMAT(F12.2,',',',',F12.2)

```

```

50  CONTINUE

```

```

      CLOSE(NDEV)

```

```

      TYPE 1055,PRFNAM

```

```

1055  FORMAT(' Plot file written: ',30A1)

```

```

      DO 70 I = 1,NDPTS

```

```

      IF(SVYNET(I) .LE. 0)SVYNET(I) = 0.01

```

```

70  CONTINUE

```

```

      OPEN(NDEV,FILE=NETFNM,STATUS='NEW')

```

```

      WRITE(NDEV,1070)EXPTID,SAMPID,EST,EEND

```

```

1070  FORMAT(' Net & profiles for ',2(1X,12A1),'; E range',F9.3,

```

```

      + ' to',F9.3' keV')

```

```

      WRITE(NDEV,1030)NDPTS

```

```

      WRITE(NDEV,1080)(XEN(I),SVYNET(I),ERYNET(I),I=1,NDPTS)

```

```

1080  FORMAT(2(F12.2,',','),F10.2)

```

```

      CLOSE(NDEV)

```

```

      TYPE 1055,NETFNM

```

```
OPEN(NDEV,FILE=REMFNM,STATUS='NEW')
WRITE(NDEV,1090)EXPTID,SAMPID,EST,EEND
1090  FORMAT(' RM for',2(1X,12A1),'; E range',F9.3,
+      ' to',F9.3,' KeV')
WRITE(NDEV,1030)NDPTS
WRITE(NDEV,1095)(RM(I),I=1,NDPTS)
1095  FORMAT(F12.2)
CLOSE(NDEV)
TYPE 1055,REMFNM

RETURN
END
```

## SUBROUTINE PRECHK

```

C
C This routine performs a preliminary check on the acquisition
C of data to determine proper sample setup and the amount of time
C required for a specified precision
C
C Modified by WDR for Pu238 system Mar-15-90
C
C     INTEGER          LID_ERR, SAM_ERR
C     INTEGER*4        ISTAT, IMASK, ISTATE
C     INTEGER*4        LID_OPEN, LID_CLOSE
C     INTEGER*4        START_ROT
C     REAL             LEPS_DT_RATIO
C     LOGICAL*1        RFLG, EFLG, PROBLEM, STATUS
C     CHARACTER        REPLY*1
C
C     INCLUDE          'CASSAY.CMN'
C     INCLUDE          'CDATE.CMN'
C     INCLUDE          'CONTRL.CMN'
C     INCLUDE          'CLUNS.CMN'
C     INCLUDE          'CVAX.CMN'
C
C     INCLUDE          'CONT90.CMN'
C     INCLUDE          'CSPECT.CMN'
C
C     INCLUDE          'ADAC.CMN'
C
C 10  CONTINUE
C     TYPE *, 'ENTERING PRECHK'
C     TYPE*
C     TYPE*, 'When sample is ready for ASSAY , '
C     type 1000, ' press RETURN ->'
C     READ(5,2000) REPLY
C 12  LID_ERR = 0
C     SAM_ERR = 0
C Test if Lid is closed
C Check if Door is closed
C     TYPE*, 'CHECK LID'
C
C     IMASK = 1
C     ISTAT = BTSTW(IMASK, LID_OPEN, C_DATA, IOSB, , ,)
C     TYPE *, 'ISTAT', LID_OPEN
C     STATUS = ERROR_CHECK(ISTAT, IOSB)
C
C     IMASK = 2
C     ISTAT = BTSTW(IMASK, LID_CLOSE, C_DATA, IOSB, , ,)
C     TYPE *, 'ISTAT', LID_CLOSE
C     STATUS = ERROR_CHECK(ISTAT, IOSB)
C
C If Lid is close (and not open) then precede with check.
C
C     IF ((LID_CLOSE.EQ.1).AND.(LID_OPEN.EQ.0)) GOTO 20
C
C     IF (LID_CLOSE.EQ.LID_OPEN)LID_ERR=1
C     IF (LID_ERR.EQ.1) THEN
C     TYPE *, 'ERROR -- Problem with LID closing or position switches'
C     GOTO 15
C     ENLIF
C 15  TYPE 1000, 'Press RETURN to retry (type A and a CLR) to ABORT'
C     READ (5,2000)REPLY

```

```
IF ((REPLY.NE.'A').AND.(REPLY.NE.'a'))GOTO 12
ABORT = .TRUE.
GOTO 900
```

```
Enable Rotation/scanning and verify.
```

```
IMASK = 512
TYPE*, 'START ROTATING'
ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
CALL ERROR_CHECK(ISTAT,IOSB)
```

```
Wait 2 sec
CALL SLEEP(2.0)
```

```
Look at sample size before checking elevator
```

```
CALL SIZESAMP
TYPE 9998,MATTYP(1)
9998 FORMAT(1X,'MATTYP ',A1)
```

```
Test if sample is rotating
```

```
TYPE*, 'CHECK ROTATING'
CALL ROTMON(RFLG,EFLG)
CALL ROTMON(RFLG,EFLG)
IF (.NOT.RFLG) THEN
  TYPE*, 'ERROR -- Unable to ROTATE Sample'
  TYPE*
  TYPE1000, 'Press RETURN to retry (type A and a <CR> to ABORT) ->'
  READ(5,2000) REPLY
  IF ((REPLY.NE.'A').AND.(REPLY.NE.'a')) GOTO 12
  ABORT = .TRUE.
  GOTO 900
```

```
ENDIF
```

```
TYPE *, 'EFLG IN PRECHK',EFLG,MATTYP(1)
IF (.NOT. EFLG .AND. MATTYP(1) .NE. 'F')THEN
TYPE *, 'ERROR -- Sample elevator not working. Check elevator !!'
TYPE*
TYPE1000, 'Press RETURN to retry (type A and a <CR> to ABORT) ->'
READ(5,2000) REPLY
IF ((REPLY.NE.'A').AND.(REPLY.NE.'a')) GOTO 12
ABORT = .TRUE.
GOTO 900
```

```
ENDIF
```

```
TYPE*, 'PLEASE WAIT 60 SECONDS WHILE ASSAY CHECK IN PROGRESS...'
```

```
Clear Time and Erase existing spectrum.
```

```
CALL CTIM(IADC,,IERR)
CALL CHK90('CTIM ',IERR)
CALL CDAT(IADC,,IERR)
CALL CHK90('CDAT ',IERR)
```

```
Set Clock Interval
```

```
Start Acquisition
```

```
PTIME = 60.0
CALL PSET(IADC,,1,1,PTIME,0.0,0,0)
```

```
Check that sample is rotating
```

```

CALL ROTMON(RFLG,EFLG)
IF (.NOT.RFLG) THEN
  TYPE*, 'ERROR -- SAMPLE NOT ROTATING. '
  TYPE1000, 'Press RETURN to retry (type A and a <CR> to ABORT) ->'
  READ(5,2000) REPLY
  IF ((REPLY.NE.'A').AND.(REPLY.NE.'a')) GOTO 20
  ABORT = .TRUE.
  GOTO 900
ENDIF

C
C Turn on ADC for collection
49  CALL COLL(IADC,,IERR)
    CALL CHK90('COLL ',IERR)

C
C Wait Until Acquisition is completed
.50  IWAIT1 = INQADC(IADC,ISTAT,MEMSIZ,IGAIN,IRANGE,IOFF,INPUTS)
     CALL CHK90('INQADC',IWAIT1)
     IF (MATTYP(1) .EQ. 'E')GOTO 51
     CALL SIZESAMP
C
TYPE 9998,MATTYP(1)

C Test if count is completed
51  IF (IWAIT1.EQ.0) GOTO 100

C
C Check that sample is rotating
52  CALL ROTMON(RFLG,EFLG)
     IF (.NOT.RFLG) THEN
       CALL ADCOFF(IADC,,IERR)
       TYPE*, 'ERROR -- SAMPLE NOT ROTATING. '
       TYPE1000, 'Press RETURN to retry (type A and a <CR> to ABORT) ->'
       READ(5,2000) REPLY
       IF ((REPLY.NE.'A').AND.(REPLY.NE.'a')) GOTO 49
       ABORT = .TRUE.
       GOTO 900
     ENDIF
     TYPE *, 'EFLG IN PRECHK',EFLG,MATTYP(1)
     IF (.NOT. EFLG .AND. MATTYP(1) .NE. 'F')THEN
       CALL ADCOFF(IADC,,IERR)
       TYPE *, 'ERROR -- Sample elevator not working. Check elevator !!!'
       TYPE*
       TYPE1000, 'Press RETURN to retry (type A and a <CR> to ABORT) ->'
       READ(5,2000) REPLY
       IF ((REPLY.NE.'A').AND.(REPLY.NE.'a')) GOTO 49
       ABORT = .TRUE.
       GOTO 900
     ENDIF

     GOTO 50

C Get Live Time (channel #0) and Clock Time (channel #1)
100  CONTINUE
     CALL FROMCA(IADC,2,0,SPCTRM,IERR,)
     CALL CHK90('FROMCA',IERR)
     TIMELV1 = SPCTRM(1)
     TIMECL1 = SPCTRM(2)
     LEPS_DT_RATIO = (TIMECL1-TIMELV1)/TIMECL1

C For assay, get 152-kev peak region and determine count rate
IF (ME .OR. MP)GOTO 900
CALL FROMCA(IADC,630,600,SPCTRM,IERR)

```

```

CALL CHK90('FROMCA',IERR)
CALL MAXVAL(2,28,JMAX,PHT,SPCTRM)
N1 = JMAX-12
IF (N1.LE.0)N1 = 1
N2 = N1 +2
M1 = JMAX + 10
M2 = M1 + 2
CALL NCNTS(N1,N2,M1,M2,ND,BGD,A152,ER152,AVBG,SPCTRM,SPCTRM)
CR152 = A152/TIMECL1
IF (CR152 .LE. 0.)THEN
  TYPE *, ' !!! No significant data collected. Check system !!! '
  TYPE *, ' Press ENTER to continue (type A and a <CR> to ABORT) ->'
  READ (5,2000) REPLY
  IF ((REPLY.EQ.'A').OR.(REPLY.EQ.'a')) THEN
    TYPE*, ' Abort assay - problem with data collection'
    ABORT = .TRUE.
    GOTO 900
  ENDIF
ENDIF
END IF
TIME = 400000./CR152
Type *, ' To count for ',time,' secs type Y'
Type *, ' To count for default time ',timepr,' secs type N'
Type 1000,' and press RETURN ->'
READ (5,2000)REPLY
IF (REPLY .EQ. 'Y')TIMEPR=TIME
RETURN
C
900 CALL AEND
RETURN
C*****
1000 FORMAT(/,1X,A,%)
2000 FORMAT(A)
END

```

```

C*****
C PRINTL
C
C Routine specific to gamma ray instruments
C Print assay
C Routine specific to the plutonium solution and solids
C gamma-ray analysis instruments
C
C SSJ 28-FEB-79
C Rev. 08-Jun-84 lc
C Rev. 06-NOV-84 WDR Revised for plutonium solution gamma-ray
C analysis
C*****
C
C SUBROUTINE PRINT1
C
C LOGICAL*1 SUMMRY
C
C The following common block files are generic for all instruments
C
C INCLUDE 'CONTRL.CMN'
C INCLUDE 'CASSAY.CMN'
C INCLUDE 'CLUNS.CMN'
C INCLUDE 'CBACKG.CMN'
C
C
C IF(LOUT .NE. 6 .AND. (.NOT. FROMDK .OR. ICYCLE .EQ. 1))
C X WRITE(LOUT,1998)
C 1998 FORMAT('1')
C SUMMRY = .FALSE.
C CALL HEADER(SUMMRY)
C
C Print final results (short printout)
C
C 200 WRITE(LOUT,1044)ASANS,ERRBAR
C 1044 FORMAT('/ Result = ',OPG13.4,' +/-',2PF6.1,' % milliwatts/gm to
C xtal Pu')
C
C
C 500 WRITE(LOUT,1998)
C CLOSE (LOUT)
C 510 RETURN
C END

```

C\*\*\*\*\*

C PUT90

C

C Routine specific to gamma ray instruments

C Get a spectrum from disk and write it to the Series 90;  
Write chunks of 256 channels at a time

C

C Calling arguments:

C Input: NADC - number of ADC, I

C IBEGIN - logical flag

C

C SSJ 28-JUN-79

C

C Rev. 13-Jun-84 lc

C Rev. 04-Apr-86 WDR

C Rev. 06-Jul-87 RDP

Modified for Micro VAX II

Modified for uVAX II w/ADAC Drivers

C

C\*\*\*\*\*

C

SUBROUTINE PUT90(NADC,IBEGIN)

C

C

C The following common block files are generic for all instruments

C

INCLUDE 'CONTRL.CMN'

INCLUDE 'CLUNS.CMN'

INCLUDE 'CBACKG.CMN'

INCLUDE 'CFILES.CMN'

INCLUDE 'CASSAY.CMN'

INCLUDE 'CDATE.CMN'

C

C The following common block files are instrument specific for  
gamma-ray instruments

C

INCLUDE 'CSPECT.CMN'

INCLUDE 'CONT90.CMN'

C

C Open disk file. Allocate channel and free when through

C

15 CONTINUE

OPEN(UNIT=DKLUN,NAME=FILNAM,ACCESS='DIRECT',TYPE='OLD',

X RECORDSIZE=128,ERR=20)

IREC = 1

GO TO 25

C

C Problem in opening the disk file

C Inform the user and return

C

20 CALL ERRMSG(1,FILNAM)

C

C Free channel

C

ABORT = .TRUE.

GO TO 500

C

C Open ADC for writing

C

25 CONTINUE

CALL OPEN80(NADC,,'OU',,IERR)

CALL CHK90('OPEN80',IERR)



ISTART = IBEGIN

```
C
C Read block in from disk
C Check for dummy variable signalling background appendage
C
)   READ(DKLUN'IREC)SPCTRM
    IF(SPCTRM(1) .EQ. -999.0)GO TO 100
C
C Write block of 128 channels to Series 90
C
    NUMCHN = 128
    CALL WRIT80(SPCTRM,NUMCHN,ISTART,IERR)
    CALL CHK90('WRIT80',IERR)
    IF(IERR .NE. 0) GO TO 100
    IREC = IREC + 1
    ISTART = ISTART + 128
    GO TO 30
C
C Close the ADC and the disk file
C
100  CONTINUE
     CALL CLOS80(IERR)
     CALL CHK90('CLOS80',IERR)
     READ(DKLUN'IREC)DUM,BK,SGBK,TIMEBK,BKDAT,BKTIM,
     X      BKOP,BKSFIL,SAMPID,OPERID,SYSID,ASDATE,ASTIME,
     X      MATTYP(1)
101  CLOSE(UNIT=DKLUN)
C
C Check if background has been written at end of data file
C
     IF(DUM .NE. -999.0) NOBACK = .TRUE.
     NBLKS = IREC-1
     TYPE 1006,NBLKS,FILNAM
1006  FORMAT('// Spectrum retrieval completed --',
     X      I6,' blocks read from ',15A)
C
C Set the live data flag to false
C
     LIVE = .FALSE.
C
500  RETURN
     END
```



SUBROUTINE RDGAMS

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C      Read previous output of GRPANL runs on various groupings to
C      reconstruct the photons/min table(s). This permits combination with
      other spectra before entering the final evaluation phase.
C
C      Last edit for RSX11 15-Jan-88 R.G.
C      Modified lines written in long printout 2-Feb-90
C
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C
      INCLUDE 'GRPLGCM.PMS/LIST'
      INCLUDE 'GRPGMS.PMS/LIST'
      INCLUDE 'GRPEVC.PMS/LIST'
C
      INCLUDE 'CONTRL.CMN/LIST'
C
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C
      BYTE LINEBF(120),SAMPID(12)
C
      MXPKS = 50
      MXFNM = 30
C
      ITEMCT = 0
      LINECT = 0
      LFG = 0
C
30      READ(INDEV,1000,END=310,ERR=300)NCHRS,(LINEBF(K),K=1,NCHRS)
1000     FORMAT(Q,120A1)
      LINECT = LINECT+1
      IF(LONGPR.AND.LINECT.LE.7)WRITE(6,1030)(LINEBF(K),K=1,NCHRS)
1030     FORMAT(120A1)
C      Look for "GRP" in beginning line
      DO 32, J=1,3
      IF(LINEBF(J) .NE. 'S')GO TO 32
      IF(LINEBF(J+1) .NE. 'p')GO TO 32
      IF(LINEBF(J+2) .EQ. 'e')GO TO 33           !Found it
32      CONTINUE
      GO TO 34
33      CONTINUE
CTYPE 1030,(LINEBF(K),K=1,NCHRS)
34      IF(NCHRS .LT. 2)GO TO 30
      IF(LFG .GT. 0)GO TO 50
C
      DO 35, I = 1,NCHRS
      IF(LINEBF(I) .EQ. '=')GO TO 85
35      CONTINUE
      DO 40, I= 1,9
      IF(LINEBF(I) .EQ. '-')GO TO 45
40      CONTINUE
      GO TO 30
C
45      IF(LINEBF(I+1) .NE. '-')GO TO 300
      LFG = 1           !Two successive dashes; switch modes.
      IF (INTFLG .NE. 2) GOTO 30
      TYPE 1045           !Running CALDET
1045     FORMAT (' Use this spectrum to calibrate the detector?',
      +         '(Y/N): ',#)
      ACCEPT 1046, M

```

```

1046  FORMAT (A1)
      IF (M .NE. 'Y') GOTO 320
      GO TO 30

C      Test for "END" flag beginning in cols 1-3
      DO 60, J=1,3
      IF(LINEBF(J) .NE. 'E')GO TO 60
      IF(LINEBF(J+1) .NE. 'N')GO TO 60
      IF(LINEBF(J+2) .EQ. 'D')GO TO 320          !End of spectrum
60     CONTINUE
      GO TO 250

85     IF(LINEBF(I+1) .NE. ' ')STOP ' SYNTAX ERR'
      M = I+2
C      ITEMCT = ITEMCT + 1
      KHKT = NCHRS - M + 1

      GO TO (100,120,130,140,150,160,170,180,185,190,200,210,
+      220)ITEMCT

100    IF(KHKT .GT. MXFNM-1)KHKT = MXFNM-1
C      DO 102 I = 1,KHKT
C102   IFNAME(I) = LINEBF(M+I-1)
      GO TO 30

120    DO 122 I=1,12
      EXPTID(I) = LINEBF(M)
      M = M + 1
122    CONTINUE
      GO TO 30

)      DO 132, I= 1,12
      SAMPID(I) = LINEBF(M+I-1)
132    CONTINUE
      GO TO 30

140    CONTINUE
C140   DO 142, I=1,6
C      DSNAME(I) = LINEBF(M+I-1)
C142   CONTINUE
      GO TO 30

150    DECODE(KHKT,1150,LINEBF(M),ERR=300)DSTTIM,ICTYR
1150   FORMAT(E21.13,I6)
      GO TO 30

160    DECODE(KHKT,1150,LINEBF(M))DTZERO,IZYR
      GO TO 30

170    DECODE(KHKT,1170,LINEBF(M))DELTAT
1170   FORMAT(E15.7)
      GO TO 30

180    DECODE(KHKT,1170,LINEBF(M))TIMEL
      GO TO 30

3      DECODE(KHKT,1170,LINEBF(M))SMPWT
      GO TO 30

190    DECODE(KHKT,1170,LINEBF(M))WT100

```

GO TO 30

200 DECODE(KHKT,1170,LINEBF(M))ANSFCT

GO TO 30

210 DECODE(KHKT,1210,LINEBF(M))PRPK(1),PRPK(2)

10 FORMAT(2E15.7)

GO TO 30

220 IF(KHKT .GT. MXFNM-1)KHKT = MXFNM-1

C DO 222 I = 1,KHKT

C222 LIBNAM(I) = LINEBF(M+I-1)

GO TO 30

C Read a line of peak data.

250 DECODE(6,1250,LINEBF)II

1250 FORMAT(I6)

IF(II .EQ. 0)GO TO 30

IG = NGAMS+1

IF(IG .GT. MXPKS)GO TO 305

DECODE(NCHRS,1252,LINEBF)I,CHANL,ENGY(IG),PKS(IG),PKE(IG),

+ COUNT,BKSMBL,GAMS(IG),FRCDER(IG)

1252 FORMAT(I7,5F13.6,A2,2E15.7)

C NGAMS = IG

GO TO 30

300 TYPE 1305,IG

1305 FORMAT(' Reached peak # ',I3)

STOP ' ERROR READING DATA FILE'

305 STOP ' Too many peaks!'

) NDFIL = 1 !EOF seen.

320 RETURN

END

SUBROUTINE RDHDR

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C      Access SAMPLE file to get sample parameters, etc.
C      Last edit for RSX-11 15-Jan-88 R.G.
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
L
      INCLUDE 'GRPLGCM.PMS'
      INCLUDE 'GRPGMS.PMS'
      INCLUDE 'CASSAY.CMN'
C
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C
      LOGICAL*1 YESNO
      INTEGER LOAD(2)
      REAL YY(2)
      EQUIVALENCE (LOAD,TEMPR)
      BYTE TDSN(6)
      EQUIVALENCE(Y(140),TDSN)
C      IF(IOFLG(1) .EQ. 0)GO TO 20
      DO 15, J=1,6
15      TDSN(J) = DSNAME(J)
C
C      Get detector-system label from spectral header
C20      DO 21, I=1,6
C21      DSNAME(I) = ' '
C      N = 0
C      DO 22, K = 1,6
C      IF (HDRBYT(90+K) .EQ. ' ')GO TO 22
C      N = N+1
C      DSNAME(N) = HDRBYT(90+K)
C22      CONTINUE
C      IF (N .EQ. 0)STOP ' BAD HEADER FILE--NO DET-SYS LABEL'
C
C      DO 30, I=1,6
C      IF(DSNAME(I) .NE. TDSN(I))GO TO 35
C30      CONTINUE
C      GO TO 40
C35      TYPE 1035,DSNAME,TDSN
C1035      FORMAT(' Warning!! DET-SYS label found in header is: ',6A1/
C      + ' which disagrees with specification:',6A1/' Continue? ',#)
C      IF(.NOT. YESNO())STOP
20      CONTINUE
C      TYPE 1036
C1036      FORMAT (' Detector/System Label (Up to 6 chars.): ',#)
C      ACCEPT 1037, DSNAME
C1037      FORMAT (6A)
C
c40      CALL GETIME (INPDEV,LIVTIM,RLTIME,DEDTIM)
      RLTIME = TIMECL
      LIVTIM = TIMELV
C      Pick up stop time for printouts
      ISTRYR = IDBUF(87)
      DSTPTM = IDBUF(88)
      DSTPTM = DSTPTM + RLHDR(45)
LLCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
      GAIN = RLHDR(38)
      ZERO = RLHDR(39)

```

```
DO 110, I=1,12
EXPTID(I) = HDRBYT(I)
SMPID(I) = HDRBYT(I+12)
CONTINUE
```

```
RLTIME = RLTIME/60.
TIMEL = LIVTIM/60.
IF(RLTIME .EQ. 0)GO TO 125
DEDTIM = 100. * (RLTIME-TIMEL)/RLTIME
ICTYR = IDBUF(83)
DSTTIM = IDBUF(84)
DSTTIM = DSTTIM + RLHDR(43)
```

```
IZYR = IDBUF(79)
DTZERO = IDBUF(80)
DTZERO = RLHDR(41) + DTZERO
IF(DTZERO .EQ. 0)GO TO 150
```

```
DELTAT = DSTTIM - DTZERO + RLTIME/2880.
IDEL = ICTYR - IZYR
IF(IDEL .EQ. 0)GO TO 150
DO 140, I=1,IDEL
DELTAT = DELTAT + 365.
J = IZYR + I - 1
IF(0.25*FLOAT(J)-FLOAT(J/4) .EQ. 0)DELTAT = DELTAT+1.
CONTINUE
```

```
IF(RLHDR(27) .NE. 0)GEOM = (RLHDR(27))
```

```
SURFC = RLHDR(10)
SMPWT = RLHDR(11)
RHO = RLHDR(12)
ANSFCT= RLHDR(13) !Normalization factor
WT100 = RLHDR(14)
VOL = RLHDR(15)
```

```
IF(SMPWT .EQ. 0)GO TO 160
IF(RHO .NE. 0)GO TO 155
IF(VOL .NE. 0)RHO = SMPWT/VOL
IF(RHO .EQ. 0)RHO = 1.
IF (SURFC .NE. 0)DEPTH=SMPWT/(SURFC*RHO)
```

```
NF1 = 0
DO 165, J=1,5
IJ = 31+3*(J-1)
MTLZ(J) = IDBUF(IJ)
LOAD(1) = IDBUF(IJ+1)
LOAD(2) = IDBUF(IJ+2)
IF(TEMPR .EQ. 0)GO TO 165
NF1=NF1+1
CMPOS(J) = TEMPR
```

```
NF2 = 0
DO 180, J=1,6
IJ = 57+3*(J-1)
NABS(J) = IDBUF(IJ)
LOAD(1) = IDBUF(IJ+1)
LOAD(2) = IDBUF(IJ+2)
IF(TEMPR .NE. 0.)NF2 = NF2+1
```

```
180  ABSRB(J) = TEMPR  
      DO 185, I = 1,32  
185  MCHBUF(I) = HDRBYT(270+I)  
      DO 190, I = 1,70  
      )  TXTBUF(I) = HDRBYT(I+192)
```

```
      RETURN  
      END
```



SUBROUTINE READ38

READIN subroutine modified for SRP Pu238 system  
Modified by WDR 30-Aug-89  
Modified for WHC by WDR 10-May-91

CC

This sub reads files required by GRPANL. The control and spectral file names are always needed, but other input may be defaulted, thus permitting batch runs of GRPANL. The control file is in the form of preloaded common blocks prepared by the separate program, EDIGRP. This file defines all peak groupings to be processed and also stores some control parameters for the analysis. Run-time options may be used with the spectral file name to control output & analysis features. Any options given are processed, the spectral file is opened and its header read before returning. RDGRP is called later from the root to read successive copies of the common block GPKS for individual processing. This file organization permits operation of the spectral analysis program GRPANL under the control of indirect files containing only source file names, if desired. This capability has been expanded to allow use of "generic" control files for which the detector-system and many other parameters are not defined until GRPANL run time.

Link with UTLIB.  
LUNs INPDEV & INDEV are defined in the root.

By J.B.Niday Last edit 5-Dec-86

CC

The following common is specific to Pu238 instrument

INCLUDE 'CFILES.CMN'  
INCLUDE 'CLUNS.CMN'  
INCLUDE 'CONTRL.CMN'  
INCLUDE 'CASSAY.CMN'  
INCLUDE 'CDATE.CMN'  
INCLUDE 'CBACKG.CMN'

The following common is specific to analysis routines

INCLUDE 'GRPLGCM.PMS'  
INCLUDE 'GRPGMS.PMS'  
  
PARAMETER MXOPTS=12  
LOGICAL\*1 CONFRM,YESNO  
BYTE FILE(18),EFILE(18),FFILE(18),ISFIL(40),CTLFIL(MXFNM)  
BYTE OPTS(MXOPTS),LOPT(MXOPTS) !Permanent & working option flag arrays  
DIMENSION IWR2(2)  
EQUIVALENCE (IWR2,LIVTIM) !Note: LIVTIM is REAL\*4  
  
BYTE DETFNM(18),SHPFNM(18),DEFNAM(18),LCLDEF(18),XTN(12)  
EQUIVALENCE (DETFNM,RM),(SHPFNM,RM(MXFNM+1))  
Save names for INITLP (files actually used for updates).  
BYTE IBUF(80)  
EQUIVALENCE (WM,IBUF) !Temp save of file input & options  
EQUIVALENCE (WM(30),LTH)

```

DATA FILE/'P','U','2','3','8','.',',','C','T','L',9*'/
DATA FFILE/'P','U','2','3','8','.',',','F','C','C',9*'/
DATA EFILE/'P','U','2','3','8','.',',','E','P','6',9*'/
DATA DETFNM/'P','U','2','3','8','.',',','D','E','T',9*0/
DATA SHPFNM/'P','U','2','3','8','.',',','S','H','P',9*0/
DATA OPTS/'I','C','S','E','P','R','G','W','F','U','L','B'/
DATA XTN/','S','H','P','.',',','D','E','T','.',',','B','K','G'/
DATA DEFNAM/'P','U','B',',',',13*'/,0/
DATA LCLDEF/17*'/,0/
DATA LOPT/'I','F',10*0/
CALL VERS(0)
ISFIL(1) = ' '

```

```

C TYPE 1001

```

```

1001 FORMAT(' GRPANL can be used for any or all of: '//
+ ' 1) Fitting peaks in peak grouping(s) '//
+ ' 2) Converting the counts to photons/min '//
+ ' 3) Storing these results in a cumulative file, and/or '//
+ ' Directly evaluating the results for probable source(s). '//
+ ' For this it requests a control file defining the peak groups',
+ ' desired' // and some other parametric data. '//
+ ' This file is prepared by the EDIGRP program ')

```

```

C Common blocks must be preloaded by EDIGRP

```

```

DO 5, I=1,LENLG
5 NWRDS(I) = 0
DO 7, I=1,LENGSP
7 IWR2(I) = 0
DO 8, I=1,2*MXNP+7
VM(I) = 0.
ERR(I) = 0.

```

```

C10 TYPE 1009

```

```

1009 FORMAT('/ Enter specs. of pre-edited GRPANL control file. '//
+ ' For fully interactive control over this run, append a "/I" '//
+ ' and a /F if a Foreign header must be processed',
+ ' (e.g., LANL files.) '//)
10 NKH = MXFNM-1
NOPT = 2

```

```

C Read default control file PU238.CTL

```

```

NKH = 0
IF(NKH .NE. 0)GO TO 30
DO 20 I = 1,18
20 CTLFIL(I) = FILE(I)
CTLFIL(19) = 0
30 OPEN(UNIT=INDEV,STATUS='OLD',READONLY,FILE=CTLFIL,
+ FORM='UNFORMATTED',ERR=35)
GO TO 40
35 TYPE 1035,(CTLFIL(I),I=1,40)
1035 FORMAT(' OPEN error on file: ',20A1)
ABORT = .TRUE.
GOTO 400
C GO TO 10
READ(INDEV,ERR=55)LNLG,LNGS,LNPK
IF(LNPK .EQ. LENPKS)GO TO 41
IF(LNPK .NE. 308)GO TO 45
41 IF(LNGS .NE. LENGSP)GO TO 45

```

```

IF(LNLG .EQ. LENLG)GO TO 50
45  TYPE 1045,LNLG,LNGS,LNPK
1045 FORMAT(' INPUT FILE BLOCKS WRONG LENGTH -USE EDIGRP',3(', ',14))
CLOSE(INDEV)
ABORT = .TRUE.
GO TO 300

50  READ(INDEV,ERR=55)(NWRDS(I),I=1,LNLG)
    READ(INDEV,ERR=55)(IWR2(I),I=1,LNGS)

C   Keep INDEV open for reading peak groups in RDGRP section.
KTGRP = 0
GO TO 60

55  TYPE 1055,INDEV,(CTLFIL(I),I=1,MXFNM-1)
    PAUSE ' READ FAILED'
1055 FORMAT(' INDEV= ',I2,', SPEC= ',29A1)
CLOSE(INDEV)
GO TO 10

60  CONTINUE

C   C   C   C   C   C   C   C   C   C
C   Open spectral file for direct access,
C   read its header, and report some contents.

90  NOPT = MXOPTS
DO 92, I = 1, MXOPTS      !Load options here in case of retrys
92  LOPT(I) = OPTS(I)

C   TYPE 1095
1095 FORMAT(/' Now enter file specs for the spectral data: ',#)

    NKH = MXFNM-1

100  OPEN(UNIT=INPDEV,NAME=filnam,STATUS='OLD',ACCESS='DIRECT',
+     READONLY,RECORDSIZE=128,ERR=105)
GO TO 110

105  CLOSE(INPDEV)
    TYPE 1105, filnam
1105 FORMAT(' Cannot open spectrum: ',A15)
ABORT = .TRUE.
RETURN

C
C ... Read trailing block of data file to get sample and operator IDs
C and measurement time and date information
110  irec=33
    READ(inpdev'irec)DUM,BK,SGBK,TIMEBK,BKDAT,BKTIM,BKOP,BKSFIL,
X   SAMPID,OPERID,SYSID,ASDATE,ASTIME,MATTYP(1)

C   TYPE *,'SDAT',ASDATE,ASTIME,SAMPID,OPERID
C
C Read default sample information file PU238.FCC
IF (ISFIL(1) .EQ. ' ') THEN
    DO 62 I = 1,18
62  ISFIL(I) = FFILE(I)
C   TYPE 9999,MATTYP(1)
C   99  FORMAT(1X,'MATTYP',A1)
    IF (MATTYP(1).EQ.'E')THEN
        DO 63 I = 1,18
63  ISFIL(I) = EFILE(I)

```

```

ENDIF
ISFIL(19) = 0
ENDIF
OPEN (UNIT=NDEV,NAME=ISFIL,TYPE='OLD',ACCESS='DIRECT',
+ READONLY,RECORDSIZE=128,ERR=64)
GOTO 66
CALL CLOSE (NDEV)

TYPE 1058,(ISFIL(I), I=1,40)
1058 FORMAT (' Cannot open SAMPLE file ',40A1)
ABORT = .TRUE.
GOTO 300
66 CALL RDBLK(256,IDBUF,1,NDEV)
CALL CLOSE (NDEV)
C C C C C C C C C
C Update pre-edited control flags vs. any option flags used.
C First check for nature of header

IF(NOPT .EQ. 0) GO TO 120

IF(RUNFLG(2) .EQ. 0)RUNFLG(2) = LOPT(9) !Test for '/F' again.

C
C
112 IF(LOPT(5) .GT. 0) IOFLG(3) = 0
IF(LOPT(5) .LT. 0) IOFLG(3) = 1
IF(LOPT(6) .GT. 0) IOFLG(4) = 0
IF(LOPT(6) .LT. 0) IOFLG(4) = 1
C IOFLG 3 & 4 are set to inhibit printing &/or residuals tables.

C Pick up nature & extent of analysis:
IF(LOPT(2) .LT. 0)RUNFLG(4) = 0
IF(LOPT(2) .GT. 0)RUNFLG(4) = 1
IF(LOPT(4) .LT. 0)RUNFLG(5) = 0
IF(LOPT(4) .GT. 0)RUNFLG(5) = 1
IF(LOPT(3) .EQ. 0)GO TO 117
IF(LOPT(3) .LT. 0)THEN
STGFIL(1) = ' '
GO TO 117
ENDIF
C
115 IF((STGFIL(1) .NE. 0).AND.(STGFIL(1) .NE. ' '))GO TO 117
TYPE 1114
1114 FORMAT(' Please enter a file name for storage of results ',#)
ACCEPT 1150,(STGFIL(I),I=1,MXFNM-1)

117 IF (LOPT(7) .EQ. 0) GO TO 118
IF(LOPT(7) .GT. 0) THEN
IOFLG(5) = 1
GO TO 120
ELSE
IF(IOFLG(5) .EQ. 1) IOFLG(5) = 0
ENDIF
118 IF((LOPT(8) .NE. 0).AND.(IOFLG(5) .NE. 1))IOFLG(5)=LOPT(8) * 2

IF(IOFLG(5) .LT. 0)IOFLG(5) = 0

-----
C Pick up spectral header info including GAIN for remaining initialization.
120 CALL RDHDR !Read LLNL-N.C.div. standard header

```

```

C      IF (GAIN.EQ.0) THEN
      TYPE 1089
1089  FORMAT(' Spectrum GAIN (keV/ch) = ',#)
      ACCEPT 1145,GAIN
      ENDIF

C
      IF((GEOM .EQ. 0) .AND. (RUNFLG(4) .NE. 0)) THEN
      TYPE 1122
1122  FORMAT(' Geometry is not given in this header but is',
+ ' required for efficiency corrections.'/' Please enter',
+ ' distance (cm.) from sample midplane to detector window: ',#)
      ACCEPT 1145, GEOM
      IF(GEOM .EQ. 0)GO TO 60          !Another try at option setting

      ENDIF

```

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C      Time to update parameters from standard files or data base,
C      unless the update option was negated (/U e.g.).
C      First priority is to look for the system default files
C      whose names are based on the DS name found in the header.
C      If no files, look for any previously stored DFILNM &/or SFILNM
C      (provided the control file was made up for the same DS).
C      Third priority is to look for locally stored default names.
C      If still no data, or if request option switch set, accept
C      entry of file names form keyboard.
C

```

```

      KK = 4 !VAX default files in N.C. Div. stored in "PUB:"
      NCDS = 0
      DO 124, I = 1,6
      IF(DSNAME(I) .EQ. ' ')GO TO 126
      NCDS = NCDS + 1
124  CONTINUE

C      Now build public & local default names for updates.
126  DO 128, I = 1,NCDS
      LCLDEF(I) = DSNAME(I)
128  DEFNAM(I+KK) = DSNAME(I)

      JDS = KK + NCDS !Remember length of default name (less XTN)

      DO 130, I = 1,4
      LCLDEF(I+NCDS) = XTN(I+4)
130  DEFNAM(I+JDS) = XTN(I+4)
      IF(RUNFLG(3) .LE. 0)GO TO 134 !Flag set by RDHDR for "new" D-S
      DFILNM(1) = 0
      SFILNM(1) = 0
      BKGDFL(1) = 0
      DO 132,I = 1,24
132  DCNST(I) = 0.

134  IF(LOPT(10) .LT. 0) THEN
      IF(RUNFLG(3) .LE. 0) GO TO 250
      O.K. to skip update when requested, provided no DS change was made.
      TYPE 1340
1340  FORMAT(' Detector-system change normally requires parameter',
+ ' updates.'/' Proceed without updating? (Y/N) ',#)

```

```

                IF(YESNO()) GO TO 250
ENDIF
145  OPEN(UNIT=NDEV,FILE=DETFNM,STATUS='OLD',READONLY,ERR=147)
      READ(NDEV,1145,END=146,ERR=147)(DCNST(I),I=1,24)
1145  FORMAT(E15.7)
146  CLOSE(NDEV)
      GO TO 165
147  CLOSE(NDEV)      !Error on OPEN/READ of DETFNM

      TYPE 1480,(DETFNM(I),I=1,MXFNM-1)
1480  FORMAT(' Cannot open/read parameter file: ',29A1)
      DO 149, I = 1,MXFNM-1
149  DETFNM(I) = ' '

      IF(RUNFLG(4) .EQ. 0) GO TO 170  !Counts only.

150  TYPE 1492
1492  FORMAT(' Please enter file specs for detector constants')
      ACCEPT 1150,(DETFNM(I),I=1,MXFNM-1)
1150  FORMAT(29A1)
      IF(DETFNM(1) .NE. ' ') GO TO 145
      IF(RUNFLG(3) .LE. 0) GO TO 165
      TYPE 1340
      IF(.NOT. YESNO()) GO TO 150
      GO TO 165

)
162  DO 162, I=1,18
      DETFNM(I) = LCLDEF(I)      !Try for a local file next
      DO 163, I = 18,MXFNM-1
163  DETFNM(I) = ' '
      KK = 0
      GO TO 145

C      Calculate approx. detector dead layer
165  CALL CALEFF(DCNST,60.,-2.813,EFF1)
      CALL CALEFF(DCNST,70.,-2.659,EFF2)
      DCNST(16) = (EFF2 - EFF1) * 1.65
      IF(DCNST(16) .LT. 0)DCNST(16) = 0.00000001
C      TYPE 1165,DCNST(16)
1165  FORMAT(' Ge dead layer calculated to be:',1PE10.2)

C      C      C      C      C      C      C      C
170  KK = 4      !Reset for public files flag.
      DO 171, I = 1,4
      LCLDEF(I+NCDS) = XTN(I)
171  DEFNAM(I+JDS) = XTN(I)  !Changes default to .SHP

C      IF (LOPT(10) .GT. 0) THEN      !/U flag entered.
C          IF(RUNFLG(1) .GT. 0) GO TO 204  !Interactive.
C          GO TO 190
      ENDIF
C180  DO 182, I = 1,JDS+4
C182  SHPFNM(I) = DEFNAM(I)
C      DO 185, I = JDS+5,MXFNM-1

```

```

C185 SHPFNM(I) = ' '
190 IF (SFILNM(1).NE.' ')THEN
DO 192, I = 1,MXFNM
192 SHPFNM(I) = SFILNM(I)
ENDIF
KK = 2

L
200 OPEN(UNIT=NDEV,FILE=SHPFNM,STATUS='OLD',READONLY,ERR=203)
I = 0
201 I = I + 1
READ(NDEV,1145,END=202,ERR=203)SHAPC(I)
IF (I .LT. 15) GO TO 201 !Protect from oversize files.
202 CLOSE(NDEV)
GO TO 250

203 TYPE 1480,(SHPFNM(I),I=1,MXFNM-1)
CLOSE(NDEV)
IF(KK .EQ. 4) THEN
IF((SFILNM(1).NE.O).AND.(SFILNM(1).NE.' '))GO TO 190
ENDIF
IF(KK .NE. 0) GO TO 220
SHPFNM(1) = 0
IF(RUNFLG(1) .LT. 0) GO TO 250

204 TYPE 1204
1204 FORMAT(' Please enter file specs for shape parameters, if any')
ACCEPT 1150,(SHPFNM(I),I=1,MXFNM-1)
IF(SHPFNM(1) .NE. ' ')GO TO 200
IF(RUNFLG(3) .EQ. 0)GO TO 250 !If = 0: not changing DS
DO 205, I = 1,15
205 SHAPC(I) = 0.
IF(KK .EQ. 0) GO TO 250

220 DO 222, I = 1,17
222 SHPFNM(I) = LCLDEF (I)
DO 223, I = 18,MXFNM-1
223 SHPFNM(I) = ' '
KK = 0
GO TO 200

C The following files could not have been read in by EDIGRP.
250 KK = 4 !Reset for "PUB:" again.

C IF(LDPT(12) .GT. 0)GO TO 260
C IF(RUNFLG(3) .EQ. 0) THEN
C IF((BKGDFL(1) .NE. 0) .AND. (BKGDFL(1) .NE. ' '))GO TO 265
C ENDIF

C252 DO 253, I = 1,4
C LCLDEF(I+NCDS) = XTN(I+8)
C253 DEFNAM(I+JDS) = XTN(I+8)
C DO 254, I = 1,17
C254 BKGDFL(I) = DEFNAM(I)
C DO 255, I = 18,MXFNM-1
C255 BKGDFL(I) = ' '
KK = 0 !Flag: Default in public has been tried.
L GO TO 265
C
C256 CLOSE(NDEV)

```

```

C      TYPE 1256,(BKGDFL(I),I=1,MXFNM-1)
C1256  FORMAT(/' Cannot read file of background peaks:',29A1)
C      IF(KK .GT. 0) GO TO 252          !Try again with default name.
C      IF(KK .EQ. 0) THEN
C      DO 257, I = 1,17
C17    BKGDFL(I) = LCLDEF(I)
C      KK = -1
C      GO TO 265
C      ENDIF
C      DO 259, I = 1,MXFNM-1
C259   BKGDFL(I) = ' '
C      IF(RUNFLG(1) .LE. 0) GO TO 270

C260   TYPE 1260
1260   FORMAT(' Please enter file specs for bkgd peak counts in',
+       ' detector, if available.'/' ',#)
C      ACCEPT 1150,(BKGDFL(I),I=1,MXFNM-1)
C      IF(BKGDFL(1) .EQ. ' ')GO TO 270

C265   OPEN(UNIT=NDEV,FILE=BKGDFL,STATUS='OLD',READONLY,ERR=256)

C      J = 1
C266   READ(NDEV,1266,ERR=256,END=268)BKPKEN(J),BKPCKT(J)
1266   FORMAT(2E15.6)
C      J = J+1
C      IF((J .LE. MXBKP) .AND. (BKPKEN(J-1) .NE. 0)) GO TO 266
C268   CLOSE(NDEV)
C      NBKPKS = J - 1

C270   IF(RUNFLG(1) .LT. 0) GO TO 277          !Not interactive
C      IF(RUNFLG(3) .EQ. 2) GO TO 273        !DS changed
C      IF((STGFIL(1) .EQ. ' ').OR.(STGFIL(1) .EQ. 0)) THEN
C          IF(RUNFLG(5) .EQ. 0) GO TO 277    !Libe not needed
C      ENDIF
C      IF((LIBNAM(1).NE. ' ').AND.(LIBNAM(1).NE.0))GO TO 277

273   TYPE 1273
1273   FORMAT(' A library will be needed for any evaluation.'/
+       15X,'Please enter file specs here :',#)
C      ACCEPT 1150,(LIBNAM(J),J=1,MXFNM-1)

277   CALL DATE(DTE)
C      CALL TIME(TBUF)

C      APRGN = GAIN
C      APRZRO= ZERO          !Preserve original values
C      Transfer temporary control parameter file name from local
C      to common for output.
C      DO 278, I=1,MXFNM
278   ICFIL(I) = CTLFIL(I)

C      IF(LONGPR)CALL INITLP

C      Read in "absorption" coefficients, when calculating photons.
C      This destroys file names of detector consts, shape parameters
C      and absorption coefficients stored by EDIBRP.

280   OPEN(UNIT=NDEV,FILE=AFILNM,STATUS='OLD',READONLY,ERR=282)

```



```

      READ(NDEV,1280,END=285,ERR=282)(NSMBL(J),EDGE(J),(COFMU(I,J),
1      I=1,8),J=1,20)
1280  FORMAT(1X,A2,F7.3,4F12.9/10X,4F12.9)
      GO TO 285
282  TYPE 1282;(AFILNM(I),I=1,MXFNM)
32  FORMAT(' Error reading: ',30A1)
      IF(NSMBL(19) .NE. 0)GO TO 285 !NSMBL(1) was used for AFILNM.
      CLOSE(NDEV)
      IF(RUNFLG(1) .LT. 0) GO TO 294
      TYPE 1283
1283  FORMAT(' Enter file specs for absorption coeffs ',*)
      ACCEPT 1150,(AFILNM(I),I=1,29)
      IF(AFILNM(1) .EQ. ' ')GO TO 294
      GO TO 280

285  CLOSE(NDEV)

C *** Adjust certain variables after printing them out in INITLP
294  IF(GEOM .NE. 0)GEOM = GEOM + DCNST(15)
      IF(NF1 .EQ. 0)GO TO 300
      DO 295 K=1,NF1
295  CMPOS(K)=0.005*CMPOS(K)*SMPWT / SURFC

300  RETURN

```

```

CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
      ENTRY RDGRP

```

```

C      This section is called repeatedly for reading successive blocks
C of group specifications (Common block /GPKS/)

```

```

      READ(INDEV,ERR=390,END=395)(NXFLG(I),I=1,LNPK)
      IF(EST .NE. 0)GO TO 400
390  TYPE 1390,KTGRP+1
1390  FORMAT('/' ERROR TRYING TO READ GROUP',I3)
395  NP = 0
      EST = 0
      CLOSE (INDEV)

400  RETURN
      END

```

```

C*****
C READB
C
C     Read background from disk
C
C     Rev. 25-Jun-84 lc
C     Rev. 04-Apr-86 WDR           Modified for Micro VAX II
C*****
C
C     SUBROUTINE READB
C
C     LOGICAL*1 FILL
C
C     The following common block files are generic for all instruments
C
C     INCLUDE 'CFILES.CMN'
C     INCLUDE 'CBACKG.CMN'
C     INCLUDE 'CLUNS.CMN'
C     INCLUDE 'CPASS.CMN'
C     INCLUDE 'CONTRL.CMN'
C
C     The following common block files are instrument specific for
C     gamma-ray instruments
C
C     INCLUDE 'CONT90.CMN'
C
C     If FROMDK, read background from spectrum file
C
C     IF(FROMDK) GO TO 10
C
C     Notify operator
C
C 1     TYPE 1003,BKFIL
1003   FORMAT(/' Reading ',A)
C
C     OPEN A DIRECT ACCESS FILE WITH NAME BKFIL
C
C     OPEN(UNIT=DKLUN,FILE=BKFIL,ACCESS='DIRECT',TYPE='OLD',
C       X      RECORDSIZE=128,ERR=999)
C
C     Read data
C
C     READ(DKLUN'1,ERR=999)BKDAT,BKTIM,BKOP,FILL,
C       X      BK,SGBK,TIMEBK,PASSWD,BKSFIL
C
C     Close file and acknowledge
C
C     CLOSE(UNIT=DKLUN)
C     IF(BEGIN)GO TO 500
C     TYPE 1000,BKFIL,BKDAT,BKTIM,BKOP
1000   FORMAT(/,1X,A,' was written ',A,2X,A,' by ',
C       X      A12)
C     GO TO 500
C
C     Read background file from spectrum file (last record)
C
C 10    TYPE 1003,FILNAM
C     OPEN(UNIT=DKLUN,NAME=FILNAM,ACCESS='DIRECT',TYPE='OLD',
C       X      RECORDSIZE=128,ERR=997)

```

```
NRECS = MEMSIZ/128
NREC = NRECS + 1
```

```
C
C Check if first entry in record is dummy
```

```
C
  READ(DKLUN,NREC,ERR=998)DUM
  IF(DUM .NE. -999.0) GO TO 998
  READ(DKLUN,NREC,ERR=998)DUM,BK,SGBK,TIMEBK,BKDAT,BKTIM,
  X      BKOP,BKSFIL
```

```
11  CLOSE(UNIT=DKLUN,ERR=11)
     TYPE 1000,FILNAM,BKDAT,BKTIM,BKOP
     GO TO 500
```

```
C
C Error in reading from spectrum file
C Do not abort, but use most recent background read
```

```
C
997  CALL ERRMSG(1,FILNAM)
     GO TO 500
998  CALL ERRMSG(2,FILNAM)
     CLOSE(UNIT=DKLUN)
     GO TO 500
```

```
C
C Error in reading
```

```
C
999  ABORT = .TRUE.
     CALL ERRMSG(1,BKSFIL)
     CALL ERRMSG(102,ARG)
```

```
C
500  RETURN
     END
```

```

C*****
C
C READC38      (READC subroutine for PU238 isotopic instrument)
C
C      Read constants for Pu-238 isotopic plutonium gamma-ray
C      analysis instrument from disk
C
C      REV. 14-JUN-84  1c
C      REV. 06-SEP-84  CMS
C      REV. 17-DEC-84  WDR      Modified for solid isotopic instrument
C      REV. 04-APR-86  WDR      Modified for Micro VAX II
C      REV. 13-JUL-87  RDP      Modified for uVAX II
C                               Convert BYTE strings to CHARACTER strings
C      Rev. 24-May-91  WDR      Modified for stabilizer settings

```

```

C*****

```

```

C      SUBROUTINE READC
C
C      LOGICAL*1 FILL,CDUM(30)
C      CHARACTER SNAM*12,SYN*3
C
C      INTEGER*2 EMC(21),EDD(41),FILLR(58)
C
C      These common block files are generic for all instruments
C
C      INCLUDE 'CFILES.CMN'
C      INCLUDE 'CONFIL.CMN'
C      INCLUDE 'CLUNS.CMN'
C      INCLUDE 'CMCONS.CMN'
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CCAL.CMN'
C      INCLUDE 'CDIAG.CMN'

```

```

C      The following common block files are instrument specific for
C      solid isotopic analysis (gamma ray) instruments
C
C      INCLUDE 'CONT90.CMN'
C      INCLUDE 'P238COM.CMN'

```

```

C      Specify names of constants files to be read in
C
C      DIMENSION DSHAPC(15)
C      DATA SNAM /'SHAPC.  /
C      DATA DSHAPC/.04,.00175,-4.,.005,3.,0.,.003,.3,.0,.0,
C      +           3.,.0,.0,.0,.0/

```

```

C      Equivalence dummy arrays to first word of common block
C      Array size = word count of common.
C
C      EQUIVALENCE (EMC(1),BITIME),(EDD(1),DIAGMX)

```

```

C      Notify operator
C
C      TYPE 1003,CONFIL,COMMNT
C      1003  FORMAT(/' Reading ',A,1X,A)
C
C      Open a direct access file with name filnam

```

```

C
C   OPEN(UNIT=DKLUN,NAME=CONFIL,ACCESS='DIRECT',TYPE='OLD',
X     RECORDSIZE=128,ERR=999)
C
C   Read data
C   Fill = Logical*1 fill to even up word boundary
C   Cdum = Dummy fill argument to blank out comment
C   Fillr = Filler for more parameter arguments
C
C     READ(DKLUN*1,ERR=999)FILDAT,FILTIM,FILOP,FILL,CDUM,
X       FILLR,EMC,EDD,PU240,PU242,NPFLG,NUFLG,TIMEPR,
X       IZPEAK,IZWNDW,ZARNG,IGPEAK,
X       IGWNDW,GARNG
C
C   Close file and acknowledge
C
C     CLOSE(UNIT=DKLUN)
C     IF(BEGIN)GO TO 505
C     TYPE 1000,CONFIL,FILDAT,FILTIM,FILOP
1000  FORMAT(/1X,A,' was written ',A,2X,A,' by ',A)
C
C   Read in constants files specific to solid isotopic plutonium
C   gamma-ray analysis
C
C   First get system id to establish file extension
C
505   DO 100 J = 1,15
      IF (CONFIL(J:J) .EQ. '.')GO TO 120
100   CONTINUE
120   K = J
      DO 150 I = 1,3
      SYN(I:I)= CONFIL(K+I:K+I)
150   CONTINUE
C
C   Read in peak shape parameters
C
C     DO 200 J = 1,3
C200  SNAM(J+6:J+6)=SYN(J:J)
C     OPEN (UNIT=DKLUN,NAME=SNAM,TYPE='OLD',ACCESS='SEQUENTIAL',ERR=995)
C     READ (DKLUN,9087,ERR=995) (OSHAPC(I), I = 1,15)
C9087  FORMAT (E15.7)
C     CALL CLOSE (DKLUN)
C     GO TO 500
C
C   Error in reading
C
995   CALL ERRMSB(1,SNAM)
      CALL ERRMSB(102,ARG)
      GO TO 500
999   CALL ERRMSB(1,CONFIL)
      CALL ERRMSB(102,ARG)
C
500   RETURN
      END

```

```

C*****
C READNM
C
C Allow operator to select or change a constants
C filename in a list of constants files
C
L LCS 27-FEB-81
C
C Rev. 16-Jun-84 lc
C Rev. 19-Feb-85 MPK Change ACCEPTS to READS with END
C and ERR exits
C Use YESNO to ask yes/no questions
C Rev. 04-Apr-86 WDR Modified for Micro VAX II
C Rev. 13-Jul-87 RDP Modified for uVAX II
C Convert BYTE strings to CHARACTER strings
C
C*****
C
C SUBROUTINE READNM
C
C LOGICAL*1 ANS,YESNO
C CHARACTER CFLNAM(5)*12,CCMMNT(5)*30
C
C The following common block files are generic for all instruments
C
C INCLUDE 'CFILES.CMN'
C INCLUDE 'CONFIL.CMN'
C INCLUDE 'CONTRL.CMN'
C INCLUDE 'CLUNS.CMN'
C
C
C NEWFIL = .TRUE.
C IF(BEGIN) GO TO 4
C
C Check if constants file read is still to be used
C
1 TYPE 1005, CONFIL,COMMNT
1005 FORMAT(/' The constants file read was: ',A15,1X,A30//
X '$Do you want to select another constants file? (Y/N) -> ')
READ (5,1050,END=1,ERR=1) ANS
IF(ANS .NE. 'Y' .AND. ANS .NE. 'N') GO TO 1
IF(ANS .NE. 'Y') NEWFIL = .FALSE.
IF(.NOT. NEWFIL) GO TO 900
C
C Open file to read names of 5 current constants files
C
4 OPEN(UNIT=DKLUN,FILE=NAMSFL,ACCESS='SEQUENTIAL',
X TYPE='OLD',ERR=5)
GO TO 80
C
C Problem in open of names file
C
5 CALL ERRMSG(1,NAMSFL)
CALL ERRMSG(102,ARG)
ABORT = .TRUE.
IF (ABORT) CALL EXIT
C
L Read names from file
C
80 DO 100 I=1,5

```

```

100      READ(DKLUN,1000) CFLNAM(I),CCMMNT(I)
1000     FORMAT(A12,A30)
C
C  If just beginning, choose file 1
C
      IF(BEGIN) JANS = 1
      IF(BEGIN) GO TO 280
C
C  Type out the current constants file names & select
C  appropriate file
C
150     TYPE 1010
1010    FORMAT(//' The current constants files are: '//)
      DO 200 I=1,5
200     TYPE 1011,I,CFLNAM(I),CCMMNT(I)
1011    FORMAT(3X,I1,3X,A12,4X,A30)
C
210     TYPE 1012
1012    FORMAT(//'#Enter constants file selected (1-5)',
      X   ' ' or 6 for new file -> ')
      READ(5,1020,END=150,ERR=150) JANS
1020    FORMAT(I1)
      IF(JANS.LT.1.OR.JANS.GT.6)GO TO 150
      IF(JANS .EQ. 6) GO TO 400
C
C  Confirm filename
C
250     TYPE 1030,CFLNAM(JANS)
1030    FORMAT(//'#Constants file = ',A12,' OK? (Y/N) -> ')
      READ(5,1050,END=250,ERR=250) ANS
      IF(ANS .EQ.'N') GO TO 210
      IF(ANS .NE.'Y') GO TO 250
C
C  Read filename & comment into common block
C
280     CONFIL=CFLNAM(JANS)
      CONFIL(11:12)=' '
350     COMMNT=CCMMNT(JANS)
      CLOSE(UNIT=DKLUN,ERR=351)
351     GO TO 700
C
C  Input new constants file name & comment
C  Attempt to open new file
C
400     TYPE 1040
1040    FORMAT(//'#Enter new constants filename -> ')
      READ (5,1050,END=400,ERR=400) CONFIL
1050    FORMAT(A30)
C
401     TYPE 1041
1041    FORMAT(//'#Enter comment (30 char) -> ')
      READ (5,1050,END=401,ERR=401) COMMNT
C
C  Check if new file is permanent & which file it replaces
C
      IF ( .NOT.YESNO('Is this file a permanent replacement')
1     ) GO TO 600
C
410     TYPE 1044
1044    FORMAT(//'#Which file does the new file replace? (1-5) -> ')

```

```
READ (5,1020,END=410,ERR=410) JANS
IF(JANS.LT.1.OR.JANS.GT.5) GO TO 410
```

```
C
420     CFLNAM(JANS)=CONFIL
421     CCOMMNT(JANS)=COMMNT
```

```
L Write new file name list to disk
```

```
C
C     REWIND DKLUN
CLOSE (UNIT=DKLUN)
OPEN(UNIT=DKLUN,NAME=NAMSFL,ACCESS='SEQUENTIAL',
X     TYPE='OLD',ERR=5)
DO 430 I=1,5
430     WRITE(DKLUN,1100) CFLNAM(I),CCMMNT(I)
1100     FORMAT(1X,A12,A30)
```

```
C Close names file
```

```
C
600     CLOSE(UNIT=DKLUN,ERR=700)
```

```
C Try opening the new constants file
```

```
C
700     OPEN(UNIT=DKLUN1,NAME=CONFIL,TYPE='OLD',ERR=999)
CLOSE(UNIT=DKLUN1,ERR=701)
701     GO TO 900
```

```
C
999     CALL ERRMSG(1,CONFIL)
GO TO 150
```

```
C
900     BEGIN = .FALSE.
RETURN
```

```
END
```



```

C*****
C  RESET38      (RESET subroutine for Pu238 isotopic instrument)
C
C  Resets default settings for Assay options
C  This subroutine is specific to the Pu238 isotopic instrument
C  CALLING ARGUMENTS:
C      IARG = 0   Default flags and parameters set
C      IARG = 1   Operator selects flags and parameters
C                  (supervisory option)
C
C  14-MAY-84 LC
C
C  Rev. 14-Jun-84 LC
C  Rev. 11-Sep-84 TES
C  Rev. 29-NOV-84 MPK      Replace list-directed with
C                          formatted I/O
C  Rev. 27-DEC-84 WDR      modified for the solid isotopic instrument
C  Rev. 19-FEB-85 MPK      Change ACCEPTS to READS with END and ERR
C                          exits
C  Rev. 01-Mar-85 WDR      modified to type values at the start
C                          and end only instead of after each entry
C  Rev. 07-MAR-85 MPK      Initialize the VAX comm channel when
C                          VAXCOM may have changed
C  Rev. 07-Apr-86 WDR      Modified for Micro VAX II
C  Rev. 13-Jul-87 RDP      Modified for uVAX II
C                          Convert BYTE strings to CHARACTER strings
C  Rev. 14-Apr-88 WDR      Removed raw data to VAX option
C  Rev. 27-Apr-90 WDR      Removed "VAX communication"
C*****

```

```

C
C  SUBROUTINE RESET(IARG)

```

```

C  LOGICAL*1 RANS

```

```

C  The following common block files are generic for all instruments

```

```

C  INCLUDE 'CFILES.CMN'
C  INCLUDE 'CONTRL.CMN'
C  INCLUDE 'CONFIL.CMN'
C  INCLUDE 'CASSAY.CMN'
C  INCLUDE 'CLUNS.CMN'
C  INCLUDE 'CCAL.CMN'

```

```

C  The following common block files are instrument specific for
C  (gamma ray) instruments and specific to the solid isotopic inst.

```

```

C  INCLUDE 'CSPECT.CMN'
C  INCLUDE 'CONT90.CMN'
C  INCLUDE 'P238COM.CMN'

```

```

C  IF(IARG .EQ. 1)GO TO 100

```

```

C  Compare name of last constants file read with the default
C  constant file name

```

```

C      IF (CONFIL.EQ.LSTFIL) GOTO 10
C      BEGIN = .TRUE.

```

```

10  CONTINUE

```

```

C If file names different, reread constants file
C (flags reset by constants file)
C
      ABORT = .FALSE.
      IF(.NOT. BEGIN) GO TO 20
      CALL READNM
      IF(ABORT) GO TO 500
      CALL READC
      IF(ABORT) GO TO 500
C
C Reread background file
C
20      CALL READB
      IF(ABORT) GO TO 500
C
C Reset parameters to default
C
      CALL INTFLG
C
C Load last file array with constants file read
C
      DO 40 I=1,15
40      LSTFIL(I:1) = CDNFIL(I:1)
C
      RETURN
C
C Supervisory option:
C Print out default options
C
*00      TYPE 1000, LONGPR, WRITNG, TIMEPR, PU240, PU242, CNDFIL
*00      FORMAT(/' 1. Long printout = ',L1/
X          /' 2. Store data with sample ID filename = ',L1/
X          /' 3. Preset assay time = ',F10.0/
X          /' 4. Pu-240 abundance (Wt.%) = ',F8.4/
X          /' 5. Pu-242 abundance (Wt.%) = ',F8.4/
X          /' 6. Constants file = ',A15)
      TYPE *,'
C
C Allow supervisor to change defaults
C
150      TYPE 1500
1500     FORMAT(/'#Enter number (1-6) to change (RETURN for no change) ->
X          )
      READ (5,1501,END=150,ERR=150) IOPT
1501     FORMAT(I1)
      IF(IOPT .EQ. 0) GO TO 500
      IF(IOPT .EQ. 1) LONGPR = .NOT. LONGPR
      IF(IOPT .EQ. 2) VAXCOM = .NOT. VAXCOM
      IF(IOPT .EQ. 2) WRITNG = .NOT. WRITNG
      IF(IOPT .EQ. 4) RAWDTA = .NOT. RAWDTA
      IF(IOPT .EQ. 3) GO TO 200
      IF(IOPT .EQ. 4) GO TO 300
      IF(IOPT .EQ. 5) GO TO 400
      IF (IOPT. EQ. 6)CALL READNM
      IF(IOPT .EQ. 6) CALL READC
      IF (IOPT.EQ.5)CALL VAXINI
      GO TO 150
C
C Get new preset time

```

```

C
200  TYPE 2000
2000  FORMAT(/'#Enter preset assay time (seconds) -> ')
      READ (5,2001,END=200,ERR=200)TIMEPR
2001  FORMAT(F10.0)
      IF(TIMEPR .LT. 10000.0) GO TO 150

C
C  Check on excessive assay time
C
2055  TYPE 2050,TIMEPR
2050  FORMAT(/'#Is ',F10.0,' seconds really correct? (Y/N) -> ')
      READ (5,1001,END=2055,ERR=2055)RANS
1001  FORMAT(15A1)
      IF(RANS .EQ. 'Y') GO TO 150
      IF(RANS .EQ. 'N') GO TO 200
      GO TO 2055

C
C  Get new value for Pu-240 abundance
C
300  TYPE 3000
3000  FORMAT(/'#Enter Pu-240 abundance (Wt.%) (0=1.97, -1=Calc.) -> ')
      READ (5,3100,END=300,ERR=300)PU240
3100  FORMAT(F8.3)
      GO TO 150

C
C  Get new value for Pu-242 abundance
C
400  TYPE 4000
4000  FORMAT(/'#Enter Pu-242 abundance (Wt.%) (0=.1) -> ')
      READ (5,3100,END=300,ERR=300)PU242
      GO TO 150

500  TYPE 1000, LONGPR, WRITNG, TIMEPR, PU240, PU242, CONFIL
      TYPE 5000
5000  FORMAT(/'#OK? (Y/N) -> ')
      READ (5,1001,END=400,ERR=400)RANS
      IF (RANS.EQ.'Y' .OR. RANS.EQ.'y')GO TO 600
      GO TO 150

C
600  RETURN
      END

```

```

C
C*****
C
C ROTMON
C
C This subroutine monitors the rotation and the translation motion
C by the sample elevator of the sample in the PU-238
C isotopic instrument through an ADAC 1616CCI interface board.
C Bit 4 in the ADAC data register will be set, if the sample stops
C rotating. If rotation is detected to have stopped, the logical
C flag "RFLG" is set FALSE. Bit 3 in the ADAC data register is set,
C if the sample elevator stops. If the elevator is detected to
C have stopped, the logical flag "EFLG" is set to FALSE.
C
C REV 08-AUG-87 RPD Using new ADLIBVMS calls.
C
C*****
C
C SUBROUTINE ROTMON(RFLG,EFLG)
C
C INTEGER*4 IMASK, STATUS, RSTATUS, ISTAT
C
C The following common block file is generic to all instruments
C
C INCLUDE 'CONTRL.CMN'
C INCLUDE 'ADAC.CMN'
C
C LOGICAL*1 RFLG,EFLG
C
C RFLG = .TRUE.
C EFLG = .TRUE.
C
C Look at Data register and see if bit 4 is set
C
C IMASK = 16
C ISTAT = BTSTW(IMASK,RSTATUS,C_DATA,IOSB,,)
C STATUS= ERROR_CHECK(ISTAT,IOSB)
C WRITE (*,*)'ROTATION',RSTATUS
C IF (RSTATUS.EQ.1) GOTO 400
C
C Rotation has stopped; set flag to suspend measurement
C
C RFLG = .FALSE.
C
C Look at Data register and see if bit 3 is set
C
C 400 IMASK = 8
C ISTAT = BTSTW(IMASK,ESTATUS,C_DATA,IOSB,,)
C STATUS= ERROR_CHECK(ISTAT,IOSB)
C WRITE (*,*)'ELEVATOR',ESTATUS
C IF (ESTATUS.EQ.1) GOTO 500
C
C Elevator has stopped moving; set flag
C
C EFLG = .FALSE.
C 500 RETURN
C END

```

C\*\*\*\*\*

C WRTLGB (WRITLG subroutine with Pu-238 isotopic abundance  
C information included in write to assay log file.)

C Write the most recent assay information into the  
C assay log file

C The format of a log entry is as follows:

Var nam	Type	Dimens	Meaning
ASDATE	L*1	9	Assay date (eg 15-JUL-85)
ASTIME	L*1	8	Assay time (eg 14:13:27)
SAMPID	C*20		Sample ID (ASCII string)
OPERID	C*8		Operator ID (ASCII string)
MATTYP	L*1	3	Material type (ASCII string)
ASTYPE	L*1	2	Assay type (assay = "AS", calib = "CA")
LGFILE	C*15		Name of raw data file (null-term ASCII string, eg SYO:FILNAM.EXT)
ICYCLE	I*2	1	Current run cycle
NUMCYC	I*2	1	Number of cycles in run
ASANS	R*4	1	Assay result
ERRBAR	R*4	1	Fractional error in result
DGFLG	L*1	10	Data diagnostic flags
BIASOK	L*1	1	MC bias run is current and in limit
PRECOK	L*1	1	MC precision run is current and in limit
BACKOK	L*1	1	Background run is current and in limit
BAD	L*1	1	Bad assay flag
CONDIT	L*1	1	Conditional assay flag
C240FG	L*1	1	Pu-240 value calculated or entered
WTPCT	R*4	6	Isotopic weight percent abundances
PCT	R*4	6	Percent errors on isotopic abundances

C Record 1 of the log file always contains two pointers:  
C LPTR is the record number of the last (most recent) entry.  
C FPTR is the record number of the first (earliest) entry.

- C SSJ 22-OCT-82
- C REV 04-FEB-83 SSJ
- C REV 25-sep-84 CMS
- C REV 20-Feb-85 CMS Added error bar and assay type to  
C log entries
- C Rev 08-Mar-85 MPK Correct spelling of DGFLG variable
- C Rev 08-Jul-85 MPK Simplify log format
- C Include more information in log
- C Handle all errors in log file
- C Rev 07-Apr-86 WDR Modified for Micro VAX II
- C Rev 13-Jul-87 RDP Modified for uVAX II
- C Convert BYTE strings to CHARACTER strings
- C Rev 20-Mar-90 WDR Modified for Pu-238 system; changed  
C LVFLAG to C240FG

C\*\*\*\*\*

C SUBROUTINE WRITLG

C The following common block files are generic for all instruments

C INCLUDE 'CLUNS.CMN'



```
C
300 IF (.NOT.WRITNG) GO TO 330
310   LGFILE=FILNAM
   GO TO 350
C
330   LGFILE=NOFILE
L
350   CONTINUE
C
C Update the pointer in record 1
C Write the most recent assay information in the next record
C Write every assay taken to the log including autocycle runs
C
   IERR=3
   WRITE (DKLUN'1,ERR=390) LPTR,FPTR
   WRITE (DKLUN'LPTR,ERR=390)
   1     ASDATE, ASTIME, SAMPID, OPERID, MATTYP, ASTYPE,
   2     LGFILE, ICYCLE, NUMCYC, ASANS, ERRBAR, DGFLG,
   3     BIASOK, PRECOK, BACKOK,   BAD, CONDIR, C240FG,
   4     (WTPCT(I),PCT(I)/100., I=1,6)
   GO TO 400
C
C Close file
C
390 CALL ERRMSG(IERR,ASLOG)
400 CLOSE (UNIT=DKLUN,ERR=500)
C
500 RETURN
   END
```

```
C
C*****
C
C  WRTDAT
C
C      Routine is used to call the subroutine that retrieves
C      the data from the hardware and writes it to disk
C      with the background data appended to the end of the
C      raw data file.  The file name is the sample id up to
C      6 characters with the instrument extension.
C
C      13-JUN-84  CMS
C*****
C
C      SUBROUTINE WRTDAT
C
C      The following common block files are instrument specific for
C      gamma-ray instruments
C
C      INCLUDE 'CONT90.CMN'
C
C      Get90 is called for the gamma instruments - a different
C      routine will need to be called for the neutron instruments
C
C      CALL GET90(IADC)           !instrument specific
C
C      RETURN
C      END
```



C\*\*\*\*\*

C WTMCLG

C Write the most recent measurement control information into the  
C measurement control log file

L The format of a log entry is as follows:

Var nam	Typ	Dim	Bias/ Prec/ backG	Meaning
ASDATE	C*9	-	B/P/G	Assay date (eg 15-JUL-85)
ASTIME	C*8	-	B/P/G	Assay time (eg 14:13:27)
SAMPID	C*1	20	B/P/G	Sample ID (ASCII string; BKG for backgr)
OPERID	C*1	8	B/P/G	Operator ID (ASCII string)
Mtyp	L*1	2	B	The string "MB"
			P	The string "MP"
			G	The string "BG"
MCRUNS	I*2	1	B	ID of calibration
			P	Number of runs in precision test
VALUE	R*4	1	B	Assay result
			P	Reduced chi squared
			G	Background value
SIGMA	R*4	1	B	Fractional error in result
			P	Measured variance
			G	S.D. of background
VAR	R*4	1	B	Standard value (for compar with result)
			P	Calculated variance
ARCHV1	R*4	1	B/P	Archival value (sent to VAX)
ARCHV2	R*4	1	B/P	Archival value (sent to VAX)
DGFLG	L*1	10	B/P/G	Data diagnostic flags
BIASOK	L*1	1	B	Bias run was OK
PRECOK	L*1	1	P	Precision run was OK
BACKOK	L*1	1	G	Background run was OK
MCERR	L*1	1	B/P/G	Error in meas control run
MCVAX	L*1	1	B/P	Entry has been sent to VAX

C Record 1 of the log file always contains two pointers:  
C LPTR is the record number of the last (most recent) entry  
C FPTR is the record number of the first (earliest) entry

C CMS 31-MAY-84

C REV 09-AUG-84

C Rev 27-Mar-85 MPK Include CRESLT common

C Rev 09-Jul-85 MPK Simplify log format

C Include more information in log

C Handle all errors in log file

C Rev 07-Apr-86 WDR Modified for Micro VAX II

C\*\*\*\*\*

C SUBROUTINE WTMCLG

C The following common block files are generic for all instruments

C INCLUDE 'CLUNS.CMN'

C INCLUDE 'CASSAY.CMN'

```
INCLUDE 'CFILES.CMN'  
INCLUDE 'CDATE.CMN'  
INCLUDE 'CONTRL.CMN'  
INCLUDE 'CDIAG.CMN'  
INCLUDE 'CMEAS.CMN'  
INCLUDE 'CMC.CMN'
```

```
INTEGER*2 LPTR,FPTR
```

```
DATA MCMAX /100/           !max number of log entries (incl pointer)  
DATA MCREC /32/           !entry length in 4-byte units
```

```
C  
C Open the existing log file
```

```
OPEN (UNIT=DKLUN,NAME=MCLOG,ACCESS='DIRECT',TYPE='OLD',  
1     RECORDSIZE=MCREC,ERR=150)  
GO TO 200
```

```
C Unsuccessful open on existing file
```

```
C Open a new file
```

```
150  IERR=1  
OPEN (UNIT=DKLUN,NAME=MCLOG,ACCESS='DIRECT',TYPE='NEW',  
1     RECORDSIZE=MCREC,MAXREC=MCMAX,ERR=390)
```

```
C Initialize pointers for the new file
```

```
LPTR=2  
FPTR=2
```

```
C Write the last record of the file to force allocation
```

```
IERR=3  
WRITE (DKLUN'MCMAX,ERR=390) 0  
GO TO 300
```

```
C Read pointer to most recent assay in record 1
```

```
200  IERR=2  
READ (DKLUN'1,ERR=390) LPTR,FPTR  
IERR=5  
IF (LPTR.LT.2 .OR. LPTR.GT.MCMAX) GO TO 390  
IF (FPTR.LT.2 .OR. FPTR.GT.MCMAX) GO TO 390
```

```
C Update pointer to latest (current) run;
```

```
C If the log has caught up with its tail, update the pointer  
C to the earliest run to drop the oldest entry
```

```
LPTR=LPTR+1  
IF (LPTR.GT.MCMAX) LPTR=2  
IF (LPTR.NE.FPTR) GO TO 300  
FPTR=FPTR+1  
IF (FPTR.GT.MCMAX) FPTR=2
```

```
C Update the pointer in record 1
```

```
C Write the most recent measurement control information in the next  
C record.
```

```
300  IERR=3
```

```
WRITE (DKLUN'1,ERR=390) LPTR,FPTR
WRITE (DKLUN'LPTR,ERR=390)
1      ASDATE, ASTIME, SAMPID, OPERID,  MTyp, MCRUNS,
2      VALUE,  SIGMA,  VAR, ARCHV1, ARCHV2,  DGFLG,
3      BIASOK, PRECOK, BACKOK,  MCERR,  MCVAX
GO TO 400
```

L

C Close file

C

```
390 CALL ERRMSG(IERR,MCLDG)
400 CLOSE (UNIT=DKLUN,ERR=500)
```

C

```
500 RETURN
END
```

```
C
C *****
C
C SIZESAMP
C
C     This subroutine checks the size of the sample in the PU-238
L     isotopic instrument through an ADAC 1616CCI interface board.
C     Bit 2 in the ADAC data register will be set, if the sample
C     is an EP60/61 container.
C
C     Version 1.0  8-May-91  WDR
C
C *****
C
C     SUBROUTINE SIZESAMP
C
C     INTEGER*4      IMASK, STATUS, SSTATUS, ISTAT
C
C     The following common block file is generic to all instruments
C
C     INCLUDE 'CONTRL.CMN'
C     INCLUDE 'CASSAY.CMN'
C     INCLUDE 'ADAC.CMN'
C
C     Look at Data register and see if bit 2 is set
C
C     IMASK = 4
C     ISTAT = BTSTW(IMASK,SSTATUS,C_DATA,IOSB,,)
C     STATUS= ERROR_CHECK(ISTAT,IOSB)
C     WRITE (*,*)'SAMPLE SIZE',SSTATUS
C     IF (SSTATUS.EQ.0) THEN
C         MATTYP(1) = 'F'
C         MATTYP(2) = 'C'
C         MATTYP(3) = 'C'
C     ELSE
C         MATTYP(1) = 'E'
C         MATTYP(2) = 'P'
C         MATTYP(3) = '6'
C     ENDIF
C
C     TYPE 9999,MATTYP(1)
C     FORMAT(1X,'MATTYP',A1)
9999  RETURN
500   END
```

```

C*****
C  SLEEP
C
C      Sleep for a given number of seconds
C
C      REV 05-JAN-83 SSJ
C      REV 08-FEB-84 MPK
C      Rev 04-Feb-85 MPK          No change in code from version
C                                used in MUADC demo.  Supersedes
C                                the version of SLEEP formerly
C                                used in C90.
C      Rev 06-Jul-87 RDP          Modified for uVAX II w/ ADAC Drivers
C      Rev 14-Jul-87 RDP          Sleeping using VAX second counter
C      Rev 19-Jan-88 WDR          Modified to use system wait function
C                                this lets system hibernate and switch to a
C                                real parameter passed to LIB$WAIT
C*****
C
C      SUBROUTINE SLEEP(time_secs)
C
C      REAL START,ELAPSE,DONE,time_secs
C
C      CALL LIB$WAIT(time_secs)
C      RETURN
C      END

```

SUBROUTINE STABSET

STABSET

=====

Lawrence Livermore National Laboratory Safeguards Technology Program
Actinide Isotopic Analysis Suite

Developed by : Applications Systems Division
and
Nuclear Chemistry Division

Lawrence Livermore National Laboratory
Livermore, CA. 94550

This work supported by LLNL Safeguards Technology Program
and performed under the auspices of the U.S. Department of
ENERGY under Contract W-7405-Eng-48.

Title : stabset
Revision : 1.00
Author : W.M. Buckley
Date : 4 June 1991
Language : FORTRAN 77 (system dependencies noted, when known)

Abstract : this routine sets the values for the control and
status reporting for the Canberra 8232 digital
spectrum stabilizer

Calling arguments: NONE

Routines called:

ci8232\_status function call status routine
ci8232\_gain\_channel sets the gain peak channel
ci8232\_gain\_width sets the gain window width
ci8232\_gain\_range sets the gain analog range
ci8232\_zero\_channel sets the zero peak channel
ci8232\_zero\_width sets the zero window width
ci8232\_zero\_range sets the zero analog range
ci8232\_report reports 8232 setup and status
ci8232\_gain\_hold sets the gain stabilization hold
ci8232\_zero\_hold sets the zero stabilization hold
ci8232\_store stores the 8232 configuration to NVRAM
SYS\$ASSIGN

Called by:
exec238 isotopic application executive routine

Revision History ( Include date, revision I.D., and details )

```

C*****
C
C      1.00      6/4/91  wmb      added to Pu-238 application for WHC
C                                  Pu-238 system, diagnostics and reporting
C                                  turned off

```

```

C*****

```

```

C      C*32      status      string corresponding to status values
C      C*32      status_reset string corresponding to status values
C      C*12      routine      string corresponding to current lib call
C
C      I*2      dsslun       logical unit for serial port connection to
C                                  CI8232
C      I*2      io_chan      I/O channel for serial port connection to
C                                  CI8232
C      I*2      stat         status return variable from 8232lib calls
C      I*2      stat_reset   status return variable from 8232lib calls
C      I*2      channel      peak channel value for gain or zero stab
C      I*2      width        window width
C      I*2      range        analog range for zero or gain stab:
C                                  0          for      1/1
C                                  1          for      1/2
C                                  2          for      1/4
C                                  3          for      1/8
C      I*2      option       stabilization mode 0/1/2 for off/hold/on

```

Definition of COMMON variables:

```

C      I*2      IZPEAK       Peak location for zero stabilization
C      I*2      IZWNDW       Peak window for zero stabilization
C      R*4      ZARNG        Analogr range for zero stabilization
C      I*2      IGPEAK       Peak location for gain stabilization
C      I*2      IGWNDW       Peak window for gain stabilization
C      R*4      GARNG        Analog range for gain stabilization

```

IMPLICIT NONE

```

X      COMMON /P238COM/PU240,PU242,
           IZPEAK,IZWNDW,ZARNG,IGPEAK,IGWNDW,GARNG

```

```

C      REAL*4      pu240,pu242,zarng,garng
C      INTEGER*2   izpeak,izwndw,igpeak,igwndw

```

```

C      CHARACTER*32 status
C      CHARACTER*32 status_reset
C      CHARACTER*12 routine

```

```

C      INTEGER*2   dsslun
C      INTEGER*2   io_chan
C      INTEGER*2   stat
C      INTEGER*2   stat_reset

```

```

C      INTEGER*2   channel
C      INTEGER*2   width
C      INTEGER*2   range
C      INTEGER*2   option

```

```

STRUCTURE /report_struct/
  INTEGER*2 gain_peak
  INTEGER*2 gain_window
  INTEGER*2 gain_range
  INTEGER*2 gain_stab
  INTEGER*2 gain_rate
  INTEGER*2 gain_corr
  INTEGER*2 zero_peak
  INTEGER*2 zero_window
  INTEGER*2 zero_range
  INTEGER*2 zero_stab
  INTEGER*2 zero_rate
  INTEGER*2 zero_corr
END STRUCTURE

```

```

RECORD /report_struct/ report

```

```

C*****
C  variable initializations
C*****

```

```

  dsslun = 8
  stat = 0
  stat_reset = 0

```

```

C ...announce function
  TYPE 2000

```

```

C ...open terminal devices

```

```

  CALL SYS$ASSIGN('TTA1:',io_chan,,)

```

```

  1 OPEN(UNIT=dsslun,FILE='TTA1:',ACCESS='SEQUENTIAL',
    FORM='FORMATTED',CARRIAGECONTROL='NONE',STATUS='OLD')

```

```

C ...generate initial report

```

```

  routine = 'REPORT'
  CALL ci8232_REPORT(dsslun,io_chan,report,stat,status)

```

```

  IF(stat.NE.0)THEN
    WRITE(6,1000)routine,stat,status
    stat = stat_reset
    status = status_reset
  ENDIF

```

```

C
C

```

```

  1 WRITE(6,1001) report.gain_peak,
  2 report.gain_window,
  report.gain_range

```

```

  1 WRITE(6,1003) report.zero_peak,
  2 report.zero_window,
  report.zero_range

```

```

  1 WRITE(6,1002) report.gain_stab,
  report.gain_rate,

```



```

2          report.gain_corr

WRITE(6,1004)  report.zero_stab,
1             report.zero_rate,
2             report.zero_corr

```

C ...set GAIN settings

```

channel = igpeak
width   = igwndw
IF(garng.EQ.1.000)THEN
  range = 0
ELSEIF(garng.EQ.0.500)THEN
  range = 1
ELSEIF(garng.EQ.0.250)THEN
  range = 2
ELSEIF(garng.EQ.0.125)THEN
  range = 3
ENDIF
option  = 2

routine = 'GAIN_PEAK'
CALL ci8232_GAIN_CHANNEL(dsslun,io_chan,channel,stat,status)

IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF

routine = 'GAIN_WINDOW'
CALL ci8232_GAIN_WIDTH(dsslun,io_chan,width,stat,status)

IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF

routine = 'GAIN_RANGE'
CALL ci8232_GAIN_RANGE(dsslun,io_chan,range,stat,status)

IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF

routine = 'GAIN_HOLD'
CALL ci8232_GAIN_HOLD(dsslun,io_chan,option,stat,status)

IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF

```

C ...set ZERO settings

```
channel = izpeak
width   = izwndw
IF(zarng.EQ.1.000)THEN
  range = 0
ELSEIF(zarng.EQ.0.500)THEN
  range = 1
ELSEIF(zarng.EQ.0.250)THEN
  range = 2
ELSEIF(zarng.EQ.0.125)THEN
  range = 3
ENDIF

routine = 'ZERO_PEAK'
CALL ci8232_ZERO_CHANNEL(dssslun,io_chan,channel,stat,status)
```

```
IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF
```

```
routine = 'ZERO_WINDOW'
CALL ci8232_ZERO_WIDTH(dssslun,io_chan,width,stat,status)
```

```
IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF
```

```
routine = 'ZERO_RANGE'
CALL ci8232_ZERO_RANGE(dssslun,io_chan,range,stat,status)
```

```
IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF
```

```
routine = 'ZERO_HOLD'
CALL ci8232_ZERO_HOLD(dssslun,io_chan,option,stat,status)
```

```
IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF
```

C ...store values in NVRAM

```
routine = 'STORE'
CALL ci8232_STORE(dssslun,io_chan,stat,status)
```

```
IF(stat.NE.0)THEN
  WRITE(6,1000)routine,stat,status
  stat = stat_reset
  status = status_reset
ENDIF
```

```

C
C Assign Data Channel for CCI
C
    DEVICE = 'ZWA1:'
    ISTAT = SYS#ASSIGN(DEVICE,C_DATA,,)
    ESTATUS = ERROR_CHECK(ISTAT,1)
    IF (.NOT.ESTATUS) THEN
        TYPE 100,DEVICE
100     FORMAT(1X,'Failure to assign channel to ',A5)
        ABORT = .TRUE.
    ENDIF

C
C Assign Data Channel for HCO
    DEVICE = 'ZWBO:'
    ISTAT = SYS#ASSIGN(DEVICE,H_DATA,,)
    ESTATUS = ERROR_CHECK(ISTAT,1)
    IF (.NOT.ESTATUS) THEN
        TYPE 200,DEVICE
200     FORMAT(1X,'Failure to assign channel to ',A5)
        ABORT = .TRUE.
    ENDIF

5     IMASK = 4096
C     Send HIGH to reset
    ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
    CALL ERROR_CHECK(ISTAT,IOSB)
C     Return to LOW
    ISTAT = BCLRW(IMASK,H_DATA,IOSB,,)
    CALL ERROR_CHECK(ISTAT,IOSB)
    IMASK = 2048
C     Send HIGH to reset
    ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
    CALL ERROR_CHECK(ISTAT,IOSB)
C     Return to LOW
    ISTAT = BCLRW(IMASK,H_DATA,IOSB,,)
    CALL ERROR_CHECK(ISTAT,IOSB)
C Turn Off Rotation
    IMASK = 1
    ISTAT = BSETW(IMASK,H_DATA,IOSB,,)
    CALL ERROR_CHECK(ISTAT,IOSB)
C Enable LOAD/UNLOAD
    IMASK = 512
    CALL = BCLRW(IMASK,H_DATA,IOSB,,)
    CALL ERROR_CHECK(ISTAT,IOSB)

C
C Initialize Series 90
C
    CALL INIT90                                !instrument specific

C
    RETURN
    END

```

```

C
C*****
C STATUS
C Routine checks the status of the hardware for failures
C prior to the beginning of an assay or responds to
C the option 'S' entered from the menu list which
C checks if hardware is busy
C
C CALLING ARGUMENTS:
C ITYP = Logical flag indicating the status of hardware
C ITYP = 0: hardware okay to run
C ITYP = 1: hardware busy checked
C
C CMS 24-JUL-84
C
C Rev. 29-Nov-84 MPK Replace list-directed with
C formatted I/O
C Rev. 07-Apr-86 WDR Modified for Micro VAX II
C Rev. 22-Sep-87 RDP Status now test to see if
C the detector is plugged in.
C*****

```

```

C
C SUBROUTINE STATUS(ITYP)
C
C INCLUDE 'ADAC.CMN'
C INTEGER*2 VALUE
C INTEGER*4 SIZE
C CHARACTER REPLY*1
C DATA SIZE/1/
C
C If argument passed is 0 - check if hardware is okay to run
C
C ABORT = .FALSE.
C IF (ITYP.EQ.1) GOTO 50
C Get Data Word from ADAC
5 ISTAT = DINW(VALUE,SIZE,C_DATA,IOSB,,)
C 5 CONTINUE
C ESTATUS= ERROR_CHECK(ISTAT,IOSB)
C TYPE*, 'Cable Value:',VALUE
C IF ((VALUE.EQ.96).OR.
C (VALUE.EQ.0).OR.(VALUE.EQ.'FFFF'X)) THEN
C-- Major Error, allow user to take corrective measure
C TYPE 500
500 FORMAT(' **** HARDWARE ERROR DETECTED ****')
C TYPE*, 'ERROR -- Problem with Cables or Connectors'
C TYPE 555, 'Press RETURN to retry (type A and a <CR> to ABORT) ->'
555 FORMAT(1X,A,#)
C READ(5,3000) REPLY
3000 FORMAT(A1)
C IF ((REPLY.NE.'A').AND.(REPLY.NE.'a')) GOTO 5
C ABORT = .TRUE.
C GOTO 20
C ENDIF
C-- Check for Idle System
6 IDLE = INQADC(1,ISTAT,MEMSIZ,IGAIN,IRANGE,IOFF,INPUTS)
C TYPE *, 'IDLE',IDLE
C IF (IDLE.NE.1)GOTO 15
14 TYPE 1500
C TYPE 555, 'Press RETURN to retry (type A and a <CR> to ABORT) ->'

```

```
READ(5,3000) REPLY
IF ((REPLY.NE.'A').AND.(REPLY.NE.'a')) GOTO 6
ABORT = .TRUE.
GOTO 20
```

```
C
15 TYPE 1000
1500 FORMAT(' **** STATUS OF HARDWARE OKAY ****')
GOTO 20
```

```
C
C If argument is 1 check if hardware is busy
```

```
C
50 IDLE = INQADC(1,ISTAT,MEMSIZ,IGAIN,IRANGE,IOFF,INPUTS)
IF (IDLE .NE. 1) GOTO 55
54 TYPE 1500
1500 FORMAT('/' **** !! SERIES 90 BUSY !! MANUAL USE?? ****',/)
GO TO 20
55 TYPE 2000
2000 FORMAT(' **** HARDWARE BUSY CHECKED ****')
20 RETURN
END
```

## SUBROUTINE SVPKS

```

C
C
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
C
      Source file name:  SVPK38.FOR
C
L
C      Sub to save results of peak fitting in a disk file for use by
C      EVAL in a separate program.  Some of the initial parameters plus
C      results for each peak are appended to any data already stored in
C      the file.  This permits optional text editing as well as combining
C      with results from other spectra of the same sample.
C      Last edit 22-Jan-85
C      Modified 4-Mar-90 for Pu-238 isotopic instrument
C
CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

      INCLUDE 'GRPLGCM.PMS'
      INCLUDE 'GRPGMS.PMS'

      INCLUDE 'CONTRL.CMN'
      INCLUDE 'CLUNS.CMN'
      INCLUDE 'CFILES.CMN'
      INCLUDE 'CASSAY.CMN'

      BYTE DEFSTG(30),EFSMBL,TMPFIL
      DATA DEFSTG /'G','R','P','O','U','T','.','T','M','P',19*' ',0/
      DATA LP/-1/      !For VAX
      DATA LP/6/
      DATA EFSMBL/' '/

      TMPFIL = 0
      IF(RUNFLG(4) .EQ. 0) EFSMBL = '*'

      IF(STGFIL(1) .EQ. ' ')STGFIL(1) = 0
      IF(STGFIL(1) .NE. 0)GO TO 5
      DO 2, I = 1,30
2      STGFIL(I) = DEFSTG(I)
      TMPFIL = 1
5      OPEN(UNIT=NDEV,NAME=STGFIL,STATUS='NEW',ERR=375)
      TYPE *,'NDEV = ',NDEV

8      WRITE(NDEV,1060)NGAMS,KTGRP,(DTE(J),J=1,9),(TBUF(I),I=1,8)
1060  +  FORMAT(// ' GRPANL results.  Photons/min for',I3,' peaks in',
      I3,' groups (calculated ',9A1,2x,8A1,')')

      LENSF = MXFNM-1
      DO 100, K = 1,MXFNM-1
      IF(FILNAM .EQ. ' ')GO TO 105
100  CONTINUE
      GO TO 110
105  LENSF = K-1
110  WRITE(NDEV,1110),FILNAM
1110  FORMAT(' Spectrum file : ',15A)
      WRITE(NDEV,1015)EXPTID
1015  FORMAT(' Expt = ',12A1)
      WRITE(NDEV,1020)SAMPID
1020  FORMAT(' Sample id = ',20A)
      WRITE(NDEV,1022)DSNAME
1022  FORMAT(' Detector-system = ',6A1)

```

```

WRITE(NDEV,1024)DSTTIM,ICTYR
1024 FORMAT(' Start time = ',F11.6,',',',I6)
WRITE(NDEV,1025)DTZERO,IZYR
1025 FORMAT(' Zero time = ',F11.6,',',',I6)
WRITE(NDEV,1026)DELTAT
36 FORMAT(' Decay time (days) = ',1PE15.7)
WRITE(NDEV,1027)TIMEL
1027 FORMAT(' Live time (mins) = ',1PE15.7)
WRITE(NDEV,1128)SMPWT
1128 FORMAT(' Sample weight = ',1PE15.7)
WRITE(NDEV,1028)WT100
1028 FORMAT(' 100% Weight = ',1PE15.7)
WRITE(NDEV,1029)ANSFCT
1029 FORMAT(' Normalization factor = '1PE15.7)
WRITE(NDEV,1030)DCNST(17),DCNST(18)
1030 FORMAT(' Pair peak constants = ',1PE15.7,',',',E15.7)
LENLB = MXFNM-1
DO 40, K = 1,MXFNM-1
IF(LIBNAM(K) .EQ. 0)GO TO 45
40 CONTINUE
GO TO 47
45 LENLB = K-1
47 WRITE(NDEV,1047),(LIBNAM(I),I=1,LENLB)
1047 FORMAT(' Library file = ',29A1)

1050 FORMAT(/' PEAK CHANNEL ENERGY PKSTART PKEND ',
+ ' COUNTS(Bkg?) GAMMAS/MIN FRCER')
1051 FORMAT(/' PEAK CHANNEL ENERGY PKSTART PKEND ',
+ ' COUNTS(Bkg?) COUNTS/MIN FRCER')

IF(RUNFLG(4) .NE. 0)THEN
WRITE(NDEV,1050)
ELSE
WRITE(NDEV,1051)
ENDIF
WRITE(NDEV,1055)
1055 FORMAT(' -----Begin table')
C IF (LONGPR) THEN
C IF(RUNFLG(4) .EQ. 0) THEN
C TYPE 1051
C ELSE
C TYPE 1050
C ENDIF
C ENDIF
DO 310, I= 1,NGAMS
C IF(LONGPR)TYPE 1200,I,CHNL(I),ENGY(I),PKS(I),PKE(I),
C + CNTS(I),BKFLGS(I),GAMS(I),FRCER(I)
WRITE(NDEV,1201)I,CHNL(I),ENGY(I),PKS(I),PKE(I),CNTS(I),
+ BKFLGS(I),GAMS(I),FRCER(I)

1200 FORMAT(I4,4(', ',F8.3),',',',F11.0,A1,',',',1PE12.4,',',',E11.2)
1201 FORMAT(I4,4(', ',F8.3),',',',F11.0,',',',A1,',',',1PE12.4,',',',E11.2)
310 CONTINUE
IF(RUNFLG(4) .EQ. 0) WRITE(NDEV,1310)
1310 FORMAT(' * Detector efficiency was assumed to be unity.')

WRITE(NDEV,1311),FILNAM
1311 FORMAT(' END output for ',15A)
CLOSE(NDEV)

```

GO TO 400

370 TYPE 1370  
1370 FORMAT(' ERROR opening default storage file')  
GO TO 420

375 TYPE 1375,(STGFIL(I),I=1,MXFNM-1)  
1375 FORMAT(' \*\*\*\*\* Error in opening disk storage file ',29A1)  
ABORT = .TRUE.  
GO TO 420

400 DO 405 I = LENSF,MXFNM-1  
405 IFNAME(I) = '  
C IF(TMPFIL .NE. 0) TYPE 1403

C IF (LONGPR)WRITE(LOUT,1400),FILNAM,(STGFIL(I),I=1,MXFNM-1)  
1400 FORMAT(' Results for ',A,'were inserted in file ',15A)  
420 RETURN  
END



C\*\*\*\*\*

C START

C Initialize Program

C SSJ 17-AUG-79

C REV. 17-APR-84 LC

C REV. 13-JUL-84 CMS

C Rev. 03-Jan-85 MPK

Go to VAX for time sync when  
setting date and time

C Rev. 30-Jan-85 MPK

Include CVAX common

C Rev. 06-Mar-85 MPK

Initialize VAX communication

C Rev. 06-Mar-85 MPK

Take out VAX Comm init and time  
synch. They are being moved to  
EXEC.

C Rev. 07-Apr-86 WDR

Modified for Micro VAX II

C Rev. 07-Aug-87 RDP

Modified for allocation of ADAC  
drivers.

C Rev. 25-Aug-87 RDP

Clear SET lights

C\*\*\*\*\*

C SUBROUTINE START

C CHARACTER DEVICE\*5  
C INTEGER\*4 ISTAT, IMASK

C Required common blocks and declaration for ADAC calls.

C INCLUDE 'ADAC.CMN'  
C INCLUDE '(\$SYSSRVNAM)'

C The following common block files are generic for all instruments

C INCLUDE 'CONTRL.CMN'  
C INCLUDE 'CLUNS.CMN'  
C INCLUDE 'CASSAY.CMN'  
C INCLUDE 'CVAX.CMN'

C The following common block files are instrument specific for  
C gamma-ray instruments

C INCLUDE 'CSPECT.CMN'

C Initialize logical flags

C CALL INTFLG

C Set up beginning

C BEGIN = .TRUE.

C Initialize hardware necessary for the instrument - checking  
C that it is operational prior to the measurement beginning

C ABORT = .FALSE.

C Allocate Channels for ADAC Control and Monitoring

C ...generate final report

```
routine = 'REPORT'  
CALL ci8232_REPORT(dssslun,io_chan,report,stat,status)
```

```
IF(stat.NE.0)THEN  
  WRITE(6,1000)routine,stat,status  
  stat = stat_reset  
  status = status_reset  
ENDIF
```

```
1 WRITE(6,1001) report.gain_peak,  
2 report.gain_window,  
report.gain_range
```

```
1 WRITE(6,1003) report.zero_peak,  
2 report.zero_window,  
report.zero_range
```

```
1 WRITE(6,1002) report.gain_stab,  
2 report.gain_rate,  
report.gain_corr
```

```
1 WRITE(6,1004) report.zero_stab,  
2 report.zero_rate,  
report.zero_corr
```

```
CLOSE(UNIT=dssslun)
```

```
C*****  
C format declarations *  
C*****
```

```
1000 FORMAT(' Error in ',A12,', ID = ',I3,', msg is ',A32)  
1001 FORMAT(' GAIN Peak = ',I5,' Window = ',I5,' Range = ',I3)  
1002 FORMAT(' GAIN Stab = ',I5,' Rate = ',I5,' % corr = ',I3)  
1003 FORMAT(' ZERO Peak = ',I5,' Window = ',I5,' Range = ',I3)  
1004 FORMAT(' ZERO Stab = ',I5,' Rate = ',I5,' % corr = ',I3)  
2000 FORMAT  
1(' 8232 Digital Stabilizer being set from values in PARMTR file')
```

```
C*****
```

```
RETURN  
END
```

```

C*****
C
C  YESNO
C      Prompt for and accept a yes-or-no answer
C
C      Calling argument:
C      PROMPT  The string to be used for prompting.  It is output,
C              with the characters " ? (Y/N) -> " appended to it.
C
C      Function value (LOGICAL*1):
C      .TRUE.  if the answer begins with the character Y
C      .FALSE. if the answer begins with the character N
C
C      If the answer does not begin "Y" or "N", the prompt
C      is repeated until it does.
C
C      Example:
C      If YESNO is called as follows
C          LOGICAL*1 CONT,YESNO
C          CONT=YESNO('Do you want to continue')
C      the following prompt is produced
C          Do you want to continue ? (Y/N) ->
C
C      04-Feb-85  MPK
C      06-Aug-87  RDP          Convert BYTE string to Character String

```

```

C*****
C
C      LOGICAL*1 FUNCTION YESNO(PROMPT)
C*****
C
C      CHARACTER PROMPT*(*),TAG*10,QUESTION*70,ANS*1
C      INTEGER    LOC
C
C      DATA TAG/'? (Y/N) ->'/
C      LOC = LEN(PROMPT)+1
C
C      QUESTION = PROMPT
C      QUESTION(LOC:LOC+10) = TAG
C
C      TYPE 1000, QUESTION
1000  FORMAT(1X,A<LOC+11>,$)
      READ (5,1002,END=10,ERR=10) ANS
1002  FORMAT (A)
      YESNO=.TRUE.
      IF (ANS.EQ.'Y') RETURN
      YESNO=.FALSE.
      IF (ANS.EQ.'N') RETURN
      GO TO 10
      END

```



C 1.02 6/5/91 wmb added QIO routine for port reads

C \*\*\*\*\*

C //  
C \\\

SUBROUTINE timed\_read(io\_chan,string,length,ier)

C \*\*\*\*\*

Title : timed\_read  
Revision : 1.00

C \*\*\*\*\*

Abstract : this routine performs a timed QIO-based read of the  
port to the Canberra 8232 stabilizer

C \*\*\*\*\*

Arguments : io\_chan I I\*2 I/O channel for serial  
communication with 8232  
ier 0 I\*2 status code(0/1 for OK/error)  
length I/O I\*2 length of read  
string 0 C\*120 character read buffer

C \*\*\*\*\*

Routines  
Called : SYS#QIOW

Libraries  
Used : none

C \*\*\*\*\*

C \*\*\*\*\*

C\*120 string calling argument  
I\*2 io\_chan calling argument  
I\*2 ier calling argument  
I\*2 length calling argument

IMPLICIT NONE

...declarations required for QIO

INCLUDE '(#IODEF)'

INTEGER\*4 SYS#QIOW  
EXTERNAL SYS#QIOW

INTEGER\*2 IDSB(4)  
INTEGER\*4 timeout  
INTEGER\*4 ierr

C .....

```
INTEGER*2      io_chan
INTEGER*2      length
INTEGER*2      ier

CHARACTER*120  string
```

```
C*****
C      variable initializations
C*****
```

```
      timeout = 30
      string = ' '
```

```
C ...perform timed read
```

```
      ierr = SYS$QIOW(,
1          %VAL(io_chan),
1          %VAL(IO$_READLBLK+IO$_M_NOECHO+IO$_M_TIMED),
1          IOSB,,,
1          %REF(string),
1          %VAL(length),
1          %VAL(timeout),,,,)

```

```
C ...obtain actual number of characters read
```

```
      length = IOSB(2)
```

```
C ...transfer return status to I*2 variable
```

```
      IF(ierr.EQ.1)THEN
          ier = 0
      ELSE
          ier = ierr
      ENDIF
```

```
C*****
```

```
      RETURN
      END
```

```
C////////////////////////////////////
C////////////////////////////////////
```

```
      SUBROUTINE ci8232_STATUS(lun,io_chan,stat,status)
```

```
C*****
```

```
      Title      :      ci8232_status
      Revision    :      1.00
```

```
C*****
```

```
      Abstract   :      this routine provides status on the previous call
                        to the Canberra 8232 stabilizer
```

```
C*****
```

```
      Arguments :      lun      I      I*2      logical unit for serial
                        communication with 8232
```

```

C          io_chan I          I*2          I/O channel for serial
C          communication with 8232
C          stat I          I*2          status code(0/1 for OK/error)
C          status 0          C*32          status string
C

```

```

C*****

```

```

L
C          Routines
C          Called      :      timed_read
C
C          Libraries
C          Used        :      none
C

```

```

C*****
C*****

```

```

C          C*32      status      calling argument
C          C*32      stat_string(6) strings corresponding to status values
C
C          C*1       esc         ASCII code for escape
C          C*1       code        1 character 8232 code for generating report
C
C          I*2       lun         calling argument
C          I*2       io_chan     calling argument
C          I*2       stat        calling argument
C
C          I*2       length      length of port read
C          I*2       clength     anticipated length of port read
C          C*120     string       character buffer for port read
C          I*2       ier         return status of timed read
C

```

IMPLICIT NONE

```

CHARACTER*32      status
CHARACTER*32      stat_string(12)

```

```

CHARACTER*1      esc
CHARACTER*1      code

```

```

INTEGER*2        lun
INTEGER*2        io_chan
INTEGER*2        stat

```

```

INTEGER*2        length
INTEGER*2        clength
INTEGER*2        ier

```

```

CHARACTER*120    string

```

```

C*****
C          variable initializations
C*****

```

```

esc = CHAR(27)
code = 'L'

```

```

stat_string(1) = 'Invalid number of characters'
stat_string(2) = 'Invalid parameter value'
stat_string(3) = 'Invalid terminator'
stat_string(4) = 'No ESC in sequence'

```





SUBROUTINE ci8232\_GAIN\_CHANNEL(lun,io\_chan,channel,stat,status)

Title : ci8232\_gain\_channel
Revision : 1.00

Abstract : this routine sets the gain channel on the Canberra 8232 stabilizer

Arguments : lun I I\*2 logical unit for serial communication with 8232
io\_chan I I\*2 I/O channel for serial communication with 8232
channel I I\*2 gain peak channel
stat 0 I\*2 status code(0/1 for OK/error)
status 0 C\*32 status string

Routines Called : ci8232\_status

Libraries Used : none

C\*1 esc ASCII code for escape
C\*1 code 1 character 8232 code for gain peak channel
C\*32 status calling argument
I\*2 lun calling argument
I\*2 io\_chan calling argument
I\*2 stat calling argument
I\*2 channel calling argument

IMPLICIT NONE

CHARACTER\*1 esc
CHARACTER\*1 code

CHARACTER\*32 status

INTEGER\*2 lun
INTEGER\*2 io\_chan
INTEGER\*2 stat
INTEGER\*2 channel

variable initializations



C\*\*\*\*\*

```

C
C      C*1      esc      ASCII code for escape
C      C*1      code     1 character 8232 code for gain window width

C      C*32     status   calling argument

C      I*2      lun      calling argument
C      I*2      io_chan  calling argument
C      I*2      stat     calling argument

C      I*2      width   calling argument

```

IMPLICIT NONE

```

CHARACTER*1      esc
CHARACTER*1      code

```

```

CHARACTER*32     status

```

```

INTEGER*2        lun
INTEGER*2        io_chan

```

```

INTEGER*2        stat
INTEGER*2        width

```

C\*\*\*\*\*

C variable initializations \*

C\*\*\*\*\*

```

      esc = CHAR(27)
      code = 'A'

```

C ...initiate gain window width set on 8232

```

      WRITE(lun,1000)esc,width,code

```

C ...read status report from 8232

```

      CALL ci8232_status(lun,io_chan,stat,status)

```

C\*\*\*\*\*

C format declarations \*

C\*\*\*\*\*

```

1000      FORMAT(1X,A1,I3.3,A1,#)

```

C\*\*\*\*\*

```

      RETURN
      END

```

```

//
//
//

```

```

SUBROUTINE ci8232_GAIN_RANGE(lun,io_chan,range,stat,status)

```

```

C*****
C
C      Title      :      ci8232_gain_range
C      Revision   :      1.00
C
C..*****
C
C      Abstract   :      this routine sets the gain analog range on the
C                       Canberra 8232 stabilizer
C
C*****
C
C      Arguments :      lun      I      I*2      logical unit for serial
C                       io_chan I      I*2      communication with 8232
C                       range  I      I*2      I/O channel for serial
C                       stat   0      I*2      communication with 8232
C                       status 0      C*32     gain analog range
C                       status 0      C*32     status code(0/1 for OK/error)
C                       status 0      C*32     status string
C
C*****
C
C      Routines
C      Called    :      ci8232_status
C
C      Libraries
C      Used      :      none
C
C*****
C*****
C
C      C*1      esc      ASCII code for escape
C      C*1      code     1 character 8232 code for gain analog range
C
C      C*32     status   calling argument
C
C      I*2      lun      calling argument
C      I*2      io_chan  calling argument
C      I*2      stat     calling argument
C
C      I*2      range   calling argument
C
C      IMPLICIT NONE
C
C      CHARACTER*1      esc
C      CHARACTER*1      code
C
C      CHARACTER*32     status
C
C      INTEGER*2        lun
C      INTEGER*2        io_chan
C
C      INTEGER*2        stat
C      INTEGER*2        range
C
C..*****
C      variable initializations
C*****

```

```
esc = CHAR(27)
code = 'B'
```

```
...initiate gain analog range set on 8232
```

```
WRITE(lun,1000)esc,range,code
```

```
C ...read status report from 8232
```

```
CALL ci8232_status(lun,io_chan,stat,status)
```

```
C*****
C format declarations *
C*****
```

```
1000 FORMAT(1X,A1,I1,A1,*)
```

```
C*****
```

```
RETURN
END
```

```
C////////////////////////////////////
C////////////////////////////////////
```

```
SUBROUTINE ci8232_ZERO_CHANNEL(lun,io_chan,channel,stat,status)
```

```
C*****
C
C Title : ci8232_zero_channel
C Revision : 1.00
```

```
C*****
C
C Abstract : this routine sets the zero channel on the
C Canberra 8232 stabilizer
```

```
C*****
C
C Arguments : lun I I*2 logical unit for serial
C communication with 8232
C io_chan I I*2 I/O channel for serial
C communication with 8232
C channel I I*2 zero peak channel
C stat O I*2 status code(0/1 for OK/error)
C status O C*32 status string
```

```
C*****
C
C Routines
C Called : ci8232_status
C
C Libraries
C Used : none
```



SUBROUTINE ci8232\_ZERO\_WIDTH(lun,io\_chan,width,stat,status)

Title : ci8232\_zero\_width
Revision : 1.00

Abstract : this routine sets the zero window width on the
Canberra 8232 stabilizer

Arguments : lun I I\*2 logical unit for serial
communication with 8232
io\_chan I I\*2 I/O channel for serial
communication with 8232
width I I\*2 zero window width
stat 0 I\*2 status code(0/1 for OK/error)
status 0 C\*32 status string

Routines
Called : ci8232\_status

Libraries
Used : none

C\*1 esc ASCII code for escape
C\*1 code 1 character 8232 code for zero window width
C\*32 status calling argument
I\*2 lun calling argument
I\*2 io\_chan calling argument
I\*2 stat calling argument
I\*2 width calling argument

IMPLICIT NONE

CHARACTER\*1 esc
CHARACTER\*1 code

CHARACTER\*32 status

INTEGER\*2 lun
INTEGER\*2 io\_chan

INTEGER\*2 stat
INTEGER\*2 width

```

C      variable initializations
C*****
      esc = CHAR(27)
      code = 'D'

C      ...initiate zero window width set on 8232

      WRITE(lun,1000)esc,width,code

C      ...read status report from 8232

      CALL ci8232_status(lun,io_chan,stat,status)

C*****
C      format declarations
C*****

1000    FORMAT(1X,A1,I3.3,A1,#)

C*****

      RETURN
      END

```

```

C////////////////////////////////////
C////////////////////////////////////

```

```

      SUBROUTINE ci8232_ZERO_RANGE(lun,io_chan,range,stat,status)

```

```

C*****
C
C      Title      :      ci8232_zero_range
C      Revision   :      1.00
C
C*****
C
C      Abstract   :      this routine sets the zero analog range on the
C                       Canberra 8232 stabilizer
C
C*****
C
C      Arguments :      lun      I      I*2      logical unit for serial
C                       io_chan I      I*2      communication with 8232
C                       range   I      I*2      I/O channel for serial
C                       stat     0      I*2      communication with 8232
C                       status  0      C*32     zero analog range
C                       status  0      I*2      status code(0/1 for OK/error)
C                       status  0      C*32     status string
C
C*****
C
C      Routines
C      Called    :      ci8232_status
C
C      Libraries

```



```

C      Used      :      none
C
C*****
C*****
C
C      C*1      esc      ASCII code for escape
C      C*1      code     1 character 8232 code for zero analog range
C
C      C*32     status   calling argument
C
C      I*2      lun      calling argument
C      I*2      io_chan  calling argument
C      I*2      stat     calling argument
C
C      I*2      range    calling argument
C
C      IMPLICIT NONE
C
C      CHARACTER*1      esc
C      CHARACTER*1      code
C
C      CHARACTER*32     status
C
C      INTEGER*2        lun
C      INTEGER*2        io_chan
C
C      INTEGER*2        stat
C      INTEGER*2        range
C*****
C      variable initializations      *
C*****
C      esc = CHAR(27)
C      code = 'E'
C
C      ...initiate zero analog range set on 8232
C
C      WRITE(lun,1000)esc,range,code
C
C      ...read status report from 8232
C
C      CALL ci8232_status(lun,io_chan,stat,status)
C*****
C      format declarations      *
C*****
C      1000  FORMAT(1X,A1,I1,A1,#)
C*****
C
C      RETURN
C      END

```

SUBROUTINE ci8232\_GAIN\_HOLD(lun,io\_chan,option,stat,status)

Title : ci8232\_gain\_hold
Revision : 1.00

Abstract : this routine sets the gain stabilization on the
Canberra 8232 stabilizer

Arguments : lun I I\*2 logical unit for serial
communication with 8232
io\_chan I I\*2 I/O channel for serial
communication with 8232
option I I\*2 0/1/2 for off/hold/on
stat 0 I\*2 status code(0/1 for OK/error)
status 0 C\*32 status string

Routines
Called : ci8232\_status

Libraries
Used : none

C\*1 esc ASCII code for escape
C\*1 code 1 character 8232 code for gain stabilization
C\*32 status calling argument
I\*2 lun calling argument
I\*2 io\_chan calling argument
I\*2 stat calling argument
I\*2 option calling argument

IMPLICIT NONE

CHARACTER\*1 esc
CHARACTER\*1 code

CHARACTER\*32 status

INTEGER\*2 lun
INTEGER\*2 io\_chan

INTEGER\*2 stat



```
C      Routines
C      Called      :      ci8232_status
C
C      Libraries
C      Used        :      none
```

```
.....
C*****
```

```
C
C      C*1      esc      ASCII code for escape
C      C*1      code     i character 8232 code for zero stabilization
C
C      C*32     status   calling argument
C
C      I*2      io_chan  calling arguemnt
C      I*2      lun      calling arguemnt
C      I*2      stat     calling argument
C
C      I*2      option   calling argument
```

```
IMPLICIT NONE
```

```
CHARACTER*1      esc
CHARACTER*1      code
```

```
CHARACTER*32     status
```

```
INTEGER*2        lun
INTEGER*2        io_chan
```

```
INTEGER*2        stat
INTEGER*2        option
```

```
C*****
C      variable initializations      *
```

```
      esc = CHAR(27)
      code = 'G'
```

```
C      ...initiate zero stabilization set on 8232
```

```
      WRITE(lun,1000)esc,option,code
```

```
C      ...read status report from 8232
```

```
      CALL ci8232_status(lun,io_chan,stat,status)
```

```
C*****
C      format declarations      *
```

```
1000      FORMAT(1X,A1,I1,A1,*)
```

```
.....
```

```
      RETURN
```

END

```
////////////////////  
C,////////////////////  
  
SUBROUTINE ci8232_STORE(lun,io_chan,stat,status)
```

```
C*****  
C  
C Title : ci8232_store  
C Revision : 1.00  
C
```

```
C*****  
C  
C Abstract : this routine stores the current settings to NVRAM  
C on the Canberra 8232 stabilizer  
C
```

```
C*****  
C  
C Arguments : lun I I*2 logical unit for serial  
C communication with 8232  
C io_chan I I*2 I/O channel for serial  
C communication with 8232  
C stat 0 I*2 status code(0/1 for OK/error)  
C status 0 C*32 status string  
C
```

```
C*****  
C  
C Routines  
C Called : ci8232_status  
C  
C Libraries  
C Used : none  
C
```

```
C*****  
C*****  
C  
C C*1 esc ASCII code for escape  
C C*1 code 1 character 8232 code for store to NVRAM  
C  
C C*32 status calling argument  
C  
C I*2 lun calling argument  
C I*2 io_chan calling argument  
C I*2 stat calling argument  
C
```

IMPLICIT NONE

```
CHARACTER*1 esc  
CHARACTER*1 code  
  
CHARACTER*32 status  
  
INTEGER*2 lun  
INTEGER*2 io_chan
```

INTEGER\*2 stat

C\*\*\*\*\*  
C variable initializations \*

esc = CHAR(27)  
code = 'M'

C ...initiate store to NVRAM on 8232

WRITE(lun,1000)esc,code

C ...read status report from 8232

CALL ci8232\_status(lun,io\_chan,stat,status)

C\*\*\*\*\*  
C format declarations \*

1000 FORMAT(1X,A1,A1,#)

C\*\*\*\*\*  
RETURN  
END

C///  
C///

SUBROUTINE ci8232\_REPORT(lun,io\_chan,report,stat,status)

C\*\*\*\*\*  
C  
C Title : ci8232\_report  
C Revision : 1.00

C\*\*\*\*\*  
C  
C Abstract : this routine retrievees setup and status information  
C on the Canberra 8232 stabalizer

C\*\*\*\*\*  
C  
C Arguments : lun I I\*2 logical unit for serial  
C communication with 8232  
C io\_chan I I\*2 I/O channel for serial  
C communication with 8232  
C report O struct status and setup report

```
C          stat      0          I*2      status code(0/1 for OK/error)
C          status    0          C*32      status string
C
```

```
C*****
```

```
C          Routines
L          Called      :      timed_read, ci8232_status
C
```

```
C          Libraries
C          Used        :      none
C
```

```
C*****
C*****
```

```
C          C*1      esc          ASCII code for escape
C          C*1      code        1 character 8232 code for store to NVRAM
C
C          C*32     status       calling argument
C
C          I*2      lun          calling arguemnt
C          I*2      io_chan      calling arguemnt
C          I*2      stat         calling argument
C
C          struct   record       calling argument
C
C          I*2      length       length of port read
C          I*2      clength      anticipated length of port read
C          C*120    string       character buffer for port read
C          I*2      ier          return status of timed read
```

```
IMPLICIT NONE
```

```
CHARACTER*1      esc
CHARACTER*1      code
```

```
CHARACTER*32     status
```

```
INTEGER*2        lun
INTEGER*2        io_chan
INTEGER*2        stat
```

```
STRUCTURE /report_struct/
  INTEGER*2 gain_peak
  INTEGER*2 gain_window
  INTEGER*2 gain_range
  INTEGER*2 gain_stab
  INTEGER*2 gain_rate
  INTEGER*2 gain_corr
  INTEGER*2 zero_peak
  INTEGER*2 zero_window
  INTEGER*2 zero_range
  INTEGER*2 zero_stab
  INTEGER*2 zero_rate
  INTEGER*2 zero_corr
END STRUCTURE
```

```
RECORD /report_struct/ report
```

```
INTEGER*2        length
```

```
INTEGER*2      clength
INTEGER*2      ier

CHARACTER*120  string
```

```
*****
L   variable initializations
C*****

      esc = CHAR(27)

C   ...initiate gain setup report on 8232

      code = 'H'
      WRITE(lun,1000)esc,code

      clength = 15
      length = clength
      CALL timed_read(io_chan,string,length,ier)
      type *,string(2:clength)
      IF((length.LT.clength-2).OR.ier.NE.0)THEN
        type *,length,clength,ier
        stat = 12
        status = 'error in QIO read'
        GOTO 500
      ENDIF
      READ(string,1001,end=400,err=450)      report.gain_peak,
1                                           report.gain_window,
2                                           report.gain_range, code

C   ...read status report from 8232
      CALL ci8232_status(lun,io_chan,stat,status)

C   ...check for error from status report
      IF(stat.NE.0)THEN
        GOTO 500
      ENDIF

C   ...check for error in response from 8232
      IF(code.NE.'H')THEN
        stat = 7
        status = 'Error in 8232 gain setup report'
        GOTO 500
      ENDIF

C   ...initiate zero setup report on 8232

      code = 'J'
      WRITE(lun,1000)esc,code

      clength = 15
      length = clength
      CALL timed_read(io_chan,string,length,ier)
      type *,string(2:clength)
      IF((length.LT.clength-2).OR.ier.NE.0)THEN
        type *,length,clength,ier
        stat = 12
        status = 'error in QIO read'
        GOTO 500
```



```

        ENDIF
        READ(string,1001,end=400,err=450)          report.zero_peak,
1                                               report.zero_window,
2                                               report.zero_range, code

...read status report from 8232
    CALL ci8232_status(lun,io_chan,stat,status)

C ...check for error from status report
    IF(stat.NE.0)THEN
        GOTO 500
    ENDIF

C ...check for error in response from 8232
    IF(code.NE.'J')THEN
        stat = 7
        status = 'Error in 8232 zero setup report'
        GOTO 500
    ENDIF

C ...initiate gain status report on 8232

    code = 'I'
    WRITE(lun,1000)esc,code

    clength = 15
    length = clength
    CALL timed_read(io_chan,string,length,ier)
    type *,string(2:clength)
    IF((length.LT.clength-2).OR.ier.NE.0)THEN
        type *,length,clength,ier
        stat = 12
        status = 'error in QIO read'
        GOTO 500
    ENDIF
    READ(string,1002,end=400,err=450)          report.gain_stab,
1                                               report.gain_rate,
2                                               report.gain_corr, code

C ...read status report from 8232
    CALL ci8232_status(lun,io_chan,stat,status)

C ...check for error from status report
    IF(stat.NE.0)THEN
        GOTO 500
    ENDIF

C ...check for error in response from 8232

    IF(code.NE.'I')THEN
        stat = 7
        status = 'Error in 8232 gain status report'
        GOTO 500
    ENDIF

L ...initiate zero status report on 8232

    code = 'K'

```

```
WRITE(lun,1000)esc,code
```

```
clength = 15
```

```
length = clength
```

```
CALL timed_read(io_chan,string,length,ier)
```

```
type *,string(1:length)
```

```
IF((length.LT.clength-2).OR.ier.NE.0)THEN
```

```
type *,length,clength,ier
```

```
stat = 12
```

```
status = 'error in QIO read'
```

```
GOTO 500
```

```
ENDIF
```

```
READ(string,1002,end=400,err=450)
```

```
report.zero_stab,
```

```
report.zero_rate,
```

```
report.zero_corr, code
```

```
1
```

```
2
```

```
C ...read status report from 8232
```

```
CALL ci8232_status(lun,io_chan,stat,status)
```

```
C ...check for error from status report
```

```
IF(stat.NE.0)THEN
```

```
GOTO 500
```

```
ENDIF
```

```
C ...check for error in response from 8232
```

```
IF(code.NE.'K')THEN
```

```
stat = 7
```

```
status = 'Error in 8232 zero status report'
```

```
GOTO 500
```

```
ENDIF
```

```
GOTO 500
```

```
400 CONTINUE
```

```
stat = 10
```

```
status = '8232 not reachable, check port'
```

```
GOTO 500
```

```
450 CONTINUE
```

```
stat = 11
```

```
status = 'Error in 8232 communication'
```

```
C ...label for error return
```

```
500 CONTINUE
```

```
C*****
```

```
C format declarations *
```

```
C*****
```

```
1000 FORMAT(1X,A1,A1,#)
```

```
1001 FORMAT(1X,I5,1X,I3,1X,I1,A1)
```

```
1002 FORMAT(1X,I1,1X,I5,1X,I3,A1)
```

```
C*****
```

```
RETURN
```

```
END
```

```
0001 SUBROUTINE AVE (NS,NE,SUM,AV,Y)
0002
0003 C
0004 C
0005 C AVE
0006 C
0007 C
0008 C
0009 C Version: 1
0010 C Author: Ray Gunnink
0011 C Date: circa 1974
0012 C
0013 C * * * * *
0014 C PURPOSE
0015 C To determine the average value of a set of data counts taken fro
0016 C channels NS to NE (inclusive) in array Y.
0017 C
0018 C ARGUMENTS
0019 C
0020 C INPUT
0021 C NS Starting channel of integration.
0022 C NE Ending channel of integration.
0023 C Y Data array
0024 C
0025 C RETURN
0026 C SUM Integrated count
0027 C AVE Average count per channel
0028 C
0029 C
0030 C
0031 C
0032 C DIMENSION Y(2)
0033 C
0034 C SUM = 0.0
0035 C SM = NE - NS + 1
0036 C DO 2 I = NS,NE
0037 C 2 SUM = SUM + Y(I)
0038 C AV = SUM / SM
0039 C RETURN
0040 C END
```



```

0058 C      SUMY = Variable used to accumulate the net count using
0059 C      the straight line for background determination.
0060 C      SUM = The final value of SUMY following the 1st pass.
0061 C      YN = Net count of a channel (temporary argument).
0062 C      YNT= Positive value only of YN (temporary argument).
0063 C      SLBG = Straight line background value at a given channel.
0064 C
0065 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0066 C
0067
0068 REAL*4 S2,AV1,AV2,Y1(2),Y2(2)
0069 INTEGER*2 NS,NE,ND
0070 C
0071
0072 DBG = AV2 - AV1           !Change in background level.
0073 ND = 0
0074 S2 = 0.0
0075 CHANS = NE - NS + 1      !Total channels in group.
0076 SLOPE = DBG / CHANS      !Average change in background per cha
0077 AVSIG = 5./CHANS*SQRT(AV2*CHANS) !Smoothing factor
0078 NN = 2
0079 IF (SLOPE .GT. 0.0) NN = 1 !One pass for positive sloping bkgr
0080
0081 DO 30 J = 1,NN
0082 SUMY = 0.0
0083 SLBG = AV1               !Initial straight line bkgrd value set a
0084
0085 DO 20 I = NS,NE
0086 SLBG = SLBG + SLOPE      !New straight line bkgrd value.
0087 YN = Y1(I) - SLBG       !Calc. net count.
0088 YNT = YN
0089 IF (YN .LT. 0.0) YNT = 0.0 !Only integrate positive values.
0090 SUMY = SUMY + YNT + AVSIG !Running sum of net counts.
0091 IF (SLOPE .GT. 0.0) GOTO 10 !Straight line only on + slope
0092 IF (J .EQ. 1) GOTO 20    !Straight line only on 1st pass
0093 YN = Y1(I) - AV1 - DBG*SUMY/SUM !Calc. smoothed-step bkgrd value
0094
0095 10 ND = ND + 1
0096 Y2(ND) = YN              !Store net channel count.
0097 S2 = S2 + YN             !Integrate net channel counts.
0098 20 CONTINUE
0099
0100 SUM = SUMY               !Store total net counts using straight line bkgr
0101 IF (SUMY .EQ. 0.)RETURN
0102 30 CONTINUE
0103
0104 RETURN
0105 END

```



```
0058 C nature of the data. ie use halves of 2 adjacen
0059 C channels.
0060 C
0061 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0062
0063 REAL*4 BGLO,BGHI,SLPLO,SLPHI,SMY,Y(2),YNET(2)
0064 INTEGER*2 NS,NE,NPTS
0065
0066 NPTS = NE-NS+1 !Number of channels in grouping.
0067 XPTS = NPTS
0068 DBG = (BGHI-BGLO) -0.5 *(SLPLO+SLPHI)*(XPTS+1.) !Level change in
0069 SLOPE = DBG / (XPTS + 1.) !Average bkgrd slope in counts/channe
0070 SLBG = BGLO !Starting (low-side) bkgrd level.
0071 SUMY1 = 0.0
0072 YSUM2 = 0.0
0073
0074 SIG = 5.*SQRT(BGLO) !Variable related to statistical level of
0075 AVSIG = SIG/XPTS
0076
0077 L = 0
0078 SLPBG = 0.0
0079
0080 DO 10 I = NS,NE !Get 1st estimate of net area of peak(s)
0081 FRAC = L / XPTS !Fractional distance through grouping.
0082 SLPBG = SLPBG+SLPLO*(1.-FRAC)+SLPHI*FRAC !Amt due to sloping bk
0083 SLBG = SLBG + SLOPE !Straight line background level.
0084 L = L + 1
0085 YN = Y(I) -SLBG -SLPBG !Net count in channel using straight lin
0086 IF(YN .LT. 0.) YN = 0. !Only accumulate positive values.
0087 10 SUMY1 = SUMY1 + YN !Approximate area using straight line.
0088 SUMY1 = SUMY1 + SIG !Add value related to statistical varian
0089 NN = 2
0090 IF (SLOPE .GE. 0.0) NN = 1 !Only 1 pass if slope is positive.
0091 DO 30 J = 1,NN !1st uses straight line. 2nd also uses step.
0092 L = 0
0093 SLBG = BGLO !Start at initial bkgrd level.
0094 YSUM1 = 0.0
0095 SLPBG = 0.0
0096 YS1 = 0.5 * BGLO !Data is histogram. Use half of 2 adjacent ch
0097 YS3 = YS1
0098 SMY = 0.0
0099
0100 DO 20 I = NS,NE
0101 SLBG = SLBG + SLOPE !Straight line background.
0102 FRAC = L / XPTS !Fractional distance through gro
0103 SLPBG = SLPBG +SLPLO *(1.-FRAC)+SLPHI *FRAC !Background slope.
0104 YN = Y(I) - SLBG - SLPBG !Net count using straight line.
0105 YNT = YN
0106 IF (YN .LT. 0.) YNT = 0. !Only integrate positive values.
0107 L = L + 1
0108
0109 IF (DBG .GE. -10.)THEN
0110 IF ((J.EQ.1) .AND.(NN.EQ.2) ) GOTO 15 !Use straight li
0111 GOTO 14
0112 ENDIF
0113 YSUM1 = YSUM1 + YNT + AVSIG !Partial area up to channel "I".
0114 YS2 = 0.5 * (BGLO + DBG * YSUM1/SUMY1) !Other half.
```

```

0115      YN = Y(I) -YS1 -YS2 -SLPBG      !Net count using straight line b
0116      IF (J .EQ. 1) GOTO 15          !Skip following on first iterati
0117
0118      YSUM2 = YSUM2+YN+AVSIG      !Partial area to channel "I" using ste
0119      YS4 = 0.5 * (BGLO + DBG * YSUM2/SUMY2)      !Other half.
0120      YN = Y(I) -YS3 -YS4 -SLPBG      !Net count resulting from refine
0121      YS3 = YS4
0122      14      YNET(L) = YN          !Transfer to YNET array on last iteratio
0123      15      SMY = SMY + YN        !Net counts (or area).
0124      YS1 = YS2
0125      20      CONTINUE
0126
0127      SUMY2 = SMY + SIG            !New value of SUMY2 for 2nd pass.
0128      30      CONTINUE
0129
0130      RETURN
0131      END

```

## PROGRAM SECTIONS

Name	Bytes	Attributes
0 #CODE	460	PIC CON REL LCL SHR EXE R
2 #LOCAL	84	PIC CON REL LCL NOSHR NDEXE R
Total Space Allocated	544	

## ENTRY POINTS

Address	Type	Name
0-00000000		BKGRD

## VARIABLES

Address	Type	Name	Address	Type	Name	Address	T
**	R*4	AVSIG	AP-00000010@	R*4	BGHI	AP-00000000@	
**	R*4	FRAC	**	I*4	I	2-00000010	
AP-00000008@	I*2	NE	2-00000018	I*4	NN	AP-00000010@	
2-00000014	R*4	SIG	**	R*4	SLBG	2-00000008	
AP-00000018@	R*4	SLPHI	AP-00000014@	R*4	SLPLO	AP-00000020@	
2-00000028	R*4	SUMY2	2-00000000	R*4	XPTS	**	
**	R*4	YS1	2-00000024	R*4	YS2	2-00000020	
**	R*4	YSUM1	2-00000010	R*4	YSUM2		



```

0001      C      FUNCTION ER
0002
0003      FUNCTION CALERR(XY,ALPHA,NTERMS,CHISQR)
0004      DOUBLE PRECISION ALPHA(3,3)
0005
0006
0007      YR=0.
0008      DO 36 I=1,NTERMS
0009      F1 = 1.
0010      IF(I.EQ.1)GOTO 30
0011      F1=XY**(I-1)
0012      30      DO 36 J=1,NTERMS
0013      F2 = 1.
0014      IF(J.EQ.1)GOTO 34
0015      F2=XY**(J-1)
0016      34      YR=YR+ALPHA(I,J)*F1*F2
0017      36      CONTINUE
0018      CALERR=SQRT(CHISQR*ABS(YR))
0019      RETURN
0020      END

```

PROGRAM SECTIONS

Name	Bytes	Attributes
0 #CODE	166	PIC CON REL LCL SHR EXE R
2 #LOCAL	24	PIC CON REL LCL NOSHR NOEXE R
Total Space Allocated	190	

ENTRY POINTS

Address	Type	Name
0-00000000	R*4	CALERR

VARIABLES

Address	Type	Name	Address	Type	Name	Address	T
AP-00000010@	R*4	CHISQR	**	R*4	F1	**	
**	I*4	J	AP-0000000C@	I*4	NTERMS	AP-00000004@	

ARRAYS

Address	Type	Name	Bytes	Dimensions
0-0000000B@	R*8	ALPHA	72	(3, 3)

```

0001 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0002 C
0003 C      Generalized GFIT
0004 C      25-MAR-80; REV 23-JUL-82 by W.D.Ruhter
0005 C      Edited 28-Feb-83 by J.B.Niday (removed index locn to caller)
0006 C      and 24-May-83(prevent crash on neg. Y values)
0007 C      Edited 25-Jun-87 by R. Gunnink (Added out-of-bound check)
0008 C
0009 C      Fit Gaussian to upper half of a peak after stripping out the
0010 C      tailing contribution. Return peak posn, height, and width. Return
0011 C      negative error flag if no valid peak in region.
0012 C      L      Approx index of peak posn in data array Y
0013 C      M      No. of points beyond L to include
0014 C      ALFA   Peak width parameter calculated (an approx received)
0015 C      EXP1   Short term tail amplitude
0016 C      EXP2   Short term tail slope
0017 C      PKHT   Peak height to be returned.
0018 C      ER     Uncertainty in peak position; unless its value is -.
0019 C             Then ER is a flag (-1= Parabola nearly degenerate.
0020 C             -2= Peak position out of bounds.)
0021 C      Y      Data array received
0022 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0023 C
0024 C
0025 C      SUBROUTINE GFIT(L,M,ALFA,EXP1,EXP2,PKPOS,PKHT,ER,Y)
0026 C      DIMENSION Y(2)
0027 C
0028 C      Initializations
0029 C      PKPOS = L      !In case of early return
0030 C      ER = 0
0031 C      SY = 0.0
0032 C      SYY=0.
0033 C      SXY=0.
0034 C      SXXY=0.
0035 C      C=M*(M+1)/3.
0036 C      Y1 = Y(L+1)
0037 C      Y2 = Y(L-1)
0038 C      TPOS = .5 * (Y2-Y1) / (Y1 - 2. * Y(L) + Y2)
0039 C      Y0 = Y(L)
0040 C
0041 C      Strip out the short term tail contribution to peak
0042 C      I=-M
0043 C      J=L+I
0044 C      XL=I
0045 C      DLX = XL - TPOS
0046 C      Y1 = Y(J)
0047 C      IF (DLX .GE. 0) GO TO 34
0048 C      Y1 =Y1-Y0 * EXP1* EXP(EXP2 * DLX) * (1. -EXP(.4*ALFA*DLX*DLX))
0049 C      34 IF(Y1 .GT. 0)GO TO 37
0050 C
0051 C      TYPE 1034,(Y(K),K=L-M,L+M)
0052 C      1034 FORMAT(' A negative data value was incountered in GFIT.'/
0053 C      + ' Its value will be set to 0.1. Data Values ='/1P(5E12.4))
0054 C
0055 C      Y1=Y(J)
0056 C      IF(Y1 .LE. 0.)Y1 = 0.1
0057 C

```

```

0058 C      Take natural log of net 'Y' data. Result should be a parabola.
0059 C      Fit to a parabola to get peak height, position, and width.
0060
0061 37      Y1=ALOG(Y1)
0062        SY=SY+Y1
0063        SYI=SYI+Y1**2
0064        SXY=SXY+XL*Y1
0065        SXXY=SXXY+(XL**2-C)*Y1
0066        I=I+1
0067        IF(I.LE.M)GOTO 30
0068
0069 C      COMPUTE THE DIAGONAL ELEMENTS OF X'X
0070        XPX11=2*M+1.
0071        XPX22=XPX11*C
0072        XPX33=XPX22*(4*M*(M+1)-3)/15.
0073 C      COMPUTE LEAST SQUARES ESTIMATES OF COEFFICIENTS
0074        A0=SY/XPX11
0075        A1=SXY/XPX22
0076        ALFA=SXXY/XPX33
0077 C      IS THE PARABOLA NEARLY DEGENERATE
0078        ER = -1.                !Default value
0079        IF(ABS(A0) .LT. 1.E-06) GOTO 50
0080
0081 C      COMPUTE THE COORDINATES OF THE VERTEX
0082        PKPOS=-0.5*A1/ALFA
0083        PKHT=A0+0.5*A1*PKPOS-ALFA*C
0084        PKPOS = PKPOS + L
0085        IF ((PKPOS .LT. (L-M)) .OR. (PKPOS .GT. (L+M))) GOTO 40 !Check b
0086        ER=SY-A0*SY-A1*SXY-ALFA*SXXY
0087        IF(M.GT.1)ER=ER/(2*(M-1))
0088        PKHT=EXP(PKHT)
0089        GOTO 60
0090
0091 40      ER = -2.                !Set flag
0092 50      PKPOS = L                !Return starting positio
0093        PKHT = Y(L)              !Return default height
0094        TYPE 9000, ER
0095        PRINT 9000, ER
0096 9000   FORMAT (' ***Error in GFIT***  Error flag = ',F3.0)
0097
0098 60      RETURN
0099        END

```

```

0001 C      GPM
0002 C
0003 C      Last edit Jun-12-1986 (J.B.N.)
0004 C      Test effect of double precision.
0005
0006 SUBROUTINE GPM (EPK,ELOG,GEOM,SURFC,DEPTH,SMU,VALUE,GAMAS,DCNST)
0007
0008 IMPLICIT DOUBLE PRECISION(A-H,O-Z)
0009 REAL*4 DCNST(2),EPK,ELOG,GEOM,SURFC,DEPTH,SMU,VALUE,GAMAS
0010
0011 EXOVR = DCNST(12)
0012 DETHT = DCNST(13)
0013 DEP = 0.5 * DEPTH
0014 RSQ = DCNST(14) * DCNST(14)
0015 XMU = -2.316 + 4.2 * EXP(-.478 * ELOG - 1.434)
0016 XMU = EXP (XMU)
0017 XM = 1.0 / XMU
0018 SEE = EXP (-DETHT * XMU)
0019 XM = XM * (1. - XMU * DETHT * SEE / (1. - SEE))
0020 GM = GEOM + XM
0021 GG = GM * GM
0022
0023 C      ROSQ = -GG + 0.5
0024 C      IF(RSQ .GT. 0.000001) THEN
0025         ROSQ = -GG + .5 * RSQ / (1. - GM/SQRT (GG+ RSQ))
0026 C      ENDIF
0027
0028 ERSQ = ROSQ * (1.1 - XM / SQRT (DETHT))
0029 GGSQ = GG + ERSQ
0030 GLOST = 1.
0031 GCORR = 1.
0032 SQGCR = 1.
0033 SATN = 1.
0034 IF (SURFC) 21,31,21
0035 21 RR = SURFC / 3.1416
0036 GM2 = GM - DEP * (1. - EXP (-SMU))
0037 GG2 = GM2 * GM2
0038 RR = RR * (1.4 - .4 / GM2)
0039 GCORR = RR / (GG2 * LOG ((GG2 + RR)/GG2))
0040 SQGCR = SQRT (GCORR)
0041 SMU = SMU * SQGCR
0042 UL = SMU * SMU
0043 TOLGA = 2. * DEP / GM
0044 ADD = 1.
0045 DFACT = 1.
0046 ASQ = GGSQ * GCORR
0047 TOTAL = (ASQ / (ASQ -DEP*DEP)) + .1 * SMU * (TOLGA **3)
0048 DO 25 I = 1,10
0049     XK = 2 * I
0050     DFACT = DFACT * XK * (XK + 1.)
0051     ADD = UL * ADD
0052     AD = ADD / DFACT
0053     IF (TOTAL - AD * 200.) 24,24,27
0054 24 TOTAL = TOTAL + AD
0055 25 CONTINUE
0056 27 TOTAL = TOTAL + AD
0057 ADD = SMU

```

```
0058          DFACT = 1.
0059          DO 28 I = 1,5
0060             XK = 2 * I - 1
0061             AD = TOLOA * ADD / ((XK + 2.) * DFACT)
0062             IF (TOTAL - AD *200.) 29,29,30
0063 29          TOTAL = TOTAL + AD
0064             ADD = ADD * UL
0065             DFACT = DFACT * (XK + 1.) * (XK + 2.)
0066 28          CONTINUE
0067 30          TOTAL = TOTAL + AD
0068             SATN = EXP (-SMU)
0069             GLOST = TOTAL * SATN
0070
0071 31          GM = GM - DEP * (1. - SATN)
0072             SGSQ = SQRT (GM * GM + ROSQ)
0073             SMU = DCNST(16) * (SGSQ - GM) / (SGSQ + GM)
0074             N = 1
0075             M = 6
0076             IF (EPK - EXOVR) 32,32,34
0077 32          N = 7
0078             M = 11
0079 34          EFLOG = DCNST(N)
0080             L = N+1
0081             DO 36 I = L,M
0082 36          EFLOG = EFLOG + DCNST(I) * ELOG ** (I-N)
0083             EFF = EXP (EFLOG)
0084             GAMAS = VALUE * GGSQ * GCORR / (EFF * GLOST)
0085             RD = .7 * DCNST(14)
0086             VALUE = (1. + RD / GG) * SQGCR
0087             RETURN
0088             END
```

```
0001          SUBROUTINE MATINV (ARRAY, MX, NORDER, DET)
0002
0003          CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0004          C      MATINV modified to carry array dimension
0005          C      Last edit 16-Jun-82
0006          CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0007
0008          DIMENSION IK(50), JK(50)
0009          DOUBLE PRECISION ARRAY (MX, MX), AMAX, SAVE
0010
0011          10      DET = 1.
0012          11      DO 100 K = 1, NORDER
0013          C
0014          C      FIND LARGEST ELEMENT ARRAY (I, J) IN REST OF MATRIX
0015          C
0016          AMAX = 0.
0017          21      DO 30 I = K, NORDER
0018          DO 30 J = K, NORDER
0019          23      IF (ABS (AMAX) - ABS (ARRAY (I, J))) 24, 24, 30
0020          24      AMAX = ARRAY (I, J)
0021          IK (K) = I
0022          JK (K) = J
0023          30      CONTINUE
0024          C
0025          C      INTERCHANGE ROWS AND COLUMNS TO PUT AMAX IN ARRAY (K, K)
0026          C
0027          31      IF (AMAX) 41, 32, 41
0028          32      DET = 0.
0029          29      GOTO 140
0030          41      I = IK (K)
0031          IF (I - K) 21, 51, 43
0032          43      DO 50 J = 1, NORDER
0033          SAVE = ARRAY (K, J)
0034          ARRAY (K, J) = ARRAY (I, J)
0035          50      ARRAY (I, J) = -SAVE
0036          51      J = JK (K)
0037          IF (J - K) 21, 61, 53
0038          53      DO 60 I = 1, NORDER
0039          SAVE = ARRAY (I, K)
0040          ARRAY (I, K) = ARRAY (I, J)
0041          60      ARRAY (I, J) = -SAVE
0042          C
0043          C      ACCUMULATE ELEMENTS OF INVERSE MATRIX
0044          C
0045          61      DO 70 I = 1, NORDER
0046          IF (I - K) 63, 70, 63
0047          63      ARRAY (I, K) = -ARRAY (I, K) / AMAX
0048          70      CONTINUE
0049          71      DO 80 I = 1, NORDER
0050          DO 80 J = 1, NORDER
0051          IF (I - K) 74, 80, 74
0052          74      IF (J - K) 75, 80, 75
0053          75      ARRAY (I, J) = ARRAY (I, J) + ARRAY (I, K) * ARRAY (K, J)
0054          80      CONTINUE
0055          81      DO 90 J = 1, NORDER
0056          IF (J - K) 83, 90, 83
0057          83      ARRAY (K, J) = ARRAY (K, J) / AMAX
```

```

0058 90      CONTINUE
0059      ARRAY(K,K) = 1. / AMAX
      50     100    CONTINUE
0061      C
0062      C      RESTORE ORDERING OF MATRIX
0063      C
0064      101    DO 130 L =1,NORDER
0065          K = NORDER - L + 1
0066          J = IK(K)
0067          IF(J-K)111,111,105
0068      105    DO 110 I =1,NORDER
0069          SAVE = ARRAY(I,K)
0070          ARRAY(I,K) = -ARRAY(I,J)
0071      110    ARRAY(I,J) = SAVE
0072      111    I = JK(K)
0073          IF(I-K)130,130,113
0074      113    DO 120 J =1,NORDER
0075          SAVE = ARRAY(K,J)
0076          ARRAY(K,J) = -ARRAY(I,J)
0077      120    ARRAY(I,J) = SAVE
0078      130    CONTINUE
0079      140    RETURN
0080          END

```

## PROGRAM SECTIONS

Name	Bytes	Attributes
0 \$CODE	693	PIC CON REL LCL SHR EXE R
2 \$LOCAL	512	PIC CON REL LCL NOSHR NOEXE R
Total Space Allocated:	1205	

## ENTRY POINTS

Address	Type	Name
0-00000000		MATINV

## VARIABLES

Address	Type	Name	Address	Type	Name	Address	T
2-00000190	R*8	AMAX	AP-00000010@	R*4	DET	**	
2-00000198	I*4	K	2-000001BC	I*4	L	AP-0000000B@	
**	R*8	SAVE					

```
0001          SUBROUTINE MAXVAL (NS,NE,JMAX,YMAX,Y)
0002
0003          C
0004          C
0005          C          MAXVAL
0006          C
0007          C
0008          C
0009          C          Version: 1
0010          C          Author: Ray Gunnink
0011          C          Date:   circa 1974
0012          C
0013          C          * * * * *
0014          C
0015          C          PURPOSE
0016          C
0017          C          To determine the position and value of the maximum value (count)
0018          C          in a set of data points from NS to NE inclusive.
0019          C
0020          C          *          *          *          *          *          *          *
0021          C
0022          C          ARGUMENTS
0023          C
0024          C          INPUT
0025          C          NS          Starting channel for search
0026          C          NE          Ending channel for search
0027          C          Y          Data array
0028          C
0029          C          RETURN
0030          C          JMAX          Position of maximum value
0031          C          YMAX          Maximum value
0032          C
0033          C
0034          C
0035          C
0036          C
0037          C          REAL*4 YMAX,Y(2)
0038          C          INTEGER*2 NS,NE,JMAX
0039          C
0040          C          YMAX = 0.0
0041          C
0042          C          DO 2 I = NS,NE
0043          C          IF (Y(I) .LE. YMAX) GOTO 2
0044          C          YMAX = Y(I)
0045          C          JMAX = I
0046          C          2 CONTINUE
0047          C
0048          C          RETURN
0049          C          END
```



```

0001
0002 C      MINVAL
0003 C      THIS ROUTINE RETURNS THE POSITION AND VALUE OF THE MINIMUM VALUE
0004 C      IN A GROUP OF DATA POINTS FROM NS TO NE INCLUSIVE.
0005      SUBROUTINE MINVAL (NS,NE,JMIN,YMIN,Y)
0006      DIMENSION Y(2)
0007      YMIN=1.0E38
0008      DO 4 I=NS,NE
0009      IF (Y(I)-YMIN)2,4,4
0010  2      JMIN=I
0011      YMIN=Y(I)
0012  4      CONTINUE
0013      RETURN
0014      END

```

PROGRAM SECTIONS

Name	Bytes	Attributes
0 #CODE	72	PIC CON REL LCL SHR EXE R
2 #LOCAL	24	PIC CON REL LCL NOSHR NOEXE R
Total Space Allocated	96	

TRY POINTS

Address	Type	Name
0-00000000		MINVAL

VARIABLES

Address	Type	Name	Address	Type	Name	Address	T
**	I*4	I	AP-00000000@	I*4	JMIN	AP-00000008@	
AP-00000010@	R*4	YMIN					

ARRAYS

Address	Type	Name	Bytes	Dimensions
AP-00000014@	R*4	Y	8	(2)

LABELS

Address	Label	Address	Label
**	2	0-00000043	4

```

0001 C NCNTS VERSION 0 19-JAN-78 AT 12:30
0002
0003 C THIS ROUTINE IS SIMILAR TO "NCTS" ROUTINE EXCEPT THAT IT RETURNS
0004 C AN AVERAGE BACKGROUND AND THE ERROR ASSOCIATED WITH THE TOTAL NET
0005 C COUNT.
0006
0007 SUBROUTINE NCNTS (N1,N2,M1,M2,ND,S1,S2,ER,AVEBG,Y1,Y2)
0008 DIMENSION Y1(2),Y2(2)
0009
0010 CALL AVE (N1,N2,S1,AV1,Y1)
0011 CALL AVE (M1,M2,S2,AV2,Y1)
0012 S1 = S1 + S2
0013 CALL BKGD (N2+1,M1-1,AV1,AV2,S2,ND,Y1,Y2)
0014 CHANS = ND
0015 XPTS = M2 + N2 - M1 - N1 + 1
0016 AVEBG = 0.5 * (AV1 + AV2)
0017 IF (S2 .LT. 0.0) S2 = 0.0
0018 ER = SQRT (S2 + AVEBG * CHANS * (1. + CHANS / XPTS))
0019 RETURN
0020 END
  
```

PROGRAM SECTIONS

Name	Bytes	Attributes
1 \$CODE	213	PIC CON REL LCL SHR EXE R
2 \$LOCAL	132	PIC CON REL LCL NOSHR NOEXE R
Total Space Allocated	345	

ENTRY POINTS

Address	Type	Name
0-00000000		NCNTS

VARIABLES

Address	Type	Name	Address	Type	Name	Address	T
2-00000000	R*4	AV1	2-00000004	R*4	AV2	AP-00000024@	
AP-00000020@	R*4	ER	AP-0000000C@	I*4	M1	AP-00000010@	
AP-00000008@	I*4	N2	AP-00000014@	I*4	ND	AP-00000018@	
**	R*4	XPTS					

```

0001          SUBROUTINE RDBLK(ISIZE,IBUFF,IRC,ICHAN)
0002
0003          CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0004          C          READ ONE BLOCK OF DIRECT ACCESS FILE
0005          CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0006
0007          INTEGER*2 IBUFF(1)
0008          READ(ICHAN,IRC,ERR=95)(IBUFF(I),I=1,ISIZE)
0009          RETURN
0010          95      TYPE 100,IRC,ISIZE,ICHAN
0011          100    FORMAT (' ERROR READING RECORD:',3I5)
0012          IF(ISIZE .LE. 256)STOP ' ATTEMPTED .GT. 256'
0013          STOP ' ATTEMPT TO PASS EOF'
0014          END

```

PROGRAM SECTIONS

Name	Bytes	Attributes
0 \$CODE	154	PIC CON REL LCL SHR EXE R
1 \$PDATA	68	PIC CON REL LCL SHR NOEXE R
2 \$LOCAL	48	PIC CON REL LCL NOSHR NOEXE R
Total Space Allocated	270	

ENTRY POINTS

Address	Type	Name
0-00000000		RDBLK

VARIABLES

Address	Type	Name	Address	Type	Name	Address	T
**	I*4	I	AP-00000010@	I*4	ICHAN	AP-00000000@	

ARRAYS

Address	Type	Name	Bytes	Dimensions
AP-00000008@	I*2	IBUFF	2	(1)

LABELS

Address	Label	Address	Label
0-00000044	95	1-00000027	100'

```

0001 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0002 C
0003 C      VERSION 0 19-JAN-78
0004 C
0005 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0006 C
0007 C      SPRAM
0008
0009 C      THIS SUBROUTINE CALCULATES SOME OF THE PARAMETERS USED TO
0010 C      CHARACTERIZE THE SHAPE OF A PEAK
0011
0012 SUBROUTINE SPRAM (EN,GAIN,SHAPC,SG,ALFA,F12M,EXP1)
0013
0014 DIMENSION SHAPC(2)
0015
0016 GSQ = GAIN * GAIN
0017 SG = (SHAPC(1) + EN * SHAPC(2)) / GSQ + 0.462
0018 ALFA = -2.7726 / SG
0019 F12M = SQRT (SG)
0020 EXP1 = EXP (EN * SHAPC(4) + SHAPC(3))
0021 RETURN
0022 END

```

PROGRAM SECTIONS

Name	Bytes	Attributes
0 \$CODE	101	PIC CON REL LCL SHR EXE R
2 \$LOCAL	20	PIC CON REL LCL NOSHR NOEXE R
Total Space Allocated	121	

ENTRY POINTS

Address	Type	Name
0-00000000		SPRAM

VARIABLES

Address	Type	Name	Address	Type	Name	Address	Type	Name
AP-00000014@	R*4	ALFA	AP-00000004@	R*4	EN	AP-00000010@		
AP-00000008@	R*4	GAIN	**	R*4	GSQ	AP-00000010@		

SUBROUTINE STRIP1 (NS,NE,PKPOS,PKHT,ALFA,TLAMP,TLSLP,WL,Y)

CC

STRIP1

Version: 1  
Author: Ray Gunnink  
Date: circa 1978

\*\*\*\*\*

PURPOSE  
To strip a peak or portion of a peak from a designated set of  
spectrum channels. In contrast to routine STRIP, this routine c  
also calculate X ray line shapes for stripping.

\* \* \* \* \*

ARGUMENTS

INPUT

NS Starting position for stripping  
NE Ending position for stripping  
PKPOS Position of peak vertex  
PKHT Height of peak to be stripped  
ALFA Width parameter of peak to be stripped  
TLAMP Amplitude of tail of peak  
TLSLP Slope of tail  
WL Full width half maximum of Lorentzian line (channels)  
WL = 0.0 for gamma rays.  
Y Data array

RETURN

Y Data array with designated peak structure stripped out.

Routines called: BWF

\* \* \* \* \*

Reference: "An Algorithm for Fitting Lorentzian-Broadened, K-Se  
X-Ray peaks of the Heavy Elements", Nucl. Inst. & M  
143, p145, (1977)

\* \* \* \* \*

REAL\*4 PKPOS,PKHT,ALFA,TLAMP,TLSLP,WL,Y(2)  
INTEGER\*2 NS,NE

IF (WL .EQ. 0.) GOTO 2

CX = SQRT (-ALFA)  
GAMA = WL \* CX  
C1 = 1. - .5642 \* GAMA

```

0058      C2 = .15915 * GAMA
0059      C3 = .25 * GAMA * GAMA
0060      C4 = .6366 * GAMA
0061      RT = C1 + C3
0062
0063      2      DO 10 I = NS,NE
0064             XI = I
0065             DX = XI - PKPOS           !Distance from peak vertex
0066             DDX = DX * DX
0067             ADX = ALFA * DDX
0068             FXX = EXP (ADX)         !Normalized Gaussian value
0069             FX = FXX
0070             IF (DX .GE. 0.) GOTO 4   !Tailing only on front side
0071             EXPN2 = EXP (TLSLP * DX) !Normalized tailing amplitude.
0072             FXX = FXX +TLAMP *EXPN2 *(1.- EXP(.4*ADX)) !Add tailing contri
0073             DX = -DX                !Make + for X ray calculations
0074             4      IF (WL .EQ. 0.) GOTO 9 !Positive value indicates a X ra
0075             CC = CX * DX
0076
0077             CALL BWF (CC) !To calculate last term in following Voigt funct
0078
0079             FXX = FXX -FX +(C1 -C2 *ADX +C3 *(1.+2.*ADX)) *FX +C4 *CC
0080             FXX = FXX / RT          !Adjust peak height for X rays.
0081             9      Y(I) = Y(I) - PKHT * FXX !Strip calc value from data arra
0082             10     CONTINUE
0083
0084             RETURN
0085             END

```

! Command file to link CHNGSI program subroutines  
!

```
$ SD [-.OBJECT]  
$ LINK/executable=[-.ASSAY]PU238_CHNG-  
    pu238_chng,-  
    chngri,-  
    dattim,-  
    errmsg,-  
    getop,-  
    juldal,-  
    nocomm,-  
    pass,-  
    readb,-  
    readc38,-  
    readrm,-  
    writb,-  
    writ38,-  
    yesno  
$ set noverify  
$sd [-.assay]
```

C\*\*\*\*\*

C PU238\_CHNG

C Change constants - [Pu-238 isotopic instrument]

C  
C

06-MAR-80 0935 SSJ

Rev. 14-Jun-84 lc

Rev. 11-Sep-84 TES

Rev. 27-Sep-84 cms

Rev. 21-DEC-84 WDR Modified for solid isotopic instrument

Rev. 03-JAN-85 MPK Restored DATTIM option argument

Rev. 04-Apr-86 WDR Modified for Micro VAX II

Rev. 14-Mar-90 WDR Modified for Pu-238 instrument

Rev. 24-May-91 WDR Modified for stabilizer control info.

C\*\*\*\*\*

C  
C

PROGRAM CHNG

LOGICAL\*1 RANS,REREAD,REWRIT

LOGICAL \* 1 VAXANS

INTEGER\*2 CHAR

C  
C

C The following common block files are generic for all instruments

C

INCLUDE 'CDATE.CMN'

INCLUDE 'CLUNS.CMN'

INCLUDE 'CFILES.CMN'

INCLUDE 'CONFIL.CMN'

INCLUDE 'CBACKG.CMN'

INCLUDE 'CMCONS.CMN'

INCLUDE 'CONTRL.CMN'

INCLUDE 'CASSAY.CMN'

INCLUDE 'CCAL.CMN'

INCLUDE 'CDIAG.CMN'

INCLUDE 'CPASS.CMN'

C

C The following common block files are instrument specific for

C Pu-238 isotopic assay (gamma ray) instrument

C

INCLUDE 'P238COM.CMN'

C

DATA LOUT/5/, DKLUN/10/, DKLUN1/11/

DATA ERRFIL/'ERRMSG.INS' //

DATA BKFIL/'BACKGR.WHB' //

DATA NAMSFL/'CONSTA.WHB' //

C

C Reset reread and rewrite flags

C Set flag for beginning constants read

C

REREAD = .FALSE.

REWRIT = .FALSE.

BEGIN = .TRUE.

CHANG = .TRUE.

FROMDK = .FALSE.

C Read background file

C Get password

C



```
CALL READB
CALL PASS(OK)
IF(.NOT. OK)GO TO 2000
```

```
C
C Get operator id
C Allow escape here
```

```
L
      ABORT = .FALSE.
      CALL GETOP
      IF(ABORT) GO TO 2000
```

```
C
C Get name of constants file to modify
C Read constants file
```

```
C
1      CALL READNM
      CALL READC
```

```
C
C Type menu
```

```
C
5      TYPE 900
900    FORMAT(/' CODE-OPTION')
6      TYPE 910
910    FORMAT(' ---- -')
7      TYPE 1000
1000   FORMAT( ' HE - Help'//
X       ' MC - Measurement Control'//
X       ' PT - Preset Time'//
X       ' PW - PassWords'//
X       ' DV - Default Values'//
X       ' RD - Read constants from Disk'//
X       ' SZ - Stabilizer settings'//
X       ' LI - LIst all constants'//
X       ' OU - Change listing device '//
X       ' EX - EXit from changes'//
X       ' Note: Enter -9999. or -9999 to input 0'//)
```

```
C
10     CALL DATTIM(1)
      TYPE 1001,NOWDAT,NOWTIM
1001   FORMAT(/1X,A9,1X,A8/'#Enter the CODE of the OPTION you want-> ')
      ACCEPT 1002,CHAR
1002   FORMAT(A2)
```

```
C
C Decide what is to be changed
```

```
C
      IF(CHAR .EQ. 'HE') GO TO 5
      IF(CHAR .EQ. 'MC') GO TO 450
      IF(CHAR .EQ. 'PT') GO TO 380
      IF(CHAR .EQ. 'DV') GO TO 500
      IF(CHAR .EQ. 'RD') GO TO 1
      IF(CHAR .EQ. 'SZ') GO TO 700
      IF(CHAR .EQ. 'LI') GO TO 820
      IF(CHAR .EQ. 'PW') GO TO 600
      IF(CHAR .EQ. 'EX') GO TO 800
      IF(CHAR .EQ. 'OU') GO TO 110
      GO TO 10
```

```
C
C Change output device for listings
```

```
L
110    TYPE 1112
1112   FORMAT(/'#Type 6 for listing to printer, 7 for terminal ->')
```

```

ACCEPT 1113,LOUT
1113  FORMAT(I)
      IF(LOUT .NE. 6 .AND. LOUT .NE. 7) GO TO 110
      GO TO 5
C
  C Preset time
L
380  CALL CHNGRL('Asay Tim',0,TIMEPR)
      GO TO 5
C
  C Change measurement control parameters
C
450  TYPE 1401
1401  FORMAT(/'#Measurement control - Bias or Precision? (B/P) -> ')
      ACCEPT 1807,RANS
      IF(RANS .EQ. 'P') GO TO 480
      IF(RANS .NE. 'B') GO TO 450
C
  C Measurement control - Bias and Background
C
      CALL CHNGRL('Biastime',0,BITIME)
      WRITE (*,*) 'Bias ref is the specific power (mW/g) on 1/1/90'
      CALL CHNGRL('Bias ref',0,BIREF)
      CALL CHNGRL('Hist std',0,HISTSD)
      GO TO 5
C
  C Measurement control - Precision
C
480  CALL CHNGRL('Prcstime',0,PRTIME)
481  CALL CHNGIN('No. runs',0,NRUNS)
      IF(NRUNS .EQ. 5 .OR. NRUNS .EQ. 15)GO TO 5
      TYPE 1450
1450  FORMAT(' *** PRECISION RUNS MUST BE 5 OR 15 ***')
      GO TO 481
C
  C Change default values for Pu-240 abundance
C
500  TYPE 1500
1500  FORMAT(' Enter Pu-240 abundance (Wt %) (0=1.97,-1=Calc.)')
      CALL CHNGRL('Pu-240 ',0,PU240)
      TYPE 1600
1600  FORMAT(' Enter Pu-242 abundance (Wt %) (0=0.1)')
      CALL CHNGRL('Pu-242 ',0,PU242)
      NUFLG = 0 !No Uranium default value
C
      TYPE 1505
1505  FORMAT(' Samples contain Uranium? (Y/N) -> ',#)
C
1507  ACCEPT 1507,RANS
      FORMAT(A1)
      IF(RANS.EQ.'Y') NUFLG = 7
      NPFLG = 0 !May be used later
      GO TO 5
C
  C Change passwords
C
600  DO 650 I=1,3
650  CALL CHNGAL('Password',1,PASSWD(I))
      GO TO 5
L
  C Get stabilizer settings for control of gain and zero
C

```

```

700    CALL CHNGIN('Zero pk ',0,IZPEAK)
      CALL CHNGIN('Zero wnd',0,IZWNDW)
      CALL CHNGRL('Zero rng',0,ZARNG)
      CALL CHNGIN('Gain pk ',0,IGPEAK)
      CALL CHNGIN('Gain wnd',0,IGWNDW)
      CALL CHNGRL('Gain rng',0,GARNG)
      GO TO 5

C
C   If changes are permanent, set rewrite flag
C
800    TYPE 1806
1806   FORMAT(/'$Do you want to write these changes to disk now?',
      X      ' (Y/N) -> ')
      ACCEPT 1807,RANS
1807   FORMAT(A1)
      IF(RANS .EQ. 'Y') REWRIT = .TRUE.
      IF(RANS .NE. 'Y' .AND. RANS .NE. 'N') GO TO 800
      GO TO 2000

C
C   List all values
C
820    WRITE(LOUT,1801)NOWDAT,NOWTIM,CONFIL,COMMNT,FILDAT,FILTIM,
      X      FILOP
1801   FORMAT('1',A9,2X,A8//1X,A,2X,A//
      X      10X,'Written ',A,2X,A,' by ',A//)

C
C   Measurement control
C
      WRITE(LOUT,1814)BITIME,PRTIME,BIREF,HISTSD,NRUNS
1814   FORMAT(///' Measurement Control'//
      X      23X,'Bias',14X,'Precision'//
      X      ' Interval = ',2(8X,G13.6)/
      X      ' Reference = ',7X,G13.6/
      X      ' Standard Dev. = ',3X,G13.6/
      X      ' Cycles = ',36X,I4)

C
C   Default values
C
      WRITE(LOUT,1817)PU240,PU242
1817   FORMAT(//,' Default analysis values'//
      1   ' Pu-240 abundance (wt. %) = ',F7.2,/,
      2   ' Pu-242 abundance (wt. %) = ',F7.2)

C
C   Stabilizer settings
C
      WRITE(LOUT,1818)
1818   FORMAT(/' Stabilizer settings for zero and gain',/)
      WRITE(LOUT,1820)
1820   FORMAT(1X,T20,'Zero settings   Gain settings')
      WRITE(LOUT,1821)IZPEAK,IGPEAK
1821   FORMAT(1X,'Peak positions',4X,I6,11X,I6)
      WRITE(LOUT,1822)IZWNDW,IGWNDW
1822   FORMAT(1X,'Peak windows ',4X,I6,11X,I6)
      WRITE(LOUT,1823)ZARNG,GARNG
1823   FORMAT(1X,'Analog ranges ',4X,G13.6,5X,G13.6)
      WRITE(LOUT,1888)
388    FORMAT('1')
      IF(LOUT .EQ. 6)CALL CLOSE(LOUT)
      GO TO 5

```

C Rewrite constants

C  
2000 IF(REWRIT) CALL WRITB  
IF(REWRIT) CALL WRITC

C  
STOP  
END

```
C*****  
C  
C NOCOMM  
C  
C      This module is used instead of COMM when no communication  
C      with the VAX is desired.  
L  
C      29-Mar-85 MPK  
C  
C*****  
C
```

```
SUBROUTINE COMM  
INCLUDE 'CVAX.CMN'  
VAXUP=.FALSE.  
RETURN  
END
```

```

C*****
C
C  GETOP
C      Get operator ID
C
C      Output argument
C
C      MTL (LOGICAL*1)          .TRUE. if material type and
C                               isotopics are needed
C
C      31-May-84 lc
C      Rev 14-Feb-85 MPK        Put in "real" VAX communications
C                               Supply end and err exits for
C                               terminal input
C      Rev 08-Oct-85 MPK        Adapted to neutron coinc counters
C                               GETOP and GETSMP combined
C      Rev 21-Oct-85 MPK        Pad "canned" sample ID properly
C      Rev 23-OCT-85 WDR        Revised for solution assay system
C      Rev 04-Apr-86 WDR        Modified for Micro VAX II
C

```

```

C*****

```

```

C
C      SUBROUTINE GETOP
C
C      The following common block files are generic for all instruments
C
C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CONTRL.CMN'
C
C      LOGICAL*1 YESNO,MTL
C
C      Get operator ID; set abort flag if it is null
C
C      10  TYPE 1000
C      1000  FORMAT (' //Enter operator ID (up to 8 chars)',
C                1  ' or RETURN to escape -> ')
C          READ (5,1002,END=10,ERR=10) OPERID
C      1002  FORMAT(A80)
C          ABORT = OPERID.EQ.' '
C          IF (ABORT) RETURN
C          RETURN
C          END

```

```

C*****
C CHNGRL
C
C     CHNGRL - Change a REAL value
C     CHNGIN - Change an INTEGER value
C     CHNGAL - Change an alphanumeric value
L
C     REV. 28-Jun-82  1500 SSJ
C*****
C
C     SUBROUTINE CHNGRL(NAME,NUMBER,VARBLE)
C
C     Change an individual real value
C
C     REAL*8 NAME
C
C     IF(NUMBER .LE. 0) TYPE 1000,NAME,VARBLE
1000  FORMAT('$',AB, '      = ',G13.6,' -> ')
C     IF(NUMBER .GT. 0) TYPE 1002,NAME,NUMBER,VARBLE
1002  FORMAT('$',AB,'(',I2,') = ',G13.6,' -> ')
C
C     ACCEPT 1001,X
1001  FORMAT(F)
C     IF(X .NE. 0.) VARBLE = X
C     IF(X .EQ. -9999.) VARBLE = 0.
C     RETURN
C     END
C     SUBROUTINE CHNGIN(NAME,NUMBER,IVARBL)
C
C     Change an individual integer value
C
C     REAL*8 NAME
C
C     IF(NUMBER .LE. 0) TYPE 1000,NAME,IVARBL
1000  FORMAT('$',AB, '      = ',I6,7X,' -> ')
C     IF(NUMBER .GT. 0) TYPE 1002,NAME,NUMBER,IVARBL
1002  FORMAT('$',AB,'(',I2,') = ',I6,7X,' -> ')
C
C     ACCEPT 1001,IX
1001  FORMAT(I)
C     IF(IX .NE. 0) IVARBL = IX
C     IF(IX .EQ. -9999) IVARBL = 0
C     RETURN
C     END
C     SUBROUTINE CHNGAL(NAME,NUMBER,ALPHVR)
C
C     Change an individual alphanumeric value
C
C     REAL*8 NAME,ALPHVR,X,BLANK
C
C     DATA BLANK/8H      /
C
C     IF(NUMBER .LE. 0) TYPE 1000,NAME,ALPHVR
1000  FORMAT('$',AB, '      = ',AB,5X,' -> ')
C     IF(NUMBER .GT. 0) TYPE 1002,NAME,NUMBER,ALPHVR
1002  FORMAT('$',AB,'(',I2,') = ',AB,5X,' -> ')
L
C     ACCEPT 1001,X
1001  FORMAT(AB)

```

```
IF(X .NE. BH          ) ALPHVR = X  
IF(X .EQ. BH#####) ALPHVR = BLANK  
RETURN  
END
```



C\*\*\*\*\*

C ERRMSG

C Read appropriate error message from the  
C error message file and output to user terminal  
C Calling arguments:

L INUM = Message calling argument  
C ARG = ASCII argument passed to be printed (15A1)

C File format:  
C INUM (Message call number) = I4  
C IARG (ASCII argument passed) = A1  
C MESSAG (ASCII error message) = 72A1

C 21-Apr-84 LC  
C REV. 12-OCT-84 MJC: CLOSED DKLUN1 IF ARG = ' '  
C Rev. 19-Apr-85 MPK Remove redundant initialization  
C of DKLUN and DKLUN1  
C Rev. 29-Apr-85 MPK Put MESSAG and ICHAR in CMSG common  
C Rev. 02-Apr-86 WDR Modified for Micro VAX II  
C Rev. 13-Jul-87 RDP Modified for uVax II  
C Convert BYTE strings to CHARACTER strings

C\*\*\*\*\*

C SUBROUTINE ERRMSG(INUM,ARG)

C CHARACTER ARG\*15,IARG\*1

C The following common block files are generic for all instruments

C INCLUDE 'CFILES.CMN'  
C INCLUDE 'CLUNS.CMN'  
C INCLUDE 'CMSG.CMN'

C Open error message file  
C Locate proper error message

C OPEN (UNIT=DKLUN1,NAME=ERRFIL,TYPE='OLD',ERR=25)  
10 READ(DKLUN1,1000,ERR=20,END=20)NUM,IARG,ICHR,  
X (MESSAG(I),I=1,ICHR)  
1000 FORMAT(I4,A1,Q,72A1)  
IF (NUM .EQ. -999)GO TO 20  
11 IF (NUM .EQ. INUM) GO TO 15  
GO TO 10

C Print error message

15 TYPE 1015, NUM, (MESSAG(I),I=1,ICHR)  
1015 FORMAT(/5X,'\*\*\*\* ERROR # ',I3/5X,72A1,' \*\*\*\*')  
IF(IARG .NE. ' ')GO TO 17  
CLOSE(UNIT=DKLUN1)  
GO TO 30  
17 TYPE 1016, ARG  
1016 FORMAT(7X,A15)  
CLOSE(UNIT=DKLUN1)  
GO TO 30

C

C Error number not found in file

C

20 TYPE 1020, INUM, ERRFIL

1020 FORMAT(/5X, ' Error # ', I4, ' not located in ', A15)

CLOSE(UNIT=DKLUN1)

GO TO 30

L

C Error opening file

C

25 TYPE 1025, ERRFIL

1025 FORMAT(/5X, '\*\*\*\* UNABLE TO OPEN ', A15, ' \*\*\*\*')

C

C

30 RETURN

END

```

C*****
C
C YESNO
C   Prompt for and accept a yes-or-no answer
C
C   Calling argument:
C   PROMPT  The string to be used for prompting.  It is output,
C           with the characters " ? (Y/N) -> " appended to it.
C
C   Function value (LOGICAL*1):
C   .TRUE.  if the answer begins with the character Y
C   .FALSE. if the answer begins with the character N
C
C   If the answer does not begin "Y" or "N", the prompt
C   is repeated until it does.
C
C   Example:
C   If YESNO is called as follows
C       LOGICAL*1 CONT,YESNO
C       CONT=YESNO('Do you want to continue')
C   the following prompt is produced
C       Do you want to continue ? (Y/N) ->
C
C   04-Feb-85  MPK
C   06-Aug-87  RDP
C
C           Convert BYTE string to Character String

```

```

C*****
C
C   LOGICAL*1 FUNCTION YESNO(PROMPT)
C*****
C
C   CHARACTER PROMPT*(*),TAG*10,QUESTION*70,ANS*1
C   INTEGER    LOC
C
C   DATA TAG/'? (Y/N) ->'/
C   LOC = LEN(PROMPT)+1
C
C   QUESTION = PROMPT
C   QUESTION(LOC:LOC+10) = TAG
C
C   TYPE 1000, QUESTION
1000  FORMAT(1X,A<LOC+11>,$)
      READ (5,1002,END=10,ERR=10) ANS
1002  FORMAT (A)
      YESNO=.TRUE.
      IF (ANS.EQ.'Y') RETURN
      YESNO=.FALSE.
      IF (ANS.EQ.'N') RETURN
      GO TO 10
      END

```

C\*\*\*\*\*

C WRITB

C

C Write background file to disk

C

Rev. 15-Jun-84 1c

L Rev. 11-Sep-84 TES

C Rev. 07-Apr-86 WDR Modified for Micro VAX II

C

C\*\*\*\*\*

C

SUBROUTINE WRITB

C

LOGICAL\*1 FILL

C

C The following common block files are generic for all instruments

C

INCLUDE 'CBACKG.CMN'

INCLUDE 'CLUNS.CMN'

INCLUDE 'CDATE.CMN'

INCLUDE 'CASSAY.CMN'

INCLUDE 'CPASS.CMN'

INCLUDE 'CONTRL.CMN'

INCLUDE 'CFILES.CMN'

INCLUDE 'CMC.CMN'

C

C

C NOTIFY OPERATOR

C

TYPE 1003,BKFIL

'003 FORMAT(/' Writing ',A)

C OPEN A DIRECT ACCESS FILE WITH NAME FILNAM

C

OPEN(UNIT=DKLUN,NAME=BKFIL,TYPE='NEW',ACCESS='DIRECT',  
X RECORDSIZE=128,MAXREC=1,ERR=999)

C

C UPDATE FILE NAME, DATE, TIME AND OPERATOR

C

IF(.NOT. BACKG) GO TO 150

BKSFIL = FILNAM

BKDAT = NOWDAT

BKTIM = NOWTIM

BKOP = OPERID

C

C WRITE DATA

C

150 WRITE(DKLUN'1,ERR=999)BKDAT,BKTIM,BKOP,FILL,  
X BK,SGBK,TIMEBK,PASSWD,BKSFIL

C

C CLOSE FILE AND ACKNOWLEDGE

C

CLOSE(UNIT=DKLUN)

TYPE 1000,BKFIL,BKDAT,BKTIM,BKOP

1000 FORMAT(/1X,A,' Written ',A,2X,A,' by ',A12)

C

GO TO 500

L

C ERROR IN WRITING

C

```
999 TYPE 1001
1001 FORMAT(/20X,'**** WRITB DISK ERROR'/
X          25X,'FILE WAS NOT WRITTEN CORRECTLY ****'/)
C
500 RETURN
END
```

C\*\*\*\*\*

C PASS

C Check the password in order to allow operator to  
C enter supervisory mode

C CALLING ARGUMENTS:

C OK = Logical flag indicating validity of operator password

C SSJ 01-Jul-82

C REV. 08-May-84 lc

C REV. 19-APR-85 MPK Remove logical\*1 OK

C (Redundant with /CPASS/

C REV. 04-APR-86 WDR Modified for Micro VAX II

C REV. 06-AUG-87 RDP Modified to accept password w/o  
C echo to screen. All character are  
C converted to UPPERCASE.

C\*\*\*\*\*

C SUBROUTINE PASS(OK)

C CHARACTER WORD\*10  
C INTEGER JSW,BIT12,NOT12,CR,LF  
C REAL\*8 PWORD

C INTEGER\*4 STATUS,VKID,LENGTH  
C INTEGER\*4 SMG\$CREATE\_VIRTUAL\_KEYBOARD  
C INTEGER\*4 SMG\$READ\_STRING  
C INTEGER\*4 SMG\$DELETE\_VIRTUAL\_KEYBOARD

C INCLUDE '(\$TRMDEF)'

C The following common block files are generic for all instruments

C INCLUDE 'CPASS.CMN'

C EQUIVALENCE (PWORD,WORD)

C DATA CR,LF /"15,"12/

C OK = .FALSE.

C Prompt operator for password

C TYPE 1000

1000 FORMAT(/' You have selected a supervisory option',  
X ' which requires a password.'/)

C Store up to 8 characters

C STATUS = SMG\$CREATE\_VIRTUAL\_KEYBOARD(VKID,'SYS\$INPUT')  
C IF (.NOT.STATUS) CALL LIB\$SIGNAL(%VAL(STATUS))

C STATUS = SMG\$READ\_STRING(VKID,  
C WORD,  
C 'Enter password ->',  
C 10,  
C TRM\$M\_TM\_NOECHO+TRM\$M\_TM\_CVTLOW)  
C IF (.NOT.STATUS) CALL LIB\$SIGNAL(%VAL(STATUS))

```

20      TYPE 1004
1004    FORMAT(1X)
C
      STATUS = SMG#DELETE_VIRTUAL_KEYBOARD(VKID)
      IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))

L  Check if operator wanted to escape
C
      IF(WORD(1:1).EQ.' ') GOTO 500
C
C  Compare with all of the passwords
C
      DO 30 I=1,3
          IF(PWORD.EQ.PASSWD(I)) GOTO 40
30      CONTINUE
C
C  Improper password
C  Inform operator
C
      TYPE 1002
1002    FORMAT(//10X,'**** IMPROPER PASSWORD - NO ENTRY',
X      ' TO SUPERVISORY MODE ****'/)
      GO TO 500
C
C  Password is OK - Enter supervisory mode
C
40      OK = .TRUE.
C
500     RETURN
      END

```

```

C*****
C READB
C
C      Read background from disk
C
C      Rev. 25-Jun-84 lc
C      Rev. 04-Apr-86 WDR          Modified for Micro VAX II
C*****
C
C      SUBROUTINE READB
C
C      LOGICAL*1 FILL
C
C      The following common block files are generic for all instruments
C
C      INCLUDE 'CFILES.CMN'
C      INCLUDE 'CBACKG.CMN'
C      INCLUDE 'CLUNS.CMN'
C      INCLUDE 'CPASS.CMN'
C      INCLUDE 'CONTRL.CMN'
C
C      The following common block files are instrument specific for
C      gamma-ray instruments
C
C      INCLUDE 'CONT90.CMN'
C
C      If FROMDK, read background from spectrum file
C
C      IF(FROMDK) GO TO 10
C
C      Notify operator
C
1      TYPE 1003,BKFIL
1003   FORMAT(/' Reading ',A)
C
C      OPEN A DIRECT ACCESS FILE WITH NAME BKFIL
C
C      OPEN(UNIT=DKLUN,FILE=BKFIL,ACCESS='DIRECT',TYPE='OLD',
X      RECORDSIZE=128,ERR=999)
C
C      Read data
C
C      READ(DKLUN/1,ERR=999)BKDAT,BKTIM,BKOP,FILL,
X      BK,SGBK,TIMEBK,PASSWD,BKSFIL
C
C      Close file and acknowledge
C
C      CLOSE(UNIT=DKLUN)
C      IF(BEGIN)GO TO 500
C      TYPE 1000,BKFIL,BKDAT,BKTIM,BKOP
1000   FORMAT(/,1X,A,' was written ',A,2X,A,' by ',
X      A12)
C      GO TO 500
C
C      Read background file from spectrum file (last record)
C
10     TYPE 1003,FILNAM
X     OPEN(UNIT=DKLUN,NAME=FILNAM,ACCESS='DIRECT',TYPE='OLD',
      RECORDSIZE=128,ERR=997)

```



NRECS = MEMSIZ/128

NREC = NRECS + 1

C  
C Check if first entry in record is dummy

READ(DKLUN/NREC,ERR=998)DUM

IF(DUM .NE. -999.0) GO TO 998

READ(DKLUN/NREC,ERR=998)DUM,BK,SGBK,TIMEBK,BKDAT,BKTIM,

X BKOP,BKSFIL

CLOSE(UNIT=DKLUN,ERR=11)

11 TYPE 1000,FILNAM,BKDAT,BKTIM,BKOP

GO TO 500

C  
C Error in reading from spectrum file  
C Do not abort, but use most recent background read

997 CALL ERRMSG(1,FILNAM)

GO TO 500

998 CALL ERRMSG(2,FILNAM)

CLOSE(UNIT=DKLUN)

GO TO 500

C  
C Error in reading

999 ABORT = .TRUE.

CALL ERRMSG(1,BKSFIL)

CALL ERRMSG(102,ARG)

500 RETURN

END

C\*\*\*\*\*

C WRIT38

C  
C Routine specific to the Pu-238 isotopic instrument  
C Write constants to disk  
C Routine needs to be modified for each instrument

L  
C REV. 14-Jun-84 LC  
C REV. 6-Sep-84 CMS  
C Last edit: 11-Sep-84 TES  
C REV. 27-DEC-84 WDR Modified for the solid isotopic  
C instrument  
C REV. 07-APR-86 WDR Modified for Micro VAX II  
C WDR Included Pu-238 common  
C Rev. 24-May-91 WDR Modified for stabilizer settings

C\*\*\*\*\*

C SUBROUTINE WRITC

C LOGICAL\*1 FILL

C INTEGER\*2 EMC(21), EDD(41), FILLR(58)

C These common block files are generic for all instruments

C INCLUDE 'CFILES.CMN'  
C INCLUDE 'CMCONS.CMN'  
C INCLUDE 'CCAL.CMN'  
C INCLUDE 'CONTRL.CMN'  
C INCLUDE 'CONFIL.CMN'  
C INCLUDE 'CLUNS.CMN'  
C INCLUDE 'CDATE.CMN'  
C INCLUDE 'CASSAY.CMN'  
C INCLUDE 'CDIAG.CMN'

C The following common block files are instrument specific for  
C Pu-238 isotopic instrument

C INCLUDE 'P238COM.CMN'  
C EQUIVALENCE (EMC(1),BITIME),(EDD(1),DIAGMX)

C Notify operator

C TYPE 1003,CONFIL  
1003 FORMAT('/ Writing ',A)

C Open a direct access file with name filnam

C OPEN(UNIT=DKLUN,NAME=CONFIL,TYPE='NEW',ACCESS='DIRECT',  
X RECORDSIZE=128,MAXREC=2,ERR=999)

C Update file date, time and operator

C FILDAT = NOWDAT  
C FILTIM = NOWTIM  
C FILOP = OPERID

```

C
C Write data
C
      WRITE(DKLUN'1,ERR=999)FILDAT,FILTIM,FILOP,FILL,COMMNT,
X      FILLR,EMC,EDD,PU240,PU242,NPFLG,NUFLG,TIMEPR,
X      IZPEAK,IZWNDW,ZARNG,IGPEAK,IGWNDW,GARNG
L
C Close file and acknowledge
C
      CLOSE(UNIT=DKLUN)
      TYPE 1000,CONFIL,FILDAT,FILTIM,FILOP
1000  FORMAT(/1X,A, ' Written ',A,2X,A, ' by ',A)
C
      GO TO 500
C
C Error in writing
C
999  TYPE 1001
1001  FORMAT(//20X,'**** WRITC DISK ERROR'/
X      25X,'FILE WAS NOT WRITTEN CORRECTLY ****'/)
C
500  RETURN
      END

```

```

C*****
C
C READC38      (READC subroutine for PU238 isotopic instrument)
C
C      Read constants for Pu-238 isotopic plutonium gamma-ray
C      analysis instrument from disk
C
C      REV. 14-JUN-84  1c
C      REV. 06-SEP-84  CMS
C      REV. 17-DEC-84  WDR      Modified for solid isotopic instrument
C      REV. 04-APR-86  WDR      Modified for Micro VAX II
C      REV. 13-JUL-87  RDP      Modified for uVAX II
C      Convert BYTE strings to CHARACTER strings
C      Rev. 24-May-91  WDR      Modified for stabilizer settings

```

```

C*****

```

```

C      SUBROUTINE READC
C
C      LOGICAL*1 FILL,CDUM(30)
C      CHARACTER SNAM*12,SYN*3
C
C      INTEGER*2 EMC(21),EDD(41),FILLR(58)
C
C      These common block files are generic for all instruments
C
C      INCLUDE 'CFILES.CMN'
C      INCLUDE 'CONFIL.CMN'
C      INCLUDE 'CLUNS.CMN'
C      INCLUDE 'CMCONS.CMN'
C      INCLUDE 'CONTRL.CMN'
C      INCLUDE 'CASSAY.CMN'
C      INCLUDE 'CCAL.CMN'
C      INCLUDE 'CDIAG.CMN'

```

```

C      The following common block files are instrument specific for
C      solid isotopic analysis (gamma ray) instruments
C
C      INCLUDE 'CONT90.CMN'
C      INCLUDE 'P238COM.CMN'

```

```

C      Specify names of constants files to be read in
C
C      DIMENSION OSHAPC(15)
C      DATA SNAM /'SHAPC.      '/
C      DATA OSHAPC/.04,.00175,-4.,.005,3.,0.,.003,.3,.0,.0,
C      +           3.,.0,.0,.0,.0/

```

```

C      Equivalence dummy arrays to first word of common block
C      Array size = word count of common.
C
C      EQUIVALENCE (EMC(1),BITIME),(EDD(1),DIAGMX)

```

```

C      Notify operator
C
C      TYPE 1003,CONFIL,COMMNT
C      1003  FORMAT(/' Reading ',A,1X,A)
C
C      Open a direct access file with name filnam

```

```

C
C OPEN(UNIT=DKLUN,NAME=CONFIL,ACCESS='DIRECT',TYPE='OLD',
X      RECORDSIZE=128,ERR=999)
C
C Read data
C Fill = Logical*1 fill to even up word boundary
C Cdum = Dummy fill argument to blank out comment
C Fillr = Filler for more parameter arguments
C
C      READ(DKLUN'1,ERR=999)FILDAT,FILTIM,FILOP,FILL,CDUM,
X      FILLR,EMC,EDD,PU240,PU242,NPFLG,NUFLG,TIMEPR,
X      IZPEAK,IZWNDW,ZARNG,IGPEAK,
X      IGWNDW,GARNG
C
C Close file and acknowledge
C
C      CLOSE(UNIT=DKLUN)
C      IF(BEGIN)GO TO 505
C      TYPE 1000,CONFIL,FILDAT,FILTIM,FILOP
1000  FORMAT(/1X,A, ' was written ',A,2X,A,' by ',A)
C
C      Read in constants files specific to solid isotopic plutonium
C      gamma-ray analysis
C
C      First get system id to establish file extension
C
505  DO 100 J = 1,15
      IF (CONFIL(J:J) .EQ. '.')GO TO 120
100  CONTINUE
120  K = J
      DO 150 I = 1,3
      SYN(I:I)= CONFIL(K+I:K+I)
150  CONTINUE
C
C      Read in peak shape parameters
C
C      DO 200 J = 1,3
C200  SNAM(J+6:J+6)=SYN(J:J)
C      OPEN (UNIT=DKLUN,NAME=SNAM,TYPE='OLD',ACCESS='SEQUENTIAL',ERR=995)
C      READ (DKLUN,9087,ERR=995) (OSHAPC(I), I = 1,15)
C9087  FORMAT (E15.7)
C      CALL CLOSE (DKLUN)
C      GO TO 500
C
C Error in reading
C
995  CALL ERRMSG(1,SNAM)
      CALL ERRMSG(102,ARG)
      GO TO 500
999  CALL ERRMSG(1,CONFIL)
      CALL ERRMSG(102,ARG)
C
500  RETURN
      END

```

C\*\*\*\*\*

C JULDAL (Livermore version of JULDAT)

C Using ASCII strings for date, time, compute Julian day  
C From Sue Johnson's Julday

C Calling arguments:  
C IDAT - Input ascii string for date 30-apr-84  
C ITIM - Input ascii string for time 13:54:17  
C XDAT - Output Julian day  
C XTIM - Output time of day in seconds

C 12-APR-83 LC  
C REV. 02-MAY-84 LC for generic package  
C REV. 31-AUG-84 CMS  
C REV. 13-DEC-84 WDR made day 0 1/1/1900  
C modified days calculation  
C Note: Day 0 1/1/1900

C Rev. 11-Oct-84 mjc: IDAT is accepted as 8 characters (1-May-84)  
C or 9 characters (10-May-84)

C Rev. 30-Oct-84 mpk Fix error in 11-Oct correction -  
C tighten code

C Rev. 31-Oct-84 mpk Tighten code some more, and handle  
C time format more flexibly. This is  
C done by copying date and time into  
C an auxiliary array, and editing the  
C array to convert all special characters  
C into commas. Then the DECODE can use  
C the comma-delimited integer field  
C feature.

C Also, correct Julian date algorithm  
C to use integer arithmetic.  
C In the correct algorithm, day 0 comes  
C out to 30-OCT-79, not 28-OCT as  
C stated previously.

C Rev. 01-Nov-84 mpk Leave spaces in date/time string  
C rather than editing to commas.

C Rev. 02-Nov-84 mpk Latest installment in space saga:  
C Ignore leading spaces and replace  
C trailing spaces with commas; this  
C makes it work right for time inputs  
C like 11:23 (seconds omitted).

C Rev. 10-DEC-85 WDR Modified to calculate number of days  
C from 1/1/1900 rather than as described  
C above.

C Rev. 3-Apr-86 WDR Modified for Micro VAX II

C\*\*\*\*\*

C SUBROUTINE JULDAT(IDAT,ITIM,XDAT,XTIM)

C CHARACTER XMONTH\*3  
C CHARACTER\*1 IDAT(9),ITIM(8), DATTIM(23)  
C LOGICAL\*1 TERMFG  
C DIMENSION NDAYS(12)

C DATA NDAYS /31,28,31,30,31,30,31,31,30,31,30,31/

C The following common block files are generic for all instruments

```

C          INCLUDE 'CMONTH.CMN'
C
C Decode ASCII strings
      DO 5 I=1,9
5      DATTIM(I)=IDAT(I)
      DO 6 I=1,8
6      DATTIM(I+9)=ITIM(I)
D      TYPE 1998,DATTIM
C
C Replace special chars ('-' or ':') and trailing spaces with commas
C
      TERMFG = .TRUE.
      DO 7 I=1,17
      IF (DATTIM(I).NE.' ') GO TO 71
      IF (TERMFG) GO TO 7
      GO TO 72
71     IF (DATTIM(I).GE.'0' .AND. DATTIM(I).LE.'9') GO TO 73
      IF (DATTIM(I).GE.'A' .AND. DATTIM(I).LE.'Z') GO TO 73
72     DATTIM(I)=','
      TERMFG = .TRUE.
      GO TO 7
73     TERMFG = .FALSE.
7     CONTINUE
D      TYPE 1998,DATTIM
D1998  FORMAT (' EDITED DATE/TIME: ',17A1)
      IDAY = 0
      XMONTH = ' '
      IYEAR = 0
      IHOUR = 0
      IMIN = 0
      ISEC = 0
      DECODE (9,1000,DATTIM,ERR=75) IDAY,XMONTH,IYEAR
1000  FORMAT (I3,A3,1X,I4)
75     DECODE (8,1001,DATTIM(10),ERR = 76) IHOUR,IMIN,ISEC
1001  FORMAT (2I5,I2)
76     CONTINUE
C
C Compare month to common
C
      IMONTH = 0
      DO 9 I=1,12
9      IF(XMONTH .EQ. AMONTH(I))IMONTH = I
D      TYPE 1999,IMONTH,IDAY,IYEAR,IHOUR,IMIN,ISEC
D1999  FORMAT(' DATE:',3I3,' TIME:',3I3)
C
C Convert time to seconds
C
      XTIM = ISEC + 60. * (IMIN + 60. * IHOUR)
C
C Convert date to days
C
      IF (IYEAR .GT. 1900) IYEAR = IYEAR - 1900
      NMONTH = IMONTH - 1
10     XDAT = AINT (365.25 * IYEAR)
C      1 + AINT (30.6001 * NMONTH)
      2 + IDAY

```

```
C      3      + XTIM/86400.      !uncomment to include fractional day
C
C      Add the number of days up to this month
C
C      IF (NMONTH.EQ.0)GO TO 16
C      DO 15 J =1,NMONTH
15      XDAT = XDAT + NDAYS(J)
16      IF(NMONTH .LT. 3)GO TO 500
C      YR = IYEAR
C
C      Include an extra day for leap year
C
C      IF(AMOD(YR,4.).EQ. 0.)XDAT = XDAT + 1.
C
C      RETURN
500      END
```



```

C*****
C  READNM
C
C    Allow operator to select or change a constants
C    filename in a list of constants files
C
C    LCS 27-FEB-81
C
C    Rev. 16-Jun-84 lc
C    Rev. 19-Feb-85 MPK      Change ACCEPTS to READS with END
C                            and ERR exits
C                            Use YESNO to ask yes/no questions
C    Rev. 04-Apr-86 WDR      Modified for Micro VAX II
C    Rev. 13-Jul-87 RDP      Modified for uVAX II
C                            Convert BYTE strings to CHARACTER strings
C*****

```

```

C
C  SUBROUTINE READNM
C
C    LOGICAL*1 ANS,YESNO
C    CHARACTER CFLNAM(5)*12,CCMMNT(5)*30

```

```

C  The following common block files are generic for all instruments
C

```

```

C    INCLUDE 'CFILES.CMN'
C    INCLUDE 'CONFIL.CMN'
C    INCLUDE 'CONTRL.CMN'
C    INCLUDE 'CLUNS.CMN'

```

```

C
C    NEWFIL = .TRUE.
C    IF(BEGIN) GO TO 4

```

```

C  Check if constants file read is still to be used
C

```

```

1    TYPE 1005, CONFIL,COMMNT
1005  FORMAT(/' The constants file read was: ',A15,1X,A30//
      X   '$Do you want to select another constants file? (Y/N) -> ')
      READ (5,1050,END=1,ERR=1) ANS
      IF(ANS .NE. 'Y' .AND. ANS .NE. 'N') GO TO 1
      IF(ANS .NE. 'Y') NEWFIL = .FALSE.
      IF(.NOT. NEWFIL) GO TO 900

```

```

C  Open file to read names of 5 current constants files
C

```

```

4    OPEN(UNIT=DKLUN,FILE=NAMSFL,ACCESS='SEQUENTIAL',
      X   TYPE='OLD',ERR=5)
      GO TO 80

```

```

C  Problem in open of names file
C

```

```

5    CALL ERRMSG(1,NAMSFL)
      CALL ERRMSG(102,ARG)
      ABORT = .TRUE.
      IF (ABORT) CALL EXIT

```

```

C  Read names from file
C

```

```

80   DO 100 I=1,5

```

```

100      READ(DKLUN,1000) CFLNAM(I),CCMMNT(I)
1000     FORMAT(A12,A30)
C
C   If just beginning, choose file 1
C
      IF(BEGIN) JANS = 1
      IF(BEGIN) GO TO 280
C
C   Type out the current constants file names & select
C   appropriate file
C
150     TYPE 1010
1010    FORMAT(//' The current constants files are: '//)
      DO 200 I=1,5
200     TYPE 1011,I,CFLNAM(I),CCMMNT(I)
1011    FORMAT(3X,I1,3X,A12,4X,A30)
C
210     TYPE 1012
1012    FORMAT(//'#Enter constants file selected (1-5)',
X      ' or 6 for new file -> ')
      READ(5,1020,END=150,ERR=150) JANS
1020    FORMAT(I1)
      IF(JANS.LT.1.OR.JANS.GT.6)GO TO 150
      IF(JANS .EQ. 6) GO TO 400
C
C   Confirm filename
C
250     TYPE 1030,CFLNAM(JANS)
1030    FORMAT(//'#Constants file = ',A12,' OK? (Y/N) -> ')
      READ(5,1050,END=250,ERR=250) ANS
      IF(ANS .EQ.'N') GO TO 210
      IF(ANS .NE.'Y') GO TO 250
C
C   Read filename & comment into common block
C
280     CONFIL=CFLNAM(JANS)
      CONFIL(11:12)=' '
350     COMMNT=CCMMNT(JANS)
      CLOSE(UNIT=DKLUN,ERR=351)
351     GO TO 700
C
C   Input new constants file name & comment
C   Attempt to open new file
C
400     TYPE 1040
1040    FORMAT(//'#Enter new constants filename -> ')
      READ (5,1050,END=400,ERR=400) CONFIL
1050    FORMAT(A30)
C
401     TYPE 1041
1041    FORMAT(//'#Enter comment (30 char) -> ')
      READ (5,1050,END=401,ERR=401) COMMNT
C
C   Check if new file is permanent & which file it replaces
C
      IF ( .NOT.YESNO('Is this file a permanent replacement')
1      ) GO TO 600
C
410     TYPE 1044
1044    FDRMAT(//'#Which file does the new file repliace? (1-5) -> ')

```

```
READ (5,1020,END=410,ERR=410) JANS
IF(JANS.LT.1.OR.JANS.GT.5) GO TO 410
```

```
C
420     CFLNAM(JANS)=CONFIL
421     CMMNT(JANS)=COMMNT
```

```
L Write new file name list to disk
```

```
C
C     REWIND DKLUN
C     CLOSE (UNIT=DKLUN)
C     OPEN(UNIT=DKLUN,NAME=NAMSFL,ACCESS='SEQUENTIAL',
X      TYPE='OLD',ERR=5)
C     DO 430 I=1,5
430     WRITE(DKLUN,1100) CFLNAM(I),CMMNT(I)
1100     FORMAT(1X,A12,A30)
```

```
C Close names file
```

```
C
600     CLOSE(UNIT=DKLUN,ERR=700)
```

```
C Try opening the new constants file
```

```
C
700     OPEN(UNIT=DKLUN1,NAME=CONFIL,TYPE='OLD',ERR=999)
C     CLOSE(UNIT=DKLUN1,ERR=701)
701     GO TO 900
```

```
C
999     CALL ERRMSG(1,CONFIL)
C     GO TO 150
```

```
C
900     BEGIN = .FALSE.
C     RETURN
```

```
END
```

C\*\*\*\*\*

C DATTIM

C  
C  
C  
C

Retrieve, verify or change the system date & time

SSJ 25-JAN-79

REV. 13-SEP-83 LC

Rev 02-Jan-85 MPK

Include option to sync to date and time received from VAX and stored in the CVAX common. Adapt to the NSR/FMF generic program. Add automatic operator input of date and time if date and time are found to be unreasonable.

Rev 14-Jul-87 RDP

Rev 18-Sep-87 RDP

Modified for uVAX II  
Remove Operator from initializing date and time.(IOPT=0)

Options:

C  
C  
C  
C  
C  
C  
C  
C  
C  
C  
C  
C

- 0 - Get initial date and time from operator input and give to monitor.
- 1 - Update date and time from monitor.
- 2 - Update and let operator verify or change.
- 3 - Sync to VAX date and time (note that these must have been set up by a call to VAX communication immediately before calling DATTIM).

C\*\*\*\*\*

C

SUBROUTINE DATTIM(OPTION)

INTEGER OPTION

C  
C  
C

The following common block files are generic for all instruments

INCLUDE 'CDATE.CMN'  
INCLUDE 'CMONTH.CMN'  
INCLUDE 'CVAX.CMN'

C

CHARACTER\*9 INPSTR  
CHARACTER\*23 DATE\_TIME

C

LOGICAL\*1 VFD

C

INTEGER\*4 STATUS, SYSTIM(2), SYS#BINTIM

C

INTEGER\*2 SYS#SETIME

C

EXTERNAL SYS#BINTIM, SYS#SETIME

C

DATA AMONTH /'JAN','FEB','MAR','APR','MAY','JUN',  
1 'JUL','AUG','SEP','OCT','NOV','DEC'/  
DATA ALOWDT /0./ !Julian date for 01-Jan-1900  
DATA AHIDAT /36525./ !Julian date for 31-Dec-99

C

Set internal option variable. If it is zero, force operator input. If it is 3, sync to date and time from VAX.

C

IOPT=OPTION

```
IF (IOPT.LE.0) GOTO 100
IF (IOPT.GE.3) GOTO 300
VFD=.FALSE.
```

```
C
C Update alphanumeric date and time.
```

```
10 CALL DATE (NOWDAT)
CALL TIME (NOWTIM)
```

```
C
C If the date and time have been operator verified, return now.
```

```
IF (VFD) GOTO 500
```

```
C
C For option 1, return if date and time from VAX are reasonable;
C otherwise, make the operator input them by setting the internal
-C option to 0.
```

```
IF (IOPT.NE.1) GOTO 40
CALL JULDAT (NOWDAT,NOWTIM,ADAT,ATIM)
VFD=.TRUE.
IF ((ADAT.GE.ALLOWDT).AND.(ADAT.LE.AHIDAT).AND.
1 (ATIM.GE.0.0).AND.(ATIM.LE.86399.0)) GOTO 500
```

```
IOPT=0
GOTO 100
40 CONTINUE
```

```
C
C Print out date and time for operator verification.
```

```
TYPE 1000,NOWDAT,NOWTIM
1000 FORMAT(/1X,A,2X,A)
VFD=.TRUE.
```

```
C Allow operator to change date.
```

```
100 WRITE(*,1003)
1003 FORMAT(' '//'$Enter Date in the form 01-FEB-85 ->')
IF (IOPT.NE.0) WRITE(*,1006)
1006 FORMAT('$ (or RETURN for no change)')
READ(5,1014,END=100,ERR=100) INPSTR
C! READ(5,1014,END=100,ERR=100) NCHAR,INPSTR
C! 1014 FORMAT(Q,A)
1014 FORMAT(A)
IF ((IOPT.NE.0).AND.(INPSTR.EQ.' ')) GOTO 200
IF (INPSTR(3:3).NE.'-') GOTO 100
D! IF (NCHAR.NE.0) GOTO 100
IF (IOPT.EQ.0) GOTO 110
GOTO 200
```

```
C
C We appear to have a date from the operator. Call JULDAT to
C parse it, and check its components for reasonableness.
```

```
110 VFD=.FALSE.
CALL JULDAT(INPSTR,NOWTIM,ADAT,ATIM)
IF (IYEAR.GE.84 .AND. IYEAR.LE.99) GOTO 120
TYPE 1024
1024 FORMAT (//20X,'**** INCORRECT YEAR NUMBER ****')
GOTO 100
120 IF (IMONTH.GE.1 .AND. IMONTH.LE.12) GOTO 130
TYPE 1025
1025 FORMAT (//20X,'**** INCORRECT MONTH NAME ****')
```

```

GOTO 100
130 IF (IDAY.GE.1 .AND. IDAY.LE.31) GOTO 170
TYPE 1020
1020 FORMAT (//20X,'**** INCORRECT DAY NUMBER ****')
GOTO 100
170 DATE_TIME(1:8)=INPSTR(1:8)
J = INDEX(DATE_TIME,'-')+1
IF (J.GT.3) GOTO 175
TYPE 8000
8000 FORMAT(//20X,' **** INCORRECT FORMAT 04-APR-86 ****')
GOTO 100
175 J = INDEX(DATE_TIME(J+1:),'-')+1+J
DATE_TIME(J:J)='1'
DATE_TIME(J+1:J+1)='9'
DATE_TIME(J+2:J+2)=INPSTR(J:J)
DATE_TIME(J+3:J+3)=INPSTR(J+1:J+1)

```

C  
C Allow operator to change time.  
C

```

200 WRITE(*,1004)
1004 FORMAT(' '//#Enter Time in the form 01:47:59 ->')
IF (IOPT.NE.0) WRITE(*,1006)
READ(5,1014,END=200,ERR=200) INPSTR
IF ((IOPT.NE.0).AND.(INPSTR.EQ.' ')) GOTO 500
IF (INPSTR(3:3).NE.':') GOTO 200
IF (IOPT.EQ.0) GOTO 210
GOTO 300

```

C  
C We appear to have a time from the operator. Call JULDAT to  
C parse it, and check its components for reasonableness.  
C

```

100 VFD=.FALSE.
CALL JULDAT (NOWDAT,INPSTR,ADAT,ATIM)
IF (IHOURL.GE.0 .AND. IHOURL.LE.23) GO TO 220
TYPE 1034
1034 FORMAT (//20X,'**** INCORRECT HOUR NUMBER ****')
GO TO 200
220 IF (IMIN.GE.0 .AND. IMIN.LE.59) GO TO 230
TYPE 1036
1036 FORMAT (//20X,'**** INCORRECT MINUTE NUMBER ****')
GO TO 200
230 IF (ISEC.GE.0 .AND. ISEC.LE.59) GO TO 270
TYPE 1038
1038 FORMAT (//20X,'**** INCORRECT SECOND NUMBER ****')
GO TO 200
270 CONTINUE

```

C  
C Load DATE\_TIME character string with time and pad with blanks

```

DATE_TIME(12:12)=' '
DATE_TIME(13:20)=INPSTR(1:8)
DATE_TIME(21:23)=' '
C
TYPE 9000,DATE_TIME
9000 FORMAT(1X,23A)
STATUS = SYS#BINTIM(DATE_TIME,SYSTEM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
STATUS = SYS#SETIME(SYSTEM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))

```

C  
C Repeat update and verification.  
C

```
IOPT=2
GOTO 10
```

```
C
C Sync to the VAX time.
C
DO CALL JULDAT(VAXDAT,VAXTIM,ADAT,ATIM)
L. Get VAX time from comm
C! CALL COMM('TS')
C Load DATE_TIME character string with time and date (modified date
C to meet Date String requirements.
DATE_TIME(1:9) = VAXDAT
DATE_TIME(10:11)= DATE_TIME(8:9)
DATE_TIME(8:9) = '19'
DATE_TIME(12:12)= ' '
DATE_TIME(13:20)= VAXTIM
DATE_TIME(21:23)= ' '
C TYPE 9000,DATE_TIME
C Update Time
STATUS = SYS#BINTIM(DATE_TIME,SYSTEM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
STATUS = SYS#SETIME(SYSTEM)
IF (.NOT.STATUS) CALL LIB#SIGNAL(%VAL(STATUS))
500 RETURN
C
C
END
```