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P.BERNARD-MOULIN

Thermo Fisher Scientific, patrick.bernard-moulin@thermofisher.com

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Abstract:

A solid dispersion from Hot Melt Extrusion gives likely questions to be answered:

- *Chemical and stage modification on final product*
- *Extrudates homogeneity*
- *Moisture content and impurities existence and quantities*
- *Interactions or modifications during processing*

Vibrational spectroscopy is a non-destructive/safe in-line measurement allowing process understanding through continuous real time and multiple measurements along the barrel

Since NIR is sensitive to influence of different process parameters it delivers valuable information for scale-up. It quantify and monitor chemically and physically all ingredients which can be detected with NIR. Raman is a complementary technique by imaging component distribution on cross-sectioned samples

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
Raman MCR Image

purple – HPMCAS, green – ibuprofen, yellow – epoxy,
red – cyanoacrylate, blue – inorganic impurity

Composant	% (Surface)	Reporté %
ibuprofen	23	25

Références

- [1] Thomas Chirkot, Sheelagh Halsey and Alan Swanborough, NIR NEWS Vol. 25 No. 2 March/April 2014
- [2] A.L. Kelly, T. Gough, R.S. Dhumal, S.A. Halsey, A. Paradkar Monitoring ibuprofen–nicotinamide cocrystal formation during solvent free continuous cocrystallization (SFCC) using near infrared spectroscopy as a PAT tool, International Journal of Pharmaceutics 426 (2012) 15– 20
- [3] Dr. Robert Heintz, Spectroscopy, November/December 2014