Project 505.6050

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0467

1987-02-12

Near Infrared Reflectance Spectroscopy. Pr.nr. 505.6050

Part of the activities of the State Institute for Quality Control of Agricultural Products consists of developing and promoting (automated) fast methods to measure quality defining parameters of foods and feeding-stuffs. Since 1980 we use within this framework near infrared reflectance instruments for compositional analysis of foods and feeding-stuffs and for research purposes. At this moment we investigate the possibility of nondestructive measurements by NIRS directly without any sample preparation. A part of this study is to make an inventory on all kinds of commercial available NIRS/NITS instruments. Therefore we ask you as well as all other producers of NIRS instruments for information about the following items:

- General information on your NIRS/NITS instruments, including prices of the different systems.
- Wavelength extension to the visible region.
- Available software.
- On-line measurements.
- Automatic sample introduction and measurements.
- Adaptable sample compartment to accomodate various forms of samples.
- Transmittance techniques (fruit-nondestructive).
- Your view on future applications, especially on NIT.
- Address of your representative in Holland or in Europe.
- References.

With kind regards,

R. Frankhuizen
NIR-Scientist

Fra/IH

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Multilab b.v.
Postbus 1180
3430 BD NIEUWEGEIN

LT Industries, Inc.
6110 Excutive Boulevard
Rockvill, Maryland 20852
United States of America

Trebor Industries Inc.
P.O. Box 2159
Gaithersburg, Maryland 2087
United States of America

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Labor MJM
P.O. Box 73
Budapest H-1450
Hongarije

Moisture Systems Corp.
120 South Street
Hopkinton, MA 1748
United States of America

USDA, Russell Research Center
Dr. G.S. Birth
P.O. Box 5677
Athens, GA 30613
United States of America

Foss Electric Ltd.
The Chantry
Bishothorpe, York, YO2 1QF
Great-Britain

Bruker Instruments Inc.
Manning Park, Billerica
Mass 01821
United States of America

Varian Benelux B.V.
Jan Rebelstraat 21
1069 CA AMSTERDAM

Multispec Inc.
23560 Lyons Avenue
Suite 209 Newhall
California 91321
United States of America

Nicolet b.v.
Meidoornkade 19
3992 AG HOUTEN

Polaron Equipment Ltd.
Greenhill Crescent Watford Business Park
Watford Hertfordshire WD1 8QS
Great Britain

LLOYD Instruments PLC
Brook Avenue
Warsash, Southampton
Hampshire SO3 6HP
Great Britain


Applied photophysics Ltd.
18-21 Gorsham St.
Londen N1 6DR
Great Britain

Mattson Instruments, Inc.
1001 Eoarier Court
Modison W1 53717
United States of America

Perkin-Elmer Nederland b.v.
Postbus 490
2800 AL GOUDA

QUANTUM 1200

BIJLA

NIR/VIS

Analyzer

Rapid Scanner For
Constituent And
Physical Properties

FOR USE IN LABS AND ON-LINE APPLICATIONS

LT Industries, Inc.

SPECIFICATIONS

Reliability: Tested for around the clock uninterrupted operation for on-line and laboratory applications.

Spectral Response: 1200-2400 nm

Spectral Expansion: 400-2400 nm

Speed of Scan: 5 complete scans per second

Modes of Measurement: Reflectance and Transmittance

Remote Measurement: Fiber Optics Extension

Spectral Resolution: 1 nm data interval

Bandwidth: adjustable

Wavelength Accuracy: ± 0.5 nm in the NIR.

Wavelength Repeatability: 0.03 nm or better

Stray Light: 0.1% or better

Detection Configuration: Optimized for minimum specular reflectance

Measurement: 16 bit data conversion corresponds to 1 part in 65,000

Software: Complete package includes:

- Data collection programs
- Data transformation programs
- Regression analysis
- Statistical interpretation
- Measured quantity derivation
- Help routines
- Diagnostics
- Extensive graphics
- Utility programs


Computer Requirements: IBM PC, XT, AT and compatibles

Auxiliary Port: Computer controlled auxiliary port for closing a process loop and controlling accessories.


Power: 120V, 3.5 Amp, 60 Hz.
or 220V, 1.8 Amp, 50 Hz.

Power Transients: Filtered


Dimensions: 12" x 15 $\frac{1}{4}$ " x 13 $\frac{1}{4}$ " (h,w,d)


● Reflectance

● Transmittance

● Fiber optics remote measurement

OPTIONS AND ACCESSORIES:

- Fiber optics scanning capabilities to allow remote on-line measurements.
- Adaptable sample compartment module to accommodate various forms of samples.

SERVICE AND WARRANTY



- One year warranty.
- Self diagnostic features to allow immediate response.
- Guaranteed rapid service and dedicated support.

SPECTRAMETRIX SOFTWARE

DATA ACQUISITION PROGRAMS—assure reliable data collection, flexibility in file and sample control, averages and ratios multiple scans.

GRAPHICS PROGRAMS—High resolution color graphics with multiple graphs display. Complete flexibility in examination of ranges, determination of points, text insertion and change of titles. Built in smoothing routines. Output compatible with graphics printers.

MATHEMATICAL PREPROCESSING PROGRAMS—Derivatives, normalization, ratioing, Kubelka-Munk, logarithms, smoothing, and more. Cascading of multiple math routines and flexible input/output file structures.

MODELING PROGRAMS—Complete package with the most advanced data analysis techniques including stepwise regression analysis, local minima multiple model generation, statistical factors and coefficient for optimization of selected model. Complete calibration and prediction results.

DIAGNOSTIC PROGRAMS—Instantaneous verification of system calibration, operation, and performance.

CUSTOM PROGRAMS—Customized programs for controlling accessories through the auxiliary port, closing a loop in a process control operation and other needs.

UTILITY PROGRAMS—Provides status and content of data files and processing history, at all levels.

```

L.T. Industries Composition Analysis
GRAPHING SUBSYSTEM (V1.28)

GRAPHING COMMANDS:      Range      Title      Y-Label      Grid
(ENTER) Draw            1200.0-2400.0nm  Samples     Reflectance  YES
(F1) Help
(F2) Samples
(F3) Range
(F4) Title
(F5) Y-Label
(F6) Grid
(F7) Exam/ins
(F8) Print
(F9) Save
(F10) QUIT

SAMPLES
1 SAMPLE 1 (211)
SAMPLE 5 (211)
3 SAMPLE 6 (211)
SAMPLE 21 (211)

GRAPHING SUBSYSTEM: Command Summary

For general operating instructions, press
(F1) again. Use the (F6) and (F8) keys on the
numeric keypad to scan the information.

(ENTER) DRAW THE GRAPH *****

Press the (Enter) key (sometimes called
(Return)) to draw the graph on the screen. Up to
four samples can be drawn. The (Enter) key is
located to the left of the numeric keypad (where the
arrows are) above the key labeled (F10).
    
```

FEATURES AND BENEFITS

- **A Turn Key System:** for instant and accurate measurement for cost savings.
- **Powerful Software:** exclusive SpectraMetrix software package.
- **Fiber Optics for Remote Measurements.**
- **Direct Modeling:** For maximum accuracy in transferring laboratory measurement to on-line.
- **Compact and Rugged Design:** Ideal for on-line process monitoring.
- **Flexible Sample Compartment Configuration:** to enable use of various sample sizes and shapes and accommodate liquid, solid, slurries, pastes and powders.
- **Wide Spectral Range:** Spectral range from VIS through NIR in reflectance and transmittance modes.
- **Computer Controlled Auxiliary Port:** for controlling accessories, or closing a loop in a process.
- **Extensive Diagnostic and Auto Calibration.**
- **Versatile Design:** Instrument can be used in any desired spatial orientation.
- **Continuous Operation and Analysis:** A unique design of the scanning mechanism provides continuous, reliable, high accuracy, and repeatability of scanning.
- **Rapid Analysis:** Scan rates 5 scans per second.
- **Superior Specifications:** high energy throughput, adjustable bandwidth, high sensitivity, and wide dynamic range.
- **Total Solution:** for your non-destructive, continuous and rapid monitoring, control and measurement needs

FEED, FORAGE AND FOOD ANALYZER

Measures Nutritional Properties of Complex Foods

TREBOR-70, THE MOST VERSATILE NEAR-IR INSTRUMENT

• Direct Measurement of Up to 12 Constituents Per Product:

- Protein
- Moisture
- Oil/Fat
- Crude Fiber
- ADF
- Heat Damaged Protein
- Ca, Mg, Ph
- Amino Acids including:
 - * Lysine, Methionine, Threonine, Tryptophan, Serine, Isoleucine, Glycine, Alanine, Valine, Proline, etc.
- Many Others

- Measures in Both Reflectance and Transmission
- Can Do Samples with High Moisture or High Oil/Fat (e.g. Liquids, Pastes)
- Can Do Dry Powders/Meals
- Can Do Whole, Unground Grains, Oilseeds, Snack Foods, Without Any Sample Preparation or Grinding
- Simple to Use, All Answers in 30 Seconds
- Delivered Factory Calibrated For Most Products (e.g., mixed feeds)
- Built-in Computer Also Allows User to Perform Calibrations

TREBOR-70 FEED, FORAGE AND FOOD ANALYZER

BACKGROUND


Near-IR instruments have become widely accepted in many industries. For example, wheat is routinely bought and sold based on its protein content, which is routinely measured by simple Near-IR instruments.

On "complex" products — e.g. mixed animal feeds, forages, prepared food products — where the recipe may change as ingredient prices and availability change, accurate measurement requires optical data at many wavelengths.

The TREBOR-70 is specifically designed to measure these complex products. It provides reflectance measurement of powdered products at every 1 nm between 1100 and 2500 nm. For high moisture products (e.g. milk, meat, silage) as well as large granular products such as whole grains and snack foods, the instrument performs transmission measurements between 600 and 1100 nm, also at every 1 nm.

This full spectrum capability is directly converted into "percent" by the TREBOR-70's powerful built-in computer.

SIMPLE TO USE

FOR NIR RESEARCHERS

The TREBOR-70 contains a precision high-energy spectrophotometer. The system provides transmission scans from 600 to 1100 nm with the data measured at every 1 nm. Reflectance scans are from 1100 to 2500 nm with data also measured at every 1 nm. Typical scan time is 20 seconds.

The TREBOR-70 can be coupled to any external computer thus allowing versatile NIR research. Research software is available for the various models of the IBM-PC and Kaypro computers.

Trebor Industries, Inc.


P.O. Box 2159, Gaithersburg, MD 20879

Telephone (301) 948-7650 TWX/TLX 710-828-9701

Model 6250 System

NEAR INFRARED SPECTROPHOTOMETER

NIR
THE NEW
SPECTROSCOPY

GARDNER
NEOTEC
Division  PACIFIC
SCIENTIFIC

THE
NEAR INFRARED
SPECTROSCOPY
SPECIALISTS

NEAR INFRARED (NIR) SPECTROSCOPY

In the laboratory, the mid-infrared (MIR) spectral region has traditionally been used for qualitative analysis of organic substances. It is in this portion of the IR spectrum that absorption bands arising from transitions involving vibrational energy states occur. These states correspond to specific bond deformations, either stretching or bending, within a molecule. In complex molecules, vibrational deformations of various functional groups—hydroxyl, methyl, carbonyl—often give rise to characteristic absorption bands whose wavelengths or frequencies usually do not vary significantly. This is due to the general insensitivity of localized functional groups to the much larger mass of the remainder of the molecule. Such characteristic absorption bands can rapidly provide direct qualitative information regarding the presence or absence of specific functional groups within a compound. The IR spectrum of a specific compound can also be considered as a "fingerprint" for that compound, and therefore, be used to either identify an unknown substance, or determine its presence in a mixture.


The widespread use of the MIR spectral region for qualitative analysis was carried over when IR techniques were later applied to quantitative analytical problems. Unfortunately, MIR analysis was somewhat of a disappointment when used quantitatively. This was primarily due to the experimental difficulties encountered in obtaining

reproducible results for most types of samples. MIR spectra of most liquids had to be recorded using extremely small pathlengths, typically much less than 0.1 mm. At such pathlengths, the tolerances required for good quantitative results are very difficult to achieve on a routine basis. It was possible to obtain transmittance spectra of solids by preparing KBr pellets or nujol mulls, or by casting thin films using a suitable solvent (if one could be found). Reflectance spectra were usually recorded using the attenuated total reflectance (ATR) technique.

The relative strength or intensity of the absorption bands in the MIR spectral region necessitated the use of such sample preparation techniques together with their subsequent problems. Absorption bands in the near-infrared region (NIR), however, are simply overtones and combinations of the bands occurring in the MIR, and as such, are much weaker than their MIR counterparts. Overtones can be best thought of as harmonics of the MIR fundamentals; they occur at near integral multiples of the fundamental absorption frequency and are often referred to as first, second, or higher order overtones. Since intensity decreases rapidly with increasing order, only the first through third overtones are generally observed


in the NIR. Combinations are bands occurring at frequencies equal to the sum or difference of two or more fundamental absorption frequencies, or an overtone and a fundamental. Due to the relative weakness of NIR absorption bands it is possible to use reasonable pathlengths (0.5mm–1.0cm) for transmittance measurements of liquids and solids, as well as the use of neat samples (often without any prior preparation) for diffuse reflectance measurements of intact solids. Simplicity of sample preparation techniques also results in good reproducibility of spectral data.

Unlike the MIR spectrum, which consists of absorption bands that are fairly well distinct from one another, the NIR spectrum of most compounds is one of overlapping bands. The overlap is compounded for spectra of complex mixtures, making it sometimes quite difficult to correlate specific bands to the components from which they arise. This is especially true when the concentration of the component to be measured is less than 10 percent. In such cases, moderate concentration variations may cause only very slight changes in the NIR spectrum.

ADVANCED MONOCHROMOTOR DESIGN

The heart of the Model 6250 System is its monochromator consisting of a cam-driven f2 concave holographic grating. This cam drive gives fast scan speeds of 2.5 scan/sec., thereby allowing averaging of multiple scans to improve precisions of both absorbance and wavelength. More importantly fast scanning allows moving the sample while measurement is being made to effectively increase sample size. The f2 grating provides high energy throughout, which greatly increases signal to noise ratios over conventional spectrophotometers.

Therefore, to effectively utilize the Near Infrared spectrum for qualitative and quantitative analysis, and take advantage of the ease of sample presentations afforded by this technique, the analyst will require: (1) a spectrophotometer—which generates spectra with extreme precision both in absorbance and wavelength, and (2) a computer—with advanced software to interpret the resulting complex spectra.

MODEL 6250 SYSTEM

The Model 6250 System incorporates the latest state-of-the-art computer hardware and software with a unique patented technology to meet the requirements of NIR quantitative analysis. This system features:

- Superior instrumental performance
- A flexible and versatile sample presentation system
- Broad software library
- True calibration transferability

SUPERIOR INSTRUMENTAL PERFORMANCE

Single Beam Operating Mode

Modern electronics achieves the spectroscopist's goal of true single beam operation which means:

■ **Low Photometric Noise**

Uniform noise levels of 25×10^{-4} log 1/R units S.D. are achieved in a total spectral scan time of 20 seconds. For those applications requiring lower noise levels, the scan time may be increased.

■ **High Wavelength Precision**

Wavelength precision is better than 0.015nm S.D. on repeat scans.


■ **Simultaneous Reflectance/Transmittance Measurement Modes**

The system operates in the reflectance mode for solid material and the transmittance mode for liquids and slurries. Lead sulfide or silicon detectors are automatically used depending on the measurement mode selected through the system's software.


■ **Fast Scan**

One complete NIR spectrum is obtained in less than 0.4 seconds. Full performance is achieved by averaging over 600 A/D conversions per wavelength in 20 seconds for a complete spectral scan.



Noise as a Function of Wavelength

Wavelength Repeatability Using Polystyrene

FLEXIBILITY / VERSATILITY

SPECTRAL REGION 3


600–1100nm is a further overtone region where absorptions are extremely weak. It is a very useful region for analyzing slurries in the transmittance mode. Pathlengths of several centimeters are possible in this region making the analysis of such difficult products as processed meat, cheese, and butter extremely easy. This region is not used for diffuse reflectance measurements because the extremely weak absorptions do not provide sufficient information in this mode.

SPECTRAL REGION 2

1100–1800nm is an overtone region whose absorptions are much less absorbing than Spectral Region 1. This region is used for transparent films or for liquid analysis using pathlengths of up to 5mm in the transmittance mode. In the diffuse reflectance mode, it is ideal for the analysis of material with coarse particles or for highly absorbing material.

SPECTRAL REGION 1

1800–2500nm is the first combination tone region containing the sharpest and most intense spectral absorptions in the NIR spectrum. It is the best region for the analysis of liquids and transparent films in the transmittance mode where pathlengths of less than 2mm can be conveniently used. In the diffuse reflectance mode, this is an ideal region for analysis of uniform powdered materials or materials with smooth surfaces.

STATIONARY SAMPLE PRESENTATION SYSTEM

The standard sample presentation system is primarily used to analyze powdered material in the reflectance mode using Spectral Regions 1 and 2. With the addition of a transmittance detection, this system can be used to analyze liquids in Spectral Regions 1 and 2, and slurries in Spectral Region 3. The first major application that demonstrated the usefulness of the NIR spectrum was the

determination of the internal quality of fruits and vegetables. These tests were performed using Spectral Region 3 in the transmittance mode. The 6250 System can be adapted to perform these same measurements with the simple addition of the large object chamber which can accommodate samples up to 6 inches in diameter

BROAD SOFTWARE LIBRARY

The system software used to interpret NIR spectra is as critical to success as the spectrophotometer itself. Pacific Scientific's Near Infrared Spectral Analysis Software (NSAS-PC) is designed to offer total calibration flexibility in addition to a comprehensive range of ancillary programs.

The American Society of Testing and Measuring (ASTM) has proposed two calibration modes for use in NIR Spectroscopy:

Closed Population: To calibrate and predict on a closed sample population.

Open Population: To calibrate on a closed sample population in order to predict an open sample population.

NSAS-PC is a comprehensive software package with a selection of programs necessary to implement both closed and open population calibrations.

Closed Population Calibration Software

- 700 Log 1/R or 1/T Wavelength Selection
- Multiterm Log 1/R, 1/T and Derivative Log 1/R, 1/T Regression
- Best Pair/Best Triplet, All Possible Combinations
- Automatic Stepwise Fast Forward Regression
- Discriminant Analysis

Open Population Calibration Software

- 1st through 4th Derivative Log 1/R and 1/T
- Pacific Scientific Stepwise Regression
- Norris Optimization Regression
- Spectral Matching



Ancillary Programs

- Sample Selection - Both Random and Greatest Spectral Variation methods are available to select the optimum calibration sample set.
- Advanced Spectral Transformations—The Log 1/R and Log 1/T spectra are quickly converted to one of several Advanced Spectral Transformations.


Reflectance R, Transmission T
 N-point Smooth Log 1/R, Log 1/T
 1st Derivative
 2nd Derivative
 3rd Derivative
 4th Derivative
 Savitsky Golay
 Kubelka Munk
 Division by Fixed Wavelength
 FTNIR

- Spectral Library—Provides a means of confirming the wavelengths chosen in calibration. User's library can be developed for raw materials and finished products.
- Statistical Quality Control—Charts are used to establish stable calibration equations and to implement SQC for routine analysis.


Master System

Correlation and Reconstructed Spectral Plots

Log 1/R and Derivative Spectral Plots


LAB vs. NIR

**MOVING
SAMPLE
PRESENTATION
SYSTEM**


(Analyze samples from less than a gram to hundreds of a gram)

The largest error source in any analysis system is always sampling error. Intact samples are not completely homogeneous. Therefore, analyzing a small portion of a sample is not completely representative of the total. The only way to improve this situation is to average multiple analyses of the same sample. This is done in a unique way with the **6250 System** using the Sample Transport Module. The **6250 System** scans the complete NIR spectrum in less than 0.4 seconds. This means that samples can be moved past the

viewport while multiple scans are averaged. The resulting average scan will be the average spectrum of the sample that was moved past the viewport. This method virtually eliminates sampling error. Stationary and flowing liquids or slurries can be analyzed without limitation to the amount of flowing sample that is analyzed. Samples with large particles which are not homogeneous can be analyzed by simply analyzing large amounts of sample.

Sample	Sheet	Slurries	Small Sample Slurries	Coarse Powders	Flowing Liquids	Static Liquids
Sample Area(cm)	18.0-73.0	18.0-73.0	4.5-38.0	17-69	6.5	3.6
Sampling Method	1	2	3	4	5	6
Reflectance	X			X		
Transmittance	X	X	X	X	X	X
Pathlength(mm)	To 1	10-15	10-15	15	5-20	1-20
Temp. Control	—	—	—	—	Yes	Yes
Range °C	—	—	—	—	25-65	25-65
Typical Products	Plastics Films Textiles Paper Coatings	Cheese Margarine Meat Salad Dressing Dough	Same as Slurries where less sample is available	Soymeal Wheat Plastic Pellets Tobacco Soaps	Polyols Oils Detergents Petroleum Products Beverages Reagents	Viscous Polyols Heavy Oils Milk Cream

■ **Calibration Transfer**—The Pacific Scientific NIR spectrophotometers enable optimum calibrations to be developed and transferred between instruments. The Calibration Transfer program matches the wavelengths and optical responses between the MASTER and SLAVE instruments and then calculates the transfer equation for the SLAVE instruments.

■ **Spectral Matching**—A fast computation of the vector dot product of a sample spectrum against a reference library of spectra for identification of the unknown or simply against a single reference spectra for qualification or fingerprinting.

3rd Party Software

Pacific Scientific is the exclusive supplier of extensive 3rd party software from:

- Honig's Scientific Ltd.
- Infrasoft International
- N.C. State University

This software includes programs for:

- Interfacing to LOTUS 1-2-3, SAS Institute, etc.
- Partial Least Square Regression
- Forage and Mixed Feed and Feed Ingredient
- Calibration and Calibration Development Software
- Fourier Transform Regression

NIR Information Network
(PSCo Bulletin Board)

The PSCo NIR Information Network (Bulletin Board) can be accessed twenty-four hours a day with the simple addition of an inexpensive MODEM to the PC. Special access codes are provided which allow either confidential communications with Pacific Scientific or general communication with all users. This Network gives the user the ability to obtain general NIR information, application information, programs, updates, etc. Pacific Scientific will answer specific questions on a confidential basis or the user can share experience with the entire Network for general comment. New applications, calibrations, software and other information such as meetings, papers related to NIR spectroscopy are continually updated on the Network. Software is also available to operate local Networks within a user's own company.

TRUE CALIBRATION TRANSFERABILITY


(To Low-cost Rugged Filter Instruments)

Prediction equations developed on the **6250 System** can be transferred to low cost, rugged filter instruments which may be located at distant experimental stations, quality control laboratories, or in industrial production environments.

Any equation defined with the **Model 6250** can be transferred directly to PSCo scanning filter instruments, including on-line instruments. The PSCo computer transfer network ensures that transfers are independent of wavelength used, transformation used, and number of terms used in the prediction equation.

Wavelengths and sensitivity at those wavelengths are matched. No slope or bias adjustment is required after transfer in this mode.

Traditional multiterm Log 1/R equations are transferred directly to fixed filter instruments, but require slope and bias adjustment for proper prediction performance.

Model 6250 System

NEAR INFRARED SPECTROPHOTOMETER

SPECIFICATIONS

Spectrophotometer

Scan Time: 0.4 seconds minimum;
20 seconds for full performance
(1750 wavelengths)

Analysis Time: 0.5 to 200 seconds

Wavelength Range: 600nm–2500nm

Wavelength Accuracy: * 0.5nm

Wavelength Repeatability: *
Less than 0.015nm S.D.

Spectral Bandwidth: *
10nm + 1nm standard
slits available in other bandwidths

Photometric Range:
3.5A (1100nm–2500nm)
6.0A (600nm–1100nm)

Linearity: * 1% of reading

Noise: * 1100nm–2500nm
 25×10^{-6} Log 1/R at 0.0 Log 1/R
units S.D.
 45×10^{-6} Log 1/R at 0.75 Log 1/R
units S.D.

Stray Light: * Less than 0.1% (at
2300nm)

Specular Reflectance: *
Infinite absorber—unmeasurable
Less than 0.7% measuring carbon
black

Detector:

Lead Sulfide, 1100nm–2500nm
Silicon, 600nm–1100nm

Power: 110 V, 60 Hz
220/240 V, 50 Hz

Sample Cells and Capacities:
2 × Power Cells, 0.5–5gm (Standard)
Granular Cell 50–200gm (Optional)
Paste Cell 50–200gm (Optional)
Liquid Cell 2–30gm (Optional)
Flow Thru Cell 1–10ml (Optional)
Textile/Films 20–2cm (Optional)
50 Plastic Liners (Optional)

Computer

IBM PC/XT or AT
Enhanced Color Monitor, 12"
Memory: 10 or 20 Mbyte
Winchester Hard Disk
512 Kbyte RAM
360 Kbyte Floppy Disk

Power: 110 V, 60 Hz

Line Printer:
80 Characters Per Line
120 Cps.
RS232C Serial Interface
32 Kbyte Buffer

Calibration Programs

1. Closed Population:
Best Pair/Best Triplet

Modified Fast Forward Stepwise
Regression
Multiple Term Linear Regression

2. Open Population:
PSCo Stepwise Forward Regression
Norris Optimization

Auxiliary Programs:

Control Charts: \bar{X} and \bar{R}

Calibration Transfer


Sample Selection:
Mahalanobis Distance or Random
Spectral Library

DIMENSIONS / WEIGHTS


Spectrophotometer
Height: 41 cm (16 in.)
Width: 69 cm (27 in.)
Depth: 48 cm (19 in.)

Computer
Height: 54 cm (21.3 in.)
Width: 65 cm (25.6 in.)
Depth: 56 cm (22.0 in.)

Weight (Total unit)
159 kg (350 lbs.) net
318 kg (700 lbs.) shipping

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4600 AG bergen op zoom telefoon 01640-36022