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
Table 4. The starch content of leaf tissue samples measured by the enzyme method, the iodine precipitation method of Pucher *et al.* (1948) and the perchloric acid/anthrone method of Clegg (1956)

Tissue ^c	Starch content (%)		
	Enzyme method ^a	Pucher <i>et al.</i> method ^b	Clegg method ^b
Soybean (<i>Glycine max</i> [L.] Merr.)	7.44 ± 0.08 ^d	7.48 (± 0.08) ^d	10.84 (± 0.12) ^d
Clover (<i>Trifolium repens</i> L.)	5.57 ± 0.06	5.42 (± 0.10)	7.25 (± 0.21)
Ryegrass (<i>Lolium perenne</i> L.)	1.44 ± 0.08	1.57 (± 0.09)	4.43 (± 0.28)
Silverbeet (<i>Beta vulgaris</i> L.)	2.96 ± 0.04	3.07 (± 0.08)	4.89 (± 0.30)

^a Calculated using starch conversion factor 0.91.

^b Calculated using starch conversion factor 0.90.

^c For all tissue the quantities analysed were approximately 10-30 mg with the enzyme method, but approximately 200 mg with the other two methods.

^d Values in parentheses indicate the spread of duplicates; values not enclosed in brackets are standard errors of the means.

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Starch

UV-method

for the determination of native starch¹ in foodstuffs and other materials


Cat. No. 207 748

Test-Combination for ca. 25 determinations

0
3
5

Principle (Ref. 1-4)

Starch is hydrolyzed to D-glucose at pH 4.6 in the presence of the enzyme amyloglucosidase (AGS) (1).

The D-glucose formed is determined with hexokinase (HK) and glucose-6-phosphate dehydrogenase (G6P-DH) at pH 7.6. D-glucose is phosphorylated to glucose-6-phosphate (G-6-P) by adenosine-5'-triphosphate (ATP) in the presence of hexokinase with the simultaneous formation of ADP (2).

In the presence of G6P-DH glucose-6-phosphate is oxidized by nicotinamide-adenine dinucleotide phosphate (NADP) to gluconate-6-phosphate with the formation of reduced nicotinamide-adenine dinucleotide phosphate (NADPH) (3).

The amount of NADPH formed in the above reaction is stoichiometric with the amount of D-glucose formed by hydrolysis of the starch. NADPH is determined by means of its absorbance at 334, 340 or 365 nm.

The Test-Combination contains

1. Bottle 1 with approx. 100 mg lyophilisate, consisting of: citric acid buffer, pH 4.6; amyloglucosidase, 84 U; stabilizers.
2. Bottle 2 with approx. 5 g powder mixture, consisting of: triethanolamine buffer, pH 7.6; NADP, 75 mg; ATP, 190 mg; magnesium sulfate; stabilizers.
3. Bottle 3 with approx. 0.7 ml enzyme suspension, consisting of: hexokinase, 200 U; glucose-6-phosphate dehydrogenase, 100 U.

Preparation of solutions

1. Dissolve contents of bottle 1 with 6.0 ml redist. water.
2. Dissolve contents of bottle 2 with 27 ml redist. water.
3. Use contents of bottle 3 undiluted.

Stability of solutions

Solution 1 is stable for 6 weeks at +4°C, or for 3 months at -20°C. Bring solution 1 to 20-25°C before use.
Solution 2 is stable for 4 weeks at +4°C, or for 2 months at -20°C. Bring solution 2 to 20-25°C before use.
The contents of bottle 3 are stable for 1 year at +4°C.

Procedure

Wavelength²: 340 nm, Hg 365 nm or Hg 334 nm
Glass cuvette³: 1 cm light path
Temperature: 55-60°C (incubation); 20-25°C (measurement)
Final volume: 2.32 ml
Read against air (without a cuvette in the light path) or against water.
Sample solution: 3-70 µg starch/cuvette⁴ (in 0.1-1.0 ml sample volume, or in 0.1-0.2 ml, when analyzing DMSO-containing solutions).

If the sample contains amounts of free D-glucose, a separate determination of D-glucose is necessary. In the presence of other oligosaccharides the alcohol extraction is recommended (see instruction pt. 1.2.).

- 1 Modified kinds of starch (phosphorylated or oxidized ones) do not react.
- 2 The absorption maximum of NADPH is at 340 nm. On spectrophotometers, measurements are taken at the absorption maximum; when spectralline photometers equipped with a mercury vapour lamp are used, measurements are taken at a wavelength of 365 nm or 334 nm.
- 3 If desired, disposable cuvettes may be used instead of glass cuvettes.
- 4 See instructions for performance of the assay.
- 5 See instructions for sample preparation, pt. 1.2 and 1.3.
- 6 Available from Boehringer Mannheim GmbH.

Not for use in *in vitro* diagnostic procedures for clinical diagnosis.

biochemical analysis food analysis

Recommendations to methods and standardized procedures see references.

Pipette into cuvettes	reagent blank	sample	sample blank ⁵
solution 1*	0.20 ml	0.20 ml	-
sample solution**	-	0.10 ml	0.10 ml
redist. water	0.10 ml	-	-
mix*, incubate for 15 min at 55-60°C (water-bath); stopper cuvettes with the lid or with Parafilm®. Addition of			
solution 2	1.00 ml	1.00 ml	1.00 ml
redist. water	1.00 ml	1.00 ml	1.20 ml
mix***. After approx. 3 min read absorbances of the solutions (A ₁). Start reaction by addition of			
suspension 3	0.02 ml	0.02 ml	0.02 ml
mix***, after completion of the reaction (ca. 10-15 min) read absorbances of the solutions (A ₂). If the reaction has not stopped after 15 min, read absorbances in 5 min intervals until the absorbance increases constantly over 5 min.			

* Pipette solution 1, the sample solution and redist. water, each, onto the bottom of the cuvette and mix by gentle swirling. When using a plastic spatula, remove it from the cuvette only directly before measuring absorbance A₁.

** Rinse the enzyme pipette or the pipette tip of the piston pipette with sample solution before dispensing the sample solution.

*** For example, with a plastic spatula or by gentle swirling after closing the cuvette with Parafilm® (registered trademark of the American Can Company, Greenwich, Ct, USA).

If the absorbances A₂ increase constantly extrapolate the absorbances to the time of addition of suspension 3. Determine the absorbance differences (A₂-A₁) for both reagent blank and sample. Subtract absorbance difference of the reagent blank from the absorbance difference of the sample.

$$\Delta A = \Delta A_{\text{sample}} - \Delta A_{\text{reagent blank}}$$

The absorbance differences measured should as a rule be at least 0.100 absorbance units to achieve sufficiently accurate results (see "Instructions for performance of assay").

Calculation

According to the general equation for calculating the concentrations:

$$c = \frac{V \times MW}{\epsilon \times d \times v \times 1000} \times \Delta A \text{ [g/l]}, \text{ where}$$

V = final volume [ml]

v = sample volume [ml]

MW = molecular weight of the substance to be assayed [g/mol]
(for starch: MW_{glucose} - MW_{water} = 162.1)

d = light path [cm]

ε = absorption coefficient of NADPH at
340 nm = 6.3 [l x mmol⁻¹ x cm⁻¹]
Hg 365 nm = 3.5 [l x mmol⁻¹ x cm⁻¹]
Hg 334 nm = 6.18 [l x mmol⁻¹ x cm⁻¹]

It follows for starch:

$$c = \frac{2.32 \times 162.1}{\epsilon \times 1 \times 0.1 \times 1000} \times \Delta A = \frac{3.761}{\epsilon} \times \Delta A \text{ [g starch/l sample solution]}$$

If the sample has been diluted during preparation the result must be multiplied by the dilution factor F.

When analyzing solid and semi-solid samples which are weighed out for sample preparation, the result is to be calculated from the amount weighed:

$$\text{content}_{\text{starch}} = \frac{c_{\text{starch}} \text{ [g/l sample solution]}}{c_{\text{sample}} \text{ [g/l sample solution]}} \times 100 \text{ [g/100 g]}$$

If the **aqueous solution** contains free D-glucose, this D-glucose is determined together with the present starch. In this case with the determination of the free D-glucose solution 1 has to be replaced by redist. water in a second assay, i.e. measurement in the absence of amyloglucosidase (sample blank). For calculation we recommend to determine the free D-glucose as glycosyl concentration (MW = 162.1) in g/l according to the above-mentioned formula. This value must be subtracted from the concentration in g/l of the first assay (measurement with solution 1, in the presence of amyloglucosidase).

Instructions for performance of assay

The starch content present in the cuvette should range between 6 µg and 70 µg (measurement at 365 nm) or 3 µg and 40 µg (measurement at 340, 334 nm), respectively. The sample solution must therefore be diluted sufficiently to yield a starch concentration between 0.06 and 0.7 g/l, or 0.03 and 0.4 g/l, respectively.

Dilution table

estimated amount of starch per liter measurements at		dilution with water	dilution factor F
340 or 334 nm	365 nm		
< 0.4 g	< 0.7 g	—	1
0.4–4.0 g	0.7–7.0 g	1 + 9	10
4.0–40 g	7.0–70 g	1 + 99	100
> 40 g	> 70 g	1 + 999	1000

If the absorbance difference measured (ΔA) is too low (e.g. < 0.100), the sample solution should be prepared anew (weigh out more sample or dilute less strongly) or the sample volume to be pipetted into the cuvette can be increased up to 1.0 ml (solution in water). The volume of water added must then be reduced so as to obtain the same final volume for the sample and blank in the cuvettes. The sample volume v must be taken into account in the calculation.

1. Instructions for sample preparation

(Starch products, flours, pastries, meat-balls and other meat products, milk products and margarine, etc.)

1.1. Solubilization

Starch or starch containing products have to be pretreated before the determination in order to convert the starch into a soluble form. Because of the simple handling solubilization with dimethylsulfoxide (DMSO) it is recommended:

Homogenize sample in a powder mill or homogenizer, and pass solid samples through a sieve of 0.2 mm pore diameter. Weigh 100 mg to 1 g of the homogenized sample accurately into a 100 ml Erlenmeyer flask, add 20 ml dimethylsulfoxide and 5 ml hydrochloric acid (8 mol/l). Concerning fat containing samples it is recommended to add first the dimethylsulfoxide whereas fat free samples or samples with a low fat content are to be treated first with hydrochloric acid. Close the Erlenmeyer flask e.g. with Parafilm[®] and incubate for 30 min (in some cases for 60 min, e.g. if bread crumbs have been used for the preparation of meat balls) at 60°C in a water-bath (shaking water-bath or heatable magnetic stirrer; care should be taken that no clotting of the sample occurs; if necessary, crush with a glass rod).

Cool quickly to room temperature, add approx. 50 ml water, and adjust to pH 4–5 with sodium hydroxide (5 mol/l) under vigorous shaking (check by pH meter). Transfer to a 100 ml volumetric flask, rinse with water, fill up to the mark with water, and filter the solution, if necessary. Use 0.1 ml (at most 0.2 ml) undiluted for the assay, immediately.

Under the aforementioned conditions, no D-glucose will be released from starch.

For simplification of the sample preparation procedure with DMSO and HCl it is recommended (ref. 13) to incubate the sample with DMSO and HCl followed by the addition of 5 ml NaOH (8 mol/l), transfer into a 100 ml volumetric flask, rinsing and filling up to the mark with citrate buffer (0.112 mol/l; pH 4). (Adjustment of the pH is not necessary in that case.)

1.2. In the presence of mono- or oligosaccharides

Sample solutions containing free D-glucose or sucrose (hydrolyzed by the DMSO/HCl pretreatment), must be determined in a second assay **without** incubation with amyloglucosidase (sample blank). Take 0.1 ml of the pretreated sample solution, add 1.2 ml water and 1.0 ml solution 2, read A_1 , then add suspension 3. Continue working as stated in the scheme.

Subtract the absorbance difference of the sample blank from the absorbance difference of the sample.

$$\Delta A = \Delta A_{\text{sample}} - \Delta A_{\text{sample blank}}$$

In the presence of maltose or its homologues the alcohol extraction is to be carried out, because maltose is also splitted by amyloglucosidase. Weigh approx. 400 mg of the homogenized and sieved sample accurately into a 100 ml centrifuge beaker (the water content of the sample is to be determined and should not exceed 20%, otherwise the sample has to be dried before). Wash the sample three times with 10 ml ethanol (40%, v/v) each and centrifuge. Filter the supernatant and discard the filtrate. Combine precipitate in the centrifuge beaker with the precipitate in the filter and transfer with 4 portions of 5 ml dimethylsulfoxide each into a 100 ml Erlenmeyer flask, add 5 ml HCl (8 mol/l), and continue working as stated above (see pt. 1.1). In this case, the determination of the sample blank can be neglected, but the reagent blank has to be determined.

1.3. Determination of starch-partial hydrolysates

a) Dextrins in beer

When determining dextrins in beer, the solubilization with dimethylsulfoxide can be neglected. The sample can be used directly for the assay with a volume of 0.1 ml (up to 1.0 ml, if necessary).

The dextrin concentration is calculated as starch after subtraction of possibly present free D-glucose.

b) "Glucose syrup" in fruit juices


Normally, fruit juices contain a high excess of D-glucose in comparison with the concentration of glucose syrup (starch). The determination of starch is carried out directly in the sample solution according to the scheme (= starch sample) without preceding solubilization with dimethylsulfoxide. The hydrolysis of the starch sample is carried out with amyloglucosidase for 30 min at 20–25°C.

In another cuvette, D-glucose is determined **without** amyloglucosidase – as stated under 1.2. – (D-glucose sample, corresponds to the sample blank when starch is determined).


Subtract the absorbance difference $\Delta A_{\text{D-glucose sample}}$ from the absorbance difference $\Delta A_{\text{starch sample}}$:

$$\Delta A^*_{\text{glucose syrup}} = \Delta A_{\text{starch sample}} - \Delta A_{\text{D-glucose sample}}$$

If the ratio of D-glucose to glucose syrup is higher than 10:1, the precision of the determination of glucose syrup is decreased. To attain a higher precision of the measurement D-glucose present in the sample solution must be removed before the determination of glucose syrup. In the presence of the enzymes glucose oxidase (GOD) and catalase D-glucose is oxidized to D-gluconate:

The hydrogen peroxide is removed by catalase:

Reagents

Glucose oxidase (GOD), Cat. No. 105 139[®]

Catalase, Cat. No. 106 810[®]

Triethanolamine hydrochloride, Cat. No. 127 426[®]

MgSO₄ · 7 H₂O

NaOH, 4 mol/l

Preparation of solutions for 10 determinations

Enzyme solution:

Dissolve 5 mg (approx. 1250 U) GOD with 0.75 ml redist. water, add 0.25 ml catalase, and mix.

Buffer solution:

Dissolve 5.6 g triethanolamine hydrochloride and 0.1 g MgSO₄ · 7 H₂O with 80 ml water, adjust to pH 7.6 with NaOH (4 mol/l) and fill up to 100 ml with water.

Stability of solutions

The enzyme solution must be prepared freshly daily before use.

The buffer solution is stable for 4 weeks at +4°C.

Performance of D-glucose oxidation

Pipette into 10 ml volumetric flask	
buffer solution	2.0 ml
sample solution (with approx. 1% D-glucose)	5.0 ml
enzyme solution	0.1 ml
Pass a current of air (O ₂) through the mixture for 1 h; during the oxidation process check the pH with indicator paper and, if necessary, neutralize the acid formed with NaOH.	

To inactivate the enzymes GOD and catalase, place the volumetric flask in a boiling water-bath for 15 min, allow to cool, and fill up to the mark with water. Mix and filter, if necessary. Use the clear solution for the determination of the glucose syrup and the remaining D-glucose.

* Glucose syrup, calculated as starch.

c) "Starch syrup" in jam with a high sucrose content

For the determination of starch syrup besides extremely high concentrations of sucrose it is recommended to prepare the homogenized sample according to the conditions of the solubilization with dimethyl-sulfoxide.

Hereby, the present sucrose in the sample is completely hydrolyzed to D-glucose and D-fructose. The excess of D-glucose formed when using this method of sample preparation should be removed as stated under 1.3.b as to attain a higher precision of the measurement.

Use up to 5 g of the sample (jam) for the assay.

1.4. Determination of modified starch

Modified starch, e.g. phosphorylated or also oxidized starch, is only determinable in a limited amount. Amyloglucosidase hydrolyzes the modified starch up to a modification grade of 1% (if 1% of the D-glucose contained in the starch molecule is modified). The released D-glucose can be determined according to the usual method, whereas the modified glucose does not react.

Modified starch with a higher grade of modification is no more recognized by amyloglucosidase as a cleavable substrate.

2. Specificity

Amyloglucosidase hydrolyzes α -1,4- and α -1,6-glucan linkages in polysaccharides, like e.g. amylose, starch, dextrin, glycogen and also of glucosyl-oligosaccharides (maltose, maltotriose etc.).

3. Further applications

The method may also be used in the examination of paper (Ref. 10.) and in research when analyzing biological samples. For details of sampling, treatment and stability of the sample see Ref. 1 and 4.

4. Detection of interferences of the test system

When the enzymatic reaction is complete after the time given in "Procedure" it can be concluded in general that the reaction is not interfered. For assurance of results a re-start of the reaction (qualitatively or quantitatively) by the addition of 'standard material' (e.g. D-glucose) can be done: a further change of absorbance proves suitability of measurements.

For the detection of gross errors when performing the assays and of interfering substances in the sample material it is recommended to analyze a sample solution in a double determination with two different sample volumes (e.g. 0.10 ml and 0.20 ml): the measured absorbance differences have to be proportional to the sample volumes.

When analyzing solid samples it is recommended to weigh in two different amounts (e.g. 1 g and 2 g) into 100 ml volumetric flasks and to perform the determinations with the same sample volume: the absorbance differences have to be proportional to the amounts weighed in.

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Enzymoplossing, 20 g/l amyloglucosidase.

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D-Glucose

UV-method

for the determination of D-glucose in foodstuffs and other materials


Cat. No. 716 251

Test-Combination for ca. 3 x 40 determinations

0
3
5

Principle (Ref. 1,2)

D-Glucose is phosphorylated to glucose-6-phosphate (G-6-P) in the presence of the enzyme hexokinase (HK) and adenosine-5'-triphosphate (ATP) with the simultaneous formation of adenosine-5'-diphosphate (ADP) (1).

In the presence of the enzyme glucose-6-phosphate dehydrogenase (G6P-DH), G-6-P is oxidized by nicotinamide-adenine dinucleotide phosphate (NADP) to gluconate-6-phosphate with the formation of reduced nicotinamide-adenine dinucleotide phosphate (NADPH) (2).

The amount of NADPH formed in this reaction is stoichiometric with the amount of D-glucose. The increase in NADPH is measured by means of its absorbance at 334, 340 or 365 nm.

The Test-Combination contains

- Three bottles 1, each with approx. 7.2 g powder mixture, consisting of: triethanolamine buffer, pH 7.6; NADP, 110 mg; ATP, 260 mg; magnesium sulfate; stabilizers.
- Three bottles 2, each with 1.1 ml enzyme suspension, consisting of: hexokinase, 320 U; glucose-6-phosphate dehydrogenase, 160 U.
- Standard solution.

Preparation of solutions

- Dissolve contents of one bottle 1 with 45 ml redist. water.
- Use contents of one bottle 2 undiluted.

Stability of solutions

Solution 1 is stable for 4 weeks at +4°C, and for 2 months at -20°C.

Bring solution 1 to 20–25°C before use.

The contents of bottle 2 are stable for 1 year at +4°C.

Procedure

Wavelength¹: 340 nm, Hg 365 nm or Hg 334 nm

Glass cuvette²: 1 cm light path

Temperature: 20–25°C

Final volume: 3.02 ml

Read against air (without a cuvette in the light path) or against water.

Sample solution: 4–100 µg D-glucose/cuvette³ (in 0.1–2.0 ml sample volume)

Pipette into cuvettes	blank	sample
solution 1	1.00 ml	1.00 ml
sample solution*	–	0.10 ml
redist. water	2.00 ml	1.90 ml
mix**, read absorbances of the solutions (A ₁) after approx. 3 min and start reaction by addition of		
suspension 2	0.02 ml	0.02 ml
mix**, wait for the end of the reaction (approx. 10–15 min), and read absorbances of the solutions (A ₂). If the reaction has not stopped after 15 min, continue to read the absorbances at 2-min intervals until the absorbances increase constantly over 2 min.***		

* Rinse the enzyme pipette or the pipette tip of the piston pipette with sample solution before dispensing the sample solution.

** For example, with a plastic spatula or by gentle swirling after closing the cuvette with Parafilm[®] (registered trademark of the American Can Company, Greenwich, Ct., USA).

*** Creep reactions occur very occasionally. They are mostly brought about by the contents of the sample solutions, such as enzymes or coloring agents. These interfering substances may be removed during sample preparation.

1 The absorption maximum of NADPH is at 340 nm. On spectrophotometers, measurements are taken at the absorption maximum; when spectralline photometers equipped with a mercury vapour lamp are used, measurements are taken at a wavelength of 365 nm or 334 nm.

2 If desired, disposable cuvettes may be used instead of glass cuvettes.

3 See instructions for performance of the assay.

Not for use in *in vitro* diagnostic procedures for clinical diagnosis

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If the absorbances A₂ increase constantly, extrapolate the absorbances to the time of addition of suspension 2.

Determine the absorbance differences (A₂–A₁) for both blank and sample. Subtract the absorbance difference of the blank from the absorbance difference of the sample, thereby obtaining ΔA_{D-glucose}.

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V = final volume [ml]

v = sample volume [ml]

MW = molecular weight of the substance to be assayed [g/mol]

d = light path [cm]

ε = absorption coefficient of NADPH at:

$$340 \text{ nm} = 6.3 \text{ [l} \times \text{mmol}^{-1} \times \text{cm}^{-1}]$$

$$\text{Hg } 365 \text{ nm} = 3.5 \text{ [l} \times \text{mmol}^{-1} \times \text{cm}^{-1}]$$

$$\text{Hg } 334 \text{ nm} = 6.18 \text{ [l} \times \text{mmol}^{-1} \times \text{cm}^{-1}]$$

It follows for D-glucose:

$$c = \frac{3.02 \times 180.16}{\epsilon \times 1 \times 0.1 \times 1000} \times \Delta A = \frac{5.441}{\epsilon} \times \Delta A \text{ [g D-glucose/l sample solution]}$$

If the sample has been diluted during preparation, the result must be multiplied by the dilution factor F.

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$$\text{content}_{\text{D-glucose}} = \frac{c_{\text{D-glucose}} \text{ [g/l sample solution]}}{c_{\text{sample}} \text{ [g/l sample solution]}} \times 100 \text{ [g/100 g]}$$

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The amount of D-glucose present in the cuvette should range between 8 µg and 100 µg (measurement at 365 nm) or 4 µg and 50 µg (measurement at 340, 334 nm), respectively. The sample solution must therefore be diluted sufficiently to yield a D-glucose concentration between 0.08 and 1.0 g/l or 0.04 and 0.5 g/l, respectively.

Dilution table

estimated amount of D-glucose per liter measurements at		dilution with water	dilution factor F
340 or 334 nm	365 nm		
< 0.5 g	< 1.0 g	–	1
0.5– 5.0 g	1.0–10.0 g	1 + 9	10
5.0–50.0 g	10.0–100 g	1 + 99	100
> 50 g	> 100 g	1 + 999	1000

If the absorbance difference measured (ΔA) is too low (e.g. < 0.100), the sample solution should be prepared anew (weigh out more sample or dilute less strongly) or the sample volume to be pipetted into the cuvette can be increased up to 2.0 ml. The volume of water added must then be reduced so as to obtain the same final volume for the sample and blank in the cuvettes. The new sample volume v must be taken into account in the calculation.

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1.1. Liquid foodstuffs

Dilute clear, nearly colorless solutions according to the dilution table and use them for the assay.

Neutralize strongly acid or strongly alkaline solutions which are to be used undiluted for the assay.

Decolorize strongly colored solutions, which are to be used undiluted for the assay, with polyamide (1 g/100 ml).

Degas carbonic acid containing solutions by whirling, by shaking or by filtration.

Filter turbid solutions (folded filter or membrane filter) or centrifuge and dilute the filtrate or the supernatant, respectively, according to the dilution table, if necessary.

Example:

Determination of D-glucose in milk

Pipette 20 ml milk into a 100 ml volumetric flask, add the following solutions and shake vigorously after each addition: 10 ml Carrez-I-solution (see 1.2.), 10 ml Carrez-II-solution (see 1.2.) and 20 ml sodium hydroxide (0.1 mol/l). Fill up to the mark with dist. water, mix and filter. Use 1.0 ml or 2.0 ml, respectively, of the filtrate for the assay.

1.2. Pasty and solid foodstuffs

Homogenize pasty samples, grind and homogenize solid foodstuffs. Afterwards extract with water, heated to 60°C, if necessary. Remove turbidity by filtration.

Concerning fat containing samples the fat is to be frozen in the cold (ice-bath or refrigerator) after the extraction and is then removed by filtration.

Deproteinize protein and enzyme containing samples with Carrez solutions (in the presence of sucrose, maltose or lactose deproteinization of the sample with acid – trichloroacetic acid or perchloric acid – is to be avoided because of the hydrolysis of these sugars and the release of D-glucose):

Weigh approx. 1 g sample accurately into a 100 ml volumetric flask, add approx. 60 ml water and keep at 70°C for 15 min. Shake from time to time. For protein precipitation add 5 ml Carrez-I-solution (3.60 g potassium hexacyanoferrate-II, $K_4[Fe(CN)_6] \cdot 3 H_2O/100$ ml), 5 ml Carrez-II-solution (7.20 g zinc sulfate, $ZnSO_4 \cdot 7 H_2O/100$ ml) and 10 ml NaOH (0.1 mol/l), shake vigorously after each addition, adjust to room temperature, fill up to the mark with water, and filter. Use the clear, possibly slightly turbid solution for the assay, dilute previously, if necessary.

2. Specificity

The method is specific for D-glucose.

3. Further applications

The method may also be used in the examination of cosmetics (Ref. 10), pharmaceuticals (Ref. 11), paper (Ref. 12), and tobacco (Ref. 13) and in research when analyzing biological samples. For sampling, treatment and stability of the sample see Ref. 1,2.

Examples:

3.1. Determination of D-glucose in urine and liquor (Ref. 2)

Dilute sample according to the dilution table and use the diluted sample for the assay (dilution factor = F).

Calculation:

$$c = \frac{5.441 \times \Delta A \times F}{\epsilon} \text{ [g D-glucose/l sample]}$$

$$c = \frac{30.2 \times \Delta A \times F}{\epsilon} \text{ [mmol D-glucose/l sample]}$$

Wavelength	Hg 365 nm	340 nm	Hg 334 nm
c [g/l]	$1.555 \times \Delta A \times F$	$0.8637 \times \Delta A \times F$	$0.8804 \times \Delta A \times F$
c [mmol/l]	$8.629 \times \Delta A \times F$	$4.794 \times \Delta A \times F$	$4.887 \times \Delta A \times F$

3.2. Determination of D-glucose in blood, serum and plasma (Ref. 2)

Mix 1.0 ml perchloric acid (0.33 mol/l) and 0.1 ml sample in a centrifuge tube, and centrifuge. Use the supernatant for the assay.

Further information see instructions of

Test-Combination D-Glucose/D-Fructose Cat. No. 139 106

Test-Combination D-Glucose/D-Fructose/
D-Sorbitol Cat. No. 724 831

Test-Combination Sucrose/D-Glucose Cat. No. 139 041

Test-Combination Sucrose/D-Glucose/
D-Fructose Cat. No. 716 260

Test-Combination D-Sorbitol/Xylitol Cat. No. 670 057

Test-Combination Starch Cat. No. 207 748

The dilution factor F (depending on sample preparation) is obtained from the sample volume (0.1 ml), the perchloric acid volume (1.0 ml), the specific gravity of the sample material (1.06 g/ml blood, 1.03 g/ml serum or plasma) and the fluid content (0.80 in case of blood and 0.92 in case of serum or plasma):

$$F_{\text{blood}} = \frac{0.1 \times 1.06 \times 0.80 + 1.0}{0.1} = 10.85$$

$$F_{\text{serum, plasma}} = \frac{0.1 \times 1.03 \times 0.92 + 1.0}{0.1} = 10.95$$

Calculation:

$$c = \frac{5.441 \times \Delta A \times F}{\epsilon} \text{ [g D-glucose/l sample]}$$

$$c = \frac{30.2 \times \Delta A \times F}{\epsilon} \text{ [mmol D-glucose/l sample]}$$

D-Glucose in blood:

Wavelength	Hg 365 nm	340 nm	Hg 334 nm
c [g/l]	$16.87 \times \Delta A$	$9.371 \times \Delta A$	$9.553 \times \Delta A$
c [mmol/l]	$93.62 \times \Delta A$	$52.01 \times \Delta A$	$53.02 \times \Delta A$

D-Glucose in serum or plasma:

Wavelength	Hg 365 nm	340 nm	Hg 334 nm
c [g/l]	$17.02 \times \Delta A$	$9.457 \times \Delta A$	$9.641 \times \Delta A$
c [mmol/l]	$94.48 \times \Delta A$	$52.49 \times \Delta A$	$53.51 \times \Delta A$

3.3. Determination of D-glucose in fermentation samples and cell culture media

Place the sample (after centrifugation, if necessary) into a water-bath at 80°C for 15 min to stop enzymatic reactions. Centrifuge and use the supernatant (diluted according to the dilution table, if necessary) for the assay. Alternatively, deproteinization can be carried out with perchloric acid, but only by absence of disaccharides, (see 1.2.) or with Carrez-solutions. See the above-mentioned examples.

Homogenize gelatinous agar media with water and treat further as described.

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D-Glucose standard solution

for the Test-Combination D-Glucose
UV-method, Cat. No. 716 251

Concentration*: see bottle label.

D-Glucose standard solution is a stabilized aqueous solution of D-glucose. It serves as standard solution for the enzymatic analysis of D-glucose in foodstuffs and other materials.

Application

1. Addition of D-glucose standard solution to the assay mixture:

Instead of the sample solution the standard solution is used for the assay.

2. Restart of the reaction, quantitatively:

After completion of the reaction with sample solution and measuring of A_2 , add 0.05 ml standard solution to the assay mixture. Read absorbance A_3 after the end of the reaction (approx. 15 min). Calculate the concentration from the difference of $(A_3 - A_2)$ according to the general equation for calculating the concentration. The altered total volume must be taken into account. Because of the dilution of the assay mixture by addition of the standard solution, the result differs insignificantly from the data stated on the bottle label.

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©1989