Afd. Koolhydraat- en Vetchemie

RAPPORT 83.48

1983-06-20

Pr.nr. 505.3000

Onderwerp: Bepaling van glucosinolaat
(isothiocyanaat (ITC) + vinyl-
thio-oxazolidon (VTO)).

Tabellen: 12.

Bijlagen: 6.

Verzendlijst: directeur, sektorhoofd (2x), direktie VKA, afd.
Koolhydraat- en Vetchemie (4x), afd. Normalisatie
(Humme), Projektbeheer, Projektleider (Muuse).

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glucosinolaat = Δ glucose x 2,03

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Revised November 1980

THE DETERMINATION OF VINYLOXAZOLIDINE-THIONE

1. Purpose and Scope

This method describes a gas-liquid chromatographic and a high performance liquid chromatographic procedure for the determination of vinyloxazolidine-thione (VOT) in animal feedingstuffs. *The lower limit of determination is of 200 mg/kg.*

2. Principle

The sample is submitted to an enzymatic treatment which liberates 2-hydroxy-3-butenyl isothiocyanate. This compound is transformed by a non enzymatic cyclization to VOT which is determined by gas-liquid chromatography or by high performance liquid chromatography.

3. Reagents

- 3.1 Light petroleum, boiling range 40-60°C and a bromine value less than 1.
- 3.2 Chloroform.
- 3.3 Citric acid solution: Dissolve 2.10 g of citric acid monohydrate in water and dilute to 100 ml.
- 3.4 Disodium hydrogen phosphate solution: dissolve 7.16 g of disodium hydrogen phosphate dodecahydrate in water and dilute to 100 ml with water.
- 3.5 Buffer solution, pH 7.0: mix 82.4 ml of disodium hydrogen phosphate solution (3.4) with 17.6 ml of citric acid solution (3.3).
- 3.6 White mustard (*Sinapis alba*). Whole seeds not more than 2 years old. Before use grind finely and remove oil and fat by cold extraction with light petroleum. The prepared ground mustard must contain active thioglucosidase (myrosinase) and not more than 500 mg/kg VOT.
- 3.7 Standard substance: pure vinyloxazolidine-thione.
- 3.8 Standard solution: weigh to the nearest 0.1 mg, 25 mg of pure vinoxazolidine-thione (3.7). Dissolve in chloroform (3.2) in a 50 ml graduated flask, make up to volume with the same solvent and mix. 1 ml of this solution contains 500 ug of VOT.

4. Apparatus

- 4.1 Laboratory mill or grinder.
- 4.2 Water bath controlled at 45°C.
- 4.3 Centrifuge tubes, capacity 100 ml, with efficient stoppers.
- 4.4 Gas chromatograph suitable for on-column injection fitted with a flame ionization detector.

The following conditions are offered for guidance:

- 4.4.1 Column: glass internal diameter 2-6 mm, and length 1.5-2.5 m.
- 4.4.2 Support: acid washed and silanized
(eg Gas Chrom Q or Diatomite CLQ). Mesh sizes of 80-100
(150-190 microns) and 100-120 (125-150 microns) are suitable.
- 4.4.3 Stationary phase: 5% OV 17
- 4.4.4 Column temperature: A temperature of 160 to 175°C is
recommended.
- 4.4.5 Injector - direct on-column injection is recommended.
- 4.4.6 Detector - this should be maintained at least 20°C above the
temperature of the column.
- 4.4.7 Carrier gas - nitrogen, flow rate of 30 ml min⁻¹.
- 4.5 High performance liquid chromatograph system
 - 4.5.1 Column: stainless steel, internal diameter 4.6 mm and
length 15 cm.
 - 4.5.2 Column packing material: Silica, average particle size 5 µm.
(eg Partisil 5, Lichrosorb SI 60 or Spherisorb 5).
 - 4.5.3 HPLC mobile phase: 2,2,4-Trimethylpentane/Ethanol 4:1,
delivered at a flow-rate of 1 ml min⁻¹.
 - 4.5.4 Detector: suitable ultraviolet absorptiometer operating in the
region of 245 nm.

5. Procedure

5.1 Preparation of the sample

Grind a sample portion so that it passes through a 1 mm sieve without residue. This operation should be carried out as rapidly as possible and precautions taken to avoid the generation of excessive heat. Analyse the ground sample without further delay.

5.2 Preparation of extracts

Weigh to the nearest 0.005 g, approximately 2.0 g of the prepared sample (5.1) and transfer to a centrifuge tube (4.3). Samples which contain more than 5% of crude oil should be defatted by adding 25 ml of light petroleum (3.1), mixing and shaking for 2 - 3 minutes. Allow to settle, centrifuge and decant the supernatant layer. Add a further portion of 25 ml of light petroleum and mix, shake and centrifuge as before. Decant the supernatant layer and remove the last traces of solvent in a current of air until the test portion is dry and friable.

Add 15.0 ml of buffer solution (3.5) and 0.20 g of white mustard (3.6) accurately weighed to the nearest 0.01 g. Stopper the centrifuge tube

Securely and shake by hand until the solid matter is dispersed. Transfer the tube to a water bath maintained at 45°C (4.2) and leave to stand for 90 minutes. Remove the tube from the bath and allow to cool. Add by pipette 40.0 ml of chloroform (3.2) and shake the stoppered tube vigorously by hand for two minutes. Centrifuge until the separated layers are clear. Remove part of the chloroform layer to a dry stoppered tube and retain for the chromatography which should be carried out without unnecessary delay.

5.3 Determination

5.3.1 By Gas Chromatography

Inject 5.0 µl of extract obtained in (5.2) into the gas chromatograph and measure the height of the VOT peak. Determine the concentration of VOT (µg per ml) in the extract by reference to the calibration curve (5.4.1).

5.3.2 By High Performance Liquid Chromatography

Inject 4.0 µl of extract obtained in (5.2) onto the HPLC column and measure the height of the VOT peak. Determine the concentration of VOT (µg per ml) in the extract by reference to the calibration curve (5.4.2).

5.4 Calibration curve

5.4.1 By Gas Chromatography

Into a series of 25 ml graduated flasks add by pipette 1.0, 2.0, 4.0, 6.0 and 8.0 ml of VOT standard solution (3.8) and make up to the mark with chloroform (3.2). These solutions contain respectively 20.0, 40.0, 80.0, 120.0 and 160.0 µg of VOT per ml.

Inject 5.0 µl of each standard into the gas chromatograph and measure the heights of the VOT peaks. Plot a graph of peak height against concentration of VOT (µg per ml) in the standard solutions.

5.4.2 By High Performance Liquid Chromatography

Prepare a series of VOT standards as in (5.4.1) and inject 4.0 µl of each solution onto the HPLC column. Measure the heights of the VOT peaks and plot a graph of peak height against concentration of VOT (µg per ml) in the standard solutions.

5.5 Blank test

Carry out a blank test using the procedure described in (5.2) and omitting only the sample. Allow for this in the calculation of the results.

6. Calculations

The VOT content of the feed in mg/kg is given by the formula:

$$\frac{C \times 40}{W}$$

where C = concentration of VOT in ug per ml as determined from the calibration graph

W = weight of sample in grams.

Notes

1. An alternative stationary phase for gas chromatography is EGSS-X. However some feeds gave interfering peaks on this column under the conditions tested.

2. Using the following gas chromatographic columns at 165°C:

column 1: 5% OV.17 on Diatomite CLQ 80-100 mesh

column 2: 1% EGSS-X on Diatomite CLQ 80-100 mesh

the retention times of VOT were found to be:

column 1: 9.0 min

column 2: 8.0 min.

3. Using the following HPLC column:

Lichrosorb SI 60 the retention time of VOT was found to be 7.5 min.

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Oilseeds and oilseed residues — Determination of isothiocyanates and vinyl thiooxazolidone

1 Scope and field of application

This International Standard specifies a method for the determination of the isothiocyanates (ITC) and vinyl thiooxazolidone (VTO) produced by the enzymic hydrolysis of glucosinolates in seeds and oilseed residues of colza and rapeseed.

Under the operating conditions described, the method does not allow determination of free isothiocyanates and free vinyl thiooxazolidone.

It also gives, in an annex, a method of determining isothiocyanates by argentimetry which may be substituted for the gas chromatographic method, by agreement between the parties. This method is not applicable, however, to oilseeds or oilseed residues having low glucosinates contents.

2 References

ISO 659, *Oilseeds — Determination of hexane extract (or light petroleum extract), called "oil content"*.

ISO 664, *Oilseeds — Reduction of contract samples to analysis samples*.

ISO 734, *Oilseed residues — Determination of hexane extract (or light petroleum extract), called "oil content"*.

3 Principle

After removal of the fats and oil (de-fatting) if necessary, and drying of the oilseed or oilseed residue, enzymatic hydrolysis of the glucosinolates, extraction of isothiocyanates with dichloromethane and of vinyl thiooxazolidone with diethyl ether.

Determination of the ITC by gas chromatography and of the VTO by ultraviolet spectrophotometry.

4 Reagents

The reagents shall be of recognized analytical quality and the water used shall be distilled water or water of at least equivalent purity.

4.1 Reagents for determination of ITC

4.1.1 Dichloromethane, or, failing this, chloroform.

4.1.2 *n*-Hexane, or, failing this, light petroleum (boiling range 40 to 60 °C).

4.1.3 Buffer solution of pH 7, commercially available, or, for example, a solution prepared as follows :

Pour 35,3 ml of 0,1 mol/l citric acid ($C_6H_8O_7 \cdot H_2O$) solution (21,01 g/l solution) into a 200 ml one-mark volumetric flask and make up to the mark with a 0,2 mol/l solution of disodium hydrogen orthophosphate dihydrate ($Na_2HPO_4 \cdot 2H_2O$) (35,61 g/l solution). Check the pH and, if necessary, adjust.

4.1.4 Enzyme source, prepared from white mustard seeds (*Sinapis alba* L.).¹⁾

Finely grind the white mustard seeds so that at least 80 % will pass through a sieve of aperture size 280 µm.

Remove the oil and fat from the grindings using cold hexane or, failing this, cold light petroleum (4.1.2), carrying out the extraction using an apparatus which will permit extraction until not more than 2 % of oil remains in the product and without the temperature exceeding 30 °C, for example a double-walled Soxhlet apparatus, or by carrying out extractions by grinding with hexane in a microgrinder cooled by running water. Remove traces of solvent at ambient temperature, preferably using a slight current of air.

Store the enzyme source thus obtained at 4 °C in a hermetically sealed glass bottle. In these conditions, it can be kept for about 6 weeks, but it is preferable to use a freshly prepared enzyme source.

1) Use seed in which more than 85 % of the grains should germinate in less than 72 h.

It is recommended that a blank test be performed to ensure that the enzyme source does not contain ITC.

4.1.5 Butyl isothiocyanate, standard solution.

Prepare a 0,500 g/l solution of butyl isothiocyanate in the dichloromethane or chloroform (4.1.1). Store at 4 °C.

This solution may be diluted if necessary, according to the expected isothiocyanates content of the sample.

4.1.6 Gases for gas chromatography.

Carrier gas : carefully dried nitrogen containing less than 10 mg of oxygen per kilogram.

Auxiliary gases : hydrogen and air.

4.2 Reagents for the determination of VTO

4.2.1 Diethyl ether, spectrophotometric quality.

4.2.2 Anti-foaming agent, for example octan-2-ol.

4.2.3 Buffer solution, of pH 7 (see 4.1.3).

4.2.4 Enzyme source (see 4.1.4).

5 Apparatus

Usual laboratory equipment and in particular

5.1 Apparatus for preparing the test sample and drying the test portion

5.1.1 Sieve, of aperture size 280 µm.

5.1.2 Apparatus for the extraction of fat and oil by the method specified in ISO 659 or ISO 734, if necessary.

5.1.3 Grinder.

5.1.4 Desiccator.

5.1.5 Analytical balance.

5.1.6 Electric oven, capable of being controlled at 103 ± 2 °C.

5.1.7 Conical flask, of capacity 25 ml.

5.2 Apparatus for the determination of ITC

5.2.1 Chromatograph, with a flame ionization detector and a recorder, comprising, for example, a column of length 2 m, external diameter 3,2 mm and internal diameter 2 mm, filled

with 10 % diethylene glycol succinate (DEGS) on GAS/CHROMOSORB P 150 to 180 µm (80 to 100 mesh) (treated with hexamethyldisilazane).

5.2.2 Stirrer/shaker system, for conical flasks, or a magnetic stirrer.

5.2.3 Microsyringe, of capacity 1 µl.

5.2.4 Pipettes, of capacities 5 ml to 10 ml.

5.3 Apparatus for the determination of VTO

5.3.1 Spectrophotometer, preferably with a recorder, suitable for measurements in the ultraviolet.

5.3.2 Stirrer/shaker system, for conical flasks, or a magnetic stirrer.

5.3.3 Conical flasks, of capacities 100 and 200 ml.

5.3.4 One-mark volumetric flasks, of capacities 25 and 100 ml.

5.3.5 Beaker, of capacity 50 ml.

5.3.6 Separating funnel, of capacity 50 ml.

5.3.7 Pipette, of capacity 1 ml.

6 Procedure

6.1 Preparation of the test sample

6.1.1 Oilseeds

After reducing the laboratory sample in accordance with ISO 664, take about 10 g of the oilseed and de-fat it using the procedure specified in ISO 659. After removing the solvent, grind if necessary, so that at least 80 % of the grindings obtained pass through the sieve (5.1.1).

6.1.2 Oilseed residues

As de-fatting is not necessary for oil contents below 5 %, use the oilseed residues as received after grinding, if necessary, so that at least 80 % of the grindings obtained pass through the sieve (5.1.1).

As results are expressed relative to the de-fatted sample, it is necessary to perform a parallel determination of the oil content of the oilseed residue, using the method described in ISO 734 or any other method.

However, if the oil content is greater than 5 % (pressed oilseed residues), carry out de-fatting of 10 g of oilseed residue using the procedure specified in ISO 734. After elimination of the solvent, grind, if necessary.

6.2 Determination of α ITC

6.2.1 Test portion

Transfer approximately 2.2 g of the test sample (6.1) to a 25 ml conical flask (5.1.7), which has been previously dried and weighed to the nearest 1 mg. Place in the oven (5.1.6), controlled at $103 \pm 2^\circ\text{C}$, for at least 8 h, and allow to cool in the desiccator (5.1.4) to room temperature.

NOTE — Simultaneously, dry the test portion to be used for the determination of VTO in the beaker if this analysis has been requested (see 6.3.1).

Weigh the conical flask to the nearest 1 mg, and determine the mass of the de-fatted, if necessary, and dried test portion (about 2 g).

6.2.2 Determination

Add to the test portion, 5 ml of the buffer solution (4.1.3), 0.25 g of the enzyme source (4.1.4) and 10 ml, using the pipette (5.2.4), of the standard butyl isothiocyanate solution (4.1.5) (using a more dilute solution, if necessary).

Shake for 2 h using the stirrer/shaker system or the magnetic stirrer (5.2.2). Allow to separate, or centrifuge if necessary.

Using the microsyringe (5.2.3), take 1 μl of the dichloromethane or chloroform phase, avoiding taking any particles in suspension, and inject this into the chromatograph (5.2.1) controlled, for example, as follows:

temperature of the injector : 135°C

temperature of the column : 110°C

temperature of the detector : 130°C

pressure or flow rate of the carrier gas : 80 kPa* or flow rate of 20 to 30 ml/min.

The order of elution is as follows:

allyl isothiocyanate,

butyl isothiocyanate,

butenyl isothiocyanate,

pentenyl isothiocyanate.

NOTES

1 For information, if using diethylene glycol succinate (DEGS) as the stationary phase, the retention times relative to butyl isothiocyanate are as follows:

allyl isothiocyanate : 0.70

butenyl isothiocyanate : 1.45

pentenyl isothiocyanate : 2.45

2 As the solution injected is very corrosive, clean the microsyringe immediately after use with a solvent.

6.3 Determination of VTO

6.3.1 Test portion

Transfer approximately 2.2 g of the test sample (6.1) to a 50 ml beaker (5.3.5), which has been previously dried and weighed to the nearest 1 mg. Place in the oven (5.1.6), controlled at $103 \pm 2^\circ\text{C}$, for at least 8 h, and allow to cool in the desiccator (5.1.4) to room temperature.

Weigh the beaker to the nearest 1 mg, and determine the mass of the de-fatted, if necessary, and dried test portion (about 2 g).

6.3.2 Determination

Quantitatively transfer the test portion into a 200 ml conical flask (5.3.3) and add 70 ml of boiling buffer solution (4.2.3), using some of which to rinse the beaker. Allow to cool to about 30°C , then add 0.50 g of the enzyme source (4.2.4) and a few drops of the anti-foaming agent (4.2.2). Shake for 2 h using the stirrer/shaker system or the magnetic stirrer (5.3.2).

Immediately transfer quantitatively to a 100 ml one mark volumetric flask (5.3.4), rinsing with water, and make up to the mark with water. Filter and collect the filtrate in a 100 ml conical flask (5.3.3). Shake gently for 30 s and, by means of the pipette (5.3.7), take 1 ml and place it in the 50 ml separating funnel (5.3.6). Carry out two extractions of the VTO using 10 ml of the diethyl ether (4.2.1) each time. Collect the ether layers in a 25 ml one-mark volumetric flask (5.3.4), and make up to the mark with diethyl ether.

Determine the absorption curve from 220 to 280 nm and subtract the absorbance read at 280 nm from that read at the maximum absorbance, which occurs at about 250 nm, in order to obtain the absorbance of the test portion.

If the absorbance values obtained are beyond the limits of the instrument, dilute as necessary with diethyl ether.

6.3.3 Blank test

Carry out a blank test under the same conditions, omitting the test portion in order to determine the absorbance due to the enzyme source.

7 Expression of results

NOTE — For animal feeding stuffs, the results are generally expressed relative to the product as received. Thus, if such a manner of expression of results is required, make the necessary corrections to take into account the water and oil contents.

7.1 ITC content

The ITC content, expressed in milligrams per gram of dry matter of the de-fatted sample, is equal to

$$\frac{m_a}{115,19 \times S_e \times m} (99,15 \times S_b + 113,18 \times S_b + 0,9 \times 127,21 \times S_p)$$

where

m is the mass, in grams, of the de-fatted and dried test portion (6.2.1);

m_a is the mass, in milligrams, of butyl isothiocyanate contained in 10 ml of solution (usually 5 mg);

S_a is the area of the peak corresponding to butyl isothiocyanate;

S_b is the area of the peak corresponding to allyl isothiocyanate;

S_b is the area of the peak corresponding to butenyl isothiocyanate;

S_p is the area of the peak corresponding to pentenyl isothiocyanate.

OTE — The value 0,9 is the response coefficient of pentenyl isothiocyanate.

7.2 VTO content

The VTO content, expressed in milligrams per gram of dry matter of the de-fatted sample, is equal to

$$(A_E - A_B) C_p \times 25 \times 100 \times \frac{1}{m} \times 10^{-3}$$

where

A_E is the absorbance of the test solution (6.3.2);

A_B is the absorbance of the blank test solution (6.3.3);

C_p is a coefficient of proportionality, equal to 8,20;

m is the mass, in grams, of the test portion (6.3.1).

A simplified formula may thus be derived

$$\frac{20,5 \times (A_E - A_B)}{m} = \frac{C}{M} \cdot K$$

Take into account any dilution of the test solution.

8 Test report

The test report shall show the method used and the results obtained. It shall also mention any operating details not specified in this International Standard, or regarded as optional, together with any incidents likely to have affected the results.

The test report shall give all the information necessary for the complete identification of the sample.

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2. Referentiemethoden

ISO 659 Oilseeds - Determination of hexane extract (or light petroleum extract), called "oil content".

ISO 664 Oilseeds - Reduction of contract samples to analysis samples.

ISO 734 Oilseed residues - Determination of hexane extract (or light petroleum extract), called "oil content".

ISO/DIS 5504 Oilseeds and oilseed residues - Determination of isothiocyanates and vinyl-thiooxazolidone.

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