provided by UKM Journal Article Repository

Volume 6, No. 1 (Supplement) June 2011 ISSN 1823-2140

National University



MEDICINE Health

The Official Journal of The Faculty of Medicine UKM



THE HbA2 LEVELS IN NORMAL, B -THALASSAEMIA AND HAEMOGLOBIN E HETEROZYGOTES AND HAEMOGLOBIN E (HbE) LEVEL IN HbE HETEROZYGOTES BY CAPILLARY ELECTROPHORESIS SYSTEM

Hafiza A¹, <u>Azma RZ¹</u>, Chong MKT¹, Irwan MA¹, Yeoh Z-N¹, Lailyvia MI¹, Nur Rabiatuladawiah MR¹, Aidilfitrina M², Malissa Y², Azlin I¹, Hamidah NH¹.

Department of ¹Pathology and ²Diagnostic Laboratory Services, Universiti Kebangsaan Malaysia Medical Centre, Kuala Lumpur, Malaysia.

Background:

Capillary electrophoresis (CE) is a technology that has recently been applied to thalassaemia screening and diagnosis. It utilizes the technique of electrokinetic separation of molecules in eight electrolyte buffer-filled capillaries. Our study aimed to establish the normal ranges of HbA $_2$ and HbF in normal population and to quantify the values of HbA $_2$ and HbE in the β -thalassaemia and haemoglobin E (HbE) heterozygous states.

Methodology:

A total of 173 normal individuals, 218 β -thalassaemia trait and 92 HbE trait cases were selected. All samples were analysed by both CE and High Performance Liquid Chromatography (HPLC) techniques.

Results:

Our results showed that the mean HbA $_2$ and HbF were 2.75% (SD 0.25%) and 0.06% (SD 0.24%) respectively, which was significantly lower than that of HPLC, 2.88%(SD 0.25%) and 0.31% (SD 0.61%) (p<0.001 for both results). For β -thalassaemia heterozygotes, the HbA $_2$ level was slightly higher than that of HPLC, where CE recorded values of 5.23% (SD 0.63%) versus 5.14% (SD 0.55%) by HPLC (p<0.001). The HbA $_2$ level for HbE heterozygotes was 3.61% (SD 0.44%), which was significantly higher than normal but lower than that of β -thalassaemia heterozygotes (P<0.00 for both). Peculiar to CE, the HbE level in HbE heterozygotes was determined as 25.41% (SD 3.38%). The correlation of HbA $_2$ levels between both methods was excellent (R 2 =0.99).

Conclusion:

The CE system was a practical option for medical laboratories for analysis of thalassaemia and haemoglobinopathies. We established the ranges of HbA2 for the normal population, β -thalassaemia and HbE heterozygotes. Unique to the CE system, the HbE level was also able to be ascertained.

Keywords:

Capillary electrophoresis, HPLC, beta thalassaemia, HbE, HbA2