

Technical University of Denmark



Improved detection of *Salmonella* through the optimized use of DNA polymerases in real-time PCR

Löfström, Charlotta; Fachmann, Mette Sofie Rousing; Josefsen, Mathilde Hasseldam; Al-Habib, Zahra Fares Sayer; Hoorfar, Jeffrey

Publication date:
2012

Document Version
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

Citation (APA):
Löfström, C., Søndergaard, M. S. R., Josefsen, M. H., Al-Habib, Z. F. S., & Hoorfar, J. (2012). Improved detection of *Salmonella* through the optimized use of DNA polymerases in real-time PCR. Poster session presented at 23rd International ICFMH Symposium, FoodMicro 2012, Istanbul, Turkey.

DTU Library
Technical Information Center of Denmark

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

Improved detection of *Salmonella* through the optimized use of DNA polymerases in real-time PCR

Charlotta Löfström, Mette S.R. Søndergaard, Mathilde H. Josefsen, Zahra Fares Al-Habib, Jeffrey Hoorfar

National Food Institute, Technical University of Denmark, Søborg.

E-mail: chalo@food.dtu.dk

Aim

To investigate the performance of different DNA polymerases in real-time PCR and their resistance towards inhibitors in a food model system. As a model detection of *Salmonella* in meat samples with a validated real-time PCR method was used.

Introduction

To improve food safety it is important to pursue fast, well performing and low-cost methods for detection of foodborne pathogens. Though real-time PCR offers several advantages compared with classical microbiology, the choice of a suitable DNA polymerase has been shown to optimize method performance considerably¹.

Conclusions

- The performances of the tested DNA polymerases varied considerably, reinforcing the importance of careful selection of an appropriate DNA polymerase for the PCR assay and sample type in question.
- Under the tested conditions for detecting *Salmonella* in minced meat samples, VeriQuest and TaKaRa Ex Taq HS were found to be the two best performing alternatives DNA polymerases.
- The experiments will be repeated using another target bacterium and sample type to assess the generality of the results.


Screening:

- Standardized *Salmonella* qPCR assay²
- Detecting a 10-fold dilution series of purified DNA
- 13 DNA polymerases and 2 DNA polymerases in kits were included

Good performance:

- AmpliTaq 360
- FastStart Taq
- HotMaster Taq
- HotStarTaq Master Mix Kit
- TaKaRa Ex Taq HS
- Tth
- VeriQuest Probe qPCR Master Mix


Figure 1. Example of PCR results using a well performing DNA polymerase (TaKaRa Ex Taq HS)



Intermediate performance:

- MyTaq
- MyTaq HS
- Phusion High-Fidelity
- PicoMaxx High Fidelity PCR System


Figure 2. Example of PCR results using an intermediate performing DNA polymerase (PicoMaxx)



Bad performance:

- Herculase II Fusion
- OneTaq
- Pfu
- TITANIUM Taq

Figure 3. Example of PCR results using a poor performing DNA polymerase (Pfu)



Further evaluation:

- 5 selected well performing DNA polymerases (Table 1)
- Standardized *Salmonella* qPCR assay²
- Minced meat samples artificially contaminated with *Salmonella* (10¹-10⁶ CFU/ml) analyzed in triplicate after DNA extraction
- Ct values, amplification efficiency (AE), linear range and linearity of standard curve (R²) were calculated

Table 1. Results from qPCR using different DNA polymerases when analysing minced meat samples enriched in BPW for 18 h followed by artificial contamination with *Salmonella*.

DNA polymerase	Ct values (mean of triplicates ^a) from qPCR for indicated no. of <i>Salmonella</i> (CFU/ml)						AE	R ²	Linear range (CFU/ml)
	10 ¹	10 ²	10 ³	10 ⁴	10 ⁵	10 ⁶			
	Tth	ND ^b	ND	35.4	32.6	27.3			
VeriQuest Probe qPCR Master Mix	ND	39.3 ^c	33.5	32.1	27.8	25.6	1.10	0.953	10 ³ -10 ⁶
FastStart Taq	ND	ND	36.3	35.7	29.5	29.3	1.45	0.874	10 ⁴ -10 ⁵
TaKaRa ExTaq HS	ND	ND	36.6	34.9	29.3	27.3	0.97	0.929	10 ³ -10 ⁶
AmpliTaq 360	ND	ND	35.4	32.0	27.3	25.5	0.94	0.842	10 ⁴ -10 ⁵

^a Standard deviation 0.1-1.9 in the trial for the different triplicate analyses

^b ND: Not Detected

^c Detection in 1/3 replicates

Materials and Methods

13 commercially available DNA polymerases and 2 polymerases in kits were included:

- AmpliTaq 360 (Applied Biosystems)
- FastStart Taq (Roche)
- Herculase II Fusion (Stratagene)
- HotMaster Taq (5 PRIME)
- HotStarTaq Master Mix Kit (QIAGEN)
- MyTaq (Bioline)
- MyTaq HS (Bioline)
- OneTaq (New England BioLabs)
- Pfu (Fermentas)
- Phusion High-Fidelity (New England BioLabs)
- PicoMaxx High Fidelity PCR System (Agilent Technologies)
- TaKaRa Ex Taq HS (TaKaRa Bio)
- TITANIUM Taq (Clontech)
- Tth (Roche)
- VeriQuest Probe qPCR Master Mix (Affymetrix)

These were evaluated on a dilution series of purified *Salmonella* ser. Typhimurium DNA analyzed by a validated standardized real-time PCR assay² using the accompanying PCR buffers for each polymerase.

The 5 best performing polymerases/kits were further evaluated using minced meat samples diluted in BPW 1:10 and enriched for 18 h at 37°C followed by artificial contaminated with *Salmonella* ser. Typhimurium. DNA extraction was performed on the samples followed by real-time PCR.

References

- Hedman et al (2012) Pre-PCR processing strategies. In: PCR Technologies, current innovations CRC Press, 3rd Ed.
- Löfström et al (2009) BMC Microbiol, 9(1):85.

