

Copyright is owned by the Author of the thesis. Permission is given for a copy to be downloaded by an individual for the purpose of research and private study only. The thesis may not be reproduced elsewhere without the permission of the Author.

IMMOBILIZATION OF LACTASE TO PERLOZA CELLULOSE RESINS

This thesis was presented in partial fulfilment of the requirements for the degree of Master of Science in Chemistry at Massey University

Min Yu

2000

AMENDMENTS

- p. ii, line 15 should read '**the** best of these results'
- p. iv, paragraph 3 – 'thank Dave Elgar for his **valuable** advice' and 'paragraph 5 – **heaps** of help?'
- p. xi 'NaOH sodium hydroxide'
- p. 2, line 1, imidazole, line 3 should be **sulphydryl**
- p. 4, Table 1.1, Streptococcus thermophilus
- p. 15, line 6, '**of** the commercial matrices available'
- p. 26, line 15 'ECH (0.75**ml**)'
- p. 28, line 28, 'the beads **were**'
- p. 30, line 24 '**rpm**' not rmp, also p. 32, line 24
- p. 31, line 14, '1 mM MgSO₄'
- p. 32, line 2 '**were** warmed to'
- p. 33, line 10, '**2.11**', not 1.11
- p. 33, line 17 should be acetonitrile and p. 34, line 2
- p. 37, line 6 '**Matsumoto**'
- p. 65 line 23, 'in other words'
- p. 72, line 8, change 'Binding Buffer **2 two** tested' to 'Binding Buffer 2 tested'
- p. 86, line 5, '**shifted**' not shift.
- p. 87, line 11, 'highest activity of immobilised lactase **being** achieved at low ECH activation levels'

ABSTRACT

A bead cellulose matrix, Perloza, was chemically modified by two attachment chemistries to prepare inexpensive resins for immobilization of lactase. A commercial product, the base-activated matrix Eupergit C was studied for comparison.

Three types of Perloza (Perloza 100 MT, Perloza 200 MT, Perloza 500 TM) were activated by epichlorohydrin (ECH) to achieve different activation levels. The best result for lactase immobilization was gained at low activation level (activated at 2% NaOH) for two attachment chemistries. The first attachment chemistry studied was that lactase immobilized directly to ECH activated Perloza. The second chemistry again used ECH activation and followed by attachment of the 6-amino caproic acid (ACA) spacer arm and then the lactase.

In the first chemistry, Perloza 100-ECH-Lactase obtained the highest activity 11.4 NLU/g (wet resin) over Perloza 200-ECH-Lactase and Perloza 500-ECH-Lactase (40 hours immobilization). In the second chemistry, Perloza 200-ECH-ACA-Lactase retained the highest activity 30.9 NLU/L (wet resin) over Perloza 100-ECH-ACA-Lactase and Perloza 500-ECH-ACA-Lactase. Overall the best results were obtained for the ECH-ACA resins. This best of these results showed about 3 times better immobilization than without ACA spacer arm.

The activity of immobilized lactase on Eupergit C obtained was 124~131.3 NLU/g (wet resin) for 24 hours immobilization. Although this result is about four times greater than Perloza, Perloza is a much cheaper matrix.

In the storage stability studies, both Perloza and Eupergit C immobilized lactase showed a sharp drop in activity initially within 1 day, then activity loss leveled out. Perloza 200-ECH-ACA-Lactase retained 82% of its original activity after 9 days storage. However, Eupergit-Lactase only retained 39% of its original activity after the same storage period.

This result indicated that Perloza 200-ECH-ACA-Lactase may possess much better storage stability than that of Eupergit-Lactase.

Studies on the inter-relationships between pH, temperature and Perloza immobilized lactase using the substrate (ONPG) indicated that maximum hydrolysis was attained at pH 6.5-7.2 and over a temperature range of 30-42°C. No shift in the pH and temperature optima in comparison to free enzyme was observed as a result of the process of immobilization of lactase on Perloza for both attachment chemistries.

The pH-activity curve of Eupergit-Lactase shifted towards more acidic pH values in the pH optimum in comparison to free lactase. The temperature optimum of Eupergit-Lactase shifted towards higher temperature compared to free lactase.

This study showed that Perloza has potential for the large scale use as a matrix of lactase immobilization.

ACKNOWLEDGEMENTS

First I would like to sincerely thank my supervisor Associate Professor David R.K. Harding for his time, advice, encouragement and especially his friendship to me throughout the two years study.

I would like to thank my supervisor Neill Haggarty for his support and chemicals supplement from NZDRI.

I would like to thank David Elgar for his valueless advice and help.

I would also like to thank Malcolm Chick for his advice on enzyme assays.

Another big thanks to the Centre for Separation Science, for those of you given me hips of help: Kate, Dick, Yishan, Claire, Satha, Rekha and Adrian.

Finally I would like to thank my family, friends for all their support.

TABLE OF CONTENTS

Abstract	ii
Acknowledgements	iv
List of Figures and Schemes	ix
List of Tables.....	x
List of Abbreviations.....	xi
CHAPTER ONE INTRODUCTION	1
1.1 The lactase problem	1
1.2 The enzyme lactase	1
1.3 Maxilact	6
1.4 Lactase application to food science.....	8
1.5 Methods for hydrolysis of lactose	10
1.6 Immobilization methods	11
1.6.1 Adsorption (and crosslinking)	11
1.6.2 Entrapment.....	12
1.6.3 Ultrafiltration membrane retainment	13
1.6.4 Covalent binding.....	13
1.6.5 Solid supports in covalent binding	15
1.7 Perloza beaded cellulose	16
1.8 Activation methods for covalent binding.....	17
1.9 Immobilized lactase systems.....	18
1.10 Activation methods studied.....	21
1.11 Analysis used in the study	22
1.12 Aim of this thesis	23
CHAPTER TWO MATERIAL AND METHODS	24
2.1 Reagents and equipment	24

2.2	Units used in this thesis.....	24
2.3	Perloza and Sepharose activation with epichlorohydrin (ECH).....	26
2.3.1	Activation method 1	26
2.3.2	Activation method 2.....	26
2.3.2	Activation method 3.....	26
2.4	Modification of Perloza 100MT with ECH and propylene oxide (PO).....	27
2.5	6-Aminocaproic acid (ACA) attachment to Perloza-ECH.....	28
2.6	Enzyme attachment	28
2.6.1	Enzyme attachment to Perloza-ECH, Sepharose-ECH and Eupergit C....	27
2.6.2	Enzyme attachment to Perloza-ECH-ACA.....	28
2.7	Titration methods	29
2.7.1	Epoxide titration.....	29
2.7.2	6-Aminocaproic acid (ACA) titration	30
2.8	Bicinchoninic acid (BCA) protein concentration determination.....	30
2.9	Lactase activity.....	31
2.9.1	Qualitative method	31
2.9.2	Quantitative method.....	31
2.9.3	Activity measurement of immobilized enzyme by "end point" method ...	32
2.10	Temperature adjustment in temperature profile study	33
2.11	pH adjustment in pH profile study	33
2.12	Purification of crude β -galactosidase (Maxilact L2000) by gel filtration.....	33

CHAPTER THREE LACTASE IMMOBILIZATION ON PERLOZA..... 35

3.1	Introduction	35
3.2	Results and discussion	39
3.2.1	Preliminary study	39
3.2.2	Further activation studies	40
3.2.3	Enzyme attachment	43
3.2.4	Lactase activity assay	54
3.2.5	Titration method.....	54

3.2.6	Determination of optimum pH for immobilized lactase	55
3.2.7	Determination of optimum temperature for immobilized lactase	58
3.2.8	Storage stability of the immobilized enzymes	60
3.2.9	Selection of storage buffer	62
3.3	Conclusions	63

CHAPTER FOUR LACTASE IMMOBILIZATION ON EUPERGIT C 65

4.1	Introduction	65
4.2.1	Structure and properties of Eupergit	65
4.1.2	Immobilization methods	66
4.2	Results and discussion	68
4.2.1	Titration of oxirane groups on Eupergit C	68
4.2.2	Optimum enzyme attachment conditions.....	68
4.2.2.1	Amount of Maxilact applied in immobilization.....	68
4.2.2.2	Optimum lactase immobilization time.....	71
4.2.2.3	Selection of pH value of the binding buffer	72
4.2.3	Lactase attachment onto Eupergit C.....	73
4.2.4	Determination of optimum pH for immobilized lactase on Eupergit C....	74
4.2.5	Determination of optimum temperature for immobilized lactase on Eupergit C	76
4.2.6	Storage stability study of immobilized lactase on Eupergit C	78
4.2.7	Comparison of stability of E-Lactase upon different immobilization times	79
4.3	Conclusions.....	81

CHAPTER FIVE SUMMARY AND FUTURE WORK..... 82

5.1	Summary	82
5.2	Future work	87

REFERENCES..... 88

LIST OF FIGURES

Figure 1.1a	Proposed mechanism of lactose hydrolysis by lactase	3
Figure 1.1b	Proposed mechanism of galactosyl transfer reaction lactase	3
Figure 3.1	The advantage of using a spacer arm	38
Figure 3.2	Activity of lactase on P-ECH and PECH-ACA	52
Figure 3.3	Lactase level on P-ECH and PECH-ACA.....	53
Figure 3.4	pH profile of free and immobilized lactase at 30°C	57
Figure 3.5	Temperature profile of free and immobilized lactase at pH 6.5.....	59
Figure 3.6	Stability study of immobilized lactase on Perloza 100	61
Figure 3.7	Stability study of immobilized lactase on Perloza 200	61
Figure 3.8	Stability study of immobilized lactase stored in different buffers	62
Figure 4.1	Preparation of Eupergit	66
Figure 4.2	Relative activity of wet catalyst on Eupergit C vs. Maxilact applied	70
Figure 4.3	Activity yield of wet catalyst on Eupergit C vs. Maxilact applied.....	70
Figure 4.4	Immobilization of Lactase/time dependency	72
Figure 4.5	pH profile of free and Eupergit C immobilized lactase at 30°C.....	75
Figure 4.6	Temperature profile of free and Eupergit C immobilized lactase at pH 6.5	77
Figure 4.7	Stability study of immobilized lactase on Eupergit C.....	78
Figure 4.8a	Activities of Eupergit-Lactase upon 24 and 96 hours immobilization vs. storage time	80
Figure 4.8b	Stability study of Eupergit-Lactase upon 24 and 96 hours immobilization	80
Figure 5.1	Comparison of storage stability of immobilize lactase on Perloza and Eupergit C	83
Figure 5.2	pH profile of free, Perloza 200-ECH-ACA and Eupergit C immobilized lactase.....	85
Figure 5.3	Temperature profile of free, Perloza 200-ECH-ACA and Eupergit C immobilized lactase.....	85

LIST OF TABLES AND SCHEMES

Table 1.1	Possible sources of lactase	4
Table 1.2	Commercial preparation of lactase.....	5
Table 1.3	Use of lactase: new product development.....	8
Table 1.4	Amino acids involved in covalent binding.....	15
Table 1.5	Selection of activation chemistries.....	18
Table 2.1	Conversions of "g dry" weight to "g" wet weight of the various resins	25
Table 3.1	Epoxide activation with different NaOH concentration.....	41
Table 3.2	Lactase attachment upon different activation levels	42
Table 3.3	Activation of three types of Perloza at NaOH concentration 2% and 3%.....	44
Scheme 3.1	Activation of Perloza-ECH-ACA and lactase attachment	47
Table 3.4	Lactase attachment upon Perloza-ECH and Perloza-ECH-ACA	51
Scheme 3.2	Titration of ACA group on Perloza-ECH-ACA.....	55
Table 4.1	Lactase immobilization with series amounts of Maxilact applied	69
Table 4.2	Activity of immobilized lactase on Eupergit C/dependency on pH of binding	73

LIST OF ABBREVIATIONS

ACA	6-aminocaproic acid
ADH	alcohol dehydrogenase
BCA	bicinchoninic acid
CDI	1,1'-carbonyldiimidazole
CMC	1-cyclohexyl-3-(2-morpholino-ethyl) carbodiimide
ECH	epichlorohydrin
EDC	1-ethyl-3-(3-dimethylaminopropyl) carbodiimide
FPLC	fast protein liquid chromatography
GRAS	generally recognized as safe
HCl	hydrochloric acid
HL	hydrolyzed lactose
K-phos	potassium phosphate buffer
NaOH	sodium hydroxyl
NaCl	sodium chloride
NHS	N-hydroxysuccinimide
ONP	o-nitrophenyl
ONPG	o-nitrophenyl- β -D-galactopyranoside
UF	ultrafiltration