



IPB  
INSTITUTE FOR BIOTECHNOLOGY AND BIOENGINEERING

ipb  
INSTITUTO POLITÉCNICO DE BRAGANÇA

LABORATÓRIO ASSOCIADO  
LSRE  
ICM

# for cationic resins in the sodium and potassium form

C. Nobre<sup>1</sup>, M. J. Santos<sup>1</sup>, A. Dominguez<sup>1</sup>, D. Torres<sup>1,2</sup>, A. M. Peres<sup>3</sup>, I. Rocha<sup>1,2</sup>, E.C. Ferreira<sup>1</sup>, L.R. Rodrigues<sup>1,2</sup>, J. A. Teixeira<sup>1</sup>

<sup>1</sup>IBB – Institute for Biotechnology and Bioengineering, Centre of Biological Engineering, University of Minho, Campus de Gualtar 4710-057 Braga Portugal

<sup>2</sup>Biotempo, Lda., AVEPARK, Zona Industrial da Gandra, Apartado 4152, 4806 - 909 Caldas das Taipas, Portugal

<sup>3</sup>LSRE – Laboratory of Separation and Reaction Engineering, Escola Superior Agrária – Instituto Politécnico de Bragança, Quinta de Santa Apolónia, Apartado 172, 5301-855 Bragança, Portugal

clarissenobre@deb.uminho.pt

## INTRODUCTION


- Separation of glucose from mixtures of fructose and sucrose in molasses is a major challenge in industrial sugar chromatographic separations.
- Sulfonated poly(styrene-co-divinylbenzene) (PS-DVB) ion exchange resins are the most frequently used for sugars separation, generally in a cationic form. The cation of the resins used in this study complex with the hydroxyl group of the sugar adsorbed leading to a selective adsorption. Thus, the conformation of the sugar determines its relative affinity for the resin and its distribution coefficient.
- The separation process is usually carried out at high temperatures. However, this implies high energy costs and high levels of hydrolysis.
- The knowledge of the adsorption isotherms of the sugars in a mixture is a very important parameter for the selection of the adsorbent to be used in the chromatographic separation.

## OBJECTIVE

- In this study, the adsorption isotherms of glucose (G), fructose (F) and sucrose (S) in mixtures of mono and multi-component at 25 and 40°C were determined for two resins of PS-DVB in the sodium and potassium forms.

## EXPERIMENTAL METHODOLOGY

### Static Method



- Sugars were quantified by High Performance Liquid Chromatography (HPLC)

- Equilibrium loading  $q$  (mg\_sugar/g\_adsorbent)

$$q = \frac{(C_i - C_f)V}{m}$$

$C_i$  - Initial concentration (mg/mL)

$C_f$  - Final concentration in bulk (mg/mL)

$V$  - Volume of sugar solution (mL)

$m$  - dry adsorbent mass (g)

- Adsorption Isotherm Model

$$q = k, C_f$$

$k$  - distribution coefficient (mL/B<sub>ads</sub>)


- Selectivity ( $\alpha$ ):

$$\alpha = \frac{k_F}{k_G} = \frac{k_G}{k_S} = \frac{k_F}{k_S}$$

- For testing the equality of the distribution coefficient a factorial design was set up as an analysis of covariance at a 5% significance level. SPSS version 17 was used.

## RESULTS AND DISCUSSION

### Adsorption Isotherms



**Table 1.** Variation of distribution coefficient ( $k$ ) from F, G, and S with resin (K/Na), component (mono/multi) and temperature (25/40°C)

	K resin				Na resin				25°C				40°C			
	25°C	40°C	Mono	Multi	25°C	40°C	Mono	Multi	Mono	Multi	Mono	Multi	Mono	Multi	Mono	Multi
F	→Multi	→Multi	25>40	25>40	→Multi	→Multi	25>40	25>40	→Multi	→Multi	25>40	25>40	→Multi	→Multi	25>40	25>40
G																
S																

#### Component

- The presence of additional sugars (multi-component) decreased the distribution coefficient, although, for the potassium resin at 40°C a synergetic effect was found.

#### Temperature


- In the mono-component mixtures the increase of temperature lead to a decrease of the distribution coefficient, therefore, the capacity of both resins decreased. Nevertheless, the capacity in multi-component mixtures was not significantly affected by the temperature in the sodium resin, and for the potassium resin a small increase in the adsorption was observed.

#### Selectivity

- The sugars that demonstrated lower selectivity values were F/G since these sugars have the same molecular weight. S/F selectivity showed a higher value as compared to the S/G one because of the complex that is formed between the fructose and the resin.
- For both resins, selectivity decreased with the increase of the temperature. Nevertheless, the selectivity for the potassium resin was higher than the selectivity of the sodium resin.

**Table 2.** Selectivity ( $\alpha$ ) between sugars:

Resin	K		Na	
T (°C)	25	40	25	40
F/G	1,28	1,18	1,07	1,03
G/S	1,56	1,49	1,04	1,10
F/S	2,00	1,76	1,11	1,13



## CONCLUSIONS

- In our case study, the potassium resin appears to be the most suitable adsorbent for an industrial purpose, being the best results achieved for an operational temperature of 25 °C.

## ACKNOWLEDGEMENTS

Agência de Inovação (AdI) – Project BIOLIFE ref. PRIME 03/347.  
The authors acknowledge to FCT - Fundação para a Ciência e a Tecnologia for the PhD grant ref. SFRH/BD/32514/2006