



International Conference on Food Contaminants 2015
Challenges in chemical mixtures

Program and Abstract Book

ICFC 2015

The logo is contained within a white, irregular hexagonal shape that is part of a larger graphic of interconnected white hexagons on a dark blue background. The logo itself consists of the text 'ICFC 2015' in blue, with a circular emblem to the right containing three icons (sun, leaf, fruit).

International Conference on Food Contaminants
Challenges in chemical mixtures

APRIL
13th and 14th

Lisbon,
Portugal

(P11) Extraction and detection of mycotoxins in medicinal and aromatic plants: a case study with *Melissa officinalis* L.

Eliana Pereira, Isabel C.F.R. Ferreira, Paula Rodrigues*

Centro de Investigação de Montanha (CIMO), ESA, Instituto Politécnico de Bragança, Portugal *(prodrigues@ipb.pt).

Plants frequently suffer contaminations by toxigenic fungi, and their mycotoxins can be produced throughout growth, harvest, drying and storage periods. The objective of this work was to validate a method for detection of toxins in medicinal and aromatic plants, through a fast and highly sensitive method, optimizing the joint co-extraction of aflatoxins (AF: AFB1, AFB2, AFG1 and AFG2) and ochratoxin A (OTA) by using *Aloysia citrodora* P. (lemon verbena) as a case study. For optimization purposes, samples were spiked (n=3) with standard solutions of a mix of the four AFs and OTA at 10 ng/g for AFB1, AFG1 and OTA, and at 6 ng/g of AFB2 and AFG2. Several extraction procedures were tested: i) ultrasound-assisted extraction in sodium chloride and methanol/water (80:20, v/v) [(OTA+AFs)1]; ii) maceration in methanol/1% NaHCO₃ (70:30, v/v) [(OTA+AFs)2]; iii) maceration in methanol/1% NaHCO₃ (70:30, v/v) (OTA1); and iv) maceration in sodium chloride and methanol/water (80:20, v/v) (AF1). AF and OTA were purified using the mycotoxin-specific immunoaffinity columns AflaTest WB and OchraTest WB (VICAM), respectively. Separation was performed with a Merck Chromolith Performance C18 column (100 x 4.6 mm) by reverse-phase HPLC coupled to a fluorescence detector (FLD) and a photochemical derivatization system (for AF). The recoveries obtained from the spiked samples showed that the single-extraction methods (OTA1 and AF1) performed better than co-extraction methods. For in-house validation of the selected methods OTA1 and AF1, recovery and precision were determined (n=6). The recovery of OTA for method OTA1 was 81%, and intermediate precision (RSD_{int}) was 1.1%. The recoveries of AFB1, AFB2, AFG1 and AFG2 ranged from 64% to 110% for method AF1, with RSD_{int} lower than 5%. Methods OTA1 and AF1 showed precision and recoveries within the legislated values and were


found to be suitable for the extraction of OTA and AF for the matrix under study.

Acknowledgments:

The authors are grateful to PRODER nº 53514, AROMAP, for financial support of the work and E. Pereira grant, to Mais Ervas Lda. for providing samples, and to FCT for financial support to CIMO (PEst-OE/AGR/UI0690/2014).

ICFC 2015

International Conference on Food Contaminants
Challenges in chemical mixtures



Certificate of Presentation

Paula Rodrigues

Presented a Poster entitled

“EXTRACTION AND DETECTION OF MYCOTOXINS IN MEDICINAL AND AROMATIC PLANTS: A
CASE STUDY WITH *Melissa officinalis* L.”

at the International Conference on Food Contaminants (ICFC 2015) held in Lisbon on the 13th
and 14th of April, INFARMED auditorium, Lisbon, Portugal.

Fernando de Almeida

Chairman of the Executive Board of the
National Institute of Health Doutor Ricardo Jorge I.P
(signature and seal)

