



# Visual Grading Analysis of image quality in pediatric abdominal images acquired by Direct Digital Radiography and Computer Radiography Systems

Poster No.:	C-2098
Congress:	ECR 2013
Туре:	Scientific Exhibit
Authors:	A. C. Oliveira <sup>1</sup> , A. P. Martins <sup>2</sup> , R. T. Avelãs <sup>3</sup> , M. S. C. D. Santos <sup>1</sup> , P. M. Martins <sup>1</sup> , S. De Francesco <sup>1</sup> , P. Sá-Couto <sup>1</sup> , <u>C. Ferreira</u> <sup>4</sup> ; <sup>1</sup> Aveiro/PT, <sup>2</sup> Viana do Castelo/PT, <sup>3</sup> Vermelha/PT, <sup>4</sup> Coimbra/PT
Keywords:	Dosimetric comparison, Comparative studies, Digital radiography, Radioprotection / Radiation dose, Pediatric, Abdomen
DOI:	10.1594/ecr2013/C-2098

Any information contained in this pdf file is automatically generated from digital material submitted to EPOS by third parties in the form of scientific presentations. References to any names, marks, products, or services of third parties or hypertext links to third-party sites or information are provided solely as a convenience to you and do not in any way constitute or imply ECR's endorsement, sponsorship or recommendation of the third party, information, product or service. ECR is not responsible for the content of these pages and does not make any representations regarding the content or accuracy of material in this file.

As per copyright regulations, any unauthorised use of the material or parts thereof as well as commercial reproduction or multiple distribution by any traditional or electronically based reproduction/publication method ist strictly prohibited.

You agree to defend, indemnify, and hold ECR harmless from and against any and all claims, damages, costs, and expenses, including attorneys' fees, arising from or related to your use of these pages.

Please note: Links to movies, ppt slideshows and any other multimedia files are not available in the pdf version of presentations.

www.myESR.org

Page 1 of 24

# Purpose

The advent of digital technology allowed for great improvements in radiology and lead the way for digital radiology, leaving behind conventional x-ray techniques, [1]. Digital post-processing of image is the main advantage of digital image systems (e.g., computed radiology and direct digital radiology) over the conventional systems. Image quality can indeed be improved avoiding the increase of patient dose and the number of unnecessary exposures [1-3].

Image quality is directly linked to the dose of radiation applied to the patient. The literature points out the need for appropriate image evaluation in order to reduce the patient dose. Optimization and practice justification are of great relevance in diagnostic radiology. In pediatric patients a good practice justification is even more important as radiation exposure on the first 10 years of life increases the possibility of negative effects on the patient's health compared to an exposure between the ages of 30 to 40 years old [4, 5].

In the past, the main concern of radiology technicians was image quality, leaving radiation dose to a second plan. Nowadays, dose reduction and the cost/benefit ratio represent the main concerns in radiology. The radiographer is responsible for applying the ALARA (As Low As Reasonable Acceptable) concept in every study involving the use of ionizing radiation [4, 6]. According to the ALARA concept, the necessary level of diagnostic image quality should be attained with the lowest patient dose possible.

In Pediatric radiology, besides the ALARA principle, the International Commission on Radiological Protection (ICRP) adds a new concept - the SMART message, also related with optimization and radiation protection when applied to pediatric radiology (Figure 1) [6].



Figure 1 - SMART message from ICRP [6]

#### Fig. 1

*References:* Figure adapted from: Alexakhin RM, Menzel H. The 2007 Recommendations of the International Commission on Radiological Protection, 2007.

There are several methods to evaluate image quality: objective (phantom based) or subjective. In the subjective methods, the evaluation is based on the observer's perception guided by some kind of image quality criteria that must be accomplished. The European Commission published a set of Guidelines on Quality Criteria in order to suggest standard protocols for a broad range of radiographic procedures responding to the diagnostic needs with acceptable dose, and establish image quality criteria that should be accomplished in the evaluation of the obtained images.

In this context, the main goal of the present study is to evaluate and compare the quality of radiographic abdominal pediatric images, acquired by direct digital radiography (DR) and computed radiography (CR) systems.

Page 3 of 24

Image quality evaluation was performed through a subjective method known as Visual Grading Analysis (VGA) [5].

#### Images for this section:



Figure 1 - SMART message from ICRP [6]

Fig. 1

Page 4 of 24

# **Methods and Materials**

#### Sample description

In this study, a sample of 44 images of pediatric abdomen was used, from children aged between 5 to 10 years. The data was collected in two institutions, referred as Institution A and B. Institution A was equipped with a DR system and Institution B with a CR system. Imaging protocols were slightly different with respect to exposure factors (Table 1).



Figure 2 - Schematic view of the study



Page 5 of 24

	Parameters	European Guidelines (Conventional Radiography)	Institution A Protocol (DR System)	Institution B Protocol (CR System)
	Patient Position	Supine or prone	Orthostatism.	Supine or orthostatism
	Radiographic device	Table, grid table	Grid table.	Table, grid table
	Nominal focal spot value	0,6 (≤ 1,3)	Thick	Thick
	Additional filtration	Up to 1 mm Al + 0,1 or 0,2 mm Cu (or equivalent)	2mm Al	2mm Al
Protocols	Anti-scatter grid	No grid for infants < 6 months of age	Grid	<ul> <li>Supine: no grid</li> <li>Orthostatism: grid</li> </ul>
	FFD	100 - 115 cm	100 cm	100 cm
	kV	65 - 85	70	67
	mAs		12,5	20
	Automatic exposure control	Chamber selected - central or both lateral	Chamber selected- central or both lateral	Chamber selected - central
	Protective shielding	Gonad capsules formale patients. Lead-rubber coverage of the thorax in the immediate proximity of the beam edge.		-

Table 1 - Protocols used in institutions participating in the study

# Table 1References: - Aveiro/PTImage quality evaluation

Image quality evaluation was performed through a subjective method known as VGA. This method can be used to assess, in a scientific way, whether or not the images comply with the established criteria [4, 7]. The criteria used in this study is in line with the one established in the European Guidelines on Quality Criteria for Diagnostic Radiographic Images in Paediatrics (EGQCDRIP), a document based on the work carried out by several renowned entities and different professionals. For more details please refer to [4].

The VGA method has two variants: the absolute, where the scores are given on an absolute scale, and the relative, where the image quality ought to be compared with a reference image. In the relative VGA, usually a five level scale (ranging from "much worse" to "much better") is used to assess whether or not the image is better than the reference image (Table 2) [5, 8, 9]. The VGA score can be calculated using the following equation (appearing as figure 8):

Page 6 of 24



Fig. 8

*References:* Tingberg A. Evaluation Of Image Quality Of Lumbar Spine Images: A Comparison Between FFE And VGA. Radiation protection dosimetry. 2005;144:53-61.

Where  $G_{i,s,o}$  is the grading (-2, -1, 0, +1, +2) for image *i*, structure *s* and observer *o*. *I* is the number of images; *S* the number of structures; *O* the number of observers [5, 9].

Evaluation	Viewed Structure
-2	Clearly lower than the reference image
-1	Slightly lower than the reference image
0	Equal to the reference image
1	Slightly better than the reference image
2	Clearly better than the reference image

Table 2 - Classification system used in VGA [5, 10]

#### Table 2 *References:* - Aveiro/PT

The images were evaluated by four observ

The images were evaluated by four observers (radiologists), with a professional experience ranging from 5 to 29 years. This assessment was based on a questionnaire grounded on the European Guidelines [4] for this specific examination (Table 3, Figure 3).

Page 7 of 24

AP/PA PEDIAT	RIC ABDOMINAL PROJECTION WITH VERTICAL/HORIZONTAL BEAM
	<ul> <li>Visualization of the abdomen, from the diaphragm to the ischial tuberosities, including the lateral abdominal walls.</li> <li>Visualization of the properitoneal fat lines consistent with age.</li> </ul>
Image criteria	<ul> <li>Visualization of the kidney outlines consistent with age and depending on bowel content.</li> <li>Visualization of the psoas outline consistent with age and depending</li> </ul>
	<ul><li>on bowel content.</li><li>Visually sharp reproduction of the bones.</li></ul>

Table 3 - Image Criteria of AP/PA Pediatric Abdominal Projection [4]

#### Table 3 *References:* - Aveiro/PT



#### Fig. 3 References: - Aveiro/PT

Since relative VGA was used in this study, an image fulfilling all the European Guidelines (EG) requirements was selected from the 44 images of the sample and was considered the reference image (Figure 4). The questionnaire was completed in accordance with a five level scale as shown in (Table 2).

Page 8 of 24



# Figure 4 - Reference image used in VGA evaluation

#### Fig. 4 *References:* - Aveiro/PT

Aiming for a more thorough description of the sample used in this study, besides the specific evaluation criteria for abdominal radiography, considered in the VGA method, a set of three items was added: positioning, centering and collimation. These items were classified as YES / NO by the observers.

Page 9 of 24

All evaluations were performed independently on the same physical space, with reduced luminance and high-resolution display systems. It is noteworthy that the observers were not allowed to manipulate the images (display window) to maintain the original contrast. They were blinded to the detector system used to acquire the images.

The Entrance Surface Dose (ESD) value, for each of the systems where the images were obtained, was estimated using MICADO software (Figure 5).

Calcul de la dos	e à la	surfac	ce d'entrée (De) à pa	rtir de					
Paramètres d'exposition			Produit Dose.Surface (PDS)						
Type d'examen	(à cho	isir da	ans l'une des deux li	stes déroulantes) aid					
Adulte : Veuillez choisir	Ċ.		ou Pédiatrie : Abdomer	n sans préparation 30 kg 💌					
			Juniant 22						
Saisie des autre	s para	mètre	es techniques	Carlité Francisa da Padiotado					
Saisie des autre Valeurs recommandées (http://www.sfmet.org)	s para dans le p	mètre uide des	es techniques procédures radiologiques de la	Société Française de Radiologie					
Saisie des autre 'Valeurs recommandées (http://www.sfrnet.org) Haute tension :	s para dans le g	mètre uide des	procédures radiologiques de la	Société Française de Radiologie					
Saisie des autre <sup>•</sup> Valeurs recommandées (http://www.sfrnet.org) Haute tension : Filtration inhérente :	dans le g	mètre uide des	s techniques procédures radiologiques de la 1 kY * [ 60.0 - 80.0 ]	Société Française de Radiologie ommandée <b>aide</b>					
Saisie des autre Valeurs recommandées (http://www.sfmet.org) Haute tension : Filtration inhérente : Filtration additionnelle :	s para dans le g 70 2,5 mm 2 mm Al	mètre uide des Al	kY * [ 60.0 - 80.0 ]	Société Française de Radiologie ommandée <b>aide</b>					
Saisie des autre Valeurs recommandées http://www.sfrnet.org) Haute tension : Filtration inhérente : Filtration additionnelle : Dist. Foyer - Peau :	70 2,5 mm 2 mm Al 100	mètre uide des Al	kV * [ 60.0 - 80.0 ] * Filtration additionnelle reco adde Dist. Foyer - Detecteur :	Société Française de Radiologie ommandée <b>aide</b> cm * ( 100.0 - 115.0 )					
Saisie des autre <sup>•</sup> Valeurs recommandées (http://www.sfrnet.org) Haute tension : Filtration inhérente : Filtration additionnelle : Dist. Foyer - Peau :	s para dans le g 70 2,5 mm 2 2 mm Al 100	mètre uide des Al	stechniques procédures radiologiques de la kV * [ 60.0 + 80.0 ] * Filtration additionnelle reco aide Dist. Foyer - Detecteur : Epaisseur patient :	Société Française de Radiologie ommandée aide cm * ( 100.0 - 115.0 ) cm aide					

Figure 5 - MICADO software interface screenshot [Available on: Institut de Radioprotection et de Sûreté Nucléaire (IRSN), website: http://nrd.irsn.fr/index.php?page=micado]

#### Fig. 5

*References:* Institut de Radioprotection et de Sûreté Nucléaire (IRSN), website: http:// nrd.irsn.fr/index.php?page=micado

#### Statistical Analysis

Statistical analysis was performed using the Statistical Package for Social Sciences (SPSS) version 19.0.

The images which at least 75% of the observers classified as having incorrect positioning were excluded from the analyses as they could introduce bias. Because of this, the sample was reduced to 39 images (20 from DR system and 19 from CR system).

Page 10 of 24

Interclass Correlation Coefficient (ICC) was used to evaluate data reliability through agreement among the four observers. The observers were considered unique and the images represent a random sample.

Student-t test was used to determine differences in the means of the two independent quantitative variables since the assumptions of normality were validated by the Kolmogorov-Smirnov test (KS).

A significance level of 0,05 (# = 0,05) was considered for all the tests performed. The final classification of each image was obtained through the average ratings given by the four observers.

#### Images for this section:



Figure 2 - Schematic view of the study

Fig. 2

Page 11 of 24

	Parameters	European Guidelines (Conventional Radiography)	Institution A Protocol (DR System)	Institution B Protocol (CR System)
	Patient Position	Supine or prone	Orthostatism.	Supine or orthostatism
	Radiographic device	Table, grid table	Grid table.	Table, grid table
	Nominal focal spot value	0,6 (≤ 1,3)	Thick	Thick
	Additional filtration	Up to 1 mm Al + 0,1 or 0,2 mm Cu (or equivalent)	2mm Al	2mm Al
Protocols	Anti-scatter grid	No grid for infants < 6 months of age	Grid	Supine: no grid     Orthostatism: grid
	FFD	100 - 115 cm	100 cm	100 cm
	kV	65 - 85	70	67
	mAs		12,5	20
	Automatic exposure control	Chamber selected - central or both lateral	Chamber selected - central or both lateral	Chamber selected - central
	Protective shielding	Gonad capsules formale patients. Lead-rubber coverage of the thorax in the immediate proximity of the beam edge.		-

Table 1 - Protocols used in institutions participating in the study

#### Table 1

Evaluation	Viewed Structure
-2	Clearly lower than the reference image
-1	Slightly lower than the reference image
0	Equal to the reference image
1	Slightly better than the reference image
2	Clearly better than the reference image

Table 2 - Classification system used in VGA [5, 10]

#### Table 2

AP/PA PEDIATI	RIC ABDOMINAL PROJECTION WITH VERTICAL/HORIZONTAL BEAM
	<ul> <li>Visualization of the abdomen, from the diaphragm to the ischial tuberosities, including the lateral abdominal walls.</li> <li>Visualization of the properitoneal fat lines consistent with age.</li> </ul>
Image criteria	<ul> <li>Visualization of the kidney outlines consistent with age and depending on bowel content.</li> </ul>
	<ul> <li>Visualization of the psoas outline consistent with age and depending on bowel content.</li> <li>Visually sharp reproduction of the bones.</li> </ul>

Table 3 - Image Criteria of AP/PA Pediatric Abdominal Projection [4]

#### Table 3



Fig. 3

Page 13 of 24





Fig. 4

Page 14 of 24

Calcul de la dos	e à la sur	face	d'entrée (De) à	a part	tir de	
Paramètres d'exposition		Produit Dose-Surface (PDS)				
Type d'examen	(à choisin	r dan	s l'une des deu	x list	es déroulantes) a	
Adulte : Veuillez choisin			ou Pédiatrie : Abd	domen s	ans préparation 30 kg 💌	
Saisie des autre	s paramè	tres	techniques	se la Soc	ciété Francaise de Radiolo	
Saisie des autre <sup>*</sup> Valeurs recommandées (http://www.sfrnet.org)	s paramè dans le guide	des pro	sulvant >> techniques occidures radiologiques d	de La Sor	ciété Française de Radiolog	
Saisie des autre <sup>•</sup> Valeurs recommandées (http://www.sfmet.org) Haute tension :	s paramè dans le guide 70	des pro	Sulvant >>> techniques océdures radiologiques d	de la Soc	ciété Française de Radiolog	
Saisie des autre <sup>•</sup> Valeurs recommandées (http://www.sfmet.org) Haute tension : Filtration inbérente :	s paramè dans le guide 70 2,5 mm Al	des pro	Sulvant >> techniques coédures radiologiques d ky *[60.0 - 80.0] * Filtration additionnelle	še la Sor e recom	ciété Française de Radiolog mandée <b>aide</b>	
Saisie des autre <sup>•</sup> Valeurs recommandées (http://www.sfrnet.org) Haute tension : Filtration inbérente : Filtration additionnelle :	s paramè dans le guide 70 2,5 mm Al 2 mm Al	des pro	Sulvant >> techniques cocédures radiologiques d kV * [ 60.0 - 80.0 ] * Filtration additionnelle aide	de La Sor e recom	ciété Française de Radiolo mandée <b>aide</b>	
Saisie des autre * Valeurs recommandées (http://www.sfrnet.org) Haute tension : Filtration inhérente : Filtration additionnelle : Dist. Foyer - Peau :	s paramè dans le guide 70 2,5 mm Al 2 mm Al 100 cm	des pro	Sulvant >> techniques cooldures radiologiques d kV * [ 60.0 - 80.0 ] * Filtration additionnelle aide Dist. Foyer - Detecteur :	de la Sor e recom	ciété Française de Radiolog mandée <b>aide</b> ] <b>cm</b> * ( 100.0 - 115.0 )	
Saisie des autre * Valeurs recommandées (http://www.sfrnet.org) Haute tension : Filtration inhérente : Filtration additionnelle : Dist. Foyer - Peau :	s paramè dans le guide 70 2,5 mm Al 2 mm Al 100 cm	des pro	Sulvant >> techniques cocidures radiologiques d kV *[60.0 - 80.0] *Filtration additionnelle aide Dist. Foyer - Detecteur : Epaisseur patient :	de la Soc e recom	clété Française de Radiolog mandée aide cm * ( 100.0 - 115.0 ) cm aide	

Figure 5 - MICADO software interface screenshot [Available on: Institut de Radioprotection et de Sûreté Nucléaire (IRSN), website: http://nrd.irsn.fr/index.php?page=micado]

Fig. 5

# Results

The results of the patient positioning evaluation, collimation and centering point are presented on Figure 6.



Figure 6 - Results from YES/NO analysis regarding positioning, collimation and centering

#### Fig. 6 References: - Aveiro/PT

The reliability study revealed an ICC value of 0,775, and that the 95% confidence interval varied between 0,633 and 0,872. As the ICC is relatively overhead, and the confidence interval has moderate amplitude, we can infer that the internal consistency degree is good. This result is statistically significant because the *p*-value obtained in the ICC (*p*<0,01) is less than 0,05 for a 95% confidence interval. Therefore we concluded that there is agreement amongst the observers. The exclusion of any of the observers did not produce an internal consistency improvement.

Based on the Levene's test for #=0,05, and because the p-value<sub>(Levene)</sub>= 0,369 is higher than #, we can conclude that the variances are similar. Consequently the *p*-value obtained in the t-Student test (*p*<0,01) is less than 0,05 which means that the mean scores obtained in the CR and DR systems are statistically different.

Page 16 of 24

			Inde	epende	nt Simp	oles Test				
		Leven for Eq Vari	e's Test uality of ances	s Test ality of nces t-test for Equality of Means						
						Sig. (2-	Mean Differenc	Std. Error Differen	95% Confide Interval Differe	of the nce Uppe
		F	Sig.	t	df	tailed)	e	ce	Lower	r
Score to tal de cada	Equal variances assumed	,827	,369	4,031	37	,000	,544053	,13495 7	,27060 3	,817 502
im a gem	Equal variances no t assum ed			4,057	35,67 2	,000	,544053	,13409 2	,27201 4	,816 091

Table 4 - Independent Samples Test

#### Table 4 References: - Aveiro/PT

The mean score found for the CR system (Mean = -0,217 and SD = 0,365) is less than the mean score found for the DR system (Mean = 0,327 and SD = 0,468), meaning that, according to the observers, the images obtained with the DR system presented a superior image quality comparing with those from the CR system.

Radiographic System	n	Mean	S.D.	Median	Minimum	Maximum
CR System	19	-0,217	0,365	-0,375	-0,667	0,542
DR System	20	0,327	0,468	0,354	-0,542	0,958

Table 5 - Statistical results

#### Table 5 *References:* - Aveiro/PT

Looking at Table 5 and the box-and-whisker plot (Figure 7) we can conclude that there are no outliers, the DR system median is higher than the one from the CR system and the variability on the CR system is lower than the one found on the results from the DR system. Both box-and-whisker plot presented symmetry.

Page 17 of 24





#### Fig. 7 References: - Aveiro/PT

Concerning the estimated ESD values we found 517  $\mu$ Gy for the DR system and 740  $\mu$ Gy for the CR System; both within the reference doses for abdominal pediatric examinations for ages between 5 and 10 years old defined by ICRP (500-800  $\mu$ Gy) [6]. These results indicate that the DR system presents a dose value relatively lower than the one from the CR System. The difference between values is around 200  $\mu$ Gy.

#### Images for this section:

Page 18 of 24





Figure 6 - Results from YES/NO analysis regarding positioning, collimation and centering

		_	Inde	epende	nt Sim	oles Test				
		Leven for Eq Vari	Levene's Test for Equality of Variances			t-test fi	or Equality	ofMeans		
							Mean	Std. Error	95% Confidenœ Interval of the Difference	
		F	Sia.	t	df	Sig. (2- tailed)	Difference	Di ffer en ce	Lower	Uppe
Score total de cada	Equal variances assum ed	,827	,369	4,031	37	,000	,544053	,13495 7	,27060 3	,817 502
im a gem	Equal variances no t assum ed			4,057	35,67 2	,000	,544053	,13409 2	,27201 4	,816 091

Table 4 - Independent Samples Test

Table 4

Fig. 6

Page 19 of 24

Radiographic System	n	Mean	S.D.	Median	Minimum	Maximum
CR System	19	-0,217	0,365	-0,375	-0,667	0,542
DR System	20	0,327	0,468	0,354	-0,542	0,958

Table 5 - Statistical results

#### Table 5





Fig. 7

Page 20 of 24

# Conclusion

In the present study, the subjective image quality, as provided by VGA scores on pediatric abdominal images, was significantly higher when using the DR system compared to the one obtained using the CR system. Moreover, based on the calculated ESD results, we can suggest that DR systems allowed a 30% dose reduction when compared with CR systems. These results indicate, that patient dose reduction can be achieved by using DR system without loss of the required diagnostic image quality [11]. This might be due to the higher quantum detection efficiency provided by DR systems and it is in line with the results reported for chest studies with a similar approach [12, 13].

It is important to notice that ESD results obtained by both systems are within the reference doses defined by the ICRP for patients aged between 5 to 10 years old, such as the patients participating in this study, although these reference values are defined for screen-film systems and not digital systems [6]. Dose reference levels for digital systems are not currently available preventing a more thorough comparison. Dose reduction can be attained if optimized protocols, devoted for digital systems, were implemented.

Although not directly related to the system where the images are acquired, general principles such as correct positioning, beam limitation and centering have a great impact on image quality and most importantly on patient dose. Unnecessary repetitions or poor image quality can be avoided if these rules are applied [4]. In this study, according to the observers, proper positioning and centering was generally achieved (80-100%) in both institutions. However, beam collimation was rated, at least by two observers, as inadequate in about 50% of the images. Beam collimation is relevant, mainly in pediatric patients due to dose concerns. Furthermore, if the field-size is too large an increment in scattering radiation leads to a decrease in contrast with a negative impact on image quality. Following best practices is of utmost importance to reduce patient doses and this matter should deserve great attention from the institutions and their professionals.

Although regulation on radiation protection is already incorporated in national legislative documents, in Portugal a real implementation of these recommendations in daily clinical practice is far from being achieved. Furthermore, typical distributions of exposure parameters in plain radiography are unknown and there are few studies about reduction/ optimization exposure doses, which become more relevant in pediatric patients [1].

The present study represents a small but important contribute to this topic since it compares image quality in CR and DR systems establishing a relationship with patient dose, in pediatric patients. Even with doses within acceptable limits, the potential for dose reduction was noted. Results should be analyzed and shared with the participating

Page 21 of 24

institutions to develop adequate imaging protocols. In the future, we recommend that other anatomical areas and direct measures of ESD are considered using a larger sample.

# References

1. Lança LJOC. A imagem radiológica em sistemas de radiografia digital: estudo do efeito dos parâmetros técnicos de exposição na qualidade diagnóstica e na optimização de dose. Tese de Doutoramento. Universidade de Aveiro. 2011.

2. Samei E, Seibert J, Andriole K, Badano A, Crawford J, Reiner B, et al. AAPM/ RSNA tutorial on equipment selection: PACS equipment overview. Radiographics. 2004;24:313-32.

3. Garmer M, Hennings S, Jager H, Schrick FDLT, Jacobs A, Hanusch A, et al. Digital radiography versus conventional radiography in chest imaging: diagnostic performance of a large-area silicon flat-panel detector in a clinical CT-controlled study. . Am J Roentgenol. 2000;174:75-80.

4. Kohn MM, Moores BM, Shibilla H, Schneider K, Stender HS, Stieve FE, et al. European Guidelines On Quality Criteria For Diagnostic Radiographic Images in Paediatrics. Brussels, Luxembourg: European Commission, 1996.

5. Tingberg A. Evaluation Of Image Quality Of Lumbar Spine Images: A Comparison Between FFE And VGA. Radiation protection dosimetry. 2005;144:53-61.

6. Alexakhin RM, Menzel H. The 2007 Recommendations of the International Commission on Radiological Protection. 2007.

7. Sund P. Comparison of visual grading analysis and determination of detective quantum efficiency for evaluating system performancein digital chest radiography. Eur Radiol. 2004.

8. Bath M. Visual grading characteristics (VGC) analysis: a non-parametric rank-invariant statistical method for image quality evaluation. British Journal of Radiology 2007;Volume 80:pp. 169-76.

Page 22 of 24

9. A. Tingberg P, C. Herrmann P, B. Lanhede M, A. Almén P, M. Sandborg P, McVey G, et al. Influence of the Characteristic Curve on the Clinical Image Quality of Lombar Spine and Chest Radiographs. British Journal of Radiology. 2004;77:204-15.

10. Tingberg A. Quantifying the Quality of Medical X-Ray Images: An Evaluation Based on Normal Anatomy Lumbar Spine and Chest Radiography: Lund University, Malmö; 2000.

11. Pereira A, Duarte A, Matela N, Pereira P. Avaliação dos Padrões de Dose em Radiologia Pediátrica e na Radiografia Simples do Abdómen. Revista Lusófona de Ciências e Tecnologias da Saúde. 2009;2:236-44.

12. Ganten M, Radeleff B, Kampschulte A, Daniels MD, Kauffmann GW, Hansmann J. Comparing Image Quality of Flat-Panel Chest Radiography with Storage Phosphor Radiography and Film-Screen Radiography. American Journal of Roentgenology. 2003;181:171-6.

13. Bacher K, Smeets P, Bonnarens H, Hauwere AD, Verstraete K, Thierens H. Dose Reduction in Patients Undergoing Chest Imaging: Digital Amorphous Silicon Flat-Panel Detector Radiography Versus Conventional Film-Screen Radiography and Phosphos-Based Computed Radiography. American Journal of Roentgenology. 2003;181:923-9.

# **Personal Information**

Oliveira, Anaïs - School of Health (ESSUA) - University of Aveiro, Aveiro, Portugal

Martins, Alexandra - School of Health (ESSUA) - University of Aveiro, Aveiro, Portugal

Avelãs, Rita - School of Health (ESSUA) - University of Aveiro, Aveiro, Portugal

Santos, Mónica - School of Health (ESSUA) - University of Aveiro, Aveiro, Portugal

Martins, Paula (BSc Radiography, MSc Speech and Hearing Sciences) - School of Health (ESSUA)/IEETA - University of Aveiro, Aveiro, Portugal

De Francesco, Sílvia (PhD Electrotecnic Engineer) - School of Health (ESSUA)/IEETA - University of Aveiro, Aveiro, Portugal

Page 23 of 24

Sá-Couto, Pedro - Department of Matematics/School of Health - University of Aveiro, Aveiro, Portugal

Ferreira, Carlos (BSc Radiography) - IBILI, University of Coimbra, Coimbra, Portugal

Page 24 of 24

European Society of Radiology | www.myESR.org