brought to you by DCORE

# Actinobaculum schaalii: identification with MALDI-TOF

#### T. Tuuminen<sup>1,2,</sup>, P. Suomala<sup>2</sup> and I. Harju<sup>3</sup>

 Department of Bacteriology and Immunology, Haartman Institute, University of Helsinki, Helsinki, 2) Eastern Finland Laboratory Centre Joint Authority Enterprise, Mikkeli District Laboratory, Finland and 3) Clinical Microbiology Laboratory, Turku University Hospital, Turku, Finland

#### Abstract

Actinobaculum schaalii is an emerging uropathogen. So far, its identification has been performed with 16S rRNA gene sequencing or PCR. The diagnosis has often been delayed due to fastidious growth and identification problems. Eleven clinical isolates of *A. schaalii* from bloodstream infections that were initially identified with 16S rRNA sequencing analysis were recovered and later identified with matrix-assisted laser desorption/ionization time-of-flight (MALDI-TOF). We present a review of bacteriological data of these patients, an algorithm for fast laboratory work-up and advocate the use of sensitized culture of urine to allow better recovery of *A. schaalii* in susceptible patients.

**Keywords:** Actinobaculum schaalii infection, bacteriological diagnostics, emerging infections, Gram-positive rod, urinary tract infection

Original Submission: 31 October 2013; Revised Submission: 2 December 2013; Accepted: 4 December 2013

Article published online: 5 February 2014

New Microbe New Infect 2014; 2: 38-41

\*Corresponding author: Tamara Tuuminen, Porrassalmenkatu 35-37, FIN-50100 Mikkeli, Finland. E-mail: tamara.tuuminen@islab.fi

#### Introduction

Actinobaculum schaalii is not yet a well-known bacterium but has been increasingly reported as an emerging pathogen causing urinary tract infections (UTIs) in patients with underlying uropathologies [1-8] or infections of other organs [9-11]. This pathogen has also been reported to cause urinary infections in children with renal or neurological disorders [12, 13]. Due to its fastidious growth in microaerophilic environments A. schaalii may be overlooked when traditional culture methods are applied. In addition, A. schaalii may be neglected and misinterpreted as a representative of normal microbiota of mucosa or skin [3, 6, 7], especially when present as a co-isolate with 'classical' pathogens of UTIs, such as *Escherichia coli*. In many reported cases diagnosis was delayed, and the final identification to the species level was achieved after 16S rRNA gene sequencing [1-13]. What makes this pathogen important is its peculiar antimicrobial resistance to ciprofloxacin and co-trimoxazole, the drugs often prescribed empirically for the treatment of UTIs.

Urine from urological or immunocompromised patients can be cultured using a sensitized procedure. This procedure includes, in addition to the conventional agar for Gram-negative bacteria, the cultivation on a chocolate agar with a 10  $\mu$ L loop to increase sensitivity.

Below we present results that advocate the use of a sensitized culture in risk-group patients. We also demonstrate the key metabolic reactions, an algorithm to increase the rate of laboratory work-up, and the suitability of matrix-assisted laser desorption/ionization time-of-flight (MALDI-TOF) for identification of A. *schaalii*.

## **Materials and Methods**

Eleven clinical isolates collected from 2008 through 2012 were from five district laboratories. All were isolated from the routine blood culture (bioMerieux, Marcy L'Etoile, France) and initially identified using 16S rRNA gene sequencing (Huslab, Helsinki, Finland) according to the previously published protocol [14]. Ethics approval was not requested because the data presented were collected only from laboratory records.

All isolates were stored at  $-70^{\circ}$ C before being retrieved and cultured on anaerobe basal agar with 10% horse blood (ABA; Oxoid Ltd, Basingstoke, UK), chocolate, and sheep blood agar at  $+35^{\circ}$ C in a CO<sub>2</sub> atmosphere. Retrospectively, phenotypical identification of all the strains was performed with Rapid ID32A; API Coryne and VITEK 2 (bioMerieux). The susceptibility testing was performed from a suspension of McFarland 6 with E-tests (bioMerieux) or with disc diffusion on ABA, 48 h.

Before MALDI-TOF analysis the isolates were sub-cultured on fastidious anaerobe agar (FAA) with 5% sheep blood and incubated for 48 h. Each isolate was initially inoculated directly from a colony onto a steel target plate on two consecutive spots. One microlitre of matrix solution ( $\alpha$ -cyano-4-hydroxycinnamic acid in 50% acetonitrile, 2.5% trifluoroacetic acid; Bruker Daltonics, Bremen, Germany) was added. The identification was performed by MALDI-TOF MS using Microflex LT

© 2014 The Authors. New Microbes and New Infections published by John Wiley & Sons Ltd on behalf of the European Society of Clinical Microbiology and Infectious Disease.

Open access under CC BY-NC-ND license

This is an open access article under the terms of the Creative Commons Attribution-NonCommercial License, which permits use,

distribution and reproduction in any medium, provided the original work is properly cited and is not used for commercial purposes.

instrument with Flex Control 3.0 software and Maldi Biotyper DB Update\_V3.2.1.1 (Bruker Daltonics) at the Laboratory of Clinical Microbiology, Turku University Hospital. If both scores from the initial run gave a value of <2.0, the bacterial material was extracted as specified by the manufacturer and re-analysed.

To test the sensitized culture midstream urine from a young healthy male was collected and artificial urine samples were prepared from three isolates with serial dilutions of  $10^3$ – $10^6$  CFU/mL. These samples and unspiked urine were cultured on ABA, chocolate agar and ChromID<sup>TM</sup>CPS3 (bioMerieux) with a 10 µL loop. The results were read 24 and 48 h later.

# **Results and Discussion**

The patient-related data and laboratory identification of recovered isolates are presented in Table I. All II strains were successfully identified to the species level (correct identification with a score value  $\geq 2.0$ ) but 10/11 strains required an extraction step with 70% formic acid for reliable identification. Only one strain (S-21528) was reliably identified to species level directly from the colony.

The simulated urine samples produced visible growth on ABA already after 24 h but on the sheep blood and chocolate agars the cultivation required 48 h. No growth was detectable on the CPS3. The number of colonies was proportional to the strength of the spiked urine. The algorithm for the identification of *A. schaalii* is presented in Fig. 1.

The retrospective analysis of the patient data indicated that urinalysis was requested from all the patients in whom A. schaalii was isolated from blood; however, aerobic mixed culture and non-significant growth of Corynebacterium species or Lactobacillus was reported in 6/11 patients. In the remainder (5/11) the urinalysis resulted in identification of bacteria associated with UTIs but as expected, A. schaalii was not recovered. In 4/11 patients co-infection with other pathogens was detected in the bloodstream (Table I). RAPID ID 32A provided the most sustained identification profile of 04XXX77717, where the X stands for any reaction (Table 1). This profile suggests Capnocytophaga spp. or Micromonas micros, but these species can easily be ruled out. This commercial test strip can also misinterpret A. schaalii for Actinomyces meyeri; however, the latter can be ruled out on the basis that (1) this pathogen is seldom recovered from blood culture, (2) the origin of infection is often related to oral cavity,



FIG. I. A proposed algorithm for identification of Actinobaculum schaalii from clinical specimens. ABA, anaerobe basal agar.

Isolate	Reported results from urine within 2 week proximity	Growth from aerobic or anaerobic bottle	Concomitant isolates in blood	Maldi Biotyper A. schaalii confirmed, score	Vitek	APICoryne	Rapid ID 32A
S-21495	Corynebacterium spp.	l anaerobic bottle	Peptoniphilus asacchari in I anaerobic bottle	2.219	<b>67</b> 64100010641	2050721	<b>04</b> 10 <b>477717</b>
S-21528	Oligella urethralis and Aerococcus spp.	l anaerobic bottle	Aerococcus species in I aerobic bottle	2.263	2601000000401	0010721	<b>04</b> 2 <b>6477717</b>
M-21846	1. No growth 2. Lactobacillus spp.	I aerobic and I anaerobic bottle each	No	2.124	6727102010661 6764102010621	2050721	<b>04</b> 3 <b>6477717</b>
I-20772	Aerobic mixed culture	L anaerobic bottle	No	2 070	6765100010441	2050721	04 6477717
P-132	Klebsiella pneumoniae	2 aerobic bottles	No	2.117	6625100010401	2050721	0436077717
P-77	Klebsiella oxytoca	l aerobic bottle	No	2.115	6625000010401	0050220	0436477717
K-26636	Aerobic mixed culture	2 anaerobic bottles	No	2.202	6724100010641	0050721	04 4477717
K-27888	Klebsiella pneumoniae	2 aerobic and 1 anaerobic bottles	Enterococcus faecalison in I anaerobic bottle Atopodium parvulum in the same anaerobic bottle as A. schaalii	2.004	2 <b>7</b> 2100000040 <b>1</b>	4050320	<b>04</b> 002 <b>71717</b>
K-27197	<ol> <li>No growth</li> <li>Mixed aerobic culture</li> </ol>	I aerobic and I anaerobic bottle	Peptosteptococcus tetradius in 1 anaerobic bottle	2.151	6725102010661	4010621	<b>04</b> 344 <b>77717</b>
T-41163	1. Mixed aerobic culture 2. Lactobacillus spp.	I anaerobic bottle	No	2.056	6764100010401 6764100010421	4150621	<b>04</b> 344 <b>77717</b>
T-42442	Enterococcus faecalis	l anaerobic bottle	No	2.183	2 <b>7</b> 21000000001	4010000	<b>04</b> 160 <b>7371</b> 4 variability in reactions

#### TABLE 1. Patient-related data, biochemical and MALDI-TOF identification of 11 Actinobaculum schaalii isolates

and (3) the bacteria are not aero-tolerant in comparison with *A. schaalii.* 

All tested isolates were susceptible to penicillin (MIC below 0.008 g/L), cephalexin (disc diffusion, mean 32, range 26–38 mm) and cefuroxime (MIC below 0.016  $\mu$ /L) being in agreement with Cattoir et al. [15]. In one study this pathogen was reported susceptible [16] but in another study [5] 28% of strains were resistant or intermediately susceptible to pivme-cillinam. Our strains had an inhibition zone to pivmecillinam in the range 28–38 mm but were unambiguously resistant to co-trimaxozole with no zone of growth inhibition whereas for ciprofloxacin a slight inhibition zone was detectable for some strains (data not shown).

In our opinion culture, phenotypical examination and MALDI-TOF identification are superior to PCR because the genomic material of *A. schaalii* can be detected with PCR in up to 13% of asymptomatic patients [17]. We do not perform microscopy of urine samples of all patients aged over 60 years, as recommended by Nielsen *et al.* [5]. We believe that screening the urine using flow cytometry sorting (FACS) [18] will help identify those samples that require special culture.

In conclusion, we advocate the sensitized culture to increase the yield of *A. schaalii* recovery from urine in risk-group patients and MALDI-TOF for fast identification.

## Acknowledgement

We thank our colleagues for sending us their A. schaalii isolates.

# **Conflict of Interest**

None declared.

#### References

- Fendukly F, Osterman B. Isolation of Actinobaculum schaalii and Actinobaculum urinale from a patient with chronic renal failure. J Clin Microbiol 2005; 43: 3567–3569.
- Sturm PD, Van Eijk J, Veltman S, Muleman E, Schüllin T. Urosepsis with Actinobaculum schaalii and Aerococcus urinae. J Clin Microbiol 2006; 44: 625–654.
- Reinhard M, Prag J, Kemp M et al. Ten cases of Actinobaculum schaalii infection: clinical relevance, bacterial identification, and antibiotic susceptibility. J Clin Microbiol 2005; 43: 5305–5308.
- Martinaud C, Gaillard T, Maslin J et al. Actinobaculum schaalii bacteria in an aged male patient. Med Mal Infect 2008; 38: 617–619.
- Nielsen H, Søby KM, Christensen JJ, Prag J. Actinobaculum schaalii: a common cause of urinary tract infection in the elderly population. Bacteriological and clinical characteristics. Scan J Infect Dis 2010; 42: 43–47.
- Gomez E, Gustafson DR, Rosenblatt JE, Patel R. Actinobaculum bacteremia: a report of 12 cases. J Clin Microbiol 2011; 49: 4311–4313.
- Beguelin C, Genne D, Vara A et al. Actinobaculum schaalii: clinical observation of 20 cases. Clin Microbiol Infect 2011; 17: 1027–1031.
- García-Bravo M, González-Fernández MB, García-Castro MA, Jaime-Muniesa ML. Urinary tract infection caused by Actinobaculum schaalii in an elderly patient. Rev Esp Quimioter 2011; 24: 52–53.
- Agudo S, Domingo D, Navarro JL, Arenal N, López-Brea M. Actinobaculum schaalii: a new cause of mastitis. Clin Microbiol Infect 2009; 15: S4.
- Hoenigl M, Leitner E, Valentin T et al. Endocarditis caused by Actinobaculum schaalii. Emerg Infect Dis 2010; 16: 1171–1173.
- 11. Vanden Bempt I, Van Trappen S, Cleenwerck I et al. Actinobaculum schaalii causing Fournier's Gangrene. J Clin Microbiol 2011; 49:2369– 2371.

- Pajkrt D, Simmons-Smit AM, Savelkoul PH, van den Hoek J, Hack WW, van Furth AM. Pyelonephritis caused by *Actinobaculum schaalii* in a child with pyeloureteral junction obstruction. *Eur J Clin Microbiol Infect Dis* 2003; 53: 679–682.
- Zimmermann P, Berlinger L, Liniger B, Grunt S, Agyeman P, Ritz N. Actinobaculum schaalii an emerging pediatric pathogen? BMC Infect Dis 2012; 12: 201.
- Jalava JM, Mäntymaa L, Ekblad U et al. Bacterial 16S rDNA polymerase chain reaction in the detection of intra-amnionitic infection. Br J Obstet Gynaecol 1996; 103: 664–669.
- Cattoir V, Varca A, Greub G, Prod'hom G, Legrand P, Lienhard R. In vitro susceptibility of Actinobaculum schaalii to 12 antimicrobial agents

and molecular analysis of fluoroquinolone resistance. J Antimicrob Chemother 2010; 65: 2514–2517.

- Andersen PK, Søby KM, Bank S, Prag J. In vitro susceptibility of Actinobaculum schaalii to mecillinam. J Antimicrob Chemother 2011; 66: 2181–2182.
- Bank S, Jensen A, Hansen TM, Søby KM, Prag J. Actinobaculum schaalii, a common uropathogen in elderly patients, Denmark. Emerg Infect Dis 2010; 16: 76–80.
- Jolkkonen S, Paattiniemi EL, Kärpänoja P, Sarkkinen H. Screening of urine samples by flow cytometry reduces the need for culture. J Clin Microbiol 2010; 48: 3117–3121.