

Technical University of Denmark



New approaches for determining solubility of organic chemicals

Birch, Heidi; Trac, Ngoc Lam; Mayer, Philipp

Publication date:
2018

Document Version
Publisher's PDF, also known as Version of record

[Link back to DTU Orbit](#)

Citation (APA):
Birch, H., Trac, N. L., & Mayer, P. (2018). New approaches for determining solubility of organic chemicals. Poster session presented at SETAC Europe 28th Annual Meeting, Rome, Italy.

DTU Library Technical Information Center of Denmark

General rights

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.



New approaches for determining solubility of volatile liquid chemicals



Heidi Birch*, Lam Ngoc Trac and Philipp Mayer
*hbir@env.dtu.dk

Introduction

Water solubility is a fundamental parameter in environmental risk assessment of chemicals, and is theoretically a simple parameter to determine.

Main challenge for liquid hydrophobic chemicals:

- avoiding formation of dispersions or micro-droplets
- while establishing equilibrium within a reasonable time frame

Two new approaches were directed at solubility determinations:

- 1) Passive dosing from a saturated silicone polymer¹
- 2) Dosing from pure phase liquid through the headspace²

Experimental setup

Passive dosing from saturated silicone

- Cleaned silicone rod (diameter: 3mm) was cut into pieces of 0.4 g (See A, bottom panel)
- Rods were loaded by immersion in the liquid test chemical
- Loaded rods were wiped with damp lint free tissue and rinsed three times with water
- Water was added to the loaded rods and equilibrated on a roller at 20 rpm (B & C)



Equilibration through headspace

- Glass wool was placed in glass inserts in 20 mL vials
- Liquid test chemical was added to the inserts
- Water was added to the vials
- Vials were shaken at 200 rpm, 10mm orbit for 23 hours



Test chemicals: dodecylbenzene, isopropyl myristate, hecylcyclohexane, + α Pinene (D)
Analysis was done using liquid/liquid extraction with iso-octane (1-5 mL iso-octane to 5-10 mL sample) and GC-MS (E). Calibration by ¹³C labeled standard for Dodecylbenzene and external standards for all other compounds.

Silicone loading kinetics

- Loading kinetics were determined for dodecylbenzene
- 30 rods were placed in dodecylbenzene and rolled at 20 rpm.
- At certain time steps, five rods were removed and wiped
- Swelling was determined by weight




Figure 1: Swelling kinetics for 3 mm silicone rod immersed in dodecylbenzene. Fit: First order biphasic model

Swelling times:

- $t_{95\%}$ of max swelling = 11 hours
- $t_{99\%}$ of max swelling = 18 hours

Cross-validation

Measured solubility





Figure 4: Solubility measured by direct passive dosing from 100 and 90% swelling of silicone as a function of solubility measured by headspace dosing. Mean and s.d., n=3

Generating saturated solutions



Legend:
 ▲ 100% Direct passive dosing ■ Headspace equilibration ● Direct passive dosing ■ Headspace
 - - - EpiSuite estimated solubility ◆ 90% Direct passive dosing - - - EpiSuite experimental database

Figure 2: Equilibrium kinetics for direct passive dosing and headspace dosing of isopropyl myristate. Mean and s.d.




Figure 3: Direct passive dosing of dodecylbenzene from silicone rods loaded at varying swelling percentages

Conclusion

- Saturation of 3 mm silicone rod by swelling can be obtained within one day
- Direct passive dosing at 90% swelling agreed with headspace equilibrated solubility measurements
- Direct passive dosing at 100% swelling produced markedly higher concentrations than headspace equilibration and direct passive dosing at 90% swelling – third phase challenges
- Headspace dosing is a promising method for solubility determinations for volatile liquid chemicals
- Direct passive dosing for solubility measurements of hydrophobic chemicals requires further work

References: ¹Stibany, F., Schmidt, S.N., Schäffer, A. and Mayer, P. 2017. Aquatic toxicity testing of liquid hydrophobic organic chemicals – passive dosing exactly at the saturation limit. Chemosphere 167: 551-558.

²Trac, L.N., Schmidt, S.N and Mayer, P. Submitted. Headspace passive dosing of volatile hydrophobic chemicals – aquatic toxicity testing exactly at the saturation level.

Acknowledgements: The authors thank Cefic LRI for financial support (Eco38), Aaron Redman and Dan Letinski for discussions and feedback and Hanne Bøggild for technical assistance.

