

PDF hosted at the Radboud Repository of the Radboud University Nijmegen

The following full text is a publisher's version.

For additional information about this publication click this link.

<http://hdl.handle.net/2066/191891>

Please be advised that this information was generated on 2019-06-01 and may be subject to change.



Clinical glycomics for the diagnosis of congenital disorders of glycosylation

Nurulamin Abu Bakar^{1,2} · Dirk J. Lefeber^{1,2} · Monique van Scherpenzeel¹

Received: 18 September 2017 / Revised: 12 January 2018 / Accepted: 17 January 2018 / Published online: 1 March 2018
© The Author(s) 2018. This article is an open access publication

Abstract

Clinical glycomics comprises a spectrum of different analytical methodologies to analyze glycan structures, which provides insights into the mechanisms of glycosylation. Within clinical diagnostics, glycomics serves as a functional readout of genetic variants, and can form a basis for therapy development, as was described for PGM1-CDG. Integration of glycomics with genomics has resulted in the elucidation of previously unknown disorders of glycosylation, namely CCDC115-CDG, TMEM199-CDG, ATP6AP1-CDG, MAN1B1-CDG, and PGM1-CDG. This review provides an introduction into protein glycosylation and presents the different glycomics methodologies ranging from gel electrophoresis to mass spectrometry (MS) and from free glycans to intact glycoproteins. The role of glycomics in the diagnosis of congenital disorders of glycosylation (CDG) is presented, including a diagnostic flow chart and an overview of glycomics data of known CDG subtypes. The review ends with some future perspectives, showing upcoming technologies as system wide mapping of the N- and O-glycoproteome, intact glycoprotein profiling and analysis of sugar metabolism. These new advances will provide additional insights and opportunities to develop personalized therapy. This is especially true for inborn errors of metabolism, which are amenable to causal therapy, because interventions through supplementation therapy can directly target the pathogenesis at the molecular level.

Abbreviations

2-AA	2-aminobenzoic acid
2-AB	2-aminobenzamide
2-AP	2-aminopyridine
AFP-L3	fucosylated serum alpha-fetoprotein
Apo C-III	apolipoprotein C-III
ARCL2	autosomal recessive cutis laxa type-2
Asn	asparagine
CDG	congenital disorders of glycosylation
CDG-I	CDG type 1
CDG-II	CDG type 2
COG	conserved oligomeric Golgi
COPD	chronic obstructive pulmonary disease (COPD)
ER	endoplasmic Reticulum

ESI	electrospray Ionization
FDA	Food and Drug Administration
FT	Fourier transform
Fuc	fucose
GA	Golgi Apparatus
GAGs	glycosaminoglycans
Gal	galactose
GalNAc	<i>N</i> -acetylgalactosamine
Glc	glucose
GlcA	glucuronic acid
GlcNAc	<i>N</i> -acetylglucosamine
HILIC	hydrophilic interaction liquid chromatography
IEF	isoelectric focusing
LC	liquid chromatography
LLO	lipid-linked oligosaccharides
m/z	mass over charge ratio
MALDI	matrix assisted laser desorption ionization
Man	mannose
MS	mass spectrometry
NGS	next-generation sequencing
PGC	porous graphitized carbon
PNGaseF	<i>N</i> -glycosidase F
QTOF	quadrupole time-of-flight
Ser	serine
Sia	sialic acid

Responsible Editor: Jaak Jaeken

✉ Monique van Scherpenzeel
Monique.vanScherpenzeel@radboudumc.nl

¹ Translational Metabolic Laboratory, Donders Institute for Brain, Cognition, and Behavior, Radboud University Medical Center, Geert Grooteplein 10, Nijmegen 6525 DA, The Netherlands

² Department of Neurology, Donders Institute for Brain, Cognition, and Behavior, Radboud University Medical Center, Nijmegen, The Netherlands

Thr	threonine
UV	ultraviolet
TOF	Time Of Flight
WES	whole exome sequencing

Introduction

Technological advances in science are creating a revolution in the world of clinical diagnostics for rare metabolic disorders. In many cases, establishing a diagnosis via the traditional care route is a complex, lengthy process involving multiple consultations by various clinical specialists. The fact that an omics technique greatly reduces this turnaround time and increases the diagnostic yield has recently been shown for a cohort of 150 patients presenting with complex neurological disorders of suspected genetic origin. They compared the number of solved cases applying either whole exome sequencing (WES) alone (29.3%) or the standard care pathway (7.3%). This improvement in diagnostic yield is significant, without increasing costs compared to the standard diagnostic trajectory (Vissers et al 2017). However, the percentage of solved cases of about 30% is still quite low, which is likely caused by either the lack of sequence coverage of the variant, by disease causes outside the coding sequences, or the presence of too many “variants of unknown significance”. In this light, the important connection of genomics with functional -omics methodologies in the diagnosis of metabolic disorders is recognized more and more. For example, a combination of WES and deep clinical phenotyping was applied to 41 patients with intellectual developmental disorder and unexplained metabolic abnormalities, which resulted in a diagnosis for 28 patients (68%), and a test for targeted intervention on 18 patients (44%) (Van Karnebeek et al 2016). For congenital disorders of glycosylation (CDG), many cases were unsolved until the inclusion of glycomics into clinical practice, to present the functional defect (Jansen et al 2016a; Jansen et al 2016b; Jansen et al 2016c; Van Damme et al 2016; Van Scherpenzeel et al 2014; Tegtmeier et al 2014; Carss et al 2013; Iqbal et al 2013). The most important advantage of integrating functional omics with genomics in the field of inherited metabolic disorders is the opportunities for therapy, which do arise from insights into functional, biochemical pathways. Initial evidence was published for PGM1-CDG (Tegtmeier et al 2014; Wong et al 2017). This review discusses glycomics, its role in CDG diagnostics, presents different methodologies, and ends with future perspectives, in which the horizon of the diagnostic laboratory needs to be broadened to functionally understand new genetic defects.

Introduction in glycosylation

Protein glycosylation is considered to be the most common post-translational modification and is ubiquitously present (Moremen et al 2012). Glycosylation is a non-template driven process involving multiple competing enzymes (e.g., glycosidases and glycosyltransferases) in the endoplasmic reticulum (ER) and the Golgi apparatus (GA) as glycoproteins traffic and mature through the secretory pathway. Therefore, glycan structures are highly diverse, with multiple possibilities for branching and linkage (micro-heterogeneity) and differences in site occupancy (macro-heterogeneity). Glycans are known to have many important biological functions, such as cell-cell, macromolecular (e.g., antibody) and pathogen interactions, protein secretion, protein signaling, and protein folding (Moremen et al 2012; Defaus et al 2014).

Glycosylation changes have been identified in various diseases, ranging from monogenetic inherited disorders such as CDG (Jaeken 2011) to multiple types of malignancies; e.g., ovarian (Leiserowitz et al 2007), colon (Sethi et al 2014) and breast cancer (Lee et al 2014), but also observed in other pathological situations, such as cirrhosis, hepatitis, and neurodegenerative diseases such as Alzheimer’s (Palmigiano et al 2016) and Parkinson disease (Russell et al 2017). Therefore, the analysis of glycosylation in complex biological matrices, which is called glycomics, has become one of the popular -omics after the era of genomics and proteomics, especially for biomarker discovery, treatment monitoring, and also to understand disease mechanisms. This has become more relevant and promising since the majority of tumor biomarkers endorsed by the Food and Drug Administration (FDA) are glycoproteins (Fuzery et al 2013). Most of the markers are only judged based on their expression level. However, for two glycoproteins, glycosylation is included as a biomarker in routine clinical diagnostics, namely carbohydrate-deficient transferrin for the detection of CDG and alcohol abuse, and fucosylated serum alpha-fetoprotein (AFP-L3) for the early diagnosis of hepatocellular carcinoma (Li et al 2001; Leerapun et al 2007). It was shown that by including glycosylation in the analysis, specificity and sensitivity were increased over the native alpha-fetoprotein, and intact transferrin mass spectrometry resulted in the direct diagnosis of several CDG-II subtypes (van Scherpenzeel et al 2015).

There are three types of protein glycosylation in mammalian cells, namely N-, O-, and C-linked glycosylation each having their own subclasses (Moremen et al 2012). The first two are the most common types of protein glycosylation in human cells (Wang et al 2014), and are therefore also the most studied and best characterized types. This review on clinical glycomics will be restricted to these two most common types. N-linked glycans are attached to the polypeptide via amide linkages to asparagine (Asn) side chains, while O-glycans are attached through glycosidic linkages to side chains of

serine (Ser) or threonine (Thr). N-glycans consist of several monosaccharides, such as *N*-acetylglucosamine (GlcNAc), mannose (Man), fucose (Fuc), galactose (Gal), glucose (Glc), and sialic acid (Sia), which are built in a specific order (Cummings 2009). Briefly, N-glycans are assembled in the ER before being further processed and modified in the GA producing three types of N-glycans, namely complex, hybrid, and high mannose (Fig. 1a).

In contrast to N-glycosylation, O-glycosylation is assembled directly onto serine or threonine residues without pre-assembly. Mucin type O-glycans are the most well-known type of O-glycans with *N*-acetylgalactosamine (GalNAc) at the reducing end. Another large group of O-glycans are the glycosaminoglycans (GAGs) on proteoglycans. GAGs are

long, unbranched carbohydrates containing repeated GalNAc or GlcNAc residues combined with glucuronic acid (GlcA) or Gal residues. The other five types of O-glycosylation are O-linked GlcNAc, -Gal, -Man, -Glc, and -Fuc. In all of these classes, one or more genetic deficiencies have been identified (Wopereis et al 2006, Hennet 2012, Freeze et al 2014).

Clinical glycomics methodologies

For protein N-glycan analysis, the sample preparation usually starts with the addition of the enzyme *N*-glycosidase F (PNGaseF), to cleave the N-glycans from the proteins. The




Fig. 1 Schematic representation of the process of protein N-glycosylation, the different types of N-glycans (a) and their analysis by mass spectrometry in CDG (b) (a) The process of N-glycosylation is started when a glycan precursor is assembled by sequential addition of monosaccharides onto the lipid anchor dolichol in the membrane of the ER. The glycan precursor containing 14 monosaccharides (Glc₃-Man₉-GlcNAc₂) is then transferred en bloc to a specific asparagine residue (N) within the consensus sequence N-X-S/T (X = any amino acid except proline, S = serine, T = threonine) in the nascent peptide chain of a protein (e.g., transferrin) being synthesized by a ribosome. After several glucose trimming by ER glucosidases, the glycoprotein is transported to the GA, where the glycans are modified in multiple steps through the action of

various glycosidases (trimming) and glycosyltransferases (prolonging). All N-glycans share the common core structure of Man₃-GlcNAc₂, and are classified into high mannose glycan (only Man residue attached to the core), complex glycan (only GlcNAc residues are attached to the core) and hybrid glycan (combination of Man and GlcNAc residues are attached to the core). (b) The two most common plasma glycomics approaches for CDG characterization are free N-glycans profiling (glycans released from whole or specific glycoproteins by PNGaseF digestion; e.g., total plasma glycoproteomic profiling) and intact protein glycoproteomic profiling (immunopurification protocol without PNGaseF digestion)

method is applicable to purified proteins, as well as on complex biological samples, such as blood, urine, and CSF. In addition to the analysis of protein-released N-glycans, advances in technology allow the analysis of intact glycoproteins, which is fast, robust, does not require enzymatic digestion, and enables high-throughput analysis (Fig. 1b).

Traditionally, derivatization is the common step after releasing N-glycans from proteins, because this type of modification increases the sensitivity of glycan detection as native N-glycans have no significant ultraviolet (UV) absorbance. One could make use of the single reactive carbonyl group at the reducing end of the glycan and perform reductive amination with UV or fluorescent tags, such as 2-aminobenzoic acid (2-AA), 2-aminobenzamide (2-AB), and 2-aminopyridine (2-AP) (Ruhaak et al 2010; Pabst et al 2009), or permethylate glycans (Ciucanu and Kerek 1984; Kang et al 2008), thereby replacing all the hydroxyl groups with methyl ethers. Recently, a method for sialic acid esterification was reported, which not only enhances the stability of acidic glycans during analysis, but also distinguishes alpha2,3 and alpha2,6 sialic acid linkage (Wheeler et al 2009; Reiding et al 2014). Disadvantages of derivatization steps are that they could be incomplete, which creates a biased analysis, as well as the need for additional purification procedures which can cause sample loss (Pabst and Altmann 2011).

Several electrophoretic and chromatographic techniques have been established to separate complex N-glycan mixtures, such as capillary electrophoresis, ion exchange chromatography, hydrophilic interaction liquid chromatography (HILIC), liquid chromatography (LC), and porous graphitized carbon (PGC) (Melmer et al 2011). Recently, PGC-LC-MS has emerged as a popular platform to efficiently separate native glycans without derivatization steps which greatly reduced turnaround time of glycomics sample preparation and enabled the separation of isomeric N-glycans (alpha and beta anomers). Native N-glycan analysis using the PGC-LC approach has been successfully applied for glycan biomarker discovery in ovarian cancer (Hua et al 2013), colorectal cancer (Sethi et al 2015), and lung cancer (Ruhaak et al 2016).

Electrospray ionization (ESI) and matrix assisted laser desorption ionization (MALDI) are the two most common ionization techniques for MS glycan analysis. Both ionization techniques are very sensitive to analyze proteins, peptides, glycans, and lipids as low as picomolar concentrations. ESI is considered as a soft ionization technique, and is able to create multiple charged ions for biomolecules, such as proteins and peptides, which enables these large molecules to be analyzed within the mass over charge ratio (m/z) ranges of the instrument. When coupled to LC and including standards for calibration, ESI is more reliable for quantitative analysis than MALDI (El-Aneed et al 2009; Sturiale et al 2011).

In comparison to ESI, MALDI has the advantages of being robust, fast, and easy to operate. Additionally, MALDI allows analyzing the acidic N-glycans in relative quantities compared to the neutral glycans when derivatization steps, such as permethylation and ethyl esterification, were performed prior to MS to neutralize and stabilize sialic acids. Fourier transform (FT-) and time of flight (TOF)-MS are currently the detectors with the highest resolution (Wuhrer et al 2004; Pabst and Altmann 2008). The second might be preferable considering its capability to analyze large molecules like intact proteins, and the high maintenance costs of FT detectors.

Unfortunately, there is no general O-glycosidase available for enzymatic release of all species of O-glycans. For O-glycan profiling, chemical release by hydrazinolysis or (reductive) beta-elimination needs to be performed, which are harsh methods that always yield some side products. The reductive beta-elimination is the most clean and therefore the most commonly used method, since the simultaneous reduction of the terminal sugar prevents the peeling of the glycan due to the alkaline conditions (degradation from its reducing end) (Thaysen-Andersen and Packer 2014). There are few methods published to simultaneously profile plasma N- and O-linked glycosylation of CDG patients (Faid et al 2007; Xia et al 2013). Because of their limited applicability so far, congenital disorders in the biosynthesis of O-glycans have been identified by genetic approaches. However, for the mucin type O-glycans, intact apolipoprotein C-III (Apo C-III) profiling is available in a diagnostic setting (Wopereis et al 2003). One way to circumvent the need to use harsh conditions to cleave O-glycans is to analyze native glycopeptides, with the additional advantage of keeping the information on the attachment site intact (Hoffmann et al 2016). System-wide mapping of the N- and O-glycoproteome is envisioned in a good review of Thaysen-Andersen and Packer in 2014 (Thaysen-Andersen and Packer 2014). This new frontier in proteomics has the merits of high-resolution MS, complementary fragmentation techniques, and bioinformatic tools.

Application of clinical glycomics for CDG diagnostics

CDG are a group of genetic defects with abnormal glycosylation of proteins, lipids or both. According to the current guidelines (Jaeken et al 2009), CDG is classified into: protein N-glycosylation defects, protein O-glycosylation defects, glycosphingolipid and glycosylphosphatidylinositol anchor glycosylation defects, and multiple glycosylation pathway defects. Currently, 105 distinct types of CDG have been reported (Jaeken and Peanne 2017; Peanne et al 2017).

However, with the emergence of next generation genomics technologies, there were also genes identified that do not directly involve the glycosylation biosynthesis pathway, for example defects in nucleotide sugar transport, defects of vesicular transport, defects in O-mannosylation, O-GlcNAcylation or defects in dolichol biosynthesis, which lead to abnormal protein glycosylation (Hennet 2012). CDG are commonly classified by the localization of the genetic defect (Jaeken et al 2009), in which CDG type 1 (CDG-I) refers to deficient synthesis of the precursor glycan in the ER, also including defects in the cytosol (e.g., PMM2- and PMI-CDG) and the transfer of the LLO to the protein, including various OST defects (e.g., DDOST, SST3A, and SST3B). CDG type 2 (CDG-2) occurs in the GA resulting in modified glycans on the glycoprotein (Lefeber et al 2011). Since the majority of these glycosylation disorders showed a defect in the biosynthetic pathway of protein N-glycosylation (Freeze et al 2014), the analysis of plasma N-glycans using MS plays a significant role in CDG research and diagnostics. The two most common clinical glycomics approaches for CDG characterization are global (N-glycans released from whole serum or plasma glycoproteins, further abbreviated as plasma glycomics) and protein-specific (e.g., intact transferrin MS & Apo C-III MS) glycoprofiling (Fig. 1).

The first application of glycomics was ESI-MS for CDG-I characterization. Transferrin profiles using ESI-MS on PMM2-CDG patients clearly showed two abnormal peaks corresponding to lack of one and both complete glycans (Wada et al 1992a; Wada et al 1992b). ESI-MS of transferrin became more ‘mature’ in 2001, when coupled to LC-MS to create a fast and high-throughput screening test for CDG (Lacey et al 2001). A decade later, the emergence of advanced quadrupole time-of-flight (QTOF) detection in combination with nanoLC-ESI-MS enabled the development of high resolution intact transferrin glycoprofiling which in turn improved the CDG diagnostics (van Scherpenzeel et al 2015). Normal transferrin IEF profiles have been observed in some CDG-I and -II cases, such as ALG14-CDG, ALG11-CDG, MOGS-CDG, SLC35A3-CDG, and SLC35C1-CDG (Lefeber 2016; Al Teneiji et al 2017), as well as in some defects in sugar metabolism like GNE-CDG (Voermans et al 2010), NANS-CDG (Van Karnebeek et al 2016), PGM3-CDG (Stray-Pedersen et al 2014), and also in a tissue-specific and GA homeostasis defect of VPS13B-CDG/Cohen syndrome (Duplomb et al 2014).

Transferrin glycopeptide analysis, which was established by Wada et al in 2004 is important to reveal the structural information on both glycans and proteins. For example, the hybrid type glycan which is diagnostic








for GA mannosidase defects of MAN1B1-CDG can only be found at Asn-432 but not on Asn-630 (Wada 2016). For O-glycan analysis, Wada et al developed plasma Apo C-III MALDI MS (Wada et al 2012).

For plasma glycomics, MALDI MS was used to analyze permethylated total plasma N-glycans (Guillard et al 2009), or combined plasma N-glycan and O-glycans MS (Xia et al 2013). MALDI MS was also suitable to analyze N-glycans released from fibroblast homogenates of SLC35A3 and healthy controls (Edvardson et al 2013), which is a useful alternative when transferrin glycosylation turned out to be normal. The patient cell homogenate showed a shift toward decreased branching of glycans compared to control, likely due to the lack of UDP-GlcNAc, which is a critical factor in the production of β 1,6-branched (tetra-antennary) structures (Sasai et al 2002).

Both total plasma and intact transferrin glycoprofiling have their own advantages in the diagnosis of CDG. Recently, intact transferrin MS has successfully identified a series of novel CDG, such as PGM1-CDG (Tegtmeyer et al 2014) and MAN1B1-CDG (van Scherpenzeel et al 2014). A unique combination of CDG-I (lack of complete glycans) and CDG-II (truncated glycans especially lacking Gal residue) in PGM1-CDG is easily seen by intact transferrin glycoprofiling. Moreover, it is also used to follow-up the biochemical improvement of several patients that underwent an oral D-galactose supplementation (Voermans et al 2017, Wong et al 2017). It has become the primary diagnostic test for CDG, also including fast identification of B4GALT1-CDG, MGAT2-CDG, SLC35A1-CDG, and SLC35A2-CDG (van Scherpenzeel et al 2015).

However, there are several CDG types which cannot be discriminated by intact transferrin glycoprofiling which required total (released) plasma N-glycans for structural insights into the complete mixture of plasma glycoproteins. For example, for diagnosis of SLC35C1-CDG (GDP-fucose transporter defect), global N-glycan mapping is more suitable, because of its higher degree of fucosylation than intact transferrin (Guillard et al 2011). It is also able to detect the N-tetrasaccharide for ALG1-CDG diagnosis and highly abundant abnormal high mannose species of 3 Man- and 4 Man-glycans in PMM2-CDG and MPI-CDG (Zhang et al 2016), and to show abnormal profiles for tissue-specific glycosylation defects for VPS13B-CDG, also called Cohen syndrome (Duplomb et al 2014). In the case of MOGS-CDG which showed a normal transferrin profile, accumulation of several high mannose species from total IgG N-glycans were detected by plasma glycomics (Sadat et al 2014). In summary, Table 1 provides a detailed overview of the glycoprofiling data of total plasma N-glycans and intact transferrin for diagnosis of different CDG subtypes.

Table 1 Overview of plasma glycosylation features from total (released) N-glycans and intact transferrin MS profiling for CDG diagnosis










Gene defects	Released plasma N-glycans MS; total or protein specific			Intact plasma transferrin MS		
	MS Glycoprofiling	MS techniques	Reference	MS Glycoprofiling	MS techniques	Reference
CDG-I (LLO assembly defects): GDP-Man synthesis defects						
<i>PMM2</i> (PMM2-CDG) & <i>MPI</i> (MPI-CDG)	<ul style="list-style-type: none"> Increase of total plasma 3 mannose glycans (Man₃-GlcNAc₂), 4 mannose glycans (Man₄-GlcNAc₂) and sialylated tetrasaccharide glycans (Sia₁-Gal₁-GlcNAc₂).  <ul style="list-style-type: none"> Low ratio of total plasma  /  could discriminate this group (PMM2-CDG and MPI-CDG) from ALG1-CDG. 	Maldi TOF MS	Zhang et al 2016	<ul style="list-style-type: none"> Increase of typical CDG-I glycoforms; lack of both and one complete biantennary glycans (Sia₂-Gal₂-Man₃-GlcNAc₄), respectively. 	ESI-LC QMS & ESI-LC MS, QTOF MS	Lacey et al 2001 & Van Scherpenzeel et al 2015
CDG-I (LLO assembly defects): mannosylation defects						
<i>ALG1</i> (ALG1-CDG)	<ul style="list-style-type: none"> Increase of total plasma sialylated tetrasaccharide glycans (Sia₁-Gal₁-GlcNAc₂), 3 mannose glycans (Man₃-GlcNAc₂) and 4 mannose glycans (Man₄-GlcNAc₂).  <ul style="list-style-type: none"> High ratio of total plasma  /  could discriminate ALG1-CDG from PMM2-CDG and MPI-CDG. 	Maldi TOF MS	Zhang et al 2016	<ul style="list-style-type: none"> Increase of two sialylated tetrasaccharide glycoforms; (Sia₁-Gal₁-GlcNAc₂) and (Sia₃-Gal₃-Man₃-GlcNAc₆). 	ESI-LC MS, QTOF MS	Bengtson et al 2016
CDG-II (glycan processing defects): ER glucosidase defects						
<i>MOGS</i> (MOGS-CDG)	<ul style="list-style-type: none"> Increase of IgG glycans of 14-saccharides (Glc₃-Man₉-GlcNAc₂), 13-saccharides (Glc₃-Man₈-GlcNAc₂) and 12-saccharides (Glc₃-Man₇-GlcNAc₂). 	Maldi TOF MS	Sadad et al 2014	<ul style="list-style-type: none"> Normal transferrin profiles. 		Lefebvre et al 2011
CDG-II (glycan processing defects): glycosyltransferase defects						
<i>MGAT2</i>	<ul style="list-style-type: none"> Most abundant of total plasma complex truncated 	Maldi TOF	Guillard et	<ul style="list-style-type: none"> Most abundant of complex truncated 	ESI-LC	Van
Gene defects						
Released plasma N-glycans MS; total or protein specific			Intact plasma transferrin MS			
Gene defects	MS Glycoprofiling	MS techniques	Reference	MS Glycoprofiling	MS techniques	Reference
(MGAT2-CDG)	<ul style="list-style-type: none"> biantennary glycan lacking one GlcNAc, one Gal and one Sia (Sia₁-Gal₁-Man₃-GlcNAc₃) with the corresponding fucosylated ones (Sia₁-Gal₁-Man₃-GlcNAc₃-Fuc₁). 	MS	al 2011	<ul style="list-style-type: none"> biantennary glycoform lacking one GlcNAc, one Gal and one Sia (Sia₂-Gal₂-Man₆-GlcNAc₆). 	MS, QTOF MS	Scherpenzeel et al 2015
<i>B4GALT1</i> (B4GALT1-CDG)	<ul style="list-style-type: none"> Most abundant of total plasma complex truncated biantennary glycan lacking all galactoses (Man₃-GlcNAc₄) with the corresponding fucosylated ones (Man₃-GlcNAc₄-Fuc₁). 	Maldi TOF MS	Guillard et al 2011	<ul style="list-style-type: none"> Most abundant of complex truncated biantennary glycoform lacking all galactoses (Man₆-GlcNAc₈). 	ESI-LC MS, QTOF MS	Van Scherpenzeel et al 2015
CDG-II (glycan processing defects): Golgi mannosidase defects						
<i>MAN1B1</i> (MAN1B1-CDG)	<ul style="list-style-type: none"> Increase of several total plasma hybrid glycans with the corresponding fucosylated ones; (Sia₁-Gal₁-Man₄-GlcNAc₃), (Sia₁-Gal₁-Man₄-GlcNAc₃-Fuc₁), (Sia₁-Gal₁-Man₅-GlcNAc₃), (Sia₁-Gal₁-Man₅-GlcNAc₃-Fuc₁), and 6 mannose glycans; (Man₆-GlcNAc₂). 	Maldi TOF MS	Rymen et al 2013	<ul style="list-style-type: none"> Increase of two hybrid glycoforms; (Sia₃-Gal₃-Man₃-GlcNAc₇) and (Sia₃-Gal₃-Man₇-GlcNAc₇). 	ESI-LC MS, QTOF MS	Van Scherpenzeel et al 2014 & Van Scherpenzeel et al 2015
CDG-II (glycan processing defects): nucleotide sugar transporter defects						
<i>SLC35C1</i> (SLC35C1-CDG)	<ul style="list-style-type: none"> Decrease of major total plasma fucosylated glycans; e.g. (Sia₂-Gal₂-Man₃-GlcNAc₄-Fuc₁), (Sia₁-Gal₂-Man₃-GlcNAc₄-Fuc₁), (Gal₂-Man₃-GlcNAc₄-Fuc₁) and (Gal₁-Man₃-GlcNAc₄-Fuc₁). 	Maldi TOF MS	Guillard et al 2011	<ul style="list-style-type: none"> Normal transferrin profiles. 		Lefebvre et al 2011
<i>SLC35A1</i> (SLC35A1-CDG)	<ul style="list-style-type: none"> Increase of total plasma complex truncated biantennary glycan lacking one Sia (Sia₁-Gal₂-Man₃-GlcNAc₄) with the corresponding fucosylated ones (Sia₁-Gal₂-Man₃- 	Maldi TOF MS	Mohamed et al 2013	<ul style="list-style-type: none"> Increase of complex truncated biantennary glycoform lacking one and two Sia, respectively; (Sia₃-Gal₄-Man₆-GlcNAc₈). 	ESI-LC MS, QTOF MS	Van Scherpenzeel et al 2015

Overview of the current CDG diagnostic workflow

As one of the most abundant glycoproteins in human plasma, transferrin has been used traditionally as a biomarker for N-

glycosylation defects. Plasma transferrin isoelectric focusing (IEF) is recognized as the classical laboratory method (Jaeken et al 1984) and widely used as a routine screening test for CDG. Human transferrin has two glycosylation sites on Asn-432 and Asn-630, carrying both a bi-antennary glycan







Table 1 (continued)

Gene defects	Released plasma N-glycans MS; total or protein specific			Intact plasma transferrin MS		
	MS Glycoprofiling	MS techniques	Reference	MS Glycoprofiling	MS techniques	Reference
	GlcNAc ₃ -Fuc ₁), (Sia ₁ -Gal ₁ -Man ₃ -GlcNAc ₄), (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄) and (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄ -Fuc ₁), and 5 mannose glycans (Man ₅ -GlcNAc ₂). 			Gal ₁ -Man ₆ -GlcNAc ₈). 		
<i>TMEM199</i> (TMEM199-CDG)	• Not reported	NA	NA	• Increase of several complex truncated glycoform lacking Gal and Sia; (Sia ₃ -Gal ₂ -Man ₆ -GlcNAc ₈), (Sia ₃ -Gal ₃ -Man ₆ -GlcNAc ₈), (Sia ₂ -Gal ₃ -Man ₆ -GlcNAc ₈), (Sia ₂ -Gal ₂ -Man ₆ -GlcNAc ₈) and (Sia ₁ -Gal ₁ -Man ₆ -GlcNAc ₈). 	ESI-LC MS, QTOF MS	Jansen et al 2016b
<i>CCDC115</i> (CCDC115-CDG)	• Increase of several total plasma complex truncated biantennary glycan lacking Gal and Sia, and fucosylated ones; (Gal ₁ -Man ₃ -GlcNAc ₄ -Fuc ₁), (Sia ₁ -Gal ₁ -Man ₃ -GlcNAc ₄), (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄) and (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄ -Fuc ₁). 	Maldi TOF MS	Jansen et al 2016c	• Increase of several complex truncated glycoform lacking Gal and Sia; (Sia ₃ -Gal ₂ -Man ₆ -GlcNAc ₈), (Sia ₃ -Gal ₃ -Man ₆ -GlcNAc ₈), (Sia ₂ -Gal ₃ -Man ₆ -GlcNAc ₈), (Sia ₂ -Gal ₂ -Man ₆ -GlcNAc ₈) and (Sia ₁ -Gal ₁ -Man ₆ -GlcNAc ₈). 	LC/MS, QTOF	Jansen et al 2016c
<i>SLC39A8</i>	• Not reported	NA	NA	• Increase of several complex truncated	ESI-LC	Park et al
Gene defects	Released plasma N-glycans MS; total or protein specific			Intact plasma transferrin MS		
	MS Glycoprofiling	MS techniques	Reference	MS Glycoprofiling	MS techniques	Reference
(<i>SLC39A8</i> -CDG)				glycoform lacking Gal and Sia; (Sia ₃ -Gal ₂ -Man ₆ -GlcNAc ₈), (Sia ₃ -Gal ₃ -Man ₆ -GlcNAc ₈), (Sia ₂ -Gal ₃ -Man ₆ -GlcNAc ₈), (Sia ₂ -Gal ₂ -Man ₆ -GlcNAc ₈), (Sia ₁ -Gal ₂ -Man ₆ -GlcNAc ₈) and (Sia ₁ -Gal ₁ -Man ₆ -GlcNAc ₈). 	TOF MS,	2015
<i>VPS13B</i> (VPS13B-CDG / Cohen syndrome)	• Increase of several total plasma complex truncated biantennary glycan lacking Gal and Sia, and fucosylated ones; (Man ₃ -GlcNAc ₄ -Fuc ₁), (Gal ₁ -Man ₃ -GlcNAc ₄), (Gal ₁ -Man ₃ -GlcNAc ₄ -Fuc ₁), (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄) and (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄ -Fuc ₁). 	Maldi TOF MS	Duplomb et al 2014	• Normal transferrin profiles.	IEF	Duplomb et al 2014
CDG-II (glycan processing defects): Golgi trafficking defects						
<i>COG1</i> (COG1-CDG)	• Increase of several total plasma complex truncated biantennary glycan lacking Gal and Sia, and fucosylated ones; (Man ₃ -GlcNAc ₄ -Fuc ₁), (Sia ₁ -Gal ₁ -Man ₃ -GlcNAc ₄), (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄) and (Sia ₁ -Gal ₂ -Man ₃ -GlcNAc ₄ -Fuc ₁), and 5 mannose glycans (Man ₅ -GlcNAc ₂). 	Maldi TOF MS	Foulquier et al 2006 & Reynders et al 2009	• Increase of complex truncated glycoform lacking one sialic acid; (Sia ₃ -Gal ₄ -Man ₆ -GlcNAc ₈) followed by minor increase of complex truncated glycoform lacking Gal and Sia; (Sia ₂ -Gal ₄ -Man ₆ -GlcNAc ₈) and (Sia ₂ -Gal ₃ -Man ₆ -GlcNAc ₈). 	ESI-LC MS, QTOF MS	Van Scherpenzeel et al 2015
<i>COG4</i> (COG4-CDG)	• Increase of several total plasma complex truncated biantennary glycan lacking Gal and Sia; (Gal ₁ -Man ₃ -GlcNAc ₄), (Sia ₁ -Gal ₁ -Man ₃ -GlcNAc ₄) and (Sia ₁ -Gal ₂ -	Maldi TOF MS	Reynders et al 2009 &	• Not reported		

transferrin isoforms with two sialic acids and no sialic acid, respectively (Fig. 2b). CDG-II IEF patterns show high variability as depicted for MGAT2-CDG (Fig. 2c), MAN1B1-

CDG (Fig. 2d), and B4GALT1-CDG (Fig. 2e), with increased trisialo-transferrin or increased asialo- and monosialo-transferrin, respectively. Untreated hereditary fructosemia and

Table 1 (continued)


Gene defects	Released plasma N-glycans MS; total or protein specific			Intact plasma transferrin MS		
	MS Glycoprofiling	MS techniques	Reference	MS Glycoprofiling	MS techniques	Reference
	 <p>Man₃-GlcNAc₄, and 5 mannose glycans (Man₅-GlcNAc₂).</p>		Xia et al 2013			
<i>COG5</i> (COG5-CDG)	<ul style="list-style-type: none"> Increase of several total plasma complex truncated biantennary glycan lacking Gal and Sia, and fucosylated ones; (Man₃-GlcNAc₄-Fuc₁), (Gal₁- Man₃-GlcNAc₄-Fuc₁), (Gal₂- Man₃-GlcNAc₄-Fuc₁), (Sia₁-Gal₂- Man₃-GlcNAc₄) and (Sia₁-Gal₂- Man₃-GlcNAc₄-Fuc₁). 	Maldi TOF MS	Palmigiano et al 2017	<ul style="list-style-type: none"> Increase of complex truncated biantennary glycoform lacking one Sia (Sia₃-Gal₄-Man₆-GlcNAc₈). 	ESI-LC MS	Paesold-Burda et al 2009
<i>COG6</i> (COG6-CDG)	<ul style="list-style-type: none"> Increase of transferrin glycans of complex truncated biantennary glycan lacking Gal and Sia; (Sia₁-Gal₁- Man₃-GlcNAc₄) and (Sia₁-Gal₂- Man₃-GlcNAc₄). 	Maldi TOF MS	Huybrechts et al 2012	<ul style="list-style-type: none"> Not reported 	NA	NA
<i>COG7</i> (COG7-CDG)	<ul style="list-style-type: none"> Increase of several total plasma complex truncated biantennary glycan lacking Gal and Sia, and fucosylated ones; (Man₃-GlcNAc₄-Fuc₁), (Gal₁- Man₃-GlcNAc₄), (Gal₁- Man₃-GlcNAc₄-Fuc₁), (Sia₁-Gal₁- Man₃-GlcNAc₄), (Sia₁-Gal₂- Man₃-GlcNAc₄) and (Sia₁-Gal₂- Man₃-GlcNAc₄-Fuc₁), and 5 mannose glycans (Man₅-GlcNAc₂). 	Maldi TOF MS	Reynders et al 2009 & Ng et al 2011 & Guillard et al 2011 & Xia et al 2013	<ul style="list-style-type: none"> Not reported 	NA	NA
<i>COG8</i> (COG8-CDG)	<ul style="list-style-type: none"> Increase of total plasma complex truncated biantennary glycan lacking Gal and Sia; (Sia₁-Gal₁- Man₃-GlcNAc₄) and (Sia₁-Gal₂- Man₃-GlcNAc₄). 	Maldi TOF MS	Foulquier et al 2007 & Kranz et al	<ul style="list-style-type: none"> Increase of complex truncated biantennary glycoform lacking one and two Sia, respectively; (Sia₃-Gal₄-Man₆-GlcNAc₈) and (Sia₂-Gal₁-Man₆-GlcNAc₈). 	ESI-LC MS	Kranz et al 2007
Gene defects	Released plasma N-glycans MS; total or protein specific			Intact plasma transferrin MS		
			2007 & Reynders et al 2009			
Mixed CDG-I & II (Mixed ER & Golgi defects)						
<i>PGM1</i> (PGM1-CDG)	<ul style="list-style-type: none"> Increase of several total plasma complex truncated biantennary glycan lacking Gal and Sia; (Man₃-GlcNAc₄), (Sia₁-Gal₁- Man₃-GlcNAc₄) and (Sia₁-Gal₂- Man₃-GlcNAc₄). 	Maldi TOF MS	Xia et al 2013	<ul style="list-style-type: none"> Increase of typical CDG-I glycoforms; lack of both and one complete biantennary glycans (Sia₂-Gal₂-Man₃-GlcNAc₄), respectively. Increase of several complex truncated glycoform lacking Gal and Sia; e.g. (Man₃-GlcNAc₄), (Sia₁-Gal₁- Man₃-GlcNAc₄), (Sia₁-Gal₁-Man₆-GlcNAc₈), (Sia₂-Gal₂-Man₆-GlcNAc₈), (Sia₃-Gal₃-Man₆-GlcNAc₈) and (Sia₃-Gal₄-Man₆-GlcNAc₈). 	LC/MS, QTOF	Tegtmeyer et al 2014 & Van Scherpenzeel et al 2015

Glycan features identified by either one of these glycomics methods were listed here. It should be noted that patient numbers are in general too small to establish statistically correct ranges of relative abundance for the different genetic defects. However, the presence of specific glycan features is often already highly characteristic for specific genetic disorders

galactosemia, alcohol abuse, hepatopathy, and bacterial sialidases are known as secondary causes of abnormal transferrin glycosylation.

The diagnostics follow-up for CDG-I profiles started with an enzyme assay in fibroblasts or leukocytes for diagnosis of PMM2-CDG and PMI-CDG. If negative, the next step was to perform lipid-linked oligosaccharides (LLO) in fibroblasts

(Gao and Lehrman 2002), or recently this has been replaced by WES using a filter for CDG-I genes and targeted sequencing of a CDG-I panel (Timal et al 2012). The recent discovery of a novel sialylated N-tetrasaccharide for ALG1-CDG and abnormal small high mannose glycan structures in PMM2-CDG and MPI-CDG has proved plasma N-glycan analysis to also be relevant for CDG-I diagnosis. Hence, plasma



glycomics is highly useful especially for ALG1-CDG, since it is challenging to identify the defect by genetics due to 14

pseudo genes that can complicate the analysis (Zhang et al 2016). Another upcoming diagnostic choice is applying

Fig. 2 Intact transferrin IEF and MS profiles, and total plasma N-glycans MS profiles in healthy control, ALG1-CDG (cdg-Ik), MGAT2-CDG (cdg-IIa), MAN1B1-CDG, and B4GALT1-CDG (cdg-IIb). Deconvoluted QTOF mass spectrum of intact transferrin showed a high intensity of the peak indicating: (a) two bi-antennary fully sialylated (complete) glycoforms (79,557 Da) corresponding with the major presence of tetrasialo-transferrin bands on IEF patterns in healthy control; (b) lack of both (75,146 Da) and one (77,351 Da) complete glycoforms (CDG-I signatures) corresponding with the major increase of asialo- and disialo-transferrin bands on IEF patterns, as well as two minor peaks indicating the N-tetrasaccharide (76,007 & 78,212 Da) in ALG1-CDG; (c) two truncated glycoforms lacking N-acetylglucosamine (78,243 & 78,900 Da) corresponding with the major increase of trisialo- and disialo-transferrin bands on IEF patterns in MGAT2-CDG; (d) two hybrid glycoforms (79,059 & 79,221 Da) corresponding with the major increase of trisialo-transferrin bands on IEF patterns in MAN1B1-CDG; (e) five truncated glycoform lacking galactose (77,541, 77,743, 77,947, 78,196, & 78,399 Da), which are mostly non-sialylated species, corresponding with the major increase of asialo- and monosialo-transferrin on IEF patterns (mimicking CDG-I profiles) in B4GALT1-CDG. Extracted compound chromatograms (ECCs) of isomeric native total plasma N-glycans profiling showed the most abundant compound peaks indicating: (f) di- and mono-sialylated bi-antennary N-glycans (5.6 & 7.2 min, respectively) in healthy control; (g) N-tetrasaccharide (6.8 & 8.2 min) in ALG1-CDG; (h) mono-sialylated N-glycans lacking N-acetylglucosamine (6.6 & 7.4 min) and its fucosylated species (7.6 & 8.3 min) in MGAT2-CDG; (i) hybrid N-glycans (6.1 & 6.7 min) and its fucosylated species (7.0 & 7.6 min) in MAN1B1-CDG; (j) two non-sialylated N-glycans lacking galactose, tri-antenna (2.9 & 3.4 min), and bi-antenna (4.2 & 4.9 min), and its fucosylated species; tri-antennary (4.1 & 4.4 min) and bi-antennary glycans (5.3 & 5.9 min) in B4GALT1-CDG

WES directly and Sanger sequencing to confirm the gene defect, but up till now it takes more time to result and is more expensive than conventional CDG screening (Van Scherpenzeel et al 2016).

The diagnostic follow-up approach for CDG-II profiles is an IEF assay of Apo C-III to profile the mucin type O-glycans (Wopereis et al 2003), to distinguish between an exclusive N-glycosylation defect and a combined disorder of N- and O-glycosylation (Wopereis et al 2007). Apo C-III has only a single O-glycan on Thr-74 which is terminally modified by up to two sialic acids to generate three main IEF isoforms: Apo C-III₀ (no sialic acid), Apo C-III₁ (one sialic acid), and Apo C-III₂ (two sialic acids). Decreased sialylation on Apo C-III profiles has been reported in conserved oligomeric Golgi (COG) defects (Spaapen et al 2005; Foulquier et al 2006; Foulquier et al 2007; Kranz et al 2007; Morava et al 2007; Wopereis et al 2007; Zeevaert et al 2008; Ng et al 2011; Palmigiano et al 2017) and autosomal recessive cutis laxa type-2 (ARCL2) due to ATP6V0A2 dysfunctions (Morava et al 2005; Kornak et al 2008). The limitation of Apo C-III IEF is that it is not able to differentiate between the three possible Apo C-III₀ isoforms; the “real unglycosylated Apo C-III” and Apo C-III with two non-sialylated monosaccharides namely Gal and GalNAc. They can easily be separated by Apo C-III MALDI MS (Wada 2016).


So for more structural insight in glycans especially in CDG-II patients, MS is always the method of choice. This rapid profiling of abnormal glycans can be linked to potential gene defects based on the knowledge of glycosylation pathways. For example (Fig. 2): Detection of N-tetrasaccharide glycans in ALG1-CDG might be explained by the modification (galactosylation and sialylation) of the chitobiose glycan core, which accumulates due to cytosolic mannosyltransferase defects; Accumulation of truncated N-glycans lacking GlcNAc in MGAT2-CDG is likely due to the N-acetylglucosaminyltransferase defects in GA; Accumulation of truncated N-glycans lacking Gal could be due to B4GALT1 gene defect, encoding for galactosyltransferase in GA. Accumulation of hybrid N-glycans in MAN1B1-CDG is due to the ER mannosyltransferase defects. The structural information is important to narrow down the number of candidate genes, thereby greatly minimizing the time to diagnosis. The complete diagnostic workflow as described above is depicted in Fig. 3.

Outlook

Based on the identification of several new genetic defects over the past four years by a combination of high-resolution intact transferrin analysis with WES, it might be concluded that clinical glycomics and NGS technology work in synergy by reducing the number of candidate genes and turn-around time to identify a CDG subtype. This approach creates opportunities to identify new genetic defects, but beyond this observation, glycomics has the benefit that the data do not only reflect the genes, but also environmental influences, which create possibilities for therapy monitoring and intervention. The first example was PGM1-CDG (Voermans et al 2017), but nowadays, several new genetic defects appear to be involved in sugar metabolism. Surprisingly, not all these genetic defects show abnormal transferrin glycosylation. Examples are NANS, GNE, PGM3, and Cohen syndrome. NANS and GNE are genetic defects in the sialic acid biosynthesis pathway, and Cohen syndrome patients have a mutation in *VPSI3B*, a protein important for proper GA function. It is known that there are tissue-specific glycosylation processes, for example gamma-glutamyl transpeptidase synthesized by the liver or the kidney (West et al 2010), and that there will be other proteins than transferrin which do show abnormal glycosylation, e.g., muscle proteins for GNE patients.

Native glycopeptide profiling would be an attractive method to obtain insight into protein-specific glycosylation. For some isolated proteins like human transferrin, human alpha-1-acid glycoprotein, influenza A virus hemagglutinin, and human IgG, this has been performed (Khatri et al 2014). If this could work for highly complex protein mixtures like serum or

Fig. 3 Clinical diagnostic flowchart for CDG. The combination of mass spectrometry with clinical exome sequencing and clinical phenotyping allows facile identification of the majority of known CDG-I and CDG-II subtypes in a diagnostic setting. The gene defects mentioned here are examples of characteristic diagnostic glycomics profiles in our laboratory. Please refer to Jaeken and Peanne (2017) and Peanne et al (2017) for clinical symptoms in CDG



plasma (Zielinska et al 2010), an enormous amount of data would show up and provide us with new mechanistic insights into tissue-specific glycosylation and likely yield several glycoprotein biomarkers for CDG and also common disorders (Yang et al 2017). Knowledge about the human glycoproteome is still very limited, but this new system-wide mapping technology will allow the study of fundamental questions in glycobiology, such as dynamics, macro- and microheterogeneity, tissue-specific glycosylation, and the function of certain glycans in specific biological contexts (Thaysen-Andersen and Packer 2014). In this context, genetic deficiencies beyond the classical CDG, which are in O-glycosylation, e.g., in O-GalNAc, O-mannose, and O-fucose glycosylation or in glycolipid or glycosaminoglycan biosynthesis (comprehensive review in Hennet 2012), have to be mentioned. There is no general diagnostic test available for these rare genetic disorders because of the structural heterogeneity of O-glycans and their tissue-specific expression. Because of the advancement in genetics and bioinformatics, we foresee glycomics becoming interlinked with metabolomics and proteomics, thereby opening research avenues to unravel protein glycosylation in a tissue- or cell-specific manner, to understand the biochemical mechanisms of glycosylation, and to ultimately develop or improve new (sugar based) therapies for this group of so far mainly untreatable disorders.

Another advance in technology is intact protein profiling. Both in the chromatographic part as well as the data analysis there are recent developments which enable the analysis of protein mixtures to have multiple biomarkers in one assay. A large number of biomarkers used in clinical laboratories are glycopeptides of which variations in glycosylation are not

taken into account, and only protein expression levels are measured. Examples are alpha-1-antitrypsin for chronic obstructive pulmonary disease (COPD), haptoglobin for gastric cancer, and human chorionic gonadotrophin for ovarian and testicular tumors (Van Scherpenzeel et al 2016). For two intact plasma glycoproteins, namely alpha-fetoprotein and transferrin, it is already known that it is essential to include glycosylation analysis for specificity and sensitivity of the marker. There is an enormous potential to improve first line diagnosis when the level of glycosylation is taken into account on top of the current protein expression levels. With the upcoming use of targeted mass spectrometry in clinical laboratories, the glycan part of the biomarker might easily be included to improve sensitivity and specificity of the marker without additional time-to-result.

With the example of effective D-galactose supplementation on the improvement of glycosylation in PGM1-CDG, which was derived from the intact transferrin glycosylation profile (Tegtmeyer et al 2014), there is emerging interest in the application of sugars as supplemental therapies for metabolic disorders or as supportive therapy to improve the mechanism of action of known therapies, such as chemotherapy, for which 2-deoxyhexose was used to enhance the therapeutic effect by inhibiting glycolysis and even induce an effective antitumor immune response (Beneteau et al 2012). The way to obtain insight into these mechanisms will be one new glycomics area, in which the building blocks of glycosylation, comprising sugar-phosphates and nucleotide sugars, are analyzed. By studying the flux through sugar metabolism, dynamic insights will be obtained, which could lead to the development of new therapies and improve current ones. Bioinformatics will

become increasingly important to extract the relevant information out of these big data sets, to visualize, and potentially integrate with other omics layers of information, for better understanding of the complex field of glycobiology.

Acknowledgements We acknowledge support of the Dutch Organization for Scientific Research (ZONMW Medium Investment Grant 40-00506-98-9001 and VIDI Grant 91713359 to DJL; VENI grant 722015012 to MvS) and the Malaysian government grant (R02087 to NAB).

Compliance with ethical standards

Conflict of interest None.

Open Access This article is distributed under the terms of the Creative Commons Attribution 4.0 International License (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution, and reproduction in any medium, provided you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons license, and indicate if changes were made.

References

- Ciucanu I, Kerek F (1984) A simple and rapid method for the permethylation of carbohydrate. *Carbohydr Res* 131:209–217
- Cummings RD (2009) The repertoire of glycan determinants in the human glycome. *Mol BioSyst* 5:1087–1104
- Defaus S, Gupta P, Andreu D, Gutierrez-Gallego R (2014) Mammalian protein glycosylation-structure versus function. *Analyst* 139:2944–2967
- Dorre K, Olczak M, Wada Y et al (2015) A new case of UDP-galactose transporter deficiency (SLC35A2-CDG): molecular basis, clinical phenotype, and therapeutic approach. *J Inher Metab Dis* 38:931–940
- Duplomb L, Duvet S, Picot D et al (2014) Cohen syndrome is associated with major glycosylation defects. *Hum Mol Genet* 23:2391–2399
- Edvardson S, Ashikov A, J alas C et al (2013) Mutations in SLC35A3 cause autism spectrum disorder, epilepsy and arthrogyriposis. *J Med Genet* 50:733–739
- El-Aneed A, Cohen A, Banoub J (2009) Mass spectrometry, review of the basics: electrospray, MALDI, and commonly used mass analyzers. *Appl Spectrosc Rev* 44:210–230
- Faid V, Chirat F, Seta N, Foulquier F, Morelle W (2007) A rapid mass spectrometric strategy for the characterization of N- and O-glycan chains in the diagnosis of defects in glycan biosynthesis. *Proteomics* 7:1800–1813
- Foulquier F, Amyere M, Jaeken J et al (2012) TMEM165 deficiency causes a congenital disorder of glycosylation. *Am J Hum Genet* 91:15–26
- Foulquier F, Ungar D, Reynders E et al (2007) A new inborn error of glycosylation due to a Cog8 deficiency reveals a critical role for the Cog1-Cog8 interaction in COG complex formation. *Hum Mol Genet* 16:717–730
- Foulquier F, Vasile E, Schollen E et al (2006) Conserved oligomeric Golgi complex subunit 1 deficiency reveals a previously uncharacterized congenital disorder of glycosylation type II. *Proc Natl Acad Sci U S A* 103:3764–3769
- Freeze HH, Chong JX, Bamshad MJ, Ng BG (2014) Solving glycosylation disorders: fundamental approaches reveal complicated pathways. *Am J Hum Genet* 94:161–175
- Fuzery AK, Levin J, Chan MM, Chan DW (2013) Translation of proteomic biomarkers into FDA approved cancer diagnostics: issues and challenges. *Clin Proteomics* 10:13
- Gao N, Lehman MA (2002) Analyses of dolichol pyrophosphate-linked oligosaccharides in cell cultures and tissues by fluorophore-assisted carbohydrate electrophoresis. *Glycobiology* 12:353–360
- Guillard M, Gloerich J, Wessels HJ, Morava E, Wevers RA, Lefeber DJ (2009) Automated measurement of permethylated serum N-glycans by MALDI-linear ion trap mass spectrometry. *Carbohydr Res* 344:1550–1557
- Guillard M, Morava E, van Delft FL et al (2011) Plasma N-glycan profiling by mass spectrometry for congenital disorders of glycosylation type II. *Clin Chem* 57:593–602
- Hennet T (2012) Diseases of glycosylation beyond classical congenital disorders of glycosylation. *Biochim Biophys Acta* 1820:1306–1317
- Hoffmann M, Marx K, Reichl U, Wuhler M, Rapp E (2016) Site specific O-glycosylation analysis of human blood plasma proteins. *Mol Cell Proteomics* 15:624–641
- Hua S, Williams CC, Dimapasoc LM et al (2013) Isomer-specific chromatographic profiling yields highly sensitive and specific potential N-glycan biomarkers for epithelial ovarian cancer. *J Chromatogr A* 1279:58–67
- Huybrechts S, De Laet C, Bontems P et al (2012) Deficiency of subunit 6 of the conserved oligomeric Golgi complex (COG6-CDG): second patient, different phenotype. *JIMD Rep* 4:103–108
- Iqbal Z, Shahzad M, Vissers LE et al (2013) A compound heterozygous mutation in DPAGT1 results in a congenital disorder of glycosylation with a relatively mild phenotype. *Eur J Hum Genet* 21:844–849
- Jaeken J (2011) Congenital disorders of glycosylation (CDG): it's (nearly) all in it! *J Inher Metab Dis* 34:853–858
- Jaeken J, Hennet T, Matthijs G, Freeze HH (2009) CDG nomenclature: time for a change! *Biochim Biophys Acta* 1792:825–826
- Jaeken J, Peanne R (2017) What is new in CDG? *J Inher Metab Dis* 40:569–586
- Jaeken J, Van Eijk HG, van der Heul C, Corbeel L, Eeckels R, Eggermont E (1984) Sialic acid-deficient serum and cerebrospinal fluid transferrin in a newly recognized genetic syndrome. *Clin Chim Acta* 144:245–247
- Jansen EJ, Timal S, Ryan M et al (2016a) ATP6AP1 deficiency causes an immunodeficiency with hepatopathy, cognitive impairment and abnormal protein glycosylation. *Nat Commun* 7:11600
- Jansen JC, Cirak S, van Scherpenzeel M et al (2016c) CCDC115 deficiency causes a disorder of Golgi homeostasis with abnormal protein glycosylation. *Am J Hum Genet* 98:310–321
- Jansen JC, Timal S, van Scherpenzeel M et al (2016b) TMEM199 deficiency is a disorder of Golgi homeostasis characterized by elevated aminotransferases, alkaline phosphatase, and cholesterol and abnormal glycosylation. *Am J Hum Genet* 98:322–330
- Kang P, Mechref Y, Novotny MV (2008) High-throughput solid-phase permethylation of glycans prior to mass spectrometry. *Rapid Commun Mass Spectrom* 22:721–734
- Khatri K, Staples GO, Leymarie N et al (2014) Confident assignment of site-specific glycosylation in complex glycoproteins in a single step. *J Proteome Res* 13:4347–4355
- Kornak U, Reynders E, Dimopoulou A et al (2008) Impaired glycosylation and cutis laxa caused by mutations in the vesicular H⁺-ATPase subunit ATP6V0A2. *Nat Genet* 40:32–34
- Kranz C, Ng BG, Sun L et al (2007) COG8 deficiency causes new congenital disorder of glycosylation type IIh. *Hum Mol Genet* 16:731–741
- Lacey JM, Bergen HR, Magera MJ, Naylor S, O'Brien JF (2001) Rapid determination of transferrin isoforms by immunoaffinity liquid chromatography and electrospray mass spectrometry. *Clin Chem* 47:513–518
- Lee LY, Thaysen-Andersen M, Baker MS, Packer NH, Hancock WS, Fanayan S (2014) Comprehensive N-glycome profiling of cultured

- human epithelial breast cells identifies unique secretome N-glycosylation signatures enabling tumorigenic subtype classification. *J Proteome Res* 13:4783–4795
- Leerapun A, Suravarapu SV, Bida JP et al (2007) The utility of Lens Culinaris agglutinin-reactive alpha-fetoprotein in the diagnosis of hepatocellular carcinoma: evaluation in a United States referral population. *Clin Gastroenterol Hepatol* 5:394–402
- Lefeber DJ (2016) Protein-specific glycoproteomics for patient diagnostics. *Clin Chem* 62:9–11
- Lefeber DJ, Morava E, Jaeken J (2011) How to find and diagnose a CDG due to defective N-glycosylation. *J Inherit Metab Dis* 34:849–852
- Leiserowitz GS, Lebrilla C, Miyamoto S et al (2007) Glycomics analysis of serum: a potential new biomarker for ovarian cancer? *Int J Gynecol Cancer* 18:470–475
- Li D, Mallory T, Satomura S (2001) AFP-L3: a new generation of tumor marker for hepatocellular carcinoma. *Clin Chim Acta* 313:15–19
- Melmer M, Stangler T, Premstaller A, Lindner W (2011) Comparison of hydrophilic-interaction, reversed-phase and porous graphitic carbon chromatography for glycan analysis. *J Chromatogr A* 1218:118–123
- Mohamed M, Ashikov A, Guillard M et al (2013) Intellectual disability and bleeding diathesis due to deficient CMP-sialic acid transport. *Neurology* 81:681–687
- Morava E, Wopereis S, Coucke P et al (2005) Defective protein glycosylation in patients with cutis laxa syndrome. *Eur J Hum Genet* 13:414–421
- Morava E, Zeevaert R, Korsch E et al (2007) A common mutation in the COG7 gene with a consistent phenotype including microcephaly, adducted thumbs, growth retardation, VSD and episodes of hyperthermia. *Eur J Hum Genet* 15:638–645
- Moremen KW, Tiemeyer M, Nairn AV (2012) Vertebrate protein glycosylation: diversity, synthesis and function. *Nat Rev Mol Cell Biol* 13:448–462
- Ng BG, Sharma V, Sun L et al (2011) Identification of the first COG-CDG patient of Indian origin. *Mol Genet Metab* 102:364–367
- Pabst M, Altmann F (2008) Influence of electrosorption, solvent, temperature, and ion polarity on the performance of LC-ESI-MS using graphitic carbon for acidic oligosaccharides. *Anal Chem* 80:7534–7542
- Pabst M, Altmann F (2011) Glycan analysis by modern instrumental methods. *Proteomics* 11:631–643
- Pabst M, Kolarich D, Polt G et al (2009) Comparison of fluorescent labels for oligosaccharides and introduction of a new postlabeling purification method. *Anal Biochem* 384:263–273
- Paesold-Burda P, Maag C, Troxler H et al (2009) Deficiency in COG5 causes a moderate form of congenital disorders of glycosylation. *Hum Mol Genet* 18:4350–4356
- Palmigiano A, Barone R, Sturiale L et al (2016) CSF N-glycoproteomics for early diagnosis in Alzheimer's disease. *J Proteome* 131:29–37
- Palmigiano A, Bua RO, Barone R et al (2017) MALDI-MS profiling of serum O- and N-glycosylation in COG5-CDG. *J Mass Spectrom* 52:372–377
- Park JH, Hogrebe M, Gruneberg M et al (2015) SLC39A8 deficiency: a disorder of manganese transport and glycosylation. *Am J Hum Genet* 97:894–903
- Peanne R, de Lonlay P, Foulquier F et al (2017) Congenital disorders of glycosylation (CDG): quo vadis? *Eur J Med Genet*. <https://doi.org/10.1016/j.ejmg.2017.10.012>
- Reiding KR, Blank D, Kuijper DM, Deelder AM, Wuhrer M (2014) High-throughput profiling of protein N-glycosylation by MALDI-TOF-MS employing linkage-specific sialic acid esterification. *Anal Chem* 86:5784–5793
- Reynders E, Foulquier F, Leao Teles E et al (2009) Golgi function and dysfunction in the first COG4-deficient CDG type II patient. *Hum Mol Genet* 18:3244–3256
- Ruhaak LR, Stroble C, Dai J et al (2016) Serum glycans as risk markers for non-small cell lung cancer. *Cancer Prev Res (Phila)* 9:317–323
- Ruhaak LR, Zauner G, Huhn C, Bruggink C, Deelder AM, Wuhrer M (2010) Glycan labeling strategies and their use in identification and quantification. *Anal Bioanal Chem* 397:3457–3481
- Russell AC, Simurina M, Garcia MT et al (2017) The N-glycosylation of immunoglobulin G as a novel biomarker of Parkinson's disease. *Glycobiology* 27:501–510
- Rymen D, Peanne R, Millon MB et al (2013) MAN1B1 deficiency: an unexpected CDG-II. *PLoS Genet* 9:e1003989
- Sadat MA, Moir S, Chun TW et al (2014) Glycosylation, hypogammaglobulinemia, and resistance to viral infections. *New Eng J Med* 370:1615–1625
- Sasai K, Ikeda Y, Fujii T, Tsuda T, Taniguchi N (2002) UDP-GlcNAc concentration is an important factor in the biosynthesis of beta1,6-branched oligosaccharides: regulation based on the kinetic properties of N-acetylglucosaminyltransferase V. *Glycobiology* 12:119–127
- Sethi MK, Kim H, Park CK et al (2015) In-depth N-glycome profiling of paired colorectal cancer and non-tumorigenic tissues reveals cancer-, stage- and EGFR-specific protein N-glycosylation. *Glycobiology* 25:1064–1078
- Sethi MK, Thaysen-Andersen M, Smith JT et al (2014) Comparative N-glycan profiling of colorectal cancer cell lines reveals unique bisecting GlcNAc and alpha-2,3-linked sialic acid determinants are associated with membrane proteins of the more metastatic/aggressive cell lines. *J Proteome Res* 13:277–288
- Spaapen LJ, Bakker JA, van der Meer SB et al (2005) Clinical and biochemical presentation of siblings with COG-7 deficiency, a lethal multiple O- and N-glycosylation disorder. *J Inherit Metab Dis* 28:707–714
- Stray-Pedersen A, Backe PH, Sorte HS et al (2014) PGM3 mutations cause a congenital disorder of glycosylation with severe immunodeficiency and skeletal dysplasia. *Am J Hum Genet* 95:96–107
- Sturiale L, Barone R, Garozzo D (2011) The impact of mass spectrometry in the diagnosis of congenital disorders of glycosylation. *J Inherit Metab Dis* 34:891–899
- Tegtmeier LC, Rust S, van Scherpenzeel M et al (2014) Multiple phenotypes in phosphoglucomutase 1 deficiency. *N Engl J Med* 370:533–542
- Thaysen-Andersen M, Packer NH (2014) Advances in LC-MS/MS-based glycoproteomics: getting closer to system-wide site-specific mapping of the N- and O-glycoproteome. *Biochim Biophys Acta* 1844:1437–1452
- Timal S, Hoischen A, Lehle L et al (2012) Gene identification in the congenital disorders of glycosylation type I by whole exome sequencing. *Hum Mol Genet* 21:4151–4161
- Van Damme T, Gardeitchik T, Mohamed M et al (2016) Mutations in ATP6V1E1 or ATP6V1A cause autosomal-recessive cutis laxa. *Am J Hum Genet* 100:216–227
- Van Karnebeek CD, Bonafe L, Wen XY et al (2016) NANS-mediated synthesis of sialic acid is required for brain and skeletal development. *Nat Genet* 48:777–784
- Van Scherpenzeel M, Steenbergen G, Morava E, Wevers RA, Lefeber DJ (2015) High-resolution mass spectrometry glycoproteomics of intact transferrin for diagnosis and subtype identification in the congenital disorders of glycosylation. *Transl Res* 166:639–649
- Van Scherpenzeel M, Timal S, Rymen D et al (2014) Diagnostic serum glycosylation profile in patients with intellectual disability as a result of MAN1B1 deficiency. *Brain* 137:1030–1038
- Van Scherpenzeel M, Willems E, Lefeber DJ (2016) Clinical diagnostics and therapy monitoring in the congenital disorders of glycosylation. *Glycon J* 33:345–358
- Vissers L, van Nimwegen KJM, Schieving JH (2017) A clinical utility study of exome sequencing versus conventional genetic testing in pediatric neurology. *Genet Med* 19:1055–1063

- Voermans NC, Guillard M, Doedee R et al (2010) Clinical features, lectin staining, and a novel *GNE* frameshift mutation in hereditary inclusion body myopathy. *Clin Neuropathol* 29:71–77
- Voermans NC, Preisler N, Madsen KL et al (2017) PGM1 deficiency: substrate use during exercise and effect of treatment with galactose. *Neuromuscul Disord* 27:370–376
- Wada Y (2016) Mass spectrometry of transferrin and apolipoprotein C-III for diagnosis and screening of congenital disorder of glycosylation. *Glycoconj J* 33:297–307
- Wada Y, Kadoya M, Okamoto N (2012) Mass spectrometry of apolipoprotein C-III, a simple analytical method for mucin-type O-glycosylation and its application to an autosomal recessive cutis laxa type-2 (ARCL2) patient. *Glycobiology* 22:1140–1144
- Wada Y, Nishikawa A, Okamoto N et al (1992b) Structure of serum transferrin in carbohydrate-deficient glycoprotein syndrome. *Biochem Biophys Res Commun* 189:832–836
- Wada Y, Tajiri M, Yoshida S (2004) Hydrophilic affinity isolation and MALDI multiple-stage tandem mass spectrometry of glycopeptides for glycoproteomics. *Anal Chem* 76:6560–6565
- Wada Y, Tamura J, Musselman BD, Kassel DB, Sakurai T, Matsuo T (1992a) Electrospray ionization mass spectra of hemoglobin and transferrin by a magnetic sector mass spectrometer. Comparison with theoretical isotopic distributions. *Rapid Commun Mass Spectrom* 6:9–13
- Wang YC, Peterson SE, Loring JF (2014) Protein post-translational modifications and regulation of pluripotency in human stem cells. *Cell Res* 24:143–160
- West MB, Segu ZM, Feasley CL et al (2010) Analysis of site-specific glycosylation of renal and hepatic gamma-glutamyl transpeptidase from normal human tissue. *J Biol Chem* 285:29511–29524
- Wheeler SF, Domann P, Harvey DJ (2009) Derivatization of sialic acids for stabilization in matrix-assisted laser desorption/ionization mass spectrometry and concomitant differentiation of alpha(2 → 3)- and alpha(2 → 6)-isomers. *Rapid Commun Mass Spectrom* 23:303–312
- Wong SY, Gadowski T, van Scherpenzeel M et al (2017) Oral D-galactose supplementation in PGM1-CDG. *Genet Med* 19:1226–1235
- Wopereis S, Grunewald S, Huijben KM et al (2007) Transferrin and apolipoprotein C-III isofocusing are complementary in the diagnosis of N- and O-glycan biosynthesis defects. *Clin Chem* 53:180–187
- Wopereis S, Grunewald S, Morava E et al (2003) Apolipoprotein C-III isofocusing in the diagnosis of genetic defects in O-glycan biosynthesis. *Clin Chem* 49:1839–1845
- Wopereis S, Lefeber DJ, Morava E, Wevers RA (2006) Mechanisms in protein O-glycan biosynthesis and clinical and molecular aspects of protein O-glycan biosynthesis defects: a review. *Clin Chem* 52:574–600
- Wuhrer M, Koeleman CA, Deelder AM, Hokke CH (2004) Normal-phase nanoscale liquid chromatography–mass spectrometry of underivatized oligosaccharides at low-femtomole sensitivity. *Anal Chem* 76:833–838
- Xia B, Zhang W, Li X et al (2013) Serum N-glycan and O-glycan analysis by mass spectrometry for diagnosis of congenital disorder of glycosylation. *Anal Biochem* 442:178–185
- Yang Y, Franc V, Heck AJR (2017) Glycoproteomics: a balance between high-throughput and in-depth analysis. *Trends Biotechnol* 35:598–609
- Zeevaert R, Foulquier F, Jaeken J, Matthijs G (2008) Deficiencies in subunits of the conserved oligomeric Golgi (COG) complex define a novel group of congenital disorders of glycosylation. *Mol Genet Metab* 93:15–21
- Zhang W, James PM, Ng BG et al (2016) A novel N-tetrasaccharide in patients with congenital disorders of glycosylation, including asparagine-linked glycosylation protein 1, phosphomannomutase 2, and mannose phosphate isomerase deficiencies. *Clin Chem* 62:208–217
- Zielinska DF, Gnad F, Wisniewski JR, Mann M (2010) Precision mapping of an in vivo N-glycoproteome reveals rigid topological and sequence constraints. *Cell* 141:897–899