

UNCLASSIFIED

# Processing Oscillatory Data with the PDV (U)



2nd Annual PDV Workshop  
August 16-17, 2007

Ted Strand, Randy Bonner,  
Daron Hester, Jerome Solberg,  
Dan White, Tony Whitworth


Lawrence Livermore National Laboratory

This work was performed under the auspices of the U.S. Department of Energy by the University of California, Lawrence Livermore National Laboratory under contract No. W-7405-Eng-48.

UNCLASSIFIED


UNCLASSIFIED

We recently fielded the PDV on a series of pulsed power tests




UNCLASSIFIED

The test fixture provides validation data for new ALE3D codes




Large magnetic pressure tends to collapse inner tube



UNCLASSIFIED


We had 4 probes mounted at the mid-plane at 90 degree angles



UNCLASSIFIED


UNCLASSIFIED

# This portable rack holds a 4-channel Photonic Doppler Velocimeter System




UNCLASSIFIED

We were told to expect low velocities  
(less than 100 m/s)



UNCLASSIFIED

The beat amplitude goes to zero at zero velocity when the moving surface changes direction




UNCLASSIFIED

UNCLASSIFIED

The spectrogram looks empty because the velocities are so low



Processed with  $1.64 \mu\text{s}$  FT windows



UNCLASSIFIED




UNCLASSIFIED

Expand the image near  $f = 0$



Processed with  $1.64 \mu\text{s}$  FT windows




UNCLASSIFIED

UNCLASSIFIED

The spectrogram shows the oscillatory behavior



Analyzing this is very cumbersome




UNCLASSIFIED

UNCLASSIFIED

This laser has a sideband at 1.5 GHz  
Expand the image around the sideband



Processed with  $1.64 \mu\text{s}$  FT windows




UNCLASSIFIED

UNCLASSIFIED

We see the sum and difference velocities  
on either side of the 1.5 GHz sideband




1165 m/s →



UNCLASSIFIED

UNCLASSIFIED


Use the analysis code to extract the desired velocities



UNCLASSIFIED

UNCLASSIFIED


Use the analysis code to extract the desired velocities



UNCLASSIFIED

UNCLASSIFIED


Extract the velocity time history  
from the spectrogram



Sideband velocity  
= 1164.62 m/s

UNCLASSIFIED


Subtract the sideband velocity and integrate to obtain displacement vs time



Velocity vs Time

integrate

Displacement vs Time



This analysis is much easier than working with the baseband solution



## Conclusions



We used the PDV to obtain data from oscillatory motions

Analyzing the baseband velocity is cumbersome

Fortunately, our laser has a 1.5 GHz sideband

Analyzing the data around the sideband is much easier

Moral: a sideband in your laser is not always bad, but we will send the laser to the shop after this series of experiments is over