

RAPID METHOD FOR CLASSIFICATION OF CHEDDAR CHEESE BASED ON FLAVOR QUALITY USING INFRARED SPECTROSCOPY

Anand Subramanian, James Harper, and Luis Rodriguez-Saona

Department of Food Science and Technology, The Ohio State University, Columbus, OH 43210.

ABSTRACT

Analysis of Cheddar cheese flavor using human taste panels is expensive and time consuming. A rapid and simple instrumental method based on Fourier transform infrared spectroscopy was developed for classifying Cheddar cheese based on flavor quality. Fifteen Cheddar cheese samples from two different commercial production plants were ground into powders using liquid nitrogen. The water-soluble compounds from the cheese powder, without interfering compounds such as complex fat and protein, were extracted using water and organic solvents. Aliquots (10 µL) of the extract were placed on the sampling accessory, dried and scanned in the spectrometer (4000 to 700 cm⁻¹). The spectra were matched with the flavor quality to build statistical classification models. The models provided 3D plots in which all the 15 cheese samples formed well separated clusters, whose orientation correlated well with their cheese flavor characteristics (fermented, unclean, low flavor, sour, good cheddar, etc.). The discrimination of the samples was mainly due to organic and fatty acids and their esters (1500 to 900 cm⁻¹), which are known to contribute significantly to cheese flavor. The total analysis time, including the sample preparation time, was less than 20 min per sample. This technique can be a rapid, inexpensive, and simple tool to the cheese industry for predicting the flavor quality of cheese.

INTRODUCTION

About 9.13 billion pounds of cheese is produced in the US every year, of which 34% is cheddar cheese. Flavor quality of cheddar cheese significantly influences its consumer acceptance, price and food processing application. Cheese flavor is directly affected by the physical, chemical and microbiological changes that occur during ripening (aging). Cheese flavor is currently determined using trained human tasting panels. This process is time consuming and very expensive. Predicting the flavor guality of cheese using instrumental methods can save time and money for the industry. However, heterogeneous composition of cheese make it difficult to develop rapid and reliable instrumental methods for flavor quality analysis.

Fourier transform Infrared (FT-IR) spectroscopy, which utilizes the light absorbing properties of chemical compounds, can be used as a rapid, inexpensive, and sensitive method to analyze cheese flavor. Unlike many chromatographic techniques, FT-IR spectroscopy provides unique overall chemical fingerprints (spectra) of cheese samples. FT-IR spectra can be analyzed through multivariate statistical techniques to rapidly discriminate Cheddar cheese based on flavor.

OBJECTIVE

The objective of this research was to develop a rapid and simple technique based on FT-IR spectroscopy to predict the flavor quality of Cheddar cheese.

A rapid instrumental method for cheese flavor quality

SIGNIFICANCE

- analysis will
- > Save time and money for the cheese industry
- Enable better quality control
- > Provide better control over the ripening process to achieve cheese of desired flavor quality
- > Initiate and accelerate further studies on the cheese ripening process using spectroscopy







 $\mathbf{0} \cdot \mathbf{S}$

& TECHNOLO

Fig 1. Typical FT-IR Spectrum of Cheddar Cheese Extract







Fig 4. Spectral Regions Responsible for Classification

DISCUSSION

RESULTS

- > The developed extraction method yielded cheese extracts whose FT-IR fingerprints were well defined, unique for each sample and very consistent (Fig 1).
- > All the 5 samples from Plant #1 formed distinct clusters in 3D space (Fig 2). Fermented and Unclean samples clustered far from the good and better samples (Good, Slight Burn and Slight Sour).
- > Samples from Plant #2 also formed clusters whose orientation correlated with their flavor. The Fermented. Sour and Low Flavor samples clustered away from the Good Cheddar (Fig 3). Samples with multiple attributes (Slight Low Flavor & Slight Sour) clustered in between Sour samples and Low Flavor samples.
- > The discrimination of the samples was due to organic and fatty acids and/or their esters (900 to 1500 cm⁻¹). which are known to contribute significantly to Cheddar cheese flavor (Fig 4).
- > The total analysis time, including the sample preparation time, was less than 20 min per sample.

- The developed sample preparation met and FT-IR technique show great poter as a rapid, inexpensive, and simple qu control tool for determination of fla quality of cheese.
- Preliminary experiments with Swiss che provided similar results. Therefore, method has potential to be applied to o types of cheese.

REFERENCES

· Rodriguez-Saona, L.E., Koca, N., Harper, W.J., and V.B. Alvarez Rapid determination of Swiss cheese composition by Fourier tra infrared/attenuated total reflectance spectroscopy. J. Dairy Sci. 89 1412.

Chen, M., Irudayaraj, J., and D.J. McMahon. 1998. Examination of and reduced fat Cheddar cheese during ripening by Fourier t infrared spectroscopy. J. Dairy Sci. 81: 2791-2797

This project was funded by DairiConcepts LP. (Springfield, MO) and M Advanced Food Manufacturing Alliance (MAFMA)

brought to

you by

R

CORE

Fig 2. 3D Classification Model for Plan