

COMB-REFERENCED COHERENT RAMAN SPECTROSCOPY ON PURE H₂

M. Lamperti, L. Rutkowski, D. Gatti, B. Alsaif, R. Gotti, G. Cerullo, D. Polli, M. Marangoni,
F. Thibault, P. Maslowski, S. Wojtewicz, P. Wcislo



POLITECNICO
MILANO 1863




THE IMPORTANCE OF H₂

- Simplest molecule and most abundant in the universe
 - Benchmark system for *ab-initio* calculations of transitions frequencies and line profiles
 - Transition frequencies:
 - Test for new physics without the need for a particle accelerator
 - Test of invariance of fundamental constants
 - Interpretation of observational data (interstellar clouds...)
- Presence of strong collisional effects
 - Low-pressure measurements are favored to extract zero-pressure transition frequencies
 - High pressure gives much better SNR and requires better collision modeling → *Raman spectroscopy*

A. Owyong, Opt. Lett. **2** (4), 1978; G. J. Rosasco and W. S. Hurst, PRA **32** (1), 1985

STIMULATED RAMAN SCATTERING

Third order nonlinear process:
stimulated emission from a virtual level



SRS SIGNAL

$$\Delta I_{P, max} \propto \chi^{(3)} L I_S I_P$$

- $\chi^{(3)}$ → H₂ Raman cross-section is high
- L → Multi-pass cell
- I_S → High-power CW laser


SRS signal estimation @ 1 atm for (1←0) Q(1) line:

- 10 mW pump at 740 nm
- 5 W Stokes at 1064 nm
- 30 m multipass cell


$$\frac{\Delta I_{P, max}}{I_P} \approx 5 \cdot 10^{-4}$$

Shot noise $\approx 3 \cdot 10^{-7}$ in 100 ms

SETUP




COMB REFERENCING




- Pump: locked
- Stokes
 - a) Lock and step the LO frequency
 - b) Scan and track the beat note

NOISE ANALYSIS

Allan deviation of the demodulated signal



DATA




SNR:

- 2000 @ 5 atm
- 200 @ 0.5 atm

Type A	< 560 KHz
Comb	40 kHz
Pump BN	300 kHz
Stokes BN tracking	400 kHz
Temperature stability	< 10 kHz
Pressure stability	< 250 kHz
Type B	< 1.53 MHz
Pressure accuracy	< 1.53 MHz
Temperature accuracy	< 90 kHz
Purity of H ₂	99.999%

RESULTS FROM RAW DATA




Literature:

$-3.4(1) \cdot 10^{-3} \text{ cm}^{-1}/\text{atm}$ [1]



$-3.14(15) \cdot 10^{-3} \text{ cm}^{-1}/\text{atm}$ [2]

[1] P. Wcislo et al., JQSRT **177**, 75 (2016)

[2] Rahn *et al.*, PRA **43** (11), 6075 (1991);



MULTISPECTRUM FIT



	γ_0 [$10^{-4} \text{ cm}^{-1}/\text{atm}$]	δ_0 [$10^{-3} \text{ cm}^{-1}/\text{atm}$]	γ_2 [$\text{cm}^{-1}/\text{atm}$]	δ_2 [$\text{cm}^{-1}/\text{atm}$]	ν_{VC} [$\text{cm}^{-1}/\text{atm}$]	η	γ_{D} [10^{-2} cm^{-1}]
Multi-line	7.3(5)	-2.83(2)	0(0.0004)	0.0018(5)	0.04499 [3]	0	1.82
Lit.	8.7(2) [1]	-3.14(15) [2]	$3.3(1)10^{-4}$ [3]	0.0023(1) [3]			
Single line	80(90)	-	0.014(13)	-0.0015(5)	0.021(12)		

[1] May *et al.*, Can. J. Phys. **39**, 1769 (1961); [2] Rahn *et al.*, PRA **43** (11), 6075 (1991); [3] Wcisło *et al.* JQRT **177**, 75 (2016)

CONCLUSIONS & OUTLOOK

- SRS + comb is a powerful approach for H₂ spectroscopy
- A better lineshape model is under preparation
- Fixing the laser intensity noise with active stabilization
< 100 kHz frequency accuracy, SNR > 10³ between 0.5 and 5 atm
- Target pure rotational lines

THANK YOU!

M. Lamperti, L. Rutkowski, D. Gatti, B. Alsaif, R. Gotti, G. Cerullo, D. Polli, M. Marangoni,
F. Thibault, P. Maslowski, S. Wojtewicz, P. Wcislo



POLITECNICO
MILANO 1863

