UNIVERSITY OF TARTU

FACULTY OF SCIENCE AND TECHNOLOGY

Quantitative analysis of Mitochondrial DNA (mtDNA) copy number in yeast S. cerevisiae strain W303 α

MASTER THESIS

Yingjian Hou

Supervisors: Juhan Sedman

1. Li	ist of abbreviations	2
2. In	troduction	2
	iterature overview	4
3.1	Basic concept of yeast	4
3.2	Structure of yeast mitochondrial DNA	4
3.3	Theory of DNA amplification	5
3.3	3.1. Principle of design primer	
3.4	Theory of Plasmid DNA	7
3.4	1.1. Principle of recombination plasmid DNA	
3.4	1.2. The type of restriction end	8
3.4	1.3. Plasmid Isolation	8
3.5	Suggestion for plasmid DNA purification	9
3.6	Suggestion for plasmid copy number calculation	10
3.7	Description of used analytical methods	10
3.7	7.1. Gel Electrophoresis	10
	7.2. UV-Vis Spectrophotometer	11
3.7	7.3. Real-time polymerase chain reaction	11
4. E	xperiment part	13
4.1	Reference plasmid preparation	13
4.2	Reference plasmid copy number calculation and uncertainty estimate	16
4.3	Total DNA preparation	19
4.4	Real-time polymerase chain reaction experiment	20
5. Re	esults and discussion	22
5.1	Absolute quantitative analysis of 15S rRNA	22
5.2	Absolute quantitative analysis of Cox3	27
5.3	Absolute quantitative analysis of COX1-exon1	28
5.4	Absolute quantitative analysis of Cox1_exon4	29
5.5	Absolute quantitative analysis of ATP9	31
5.6	Absolute quantitative analysis of Cob4	32
5.7	Absolute quantitative analysis of SceI	33
5.8	Absolute quantitative analysis of FIT2	34

6. Summary	<i>37</i>
7. Summary in Estonian	38
8. Acknowledgments	38
9. Reference	38
10. Appendices	
Appendice1. ANOVA statistical analysis result	41
Appendice2. qPCR amplify curve and melting curve_	44
Appendice3. DNA extraction protocol	61
Appendice4. Calculation example	62
Appendice5. Determination uncertainty of Pipette	63

1. List of abbreviations

mtDNA - Mitochonarial DNA

qPCR - real-time polymerase chain reaction

S.cervisiae - Saccharomyces cerevisiae

C. albicans - Candida albicans

CsCl - cesium chloride

(EtBr) - ethidium bromide

Ct – threshold cycles

CN – copy number

ddH₂O - Deionized distilled water

EtOH - ethanol

KoAC - Potassium acetate

MgCl2 - Magnesium chloride

UV-VIS – Ultraviolet-visible spectroscopy

i-PrOH - iso-propanol

2. Introduction

Mitochondrion is a double phospholipids bilayer membrane-enclosed organelle

found in most eukaryotic cells, contains DNA, and carries out oxidative phosphorylation, The mitochondria are the sites for much of the metabolism necessary for the production of enrgy in the form of ATP^{1,2}.

Mitochondrial DNA (mtDNA) copy number (molecules or genomes per cell) is believed to be relatively constant. Early reports demonstrated that the control of mtDNA synthesis is independent of the mechanisms controlling mtDNA replication during S phase of the cell cycle. Mitochondrial DNA encodes for a number of respiratory chain peptides and RNA molecules for protein biosynthesis in mitochondria. Therefore, stability of mtDNA copy number is essential to produce fully functional cells.

The regulation of mtDNA synthesis, the defects in oxdative phosphorylation enzymatic machinery during human disease and different fundamentals questions concerning DNA cis acting elements can only be understood when the data on correct copy number of mtDNA will be available³.

mtDNA copy number has only been measured in a population of cells and quantitative mtDNA copy number in single cell had not been published. The budding yeast *Saccharomyces cerevisiae* which is studied as organism model in our lab has been widely used to analyze the mechanisms of mtDNA replication and stability, There are few reports on proteins that could play a role in mtDNA copy number control, mtDNA binding protein Abf2p has been proposed to be a key character in these processes (Zelenaya-Troitskaya et al., 1998); Pif1p is a conserved DNA helicase which occurs both in the nucleus and mitochondria in S.cerevisiae (Foury and Lahaye, 1987; Schulz and Zakian, 1994), and previous research postulated that Pif1p plays a role in the repair or tolerance of oxidative mtDNA damage, perhaps by governing the rate of mtDNA replication or regulating mtDNA copy number⁴. However, the fundament of mechanism for maintaining mtDNA copy number has remained a mystery.

The main goal for this thesis was to develop reliable quantitative analysis method to measure mtDNA copy number in yeast *S. cerevisae* cells. qPCR based system was next used to ask questions that have biological relevance. Gene copy number for

different mitochondrial genomic loci was analyzed. We also discuss the possibilities to use the method for single cell analysis.

3. Literature overview

3.1 Basic concept of yeast

Yeast Serves as a Minimal Model Eukaryote

In order to analyze the internal workings of the eukaryotic cell without the additional problems of multi-cellular development, we usually use a species that is unicellular and as simple as possible. Therefore, the yeast species is a popular choice for this role. *S. cerevisiae* is a small, single-celled member of the kingdom of fungi (Figue 3.1) and thus, according to modern views, at least as closely related to animals as it is to plants⁵.




Figure 3.1 The yeast *Saccharomyces cerevisiae*. (A) A scanning electron micrograph of a cluster of the cells (B) A transmission electron micrograph of a cross section of a yeast cell, showing its nucleus, mitochondrion and thick cell wall. (Alberts. Molecular Biology Of The Cell.5th.Ed)

3.2 Structure of yeast mitochondrial DNA


The mtDNA of strain FY1679, an isogenic derivative S288C, is a 85779 kb double-stranded molecule, which contains the following genes: 2 rRNA genes (15s rRNA and 21s rRNA), cytochrome c oxidase subunit I, II and III (cox1, cox2 and

cox3), ATP synthase subunits 6, 8 and 9 (atp6, atp8, and atp9), apocytochrome b (cytb), a ribosomal protein (var1) and several intron-related open reading frames (ORFs)^{6,7,8}.

Yeast mtDNA has often been considered as a circle form during past 40 years. However, this view is probably not correct and instead of circular molecules the mtDNA consists of a phage T4 like network, our previous research found that the *C. albicans* mtDNA forms a complex branched network that does not contain detectable amounts of circular molecules by using conventional and pulsed-field gel electrophoresis, two-dimensional gel electrophoresis and DNA sequence analysis⁹. The aerobic *C. albicans* is a suitable model organism to study mtDNA maintenance in yeast, Unfortunately there are currently no sufficientdata on mtDNA isoforms in *S. cerevisiae* wild type cells

3.3 Theory of DNA amplification


Genes can be selectively amplified by polymerase chain reaction (PCR). Conventionally, the PCR amplification reaction is described as at first stochastic, then exponential and finally stagnant¹⁰. The PCR method is extremely sensitive, and even can detect the single DNA molecule.



Scheme 3.2 First cycle of PCR amplification process(Molecular Biology Of The Cell.5th.Ed)

In the first cycle, double-stranded DNA is separated to single strands by heating. Then the pair of DNA oligonucleotides which are chosen to flank the desired nucleotide sequence of the gene are used to prime DNA synthesis on single strands.

The newly synthesize double-stranded DNA is produced in a reaction catalyzed in vitro by a purified DNA polymerase⁵.



Scheme 3.3 PCR amplification overview (Molecular Biology Of The Cell.5th)

In principle, every cycle doubles the amount of DNA synthesized in the previous cycle. After strand separation, the two primer DNA oligonucleotides allowed to hybridize to complementary sequences in the two DNA strands by cooling the reaction mixture. The synthesize direction is from 3'-OH terminus to 5'terminus on template strand DNA.

3.3.1. Principle of design primer

Good primer design is essential for successful PCR reactions. The important design has to follow several key principles which are described below.

Primer length: the common length of PCR primers is 18-27bp. It is not suitable to use primers longer than 38bp. Too long primers can lead annealing temperatures higher than 74°C and this is not suitable for TAQ DNA polymerase reaction.

Primer annealing temperature: The primer melting temperature is the estimate of DNA hybrid stability and is a critical parameter of the system. Too high annealing temperature will produce insufficient primer template hybridization, leading low PCR product yield. Too low temperature may lead to non-specific products. Therefore, it is necessary to run melting curve within qPCR experiment and afterwards possibly to check the product by gel electrophoresis.

GC content: the optimal content should be 40% - 60%. Furthermore, there should be no large melting temperature differences between the upstream and the downstream primer.

Avoid non-specificity amplification: It is necessary to avoid regions of homology, in order to improve specificity of primers. Therefore primers designed for a particular sequence must not amplify other genes in the target DNA sequence. Commonly, we tested the primers by BLAST on Saccharomyces GENOME DATABASE (http://www.yeastgenome.org/).

3.4 Theory of Plasmid DNA

3.4.1. Principle of recombination plasmid DNA

Plasmids are widely used as vectors to carry the genes into bacterial cells, Plasimds are small circular molecules of double-stranded DNA. In this research work, the pGEM-7Zf which produced by Promega Corporation (USA) has been used as cloning vectors to insert the standard gene fragment that used as reference material for create of the calibration curve. To use as a vector, the purified plasmid DNA has to be cleaved with restriction nuclease to create a liner DNA molecule. The insert gene fragment of interest was amplified with PCR and primers with added the same restriction nuclease sites and after cleavage with restriction enzymes the fragments of the gene of the interest to be cloned are then added into the liner plasmid DNA by annealingvia their cohesive ends to form recombinant DNA circles. These recombinant DNA containing interest fragment inserts are then covalently linked with enzyme DNA ligase⁵. The Figure 3.4 demonstrates the process of introducing foreign DNA fragment into plasmid DNA.




Figure 3.4 the process of introducing foreign DNA fragment into plasmid DNA (Molecular Biology Of The Cell.5th)

3.4.2. The type of restriction end

The recombinant plasmid DNA was transformed into bacterial cells (E.coli), in this work competent strain DH5 alpha competent cell were used as transform host cells to obtain colonies for screening the .the recombinant plasmids and for large-scale plasmid purification. Traditional methods to introduce DNA into bacterial cells involve incubating bacterial cells in calcium salt solution that makes the cell membrane leaky, permeable to the plasmid DNA

3.4.3. Plasmid Isolation

To use plasmid DNA which contains a target gene as a reference material, high purity plasmid must be obtained. Generally, in order to obtain plasmid DNA, the cells are harvested, the bacteria are lysed, and plasmid DNA recovered by alcohol precipitation. The best method for further DNA purificationis cesium chloride (CsCl)/ethidium bromide (EtBr) centrifugation¹¹. Although many techniques have been invented for plasmid isolation, Cesium Chloride (CsCl) Density Gradient Centrifugation technique has always being considered as the best method to obtain the highest purity plasmid¹². Therefore Cesium Chloride (CsCl) Density Gradient

Centrifugation technique was applied in this work.

CsCl is highly soluble in water, under high centrifugal force, CsCl molecules will dissociate and the heavy Cs⁺ atoms will be forced towards the outer end of the tube, thus forming a shallow density gradient. DNA molecules placed in this gradient will migrate to the point where they have the same density as the gradient (the neutral buoyancy or isopycnic point). The gradient is sufficient to separate different types of DNA with slight difference in density due to differing G+C content, different isotopic content or due to different EtBr binding. The latter is the basis for separation of plasmid and chromosomal DNA, or physical form. In this work, after 15hours centrifugagation, plasmid DNA separated from chromosomal DNA, was detected under the UV light as a separate heavier band in the CsCl density gradient.




Figure 3.8 The procedure of CsCl Density Gradient Centrifugation. (A) position comparison between initial and final procedure; (B) the position of plasmid DNA after centrifugation.

3.5 Suggestion for plasmid DNA purification

Impure plasmid DNA may impact qPCR amplification, leading inaccurate initial copy number estimation in the reference plasmid material. Hence,in this work the plasmids were—treated with 2mg/mL RNase and purified on CsCl Density Gradient Centrifugation to remove large size—RNA. However, there was trace amount of RNA in plasmid DNA preparations. Therefore urther purification of DNA was performed with gel filtration.

The technique of gel filtration chromatography is used to separate molecules of different sizes and shapes¹³. The process of gel filtration involves letting the mixture of RNA and DNA be carried by an elution buffer through a filtration matrix contained in a glass column. The matrix consists of microscopic beads which are porous¹³. Compact molecules which used in this work were Sephacryl S-400 HR tend to get trapped in the pores for awhile before passing through the matrix¹³. Therefore, a heavier molecule will pass through the matrix faster than a lighter one, meaning the plasmid DNA should pass through the matrix faster than RNA. Hence, plasmid DNA allows to be collect and RNA stays remain in matrix after forcing brief micro-centrifugation.

3.6 Suggestion for plasmid copy number calculation

As the molecular weight of the plasmid and insert gene are know, it is possible to calculate the copy number as follow¹⁴:

$$copy number = \frac{6.02 \times 10^{23} (copy/mol) \times DNA \text{ amount(g)}}{DNA \text{ length(bp)} \times 660 (g/mol/bp)}$$
(1)

(where: bp = base pairs, weight in Daltons(g/mol)=(bp size of dsDNA product)(330 Da \times 2nt/bp, ds = double-stranded, nt = nucleotides)

3.7 Description of used analytical methods

3.7.1.Gel Electrophoresis

Gel electrophoresis separates molecules according to their size and shape, but also according to their charge (JAC Biology NYA Laboratory Outlines, 2001). The samples to be analyzed are placed into a gel bed and when a voltage is applied across the gel bed, the molecules will migrate towards the pole they are attracted to according to their charge. In the case of DNA and RNA the charge-to-mass ratio on the two types of nucleic acids is not a factor, because there are always two negatively charged phosphate groups for each base pair charge (JAC Biology NYA Laboratory Outlines, 2001). Both types of nucleic acid will migrate towards the positive pole, however the RNA will move faster through the gel bed because it is smaller¹⁵.In this

work, gel electrophoresis method mainly used to identify specific DNA fragments.

3.7.2. UV-Vis Spectrophotometer

Spectrophotometer is the instrument for quantitative measurement of the reflection or transmission properties of a material as a function of wavelength¹⁶. The most common spectrophotometers are used in the UV and visible region of the spectrum.

Analysis of nucleic acids is regularly performed to determine the average concentrations of DNA or RNA present in a mixture, as well as their purity. The Beer-Lambert law is used to determine unknown concentrations without the need for standard curves. In essence, the Beer Lambert Law makes it possible to relate the amount of light absorbed to the concentration of the absorbing molecule.

Usually A260 (ultraviolet light absorption at the wavelength of 260 nm) is used as a quantitative measure for nucleic acids.

The secondary benefit of using spectrophotometric analysis for nucleic acid quantitation is the possibility to determine sample purity using the 260 nm/280 nm ratio calculation. A pure DNA sample will yield an A260/280 of approximately 1.8. The pure RNA should have an A260/280 ratio of approximately 2.0. These ratios are commonly used to assess the amount of protein contamination that is left from the nucleic acid isolation process since proteins absorb at 280 nm¹⁷.

The spectrophotometer which used in this study was "NanoDrop 2000" produced by Thermo scientific company.

3.7.3. Real-time polymerase chain reaction

Real-time polymerase chain reaction also called real-time quantitative PCR (qPCR) offers fast and reliable quantification of any target sequence in a DNA sample. The quantity can be either an absolute copy number or relative amount of target DNA when the housekeeper gene is being used as a reference in mRNA quantifications performed on cDNA¹⁸. Nowadays, many qPCR assays have been widely used in several applications such as environmental, clinical, forensic analysis. Real-time qPCR assays are being used in those applications due to high levels of sensitivity, specificity and precision. The technique allows for the continuous collection

fluorescence signal of PCR products during the reaction process. Ideally, these products will be exponentially amplidied, which means their quantities double with each thermal reactions cycles. Hence, qPCR can be applied to determine a threshold where the accumulation of amplified product is first significantly visible in the data. The fractional cycle number where PCR product accumulation passes this fixed threshold is called the threshold cycles (Ct)¹⁹.A basic PCR run can be broken up into three phases: Exponential, Linear and Plateau.





Figure 3.8 PCR phase (source: http://www.invitrogen.com)

At the linear phase, some of reagents are being consumed up as a result of amplification. The reactions start to slow down and the PCR product is no longer being doubled at each cycle, which means this phase generates high variable copy number. By contrast, the reaction in the exponential phase is very specific and precise, because all of the reagents are fresh and available, the kinetics of the reaction lead the reaction to favor doubling of the amplicon. Therefore, the noise baseline has to be set up as lower as possible, the threshold line should cross the exponential phase.



In an Absolute Quantification analysis, a standard curve is used to determine the concentration of unknown samples. In a standard curve, the concentrations of standard samples are plotted against the Ct of the samples. There is a log-linear relationship between the initial amount of DNA target in the reaction and the Ct value that is obtained.

Real-time PCR employs fluorescent dyes or probes that interact with the PCR products. The two primary types of fluorescent detection are DNA binding dyes, such as SYBR Green, or fluorescently tagged sequence specific probes, such as TaqMan or Molecular Beacon probes²⁰. SYBR Green dye was used in this study. SYBR Green is a dye that binds the Minor Groove of double stranded DNA. When SYBR Green dye binds to double stranded DNA, the intensity of the fluorescent emissions increases. As more double stranded amplicons are produced, SYBR Green dye signal will increase. SYBR Green dye will bind to any double stranded DNA molecule, Since detection of fluorescent signal from these dyes is not sequence specific, melting temperature analysis must be performed to ensure the production of a single PCR product²¹.

4. Experiment part

4.1 Reference plasmid preparation

Plasmid which contains specific sequence of analyte was used in qPCR

measurement to set up the calibration curve to calculate sample's copy number.

pGEM-7Zf plasmid vector (Promega, Madison, Wisconsin, USA) was used in this study. Specific fragments of yeast mitochondrial DNA were amplified with PCR. The oligonucleotides used for amplification were also used to introduce restriction sites for cloing at the fragment end. The fragments were purified from agarose gel by UltraClean®15 DNA purification kit (MO BIO, Carlsbad, CA USA). Fragments of 15S rRNA, COX3, COX1-exon1, COX1-exon4, SceI and ATP9 genes were inserted between the restriction sites XbaI and EcoRI. A fragments of COB4 was inserted between XbaI and SmaI, SmaI was used to create a blunt ends, because the inserter COB4 fragment sequence contains a restriction site EcoRI (comfirmed by WebCutter program). Therefore COB4 could not be cloned by using EcoRI.

The list of oligonucleotides used for amplification of mitochondrial DNA fragments is presented in the Table 4.1.

	Table 4.1 Oligonucleotides for amplification of		
	mitochondrial DNA fragments		
		Product	Unavaila
Name	Sequences $(5' \rightarrow 3')$	size	ble
		(bp)	enzyme
15S rRNA	L-Xba: CTATCTAGA TCAAGCCAATAATGGTTTAGGT	169	
135 IKNA	R-EcoRI: CTGAATTC AACTGGATCAATCTTTCGATCA	109	
COX3	L-Xba: CTATCTAGA TCCATTCAGCTATGAGTCCTGA	- 281	
COAS	R-EcoRI: CTGAATTC TGAACCATAAACACCATCAGAGA	201	
COV11	L-Xba: CTATCTAGA TGATTATATTCAACAAATGCAAAAGA	225	
COX1-exon1	R-EcoRI: CTGAATTC CCATGTAACACTTAACTCCACCT	325	
GOV4 4	L-Xba: CTATCTAGA CTACACGTGTTGCACCCATT	266	
COX1-exon4	R-EcoRI: CTGAATTC GGTCCTGAATGTGCCTGAAT	266	
A TEDO	L-Xba: CTATCTAGA GCAGCTAAATATATTGGAGCAGGTA	102	
ATP9	R-EcoRI: CTGAATTC TGAAACCATTAAACAGAATAAACCTG	192	
СОВ4	L-Xba: CTATCTAGA AAATAGGATTAAGATATAGTCCGAACA	246	EcoRI

	R-EcoRI: CTGAATTC TGAATGCATTGGAATTCTATCTAAA		
	L-Xba: CTATCTAGA ATTTACCCCCTTGTCCCATT		
SceI	R-EcoRI: CTGAATTC CCTCCATCATCTATAAATCAATATGC	736	HindIII

The list of constructed plasmid is presented in Table 4.2.

	Table 4.2 List of constructed plasmid											
Name	Upstream site	Downstream site	Plasmid size(bp)									
15S rRNA	Xbal	EcoRI	3166									
COX3	Xbal	EcoRI	3278									
COX1-exon1	Xbal	EcoRI	3322									
COX1-exon4	Xbal	EcoRI	3263									
АТР9	Xbal	EcoRI	3189									
COB4	Xbal	SmaI	3230									
SceI	Xbal	EcoRI	3733									
	(Plasmid size = pGEM7/Zf size + inserter size)											

The reference plasmids were purified from 100-200ml bacterial cultures by alkaline lysis, Ribonouclease A treatment (20 mkg RNAse A in 1 ml, 56 °C, for 30 min) and CsCl Density Gradient Centrifugation. The traces of low-molecular RNA were removed by additional purification step on Sephacryl S-400 HR gel filtration column. The plasmid preparations were inspected by agarose gel electrophoreses.




Figure 4.1 Inspection of plasmid purification of plasimd by agarose gel electrophoreses. Lane 1,3,5,7,9,11,13 corresponds SceI,cox3,cox1-1,ATP9,Cob4,15s,cox1-4 after gel filtration treatment, respectively. Column 2,4,6,8,10,12,14 corresponds SceI,cox3,cox1-1,ATP9,Cob4,15s,cox1-4 before gel filtration, respectively. The result demonstrated that there was no rsidual RNA remains in the Plasmid preparations.


Chromosome reference gene was FIT2. FIT2 was amplified by PCR, and then fragment purified from agarose gel by UltraClean®15 DNA purification kit (MO BIO, Carlsbad, CA USA).

The stocks of reference plasmid were stored in $T_{10}E_{0.1}$ at +4 °C.

4.2 Reference plasmid copy number calculation and uncertainty estimate

The copy number calculation formula was given by (1). In order to calculate the copy number of plasmid, concentration of plasmid was measured by Nanodrop 2000 spectrophotometer by 5 replications. The uncertainty source flow is presented in Figure 4.2. If the DNA is not of sufficient quality, the downstream PCR efficiency may be affected. The quality of DNA can be affected by its size, stability, level of degradation, structural integrity, and presence of inhibitors and modifications, which may be inherent in the sample matrix²²

Fig4.2 uncertainty source flow



It is important to ensure that the sample was mixed well before quantification, and that the spectrophotometer has been cleaned to ensure no sample carryover between redings. Based on manufacturer's guidelines, a $2\mu l$ drop is used for t measurements. The blank sample was $T_{10}E_{0.1}$ which was also used to dissolve plasmid DNA.

The plasmid DNA was diluted to a proper concentration, keep in $T_{10}E_{0.1}$ at +4 °C for routine use.

The series dilution process is presented in Table 4.3

	Table 4.3 Series dilution process														
Name	Initial		Dilution No.1			Dilution	No.2	Final							
	concentration	on	con	centration		concentr	ation	concentration							
158	147.82ng/μl	5μl	68.91µlTE	10ng/μl	<i>Σ</i> μΙ	45µlTE		1ng/μl							
COX1-1	970.18 ng/μl	Δ 5μl	43.509μlTE	100 ng/μl	5μl	45μlTE	5μl	45μlTE	lng/μl						

			Z		\mathcal{A}	Z	\mathcal{A}	<i>2</i> /	
COX1-4	2394.16 ng/μl	Δ 5μl	114.71µlTE	100 ng/μl	Δ 5μl	45μlTE	- Δ 5μl	45μlTE	1ng/μl
		7							
COX3	1705.12 ng/μl	5μl	80.25µlTE	100 ng/μl	5μl	45μlTE	- 5μl	45µlTE	1ng/μl
		7	\mathcal{Y}		(1	\mathcal{Y}		\mathcal{I}	
ATP9	1474.4 ng/μl	5μl	68.72μlTE	100 ng/μl	5μl	45μlTE	5μl	45μlTE	1ng/μl
		(1			(1		(1		
СОВ4	2473.5 ng/μl	<i>Δ</i> 5μl	118.68µlTE	100 ng/μl	 5μl	45μlTE	 5μl	45μlTE	1ng/μl
								7	
SCEI	2173.14 ng/μl	5μl	103.66µlTE	100 ng/μl	5μl	45μlTE	5μl	45µlTE	1ng/μl
			<i></i>			7			
FIT2	18.52 ng/μl	10µl	8.52μlTE	10 ng/μl	Δ 5μl	45µlTE			1ng/μl

The pipette was set at the proper volume to calculate dilution uncertainty by 5 replications. The measurement data is presented in Appendices Table 4.

Uncertainty is associated with dilution procedures prior the calibrator using in the work. The exactly uncertainty of the dilution factor depends on the series dilution process. The numerator and denominator in the formula for calculating the dilution factor are correlated (same a). Therefore, it is not completely correct to apply the rules of combining independent components of uncertainty when estimating the uncertainty of a dilution factor²³. The uncertainty variance of a dilution step is obtained from:

$$u_f^2 = \frac{a^2 u_b^2 + b^2 u_a^2}{a^4} = \frac{u_b^2}{a^2} + \frac{b^2 u_a^2}{a^4}$$
 (2)

And the RSD squared from:

$$w_f^2 = \frac{1}{(a+b)^2} \left[\frac{a^2 u_b^2 + b^2 u_a^2}{a^2} \right] = \frac{b^2}{(a+b)^2} (w_b^2 + w_a^2)$$
 (3)

a = suspension transfer volume

b = diluton blank volume

 u_a = standard uncertainty of a

 u_b = standard uncertainty of b

 w_a = relative standard uncertainty of α

 w_b = relative standard uncertainty of b

If the total dilution contains of k similar steps the combined RSD squared of the dilution factor is calculated from:

$$w_F^2 = k w_f^2 \tag{4}$$


If the total dilution contains several different steps, like each time different volume to perform dilution, the relative uncertainty of each step should be separately estimated, the combined RSD squared of the dilution factor is calculated from:

$$w_F^2 = w_{f_1}^2 + w_{f_2}^2 + \dots + w_{f_k}^2 \tag{5}$$

4.3 Total DNA preparation

The cultures of yeast *S.cerevisiae* strain W303 α were grown in 12ml YPG at 25 °C to OD600 = 0.70 (logarithmic growth phase), then transformed to 70ml YPG and grown to OD600 = 1.02, 10ml of the culture was spun down at that moment and used for DNA preparation.YPG contains yeast extract, peptone and glycerol(C₃H₈O₃).

DNA was extracted with the following protocol that minimizes damage to high molecular weight molecules. And then keep DNA at -80 °C in small alignotes. Once DNA was resuspended in the $T_{10}E_{0.1}$ (PH7.5), tips with cut end were used to handle DNA solution. The concentration of total DNA was 147.8ng/µl. Then DNA was diluted to 1000 and 10,000 times, separated in 6 tubes (15µl per tube). The protocol details are presented in Appendices 3.



4.4 Real-time polymerase chain reaction experiment

The calibration curves consisted each 5 calibration points: 10⁵, 10⁴, 10³, 10², 10¹ copies of target gene per PCR assay. The copies of calibration plasmid were calculated based on formula (1). qPCR measurements were performed in MicroAmp® Optical 96-Well Reaction plates (Applied Biosystems, Lennik, Belgium) on a Roche LightCycler® 480 Real-Time PCR System using DNA binding dyes SYBR Green PCR MasterMix. Briefly, PCR cycling conditions were as follows: activation of DNA polymerase and initial denaturation for 15 mins at 95 °C; 50-60 cycles of denaturation for 10 s at 95 °C, annealing for 10s at 53 °C and extension for 15s at 72 °C, followed by a melting curve program, 1s quick denaturation at 95 °C, cooling down to 53 °C for 30s, followed by continued detection at 95 °C. All PCR assays were performed in triplicate using 10μl reaction system.

	1	2	3	4	5	6	7	8	9	10	11	12
A	A1	A1	A1	B1	B1	B1	C1	C1	C1	D1	D1	D1
В	A2	A2	A2	B2	B2	B2	C2	C2	C2	D2	D2	D2
C	A3	A3	A3	В3	В3	В3	C3	C3	C3	D3	D3	D3
D	A4	A4	A4	B4	B4	B4	C4	C4	C4	D4	D4	D4
E	A5	A5	A5	B5	B5	B5	C5	C5	C5	D5	D5	D5
F	X	X	X	X	X	X	X	X	X	X	X	X
G	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
H	N	N	N	N	N	N	N	N	N	N	N	N

Fig. 4.4 Schematic representation of the PCR plate design used in this study. The wells are marked with the following numeric coding for the respective series dilution of plasmid DNA (A = plasmid A, B = plasmid B, C = plasmid C, D = plasmid D). The rectangles represent triplicates for each point of the different calibration curves. The wells labelled X contain the genomic DNA samples, Y was stochastic variables, ie, 10x dilution of previously loaded genomic DNA, spiked plasmid DNA etc. The wells were marked N is negative control point (there was no DNA in it).


	qPCR reaction components	
Chemical	volume per reaction (μl)	final concentration
10X homemade PCl	R 1	1x
buffer		
2mM dNTPs	1	0.2mM
20μM primers	0.3 + 0.3	600nM
20mM MgCl ₂	0/0.25-0.5	3-4mM
1:10,000 SYBR Green	1	1/10 of reaction volume
dye		
Genomic DNA	1	
Hot TAQ	0.1	1/10 of reaction volume
Deionized distilled wat	er 4.8-5.3	

- ◆ 10x Homemade Core PCR Buffer consists 30mM MgCl₂, 200mM KCl,100mM Tris HCl (pH8.5), 1% Triton X-100
 - ◆ MgCl₂ concentration was fluctuated for different genes, decision made by

reaction efficiency, As a rule, higher MgCl₂ concentration can significantly raise efficiency, but can also lead to non-specificity amplification²⁴.

10 fold series plasmid dilution made from $1 \text{ng/}\mu\text{l}$ stock solution to $0.0001 \text{pg/}\mu\text{l}$ by a standard scheme (5 μ l DNA solution mix with 45 μ l ddwater).

Fig. 4.5 A example for mixing reaction components



★ The loading sample for each individual well was 1µl.

Suggestion for avoiding contamination

- \triangle Keep your stock solution in the safe place
- \triangle Use new tips for each experiment
- △ Prepare several gloves for experiment, change them frequently
- △ Keep working environment sterile

5. Results and discussion

5.1 Absolute quantitative analysis of 15S rRNA

Uncertainty estimate

15S rRNA copy number was measured 4 times at 2 different days. The concentration of stock plasmid calibrator was 147.82ng/µl with 0.25% relative

uncertainty. According to dilution process, following the formula (1), (3) and (5), the 5 series 10 fold calibration points had 1.24%, 1.38%, 1.51%, 1.63% and 1.74% relative uncertainty, respectively.

However, the concentration or copy number of the calibrator does not correspond to the Ct value, but the log₁₀(copy number) was used to establish the calibration curve with Ct. The uncertainty of log₁₀(copy number) was calculated from formula:

$$f_{(x)} = \log_{10}(CN) \tag{6}$$

$$u_{f_{(x)}} = \sqrt{\left(\frac{\partial_{f_{(x)}}}{\partial_{CN}} \times u_{CN}\right)^2} \tag{7}$$

$$\frac{\partial_{f(x)}}{\partial_{CN}} = \frac{1}{x \times In10} \tag{8}$$

CN = copy number

 $u_{f_{(x)}} = \text{uncertainty of } \log_{10} \text{copy number}$

Hence, the uncertainty of log₁₀copy number of plasmid standard DNA was 0.0054, 0.0060, 0.0066, 0.0071, and 0.0076, respectively.

The uncertainty of the Ct was observed by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in logarithmic form; uncertainty of this result was estimated by Kragten approach. The formula is presented as following:

$$\log_{10}(CN) = \frac{Ct - b}{a} \tag{9}$$

CN = copy number

Ct = threshold cross value

b = the intercept of the calibration function

a = the slope of the calibration function

The sample initial copy number and its uncertainty were calculated by following formula:

$$CN_{IN-run} = f_{(CI)} = 10^{\frac{CI-b}{a}}$$

$$CN_{IN-run} = 10^{\log_{10}(CN)}$$

$$STDVEV_{In-run} = \sqrt{\left(\frac{\partial f_{(y)}}{\partial_{\log_{10}(CN)}} \times STDVEV_{\log_{10}(CN)}\right)^{2}}$$

$$\frac{\partial f_{(y)}}{\partial_{\log_{10}(CN)}} = In 10 \times 10^{\log_{10}(CN)}$$

$$(10)$$

The recovery factor was measured by spiking same amount of plasmid DNA. The result was 1.03, which was close to 1, and considering this it is not necessary to correct the result.

The the copy number was converted to its undiluted a smple copy number by correcting with the dilution factor:

$$CN_{undiluted} = CN_{IN-run} \times d$$

$$u_{CN_{undiluted}} = CN_{undiluted} \times \sqrt{\left(\frac{u_{CN_{IN-run}}}{CN_{IN-run}}\right)^2 + \left(\frac{u_d}{d}\right)^2}$$

$$d = dilution \ factor$$

$$CN_{IN-run} = copy \ number \ of \ loading \ sample$$

$$CN_{undiluted} = copy \ number \ of \ undiluted \ sample$$

$$CN_{undiluted} = copy \ number \ of \ undiluted \ sample$$

The result are presented in following table:

Gene ID	No.	Data	Copy_undiluted	u_rel	R^2	E					
	1	1-1	241851756	11.64%	0.9872	99.51%					
15s RNA	2	1-2	229206668	13.17%	0.9832	96.89%					
155 KNA	3	2	260983230	12.47%	0.9961	97.73%					
	4	3	208288676	27.74%	0.9945	99.21%					
E: efficie	E: efficiency of amplification; $E = 10^{(-1/slope)} - 1$										

Linearity




Fig.5.1 15S rRNA calibration curve comparison in different days. R² was less than 0.99, the linearity should be caution.

For each of the 4 platforms independently, a two-way analysis of variance (ANOVA, "GraphPad Prism 5") was used to test for significant differences. The results from both platforms showed that at the 95% level of confidence, there was a non-significant interaction (P>0.005) between the each calibration samples points and Ct value within different platforms. However, the P value of Ct less than 0.05. Considering it is significant. The calibration curve also was compared by 95% confidence interval as following Fig.5.1.2. There are no extreme result which is out of 95% CI occurred. (The statistical result form is presented in Appendices)




Fig.5.1.2 95% confidence interval test

Specificity

For each run the melting curve detection procedure was applied. The specific gene amplification products are homogeneous molecules that have identical melting temperature, therefore if the melting curve demonstrates only one peak, the specificity can be concluded²⁵.




Fig.5.1.3 The melting curve at day 3. Because the running cycle was set up at 60 cycles, the single gene copy amplification should end before the 40th cycle by using those primers. The negative control samples present specific amplification from contaminating mtDNA traces and the singal. In these samples demonstrates the extreme sensitivity of the method.

5.2 Absolute quantitative analysis of Cox3

Cox3 copy number was measured 2 times at 2 different days. The concentration of stock plasmid calibrator was 1707.12ng/µl with 0.45% relative uncertainty. According to dilution process, the uncertainty of log₁₀copy number of five 10-fold plasmid standard DNA was 0.0062, 0.0068, 0.0073, 0.0077 and 0.0082, respectively.

Uncertainty of the Ct was derived by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in logarithmic form; uncertainty of this result was estimated by Kragten approach, and then transformed to absolute copy number by formula (10), (11), (12).

The recovery factor was measured by spiking same amount of plasmid DNA. The result was 1.01, which was close to 1, and considering this it is not necessary to correct the result.

The result are presented in the following table:

Gene ID	No.	Data	Copy_undiluted	u_rel	R ²	Е
Cav2	1	2	108060821	29.49%	0.9987	89.34%
Cox3	2	5	147951996	8.04%	0.9964	93.26%
E: efficie	ency of a	mplifica	ation; $E = 10^{(-1/slope)}$) -1		




Fig.2 cox3 calibration curve comparison in different days. R² was greater than 0.99, the linearity was acceptable. The graph demonstrated the slopes are not significantly different. If the overall elevations were identical, intercepts are not significantly different. There was a non-significant interaction (P>0.005) between the each calibration samples points and Ct value within different platforms. However, the P value of Ct less than 0.05 and F-ratio was greater. Considering it is significant.

5.3 Absolute quantitative analysis of COX1-exon1

Cox1-exon1 copy number was measured 3 times at 3 different days. The concentration of stock plasmid calibrator was 970.18ng/µl with 0.44% relative uncertainty. According to dilution process, the uncertainty of log₁₀copy number of five 10-fold plasmid standard DNA was 0.0056, 0.0062, 0.0068, 0.0073 and 0.0077, respectively.

Uncertainty of the Ct was derived by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in

logarithmic form; uncertainty of this result was estimated by Kragten approach, and then transformed to absolute copy number by formula (10), (11), (12).

The recovery factor was measured by spiking same amount of plasmid DNA. The result was 1.03, which was close to 1, and considering this it is not necessary to correct the result.

The result are presented in the following table:

Gene ID	No.	Data	Copy_undiluted	u_rel	R^2	E
	1	2	159588837	19.33%	0.9901	94.28%
Cox1_exon1	2	5	143168448	23.60%	0.9982	73.97%
	3	6	85824549	12.27%	0.9906	76.17%
E: efficie	ency of a	mplifica	ation: $E = 10^{(-1/slope)}$) -1		




Fig.5.3 cox1_1 calibration curve comparison in different days. R² was greater than 0.99, the linearity was acceptable. There wasn't a non-significant interaction (P<0.005) between the each calibration samples points and Ct value within different platforms. The P value of Ct less than 0.05 and F-ratio was greater. Considering it is significant.

5.4 Absolute quantitative analysis of Cox1 exon4

Cox1 4 copy number was measured 3 times at 3 different days. The

concentration of stock plasmid calibrator was 2394.2ng/µl with 0.56% relative uncertainty. According to dilution process, the uncertainty of log₁₀copy number of five 10-fold plasmid standard DNA was 0.0064, 0.0070, 0.0074, 0.0079 and 0.0083, respectively.

Uncertainty of the Ct was derived by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in logarithmic form; uncertainty of this result was estimated by Kragten approach, and then transformed to absolute copy number by formula (10), (11), (12).

The result is presented in following table:

Gene ID	No.	Data	Copy_undiluted	u_rel	R²	E	
	1	2	155989983	20.39%	0.9826	102.23%	
Cox1_exon4	2	5	140773979	18.28%	0.9986	100.65%	
	3	6	96891310	8.31%	0.9978	90.10%	
E: efficiency of amplification; $E = 10^{(-1/slope)} - 1$							




Fig.5.4 cox1_4 calibration curve comparison in different days. R² was greater than 0.99, the linearity was acceptable. There wasn't a non-significant interaction (P<0.005) between the each calibration samples points and Ct value within different platforms. The P value of Ct less than 0.05 and F-ratio was greater. Considering it is significant.

5.5 Absolute quantitative analysis of ATP9

ATP9 copy number was measured 3 times at 2 different days. The concentration of stock plasmid calibrator was 1474.4ng/μl with 0.27% relative uncertainty. According to dilution process, the uncertainty of log₁₀copy number of five 10-fold plasmid standard DNA was 0.0060, 0.0066, 0.0071, 0.0076 and 0.0080, respectively.

Uncertainty of the Ct was derived by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in logarithmic form; uncertainty of this result was estimated by Kragten approach, and then transformed to absolute copy number by formula (10), (11), (12).

The result are presented in the following table:

Gene ID	No.	Data	Copy_undiluted	u_rel	R^2	E	
	1	2	171712693	29.29%	0.9807	89.38%	
ATP9	2	7	119323099	43.41%	0.9812	83.36%	
	3	7	134358399	30.90%	0.9816	83.20%	
E: efficiency of amplification; $E = 10^{(-1/slope)} - 1$							




Fig. 5.5 ATP9 calibration curve comparison in different days. R^2 was greater than 0.98, the linearity was acceptable. There was a non-significant interaction (P > 0.005) between the each calibration samples points and Ct value within different platforms. The P value of Ct less than 0.05 and F-ratio was greater. Considering it is significant.

5.6 Absolute quantitative analysis of Cob4

Cob4 copy number was measured 3 times at 2 different days. The concentration of stock plasmid calibrator was 2473.5ng/µl with 0.55% relative uncertainty. According to dilution process, the uncertainty of log₁₀copy number of five 10-fold plasmid standard DNA was 0.0063, 0.0069, 0.0074, 0.0070 and 0.0083, respectively.

Uncertainty of the Ct was derived by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in logarithmic form; uncertainty of this result was estimated by Kragten approach, and then transformed to absolute copy number by formula (10), (11), (12).

The result are presented in the following table:

Gene ID	No.	Data	Copy_undiluted	u_rel	R ²	E
	1	2	103537584	23.55%	0.9918	96.50%
Cob4	2	7	178176340	24.42%	0.9983	84.48%
	3	7	123255941	31.56%	0.09981	83.47%
T 10(-1/slang) 1						

E: efficiency of amplification; $E = 10^{(-1/slope)} - 1$




Fig.5.6 cob4 calibration curve comparison in different days. R² was greater than 0.98, the linearity was acceptable. There wasn't a non-significant interaction (P<0.005) between the each calibration samples points and Ct value within different platforms. The P value of Ct less than 0.05 and F-ratio was greater. Considering it is significant.

5.7 Absolute quantitative analysis of SceI

SceI copy number was measured 4 times at 3 different days. The concentration of stock plasmid calibrator was 2173.14ng/µl with 0.28% relative uncertainty. According to dilution process, the uncertainty of log₁₀copy number of five 10-fold plasmid standard DNA was 0.0060, 0.0066, 0.0071, 0.0076 and 0.0080, respectively.

Uncertainty of the Ct was derived by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in logarithmic form; uncertainty of this result was estimated by Kragten approach, and then transformed to absolute copy number by formula (10), (11), (12).

The result are presented in the following table:

Gene ID	No.	Data	Copy_undiluted	u_rel	R ²	E
SceI	1	8	117698439	12.71%	0.9985	60.67%

	2	5	100346614	27.30%	0.9661	69.31%		
	3	7	111995901	9.40%	0.9961	62.50%		
	4	7	96504569	14.43%	0.9976	63.05%		
E: efficie	E: efficiency of amplification; $E = 10^{(-1/slope)} - 1$							




Fig.5.7 Scel calibration curve comparison in different days, At day 5, the R² was 0.96 and associated with higher relative uncertainty, therefore the result was eliminated from final comparison.

5.8 Absolute quantitative analysis of FIT2

FIT2 copy number was measured 5 times at 4 different days. The concentration of stock plasmid calibrator was $18.5 \text{ng/}\mu\text{l}$ with 0.47% relative uncertainty. According to dilution process, the uncertainty of $\log_{10} \text{copy}$ number of five 10-fold plasmid standard DNA was 0.0062, 0.0068, 0.0073, 0.0074 and 0.0082, respectively.

Uncertainty of the Ct was derived by calculating the standard deviation of mean of triplicates reading.

The directive result from the calibration curve was sample initial copy number in logarithmic form; uncertainty of this result was estimated by Kragten approach, and

then transformed to absolute copy number by formula (10), (11), (12).

The result is	presented in	n follo	owing	table:
The result is	presented if	11011	J W 1115	table.

Gene ID	No.	Data	Copy_undiluted	u_rel	R ²	E
	1	2	24360613	21.82%	0.9989	78.21%
FIT2	2	3	21578152	18.50%	0.9923	76.04%
F112	3	5	27712454	17.90%	0.9951	77.26%
	4	8	45093626	9.56%	0.9966	74.00%
	5	8	50287795	11.61%	0.9968	73.67%
Fractionary of amplifications $E = 10^{(-1/slope)} - 1$						




Fig.5.8 cob4 calibration curve comparison in different days. R^2 was greater than 0.99, the linearity was acceptable. There was a non-significant interaction (P > 0.005) between the each calibration samples points and Ct value within different platforms. The P value of Ct less than 0.05 and F-ratio was greater. Considering it is significant.

5.9 mtDNA copy number in single cell

The chromosome gene FIT2 is single copy chromosomal gene, that is present as a single copy in haploid cells during most of the cell cycle. Therefore the ratio of mtDNA gene to FIT2 is approximates to the mitochondrial genome copy number in a single cell.

Each gene was observed several independent copy number result from different days. The result was treated with box and whisker plot (presented below), using the description outlined above where the whiskers encompass 5% to 95% of the range of results. An outlier can be defined as a data point that does not follow the typical distribution of the rest of the data set and can be regarded as an irregular observation²⁷. There was no significant outlier, hence the mean value were able to be a true value as used to compare with FIT2 copy number. The pooled standard deviation was used to calculate uncertainty in this case.

Box and Whiskers Plot





Table 5.9 mtDNA genes ratio				
Gene	Ratio		Relative	
			Uncertainty	
15s	8. 81	土	29. 48%	
cox1-1	4. 67	土	29. 76%	
cox1-4	4. 73	土	28. 63%	
cox3	4. 62	±	27. 62%	
ATP9	5. 12	±	40. 96%	
Cob4	4. 45	±	37. 49%	



6. Summary

The main goal of this work was developing the qPCR method for yeast mitochondrial DNA analysis, and to estimate the ratio of between mitochondrial and nuclear genes.

This study revealed that the mtDNA copy number in yeast S. cerevisiae stain w303 α does not have a constant value—for each mitochondrial genetic elements. The analyzed protein encoding genes, cytochrome c oxidase subunit I, III (cox1_exon1, cox1-exon4 and cox3), had roughly the same copy number. However, the copy number of the mitochondrial RNA gene, 15S, is approximately twofold higher. Hence, based on this study, we can speculate that the mitochondrial gene copy number depends on the gene function and obviously can be regulated independently, according to the metabolic needs of the organelle. It will be therefore interesting to test how the changes of copy number of individual mitochondrial genes during different growth conditions, for example on different carbon sources. It will also be interesting to analyze other yeast strains to confirm the detected copy number difference in other genetic

background.

7. Summary in Estonian

Töö peamiseks eesmärgiks oli qPCR-l baseeruvate mõõtmismeetodite väljaarendamine, mille abil oleks võimalik kvantitatiivselt iseloomustada mitokondriaalse DNA koopiaarvu pärmi *S. cerevisiae* rakkudes.

Töö käigus läbiviidud analüüs näitas, et *S. cerevisiae* tüves w303 on erinevate mitokondriaalsete geenide koopiaarv erinev. Valke kodeerivad tsütokroomi oksidaasi subühiku I (eksonid 1 ja 4) ja III geenide koopiaarv on ligikaudu võrdne, ent mitokondriaalse ribosoomi 15S rRNA geeni koopiaarv on kaks korda kõrgem. Seega on mitokondriaalse DNA koopiaarv sõltuv konkreetselt analüüsitavast geenist ja võib olla individuaalselt reguleeritud vastavalt metabolismi vajadustele. Huvitav oleks analüüsida, kuidas muutub koopiaarv erinevates kasvutingimustes, näiteks erinevaid süsinikuallikaid kasutades ning kuidas sõltub mitokondriaalse DNA koopiaarv tüve tuuma genoomi erisustest.

8. Acknowledgments

Hereby, I would like to express my sincere gratitude to all the persons who have helped during my studies in Estonia.

I want to thank Professor Juhan Sedman for accepting me as a research student in his group and for guidance during these years.

I would like to thank very much PhD Tiina Sedman – for support and limitless help during practical work in the laboratory.

Juhan Sedman is the best supervisor what I ever had.

Also I would like to thank all my collages and friends from biochemistry laboratory.

I want to thank all my classmates for great time tougher during study and leisure time.

Last, but not the least, I would like to thank my family for support, understanding and true motivation during hard time.

9. Reference

¹ Henze K, Martin W; Martin, William (2003). "Evolutionary biology: essence of mitochondria". Nature 426 (6963): 127–8.

Lodish H, Baltimore D, Berk A, Zipursky S.L, Matsudaira P, Darnell J."MOLECULAR CELL

BIOLOGY"Third Edition.

- ³ Copeland, W.C. (2008). Annu. Rev. Med. 59, 131–146.
- ⁴ Taylor S.D, Zhang H, Eaton J.S, Rodeheffer M.S, Lebedeva M.A, O'Rourke T.W, Siede W, and Shadel G.S "The Conserved Mec1/Rad53 Nuclear Checkpoint Pathway Regulates Mitochondrial DNA Copy Number in Saccharomyces cerevisiae". Molecular Biology of the Cell Vol. 16, 3010–3018, June 2005
- ⁵ Alberts B, Johnson A, Lewis J, Raff M, Roberts K, Walter P. "Molecular Biology of The Cell" Fifth Edition.
- ⁶ Tzagoloj, A. and Myers, A.M. (1986) Annu. Rev. Biochem. 55, 249-285
- ⁷ de Zamaroczy, M. and Bernardi, G. (1986) Gene 47, 155-177.
- ⁸ Foury F, Roganti T, Lecrenier N, and Purnelle B "The complete sequence of the mitochonarial genome of Saccharomyces cerevisiae". FEBS Letters 440 (1998) 325-331
- ⁹ Gerhold, J.M., Aun, A., Sedman, T., Joers, P., and Sedman, J. (2010). Mol. Cell 39, this issue, 851–861.
- Wittwer CTm herrmann MG, Moss AA, Rasmussen, RP: Continuous fluorescence monitoring frapid cycle DNA amplification. Biotechniques 1997, 22:130-138.
- Sambrook, Joseph; David Russell (2001). *Molecular Cloning: A Laboratory Manual.* New York: Cold Spring Harbor Laboratory Press.
- 12 http://en.wikipedia.org/wiki/Sticky_and_blunt_ends
- ¹³ Watson, J.D., Gilman, M., Witkovski., Zoller, M. Recombinant DNA, 78
- Weaver R.F "Molecular Biology", Fourth Edition, Chapter 4, Molecular Cloning Methods.
- ¹⁵ Radloff, R., Bauer, W. and Vinograd, J. 1967. Proc. Nat. Acad. Sci.USA 57: 1514-1521
- ¹⁶ Michael R. Green, *Howard Hughes Medical Institute, University of Massachusetts Medical School;* Joseph Sambrook, *Peter MacCallum Cancer Institute, Melbourne, Australia* Molecular Cloning: A Laboratory Manual (Fourth Edition).
- ¹⁷ S. M. Carr and O. M. Griffith (1987), "Biochemical Genetics," 25:385-390
- ¹⁸ JAC Biology NYA Laboratory Outlines, 2001.
- Whelan, J.A., Russel, N.B., Whelan, M.A., 2003. A method for the absolute quantification of cDNA using real time PCR. J. Immunol. Meth. 278, 261–269.
- ²⁰ Kryndushkin DS, Alexandrov IM, Ter-Avanesyan MD, Kushnirov VV (2003). "Yeast [PSI+] prion aggregates are formed by small Sup35 polymers fragmented by Hsp104". Journal of Biological Chemistry 278 (49): 49636–43.
- Allen, D., Cooksey, C., & Tsai, B. (2010, October 5). Spectrophotometry.
- Glasel, J.A. (1995) Validity of Nucleic Acid Purities Monitored by A260/A280 Absorbance Ratios, Biotechniques 18:62-63.
- ²³ Burgos, J.S., Ramirez, C., Tenorio, R., Sastre, I., Bullido, M.J., 2002. Influence of reagents formulation on real-time PCR parameters. Mol. Cell. Probe 16, 257–260.
- Higuchi R, Fockler C, Dollinger G, Watson R: Kinetic PCR analysis: real-time monitoring of DNA amplification reactions. Biotechnology 1993, 11:1026-1030.
- ²⁵ Schneeberger C, Speiser P, Kury F, Zeillinger R. Quantitative detection of reverse transcriptase-PCR products by means of a novel and sensitive DNA stain. PCR Methods Appl. 1995;4:234-8.
- Ririe K, Rasmussen R, Wittwer C. Product differentiation by analysis of DNA melting curves during the polymerase chain reaction. Anal Biochem. 1997;245:154-60.
- ²⁷ Terry C, Harris N, Parkes H (2002) J AOAC Int 85:768–774.
- ²⁸ Seppo I. Niemelä. Uncertainty of quantitative determinations derived by cultivation of microorganisms.
- ²⁹ Stephen A. BUSTIN. A-Z of Quantitative PCR: 414
- Rasha H. Soliman, Ahmad A. Othman, Evaluation of DNA Melting Curve Analysis Real-Time

PCR for Detection and Differentiation of Cryptosporidium Species. Parasitologists United Journal (PUJ) Vol. 2, No. 1 , 2009.

Malcolm J Burns*, Gavin J Nixon, Carole A Foy and Neil Harris. Standardisation of data from real-time quantitative PCR methods - evaluation of outliers and comparison of calibration curves.

10. Appendices

Appendice1. ANOVA statistical analysis result

• 15s

Two-way RM ANOVA	Matching by rows	
Source of Variation	% of total variation	P value
Interaction	0.44	0.2832
Plate	97.56	< 0.0001
Ct	0.50	0.0032

Source of Variation	Df	Sum-of-squares	Mean square	F
Interaction	12	6.220	0.5183	1.275
Plate	4	1365	341.3	391.0
Ct	3	6.972	2.324	5.715
Subjects (matching)	10	8.729	0.8729	2.146
Residual	30	12.20	0.4067	

• Cox3

Two-way RM ANOVA	Matching by rows	
Source of Variation	% of total variation	P value
Interaction	0.11	0.8094
Ct	92.26	< 0.0001
plate with different calibration opint	5.98	< 0.0001
Subjects (matching)	0.9286	0.3473

Source of Variation	Df	Sum-of-squares	Mean square	F
Interaction	4	0.9281	0.2320	0.3927
Ct	4	757.6	189.4	248.4
plate with different calibration opint	1	49.10	49.10	83.10
Subjects (matching)	10	7.625	0.7625	1.290
Residual	10	5.909	0.5909	

• COX1-excon1

Two-way RM ANOVA	Matching by rows	
Source of Variation	% of total variation	P value
Interaction	1.08	0.0013
Plate	91.64	< 0.0001
Ct	6.46	< 0.0001
Subjects (matching)	0.3037	0.3610

Source of Variation	Df	Sum-of-squares	Mean square	F
Interaction	8	16.06	2.008	5.210
Plate	4	1369	342.2	754.5
Ct	2	96.52	48.26	125.2
Subjects (matching)	10	4.535	0.4535	1.177
Residual	20	7.708	0.3854	

• COX1-excon4

Two-way RM ANOVA	Matching by rows	
Source of Variation	% of total variation	P value
Interaction	0.55	0.0088
Ct	97.44	< 0.0001
plate with different calibration points	1.47	< 0.0001
Subjects (matching)	0.1725	0.5333

Source of Variation	Df	Sum-of-squares	Mean square	F
Interaction	8	5.819	0.7274	3.659
Ct	4	1035	258.8	1412
plate with different calibration points	2	15.58	7.788	39.18
Subjects (matching)	10	1.833	0.1833	0.9220
Residual	20	3.976	0.1988	

• SceI

Two-way RM ANOVA	Matching by rows	
Source of Variation	% of total variation	P value
Interaction	1.04	0.2217
Ct_secl	93.06	< 0.0001
plates with differents calibrator points	3.28	< 0.0001
Subjects (matching)	0.7640	0.3110

Source of Variation	Df	Sum-of-squares	Mean square	F
Interaction	12	27.80	2.316	1.396
Ct_secl	4	2492	622.9	304.5
plates with differents calibrator points	3	87.69	29.23	17.61
Subjects (matching)	10	20.46	2.046	1.233
Residual	30	49.78	1.659	

Two-way RM ANOVA	Matching by rows	
Source of Variation	% of total variation	P value
Interaction	0.55	0.5690
Ct_ATP9	89.17	< 0.0001
Plates with different calibration points	5.75	< 0.0001
Subjects (matching)	2.9094	0.0071

Source of Variation	Df	Sum-of-squares	Mean square	F
Interaction	8	7.881	0.9851	0.8536
Ct_ATP9	4	1273	318.2	76.62
Plates with different calibration points	2	82.10	41.05	35.57
Subjects (matching)	10	41.53	4.153	3.599
Residual	20	23.08	1.154	

• Cob4

Source of Variation	% of total variation	P value
Interaction	0.40	0.0101
Ct_cbo4	89.60	< 0.0001
Plates with different calibration points	9.10	< 0.0001
Subjects (matching)	0.6092	0.0027

Source of Variation	Df	Sum-of-squares	Mean square	F
Interaction	8	5.419	0.6774	3.557
Ct_cbo4	4	1204	300.9	367.7
Plates with different calibration points	2	122.3	61.14	321.0
Subjects (matching)	10	8.184	0.8184	4.297
Residual	20	3.809	0.1904	
				1

• FIT2

Two-way RM ANOVA	Matching by rows	
Source of Variation	% of total variation	P value
Interaction	0.10	0.9626
Ct_FIT2	98.03	< 0.0001
Plates with different calibration points	1.08	< 0.0001
Subjects (matching)	0.1522	0.7033
	_	


Source of Variation D Interaction 1:		Sum-of-squares 1.958	Mean square	
Interaction 1:	2	1.050	0.4000	
		1.900	0.1632	0.3747
Ct_FIT2 4	ļ	2007	501.7	1610
Plates with different calibration points 3	3	22.20	7.399	16.99
Subjects (matching) 10	10	3.116	0.3116	0.7156
Residual 3	30	13.06	0.4355	

Appendice2. qPCR amplify curve and melting curve

15s










Amplification Curves

A7: 1.15a 1 pg	- AB: 1.159 1pg	B7: 2:159	— Bit 2.15g	B9: 2.152	C7: 3.159
- CB: 3.15s	- C9: 3.15s	- D7: 4.15s	- DB: 4.15s	- D9: 4.15s	- B8: 5.15s
£9: 5.15s	67: 10000: 15s	- Fo: 1000X 15s	F9:1000815s	H7: N.156	 Htt. N.15s
- H9: N.15:	 Noisehand 				





Amplification Curves






Helting Curves






Melting Peaks

1	A7: 1.152 1pg	AB 1.1591pg	A9: 1.152 1pg	B7: 2.152	BB: 2.150	B9: 2.15s
	C7: 3.15a	C8: 3.15s	C9: 3.15s	D7: 4.15s	— D8: 4.15≈	D9: 4.15s
	E7: 5.15s	- Ett 5.15s	£8: 5.15s	P7: 1000X 158	- Fo: 1000X 15s	FR:1000X 15s
	- 37: 1000000 15s	- GR 10000X 155	09: 1000000 15s	H7: N.15s	HB: N:15s	





Cox1_excon1

Amplification Curves





Amplification Curves






Amplification Curves










Cox1_excon4


Amplification Curves





Amplification Curves







Settings

Channel	465-510				
Color Compensation	Off				
Standard Curve	In-run				
Program	Program			L	Jnits
First Cycle	1	Last Cycle	55	Background	2-6
Noiseband Method	Noiseband (Fluor)	Noiseband	0.2542	STD Dev Multiplier	3.6906
# of Fit Points	2	Threshold Method	Auto	Threshold	0.3050


Subset Name cox1-4 22









Melting Peaks






Cox3



18 19 20 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 38 39 40 41 42 43 44 45 46 47 43 49 50 Cyclox


Melting Curves







Melting Peaks






Cob4






ATP9




Abs Quant/Fit Points for ATP9 (Abs Quant/Fit Points) Settings


olor Compensation	oe.				
olor Compensation	Oli				
Standard Curve	In-run				
Program	Program	Program			Units
First Cycle	1	Last Cycle	50	Backgroun	2-6
Noiseband Method	Noiseband (Fluor)	Noiseband	0.3675	STD Dev Multiplie	4.6816
# of Fit Points	2	Threshold Method	Auto	Threshol	0.5566



Helting Peaks



SceI




Abs Quant/Fit Points for sce (Abs Quant/Fit Points) Settings

Channel	465-510				
Color Compensation	Off				
Standard Curve	In-run				
Program	Program	ogram			Jnits
First Cycle	1	Last Cycle	50	Background	2-6
Noiseband Method	Noiseband (Fluor)	Noiseband	0.2407	STD Dev Multiplier	2.5602
			Auto		


Helting Curves






Melting Peaks


- A4 1.sec 1pg	A5: 1.sec 1pg	- B4 2.sec	B5: 2.990	B6: 2.sec	- C4: 3.sec	_
— C5: 3.sec	— СВ: 3.sec	— D4 4.sec	D5: 4.sec	D6: 4.sec	E4 5.sec	
ES 5.xec	— (ft. 5)sec	- G4: 200pg sce	- 95, 200pg see	- G6: 200pg sce	— H4: N.sce	
- HS: Nusce	- HR Nisce					



FIT2










Amplification Curves





Amplification Curves






Abs Quant/Fit Points for FIT 22 (Abs Quant/Fit Points) Settings

Channel	465-510				
Color Compensation	Off				
Standard Curve	In-run				
Program	Program	ogram			Units
First Cycle	1	Last Cycle	55	Background	2-6
Noiseband Method	Noiseband (Fluor)	Noiseband	0.2914	STD Dev Multiplie	3.7402
	2	Threshold Method		Threshok	


Helting Curves





Melting Peaks

		A12: 1.FT2-0.1pg	
	B12: 2 FIT2		C11: 3.FiT2
	- D10: 4.F1T2	D11: 4.FTZ	- D12: 4,7112
E10: 5.FIT2		E12: 5.FT2	
— F11: W8 09.0point.FIT2			
- HIZ NETZ			



Appendice3. DNA extraction protocol

- 1. Pour out 10ml YPG culture ,spin 5min at 3500rpm (+4°C)
- 2. Discard the liquid solution, resuspend pellet into 0.5 ml sorbitol buffer(1M sorbitol; 10mM Tris Ph8.0; 1mM EDTA Ph8)
 - 3. weight empty epp-tube
 - 4. tranfer resuspend solution into epp-tube, spin 1min at 5000 rpm
 - 5. Discard solution ,weight epp-tube again, get the weight of pellet(cells)
- 6. Add Zymolyase 20T(4mg/per 1g of cell),add o.5ml sorbitol buffer, incubate 30mins at 37°C (shake the tube during this time)
 - 7. spin 1min at 5000rpm, dicard suspention
 - 8. resuspend pellet in 0.5ml of $T_{50}E_{20}$
 - 9. Add 50µl of 10%SDS,mix it well
 - 10. Incubate 20mins at 65 °C
 - 11. transfer tubers on ice
 - 12. Add 200µl 5M KoAC, mix it, incubate on ice for 40mins
 - 13. spin 5mins at the max speed ,+4°C
 - 14. Transfer suspention to the fresh tube
 - 15. Add 1vol of the isopropanaol
 - 16. Mix it, let it stay 5mins at the room temperature
 - 17. Spin 1 min at max speed, dicard sup.
 - 18. wash pellet with 80% EtoH
 - 19. Keep it in the $T_{10}E_{0.1}$
 - 20. storage in -80°C

Appendice4. Calculation example

15s


Estimate copy number from plasmid concentration

		7/15x:	3166	bp			
		mumber of Nt:	6314				
		copics in	1	PE		288098, 9299	сра
nanodrop concent	ration						co
No.		concentration (ng/u1)	copies				co
	1	149.1	42955550451		107715811.4		
	2	148. 2	42696261414		107065615.3		
	3	147.4	42465782270		106487663.3		
	4	147.1	42379352591		106270931. 3		
	5	147. 3	42436972377		106415419.3		
mean		147.82	42586783821		106791088.1		
at. dv		0.828854631	238792132.4		598797.7808		
st.dv_mean		0.37067506	106791088.1		267790, 5085		
rt.dev_re1%		0.250761101	0.250761101		0.250761101		

Uncertaity of Dilution factor:

					dilution				1	
					volume(u1)	u (u1)	u rel			
					68. 9		0.001896182			
					45		0.001396132			
					40	0.100985892	0.002377464			
Spette:	0-10			stock soluti	ion					
et volume:		8.9		C. (ng/u1)	u_rel	DF	u_DF_re1	cumulative DF	u_DF_cumulative	u comp rel
emperature of the water				147.82						
No.	weight	volume		10	0.002507611	14.78	3.78952E-05			
1	0.0088	8.81767503		1	0.002507611	10	0.006080937	147.8	0.006081055	0.006577792
2	0.0087	8.717474177		0.1	0.002507611	10	0.006080937	1478		0.008957965
3	0.0088	8.81767503			0.002507611	10	0.006080937	14780		0.010826954
4	0.0089	8.917875882	1		0.002507611	10	0.006090937	147800		0.012417759
5	0.0089	8.917875882	2		0.002507511	10	0.006080937	1478000		0.013826732
	Mean:	8.8377152	3		0.002507611	10	0.006080937	14780000		0.015104844
	Stdev	0.083834048	4		0.002507611	10	0.006090937	147800000		0.01628294
Stdev of mean		0.037492	5		0.002507611	10	0.006080937	1478000000		0.017381367
Pipette:	0-10			Pipette:	0-10			Pipette:	0-10	
et volume:		5		Set volume:		10		Set volume:		1
emperature of the water				Temperature of				Temperature of the		
No.	weight	volume		No.	weight	volume		No.	weight	volume
1	0.0049	4.909841778		1	0.0101	10.12028611		1	0.0011	1.102209379
2	0.005	5.01004263		2	0.01	10.02008526		2	0.0012	1.202410231
3	0.005	5.01004263		3	0.0101	10.12028611		3	0.0009	0.901807673
4	0.005	5.01004263		4	0.0102	10.22048697		4	0.0009	0.901807673
5	0.0051	5.110243483		5	0.0099	9.919884408		5	0.0012	1.202410231
	Mean:	5.01004263			Mean:	10.08020577			Mean:	1.062129038
	Stdev	0.070852702			Stdev	0.11424655			Stdev	0.151962117
Stdev of mean		0.031686		Stdev of mea	n	0.051093		Stdev of mean		0.0679595
st.d_rel		0.006324497								

Uncertainty estimate by Kragten method:



No.	alue(logcopies	u_log												
1	5. 459541645	5. 39E-03		5. 464934609	5. 459541645	5. 459541645	5. 459541645	5. 459541645	5. 459541645	5. 459541645	5. 459541645	5. 459541645	5. 4595416	5. 4595416
2	4. 459541645	6. 00E-03		4. 459541645	4. 465546519	4. 459541645	4. 459541645	4. 459541645	4. 459541645	4. 459541645	4. 459541645	4. 459541645	4. 4595416	4. 4595416
3	3. 459541645	6. 56E-03		3. 459541645	3. 459541645	3. 466101596	3. 459541645	3. 459541645	3. 459541645	3. 459541645	3. 459541645	3. 459541645		3. 45954164
4	2. 459541645	7. 07E-03		2. 459541645	2. 459541645	2. 459541645	2. 466613236	2. 459541645	2. 459541645	2. 459541645	2. 459541645	2. 459541645	2. 4595416	2. 45954164
5	1.459541645	7. 55E-03		1. 459541645	1. 459541645	1. 459541645	1. 459541645	1.467090277	1.459541645	1. 459541645	1. 459541645	1. 459541645	1. 4595416	1.45954164
A1	19.07	0.56		19. 07333333	19. 07333333	19. 07333333	19. 07333333	19.07333333	19.64	19. 07333333	19. 07333333	19. 07333333		19. 0733333
A2	22. 03	0.19		22. 03333333	22. 03333333	22. 03333333	22. 03333333	22. 03333333	22. 03333333	22. 23	22. 03333333	22. 03333333	22. 033333	22. 0333333
A3	24. 80	0.22		24.8	24.8	24. 8	24. 8	24.8	24.8	24.8	25. 02	24.8	24. 8	24.
A4	29. 08	0.07		29. 07666667	29. 07666667	29. 07666667	29. 07666667	29. 07666667	29. 07666667	29. 07666667	29. 07666667		29. 076667	29. 0766666
A5	32. 26	0.14		32. 25666667	32, 25666667	32. 25666667	32. 25666667	32. 25666667	32. 25666667	32. 25666667	32, 25666667	32, 25666667	32. 39	32, 2566666
A_Sample	19. 23666667	0. 21		19. 23666667	19, 23666667	19. 23666667	19. 23666667	19. 23666667	19. 23666667	19. 23666667	19, 23666667	19. 23666667	19. 236667	19.44
				5. 321843885	5. 32096309	5. 319747538	5. 318884999	5. 31742607	5. 418392923	5. 341117119	5. 331948261	5. 318978933	5. 3117623	5. 25680793
			Components	0.003178225	0.00229743	0.001081878	0.000219339	-0.00123959	0.099727263	0.022451459	0.013282601	0.000313273		-0.06185772
			Compoents^2	1. 01011E-05	5. 27818E-06	1. 17046E-06	4. 81098E-08	1.53658E-06	0.009945527	0.000504068	0.000176428	9.81398E-08	4. 766E-05	0.00382637
			sum(Com^2)	0.01451829										
			u_c_logCN.	0.120491868										
			indexes	0.07%	0.04%	0.01%	0.00%	0.01%	68.50%	3.47%	1. 22%	0.00%	0. 33%	26.36

Result was correct by dilution factor:

							u_rel_d=	1. 22%	12. 16187416		
							u_rel_d= u_rel_C=	27.74%	57788. 18935		
	C_Sample_uncorrect=	5. 31866566	log(copies)								
lesult											
	C_Sample_uncorrect=(208289	±	115576. 3787)copies/µ1, k=2	1	C_undilute=	208288676.3	±	57788744. 57	27.749

Appendice5. Determination uncertainty of Pipette

Pipette:	0-10		1 1	Pipette:	0-10		Pipette:	0-10	
Set volume:				Set volume:		10	Set volume:		8.9
Temperatu re of the water:				Temperatu re of the water:			Tempera re of the water:	tu	
No.	weight	volume		No.	weight	volume	No.	weight	volume
1	0.0049	4.909841778		1	0.0101	10.12028611	1	0.0088	8.81767503
2	0.005	5.01004263		2	0.01	10.02008526	2	0.0087	8.717474177
3	0.005	5.01004263		3	0.0101	10.12028611	3	0.0088	8.81767503
4	0.005	5.01004263		4	0.0102	10.22048697	4	0.0089	8.917875882
5	0.0051	5.110243483		5	0.0099	9.919884408	5	0.0089	8.917875882
	Mean:	5.01004263			Mean:	10.08020577		Mean:	8.8377152
	Stdev	0.070852702			Stdev	0.11424655		Stdev	0.083834048

Pipette:	0-10		Pipe	ette:	0-10			Pipette:	0-10		Pipette:	0-10	
Set		8.52	Set			3.509	9	Set		6.6	Set		
volume:			volu	ume:			١	volume:			volume:		
Temperatu			Ten	nperatu			· ·	Temperatu			Temperatu		
re of the			re o	fthe				re of the			re of the		
water:			wat	ter:			١	water:			water:		
No.	weight	volume		No.	weight	volume		No.	weight	volume	No.	weight	volume
1	0.0083	8.316670767		1	0.0034	3.406828989		1	0.0065	6.51305542	1	0.0011	1.10220937
2	0.0083	8.316670767		2	0.0035	3.507029841		2	0.0066	6.613256272	2	0.0012	1.20241023
В	0.0086	8.617273324		3	0.0036	3.607230694		2	0.0066	6.613256272	3	0.0009	0.90180767
4	0.0085	8.517072472		4	0.0035	3.507029841		4	0.0065	6.51305542	4	0.0009	0.90180767
5	0.0086	8.617273324		5	0.0035	3.507029841		5	0.0065	6.51305542	5	0.0012	1.20241023
	Mean:	8.476992131			Mean:	3.507029841			Mean:	6.553135761		Mean:	1.06212903
	Stdev	0.151962117			Stdev	0.070852702			Stdev	0.054882267		Stdev	0.15196211
Stdev of	mean	0.06796	Sto	dev of	mean	0.031686		Stdev of	mean	0.024544	Stdev of	mean	0.0679

Pipette:	0-10		Pipette:	0-200		Pipette:	0-200	
Set		7.02	Set		68.9	Set		80.256
volume:			volume:			volume:		
Temperatu			Temperatu			Temperatu		
re of the	1	- 1	re of the			re of the		
water:			water:			water:		
No.	weight	volume	No.	weight	volume	No.	weight	volume
1	0.007	7.014059683	1	0.0693	69.43919086	1	0.0798	79.96028038
2	0.0071	7.114260535	2	0.0685	68.63758404	2	0.0799	80.06048123
3	0.007	7.014059683	3	0.0691	69.23878915	3	0.0798	79.96028038
4	0.007	7.014059683	4	0.0692	69.33899001	4	0.0797	79.86007953
5	0.007	7.014059683	5	0.0689	69.03838745	5	0.0796	79.75987868
	Mean:	7.034099853		Mean:	69.1385883		Mean:	79.92020004
	Stdev	0.044811184		Stdev	0.316862918		Stdev	0.11424655
		0.02004	Stdev of		0.141705	Stdev of		0.051093
	mean 0-200		Pipette:	0-200		Stdev of	0-200	0.051093
Pipette: Set		114.708	Pipette: Set		0.141705	Pipette: Set		0.051093
Pipette: Set			Pipette: Set volume:			Pipette: Set volume:		
Pipette: Set volume: Temperatu			Pipette: Set volume: Temperatu			Pipette: Set volume: Temperatu		
Pipette: Set volume: Temperatu re of the			Pipette: Set volume: Temperatu re of the			Pipette: Set volume: Temperatu re of the		
Pipette: Set volume: Temperatu re of the water:	0-200	114.708	Pipette: Set volume: Temperatu re of the water:	0-200	118.675	Pipette: Set volume: Temperatu re of the water:	0-200	103.6
Pipette: Set volume: Temperatu re of the water: No.	0-200 weight	114.708 volume	Pipette: Set volume: Temperatu re of the water: No.	0-200 weight	118.675 volume	Pipette: Set volume: Temperatu re of the water: No.	0-200 weight	103.6
Pipette: Set volume: Temperatu re of the water: No.	0-200 weight 0.1139	114.708 volume 114.1287711	Pipette: Set volume: Temperatu re of the water: No.	0-200 weight 0.1174	118.675 volume 117.635801	Pipette: Set volume: Temperatu re of the water: No.	0-200 weight 0.1021	103.6 volume 102.3050705
Pipette: Set volume: Temperatu re of the water: No. 1	0-200 weight 0.1139 0.1133	volume 114.1287711 113.527566	Pipette: Set volume: Temperatu re of the water: No. 1	0-200 weight 0.1174 0.1173	volume 117.635801 117.5356001	Pipette: Set volume: Temperatu re of the water: No. 1	0-200 weight 0.1021 0.1027	volume 102.3050705 102.9062756
Pipette: Set volume: Temperatu re of the water: No. 1 2 3	0-200 weight 0.1139 0.1133 0.1148	volume 114.1287711 113.527566 115.0305788	Pipette: Set volume: Temperatu re of the water: No. 1 2 3	0-200 weight 0.1174 0.1173 0.1174	volume 117.635801 117.635801 117.635801	Pipette: Set volume: Temperatu re of the water: No. 1 2 3	0-200 weight 0.1021 0.1027 0.1026	volume 102.3050705 102.9062756 102.8060748
Pipette: Set volume: Temperatu re of the water: No. 1 2 3	0-200 weight 0.1139 0.1133 0.1148	volume 114.1287711 113.527566 115.0305788 114.8301771	Pipette: Set volume: Temperatu re of the water: No. 1 2 3 4	0-200 weight 0.1174 0.1173 0.1174 0.118	volume 117.635801 117.5356001 117.635801 118.2370061	Pipette: Set volume: Temperatu re of the water: No. 1 2 3 4	0-200 weight 0.1021 0.1027 0.1026 0.1029	volume 102.3050705 102.9062756 102.8060748 103.1066773
Pipette: Set volume: Temperatu re of the water: No. 1 2 3	0-200 weight 0.1139 0.1133 0.1148 0.1146	volume 114.1287711 113.527566 115.0305788 114.8301771 114.7299762	Pipette: Set volume: Temperatu re of the water: No. 1 2 3	0-200 weight 0.1174 0.1173 0.1174 0.118	volume 117.635801 117.5356001 117.635801 118.2370061 118.4374078	Pipette: Set volume: Temperatu re of the water: No. 1 2 3	0-200 weight 0.1021 0.1027 0.1026 0.1029	volume 102.3050705 102.9062756 102.8060748 103.1066773 103.2068782
1 2 3 4	0-200 weight 0.1139 0.1133 0.1148	volume 114.1287711 113.527566 115.0305788 114.8301771	Pipette: Set volume: Temperatu re of the water: No. 1 2 3 4	0-200 weight 0.1174 0.1173 0.1174 0.118	volume 117.635801 117.5356001 117.635801 118.2370061	Pipette: Set volume: Temperatu re of the water: No. 1 2 3 4	0-200 weight 0.1021 0.1027 0.1026 0.1029	volume 102.3050705 102.9062756 102.8060748 103.1066773