



Intensification of bioethanol production

Simultaneous saccharification and fermentation in the oscillatory baffled reactor (OBR)

Joseph Ikwebe, Adam P. Harvey

**Bioenergy III:
Present and New Perspectives on Biorefineries**
May 22-27, 2011
Lanzarote, Canary Islands, Spain

Outline


- Introduction
- Motivation
- Oscillatory baffled reactor (OBR)
- Simultaneous saccharification and fermentation (SSF)
- Results
- Conclusion

Introduction

- Energy need has been met principally by the use of fossil fuel resources.
- Combustion of fossil fuels inevitably contributes significantly to global warming.
- Global crude oil reserves are finite.
- Society recognises the need to find renewable alternative energy sources.


Bioethanol

- Currently the dominant global renewable transport fuel.
- Produced primarily by the microbial fermentation of simple sugars, mostly hexoses from starch or cellulosic materials.
- It offers greenhouse gas savings of up to 80% over conventional fossil fuels.






State of technology estimates

source: NREL 2009




Oscillatory Baffled Reactor (OBR)






Flow patterns in un baffled and OBR cells

Current technology using batch processing



Process Intensification using the OBR



Simultaneous saccharification and fermentation (SSF)

- SSF is a process in which the saccharification and fermentation are combined into one step.
- Better than SHF.



Why?

1. Glucose is immediately consumed.
2. Capital cost is reduced, as two process steps in one vessel.


Simultaneous saccharification and fermentation (SSF)

- SSF was carried out in both OBR and shake flasks under the following conditions:
 - 38 °C, pH 4.8
 - Microcrystalline cellulose 2.5%
 - Cellulase 40 FPU/g cellulose
 - B- glucosidase 10% cellulase
 - *Saccharomyces cerevisiae* (10% inoculum)
 - YP medium (10 g/L yeast extract and 20 g/L peptone)


Ethanol fermentation profile in shake flask and OBR




OBR vs Shake flask ethanol profile



OBR vs STR saccharification of 2.5% cellulose at 50 °C and at 120 Wm⁻³




A



HV	Mag	Spot	WD	Det
20 kV	1500 x	3	8.1 mm	SE


20 µm

B




HV	Mag	Spot	WD	Det
20 kV	1500 x	3	8.2 mm	SE

20 µm



HV	Mag	Spot	WD	Det
20 kV	1500 x	3	8.2 mm	SE

20 µm



HV	Mag	Spot	WD	Det
20 kV	1500 x	3	8.2 mm	SE

20 µm

OBR vs STR vs Shake flask

	conversion of 2.5% cellulose %	Ethanol g/L	Ethanol yield %	Power density Wm^{-3}
OBR	78	8.0	89.8	120
STR	66	-	-	120
Shake flask	56.2	7.7	81.29	-

Conclusion

- Although the concentration of ethanol from the two systems were similar (8 g/L OBR vs. 7.7 g/L shake flask), the OBR attained this concentration 25 h earlier than the shake flask.
- For the same size of reactor this would lead to a ~30% higher productivity.
- At the same power density, 12% more glucose (g/L) will be available to the yeast in the OBR than in the STR.
- A better mixed SSF environment in the OBR which promoted better mass transfer characteristics was responsible for the increased yield and reduced process time.
- The OBR is linearly scaleable but in conventional STRs reaction times increase with scale, as good mixing becomes increasingly difficult to achieve.
- Minimum shear stress compared to STR, hence suitable for bioprocesses.
- OBR has potential to reduce process times and cost by 50%.

Acknowledgement

- Petroleum Technology and Development Fund (PTDF) & the Nigerian government for funding this work.
- Novozymes A/S, Bagsvaerd (Denmark) for kindly supplying the enzymes.



Thanking you