

Spring 5-23-2012

Automation and multiplexing of immunoassays: Improving precision and throughput

Ilia Tikhonov

PPD

Follow this and additional works at: http://dc.engconfintl.org/vaccine_iv



Part of the [Biomedical Engineering and Bioengineering Commons](#)

Recommended Citation

Ilia Tikhonov, "Automation and multiplexing of immunoassays: Improving precision and throughput" in "Vaccine Technology IV", B. Buckland, University College London, UK; J. Aunins, Janis Biologics, LLC; P. Alves, ITQB/IBET; K. Jansen, Wyeth Vaccine Research Eds, ECI Symposium Series, (2013). http://dc.engconfintl.org/vaccine_iv/35

This Conference Proceeding is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Vaccine Technology IV by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.



Automation of Assays: Improving Precision and Throughput

Ilia I. Tikhonov, Director of Vaccines & Biologics Lab, PPD

23 May 2012

Introduction

- Manual immunoassays, especially functional assays involving living organisms are inherently variable.
 - Assay variability is a consequence of the assay complexity and high sensitivity of live organisms to external factors.
 - Multiple external factors affecting the assay, such as variability of individual techniques between the analysts, are often hard to control.
 - Manual assays usually require extensive technical training of analysts.



Potential Benefits of Assay Automation

- Numerous potential benefits are considered when the possibility of automation in the laboratory are evaluated, including:
 - improvement of the quality and reliability of test results;
 - increased workplace ergonomics, employee productivity and safety;
 - enhancement of the assay efficiency/throughput;
 - standardization of lab techniques;
- Complete automation presents operational and financial challenges and is not always economically justifiable; partial automation of the process workflow may be beneficial.


Example 1.
Automation of the Meningococcal Serum Bactericidal Assay
(Mn SBA) for Serogroups A, C, Y and W135.

Principle of the Serum Bactericidal Assay


Step 1: Mixing bacteria with test serum



Step 2: Adding complement (●)



Step 3: Complement-mediated killing (defined as reduction in colony counts)



Manual Mn SBA ACYW135 assay

- **Classic Meningococcal Serum Bactericidal Assay.**
 - Four individual manual SBAs are performed to measure bactericidal antibody titers to *Neisseria meningitidis* serogroup A, C, Y and W135.
 - The SBA is a complex assay involving multiple manual dilutions performed by analysts.
- **Qualification studies were conducted to objectively evaluate the assay parameters and performance.**
 - The statistical analysis revealed that total variability was higher than expected.
- **Further analysis demonstrated that among the ruggedness factors, the variability associated with the analyst showed the highest impact on the assay.**

Manual Mn SBA ACYW135 Assay Variability

Serogroup	Variance estimates (%RSD)		
	Analyst to Analyst	Run within Analyst	Total
A	91%	90%	175%
C	15%	74%	100%
W135	29%	39%	65%
Y	19%	34%	74%

- The analysts were producing substantially different titers
 - Within analyst variability changed from analyst to analyst.
 - High variability was observed for both assay controls and test samples.


Automation of the sample processing for the Mn SBA ACYW135: Hamilton Microlab STAR



Summary of The Mn SBA ACYW135 Automation Effort

- Assay scripts were developed and validated. Scripts were created based on the manual process described in the SOP.
- The control and sample processing was performed by Hamilton.
- The use of the sample master dilution plate enabled simultaneous testing of the four serogroups improving assay throughput.


Standardize Sample Processing: Sample Assay Plate Transfer



Deck is reloaded 1-2 times dependent on the number of Master plates.

Improvement of Overall Assay Precision After the Implementation of the Automation Step, per Serogroup


Serogroup	Total Assay Variability (%RSD)	
	Before Automation	After Automation
A	175%	41%
C	100%	36%
Y	65%	27%
W135	74%	35%




Example 2.
Automation of Seasonal Influenza HAI

Principle of the Hemagglutination Inhibition Assay

Hemagglutination




Inhibition of hemagglutination




Seasonal Flu HAI

- The HAI assay is used to test various seasonal strains of Influenza virus
- The assay procedure consists of two analyst serially titrating a sample **replicate** with a 2-fold serial dilution for up to 12 microtiter assay plates per analyst, per strain
- Variability is increased because two analysts are generating replicate plates (and independent results)

Analyst 1: Sample Replicate 1



Analyst 2: Sample Replicate 2




Seasonal FLU HAI: Sample Processing Applications

- Similar to the MnSBA sample processing standardization
- Hamilton automation methods were designed to facilitate the testing of multiple virus strains against the same sample dilution and reduce sample repeat rates



Seasonal Flu HAI: Replicate Repeat Rate Results




Improved Biofunctional Assay Variability Conclusions

- Utilizing automation for high throughput sample processing applications has improved two of our biofunctional assays by reducing:
 - Variability and sample retest rate for Mn SBA ACYW135
 - Variability and sample retest rate for Seasonal Flu HAI

Example 3. Multiplexing and Automation of PCR Assays.

General PCR Workflow



Each of these steps can be automated

Nucleic Acid Extraction: QIAGEN BioRobot MDx (DNA) and BioMérieux NucliSens easyMAG (RNA)

- The Qiagen BioRobot MDx Workstation
 - Fully automated
 - High-throughput
 - Accurate liquid handling
 - filtered tips
 - On-board vacuum processing
 - Locked down once program begins


- The NucliSens easyMAG
 - Automated system
 - Magnetic silica beads
 - On-board magnet
 - Requires little hands-on interaction during a run



Increased Efficiencies – DNA and RNA Extraction

- Manual DNA Extraction (QIAGEN Blood Kit)

- 6 hours touch time = 384 samples



- Automated Extraction (Qiagen BioRobot MDx)

- 3 hours touch time = **1152** samples

- QIAGEN Manual Method
- 4.5 hours touch time = 72 samples

- **OR**

- Boom Silica Manual Method
- 5 hours touch time = 48 samples




- BioMérieux NucliSens easyMAG


- 3 hours touch time = **144** samples

Automation of Sample Processing by TECAN Freedom Evo

- The TECAN Freedom Evo
 - Automated liquid handling system
 - Highly accurate pipetting
 - Decreased cross-contamination
 - Can be customized
 - Different arm configurations
 - Integration of other lab equipment.
- The Molecular Testing Laboratory at PPD VBL utilizes two separate Freedom Evo configurations to maximize the sample preparation potential of our high throughput PCR process.
- TECAN Freedom Evo can process 50 assay plates/FTE/Day with 1.5-2 h of touch time, resulting in a throughput of
 - Up to 200 plates per day
 - Or approximately 6 times the throughput of a single analyst




Principle of a Singleplex PCR assay



Throughput calculation,
Example for 90 samples:

- 1 plate
- 15 different strains/organisms
- 3 Open Reading Frames
- 45 Assay Plates
- **1.5 FTE**

Principle of a Multiplex PCR assay




Throughput calculation

Example for 90 samples:

- 1 plate
- 15 different strains/organisms
- 3 Open Reading Frames measured simultaneously
- 15 Assay Plates
- **0.5 FTE**
- Note: Target genes may come from the same or a different organism.

Reaction Efficiency and Linearity of a Multiplexed PCR Assay



- Typical characteristics of the assay:
 - The assay is linear over 6 logarithms of concentration
 - Each target is amplified equivalently
 - Each target has similar sensitivity

Benefits of Automating and Multiplexing of the Assays

- Increases quality of the data compared to manual assays (reduces error rate, including sample dilutions, placement, cross-contamination, etc.)
- Increases assay throughput
- Increases the assay development flexibility
 - Combination of different targets (i.e. antigens) can be tested simultaneously
- Improves assay parameters
 - Simultaneous detection of different targets to enhance specificity and sensitivity of detection
 - Reduced variability



PPD[®]

- BACKUP

Considerations to be taken into account when automating an assay



- **Why is automation required?** What is the reason for automating a process? What is this going to achieve for the laboratory?
 - What is the testing volume?
 - What are the potential benefits (throughput, retest rate, precision)?
- **How much automation is required?** Is a fully automated system needed or a semi-automated system?
- **What systems are available?**

Potential drawbacks of automation


- Once that step is taken it is very hard to remove automation from a laboratory. It becomes a part of the validated process.
- Logistical issues. In partial automation (automation of individual steps) balancing the throughput of different steps and eliminating bottlenecks may be challenging.
- Re-training of personnel is needed.
- Significant upfront investment may be required.

Principle of the Serum Bactericidal Assay


Step 1: Mixing bacteria with test serum



Step 2: Adding complement (●)




Step 3: Complement-mediated killing (defined as reduction in colony counts)




Variability: Summary for SBAs

		Variance Estimates (% RSD)					
		Before automation			After automation		
Sero group	Analyst	Run Within Analyst	Plate Within Run	Total	Run Within Analyst	Plate Within Run	Total
A	1	10%	0%	44%	0%	5%	31%
C		86%	19%	103%	18%	20%	44%
W		0%	18%	35%	2%	4%	19%
Y		14%	13%	34%	0%	16%	23%
A	2	256%	20%	337%	49%	0%	59%
C		98%	29%	140%	17%	12%	33%
W		78%	9%	97%	0%	11%	20%
Y		41%	28%	83%	0%	0%	16%
A	3	20%	0%	74%	2%	10%	20%
C		6%	19%	41%	9%	13%	29%
W		26%	13%	41%	44%	8%	52%
Y		38%	65%	145%	0%	0%	23%

Standardize Sample Processing: Sample Master Dilution Plates




Standardize Sample Processing: Sample Assay Plate Transfer




Deck is reloaded 1-2 times dependent on the number of Master plates.

Improvement of Overall Assay Precision After the Implementation of the Automation step (Serogroups C and Y).



Improvement of Overall Assay Precision After the Implementation of the Automation Step



Applied Biosystem 3730xl DNA Analyzer-Sequencing

- The ABI 3730xl DNA Analyzer
 - Fully automated sequence analyzer
 - 16-plate stacker
 - Monitoring software
 - will alert an analyst of any errors



Increased Efficiencies - Sequencing

- GE MegaBACE 1000
 - Assay 1 = 3.0 hour run / 96-well plate
 - Plus 30 minutes touch time / plate
 - Max three 96-well plates / day
 - 288 samples
 - 1 FTE, 1.5 hours touch time
- 1.5 hours touch time = 288 samples
 - Assay 2 = 5.0 hour run / 96-well plate
 - Plus 30 minutes touch time / plate
 - Max two 96-well plates / day
 - 192 samples
 - 1 FTE, 1.0 hours touch time
- 1.0 hours touch time = 192 samples
- ABI 3730xl
 - Weekday
 - Assay 1 and 2 = 2 hour run / 96-well plate
 - Miniscule amount of touch time to load run settings
 - Max 12 plates in 24 hours
 - 1, 152 samples
 - 1 FTE, 30 mins touch time
- 30 minutes touch time = 1, 152 samples
 - Friday Night
 - Max 16 plates / stacker
 - 1, 536 samples
 - 1 FTE, 30 mins touch time
- 30 minutes touch time = 1, 536 samples

Improvement of Analyst-related Variability After Implementation of the Automation Step

Serogroup	Variability (%RSD) Attributed to Analyst	
	Before Automation	After Automation
A	91%	5%
C	15%	0%
Y	29%	0%
W135	19%	0%

Serogroup	Fold-Difference in Titer Between Analysts	
	Before Automation	After Automation
A	5.7	1.4
C	2.1	1.3
Y	2.4	1.5
W135	1.5	1.1

Benefits of Automating Immunoassays

- A contemporary immunoassay used for clinical trials is typically a complex multistep process which is qualified and validated.
- Automation of the most labor intensive and manual steps of immunoassays often improves assay precision and reduces the possibility of costly human errors.
- Automated liquid handling technology should be incorporated in the assay workflow as early as possible, ideally, at the development stage.