

6-13-2016

# Controlled, pulsatile release of thermostabilized inactivated polio vaccine from PLGA-based microspheres

Stephany Y. Tzeng

*Massachusetts Institute of Technology, sytzeng@mit.edu*

Rohiverth Guarecuco

*Massachusetts Institute of Technology*

Kevin J. McHugh

*Massachusetts Institute of Technology*

Evan M. Rosenberg

*Massachusetts Institute of Technology*

Yingying Zeng

*Massachusetts Institute of Technology*

*See next page for additional authors*

Follow this and additional works at: [http://dc.engconfintl.org/vaccine\\_vi](http://dc.engconfintl.org/vaccine_vi)

 Part of the [Engineering Commons](#)

---

## Recommended Citation

Stephany Y. Tzeng, Rohiverth Guarecuco, Kevin J. McHugh, Evan M. Rosenberg, Yingying Zeng, Robert Langer, and Ana Jaklenec, "Controlled, pulsatile release of thermostabilized inactivated polio vaccine from PLGA-based microspheres" in "Vaccine Technology VI", Laura Palomares, UNAM, Mexico Manon Cox, Protein Sciences Corporation, USA Tarit Mukhopadhyay, University College London, UK Nathalie Garçon, BIOASTER Technology Research Institute, FR Eds, ECI Symposium Series, (2016).  
[http://dc.engconfintl.org/vaccine\\_vi/15](http://dc.engconfintl.org/vaccine_vi/15)

This Abstract and Presentation is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Vaccine Technology VI by an authorized administrator of ECI Digital Archives. For more information, please contact [franco@bepress.com](mailto:franco@bepress.com).

---

**Authors**

Stephany Y. Tzeng, Rohiverth Guarecuco, Kevin J. McHugh, Evan M. Rosenberg, Yingying Zeng, Robert Langer, and Ana Jaklenec



# Controlled, pulsatile release of stable inactivated polio vaccine from PLGA-based microspheres

Stephany Y. Tzeng


Jaklenec group, Langer Lab  
Massachusetts Institute of Technology

June 13, 2016


BILL & MELINDA  
GATES *foundation*

# Inactivated polio vaccine (IPV): background

- Poliomyelitis is a potentially fatal infectious disease that can be prevented by vaccination
- For vaccine efficacy, IPV must be administered in 2-3 bolus injections (IM) spread over weeks or months




- Vaccine coverage is hampered in developing countries by difficulties of patient access in at-risk populations



# Poly(lactic-co-glycolic acid) (PLGA) microparticles for controlled release


- Biodegradable polyester

- Hydrolytically cleaved at ester linkages
- Bulk-degrading
- Degradation rate depends on lactide-to-glycolide ratio, molecular weight, and formulation parameters




- Polymer degradation leads to protein release


*IPV encapsulated in PLGA matrix*




*Burst #1: IPV near surface (diffusion)*



*Burst #2: PLGA degrades and forms pores that allow IPV release*




Pulsatile release at relevant time points can be measured from microspheres loaded with **stable antigens**



# Poly(lactic-co-glycolic acid) (PLGA) microparticles for controlled release


- Biodegradable polyester

- Hydrolytically cleaved at ester linkages
- Bulk-degrading
- Degradation rate depends on lactide-to-glycolide ratio, molecular weight, and formulation parameters



- Polymer degradation leads to protein release


- All 3 serotypes of IPV must survive encapsulation, incubation, and release without denaturation



# Poly(lactic-co-glycolic acid) (PLGA) microparticles for controlled release


- Biodegradable polyester

- Hydrolytically cleaved at ester linkages
- Bulk-degrading
- Degradation rate depends on lactide-to-glycolide ratio, molecular weight, and formulation parameters




- Polymer degradation leads to protein release

- All 3 serotypes of IPV must survive encapsulation, incubation, and release without denaturation



# IPV stability: thermostability

## Stability at 37°C



- IPV consists of three serotypes (types 1, 2, and 3), which have different stability properties
- Recovery was calculated from the D-antigenicity as measured by ELISA
- IPV was incubated at 37°C in PBS (pH 7.4) buffer
- IPV is unstable at elevated temperatures, with >50% of serotype 1 denatured after 1 week at 37°C.
- All initial thermostability studies were done with IPV in aqueous solution (not in microspheres)





# Types of excipients

<b>Sugar</b>	Protein	Amino acid(s)	Ions	pH modulator	Other
<b>Sorbitol</b>	<b>Gelatin</b>	<b>MSG (glutamate)</b>	<b>MgCl<sub>2</sub></b>	Mg(OH) <sub>2</sub>	D <sub>2</sub> O
<b>Sucrose</b>	BSA	Arginine		Mg(CO) <sub>3</sub>	
<b>Trehalose</b>	Casein	Lysine		Arginine	
<b>Maltodextrin</b>		Methionine		EPO	
Inulin		Aspartate			
Chitin/chitosan		Glutathione			

Vaccine formulations often contain a mixture of multiple excipients. Some of these are remnants of the cell culture used for vaccine production, some are added deliberately to preserve stability, and some are a combination of the above.



# IPV stability: incubation at 37°C



**Sugars** on their own did not confer improved stability on IPV during incubation at 37°C.

**Sugars** performed better when combined with **MSG** and **MgCl<sub>2</sub>** and also performed better at higher concentrations. Of the sugars tested, **maltodextrin** and **sucrose** were the best.



# IPV stability: incubation at 37°C



**Sugars** on their own did not confer improved stability on IPV during incubation at 37°C.

**Sugars** performed better when combined with **MSG** and **MgCl<sub>2</sub>** and also performed better at higher concentrations. Of the sugars tested, **maltodextrin** and **sucrose** were the best.



# IPV stability: incubation at 37°C



**Sugars** on their own did not confer improved stability on IPV during incubation at 37°C.

**Sugars** performed better when combined with **MSG** and **MgCl<sub>2</sub>** and also performed better at higher concentrations. Of the sugars tested, **maltodextrin** and **sucrose** were the best.


# IPV stability: incubation at 37°C



**Sugars** on their own did not confer improved stability on IPV during incubation at 37°C.

**Sugars** performed better when combined with **MSG** and **MgCl<sub>2</sub>** and also performed better at higher concentrations. Of the sugars tested, **maltodextrin** and **sucrose** were the best.

# IPV stability: incubation at 37°C




After two months of incubation at 37°C, IPV formulation with maltodextrin + MSG + MgCl<sub>2</sub> saw good recovery of all three serotypes, with 30-60% recovery of type 1 and 40-70% recovery of types 2 and 3.

# Poly(lactic-co-glycolic acid) (PLGA) microparticles for controlled release


- Biodegradable polyester

- Hydrolytically cleaved at ester linkages
- Bulk-degrading
- Degradation rate depends on lactide-to-glycolide ratio, molecular weight, and formulation parameters




- Polymer degradation leads to protein release

- All 3 serotypes of IPV must survive encapsulation, incubation, and release without denaturation



# IPV stability: Vacuum drying



An aqueous solution of IPV was dried under vacuum for 1 hr at room temperature.


IPV alone (without added excipients) was denatured by drying.

Excipients and excipient combinations were added to improve stability.



# IPV stability: Vacuum drying

Stability after drying




**Gelatin** and **sugars** both showed good results after drying.

**Sugars** generally showed increased protective ability in combination with amino acids (glutamate, MSG) and/or ions (MgCl<sub>2</sub>).

# IPV stability: Vacuum drying

Stability after drying




**Gelatin** and **sugars** both showed good results after drying.

**Sugars** generally showed increased protective ability in combination with amino acids (glutamate, MSG) and/or ions (MgCl<sub>2</sub>).

# IPV stability: Vacuum drying

Stability after drying





**Gelatin** and **sugars** both showed good results after drying.


**Sugars** generally showed increased protective ability in combination with amino acids (glutamate, MSG) and/or ions (MgCl<sub>2</sub>).

# Poly(lactic-co-glycolic acid) (PLGA) microparticles for controlled release

- Biodegradable polyester
  - Hydrolytically cleaved at ester linkages
  - Bulk-degrading
  - Degradation rate depends on lactide-to-glycolide ratio, molecular weight, and formulation parameters
- Polymer degradation leads to protein release
- All 3 serotypes of IPV must survive encapsulation, incubation, and release without denaturation



# IPV stability: Emulsification (sonication)




**Gelatin** and showed the best protection of IPV during the sonication step.

**Sugars and salts** generally had minimal effects on the stability during the emulsification step.

# IPV stability: Emulsification (sonication)


- Incubation at 37°C:
  - Sucrose/MSG/MgCl<sub>2</sub>
  - Maltodextrin/MSG/MgCl<sub>2</sub>

Stability after 28 days  
of incubation at 37°C: Type 1 only




- Drying:
  - Sorbitol/MSG/MgCl<sub>2</sub>
  - Maltodextrin/MSG/MgCl<sub>2</sub>

Stability after drying




Stability after sonication




- Emulsification:
  - Gelatin

# IPV release from PLGA particles

Pulsatile release









Cumulative release



*In vitro* cumulative release results are plotted as percentages of the human dose (40 DU type 1, 8 DU type 2, 32 DU type 3)

# IPV release from PLGA particles





- Sugar-based excipients co-encapsulated in PLGA particles promote IPV release in two bursts; however, very little later release is seen for types 1 and 3.



# PLGA microspheres: pH



- PLGA degrades into acidic molecules
- PLGA degrades by bulk erosion
  - Acid buildup inside microspheres
- pH of release medium (outside particles) over time suggests acid is also building up inside particles over time



# PLGA microspheres:

## Insoluble excipients buffer acidic products



- Water-soluble bases raise the pH of the IPV environment and cause denaturation
- Small molecule bases diffuse out of the microsphere faster than the antigen, depleting the basic or buffering component before PLGA degradation is completed
- $\text{Mg}(\text{OH})_2$  is an insoluble base that can be dispersed in the polymer matrix




# PLGA microspheres:

## Insoluble excipients buffer acidic products

- Cationic polymer Eudragit E PO is insoluble in water at neutral pH but soluble at  $< 3$  pH
  - Basic functional groups raise the local pH of PLGA
  - The local basic pH accelerates PLGA degradation, forming acid components that are buffered by the basic Eudragit E




# PLGA microspheres: cationic polymers as pH-modulating dopants




# Human doses released/50 mg

Burst Time (d)	Type 1	Type 2	Type 3
0-4	0.92	1.32	0.75
7-14	0.25	0.62	0.44



# Human doses released/50 mg

Burst Time (d)	Type 1	Type 2	Type 3
0-4	1.41	1.09	1.18
11-21	0.47	0.74	0.55




# Human doses released/50 mg

Burst Time (d)	Type 1	Type 2	Type 3
0-4	1.04	0.80	0.87
18-28	0.32	0.97	0.52

- The amount of Eudragit E in the particles controls the spacing of burst release
- The formulations shown above release a total of ~2 human doses spread between 2 distinct bursts, mimicking multiple bolus injections


# IPV Type 1 Neutralizing Antibodies in Serum



All groups:  
first injection  
at t=0

One bolus administration of IPV elicits no detectable neutralizing antibodies against serotype 1.

## IPV Type 1 Neutralizing Antibodies in Serum




All groups:  
first injection  
at t=0

Boost:  
t=4 wk

Two boluses administered 4 weeks apart are able to elicit neutralizing antibodies against type 1 poliovirus, although the proportion of animals protected against the virus decreased over time.

## IPV Type 1 Neutralizing Antibodies in Serum




All groups:  
first injection  
at t=0

Boost:  
t=4 wk

A single injection of IPV-containing microspheres (formulation F1) elicits a strong neutralizing response immediately after injection, without needing a second administration.

A single injection of microspheres is equal or superior to two bolus injections and appears to have even longer duration of immunity.

## IPV Type 2 Neutralizing Antibodies in Serum




All groups:  
first injection  
at t=0

Boost:  
t=4 wk

Immunity against type 2, the most stable of the antigens, is elicited well by the microsphere formulation (red line), matching or exceeding the multiple-bolus controls.



## IPV Type 3 Neutralizing Antibodies in Serum



All groups:  
first injection  
at t=0


Boost:  
t=4 wk

Interestingly, type 3 stability in the microspheres was poorer *in vivo* than in the initial *in vitro* experiments.


Further stability studies were conducted to optimize type 3 stability.

# PLGA microspheres: hydrophilic polycations for improved type 3 stability

IPV Stability



Maltodextrin, MSG, MgCl<sub>2</sub>  
+ 3% b-PEI(1.8k) + 0.008% Poly(lysine)




# Human doses released/50 mg

Burst Time (d)	Type 1	Type 2	Type 3
0-4	3.2	9.1	3.2
18-28	2.2	7.7	2.8

- An optimized concentration of poly(L-lysine) (PLL) improves IPV stability, particularly for types 1 and 3.
- Higher overall release and higher late type 3 release are observed from microspheres that co-encapsulate polycations, such as PLL and low MW polyethylenimine (PEI).

# IPV *in vivo* experiment (next iteration): Neutralizing antibodies (2 weeks)



Type 1: F1 microspheres (red) elicited neutralizing antibodies within 2 weeks, while none of the bolus controls did. The new formulation F5 (purple) is superior to F1 at the early time point.


**Type 3: The new formulation F5 (purple) is superior to the previous formulation F1 and, importantly, superior to the bolus controls. Importantly, the low dose of F5, in which the dosage was calculated by theoretical loading (i.e. assuming no loss of IPV doses due to processing or stability) is still superior to all controls and previous formulations.**

# Summary


- Excipients can stabilize IPV against thermal and physical stresses over time
- IPV co-encapsulated with stabilizers and pH-modulators can be released in a pulsatile manner from PLGA particles
- Encapsulated IPV elicits a more potent neutralizing antibody response *in vivo* than free IPV injected as a bolus

## Future Directions

- Expand pulsatile or continuous release platform to other vaccines
  - Sabin IPV (sIPV), experimental HIV vaccines, etc.
- Investigate the effect of types of release kinetics on immune response
- Increase the number of pulses that are released *in vivo*
- Pulsatile release of vaccines from core-shell particles



Poster #23



# Acknowledgements

- MIT: Langer lab  
**Robert Langer, ScD**  
**Ana Jaklenec, PhD**

**Rohiverth Guarecuco**  
**Jennifer Lu**  
**David Yang**  
**James Sugarman**  
**Allison Linehan**  
**Sviatlana Rose**  
**Thanh Nguyen, PhD**  
**Matthias Oberli, PhD**  
**Kevin McHugh, PhD**  
**Aaron Anselmo, PhD**  
**Adam Behrens, PhD**  
**Shiran Ferber, PhD**  
**Corey Bishop, PhD**  
**James Norman, PhD**



- Bill and Melinda Gates Foundation

