

Engineering Conferences International ECI Digital Archives

Cell Culture Engineering XV

Proceedings

Spring 5-10-2016

Investigating the impact of process optimization on productivity, product quality, cell metabolism, and intracellular environment

Shailendra Singh

MedImmune/AstraZeneca, singhsh@medimmune.com

Follow this and additional works at: http://dc.engconfintl.org/cellculture_xv

 Part of the [Biomedical Engineering and Bioengineering Commons](#)

Recommended Citation

Shailendra Singh, "Investigating the impact of process optimization on productivity, product quality, cell metabolism, and intracellular environment" in "Cell Culture Engineering XV", Robert Kiss, Genentech Sarah Harcum, Clemson University Jeff Chalmers, Ohio State University Eds, ECI Symposium Series, (2016). http://dc.engconfintl.org/cellculture_xv/93

This Abstract is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Cell Culture Engineering XV by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

INVESTIGATING THE IMPACT OF PROCESS OPTIMIZATION ON PRODUCTIVITY, PRODUCT QUALITY, CELL METABOLISM, AND INTRACELLULAR ENVIRONMENT

Shailendra Singh, MedImmune/AstraZeneca
singhsh@medimmune.com

Allison Lee, MedImmune/AstraZeneca
Michael Handlogten, MedImmune/AstraZeneca
Gargi Roy, MedImmune/AstraZeneca
Sophia Levitskaya, MedImmune/AstraZeneca
Sanjeev Ahuja, MedImmune/AstraZeneca
Raghavan Venkat, MedImmune/AstraZeneca

Key words: Antibody Production, Product Quality, Metabolism, Oxidative Stress, UPR

One of the key goals in process development for monoclonal antibodies is to improve productivity and product quality as needed. The early stage cell culture process developed for an antibody had titers averaging 4 g/L and variable aggregate levels within cell culture. Through process optimization work, involving changes in media formulations, feeding strategy, and process parameters, the final optimized process achieved industry leading titers (greater than 10 g/L) with consistently lower aggregate levels. To understand the impact of process optimization on the CHO cell metabolism and intracellular environment we evaluated 4 conditions: Early stage cell culture process, 2 intermediate processes, and the final optimized process. In this study we monitored multiple process performance attributes, cellular metabolites and indicators of intracellular health. The higher productivity processes had increased nutrient consumption, moderate oxidative stress and a more active unfolded protein response when compared with the early stage process. Additionally, aggregate levels were impacted by many process parameters. These findings shed important understanding into how cell culture process changes affect productivity and internal cell metabolism and are applicable to process optimization efforts to enhance productivity and product quality.