

Engineering Conferences International ECI Digital Archives

Separations Technology IX: New Frontiers in
Media, Techniques, and Technologies

Proceedings

3-7-2017

New adsorbers for the removal of genotoxic impurities from active pharmaceutical ingredients

Teresa Esteves

iBB-Institute for Bioengineering and Biosciences, Department of Bioengineering, Instituto Superior Técnico, Universidade de Lisboa, Avenida Rovisco Pais, 1049-001 Lisbon, Portugal, teresa.esteves@tecnico.ulisboa.pt

Frederico C. Ferreira

iBB-Institute for Bioengineering and Biosciences, Department of Bioengineering, Instituto Superior Técnico, Universidade de Lisboa, Avenida Rovisco Pais, 1049-001 Lisbon, Portugal

Ana I. Vicente

Research Institute for Medicine (iMED, Lisboa); Faculty of Pharmacy, Universidade de Lisboa; Av. Prof. Gama Pinto, 1649-009 Lisboa, Portugal.

Carlos A. M. Afonso

Research Institute for Medicine (iMED, Lisboa); Faculty of Pharmacy, Universidade de Lisboa; Av. Prof. Gama Pinto, 1649-009 Lisboa, Portugal.

Follow this and additional works at: http://dc.engconfintl.org/separations_technology_ix

Recommended Citation

Teresa Esteves, Frederico C. Ferreira, Ana I. Vicente, and Carlos A. M. Afonso, "New adsorbers for the removal of genotoxic impurities from active pharmaceutical ingredients" in "Separations Technology IX: New Frontiers in Media, Techniques, and Technologies", Kamallesh K. Sirkar, New Jersey Institute of Technology, USA Steven M. Crame, Rensselaer Polytechnic Institute, USA João G. Crespo, LAQV-Requimte, FCT-Universidade Nova de Lisboa, Caparica, Portugal Marco Mazzotti, ETH Zurich, Switzerland Eds, ECI Symposium Series, (2017). http://dc.engconfintl.org/separations_technology_ix/48

This Abstract and Presentation is brought to you for free and open access by the Proceedings at ECI Digital Archives. It has been accepted for inclusion in Separations Technology IX: New Frontiers in Media, Techniques, and Technologies by an authorized administrator of ECI Digital Archives. For more information, please contact franco@bepress.com.

New Adsorbers for the Removal of Genotoxic Impurities from Active Pharmaceutical Ingredients

T. Esteves, A. I. Vicente, C. A. M. Afonso, F. C. Ferreira


Institute for Bioengineering and Biosciences, Department of Bioengineering, Instituto Superior Técnico, Portugal
Instituto de Investigação do Medicamento (iMed.Ulisboa), Faculdade de Farmácia da Universidade de Lisboa, Portugal

Outline

- 1. Introduction**
- 2. Main Goal**
- 3. Results and Discussion**
- 4. Concluding Remarks**
- 5. Acknowledgements**

1. Introduction

□ Sources of GTIs in API production:




1. Introduction

□ GTIs:

- Broad range of chemical families (structural alerts)
- Electrophilic species
- React with DNA; can lead to strand breaks
- Associated carcinogenic risk

 Targeted nucleophilic sites of the DNA bases.



1. Introduction

❑ Removal of GTIs from APIs is of major importance:

- Strict regulations (FDA, EMEA) defined a **Threshold of Toxicological Concern (TTC)** :

$$\text{TTC} = 1.5 \mu\text{g/day}$$

- **TTC** corresponds to the probability of one patient in 1,000,000 to manifest the risk of having cancer.
- Below the **TTC** there is no appreciable risk to human health.
- GTI limits in APIs are calculated by dividing the **TTC** value by the maximum daily dose (g/day).

Székely G. et al., Sep. Purif. Technol., 2012, 86, 190-198.

EMA Guidelines on the "Limits on Genotoxic Impurities", EMA/CHMP/QWP/251344/2006, 2006.

FDA Guidance for Industry Genotoxic and Carcinogenic Impurities in Drug Substances and Products: Recommended Approaches; U.S. Department of Health and Human Services, 2008.

1. Introduction

❑ **Several purification stages**

to remove GTIs from APIs

❑ **Conventional techniques**

- Preparative column chromatography
- Recrystallization
- Phase exchange
- Resins
- Distillation...

**Undesirable
API losses**


❑ To achieve **GTI** content in the API at **low levels** is in many cases **extremely difficult**.

❑ Production of APIs with low GTI contents is a **major** concern for API-manufacturing companies.


2. Main Goal

To design new adsorbers as GTI scavengers.
Advanced purification technique - organic solvent compatible
platforms.

□ MIP (Molecularly Imprinted Polymer)



□ DNA base containing polymer



3. Results and Discussion

3.1. MIP

- Concept
- Model compounds
- Synthesis
- Adsorber selection
- Experimental Design - GTI Binding
- Scale-up (SPE)
- Process Design

3.2. PBI-adenine


- Concept
- Model compounds
- Synthesis
- Adsorber selection
- Characterization
- Adsorber versatility
- Process Design

3.1. MIP: Concept

□ Molecularly Imprinted Polymer (MIP)

- polymerization in the presence of a **template** molecule
- after polymerization, the template is removed and a **cavity** remains
- the polymer binds **specifically** target analytes, providing an accurate mechanism of **recognition**

- Drug delivery systems
- Sensors
- Solid phase extraction
- Chromatography




Sellergren B. et al., WO2012172075 A1, 2012.

Székely G. et al., Sep.Purif. Technol., 2012, 86, 190-198.

Kupai et al., ACS Appl. Mater. Interfaces, 2015, 7, 9516-9525.

3.1. MIP: Model compounds

GTI: Methanesulfonyl chloride
(MsCl, 114.54 Da)



GTI: 4-Dimethylaminopyridine
(DMAP, 122.17 Da)


API: Mometasone furoate
(META, 521.43 Da)

API: glucocorticoid steroid used topically to reduce inflammation of the skin or in the airways:



- treatment of inflammatory skin disorders (such as eczema and psoriasis)
- allergic rhinitis (such as hay fever)
- asthma for patients unresponsive to less potent corticosteroids

3.1. MIP: Synthesis

□ DMAP (GTI) genotoxicity:




- two structural alerting groups: aromatic and alkyl amine
- aromatic amine *in vivo* decomposition leads to electrophilic reactive species:
 - attack nucleophilic centre(s) of DNA
 - associated carcinogenic risk




3.1. MIP: Adsorber Selection

- Porogen: DCM
- T: template
- MAA: functional monomer
- EGDMA: cross-linker
- Method 1: 16h at 40 °C + 4 h at 65°C
- Method 2: 16 h at 65°C

	Stoichiometry (mmol)			Method
	T	MAA	EGDMA	
MIP1	0.1	0.4	1	1
MIP2	0.1	0.4	2	1
MIP3	0.4	0.4	4	1
MIP4	0.4	0.4	4	2




3.1. MIP: Adsorber selection



□ Kinetic experiments: maximum 93% of GTI binding reached in less than 5 min.


□ Physical properties of **MIP2** and **NIP2** by multipoint BET method.

	Surface area (m ² ·g ⁻¹)	Pore volume (cm ³ ·g ⁻¹)	Average pore Diameter (nm)
MIP2	207	0,024	5.6
NIP2	242	0,025	7.3



3.1. MIP: Experimental Design - GTI Binding

□ Model:



3.1. MIP: Experimental Design - GTI Binding

Multifactorial design:

x_1 – DMAP concentration in ppm

x_2 – **MIP2** in mg

x_3 – Solution volume in mL

Factor	Low level (-1)	Central point (0)	High level (+1)
x_1 (ppm)	7	100	600
x_2 (mg)	37.5	75	100
x_3 (mL)	1.5	3	5

Univariable design:

Bind. – binding percentage (DMAP removal from solution)

A two level face-centered design was performed in order to optimize the 3 factors.


3.1. MIP: Experimental Design - GTI Binding

	SS	DF	MS	F-value	p-value
DMAP	33.61	1	33.61	442.00	0.0302
DMAP ²	3021.29				
MIP	830.7				
Volume	842.1				
DMAP x MIP	485.4				
DMAP x Volume	909.6				
Lack of fit	96.7				
Pure Error	0.07				
Total SS	6169.7				


Model predictions:
250 < DMAP (ppm) < 350
75 < MIP (mg) < 100
Volume = 1.5 mL

ANOVA performed to the model with only statistically significant terms considered ($p < 0.05$).

(A)



(B)



Response surface plots. Effect of:
A) DMAP concentration and MIP quantity;
B) DMAP concentration and solution volume on the binding.

3.1. MIP: Scale-up (SPE)

	Batch	SPE 1 pass through
DMAP binding	93%	88%
Meta binding	12%	4%
mg DMAP / g Meta	0.79	1.25

T = RT

MIP2 = 50 mg

Detection: RP-HPLC

Solvent: DCM

flow: 0.15 – 0.42 mL/min


DMAP: 100 ppm

Meta: 10,000 ppm


3.1. MIP: Process Design

□ SPE:

Application	Nasal Spray	Cream
Maximum daily dose	200 µg	2 mg
Nº Steps	4	5
API loss	16%	20%
GTI removal	98%	> 99%
Target (mg DMAP / g Meta)	7.50	0.75
Obtained (mg DMAP / g Meta)	1.65	< 0.37




3.1. MIP: Process Design



	Rejection (%) at 10 bar
Meta (API)	99.0 ± 0.1
DMAP (GTI)	15.1 ± 0.3

□ Hybrid process:




Application	Nasal Spray	Cream
Nº Steps	1	2
API loss	4%	8%
Target (mg DMAP / g Meta)	7.50	0.75
Obtained (mg DMAP / g Meta)	3.12	0.38

3.1. MIP: Process Design

Application	Nasal Spray		Cream	
Maximum daily dose	200 µg		2 mg	
Method	SPE	OSN	SPE	OSN + SPE
Nº Steps	4	1	5	2
API loss	16%	4%	20%	8%
GTI removal (%)	98%	97%	> 99%	> 99%
Target (mg DMAP / g Meta)	7.50		0.75	
Obtained (mg DMAP / g Meta)	1.65	3.12	< 0.37	0.38



3.2. PBI-adenine: Concept

□ Polymer: PBI



(polybenzimidazole)




□ DNA base: adenine





3.2. PBI-adenine: Model Compounds

□ GTIs: DNA alkylating agents


□ APIs




- Alkyl tosylate
- Alkyl mesylate
- Dihaloalkane



- Epoxide
- Dimethyl sulfate




- Meta





- Beta

3.2. PBI-adenine: Synthesis

□ Synthetic strategy for alkylation of adenine and the PBI-A_{x%} polymer:




3.2. PBI-adenine: Adsorber selection






t = 24 h
T = RT
Solvent: DCM

rpm = 200
V = 1 mL
Detection: HPLC



3.2. PBI-adenine: Characterization

□ ¹H NMR studies: DMSO-d₆



$$Eq_{mol} = 2 \frac{A_{8.7 \text{ ppm}}}{A_{8.7 \text{ ppm}} + A_{9.1 \text{ ppm}}}$$

4. Concluding Remarks: MIP

- ❑ A **MIP** designed to target DMAP GTI was successfully synthesized.

- ❑ Effective removal of **93%** of DMAP (2 mg DMAP / g MIP) with a loss of **12%** of Meta (API) due to non-specific binding in batch experiments.

- ❑ Two step **hybrid process** design:
 - **OSN + SPE (MIP)**
 - **92% API recovery**
 - **> 99% GTI removal**

4. Concluding Remarks: PBI-adenine

- ❑ An **adsorber** was successfully designed to target **DNA alkylating agents**.

- ❑ Process design:
 - Effective removal of **96%** of **MPTS** (2 mg MPTS / g PBI-A_{12%})

5. Acknowledgements



Institute for Bioengineering and Biosciences

Frederico Ferreira
Flávio Ferreira
Filipa Ribeiro



Carlos Afonso
Ana Vicente



João Bandarra
William Heggie



Teresa Casimiro
Raquel Viveiros
José Pinto



Fundação para a Ciência e a Tecnologia
MINISTÉRIO DA EDUCAÇÃO E CIÊNCIA

PTDC/QEQ-PRS/2757/2012 (Genosep)

PTDC/QEQ-PRS/4157/2014 (SelectHost)