IEEE TRANSACTIONS ON NUCLEAR SCIENCE, VOL. 53, NO. 5, OCTOBER 2006

An Advanced Preamplifier for Highly Segmented Germanium Detectors

A. Pullia, F. Zocca, and G. Pascovici

Abstract—We present a fast low-noise hybrid charge-sensitive preamplifier for germanium position-sensitive γ -ray detectors. In conjunction with a bulky 36-fold segmented detector it provided an excellent resolution of 1.71/1.77 keV fwhm on the 1.17/1.33 MeV⁶⁰ Co lines. The preamplifier rise time, as measured at the test bench, is as fast as 7.5 ns, with a detector capacitance of 21 pF and with a 5 m 50 Ω twisted-pair cable connected at its output. The dynamic range of the preamplifier input stage is as large as 92 dB, ranging from 0.275 fC to 9.9 pC, i.e., from 5 keV to 180 MeV in terms of photon energy. On signals larger than 2 MeV a fast reset (~ 10 MeV/ μ s) is enforced to reduce the system dead time. An estimate of the amplitude of such large signals is derived from the reset time, still obtaining a high resolution. Using this technique we achieved an energy resolution of 0.3% at 16.7 MeV.

Index Terms—Charge-sensitive preamplifier, gamma ray spectroscopy, germanium detector.

I. INTRODUCTION

NEW generation of nuclear-physics experiments with radioactive ion beams is foreseen in the near future. In these experiments 4π arrays of monolithic multi-electrode positionsensitive germanium detectors will be used for spectroscopy and tracking of γ rays [1], [2]. In these experiments besides having a low noise, a fast rise time is required because the shape of the leading edge of the preamplifier signal is a key signature of the physical positions of the γ -photon interaction points inside the detector [3]. A large dynamic range is also required for a class of experiments where giant-dipole resonances can produce γ photons in the 20 MeV range [4]. Another key requirement is a reduced dead time in the presence of a background of charged particles (pions, kaons), in the 10 to 100 MeV energy range.

We have developed a dedicated preamplifier, optimized in such a way to feature low noise, a fast rise time, and a large dynamic range. In addition the preamplifier is equipped with a time-variant device which quickly resets the circuitry on large signals in excess of a few MeV. This device substantially reduces the dead time caused by large signals or random walks that could saturate the output buffer or the ADC voltage range.

Two versions have been realized. One is optimized for cathode detector electrodes. One other for anode detector electrodes. Both include a fast differential buffer for low-noise signal transmission to a remote ADC.

Manuscript received December 4, 2005; revised June 23, 2006.

A. Pullia and F. Zocca are with the University of Milan, Department of Physics and Istituto Nazionale di Fisica Nucleare, I-20133 Milano, Italy. G. Pascovici is with the University of Koeln, Institüt für Kernphysik, D-50937

Koeln, Germany.

Digital Object Identifier 10.1109/TNS.2006.880402

 $\begin{array}{c} \text{detector} \\ \text{FET} \\ \hline \end{array} \\ \hline \end{array} \\ \hline \\ \text{Fig. 1. Structure of the charge-sensing stage. The chosen values for the feed-$

II. PREAMPLIFIER STRUCTURE

back devices are $C_F = 1 \text{ pF}$ and $R_F = 1 \text{ G}\Omega$.

The charge-sensing stage of the preamplifier comprises (a)an input Field-Effect Transistor and a feedback passive network mounted in the cryostat as close as possible to the detector electrode, and (b) a Transimpedance Amplifier (TA) operated at room temperature outside the cryostat, as shown in Fig. 1. In Fig. 2 the simplified schematic diagrams of the realized transimpedance amplifiers are shown, for the cases of cathode and anode detector signals. In both cases a 'folded cascode' structure is adopted for the gain stage, with the unusual solution of using two cascaded bipolar transistors in the common-base path (T1/T2). This solution has these advantages: (i) it yields a substantial increase of the loop gain and (ii) it permits to maximize the negative/positive output-voltage range for either cathode [Fig. 2(a)] or anode [Fig. 2(b)] detector signals. The drawback of the double-transistor cascode solution is an increase of the electronic noise, which, however, may be made negligible by properly biasing the circuit structure. The output voltage range of the TA is about 10 V for either cathode or anode signal versions. As the closed loop gain is \sim 55 mV/MeV $(C_F = 1 \text{ pF})$ we obtain an energy range of ~ 180 MeV. From experimental tests we found that the low-level discrimination threshold can be set down to ~ 5 keV. So we get a dynamic range of as much as 92 dB for the integration stage.

To get a fast rise time we have used a fast input FET, mod. BF862, and broadband bipolar transistors, mod. BFT92 and BFR92, with a cut-off frequency in excess of 1 GHz at the used collector bias currents. We have carefully chosen the collector current I_C of the cascode transistors in order to guarantee both a low noise and a wide bandwidth, finding an optimum value for I_C of 380 μ A. A further increase of I_C yields a higher noise and no significant improvement of the preamplifier bandwidth.

In Fig. 3 the simplified schematic diagram of the entire preamplifier (version for anode signals) is shown. The



2869



Fig. 2. Simplified schematic diagram of the realized transimpedance amplifier. (a) Case for cathode detector signals, or negative output voltage swing. (b) Case for anode detector signals, or positive output voltage swing. Capacitance C_C is used to stabilize the feedback loop.



Fig. 3. Simplified schematic diagram of the preamplifier. It consists of an integration stage, a Pole/Zero compensation network, and a fast-reset circuitry which automatically switches on when the output voltage overcomes a preset threshold.

charge-sensitive stage of the preamplifier is followed by a Pole/Zero stage and a differential line driver. The differential signal transmission mode is useful to enhance the rejection to common-mode noise and disturbances picked up along the output cable. The signal can so be transmitted to a remote ADC module through a cable of several meters preserving its original quality, which makes the preamplifier suitable for noisy experimental environments. The fast-reset circuitry is put in evidence by a dashed line in the lower part of the schematic diagram. It consists of a Schmitt trigger comparator which switches on/off a temperature-compensated current sink. When a large signal arrives the current sink is switched on, which discharges capacitance C and quickly restores voltage y(t) to 0. As soon as y(t) (as well as the output voltage) reaches 0 V, the current sink is switched off, and normal operation is restored [5].

III. FEEDBACK LOOP ANALYSIS

Maximization of the preamplifier loop gain is of primary importance in gamma-ray spectroscopy with germanium detectors. In fact the energy resolution obtainable with these detectors is as high as 0.2% at 1 MeV. In order to achieve and maintain such a good resolution an outstanding long-term gain stability of less than 0.2% is required. Provided that a high-quality feedback capacitance is used, the larger the loop gain and the better the preamplifier long-term gain stability.

The double-transistor cascode structure yields (i) a larger overall load resistance at the cascode output, and (ii) a reduced Miller effect on the gate-to-drain capacitance of the input FET. Both of these features contribute to maximize the open-loop gain as well as the loop gain of the circuit, as is described next.



Fig. 4. Open loop gain and loop gain of the integration stage versus frequency. Three different cases are shown, which subsequently yield an improvement of the gain value at low frequency.

A. Open-Loop Gain

Let us consider the "gain node" of the circuit, i.e., respectively point P or Q in Figs. 2(a) and (b). At low frequency it can be simply modeled by a current sink and an equivalent load resistance using the Norton theorem. The current comes from the input FET, which acts as a transconductor transforming voltage signals upon its gate into current signals through its drain. Let us call "g_m" its transconductance. The drain current signal flows mostly through the common base path of the cascode structure and reaches thereby the load resistance R_L . The open loop gain A_0 is given by

$$A_0 = \alpha \quad g_m \quad R_L \tag{1}$$

where α is a loss factor which takes into account that a fraction of the FET current does not enter T1's emitter but flows in part through R_D and in part through the Early resistance of the FET. In this case $g_m \approx 38$ mA/V and α is found to be in the range from 0.85 and 0.95 for the used bias points.

The load resistance is given by

$$R_{\rm L} = R_{\rm B3} ||R_{\rm I1}||R_{\rm C2} \tag{2}$$

where symbol "||" stands for "parallel connection of". R_{B3} is the dynamic resistance looking into the base of transistor T_3 , R_{I1} is the dynamic resistance of active load I_1 , R_{C2} is the dynamic resistance looking into the collector of transistor T_2 .

We are to maximize R_L and so also A_0 . For this purpose let us consider three cases: (a) single-transistor cascode and resistive load in place of I_2 , (b) single-transistor cascode and active load I_2 , (c) two-transistor cascode and active load I_2 . As can be seen

TABLE I RELEVANT RESISTANCE VALUES

	R _{B3}	R _{II}	R _{C2}	R _L	A_0
(a)	312 kΩ	4.92 MΩ	487 kΩ	183 kΩ	6.22×10^3
(b)	3.4 MΩ	4.92 MΩ	$487 \ k\Omega$	388 kΩ	1.32×10^4
(c)	3.4 MΩ	4.92 MΩ	1.68 MΩ	917 kΩ	3.03×10^{4}



Fig. 5. Open-loop equivalent circuit used for the calculation of the loop gain. The Transimpedance Amplifier structure is shown in Figs. 2(a) and Fig. 2(b). The capacitance values are shown in Table II.

in Table I, going from (a) to (b) yields an increased value of $R_{\rm B3}.$ Going from (b) to (c) yields an increased value of $R_{\rm C2}.$ So A_0 is improved step by step. In Table I the values of $R_{\rm B3}, R_{\rm I1}, R_{\rm C2}$ and A_0 as obtained from a computer simulation are shown. The open-loop gain A_0 is found to be in good agreement with (1). The open-loop gain vs. frequency plot is shown in Fig. 4 for cases (a), (b) and (c). The dominant pole is due to capacitance C_C and can be adjusted to ensure circuit stability.



Fig. 6. Noise of preamplifier. A quasi-Gaussian shaper with 6 μ s shaping time was used to filter the signal. The noise decreases when mounting the FET on a ptfe board.

TABLE II Relevant Capacitance Values						
	C _F	CD	C _{GS}	C _{GD}		
1 pF		21 pF	4.4 pF	3.4 pF		
TABLE III RELEVANT PARAMETERS						
(a)	9	30.7 pF	1.1×10^{2}	57		
(b)	15	51.1 pF	1.7×10^2	78		
(c)	2.5	8.5 pF	0.84×10^{3}	36		

B. Loop Gain and Miller Effect on C_{GD}

In Fig. 5 the equivalent circuit used to calculate the loop gain G_0 is shown. Taking into account the Miller effect acting on C_{GD} , G_0 is easily obtained, or

$$G_0 = \frac{\tilde{V}_o}{V_T}$$

= $-\frac{C_F}{C_F + C_D + C_{\rm GS} + C_{\rm GD}(1+K)} A_0$ (3)

where C_F and C_D are the feedback and detector capacitances, C_{GS} and C_{GD} are the gate-to-source and gate-to-drain capacitances of the input FET, and K is the gate-to-drain voltage gain responsible for the Miller effect on C_{GD} . It can be seen that the loss term brought about by the capacitance divider decreases as K is increased. In Tables II and III the relevant parameters as obtained from a computer simulation are shown.

Note that going from (a) to (b) the gate-to-drain voltage gain K nearly doubles. Here is why: the increased value of the load resistance yields an increased value of the resistance looking



Fig. 7. Electronic noise versus shaping time.

into T_1 's emitter, whose value ranges from $1/g_m$ to β/g_m when a high-value resistance is connected from its collector to ground [6]. The increased value of the resistance looking into T_1 's emitter yields in turn an increased value of K and a stronger Miller effect on C_{GD} . In order to reduce K the double-transistor cascode solution is introduced, which strongly reduces the value of the resistance connected from T_1 's collector to ground. The beneficial effect of such a circuit solution can be seen in Table II going from (b) to (c), where K drops substantially and G_0 subsequently grows by almost one order of magnitude. The loop gain vs. frequency plot is shown in Fig. 4 for cases (a), (b) and (c).

IV. EXPERIMENTAL RESULTS

We fully characterized the preamplifier at the test bench, using a capacitor to simulate the impedance of the detector electrode and a fast pulser to simulate the events. In Fig. 6 the noise of the preamplifier is shown for 36 different preamplifier channels, as measured operating the FET at room temperature and using a 21 pF detector capacitance. The noise as a function



Fig. 8. Rise time of preamplifier. Results obtained at room temperature with $C_{\rm det}$ = 21 pF.

TABLE IV PREAMPLIFIER SPECIFICATIONS

Conversion gain – Ge det. (terminated)	100 mV/MeV +/- 10 mV/MeV	
Rise time	7.5 ns (21 pF) differential signal, 5m terminated cable, warm FET (mod. BF862)	
Slope	< 0.4 ns / pF (0 to 79 pF range)	
Fall time	50 µs +/- 5 %	
	0.80 keV fwhm (0 pF) @ 3 µs shap time [*]	
Noise	0.77 keV fwhm (0 pF) @ 6 μs shap time ້	
	*differential signal, warm FET (mod. BF862)	
Noise slope	7.7 eV / pF @ 3 μs shap time 5.0 eV / pF @ 6 μs shap time (in 0 to 68 pF range)	
Power required (including FET)	+12 V ; 8.5 mA (I _D =7.8 mA, V _D =2.2 V) -12 V ; 3.7 mA +6 V ; 16 mA -6V ; 6.7 mA 283 mW	
Integral nonlinearity	0.04 % for a full scale range of +/- 3.5 V (unterminated)	
Mechanical dimension (for a triple preamplifier)	57x56x5 mm (including connectors)	

of the shaping time is shown in Fig. 7. It reaches a minimum of about 140 electrons r.m.s. at a shaping time of 6 μ s.

In the AGATA (Advanced GAmma ray Tracking Array) [1] setup the preamplifiers will be connected to sampling modules through long 50 Ω -terminated cables. In order to (*a*) maintain a wide bandwidth while driving such a low impedance load, and (*b*) guarantee an outstanding gain stability, use of a high-quality output stage as operated in a closed-loop configuration is mandatory. The adopted output stage consists of a dual operational amplifier configured as a differential buffer, which helps reducing the common-mode disturbances picked up along the cable. The results are shown in Fig. 8 observing the signal after a 5m 50 Ω -terminated twisted-pair cable. Two different dual operational amplifiers have been tested, both quoted for low



Fig. 9. Zoom of the 1.17/1.33 MeV lines of 60 Co spectrum.



Fig. 10. (a) Large signal (H) upon a tail (h) at the PZ stage output. (b) Schmitt-trigger comparator signal.

noise, low power and wide bandwidth. A rise time of 13 ns was obtained with the LM6172 and of as little as 7.5 ns with the AD8012. It is worth observing that the full bandwidth of the charge sensing stage can be exploited only if the output buffer maintains its performance when driving a low-impedance load.

If the detector capacitance is decreased from 21 pF to 5–10 pF as is typical for large multichannel systems, the loop gain (3) grows and the phase margin decreases. A ringing appears in the preamplifier response. However it can be easily removed stabilizing the feedback loop by increasing capacitance C_C of Figs. 2(a) and (b), and so inherently reducing the gain-bandwidth product of the loop gain. In Table IV the specifications of the preamplifier are summarized.

Finally we connected the preamplifier to a cathode segment of an HPGe (High Purity Germanium) detector for the AGATA demonstrator [1] and collected the⁶⁰ Co spectrum shown in Fig. 9. An excellent resolution was achieved of 1.71/1.77 keV fwhm on the 1.17/1.33 MeV lines. A resolution of 770 eV fwhm was obtained on the 59 keV²⁴¹ Am line, which is a good result for γ -ray spectroscopy systems. Of course better resolutions could be obtained at 59 keV with



Fig. 11. Spectrum of 152 Eu with a distribution of high-energy pulser lines. Note the scale break along the x axis. The average event rate is 6 kHz.

preamplifiers designed mostly for X-ray spectroscopy, at the expense of a much narrower available energy range.

In the AGATA setup the power dissipation of the input FET, operated at cryogenic temperature, is as low as 17.1 mW (drain voltage = 2.2 V, drain current = 7.8 mA) which yields an overall power dissipation of 630 mW for each 37-electrode detector crystal. No micro-boiling of the cryogen has been observed. The inter-channel cross-talk has been reduced to less than 0.1% by carefully shielding the region where the detector electrode and the input FET are electrically connected.

V. RECONSTRUCTION OF LARGE SIGNALS FROM RESET TIME

The typical energy range of preamplifiers for γ -ray spectroscopy is from a few keV to a few MeV. However in the next generation of nuclear physics experiments with exotic beams, a background of highly energetic charged particles (up to 100 MeV) is foreseen. Moreover, at high count rates, bursts of piled-up signals can easily exceed the typical energy range. These events cause ADC saturation and introduce a significant system dead time, which could even lead to a paralysis of the acquisition system.

The fast-reset device (see Fig. 3) has been implemented to address this issue. It yields a full reset of the preamplifier whenever it receives a large signal, i.e., a signal which overcomes a preset threshold, corresponding to the ADC saturation point. The reset time ranges from 0.5 to 150 μ s, for signals in the 5 to 150 MeV range. For signals below the threshold a normal resistive continuous reset is performed. As compared to pure continuous or pulsed reset solutions [7]–[9], such operation mode provides an improved performance: it maximizes the live time in the typical energy range and minimizes the dead time caused by amplitude saturation. Moreover it allows to measure the amplitude of large individual signals from the reset time, when the ADC overflow condition would give no possibility to make a direct measurement of the pulse height. As a result, the dynamic range is substantially increased.

In Fig. 10(a) we show the typical reset transient seen at the P/Z amplifier output (see Fig. 3), when a large individual signal (H) occurs upon a tail of a regular signal (h). Fig. 10(b) shows the corresponding waveform seen at the Schmitt-trigger comparator output, which delivers a rectangular pulse of width T, corresponding to the fast-reset time. The regular signal which provides the tail h has the typical exponential decay time of about 50 μ s. The large signal H is instead over threshold and puts the preamplifier into the fast reset mode. By accurately calculating the expression of the reset transient at the P/Z stage output, we found an exponential relation between the reset time T and the sum of the energy E_S of the large event and the energy E_C of the tail, or

$$E_S + E_C = A \left(e^{T/\tau} - 1 \right), \tag{4}$$

where A and τ are coefficients depending on circuit parameters [5]. Because the reset time T is substantially shorter than the time constant τ , the relation can be linearized with little loss of accuracy by expanding the exponential term. We can then express the large signal energy E_S as a function of the reset time T and the energy of the tail E_C , obtaining

$$E_{S} = b_{1}T + b_{2}T^{2} - E_{C},$$
 (5)



Fig. 12. Energy resolution obtained from the reset time at high energy for two different values of ¹⁵² Eu "background" event rate.

where b_1 and b_2 are the the first- and second-order coefficients of the Taylor expansion. The reset time T can be obtained through a time over threshold measurement performed on the waveform acquired at the comparator output [see Fig. 10(b)], while the tail energy E_C can be obtained from the difference of the pre- and post-large-event voltage levels evaluated on the output signal of the preamplifier. So (5) can be written as

$$E_{\rm S} = b_1 T + b_2 T^2 - (V_1 - V_2)/G + E_0, \qquad (6)$$

where (V1 - V2) is the voltage amplitude of the tail, G is the energy-to-voltage gain of the preamplifier and E_0 is an offset term. We can conclude that (6) permits to indirectly derive the large signal energy E_S . In this relation T, V_1 and V_2 must be evaluated every time an over-threshold signal occurs, while b_1 , b_2 , G and E_0 are parameters obtained through preliminary calibrations.

We tested this operation mode using an ORTEC coaxial HPGe detector [10], a calibration radioactive source and a pulser to simulate the arrival of energetic events. Fig. 11 shows a spectrum of¹⁵² Eu obtained while also injecting large pulser signals in the 20 MeV range through a calibrated test capacitance. The¹⁵² Eu lines have been histogrammed using a quasi Gaussian shaper amplifier and a Pulse Height Analyzer, obtaining an energy resolution of 1.25 keV fwhm on the 122 keV line and 2 keV fwhm on the 1.4 MeV line. The large pulser signals, instead, have been histogrammed using the reset times and a proper calibration algorithm, according to (6). As can be seen, we could histogram very large signals of up to 27 MeV, enhancing the available preamplifier dynamic range by more than one order of magnitude as compared to the typical case. In Fig. 12 the energy resolutions obtained at high energy are shown for two different values of ¹⁵² Eu "background" events.

Note that such a large dynamic range was obtained with no sacrifice in terms of energy resolution in the typical energy range for gamma spectroscopy. Namely we achieved a good energy resolution of about two per thousand in the 1 MeV energy range and still a good energy resolution of about three per thousand in the 20 MeV energy range. Such a large energy range, from \sim 5 keV up to \sim 20 or 30 MeV, can be of great interest for those experiments in which giant-dipole resonances up to 20 MeV are to be measured.

VI. CONCLUSION

A low-noise fast charge preamplifier purposely conceived for position-sensitive germanium detectors has been designed, constructed, and characterized. It comprises a fast-reset circuitry to reduce the dead time in the presence of a background of energetic events or large random walks. We could estimate the amplitude of the large events from the measurement of the reset time. The circuit structure of the preamplifier has been presented and discussed. We also carried out and presented an experimental characterization of the circuit as connected to a state-ofthe-art HPGe detector.

ACKNOWLEDGMENT

The authors would like to thank D. Bazzacco, R. Isocrate, C. Boiano, R. Bassini, and J. Eberth for their suggestions and stimulating discussions.

REFERENCES

- D. Bazzacco, B. Cederwall, J. Cresswell, G. Duchene, J. Eberth, and W. Gast, AGATA Tech. Proposal [Online]. Available: http://agata.pd. infn.it/documents/Agata-proposal.pdf
- [2] I. Y. Lee, "Gamma-ray tracking detectors," Nucl. Instrum. Methods Phys. Res. A, vol. A422, pp. 195–200, 1999.
- [3] Th. Kröll and D. Bazzacco, "Simulation and analysis of pulse shapes from highly segmented HPGe detectors for the gamma-ray tracking array MARS," *Nucl. Instrum. Methods Phys. Res. A*, vol. A463, pp. 227–249, 2001.
- [4] P. F. Bortignon, A. Bracco, and A. R. Broglia, *Giant Resonances-Nuclear Structure at Finite Temperature*. Australia: Harwood Academic Publishers, 1998.
- [5] A. Pullia and F. Zocca, "A low-noise preamplifier for γ -ray sensors with add-on device for large-signal management," *Nucl. Instrum. Methods Phys. Res. A*, vol. A545, no. 3, pp. 784–92, 2005.
- [6] J. Millman and A. Grabel, *Microelectronics*, 2nd ed. New York: Mc-Graw-Hill, 1987, ch. 10.
- [7] D. A. Landis, N. W. Madden, and F. S. Goulding, "Bipolar pulsed reset for AC coupled charge-sensitive preamplifiers," in *Proc. IEEE Nuclear Science Symp. Conf. Rec.*, Albuquerque, NM, Nov. 9-15, 1997, vol. 1, pp. 887–891.
- [8] T. Lakatos, G. Hegyesi, and G. Kalinka, "A simple pulsed drain feedback preamplifier for high resolution high rate nuclear spectroscopy," *Nucl. Instrum. Methods Phys. Res. A*, vol. 378, no. 3, pp. 583–588, 1996.
- [9] D. A. Landis, C. P. Cork, N. W. Madden, and F. S. Goulding, "Transistor reset preamplifier for high resolution spectroscopy," *IEEE Trans. Nucl. Sci.*, vol. NS-29, no. 1, pp. 619–624, 1982.
- [10] Th. Kröll, D. Bazzacco, R. Venturelli, B. Quintana Arnes, and C. A. Ur, "Gamma-ray tracking with the MARS detector," *Eur. J. Phys.*, vol. 20, no. 1, pp. 205–206, 2003.