

School of Biomedical Engineering, Science and Health Systems

Biomedical Technology Showcase, 2006



Drexel E-Repository and Archive (iDEA)

<http://idea.library.drexel.edu/>

Drexel University Libraries

www.library.drexel.edu

The following item is made available as a courtesy to scholars by the author(s) and Drexel University Library and may contain materials and content, including computer code and tags, artwork, text, graphics, images, and illustrations (Material) which may be protected by copyright law. Unless otherwise noted, the Material is made available for non profit and educational purposes, such as research, teaching and private study. For these limited purposes, you may reproduce (print, download or make copies) the Material without prior permission. All copies must include any copyright notice originally included with the Material. **You must seek permission from the authors or copyright owners for all uses that are not allowed by fair use and other provisions of the U.S. Copyright Law.** The responsibility for making an independent legal assessment and securing any necessary permission rests with persons desiring to reproduce or use the Material.

Please direct questions to archives@drexel.edu




Improving the efficacy of cellular therapy by magnetic cell targeting

¹Zachary G. Forbes, PhD, ²Kenneth A. Barbee, PhD, ¹Frederick Stoddard, MD, ¹Ari D. Brooks, MD, ³Doris A. Morgan, PhD, ⁴Benjamin B. Yellen, PhD, ⁵Gary Friedman, PhD.






¹Drexel University College of Medicine, Dept. of Surgery Philadelphia, PA, ²Drexel University School of Biomedical Engineering, Philadelphia, PA, ³Drexel University College of Medicine, Div. of Hematology/Oncology, Philadelphia, PA, ⁴Duke University Dept of Mechanical Engineering and Materials Science, Durham, NC, ⁵Drexel University Dept of Electrical and Computer Engineering, Philadelphia, PA



Background





Objectives <ul style="list-style-type: none">To extend our concurrent efforts in magnetic drug delivery to coronary stents using magnetic nanoparticles to the delivery of cells and biological agentsTo accomplish this <i>in vivo</i> by applying high magnetic field gradients within the body to an injected suspension of magnetized cells, without relying on large external magnetic gradients to apply forces	Why Target Cells with Magnets? <ul style="list-style-type: none">Utilizes our effective drug targeting techniqueCapabilities:<ul style="list-style-type: none">Efficient localization of cellsManage small resource supplyMinimally invasiveImprove quality of treatmentAutologous or Non-immunogenic Cell SourceApplications:<ul style="list-style-type: none">Implant endothelializationVascular tissue regenerationCells as a drug delivery vehicle 	High Gradient Magnetic Delivery <ul style="list-style-type: none">The Magnetic Implant Based Cell & Drug Delivery System<ul style="list-style-type: none">Inject magnetically loaded cells with the aid of a modest external field, allowing highly efficient targeting of cells or drugsThe magnetic forces of the implant, when exposed to a uniform magnetic field, can apply the forces needed to capture susceptible agents from blood flow 	Our method: Two magnetic sources  <ol style="list-style-type: none">One source: Maximize the carrier's moment $\vec{m} = \chi V \vec{H}$Another source: Maximize the field gradient $\vec{F} = \mu_0 (\vec{m} \cdot \nabla) \vec{H}$ <p><i>Independently controlling these variables can optimize the forces which can be applied on injected drug carriers.</i></p>
---	---	---	---

Experiments





Cell Culture Models <ul style="list-style-type: none">Objectives<ol style="list-style-type: none">Develop a reproducible method for quantifiable loading of endothelial cells with superparamagnetic particlesParticle uptake in suspension or adherent culture, and subsequent biological responseMagnetic capture of endothelial cells with internalized nanoparticles in parallel plate flow and in vivo	Model Implant: Wire Mesh  <ul style="list-style-type: none">140 µm wire400 µm spacing200 L, 10100cm <p>Stainless Steel Co-Ni Magnetic Coating</p>	Magnetic Measurements  <p>Determination of a 0.02 Tesla saturation for the applied magnetic field by Alternating Gradient Magnetometer measurement of a single magnetically plated wire</p>	Cell Loading Groups <ul style="list-style-type: none">Bovine Aortic Endothelial Cells<ul style="list-style-type: none">350nm diameter, rhodamine-stained polystyrene magnetic particles, carboxyl terminusPreliminary non-degradable particle modelHuman Umbilical Vein Endothelial Cells<ul style="list-style-type: none">130nm diameter, dextran magnetic particles, NH₂ terminusDegradable particles suitable for in vivo applications
Particle loading to ECs <ul style="list-style-type: none">Particle suspension vortexedParticles seeded to adherent cultures in BAECs, in suspension to HUVECsNo magnetic field source usedCells with endocytosed particles<ul style="list-style-type: none">Magnetometer analysis concludes approximately 10 particles per BAEC at saturationNo saturation data for 130nm particles, but contribute load at ~10³ particles per cell 	In Vitro Experimental Setup  <p>Flow Chamber FLOW IN FLOW OUT Flow Chamber</p> <p>Applied magnetic field = 0.05 Tesla</p> <p>Flow velocity = 15 cm/s</p> <ul style="list-style-type: none">Magnetically endothelialized "wired" stent encased in a shaped stentStent is stored in a vial and collagen coated to facilitate tissue integration	Current In Vivo Studies (Unloaded Magnetic Particles) <ul style="list-style-type: none">Sprague Dawley Rats – non survival130nm magnetic particlesDextran encapsulatedTc99m linked to PEG-NH₂1mm diameter, 5mm long stainless steel stent implanted in common iliac arteryDeliver dose intra-venously under 0.05 Tesla applied uniform magnetic field for varied time stepsHarvest device, organs, for biodistribution counts	Surgery Images  <p>A rat receives a magnetic particle dose intravenously under exposure to a uniform static magnetic field at 0.05 Tesla</p>

Results


In vitro Loading Studies and Parallel Plate Flow Chamber Analysis

Magnetic Endothelial Cells  <p>(TOP) Fluorescent images of Bovine Aortic Endothelial Cells loaded with ~10³ 350nm diameter magnetic particles. (BOTTOM) Confocal images of GFP stained BAECs loaded with 130nm magnetic particles</p>	Capture Results – 30 minutes  <p>5X 10X</p> <p>Dense capture of magnetically loaded BAECs is seen along the wires of the magnetic stent, with virtually no cells in between the struts.</p>	Capture Results – 6 Hours  <p>5X 10X</p> <p>After 24 hours in culture at 37°C 5% CO₂</p> <p>Magnetically loaded BAECs are beginning to migrate off the intersections and attach elsewhere</p>	Capture Results – 24 hours  <p>5X 10X</p> <p>After 24 hours in culture at 37°C 5% CO₂</p> <p>At 24 hours, BAECs have continued to migrate and divide across the stent, indicating the ability of mature cells loaded with magnetic particles to withstand the trauma of applied magnetic forces</p>
--	---	---	--

In Vivo Magnetic Particle Delivery to Stents in a Rat Common Iliac Artery

Qualitative Results  <p>Image of implanted stent (contrast added) and gamma image of isolated stent.</p> <p>The gamma images clearly indicate the high activity of the stent implant site (this activity is 40 times as high as the rest of the body) and the low level of activity elsewhere.</p>	Dose Dependent Delivery  <p>Dose dependent capture was seen over a range of doses, beginning to level off at 1 mT. This data reflects only one field strength and exposure time.</p>	Biodistribution of Uncaptured Particles  <p>Concurrent relationships in biodistribution shows, encouraging in the capacity to reduce uptake in the lungs and liver while maintaining robust stent capture</p>	Time Varied Magnetic Dosing  <p>A 25% increase in particle capture was seen over a 15 minute increase in field exposure before sacrifice and organ harvest</p>
--	--	--	--

Conclusions

Proof of Principle <ul style="list-style-type: none">Functionality of the our method<ul style="list-style-type: none">Sub-micron and nano-scale magnetic particles can be captured, by controllable means, to the surface of a magnetizable implantFlexibility in application<ol style="list-style-type: none">A microlayer of soft magnetic deposition on a stent or other non-vascular implantsVarying the alloy content of a stent or other non-vascular implantsImpregnating magnetic crystals into a polymer scaffold or woven graft	Transitions to Cellular Therapy <ul style="list-style-type: none">Biocompatibility<ul style="list-style-type: none">Vascular endothelial cells have been shown to grow upon magnetized surfaces, uptake magnetic nanoparticles, maintaining normal growth rates, morphology, and behavior.These endothelial cells can be delivered magnetically to the surface of a stent and withstand potentially traumatic forces	Implementation <p>From particles, to mature cells, to stem cells</p> <ol style="list-style-type: none">As we have recently with lab grown biofilms comprised of mature cells, in many cases it is not necessary to involve stem cells when adequate endothelial structures of mature cells are availableThis project is following a logical path of adapting to mature vascular endothelial cells, while beginning to explore efforts to an off-targeted in cord blood stem cells loaded with magnetic nanoparticlesAs we have demonstrated that cells can be loaded in suspension in controlled fashion, in the clinic this may limit procedural steps and risks of contaminationEndothelial cell suspensions loaded with nanoparticles will comprise a magnetic susceptibility factors of 10 ahead of individual particles, already shown to be accurately targeted on their own in vivo, for addressing further reduced amounts of magnetic material needed, and in turn a more efficient method of treatment	Future Work <ul style="list-style-type: none">NZW Rabbit in vivo model<ul style="list-style-type: none">Pharmacokinetics of magnetic ECsToxicity / RiskTissue integrationProduct prototyping (magnetic stent, magnetic collagen scaffold)Expanding applications<ul style="list-style-type: none">Card blood stem cellsMagnetic bioabsorbable staples, suturesWound careYour idea here! 
---	--	--	--