Quantitative analysis of hydroxyapatite-binding plasma proteins in genotyped individuals with late-stage age-related macular degeneration

Swati Arya ^a, Eszter Emri ^b, Silvia A. Synowsky ^c, Sally L. Shirran ^c, Neda Barzegar-Befroei^d, Tunde Peto^{b,e}, Catherine H. Botting ^c, Imre Lengyel ^{b,d}, Alan J. Stewart ^{a,*}

Running title: Hydroxyapatite-binding plasma proteins in AMD

Key words: drusen, sub-retinal pigment epithelial deposits, mineral-protein interactions, retinal disease, quantitative proteomics

Abbreviations: AMD, age-related macular degeneration; APOC4, apolipoprotein C4; BCA, bicinchoninic acid; C4A, complement factor C4A; C4B, complement factor C4B; C8A, complement factor C8A; C8B, complement factor C8B; CFH, complement factor H; ELISA, enzyme-linked immunosorbent assay; FDR, false discovery rate; FHR1, factor H-related protein 1; FHR3, factor H-related protein 3; HAP, hydroxyapatite; IC1, plasma protease C1 inhibitor; iRT, indexed retention time; PZP, pregnancy zone protein; RPE, retinal pigment epithelium; SNP, single nucleotide polymorphism; SWATH-MS, sequential window acquisition of all theoretical fragment-ion spectra by mass spectrometry.

^a School of Medicine, University of St Andrews, St Andrews, KY16 9TF, UK.

^b Centre for Experimental Medicine, Queen's University Belfast, Belfast, BT9 7JL, UK.

^c Biomedical Sciences Research Complex, University of St Andrews, St Andrews, KY16 9ST, UK.

^d UCL Institute of Ophthalmology, University College London, London EC1V 9EL

^e Moorfields Eye Hospital, NHS Foundation Trust, London, EC1V 2PD

^{*}To whom correspondence should be addressed: Dr Alan J. Stewart, Tel: +44 (0)1334 463627; E-mail: ajs21@st-andrews.ac.uk

Abstract

Age-related macular degeneration (AMD) is associated with the formation of sub-retinal pigment epithelial (RPE) deposits that block circulatory exchange with the retina. The factors that contribute to deposit formation are not well understood. Recently, we identified the presence of spherular hydroxyapatite (HAP) structures within sub-RPE deposits to which several AMD-associated proteins were bound. This suggested that protein binding to HAP represents a potential mechanism for the retention of proteins in the sub-RPE space. Here we performed quantitative proteomics using Sequential Window Acquisition of all THeoretical fragment-ion spectra-Mass Spectrometry (SWATH-MS) on plasma samples from 23 patients with late-stage neovascular AMD following HAP-binding. Individuals were genotyped for the high risk CFH variant (T1277C) and binding to HAP was compared between wild type and risk variants. From a library of 242 HAP binding plasma proteins (1% false discovery rate), SWATH-MS revealed significant quantitative differences in the abundance of 32 HAPbinding proteins (p<0.05) between the two homozygous groups. The concentrations of six proteins (FHR1, FHR3, APOC4, C4A, C4B and PZP) in the HAP eluted fractions and whole plasma were further analysed using ELISA and their presence in sections from human cadaver eyes was examined using immunofluorescence. All six proteins were found to be present in the RPE/choroid interface, and four of these (FHR1, FHR3, APOC4 and PZP) were associated with spherules in sub-RPE space. This study provides qualitative and quantitative information relating to the degree by which plasma proteins may contribute to sub-RPE deposit formation through binding to HAP spherules and how genetic differences might contribute to deposit formation.

1. Introduction

Age-related macular degeneration (AMD) is the leading cause of visual impairment in the elderly in the developed world (Congdon *et al.*, 2003). One of the hallmarks of AMD is the formation and growth of sub-retinal pigment epithelial (RPE) deposits in the macula that can be focal (drusen) or diffuse (basal linear laminar deposits; Bird *et al.*, 2010; Sarks, 1976) There are several potential sources of accumulating proteins in the sub-RPE space. These include circulatory proteins from the vasculature and RPE secreted proteins. Curcio *et al.* have drawn an analogy between drusen formation in the retina and plaque formation in arterial walls (Curcio *et al.*, 2001). It was also suggested that drusen originates primarily from incomplete digestion of photoreceptor outer segments (Farkas *et al.*, 1971). In fact, the most likely scenario is that sub-RPE deposits originate from a range of sources, and Crabb *et al.* have identified constituents of drusen consistent with this (Crabb *et al.*, 2002). It was, however, not clear how these proteins are retained in the sub-RPE space until it was demonstrated that all sub-RPE deposits appear to contain hydroxyapatite (HAP) spherules that, via their binding capacity to proteins, can serve as a retainer (Thompson *et al.*, 2015)

Based on immunohistochemical analysis of sub-RPE deposits, it was proposed that AMD is associated with complement attack (Hageman *et al.*, 2001). Further support for the involvement of the complement system in AMD was provided by genetic studies that identified a mutation in the complement factor H (CFH) gene that in turn increased the risk to develop AMD by several folds (Edwards *et al.*, 2005; Haines *et al.* 2005; Klein *et al.*, 2005; Toomey *et al.*, 2015). Since then, further genetic associations between the complement system and AMD have been uncovered. However, the mechanism(s) with which CFH or its genetic polymorphisms contribute to AMD is not fully understood although information is emerging (van Asten *et al.*, 2017).

CFH is a component of sub-RPE deposits (Klein *et al.*, 2005), and is associated with HAP spherules (Thompson *et al.*, 2015). It can form multimeric complexes (Perkins *et al.*, 2012) or can undergo zinc-induced self-association with complement C3b (Nan *et al.*, 2013), all of which directly influence the regulation of the complement cascade. As the RPE-expressed CFH appears to be preferentially secreted apically towards to the photoreceptor (Kim *et al.*, 2009) the sub-RPE CFH could be derived from the blood circulation. As other proteins from the circulation appear to contribute to sub-RPE deposit formation (Crabb *et al.*, 2002), it is important to define which of these could be retained in the sub-RPE space due to binding to

HAP. Changes in plasma composition have been associated with early onset drusen formation (Kobayashi *et al.*, 2014), it is thus important to determine whether CFH genotype affects the plasma proteome and thus the HAP interactome.

Here, we quantitatively identify HAP-binding proteins in the plasma of individuals genotyped for the most common AMD-associated CFH polymorphism (T1277C) in late-stage AMD patients by Sequential Window Acquisition of all THeoretical Mass Spectra (SWATH-MS). The concentrations of several proteins that exhibited quantitative differences in HAP eluted fractions were further examined using ELISA and their presence in RPE/choroid interface in sections from human cadaver eyes assessed using immunofluorescence. This study explores the influence of an AMD-associated CFH polymorphism on the plasma abundance of HAP-binding proteins. We believe the presented findings will be useful in determining how genetic factors may contribute to the formation of sub-RPE deposits.

2. Materials and Methods

2.1. Patient information and sample collection

Plasma samples were taken from 23 individuals aged 65-90 with late stage AMD, each displaying drusen and choroidal neovascularisation in clinical images and attending the anti-VEGF injection clinic at Moorfields Eye Hospital, London. Following consent, the haematocrit values were determined before blood collection into K₂EDTA bags (Macopharma, Twickenham, UK) for plasma separation by centrifugation (10 minutes at 1300 RCF at room temperature). The plasma was then aliquoted and frozen immediately at -80C. Red-cupped blood collection tubes were used for DNA isolation (Becton and Dickinson, Plymouth, UK). To determine the genotype of the blood samples PCR (with GoTaq DNA polymerase, Promega, Southampton, UK) and Sanger sequencing (with Big Dye v3.1; Life Technologies, Paisley, UK) was carried out according to manufacturer's instructions. The samples were run on ABI3730 DNA Sequencer. For the proteomic analysis, samples were run with technical duplicates and we included no fewer than 6 biological samples from each of the CFH genotypic groups: homozygous T1277, homozygous C1277 and heterozygous T/C1277. There was no statistically significant difference in age between the participants assigned to each of the three genotypic groups (one-way ANOVA; Kruskal-Wallis test, p=0.1373). The study was approved by Bromley Local Research Ethics Committee (REC reference number: 08/H0805/6) and adhered to the conventions of the Helsinki Declaration.

2.2. Sample preparation for mass spectroscopy

Plasma samples were filtered using 0.45 µm syringe filter (precautionary stto remove particles that could block the columns) and depleted of albumin and IgG with single-use spin columns using the Albumin and IgG Depletion SpinTrap kit (GE Healthcare, Buckinghamshire, UK The resulting plasma was quantitated using BCA/bicinchoninic acid assay (Thermo Fisher Scientific, Paisley, UK) and 600 µg of plasma proteins from each sample werebound to 10 mg of CHT ceramic HAP type-I beads (Bio-Rad, Hertfordshire, UK), in individual microcentrifuge tubes The HAP beads were used once (and not recycled) in order to avoid any cross-contamination. These spherical beads are coated with HAP at their surface providing a similar structure and chemical surface to that present *in vivo* in the human eye (Thompson *et al.*, 2015). They therefore represent an ideal medium for the *in vitro* identification of HAP binding proteins. Plasma protein binding was carried out in a binding buffer composed of 10 mM sodium phosphate, 75 mM sodium chloride, pH 7.4. Bound proteins were eluted in 250 mM sodium phosphate, 10 mM Tris-Cl, 0.5 M sodium chloride, pH 7.4. The resulting proteins were quantified using a bicinchoninic acid assay.

All samples for mass spectrometric analysis were processed as follows: Samples containing 30 µg of protein were denatured in 6 M urea in 50 mM ammonium bicarbonate and were then reduced and alkylated with 1 mM tris (2-carboxyethyl) phosphine and 5 mM iodoacetamide. The reaction was quenched with 10 mM dithiothreitol. Samples were diluted with 50 mM ammonium bicarbonate to a final urea concentration of 1.5 M. The resulting samples were then digested with trypsin (1:50 ratio (w/w), 0.2 µg/µl trypsin; Promega, Southampton, UK), overnight at 30°C. The digestion was stopped upon addition of 0.5% (v/v) trifluoroacetic acid. Peptides were desalted using a C18 SepPak cartridge (Thermo Fisher Scientific) and dried under vacuum.

2.3 Protein identification after HAP binding by LC-ESI-MSMS

Albumin/IgG-depleted plasma proteins were analysed after HAP binding. For this, 1 µg of the respective sample peptides were injected on an AB Sciex TripleTOF 5600+ system mass spectrometer (Sciex, Framingham, MA, USA) coupled to an Eksigent nanoLC AS-2/2Dplus system (Sciex). Prior to mass spectrometric analysis, iRT reference peptides (Biognosys, Schlieren, Switzerland) were added to each sample according to the manufacturer's

specifications. The samples were prepared in loading buffer (2% acetonitrile, 0.05% trifluoroacetic acid), bound to an Acclaim Pepmap 100 µm × 2 cm trap (Thermo Fisher Scientific) and washed for 10 min to waste after which the trap was turned in-line with the analytical column (Acclaim Pepmap RSLC 75 µm × 15 cm). The analytical solvent system consisted of buffer A (2% acetonitrile, 0.1% formic acid in water) and buffer B (2% water, 0.1% formic acid in acetonitrile) at a flow rate of 0.3 µl/min with the following gradient: linear 1-20% of buffer B over 90 min, linear 20-40% of buffer B for 30 min, linear 40-99% of buffer B for 10 min, isocratic 99% of buffer B for 5 min, linear 99-1% of buffer B for 2.5 min and isocratic 1% solvent buffer B for 12.5 min. The mass spectrometer was operated in DDA top 20 positive ion mode, with 250 and 150 ms acquisition time for the MS1 (m/z 400-1200) and MS2 (m/z 230-1800) scans respectively, and 15 s dynamic exclusion. Rolling collision energy with a collision energy spread of 5 eV was used for fragmentation. The data files were searched using Protein Pilot v5.0.1 (Sciex) with the following search parameters: urea denaturation as special factors, trypsin as the cleavage enzyme and carbamidomethylation as a fixed modification of cysteines. The search was carried out in "rapid ID" mode with a detected protein threshold of 10% plus false discovery rate analysis against the Swissprot database downloaded June 2015, restricted only to proteins from humans. Note that the iRT peptides were added to this database.

2.4. Sequential window acquisition of all theoretical fragment-ion spectra-mass spectrometry (SWATH-MS)

2.4.1 LC-ESI-MS/MS analysis for spectral library generation

For spectral library generation, 1 µg of peptides from each individually digested HAP eluate sample from each patient were combined to be a representative sample. This sample was analysed on an AB Sciex TripleTOF 5600+ system mass spectrometer coupled to an Eksigent nanoLC AS-2/2Dplus system as described in Section 2.3. The data was searched against three databases resulting in three different libraries. The first was the original Swissprot database (downloaded June 2015), restricted only to proteins from humans. This database contains FHR1 in its original form (with the amino acid sequence, STDTSCVNPPTVQNAHILSR at position 142 to 161; FHR1a). The second database is identical to this original Swissprot database but contains FHR1b instead of FHR1a (with the amino acid sequence, STDTSCVNPPTVQNAYIVSR at position 142 to 161). The third database was the original Swissprot database but the C-terminus but the protein sequence for FHR1a was extended to

include the following sequence: STDTSCVNPPTVQNAYIVSR (taken from FHR1b). The sequences of the iRT peptides were added to each database.

2.4.2 SWATH-MS data acquisition and analysis

For SWATH-MS data acquisition, the same mass spectrometer and LC-MS/MS setup as described above was used, but operated in SWATH mode. The method uses 50 windows of variable Da effective isolation width with a 1 Da overlap using AB Sciex Variable Window Calculator tool. Each window has a dwell time of 150 ms to cover the mass range of 400-1250 m/z in TOF-MS mode and MS/MS data is acquired over a range of 230-1800 m/z with high sensitivity setting and a dwell time of 70 ms, resulting in a cycle time of 3.6 s. The collision energy for each window was set using the collision energy of a 2⁺ ion centred in the middle of the window with a spread of 5 eV. During data processing identified proteins with a 1% false detection rate (FDR) at the global protein level (242 proteins) were imported into PeakView 2.2 (Sciex) for spectral library generation. SWATH-MS results were analysed using SWATH micro App (v2.0) embedded in PeakView. Sample peptide retention times were calibrated against the iRT peptide retention times. The data was then further processed with the following settings: 10 peptides/protein, 5 transitions/peptide, peptide confidence threshold of 95%, FDR threshold of 1%. Modified and shared peptides were excluded. The XIC extraction window was set to 6 min and 10 ppm. The quantified protein peak areas were then exported into Markerview 1.2.1.1 (Sciex) for data normalisation. Relative protein quantitation and statistical tests were performed using Microsoft Excel. Basic statistical calculations were performed using Markerview software. The p-values were calculated using students t-test and the q-values were calculated using the Benjamini Hochberg method. The mass spectrometry proteomics data is available at the ProteomeXchange Consortium, PRIDE with the dataset identifier PXD007944.

2.5. STRING analysis

The protein interaction network for significantly changing HAP binding plasma proteins (between the homozygous CFH genotyped groups) was analysed by STRING software (https://string-db.org/). As a reference dataset we used the human plasma proteome from the Plasma Proteome Database's official website (http://plasmaproteomedatabase.org/).

2.6. Enzyme-linked immunosorbent assays (ELISAs)

ELISAs were performed on patient plasma prior to and following HAP bead binding and elution. Patient plasma was filtered through 0.45 μm filter and incubated with and eluted from 10 mg of CHT ceramic HAP beads as described above for SWATH-MS preparation. The ELISAs were performed using commercially available kits (all from Cloud-Clone Corp., Buckingham, UK) for pregnancy zone protein (PZP, #SEG324Hu), complement factor H-related protein 1(FHR1, #SEL327HU), complement factor H-related protein 3 (FHR3, #SEL329Hu), apolipoprotein C4 (APOC4, #SEB828Hu), complement factor 4A (C4A, #SEA389Hu) and complement factor 4B (C4B, #SEB305Hu) in accordance with the manufacturer's instructions. Statistical comparisons (p-values) were done by one-way ANOVA using Dunnet's multiple comparison test.

2.7. Immunofluorescence imaging

Tissues were obtained from the Eye Depository at the UCL Institute of Ophthalmology, London. Approval for the use of these tissues was obtained from the UCL Institute of Ophthalmology Local Research Ethics Committee (REC reference number: 10/H0106/57-2012ETR26). The donor eyes were fixed <24 hours after death. Samples were selected based on the presence of drusen from male and female donors, aged between 76 and 94 years, without clinical and genotype information. Four micron sections containing several drusen were selected after microscopic evaluation of formalin-fixed paraffin embedded tissues following hematoxilin-eosin staining. In addition, adjacent sections were labelled with von Kossa (Puchtler and Meloen, 1978) and OsteoSense 680EX (Figueiredo *et al.*, 2008) to verify the presence of mineral in the deposits.

Antigen retrieval was carried out in 10 mM citrate buffer (pH 6). The sections were heated for 9 min in a microwave following de-paraffinisation. The samples were blocked with 20% goat sera for 1 hr at room temperature. Afterwards, the sections from 3-5 donors were incubated with primary antibodies overnight at 4°C; C4A (ab170942, Abcam, Cambridge, UK, 1:50 dilution); C4B (ab181241, Abcam, 1:250 dilution); FHR1 (ab182652, Abcam, 1:50 dilution); APOC4 (ab199163, Abcam, 1:50 dilution); PZP (Abcam: ab122718, 1:20 dilution); FHR3 (16583-1-AP, Proteintech, Manchester, UK, 1:100 dilution); amyloid-β 1-16 (SIG-39320, Covance, London, UK, 1:100) followed by incubation with goat anti-rabbit or goat anti-mouse secondary antibody conjugated to alexa-fluor 488 or 568, respectively (Life Technologies, Paisley, UK), at 1:200 dilution for 1 hr at room temperature. Slides were counterstained for nuclei with DAPI (Thermo Fisher Scientific, 1:1000 dilution) for 15 min at

room temperature and mounted with antifade mounting media (Vector Laboratories, Peterborough, UK) using glass coverslips (Thermo Fisher Scientific). Optical sections were generated using a Leica TCS SP8 confocal laser scanning microscope (Leica Microsystems, Milton Keynes, UK). The images presented are in all cases representative.

3. Results

3.1. Quantitative analysis of HAP-binding plasma proteins in individuals with late-stage AMD

For SWATH-MS quantitation, a sample specific reference spectral library was generated (for use as a reference in generating peptide query parameters for the peptide centric analysis of the data) by data-dependent acquisition analysis of the HAP-binding plasma proteins. This resulted in a spectral library consisting of 242 human plasma proteins with a FDR of 1% on the protein level. To determine the differences in levels of HAP-binding proteins present across the different CFH genotypic groups, we measured quantitative proteome profiles of HAP-eluted plasma proteins in each individual sample by SWATH-MS. Comparisons of the obtained quantitative data identified 210 proteins (1% FDR) that were consistently detected in all samples. The quantitative information for the 174 proteins that were quantified based on 2 or more peptides is provided in Table S2. To account for the probability of minor sample variability due to multiple steps in sample processing, all the samples were prepared in duplicates and ran as technical replicates. A principal component analysis of the proteomic data obtained from individual samples (PC1=13.6%, PC2=7.9%) revealed that the samples did not vary greatly as a collective, with technical replicates clustering together closely (indicating strong reproducibility) as shown in Figure S1. A volcano plot showing a comparison of relative abundance of all HAP-binding plasma proteins between T/T1277 and C/C1277 groups is presented in Figure 1A. Of the proteins consistently quantified, the relative abundance of a total of 32 unique proteins (14%) was statistically significantly altered (p-value ≤0.05) between T/T1277 and C/C1277 groups. Multiple hypothesis testing was also performed with the respective q-values calculated. However, these were too stringent for this study and as such were not used in determining significance. An analysis of the STRING database for these 32 proteins (shown in Figure 1B) indicates that they have more interactions among themselves than would be expected for a random set of 32 proteins drawn from the proteome. Such enrichment indicates that the proteins might be, at least partially, biologically connected as a group. The quantitative data for all significantly altered proteins is presented in Figure 1C.

Twenty-one plasma proteins were found to exhibit a statistically significant increase ($p \le 0.05$) in mean relative abundance (represented as average peak intensity in Table S2) in the C/C1277 group (compared to the T/T1277 group). Of those displaying the greatest increase were IgG delta chain C region (~2.8-fold change; p=0.04), thrombospondin-I (~2.2-fold change; p=0.0342), FHR1 (~2.1-fold change; p=6.58 \times 10⁻⁸), platelet basic protein (~2-fold change; p=0.04). The observed increased abundance of FHR1 protein in the C/C1277 group was interesting as a variant (FHR1a) of the corresponding gene is associated with AMD, and is genetically linked to the C1277 CFH allele (Timmann et al., 1991; Martínez-Barricarte et al., 2012). FHR1a encodes a protein, which is the same length as (the other variant) FHR1b but differs in amino acid sequence at 3 residues (H157Y, L159V and E175Q). Examination of FHR1-derived peptides in the SWATH-MS spectra revealed that FHR1a and FHR1b variants could be identified within the samples based upon the presence of a peptide unique to each form. This was achieved by generating separate peptide libraries containing peptides derived from either FHR1a or FHR1b variant proteins. The presence of FHR1a and FHR1b variants within the C/C1277 and T/T1277 groups were analysed by quantitative mass spectrometric analysis as shown in Figure 2. The p-values were calculated by one-way ANOVA using Dunnet's multiple comparison test. The results confirmed strong associations between the C1277 CFH allele and FHR1a and the T1277 and FHR1b, respectively.

Eleven plasma proteins found to exhibit a statistically significant increase (p≤0.05) in mean relative abundance in the T/T1277 group (compared to the C/C1277 group). Of those displaying the greatest increase were PZP (~2-fold change; p=0.01), APOC4 (~1.7-fold change; p=0.0118), FHR3 (~1.7-fold change; p=0.001), C4B (~1.7-fold change; p=0.01), vitamin K-dependent protein Z (~1.6-fold change; p=0.03) and C4A (~1.6-fold change; p=0.01). Although PZP exhibited the greatest increase in mean relative abundance in the T/T1277 group, the protein has been shown previously to be expressed at higher levels in females when compared to levels in males (Folkersen *et al.*, 1981). To determine whether this observation may have been due to sex differences in individuals within the groups, we carried out a separate analysis of the SWATH-MS data where samples from females were compared with those from males (as shown in Table S3). This analysis suggested that PZP may only be increased in the T/T1277 group as it is more prevalent in females as there was a bias toward females in this genetic group in our study. The abundance of six other proteins (protein S100-A9, Ig alpha-2 chain C region, cystatin C, Ig kappa chain V-III region POM, complement

factor D and carboxypeptidase B2) was also found to significantly differ between male and female participants (in all these cases levels were higher in males) but did not vary significantly across the three CFH genotypic groups.

3.2. Determination of concentrations of plasma proteins exhibiting quantitative differences between genotypic groups

Several proteins that exhibited differing mean relative abundances in the different CFH genotypic groups were chosen for further analysis – these were PZP, FHR1, FHR3, APOC4, C4A and C4B. Although APOC4 was quantified based on single peptide we thought it also merited further analysis given that it displayed a statistically significant p-value (0.01) and relatively high fold-change (1.7), in the T/T1277 vs C/C1277 comparison and a confident signal in SWATH-MS. To validate the findings from the SWATH-MS analysis and to determine whether differences in relative abundance of these proteins between samples from different genotypic groups were due to affinity/competition effects of HAP-binding or their plasma concentration prior to binding, ELISA assays were carried out both on untreated plasma samples and the HAP-eluted samples subjected to the SWATH-MS analysis. The resultant data is shown in Figure 3. Statistically significant differences in the concentrations of each of the six proteins were observed in both the native plasma and processed samples, suggesting that the differences in abundance of these proteins between groups were most likely due to initial differences in the plasma concentrations of these proteins prior to processing.

3.3. Immunohistochemistry of identified plasma proteins in cadaveric retinal tissue

Human cadaveric eye sections pre-screened according to Thompson *et al.* (2015) contained several drusen and HAP spherules were immunolabelled for the detection of C4A, C4B, APOC4, PZP, FHR3 and FHR1 (Figure 4). For visualisation of HAP spherules, we used amyloid- β immunolabelling as the acidity of the citrate buffer used for antigen retrieval inhibited HAP labelling. The amyloid- β antibody 6E10 reliably co-localises with some of the HAP spherules (Dentchev *et al.*, 2003; Anderson *et al.*, 2004; Luibl *et al.*, 2006; Thompson *et al.*, 2015). FHR1, FHR3, PZP and APOC4 showed distinct spheroid shape immunostaining that co-localised with amyloid- β (Figure 4A-F). C4A and C4B immunolabelling was localised in most cases to choriocapillary pillars (Figure 4G-I). Complement C4A was also observed in drusen, but did not appear in association with amyloid- β (Fig 4G). APOC4, PZP and FHR3 were also detected in RPE cells (Figure 4B-F) and FHR1 and PZP were present in

the Bruch's membrane (Figure 4A and E). In addition, FHR1 immunolabelling was observed in the choriocapillary pillars (Figure 4A).

4. Discussion

Several reports have described circulatory plasma biomarkers for AMD (Gu et al., 2010; Lambert et al., 2016; Kersten et al., 2017). However, how these contribute to the development and progression of the AMD is not well understood. Plasma proteins appear to contribute to the development and growth of sub-RPE deposits (Crabb et al., 2002) and as sub-RPE deposit accumulation is a hallmark of AMD (Bird, 2003), understanding how plasma proteins may contribute to deposit formation could be important to identify early events in AMD pathogenesis. In this study, we examined the relationship between genotype and abundance of HAP-binding plasma proteins that are potentially involved in deposit formation. We made two new and significant observations: 1) at least 242 plasma proteins can readily bind HAP (either directly or, in some cases, via association with another HAP-binding protein); 2) genetic variations in CFH can influence the plasma concentration of some of these proteins.

In our SWATH-MS analysis 14% of quantified plasma proteins showed statistically significant difference when CFH genotyped samples (T/T1277 and C/C1277) were compared (32 out of 242). A STRING analysis of the 32 proteins revealed that many of these proteins are biologically connected: for example, we found 8 complement or complement-related proteins as well as 2 apolipoproteins (APOC4 was not included) amongst these proteins. Such connections might be relevant because they align with the concept that complement dysregulation and changes in lipid homeostasis are mechanistically involved in AMD pathogenesis. In addition, protein-protein interactions could lead to the growth of sub-RPE deposits.

The HAP-binding protein exhibiting the most significant relative decrease in samples from the T/T1277 compared to C/C1277 group, was FHR1 (fold–change of 0.47; p-value ≥6.58×10⁻⁸). An ELISA-based analysis confirmed this finding in both whole plasma samples and those eluted from HAP beads. FHR1 is a CFH-related protein that exists in two polymorphic forms, FHR1a and FHR1b derived from distinct FHR1 alleles. Further inspection of the SWATH-MS data revealed that a peptide specific to the FHR1a form was found in all samples derived from C/C1277 subjects and absent in all but one sample derived

from the T/T1277 subjects. Similarly, a peptide specific to the FHR1b protein was found in all T/T1277 derived samples and absent in all but one sample derived from the C/C1277 subjects. This is the first time that such an association has been shown at the protein level, supporting a previous study suggesting a link between the corresponding CFH and FHR1 genotypes (Martinez-Barricante *et al.*, 2012). Interestingly, deletion of the gene encoding FHR1 together with FHR3, which we conversely found to be more abundant in the T/T1277 subjects, is associated with a reduced risk of developing AMD. The specific role(s) these proteins play in the aetiology of AMD is unclear, but it is thought that they may function to modulate CFH activity in vivo (Fritsche *et al.*, 2010). In addition, our study demonstrates that both FHR1 and FHR3 could bind to HAP and associate with spherules at the RPE/choroid interface.

In our analysis of the HAP-eluted plasma proteins, significantly higher plasma levels of both C4A and C4B (1.68 and 1.57-fold, respectively) were observed in individuals with the T/T1277 CFH genotype than those with the C/C1277 genotype, with this trend confirmed by ELISA both after HAP-binding and in whole plasma. This is an interesting finding as it suggests that an allelic difference in CFH, which encodes a protein associated with the alternative complement pathway, influences the plasma level of C4, which is essentially a classical pathway component. In eye sections both C4A and C4B were found to be present in the choriocapilliary pillars in every section we examined, and C4A occasionally associated with drusen but neither were observed to be associated with HAP spherules in the sections we examined. Examining further sections may clarify whether HAP-association exists. In a previous quantitative proteomics analysis of Bruch's membrane/choroid complex from patients with AMD, it was revealed that C4A and C4B are present in this complex and that C4A was statistically significantly more prevalent in the AMD group compared to controls (Yuan et al., 2010).

Another protein that exhibited higher plasma levels in T/T1277 individuals was plasma protease C1 inhibitor (IC1; 1.4-fold), a negative regulator of the classical pathway. IC1 is a glycoprotein inhibitor of activated complement factor 1r and 1s proteins, which form part of the first complement component (C1). Deficiency of IC1 is associated with hereditary angioneurotic oedema (Cugno *et al.*, 2009). It is possible that the relatively lower abundance of either IC1 or C4 proteins (in C/C1277 group relative to T/T1277 group) could be interdependent and that the change in concentration of either occurs to maintain balance of

the classical pathway. Other components of the complement cascade, C8 (both C8A and C8B proteins) and complement factor I (negative regulator of alternative pathway), were also found to be higher in the T/T1277 group, albeit by less than 1.3-fold. We hasten to add that though these results are interesting and may have functional implications, it is not possible based on this data alone to relate these observations directly to complement dysregulation in AMD

Defective lipid metabolism and genotypic differences in genes encoding proteins involved in lipid transport have also been associated with the risk of developing AMD (Cheung *et al.*, 2017). Such genes include those encoding apolipoproteins such as ApoE and ApoB (Klaver *et al.*, 1998, Baird *et al.*, 2006, Curcio *et al.*, 2010). In our SWATH-MS analysis we observed a higher abundance of APOC4 and ApoE in the HAP-eluted protein fractions derived from the T/T1277 group, relative to the C/C1277 group (1.78-fold and 1.24-fold, respectively). The difference in APOC4 levels (in the HAP-eluted fractions) between T/T1277 and C/C1277 groups were also confirmed by ELISA, with levels of APOC4 in whole plasma also being significantly higher in the T/T1277 group. In retinal sections APOC4 was found to colocalise with spherular structures and was also found to be expressed in RPE cells. Associations between APOC4 and AMD have been reported previously based on genomewide candidate gene association studies (Pennington and DeAngelis, 2016). It is therefore possible that this protein may be involved in AMD pathogenesis by binding to HAP in the sub-RPE space.

We also identified new sub-RPE deposit-associated proteins that might be relevant for AMD progression. The protein that exhibited the greatest increase in relative abundance in the T/T1277 group compared to the C/C1277 group was PZP, an α2-macroglobulin-like endopeptidase inhibitory protein recently implicated in Alzheimer's disease and identified within amyloid plaques in cadaveric brain tissue (Ijsselstijn *et al.*, 2011; Nijholt *et al.*, 2015). However, PZP was found to be higher in samples taken from females rather than males, which could account for this difference (as the T/T1277 group consisted exclusively of female subjects), although it should be noted that females are more predisposed to developing AMD than males (Owen *et al.*, 2012; Rudnicka *et al.*, 2012). PZP has not previously been associated with sub-RPE deposit formation or AMD. Here we found that the protein is present in drusen and associated with spherular structures. PZP was also present in the Bruch's membrane and in the RPE. These observations support the hypothesis that PZP may

contribute to the growth of deposits, thickening of Bruch's membrane and thus the progression of AMD.

A number of proteins present in plasma, including complement factors, have already been identified in drusen (Crabb *et al.*, 2002; Hollyfield *et al.*, 2003), with some of these, such as haemoglobin, specific to the blood (Crabb *et al.*, 2002; Beattie *et al.*, 2010). These data together with our findings collectively suggest that plasma proteins are likely to contribute to sub-RPE deposit formation. The differences in concentrations of potential deposit-forming plasma proteins in AMD-affected individuals of different genotypes could points towards different aetiologies of disease. For example, based on our data, FHR1 may play a more prominent role in sub-RPE deposit formation in individuals with the C/C1277 CFH genotype than those with the T/T1277 genotype.

In conclusion, this study has utilised a new quantitative proteomic approach to understand the potential molecular events that give rise to sub-RPE deposit formation, a feature central to the development of AMD. More specifically, we have identified a number of a plasma proteins that can bind HAP and are differentially expressed in plasmas from individuals belonging to different genotypic groups. These observations point towards different aetiologies of sub-RPE deposit formation and heterogeneity and why people from different genetic backgrounds develop AMD. Given that the retention of proteins in the sub-RPE space occurs at a very early stage of the disease process, targeting the binding of key proteins to HAP may provide a robust therapeutic strategy that could slow or even prevent drusen formation.

Contributors

Study concept and design: S.A., E.E., C.H.B., I.L., A.J.S.; collection of blood samples: N.B.-B., T.P.; acquisition of data: S.A., E.E., S.A.S., S.L.S.; analysis and interpretation of data: S.A., E.E., S.A.S., C.H.B., I.L., A.J.S.; drafting of the manuscript: S.A., E.E., S.A.S., I.L., A.J.S.; and all authors critically reviewed the manuscript and approved the final version.

Competing Interests

The authors have no competing interests to declare.

Acknowledgements

The investigators would like to thank the patients for their donation of blood samples for the study. Yash Bansod and Lasse Olesen are acknowledged for their technical help with the immunohistochemistry examination. We are also grateful to Fight for Sight for financial support (project grant to A.J.S and I.L.; grant ref.: 1586/1587). This research was also supported by "Eye-Risk" a European Union's Horizon 2020 research and innovation program (grant ref.: 634479) Bill Brown Charitable Trust, MEH Special Trustees and Mercer Fund (I.L). This work was also supported by the Wellcome Trust (grant ref.: 094476/Z/10/Z) for funding the purchase of the TripleTOF 5600 mass spectrometer at the BSRC Mass Spectrometry and Proteomics Facility, University of St Andrews.

References

- Anderson, D.H., Talaga, K.C., Rivest, A.J., Barron, E., Hageman, G., Johnson, L.V., 2004. Characterization of β amyloid assemblies in drusen: the deposits associated with aging and age-related macular degeneration. Exp. Eye Res. 78, 243-256.
- Baird, P.N., Richardson, A.J., Robman, L.D., Dimitrov, P.N., Tikellis, G., McCarty, C.A., Guymer, R.H., 2006. Apolipoprotein (APOE) gene is associated with progression of age-related macular degeneration (AMD). Hum. Mutat. 27, 337-342.
- Beattie, J.R., Pawlak, A.M., Boulton, M.E., Zhang, J., Monnier, V.M, McGarvey, J.J., Stitt, A.W., 2010. Multiplex analysis of age-related protein and lipid modifications in human Bruch's membrane. FASEB J. 24, 4816-4824.
- Bird, A.C., 2003. The Bowman lecture. Towards an understanding of age-related macular disease. Eye 17, 457-466.
- Bird, A.C., 2010. Therapeutic targets in age-related macular disease. J. Clin. Invest. 120, 3033–3041.
- Cheung, C.M.G., Gan, A., Fan, Q., Chee, M.L., Apte, R.S., Khor, C.C., Yeo, I., Mathur, R., Cheng, C.Y., Wong, T.Y., 2017. Plasma lipoprotein subfraction concentrations are associated with lipid metabolism and age-related macular degeneration. J. Lipid Res. 58, 1785-1796.
- Congdon, N.G., Friedman, D.S. and Lietman, T., 2003. Important causes of visual impairment in the world today. JAMA 290, 2057-2060
- Crabb, J.W., Miyagi, M., Gu, X., Shadrach, K., West, K.A., Sakaguchi, H., Kamei, M., Hasan, A., Yan, L., Rayborn, M.E., Salomon, R.G., Hollyfield, J.G., 2002. Drusen proteome analysis: an approach to the etiology of age-related macular degeneration. Proc. Natl. Acad. Sci. USA. 99, 14682-14687.

- Cugno, M., Zanichelli, A., Foieni, F., Caccia, S., Cicardi, M., 2009. C1-inhibitor deficiency and angioedema: molecular mechanisms and clinical progress. Trends Mol. Med. 15, 69-78.
- Curcio, C.A., Johnson, M., Huang, J.D., Rudolf, M., 2010. Apolipoprotein B-containing lipoproteins in retinal aging and age-related macular degeneration. J. Lipid Res. 51, 451–467.
- Curcio, C.A., Millican, C.L., Bailey, T., Kruth, H.S., 2001. Accumulation of cholesterol with age in human Bruch's membrane. Invest. Ophthalmol. Vis. Sci. 42, 265-274.
- Dentchev, T., Milam, A.H., Lee, V.M., Trojanowski, J.Q., Dunaief, J.L., 2003. Amyloid-beta is found in drusen from some age-related macular degeneration retinas, but not in drusen from normal retinas. Mol. Vis. 9, 184-190.
- Edwards, O.A., Ritter, R., Abel, K.J., Manning, A., Panhuysen, C., Farrer, L.A., 2005. Complement factor H polymorphism and age-related macular degeneration. Science 308, 421-424.
- Farkas, T.G., Sylvester, V., Archer, D., 1971. The ultrastructure of drusen. Am. J. Ophthalmol. 71, 1196-1205.
- Figueiredo, J.L., Passerotti, C.C., Sponholtz, T., Nguyen, H.T., Weissleder R., 2008. A novel method of imaging calcium urolithiasis using fluorescence. J. Urol. 179, 1610-1614.
- Folkersen, J., Teisner, B., Grunnet, N., Grudzinskas, J.G., Westergaard, J.G., Hindersson, P., 1981. Circulating levels of pregnancy zone protein: normal range and the influence of age and gender. Clin. Chim. Acta. 110, 139-145.
- Fritsche, L.G., Lauer, N., Hartmann, A., Stippa, S., Keilhauer, C.N., Oppermann, M., Pandey, M.K., Köhl, J., Zipfel, P.F., Weber, B.H., Skerka, C., 2010. An imbalance of human complement regulatory proteins CFHR1, CFHR3 and factor H influences risk for agerelated macular degeneration (AMD). Hum. Mol. Genet. 19, 4694-4704.
- Gu, J., Pauer, G.J., Yue, X., Narendra, U., Sturgill, G.M., Bena, J., Gu, X., Peachey, N.S., Salomon, R.G., Hagstrom, S.A., Crabb, J.W., Clinical Genomic and Proteomic AMD Study Group., 2010. Proteomic and genomic biomarkers for age-related macular degeneration. Adv. Exp. Med. Biol. 664, 411-417.
- Hageman, G.S., Luthert, P.J., Victor, C.N.H., Johnson, L.V., Anderson, D.H., Mullins, R.F., 2001. An integrated hypothesis that considers drusen as biomarkers of Immune-mediated processes at the RPE-Bruch's membrane interface in aging and age related macular degeneration. Prog. Retin. Eye Res. 20, 705-732.

- Haines, J.L., Hauser, M.A., Schmidt, S., Scott, W.K., Olson, L.M., Gallins, P., Spencer, K.L.,
 Kwan, S.Y., Noureddine, M., Gilbert, J.R., Schnetz-Boutaud, N., Agarwal, A., Postel,
 E.A., Pericak-Vance, M.A., 2005. Complement factor H variant increases the risk of
 age-related macular degeneration. Science 308, 419-421
- Hollyfield, J.G., Salomon, R.G., Crabb, J.W., 2003. Proteomic approaches to understanding age-related macular degeneration. Adv. Exp. Med. Biol. 533, 83-89.
- Ijsselstijn, L., Dekker, L.J., Stingl, C., van der Weiden, M.M., Hofman. A., Kros, J.M., Koudstaal, P.J., Sillevis Smitt, P.A., Ikram, M.A., Breteler, M.M., Luider, T.M., 2011. Serum levels of pregnancy zone protein are elevated in presymptomatic Alzheimer's disease. J. Proteome Res. 10, 4902-4910.
- Kersten, E., Paun, C.C., Schellevis, R.L., Hoyng, C.B., Delcourt, C., Lengyel, I., Peto, T., Ueffing, M., Klaver, C.C.W., Dammeier, S., den Hollander, A., de Jong, E.K., 2017. Systemic and ocular fluid compounds as potential biomarkers in age-related macular degeneration. Surv. Ophthalmol. 6257, 30262-30264.
- Kim, Y.H., He, S., Kase, S., Kitamura, M., Ryan, S.J., Hinton, D.R., 2009. Regulated secretion of complement factor H by RPE and its role in RPE migration. Graefes Arch. Clin. Exp. Ophthalmol. 247, 651-659.
- Klaver, C.C., Kliffen, M., van Duijn, C.M., Hofman, A., Cruts, M., Grobbee, D.E., van Broeckhoven, C., de Jong, P.T., 1998. Genetic association of apolipoprotein E with age-related macular degeneration. Am. J. Hum. Genet. 63, 200-206.
- Klein, R.J., Zeiss, C., Chew, E.Y., Tsai, J.Y., Sackler, R.S., Haynes, C., Henning, A. K., SanGiovanni, P., Mane, S.M., Mayne, S.T., Bracken, M.B., Ferris, F.L., Ott, J., Barnstable, C., Hoh, J., 2005. Complement factor H polymorphism in age-related macular degeneration. Science 308, 385–389.
- Kobayashi, H., Okamoto, H., Murakami, A., Iwata, T., 2014. Plasma proteome analysis on cynomolgus monkey (*Macaca fascicularis*) pedigrees with early onset drusen formation. Exp. Anim. 63, 305-310.
- Lambert, N.G., El Shelmani, H., Singh, M.K., Mansergh, F.C., Wride, M.A., Padilla, M., Keegan, D., Hogg, R.E., Ambati, B.K., 2016. Risk factors and biomarkers of agerelated macular degeneration. Prog. Retin. Eye Res. 54, 64-102.
- Luibl, V., Isas, J.M., Kayed, R., Glabe, C.G., Langen, R., Chen, J., 2006. Drusen deposits associated with aging and age-related macular degeneration contain nonfibrillar amyloid oligomers. J. Clin. Invest. 116, 378-385.

- Martínez-Barricarte, R., Recalde, S., Fernández-Robredo, P., Olavarrieta, I.M., Viñuela, A., Pérez-Pérez, J., Layana, A.G., Rodríguez de Córdoba, S., Spanish Multicenter Group on AMD, 2012. Relevance of complement factor H–related 1 (CFHR1) genotypes in age-related macular degeneration. Invest. Ophthalmol. Vis. Sci. 53, 1087-1094.
- Nan, R., Tetchner, S., Rodriguez, E., Pao, P.J., Gor, J., Lengyel, I., Perkins, S.J., 2013. Zinc-induced self-association of complement C3b and factor H: implications for inflammation and age-related macular degeneration. J. Biol. Chem. 288, 19197-19210.
- Nijholt, D.A., Ijsselstijn, L., van der Weiden, M.M., Zheng, P.P., Sillevis, Smitt, P.A., Koudstaal, P.J., Luider, T.M., Kros, J.M., 2015. Pregnancy zone protein is increased in the Alzheimer's disease brain and associates with senile plaques. J. Alzheimers Dis. 46, 227-238.
- Owen, C.G., Jarrar, Z., Wormald, R., Cook, D.G., Fletcher, A.E., Rudnicka, A.R., 2012. The estimated prevalence and incidence of late stage age related macular degeneration in the UK. Br. J. Ophthalmol. 96, 752-756.
- Pennington, K.L., DeAngelis, M.M., 2016. Epidemiology of age-related macular degeneration (AMD): associations with cardiovascular disease phenotypes and lipid factors. Eye Vis. 3, 34-54.
- Perkins, S.J., Nan, R., Li, K., Khan, S., Miller, A., 2012. Complement factor H-ligand interactions: self-association, multivalency and dissociation constants. Immunobiology 217, 281-297.
- Puchtler, H., Meloen, S.N., 1978. Demonstration of phosphates in calcium deposits: a modification of von Kossa's reaction. Histochemistry 56, 177-185.
- Rudnicka, A.R., Jarrar, Z., Wormald, R., Cook, D.G., Fletcher, A., Owen, C.G., 2012. Age and gender variations in age-related macular degeneration prevalence in populations of European ancestry: a meta-analysis. Ophthalmology 119, 571-580.
- Sarks, S.H., 1976. Ageing and degeneration in the macular region: A clinico-pathological study. Br. J. Ophthalmol. 60, 324–341.
- Timmann, C., Leippe, M., Horstmann, R. D., 1991. Two major serum components antigenically related to complement factor H are different glycosylation forms of a single protein with no factor H-like complement regulatory functions. J. Immunol. 146, 1265–1270.
- Thompson, R.B., Reffatto, V., Bundy, J.G., Kortvely, E., Flinn, J.M., Lanzirotti, A., Jones, E.A., McPhail, D.S., Fearn, S., Boldt, K., Ueffing, M., Ratu, S.G., Pauleikhoff, L., Bird, A.C., Lengyel, I., 2015. Identification of hydroxyapatite spherules provides new insight

- into subretinal pigment epithelial deposit formation in the aging eye. Proc. Natl. Acad. Sci. USA. 112, 1565-1570.
- Toomey, C.B., Kelly, U., Saban, D.R., Bowes Rickman, C., 2015. Regulation of age-related macular degeneration-like pathology by complement factor H. Proc. Natl. Acad. Sci. USA. 112, 3040-3049.
- van Asten, F., Simmons, M., Singhal, A., Keenan, T.D., Ratnapriya, R., Agron, E., Clemons, T.E., Swaroop, A., Lu, Z., Chew, E.Y., Age-Related Eye Disease Study 2 Research Group. A Deep Phenotype Association Study Reveals Specific Phenotype Associations with Genetic Variants in Age-related Macular Degeneration: Age-Related Eye Disease Study 2 (AREDS2) Report No. 14. Ophthamology, doi: 10.1016/j.ophtha.2017.09.023.
- Yuan, X., Gu, X., Crabb, J.S., Yue, X., Shadrach, K, Hollyfield, J.G., Crabb, J.W., 2010. Quantitative proteomics: comparison of the macular Bruch membrane/choroid complex from age-related macular degeneration and normal eyes. Mol. Cell. Proteomics 9, 1031-1046.

Figure Legends


Figure 1. Quantitative analysis of HAP-binding plasma proteins examined by SWATH-

MS. (A) Volcano plot of all quantified proteins displaying differences in relative abundance between CFH T/T1277 (TT) and C/C1277 (CC) groups. The plot represents the negative log₁₀ of the p-value against the log₂ of the fold change (FC). FC was calculated based on the relative abundance in TT/CC. Proteins shown by the open circles exhibited statistically significant differences between the homozygous groups (p<0.05) following binding and elution from HAP beads. The six proteins selected for further investigation are shown in bold. (B) The protein interaction network for 32 significantly changing HAP binding plasma proteins (between the homozygous CFH genotyped groups) as analysed by STRING software (https://string-db.org/). A network was drawn with seven differently coloured lines that represent different types of associations. A green line indicates neighbourhood evidence; red line indicates the presence of gene fusion evidence; a blue line indicates co-occurrence evidence; a purple line indicates experimental evidence; a yellow line indicates text-mining evidence; a light blue line indicates database evidence; and a black line indicates coexpression evidence. (C) HAP-binding plasma proteins displaying a significant difference in abundance between CFH T/T1277 (TT) and C/C1277 (CC) groups. Proteins with ≥1.5 relative fold-change (both CC vs TT and TT vs CC) are shaded. FC was calculated based on the relative abundance in TT/CC. The six proteins selected for further investigation are shown in bold.


Figure 2. Detection of FHR1a and FHR1b-specific peptides in HAP-eluted fractions derived from CFH T/T1277 (TT), C/T1277 (CT) and C/C1277 (CC) plasma samples by SWATH-MS. (A) Amino acid sequence alignment of FHR1a and FHR1b proteins. Alignment was performed using ClustalW ("*" denote fully conserved residues; ":" denotes semi-conserved residues). Location of specific peptides in the sequences are shown by the grey-shaded box. (B) Plot of average intensities of the FHR1a-specific peptide in T/T1277 (TT), T/C1277 (CT) and C/C1277 (CC) groups. (C) Plot of average intensities of the FHR1b-specific peptide in T/T1277 (TT), T/C1277 (CT) and C/C1277 (CC) groups. Data in (B) and (C) are based on SWATH-MS data from HAP-eluted plasma-derived samples taken from the genotyped subjects. Statistical comparisons to the TT group were carried out by one-way ANOVA using Dunnet's multiple comparison test (ns, not significant; *, p<0.05; ***, p<0.001; ****, p<0.0001).

Figure 3. Quantitative analysis of six HAP-binding plasma proteins examined by enzyme-linked immunosorbent assay (ELISA) in HAP-eluted fractions and whole plasma. Relative concentrations of FHR1, FHR3, PZP, APOC4, C4A and C4B protein in HAP-eluted fractions (A) and whole plasma samples (B) from CFH T/T1277 (TT), C/T1277 (CT) and C/C1277 (CC) groups. Statistical comparisons to the TT group were carried out by one-way ANOVA using Dunnet's multiple comparison test (ns, not significant; *, p<0.05; **, p<0.01; ***, p<0.001; ****, p<0.001).


Figure 4. Immunolocalisation of HAP-binding proteins in the human RPE-choroid interface. Immunostaining was performed on sections from donated human eyes. RPE pigmentation shows natural auto-fluorescence. Most spherules (yellow arrows) in deposits (yellow star) are visualised by amyloid-β immunolabelling using Alexa Fluor 568 (red). Localisation of complement factor H-related protein 1 (A), Complement factor H-related protein 3 (B-C), pregnancy zone protein (D-E), apolipoprotein C4 (F), complement 4A (G-H) and complement 4B (I) are visualised by immunolabelling using Alexa Fluor 488 (green and yellow in case of co-labelling with amyloid-β; white arrows). The presented observations were found in all samples analysed, therefore the images presented are representative. Scale bar is 10 μm.




















Supplementary Data

Table S1: Participant information showing age, sex and CFH genotype (either T/T1277 (TT), T/C1277 (CT) or C/C1277 (CC)) of the 23 individuals in this study.

	1	1	1
Participant	Genotype	Sex	Age
Number			
1	СТ	F	77
2	СТ	M	81
3	CC	M	79
4	СТ	M	75
5	СТ	F	86
6	CC	M	65
7	СТ	M	85
8	TT	F	90
9	СТ	M	80
10	CC	F	64
11	TT	F	77
12	TT	F	84
13	CC	M	80
14	CC	F	84
15	CC	M	68
16	TT	F	86
17	СТ	M	84
18	TT	F	81
19	CC	M	80
20	CC	F	87
21	TT	F	82
22	CC	F	69
23	СТ	M	78

Table S2: List of HAP-binding plasma proteins detected and quantified (based on two or more peptides) in all 23 participant samples. For each protein the Uniprot ID and name, number of peptides, and average peak intensities for C/T1277 (CT), C/C1277 (CC) or T/T1277 (TT) groups are provided. The peptide number indicates the total number of peptides used for protein quantitation that were consistently detected across all samples. The quantitative value for each unique peptide originates from summing the integrated area of the selected b and y-ions for this peptide and is an average value for each genotypic group (indicated as average intensity). 'σ' represents the standard deviation

				Average	Average	Average			
Ind	Uniprot	Unipro	Pepti	peak	peak	peak	CT.	66	
ex	number	t name	des	intensity	intensity	intensity	σCT	σ СС	σTT
				CT	cc	TT			
1	P04114	APOB	10	1.20E+07	9.99E+06	1.09E+07	2.86E+06	2.27E+06	3.50E+06
2	P01024	CO3	10	2.92E+07	2.92E+07	2.84E+07	9.29E+06	7.98E+06	8.59E+06
3	P08603	CFAH	10	1.26E+07	1.15E+07	1.09E+07	8.65E+05	1.59E+06	2.81E+06
4	P01023	A2MG	10	1.21E+07	1.50E+07	1.41E+07	2.58E+06	3.63E+06	2.15E+06
5	P02768	ALBU	10	1.43E+07	1.51E+07	1.31E+07	1.12E+07	1.42E+07	1.79E+07
6	P00450	CERU	10	1.52E+07	1.41E+07	1.50E+07	2.81E+06	3.07E+06	1.91E+06
7	P02787	TRFE	10	1.41E+07	1.70E+07	1.36E+07	1.23E+06	2.94E+06	2.64E+06
8	P01031	CO5	10	3.01E+06	2.81E+06	3.06E+06	8.67E+05	1.07E+06	1.20E+06
9	P02671	FIBA	10	1.33E+07	1.31E+07	1.57E+07	3.11E+06	3.38E+06	6.05E+06
10	P00751	CFAB	10	1.09E+07	9.62E+06	9.92E+06	1.66E+06	2.76E+06	2.67E+06
11	P02675	FIBB	10	1.11E+07	1.09E+07	1.33E+07	2.26E+06	3.95E+06	5.49E+06
12	P02751	FINC	10	2.91E+06	3.08E+06	2.26E+06	2.08E+06	2.36E+06	2.28E+06
13	P00747	PLMN	10	1.96E+06	1.94E+06	2.07E+06	5.63E+05	3.95E+05	5.86E+05
14	P02679	FIBG	10	9.67E+06	9.71E+06	1.21E+07	2.15E+06	3.12E+06	4.81E+06
15	P04003	C4BPA	10	1.23E+07	1.19E+07	1.16E+07	1.88E+06	2.81E+06	1.65E+06
16	P19823	ITIH2	10	9.29E+06	9.42E+06	1.05E+07	1.28E+06	2.78E+06	3.58E+06
17	P00734	THRB	10	1.32E+07	1.03E+07	1.19E+07	4.41E+06	4.32E+06	4.35E+06
18	P02748	CO9	10	7.56E+06	7.62E+06	7.71E+06	2.02E+06	2.15E+06	1.98E+06
19	P06396	GELS	10	2.61E+06	2.32E+06	2.53E+06	4.44E+05	4.60E+05	3.28E+05
20	P02647	APOA1	10	1.06E+07	1.39E+07	1.31E+07	1.63E+06	1.48E+06	1.43E+06
21	P10643	CO7	10	1.62E+06	1.47E+06	1.58E+06	2.94E+05	3.95E+05	3.53E+05
22	P00738	HPT	10	1.24E+07	1.09E+07	1.13E+07	2.35E+06	4.73E+06	2.80E+06
23	P01871	IGHM	10	7.89E+06	6.45E+06	7.04E+06	2.89E+06	2.65E+06	4.02E+06
24	P19827	ITIH1	10	7.04E+06	7.44E+06	7.99E+06	1.05E+06	2.25E+06	2.85E+06
25	P01009	A1AT	10	7.58E+06	8.97E+06	8.13E+06	9.28E+05	2.11E+06	1.54E+06
26	Q14624	ITIH4	10	1.48E+06	2.05E+06	1.80E+06	2.36E+05	4.35E+05	5.28E+05
27	P06727	APOA4	10	1.14E+06	1.23E+06	1.41E+06	4.74E+05	3.89E+05	2.42E+05
28	P07225	PROS	10	4.14E+06	3.74E+06	3.73E+06	1.06E+06	9.29E+05	6.66E+05
29	P13671	CO6	10	1.69E+06	1.37E+06	1.60E+06	4.34E+05	2.86E+05	5.30E+05
30	P02774	VTDB	10	2.94E+06	3.65E+06	3.11E+06	4.64E+05	8.26E+05	7.87E+05
31	P01042	KNG1	10	4.92E+06	5.51E+06	5.63E+06	1.93E+06	1.63E+06	1.62E+06
32	P02790	HEMO	10	1.04E+07	1.36E+07	1.05E+07	1.30E+06	1.92E+06	2.02E+06
33	P07358	CO8B	10	2.57E+06	2.09E+06	2.48E+06	3.64E+05	2.98E+05	4.01E+05
34	P04196	HRG	10	6.26E+06	5.78E+06	6.94E+06	1.59E+06	2.81E+06	1.53E+06
35	P02749	APOH	10	2.20E+07	1.77E+07	1.86E+07	3.59E+06	4.00E+06	3.72E+06
36	P03952	KLKB1	10	2.12E+06	1.97E+06	2.17E+06	5.85E+05	5.50E+05	3.84E+05
37	P02649	APOE	10	6.95E+06	5.04E+06	6.31E+06	2.80E+06	8.61E+05	2.18E+06

P01011 AACT 10 5.89E+06 7.41E+06 6.83E+06 8.50E+05 1.47E+06 2.08E+06 40 P36955 P6DF 10 4.31E+06 3.20E+06 3.19E+06 3.46E+05 6.93E+05 4.87E+05 41 P01008 ANT3 10 3.41E+06 4.15E+06 3.38E+06 3.46E+05 6.93E+05 4.87E+05 42 P04217 A1BG 10 3.29E+06 4.31E+06 3.60E+06 3.47E+05 6.01E+05 7.21E+05 42 P04217 A1BG 10 3.29E+06 4.31E+06 3.60E+06 3.47E+05 6.01E+05 7.21E+05 42 P04217 A1BG 10 3.29E+06 4.76E+05 4.56E+05 3.34E+05 6.01E+05 7.21E+05 42 P04217 A1BG 10 5.15E+05 4.76E+05 5.15E+05 4.33E+05 1.54E+06 1.64E+06 4.49E+07 4.49E+07	38	P09871	C1S	10	1.39E+06	1.38E+06	1.52E+06	3.31E+05	3.03E+05	3.51E+05
A0			ļ							
P1	-									
P04217 A1BG 10 3.29E+06 4.31E+06 3.60E+06 3.47E+05 6.01E+05 7.21E+05 4.3 P12259 FAS 10 5.15E+05 4.75E+05 4.5E+06 7.3E+05 2.28E+05 1.22E+05 4.22E+05	-		1							-
43										
44 P00748 FA12 10										
45 P00736 C1R 10 1.59E+06 1.54E+06 1.66E+06 4.10E+05 4.22E+05 3.22E+05 46 P05155 IC1 10 5.30E+06 6.16E+06 8.65E+06 1.42E+06 2.90E+06 3.30E+06 47 P08519 APOA 10 6.71E+05 7.63E+05 7.35E+05 1.00E+06 9.08E+05 6.65E+05 48 P07357 CO8A 10 2.74E+06 2.41E+06 2.72E+06 3.74E+05 3.74E+05 3.74E+05 3.74E+05 3.74E+05 3.74E+05 3.74E+05 3.26E+06 3.32E+05 7.06E+05 5.71E+05 50 P06681 CO2 10 6.21E+05 6.34E+05 7.08E+05 1.24E+05 2.02E+05 1.94E+06 51 P02760 AMBP 10 1.24E+06 1.13E+06 1.45E+06 2.51E+05 2.25E+05 3.88E+05 53 P35858 ALS 10 1.10E+06 1.01E+06 1.08E+06 3.96E+05 2.21E+05 2.36E+05	-									
46 P05155 IC1 10 5.30E+06 6.16E+06 8.65E+06 1.42E+06 2.90E+06 3.30E+06 47 P08519 APOA 10 6.71E+05 7.63E+05 7.35E+05 1.00E+06 9.08E+05 6.65E+05 48 P07357 CO8A 10 2.74E+06 2.41E+06 2.72E+06 3.74E+05 1.43E+05 3.71E+05 49 P10909 CLUS 10 2.99E+06 3.54E+06 3.26E+06 2.39E+05 7.06E+05 5.71E+05 50 P06681 CO2 10 6.21E+05 6.34E+05 7.08E+05 1.24E+05 2.02E+05 1.94E+05 51 P0760 AMP 10 2.96E+06 3.26E+06 3.67E+06 7.42E+05 1.26E+06 1.21E+05 51 P02760 AMP 10 2.96E+06 3.26E+06 3.67E+06 7.42E+05 1.26E+06 1.21E+05 51 P036 A986 CDAI 10 1.40E+06 1.01E+06 1.08E+06 3.02E+06										
47 P08519 APOA 10 6.71E+05 7.63E+05 7.35E+05 1.00E+06 9.08E+05 6.65E+05 48 P07357 COBA 10 2.74E+06 2.41E+06 2.72E+06 3.74E+05 1.43E+05 3.71E+05 49 P10909 CLUS 10 2.99E+06 3.54E+06 3.26E+06 2.39E+05 7.06E+05 5.71E+05 50 P06681 CO2 10 6.21E+05 6.34E+05 7.08E+05 1.24E+05 2.02E+05 1.94E+05 51 P02760 AMBP 10 2.96E+06 3.26E+06 3.67E+06 7.42E+05 1.26E+06 1.21E+06 52 P05156 CFAI 10 1.24E+06 1.13E+06 1.45E+06 2.51E+05 2.25E+05 3.88E+05 53 P35858 ALS 10 1.10E+06 1.01E+06 1.08E+06 3.96E+05 2.21E+05 2.36E+05 54 P43652 AFAM 10 5.12E+05 5.08E+05 4.75E+05 2.08E+05 1.22E+06										
48 P07357 CO8A 10 2.74E+06 2.41E+06 2.72E+06 3.74E+05 1.43E+05 3.71E+05 49 P10909 CLUS 10 2.99E+06 3.54E+06 3.26E+06 2.39E+05 7.06E+05 5.71E+05 50 P06681 CO2 10 6.21E+05 6.34E+05 7.08E+05 1.24E+05 2.02E+05 1.94E+05 51 P02760 AMBP 10 2.96E+06 3.26E+06 3.67E+06 7.42E+05 1.26E+06 1.21E+06 52 P05156 CFAI 10 1.24E+06 1.13E+06 1.45E+06 2.51E+05 2.25E+05 3.88E+05 53 P35858 ALS 10 1.10E+06 1.01E+06 1.08E+06 3.96E+05 2.21E+05 2.36E+05 5.3E+05 5.478E+05 3.68E+05 1.27E+05 2.96E+05 1.63E+06 1.22E+06 2.07E+06 5.0E+04 1.3E+05 1.29E+05 6.41E+04 5.72E+05 3.68E+05 1.29E+05 1.29E+05 6.41E+04 1.25E+05 3.28E+05 <td>-</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td>	-									
49 P10909 CLUS 10 2.99E+06 3.54E+06 3.26E+06 2.39E+05 7.06E+05 5.71E+05 50 P06681 CO2 10 6.21E+05 6.34E+05 7.08E+05 1.24E+06 2.02E+05 1.94E+05 51 P0760 AMBP 10 2.96E+06 3.26E+06 3.67E+06 7.42E+05 1.26E+06 1.21E+06 52 P05156 CFAI 10 1.24E+06 1.08E+06 3.96E+05 2.25E+05 3.88E+05 53 P35858 ALS 10 1.10E+06 1.01E+06 1.08E+06 3.96E+05 2.21E+05 2.36E+05 54 P43652 AFAM 10 5.20E+05 5.94E+05 4.65E+05 1.18E+05 1.27E+05 5.91E+04 55 O43866 CDSL 10 3.41E+06 2.67E+06 3.02E+06 1.63E+06 1.22E+05 5.91E+04 56 Q92954 PRG4 10 5.12E+05 5.08E+05 4.75E+05 2.09E+05 1.22E+06 3.68E+05<										
50 P06681 CO2 10 6.21E+05 6.34E+05 7.08E+05 1.24E+05 2.02E+05 1.94E+05 51 P02760 AMBP 10 2.96E+06 3.26E+06 7.42E+05 1.26E+06 1.21E+06 52 P05156 CFAI 10 1.24E+06 1.13E+06 1.45E+06 2.51E+05 2.25E+05 3.88E+05 53 P35858 ALS 10 1.00E+06 1.00E+06 1.08E+06 3.96E+05 2.21E+05 2.36E+05 54 P43652 AFAM 10 5.20E+05 5.94E+05 4.65E+05 1.18E+05 1.27E+05 5.91E+04 55 O43866 CD5L 10 3.41E+06 2.67E+06 3.02E+06 1.63E+06 1.22E+06 2.07E+06 56 Q92954 PR64 10 5.12E+05 5.08E+05 4.20E+05 2.09E+05 1.29E+05 6.41E+04 57 P02765 FETUA 10 2.60E+06 3.52E+06 2.0E+05 8.06E+04 1.25E+05 7.49E+04										
51 P02760 AMBP 10 2.96E+06 3.26E+06 3.67E+06 7.42E+05 1.26E+06 1.21E+06 52 P05156 CFAI 10 1.24E+06 1.13E+06 1.45E+06 2.51E+05 2.25E+05 3.88E+05 53 P35858 ALS 10 1.10E+06 1.01E+06 1.08E+06 2.51E+05 2.21E+05 2.36E+05 54 P43652 AFAM 10 5.20E+05 5.94E+05 4.65E+05 1.18E+05 1.27E+05 5.91E+04 55 O43866 CD51 10 3.41E+06 2.67E+06 3.02E+06 1.63E+06 1.22E+05 5.91E+04 56 Q92954 PRG4 10 5.12E+05 5.08E+05 4.75E+05 2.09E+05 1.29E+05 6.41E+04 57 P02765 FETUA 10 2.60E+06 3.52E+06 2.64E+06 7.16E+05 9.57E+05 3.68E+05 58 P08690 PAGP 10 3.77E+05 4.78E+05 3.64E+05 7.18E+05 7.38E+	-									
52 P05156 CFAI 10 1.24E+06 1.13E+06 1.45E+06 2.51E+05 2.25E+05 3.88E+05 53 P35858 ALS 10 1.10E+06 1.01E+06 1.08E+06 3.96E+05 2.21E+05 2.36E+05 54 P43652 AFAM 10 5.20E+05 5.94E+05 1.18E+05 1.27E+05 5.91E+04 55 O43866 CD5L 10 3.41E+06 2.67E+06 3.02E+06 1.63E+06 1.22E+05 5.91E+04 56 O92954 PR64 10 5.12E+05 5.08E+05 4.75E+05 2.09E+05 1.29E+05 6.41E+04 57 P02765 FETUA 10 2.60E+06 3.52E+06 2.64E+06 7.16E+05 9.57E+05 3.68E+05 58 P08697 A2AP 10 4.12E+05 5.05E+05 4.20E+05 8.06E+04 1.25E+05 7.49E+04 59 Q96PD5 PGRP2 10 3.77E+05 4.78E+05 5.56E+04 7.14E+04 6.31E+04 <tr< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></tr<>										
53 P35858 ALS 10 1.10E+06 1.01E+06 1.08E+06 3.96E+05 2.21E+05 2.36E+05 54 P43652 AFAM 10 5.20E+05 5.94E+05 4.65E+05 1.18E+05 1.27E+05 5.91E+04 55 O43866 CDSL 10 3.41E+06 2.67E+06 3.02E+06 1.63E+06 1.22E+06 2.07E+06 56 Q92954 PRG4 10 5.12E+05 5.08E+05 4.75E+05 2.09E+05 1.29E+05 6.1E+04 57 P02765 FETUA 10 2.60E+06 3.52E+06 2.64E+06 7.16E+05 9.57E+05 3.68E+05 58 P08697 A2AP 10 4.12E+05 5.05E+05 4.20E+05 8.06E+04 1.25E+05 7.49E+04 69 Q96PD5 PGRP2 10 3.77E+05 4.78E+05 3.54E+05 5.56E+04 7.14E+04 6.31E+04 60 P29622 KAIN 10 8.10E+05 7.21E+05 7.38E+05 1.39E+05 2.13E+			ļ							
54 P43652 AFAM 10 5.20E+05 5.94E+05 4.65E+05 1.18E+05 1.27E+05 5.91E+04 55 O43866 CDSL 10 3.41E+06 2.67E+06 3.02E+06 1.63E+06 1.22E+06 2.07E+06 56 Q92954 PRG4 10 5.12E+05 5.08E+05 4.75E+05 2.09E+05 1.29E+05 6.41E+04 57 P02765 FETUA 10 2.60E+06 3.52E+06 2.64E+06 7.16E+05 9.57E+05 3.68E+05 58 P08697 A2AP 10 4.12E+05 5.05E+05 4.20E+05 8.06E+04 1.25E+05 7.49E+04 59 Q96P05 PGRP2 10 3.77E+05 4.78E+05 3.54E+05 5.56E+04 7.14E+04 6.31E+04 60 P29622 KAIN 10 8.10E+05 7.21E+05 7.33E+05 1.39E+05 2.13E+05 2.00E+05 61 P03951 FA11 10 6.92E+05 6.94E+05 6.98E+05 1.70E+05 2.88										
55 O43866 CD5L 10 3.41E+06 2.67E+06 3.02E+06 1.63E+06 1.22E+06 2.07E+06 56 Q92954 PRG4 10 5.12E+05 5.08E+05 4.75E+05 2.09E+05 1.29E+05 6.41E+04 57 P02765 FETUA 10 2.60E+06 3.52E+06 2.64E+06 7.16E+05 9.57E+05 3.68E+05 58 P08697 A2AP 10 4.12E+05 5.05E+05 4.20E+05 8.06E+04 1.25E+05 7.49E+04 59 Q96PD5 PGRP2 10 3.77E+05 4.78E+05 3.54E+05 5.56E+04 7.14E+04 6.31E+04 60 P29622 KAIN 10 8.10E+05 7.21E+05 7.83E+05 1.39E+05 2.13E+05 2.00E+05 61 P03951 FA11 10 6.92E+05 6.98E+05 1.70E+05 2.88E+05 2.01E+05 62 P05452 TETN 10 9.76E+05 1.14E+06 9.22E+05 2.32E+05 2.35E+05 2.31										
56 Q92954 PRG4 10 5.12E+05 5.08E+05 4.75E+05 2.09E+05 1.29E+05 6.41E+04 57 P02765 FETUA 10 2.60E+06 3.52E+06 2.64E+06 7.16E+05 9.57E+05 3.68E+05 58 P08697 A2AP 10 4.12E+05 5.05E+05 4.20E+05 8.06E+04 1.25E+05 7.49E+04 59 Q96PD5 PGRP2 10 3.77E+05 4.78E+05 3.54E+05 5.56E+04 7.14E+04 6.31E+04 60 P29622 KAIN 10 8.10E+05 7.21E+05 7.83E+05 1.39E+05 2.13E+05 2.00E+05 61 P03951 FA11 10 6.92E+05 6.98E+05 1.70E+05 2.18E+05 2.14E+05 62 P05452 TETN 10 9.76E+05 1.14E+06 9.22E+05 2.32E+05 2.47E+05 2.14E+05 63 P27169 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.35E+05 2.3E										
57 P02765 FETUA 10 2.60E+06 3.52E+06 2.64E+06 7.16E+05 9.57E+05 3.68E+05 58 P08697 A2AP 10 4.12E+05 5.05E+05 4.20E+05 8.06E+04 1.25E+05 7.49E+04 59 Q96PD5 PGRP2 10 3.77E+05 4.78E+05 3.54E+05 5.56E+04 7.14E+04 6.31E+04 60 P29622 KAIN 10 8.10E+05 7.21E+05 7.83E+05 1.39E+05 2.13E+05 2.00E+05 61 P03951 FA11 10 6.92E+05 6.98E+05 1.70E+05 2.88E+05 2.01E+05 62 P05452 TETN 10 9.76E+05 5.92E+05 6.98E+05 1.70E+05 2.14E+05 63 P27169 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.33E+05 2.14E+05 64 P05546 HEP2 10 1.19E+06 1.22E+06 1.13E+06 3.05E+05 2.38E+05 2.51E+05 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>										
58 P08697 A2AP 10 4.12E+05 5.05E+05 4.20E+05 8.06E+04 1.25E+05 7.49E+04 59 Q96PD5 PGRP2 10 3.77E+05 4.78E+05 3.54E+05 5.56E+04 7.14E+04 6.31E+04 60 P29622 KAIN 10 8.10E+05 7.21E+05 7.83E+05 1.39E+05 2.13E+05 2.00E+05 61 P03951 FA11 10 6.92E+05 6.94E+05 6.98E+05 1.70E+05 2.88E+05 2.01E+05 62 P05452 TETN 10 9.76E+05 1.14E+06 9.22E+05 2.32E+05 2.47E+05 2.14E+05 63 P27169 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.35E+05 2.96E+05 64 P05546 HEP2 10 1.19E+06 1.22E+06 1.33E+06 3.05E+05 2.35E+05 2.51E+05 65 P25311 ZA2G 10 1.12E+06 8.30E+05 9.5E+05 3.86E+04 1.59E+		-								
59 Q96PD5 PGRP2 10 3.77E+05 4.78E+05 3.54E+05 5.56E+04 7.14E+04 6.31E+04 60 P29622 KAIN 10 8.10E+05 7.21E+05 7.83E+05 1.39E+05 2.13E+05 2.00E+05 61 P03951 FA11 10 6.92E+05 6.94E+05 6.98E+05 1.70E+05 2.88E+05 2.01E+05 62 P05452 TETN 10 9.76E+05 1.14E+06 9.22E+05 2.32E+05 2.47E+05 2.14E+05 63 P27369 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.53E+05 2.96E+05 64 P05546 HEP2 10 1.19E+06 1.22E+06 1.13E+06 3.05E+05 2.35E+05 2.51E+05 65 P25311 ZA2G 10 3.36E+05 3.96E+05 4.04E+05 9.61E+04 1.59E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E										
60 P29622 KAIN 10 8.10E+05 7.21E+05 7.83E+05 1.39E+05 2.13E+05 2.00E+05 61 P03951 FA11 10 6.92E+05 6.94E+05 6.98E+05 1.70E+05 2.88E+05 2.01E+05 62 P05452 TETN 10 9.76E+05 1.14E+06 9.22E+05 2.32E+05 2.47E+05 2.14E+05 63 P27169 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.53E+05 2.96E+05 64 P05546 HEP2 10 1.19E+06 1.22E+06 1.13E+06 3.05E+05 2.35E+05 2.51E+05 65 P25311 ZA2G 10 3.36E+05 3.96E+05 4.04E+05 9.61E+04 1.59E+05 6.55E+04 66 P00742 FA10 10 1.12E+06 8.30E+05 9.95E+05 5.05E+05 3.86E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.03E+06 4.61E+05 1.51E+										
61 P03951 FA11 10 6.92E+05 6.94E+05 6.98E+05 1.70E+05 2.88E+05 2.01E+05 62 P05452 TETN 10 9.76E+05 1.14E+06 9.22E+05 2.32E+05 2.47E+05 2.14E+05 63 P27169 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.53E+05 2.96E+05 64 P05546 HEP2 10 1.19E+06 1.22E+06 1.13E+06 3.05E+05 2.35E+05 2.51E+05 65 P25311 ZA2G 10 3.36E+05 3.96E+05 4.04E+05 9.61E+04 1.59E+05 6.55E+04 66 P00742 FA10 10 1.12E+06 8.30E+05 9.95E+05 5.05E+05 3.86E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E+06 5.25E+05 68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+										
62 P05452 TETN 10 9.76E+05 1.14E+06 9.22E+05 2.32E+05 2.47E+05 2.14E+05 63 P27169 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.53E+05 2.96E+05 64 P05546 HEP2 10 1.19E+06 1.22E+06 1.13E+06 3.05E+05 2.35E+05 2.51E+05 65 P25311 ZA2G 10 3.36E+05 3.96E+05 4.04E+05 9.61E+04 1.59E+05 6.55E+04 66 P00742 FA10 10 1.12E+06 8.30E+05 9.95E+05 5.05E+05 3.86E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E+06 5.25E+05 68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+06 3.20E+06 69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E			+							
63 P27169 PON1 10 6.67E+05 5.92E+05 6.90E+05 1.45E+05 2.53E+05 2.96E+05 64 P05546 HEP2 10 1.19E+06 1.22E+06 1.13E+06 3.05E+05 2.35E+05 2.51E+05 65 P25311 ZA2G 10 3.36E+05 3.96E+05 4.04E+05 9.61E+04 1.59E+05 6.55E+04 66 P00742 FA10 10 1.12E+06 8.30E+05 9.95E+05 5.05E+05 3.86E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E+06 5.25E+05 68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+06 3.20E+06 69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E+05 2.13E+05 70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+	-		+							
64 P05546 HEP2 10 1.19E+06 1.22E+06 1.13E+06 3.05E+05 2.35E+05 2.51E+05 65 P25311 ZA2G 10 3.36E+05 3.96E+05 4.04E+05 9.61E+04 1.59E+05 6.55E+04 66 P00742 FA10 10 1.12E+06 8.30E+05 9.95E+05 5.05E+05 3.86E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E+06 5.25E+05 68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+06 3.20E+06 69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E+05 2.13E+05 70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+05 1.29E+05 71 P05090 APOD 10 1.13E+06 1.23E+06 1.23E+05 2.68E+05 3.69E+										
65 P25311 ZA2G 10 3.36E+05 3.96E+05 4.04E+05 9.61E+04 1.59E+05 6.55E+04 66 P00742 FA10 10 1.12E+06 8.30E+05 9.95E+05 5.05E+05 3.86E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E+06 5.25E+05 68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+06 3.20E+06 69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E+05 2.13E+05 70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+05 1.24E+05 71 P05090 APOD 10 1.13E+06 1.23E+06 1.23E+06 2.68E+05 3.69E+05 5.97E+05 72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+										
66 P00742 FA10 10 1.12E+06 8.30E+05 9.95E+05 5.05E+05 3.86E+05 5.28E+05 67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E+06 5.25E+05 68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+06 3.20E+06 69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E+05 2.13E+05 70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+05 1.24E+05 71 P05090 APOD 10 1.13E+06 1.23E+06 1.23E+06 2.68E+05 3.69E+05 5.97E+05 72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+05 3.40E+05 73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 8.75E	-									-
67 P02766 TTHY 10 3.32E+06 4.19E+06 3.05E+06 4.61E+05 1.61E+06 5.25E+05 68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+06 3.20E+06 69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E+05 2.13E+05 70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+05 1.24E+05 71 P05090 APOD 10 1.13E+06 1.23E+06 2.68E+05 3.69E+05 5.97E+05 72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+05 3.40E+05 73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 2.16E+05 1.29E+05 74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E			1							
68 P01834 IGKC 10 8.52E+06 7.40E+06 7.59E+06 3.03E+06 1.94E+06 3.20E+06 69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E+05 2.13E+05 70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+05 1.24E+05 71 P05090 APOD 10 1.13E+06 1.23E+06 2.68E+05 3.69E+05 5.97E+05 72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+05 3.40E+05 73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 2.16E+05 1.29E+05 74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E+05 75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E			+							
69 Q14520 HABP2 10 4.81E+05 4.47E+05 5.49E+05 1.64E+05 1.59E+05 2.13E+05 70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+05 1.24E+05 71 P05090 APOD 10 1.13E+06 1.23E+06 2.68E+05 3.69E+05 5.97E+05 72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+05 3.40E+05 73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 2.16E+05 1.29E+05 74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E+05 75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E+05 76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E	-									
70 P51884 LUM 10 6.00E+05 7.61E+05 6.52E+05 1.74E+05 2.68E+05 1.24E+05 71 P05090 APOD 10 1.13E+06 1.23E+06 2.68E+05 3.69E+05 5.97E+05 72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+05 3.40E+05 73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 2.16E+05 1.29E+05 74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E+05 75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E+05 76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E+05 78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.40E+05 5.08E										
71 P05090 APOD 10 1.13E+06 1.23E+06 2.68E+05 3.69E+05 5.97E+05 72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+05 3.40E+05 73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 2.16E+05 1.29E+05 74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E+05 75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E+05 76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E+05 77 P02750 A2GL 10 6.15E+05 8.40E+05 9.34E+05 1.84E+05 2.33E+05 5.05E+05 78 P20851 C4BPB 10 2.19E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 <t< td=""><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td><td></td></t<>										
72 P02753 RET4 10 9.65E+05 1.08E+06 8.54E+05 2.82E+05 2.98E+05 3.40E+05 73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 2.16E+05 1.29E+05 74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E+05 75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E+05 76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E+05 77 P02750 A2GL 10 6.15E+05 8.40E+05 9.34E+05 1.84E+05 2.33E+05 5.05E+05 78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E	-									
73 Q96IY4 CBPB2 10 7.73E+05 6.92E+05 6.45E+05 1.36E+05 2.16E+05 1.29E+05 74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E+05 75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E+05 76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E+05 77 P02750 A2GL 10 6.15E+05 8.40E+05 9.34E+05 1.84E+05 2.33E+05 5.05E+05 78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E+05 1.50E+06 80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E			1							
74 P00739 HPTR 10 1.36E+06 1.50E+06 1.34E+06 4.69E+05 8.75E+05 7.95E+05 75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E+05 76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E+05 77 P02750 A2GL 10 6.15E+05 8.40E+05 9.34E+05 1.84E+05 2.33E+05 5.05E+05 78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E+05 1.50E+06 80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E+06 2.65E+06 81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+0	-									
75 P01019 ANGT 10 4.90E+05 6.94E+05 7.47E+05 1.44E+05 9.97E+04 3.35E+05 76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E+05 77 P02750 A2GL 10 6.15E+05 8.40E+05 9.34E+05 1.84E+05 2.33E+05 5.05E+05 78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E+05 1.50E+06 80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E+06 2.65E+06 81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+05 5.64E+05			1	10	7.73E+05	6.92E+05	6.45E+05	1.36E+05		
76 P02743 SAMP 10 2.62E+06 1.80E+06 9.84E+05 1.90E+06 1.53E+06 6.13E+05 77 P02750 A2GL 10 6.15E+05 8.40E+05 9.34E+05 1.84E+05 2.33E+05 5.05E+05 78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E+05 1.50E+06 80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E+06 2.65E+06 81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+05 5.64E+05				10						
77 P02750 A2GL 10 6.15E+05 8.40E+05 9.34E+05 1.84E+05 2.33E+05 5.05E+05 78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E+05 1.50E+06 80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E+06 2.65E+06 81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+05 5.64E+05										
78 P20851 C4BPB 10 2.19E+06 2.10E+06 2.23E+06 5.28E+05 7.08E+05 5.08E+05 79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E+05 1.50E+06 80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E+06 2.65E+06 81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+05 5.64E+05	76	P02743	SAMP	10		1.80E+06	9.84E+05	1.90E+06	1.53E+06	6.13E+05
79 P01591 IGJ 10 2.84E+06 2.53E+06 2.45E+06 1.23E+06 7.40E+05 1.50E+06 80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E+06 2.65E+06 81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+05 5.64E+05	77	P02750	A2GL	10	6.15E+05	8.40E+05	9.34E+05	1.84E+05	2.33E+05	5.05E+05
80 P01876 IGHA1 9 7.88E+06 9.82E+06 8.60E+06 2.52E+06 4.66E+06 2.65E+06 81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+05 5.64E+05	78			10	2.19E+06		2.23E+06		7.08E+05	5.08E+05
81 P04004 VTNC 9 3.44E+06 3.96E+06 3.43E+06 1.56E+06 8.31E+05 5.64E+05	79	P01591	IGJ	10	2.84E+06	2.53E+06	2.45E+06	1.23E+06	7.40E+05	1.50E+06
	80	P01876	IGHA1	9	7.88E+06	9.82E+06		2.52E+06	4.66E+06	2.65E+06
93 037019 0000 0 4 605.05 4 305.05 3 705.05 0 305.04 4 305.05 4 435.05	81	P04004	VTNC	9	3.44E+06	3.96E+06	3.43E+06	1.56E+06	8.31E+05	5.64E+05
02 FZ/318 FROP 3 4.09E+U5 4.39E+U5 3.79E+U5 9.39E+U4 1.38E+U5 1.12E+U5	82	P27918	PROP	9	4.69E+05	4.39E+05	3.79E+05	9.39E+04	1.38E+05	1.12E+05
83 P00740 FA9 9 4.52E+05 4.55E+05 5.50E+05 3.10E+05 2.62E+05 2.91E+05	83	P00740	FA9	9	4.52E+05	4.55E+05	5.50E+05	3.10E+05	2.62E+05	2.91E+05
84 P07360 CO8G 9 2.38E+06 2.08E+06 2.05E+06 4.15E+05 3.26E+05 3.96E+05	84	P07360	CO8G	9	2.38E+06	2.08E+06	2.05E+06	4.15E+05	3.26E+05	3.96E+05
85 P68871 HBB 8 8.85E+05 1.04E+06 6.79E+05 3.32E+05 4.11E+05 2.41E+05	85	P68871	HBB	8	8.85E+05	1.04E+06	6.79E+05	3.32E+05	4.11E+05	2.41E+05
86 P02763 A1AG1 8 3.44E+06 4.38E+06 3.40E+06 3.97E+05 1.17E+06 1.12E+06	86	P02763	A1AG1	8	3.44E+06	4.38E+06	3.40E+06	3.97E+05	1.17E+06	1.12E+06
87 Q9BXR6 FHR5 8 2.95E+05 2.96E+05 2.91E+05 6.75E+04 1.37E+05 1.04E+05	87	Q9BXR6	FHR5	8	2.95E+05	2.96E+05	2.91E+05	6.75E+04	1.37E+05	1.04E+05

88 P02747 C1QC 8 1.23E+06 1.24E+05 1.35E+06 4.48E+05 6.70E+05 6.72E+01 89 P19652 AFAG2 8 5.84E+05 6.08E+05 5.14E+05 1.14E+05 1.85E+05 1.19E+01 90 P02652 APOA2 7 4.73E+05 5.65E+05 4.92E+05 9.35E+04 9.81E+04 6.75E+00 91 P02746 C1QB 7 9.06E+05 9.57E+05 1.03E+06 3.20E+05 5.34E+05 4.09E+05 92 P69905 HBA 7 1.02E+06 1.12E+06 3.07E+04 2.71E+04 3.35E+05 93 0.06033 ITHB 7 8.00E+05 1.37E+05 2.31E+05 4.62E+04 9.73E+04 7.21E+04 3.35E+06 95 P01857 IGHG1 6 2.95E+06 2.12E+06 1.43E+05 3.36E+06 3.18E+06 1.81E+09 96 B9A064 IGLLS 6 7.3E+05 8.32E+05 5.06E+05 1.36E+05 1.35E+05
90 P02652 APOA2 7 4.73E+05 5.65E+05 4.92E+05 9.35E+04 9.81E+04 6.75E+04 91 P02746 C1QB 7 9.06E+05 9.57E+05 1.03E+06 3.20E+05 5.34E+05 4.90E+05 92 P69905 HBA 7 1.02E+06 1.22E+06 8.11E+05 4.09E+05 4.48E+05 2.43E+01 93 Q06033 ITIH3 7 8.00E+04 1.13E+05 2.07E+04 2.71E+04 3.35E+06 94 POCOL5 CO4B 6 1.36E+05 1.37E+05 2.31E+05 4.62E+04 9.73E+04 7.21E+04 95 P01857 IGHG1 6 2.95E+06 2.12E+06 1.43E+06 3.36E+06 3.38E+06 1.81E+06 96 B9A064 IGLL5 6 7.34E+05 8.37E+05 5.81E+05 4.21E+05 5.36E+05 1.77E+01 98 O95445 APOM 6 1.50E+05 1.63E+05 1.35E+05 1.55E+05 100 <
P02746
92 P69905 HBA 7 1.02E+06 1.22E+06 8.11E+05 4.09E+05 4.48E+05 2.43E+09 93 Q06033 ITIH3 7 8.00E+04 1.13E+05 1.13E+05 2.07E+04 2.71E+04 3.35E+04 95 P01857 IGHG1 6 2.95E+06 2.12E+06 1.43E+06 3.36E+06 3.18E+06 3.18E+06 1.81E+06 95 P01857 IGHG1 6 2.95E+06 2.12E+06 1.43E+06 3.36E+06 3.18E+06 3.18E+06 96 B9A064 IGLL5 6 7.34E+05 8.37E+05 5.81E+05 4.21E+05 5.36E+05 1.77E+09 97 P00746 CFAD 6 6.39E+05 5.82E+05 5.06E+05 1.36E+05 1.35E+05 1.75E+09 97 P00746 APOM 6 1.50E+05 1.61E+05 1.46E+05 3.95E+04 4.45E+04 5.96E+04 99 P22792 CPN2 6 1.17E+05 1.45E+05 1.63E+05 3.32E+04 4.45E+04 5.96E+05 1.00 P05543 THBG 6 1.49E+05 1.82E+05 1.63E+05 3.32E+04 6.09E+04 4.89E+04 1.01 0.14791 APOL1 6 2.31E+05 2.30E+05 1.80E+05 7.55E+04 8.08E+04 7.55E+00 1.02 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 7.03E+05 3.96E+04 1.00E+05 1.03 P01859 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 7.35E+04 1.00E+05 1.03 P01859 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 1.14E+06 1.85E+01 1.05 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+05 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 1.00E+05 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.37E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+00 1.07 P49908 SEPP1 5 3.61E+05 3.39E+05 3.
93 Q06033 ITIH3 7 8.00E+04 1.13E+05 2.07E+04 2.71E+04 3.35E+04 94 POCOL5 COAB 6 1.36E+05 1.37E+05 2.31E+05 4.62E+04 9.73E+04 7.21E+00 95 PO1857 IGHGI 6 2.95E+06 2.12E+06 1.43E+06 3.36E+06 3.18E+06 1.81E+06 96 B9A064 IGLIS 6 7.34E+05 8.37E+05 5.06E+05 1.35E+05 1.35E+05 1.77E+06 97 P00746 CFAD 6 6.39E+05 5.86E+05 3.06E+05 1.35E+05 1.55E+04 4.45E+04 5.96E+04 98 O95445 APOM 6 1.50E+05 1.61E+05 1.46E+05 3.95E+04 4.45E+04 5.96E+04 99 P22792 CPN2 6 1.17E+05 1.45E+05 1.68E+05 3.32E+04 4.56E+04 5.87E+00 101 Ol4791 APOL1 6 2.31E+05 2.30E+05 1.95E+05 3.17E+04 3.96E+04<
94 POCOLS CO4B 6 1.36E+05 1.37E+05 2.31E+05 4.62E+04 9.73E+04 7.21E+04 95 PO1857 IGHG1 6 2.95E+06 2.12E+06 1.43E+06 3.36E+06 3.18E+05 1.81E+00 96 B9A064 IGLIS 6 7.34E+05 8.37E+05 5.81E+05 4.21E+05 5.36E+05 1.77E+00 97 P00746 CFAD 6 6.39E+05 5.82E+05 5.06E+05 1.35E+05 1.35E+05 1.17E+05 98 O95445 APOM 6 1.50E+05 1.61E+05 1.46E+05 3.95E+04 4.45E+04 5.96E+00 100 P05543 THBG 6 1.49E+05 1.45E+05 1.63E+05 3.32E+04 6.09E+04 4.89E+00 101 O14791 APOL1 6 2.31E+05 1.26E+05 3.32E+04 6.09E+04 4.89E+00 102 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 4.09E+04 7.51E+05<
95 P01857 IGHG1 6 2.95E+06 2.12E+06 1.43E+06 3.36E+06 3.18E+06 1.81E+06 96 B9A064 IGLL5 6 7.34E+05 8.37E+05 5.81E+05 4.21E+05 5.36E+05 1.77E+09 97 P00746 CFAD 6 6.39E+05 5.82E+05 5.06E+05 1.36E+05 1.35E+05 1.55E+04 4.95E+04 5.96E+04 98 O95445 APOM 6 1.50E+05 1.61E+05 1.46E+05 3.95E+04 4.45E+04 5.96E+04 5.96E+04 99 P22792 CPN2 6 1.17E+05 1.45E+05 1.63E+05 3.49E+04 4.56E+04 5.96E+04 1.00E+04 4.89E+04 4.56E+04 4.89E+04 1.00E+01 1.00E+03 1.45E+04 4.95E+04 4.89E+04 4.56E+04 1.00E+03 3.3EE+04 4.09E+04 4.89E+04 4.89E+04 4.89E+04 4.95E+04 4.95E+04 4.95E+04 4.95E+04 4.95E+04 4.95E+04 4.95E+04 4.95E+04 1.00E+03 3.80E+04
96 B9A064 IGLLS 6 7.34E+05 8.37E+05 5.81E+05 4.21E+05 5.36E+05 1.77E+03 97 P00746 CFAD 6 6.39E+05 5.82E+05 5.06E+05 1.36E+05 1.35E+05 1.15E+03 98 O95445 APOM 6 1.50E+05 1.61E+05 1.46E+05 3.95E+04 4.45E+04 5.96E+04 99 P22792 CPN2 6 1.17E+05 1.45E+05 1.63E+05 3.95E+04 4.45E+04 5.96E+04 100 P05543 THBG 6 1.49E+05 1.82E+05 1.68E+05 3.32E+04 6.09E+04 4.89E+04 101 O14791 APOL1 6 2.31E+05 2.30E+05 1.95E+05 7.55E+04 8.08E+04 7.55E+04 102 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 3.96E+04 1.00E+05 103 P01889 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05
97 P00746 CFAD 6 6.39E+05 5.82E+05 5.06E+05 1.36E+05 1.35E+05 1.15E+03 98 O95445 APOM 6 1.50E+05 1.61E+05 1.46E+05 3.95E+04 4.45E+04 5.96E+04 99 P22792 CPN2 6 1.17E+05 1.45E+05 1.63E+05 3.49E+04 4.56E+04 5.87E+04 100 P05543 THBG 6 1.49E+05 1.80E+05 3.32E+04 4.09E+04 4.89E+04 101 O14791 APOL1 6 2.31E+05 2.30E+05 1.95E+05 7.55E+04 8.08E+04 7.55E+04 102 P18428 IBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 3.96E+04 1.00E+05 103 P01859 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05 7.71E+05 104 P18136 KV313 5 2.80E+06 2.66E+06 2.48E+06 8.35E+05 1.44E+06 1.85E+0
98 O95445 APOM 6 1.50E+05 1.61E+05 1.46E+05 3.95E+04 4.45E+04 5.96E+04 99 P22792 CPN2 6 1.17E+05 1.45E+05 1.63E+05 4.49E+04 4.56E+04 5.87E+04 100 P05543 THBG 6 1.49E+05 1.82E+05 1.68E+05 3.3ZE+04 6.09E+04 4.89E+0-10 101 O14791 APOL1 6 2.31E+05 2.30E+05 1.95E+05 7.55E+04 8.08E+04 7.55E+01 102 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 3.96E+04 1.00E+03 103 P01880 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05 7.71E+06 105 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+01 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.14E+05 2.64
99 P22792 CPN2 6 1.17E+05 1.45E+05 1.63E+05 4.49E+04 4.56E+04 5.87E+04 100 P05543 THBG 6 1.49E+05 1.82E+05 1.68E+05 3.32E+04 6.09E+04 4.89E+06 101 O14791 APOL1 6 2.31E+05 2.30E+05 1.95E+05 7.55E+04 8.08E+04 7.55E+04 102 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 3.96E+04 1.00E+05 103 P01859 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05 7.71E+06 104 P18136 KV313 5 2.80E+06 2.66E+06 2.48E+06 8.35E+05 1.14E+06 1.85E+06 105 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+06 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.14E+05 2.64
100 P05543 THBG 6 1.49E+05 1.82E+05 1.68E+05 3.32E+04 6.09E+04 4.89E+04 101 O14791 APOL1 6 2.31E+05 2.30E+05 1.95E+05 7.55E+04 8.08E+04 7.55E+04 102 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 3.96E+04 1.00E+05 103 P01889 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05 7.71E+06 105 P01880 IGHD 5 2.06E+06 2.48E+06 8.35E+05 1.14E+06 1.85E+00 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 6.09E+05 4.39E+05 1.58E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.
101 014791 APOL1 6 2.31E+05 2.30E+05 1.95E+05 7.55E+04 8.08E+04 7.55E+04 102 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 3.96E+04 1.00E+05 103 P01859 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05 7.71E+05 104 P18136 KV313 5 2.80E+06 2.66E+06 2.48E+06 8.35E+05 1.14E+06 1.85E+06 105 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+05 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.49E+05 4.39E+05 1.58E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 108 P61626 LYSC 5 2.03E+05 1.38E+05 1.58E+05 8.07E+04 6.
102 P18428 LBP 6 1.35E+05 1.42E+05 1.80E+05 3.17E+04 3.96E+04 1.00E+05 103 P01859 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05 7.71E+05 104 P18136 KV313 5 2.80E+06 2.66E+06 2.48E+06 8.35E+05 1.14E+06 1.85E+06 105 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+05 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.14E+05 2.64E+05 1.20E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 108 P61626 LYSC 5 2.03E+05 1.38E+05 1.58E+05 6.94E+04 6.12E+04 7.37E+04 109 O75636 FCN3 5 2.82E+05 2.36E+05 3.72E+05 8.07E+04 6.1
103 P01859 IGHG2 5 7.71E+05 5.27E+05 6.64E+05 7.03E+05 8.26E+05 7.71E+05 104 P18136 KV313 5 2.80E+06 2.66E+06 2.48E+06 8.35E+05 1.14E+06 1.85E+06 105 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+05 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.14E+05 2.64E+05 1.20E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 108 P61626 LYSC 5 2.03E+05 1.38E+05 1.58E+04 6.12E+04 7.37E+04 109 O75636 FCN3 5 2.82E+05 2.36E+05 2.59E+05 8.07E+04 6.11E+04 7.89E+04 110 P17936 IBP3 5 3.99E+05 3.90E+05 3.72E+05 1.45E+05 9.77E+04 7.
104 P18136 KV313 5 2.80E+06 2.66E+06 2.48E+06 8.35E+05 1.14E+06 1.85E+06 105 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+05 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.14E+05 2.64E+05 1.20E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 108 P61626 LYSC 5 2.03E+05 1.38E+05 1.58E+05 6.94E+04 6.12E+04 7.37E+04 109 O75636 FCN3 5 2.82E+05 2.36E+05 2.59E+05 8.07E+04 6.11E+04 7.89E+04 110 P17936 IBP3 5 3.99E+05 3.90E+05 3.72E+05 1.45E+05 9.77E+04 7.90E+04 111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.
105 P01880 IGHD 5 5.06E+05 4.22E+05 1.46E+05 6.09E+05 4.39E+05 1.58E+05 106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.14E+05 2.64E+05 1.20E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 108 P61626 LYSC 5 2.03E+05 1.38E+05 6.94E+04 6.12E+04 7.37E+04 109 O75636 FCN3 5 2.82E+05 2.36E+05 2.59E+05 8.07E+04 6.11E+04 7.89E+04 110 P17936 IBP3 5 3.99E+05 3.90E+05 3.72E+05 1.45E+05 9.77E+04 7.90E+04 111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.17E+05 2.15E+05 112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.
106 P02655 APOC2 5 1.03E+06 7.67E+05 7.46E+05 7.14E+05 2.64E+05 1.20E+05 107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 108 P61626 LYSC 5 2.03E+05 1.38E+05 1.58E+05 6.94E+04 6.12E+04 7.37E+04 109 O75636 FCN3 5 2.82E+05 2.36E+05 2.59E+05 8.07E+04 6.11E+04 7.89E+04 110 P17936 IBP3 5 3.99E+05 3.90E+05 3.72E+05 1.45E+05 9.77E+04 7.90E+04 111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.17E+05 2.15E+05 112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.15E+05 113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.
107 P49908 SEPP1 5 3.61E+05 3.78E+05 4.09E+05 1.38E+05 8.65E+04 1.65E+05 108 P61626 LYSC 5 2.03E+05 1.38E+05 1.58E+05 6.94E+04 6.12E+04 7.37E+04 109 O75636 FCN3 5 2.82E+05 2.36E+05 2.59E+05 8.07E+04 6.11E+04 7.89E+04 110 P17936 IBP3 5 3.99E+05 3.90E+05 3.72E+05 1.45E+05 9.77E+04 7.90E+04 111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.17E+05 2.15E+05 112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.15E+05 113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.35E+05 1.21E+05 114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.27E+04 4.4
108 P61626 LYSC 5 2.03E+05 1.38E+05 1.58E+05 6.94E+04 6.12E+04 7.37E+04 109 O75636 FCN3 5 2.82E+05 2.36E+05 2.59E+05 8.07E+04 6.11E+04 7.89E+04 110 P17936 IBP3 5 3.99E+05 3.90E+05 1.45E+05 9.77E+04 7.90E+04 111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.17E+05 2.15E+05 112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.66E+05 113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.35E+05 1.21E+05 114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.32E+04 5.02E+04 4.49E+04 115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.9
109 O75636 FCN3 5 2.82E+05 2.36E+05 2.59E+05 8.07E+04 6.11E+04 7.89E+04 110 P17936 IBP3 5 3.99E+05 3.90E+05 3.72E+05 1.45E+05 9.77E+04 7.90E+04 111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.17E+05 2.15E+05 112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.66E+05 113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.35E+05 1.21E+05 114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.32E+04 5.02E+04 4.49E+04 115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.99E+04 116 P0CG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.1
110 P17936 IBP3 5 3.99E+05 3.90E+05 3.72E+05 1.45E+05 9.77E+04 7.90E+04 111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.17E+05 2.15E+05 112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.66E+05 113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.35E+05 1.21E+05 114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.32E+04 5.02E+04 4.49E+04 115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.99E+04 116 P0CG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.18E+04 5.21E+04 117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.7
111 P02656 APOC3 5 2.52E+06 1.73E+06 1.79E+06 1.16E+06 6.17E+05 2.15E+05 112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.66E+05 113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.35E+05 1.21E+05 114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.32E+04 5.02E+04 4.49E+04 115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.99E+04 116 P0CG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.18E+04 5.21E+04 117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.73E+04 1.10E+05 118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57
112 P01877 IGHA2 5 9.20E+05 6.94E+05 6.12E+05 7.20E+05 4.03E+05 2.66E+05 113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.35E+05 1.21E+05 114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.32E+04 5.02E+04 4.49E+04 115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.99E+04 116 P0CG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.18E+04 5.21E+04 117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.73E+04 1.10E+05 118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57E+05 2.12E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68
113 Q04756 HGFA 5 1.39E+05 1.59E+05 1.47E+05 9.84E+04 1.35E+05 1.21E+05 114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.32E+04 5.02E+04 4.49E+04 115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.99E+04 116 P0CG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.18E+04 5.21E+04 117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.73E+04 1.10E+05 118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57E+05 2.12E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E
114 P04070 PROC 5 1.17E+05 8.53E+04 1.03E+05 4.32E+04 5.02E+04 4.49E+04 115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.99E+04 116 P0CG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.18E+04 5.21E+04 117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.73E+04 1.10E+05 118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57E+05 2.12E+05 119 P04220 MUCB 4 4.45E+05 3.34E+05 4.17E+05 1.89E+05 2.10E+05 2.85E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.71E+05 1.94E+05 6.83E+04 1.06E
115 Q9NQ79 CRAC1 5 1.29E+05 1.25E+05 1.84E+05 4.27E+04 5.78E+04 4.99E+04 116 POCG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.18E+04 5.21E+04 117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.73E+04 1.10E+05 118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57E+05 2.12E+05 119 P04220 MUCB 4 4.45E+05 3.34E+05 4.17E+05 1.89E+05 2.10E+05 2.85E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E+04 3.38E+04 122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E
116 POCG05 LAC2 4 8.52E+04 7.31E+04 1.06E+05 3.37E+04 4.18E+04 5.21E+04 117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.73E+04 1.10E+05 118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57E+05 2.12E+05 119 P04220 MUCB 4 4.45E+05 3.34E+05 4.17E+05 1.89E+05 2.10E+05 2.85E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E+04 3.38E+04 122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E+05 8.53E+04 123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04<
117 P36980 FHR2 4 2.95E+05 3.77E+05 3.00E+05 5.98E+04 5.73E+04 1.10E+05 118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57E+05 2.12E+05 119 P04220 MUCB 4 4.45E+05 3.34E+05 4.17E+05 1.89E+05 2.10E+05 2.85E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E+04 3.38E+04 122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E+05 8.53E+04 123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+
118 P04275 VWF 4 6.28E+04 1.10E+05 1.43E+05 6.19E+04 1.57E+05 2.12E+05 119 P04220 MUCB 4 4.45E+05 3.34E+05 4.17E+05 1.89E+05 2.10E+05 2.85E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E+04 3.38E+04 122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E+05 8.53E+04 123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+04 1.68E+04 126 P0C0L4 C04A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+
119 P04220 MUCB 4 4.45E+05 3.34E+05 4.17E+05 1.89E+05 2.10E+05 2.85E+05 120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E+04 3.38E+04 122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E+05 8.53E+04 123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+04 1.68E+04 126 P0C0L4 C04A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+05 1.11E+05
120 P01602 KV110 4 3.05E+05 2.26E+05 4.14E+05 1.20E+05 7.68E+04 4.50E+05 121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E+04 3.38E+04 122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E+05 8.53E+04 123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+04 1.68E+04 126 P0C0L4 CO4A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+05 1.11E+05
121 Q16610 ECM1 4 1.05E+05 1.43E+05 1.17E+05 4.26E+04 3.80E+04 3.38E+04 122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E+05 8.53E+04 123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+04 1.68E+04 126 P0C0L4 CO4A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+05 1.11E+05
122 P02745 C1QA 4 1.65E+05 1.71E+05 1.94E+05 6.83E+04 1.06E+05 8.53E+04 123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+04 1.68E+04 126 P0C0L4 CO4A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+05 1.11E+05
123 P02775 CXCL7 4 2.46E+05 2.20E+05 1.12E+05 2.28E+05 1.73E+05 4.42E+04 124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+04 1.68E+04 126 P0C0L4 CO4A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+05 1.11E+05
124 Q9UK55 ZPI 4 8.22E+04 6.78E+04 5.40E+04 7.20E+04 5.33E+04 1.68E+04 126 P0C0L4 CO4A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+05 1.11E+05
126 POCOL4 CO4A 4 2.29E+05 1.84E+05 2.89E+05 9.87E+04 1.06E+05 1.11E+05
127 P01861 IGHG4 4 2.79E+05 3.98E+04 3.42E+04 6.16E+05 5.25E+04 3.58E+04
128 P20742 PZP 4 2.56E+04 6.91E+04 1.42E+05 3.45E+04 5.55E+04 1.04E+05
129 P01860 IGHG3 4 2.98E+05 3.08E+05 1.75E+05 2.85E+05 5.59E+05 2.47E+05
130 P22891 PROZ 4 5.53E+04 5.16E+04 8.28E+04 4.03E+04 2.68E+04 5.07E+04
131 P15169 CBPN 4 9.87E+04 1.28E+05 1.11E+05 2.21E+04 5.17E+04 3.88E+04
125 Q03591 FHR1 3 6.08E+05 7.50E+05 3.55E+05 1.14E+05 1.60E+05 1.21E+05
132 P01764 HV303 3 1.37E+05 1.02E+05 1.81E+05 9.61E+04 4.48E+04 2.35E+05
133 P35542 SAA4 3 7.16E+05 6.77E+05 5.47E+05 2.14E+05 1.93E+05 1.26E+05
134 P01781 HV320 3 1.51E+05 7.72E+04 6.45E+04 1.80E+05 4.65E+04 3.17E+04
135 P06310 KV206 3 2.26E+05 2.11E+05 1.57E+05 8.52E+04 1.01E+05 5.30E+04
136 P11226 MBL2 3 1.12E+05 8.37E+04 1.28E+05 7.31E+04 5.64E+04 1.04E+05
137 PODJI8 SAA1 3 9.76E+04 1.49E+05 2.14E+05 8.56E+04 1.33E+05 2.75E+05

139 O75882 ATRN 3 1.29E+04 1.66E+04 1.28E+05 2.21E+03 6.71E+03 2.99E+03 140 P02654 APOC1 3 7.28E+05 5.97E+05 5.05E+05 4.00E+05 3.29E+05 2.61E+05 141 Q151313 PCOC1 3 1.90E+05 2.04E+05 5.97E+05 9.37E+04 1.08E+05 2.99E+04 142 Q15582 BGH3 3 2.88E+04 2.98E+04 5.76E+04 1.08E+03 3.54E+03 143 Q02985 FHR3 3 8.43E+04 6.06E+04 1.07E+05 2.76E+04 1.76E+04 4.99E+04 144 P80108 PHLD 3 5.60E+04 6.67E+04 4.28E+04 1.08E+04 3.33E+04 1.90E+04 144 P80108 PHLD 3 5.60E+04 1.07E+05 1.08E+04 1.94E+04 1.37E+04 145 P08315 CBG 3 6.95E+04 1.07E+05 1.08E+04 1.57E+04 1.94E+04 3.36E+04 1.2E			1		1	1				
140 P02654 APOC1 3 7.28E+05 5.97E+05 5.05E+05 4.00E+05 3.29E+05 2.61E+05 141 Q15513 PCOC1 3 1.90E+05 2.04E+05 1.59E+05 9.37E+04 1.08E+05 2.99E+04 142 Q15582 BGH3 3 2.88E+04 2.98E+04 2.58E+04 1.07E+05 2.76E+04 1.66E+04 4.48E+04 144 P80108 PHID 3 5.60E+04 6.67E+04 4.28E+04 1.08E+04 3.33E+04 1.90E+04 144 P80108 PHID 3 5.60E+04 6.67E+04 4.28E+04 1.08E+04 3.33E+04 1.90E+04 144 P803185 CBG 3 6.95E+04 1.07E+05 1.10E+05 2.64E+04 3.94E+04 3.80E+04 147 O00391 GSOX1 3 1.14E+05 1.10E+05 1.34E+04 4.31E+03 3.42E+04 148 P23142 FBIN1 3 1.43E+04 2.02E+04 4.38E+04 4.31E+03 8.	138	P06313	KV403	3	1.19E+06	1.05E+06	1.43E+06	4.07E+05	3.99E+05	7.06E+05
141 Q15113 PCOC1 3 1.90E+05 2.04E+05 1.59E+05 9.37E+04 1.08E+05 2.99E+04 142 Q15582 BGH3 3 2.88E+04 2.98E+04 2.58E+04 5.76E+03 6.69E+03 3.54E+03 143 Q02985 FHR3 3 8.43E+04 6.06E+04 1.07E+05 2.76E+04 1.76E+04 4.99E+04 144 P80108 PHLD 3 5.60E+04 6.67E+04 4.28E+04 1.08E+04 3.33E+04 1.9E+04 1.37E+04 145 P08155 CBG 3 6.95E+04 1.07E+05 1.10E+05 2.64E+04 3.94E+04 1.37E+04 146 P08185 CBG 3 6.95E+04 1.07E+05 1.10E+05 2.64E+04 3.94E+04 3.80E+04 147 P00391 QSOX1 3 1.43E+04 2.02E+04 1.88E+04 4.31E+03 8.60E+03 1.22E+04 148 P23142 FBLM1 3 1.43E+04 2.02E+04 1.88E+04 4.31E+	139	075882	ATRN	3	1.29E+04		1.28E+04	2.21E+03	6.71E+03	2.99E+03
142 Q15582 BGH3 3 2.88E+04 2.98E+04 2.58E+04 5.76E+03 6.69E+03 3.54E+03 143 Q02985 FHR3 3 8.43E+04 6.06E+04 1.07E+05 2.76E+04 1.76E+04 4.99E+04 144 P80108 PHLD 3 5.60E+04 6.67E+04 4.28E+04 1.08E+04 3.33E+04 1.90E+04 145 P08571 CD14 3 4.64E+04 5.38E+04 1.157E+04 1.94E+04 1.37E+04 146 P08185 CBG 3 6.695E+04 1.07E+05 1.08E+05 2.64E+04 3.94E+04 1.37E+04 147 700391 QSOX1 3 1.14E+05 1.0E+05 1.08E+05 1.34E+04 8.43E+03 1.22E+04 148 P23142 FBLN1 3 1.43E+04 2.02E+04 1.89E+04 4.31E+03 8.60E+03 4.02E+03 149 P63261 ACTG 3 3.68E+04 5.10E+04 2.53E+04 4.56E+04 4.2E+03 <tr< td=""><td>140</td><td>P02654</td><td>APOC1</td><td></td><td>7.28E+05</td><td>5.97E+05</td><td>5.05E+05</td><td>4.00E+05</td><td>3.29E+05</td><td>2.61E+05</td></tr<>	140	P02654	APOC1		7.28E+05	5.97E+05	5.05E+05	4.00E+05	3.29E+05	2.61E+05
143 Q02985 FHR3 3 8.43E+04 6.06E+04 1.07E+05 2.76E+04 1.76E+04 4.99E+04 144 P80108 PHLD 3 5.60E+04 6.67E+04 4.28E+04 1.08E+04 3.33E+04 1.90E+04 145 P08571 CD14 3 4.64E+04 5.38E+04 4.82E+04 1.57E+04 1.94E+04 1.37E+04 146 P08185 CBG 3 6.95E+04 1.07E+05 1.10E+05 2.64E+04 3.94E+04 3.80E+04 147 O00391 QSOX1 3 1.43E+04 2.02E+04 1.89E+04 4.31E+03 8.60E+03 1.22E+04 148 P23142 FBLN1 3 3.68E+04 5.10E+04 3.23E+04 9.76E+03 1.30E+04 7.00E+03 150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 3.56E+04 4.66E+04 1.47E+04 151 P04278 SHBG 2 2.29E+04 3.65E+04 4.25E+04 9.76E+03 2.42E	141	Q15113	PCOC1	3	1.90E+05	2.04E+05	1.59E+05	9.37E+04	1.08E+05	2.99E+04
144 P80108 PHLD 3 5.60E+04 6.67E+04 4.28E+04 1.08E+04 3.33E+04 1.90E+04 145 P08571 CD14 3 4.64E+04 5.38E+04 4.82E+04 1.57E+04 1.94E+04 1.37E+04 146 P08185 CBG 3 6.95E+04 1.07E+05 1.10E+05 2.64E+04 3.94E+04 3.80E+04 147 O00391 QSOX1 3 1.14E+05 1.10E+05 1.08E+05 1.34E+04 8.43E+03 1.22E+04 148 P23142 FBLN1 3 1.43E+04 2.02E+04 1.89E+04 4.31E+03 8.60E+03 4.02E+04 149 P63261 ACTG 3 3.68E+04 2.02E+04 3.23E+04 9.76E+03 1.30E+04 7.00E+03 150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 3.56E+04 4.6E+04 1.4TE+04 151 P04278 SHBG 2 2.29E+04 3.65E+04 4.25E+04 9.76E+03 2.42E+	142	Q15582	BGH3	3	2.88E+04	2.98E+04	2.58E+04	5.76E+03	6.69E+03	3.54E+03
145 P08571 CD14 3 4.64E+04 5.38E+04 4.82E+04 1.57E+04 1.94E+04 3.37E+04 146 P08185 CBG 3 6.95E+04 1.07E+05 1.10E+05 2.64E+04 3.94E+04 3.80E+04 147 O00391 QSOX1 3 1.14E+05 1.10E+05 1.08E+05 1.34E+04 8.43E+03 1.22E+04 148 P23142 FBLN1 3 1.43E+04 2.02E+04 1.89E+04 4.31E+03 8.60E+03 4.02E+03 149 P63261 ACTG 3 3.68E+04 5.10E+04 2.53E+04 9.76E+03 1.30E+04 7.00E+03 150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 9.76E+03 1.30E+04 7.49E+04 151 P04278 SHBG 2 2.29E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 152 P06316 LV107 2 9.37E+04 7.39E+04 4.50E+02 5.47E+04 8.28	143	Q02985	FHR3	3	8.43E+04	6.06E+04	1.07E+05	2.76E+04	1.76E+04	4.99E+04
146 P08185 CBG 3 6.95E+04 1.07E+05 1.10E+05 2.64E+04 3.94E+04 3.80E+04 147 O00391 QSOX1 3 1.14E+05 1.10E+05 1.08E+05 1.34E+04 8.43E+03 1.22E+04 148 P23142 FBLN1 3 1.43E+04 2.02E+04 1.89E+04 4.31E+03 8.60E+03 4.02E+03 149 P63261 ACTG 3 3.68E+04 5.10E+04 3.23E+04 9.76E+03 1.30E+04 7.00E+03 150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 3.56E+04 4.66E+04 1.47E+04 151 P04278 SHBG 2 2.29E+04 3.65E+04 4.25E+04 9.76E+03 2.42E+04 2.86E+04 152 P06316 LV107 2 9.37E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 153 P02741 CRP 2 5.13E+04 7.99E+02 4.50E+02 5.47E+04 8.28E	144	P80108	PHLD	3	5.60E+04	6.67E+04	4.28E+04	1.08E+04	3.33E+04	1.90E+04
147 O00391 QSOX1 3 1.14E+05 1.10E+05 1.08E+05 1.34E+04 8.43E+03 1.22E+04 148 P23142 FBLN1 3 1.43E+04 2.02E+04 1.89E+04 4.31E+03 8.60E+03 4.02E+03 149 P63261 ACTG 3 3.68E+04 5.10E+04 3.23E+04 9.76E+03 1.30E+04 7.00E+03 150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 3.56E+04 4.66E+04 1.47E+04 151 P04278 SHBG 2 2.29E+04 7.30E+04 4.25E+04 9.76E+03 2.42E+04 2.86E+04 152 P06316 LV107 2 9.37E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 153 P02741 CRP 2 5.13E+04 6.67E+04 1.12E+05 4.09E+04 5.64E+04 1.24E+05 154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.2	145	P08571	CD14	3	4.64E+04	5.38E+04	4.82E+04	1.57E+04	1.94E+04	1.37E+04
148 P23142 FBLN1 3 1.43E+04 2.02E+04 1.89E+04 4.31E+03 8.60E+03 4.02E+03 149 P63261 ACTG 3 3.68E+04 5.10E+04 3.23E+04 9.76E+03 1.30E+04 7.00E+03 150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 3.56E+04 4.66E+04 1.47E+04 151 P04278 SHBG 2 2.29E+04 3.65E+04 4.25E+04 9.76E+03 2.42E+04 2.86E+04 152 P06316 LV107 2 9.37E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 153 P02741 CRP 2 5.13E+04 6.67E+04 1.12E+05 4.09E+04 5.64E+04 1.24E+05 154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.28E+02 2.43E+02 155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.2	146	P08185	CBG	3	6.95E+04	1.07E+05	1.10E+05	2.64E+04	3.94E+04	3.80E+04
149 P63261 ACTG 3 3.68E+04 5.10E+04 3.23E+04 9.76E+03 1.30E+04 7.00E+03 150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 3.56E+04 4.66E+04 1.47E+04 151 P04278 SHBG 2 2.29E+04 3.65E+04 4.25E+04 9.76E+03 2.42E+04 2.86E+04 152 P06316 LV107 2 9.37E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 153 P02741 CRP 2 5.13E+04 6.67E+04 1.12E+05 4.09E+04 5.64E+04 1.24E+05 154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.28E+02 2.43E+05 155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.25E+04 3.49E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.7	147	000391	QSOX1	3	1.14E+05	1.10E+05	1.08E+05	1.34E+04	8.43E+03	1.22E+04
150 P07996 TSP1 3 5.07E+04 5.61E+04 2.53E+04 3.56E+04 4.66E+04 1.47E+04 151 P04278 SHBG 2 2.29E+04 3.65E+04 4.25E+04 9.76E+03 2.42E+04 2.86E+04 152 P06316 LV107 2 9.37E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 153 P02741 CRP 2 5.13E+04 6.67E+04 1.12E+05 4.09E+04 5.64E+04 1.24E+05 154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.28E+02 2.43E+02 155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.25E+04 3.49E+04 156 P05160 F13B 2 5.30E+04 5.99E+04 6.59E+04 1.53E+04 2.01E+04 2.31E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.7	148	P23142	FBLN1	3	1.43E+04	2.02E+04	1.89E+04	4.31E+03	8.60E+03	4.02E+03
151 P04278 SHBG 2 2.29E+04 3.65E+04 4.25E+04 9.76E+03 2.42E+04 2.86E+04 152 P06316 LV107 2 9.37E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 153 P02741 CRP 2 5.13E+04 6.67E+04 1.12E+05 4.09E+04 5.64E+04 1.24E+05 154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.28E+02 2.43E+02 155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.25E+04 3.49E+04 156 P05160 F13B 2 5.30E+04 5.99E+04 6.59E+04 1.53E+04 2.01E+04 2.31E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.	149	P63261	ACTG	3	3.68E+04	5.10E+04	3.23E+04	9.76E+03	1.30E+04	7.00E+03
152 P06316 LV107 2 9.37E+04 7.30E+04 7.95E+04 5.43E+04 5.58E+04 7.73E+04 153 P02741 CRP 2 5.13E+04 6.67E+04 1.12E+05 4.09E+04 5.64E+04 1.24E+05 154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.28E+02 2.43E+02 155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.25E+04 3.49E+04 156 P05160 F13B 2 5.30E+04 5.99E+04 6.59E+04 1.53E+04 2.01E+04 2.31E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+03 158 P0CG06 LAC3 2 1.66E+04 1.28E+04 1.47E+04 8.48E+03 1.	150	P07996	TSP1	3	5.07E+04	5.61E+04	2.53E+04	3.56E+04	4.66E+04	1.47E+04
153 P02741 CRP 2 5.13E+04 6.67E+04 1.12E+05 4.09E+04 5.64E+04 1.24E+05 154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.28E+02 2.43E+02 155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.25E+04 3.49E+04 156 P05160 F13B 2 5.30E+04 5.99E+04 6.59E+04 1.53E+04 2.01E+04 2.31E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+04 158 P0CG06 LAC3 2 1.66E+04 1.28E+04 1.47E+04 8.48E+03 1.09E+04 8.57E+03 159 P01617 KV204 2 4.81E+05 4.68E+05 3.10E+05 1.79E+05 1.94E+05 1.18E+05 160 P01593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.	151	P04278	SHBG	2	2.29E+04	3.65E+04	4.25E+04	9.76E+03	2.42E+04	2.86E+04
154 P01603 KV111 2 2.02E+04 7.99E+02 4.50E+02 5.47E+04 8.28E+02 2.43E+02 155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.25E+04 3.49E+04 156 P05160 F13B 2 5.30E+04 5.99E+04 6.59E+04 1.53E+04 2.01E+04 2.31E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+03 158 P0CG06 LAC3 2 1.66E+04 1.28E+04 1.47E+04 8.48E+03 1.09E+04 8.57E+03 159 P01617 KV204 2 4.81E+05 4.68E+05 3.10E+05 1.79E+05 1.94E+05 1.18E+05 160 P01593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.37E+04 3.03E+04 161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05	152	P06316	LV107	2	9.37E+04	7.30E+04	7.95E+04	5.43E+04	5.58E+04	7.73E+04
155 P04208 LV106 2 1.73E+05 1.66E+05 1.34E+05 8.76E+04 5.25E+04 3.49E+04 156 P05160 F13B 2 5.30E+04 5.99E+04 6.59E+04 1.53E+04 2.01E+04 2.31E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+03 158 POCG06 LAC3 2 1.66E+04 1.28E+04 1.47E+04 8.48E+03 1.09E+04 8.57E+03 159 P01617 KV204 2 4.81E+05 4.68E+05 3.10E+05 1.79E+05 1.94E+05 1.18E+05 160 P01593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.37E+04 3.03E+04 161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05 3.23E+05 1.62E+05 162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03	153	P02741	CRP	2	5.13E+04	6.67E+04	1.12E+05	4.09E+04	5.64E+04	1.24E+05
156 P05160 F13B 2 5.30E+04 5.99E+04 6.59E+04 1.53E+04 2.01E+04 2.31E+04 157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+03 158 POCG06 LAC3 2 1.66E+04 1.28E+04 1.47E+04 8.48E+03 1.09E+04 8.57E+03 159 P01617 KV204 2 4.81E+05 4.68E+05 3.10E+05 1.79E+05 1.94E+05 1.18E+05 160 P01593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.37E+04 3.03E+04 161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05 3.23E+05 1.62E+05 162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03 1.16E+04 9.50E+03 163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4	154	P01603	KV111	2	2.02E+04	7.99E+02	4.50E+02	5.47E+04	8.28E+02	2.43E+02
157 Q8NF91 SYNE1 2 7.91E+03 4.85E+03 3.97E+03 8.15E+03 3.78E+03 2.13E+03 158 POCG06 LAC3 2 1.66E+04 1.28E+04 1.47E+04 8.48E+03 1.09E+04 8.57E+03 159 PO1617 KV204 2 4.81E+05 4.68E+05 3.10E+05 1.79E+05 1.94E+05 1.18E+05 160 PO1593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.37E+04 3.03E+04 161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05 3.23E+05 1.62E+05 162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03 1.16E+04 9.50E+03 163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4.41E+04 2.76E+04 164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1	155	P04208	LV106	2	1.73E+05	1.66E+05	1.34E+05	8.76E+04	5.25E+04	3.49E+04
158 POCG06 LAC3 2 1.66E+04 1.28E+04 1.47E+04 8.48E+03 1.09E+04 8.57E+03 159 PO1617 KV204 2 4.81E+05 4.68E+05 3.10E+05 1.79E+05 1.94E+05 1.18E+05 160 PO1593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.37E+04 3.03E+04 161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05 3.23E+05 1.62E+05 162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03 1.16E+04 9.50E+03 163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4.41E+04 2.76E+04 164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1.01E+04 4.86E+03 165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1	156	P05160	F13B	2	5.30E+04	5.99E+04	6.59E+04	1.53E+04	2.01E+04	2.31E+04
159 P01617 KV204 2 4.81E+05 4.68E+05 3.10E+05 1.79E+05 1.94E+05 1.18E+05 160 P01593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.37E+04 3.03E+04 161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05 3.23E+05 1.62E+05 162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03 1.16E+04 9.50E+03 163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4.41E+04 2.76E+04 164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1.01E+04 4.86E+03 165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1.21E+05 1.70E+05 166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.	157	Q8NF91	SYNE1	2	7.91E+03	4.85E+03	3.97E+03	8.15E+03	3.78E+03	2.13E+03
160 P01593 KV101 2 7.70E+04 5.68E+04 5.30E+04 1.81E+04 3.37E+04 3.03E+04 161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05 3.23E+05 1.62E+05 162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03 1.16E+04 9.50E+03 163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4.41E+04 2.76E+04 164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1.01E+04 4.86E+03 165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1.21E+05 1.70E+05 166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 2.42E+04 2.30E+04 9.33E+04 167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1	158	POCG06	LAC3	2	1.66E+04	1.28E+04	1.47E+04	8.48E+03	1.09E+04	8.57E+03
161 P80748 LV302 2 6.18E+05 6.69E+05 6.81E+05 2.26E+05 3.23E+05 1.62E+05 162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03 1.16E+04 9.50E+03 163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4.41E+04 2.76E+04 164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1.01E+04 4.86E+03 165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1.21E+05 1.70E+05 166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 2.42E+04 2.30E+04 9.33E+04 167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1	159	P01617	KV204	2	4.81E+05	4.68E+05	3.10E+05	1.79E+05	1.94E+05	1.18E+05
162 P14151 LYAM1 2 2.70E+04 2.71E+04 2.70E+04 9.98E+03 1.16E+04 9.50E+03 163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4.41E+04 2.76E+04 164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1.01E+04 4.86E+03 165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1.21E+05 1.70E+05 166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 2.42E+04 2.30E+04 9.33E+04 167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1.84E+03 2.44E+03 169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3	160	P01593	KV101	2	7.70E+04	5.68E+04	5.30E+04	1.81E+04	3.37E+04	3.03E+04
163 P10720 PF4V 2 7.14E+04 6.52E+04 4.09E+04 6.72E+04 4.41E+04 2.76E+04 164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1.01E+04 4.86E+03 165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1.21E+05 1.70E+05 166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 2.42E+04 2.30E+04 9.33E+04 167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1.84E+03 2.44E+03 169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3.50E+05 8.95E+05	161	P80748	LV302	2	6.18E+05	6.69E+05	6.81E+05	2.26E+05	3.23E+05	1.62E+05
164 P08253 MMP2 2 1.68E+04 2.05E+04 2.06E+04 4.75E+03 1.01E+04 4.86E+03 165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1.21E+05 1.70E+05 166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 2.42E+04 2.30E+04 9.33E+04 167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1.84E+03 2.44E+03 169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3.50E+05 8.95E+05	162	P14151	LYAM1	2	2.70E+04	2.71E+04	2.70E+04	9.98E+03	1.16E+04	9.50E+03
165 P06331 HV209 2 2.79E+05 2.60E+05 2.29E+05 8.18E+04 1.21E+05 1.70E+05 166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 2.42E+04 2.30E+04 9.33E+04 167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1.84E+03 2.44E+03 169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3.50E+05 8.95E+05	163	P10720	PF4V	2	7.14E+04	6.52E+04	4.09E+04	6.72E+04	4.41E+04	2.76E+04
166 P04433 KV309 2 7.68E+04 7.35E+04 8.07E+04 2.42E+04 2.30E+04 9.33E+04 167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1.84E+03 2.44E+03 169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3.50E+05 8.95E+05	164	P08253	MMP2	2	1.68E+04	2.05E+04	2.06E+04	4.75E+03	1.01E+04	4.86E+03
167 P01833 PIGR 2 2.64E+04 2.50E+04 3.19E+04 1.29E+04 1.29E+04 2.12E+04 168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1.84E+03 2.44E+03 169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3.50E+05 8.95E+05	165	P06331	HV209	2	2.79E+05	2.60E+05	2.29E+05	8.18E+04	1.21E+05	1.70E+05
168 P01612 KV120 2 3.50E+05 3.85E+03 4.87E+03 9.69E+05 1.84E+03 2.44E+03 169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3.50E+05 8.95E+05	166	P04433	KV309	2	7.68E+04	7.35E+04	8.07E+04	2.42E+04	2.30E+04	9.33E+04
169 P01779 HV318 2 8.68E+05 6.51E+05 9.21E+05 4.41E+05 3.50E+05 8.95E+05	167	P01833	PIGR	2	2.64E+04	2.50E+04	3.19E+04	1.29E+04	1.29E+04	2.12E+04
	168	P01612	KV120	2	3.50E+05	3.85E+03	4.87E+03	9.69E+05	1.84E+03	2.44E+03
170 P22105 TENX 2 1.52F+04 2.60F+04 2.43F+04 6.32F+03 1.22F+04 1.18F+04	169	P01779	HV318	2	8.68E+05	6.51E+05	9.21E+05	4.41E+05	3.50E+05	8.95E+05
170 122103 12107 2 1322101 21322101 3322103 21222101 21322101	170	P22105	TENX	2	1.52E+04	2.60E+04	2.43E+04	6.32E+03	1.22E+04	1.18E+04
171 Q12805 FBLN3 2 1.22E+04 1.52E+04 1.09E+04 1.98E+03 5.36E+03 4.68E+03	171	Q12805	FBLN3	2	1.22E+04	1.52E+04	1.09E+04	1.98E+03	5.36E+03	4.68E+03
172 P22352 GPX3 2 6.27E+04 1.02E+05 7.60E+04 1.53E+04 2.86E+04 2.71E+04	172	P22352	GPX3	2	6.27E+04	1.02E+05	7.60E+04	1.53E+04	2.86E+04	2.71E+04
173 P03950 ANGI 2 1.31E+04 1.10E+04 8.34E+03 3.53E+03 6.07E+03 1.51E+03	173	P03950	ANGI	2	1.31E+04	1.10E+04	8.34E+03	3.53E+03	6.07E+03	1.51E+03
174 Q13790 APOF 2 1.96E+04 2.63E+04 2.16E+04 8.19E+03 6.30E+03 7.20E+03	174	Q13790	APOF	2	1.96E+04	2.63E+04	2.16E+04	8.19E+03	6.30E+03	7.20E+03

Table S3: List of HAP binding plasma proteins displaying a statistically significant difference (p-value <0.05) in abundance between male and female participants.

No.	Uniprot ID	Uniprot	Group	Fold change	p-value
		name		Male/Female	
1	P20742	PZP	Pregnancy zone protein	0.25	3.88E-03
2	P06702	S10A9	Protein S100-A9	2.95	1.70E-02
3	P01877	IGHA2	Ig alpha-2 chain C region	1.94	2.10E-02
4	P01034	CYTC	Cystatin-C	1.40	2.24E-02
5	P01624	KV306	Ig kappa chain V-III region POM	1.92	2.36E-02
6	P00746	CFAD	Complement factor D	1.29	3.67E-02
7	Q96IY4	CBPB2	Carboxypeptidase B2	1.24	4.81E-02

Figure S1: Principal component analysis of SWATH-MS data from technical duplicates of HAP-eluted plasma samples from all participants. Points are labelled with the participant number (as given in Table S1). Samples were also designated as replicate "a" or "b" and assigned to their respective CFH genotypic group; T/T1277 (TT), T/C1277 (CT) or C/C1277 (CC).

