



THE UNIVERSITY *of* EDINBURGH

Edinburgh Research Explorer

Notochord Injury Assays that Stimulate Transcriptional Responses in Zebrafish Larvae

Citation for published version:

Zeng, Z, Lopez-Baez, JC, Lleras Forero, L, Brunsdon, H, Wyatt, C, Rybski, W, Hastie, N, Schulte-Merker, S & Patton, E 2018, 'Notochord Injury Assays that Stimulate Transcriptional Responses in Zebrafish Larvae' Bio-protocol, vol. 8, no. 23, pp. e3100. DOI: 10.21769/BioProtoc.3100

Digital Object Identifier (DOI):

[10.21769/BioProtoc.3100](https://doi.org/10.21769/BioProtoc.3100)

Link:

[Link to publication record in Edinburgh Research Explorer](#)

Document Version:

Peer reviewed version

Published In:

Bio-protocol

General rights

Copyright for the publications made accessible via the Edinburgh Research Explorer is retained by the author(s) and / or other copyright owners and it is a condition of accessing these publications that users recognise and abide by the legal requirements associated with these rights.

Take down policy

The University of Edinburgh has made every reasonable effort to ensure that Edinburgh Research Explorer content complies with UK legislation. If you believe that the public display of this file breaches copyright please contact openaccess@ed.ac.uk providing details, and we will remove access to the work immediately and investigate your claim.



Notochord Injury Assays that Stimulate Transcriptional Responses in Zebrafish Larvae

Zhiqiang Zeng^{1,2}, Juan C. Lopez-Baez^{1,2}, Laura Lleras-Forero^{3,4,5}, Hannah Brunson^{1,2}, Cameron Wyatt¹, Witold Rybski¹, Nicholas D Hastie¹, Stefan Schulte-Merker^{3,4,5} and E. Elizabeth Patton^{1,2}. *

¹MRC Human Genetics Unit, MRC Institute of Genetics and Molecular Medicine, University of Edinburgh, Crewe Road, Edinburgh, EH4 2XU, The United Kingdom; ²CRUK Edinburgh Centre, MRC Institute of Genetics and Molecular Medicine, University of Edinburgh, Crewe Road, Edinburgh, EH4 2XU, The United Kingdom; ³Hubrecht Institute - KNAW & UMC Utrecht, Utrecht, Netherlands; ⁴Faculty of Medicine, Institute for Cardiovascular Organogenesis and Regeneration, WWU Münster, Münster, Germany; ⁵CiM Cluster of Excellence, Münster, Germany

*For correspondence: e.patton@igmm.ed.ac.uk

[Abstract] Zebrafish have become an increasingly important model organism in the field of wound healing and regenerative medicine, due to their high regenerative capacity coupled with high-resolution imaging in living animals. In a recent study, we described multiple physical and chemical methods to induce notochord injury that led to highly specific transcriptional responses in notochord cellular subpopulations. The notochord is a critical embryonic structure that functions to shape and pattern the vertebrae and spinal column. Here, we describe precision needle injury, tail-notochord amputation, and chemical inhibition of caveolin that trigger a wound-specific *wt1b* expression response in the notochord sheath cell subpopulation. We propose that these procedures can be used to study distinct cell populations that make up the cellular processes of notochord repair.


Keywords: Zebrafish, Notochord, Tail fin, Injury, Amputation, Tissue repair, Tungsten wire, Nystatin

[Background] The notochord is a transient embryonic structure that provides axial support and signalling information to the developing embryo (Ellis *et al.*, 2013). It is comprised of two structurally distinct cell populations: the inner vacuolated cells that provide embryo support and structure, and the outer sheath cells that maintain turgor pressure for the vacuolated cells as well as patterning the developing vertebrate spine (Wopat *et al.*, 2018; Lleras Forero *et al.*, 2018; Figure 1A). We have recently discovered that *wilms' tumour 1b* (*wt1b*) is specifically expressed in a notochord sheath cell subpopulation that emerges at the site of damage and is maintained throughout repair and formation of adult vertebra structure in zebrafish (Figures 1B and 1C; Lopez-Baez *et al.*, 2018). WT1 is a zinc-finger transcription factor involved in mesodermal tissue development, adult tissue homeostasis, and becomes reactivated during epicardial tissue damage (Hastie *et al.*, 2017). Our discovery that *wt1b* becomes expressed at the notochord wound may have important implications for the development of therapies for vertebrae spinal injuries or degenerative processes.

In zebrafish wounding and regeneration models, injury is induced by a variety of methods such as amputation, surgical resection, irradiation, laser ablation and genetic ablation (Gemberling *et al.*, 2013).

40 For example, in larval zebrafish, syringe needles of various sizes have been used for tail fin amputation
41 and spinal cord injury experiments (Lisse *et al.*, 2015; Wehner *et al.*, 2017).

42 We have conducted notochord injury assays in zebrafish larvae using physical and chemical
43 approaches (Lopez-Baez *et al.*, 2018). Electrolysis-sharpened tungsten wire and insect pins described
44 here and in our recent paper induce precise, localized injury and trigger wound-specific *wt1b* expression
45 (Figures 1B and 1C). The structural integrity of the notochord can also be disrupted chemically by
46 treating embryos with nystatin, a small molecule which binds sterols and disassembles caveolae
47 (Rothberg *et al.*, 1992), which are particularly abundant in the notochord (Lim *et al.*, 2017). We detected
48 increased *wt1b* expression in nystatin treated notochords suggesting changes in caveolae caused by
49 non-physical damage and stress may also induce *wt1b* expression.



50
51

52 **Figure 1. Cell populations of the notochord and the *wt1b* notochord wound response.**

53 A. Schematic of the cell populations of the notochord. The notochord is comprised of two physically
54 distinct cell populations: an epithelial-like notochord sheath cell population (outer cells; red) and a large
55 vacuolated notochord cell population (inner cells, green), which are tightly wrapped by a thick, elastic
56 extracellular basement membrane (peri-notochordal sheath). B. Schematic of the zebrafish embryo and
57 the site of the notochord wound at the [end of the yolk sac \(YS\)](#). C. Needle injury triggers localised
58 *wt1b:gfp* expression in the notochord at the site of damage by 24 h post injury (hpi; arrow). Scale bars
59 = 100 μ m in Panel C.

60

61 **Materials and Reagents**

62


- 63 1. \varnothing 0.25 mm Tungsten wire (Alfa Aesar, catalog number: 010073.G2)
- 64 2. Metal needle holder (VWR International Ltd. UK, catalog number: MURRL110/01)
- 65 3. Sterile scalpel blade (Swann Morton UK, catalog number: 11708353)
- 66 4. 0.10 mm Austerlitz insect pins, [stainless steel \(Fine Science Tools, catalog number: 26002-10\)](#)
- 67 5. Glass Pasteur Pipettes length: 145 mm (Brand, catalog number: 7477 15)
- 68 6. Petri dish (Thermo Scientific, catalog number: 15370366)

- 69 7. Zebrafish larvae (3-7days post fertilization)
70 Non-pigmented *mitfa* mutant (*nacre* allele) (Lister *et al.*, 1999) may be preferable for ease of
71 imaging. For our experiments, we used the Tg(R2col2a1a:mCherry) transgenic line to visualise
72 notochord sheath cells (Dale and Topczewski, 2011), and the Tg(*wt1b:gfp*) line to study the
73 wound response in the notochord (Perner *et al.*, 2007).
74 *Important: Zebrafish older than 5 days post fertilization are protected animals by UK and EU*
75 *law, and require proper animal procedure licenses and approval from institutional ethics*
76 *committees.*
77 8. Dimethylsulfoxide (DMSO) (Sigma-Aldrich, catalog number: D2650-100ML)
78 9. Nystatin (Sigma-Aldrich, catalog number: N6261-500KU)
79 10. Agarose (Invitrogen, catalog number: 15510-027)
80 11. Tricaine (MS-222, 3-amino benzoic acid ethyl ester, Sigma, catalog number: A-5040)
81 12. NaOH pellets (Sigma, catalog number: S8045-500G)
82 13. Sodium chloride (NaCl) (Sigma-Aldrich, catalog number: S7653-1KG)
83 14. Potassium chloride (KCl) (Sigma-Aldrich, catalog number: P9333-500G)
84 15. Calcium chloride dihydrate (CaCl₂·2H₂O) (Sigma-Aldrich, catalog number: 223506-500G)
85 16. Magnesium chloride hexahydrate (MgCl₂·6H₂O) (Sigma-Aldrich, catalog number: M2393-500G)
86 17. E3 embryo medium (60x stock solution) (see Recipes)
87 18. Tricaine (MS-222) (1x working solution) (see Recipes)
88 19. NaOH (5 M) (see Recipes)

89
90 **Equipment**

- 91
92 1. Microscopes:
93 Upright Stereomicroscope (e.g., Nikon, model: SM2645)
94 Light microscope (e.g., Olympus, model: SZX16)
95 Confocal microscope (e.g., Nikon, model: A1R)
96 2. Electrolysis device for tungsten wire sharpening (the device shown in this protocol is custom-
97 made and no longer in production; Figure 2A)
98 However, users can set up their own equipment.
99 1) Parts required are: a). DC power supply (3-20V, such as Bosch C3 smart car battery
100 charger, Part No: 0092C35000), b). carbon electrode rod (Lasec, Catalog number:
101 ERDI9470), c). crocodile clips (optional, DC power supply may come with leads with
102 crocodile clips), d) a 200ml glass container (jar) with lid.
103 2) Assemble the equipment: Plug the leads into the charger. Using crocodile clips,
104 connect the carbon electrode rod to the negative terminal (black lead), connect the
105 metal holder to the positive terminal (red lead). Set the voltage to 6V and fill up the jar
106 with 100ml of NaOH solution. The equipment is now ready for use.
107 3. Bunsen burner

108



109

110

111

112

113

114

115

116


117

118


Figure 2. Preparation of electrolysis-sharpened needle. A. Electrolysis-based apparatus, consisting of a rectified DC transformer, an anode with crocodile clip (red), a cathode with carbon electrode (yellow), electrolysis chamber with 5M NaOH electrolyte. B-C. Lateral and dorsal views of a tungsten wire needle being sharpened. With the power on (6 V), the mounted needle is held vertically and dipped into and out of the electrolyte steadily and slowly until the desired tip is achieved. D. A finished electrolysis-sharpened tungsten wire needle.

Procedure

- 119 A. Notochord needle injury of zebrafish larvae
- 120 1. Prepare electrolysis-sharpened tungsten wire (adapted from Brady, 1965) (Figure 2).
- 121 a. Cut off 3 cm of tungsten wire and mount it into a needle holder
- 122 b. Connect the metal handle of the needle holder to the (+) terminal of the transformer
- 123 using a crocodile clip, connect the cathode with carbon plate to the (-) terminal.
- 124 c. Place the carbon electrode in the glass chamber with 100 ml of 5 M NaOH solution and
- 125 switch on the transformer, set the output voltage to 6 V. With the mounted needle held
- 126 vertically, dip the needle into and out of NaOH, slowly and steadily until desired tip is
- 127 produced. Faster movement = Longer slope on needle, Slower movement = Shorter
- 128 tip with more angled slope. The dial gauge of current reads between 0 and 1 A when
- 129 the needle moves up and down. It takes about 2.5 min to sharpen a needle from 0.25
- 130 mm to 0.02 mm in diameter.




- 131
- 132 2. Alternatively, 0.1 mm insect pins can also be used to injure the notochord
- 133 a. Take a clean glass pipette and using a Bunsen burner, bend the thin side in the middle in
- 134 order to create a 45 degree angle. This will help with the injury manipulation procedure.
- 135 b. Close the hole of the thin side by about three-quarters using the Bunsen burner. This is
- 136 done by placing the tip in the strongest part of the flame and rotating it in a circular manner.
- 137 c. Take the insect pin with forceps and place it in the hole, taking care that the sharp side is
- 138 facing the outside.
- 139 d. Carefully continue burning the tip of the glass pipette in the weakest part of the flame, until
- 140 the hole is closed and the insect pin secure. The insect pin will burn if exposed to too much
- 141 heat. Try to close the hole as much as possible before inserting the pin.



142
143
144

Figure 3. Preparation of insect pin. A. The required equipment. B. The finished instrument. C. A close up image of the tip.



145
146
147
148
149
150
151
152
153
154
155
156
157
158

3. Prepare 1.5% agarose: Weight 1.5g agarose and melt in 100ml E3 embryo medium in a 250ml glass bottle using a microwave oven, pour a thin layer into a ϕ 90mm Petri dish. About 20ml of 1.5% agarose is used. Let the agarose solidify.
4. Anaesthetize larvae in tricaine solution: Prepare 1:10,000 tricaine solution using E3 embryo medium (please see Recipes), and pour 30ml in a separate ϕ 90mm Petri dish. Transfer one larva into the tricaine solution and wait until it is anesthetized.
5. Under a stereomicroscope, place one larva on it's side on a Petri dish coated with agarose so that the lateral side can be accessed with needle from above. Remove as much liquid as possible so that the surface tension adheres the larvae to the dish and prevents it from slipping. Gently insert the tip of the tungsten wire into the notochord vertically at the level of the end of the yolk sac (Figure 1B), then withdraw the wire.

- 159 6. Transfer injured larvae to a petri dish with fresh E3 medium to recover and place the dish at
160 28.5 °C to grow the larvae to the desired stages. Keep uninjured age-matched larvae as non-
161 injured controls.

Video3: Notochord Needle Injury


162

163 B. Chemically-induced disruption of notochord

- 164 1. Cross fish carrying the *Tg(wt1b:GFP)* and the notochord-marking *Tg(R2-col2a1a:mCherry)*
165 transgenes in an unpigmented *nacre*^{-/-} background, to obtain *Tg(wt1b:GFP;R2-*
166 *col2a1a:mCherry);nacre*^{-/-} embryos.
- 167 2. Prepare fresh 5 mg/ml nystatin stock solution (5.4mM) before each use by dissolving in DMSO.
- 168 3. Dilute nystatin stock solution in E3 to obtain 20 µM final working concentration. Add this to
169 dechorionated 48 hpf embryos in a 6-well plate. Add 0.4% DMSO to control embryos.
- 170 4. Incubate embryos at 28.5 °C for up to 48 h. After 24 hours of nystatin treatment, lesions appear
171 along the length of the notochord. They tend to appear first in regions that are naturally
172 compressed as the embryo moves, and then spread along the length of the notochord. The
173 majority of embryos acquire notochord lesions, however their size and severity can be variable.
174 Therefore, regular screening for lesions and/or the onset of *wt1b:GFP* expression is
175 recommended in order to identify embryos with the desired level of notochord damage.
- 176 5. For imaging, anaesthetise embryos in tricaine (1:10,000), and mount sagittally in 1% low-melt
177 agarose. Brightfield images are taken using a light microscope (Figures 4A and 4B). Expression
178 of the *R2-col2a1a:mCherry* transgene, which marks the notochord, and the induction of the
179 *wt1b:GFP* transgene at sites of notochord damage is visualised using confocal microscopy
180 (Figures 4A' and 4B").

181

182




183
184
185
186
187
188
189
190
191
192
193
194
195

Figure 4. Disrupting notochord structure using nystatin (modified from Lopez-Baez, 2018). Nystatin is a small molecule which binds sterols and leads to the disassembly of caveolae, a component abundant in the notochord (Lim *et al.*, 2017). *Tg(wt1b:GFP;R2-col2a1a:mCherry);nacre^{-/-}* zebrafish embryos are treated with either DMSO or 20µM Nystatin from 48 hpf to 72 hpf. When observed under a light microscope, the notochord structure of (A) DMSO-treated embryos appears normal, however lesions can be observed in (B) nystatin-treated embryos. (A' and B') *wt1b:GFP* expression is induced at lesion sites, but not in control notochords. *R2-col2a1a:mCherry* expression in notochord sheath cells also shows increased cellularity at (B'') lesion sites of nystatin-treated embryos compared to (A'') DMSO controls. Scale bars are 50 µm.

196 C. Tail amputation

- 197 1. Prepare 1.5% agarose using E3 embryo medium and pour a thin layer into a Petri dish. Let the
198 agarose solidify.
- 199 2. Anaesthetize larvae in tricaine solution.
- 200 3. Under a stereomicroscope, place one larva on it's side onto the solidified agarose. Remove as
201 much as liquid as possible so the surface tension adheres the larvae to the dish and prevents it
202 from slipping, then amputate the tail with a sterile scalpel blade with slight pressure. Amputation
203 sites are dependent on experiments being performed (Figure 5). Amputations at the **tail fin and**
204 **tip of notochord** site do not include notochord tissue, and do not stimulate a *wt1b:gfp* notochord
205 injury response. Amputations beyond the **tail fin** and into the notochord (**before caudal vein, past**
206 **caudal vein**) stimulate a *wt1b:gfp* expression.
- 207 4. Transfer injured larvae to a Petri dish with fresh E3 medium to recover, and place the dish at
208 28.5 °C to grow the larvae to the desired stages. Keep uninjured age-matched larvae as non-
209 injured controls.
- 210



211 **Figure 5. Selected tail amputations uncover the notochord specificity of the response**
212 (modified from Lopez-Baez, 2018). (A) Illustration of tail amputations at different tail sites and
213 time points when images are taken. Tail fin (TF), tip of the notochord (TN), before caudal vein
214 (BCV), past caudal vein (PCV), somite (S), notochord (N), caudal vein (CV). (B) TF and TN
215 amputated larvae showed no GFP upregulation in their notochord after the injury, but show
216 marked fin regeneration (arrow head). BCV and PCV amputated groups both showed strong
217 GFP upregulations by 72hpa (arrows), with PCV amputated larvae showing an overall stronger
218 and faster upregulation than BCV amputated larvae.

221

222 **Notes**

223

224 1. In the UK and EU, all animal procedures need to be approved by the Home Office (UK) or its
225 equivalent. Appropriate Personal Project License (PPL) and Personal individual License (PIL)
226 are required.

227 2. The procedure of notochord needle injury requires a fair amount of practice, and care should
228 be taken to not cause injury outside the needle injury site. The appearance of a small bulge
229 structure at the site of injury within the notochord about 5 min post-surgery indicates a
230 successful operation. It is achievable to injure 30 larvae during a period of an hour.

231 3. Optimising nystatin dosage for the first-time use is recommended, as there is batch-to-batch
232 variation. A longer nystatin incubation period can be attempted however the adverse off-target
233 effects of nystatin cause gross developmental abnormalities and embryos do not survive long-
234 term. *In our hands, treatment with 20µM nystatin from 48hpf for 24 hours gives the most
235 consistent results. Treatment with nystatin before 48hpf is possible, however due to its off-target
236 effects, more toxicity is seen. Embryos tolerate later nystatin treatment from 72 and 96 hpf much
237 better, with 20µM nystatin producing notochord lesions in 60-80% of embryos after 24 hours,
238 although these lesions are smaller in size and fewer per embryo compared those shown in
239 Figure 4.*

240

241 **Recipes**

242

243 1. E3 embryo medium (60x stock solution)

244 17.4 g NaCl

245 0.8 g KCl

246 2.9 g CaCl₂·2H₂O

247 4.89 g MgCl₂·6H₂O

248 Dissolve in 1 L H₂O

249 2. Tricaine (MS-222) (1x working solution)

250 Dissolve 0.1g of Tricaine powder in 1 L of 1x E3 medium, adjust pH to 7.0

251 3. NaOH (5 M)

252 Dissolve 200 g NaOH pellets in 1 L H₂O

253

254 **Acknowledgments**

255

256 Funding sources: MRC Human Genetics Unit Programme (MC_PC_U127585840,
257 MC_PC_U127527180), European Research Council (ZF-MEL-CHEMBIO-648489) and L'Oreal-
258 Melanoma Research Alliance (401181), Cells in Motion - Cluster of Excellence (EXC 1003-CiM).

259 These procedures are adapted from our recent eLife paper (Lopez-Baez *et al.*, 2018)

260
261
262
263
264
265
266
267
268
269
270
271
272
273
274
275
276
277
278
279
280
281
282
283
284
285
286
287
288
289
290
291
292
293
294
295
296
297
298

Competing interests

The authors have no competing interests.

Ethics

Procedures presented here were approved by the University of Edinburgh Ethics Committee, and performed under the Home Office Project License 70/8000 to EEP at the University of Edinburgh, United Kingdom; and by the Animal Experimentation Committee (DEC) of the Royal Netherlands Academy of Arts and Sciences to SSM at the Hubrecht Institute and the Institute of Cardiovascular Organogenesis and Regeneration WWU Münster, Germany.

References

1. Brady, J. (1965). [A simple technique for making very fine, durable dissecting needles by sharpening tungsten wire electrolytically.](#) *Bull World Health Organ* 32(1): 143-144.
2. Dale, R. M. and Topczewski, J. (2011). [Identification of an evolutionarily conserved regulatory element of the zebrafish *col2a1a* gene.](#) *Dev Biol* 357(2): 518-531.
3. Ellis, K., Hoffman, B. D. and Bagnat, M. (2013). [The vacuole within: how cellular organization dictates notochord function.](#) *Bioarchitecture* 3(3): 64-68.
4. Gemberling, M., Bailey, T. J., Hyde, D. R. and Poss, K. D. (2013). [The zebrafish as a model for complex tissue regeneration.](#) *Trends Genet* 29(11): 611-620.
5. Hastie, N. D. (2017). [Wilms' tumour 1 \(WT1\) in development, homeostasis and disease.](#) *Development* 144(16): 2862-2872.
6. Lim, Y. W., Lo, H. P., Ferguson, C., Martel, N., Giacomotto, J., Gomez, G. A., Yap, A. S., Hall, T. E. and Parton, R. G. (2017). [Caveolae protect notochord cells against catastrophic mechanical failure during development.](#) *Curr Biol* 27(13): 1968-1981 e1967.
7. Lister, J. A., Robertson, C. P., Lepage, T., Johnson, S. L. and Raible, D. W. (1999). [Nacre encodes a zebrafish microphthalmia-related protein that regulates neural-crest-derived pigment cell fate.](#) *Development* 126(17): 3757-3767.
8. Lisse, T. S., Brochu, E. A. and Rieger, S. (2015). [Capturing tissue repair in zebrafish larvae with time-lapse brightfield stereomicroscopy.](#) *J Vis Exp*(95).
9. Lleras Forero, L., Narayanan, R., Huitema, L. F., VanBergen, M., Apschner, A., Peterson-Maduro, J., Logister, I., Valentin, G., Morelli, L. G., Oates, A. C. and Schulte-Merker, S. (2018). [Segmentation of the zebrafish axial skeleton relies on notochord sheath cells and not on the segmentation clock.](#) *Elife* 7: e33843.
10. Lopez-Baez, J. C., Simpson, D. J., L, L. L. F., Zeng, Z., Brunsdon, H., Salzano, A., Brombin, A., Wyatt, C., Rybski, W., Huitema, L. F. A., Dale, R. M., Kawakami, K., Englert, C., Chandra, T.,

- 299 Schulte-Merker, S., Hastie, N. D. and Patton, E. E. (2018). [Wilms Tumor 1b defines a wound-](#)
300 [specific sheath cell subpopulation associated with notochord repair.](#) *Elife* 7: e30657.
- 301 11. Perner, B., Englert, C. and Bollig, F. (2007). [The Wilms tumor genes wt1a and wt1b control](#)
302 [different steps during formation of the zebrafish pronephros.](#) *Dev Biol* 309(1): 87-96.
- 303 12. Rothberg, K. G., Heuser, J. E., Donzell, W. C., Ying, Y. S., Glenney, J. R. and Anderson, R. G.
304 (1992). [Caveolin, a protein component of caveolae membrane coats.](#) *Cell* 68(4): 673-682.
- 305 13. Wehner, D., Tsarouchas, T. M., Michael, A., Haase, C., Weidinger, G., Reimer, M. M., Becker,
306 T. and Becker, C. G. (2017). [Wnt signaling controls pro-regenerative Collagen XII in functional](#)
307 [spinal cord regeneration in zebrafish.](#) *Nat Commun* 8(1): 126.
- 308 14. Wopat, S., Bagwell, J., Sumigray, K. D., Dickson, A. L., Huitema, L. F. A., Poss, K. D., Schulte-
309 Merker, S. and Bagnat, M. (2018). [Spine patterning is guided by segmentation of the notochord](#)
310 [sheath.](#) *Cell Rep* 22(8): 2026-2038.
- 311
- 312