

How to Assign a (3+1)D Superspace Group to an Incommensurately Modulated Biological Macromolecular Crystal

Gloria Borgstahl

University of Nebraska Medical Center, gborgstahl@unmc.edu

Jason Porta

University of Nebraska Medical Center

Jeff Lovelace

University of Nebraska Medical Center

Follow this and additional works at: <https://lib.dr.iastate.edu/aperiodic2018>



Part of the [Chemistry Commons](#), and the [Materials Science and Engineering Commons](#)

Borgstahl, Gloria; Porta, Jason; and Lovelace, Jeff, "How to Assign a (3+1)D Superspace Group to an Incommensurately Modulated Biological Macromolecular Crystal" (2018). *Aperiodic 2018 ("9th Conference on Aperiodic Crystals")*. 41.

<https://lib.dr.iastate.edu/aperiodic2018/2018/abstracts/41>

This Oral Presentation is brought to you for free and open access by the Conferences and Symposia at Iowa State University Digital Repository. It has been accepted for inclusion in Aperiodic 2018 ("9th Conference on Aperiodic Crystals") by an authorized administrator of Iowa State University Digital Repository. For more information, please contact digirep@iastate.edu.

How to Assign a (3+1)D Superspace Group to an Incommensurately Modulated Biological Macromolecular Crystal

Gloria E. O. Borgstahl^{a,b*}, Jason Porta^a, and Jeff Lovelace^b

^aThe Eppley Institute for Research in Cancer and Allied Diseases and Department of Biochemistry and Molecular Biology, University of Nebraska Medical Center, 987696 Nebraska Medical Center, Omaha, NE, 68198, USA

^bThe Eppley Institute for Research in Cancer and Allied Diseases, University of Nebraska Medical Center, 987696 Nebraska Medical Center, Omaha, NE, 68198, USA

*gborgstahl@unmc.edu

Periodic crystal diffraction is described using a three dimensional (3D) unit cell and 3D space group symmetry. Incommensurately modulated crystals are a subset of aperiodic crystals that need four to six dimensions to describe the observed diffraction pattern and they have characteristic satellite reflections that are offset from the main reflections. These satellites have a non-integral relationship to the primary lattice and require a \mathbf{q} vector(s) for processing. Incommensurately modulated biological macromolecular crystals have been frequently observed but so far not solved. We have been spearheading an initiative to solve this type of crystal. The first step toward structure solution is to collect the diffraction data making sure that the satellite reflections are well separated from the main reflections. Once collected they can be integrated with *Eval15* and then scaled with *SADABS* software. Then the assignment of the superspace group is needed. The most common form of modulation is in only one extra direction and can be described with a (3+1)D superspace group. The (3+1)D superspace groups for chemical crystallographers are fully described in Volume C of the International Tables for Crystallography in Chapter 9.8. This text includes all types of crystallographic symmetry elements found in small molecule crystals and can be difficult for structural biologists to understand and apply to their crystals. Here we provide an explanation for structural biologists that includes only the subset of biological symmetry elements and demonstrate the application to a real-life example of an incommensurately modulated protein crystal.