

THE
UNIVERSITY
OF RHODE ISLAND

University of Rhode Island
DigitalCommons@URI

Department of Electrical, Computer, and
Biomedical Engineering Faculty Publications

Department of Electrical, Computer, and
Biomedical Engineering

2016

Smartwatch-driven multisensory recorder: Design and testing of a smartwatch-based framework to support psychiatric disorders

Matthew Constant

Leslie Brick

See next page for additional authors

Follow this and additional works at: https://digitalcommons.uri.edu/ele_facpubs

**The University of Rhode Island Faculty have made this article openly available.
Please let us know how Open Access to this research benefits you.**

Terms of Use

This article is made available under the terms and conditions applicable towards Open Access Policy Articles, as set forth in our [Terms of Use](#).

Authors

Matthew Constant, Leslie Brick, Nicole Nugent, Michael Armey, and Kunal Mankodiya

Smartwatch-driven Multisensory Recorder

Design and testing of a smartwatch-based framework to support psychiatric disorders

Matthew Constant¹, Leslie Brick², Nicole Nugent², Michael Armey³, Kunal Mankodiya¹

¹Department of Electrical, Computer, and Biomedical Engineering, University of Rhode Island, RI, USA

²Department of Psychiatry, Rhode Island Hospital, Providence, RI, USA.

³Butler Hospital & Brown University School of Medicine, Providence, RI, USA

Abstract—The advancement of technologies such as Bluetooth Low Energy (BLE) and the availability of low cost sensors has allowed the Internet-of-Things (IoT) industry to develop rapidly. The Android application described in this paper takes advantage of the opportunities provided by these technologies. This application can simultaneously collect data from all sensors available on the Microsoft Band 2. It then stores this data to be later viewed and analyzed. This method of data collection could help in many cases, including psychiatric disorders. For example, our app is aimed at monitoring physiological data (heart rate and galvanic skin response) and ambient environment data (sound and lighting conditions) of trauma exposed individuals with post-traumatic stress disorder (PTSD). The preliminary data collection study on healthy individuals demonstrate that our framework efficiently handles the continuous recording of multiple sensors.

Keywords—Smart Textile; Motion Sensing, Parkinson's Disease; Wearable Sensors; Internet-of-Things.

I. INTRODUCTION

Today, an overwhelming number of citizens face mental health challenges. For example, post traumatic stress disorder (PTSD) that is a kind of psychiatric disorder manifesting various emotional and behavioral symptoms such as anxiety, mood swings, depression, anger, impulsivity, and many other issues is prevalent in adult Americans at the rate of 6.8% [1]. These individuals with PTSD receive various behavioral therapies. It is very difficult to personalize the therapy based on the manifested symptoms. Therefore, it is very important for clinicians to monitor their patients with PTSD when they are at home. Our paper presents a smartwatch-based data collection framework, called “Android Electronically Activated Recorder (anEAR).

The anEAR framework uses the multi-dimensional sensors built into the Microsoft Band 2, alongside a Bluetooth capable Android device. The framework includes a custom application that uses the smartwatch from which the Android device receives sensor data. anEAR is able to simultaneously collect sensor data from several sensors while also periodically recording audio using the device’s microphone. This allows users to view various individual’s biometrics while also being able to pinpoint the cause of irregularities by matching the audio recordings to the sensor data.




Figure 1: anEAR Smartwatch Framework

II. BACKGROUND

A. Android App Development

Android is an open-source operating system based on the Linux kernel and sponsored by Google. It runs on many different systems including smartphones, tablets, smartwatches, and more. Android’s Software Development Kit (SDK) provides many different Application Programming Interfaces (APIs) that allow developers to take full advantage of all of the technologies available on the Android device. Due to its ability to run on so many diverse devices along with all of the APIs available, Android is the ideal operating system upon which to build the anEAR framework.

B. Smartwatch - Microsoft Band 2

The Microsoft Band 2 is an activity tracker developed by Microsoft which features continuous streaming of heart rate data along with ten other sensors. The eleven sensors are: optical heart rate sensor, three axis accelerometer, gyrometer, Global Positioning System (GPS), ambient light sensor, skin temperature, ultra-violet sensor, capacitive sensor, Galvanic Skin Response (GRS), microphone, and barometer. A description of each sensor can be seen in Table 1. The Band can connect to any iPhone or smartphone running Android or Windows operating system. The band connects using its own application which must be installed on the device at all times. Furthermore, the band also provides its own SDK for iOS, Android, and Windows that allows developers to access the band’s sensors at any time and to display notifications on the

band's screen. package provided by the band was the driving factor in our decision to use it in the anEAR framework.

Table 1: List of built-in sensors of Microsoft Band 2.0 [2]

Sensor Name	Description	Sampling Rate
Accelerometer	Provides X, Y, and Z acceleration in g units. 1 g = 9.81 meters per second squared (m/s ²).	62/31/8 Hz
Gyroscope	Provides X, Y, and Z angular velocity in degrees per second (°/sec) units.	62/31/8 Hz
Distance	Provides the total distance in centimeters, current speed in centimeters per second (cm/s), current pace in milliseconds per meter (ms/m), and the current pedometer mode (such as walking or running)	1 Hz
Heart Rate	Provides the number of beats per minute; also indicates if the heart rate sensor is fully locked on to the wearer's heart rate.	1 Hz
Pedometer	Provides the total number of steps the wearer has taken since the Band was last factory-reset.	Value Change
Skin Temperature	Provides the current skin temperature of the wearer in degrees Celsius.	1 Hz
UV	Provides the current ultraviolet radiation exposure intensity.	1 Hz
Band Contact	Provides the current state of the Band as being worn/not worn.	Value Change
Calories	Provides the total number of calories the wearer has burned since the Band was last factory-reset.	Value Change
GSR	Provides the current skin resistance of the wearer in kohms.	0.2/5 Hz
RR Interval	Provides the interval in seconds between the last two continuous heart beats.	Value Change
Ambient Light	Provides the current light intensity (illuminance) in lux (Lumes per sq. meter).	2 Hz
Barometer	Provides the current raw air pressure in hPa (hectopascals) and raw temperature in degrees Celsius.	1 Hz
Altimeter	Provides current elevation data like total gain/loss, steps ascended/descended, flights ascended/descended, and elevation rate.	1 Hz

III. APP DEVELOPMENT & PROGRAMMING

The Android application for anEAR was designed to have a straight-forward User Interface (UI) while handling as many errors as possible without stopping the entire application. In order to do this, the application follows several fundamental design guidelines created specifically for anEAR.

A. anEAR Activities

In Android, *Activities* are a main component of any application since they define the UI and handle any events corresponding to the UI such as button clicks. In anEAR, every Activity provides a way to graphically show messages to the User and contains static methods that handle all responsibilities related to the Activity. For example, the Settings Activity provides static methods to access and manipulate the settings that anEAR uses. This allows for each Activity to be completely independent of each other while providing a layer of abstraction between the client, in this case the Activity accessing the settings, and the provider, in this case the Settings Activity. The reason for this is to make further development on anEAR as efficient as possible and this allows new Activities to be easily implemented and existing Activities to be easily updated. For example, if the settings on anEAR needed to be moved from residing in a file to being stored in a database, only the Settings Activity would have to be modified and the rest of the Activities would not even be aware of the change. anEAR consists of four Activities: Settings Activity, Lock Screen Activity, Blackout Settings Activity, and Main Activity.




Figure 2: A workflow diagram describing the process to begin sensor recordings.

Main Activity is the first activity that is launched when the User starts the application. First, this Activity checks that the device has all the correct permissions and contains all the necessary hardware. If not, it displays the appropriate error message to the User. Next, it determines whether or not the application has a pattern lock set and starts the Lock Screen Activity, passing this information to it. anEAR is meant to always be locked with a pattern lock, so if Main Activity does not find a pattern lock saved, it prompts the User to create one by starting the Lock Screen Activity to create a pattern lock.




Figure 3: App screenshots of AnEAR. Pattern lock restricts the user to change the settings. Settings are for clinicians to provide therapies who can select sensors and their sampling periods. Sensor data are logged onto the phone memory.

Lock Screen Activity is responsible for creating and changing the pattern lock for anEAR as well as making the User enter the correct pattern lock to continue in the application. anEAR uses a third party library [3], to display the pattern lock as well as providing information about it, such as what the pattern is.

Settings Activity is responsible for displaying all of the current settings for anEAR as well as providing the User with a way to manipulate each setting. The User can edit the duration of an audio or sensor recording, the delay between audio or sensor recordings, the sensors to be recorded, the sensors to trigger automatic recordings, and the name of the root folder where all of the data is saved.

Blackout Settings Activity is responsible for creating, updating, and deleting blackout periods. Blackout periods are periods of time where anEAR will not record and periodic audio recordings. This does not affect sensor recordings, which are toggled in a way that will be explained later. anEAR allows the User to specify a start time, an end time, and whether this should take effect on the weekends or weekdays. The values are then encoded into binary strings and saved into a SQLite database on the device.

B. anEAR Services

Services are another main component in Android. *Services* are very similar to *Activities* except they do not have any UI and can be run on a separate thread than the application's main, or UI, thread. In anEAR, *Services* are used to connect and collect sensor data from the Microsoft Band as well as saving this data on the device. anEAR only connects to the Band, and therefore only collects data, via *Collection Services*. These *Services* are made specifically for anEAR collecting data from the Band. They automatically connect to the Band when they are created, start a timer once connected, and disconnect from the Band once the timer is finished. The duration of the timer is based on the amount of time the User specifies for recording sensor data. By using this method to connect to the Band, cases of unnecessary or unclosed connections are significantly

decreased which significantly increases battery life of both the device and the Band. Each sensor that is to be streamed from the device creates a new *Collection Service*. By using individual *Services* for each sensor, anEAR is able to gracefully handle any error related to one sensor while streaming all other sensors without issue.

anEAR uses one main *Service*, named *Record Manager Service*, to coordinate all *Collection Services*. *Record Manager Service* is triggered anytime either a sensor or audio recording session is scheduled to begin or should be canceled. The *Service* then checks that it is able to and allowed to carry out whatever action it is supposed to carry out and executes it out if it is. For example, if *Record Manager Service* is started with the intention of beginning one of the periodic audio recordings, it will check that it is not in a blackout period and then start the *Audio Collection Service*.




Figure 4: Diagram demonstrating how Record Manager Service uses separate Collection Services for each sensor, allowing for each sensor and its errors to be handled separately.

IV. EXPERIMENTS AND RESULTS

anEAR produces Comma Separated Values (CSV) files for each sensor's recording session. This allows for Users to view to data in any variety of ways including spreadsheet tools such as Microsoft Excel or a computing environment such as Matlab. These files take up very little memory, for example a one minute recording of heart rate data creates a CSV file with

a size of around 1.90 kilobytes. By using the CSV file format, the sensor data that is stored by anEAR can be very easily viewed but also gives the User freedom to use the data however he or she chooses.

anEAR stores the audio recordings in a Waveform Audio File (WAV) format. WAV files contain a 44 byte header followed by a body of uncompressed audio data. Although using the WAV format sacrifices some memory space due to the fact that the audio is not compressed, it allows for a higher quality sound which could be very important for many of anEAR's various use-cases.

In a further attempt to reduce errors, anEAR also contains a CSV file named Audio Record Log. This file contains the start time and date, end time and date, along with other useful information about every audio recording session. This allows Users to go back and ensure that anEAR worked successfully throughout the time it was used and also to see if any audio recordings failed.

Table 2: Resulting data from the data collection study.

Overall Size of Day's Folder	49.95 mB
Number of Audio Recordings	57 WAV files
Average Size of Audio File	465.61 kB/file
Average Size of Accelerometer files	96.18 kB/file (225 files, 21.64 mB)
Average Size of Band Contact files	0.06025 kB/file (319 files, 19.22 kB)
Battery Life of Microsoft Band	~25 hours




Figure 5: Time series graph created using heart rate data recorded by anEAR using the Microsoft Band 2's heart rate sensor.




Figure 6: Time series graph created using GSR data recorded by anEAR using the Microsoft Band 2's GSR sensor.




Figure 7: Time series graph created using RR Interval data recorded by anEAR using the Microsoft Band 2's RR Interval sensor.




Figure 8: Time series graph created using the skin temperature data recorded by anEAR using the Microsoft Band 2's skin temperature sensor.

V. CONCLUSIONS

We developed a smartwatch-based framework called anEAR to record continuous physiological and environmental sensor data, providing insight of how the wearer is interacting in social and personal contexts. anEAR is particularly aimed at tracking emotional arousals in daily lives of individuals who develop post traumatic stress disorder (PTSD). The paper presents the architecture of the system. We successfully demonstrated anEAR in operation to record the multisensor data from smartwatch. Our future plans include the deployment of this app in study participants with PTSD.

ACKNOWLEDGMENT

This work is supported by NIH Grants R01MH108641 and R01MH105379. We are thankful to Cody Goldberg for his help in developing smartwatch apps.

REFERENCES

- [1] Kessler, R.C., Berglund, P., Delmer, O., Jin, R., Merikangas, K.R., & Walters, E.E. (2005). Lifetime prevalence and age-of-onset distributions of DSM-IV disorders in the National Comorbidity Survey Replication. *Archives of General Psychiatry*, 62(6): 593-602.
- [2] Microsoft Band 2 Sensors Table (<http://developer.microsoftband.com/Content/docs/Microsoft%20Band%20SDK.pdf>)
- [3] Pattern Lock Library (<https://github.com/DreaminginCodeZH>)