

Highly Pure, Multi-Epitopic Lipopeptide Vaccine Delivery System: Synthesis and Investigation

Peter M Moyle,^a Colleen Olive,^b Aniko Horváth,^a Levente Karpati,^a
Mei-Fong Ho,^b Melinda Burgess,^b Michael F Good,^b Istvan Toth^{a,*}

^a School of Pharmacy and School of Molecular and Microbial Sciences (SMMS),
the University of Queensland, St Lucia 4072, QLD, AUSTRALIA

^b Queensland Institute of Medical Research (QIMR), Herston, QLD, AUSTRALIA



Queensland Institute of
Medical Research



THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA


www.uq.edu.au

Introduction

Vaccination is the most effective/cost-effective public health intervention

- Disease prevention
- Reduces health care costs
- Reduces lost work time due to sickness

Adapted from: Med Res Rev 1997;17:277



THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

Introduction

Subunit Vaccines

Contain the minimal microbial components necessary to stimulate an appropriate immune response

- Vaccines are administered to healthy individuals (normally children).
- These people are being asked to take a medication when they are well.
- Therefore adverse effects must be minimal.

Advantage

- Removing unnecessary components, reduces the risk of auto-immune diseases and adverse effects.
- Not infectious; No reversion to virulence.
- Can customise the vaccine components to tailor an appropriate immune response.

Problem

- Removing unnecessary components often removes danger signals.
- Need strong adjuvant ('immune stimulating agent').
- In the case of peptides:
 - Small molecular weight limits their capacity to elicit immune responses.
 - Peptides lack the T-helper epitopes required for efficacy in an outbred population.

} **Carrier Molecule**




**THE UNIVERSITY
OF QUEENSLAND**
AUSTRALIA

The Lipid-Core Peptide (LCP) System

Poly-lysine
Multiple Antigen Peptide (MAP) System
(Carrier)
 PNAS 1988;85:5409


Tetrahedron Lett 1993;34:3925

Lipoamino acid
 Liebigs Ann Chem 1990;(12):1175



Peptides
(Antigen)

Lipid Core
(Adjuvant)
 * Mimics Pam₃Cys




THE UNIVERSITY
OF QUEENSLAND
 AUSTRALIA

LCP Synthesis

1. Synthesize LCP Lipid Core

- Using stepwise solid-phase peptide synthesis



C12: 2-amino-D,L-dodecanoic acid (Liebigs Ann Chem 1990;(12):1175)




THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

LCP Synthesis

2. Synthesize Tetravalent MAP


- Using stepwise solid-phase peptide synthesis



LCP Synthesis

3. Synthesize Peptide Antigens

- Using stepwise solid-phase peptide synthesis




THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

LCP Synthesis


3. Cleave Peptide From Resin and Purify

- Cleave using hydrogen fluoride
- Purify by gel filtration/HPLC


Peptide Antigen



Peptide Antigen



Peptide Antigen




Peptide Antigen



THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA


LCP-88/30-J8

$C_{536}H_{918}N_{164}O_{166}S_2$
12380.16g/mol



J8: QAEDK VKQSR EAKKQ VEKAL KQLED KVVQ (28mer)

88/30: DNGKA IYERA RERAL QELGP C (21mer)



A/ 0.1% TFA/H₂O
B/ 90% IPA/0.1% TFA/H₂O

Gradient: 0-100%B over 30min
Flowrate: 1mL/min
Detection: 214nm
Column: Vydac 214TP54 (5μm; 250 × 4.6mm)
 t_R : 19.625 min




THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

LCP-88/30-J8

Systemic IgG Antibody Titers (ELISA)

(Day 56)




Subcutaneous Immunization

B10.BR (H-2^K) mice (n=10) 4-6 week old female
Prime: 30µg LCP-88/30-J8 either 1:1 in CFA or in 50µL PBS
Boost: 3µg in PBS, days 21, 28, 35, 42, & 49

Intra-peritoneal Challenge

400µL (1 × 10⁵ CFU/mL 88/30 GAS)




*, P 0.0001; **, P 0.0001; ***, P 0.0001; ****, P 0.0001; *****, P 0.0001; *****, P 0.002; *****, P 0.001.




LCP-88/30-J8

Systemic IgG Antibody Titers (Intranasal)




Systemic IgG Antibody Titers (Oral)



Intra-peritoneal Challenge

400µL (1-2 × 10⁵ CFU/mL 88/30 GAS)



Intranasal Immunization

B10.BR (H-2^K) mice 4-6 week old female

Prime: 30-60µg LCP-88/30-J8

±10µg cholera toxin B-subunit in 20µL PBS

Boost: days 7, 14, 21, 28, 35

Oral Immunization

Prime: 100µg LCP-88/30-J8

±10µg cholera toxin in 0.2M sodium bicarbonate

Boost: weekly intervals; 5 (with CT) or 7 (without CT)



**THE UNIVERSITY
OF QUEENSLAND**
AUSTRALIA

Project Aims

Lipid Core Peptide System

Advantages:

- High antibody (IgG) titers against attached peptides
- Comparable with the highly toxic adjuvant complete Freund's adjuvant (CFA)
- Potentially safe (non-toxic) for use in humans

Disadvantages:

- Difficult to purify
- ∴ not suitable for use in human clinical trials


Project Aim:

- To develop a method to enable the synthesis of highly pure, easily characterized analogues of the lipid core peptide system
- Techniques to be assessed:
 - Solution- and solid-phase native chemical ligation
 - Fragment condensation



Native Chemical Ligation (NCL)

- Formation of "Native" peptide bond



Curr Opin Biotech 1998;9:412

C-terminal Peptide

- Contains **N-terminal Cysteine**

N-terminal Peptide

- Contains **C-terminal Thioester**

- Aqueous denaturing conditions

- 6M Gdn.HCl

- Urea

- Phosphate buffer


- Performed at pH 7-8

- Minimal side reactions



THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

Poor solubility in aqueous buffers



Based on:

J Org Chem 2000;65(12):3829

Resin:

Diaminodipropylamine (DAPDA) derivatized 4% crosslinked agarose beads (16 μ mol NH₂/mL; Pierce Biotechnology, Rockford IL)

SCAL Linker:

Tetrahedron Lett 1991;32(31):3891

Boc-safety-catch acid labile linker (CSPS pharmaceuticals, San Diego CA)

In oxidised (SO) form:

Stable to TFA, HF, 50% piperidine, Pd(0)

In reduced (S) form:

Cleaved by 50% TFA

Reducing agent:

SiCl₄



THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

Problems with Solid-Phase NCL

Problems:

- Poor solubility of lipid adjuvant in aqueous buffers
- Addition of organic solvents (e.g. TFE, MeCN, DMF, dioxane)
 - Solubilizes lipidic adjuvant
 - Ligation does not occur
- Need excess of thioester peptide to push ligation to completion (wasteful)
- Monitoring of ligation reactions and protecting group removals difficult
 - RP-HPLC provides some quantitative data
- Cleavage of product from the resin is problematic
- The resin is not completely stable to the conditions used for ligation and protecting group removal

Possible Solution to Solubility Issue:

- Use fragment condensation to couple lipid adjuvant to resin, then use NCL to ligate immunogenic peptides



Fragment Condensation

Fragment Condensation:


- Solubility in DMF of lipid adjuvant good.
- Coupling took over 24 hours

Native Chemical Ligation:

- Only 33% complete despite using 2eq thioester peptide
- Subsequent ligations worse


Conclusion:

- Difficult, expensive, and wasteful
- Use solution phase ligation



THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

Solution-Phase NCL



Problems:


Poor solubility of lipid adjuvant in aqueous buffers

Ligation does not occur

Addition of organic solvents (e.g. TFE, MeCN, DMF, dioxane)


- Solubilizes lipidic adjuvant
- Ligation does not occur

NCL + SDS



THE UNIVERSITY
OF QUEENSLAND
AUSTRALIA

Synthesis of a Highly Pure LCP-Analogue




RP-HPLC: Solvent A: 0.1% TFA/H₂O; Solvent B: 90% ACN/0.1% TFA/H₂O; Flowrate: 1mL/min; Column: Vydac C4 (214TP54; 300Å; 5µm, 4.6 x 250mm); Gradient: 0% to 100% B over 30 min; Detection: 214nm; t_R : 18.733 min; Purity: 97.7%.




THE UNIVERSITY OF QUEENSLAND AUSTRALIA

LCP-analogue 1




J Med Chem 2006;49(21):6364

LCP-analogue 2



J Org Chem 2006;71(18):6846

LCP-system







i) THIOL EXCHANGE




MESNA
0.1M Phosphate Buffer pH 7.6

ii) LIGATION

1% (w/v) SDS
TCEP
0.1M Phosphate Buffer pH 7.6



Subcutaneous Immunisation



Mice: 4-6 week old ♀ B10.BR (H-2^k)

Immunised at the tail base

1^o: 30 μ g in 50 μ L PBS or 1:1 CFA

Boosts: 3 μ g in PBS (days 21, 28, 35, 42, & 49)

Conclusions



Photo of Brisbane River Skyline

- Demonstrated a method for the synthesis of highly pure, multi-epitopic, self-adjuvanting lipopeptide vaccines.
 - Required the use of SDS
- May prove useful for the synthesis of multi-epitopic vaccines against diseases caused by other microorganisms.

Acknowledgements



Queensland Institute of
Medical Research

Dr Colleen Olive
Prof Michael Good.
Immunology experiments.



THE UNIVERSITY
OF QUEENSLAND

Prof Istvan Toth
Dr Joanne Blanchfield
Toth group



Australian Government
National Health and Medical Research Council

Research Funding \$\$\$



*The Prince Charles
Hospital Foundation*

Research Funding \$\$\$



**Heart
Foundation**

Research Funding & Travel Grant \$\$\$



**Queensland
Government**

Research Funding \$\$\$

This research is proudly supported by the Queensland Government's Growing the Smart State PhD Funding Program. The State of Queensland accepts no responsibility for decisions or actions resulting from any information supplied. The views and information contained in the research do not necessarily represent the views or opinions of the Queensland Government and carry no endorsement by the Queensland Government.



**THE UNIVERSITY
OF QUEENSLAND**
AUSTRALIA