



# OSTEOGENESIS OF DENTAL STEM CELLS

NOVA SOUTHEASTERN UNIVERSITY COLLEGE OF DENTAL MEDICINE


BIOL4990-INDEPENDENT STUDY

DR. UMADEVI KANDALAM

JERRY ENNOLIKARA

# RESEARCH LABORATORY OVERVIEW

- This research lab is dedicated to the study of bioengineering bone tissue, an application of biology research within medicine known as “regenerative medicine”, through a variety of in-vivo and in-vitro projects.
- The primary focus of our lab is to replace a lifetime of bone graft surgeries associated with cleft palate, a congenital disorder where the palatal process of the maxilla fails to fuse during development, with a noninvasive injectable treatment
- The lab upholds the central motif of tissue engineering of “cells, scaffolds, and growth factors” with our research utilizing periodontal ligamental stem cells/human gingival mesenchymal stem cells; 3D hydrogel scaffolds such as Puramatrix™ and Viagraft™ or bone powder; and osteogenic growth factors such as BMP-2 and more recently, FGF and VEGF in our newer studies
- The research is focused in a dental maxillofacial model yet also has many implications for medicine within the orthopedic field.




# STUDENT RESEARCHER EXPERIENCES

- The lab is located primarily within the College of Dental Medicine, yet students ranging from the pre-dental to the pre-medical track work in the lab.
- Student researchers have the opportunity to work and learn in many cross-curricular laboratory processes/techniques and connect the skills learned from the lab sciences to the real dental and medical field, such as:
  - Genetics: PCR, Real Time PCR
  - Biochemistry: Western blotting, Enzymatic assaying
  - Microbiology: Cell culturing, Cell assaying, Cell staining
  - Histology: Microscopy, Slide preparation, Tissue analysis
- Most students researcher participate in research-for-credit by setting up and registering for a BIOL4990 Independent Study course section (1-3 credits) with the Biology department
- Each credit is valued to three hours of lab work for every week, or else otherwise as coordinated with Dr. Kandalam (or whoever your research professor is; this process is general for most research-for-credit courses)

# Isolation and Osteogenic Differentiation of Dental Stem Cells in a Fibrous Hydrogel

Dr. Arshad Ali & Jerry Enslin  
Dr. Umadevi Kishan  
Southwestern University



## Cleft Palate



- Most common birth defect
- One in 700 births in the US
- Incomplete fusion of two palatal shelves
- Only current treatment: multiple surgeries within first two years, multiple follow up surgeries throughout life

## Methodology

- PDLSCs were isolated and tested for surface markers.
- Cell proliferation was measured by WST assay over intervals of time (1, 2, 3 days)
- Viability was measured with Live/Dead cell assays.
- Gene expression of alkaline phosphatase (ALP) was measured by qPCR
- Differentiation was observed by measuring ALP enzyme activity.

## Methodology

- PDLSCs were isolated and tested for surface markers.
- Cell proliferation was measured by WST assay over intervals of time (1, 2, 3 days)
- Viability was measured with Live/Dead cell assays.
- Gene expression of alkaline phosphatase (ALP) was measured by qPCR
- Differentiation was observed by measuring ALP enzyme activity.



## Detection of Surface Markers: Flow Cytometry




- International Society for Cellular Therapy (ISCT) sets standard for a cell population to be MSCs
- Surface markers CD105+, CD73+, and CD90+ must be identified from flow cytometry

## Cell Viability: Live/Dead Cell Assay

- Images captured highlight cell viability/death
- Light sensitive assay utilizing calcein-AM and ethidium homodimer (EthD) to dye live and dead cells, respectively.
- Composite images show only live cells

## Methodology

- PDLSCs were isolated and tested for surface markers.
- Cell proliferation was measured by WST assay over intervals of time (1, 2, 3 days)
- Viability was measured with Live/Dead cell assays.
- Gene expression of alkaline phosphatase (ALP) was measured by qPCR
- Differentiation was observed by measuring ALP enzyme activity.



## Osteogenic Differentiation: ALP Enzyme Activity

- Alkaline phosphatase (ALP) enzyme activity indicates differentiation
- pNPP assay indicates amount of ALP enzyme activity
- Cells treated for 1 week, ALP activity was measured on day 7

## Long Term Goal

Repair and regenerate bone in the cleft palate region.

## Study Objective

Develop injectable scaffold systems with stem cells to regenerate bone tissue.

## Methodology

- PDLSCs were isolated and tested for surface markers.
- Cell proliferation was measured by WST assay over intervals of time (1, 2, 3 days)
- Viability was measured with Live/Dead cell assays.
- Gene expression of alkaline phosphatase (ALP) was measured by qPCR
- Differentiation was observed by measuring ALP enzyme activity.

## Methodology

- PDLSCs were isolated and tested for surface markers.
- Cell proliferation was measured by WST assay over intervals of time (1, 2, 3 days)
- Viability was measured with Live/Dead cell assays.
- Gene expression of alkaline phosphatase (ALP) was measured by qPCR
- Differentiation was observed by measuring ALP enzyme activity.

## P-Nitrophenylphosphate (pNPP) Assay

- Hydrolysis of p-Nitrophenylphosphate (pNPP) by alkaline phosphatase creates p-NITROPHENOL
- Cells grown in complete medium (CM) and osteogenic medium (OM) treated with pNPP
- Absorbance results indicate ALP enzyme activity

## Cells: Periodontal Ligament Stem Cells

- Mesenchymal stem cells (MSCs) associated with periodontal regeneration
- Basic cell type used throughout project

## Cell Proliferation: Encapsulation in Hydrogel

- 10,000x, 50,000x, and 300,000x PDLSCs were added to Puramatrix™ solution
- Culture medium added for peptide units to assemble into 3D matrix
- 3-D encapsulated cells were grown and used for WST assay

## Differentiation: Complete vs. Osteogenic Medium

- Two million cells per gel were encapsulated.
- Cells grown in complete medium (CM) and osteogenic medium (CM with Ascorbic acid,  $\beta$ -Glycerolphosphate, Dexamethasone)
- Samples were tested for differentiation markers.



## Scaffold: Hydrogel

Puramatrix™ (PM) is a commercially available self-assembling peptide units form nanofiber structures (~100nm in diameter) similar to the natural extracellular matrix.

## WST Assay

- Colorimetric assay performed after 1, 2, and 3 days
- Dehydrogenase enzymes reduce WST-1, a tetrazolium salt
- Depth of orange color from metabolite indicates amount of proliferation

## Osteogenic Differentiation

- Mesenchymal stem cells differentiate into mature bone cells. This leads to bone formation.
- Differentiation begins

## Conclusion

- Puramatrix™ scaffolds promote cell growth and differentiation
- This tissue engineering approach shows the potential for bone regeneration

# HOW TO GET INVOLVED

- The lab is open to any undergraduate students who are passionate about the lab sciences.
- To get involved, please email Dr. Umadevi Kandalam at [umadevi@nova.edu](mailto:umadevi@nova.edu) with your resume and a letter of interest. From there, she or someone else will get back to you to potentially set up a date to interview with Dr. Kandalam herself and see the lab!
- If you have any other questions, please email me (Jerry Ennolikara) at [je788@mynsu.nova.edu](mailto:je788@mynsu.nova.edu)
- Additional contacts:
  - Nicole DeLorenzo: [ndelorenzo@mynsu.nova.edu](mailto:ndelorenzo@mynsu.nova.edu) (Lab Manager)
  - Shreya Patel: [sp1743@mynsu.nova.edu](mailto:sp1743@mynsu.nova.edu) (Student Researcher)