

PROCEEDINGS OF SPIE

SPIDigitalLibrary.org/conference-proceedings-of-spie

Extending the wavelength range of multi-spectral microscope systems with Fourier ptychography

Antony Chi Shing Chan, Cheng Shen, Elliott Williams, Xiaoyu Lyu, Hangwen Lu, et al.

Antony Chi Shing Chan, Cheng Shen, Elliott Williams, Xiaoyu Lyu, Hangwen Lu, Craig Ives, Ali Hajimiri, Changhui Yang, "Extending the wavelength range of multi-spectral microscope systems with Fourier ptychography," Proc. SPIE 10890, Label-free Biomedical Imaging and Sensing (LBIS) 2019, 108902O (25 March 2019); doi: 10.1117/12.2510875

SPIE.

Event: SPIE BiOS, 2019, San Francisco, California, United States

Extending the wavelength range of multi-spectral microscope systems with Fourier ptychography

SPIE.

Antony Chi Shing Chan¹, Cheng Shen¹, Elliott Williams¹, Xiaoyu Lyu², Hangwen Lu³, Craig Ives¹, Ali Hajimiri¹ and Changhui Yang^{1*}

Caltech

¹Department of Electrical Engineering, Caltech (United States) ²Academy of Opto-Electronics, CAS (China) ³Google LLC
*chyang@caltech.edu

Paper No.: 10890-96

Introduction

Limited by manufacturing and optical materials, even the best designed multi-spectral objective will exhibit residual monochromatic aberration and chromatic aberration. Here we propose a multi-spectral imaging scheme based on the Fourier Ptychographic Microscopy (FPM), which can render in-focus images by digital aberration correction. Also, the wavelength extension to Near Infrared (NIR) regime broadens the application of our system, like imaging through bonded silicon wafer devices for locating defects.




Fig. 1 Schematic of multi-spectral Fourier ptychographic microscope, including transmission and reflection mode.

Aberration Characterization

Principle

$$\arg \min_w L(\tilde{B}_{\text{on-axis}}(u, v) e^{i2\pi z(u, v)} - \tilde{B}_{\text{off-axis}}(u, v))$$

$$z(u, v) = \sum_{i=1}^8 w_i Z_i(u, v)$$




Fig. 2 Digital aberration calibration based on Zernike mode decomposition, taking the complex amplitude of an on-axis microbead as reference.

Aberration Characterization

Result

For an objective with a low NA but a large field-of-view (FOV), its aberration is spatially variant. Usually, aberrations increase as a function of distance from the image center.

$$w_i(x, y), i = 1, \dots, 8$$




Fig. 3 Spatially varying aberration characterization for the whole FOV under the wavelength of 532nm. (a)-(f) fitted 2D surface for each type of aberration (first order Zernike modes are normally not considered as aberrations, thus not shown).

Digital Aberration Correction

Since the location-dependent aberration composed of different Zernike modes has been calibrated first, digital aberration correction can be done for other samples to obtain the clear and undistorted images.

1. A large FOV (~1.92x1.53mm²)
2. A better way of pixel super-resolution




Fig. 4 Demonstration of digital aberration characterization under the wavelength of 1552nm. (a) the whole FOV raw data. (b) Comparison between the pixel upsampling result and FPM reconstruction without/with aberration correction.

Imaging of Silicon Chip

With our system, the non-invasive detection on the silicon wafer devices, like flip-chips, is possible. Usually, the front-side imaging could directly visualize the top metallic layers, but the ones deep inside can only be detected from the backside with the aid of NIR light thanks to the fairly low absorption of silicon to this regime. Besides, FPM reconstruction can further enhance the imaging resolution and digitally refocus to different components beneath.




Fig. 5 Imaging of a flip-chip under reflection mode. (a) the schematic of the flip-chip. (b)-(c) the front and back side imaging of the flip-chip. The bottom metal layers could only be visible under NIR light. Different components are enlarged in the enlargements.

Conclusions

Our multi-spectral FPM system enables:

1. Digital aberration correction
2. Multi-spectral high-resolution imaging with large FOV
3. Digital refocusing
4. Non-invasive defect detection of silicon wafer devices

References

1. Zheng, G., Ou, X., Horstmeyer, R., & Yang, C. (2013). Characterization of spatially varying aberrations for wide field-of-view microscopy. *Optics express*, 21(13), 15131-15143.
2. Zheng, G., Horstmeyer, R., & Yang, C. (2013). Wide-field, high-resolution Fourier ptychographic microscopy. *Nature photonics*, 7(9), 739-745.
3. Ou, X., Chung, J., Horstmeyer, R., & Yang, C. (2016). Aperture scanning Fourier ptychographic microscopy. *Biomedical optics express*, 7(8), 3140-3150.