

LA--10544

LA-10544

DE86 005840

UC-10

Issued: November 1985

## Determining Yttrium in Plutonium by Anion-Exchange X-Ray Fluorescence

Calvin J. Martell  
James M. Hansel

### DISCLAIMER

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

Los Alamos Los Alamos National Laboratory  
Los Alamos, New Mexico 87545

  
DISTRIBUTION OF THIS DOCUMENT IS UNLIMITED

# DETERMINING YTTRIUM IN PLUTONIUM BY ANION-EXCHANGE X-RAY FLUORESCENCE

by

Calvin J. Martell and James M. Hansel

## ABSTRACT

This report describes a method for determining yttrium in plutonium using an anion-exchange separation and x-ray fluorescence. We add zirconium to the plutonium solution as an internal standard. We oxidize the plutonium to Pu+4 and pass the solution through an anion-exchange column with 8M HCl. The Pu+4 sorbs to the resin and the yttrium and zirconium pass through completely. We evaporate the eluate solution containing the yttrium and zirconium and transfer it to a 10-ml volumetric flask. We add a portion of this solution to an x-ray cell and measure the  $K\alpha$  x-ray line for both yttrium and zirconium. The ratio of yttrium to zirconium is then compared with standards. This method has a precision of 0.84% relative standard deviation for yttrium over a concentration range of 0.5 to 3.5 mg in a 10-ml volume.

## INTRODUCTION

Yttrium, a detector element added to plutonium, is used to measure neutron energies of underground nuclear explosions. To properly determine the strength of the nuclear device, we must know the concentration of the yttrium in plutonium.

It is difficult to accurately determine a small amount of yttrium in several hundred milligrams of plutonium. Therefore, we decided to separate the yttrium from the plutonium. The sample can then be handled more easily and safely with the plutonium removed.

To achieve the best precision possible, we add zirconium as an internal standard at the beginning of the analysis. Because it behaves like yttrium, zirconium is able to compensate for any changes that may occur in the sample during processing.

We transfer the eluate solution containing both yttrium and zirconium to a 10-ml volumetric flask. We add a portion of this solution to an x-ray cell for the analysis. By working with solutions, we reduce the analysis time compared with the time needed to evaporate the solution to a small volume and then pipette it onto a resin paper disc.

## Equipment and Reagents

In this work, we used the following equipment:

- Aluminum cell holder, 50-mm o.d.
- Columns, glass, ion exchange, 32 by 200 mm
- Glass rod, bent and flattened on one end
- Heat lamps, infrared, red end, 250 W
- Hot plate, 150 by 180 mm
- Mylar film, 6- $\mu$ m thick
- Sample cells, x-ray, Chemplex, plastic, series 1400, 32-mm o.d.
- X-ray spectrograph; vacuum, Philips Electronics, Inc., LiF crystal; tungsten-target x-ray tube; and scintillation detector with associated electronic counting circuitry

We also used the following reagents:

- Hydrochloric acid, 12M, 8M, and 0.1M
- Hydrofluoric acid, 1.3M
- Ion-exchange resin, analytical grade anion-exchange resin (AG1X4), 50 to 100 mesh, chloride form, Bio-Rad Laboratories

- Nitric acid, 15.7M
- Yttrium solution, 0.500 mg/ml
- Zirconium solution, 0.500 mg/ml

## SEPARATION OF YTTRIUM

To dissolve the metal, we place a 1-g sample of plutonium in a covered 50-ml beaker and add 5 to 10 ml of 8M HCl. (If other plutonium compounds need to be dissolved, see our previous report.<sup>1</sup>) We transfer the dissolved plutonium to an appropriate size volumetric flask and add 5 ml of 15.7M HNO<sub>3</sub> to the flask to help oxidize the Pu+3 to Pu+4. Then we add water to bring the plutonium solution to the flask mark.

The weight of the plutonium sample, the size of the volumetric flask, and the size of the aliquot all depend on the amount of yttrium in the plutonium. If we know the approximate yttrium concentration, we can arrange these three variables so that the amount of plutonium will provide sufficient yttrium to fall within the range of our standards. If the approximate yttrium concentration is not known, we take one small aliquot in case the concentration is high, and one large aliquot in case the yttrium concentration is low.

We add 2 ml of the 0.500-mg/ml zirconium solution to a 20-ml beaker. After having decided on the appropriate size of the aliquot, we add it to the 20-ml beaker along with the internal standard, zirconium. To this solution, we add 0.5 ml of 15.7M HNO<sub>3</sub> to ensure that all the Pu+3 has been converted to Pu+4. We heat the covered 20-ml beaker with a heat lamp. The beaker is covered because there is effervescence at this point. (The effervescence can be greatly reduced by allowing the beaker to sit uncovered overnight.) When the effervescence has stopped, we rinse off the cover into the beaker with water and continue heating to incipient dryness. We add 2-3 ml of 8M HCl to the beaker and dissolve the sample, take it to incipient dryness once more, and redissolve the sample in 5-10 ml of 8M HCl.

We fill the anion-exchange column with resin to within 2-3 cm of the top and condition the column with 8-10 ml of 8M HCl four times. We place a 20-ml beaker under the column to catch the solution that will pass through and then add the plutonium solution to the top of the column by pouring the solution down a bent and flattened glass rod. After the plutonium solution has stopped flowing, we rinse the 20-ml sample beaker with 2-3 ml of 8M HCl and a few drops of 1.3M HF and add this to the column. We continue to rinse with the 8M HCl until 15-20 ml of solution has collected in the beaker under the column.

Passing the sample through the column allows the Pu+4 to sorb to the resin and allows the yttrium and zirconium to pass completely through. The flattened

glass rod helps keep the plutonium near the top of the column so more resin is available to sorb the plutonium. This process achieves our separation. Then we need only start the evaporation of the yttrium and zirconium solution and recover the plutonium.

We evaporate our solution so that it will go into a 10-ml volumetric flask. We place an uncontaminated paper towel on the floor of the plutonium open-front box. Then working "hot" with one glove-covered hand, we pour the yttrium and zirconium solution into an x-ray cell held by our other glove-covered hand, which touches only uncontaminated objects. In this way, we transfer the essentially Pu-free but still contaminated solution to an x-ray cell, cover the cell with Mylar film, and then place it into our x-ray spectrograph with very little likelihood of contaminating our equipment.

We place the covered x-ray cell containing our sample solution in an aluminum cell holder, which also has a Mylar cover. The sample solution is excited by the tungsten-target x-ray tube. We measure, for 60 s each, the following 2θ settings:

2θ	Measurement
22.00	Background 1
22.48	Zr Ka
23.00	Background 2
23.72	Y Ka
24.50	Background 3

Figure 1 shows a plot of intensity vs 2θ for yttrium and zirconium after separation from a 250-g sample of plutonium by anion exchange. The plot shows 1 mg of zirconium and 2 mg of yttrium. The intensity of yttrium is approximately 9000 c/s.




Fig. 1. Intensity vs 2θ for 1 mg zirconium and 2 mg yttrium after their separation from a 250-g sample of plutonium.

We analyze the samples with the spectrographic conditions shown in the following table:

X-Ray Operating Parameters	
X-ray tube	Tungsten
Potential	50 kV
Current	30 mA
Counter	Scintillation 960 V
Analyzing crystal	LiF
Counting time	60 s

### STANDARDS

We prepared standards containing known amounts of yttrium and zirconium as listed below.

Standard No.	Yttrium (mg)	Zirconium (mg)	Plutonium (mg)
Y-1	0.0	1.00	100
Y-2	0.0	1.00	100
Y-3	1.0	1.00	100
Y-4	1.0	1.00	100
Y-5	3.0	1.00	100
Y-6	3.0	1.00 </td <td>100</td>	100

We place 10 g of high-purity plutonium in a covered 250-ml beaker and dissolve with 8M HCl and a small amount of 1.3M HF. Next, we transfer the plutonium solution to a 100-ml volumetric flask and add 5 ml of 15.7M HNO<sub>3</sub> to the flask to oxidize the plutonium to Pu<sup>+4</sup>.

Next, we pipette the appropriate amounts of the yttrium solution into each of six 20-ml beakers. Then we add 1.0 mg of the internal standard, zirconium, to each of the six beakers and pipette an appropriate amount of the high-purity plutonium solution, that is, 100 mg, into each of the beakers.

We analyze these standards as described in the "Separation of Yttrium" section. Figure 2 shows a plot of intensity ratios (Y/Zr) vs concentration of yttrium. These data have a correlation coefficient for fit = 0.9998.

### CALCULATIONS

To do the background corrections, we convert the intensities for all the readings to counts per second (c/s)




Fig. 2. Intensity vs concentration of yttrium.

and average the c/s for background 1 and background 2 and subtract this average from the zirconium c/s. Then we average background 2 and background 3 and subtract this average from the yttrium c/s. To obtain the Y/Zr ratio, we divide the yttrium c/s by the zirconium c/s.

We do a least squares fit for the standards relating Y/Zr ratios vs the yttrium concentrations. We use the coefficients for the equation representing this standard curve to calculate the concentration of yttrium in the analyzed samples.

### PRECISION

To establish the precision of the method, we analyze samples repeatedly. Some of the samples were pure plutonium with known amounts of yttrium added. The results are shown in Table I.

The average relative standard deviation for 0.5 to 3.5 mg yttrium is 0.84.

TABLE I. Precision of Method

Yttrium (mg)	Number of Determinations	Relative Standard Deviation (%)
0.5	6	1.13
1.0	6	0.90
1.5	6	0.67
3.5	6	0.64

## ACCURACY

An official yttrium reference solution or compound was not available. We had  $Y_2O_3$  assayed and analyzed. Using this well characterized material, we prepare a concentrated solution of yttrium (12.5 mg/ml) and use this to prepare the yttrium solution (0.5 mg/ml) used in this work.

We add pure plutonium and zirconium to known amounts of yttrium at three levels and treat these solutions as described in the "Separation of Yttrium" section. The results are given in Table II.

The average recovery for the range 0.500 to 3.000 mg yttrium is 100.1.

Yttrium Added (mg)	Yttrium Found (mg)	Yttrium Recovery (%)
0.500	0.504	100.8
1.500	1.499	99.9
3.000	2.991	99.7

## DISCUSSION

The precision of the Y/Zr ratio is very good while that for just yttrium is relatively poor. The zirconium behaves so much like yttrium that whatever causes the variations for yttrium also causes them for zirconium.

The amount of plutonium taken as a sample is one variable in this method. The relatively large ion-ex-

change columns (32 by 200 mm) in use permit as much as 500 mg of plutonium to sorb on to the resin. At this level, approximately 1% of the plutonium passes through the resin along with the yttrium and the zirconium. At about 250 mg of plutonium, somewhat less than 1% passes through the column.

When the sample size is between 200 and 10 mg plutonium, the retention of plutonium is nearly 100%. This figure assumes filling the column with resin to within 2-3 cm from the top and using the bent and flattened glass rod to keep the plutonium at the top of the column of resin.

Use of large sample sizes (250 to 500 mg) causes a deterioration of precision to approximately 5%. The reason for this is not known, but we speculate that the plutonium on the column upsets the relationship between yttrium and zirconium so that they no longer behave identically.

When a relatively small sample size is used and nearly 100% of the plutonium has been removed, the standards and samples may be read on an energy dispersive x-ray unit. The x-ray tube target should be other than molybdenum (perhaps rhodium or silver) because the molybdenum K $\alpha$  line will excite the zirconium only half as well as the yttrium line. This is because of the position of the absorption edges of these two elements.

## REFERENCE

1. Calvin J. Martell and James M. Hansel, "Anion Exchange X-Ray Fluorescence Determination of Uranium in Plutonium," Los Alamos National Laboratory report LA-9552 (February 1983).