

LBL-36511
UC-401

Hydrogen-Ion Titrations of Amino Acids and Proteins in Solutions Containing Concentrated Electrolyte

Florian Fergg,¹ Daniel E. Kuehner, Harvey W. Blanch, and John M. Prausnitz*

Department of Chemical Engineering, University of California
and
Chemical Sciences Division, Lawrence Berkeley Laboratory
University of California, Berkeley, California 94720

¹Visiting Researcher, Technische Universität München,
Lehrstühle B für Verfahrenstechnik

*Author to whom correspondence should be addressed

December 1994

MASTER

This work was supported by the Director, Office of Energy Research, Office of Basic Energy Sciences,
Chemical Sciences Division, of the U.S. Department of Energy under Contract No. DE-AC03-76SF00098.

DISTRIBUTION OF THIS DOCUMENT IS UNLIMITED

fr

DISCLAIMER

This report was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor any agency thereof, nor any of their employees, make any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof.

DISCLAIMER

Portions of this document may be illegible in electronic image products. Images are produced from the best available original document.

Contents

Abstract.....	3
1 Introduction	4
2 Theory.....	6
3 Experimental Methods.....	12
3.1 Titration Apparatus	12
3.2 Software Description	15
3.2.1 CAL-TITR.....	15
3.2.2 CAL-ANAL.....	17
3.2.3 TITRATE.....	17
3.2.4 SETUPBTA.....	18
3.2.5 BETA90.....	18
3.3 Experimental Procedure.....	19
4 Results/Discussion	20
5 Conclusion	34
Appendix	
A Amino-Acid Sequences of Lysozyme and α-Chymotrypsin.....	35
B Electrode Calibration	37
C Details on Dosimat Pump Configuration and Operation	39
D Lid Construction.....	40
E Details of Apparatus Troubleshooting.....	41
F Program Flow Sheet for CAL-TITR and TITRATE.....	42

G	Setup and Operation of BETA90.....	44
	G.1 SETUPBTA Inputs.....	44
	G.2 Refining Amino-Acid pKa's with BETA90.....	45
H	Experimental Procedure for Titration	47
I	Experimental Data for all Measured p^oK_a.....	53
K	Sample Titration-Data Recording.....	56
	Nomenclature	59
	References	60

Abstract

This report describes a first attempt to quantify the net charge as a function of solution pH for lysozyme and α -chymotrypsin at 0.1 M, 1.0 M and 3.0 M ionic strength, (IS). The calculations are based on the residue (titratable group) pK_a 's in the amino-acid sequence of the protein. To determine these pK_a 's, a simple theory was used which assumes that the pK_a 's are independent from each other in the protein and are equal to their pK_a values in free amino-acid solution (Independent-Site Theory, IST).

Residue pK_a 's were obtained from amino-acid hydrogen-ion titrations at three different KCl concentrations corresponding to 0.1M, 1.0M and 3.0M ionic strength. After construction of a suitable apparatus, the experimental procedure and data reduction were computerized to perform a large number of titrations. Most measured pK_a 's showed high reproducibility (the difference of pK_a values observed between two experiments was less than 0.05). For IS = 0.1M, observed pK_a 's agreed with literature values to within a few hundredths of a pH unit. Furthermore, the ionic-strength dependence of the pK_a 's followed the trends reported in the literature, viz. pK_a values decrease with increasing ionic strength until they reach a minimum at about IS = 0.5M. At still higher IS, pK_a 's increase as the ionic strength rises to 3M.

The known pK_a 's of all titratable groups in a protein were used with the IST to give a first approximation of how the protein net charge varies with pH at high ionic strength. A comparison of the titration curves based on the IST with experimental lysozyme and α -chymotrypsin titration data indicates acceptable agreement at IS = 0.1M. However, comparison of measured and calculated titration curves at IS = 1M and IS = 3M indicates only qualitative agreement.

1 Introduction

Protein precipitation is a convenient method for concentrating protein from dilute solutions. Salting-out, the precipitation of protein by high electrolyte concentration, is a commonly used technique. Experimental work and theoretical modeling are conducted to gain a quantitative understanding of the nature of protein interactions in concentrated salt solutions [1-3]. To support modeling efforts, it is important to quantify the dependence of protein net charge on solution pH at high salt ionic strength, e.g. 3M. The current literature provides charge/pH information for some proteins only at low ionic strengths, e.g. 0.01 - 0.2M [4].

Hydrogen-ion titration experiments yield the relative net charge of a protein as a function of pH. To determine absolute net charge, it is necessary to know one charge/pH point exactly. An experimental determination of the pH point at which the protein has zero net charge (the isoelectric point, pI) can be made by isoelectric focusing (IEF). However, in high ionic-strength solutions, conventional IEF is not possible. It is necessary to have an independent method for measuring or calculating the pI at high salt concentrations.

In an attempt to calculate the pI's of chymotrypsin and lysozyme in high IS solutions, we used a simple theory for calculating the net charge as a function of pH. We assume that the overall charge of the protein is equal to the sum of the charges of all the protonatable groups in amino-acid sequence of this protein. With knowledge of all the individual hydrogen ion dissociation equilibrium constants, the K_a 's, of the protonatable sites in a protein, it is simple to calculate the absolute net charge of a protein as a function of pH. Determination of individual amino-acid K_a 's in proteins is difficult and is a topic of much research, especially for residues buried in the interior of proteins [5]. Our independent-site theory (IST) makes the (large) assumption that the K_a 's of amino acids in the protein are equal to their K_a 's in free solution. These are measurable in high ionic strength by performing hydrogen-ion titrations of amino acids in solutions containing concentrated electrolyte. The amino-acid sequence of the protein of interest determines which amino acids need to be titrated to find their K_a 's for theoretical net-charge calculations.

In this work we present a detailed description of an apparatus that we constructed to perform high-precision hydrogen-ion titrations. We present and discuss the result of

amino-acid titrations in solutions of KCl at 0.1M, 1.0M and 3.0M ionic strength. Also presented are preliminary results for titrations of lysozyme and α -chymotrypsin in a few different solutions containing KCl at high ionic strength. A comparison is given for experimental and theoretical protein net-charge dependence on solution pH. Results and applicability of the IST are discussed.

2 Theory

Each protein molecule has in its amino-acid sequence many different residues with potential hydrogen-ion equilibria. At any given pH, these residues will be in various protonation states, depending on their individual hydrogen-ion dissociation constant, or pK_a 's. Thus the protein will have a net charge which varies with pH. The dissociation of hydrogen ions will also be affected by the presence of added electrolyte, and therefore the net charge/pH behavior of a protein will be a function of solution ionic strength.

Protein hydrogen-ion titration experiments yield relative net charge as a function of pH. To determine absolute net charge, it is required to know one charge/pH point exactly. In solution with low concentration of added electrolyte an experimental determination of the pH point (the protein has zero net charge) can be made by isoelectric focusing (IEF) [6]. However, in high ionic-strength solutions, conventional IEF is not possible. An independent method is necessary for either measuring or calculating the isoelectric point (pI) of the protein .

With knowledge of all the individual pK_a 's of all the possible hydrogen-ion equilibria in a protein, it is possible to calculate the absolute net charge of a protein as a function of pH. Experimental determination of the pK_a 's of equilibria in the interior of proteins is difficult and is a topic of much research [5]. In an attempt to calculate the pI's of chymotrypsin and lysozyme in high ionic-strength solutions, we used an independent-site theory (IST). Our independent-site theory makes the (large) assumption that the pK_a 's of amino acids in the protein are equal to their free-solution values, and that each titratable site is independent of all others.

Independent-Site Theory

In every protein molecule there are n_i acid/base equilibria of type i , each described by a hydrogen-ion dissociation constant K_{ai} such that, in dilute solutions where all activity coefficients are unity [7]:



$$K_{ai} = \frac{[H^+] \cdot [A^-]_i}{[HA]_i} \cdot \frac{\gamma_{H^+} \cdot \gamma_{A^-}}{\gamma_{HA_i}} \quad (2)$$

and

$$K_{ai}^o = \frac{[H^+] \cdot [A^-]_i}{[HA]_i} = \frac{K_{ai}}{K_\gamma} \quad (3)$$

The individual n_i 's are known from the amino-acid sequence of a given protein. These equilibria occur on the side chains of acidic or basic residues and at the amino and carboxyl termini of the amino-acid chain in the protein. Taking the logarithm of equation (3) yields the well-known Henderson-Hasselbach equation [7]:

$$p^o H - p^o K_{ai} = \log \frac{[A^-]_i}{[HA]_i} \quad (4)$$

where

$$p^o H = -\log[H^+] \quad (5)$$

$$p^o K_{ai} = -\log K_{ai}^o \quad (6)$$

A mass balance on titratable sites of type i dictates that the total concentration of sites (in moles of sites per liter of protein solution) must equal the sum of the concentrations of unprotonated and protonated sites of type i :

$$[S]_i = [A^-]_i + [HA]_i \quad (7)$$

Here $[S]_i = n_i \cdot [P]$, where $[P]$ is the total protein concentration, in moles/liter. Combining equations (4) and (7) yields expressions for the molar concentrations of protonated and unprotonated sites of type i in terms of the quantities pH, pK_{ai} , and n_i :

$$[A^-]_i = [S]_i - \left(\frac{[A^-]_i}{10^{(p^o H - p^o K_{ai})}} \right) = \frac{[S]_i \cdot 10^{(p^o H - p^o K_{ai})}}{(10^{(p^o H - p^o K_{ai})} + 1)} \quad (8)$$


and

$$[HA]_i = [S]_i \cdot \left(1 - \frac{10^{(p^o H - p^o K_{ai})}}{(10^{(p^o H - p^o K_{ai})} + 1)} \right) \quad (9)$$


Note that species A_i^- is not necessarily negatively charged; it may be either negative or neutral, depending on the particular dissociation equilibrium. Therefore, we define a function $f(x)$, which accounts for the valence of each species, as:

$$f(x) = \text{charge on species } x.$$

For an acidic side chain, e.g. aspartic acid,



and for a basic side chain, e.g. lysine,



Then the net charge on a single protein molecule is:

$$z_{net} = \sum_i z_i = \frac{1}{[P]} \sum_i \left(f(A_i^-) \cdot [A^-]_i + f(HA_i) \cdot [HA]_i \right) \quad (10)$$

or

$$z_{net} = \sum_i \left[f(A_i^-) \cdot \left\{ \frac{n_i 10^{(p^{\circ}H - p^{\circ}K_{ai})}}{10^{(p^{\circ}H - p^{\circ}K_{ai})} + 1} \right\} + f(HA_i) \cdot n_i \cdot \left\{ 1 - \frac{10^{(p^{\circ}H - p^{\circ}K_{ai})}}{10^{(p^{\circ}H - p^{\circ}K_{ai})} + 1} \right\} \right] \quad (11)$$

This expression allows us to calculate the net charge on a protein as a function of pH and $\{pK_{ai}, n_i\}$ within the abovementioned assumptions of this independent-site theory.

Of the 20 natural amino acids found in most proteins, there are seven which have side chains with hydrogen-ion equilibria [7], namely:

arginine	cysteine	histidine	tyrosine
aspartic acid	glutamic acid	lysine	

All of these side-chain equilibria, except the cysteine sulfhydryl equilibrium, are found in lysozyme and α -chymotrypsin. Also, the amino-acid sequences of lysozyme and

α -chymotrypsin (see Appendix A) give the identity of titratable carboxyl and amino end groups on all the protein chains. Note α -chymotrypsin has three chains linked via cysteine disulfide bridges, and thus six titratable end groups, while lysozyme has only one unbroken chain and thus 2 titratable end groups. The relevant end-group equilibria are:

alanine (amino)	isoleucine (amino)	tyrosine (carboxyl)
asparagine (carboxyl)	leucine (carboxyl)	
cysteine (amino)	lysine (amino)	

To calculate protein net charge as a function of pH, it is necessary to measure the pK_a 's of all the relevant equilibria mentioned above via amino-acid hydrogen-ion titration. These titrations must be carried out in amino-acid solutions with ionic strengths corresponding to those found in our theoretical model. Also, protein titrations must be performed at solution conditions identical to those for amino-acid titrations to compare experimental protein-charge vs. pH curves to curves calculated from amino-acid pK_a 's via IST.

Effect of Ionic Strength on Hydrogen-Ion Equilibria: Calibration Titration

In solutions of high ionic strength, care must be exercised in interpretation of experimental pH data. The pH of a solution, is strictly defined as [8]:

$$pH = -\log(a_{H^+}) \quad (12)$$

Hydrogen-ion activity, a_{H^+} , is the quantity observed in a pH measurement [8]. In dilute solutions, it is commonly assumed that the activity of the hydrogen ion is equal to its concentration [9]. In concentrated salt solutions, however, the hydrogen-ion concentration and activity are not equal; they are related by the activity coefficient, γ_{H^+} . Hence pH in concentrated salt solutions is:

$$pH = -\log(a_{H^+}) = -\log(c_{H^+} \cdot \gamma_{H^+}) \quad (13)$$

Recall from equation 5 that quantity p^0H is equal to the pH only in the infinite-dilution limit, when the concentration and activity of hydrogen ion are equal:

$$p^{\circ}H = -\log(c_{H^+}) \quad (14)$$

In general, and particularly for hydrogen ions in a concentrated salt solution:

$$pH = p^{\circ}H - \log(\gamma_{H^+}) \quad (15)$$

All theoretical calculations outlined above require c_{H^+} to determine net charge on a protein as a function of $p^{\circ}H$. Thus, it is necessary to relate titration data, which are based on activity measurements, to hydrogen-ion concentrations. Rather than attempting to estimate γ_{H^+} as a function of solution conditions, it is possible to account for ionic strength effects directly in each ligand titration experiment. This is done by first performing a separate pH electrode calibration titration with strong acid/strong base titrants of exactly known c_{H^+} and c_{OH^-} , respectively [4, 10]. These titrants contain sufficient supporting electrolyte to set their ionic strength equal to the desired value for the subsequent ligand titration. This calibration titration spans a broad range of hydrogen-ion activities (corresponding to $2 < p^{\circ}H < 12$) and is usually performed by titrating acid with base (increasing $p^{\circ}H$). At each step in the titration, c_{H^+} is known from the acid and basic titrant concentrations and from the total volume and volume basic titrant added. Also known at each step is a_{H^+} from the electrode response, expressed in mV. Data from a calibration titration will therefore consist of a set of (c_{H^+}, a_{H^+}) points.

The operating relation for a pH electrode is a modified form of the Nernst equation [11]:

$$E = E_0' + \frac{2.303 \cdot R \cdot T}{F} \cdot \log(c_{H^+} \cdot \gamma_{H^+}) \quad (16)$$

where E is the measured hydrogen-ion activity (mV), E_0' is the electrode reference potential (mV), R is the gas constant ($8.314 \text{ J/mol} \cdot \text{K}$), T is the system temperature (K), and F is the Faraday constant (96.485 C/mol). This is equivalent to:

$$E = E_0' + S \cdot \log(c_{H^+} \cdot \gamma_{H^+}) \quad (17)$$

where $S = 2.303 \cdot R \cdot T / F$.

Here it may be asserted that γ_{H^+} is essentially a constant over the range of the calibration titration since the titrants both have the same concentrations of supporting electrolyte. However it is expected that γ_{H^+} will be different in different ionic-strength calibration titrations, due to the effect of electrolyte concentration on a_{H^+} [10]. Rearranging equation (17) allows us to relate hydrogen-ion activity (E) to hydrogen-ion concentration:

$$E = E_0 + S \cdot \log(c_{H^+}) = E_0 - S \cdot (p^{\circ}H) \quad (18)$$

where $E_0 = E_0' + S \cdot \log(\gamma_{H^+})$ and is constant over a given calibration titration. From equation (18) it is apparent that E_0 and S must be known to convert measured hydrogen-ion activity (E) into concentration ($p^{\circ}H$). A Marquardt least-squares algorithm was used to fit E_0 and S from calibration titration data [10]; details of this procedure are in Appendix B.

This calibration-titration process is performed prior to each ligand titration. The regressed electrode parameters E_0 and S are used to convert measured hydrogen-ion activities to concentrations in the ligand titrations, to support the theoretical net-charge calculations described above.

3 Experimental Methods

In each titration experiment, several steps occur before the titration of the ligand of interest (in our case amino acid or protein). Figure 2a illustrates the order of a typical titration experiment. Step (1)* is a titration of stock solutions of HCl and KOH (our acidic and basic titrants in all experiments) against potassium hydrogen phthalate and phenolphthalein to determine their exact molar concentrations of H^+ and OH^- , respectively. This is important information for step (2), a titration of HCl with KOH from pH 2-12, for the purpose of performing a broad range calibration of the pH electrode. Since many hydrogen-ion equilibria in amino acids and proteins have extreme pK_a 's (e.g. aspartic acid carboxyl $pK_a \sim 2$, and arginine ϵ -imino $pK_a \sim 12$), careful calibrations are vital to measure these pK_a 's accurately. With a properly calibrated electrode, it is then possible to gather reliable titration data for amino acids and proteins over the same pH range as the calibration (step (3)). Steps (2) and (3) are automated with software described below.

Regression of amino-acid p^0K_a 's from experimental titration data is carried out by sophisticated software (program BETA90, described in §3.2.5 below) developed in the College of Chemistry at the University of California at Berkeley [10]. Once the relevant amino-acid p^0K_a 's have been determined at a certain solution ionic strength, it is possible to calculate via IST the p^0H dependence of protein net charge at the same ionic strength. Comparison of calculated protein titration curves with experimental data is discussed in §4.

3.1 Titration Apparatus

The apparatus used to conduct the experiments is shown in Figure 1. Hydrogen-ion titrations were carried out in a 500-mL jacketed reaction vessel (5). Prepared solutions of ligand (either amino acid or protein) were loaded into the vessel and agitated by a magnetic stirrer. Small aliquots of titrant, either acid or base, were injected through a fritted line (6) by a Dosimat 665 (2) (Metrohm, Westbury, NY) precision pump (see Appendix C for details on pump configuration and operation). pH was measured by a Ross combination

* Step (1) need be performed only when new reservoir volumes (~1.5 L) of titrant are prepared; a typical reservoir is sufficient for approximately 15 ligand titrations. Since the reservoirs are blanketed with argon at a small positive pressure, carbon dioxide dissolution into the titrant is considered to be negligible.

semimicro electrode (7) (Orion, Boston, MA) [12] and a Chemtrix 60A pH meter (8). The vessel was sealed with a specially-designed plexiglas lid, allowing experiments to be conducted in a CO₂-free environment. The contents of the vessel were maintained at 25°C by recirculation of water from a temperature-controlled bath (9) through the jacket.

An IBM PC-286 (1) was used both to control the Dosimat titrant injection pump and to acquire pH and volumetric data for each titration. With an analog-to-digital signal conversion board and the appropriate software, it was possible to specify aliquot volume and injection frequency and to control sampling of the mV output from the pH meter. For each injection of titrant, the aliquot volume and the resulting pH were recorded in a data file. Results from a typical titration experiment consisted of a data file containing approximately 200 of these (volume,pH) points.

The presence of CO₂ in the headspace of the titration vessel and the titrant reservoir bottles would cause a change in pH due to CO₂ partitioning into the solution [10]. In order to eliminate this effect, the titration vessel and reservoirs were blanketed with argon at a small positive pressure. Argon from a high pressure cylinder passed through a column (13) containing Drierite and acid-gas-specific molecular sieves to remove any residual water and CO₂, before passing to a distribution manifold (14) connected to the titration vessel and the titrant reservoir bottles (3).

A temperature bath equipped with a thermostatted heater/mixer unit (10) was used to circulate water at 25°C through the titration vessel jacket. The bath was also equipped with a coil for house cold water (11) to provide cooling in the case of ambient temperatures greater than 25°C.

The Plexiglas lid for the titration vessel was designed and constructed in-house. It is shown in Appendix D. The holes in the lid were specially-machined to fit the various dimensions of the titrant injection line, the pH probe and the argon lines. During a titration, the lid rested on an O-ring mounted in a groove on the top lip of the titration vessel. Since the operating pressure was only slightly greater than atmospheric, the lid was secured to the titration vessel with Parafilm. However, in the future a compression clamp would be preferable.

Appendix E describes some important points that must be considered to ensure that the apparatus is properly configured for gathering titration data with minimal experimental error.




Figure 1: Titration Hardware

- | | |
|--|--------------------------------------|
| (1) IBM PC-286 titration controller | (8) Chemtrix 60A pH meter |
| (2) Dosimat 665 injection pump | (9) Water bath |
| (3) Titrant reservoir bottle | (10) Thermostatted heater/mixer |
| (4) Injection syringe | (11) Cooling water coil |
| (5) Jacketed reaction vessel | (12) Argon supply |
| (6) Fritted titrant injection line | (13) Drierite/molecular sieve column |
| (7) Ross combination semimicro electrode | (14) Gas manifold |

3.2 Software Description

The main function of the control software was twofold: to provide automation of the titration equipment; and to gather and record experimental data. Figure 2b illustrates the interrelation and succession of all the titration programs and provides a schematic representation of the order and flow of information in each titration experiment. The following sections describe the function of each program and how the output of each program is related to the subsequent program step.

3.2.1 CAL-TITR

This program controls the HCl/KOH titration necessary to calibrate the pH electrode before every ligand titration. To perform these calibrations, HCl was titrated with KOH; initial pH was always near 2, and final pH was always near 12. Results of the calibration titration were stored in a data file named CALIB.DAT, in the form of volume KOH added vs. electrode mV response. The required inputs are the following:

1) Total volume of titrant to be added. The program which calculates the electrode calibration constants (CAL-ANAL) requires 30 volume-added vs. mV-response data points from the HCl/KOH titration. Based on this number of points, the pH range over which it is desired to calibrate the electrode (2-12), the concentration of HCl and KOH, and the volume of the titration vessel, it was possible to calculate the total volume of titrant (KOH) required to perform a calibration titration. In all cases this volume was 21 mL since, in all ligand titrations, the parameters mentioned above were held constant. From this it was possible to calculate the titrant aliquot volume per injection to be 0.7 mL

$$(21\text{mL}/30\text{pts} = 0.7\text{ mL/pt}).$$

2) The delay time between titrant injections. This quantity was chosen to allow enough time for all hydrogen-ion reactions to equilibrate after injection of titrant aliquot. A 40-second delay was assumed to be sufficient, for the amino acid and protein ligand titrations, in all cases, as well as for the HCl/KOH calibration titration. This was justified by the observation that mV response after each injection of titrant stopped changing well before 40 seconds had elapsed.




Figure 2a: Experimental Titration outline




Figure 2b: Titration software flow sheet (* = user-defined file name)

It is important to realize that in the calibration titration, the measured quantity is electrode mV response to a change in H^+ concentration, not pH. Before mV readings could be converted to p^oH , the electrode slope and intercept were calculated by the program CAL-ANAL, using the calibration titration results stored by CAL-TITR in the file CALIB.DAT. Further details on the operation of CAL-TITR are provided in Appendix F.

3.2.2 CAL-ANAL

This program calculates the pH electrode slope (mV/ p^oH unit) and intercept (E_0) based on the results of the HCl/KOH calibration titration [10]. It performs a non-linear least squares fit of the calibration titration data to the simplified Nernst equation given in §2. Required inputs for CAL-ANAL are the following:

- 1) The volume titrant added vs. mV response data from the HCl/KOH titration. These are read into the program from the file CALIB.DAT created by CAL-TITR.
- 2) The exact molar concentrations of HCl/KOH, given by the potassium hydrogen phthalate/phenolphthalein titration.
- 3) The initial volume and concentration of HCl added to the vessel before titration with KOH.

With this information CAL-ANAL converts the volume added data to p^oH and then performs the p^oH/mV fit to obtain electrode slope and intercept. These are stored in a file named CALIB.PAR, which is needed by the ligand-titration control program, TITRATE.

3.2.3 TITRATE

This program controls the titration of ligand (amino acid or protein) with either strong acid or strong base titrant. TITRATE is similar to CAL-TITR described above (see Appendix F). The required inputs are the following:

- 1) Whether the titration is acidic (decreasing pH) or basic (increasing pH).
- 2) Endpoint pH of the titration.

- 3) Aliquot volume for titrant injection. This was chosen small enough to capture the important features of the titration curve, in hopes of optimizing the accuracy of the $p^{\circ}K_a$ fitting. Typical ΔV 's were on the order of 0.4 mL, which yielded between 200 and 300 data points for the whole titration, depending on the ligand.
- 4) Data file CALIB.PAR. This file contains the slope and intercept for the pH electrode, calculated by CAL-ANAL. These parameters are necessary to convert the mV response to addition of a titrant aliquot to $p^{\circ}H$ change for each data point in a ligand titration. This file is automatically supplied by CAL-ANAL following the electrode calibration titration described above.
- 5) Delay time between injections of aliquots. As previously described a delay time of 40 seconds was found to be sufficient for hydrogen-ion reactions to equilibrate in all our ligand titrations.
- 6) Output data file name *.DAT. Results of the ligand titration are written to a data file with this name (where * is specified by the operator). In this case, the form of saved data for each point is volume titrant added vs. $p^{\circ}H$. This data file contains all the experimental information necessary to calculate amino acid $p^{\circ}K_a$'s.
- 7) Operating mode of Dosimat pump.

3.2.4 SETUPBTA

This program creates an input file for the $p^{\circ}K_a$ fitting program (BETA90) [10]; the input file name is *.INP. This file contains all the experimental ligand-titration data points in *.DAT, plus further information about the titration conditions and ligand chemistry required for BETA90 to regress ligand $p^{\circ}K_a$'s and protonation states. For further details on the use of SETUPBTA and the contents and format of *.INP, see Appendix G.1.

3.2.5 BETA90

This program regresses amino acid $p^{\circ}K_a$'s from experimental titration data [10]. Based on the data contained in *.INP and initial approximations of the ligand $p^{\circ}K_a$'s, this program iteratively solves for the experimental $p^{\circ}K_a$'s. It fits the ligand $p^{\circ}K_a$'s by minimizing differences between the experimental titration curve and the curve calculated with current values of the $p^{\circ}K_a$'s. It also calculates $p^{\circ}I$ of the ligand, based on titration data

and information contained in *.INP. When the best set of $p^{\circ}K_a$'s is found, the program terminates and writes the fit values of $p^{\circ}K_a$'s, the calculated charge vs. $p^{\circ}H$ curve calculated from the fit $p^{\circ}K_a$'s and other summary information (quality of fit, # iterations, etc.) to the output files PLOT.O&C and *.SUM.

It should be noted that BETA90 in its current form has the ability to treat ligands having 20 or fewer hydrogen ion equilibria [10]. It handles amino acid titrations (3 $p^{\circ}K_a$'s) easily, but cannot be applied to our protein titration data (lysozyme = 32 $p^{\circ}K_a$'s; chymotrypsin = 43 $p^{\circ}K_a$'s). Also BETA90 has the ability to fit $p^{\circ}K_a$'s outside the $p^{\circ}H$ range of the experimental titration data, given a reasonable initial estimate. Initial estimates were taken from tabulated values for low ionic strength conditions, found in the literature [7]. Appendix G.2 gives further detail on the operation of BETA90.

3.3 Experimental Procedure

Eleven amino acids having thirteen relevant hydrogen-ion equilibria were titrated to determine their $p^{\circ}K_a$'s. Also, lysozyme and α -chymotrypsin were titrated to determine their relative charges as a function of $p^{\circ}H$ and for comparison to IST calculations. The titrations were conducted in KCl solutions at three different ionic strengths: 0.1, 1.0 and 3.0 M. Instructions for every step in the experimental titration procedure are found in Appendix H.

To check the amino-acid and protein-titration results for reproducibility, every titration experiment was performed at least twice under the same experimental conditions. If significant differences in amino acid $p^{\circ}K_a$'s (0.05 pH units) were observed over two titrations, a third titration was performed.

4 Results/Discussion

Amino-Acid Hydrogen-Ion Titration

Hydrogen-ion titrations were performed over the $p^{\circ}H$ range 2 - 11 for the 11 amino acids listed in § 2. Of the 29 hydrogen-ion equilibria probed by conducting titrations, 13 were relevant for protein net-charge calculations using the independent-site theory described previously. The $p^{\circ}K_a$'s for these equilibria are listed in Table 1; they are discussed below. The remaining 16 $p^{\circ}K_a$'s are tabulated in Appendix I. As indicated in § 2, all $p^{\circ}K_a$'s are expressed in term of concentrations rather than activities.

Figure 3 shows $p^{\circ}H$ vs. net charge for a titration of aspartic acid in 3.0 M KCl. The open circles represent experimental titration points (~200) converted from volumetric data to net charge, while the solid line is the fitted curve calculated by BETA90. Generally, the difference between the two curves is much less than 0.05 pH units, as indicated by the dashed line. Near the water-equivalence point, where data are sparse because the $p^{\circ}H$ changes rapidly with addition of small titrant aliquots, the error is slightly greater. The $p^{\circ}K_a$'s correspond to points where the $p^{\circ}H$ vs. charge curve has an inflection point of minimum slope [7]. Note from Figure 3 that BETA90 has the ability to determine $p^{\circ}K_a$'s which may not be readily observable by examining the titration curve, as in the case of the γ -carboxyl substituent of aspartic acid ($p^{\circ}K_a = 3.95$). (Refer to § 3.2.5 above and Appendix G.2 for further details on the BETA90 curve-fitting algorithm and $p^{\circ}K_a$ refinement)

Table 1 shows the experimental $p^{\circ}K_a$ values of interest in IST net-charge calculations for lysozyme and α -chymotrypsin. Results are shown for three different values of KCl ionic strength: 0.1 M, 1.0 M and 3.0 M. At each ionic strength each amino acid was titrated at least twice for reproducibility. A third titration was performed in cases where the first two $p^{\circ}K_a$'s measured were different by more than 0.05 pH units. In Table 1, all measured $p^{\circ}K_a$'s are shown for each amino acid, for completeness. The boldfaced numbers are the values of $p^{\circ}K_a$ chosen for use in the theoretical calculations of protein net charge. For titrations where the first two $p^{\circ}K_a$'s measured differ by less than 0.05 pH units, the chosen value was the average of the two. However, in the case where a third titration was necessary, an appropriate average was taken, depending on the resulting $p^{\circ}K_a$'s and on an estimation of experimental error (quality of fit).

AMINO ACID (titrated moiety)	IS = 0.1M			IS = 1.0M		IS = 3.0M	
	Experiment		Lit.				
Alanine (amino)*	9.85			9.74		10.06	
	9.88	9.86	9.69	9.78	9.76	10.13	10.1
Arginine (substituent)	13.42			12.76		12.82	
	12.68	12.7		12.44	12.9	12.85	12.83
	12.87			13.40			
Asparagine (carboxyl)	2.18			2.35		2.45	
	2.19	2.18	2.14	2.28	2.3	2.47	2.46
				2.27			
Aspartic Acid (substituent)*	3.73			3.65		3.95	
	3.72	3.72	3.70	3.70	3.67	3.94	3.95
Cysteine (amino)*	10.25			10.22		10.37	
	10.44	10.42	10.29	10.29	10.26	10.28	10.4
	10.40					10.51	
Glutamic Acid (substituent)*	4.18			4.06		4.44	
	4.16	4.17	4.20	4.15	4.14	4.43	4.43
				4.14			
Histidine (substituent)	6.06			6.28		6.62	
	6.06	6.06	6.02	6.23	6.26	6.58	6.60
Isoleucine (amino)	9.63			9.76		10.06	
	9.65	9.64	9.62	9.68	9.71	10.03	10.04
				9.73			
Leucine (carboxyl)	2.33			2.44		2.70	
	2.32	2.33	2.35	2.43	2.43	2.70	2.70
Lysine (amino)	9.22			9.37		9.71	
	9.21	9.21	9.08	9.40	9.38	9.74	9.72
(substituent)	10.83			10.85		11.04	
	10.79	10.81	10.69	10.84	10.85	10.95	11.04
						11.05	
Tyrosine (carboxyl)	2.34			2.29		2.59	
	2.10	2.3	2.17	2.38	2.3	2.64	2.62
(substituent)*	10.18			10.02		10.21	
	10.15	10.17	10.14	10.02	10.02	10.15	10.18

Table 1: Summary of amino-acid titration results, showing experimental $p^{\circ}K_a$'s as a function of KCl concentration. Literature values for IS=0.1 M are shown for comparison [4]. The $p^{\circ}K_a$'s shown here are the values necessary for IST calculations of protein net charge. More complete data are found in Appendix I. The amino-acid concentrations in the hydrogen-ion titrations are 10 mM, except Tyrosine 2 mM and Alanine, Isoleucine and Leucine 20 mM

In most titrations, reproducibility of $p^{\circ}K_a$'s was high. However, two special cases are worthy of distinction. For the arginine ϵ -imino group, the $p^{\circ}K_a$ is near 13. This was the most basic $p^{\circ}K_a$ probed in all our measurements and was well outside the $p^{\circ}H$ range of our titrations. Therefore the uncertainty in BETA90 $p^{\circ}K_a$ fitting for this moiety was much higher than for equilibria with less extreme $p^{\circ}K_a$'s, especially in the 0.1 M and 1.0 M cases. Also, it proved very difficult for BETA90 to perform $p^{\circ}K_a$ fitting for tyrosine in 0.1 M KCl. Of the 5 tyrosine titrations conducted at this ionic strength, only two yielded values for $p^{\circ}K_a$'s of the carboxyl-endgroup and phenolic-substituent equilibria. Consequently, great variation was seen in the carboxyl $p^{\circ}K_a$; the value 2.3 was chosen since it had a lower G.O.F.. This abnormal behavior may be a consequence of the very low solubility of tyrosine (0.5 mM).

Table 1 also shows comparison of $p^{\circ}K_a$'s measured in 0.1 M KCl at 25°C to values found in the literature [4]. In many cases the measured $p^{\circ}K_a$'s agree with the literature values to within a few hundredths of a pH unit, while the remaining $p^{\circ}K_a$'s differ by at most 0.15 pH units. Based on this substantial agreement between our titration and fitting results and the literature values, we are confident that the $p^{\circ}K_a$ values shown in Table 1 for amino acids in 1.0 M and 3.0 M KCl solutions are accurate, as well.

We observed that several amino acid hydrogen ion equilibria showed an initial decrease in $p^{\circ}K_a$'s as ionic strength rises from 0.1 M to 1.0 M, followed by an increase as the ionic strength rises to 3.0 M. The $p^{\circ}K_a$ at 3.0 M is generally higher than the value at 0.1 M. It is important to note that not all equilibria exhibit this effect, many demonstrate a simple steady increase in $p^{\circ}K_a$ with increasing ionic strength. There seems to be no generalization in the chemistry of the amino acids with which to correlate this behavior. The equilibria exhibiting this effect are denoted with an asterisk on Table 1.

Martell and Smith noted similar trends for certain amino-acid hydrogen-ion equilibria in their review of the amino-acid titration literature [4]. They also found that, for many simple amino acids, certain $p^{\circ}K_a$'s decrease with increasing ionic strength and reach a minimum near $IS=0.5$ M, followed by a steady increase as ionic strength goes to 3.0 M. Literature data is extremely rare for $p^{\circ}K_a$'s in solutions of $IS > 3.0$ M [4]. They show that, in these cases, $p^{\circ}K_a$ values are often nearly equal at $IS=0.1$ M and at $IS=1.0$ M. It may be possible to explain this behavior by considering the effect of electrolyte concentration on the activity of the hydrogen ion. Recall that in § 2, we asserted that γ_{H^+} was constant over a given titration. However, γ_{H^+} is a function of ionic strength and differs for titrations of

different IS. If we define our hydrogen-ion equilibrium constant, K_a , for a specific moiety in a given titration as [9]:

$$K_a = \frac{c_{A^-} \cdot c_{H^+}}{c_{HA}} \cdot \frac{\gamma_{A^-} \cdot \gamma_{H^+}}{\gamma_{HA}} = K_a^o \cdot \frac{\gamma_{A^-} \cdot \gamma_{H^+}}{\gamma_{HA}} \quad (19)$$

and make the assumption that ligand is very dilute with respect to electrolyte (typical ligand concentration were ~0.5 - 5 mM), then it is appropriate to assume ligand activity coefficients are equal to unity, and thus equation (19) reduces to:

$$K_a = K_a^o \cdot \gamma_{H^+} \quad (20)$$

Since K_a is based on activity and therefore a function of temperature only, we see that changes in K_a^o with changing electrolyte concentration must be counterbalanced by appropriate changes in γ_{H^+} [7]. Under these assumptions, it is possible to quantify approximate changes in γ_{H^+} as ionic strength increases, in the following manner:

$$\Delta p^o K_a = (p^o K_a)_{IS2} - (p^o K_a)_{IS1} = -\log \left[\frac{(\gamma_{H^+})_{IS1}}{(\gamma_{H^+})_{IS2}} \right] \quad (21)$$

From the data in Table 1 we estimate that γ_{H^+} decreases by roughly 20% as ionic strength goes from 0.1 M to 1.0 M and increases by roughly 50% as ionic strength increases from 1.0 M to 3.0 M depending on the equilibrium examined. In general, these activity arguments should be valid for all hydrogen-ion equilibria. However, as previously mentioned, it is not presently known why some equilibria show this trend in $p^o K_a$, while others demonstrate a steady increase with ionic strength instead.

Independent-Site Theory

After we determined the relevant $p^o K_a$'s shown in Table 1, it was possible to calculate the protein net charge as a function of $p^o H$ using the independent-site theory. Figures 4a and 4b show these calculated hydrogen-ion titration curves of α -chymotrypsin and lysozyme for three different solutions of KCl at 0.1 M, 1.0 M and 3.0 M ionic strength.

In both diagrams the titration curves for 0.1 M and 1.0 M ionic strength overlap almost completely, because the amino-acid $p^{\circ}K_a$'s are relatively the same for both ionic strength solutions. Note that the $p^{\circ}K_a$'s determine the shape of the curves predicted by the IST as discussed in § 2. At the pH extremes, the protein is either totally protonated or unprotonated because the calculated titration curve covers the full range of $p^{\circ}K_a$'s. The three different curves of Figures 4a and 4b meet at the endpoints where the protein has either its maximum (low $p^{\circ}H$) or minimum (high $p^{\circ}H$) charge.

The 3.0M-ionic-strength titration curves differ from those curves of the lower ionic strengths in the regions of $p^{\circ}H$ 2 to 5 and 9 to 11.5, where most of the $p^{\circ}K_a$ values are found. The accumulation of $p^{\circ}K_a$'s in these two regions will magnify the net charge difference. Since the $p^{\circ}K_a$'s are the highest in the 3.0 M ionic-strength salt solution, the respective curve for IS=3.0 M will be shifted to the right of the lower ionic strength curves.

A special point of interest in every protein titration is the isoelectric point, pI , defined as the pH where the net charge is equal to zero. With the IST calculations the $p^{\circ}I$ of lysozyme and α -chymotrypsin can be determined even for high ionic strength, where conventional isoelectric focusing fails. The results are collected in Table 2. Since the IST calculates the $p^{\circ}I$ from free amino-acid $p^{\circ}K_a$'s, the $p^{\circ}I$ values are the same for IS=0.1 M and 1.0 M, whereas the $p^{\circ}I$ is higher for the 3.0M ionic strength.

	IS=0.1M	IS=1.0M	IS=3.0M
Lysozyme	11.131	11.159	11.311
α -Chymotrypsin	9.676	9.625	9.831


Table 2: Calculated $p^{\circ}I$'s for lysozyme and α -chymotrypsin using IST

Comparison to Protein Hydrogen-Ion Titrations

Protein hydrogen-ion titrations yield data points consisting of the volume of titrant added and the corresponding $p^{\circ}H$ response. In order to calculate protein net charge as a

function of $p^{\circ}H$ for the purpose of comparing theory to experiment, it was necessary to convert volumetric data to charge data. This was done simply as follows.

In a basic protein titration every aliquot of KOH titrant added contains a discrete number of OH^- -ions, n_{OH} . Some of these injected OH^- -ions, $n_{\Delta pH}$, will raise the $p^{\circ}H$ of the solution by the observed value $\Delta p^{\circ}H$. The remaining number of OH^- -ions, $n_{\Delta Z}$, will react with the titratable groups of the proteins to change their net charge by ΔZ , as illustrated below:



$$n_{OH} = n_{\Delta Z} + n_{\Delta pH} \quad (22)$$

For a known amount of protein in the solution, n_{prot} , the change in net charge of one protein is then:

$$\Delta Z = \frac{n_{\Delta Z}}{n_{prot}} = \frac{n_{OH\text{added}} - n_{\Delta pH}}{n_{prot}} \quad (23)$$

The change of the net charge for a protein molecule titrated from starting point '0' ($p^{\circ}H_0; mL_0$) to any point 'i' ($p^{\circ}H_i; mL_i$) can be calculated as follows.

The amount of OH^- -ions added in one or more aliquots of basic titrant (KOH) till point 'i' is reached:

$$n_{OH\text{added}} = mL_{\text{added}} \cdot c_{KOH} \quad (24)$$

The amount of OH^- -ions needed to raise the solution $p^{\circ}H$ from $p^{\circ}H_0$ to $p^{\circ}H_i$:

$$n_{\Delta pH} = n_{OH_i} - n_{OH_0} = V_i \cdot c_{OH_i} - V_0 \cdot c_{OH_0} \quad (25)$$

with

$$K^{\circ}_w = [H^+] \cdot [OH^-] \text{ and } V_i = V_0 + mL_{\text{added}}$$

follows that

$$n_{\Delta pH} = (V_0 + mL_{\text{added}}) \frac{K^{\circ}_w}{10^{-p^{\circ}H_i}} - V_0 \frac{K^{\circ}_w}{10^{-p^{\circ}H_0}} \quad (26)$$

and therefore the change of the net charge for one titrated protein molecule is:

$$\Delta z = \frac{\Delta mL_{added} \cdot c_{KOH} - \left[(V_0 + mL_{added}) \frac{K_w^o}{10^{-p^oH_i}} - V_0 \frac{K_w^o}{10^{-p^oH_0}} \right]}{n_{prot}} \quad (27)$$

Using this equation, only the net charge differences Δz , relative to the starting point of the titration, could be calculated; in a basic titration Δz is always negative. These $\Delta z/p^oH$ titration data can only be converted into a absolute net-charge/ p^oH titration curve if one absolute net charge is known for a certain p^oH . In low ionic-strength solutions, it is possible to measure p^oI by isoelectric-focusing electrophoresis. However, this is not possible in high ionic-strength solutions. Therefore, for the purpose of comparing experimental protein titration curves with our IST, the experimental p^oI 's of the proteins were considered as a first approximation to be equal to the ones calculated from the IST. The calculated p^oI 's for lysozyme and α -chymotrypsin are shown in Table 2.

With this knowledge, equation (14) can now be extended by an intercept t , allowing to calculate absolute net charge z_{net} as a function of p^oH :

$$z_{net} = \frac{\Delta mL_{added} \cdot c_{KOH} - \left[(V_0 + mL_{added}) \frac{K_w^o}{10^{-p^oH_i}} - V_0 \frac{K_w^o}{10^{-p^oH_0}} \right]}{n_{prot}} + t \quad (28)$$

The intercept t is chosen so that the absolute charge z_{net} is zero when p^oH equals the p^oI calculated via IST.

The experimental data of a protein hydrogen-ion titration can be compared to the calculated titration curve based on the IST by plotting these two curves in the same diagram fitted together at the calculated p^oI . Figures 5a and 5b show experimental and calculated titration results for lysozyme in KCl solutions of IS=0.1 M and IS=1.0 M. The IST results fit the experimental data well in the lower ionic strength (0.1 M). The IST seems to be more sensitive to p^oH changes than the experimental protein titration, which is depicted by the sharper bends in the IST graph. Furthermore literature data [13] of IS=0.15 added to the graphs show that the general tendency of the IST is correct. In this low ionic strength, the calculation based on free amino-acid p^oK_a 's can be seen as a success for the simplicity and

large assumptions of our model. In IS=1.0 M the calculated and measured curve still conform in the region between p^oH 4 and 11.3. The deviation in this section can be explained by the p^oI being fixed in the steep part of the curve so that small experimental inaccuracies would cause a discernible vertical shift. The acidic and basic ends of the measured protein titration curve, however, differ by larger amounts. This difference is even higher in IS=3.0, making a graphic illustration for this ionic strength useless.

Figures 6a and 6b show measured and calculated titration curves of chymotrypsin in two higher ionic strengths (1.0M and 3.0M). Both diagrams display similar characteristics. Again the IST fulfills the general titration-curve tendencies especially in the p^oH range from 4 to 11. For this protein the deviation of the IST seems to be independent of the increasing ionic strength. There is a drift in the acidic region of the titration curve similar to that of lysozyme in IS=1.0 M KCl. The experimental titration is more likely to be at fault. For example at a p^oH below 3.5 for chymotrypsin an absolute net charge of 30 is observed. This is impossible since chymotrypsin contains only 22 titratable groups in its amino-acid sequence which can be positively charged. The higher value of the net charge stretches the protein-titration curve leading to a deviation that makes a comparison to the calculated curve impossible. These errors can be explained by the special protein peculiarities, e.g. the folded structure [6]. Taking this into account goes beyond the scope of this work. Thus, in high ionic strength an assessment of the IST fails because of the incorrect protein titration curves.

Figure 7 illustrates the reproducibility of two α -chymotrypsin experiments in the same 1.0M ionic strength. Because we know the plots are wrong, yet reproducible, proteins can not be titrated with the simple procedure used for the amino-acid titrations. Therefore further development of protein titration methods using the available equipment is necessary.

The IST has been shown to be quantitative in low ionic strength, but for high ionic strength, without an accurate experimental titration curve, it is hard to determine the IST's effectiveness.

Comparison of Basic and Acidic Protein Titration

Even if the experimental titration curves show consistent error, it is still possible to observe basic trends of the hydrogen-ion titration behavior. This is reasonable, because in most of the protein titrations performed, the reproducibility was high as mentioned above.

Figure 8 shows a basic (increasing $p^{\circ}H$) and subsequent acidic (decreasing $p^{\circ}H$) titration of an α -chymotrypsin protein solution in 1.0 M ionic strength. For illustration, the basic titration curve is again fixed on the calculated $p^{\circ}I$ of the IST. The starting point of the acidic titration curve was fitted to the end point of the previously conducted basic titration because the net charge of the protein does not change between the end of the basic titration and the beginning of the acidic titration. We observed that the basic and acidic titration curves differ for most of the titration range, but converge again in the acidic extreme of the titration when all titratable groups are protonated and the maximum net charge is reached.

The deviation between the curves might be due to hysteresis effects, caused by reversible changes in the protein molecule configurations [6]. It is also possible that irreversible configuration changes happen in the basic extreme of the titration.

Aspartic Acid Titration Ionic Strength 3.0 M KCl




Figure 3: Obtained and calculated titration curve for aspartic acid in 3.0 M KCl

Calculated IST Titration Curves





Figure 4: Calculated IST hydrogen-ion titration curves in ionic strength 0.1M, 1.0M and 3.0 M KCl for a) a-chymotrypsin and b) lysozyme

Lysozyme Hydrogen-Ion Titration

a) Ionic Strength = 0.1 M KCl



b) Ionic Strength = 1.0 M KCl





Figure 5: Lysozyme hydrogen-ion titration comparison in a) 0.1 M IS and b) 1.0 M IS

α -Chymotrypsin Hydrogen-Ion Titration

a) Ionic Strength = 1.0 M KCl



b) Ionic Strength = 3.0 M KCl




Figure 6: α -Chymotrypsin hydrogen-ion titration comparison in
a) 1.0 M IS and b) 3.0 M IS

**Reproducibility of
Chymotrypsin Hydrogen-Ion Titrations**
Ionic Strength = 1.0M KCl




Figure 7: Reproducibility of chymotrypsin hydrogen-ion titration in IS=1.0 M KCl

Basic and Acidic Chymotrypsin Titration

Ionic Strength = 1.0M KCl




Figure 8: Basic and acidic titration of chymotrypsin in IS=1.0 M KCl

5 Conclusion

In this work, an apparatus was constructed and a method developed to perform high precision titrations that measure the concentration of the hydrogen ions instead of their activity. The ligands were titrated efficiently by computerizing the experiment.

Titration results were obtained for 11 different amino acids dissolved in KCl at three separate ionic strengths (0.1M, 1.0M and 3.0M). Table 1 shows the $p^{\circ}K_a$ values of the 13 titratable groups which are subsequently used for net charge vs. pH calculation for lysozyme and α -chymotrypsin by the independent-site theory. In the majority of the titrations, the $p^{\circ}K_a$'s were reproducible to within 0.05 for each experiment. For the 0.1 M KCl experiments, most of the $p^{\circ}K_a$'s differ from literature values by less than a few hundredths of a pH unit. Based on these arguments, we are confident that the experimental $p^{\circ}K_a$'s determined for the amino acid in IS = 1.0M and IS = 3.0M KCl are also accurate. The literature $p^{\circ}K_a$'s show a uniform trend with ionic strength that is matched by the experimentally determined $p^{\circ}K_a$'s. The values decrease initially with increasing ionic strength and reach a minimum near IS = 0.5M, followed by a steady increase as ionic strength goes to 3.0M. Because of this dependence, the $p^{\circ}K_a$ values are often nearly equal at IS = 0.1M and IS = 1.0M. Appendix I contains all the measured $p^{\circ}K_a$'s of the eleven amino acids, including the values for titratable groups that are not relevant for the protein net-charge calculations. These data can be used for further determination of the amino acid's $p^{\circ}K_a$ dependence on ionic strength.

The net charges of lysozyme and α -chymotrypsin were calculated as a function of pH for each ionic strength using the experimental amino-acid $p^{\circ}K_a$'s and the IST. The IST assumes that all of the titratable groups of the protein are independent of each other and that their $p^{\circ}K_a$'s are equal to those in free amino-acid solution. To compare this theory with experimental results, we also performed hydrogen-ion titrations with lysozyme and α -chymotrypsin at all three ionic strengths. For the low ionic-strength case (0.1 M), the protein-titration data and the results obtained from the IST calculation agree reasonably well considering the large IST assumptions. Although the experimental and calculated titration curves show the same trends in high ionic strength (1.0 M, 3.0M), an absolute comparison fails because of the uncertainties in the protein-titration data.

The IST is a first approximation for the calculation of protein net charge at high ionic strength. Improvements in the experimental titrations of proteins are required to assess

accurately the validity of the IST. The amino-acid $p^{\circ}K_a$ obtained here give a basis for further development of this theory.

A Amino-Acid Sequences of Lysozyme and α -Chymotrypsin

The amino acid sequences of lysozyme and α -chymotrypsin are shown in Table 3 and Figure 9 respectively. These illustrations provide the type and number of all titratable groups contributing to the net charge calculation for each protein. These data are summarized in Table 4.

	Amino acid		Amino acid		Amino acid
1	Lys	45	Arg	90	Ala
2	Val	46	Asn	91	Ser
3	Phe	47	Thr	92	Val
4	Gly	48	Asp	93	Asn
5	Arg	49	Gly	94	Cys
6	Cys	50	Ser	95	Ala
7	Glu	51	Thr	96	Lys
8	Leu	52	Asp	97	Lys
9	Ala	53	Tyr	98	Ile
10	Ala	54	Gly	99	Val
11	Ala	55	Ile	100	Ser
12	Met	56	Leu	101	Asp
13	Lys	57	Gln	102	Gly
14	Arg	58	Ile	103	Asp
15	His	59	Asn	104	Gly
16	Gly	60	Ser	105	Met
17	Leu	61	Arg	106	Asn
18	Asp	62	Trp	107	Ala
19	Asn	63	Trp	108	Trp
20	Tyr	64	Cys	109	Val
21	Arg	65	Asn	110	Ala
22	Gly	66	Asp	111	Trp
23	Tyr	67	Gly	112	Arg
24	Ser	68	Arg	113	Asn
25	Leu	69	Thr	114	Arg
26	Gly	70	Pro	115	Cys
27	Asn	71	Gly	116	Lys
28	Trp	72	Ser	117	Gly
29	Val	73	Arg	118	Thr
30	Cys	74	Asn	119	Asp
31	Ala	75	Leu	120	Val
32	Ala	76	Cys	121	Gln
33	Lys	77	Asn	122	Ala
34	Phe	78	Ile	123	Trp
35	Glu	79	Pro	124	Ile
36	Ser	80	Cys	125	Arg
37	Asn	81	Ser	126	Gly
38	Phe	82	Ala	127	Cys
39	Asn	83	Leu	128	Arg
40	Thr	84	Leu	129	Leu
41	Gln	85	Ser		
42	Ala	86	Ser		
43	Thr	87	Asp		
44	Asn	88	Ile		
		89	Thr		

Table 3: Amino-acid sequence of lysozyme [14]




Figure 9: Amino-acid sequence of chymotrypsin [15]

	Amino Acid	Chymotrypsin	Lysozyme
Substituent	Arginine	3	11
	Aspartic Acid	9	7
	Glutamic Acid	5	2
	Histidine	2	1
	Lysine	14	6
	Tyrosine	4	3
Carboxylic endgroup	Asparagine	1	
	Leucine	1	1
	Tyrosine	1	
Amino Endgroups	Alanine	1	
	Cysteine	1	
	Isoleucine	1	
	Lysine		1

Table 4: Titratable groups in lysozyme and chymotrypsin

B Electrode Calibration

As derived in §2, for each ligand titration a calibration titration and subsequent regression of the electrode intercept and slope were performed in order to convert measured hydrogen ion activity (E) into concentration ($p^{\circ}H$). The parameters E_0 and S are determined by fitting a line through the calibration-titration data points using a least squares Marquardt algorithm. Therefore these data must consist of a set of (c_{H^+} , a_{H^+}) points. The calibration titration yields only ($mL_{\text{base added}}$; a_{H^+}) points. The hydrogen-ion concentration $p^{\circ}H$ can be calculated from the total solution volume, the volume of acidic titrant added before the start of the titration, the volume of basic titrant added ($mL_{\text{base added}}$) and the exact concentrations of acidic and basic titrant.

The calculation follows a strong acid/strong base titration in which a given amount of added OH^- ions will consume a stoichiometric amount of H^+ ions present in solution [7]. Equation 29 and 30 describe the calculation for $p^{\circ}H$ below

$$p^{\circ}H = -\log\left(c_{\text{acid}} \cdot \frac{V_{\text{acid}}}{V_{\text{total}}} - c_{\text{base}} \cdot \frac{V_{\text{base}}}{V_{\text{total}}}\right) \quad (29)$$

and $p^{\circ}H$ above the equivalence point

$$p^{\circ}H = p^{\circ}K_W + \log\left(c_{\text{base}} \cdot \frac{V_{\text{base}}}{V_{\text{total}}} - c_{\text{acid}} \cdot \frac{V_{\text{acid}}}{V_{\text{total}}}\right) \quad (30)$$

From equation 30, it follows that a $p^{\circ}K_W$ is necessary to determine the $p^{\circ}H$ in the basic range. This $p^{\circ}K_W$ is based on a K°_W defined as:

$$K^{\circ}_W = \frac{c_{H^+} \cdot c_{OH^-}}{c_{H_2O}} = \frac{K_W}{\frac{\gamma_{H^+} \cdot \gamma_{OH^-}}{\gamma_{H_2O}}} \quad (31)$$

K°_W is a constant for a certain ionic strength and temperature, as the water-dissociation constant K_W depends only on the temperature and all activity coefficients are constant at a given high ionic strength of supporting electrolyte (see § 2).

Although K_w^0 is constant for a certain ionic strength its value differs from the known K_w due to activity effects, and it must be fitted, as well. However, the p^0H in the basic range will change with a changing p^0K_w . The Marquadt algorithm varies p^0K_w based on an initial approximation until the best fit is reached [10].

When performing the regression of S and E_0 , the refinement program discards data points that differ more than a certain amount from the fitted function. Thus, a good fit quality of a calibration titration is represented by a low sum of squared residuals (SSR) and a high number of points used for the fit.

Figure 10 shows a calibration titration result in $IS=0.1$ M containing the fitted parameters E_0 , slope S and p^0K_w as well as the selected data points.

```

MARQUARDT FIT of CALIB.DAT on 05-16-1994

  [Base] = 0.09940   [Acid] = 0.10040
  10.000 mL Titrand, 250.00 mL Electrolyte

E(0) = 400.018
slope = 58.902 (Eff. = 99.56%)
Gamma Titrant = 1.003
pKw = 13.822
SSR = 0.029479 using 18 points
ITERATIONS = 1 ( 3 Passes)

  Volume (mL)   mv (meas.)   mv (calc.)   DIFFERENCE   p[H]obs   p[H]calc
  *****
    0.000       257.810       257.875      -0.065       2.414     2.413
    0.700       255.940       255.963      -0.023       2.446     2.446
    1.400       253.910       253.909      +0.001       2.481     2.481
    2.100       251.710       251.687      +0.023       2.518     2.518
    2.800       249.290       249.268      +0.022       2.559     2.559
    3.500       246.650       246.611      +0.039       2.604     2.604
    4.200       243.710       243.662      +0.048       2.654     2.655
    4.900       240.370       240.347      +0.023       2.710     2.711
    5.600       236.570       236.560      +0.010       2.775     2.775
    6.300       232.060       232.138      -0.078       2.852     2.850
   13.300      -241.550      -241.620      +0.070      10.892     10.893
   14.000      -246.590      -246.575      -0.015      10.978     10.978
   14.700      -250.760      -250.704      -0.056      11.049     11.048
   15.400      -254.280      -254.242      -0.038      11.108     11.108
   16.100      -257.340      -257.335      -0.005      11.160     11.160
   16.800      -260.100      -260.080      -0.020      11.207     11.207
   17.500      -262.540      -262.548      +0.008      11.249     11.249
   18.200      -264.730      -264.787      +0.057      11.286     11.287
  *****
  
```

Figure 10: Refinement Printout of Calibration Titration

C Details on Dosimat Pump Configuration and Operation

Operating units

Injection switch: Controls the plunger in the syringe. Lifting the switch up injects a certain volume depending on working mode. Pushing down refills the syringe with titrant from the reservoir bottle.

Memory/Re-zero switch: Lifting the switch up stores the actual value of the digital counter in the memory. Pushing switch down re-zeros the digital counter.

Pump speed dial: Controls the injection speed of titrant (scale from 1 to 10).

Mode-selection switch: Controls the operating mode.

Digital meter: Displays injected volume, accurate to 0.01 mL.

Operating modes (explained in symbols on top of the Dosimat)

Mode 1: Titrant will be added at determined pump speed if the injection switch is lifted.

Mode 2: Every 5 lifts of the injection switch cause 0.01 mL titrant to be added.

Mode 3: Every lift of the injection switch will add the titrant aliquot stored in the memory.

Mode 4: Same as mode 3 except syringe refills after each injection.

Mode 7: Allows entry of desired aliquot volume into resident memory by manipulation of the injection switch.




Figure 11: Dosimat Pump

D Lid Construction




Figure 12: Lid-construction drawing

E Details of Apparatus Troubleshooting

In the process of constructing and automating the titration apparatus, several problems were encountered along the way. The following is a list of these problems and the changes that were made to correct them. These are important for future operators in order to minimize experimental errors and time spent trouble shooting.

- (1) For mV inputs to the A-to-D board greater than ± 200 mV, the output from the board lost correlation with the mV inputs and the response became random. Proper grounding of the A-to-D board remedied this problem, allowing signals of greater than ± 200 mV to be measured accurately [16].
- (2) The 12-bit resolution of the A-to-D board made single measurements of mV response from the pH electrode inaccurate [16]. Therefore, it was decided to make each p^oH point measured consist of an average of 300 separate mV readings, collected over 3 seconds with a frequency of 100 readings per second.
- (3) An external digital voltmeter was installed in parallel with the A-to-D board to provide confirmation of the stability of the mV output from the pH electrode and an independent check on the performance of the A-to-D board.
- (4) It was found that a certain type of Dosimat syringe was constructed in a manner that allowed a potential gradient to be established between the Dosimat itself and the pH electrode though the titrant solution. This introduced an error of up to ± 3 mV in our p^oH readings corresponding to a Δ p^oH of 0.05. A different syringe with more extensive insulation was then employed, eliminating the external voltage drop.
- (5) Originally it had been planned to conduct amino-acid and protein titrations in (NH₄)₂SO₄ solutions. However, the pH electrode calibrating must be conducted in solutions of strong acid and strong base which have no secondary hydrogen-ion equilibria in order for the program CAL-ANAL to determine electrode slope and intercept accurately. Therefore it was decided to study amino acid and proteins in solutions of KCl instead.

F Program Flow Sheet for CAL-TITR and TITRATE

Figure 13 is a flow sheet for the programs CAL-TITR and TITRATE. Before any titrant is added, the initial p^oH of the titrant solution is measured. This is done by averaging 300 separate mV readings taken at a frequency of 100 readings per second. The mV readings are converted to p^oH using the appropriate values of electrode slope and intercept from CALIB.PAR. The electrode slope and intercept will be the generic default values in the case of CAL-TITR (electrode-calibration titration; slope = 60 mV/pH, E₀ = 400 mV). In the case of TITRATE (titration of ligand after electrode has been calibrated), they will be equal to the fit values calculated by CAL-ANAL. Four of these converted p^oH measurements (each composed of the average of 300 mV readings) are made. If the average of the first three is within 0.01 pH units of the last, the p^oH is considered to be in equilibrium and the p^oH point is written to a data file along with the volume of titrant added (zero in the case of the initial p^oH measurement when no titrant has been added yet). However, if this convergence criteria is not met, the p^oH is considered to be not yet equilibrated, and the program collects another p^oH measurement (again composed of the average of 300 mV readings taken over 3 seconds) and tests for convergence with the four most recent p^oH readings. The process of convergence testing proceeds for up to 9 separate p^oH measurements, after which if convergence has not been archived, the average of the last 4 p^oH measurements is taken as the value of p^oH and is written to the data file*. After the most recent p^oH and volume-added titration point is written to the file, the program checks to see if the titration endpoint has been reached (final volume for CAL-TITR; final p^oH for TITRATE). If not, then an aliquot of titrant is injected, the defined delay time elapses, and the program again begins measuring and averaging the new p^oH. When the titration endpoint is reached, the program will terminate.

*In all our calibration and ligand titrations, the hydrogen-ion dissociation reactants equilibrated very quickly, and convergence was usually achieved with the first four p^oH measurements.




Figure 13: Program Flow Sheet for CAL-TITR and TITRATE

G Setup and Operation of BETA90

The potentiometric data refinement described in this chapter can only be used for amino-acid titration data, as the maximum number of pK_a 's in BETA90 is limited to 20. The refinement of amino-acid pK_a 's takes two steps: the creating of an BETA90 input file with SETUPBTA and the modification of some values within this file to improve the goodness of fit (G.O.F.).

G.1 SETUPBTA Inputs

The program BETA90 regresses amino-acid pK_a 's from experimental titration data. BETA90 requires input of all the relevant information to perform pK_a fitting in very a specific format. The program SETUPBTA is executed prior to BETA90 [10]; it assembles all the relevant information and creates a file called *.INP with the proper format for BETA90. An exact replica of a SETUPBTA input screen is shown below, with explanation of all information prompts.

<u>Input Prompt</u>	<u>Correct Response</u>	<u>Explanation</u>
• Input titration data file name	*.DAT	
• File name for BETA input file	*.INP	
• Run descriptor	*.INP	header of summary
• # of cycles (NC)	9	
• # of formation constants to be refined (NV)	2 or 3	number of pK_a 's to be re-gressed (should do all pK_a 's)
• # of reacting components (NX)	2	amino acid and hydrogen ion
• Analytical derivatives	enter	use yes default
• Computer weighting scheme	enter	use yes default
• Total # of formation constants	2 or 3	total number of pK_a 's in amino acid

- Enter key integer

guesses for pKa's in an additive form starting with highest one

Example for an amino acid with 3 pKa's

1 , 10 , 1 , 1
1 , 14 , 1 , 2
1 , 16 , 1 , 3

opt 0: don't regress number of ligands (always 1) total number of hydrogen ions titrated
opt 1: regress

- intercept enter use 0 default
- slope enter use 1 default
- pKw p^oKw from calibration titration

At this point, SETUPBTA reads in the experimental data from *.DAT file

- enter # of protonated sites (NH) 1 or 2 # of protonated sites on amino acid upon dissolving in KCl; e.g. glu=2, his=1
- analytical millimoles of ligand 2.5 *the exact values
- concentration of titre (base is -) -0.1 *from the
- millimoles added acid at beginning 2.000 *titration sheet
- initial volume at beginning of titration (mL) 270 250 + Volume titrant added to reach starting p^oH
- initial guess for pL enter use default
- sigma of pH meter enter use default
- sigma of buret enter use default

OK. Your BETA input file is setup. It is named *.INP

G.2 Refining Amino-Acid pK_a's with BETA90

With the properly-configured *.INP file as input, BETA90 will fit p^oK_a's from the experimental titration data. After the program has terminated and written the results to the files *.SUM and PLOT.O&C, it is important to determine the quality of the p^oK_a fits performed. This is done by checking several quantities.

A first step is to issue the command "see". This creates a plot of measured p^oH, calculated p^oH (based on p^oK_a's found by BETA90) and the difference between measured and calculated p^oH's against the net charge of the amino acid (which is also calculated from the fit p^oK_a's). It is generally observed that the Δp^oH has a spike near the equivalence point in the ligand titration curve. This, is the point where the p^oH changes most rapidly with addition of titrant. Thus there are fewer data at this point in the titration than at the p^oH extremes, decreasing the sensitivity of the fit performed by BETA90 near the equivalence point; hence, we see a maximum in the Δp^oH curve here. A smaller Δp^oH spike indicates a better fit.

Another equally important indicator of the p^oK_a fit quality is the G.O.F. parameter contained in the *.SUM file. This parameter is related to the overall error in the p^oK_a fitting procedure; a smaller G.O.F. indicates a more accurate fit. A G.O.F. of 6 is considered the maximum allowable for a good fit.

It is possible to minimize the G.O.F. and the Δp^oH spike for a given titration by adjusting very slightly certain of the parameters contained in the *.INP file [10]. This is done by invoking the screen editor, taking care to preserve the exact input format required for BETA90 to recognize all the contents of the *.INP file.

Experience has shown that manual variation of 3 parameters - mass of ligand, amount of acid initially added to the titration vessel, and fit p^oK_w - by less than 1% each can improve the G.O.F. dramatically for a given set of titration data [10]. Variation in the amount of ligand results in general improvement of G.O.F.. Variation of the amount of initial acid reduces the Δp^oH spike near the equivalence point. Deviations in this curve in basic regime are diminished by varying Δp^oH. It is important to check the G.O.F. and view the contents of the PLOT.O&C file via "see" after making each change in order to determine whether overall improvement has resulted.

While this parameter variation may often result in improvement of fit quality for the calculated curve, it does not produce more than 0.2 pH units change in fit p^oK_a's.

H Experimental Procedure for Titration

The aim of this chapter is to give guidance for conducting titration experiments, collecting all necessary experimental data and especially creating the $\text{mL}_{\text{added}}/\text{p}^{\circ}\text{H}$ -output file describing this titration run. This *.DAT file is the basis for discussion of the results. In the case of protein experiments it can be converted to a net-charge vs. p°H -file (compare § 4 Discussion) or in the case of amino-acid experiments it is needed for the subsequent $\text{p}^{\circ}\text{K}_a$ -refinement.

1) Solution Preparations

Acid/Base Titrant: In order to guarantee a sufficient number of titration data points to perform accurate $\text{p}^{\circ}\text{K}_a$ fitting, it was decided to set the HCl/KOH concentration in the acidic/basic titrant to 0.1 N. Therefore, in our high IS titrations, the balance of the ionic strength was made up by adding of KCl to the titrant solutions.

- Weigh the amount of KCl necessary to reach the desired ionic strength
- Dissolve salt in ~250 mL nanopure water in a clean 2 L volumetric flask
- Add the volume of 1.00 N HCl or KOH stock-solution (Fisher Chemical SA48-4 and SP208-500) necessary to create a 0.1 N solution
- Top off the volumetric flask with nanopure water, then mix thoroughly with magnetic stirrer
- Disperse titrant into appropriate reservoir bottle; seal bottle and blanket with argon

Potassium Hydrogen Phthalate (KHP) Solution:

- Weigh 0.2 mg KHP in weighing dish on analytical balance; record exact weight to 5 decimal places
- Dissolve KHP in ~100 mL nanopure water in a clean 250 mL volumetric flask, extreme care must be exercised to make sure all KHP is transferred into the volumetric flask. Accuracy of the HCl/KOH concentration measurement depends strongly on this step

- Weigh the amount of KCl necessary to achieve the desired solution ionic strength •
Dissolve KCl in KHP solution, top off the volumetric flask with nanopure water

KCl Solution: (used in HCl/KOH titrations for pH electrode calibration)

- Weigh the amount of KCl necessary to achieve the desired ionic strength
- In a clean 1 L volumetric flask dissolve KCl in nanopure water, top off volumetric flask

Titration Solution with Amino Acid as Ligand:

- Weigh the appropriate amount of amino acid determined by solubility considerations directly into a clean, dry 250 mL volumetric flask on a analytical balance; record exact weigh to 4 decimal places
- Weigh the amount of KCl necessary to achieve the desired ionic strength into a separate weigh dish
- Rinse KCl into the 250 mL volumetric flask containing amino acid; dissolve contents of flask in nanopure water; top off volumetric flask

Titration Solution with Protein as Ligand:

- Weigh the determined amount of protein into a beaker on a tared analytical balance
- Add enough KCl solution of the desired ionic strength described above to the beaker to dissolve the powdered protein
- Transfer protein solution to a clean and dry 250 mL volumetric flask; top off volumetric flask with KCl solution; mix thoroughly
- Reserve 10 mL of the protein solution for later UV spectrophotometric measurements of protein concentration

2) Titrant Standardization

The following is a procedural outline for the measurement of exact concentrations of HCl and KOH in the acidic and basic titrant, respectively. Phenolphthalein is used as a pH indicator.

Base Titrant Standardization (KOH):

- Activate argon supply to titrant-reservoir bottles
- Turn on the temperature bath and wait until the desired temperature is reached
- Switch Dosimat to 'mode 1' to allow manual titrant dispensing (see Appendix C)
- Place syringe unit with basic reservoir bottle in the Dosimat and flush any air bubbles from the fritted injection line by dispensing about 3 mL titrant into a beaker on pump speed 6 or greater.
- Place the potassium hydrogen phthalate (KHP) KCl solution into the titration vessel
- Add 10 drops of a 0.15% phenolphthalein solution in ethanol; stir magnetically
- Rinse the basic titrant fritted injection line and place it in the vessel, securing it with a clamp
- Add carefully, beginning on pump speed 5 and slowing down to speed 1, until the pH endpoint is reached. The endpoint is defined as the first pink color that does not disappear after one minute
- Record the volume of basic titrant required to reach the pH endpoint. This volume is necessary in order to calculate the KOH concentration, using the following equation:

$$c_{base} = \frac{m_{KHP}}{M_{KHP}} / V_{base}$$

where m_{KHP} = g KHP

M_{KHP} = molecular weight of KHP, =204.22 g/mol

V_{base} = measured volume of basic titrant added

- Discard titration solution, rinse vessel and refill Dosimat dispenser syringe
- Repeat the above steps at least 2 times or until a reproducible KOH concentration is measured

Acid Standardization (HCl):

- Place 250 mL KCl solution of desired ionic strength in titration vessel; add about 10 drops of a 0.15% phenolphthalein solution in ethanol; stir magnetically
- Place syringe unit with acid reservoir bottle in the Dosimat and flush all air bubbles from the fritted injection line
- Rinse the fritted injection line and place it in the vessel, securing it with a clamp
- Add exact 15 mL of acid with the pump in programming 'mode 4' (see Appendix C)
- Place syringe unit with basic reservoir bottle in the Dosimat and flush all air bubbles from the fritted injection line
- Perform a neutralization titration as described in the basic titrant-standardization section above
- Record the volume of basic titrant required to reach the pH endpoint. Calculate the acid concentration using the following equation:

$$c_{acid} = \frac{(c_{base} \cdot V_{base})}{V_{acid}}$$

- Repeat the above steps at least 2 times or until reproducibility in c_{acid} is demonstrated

3) pH Electrode Calibration

- Activate argon supply to the titrant-reservoir bottles
- Turn on the temperature bath and wait until the desired temperature is reached
- Carefully rinse the titration vessel and stir bar with nanopure water
- Place the acid titrant syringe unit into the Dosimat and flush all air bubbles from the injection line using pump 'mode 1'
- Place 250 mL of KCl solution, of the same ionic strength as the subsequent potentiometric titration, in the titration vessel; stir magnetically
- Seal vessel with lid and blanket with argon
- Insert pH electrode into the titration vessel through the proper hole in the lid. Make sure rubber ring around the top of electrode makes a tight seal with the lid. Switch

pH meter to mV mode. Note: always switch the pH meter to 'standby' when removing electrode from solution

- Dispense 10 mL of acid titrant into the vessel using the Dosimat in 'mode 4'
- Exchange the acidic titrant syringe unit with the basic titrant syringe unit. When the basic titrant syringe unit is secure, flush all air bubbles from the injection line (mode1)
- Allow 5 minutes for solution temperature and electrode response to equilibrate
- Program Dosimat to add basic titrant aliquots of 0.7 mL in 'mode 7', then switch to 'mode 3' (see Appendix C)
- Run calibration-titration program CAL-TITR
- When calibration titration is complete (~20 minutes) place pH meter on 'standby' and remove electrode from the vessel immediately, rinse the electrode and place it in the storage buffer
- Run program CAL-ANAL to analyze electrode-calibration data. Record calibration parameters in the electrode logbook. Attach printouts of calibration data and results for record keeping (see figure 10 in Appendix B)

4 Basic Ligand Titration (KOH)

- Rinse and dry the titration vessel
- Place 250 mL of amino acid or protein solution in the vessel
- Flush all titrant injection lines
- Seal vessel with lid and blanket with argon
- Rinse pH electrode and place it into vessel as described above. Set pH meter on mV mode
- Begin gentle stirring. Allow 5 minutes for the solution temperature and electrode response to equilibrate and for argon blanket to establish
- Add the determined volume of acid using 'mode 1' or '4' until p^oH equals 2.5
- Program the basic titrant increment to be added in 'mode 7'; switch to 'mode 3'
- Run titration program TITRATE
- Record all relevant information on the titration information sheet
- When titration is complete, set pH meter to 'standby', remove and rinse pH electrode, close electrode reference solution fill-cap and place in storage solution.
- Make a back-up copy of *.DAT on floppy for archiving and transfer to computer with program BETA90

For protein titrations it is possible to reduce the cost of a titration by preparing 50 mL of protein solution instead of 250 mL. After the electrode-calibration titration, place this 50 L of protein solution in a 150 mL beaker on a platform inside the titration vessel. For temperature control, fill the titration vessel with water so that the sides of the beaker containing the protein solution are covered.

For best pH electrode response and accurate measurements it is recommended to alternate the two pH electrodes in use. Furthermore, an exchange of the pH electrode filling solution every 5 titration experiments contributes to a faster recover of the pH electrode after a ligand titration.

I Experimental Data for all Measured pKa

Ionic Strength 0.1 M KCl

Aminoacid Carbox.acid Ammon.group Substituent G.O.F.

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Ala	2.361	9.849		3.238
	2.353	9.875		3.858

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Arg	2.146	9.163	13.424	1.935
	2.058	9.072	12.676	3.050
	2.047	9.042	12.872	1.977

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Asn	2.175	8.731		0.878
	2.191	8.729		0.886

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Asp	2.000	9.758	3.727	1.781
	1.992	9.685	3.719	1.884

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Cys	2.001	10.247	8.194	1.823
	1.924	10.440	8.254	1.218
	1.938	10.404	8.233	1.770

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Glu	2.238	9.647	4.184	2.236
	2.206	9.638	4.161	2.070

Aminoacid Carbox.acid Ammon.group Substituent G.O.F.

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
His	1.731	9.160	6.062	2.99
	1.698	9.143	6.059	2.436

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Ile	2.330	9.632		2.492
	2.328	9.645		1.553

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Leu	2.332	9.651		1.764
	2.318	9.631		1.303

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Lys	2.106	9.216	10.826	5.298
	2.089	9.205	10.788	1.677

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Tyr	2.343	9.068	10.181	0.977
	2.104	9.004	10.152	6.886

Aminoacid Carbox.acid Ammon.group Substituent G.O.F.

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Ala	2.442	9.735		4.778
	2.449	9.783		3.835

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Arg	2.285	9.264	12.758	3.807
	2.275	9.228	12.437	4.958
	2.264	9.339	13.402	2.250

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Asn	2.346	8.787		2.591
	2.223	8.892		3.666
	2.273	8.816		2.205

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Asp	2.051	9.512	3.652	1.981
	2.040	9.649	3.699	3.236

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Cys	2.013	10.224	8.200	0.968
	2.009	10.287	8.273	1.817
	(1.998)	(10.376)	(8.304)	2.540

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Glu	2.066	9.442	4.062	2.163
	2.258	9.510	4.145	5.357
	2.230	9.505	4.135	2.477

Aminoacid Carbox.acid Ammon.group Substituent G.O.F.

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
His	1.960	9.206	6.284	2.373
	1.921	9.146	6.232	1.577

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Ile	2.439	9.758		3.635
	2.426	9.679		2.944
	2.431	9.728		4.407

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Leu	2.437	9.793		3.962
	2.426	9.702		3.337

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Lys	2.272	9.374	10.852	1.586
	2.313	9.400	10.842	5.097

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Tyr	2.288	9.083	10.022	1.921
	2.376	9.105	10.019	2.745
	(2.457)	(9.103)	(9.974)	5.551

Ionic Strength 1.0 M KCl

Aminoacid Carbox.acid Ammon.group Substituent G.O.F.

Aminoacid	Carbox.acid	Ammon.group	Substituent	G.O.F.
Ala	2.713	10.059		6.767
	2.727	10.132		6.066

Arg	2.549	9.633	12.824	4.839
	2.558	9.636	12.846	5.282

Asn	2.452	9.094		4.074
	2.471	9.172		2.637

Asp	2.244	9.925	3.954	3.988
	2.273	9.856	3.938	4.691

Cys	2.269	10.372	8.433	3.918
	2.277	10.280	8.356	4.184
	2.268	10.513	8.517	2.224

Glu	2.493	9.766	4.437	4.142
	2.507	9.711	4.428	5.979

Aminoacid Carbox.acid Ammon.group Substituent G.O.F.

His	2.161	9.465	6.617	4.548
	2.140	9.435	6.581	4.339

Ile	2.718	10.055		5.963
	2.711	10.030		6.876

Leu	2.700	9.974		5.505
	2.701	10.044		5.724

Lys	2.603	9.711	11.043	6.122
	2.641	9.659	10.949	5.254
	2.612	9.738	11.045	3.155

Tyr	2.590	9.380	10.210	2.558
	2.641	9.346	10.151	2.120

Ionic Strength 3.0 M KCl

K Sample Titration-Data Recording

Amino Acid - Titrations		Date : 7 / 20 / 94	
Ligand :	<u>Glutamic Acid</u>	IS = <u>3.0 M</u>	
File :	<u>Gluti301.DAT</u>	Electr.Nr. <u>3</u>	
pKa :	Carbox.acid	Ammon.group	Substituent
	<u>2.493</u>	<u>9.766</u>	<u>4.437</u>
GOF :	<u>4.142</u>		
Conc. Ligand : n = 2.6173 mmol m = 0.385 g			
Acid/Base-Pair : KOH c = 0.0990 n			
HCl c = 0.0996 n			
V add : 20 ml HCl -> pH = 2.513			
dV : 0.4 ml KOH -> pH = 11			
pKw: 14.086			

Figure 14: Example for fillout titration sheet for record keeping

LEAST-SQUARES REFINEMENT SUMMARY 7-18-1994
Program: BETA90 (1-14-91) 9:29:21

GLUT301.INP

pH-OBS CORRECTED WITH INTERCEPT .0000 AND SLOPE 1.0000
pK_w = 14.00400
NUMBER OF OBSERVATIONS IS 188 pH RANGE: 2.513 TO 11.018
TOTAL VOLUME OF TITRANT ADDED: 74.800 mL
TOTAL NO OF ACID/BASE SITES 2
NUMBER OF DEPROTONATED(SALT) SITES 0
NUMBER OF PROTONATED SITES 2

MILLIMOLES OF LIGAND 2.61600000
MOLARITY OF TITRANT -.099000
VOLUME OF ADDED ACID .000000
MOLARITY OF ADDED ACID .000000
INITIAL SAMPLE VOLUME 270.000000
MILLIMOLES OF STOCK ACID 1.96100000
Sigma(meter) = .003 Sigma(buret) = .002

INPUT DATA

#	BETA	DAMPING	DP(J)	INDICES
1	10.00	1	.000E+00	1 1
2	14.40	1	.000E+00	1 2
3	16.60	1	.000E+00	1 3

PARAMETERS AFTER FINAL REFINEMENT

	OLD	CHANGE	NEW	ERROR	SHIFT/ERROR
1	9.7656	.0000	9.7656	.0017	.0026
2	14.2026	.0000	14.2026	.0024	.0104
3	16.6961	-.0004	16.6957	.0038	.1014

Value of convergence function DELGOF= .0008

FINAL AGREEMENT FACTORS:

Sum|O-C| IS 1.646E+00 Sum[W*(O-C)**2] IS 3.174E+03
Sum[Obs] IS 1.152E+03 Sum[W*(Obs)**2] IS 9.580E+08
Sqrt[Sum(W*(O-C)**2)/(No-Nv)] = 4.1421 (G.O.F.)
R = .00143 Rw = .00182

PROGRAM TERMINATED FOR REASON 3

CORRELATION MATRIX:

1	1.0000	.6902	.4368
2	.0000	1.0000	.4779
3	.0000	.0000	1.0000

Figure 15: Printout of summary file *.SUM containing all titration run information

Glutamic Acid Titration #301
July 20, 1994




Figure 16: Plot of titration curves from PLOT.O&C file

Nomenclature

a	-	activity
c	mol/l	concentration
E	mV	measured hydrogen-ion activity
E ₀ '	mV	electrode reference potential
E ₀	mV	electrode intercept
K _a	-	hydrogen-ion dissociation constant
m	g	mass
M	g/mol	molecular weight
n	-	number of site
P	-	protein
R	J/mol·K	gas constant
S	mV/pH	slope
T	K	temperature
V	mL	volume
z	-	charge per amino acid or protein

Δ		difference
γ	-	activity coefficient

Subscript:

i	species i
net	net charge
prot	protein

Superscript:

o	calculation based on concentration
---	------------------------------------

References

1. Shih, Y.-C., J.M. Prausnitz, and H.W. Blanch, *Some Characteristics of Protein Precipitation by Salt*. Biotech. Bioeng., 1992. **40**: p. 1155.
2. Vlachy, V., H.W. Blanch, and J.M. Prausnitz, *Liquid-Liquid Phase Separation in Aqueous Solutions of Globular Proteins*. AIChE J., 1993. **39**: p. 215.
3. Coen, C.J., H.W. Blanch, and J.M. Prausnitz, *Protein in Aqueous Electrolyte Solutions: Measurement of Salting-Out Phase Equilibria and Comments on Protein Intermolecular Potentials*. 1992, University of California, Berkeley.
4. Martell, A.E. and R.M. Smith, *Critical Stability Constants*. Vol. Volume 1: Amino Acids. 1974, New York: Plenum Press.
5. Bashford, D. and M. Karplus, *pKa's of Ionizable Groups in Proteins: Atomic Detail from a Continuum Electrostatic Model*. Biochemistry, 90. **29**(44): p. 10219.
6. Tanford, C., *Physical Chemistry of Macromolecules*. 1961, New York: John Wiley & Sons.
7. Harris, D.C., *Quantitative Chemical Analysis*. 3rd ed. 1991, New York: W. H. Freeman and Company.
8. Bates, R.G., *Determination of pH Theory and Practice*. 2nd ed. 1973, New York: John Wiley & Sons.
9. Zemaitis, J.F., et al., *Handbook of Aqueous Electrolyte Thermodynamics*. 1986, New York: AIChE, DIPPR.
10. Franczyk, T.S., *Potentiometric Titrator Software Guide*. 1990, University of California, Berkeley.
11. Tinoco, S., Wang, *Physical Chemistry*. 2nd ed. 1985, Englewood Cliffs, New Jersey: Prentice-Hall, Inc.
12. Orion, *ROSS pH Electrode Instruction Manual*. 1990.
13. Tanford, C. and M.L. Wagner, *Hydrogen Ion Equilibria of Lysozyme*. JACS, 1954. **76**: p. 3332.
14. Boyer, P.D., *The Enzymes*. Vol. Volume 7. 1972, New York: Academic Press.
15. Fersht, A., *Enzyme Structure and Mechanism*. 2 ed. 1985, New York: W. H. Freeman and Company.
16. Keithley, *User Guide for Data Acquisition Boards*. 1993.