


UV-Vis Characterization of Aromatic Content in Bio-oil

Emily Bladorn¹, Asanga Padmaperuma²

(1) STAR (2) Pacific Northwest National Laboratory




Background:

- Bio-oil is an alternative source currently being studied for its use as an alternative and renewable transportation fuel.
- Lignin present in the biomass feedstock is converted to aromatics during the oil production.
- Absorbance of these functionalities can be used to quantify the aromatic content in bio-oils.
- UV-Vis is an inexpensive method by which specific target compounds can be quantified within a mixture.
- Beer-Lambert law can be used to quantify the amount of specific molecules within a sample based on absorption.
- The relationship between concentration of bio-oil mixture and the absorption of UV light is shown in Figure 1.
 - Linear response indicates that the unknown will obey the Beer-Lambert law.
 - Different oil types have different slopes.
- This plot illustrates why we can use colorimetric techniques to analyze bio-oil.
- We used the Standard Addition method for the analyses, which adds a known amount of a standard to a unknown mixture. The mixture contains the standard, so when plotted, the amount of target compound within the sample can be extrapolated.



$$A = \epsilon bc$$

A = Absorption
 ϵ = Molar Absorbivity
 b = Path Length (cm)
 c = Concentration (mol/L)



Sample Preparation:

- Unknown oils were prepared using fast pyrolysis (A), catalytic fast pyrolysis (B and C). For the samples B and C were prepared using different catalysts.
- Solutions were created by using a fixed amount of bio-oil and an increasing amount of Eugenol Standard. All samples were prepared by weight.


Results and Future Work:

- The calculated amount of aromatic groups in each bio-oil sample is summarized in the Table.
- It is possible with this method to quantify certain compounds within bio-oil by using the Beer-Lambert law.
- Further testing and repeated trials with each type of bio-oil is of importance for verifying the reproducibility of these results.
- As this method was tested with one standard, testing different compounds would be of interest. Bio-oil contains a vast number of compounds hence creating a need for development of a method to streamline the sample preparation and data collection.

Colorimetric Analysis:

Each sample was run on a SpectraMax M5 UV-Vis by Molecular Designs. Absorbance spectra were collected between 200-400 nm.

Sample ID	Aromatic content (Eugenol (mol/L))
A1	4.32×10^{-5}
A2	4.49×10^{-5}
A3	3.93×10^{-5}
A4	4.46×10^{-5}
B1	1.17×10^{-4}
C1	6.99×10^{-5}



Data Analysis:

- Data was plotted as Absorbance vs. Concentration of Eugenol. The trend line was extrapolated to the x axis in order to calculate concentration of the standard in the bio-oil.
- The aromatic content is reported as moles/L of Eugenol.

The 2018 STEM Teacher and Researcher Program and this project have been made possible through support from Chevron (www.chevron.com), the National Marine Sanctuary Foundation (www.marinesanctuary.org), the National Science Foundation through the Robert Noyce Program under Grant #1836335 and 1340110, the California State University Office of the Chancellor, and California Polytechnic State University in partnership with PNNL. Any opinions, findings, and conclusions or recommendations expressed in this material are those of the author(s) and do not necessarily reflect the views of the funders.