

CEJ Ankle Support Report

BMED 456

Christine Prothe, Erik Espinoza, James Baldwin

Table of Contents	Pages
1.0 Executive Summary	2
2.0 Introduction and Background	2
3.0 Customer Requirements and Design Specifications	2
3.1 IFU	2
3.2 Product Design Specifications	2
3.3 House of Quality	4
4.0 Stage Gate Process	5
4.1 Concept Review	5
4.2 Design Freeze	7
4.3 Design Review	8
5.0 Description of Final Prototype Design	8
5.1 Overview	8
5.2 Design Justification	8
5.3 Analysis	8
5.4 Cost Breakdown	8
5.5 Safety Considerations	8
6.0 Prototype Development	8
6.1 Model Analyses	8
6.2 Evolution of Prototypes	9
6.3 Manufacturing Process	10
6.4 Divergence Between Final Design and Final Function Prototype	12
7.0 IQ/OQ/PQ	12
7.1 DOE	12
7.2 Verification and Validation	13
8.0 Conclusions and Recommendations	15
8.1 Recommendations	15
8.2 Conclusions	15
9.0 Acknowledgments	15
10.0 Appendices	15
10.1 Appendix A: References	15
10.2 Appendix B: Project Plan (PERT Chart)	16
10.3 Appendix C: CAD Drawings	16
10.4 Appendix D: FMEA, Hazard & Risk Assessment	17
10.5 Appendix E: Pugh Chart	19
10.6 Appendix F: Vendor Information, Specifications, and Data Sheets	20
20	
10.7 Appendix G: Budget	20
10.8 Appendix H: DHF	22

1. Executive Summary

This document provides information about an ankle support device designed to allow for ankle dorsiflexion and plantar flexion. In this document the product's specifications, objectives and expectations are outlined. Included in this document is background information, stage gate process, prototype design development and considerations, and qualification requirements.

2. Introduction and Background

The CEJ is an ankle foot orthosis (AFO). AFOs are intended for users with multiple sclerosis (MS), or other diseases or disorders, who suffer from foot drop. Foot drop is caused by weakness or paralysis of the muscles involved in lifting the front part of the foot, and can be the result of MS, nerve injuries, diabetes, or spinal cord and brain disorders [1]. Current AFO devices only allow for the users foot to be held in dorsiflexion, preventing the front part of the foot from dropping during the swing phase of the gait cycle. While this is an essential aspect of the device, sitting for prolonged periods of time with the foot being held in dorsiflexion can become uncomfortable and can be tiring on the muscles of the lower leg. The CEJ AFO will hold the users foot in dorsiflexion during gait, but also have a buckle release to allow the user to release their foot into plantar flexion while seated.

3. Customer Requirements and Design Specifications

3.1. IFU

The CEJ Ankle Support is an ankle-foot orthosis (AFO) designed to allow the user to extend their ankle while seated as well as support the ankle in flexion during gait. This AFO has a quick release to allow foot drop. Traditional AFOs do not allow foot drop which causes strain on the ankle when sitting/resting due to the constant ankle flexion. This device is designed to be lightweight, affordable, and safe for the user. The ankle support is aimed to improve mobility and comfort of patients with multiple sclerosis (MS). The device is for use for patients less than 200 lbs.

3.2. Product Design Specifications

Table 1. Product Specification Matrix

	Customer Requirement	Engineering Specifications	Rationale	Testing Protocol
1	Functional	Allows ankle to extend with a quick release mechanism The brace's material does not add more than 1 inch to the circumference of the ankle	Device's function is to support ankle in flexion but should allow for ankle extension when desired Able to wear under clothes	Functional Week Long Testing to determine if releasing mechanism is effective
2	Lightweight	Weighs no more than 20 ounces	The upper limit weight of other ankle support devices is 15 ounces	Weigh on a scale
3	Safe	The ankle support will have no sharp or abrasive parts. The ankle support can statically hold 200 lbs and can withstand the impact force associated with a 200 lb user walking. The ankle support will prevent hyperextension of the knee.	The user should not be harmed when using the product. The user should be able to move without the product breaking and causing harm to the user.	Functional Week Long Testing to determine if any irritation or worse harm is caused to the user
4	Affordable	Cost does not exceed \$300.00	Similar ankle support devices have a cost of about \$160.00	Calculate the cost of manufacturing, materials, and distribution
5	Durable	Able to withstand pressures of foot when in flexion and extension up to 250 lbs Temperature resistant (0-120 F) Humidity Resistance (90%)	Needs to be able to support weight of a user in multiple environments.	Fatigue testing under maximum weight (250 lbs)

The product specification matrix, as seen in Table 1, as well as other specifications were used to conduct a conjoint analysis to determine the most important customer requirements.

Table 2. Product Design Factors and Levels.

Factor	Level 1 (1)	Level 2 (0)
Cost	\$100	\$115
Weight	10oz	12oz
flexibility	keeps foot in dorsiflexion	allows plantar flexion
rigidity	no flex with movement of the leg	allows some flex with movement of the leg
strap	velcro	strap with a snap
locking mechanism	three click release	button release
height	12 in	14 in

A regression analysis of variance was conducted on eight different conjoint cards of the different levels for the design to determine the most important factors. The conjoint cards were then voted on from best to worst by classmates. Based on the data obtained from our classmates, the cost is the most important factor (p-value = 1.06E-05) and rigidity (p-value = 0.043) was another. Rigidity was a factor we originally didn't think would be too significant, so we were surprised when it was one of the lower p-values. After talking to our sponsor, we found that it is actually one of the more important factors. This is due to the fact that it can prevent knee hyperextension. Cost was a factor we expected to be important, and we are planning on keeping the price low. The regression analysis of variance also revealed which level was most important. For cost, the conjoint analysis revealed that a cost of \$115 is preferred to \$100. As for rigidity, the conjoint analysis revealed that a some flex is important to customers but we have to be sure not to compromise the prevention of hyperextension. A factor that was not found to be important but is, is the ability to allow the foot to go into plantar flexion, not just being kept in dorsiflexion at all times.

3.3. House of Quality

Table 3. House of Quality (Rooms 1, 2, 4, 5, 6)

		Engineering Characteristics												
Improvement Directions		↓	↓		↑	↓	↑	↑	↑	↓	↑			
Units		lbs _f	%	Degrees	Degrees	Oz	n/a	psi	psi	lbs _f	n/a	Competitor Rankings		
Customer Requirements	Importance Weight Factor	Quick release	Jamming probability	Angle support in flexion	Angle at rest	Weight	Weather resistance	Material strength	Material stiffness	Straps release ease	Comfortability	OSSUR AFO Leaf Spring Foot Drop Brace	Swedish Ankle Foot Orthosis	OSSUR AFO Dynamic Drop Orthosis
Allows Plantar Flexion	5	9	9		9							1	1	1
Lightweight	4					9		3				4	4	4
Safe: Prevents hyperextension	5							9				2	2	2
Affordable	3	3						3	3	3		5	5	2
Durable	4						9	9	9			3	4	4
Raw Score (372)		54	45	0	45	36	36	102	45	9	0			
Relative Weight %		14.5	12.1	0	12.1	9.7	9.7	27.4	12.1	2.4	0			
Rank Order		2	3	9	3	6	6	1	3	8	9			

The house of quality, shown in Table 3, was used to determine the most important engineering characteristic, based on our customer requirements. It was determined that the material strength was the

most important engineering characteristic, as it had influences on the weight, the safety, the price, and the durability of the product.

4. Stage Gate Process

4.1. Concept Review

On February 20th, 2019, group members presented a concept review of the CEJ Ankle Support. The presentation included the following: background information, Indication For Use (IFU), Total Available Market (TAM), Regulatory Plan, Modified Budget, Project Plan, House of Quality, Failure Mode and Effects Analysis (FMEA), Hazard and Risk considerations, Potential Concepts Considered, Pugh Chart, and Front Runner Design description. A few of the presented topics can be found throughout the report. Below are the topics not already previously covered in the report.

Total Available Market (TAM)

The TAM was defined by identifying how many individuals have multiple sclerosis in the U.S. There are at least 400,000 individuals in the U.S with MS, with 200 new cases being diagnosed each week [3]. Foot drop is a common symptom of MS along with nerve injuries, diabetes, and spinal cord and brain disorders [1]. Foot drop is caused by weakness or paralysis of the muscles involved in lifting the front end of the foot. Because foot-droop is a common symptom in patients with MS, it is estimated that each individual can make use of the CEJ Ankle Support. If our device is roughly \$100, the TAM is approximately \$40,000,000.

Regulatory Plan

The CEJ is a Class 1 Device and is 510K exempt.

Potential Concepts Considered




Figure 1. Concept 1: Button Release

Concept 1, figure 1: A button release to allow the foot to be released from dorsiflexion to plantar flexion. The male head will be attached to the rubber joint that's attached on the lower half of the ankle support. The female back attachment will be attached to the calf support/upper half of the ankle support. When clipped the ankle support will hold the foot in a dorsiflexed position. When unclipped the ankle will have the ability to go into plantar flexion.




Figure 2. Concept 2: Rear Bar Extension

Concept 2, figure 2: A bar extension on the back of the brace, that when fully down holds the ankle in dorsiflexion, and when fully up allows the ankle to fall into plantar flexion. Both pieces of the ankle support brace will be very similar to the existing brace our patient currently wears, however, the rear adjustment will replace the stopper currently in place. The bar will push down on the heel cup and raise the ankle into dorsiflexion. When the bar is up it allows the heel cup and ankle to be in a plantar flexed position.




Figure 3. Concept 3: Double Sided Buckle Release

Concept 3, figure 3: A double side release buckle that when clicked into the buckle holds the foot in dorsiflexion, and when unclicked from the buckle allows the joint to release and get longer, and allow the foot to fall into plantar flexion. This brace would use a similar joint mechanism to the brace currently being used, but the rubber joint on this brace could be released by the buckle.

Our third concept, the buckle release, was determined to be the most beneficial design concept. This concept was compared to the other two concepts we've developed and the "gold standard" brace, the Swedish AFO. The design we selected was found to be similar in price to the Swedish AFO and cheaper to manufacture than the other two designs. This brace was also similar to the gold standard or the others,

or better than the standard in the other preferred areas. The difference in cost than the other two was the primary determination in selecting the buckle design over the button and the posterior extension bar. The selected design had similar pluses that applied to the other two as well, but this design had no negatives.

4.2 Design Freeze




Figure 4. Mock-up of final design

4.3 Design Review

The use of a clip will allow for extension of the ankle. Just connecting the Tamarack Joint directly to the buckle could allow the possibility of rotation of the joint and the buckle at the connection point to the brace. To reduce the possibility of rotation, we secured the buckle with sheet metal, providing 2 connection points at the brace, and a 3D printed joint attachment to the buckle, also to give 2 connection points. These adjustments allow for extra strength and stability of the joint.

5. **Description of Final Prototype Design**

5.1 Overview

The final prototype design utilizes COBRA buckles along with some aluminum sheet metal and a 3D printed Tamarack Joint attachment to allow for the option of extension.

5.2 Design Justification

This design fulfills the customer requirements of providing flexion of the ankle during the swing phase of the gait cycle, while allowing extension while the user is seated. The clip allows for release of the joint, while the sheet metal and the 3D printed Tamarack Joint attachment provide stability of the joint and a connection to the clip.

5.3 Analysis

The final prototype was able to prevent foot-drop by holding the foot at about 83° in dorsiflexion. The release also allowed for a more comfortable seated position of the ankle in plantar flexion. While these metrics were met, the prevention of hyperextension of the knee was not accomplished.

5.4 Cost Breakdown

The cost breakdown can be seen in the budget, Table 8, in Appendix G. The breakdown only accounts for the materials used to make the brace. There the manufacturing and testing machines were not allotted for because they were available at no cost in labs on campus.

5.5 Safety Considerations

For our safety, we will be using heat resistant gloves to handle the hot polypropylene when modeling the device around the mold of the users leg. For the users safety, we will test each material extensively to ensure the durability of each piece to be sure it will not fail during use.

6. **Prototype Development**

6.1. Model Analyses




Figure 5. 3D scan of the patient's leg to create a positive mold from.

Using a 3D scan of the patient's leg to mold the polypropylene of our AFO. The foot mold will allow to create a comfortable mold around her ankle for optimal fit and comfort around a rigid part of the body.




Figure 6. The stacked cardboard slices from the 3D leg model, and the model covered in Bondo

6.2. Evolution of Prototypes




Figure 7. Front and side views of preliminary prototype



Figure 8. Front and side views of final prototype

The prototypes will be tested in tension and compressive forces to determine if the side-release buckles will be strong enough to withstand the weight transfer during walking. If the tensile force of the side-release buckle is stronger than the Tamarack side joints then our only worries is fatigue and impact causing the side-release buckle to fail.

6.3 Manufacturing Process

1. Scan foot and ankle using a camera (XBOX 360 Kinect) to produce 3D Model.

1.1 ReconstructMe and an Xbox Kinect camera were used to scan and create a 3D computer model of the lower leg and foot.

2. Produce a smooth mesh of the 3D model, and create slices of the model.

2.1 MeshLab and Autodesk Meshmixer were used to crop out the unwanted excess pieces of the 3D computer model from background objects, and to create a smooth mesh model of the lower leg.

2.2 Autodesk Fusion 360 and Slicer for Fusion 360 were used to break the leg model into slices that can be cut out and stacked to make a physical model

3. Use a laser cutter to cut out the slices from step 2, and stack the slices to create a physical model.

4. Cover model with Bondo.

5. Vacuum form polypropylene to mold.

5.1 Heat polypropylene to 380 F

5.2 Move oven back, so it is no longer over the polypropylene sheet

5.3 Raise the platform with the mold to the heated polypropylene sheet

5.4 Turn on the vacuum to allow the sheet to form around the mold

6. Add attachments.

6.1 Align the Tamarack Joint so that it sits right on the lateral malleolus and the medial malleolus and mark the corresponding point on the foot piece where the attachment screw should be placed

6.2 Mark on the calf piece where the top part of the Tamarack Joint lines up

6.3 Measure up from that point about 4 inches and drill 2 holes about 1.5 inches apart for the aluminum connection

6.4 Loop the aluminum through both ends of the COBRA buckle, and attach the top end to the 2 drilled holes, and attach the bottom part to the 3D printed joint attachment

6.5 Attach the Tamarack Joint to the 3D printed attachment and attach the bottom part of the joint to the foot piece

6.6 Place padding over the screw heads on the interior of the brace

6.4 Divergence Between Final Design and Final Functional Prototype

The final prototype does not deviate from the final design, but it does from the preliminary functional prototype. The preliminary functional prototype did not account for the curvature of the calf piece of the brace, so with the modified joint pieces, the calf and the foot pieces were slightly unaligned. For our final prototype, we flattened the sides to allow for a straight connection.

7. **IQ/OQ/PQ**

7.1. DOE

Table 4. Design of Experiments

Engineering Metric	Specification	Test Method	Test Apparatus Location	Apparatus Experience / Training	Sample Size	Power
Costs less than \$300	The budget give to us by the BMED department, and external funding	Make sure all modifications to brace cost < \$300	Anywhere, computer tracking receipts	None	1	95%
Fits as well as current brace	Want to make sure the support fits on Denise as well as	Put ankle support on Denise to verify the fit.	Anywhere Denise and Kim are.	Nurse/ Physician Training (Kim)	1	100%
Tamarack joint - Ultimate Strength*	The hinge needs to be strong enough to withstand forces associated with the ankle	Tension/ Compression test to find the Ultimate Strength of the hinge	192-135	Bluehill software manual needs to be followed. Minimal Experience	10	95%

Tamarack joint - Fatigue Testing*	The hinge needs to be durable enough to withstand the many times denise stands and sits throughout the day	Cyclic/ Fatigue Testing	192-135	Fatigue testing training	5	95%
Polypropylene- Fatigue Testing*	The support structure needs to be able to support the patient's weight	Tension/ Compression test to find the Ultimate Strength	192-135	Bluehill software manual needs to be followed. Minimal Experience	5	95%
Tamarack Connector pull out strength	The connector needs to be able to withstand the static and dynamic forces applied when in use. (200 lbf / 2 sides) * 3(F.O.S) = 300lbf	Tension test to find the Ultimate Strength	192-135	Bluehill software manual needs to be followed. Minimal Experience	5	95%
Sheet Metal pull out strength	The sheet metal connection needs to be able to withstand the static and dynamic forces applied when in use. (200 lbf / 2 sides) * 3(F.O.S) = 300lbf	Tension test to find the Ultimate Strength	192-135	Bluehill software manual needs to be followed. Minimal Experience	5	95%

* These engineering metrics will not be tested as they are not products we have designed and have been tested to work in our application.

The design of experiments, Table 4, shows the engineering metrics important to this device and the tests necessary to determine the individual pieces reliability. These testing protocols were necessary to ensure the overall reliability and effectiveness of the device.

7.2 Verification and Validation

The brace could withstand a tensile force up to 277.9lbs, shown in figure 9. In compression, the brace folded in half at the joint under 200lbs of compression and did not break at any point. For the specifications, it met most requirements, Table 5, but was not able to prevent hyperextension of the knee.




Figure 9. Graph from tensile test results

Table 5. Verification and validation of device specifications

Specifications	Engineering Metric	Results	Requirement Met? (Y/N)
Allows ankle to extend with a quick release mechanism	Measure angle ankle when released from dorsiflexion, minimum 30°	Ankle angle increases 36.845°	Y
The brace's material does not excessively add to the circumference of the ankle	Measure circumference of device, ensure it is less than 2.5 inch larger than ankle circumference	11 11/16 in -9 5/16 in=2 3/8 in	Y
Weight	Weighs no more than 20 ounces	17.45 oz	Y
The ankle support will have no sharp or abrasive parts.	Sand all parts	Smooth edges	Y

Withstand static and impact forces associated with user walking	200 lb tension/compression test for Cobra Buckle, Tamarack connector, and sheet metal assembly	Withstood 277.9 lbs in tension, 200 lbs in compression. Polypropylene used for current application, impact forces not tested.	Y
Cost	Less than \$300	\$176.58* *Does not include shipping and handling	Y
Attachments withstand pressures of foot when in dorsiflexion and plantar flexion.	200 lb tension/compression test for Cobra Buckle, Tamarack connector, and sheet metal assembly	Withstood 277.9 lbs in tension, 200 lbs in compression.	Y
Prevent hyperextension of the knee.	Measure angle of knee between mid-stance and toe off of gait cycle. Ensure less than 23°	>23°	N

8. Conclusions and Recommendations

8.1 Recommendations

Some recommendations for future models of this brace would be the integration of a guide wire for the COBRA buckle, the use of buckles with a shorter male head, rounding out the body of the brace, and setting the buckle into the polypropylene. The addition of the guide wire and the shorter male head would allow for more easy reattachment of the joint to regain the flexion held position. The buckles being set into the polypropylene and the rounding out of the body of the brace would reduce the bulk of the brace.

8.2 Conclusions

We found that the brace was able to prevent foot drop by holding the foot in dorsiflexion during gait, but it wasn't able to prevent hyperextension of the knee. We also came across some issues with the molding process, but through trouble-shooting the process, we were able to achieve a polypropylene mold from the Bondo/ cardboard mold of the user's calf. We also found that her current brace was not a perfect mold of her lower leg, and was slightly bigger than the mold we created.

9. Acknowledgements

The development of the CEJ would not have been possible without the support of our sponsors Kim Chartrand, Denise and David West, Senior Project faculty advisors John Gerrity


and Dr. Whitt, our financial sponsor Dr. Christopher Porterfield and the Engineering Department at California Polytechnic State University in San Luis Obispo.

10. Appendices

10.1. Appendix A: References

[1] Foot Drop Symptoms and Causes. Mayo Clinic.

10.2. Appendix B: Project Plan (PERT Chart)



		Task Mode ▾	Task Name ▾	Duration ▾	Start ▾	Finish ▾
1			Project Specifications	13 days	Mon 1/14/19	Wed 1/30/19
2			Overall product specifications	1 wk	Mon 1/14/19	Fri 1/18/19
3			Hardware Specifications	4 days	Mon 1/21/19	Thu 1/24/19
4			Market Research	4 days	Fri 1/25/19	Wed 1/30/19
5			Supplier Specifications	15 days	Fri 1/25/19	Thu 2/14/19
6			Hardware	4 days	Fri 1/25/19	Wed 1/30/19
7			Market Research	4 days	Thu 1/31/19	Tue 2/5/19
8			Concept Review	7 days	Wed 2/6/19	Thu 2/14/19
9			Product Design	15 days	Fri 2/15/19	Thu 3/7/19
10			Hinge	1 wk	Fri 2/15/19	Thu 2/21/19
11			Outer Cover	1 wk	Fri 2/22/19	Thu 2/28/19
12			Functionality	1 wk	Fri 3/1/19	Thu 3/7/19
13			Straps	3 days	Fri 2/22/19	Tue 2/26/19
14			Project Integration	18 days	Wed 2/27/19	Fri 3/22/19
15			Prototypes Designed	8 days	Wed 2/27/19	Fri 3/8/19
16			Design Freeze	4 days	Mon 3/11/19	Thu 3/14/19
17			Preliminary Prototype Plan	6 days	Fri 3/15/19	Fri 3/22/19
18			Next Quarter	55 days	Mon 4/1/19	Fri 6/14/19
19			Test and Manufacturing Plan Presentation	4.2 wks	Mon 3/25/19	Mon 4/22/19
22			Functional Prototype Presentation	31 days	Mon 5/6/19	Mon 6/17/19
23			Project Report Draft Due	17 days	Mon 5/6/19	Tue 5/28/19
24			Poster Due	17 days	Mon 5/6/19	Tue 5/28/19
25			Project Expo	4 days	Tue 5/28/19	Fri 5/31/19
26			Final Design and Prototype Presentation	21 days	Mon 5/6/19	Mon 6/3/19
27			Final Report	26 days	Wed 5/8/19	Wed 6/12/19
28			Give Denise Brace	1 day	Fri 6/14/19	Fri 6/14/19
20			Tensile Testing At MATE	22 days	Tue 4/23/19	Wed 5/22/19
21			Assemble Prototype	10 days	Tue 4/23/19	Mon 5/6/19

Figure 10. PERT Chart through spring quarter

10.3. Appendix C: CAD Drawings




Figure 11. CAD drawing of 3D printed joint attachment.

10.4. Appendix D: FMEA, Hazard & Risk Assessment

Table 6. Failure Mode and Effects Analysis (FMEA)

Component Name	Possible Failure Mode	Type	Cause of Failure	O C C	D E T	S E V	RPN	Effect of Failure on System	Failure Improvement Alternative Actions (actions to fix the problem...)
COBRA Buckle	broken	C	improper use of the buckle, or being crushed under high pressure could cause the buckle to break	4	6	6	144	would keep ankle in plantar flexion, wouldn't provide dorsiflexion support	perform stress tests on the material used to make the buckle to make sure it won't break under slightly higher pressures than would normally be applied
Velcro	worn	C	inadequate use, user applies too much wear and tear, dirt and water	5	5	5	125	could impact the device being held in place	testing the velcro under unusual circumstances to ensure reliability
Polyester strap	worn	C	extra strain applied to it due to improper use	2	3	5	30	wouldn't hold the leg piece and the foot piece together when the buckle is released to allow plantar flexion	repeated stress tests to be sure that the strap can handle repeated use and stresses
Tamarack Joint (M-95)	worn	M	Manufactured incorrectly	2	2	6	24	would prevent the ankle from being held in dorsiflexion	perform repeated stress tests to be sure the joint won't fail after repeated use
Polypropylene	cracked	W	improper manufacturing of material, wears out too easily	3	1	5	15	if the crack grows, device could be unusable	stress tests to be sure the polypropylene will withstand the basic pressures applied by the user

Metal strap loop	bent	C	improper use, being used for more than the design allows	3	1	5	15	would impact the strap, could make uncomfortable or difficult to use strap	pressure test to be sure the loop won't bend under unusual pressure
Rivet	bent	C	improper use, user applies too much pressure, using for improper activities	3	1	5	15	would impact the release from dorsiflexion	stress test to be sure that repeated use and release of the joint won't affect the piece

The FMEA, Table 6, shows the components of the device that have the highest potential of failure. We determined that the double sided buckle was the most significant, due to the fact that if it fails, the device will no longer be able to maintain the support in dorsiflexion. This is crucial, as the dorsiflexed position is the primary goal of this device.

Table 7. Hazard and Risk Assessment

Description of Hazard	Planned Corrective Action	Planned Date	Actual Date
User of ankle support can venture into extreme environmental conditions such as fog, humidity, cold, and hot temperatures.	Ankle support will be water resistant and made from high and low temperature resistant materials.	1/27/2019	1/30/19
Burns from molding Polypropylene to positive model of leg.	Heat resistant gloves and careful handling of heat remoldable polypropylene.	2/16/2019	6/4/19

The hazards and risks, are precautions, possible from the manufacturing and use of this device are shown in Table 7. We will follow safety regulations along with taking extra precautions to prevent burns, and we will test the materials to prevent injury to the user.

10.5. Appendix E: Pugh Chart

Table 8. Pugh Chart

		Concepts		
Criteria	Swedish AFO	1	2	3
Cost	D A T U M	-	-	S
Weight		S	S	S
Allows Plantar Flexion		+	+	+
Prevents Knee Hyperextension		+	+	+
Durable		S	S	S
# of Pluses		2	2	2
# of Minuses		1	1	0

The Pugh Chart, Table 8, shows our three proposed concepts compared to a current device available on the market that is similar to the custom device the user is currently using. We found that the double sided buckle release would be the best of our three concepts for accomplishing the user's desired specifications.

10.6. Appendix F: Vendor Information, Specifications, and Data Sheets

McMaster-Carr: <https://www.mcmaster.com/>

Kingsley Orthotic and Prosthetic Supply:
<http://www.kingsleymfg.com/KMFGStore/Default.asp>

Fastenal: <https://www.fastenal.com/>

10.7. Appendix G: Budget

Table 9. Budget

Item Description	Manufacturer/ Distributor	Product Number	Purpose	Associated Task	Planned			
					Unit	Quantity	Cost/Unit	Total Cost
Light-Duty Cam Buckles	McMaster Carr	29705T31	To fasten strap	Keep strap locked	1	5	\$15.84	\$15.84
Low-Stretch Polyester Webbing	McMaster Carr	3444T21	Strap to attach support to ankle	Keep support on ankle	5	1	\$0.67	\$6.70
Polypropylene Sheet (3/16")	McMaster Carr	2898K42	Distribute load on ankle	Support the ankle	1	1	\$26.81	\$26.81
Tamarack Joint*	Kingsley	B742-85-M	Provide semi-rigid support	Allow attachment of supports, allow smooth gait pattern	2	1	\$33.59	\$33.59
Sheet Metal (Aluminum)	McMaster Carr	9015T131	Provide attachment from buckle to brace	Keeps buckle in a rigid position	1	1	\$14.14	\$14.14
GT COBRA Quick Release Buckle	AustriAlpin	PLH-SRB-CO B-000	To allow joint release	Allows for release from flexion to extension	2	2	\$10.00	\$20.00
Resilient Polyurethane Foam Sheet	McMaster Carr	86375K111	Provide comfort and support for the foot in the brace	Allow for non direct skin to polypropylene contact	1	1	\$11.18	\$11.18

Metric 18-8 Stainless Steel Hex Drive Flat head Screws	Fastenal	11547539	Provide connection between COBRA buckle, sheet metal, and brace	Allow for the various pieces making up the releasable joint to connect to each other and the brace	1	10	\$23.20	\$23.20
Metric Medium-Strength Steel Hex Nuts	McMaster Carr	90592A090	Screw joint to attachment	Holds the screw in place	1	100	\$1.32	\$1.32
Heat-Set Inserts for Plastic	McMaster Carr	94459A150	Provide a strong hold for the bolt in the plastic	Allow for the bolt to thread into the polypropylene	1	50	\$9.37	\$9.37
Washers	McMaster Carr	93475A230	Provide distribution of the pressure on the nut evenly over the surface	Allow for distribution of pressure, and provide a flat surface for the nut to sit against	1	100	\$1.86	\$1.86
Bondo	McMaster Carr	92930A420	To cover the mold of the patients lower leg	To allow for the part to be molded around the model	1	32oz	\$12.57	\$12.57
							Total \$176.58*	

Our current projected budget is shown in Table 9. Our total allocated budget is \$200, so we are currently about \$25 below our total amount. Having this buffer room will be beneficial for the possibility of prototypes that don't work as well as planned. *The total doesn't account for shipping and handling.

10.8. Appendix H: DHF

Engineering Specifications and Product Specifications

Table 10. Product Specification Matrix

	Customer Requirement	Engineering Specifications	Rationale	Testing Protocol
1	Functional	Allows ankle to extend with a quick release mechanism The brace's material does not add more than 2 inches to the circumference of the ankle	Device's function is to support ankle in flexion but should allow for ankle extension when desired Able to wear under clothes	Functional Week Long Testing to determine if releasing mechanism is effective
2	Lightweight	Weighs no more than 10 ounces	The upper limit weight of other ankle support devices is 11.5 ounces	Weigh on a scale
3	Safe	The ankle support will have no sharp or abrasive parts. The ankle support can statically hold 200 lbs and can withstand the impact force associated with a 200 lb user walking. The ankle support will prevent hyperextension of the knee.	The user should not be harmed when using the product. The user should be able to move without the product breaking and causing harm to the user.	Functional Week Long Testing to determine if any irritation or worse harm is caused to the user
4	Affordable	Cost does not exceed \$300.00	Similar ankle support devices have a cost of about \$160.00	Calculate the cost of manufacturing, materials, and distribution
5	Durable	Able to withstand pressures of foot when in flexion and extension up to 250 lbs Temperature resistant (0-120 F) Humidity Resistance (90%)	Needs to be able to support weight of a user in multiple environments.	Fatigue testing under maximum weight (250 lbs)

Design History Record (DHR)

Table 11. Design History Record

Process	Completed by	Expected Date	Actual Date
3D scan of leg	James Baldwin, Erik Espinoza	2/25/19	2/25/19
Laser cut cardboard slices	James Baldwin, Erik Espinoza, Christine Prothe	4/10/19	4/17/19
Stack cardboard slices to create mold	James Baldwin, Erik Espinoza, Christine Prothe	4/10/19	4/17/19
Cover mold in Bondo	Christine Prothe	4/30/19	5/1/19
Create joint piece	James Baldwin, Erik Espinoza	5/1/19	5/1/19
Test Prototype	James Baldwin, Erik Espinoza, Christine Prothe	5/22/19	5/24/19
Mold final brace piece	James Baldwin, Erik Espinoza, Christine Prothe	5/29/19	6/4/19
Create final prototype	James Baldwin, Erik Espinoza	5/29/19	6/5/19

Manufacturing Process Instructions (MPI)

1. Scan foot and ankle using a camera (XBOX 360 Kinect) to produce 3D Model.

1.1 ReconstructMe and an Xbox Kinect camera were used to scan and create a 3D computer model of the lower leg and foot.

2. Produce a smooth mesh of the 3D model, and create slices of the model.

2.1 MeshLab and Autodesk Meshmixer were used to crop out the unwanted excess pieces of the 3D computer model from background objects, and to create a smooth mesh model of the lower leg.

2.2 Autodesk Fusion 360 and Slicer for Fusion 360 were used to break the leg model into slices that can be cut out and stacked to make a physical model

3. Use a laser cutter to cut out the slices from step 2, and stack the slices to create a physical model.
4. Cover model with Bondo.
5. Vacuum form polypropylene (4) to mold.
 - 5.1 Heat polypropylene to 380 F
 - 5.2 Move the oven back from the polypropylene sheet
 - 5.3 Raise the mold to the heated polypropylene sheet
 - 5.4 Turn on the vacuum to allow the sheet to form around the mold
6. Add attachments.
 - 6.1 Align the Tamarack Joint (1) so that it sits right on the lateral malleolus and the medial malleolus and mark the corresponding point on the foot piece where the attachment screw should be placed
 - 6.2 Mark on the calf piece where the top part of the Tamarack Joint lines up
 - 6.3 Measure up from that point about 4 inches and drill 2 holes about 1.5 inches apart for the aluminum connection
 - 6.4 Loop the aluminum (5) through both ends of the COBRA buckle (6), and attach the top end to the 2 drilled holes, and attach the bottom part to the 3D printed joint attachment (12)
 - 6.5 Attach the Tamarack Joint to the 3D printed attachment and attach the bottom part of the joint to the foot piece
 - 6.6 Place padding (7) over the screw heads on the interior of the brace

Installation Qualification (IQ) and Operations Qualifications (OQ)

Table 12. Installation Qualification

Item #	Part Number	Part Name	Supplier/ Manufacturer
1	B742-85-M	Tamarack Joint	Kingsley
2	29705T31	Light-Duty Cam Buckles	McMaster Carr
3	3444T21	Low Stretch Polyester Webbing	McMaster Carr
4	2898K42	Polypropylene Sheet (3/16")	McMaster Carr
5	9015T131	Sheet Metal (Aluminum)	McMaster Carr
6	PLH-SRB-COB-000	GT COBRA Quick Release Buckle	AustriAlpin
7	86375K111	Resilient Polyurethane Foam Sheet	McMaster Carr
8	94459A150	Heat-Set Inserts for Plastic	McMaster Carr
9	11547539	Metric 18-8 Stainless Steel Hex Drive Flat Head Screws	Fastenal
10	90592A090	Metric Medium-Strength Steel Hex Nuts	McMaster Carr
11	93475A230	18-8 stainless steel washers	McMaster Carr
12	N/A	Tamarack Joint connector	3D printed - Innovation Sandbox (197-205)

Table 13. Operations Qualifications

Specifications	Engineering Metric	Results	Requirement Met? (Y/N)
Allows ankle to extend with a quick release mechanism	Measure angle ankle when released from dorsiflexion, minimum 30°	Ankle angle increases 36.845°	Y
The brace's material does not excessively add to the circumference of the ankle	Measure circumference of device, ensure it is less than 2.5 inch larger than ankle circumference	11 11/16 in -9 5/16 in=2 3/8 in	Y
Weight	Weighs no more than 20 ounces	17.45 oz	Y
The ankle support will have no sharp or abrasive parts.	Sand all parts	Smooth edges	Y
Withstand static and impact forces associated with user walking	200 lb tension/compression test for Cobra Buckle, Tamarack connector, and sheet metal assembly	Withstood 277.9 lbs in tension, 200 lbs in compression. Polypropylene used for current application, impact forces not tested.	Y
Cost	Less than \$300	\$176.58* *Does not include shipping and handling	Y
Attachments withstand pressures of foot when in dorsiflexion and plantar flexion.	200 lb tension/compression test for Cobra Buckle, Tamarack connector, and sheet metal assembly	Withstood 277.9 lbs in tension, 200 lbs in compression.	Y
Prevent hyperextension of the knee.	Measure angle of knee between mid-stance and toe off of gait cycle. Ensure less than 23°	>23°	N

Bill of Materials

Table 14. Bill of materials

Item #	Part Number	Part Name	Description	Quantity	Units	Supplier/ Manufacturer
1	B742-85-M	Tamarack Joint	Keeps foot in flexion during gait	2	1	Kingsley
2	29705T31	Light-Duty Cam Buckles	Strap adjustment	2	1	McMaster Carr
3	3444T21	Low Stretch Polyester Webbing	Strap to hold the brace on the leg	2	1	McMaster Carr
4	2898K42	Polypropylene Sheet (3/16")	To make the brace	1	1	McMaster Carr
5	9015T131	Sheet Metal (Aluminum)	Rigid attachment to the brace for the buckles	1	1	McMaster Carr
6	PLH-SRB-COB-000	GT COBRA Quick Release Buckle	Allows for release from flexion to extension	2	2	AustriAlpin
7	86375K111	Resilient Polyurethane Foam Sheet	Creates a cushion between the user and the polyethylene	1	1	McMaster Carr
8	94459A150	Heat-Set Inserts for Plastic	Allows the screw to connect to the polypropylene	8	1	McMaster Carr
9	11547539	Metric 18-8 Stainless Steel Hex Drive Flat Head Screws	Provide connection between sheet metal and brace	8	1	Fastenal

10	90592A090	Metric Medium-Strength Steel Hex Nuts	Provide the screw stability for the joint metal attachment	4	1	McMaster Carr
11	93475A230	18-8 stainless steel washers	Provide a strong hold for the bolt in the plastic	4	1	McMaster Carr

Failure Mode and Effect Analysis (FMEA)

Table 15. Failure Mode and Effects Analysis (FMEA)

Component Name	Possible Failure Mode	Type	Cause of Failure	O C C	D E T	S E V	RPN	Effect of Failure on System	Failure Improvement Alternative Actions (actions to fix the problem...)
Double Sided Buckle	broken	C	improper use of the buckle, or being crushed under high pressure could cause the buckle to break	4	6	6	144	would keep ankle in plantar flexion, wouldn't provide dorsiflexion support	perform stress tests on the material used to make the buckle to make sure it won't break under slightly higher pressures than would normally be applied
Velcro	worn	C	inadequate use, user applies too much wear and tear, dirt and water	5	5	5	125	could impact the device being held in place	testing the velcro under unusual circumstances to ensure reliability
Polyester strap	worn	C	extra strain applied to it due to improper use	2	3	5	30	wouldn't hold the leg piece and the foot piece together when the buckle is released to allow plantar flexion	repeated stress tests to be sure that the strap can handle repeated use and stresses

Tamarack Joint	worn	M	Manufactured incorrectly	2	2	6	24	would prevent the ankle from being held in dorsiflexion	perform repeated stress tests to be sure the joint won't fail after repeated use
Polypropylene	cracked	W	improper manufacturing of material, wears out too easily	3	1	5	15	if the crack grows, device could be unusable	stress tests to be sure the polypropylene will withstand the basic pressures applied by the user
Metal strap loop	bent	C	improper use, being used for more than the design allows	3	1	5	15	would impact the strap, could make uncomfortable or difficult to use strap	pressure test to be sure the loop won't bend under unusual pressure
Rivet	bent	C	improper use, user applies too much pressure, using for improper activities	3	1	5	15	would impact the release from dorsiflexion	stress test to be sure that repeated use and release of the joint won't affect the piece



Figure 12. Photo of senior project team and sponsors.

From left to right: Kim (sponsor), Denise (patient), David West, James, Erik, and Christine.