

NEAR-INFRARED MOLECULAR SPECTROSCOPY USING NICE-OHMS WITH HIGH FINESSE CAVITY

TZULING CHEN, *Division of Chemistry and Chemical Engineering, California Institute of Technology, Pasadena, CA, USA*; YIWEI LIU, *Department of Physics, National Tsing Hua University, Hsinchu, Taiwan*.

Noise-immune cavity-enhanced optical heterodyne molecular spectroscopy (NICE-OHMS), taking advantage of combining the cavity enhancement and frequency modulation techniques, provides an excellent method to enable ultra sensitive detection. The advantage of noise immunity enables for shot-noise limit detection, particularly for those molecules with ultra-small dipole moments, such as overtone transitions and symmetrical molecules, NICE-OHMS provides a superior strategy to achieve sub-Doppler saturation spectroscopy. Using an optical cavity with a high finesse $>100,000$, in our previous work, we reported the sub-Doppler saturation NICE-OHMS spectroscopy for nitrous oxide (N_2O) overtone transitions using the quantum-dot (QD) laser developed at $1.28\ \mu\text{m}$. At a pressure of several mTorr, the saturation dip is observed with a full width at half-maximum of about 2 MHz. The noise equivalent bandwidth-reduced sensitivity is $1.6 \times 10^{-11} \text{cm}^{-1} \text{Hz}^{-1/2}$. The QD laser is then locked to this dispersion signal with a stability of 15 kHz at 1 s integration time. We demonstrate the potential of the (N_2O) as a marker because of its particularly rich spectrum in the vicinity of $1.28\ \mu\text{m}$, where there are several important forbidden transitions of atomic parity violation measurements. In current work, we have used a new QD laser system coupled to a high finesse cavity for NICE-OHMS of the H_2 overtone transition S1, where the dipole moment is only 30 mD. The QD laser developed at $1.16\ \mu\text{m}$ is gain-chip based and mounted in an integrated mechanism system to reduce the passive laser linewidth. The cavity finesse was measured to be 234,000 (8000) by using laser-swept cavity ring-down time measurements. After the laser locking, the laser power will be amplified by a fiber amplifier to satisfy saturation condition to achieve sub-Doppler NICE-OHMS spectroscopy.