


Third International Conference on Radiation and Applications in Various Fields of Research

> June 8 - 12 | 2015 Slovenska Plaža Budva | Montenegro www.rad-conference.org

PROCEEDINGS



PUBLISHER: RAD Association, Niš, Serbia

www.rad-association.org

FOR THE PUBLISHER: Prof. Dr. Goran Ristić

EDITOR: Prof. Dr. Goran Ristić

COVER DESIGN: Vladan Nikolić, M.Sc.

TECHNICAL EDITING: Sasa Trenčić and Vladan Nikolić

PROOF-READING: Saša Trenčić, MA and Mila Aleksov, BA

ISBN: 978-86-80300-01-6

RAD RAD RAD

СІР - Каталогизација у публикацији - Народна библиотека Србије, Београд

539.16(082)(0.034.2)

INTERNATIONAL Conference on Radiation and Applications in Various Fields of Research (3rd; 2015; Budva)

Proceedings [Flektronski izworl / Third International Conference on

Proceedings [Elektronski izvor] / Third International Conference on Radiation and Applications in Various Fields of Research, RAD 2015, June 8-12, 2015, Budva, Montenegro; [editor Goran Ristić]. - Niš: RAD Association, 2015 (Niš: RAD Association). - 1 elektronski optički disk (CD-ROM); 12 cm

Sistemski zahtevi: Nisu navedeni. - Nasl. sa naslovne strane dokumenta. - Tiraž 400. - Bibliografija uz svaki rad.

ISBN 978-86-80300-01-6

а) Јонизујуће зрачење - Дозиметрија - Зборници

COBISS.SR-ID 215655436



ASSESSMENT OF OCCUPATIONAL DOSE IN FLUOROSCOPY PROCEDURES WHEN INDIVIDUAL MONITORING IS NOT UTILIZED

Olivera Ciraj-Bjelac¹, Danijela Arandjic¹, Predrag Bozovic¹, Sandra Ceklic², Jelena Stankovic¹, Djordje Lazarevic²

¹ School of Electrical Engineering and Vinca Institute of Nuclear Sciences, University of Belgrade, Belgrade, Serbia ²Vinca Institute of Nuclear Sciences, University of Belgrade, Belgrade, Serbia

Fluoroscopy procedures may lead to increased radiation exposure of radiologists and other staff members. The objective of the study is to assess whole body radiation doses and doses to the eye in fluoroscopy procedures, based on measurements that allow for estimates of occupational doses when personal dosimeters have not been used. Four geometrical configurations were considered: overcouch x-ray tube with and without table mounted lead rubber curtains in place, undercouch x-ray tube with horizontal x-ray beam and undercouch x-ray tube with vertical x-ray beam. The doses were estimated using distribution of the scattered radiation and typical workload, as well as the pattern of used of protective tools. Estimated effective dose was in the range from few to 60 μ Sv per procedure for radiologist, and from few to 20 μ Sv per procedure for radiographer, depending on the geometrical configuration and level of personal protection. Corresponding unprotected eye doses were estimated to be in the range 0.03 – 2.8 mSv per procedure for radiologist and 0.02-0.78 mSv for radiographer. The presented results allowed for realistic estimations of the occupational whole body dose and dose to the eyes from the workload of the staff members and from the level of use of radiation protection tools when personal dosimeters have not been regularly used.

Key words: fluoroscopy, radiation dose, eye lens dose, scatter radiation, x-rays

1. Introduction

The use of x-rays for diagnostic or interventional procedures has increased in recent years. Fluoroscopy unit are of particular concern because of potential for delivering relatively high exposures to staff and patients. Among other categories of professional exposure, occupational exposure in fluoroscopy in one of the area in which increased exposure is likely to occur, in particular for certain organs as extremities including eye lens [1-4]. As the eye lens is one of the most radiosensitive tissues in the human body, occupational exposure to radiation has become a real concern for many medical professions, in particular for interventional medical staff since it may lead to the development of lens opacities [5-7].

During these procedures, if staff remain near the xray source and within a high scatter radiation fields for several hours per day, and do not use radiological protection tools and methods, the risk may become substantial [5]. The primary operators who usually stay longer and stand closer to the patient's right side are most likely to receive the highest exposure among all staff involved in interventional procedures [8]. In addition, there have been recent reports of opacities detected in the lens of the eye among some groups of healthcare workers using fluoroscopy in interventional radiology and cardiology [7,9,10]. If these effects are seen in interventional radiology and cardiology, the potential for such effects exists for other uses of fluoroscopy outside imaging departments. Apparently, this appears to be the only circumstance where occupational exposures to ionizing radiation may be routinely resulting in clinically observable tissue reactions [5].

In most clinical situations individual monitoring is performed using a single dosimeter aimed for whole body dose assessment or using two dosemeters positioned under and above the lead apron (double dosimetry). Although double dosimetry is well interventional procedures established in recommended by numerous professional organizations, practical utilization of this approach is not well known mainly in terms of availability and correctness in use of two dosimeters among fluoroscopy staff. In addition, many studies reported irregular use of personal dosimeters [11,12], and even if used, it position on the body remains uncertain. Contrary to whole body dosimetry, eye dosimetry is not well established. Numerous studies have highlighted many drawbacks in monitoring occupational exposure in fluoroscopy including the lack of comprehensive personal dose records [11,13].

Many professionals are still working without using (or not using regularly) any existing simple radiation protection tools (e.g., ceiling suspended screen, glasses) or/and without personal dosemeters [11,13]. In the cases where personal dosimeters have not been used regularly, the approaches taken to estimate retrospectively the dose to staff members often use approach based on the availability of values from experimental measurements made at certain positions at the fluoroscopy room, representative of the typical position of the particular staff members. To improve the accuracy, these measurements are usually made without the protection of the ceiling suspended screen, and the obtained scatter dose values can be corrected by the attenuation factor of the protection tools used [11].

The objective of this study is to assess whole body radiation doses and doses to the eye in fluoroscopy



rooms, based on measurements that allow for estimates of occupational doses when personal dosimeters have not been used.

2. METHODS

Four geometrical configurations were considered: overcouch x-ray tube with and without table mounted lead rubber curtains in place, undercouch x-ray tube with horizontal x-ray beam and undercouch x-ray tube with vertical x-ray beam (Figure 1). Scatter radiation dose rate was measured using a calibrated dosimeter in terms of operational dosimetric quantity H*(10), in the presence of the water phantom 25 cm x 25 cm x 15 cm as a scattering medium using a calibrated pressurized ionizing chamber Victoreen 451 P (Fluke Biomedical, USA). The dose rate measurements was performed in steps of 20 cm from 40 to 200 cm from the floor in vertical direction and from 20 to 100 cm from patient table side in horizontal direction [14]. Therefore, scatter radiation dose was measured in the plane perpendicular and plane parallel to the x-ray bema axis in steps of 20 cm. The measured values were further interpolated to obtain isodose curves.

Two models of conventional x-ray units were used for measurements: Superix 1250, Telestatix (EI Nis) with overcouch x-ray tube geometry and Superix 1250, Undistat (EI Nis) with undercouch x-ray tube geometry.




Fig. 1. Irradiation geometry: undercouch, vertical beam (a), undercouch tube, horizontal beam (b), overcouch tube, no protection (c), overcouch tube with protection (d)

These results were used for assessment of weighted surface dose, i.e. surface dose at position of a particular organ contribution to the effective dose, as this quantity could be used as a conservative estimation of the effective dose. The position of organs was deduced from the MIRD5 mathematical phantom [15, 16]. Information on typical exposure parameters was taken from clinical practice (average:95 kV, 3 mA, 10 min of fluorospy time and 6 exposures/examination of 82 mAs each) and was further used to assess effective dose and eye dose for radiologist and radiographer performing a fluoroscopy procedure. The weighted surface dose was calculated using a H*(10) value at a level of each particular organ and ICRP tissue

weighting factor [15]. This quantity was calculated for a typical position of radiologist and radiographer for each irradiation geometry. Effective dose was further estimated using following equation:

$$\widetilde{E}_k = I \cdot t \cdot H_{k,p} + n_g \cdot mAs \cdot \dot{H}_{k,g} \quad (1)$$

where, t is fluoroscopy time, I is tube current in fluoroscopy mode, n_g is number of images, mAs tube loading per image and $\dot{H}_{k,p}$, $\dot{H}_{k,g}$ are wiighted surface doses for fluorosopy and radiography in $\mu Sv/mAs$, respectively. Typical exposure paarmeter for fluorosopy rpocedures were used to assess the efective dose. Furthermore, different arrangements of personal protective tools were considered, as aprons of 0.25 mm, 0.35 mm and 0.5 mm lead equivalence and combination of lead apron and collar shielding of 0.25 mm lead equivalence. Eye lens dose assessment was assessed based on the information on the typical exposure parameters and scatter dose rate at the level of 160 cm at the position of radiologist and radiographer.

3. RESULTS AND DISCUSSION

The radiation profiles are presented in Figs. 2-5. The essential difference in radiation safety between two types of fluoroscopic unit configurations lies in scattered radiation aspects which has a huge impact on the occupational exposure.




Fig. 2. Vertical scatter radiation profiles for 80~kV, for an overcouch X-ray tube fluoroscopic installation. Pb indicates the presence of protective drapes/curtains




Fig. 3. Vertical scatter radiation profiles for 80 kV, for an undercouch X-ray tube fluoroscopic installation, vertical x-ray beam. Pb indicates the presence of protective drapes/curtains




Fig. 4. Vertical scatter radiation profiles in a vertical plane for an undercouch X-ray tube fluoroscopic installation, horizontal x-ray beam.

Presented figures demonstrate that the protective tools have a considerable effect; with the screens removed the dose profiles are flattened in shape. On removal of these screens the dose rate is increased by more than factor of ten.

The resulting body exposure is extremely nonuniform. Under such conditions the partially unshielded organs in the trunk, together with tissues and organs in the head and neck region influence the effective dose. The eye lens dose is of particular concern. The assessed whole body and eye doses for different x-ray unit geometries and protective tool arrangements are presented in Table 1. Estimated effective dose was in the range from less than 0.01 to 0.06 mSv per procedure for radiologist, and from less than 0.01 to 0.02 mSv per procedure for radiographer, depending on the geometrical configuration and level of personal protection. Corresponding unprotected eye doses were estimated to be in the range 0.03 - 2.8 mSv per procedure for radiologist and 0.02-0.78 mSv for radiographer. The highest exposures are associated with overcouch tube geometry and absence of table mounted protective tools.




Fig. 5. Vertical scatter radiation profiles in a vertical plane for an undercouch X-ray tube fluoroscopic installation, couch in vertical position

Major problems and sources of uncertainty in occupational dose assessment are related to the variations in staff positions and variable technical factors used during fluoroscopy, in particular fluoroscopy time and number of images or series.

It was assumed that the staff wears lead-rubber protective devices of different lead equivalent and transmission factor. Table 1 provides the weighted surface doses for the cases in with no protective clothing and for the for different types of personal protective devices. When applicable, a factor to take account of dose reduction achieved by the lead apron a partial body shielding was applied. In that manner, weighted surface dose contains contributions from organs lying outside the protective apron. The calculations were performed for protective aprons of lead equivalent of 0.25, 0.35 and 0.5 mm and for the combination of apron of lead equivalent of 0.35 mm and thyroid shield of lead equivalent of 0.25 mm.

Similar calculations were performed to estimate eye doses and these values are also given in Table 1.

Table 1 Estimated doses to staff during fluoroscopic procedures

| | Dose per procedure [mSv] | | | | | |
|-------------------------------|--------------------------|---------------------------------|----------------|----------------|----------------|---------------|
| Category of staff | No apron | With apron (mm lead equivalent) | | | | Eye dose |
| | | 0.25 | 0.35 | 0.5 | 0.35 +0.25* | [mSv] |
| Radiologist | | | | | | |
| overcouch X-ray tube unit | 0.28 | 0.06 | 0.05 | 0.04 | 0.02 | 2.84 |
| (without protective curtains) | (0.08-0.67) | (0.02 - 0.13) | (0.02 - 0.11) | (0.01-0.10) | (0.006-0.05) | (0.85-6.67) |
| overcouch X-ray tube unit | 0.08 | 0.02 | 0.01 | 0.01 | <0.01 | 0.20 |
| (with protective curtains) | (0.02 - 0.18) | (0.01-0.05) | (0.003 - 0.02) | (0.002 - 0.02) | | 0.06-0.47) |
| undercouch X-ray tube unit | 0.05 | 0.01 | 0.01 | 0.006 | <0.01 | 0.02 |
| (horizontal) | (0.02 - 0.12) | (0.002 - 0.02) | (0.003 - 0.02) | (0.002 - 0.01) | | (0.005-0.04) |
| undercouch X-ray tube unit | 0.01 | <0.01 | <0.01 | <0.001 | <0.01 | 0.03 |
| (vertical) | (0.002 - 0.02) | | | | | (0.01-0.07) |
| Radiographer | | | | | | |
| overcouch X-ray tube unit | 0.16 | 0.02 | 0.02 | 0.02 | 0.01 | 0.78 |
| (without protective curtains) | (0.05-0.37) | (0.01-0.05) | (0.01-0.04) | (0.01-0.04) | (0.003 - 0.02) | (0.23-1.82) |
| overcouch X-ray tube unit | 0.04 | 0.01 | 0.01 | <0.01 | <0.01 | 0.06 |
| (with protective curtains) | (0.01-0.08) | (0.002 - 0.02) | (0.002 - 0.02) | | | (0.02 - 0.15) |
| undercouch X-ray tube unit | 0.04 | < 0.01 | <0.01 | <0.01 | <0.01 | 0.04 |
| (horizontal) | (0.01-0.10) | \0.01 | \0.01 | | | (0.01-0.10) |
| undercouch X-ray tube unit | <0.01 | <0.01 | <0.01 | <0.01 | <0.01 | 0.02 |
| (vertical) | | | | | | (0.01-0.05) |

The data in the Table 1 emphasize the importance of wearing a lead apron. It is worth mentioning that increasing lead thickness above 0.35 mm dose not deliver much dose reduction, but substantially increase the apron weight. Additional thyroid shielding results in more completed shielding of the body, so the combination of apron of lead equivalent of 0.35 mm and thyroid shielding thickness of 0.25 mm Pb would be very efficient in reduction the effective dose.

Assuming a radiologist performs five procedures daily, five-day week, 48 week-year, a maximum of about 1200 examinations could be performed annually. In practice, in most cases annual workload is only a fraction of this number.

3. Conclusions

Scattered radiation rates in the vicinity of fluoroscopy equipment can be quite high, depending on type of fluoroscopic procedure and equipment configuration, as well as on the use of protective tools. In overcouch-tube fluoroscopy with the radiologist at the couch side, a significant dose to the eye may result from modest number of procedures. Under these circumstances, protective actions need to be taken.

In undercouch-tube fluoroscopy, the estimated dose is well within dose limits. Dose to eyes is, however, limiting factor. Consequently, wearing eye protective glasses should be mandatory.

In addition, it is important to consider personal monitoring arrangements for each type of procedure and equipment type, as the use of individual monitoring is the best way to track occupational exposure of an individual.

Presented methodology could be used as an option when individual monitoring is not available or regular. Using presented data, extrapolation to different diagnostics techniques can be applied, in respect to differences in workload, fluoroscopy time, exposure factors, field sizes and use of protective tools.

The presented results allowed for realistic estimations of the occupational whole body dose and dose to the eyes from the workload of the staff members and from the level of use of radiation protection tools when personal dosimeters have not been regularly used.

Acknowledgement: This work was supported by the Ministry of Education and Science of the Republic of Serbia (grant agreement 43009 and 171007).

REFERENCES

- International Commission on Radiological Protection. Avoidance of radiation injuries from medical interventional procedures. ICRP Publication 85; Annals of the ICRP 30; 2000.
- Faulkner, K., Moores, B. An Assessment of the Radiation Dose Received by Staff using Fluoroscopic Equipment. Br. J. Radiol. 55(1982)272-276.
- 3. Vano E, Kleiman NJ, Duran A, et al. Radiationassociated lens opacities in catheterization personnel: results of a survey and direct assessments. J Vasc Interv Radiol. 2013 Feb;24(2):197-204.
- 4. Vano E, Kleiman NJ, Duran A, et al. Radiation cataract risk in interventional cardiology personnel. Radiat Res. 2010 Oct;174(4):490-5.
- 5. Rehani MM, Ciraj-Bjelac O, Vano E, eta l.. ICRP Publication 117. Radiological protection in fluoroscopically guided procedures performed outside the imaging department. Ann ICRP. 2010 Dec;40(6):1-102
- International Commission on Radiological Protection. Statement on tissue reactions. Approved by the Commission on April 21, 2011. Available at www.icrp.org/docs/ICRP Statement on Tissue Reactions.pdf.
- Ciraj-Bjelac, O, Rehani M, Minamoto A, et al. Radiation-induced eye lens changes and risk for cataract in interventional cardiology. Cardiology. 2012;123(3):168-71.
- Kong I Y, Struelens L, Vanhavere F, et al. Influence of standing positions and beam projections on effective dose and eye lens dose of anaesthetists in interventional procedures. Radiat Prot Dosimetry (2015) 163 (2): 181-187.
- 9. Vano, E, Gonzalez L, Fernandez JM, et al. Eye Lens Exposure to Radiation in Interventional Suites: Caution Is Warranted .Radiology: Volume 248: Number 3—September 2008, 945-53
- Vano E, Gonzalez L, Guibelalde E et alRadiation exposure to medical staff in interventional and cardiac radiology Br J Radiol. 1998 Sep;71(849):954-60
- 11. Vano E, Fernandez JM, Sanchez R, eta l. Realistic approach to estimate lens doses and cataract radiation risk in cardiology when personal dosimeters have not been regularly used. Health Phys. 105(4):330Y339; 2013
- Martin, C J, A Review of Radiology Staff Doses and Dose Monitoring Requirements "Radiation Protection Dosimetry (2009), Vol. 136, No. 3, pp. 140–157
- Ciraj-Bjelac O, Rehani M.M.Eye Dosimetry in Interventional Radiology and Cardiology: Current Challenges and Practical Considerations. Radiat Prot Dosimetry. 2013 Nov 20.
- 14. IEC. Medical electrical equipment. Part 1: General requirements for safety. 3. Collateral standard: General requirements for radiation protection in diagnostic X-ray equipment. IEC 601-1-3 (1994).
- International Commission on Radiological Protection. The 2007 recommendations of the International Commission on Radiological Protection. ICRP publication 103; Ann ICRP 37, 2007.
- 16. Handbook of Anatomical Models for Radiation Dosimetry (Series in Medical Physics and Biomedical Engineering) ed. Xie George Xu and Keith F. Eckerman, Taylor and Francis, 2010.