

**Vaccination with Human Hookworm Vaccine**  
**“*Necator americanus* Aspartic Protease-1 M74”**  
**Generates Neutralizing Antibodies and a**  
**Potent Immune Response in BALB/c Mice**

Amar R. Jariwala, MD, MS

Assistant Research Professor

The George Washington University Medical Center


---

**THE GEORGE  
WASHINGTON  
UNIVERSITY**


---

WASHINGTON, DC


# Life Cycle of *Necator americanus* Human Hookworm



# *Necator americanus* Degradation of Host Hemoglobin and Vaccine Targets




# Method<sup>†</sup> (MOD=Fluorescence)



<sup>†</sup>Reaction Buffer - 50 mM Sodium Acetate (pH = 6.0).

## Steps for Generation and Purification of Polyclonal Mouse IgG<sup>+</sup>



<sup>†</sup>50 BALB/c mice vaccinated twice with 9.33µg Na-APR-1 M74 + 74.64µg Alhydrogel®.


**Neutralizing Capacity of purified IgG with 250 ng of *Na*-APR-1wt and 1μM Cathepsin-D Substrate in 50mM Sodium Acetate, pH 6.0.**

	IgG (μg)	Operators	% Inhibition*	Intra-plate %CV	Inter Operator %CV	All Runs %CV
Run 1	5	1	48.93	4.86	9.81	14.57
		2	42.58	6.98		
Run 2		1	60.44	19.39	10.22	
		2	52.29	7.68		

\*RFU of the purified IgG compared to RFU of Negative IgG (commercial mouse IgG).

Note: All the components were incubated at 37°C for 30 Min. RFU=Relative fluorescence Units.



## Neutralizing Capacity of 8, 4, 2 and 1 $\mu\text{g}$ of IgG with 250ng of *Na*-APR-1wt




Note: All the components were incubated at 37°C for 30 Min.  
RFU of Positive IgG compared to RFU of Negative IgG (commercial mouse IgG).  
RFU=Relative fluorescence Units.




# Standard Calibration Curves (22 curves from 22 Anti *Na*-APR-1 M74 IgG ELISA Plates)



## LOQ and RT



## Standard Reference Serum




9.33 $\mu$ g *Na*-APR-1 M74 + 74.64 $\mu$ g Alhydrogel<sup>®</sup>



# Responders on Study Day 28 Post-Prime Vaccination

Fractional Doses ( $\mu\text{g}$ )		Responders	
		Month	
<i>Na</i> -APR-1 M74	Alhydrogel®	1	7
N/A	400	0	0
50.00	N/A	0	0
50.00	400.00	9	10
28.57	228.56	10	6
16.33	130.64	6	7
9.33	74.64	3	8
5.33	42.64	0	1
3.05	24.40	0	0
1.74	13.92	0	0
0.99	7.92	0	0
ED50 (Effective Dose 50) ( $\mu\text{g}$ )		14.15	11.46
95% Fiducial Limits ( $\mu\text{g}$ )		10.47 -- 18.93	4.86 -- 27.42

# ED50 on Study Day 28 Post-Prime Vaccination



ED50= Antilog(1.06)= 11.46 µg

ED50= Antilog(1.15)= 14.15 µg

Log (Dilution of Vaccine)

Na-APR-1(µg)

Alhydrogel® (µg)

0.99

1.74

3.05

5.33

9.33

16.33

28.57

50

7.96

13.93

24.37

42.65


74.64

130.61

228.57

400

# Relative Potency (RP) of *Na*-APR-1 M74 Vaccine



# Summary

- Five microgram of purified IgG from BALB/c mice vaccinated with 9.33  $\mu\text{g}$  *Na*-APR-1 M74 and 76.64  $\mu\text{g}$  Alhydrogel® neutralized 51.06% of the enzymatic activity of 250 ng of *Na*-APR-1wt
- An excellent dose response (% inhibition vs IgG) was observed
- The standard reference serum generated an excellent standard calibration curves (SCCs) as well as an excellent global standard calibration curve (GSCC) using an Anti-*Na*-APR-1 M74 IgG ELISA
- *Na*-APR-1 M74 vaccine generated a potent immune response in Balb/c as evident by the generation of ED50 at time 1 and 7 month post-manufacture
- *Na*-APR-1 M74 vaccine became 1.23 times more potent at time 7 month when compared to its potency at time 1 month post-manufacture

# Acknowledgments

- This project is supported by the Sabin Vaccine Institute through funding from the Bill and Melinda Gates Foundation and the Dutch Government