

Rochester Institute of Technology

RIT Scholar Works

Theses

3-25-1982

Microprocessor-based in-process infrared densitometer

Steven Cox

Follow this and additional works at: <https://scholarworks.rit.edu/theses>

Recommended Citation

Cox, Steven, "Microprocessor-based in-process infrared densitometer" (1982). Thesis. Rochester Institute of Technology. Accessed from

This Thesis is brought to you for free and open access by RIT Scholar Works. It has been accepted for inclusion in Theses by an authorized administrator of RIT Scholar Works. For more information, please contact ritscholarworks@rit.edu.

MICROPROCESSOR-BASED IN-PROCESS

INFRARED DENSITOMETER

by

Steven P. Cox

B.S.E.E. California State University, Northridge

(1978)

A thesis submitted in partial fulfillment
of the requirements for the degree of
Master of Science in the School of
Photographic Arts and Sciences in the
College of Graphic Arts and Photography
of the Rochester Institute of Technology

March, 1982

Steven P. Cox

Signature of the Author.....

Photographic Science
and Instrumentation

Ronald Francis

Accepted by.....

Coordinator, Graduate Program

School of Photographic Arts and Sciences
Rochester Institute of Technology
Rochester, New York

CERTIFICATE OF APPROVAL

MASTER'S THESIS

The Master's Thesis of Steven P. Cox
has been examined and approved
by the thesis committee as satisfactory
for the thesis requirement for the
Master of Science degree

John F. Carson

.....
Professor John F. Carson, Thesis Advisor

J. S. Wirtz

.....
John S. Wirtz

Burt H. Carroll

.....
Dr. Burt H. Carroll

..25 March 1982.....

Date

MICROPROCESSOR-BASED IN-PROCESS
INFRARED DENSITOMETER

by

Steven P. Cox

Submitted to the Photographic Science and
Instrumentation Division in partial fulfillment
of the requirements for the Master of Science
degree at the Rochester Institute of Technology

ABSTRACT

By the use of high performance solid-state infrared emitting diodes, matched infrared detectors, and a single chip 8-bit microprocessor, an in-process infrared densitometer has been designed and constructed. The device is capable of recording the build-up of optical transmission density on an exposed film sample as it develops, with no effect on the normal development process. The most novel feature of this system is a special development chamber with eleven built-in infrared densitometers. These densitometers are located as to read every other step of a standard Kodak #2 step-tablet exposure, and to allow the passage of developer across the film sample.

As an example, every other step of a 21-step exposure can be measured to produce a series of eleven density

readings in approximately 22 milliseconds. These measurements can be repeated from once every second to once per minute at any time during the development process. The microprocessor that controls the infrared densitometer has memory capacity to store the data required to generate 187 complete D-Log exposure curves. At the completion of processing, the density readings recorded earlier are available for user viewing on a large seven segment digital display system, or can be reproduced on a line printer.

D-Log exposure curves obtained from the device are in good agreement with curves obtained by conventional means. The device has a useful range of 0 to 3.00 density units with an accuracy of $\pm .02$ density units.

ACKNOWLEDGEMENTS

I would like to thank the United States Central Intelligence Agency for their grant that funded the majority of this project. I would also like to thank the Polaroid Corporation for their grant that funded the construction of the development chamber. The contributions to this project by my Department Chairman, Dr. Ronald F. Francis, have not gone unnoticed, and are appreciated also. I extend my sincere appreciation to my thesis advisor John F. Carson and to my committee members John S. Wirtz and Dr. Burt H. Carroll for their time and patience. Lastly, I sincerely thank the entire Engineering Department at the Graphic Equipment Division of the Itek Corporation for their interest and invaluable technical contributions that really made this project possible.

Table of Contents

Introduction.....	1
Prior Work.....	3
General Design Considerations.....	9
IR Densitometer System Description.....	11
Data Acquisition System Analysis.....	18
Experimental Verification.....	29
Experimental Conclusions and Recommendations.....	60
In-Process IR Densitometer Operating Instructions.....	63
Emitter Calibration.....	65
Software Calibration.....	69
Film Sample Preparation.....	74
Data Entry.....	76
Fluid Transport System (FTS).....	79
Prompt Description and Error Messages.....	85
Printer Maintenance.....	93
References.....	96
Appendix 1 : Input/Output Port And Peripheral Addressing.....	98
Appendix 2 : Operating System Program Listings.....	108
Appendix 3 : Raw Data.....	177
Appendix 4 : Infrared Emitting Diode LED55C Specifications.....	190
Appendix 5 : UDT-450 Specifications.....	192

Appendix 6 : Model 757N Logarithmic Ratio Amplifier.....	194
Appendix 7 : SDK-85 Specifications.....	195
Appendix 8 : System Drawings.....	197
Vita.....	226

List of Tables

Table 1 :	Repetition data for Fine Grain Release Positive, type 5302, developed in D-19.....	47
Table 2 :	Repetition data for Fine Grain Release Positive, type 5302, developed in D-76.....	49
Table 3 :	Repetition data for Commercial Film, type 6127, developed in D-76.....	51
Table 4 :	Process data for Fine Grain Release Positive, type 5302, developed in D-19 demonstrating effect of developer heating.....	53
Table 5 :	Process data for Fine Grain Release Positive, type 5302, developed in D-76 demonstrating effect of developer heating.....	54
Table 6 :	Log exposure values used for all experiments.....	62

List of Figures

Figure 1:	Functional block diagram.....	12
Figure 2:	Oblique sketch of the development chamber.....	13
Figure 3:	System block diagram.....	16
Figure 4:	Exploded view of a single densitometer.....	19
Figure 5:	Ideal analog system configuration.....	20
Figure 6:	Non-ideal analog system configuration.....	24
Figure 7:	Output density vs. input density for different levels of amplifier drift.....	27
Figure 8:	Density as a function of Log exposure for Eastman Kodak Fine Grain Release Positive type 5302, developed in D-19.....	31
Figure 9:	Density as a function of Log exposure for Eastman Kodak Fine Grain Release Positive type 5302, developed in D-76.....	32
Figure 10:	Density as a function of Log exposure for Eastman Kodak Commercial Film, type 6127, developed in D-76.....	33
Figure 11:	Unfixed, wet IR density vs. fixed-out, dry diffuse density for Fine Grain Release Positive, type 5302, in D-19.....	36
Figure 12:	Unfixed, wet IR density vs. fixed-out, dry diffuse density for Fine Grain Release Positive, type 5302, in D-76.....	37
Figure 13:	Unfixed, wet IR density vs. fixed-out, dry diffuse density for Commercial Film, type 6127, in D-76.....	38
Figure 14:	Comparison of D-Log H curves for wet unfixed, and dry fixed-out Fine Grain Release Positive, type 5302, in D-19.....	42

Figure 15: Comparison of D-Log H curves for wet Unfixed, and dry fixed-out Fine Grain Release Positive, type 5302, in D-76.....	43
Figure 16: Comparison of D-Log H curves for wet unfixed, and dry fixed-out Commercial Film, type 6127, in D-76.....	44
Figure 17: Density as a function of position and flowrate for a sample of Pan-X Recording Film, SO-164, exposed to room lights for 5 minutes and developed in D-76.....	58
Figure 18: Density as a function of position and flow- rate for a uniformly exposed sample of Pan-X Recording Film, SO-164, developed in D-76.....	59
Figure 19: Software calibration plot.....	71
Figure 20: Flowrate as a Function of Pump Setting.....	80
Figure 21: Fluid Transport System Hydraulic Schematic....	83
Figure 22: Power Output vs. Input Current for LED55C....	191
Figure 23: Typical Radiation Pattern for LED55C.....	191
Figure 24: Spectral Response of UDT-450.....	193
Figure 25: Output voltage of UDT-450 as a function of incident energy.....	193
Figure 26: Logarithmic Amplifier Schematic.....	197
Figure 27: Keyboard Schematic.....	198
Figure 28: FTS Manual Override Schematic.....	199
Figure 29: Analog to Digital Converter Assembly.....	200
Figure 30: Analog to Digital Converter Schematic.....	201
Figure 31: 12 Character Display Assembly.....	202
Figure 32: 12 Character Display Schematic.....	203
Figure 33: Opto Isolated I/O Assembly.....	204
Figure 34: Opto Isolated I/O Schematic.....	205

Figure 35: 33 Digit Display Assembly.....	206
Figure 36: 33 Digit Display Schematic.....	207
Figure 37: IRED Driver Board Assembly.....	208
Figure 38: IRED Driver Board Schematic.....	209
Figure 39: Driver Board Assembly.....	210
Figure 40: Driver Board Schematic.....	211
Figure 41: Chamber Front.....	212
Figure 42: Chamber Back.....	213
Figure 43: Lock.....	214
Figure 44: Hinge.....	215
Figure 45: Chamber Mount.....	216
Figure 46: Cover Plate.....	217
Figure 47: Lens Tube.....	218
Figure 48: Cable Support.....	219
Figure 49: Cable Support.....	220
Figure 50: Lock Pin.....	221
Figure 51: Input Pipe.....	222
Figure 52: Output Pipe.....	223
Figure 53: Hinge Shaft.....	224
Figure 54: Registration Pin.....	225

INTRODUCTION

The design and construction of a device that would allow for the procurement of densitometric data from a sample of photographic film as it developed has been the subject of numerous thesis projects at the Rochester Institute of Technology. Some of these devices have been constructed, and feasibility proven, but their usefulness was extremely limited due to their primitive nature. The purpose of this thesis is to tie all previous efforts together in the design, construction, and testing of a completely new, workable, and reliable in-process IR densitometer.

The conceptual requirements for IR densitometry have been well known for years, but prior workers have either lacked the technology, or the expertise to construct a truly significant device. This project makes use of state-of-the-art technology along with industry standard development and assembly techniques in its design and construction. Therefore, the general aim of this project is to produce a research instrument of the highest possible caliber to allow those qualified to study photographic processing.

Of the many possible uses for this device, induction period studies have received foremost consideration in its design. The ability to monitor changes in the rate of development and to know exactly when the developer makes actual contact with the emulsion are absolute necessities when attempting to determine the length of the induction period. This is where the conventional method of arresting development fails. The

procedure becomes inaccurate at short development times. The IR densitometer, however, has been designed to monitor density build up on a second-by-second basis, thus providing an accurate record of exactly how the film sample developed.

Prior Work:

Conceptually, the design of an in-process IR densitometer is quite simple. All that is really needed is:

- 1.) A source of infrared radiation to which the film being investigated is insensitive.
- 2.) Some method of directing this radiation through the developing film sample.
- 3.) The collection and detection of the radiation after it has passed through the sample.
- 4.) A method of relating the detected signal to some type of optical transmission density.
- 5.) An output system which will deliver the desired sensitometric data to the operator in some convenient form.

However simple this may sound, it has taken almost fifteen years and the contributions of many individuals to reach the point where all prior efforts could be concentrated into the construction of a truly significant and useful in-process IR densitometer.

One of the earliest uses of infrared densitometry was made in 1953 by Fortmiller and James¹. They used a device that would be the predecessor for future designs at the Rochester Institute of Technology. Their apparatus was used to study the kinetics of development of a fine grain motion picture positive emulsion by vanadous ion.

The first attempt to actually build such a device at R.I.T was made by Hughes in 1964². As it was only a B.S. thesis, the project was somewhat limited in scope, but the conceptual feasibility of in-process IR densitometry was shown.

Hughes IR source consisted of a tungsten display case lamp filtered with a Wratten No. 70 filter. The radiation was directed onto a developing sensitometric exposure fastened to the bottom of a plastic tray. From the other side of the tray, 84 plastic rods piped the IR radiation to a large collar with a single L-shaped light pipe mounted in the center. This L-shaped rod then had a silicon photovoltaic cell located on axis so as the rod was rotated, it could scan all 84 pipes and direct their output onto the single photocell. The pipes were scanned by the rotating rod at such a rate as to produce 40 characteristic curves per second. The output of this cell was then logarithmically amplified and displayed on an oscilloscope screen. The screen was photographed by a single frame motion picture camera and the resulting images could be projected onto an empirically determined transfer curve relating wet IR density to dry white light density.

The problems encountered in this system were numerous, the major ones being non-uniformity in transmission between the light pipes and poor photocell sensitivity at high densities. However, the five necessary sub-systems called for in the earlier paragraph concerning conceptual simplicity were present.

- 1.) IR source: tungsten lamp and Wratten filter.
- 2.) Radiation through-put: plastic light pipes.
- 3.) Detection: rotating silicon photocell.
- 4.) Relation of system output to density: log amp and transfer curve.
- 5.) Display: oscilloscope, movie camera, and projector combination.

In 1970, a second endeavor was made to construct an IR densitometer by Hisler and Casinelli³. The five major sub-systems of their design can be outlined as follows:

- 1.) IR source: one six-volt, fifty watt tungsten lamp filtered by three Kodak No. 29 Wratten filters.
- 2.) Radiation through-put: single emitter/detector combination. Developing sensi-strip was mechanically moved through IR beam by being mounted inside a revolving plastic cylindrical tank. Method of rotation was a phonographic turntable which would produce one complete scan of the sensi-strip every .77 second.
- 3.) Detection: single photomultiplier tube from a Macbeth TD-102 densitometer with an S-4 response.
- 4.) Relation of system output to density: output of TD-102 fed directly to oscilloscope, logging and scaling performed by TD-102 electronics.
- 5.) Display: oscilloscope screen with copies being made by a Polaroid scope camera. Type 107 film was used.

The major problem with this system was the mis-match in peak spectral responses of the source/detector combination. Add to this the necessity of three Wratten No. 29 filters to prevent film fogging and it is understandable why no useable output was obtained.

A third series of design improvements were proposed and implemented in 1971 by Beaupre and Jasper⁴. Their system consisted of:

- 1.) IR source: a single 120-volt, 500 watt Sylvania EHA tungsten-halogen projector lamp filtered by Kodak Wratten filters Nos. 87 and 87C.
- 2.) Radiation through-put: same rotating turntable arrangement used by Hisler and Casinelli.
- 3.) Detection: Hewlett-Packard PIN photodiode.
- 4.) Density transformation was performed by a solid state logarithmic amplifier.
- 5.) Display was same as that used by Hisler and Casinelli.

that is, an oscilloscope and Polaroid scope camera.

This device was tested with Eastman Kodak Fine Grain Release Positive Film, type 5302, developed in DK-50 and D-76. The results of their experiment demonstrated good agreement between conventionally obtained characteristic curves and their IR obtained ones. However, beyond a density of 2.5, curve deviation due to non-linearities in their amplifier became significant, thereby limiting the useful range of the device.

The first attempt to use the above device as an actual research tool, and not a single project in itself, was made by Turbide and Williams in 1972⁵. They investigated the phenomenon of lith or infectious development. Film type used was again 5302 and ortho Kodalith type 3, 2556. No modifications were effected upon the device and the only significant contribution of this endeavor was to simply repeat those results obtained by Beaupre and Jasper.

In an attempt to increase the useful density range of the instrument, Turbide and Williams tried removing the anti-halation backing of their sensi-strips before development. But even with this, the device gave linear results only out to a density of about 2.4.

A more comprehensive write-up and representative curves are contained in an article written by B.H. Carroll⁶ summarizing the work of Jasper, Turbide, and Williams. Technically, the paper deals little with the construction and operation of the densitometer, but rather with the subject for which the project was performed; photographic chemistry.

Industry references to IR densitometry are scarce, but two major ones will be mentioned here in the hopes of demonstrating that the concept is not merely an academic

curiosity. In 1966, Spitzak⁷ made use of solid state IR emitters and silicon photodetectors in a film scanning system. The ultimate purpose of the system was to scan photographs of planets obtained by spacecraft, digitize the signal, and transmit them to earth. As a side note, the author mentions the possibility of using such a system to perform density measurements before fixing the film to determine the degree of development.

More in keeping with the scope of this project, Keemink and Van der Wildt⁸ designed and built a device which they called a gammascope. The main intention of this project was to monitor the density build-up of a piece of film while it was developing, and then to arrest the process when a certain value of gamma was reached.

The conceptual layout of this system is very similar to that under consideration in this paper. An array of discrete IREDs (infrared emitting diodes) with matching discrete photocells situated such that twelve positions of a continuous density wedge (maximum density of 2.00), could be scanned during development. The output of the photocells is then electronically transformed into density values by a logarithmic amplifier. From here, the materials characteristic curve is displayed upon an oscilloscope screen as a series of discrete "dots". Gamma could then be determined by fitting a line through the linear portion of the curve and calculating its slope.

It is within this body of work where the first mention of problems involving the Herschel effect is encountered. Decreases in density were stated to occur when repeated scans were made in the earliest states of development. It is at this point when latent image centers are not yet rendered fully developable and are most susceptible to deterioration by IR radiation. This condition can be further aggravated by the use of very dilute developer solutions. The authors suggest that

this problem can be avoided by performing no scans during the first minutes of development. However, one of the most important design aspects of an in-process IR densitometer is to gather sensiometric data during precisely the first moments of development in order to study induction effects.

The most recent work in the field of IR densitometry was performed in 1976 by Piskacek⁹. His M.S. thesis was a natural progression of previous work leading up to a final proof of in-process IR densitometry feasibility.

The body of his project dealt with the design of a special chamber which would allow a sensiometrically exposed piece of 35mm film to be scanned by IR radiation while it developed. His specifications called for an array of eleven solid state IR emitter/detector combinations to be mounted within the chamber to facilitate density measurements being taken from a 21-step sensiometric exposure.

Not only was the chamber designed, but the materials called out by his plans were tested for their suitability in a developer environment (e.g. extreme excursions in pH). Also, any material which was to be used in the transmission of IR was tested for its optical properties and possible attenuation in the IR.

Although Piskacek was unable to actually build his chamber, he was able to bench test a small portion of the electro-optical system. Using an IR emitting diode, (Monsanto ME5) and a solid state detector, (UDT-500, effective area .05 cm²), he was able to obtain a linear correlation between unfixed IR density and fixed out white light density. Agreement was good out to a maximum density of 2.00.

General Design Considerations:

Piskacek's work in 1976 was one of the last design steps necessary for the construction of a truly significant and useful in-process IR densitometer. Very little was said, however, about the design of a working densitometer from a systems approach. Piskacek used a simple block diagram point of view for detailing how he thought the remainder of the densitometer should be configured. In all actuality, he really intended for the remaining work to be carried out by someone else, and purposely left any electronic designs vague.

Picking up where Piskacek left off, a complete IR densitometer system can be generalized. His selection of IR emitters and detectors mounted in a specially designed development chamber would produce eleven discrete optical transmission signals. Some of the most important design considerations for this system are that these signals be generated quickly (on the order of milliseconds), accurately (within $\pm .02$ density units), and over a density range of 0 to 3.00. Another important consideration is an efficient and convenient method of storage and display of density data. The old method of photographing an oscilloscope screen would not be suitable for the system being considered here.

To meet the above generalized requirements, a single chip 8-bit microprocessor has been selected as the system controller. Microprocessors are ideally suited for the rapid control, storage, and display of data. For this particular system, the microprocessor will control the sequential pulsing of the eleven emitter/detector pairs to generate eleven optical transmission signals. By the use of an analog multiplexer, it will direct these signals to a high-performance logarithmic amplifier, and then put them into digital form for input to the microprocessor. The output system will consist of a light emitting diode display

of the actual density values, or permanent records of density information can be generated on a line printer.

The overall design philosophy is to construct a device that can be used as a research instrument by scientific personnel.

IR Densitometer System Description:

The design of the in-process IR densitometer can be broken up into three major functional groups. These three groups are: data acquisition, data management, and data output. Although all groups are under total control of the microprocessor, this functional grouping will aid in understanding the system as a whole (see figure 1).

Data Acquisition:

The major component of this sub-system is the development chamber. The chamber is a precisely machined block of stainless steel designed to hold a strip of 35mm film in registration against eleven solid state IR densitometers while developer is pumped across the sample. The chamber is made up of two matched sides that mate to form the development cavity. It is connected at the bottom by a hinge pin (see figure 2) and sealed at the top by a locking pin. The top section contains the eleven IR emitting diodes (General Electric LED55C, see appendix 4 for complete specifications). Each diode is mounted in a small tube with a lens at the end ($f=5\text{mm}$). This assembly is mounted in a guide hole behind a plastic window and focused onto its companion detector below. The bottom section contains the eleven IR detectors. They too are mounted behind a plastic window and when the chamber is sealed, are in perfect registration with the emitters on top and create eleven discrete, solid state mini IR densitometers. The detectors (UDT-450's, a photodiode-operational amplifier combination, see appendix 5) produce an optical transmission signal directly. These eleven outputs are fed into an analog 16-channel multiplexer. A multiplexer is simply a 16-position switch with one common output. The switch selected, however, depends on the state of four address lines into the multiplexer. In this way, the output of an individual detector




Figure 1: Functional block diagram




Figure 2: Oblique sketch of the development chamber

can be selected and processed within the array of eleven. Due to this one-at-a-time selection of the detectors, the term "array scan" becomes useful in describing how the eleven emitter/detector pairs are sequentially pulsed to make a complete reading of the optical transmission density (commonly called "density") of the developing film strip. To take a density measurement of the film sample, the first emitter/detector pair is switched on, a reading made, switched off, the second pair turned on, etc. Hereafter, a complete set of density measurements for a film sample will be called a "scan". A scan produces eleven density readings and takes approximately twenty-two milliseconds to perform. Keep in mind that all scan timing and signal processing is performed under the direction of the microprocessor.

From the multiplexer, the analog transmission signal is fed into a logarithmic amplifier (see appendix 6). Here the log of the input voltage is taken to produce an analog transmission density signal. This log amplifier has been scaled to produce a two volt change in output for a ten times change on the input. In other words, a sample density of 1.00 produces an output voltage of 2.00 volts, a density of 2.00 gives an output of 4.00 volts, etc. This signal becomes the input to a 10-bit analog to digital converter. The output of this device is a 10-bit digital density signal, which is in a form ready for input to the microprocessor.

These components make up the data acquisition group. Note that any component that operates on analog signals is considered part of the data acquisition system. Once the signal becomes a 10-bit digital word, it enters the digital domain and remains there until final output.

Data Management:

The data management sub-system is made up entirely by the Intel SDK-85 single board computer (see appendix 7). This board contains the Central Processing Unit (CPU) or what is better known as the microprocessor. It also contains all memory elements, input/output devices, timebases, and display controllers. A complete description of how this computer works is contained in the SDK-85 System Design Kit User's Manual #9800451B.

The reason the IR densitometer was designed to be controlled by a computer is twofold. Of the many problems encountered in earlier IR densitometer designs, accurate timing of the density scans and a convenient means of data storage and display were two of the most difficult to solve. These are areas, however, where microprocessors excel. They operate on a very accurate crystal time base, and are able to store large amounts of data very quickly.

Once the 10-bit digital density word leaves the A/D converter, it is brought into the SDK-85 bus structure through two 8-bit I/O ports (see figure 3). The system bus is the means by which data is transported within the computer itself. This data is then immediately stored in random access memory (RAM) for later processing. As each scan produces a burst of eleven density readings, this data is stored as eleven 10-bit word blocks. From RAM the data is processed further (i.e. scaled in accordance with certain calibration procedures) and becomes ready for output.

Data Output:

For computer data to be of any use, there must be some convenient method to present this data to the user in the real world. For the IR densitometer, the user has three methods of




Figure 3: System block diagram

data output. The first output system is for the immediate viewing of density data. The computer can be directed to display any density scan on a 33-digit light emitting diode (LED) display. The system will display all eleven density readings from a single scan. The format for a single reading is X.YZ. That is, each reading is accurate to three significant figures, i.e., with a resolution of $\pm .01$ optical density units.

The second method of data output is through a 40-column line printer. Any set of density readings on the LED display can be reproduced on the line printer. This is to allow the user to generate a permanent record of any particular density scan. For a complete description of these two output modes, please consult the In-Process IR Densitometer User's Manual*.

A third method of data output has been partially provided. Connectors have been installed to allow for the connection of this computer to a higher level computer. With this type of inter-connection, data could be transferred to the higher level system for even more processing of density data. As an example, computer graphics could be used to create D-Log E curves on a CRT screen. Completion of this interface depends on what type of higher level computer is selected, and therefore has been left as a separate project.

* The IR densitometer user's manual begins on page 63.

Data Acquisition System Analysis:

Information within the IR densitometer can be either in analog form or digital form. An analog signal is a voltage that can vary infinitely over a given range. A digital signal is a voltage level that can vary only in discrete steps. This particular digital system is binary; the voltage signal has only two levels, low and high ("0" and "1"). When working with digital systems, very strict and well known rules of logic apply. For analog systems, however, events are not so well defined. It is the intent of this section to provide a set of equations that will describe the operation of the analog system within the IR densitometer.

Figure 4 shows the physical relationship of the components that make up a single densitometer. Figure 5 gives a schematic representation of the same densitometer, along with the other components that make up the analog system. The figure presents an ideal system configuration where all devices are considered to be perfect and introduce no errors.

The IR emitting section consists of the IRED, lens tube, and lens which is used to focus the output of the IRED onto the detector surface. The power per unit area incident on this surface is the irradiance and is expressed as H watt/cm². The detector surface has an active area of A cm² and a responsivity of R amp/watt. The peak output of the IRED and peak responsivity of the detector are so well matched that no wavelength term need be included in the calculations that follow.

With these terms now defined, the current generated by the photo-diode is given by:

$$I_T = (T)(H \text{ watt/cm}^2)(A \text{ cm}^2)(R \text{ amp/watt})$$




FIGURE 4: EXPLODED VIEW OF A SINGLE DENSITOMETER




FIGURE 5: IDEAL ANALOG SYSTEM CONFIGURATION

The output voltage of the operational amplifier is:

$$V_T = R_F \cdot I_T$$

where R_F is the feedback resistance in ohms.

Now the sample transmission T , can also be expressed as a density:

$$T = 10^{-D}$$

Substituting:

$$V_T = H \cdot R \cdot A \cdot R_F \cdot 10^{-D}$$

The irradiance, responsivity, active area, and feedback resistor all remain constants, and can be collected into the term C .

$$V_T = C \cdot 10^{-D}$$

Neglecting any effects from the analog multiplexer, this V_T is input to the logarithmic amplifier. The transfer equation for the log amplifier is:

$$V_D = -K \cdot \text{Log}(I_S/I_R)$$

where I_S is the signal current and I_R is a reference current. To place the expression in voltage terms:

$$I_S = V_T/R_S \quad \text{and} \quad I_R = V_R/R_R$$

Substituting:

$$V_D = -K \cdot \text{Log}(V_T/V_R \cdot R_R/R_S)$$

V_T is known and can be substituted:

$$V_D = -K \cdot \text{Log}(C \cdot 10^{-D} / V_R \cdot R_R / R_S)$$

V_R is a precision reference, and therefore a constant.
 R_R and R_S are resistors; also constants and can be removed.

Let:

$$M = C \cdot R_R / V_R \cdot R_S$$

And:

$$V_D = -K \cdot \text{Log}(M \cdot 10^{-D})$$

$$V_D = K \cdot D - K \cdot \text{Log}(M)$$

Notice that $-K \cdot \text{Log}(M)$ is just another constant, therefore combine it form N .

Lastly:

$$V_D = K \cdot D - N$$

This analysis shows that the output voltage of the log amplifier is simply a scaled version of the film sample density, plus some additive constant, N . For this particular system, K was made to equal 2.0. The constant N presents no problem as it can be nulled nearly to zero by an external offset adjust on the Log amp. Any remaining offset error can then be removed by software calibration routines. Therefore:

$$V_D = 2 \cdot D - N$$

The log amplifier is the last element in the analog system. From here, V_D is converted to a 10-bit digital word and then enters the digital system.

In reality, all the components of the analog system are non-ideal and each has a certain amount of error associated with it. The most critical types of errors are voltage offsets and non-linearities. Within the entire system, only two devices have these problems to any appreciable degree. They are the operational amplifiers and the logarithmic amplifier. The remaining devices introduce errors so small as to be negligible, or they cannot be characterized. Figure 6 again shows the analog system configuration, but this time the two major sources of error are shown as offsets D_1 and D_2 .

An ideal op-amp would have an output voltage of zero volts when the photo-diode connected to it was in complete darkness. This, however, is not the case. Due to temperature drift, noise pick-up and other effects, a trim pot has to be connected to each amplifier to zero out these offsets. No offset adjustment is perfect, and the remaining offset is called D_1 .

The output of the op-amp now becomes:

$$V_T = R_F \cdot I_T + D_1$$

The operation of the logarithmic amplifier is also non-ideal and has two types of errors associated with it. The first is an offset D_2 which is present for similar reasons as in the op-amp. The other type of error is something called Log Conformity Error (LCE). LCE is the error expressed in the output's ability to perform an accurate logging function of the input. Hence the LCE is always expressed Relative To Input (RTI). For the Analog Devices log ratio amplifier:




FIGURE 6 : NON-IDEAL ANALOG SYSTEM CONFIGURATION

$$\text{LCE} = \pm 1\% \text{ RTI}$$

These errors manifest themselves in the transfer equation for the device as:

$$V_D = -K \cdot \text{Log}(1.01 \cdot (V_T/V_R \cdot R_R/R_S)) + D_2$$

Substituting in for V_T :

$$\begin{aligned} V_D &= -K \cdot \text{Log}(1.01 \cdot ((C \cdot 10^{-D} + D_1)/V_R \cdot R_R/R_S)) + D_2 \\ &= -K \cdot \text{Log}(1.01 \cdot C \cdot R_R \cdot 10^{-D}/V_R \cdot R_S + 1.01 \cdot R_R \cdot D_1/V_R \cdot R_S) + D_2 \\ &= -K \cdot \text{Log}(C \cdot R_R \cdot 10^{-D}/V_R \cdot R_S \cdot (1.01 + 1.01 \cdot D_1/C \cdot 10^{-D})) + D_2 \end{aligned}$$

Since:

$$M = C \cdot R_R/V_R \cdot R_S$$

Then:

$$\begin{aligned} V_D &= -K \cdot \text{Log}(M \cdot 10^{-D} \cdot (1.01 + 1.01 \cdot D_1/C \cdot 10^{-D})) + D_2 \\ &= -K \cdot \text{Log}(M \cdot 10^{-D}) - K \cdot \text{Log}(1.01 + 1.01 \cdot D_1/C \cdot 10^{-D}) + D_2 \end{aligned}$$

From ideal case :

$$V_D = -K \cdot \text{Log}(M \cdot 10^{-D})$$

$$= K \cdot D - N$$

Therefore, in the general case:

$$V_D = K \cdot D - K \cdot \text{Log}(1.01 + 1.01 \cdot D_1 / C \cdot 10^{-D}) + D_2 - N$$

Ideal
Term

Error
Term

The output voltage of the log amplifier can now be expressed as the sum of an ideal term and an error term. Note that the error term is not a constant but rather a function of the film sample density. Note also that if the LCE were 1.00 (ideal) and the op-amp offset D_1 were zero, the error term would vanish and only D_2 would remain.

This error analysis details which sources of error are important and which ones are not. The offset errors D_2 and N are simple additive constants and can be easily subtracted out later by a computer software routine. The LCE and the offset D_1 , however, are not insignificant. They produce a logarithmic error which is a function of density. It is extremely difficult to perform logarithmic functions with assembly level software, therefore, these errors cannot be removed by any calibration routine. What must be done is to minimize the effects of these errors. Log Conformity Error is fixed. It will always be present and cannot be removed. On the other hand, D_1 can be made to be very close to zero. Any error introduced by the offset D_1 becomes appreciable only at higher densities. For this particular system, D_1 can be kept to about + 3 millivolts. This along with the LCE translates into an error of about -.06 density units at a sample density of 3.00. This is indeed acceptable, and the error is almost non-existent at lower densities.

Figure 7 is a graphical representation of how detector offset errors affect density measurements. Input density is to be considered actual sample density, while output density is the




Figure 7 : Output density vs. input density for different levels of amplifier drift.

density reported by the system with the above transfer equation. The numbers used in that equation for the generation of figure 7 are based on actual constants used in the physical system. Only positive offset errors are shown in the figure as negative offsets have a completely different effect on the system. This is because the physical log amplifier produces undefined output voltages for negative inputs. What really tends to happen is the amplifier output will saturate (go to maximum output voltage, which is about 14.3 volts). If the analog to digital converter were able to resolve this voltage (10.00 volts is it's maximum), the computer would think a sample density of 7.0 were being read. Therefore, to prevent such events from happening, the offset null trim pots for the eleven detectors are purposely adjusted to produce a positive offset of a few millivolts at all times. This precludes the possibility that the detectors might drift and produce a negative output. This method, of course, introduces some slight error, but to have a density reading that is somewhat lower than reality, is considered far better than a reading of 3.0 being reported as 7.0.

Experimental Verification:

The purpose of this section is to demonstrate two of the most important operational capabilities of the IR densitometer. They are:

- 1.) The ability of the device to relate wet, un-fixed IR density to dry, fixed-out diffuse density in a linear fashion. Also that for the same development time, wet, un-fixed IR D-Log exposure curves are in good agreement with dry, fixed-out diffuse D-Log exposure curves.
- 2.) Repeatability. Most data has been replicated three times to demonstrate the repeatability of the IR densitometer. This data has been tabulated for ease of viewing (see tables 1,2, and 3).

Three series of experiments have been run to obtain sufficient data for conclusions to be drawn about the above two statements. The two film types used were selected because of their relatively simple emulsions. It was also desired to run the experiments with a fairly active developer, and a much milder developer, hence D-19 and D-76 were selected. Details on how these experiments were run are listed below.

Series 1:

Eastman Kodak Fine Grain Release Positive, type 5302, was processed in D-19 at 70° F for 5 minutes. Sensitometry consisted of film samples being exposed in a Kodak model 101 sensitometer. A Kodak #5 step tablet was used along with a 1.20 Inconel N.D. filter. The sensitometer provided 340 lux·seconds at the wedge. Three runs were made with fresh developer and rinse water. To prevent contamination of the Fluid Transport System, no other chemistry was introduced into the system. Scans

were taken once per second for the first minute, then once every 10 seconds for the remaining 4 minutes. No scans were taken during the one minute water rinse that followed development.

All processing was done under Kodak 1A safelight conditions. After the water wash, film samples were tray-fixed for five minutes then washed for twenty minutes and dried.

Results from the first run of this series are shown in figure 8. Wet, un-fixed IR density is plotted as a function of step # (actual Log exposure values can be found at the end of this section). Curves obtained after 30 seconds, 2 minutes and 5 minutes of processing are shown, with the 5 minute curve being the last one obtained before the water wash.

The curves show just about what one would expect for this type of developer/film combination. After about 2 minutes of development, there is no further increase in gamma. At about the same time, toe densities, as well as the over-all sample density, begin to pick-up dramatically. For this sample as with all the others, the IR densitometer would report an initial density of approximately .18 a few seconds after developer contact. This increase in density is due to the fact that the specular optical density of AgBr microcrystals embedded in gelatin is strongly dependent on the water content of the emulsion¹⁰. Also, the optical density does not depend directly on the thickness of the gelatin, but rather on its index of refraction. As water is absorbed by the emulsion, the index of refraction will increase, causing an increase in scattering and therefore, and overall increase in density will be observed.

Figure 11 shows unfixed, wet IR density plotted against fixed-out dry visual diffuse density. All dry density readings are made on a Macbeth TD-504 digital densitometer. Some of the data points are rather far apart because of the high contrast for




Figure 8 : Density as a function of Log exposure for Eastman Kodak Fine Grain Release Positive, type 5302, developed in D-19.




Figure 9 : Density as a function of Log exposure for Eastman Kodak Fine Grain Release Positive, type 5302, developed in D-76.




Figure 10 : Density as a function of Log exposure for Eastman Kodak Commercial Film, type 6127, developed in D-76.

this process, but the linearity between the two types of densities is clearly shown.

The last graph for this series, figure 14, shows both wet and dry D-Log exposure curves for a single film sample after 5 minutes of development. Agreement between the two curves is good with the largest error occurring in the base plus fog portion of the curves. Whenever wet IR density and dry diffuse density are plotted together like this, the wet IR density will always be somewhat higher in the toe region. This again is probably due to the increase in scattering from water absorption combined with the fact that the device is calibrated with a dry film base sample (refer to the Software Calibration section contained in the operating manual for details). Another interesting point is that the two curves cross one another close to a density of .70. This is the density of the mid-point calibration sample used at the start of the process. At and around densities of .70, the calibration routine used by the IR densitometer forces differences nearly to zero.

Series 2:

Eastman Kodak Fine Grain Release Positive, type 5302, was processed in D-76 at 70° F for 5 minutes. Sensitometry consisted of film samples being exposed in a Kodak model 101 sensitometer. A Kodak #5 step tablet was used along with a 0.78 Inconel N.D. filter. The sensitometer provided 340 lux·seconds at the wedge. Three runs were made with fresh developer and rinse water. Scans were taken once per second for the first minute, then once every 10 seconds for the remaining 4 minutes. No scans were taken during the one minute rinse that followed development.

All processing was done under Kodak 1A safelight conditions. After the water wash, film samples were tray-fixed

for five minutes then washed for twenty minutes and dried.

Results from the first run of this second series are shown in figure 9. Wet, un-fixed IR density is plotted as a function of step # (actual Log exposure values can be found at the end of this section). Curves after 30 seconds, 1 minute 30 seconds, and 5 minutes of processing are shown, with the 5 minute curve being the last one obtained before the water wash.

The IR densitometer reported an initial density of approximately .18 due to the usual emulsion swelling from developer absorption. This is not considered to be significantly different from what was observed in series 1 where D-19 was used. These curves appear typical, but once again something unusual is happening in the toe region. The curves show a decrease in base plus fog density with increased development time. This decrease in toe densities is probably due to the slight fixation of the film sample by the developer. D-76 has a relatively high sulfite content (100 g/l)¹¹, and sulfite is a well-known silver-halide solvent.

Figure 12 shows un-fixed, wet IR density plotted against fixed-out dry diffuse density. These data were obtained after five minutes of development, which is always the last scan taken. More data points are present in this plot because the contrast for this process was lower. There also seems to be a slight S-shape to this plot, but a straight line fit still works well if errors less than $\pm .03$ D are acceptable. Straight lines have been fitted to the three graphs of this type. If the plot is extended, it will intercept the x-axis at approximately .1 dry diffuse density units. The reason for this probably lies with the calibration routine used, and the fact that zero density is derived from a dry film sample.




Figure 11 : Unfixed, wet IR density vs. fixed out, dry diffuse density for Fine Grain Release Positive in D-19.




Figure 12 : Unfixed, wet IR density vs. fixed out, dry diffuse density for Fine Grain Release Positive in D-76.




Figure 13 : Unfixed, wet IR density vs. fixed out, dry diffuse density for Commercial Film, type 6127, in D-76.

The last graph for this series, figure 15, shows both wet and dry D-Log exposure curves for the same film sample after 5 minutes of development. This time the curves cross around .80, but the reason is still probably due to the calibration routine. For the toe region of the curve, the comments made about figure 14 apply here also.

Series 3:

Eastman Kodak Commercial Film, type 6127, was processed in D-76 at 70° F for five minutes. Sensitometry consisted of film samples being exposed in a Kodak model 101 sensitometer. A Kodak #5 step tablet was used along with a 1.20 Inconel N.D. filter. The sensitometer provided 340 lux·seconds at the wedge. Three runs were made with fresh developer and rinse. Scans were programmed just as they were in the first two series. Darkroom conditions and fixing procedures were the same as used in the first two series.

Results from the first run of this third series are shown in figure 10. Wet, un-fixed IR density is plotted as a function of step number (actual Log exposure values for each step number can be found at the end of this section). Curves after 50 seconds, 2 minutes, and 5 minutes of processing are shown, with the 5 minute curve being the last one obtained before the water wash.

The IR densitometer reported an initial density of approximately .19 due to the usual absorption of developer. This "wet emulsion density" was noted to be just about the same for all three series run. The amount of time after the start of processing it took to reach this density was about the same for all three series also (5 to 9 seconds).

The curves in figure 10 show very little increase or

decrease in base plus fog density with increased development time. It might help to explain this by breaking the development of the toe region into two competing processes. The first process is the reduction of silver-halide to silver metal. This, of course, produces density, and may be brought about by the development of latent image centers, or the random development of silver halide crystals to produce fog. The second process is the slight fixation of the film by the high sulfite developer. This would cause a reduction in density as the effect of scattering would decrease.

For the film/developer combinations considered here, these two competing processes have various outcomes. The first was Fine Grain Release Positive developed in D-19. Toe densities increased with time because D-19 is an active developer and fog production outweighed fixation. The second combination was Fine Grain Release Positive in D-76. D-76 is a milder developer than D-19 and fixation proceeded faster than fog and an overall decrease in toe densities was observed. The last combination was Commercial Film developed in D-76. Toe densities were noted to remain relatively constant for this series. A probable reason is Commercial Film is coarse grained compared to Fine Grain Release Positive and has almost twice as much silver¹². Perhaps this extra amount of silver halide (coated in two coats), along with the use of the milder D-76 produced a stalemate between the two processes described above. That is, perhaps the rate of density decrease due to fixation was equalled by the rate of increase of density by fog production and therefore little or no net change in toe density resulted.

Figure 13 shows un-fixed, wet IR density plotted against fixed-out, dry diffuse density. Again a straight line can be fit to the data, but the s-shape is more pronounced in this plot than any of the others. It is felt that a linear model is still a good choice to describe the relationship, but further study is

indicated. It is interesting to note that this same type of S-shaped relationship appeared in early studies by M. Piskacek⁸ in 1976.

The last graph of this series, figure 16, shows both wet and dry D Log exposure curves for the same film sample after 5 minutes of development. The slopes of the two curves are roughly parallel and therefore differ mostly by some additive constant. This constant is probably due to the base plus fog density remaining almost constant through out the process as was explained earlier. This is only supposition, and further study is indicated.




Figure 14 : Comparison of D-Log H curves for wet unfixed, and dry fixed out Fine Grain Release Positive in D-19.




Figure 15 : Comparison of D-Log H curves for wet unfixed, and dry fixed out Fine Grain Release Positive in D-76.




Figure 16 : Comparison of D-Log H curves for wet unfixed, and dry fixed out Commercial film, type 6127, in D-76.

Repeatability Study:

From each of the three different film/developer combinations, three sets, or replications, were obtained. Recall from the scan programming used, 100 scans were made during each process run. As there were 3 replications of each of the 3 series, this makes for a total of 900 scans being taken! Of the 100 scans taken for each process, 8 were selected as being representative, and are listed as raw data in the appendix (for a total of 72 scans for all 3 series). To be able to make statements about the repeatability of the device, the 8 scans recorded were taken at the same time within each process. Of these 8 scans taken, 3 from each process run are used in the following analysis. The three elapsed times selected correspond to the ones used for plotting the family of curves shown in figs. 8, 9, and 10.

Table 1 is data from series 1. The data are grouped as to when during development they were acquired. The first column on the left is the step number. As the step numbers increase, so does the amount of exposure that particular step received, and hence an increase in density is noted. The next three columns are the replicated data sets for the development time under study. The next \bar{D} gives the average density of each step for all three replicates, followed by the sample standard deviation, S_D , for each step. The last column is simply $D_{\max} - D_{\min}$ for the three readings taken at each step just to give some idea of what the data spread is. Table 2 is for data taken during series 2 and table 3 is for data taken during series 3.

Tables 4 and 5 show data taken from a fourth run not mentioned earlier. This data was taken as a fourth replicate of experimental series 1 and 2, but with one major difference; these samples were processed without a fresh change of developer. Also, the tubes leading to the development chamber were not given

a chance to drain, as had been the case between all other process runs. The purpose of this is to demonstrate what effect on data repeatability using developer warmed by an earlier process run can have. It was noted that after a process run, the developer temperature would be 2 or 3 degrees higher than at the start of processing. Therefore, an increase in density due to an increase in temperature can be expected. The data in these two tables support this statement. Average density, \bar{D} , and the sample standard deviation, S_D , are included to show that the increase is a real one and not simple variability within the process.

Repeatability study: Series 1Table 1 : Repetition data for Fine Grain Release Positive
developed in D-19.5 min 0 sec

<u>Step #</u>	<u>Rep 1</u>	<u>Rep 2</u>	<u>Rep 3</u>	<u>\bar{D}</u>	<u>S_D</u>	<u>Spread</u>
1	.32	.30	.33	.32	.015	.03
2	.33	.30	.33	.32	.017	.03
3	.33	.31	.33	.32	.012	.02
4	.36	.34	.36	.35	.012	.02
5	.40	.36	.40	.39	.023	.04
6	.51	.51	.51	.51	0	0
7	.84	.80	.87	.84	.035	.07
8	1.49	1.43	1.50	1.47	.038	.07
9	2.27	2.24	2.31	2.27	.035	.07
10	2.94	2.90	2.93	2.92	.020	.04
11	3.08	3.06	3.08	3.07	.012	.02

2 min 0 sec

1	.19	.18	.19	.19	.006	.01
2	.20	.17	.19	.19	.015	.03
3	.20	.19	.19	.19	.006	.01
4	.22	.21	.22	.22	.006	.01
5	.24	.22	.23	.23	.010	.02
6	.29	.31	.30	.30	.010	.02
7	.46	.45	.50	.47	.026	.04
8	.87	.86	.89	.87	.015	.03
9	1.47	1.42	1.48	1.46	.032	.06
10	2.15	2.09	2.17	2.14	.042	.06
11	2.60	2.56	2.61	2.59	.026	.05

Series 1 (con't)

30 sec

<u>Step #</u>	<u>Rep 1</u>	<u>Rep 2</u>	<u>Rep 3</u>	<u>\bar{D}</u>	<u>SD</u>	<u>Spread</u>
1	.17	.17	.17	.17	0	0
2	.17	.16	.18	.16	.010	.02
3	.18	.18	.17	.18	.006	.01
4	.19	.19	.19	.19	0	0
5	.20	.18	.19	.19	.010	.02
6	.20	.23	.21	.21	.015	.03
7	.24	.25	.28	.26	.021	.04
8	.38	.39	.41	.39	.015	.03
9	.64	.59	.60	.61	.026	.06
10	.91	.86	.89	.89	.025	.05
11	1.17	1.10	1.17	1.15	.040	.07

Repeatability study:

Series 2

Table 2: Repetition data Fine Grain Release Positive developed in D-76.

5 min 0 sec

<u>Step #</u>	<u>Rep 1</u>	<u>Rep 2</u>	<u>Rep 3</u>	<u>\bar{D}</u>	<u>S_D</u>	<u>Spread</u>
1	.10	.08	.08	.09	.012	.02
2	.10	.10	.08	.09	.012	.02
3	.10	.10	.09	.10	.006	.01
4	.16	.16	.15	.16	.006	.01
5	.34	.33	.33	.33	.006	.01
6	.65	.64	.65	.65	.006	.01
7	1.05	1.07	1.07	1.06	.012	.02
8	1.48	1.50	1.50	1.49	.012	.02
9	1.97	1.99	1.98	1.98	.010	.02
10	2.42	2.47	2.46	2.45	.026	.03
11	2.72	2.74	2.80	2.75	.042	.08

1 min 30 sec

1	.16	.15	.14	.15	.010	.02
2	.17	.17	.13	.16	.023	.04
3	.15	.16	.13	.15	.012	.03
4	.18	.20	.16	.18	.020	.04
5	.25	.23	.24	.24	.010	.02
6	.39	.38	.40	.39	.010	.02
7	.59	.59	.60	.59	.006	.01
8	.81	.83	.82	.82	.010	.02
9	1.10	1.12	1.10	1.11	.012	.02
10	1.37	1.37	1.37	1.37	0	0
11	1.54	1.57	1.56	1.56	.015	.03

Series 2 (con't)

30 sec

<u>Step #</u>	<u>Rep 1</u>	<u>Rep 2</u>	<u>Rep 3</u>	<u>\bar{D}</u>	<u>S_D</u>	<u>Spread</u>
1	.18	.17	.17	.17	.006	.01
2	.19	.20	.16	.18	.020	.04
3	.16	.19	.16	.17	.017	.03
4	.18	.20	.17	.18	.015	.03
5	.17	.17	.16	.17	.006	.01
6	.21	.19	.22	.21	.016	.03
7	.25	.23	.25	.24	.012	.02
8	.28	.32	.29	.30	.021	.04
9	.37	.39	.38	.38	.010	.02
10	.47	.44	.47	.46	.017	.03
11	.50	.54	.51	.52	.021	.04

Repeatability study:

Series 3

Table 3: Repetition data for Commercial Film developed in D-76.

5 min 0 sec

<u>Step #</u>	<u>Rep 1</u>	<u>Rep 2</u>	<u>Rep 3</u>	<u>\bar{D}</u>	<u>S_D</u>	<u>Spread</u>
1	.21	.20	.21	.21	.006	.01
2	.26	.26	.26	.26	0	0
3	.41	.41	.41	.41	0	0
4	.69	.68	.69	.69	.006	.01
5	.99	.98	.99	.99	.006	.01
6	1.28	1.27	1.29	1.28	.010	.02
7	1.58	1.55	1.58	1.57	.018	.03
8	1.88	1.87	1.89	1.88	.010	.02
9	2.31	2.30	2.33	2.31	.015	.03
10	3.32	3.75	3.49	-	-	-
11	3.76	3.76	3.76	-	-	-

2 min 0 sec

1	.19	.18	.19	.19	.006	.01
2	.21	.19	.21	.20	.012	.02
3	.26	.26	.26	.26	0	0
4	.39	.37	.38	.38	.010	.02
5	.53	.52	.54	.53	.006	.02
6	.71	.70	.73	.71	.015	.03
7	.90	.89	.93	.91	.020	.04
8	1.13	1.13	1.16	1.14	.017	.03
9	1.38	1.38	1.42	1.39	.023	.04
10	1.68	1.69	1.73	1.70	.026	.05
11	1.96	1.97	2.03	1.99	.038	.07

Series 3 (con't)

50 sec

<u>Step #</u>	<u>Rep 1</u>	<u>Rep 2</u>	<u>Rep 3</u>	<u>\bar{D}</u>	<u>S_D</u>	<u>Spread</u>
1	.19	.18	.19	.19	.006	.01
2	.19	.18	.19	.19	.006	.01
3	.19	.19	.19	.19	0	0
4	.22	.20	.21	.21	.010	.02
5	.23	.22	.24	.23	.010	.02
6	.28	.26	.29	.28	.015	.03
7	.33	.33	.35	.34	.012	.02
8	.41	.41	.44	.42	.017	.03
9	.51	.52	.55	.53	.021	.04
10	.64	.67	.69	.67	.025	.05
11	.75	.77	.81	.78	.031	.06

Repeatability study:

Series 1

Table 4 : Process data for Fine Grain Release Positive developed in D-19 demonstrating effect of developer heating.

5 min 0 sec

<u>Step #</u>	<u>\bar{D}</u>	<u>S_D</u>	<u>Spread</u>	<u>Run#4</u>	<u>$D_i - \bar{D}$</u>
1	.32	.015	.03	.36	.04
2	.32	.017	.03	.36	.04
3	.32	.012	.02	.38	.06
4	.35	.012	.02	.38	.03
5	.39	.023	.04	.43	.04
6	.51	0	0	.56	.05
7	.84	.035	.07	.94	.10
8	1.47	.038	.07	1.61	.14
9	2.27	.035	.07	2.40	.13
10	2.92	.020	.04	3.07	.15
11	3.07	.012	.02	3.12	.05

2 min 0 sec

1	.19	.006	.01	.22	.03
2	.19	.015	.03	.20	.01
3	.19	.006	.01	.22	.03
4	.22	.006	.01	.23	.01
5	.23	.010	.02	.23	0
6	.30	.010	.02	.32	.02
7	.47	.026	.04	.54	.07
8	.87	.015	.03	.96	.09
9	1.46	.032	.06	1.59	.13
10	2.14	.042	.06	2.29	.16
11	2.59	.026	.05	2.75	.16

Repeatability study:

Series 2

Table 5: Process data for Fine Grain Release Positive developed in D-76 demonstrating the effect of developer heating.

5 min 0 sec

<u>Step #</u>	<u>\bar{D}</u>	<u>S_D</u>	<u>Spread</u>	<u>Run#4</u>	<u>$D_i - \bar{D}$</u>
1	.09	.012	.02	.09	0
2	.09	.012	.02	.10	.01
3	.10	.006	.01	.10	0
4	.16	.006	.01	.17	.01
5	.33	.006	.01	.36	.03
6	.65	.006	.01	.68	.03
7	1.06	.012	.02	1.11	.05
8	1.49	.012	.02	1.54	.06
9	1.98	.010	.02	2.02	.04
10	2.45	.026	.03	2.50	.05
11	2.75	.042	.08	2.85	.10

1 min 30 sec

1	.15	.010	.02	.15	0
2	.16	.023	.04	.16	0
3	.15	.012	.03	.14	-.01
4	.18	.020	.04	.19	.01
5	.24	.010	.02	.26	.02
6	.39	.010	.02	.42	.03
7	.59	.006	.01	.64	.05
8	.82	.010	.02	.88	.06
9	1.11	.012	.02	1.16	.05
10	1.37	0	0	1.43	.06
11	1.56	.015	.03	1.64	.08

Uniformity of Development:

The Fluid Transport System (FTS), which includes the chamber, pump, tubing and valves, presented some of the most complex design problems encountered in the project. The most important question was; would the flow of developer through the development chamber give uniform development, or would there be some type of directional effect. The problem is really one of hydraulics, and therefore, design procedures are not always as straightforward as one would like. The main design criterion was to produce a laminar flow of fluid through the chamber. Laminar flow is the smooth, uniform movement of fluid that is free of turbulence and bubbles.

To determine what level of development uniformity could be obtained once the FTS was actually constructed, the following experiment was performed.

- 1.) Samples of Eastman Kodak Pan-X Recording Film, Type S0-164, were given two levels of uniform exposure. The first set were exposed to room lighting for five minutes to give a developed density of 2.4. The second set received their exposure from an open-gated Durst enlarger. The level of exposure was enough to produce a developed density of approximately 1.5.
- 2.) All samples were developed in D-76 at 68° F for 5 minutes. Development was followed by a one minute rinse, then by a two minute fix. Samples were then washed and dried in the normal manner.
- 3.) The pump used by this system has 7 different speed settings (see Figure 20 for flowrate characteristics). This is about the only variable over which the user of the IR densitometer has control when it comes to flow

studies. Therefore, the uniformly exposed samples from above were each developed at a different flowrate.

The results of this experiment are shown in figures 17 and 18. Figure 17 shows dry, fixed-out diffuse density as a function of position for the sample that received the room light exposure. The position marked "IN" is where the developer first enters the chamber and strikes the film sample. The first density reading is taken there. The next four readings are taken at equally spaced intervals along the detector region, i.e. the region of the sample where the eleven IR densitometers are located. The last reading called "OUT" is from the area just below the chamber exit port. Figure 18 shows the same type of graph, but for the samples exposed for a developed density of approximately 1.5.

The data presented here is not exhaustive, but general statements about uniformity of development can be made. The data suggest that as the flowrate of developer is increased through the chamber, directional effects begin to decrease. However, this same increase in developer flowrate also causes an increase in the overall density of the film sample. This all sounds reasonable as one would expect an increase in developer flow to bring in fresh developer faster and remove development by-products that could retard development. Another general effect is that densities always are somewhat higher where the developer enters and somewhat lower where it exits. This is really not of much concern as the graphs tell us that if the pump flowrate is made high enough, very uniform development down the length of the film where the eleven densitometers read can be obtained. Also, figure 18 seems to show that even better uniformity can be achieved as the overall exposure level is reduced. This again sounds reasonable because as the exposure level is decreased, one would expect the amount of developer by-products to decrease and therefore directional effects should lessen.

Lastly, if uniformity is critical, users may wish to perform more exacting uniformity experiments for their particular film/developer combination. They could then possibly generate a methodology for correcting their raw data for directional effects.




Figure 17 : Density as a function of position and flow rate for a sample of Pan-X Recording Film, S0-164, exposed to room lights for 5 minutes and developed in D-76.




Figure 18 : Density as a function of position and flow rate for a uniformly exposed sample of Pan-X Recording Film developed in D-76.

Experimental Conclusions and Recommendations:

As should be obvious by now, the amount of additional research that could be performed on the IR densitometer is tremendous. However, it is felt that enough data has been collected to allow some generalized statements about it's performance to be made.

First, the assumption that an approximately linear relationship exists between wet, un-fixed IR density and dry, fixed-out diffuse density seems reasonable. More study is indicated in this area, with particular focus on what is happening optically to a wet, developing emulsion. Perhaps experiments should be conducted with software calibrations being performed when film samples are wet instead of dry, as is now currently recommended. Also, all experiments were performed with the anti-halation backing untouched. Although the backside of the film is never wetted during processing, perhaps the infrared radiation has some effect on the backing that could degrade density readings later on in the process.

Second, the run-to-run repeatability of the device is considered to be very good. The data from tables 1, 2, and 3 seem to indicate that the best repeatability is obtained at lower densities. As an example, the device is able to hold an average repeatability of $\pm .02$ out to a density of .50 for all three series, or a variability of approximately 5%. This value becomes only $\pm .03$ out to a density of 2.5. Also observed, was that the more active the developer, the poorer the run-to-run repeatability. Other important results from the experiments are the effect on repeatability due to developer temperature variations. It was noted that a five minute development run would raise the temperature of the developer by 2 or 3 degrees. This is due to the heat generated by the pump and valves, and general frictional effects. This rise in temperature is enough to cause a

significant increase in density if another process run is made using the same tank of developer. Also, if the tubes leading to the development chamber are not drained between runs, repeatable results are difficult to obtain due to solution carry overs. The full extent of what solution carry over will do, however, is not known and further study is indicated.

Lastly, uniformity of development across the region where the IR densitometers are located seems to be very good and not a source of great concern. There are major directional effects near the edges of the film sample where the O-ring seal is located, but these areas are not seen by the IR densitometers. As was stated earlier, if extreme precision is required, specific studies can be made to better characterize what little non-uniformity there is, and correction procedures can be implemented. Suggested improvements might be a temperature control for the FTS and a pump of higher and more consistent flowrate.

Table 6: Log exposure values used for all experiments.

Sensitometer : Kodak Model 101.

Wedge : Kodak Sensitometric Step Tablet #5, S/N R903-12-2.

Open-gate exposure : 340 lux·seconds.

Neutral density filter type : Inconel

<u>Step</u>	<u>Step Density</u>	<u>Log E w/.78 ND</u>	<u>Log E w/1.20 ND</u>
1	3.07	-1.32	-1.74
2	2.77	-1.02	-1.44
3	2.46	-.71	-1.13
4	2.16	-.41	-.81
5	1.86	-.11	-.53
6	1.54	.21	-.21
7	1.24	.51	.09
8	.94	.81	.39
9	.64	1.11	.69
10	.34	1.41	.99
11	.04	1.71	1.29

In-process IR Densitometer Operating Instructions:

For best results, the user of this device should have a complete understanding of its operating principles. Not only does this involve just knowing how to make the device run, but also a thorough knowledge of how the densitometer works.

This operations manual has been broken up into six main function groups. They are:

- A.) Emitter calibration.
- B.) Software calibration.
- C.) Film sample preparation.
- D.) Data entry.
- E.) Fluid Transport System (FTS).
- F.) Prompt descriptions, and error messages.
- G.) Printer Maintenance.

The user should read and understand these instructions and descriptions before attempting to use the densitometer. In this way, the user can be assured the best possible results, and wear and tear on the device can be minimized.

Equipment List:

Below are listed the various components needed to operate the In-process IR Densitometer. Some items are kept with the densitometer and will be found on the lower shelf of the rolling service cart. Other items are of general purpose and will have to be located before operating the device.

Items to be kept with the densitometer:

1. Special two-hole registration punch.
2. Modified head for the Kodak Model 101 sensitometer.
3. Calibration samples.
4. Chamber locking pin.

User supplied items:

1. Voltmeter. Range of 20 volts with accuracy of $\pm .02$ volt.
2. BNC-to-banana-jack adaptor, or some type of adaptor to fit a female BNC connector to voltmeter being used.
3. Test lead. A clip lead of some sort is necessary to effect a ground connection between the voltmeter and the chamber support frame.
4. Small screwdriver.
5. Darkroom. The densitometer must be loaded and used in total darkness.
6. Source of fresh water to clean the device when processing is completed.

Procedural Descriptions:

A.) Emitter Calibration

This densitometer system uses gallium arsenide infrared emitting diodes (IREDS) as its primary source of radiation for the measurment of optical transmission, hereafter referred to as "density". As these sources are diodes, their output intensity is a function of the forward current passing through the device, rather than the voltage across the diode.

Two other important considerations must be pointed out before detailing the calibration procedure. The first is the type of detector used with the above emitters. Due to their mode of operation, they can only detect a little more than a three decade change in light with an acceptable signal to noise ratio. The second consideration is the operating principle of the entire IR densitometer itself: the densitometer is to measure density of a piece of film as it develops. This implies that the film sample may still have an anti-halation backing on it while readings are being made. So not only will the density of the developing silver be measured, the infrared density of the wet anti-halation backing will also be measured.

The point to all of this is that the user would like to minimize the effect of the anti-halation backing on the density measurments. At the same time, the user would also like to maximize the dynamic range of the densitometer to allow for density readings of at least 3.00 and still maintain some level of accuracy. Therefore, it is necessary to adjust the output of the IR emitters to compensate for the various types of anti-halation backings likely to be processed by this device.

This is done by placing a sample of the particular film under study, dry and unexposed, into the optical path of each of

the eleven emitter/detector pairs. The current passing through each IR diode is then adjusted until the output voltage of the companion detector just begins to saturate. Saturation means the maximum output voltage the device is capable of producing while remaining linear. The exact procedure is outlined below, but in this way each emitter/detector pair is optimized for maximum dynamic range in density.

It should be noted that this method of emitter calibration defines what will be considered zero density. That is, a dry, unexposed piece of raw film stock will have, by definition, zero density. All density measurements are therefore relative to dry unexposed film. However, note that once a film sample becomes wet, it's IR density will increase due to emulsion swelling.

Procedure:

- a). With the main power switch in the off position, disconnect the Analog Transmission Signal cable from the back of the control cabinet. This is the black coaxial cable with the large female BNC connector on the end. The control cabinet is thoroughly labeled and will direct the user to this connector.
- b.) This part of the calibration can be performed with just about any type of voltmeter. At the time of this writing, however, a Hewlett-Packard digital voltmeter, (H/P DVM) was available, and thus the procedure will be detailed for this particular type of voltmeter.

Insert the banana-plug-to-BNC adaptor into the H/P DVM. The adaptor should be orientated such that the banana plug marked "COM" fits into the red ground socket on the front panel of the voltmeter. The other plug fits into the red socket marked "VOLTS". DO NOT INSERT BACKWARDS!! To do so will ground the output of

the analog multiplexer and destroy it.

- c). Connect the Analog Transmission Signal cable to the other side of the BNC-to-banana adaptor. Insert a banana jack test lead into the side of the adaptor marked "COM". Attach the other end to the metal chamber support structure. In this way, an electrical ground path is established between the voltmeter and densitometer.
- d). Select the D.C. volts button on the front panel of the DVM. Select the 20 volt range button. Turn the DVM on. Now turn on the main power to the IR densitometer and wait 10 minutes for it to warm up. The DVM should now be reading +15 volts +/- .1. This is a test voltage upon initial power-up. The output of these detectors is directly proportional to the amount of incident flux striking them, therefore, the output is considered a transmission signal.
- e). In complete darkness, or under safelight conditions if the film will allow, punch a set of registration holes in a sample of the film type under study. Cut it long enough to cover all eleven detectors, about 8 inches. This is then mounted in the chamber, placing the strip such that the registration holes fit over the two pins at the rear of the chamber. It is most important that the film be mounted emulsion side up, toward the developer channel. Lastly, close the chamber and seal it with the locking pin. The device is now ready to be calibrated.
- f). When the main power was first applied, the system came up and prompted: "EMITTER CAL?". If the user is going to process a film type different from that for which the device was last calibrated, an emitter calibration will be necessary. To select the calibrate mode, enter any number from 1-9 , followed by an "E" depression. The system will respond with the prompt: "EMITTER 1". This

means that the emitter/detector pair at the top of the chamber (#1) has been activated.

- g). The DVM should now be reading some voltage. Adjust the grey trimpot marked #1 located behind the chamber with a small screwdriver until the meter reads approximately 13.0 volts +/- .2 volt. Turn the pot clockwise to increase the voltage. Turn it counter clockwise to decrease the voltage.
- h). Press the "E" key and the display will change to: "EMITTER 2". Repeat as above for all eleven emitters. When complete, the system will again prompt: "EMITTER CAL?". To leave this calibration mode, enter zero, followed by "E". This will jump the user to the next prompt.
- i). When all eleven emitters have been calibrated, turn the machine off and reconnect the AnalogTransmissionSignal cable to the control cabinet. Remove the film base sample from the chamber and continue with the software calibration routines.

Note that this emitter calibration will be valid as long as the same type of film is to be used. Turning off the power will not change the settings of the trim pots, therefore, this calibration does not have to be performed each time the device is powered up. The only time the emitters have to be calibrated is when processing a film type different from that for which the device was last calibrated.

B.) Software Calibration:

The basic premise behind the software calibration is that the relationship between wet unfixed IR density and dry, fixed-out diffuse density is a linear one. This important fact was amply demonstrated by M. Piskacek in 1977.

The operation of the software calibration is very straight forward. Certain conditions are established within the development chamber, (detailed below) and a single complete density scan is made. The eleven readings that result are considered to be zero density, and stored in memory as MD0 (measured zero density). New conditions are again established within the chamber, this time with a film sample of uniform density large enough to cover all eleven emitter/detector pairs. The density of this sample should be near the mid-point of the maximum density the user expects from the particular film/developer combination under study. For example, if a maximum density of 2.00 is expected, the calibration standard would have a density of 1.00. The density of this sample is known, and entered into the computer as ADMAX (actual maximum density). ADMAX is only a computer variable name, and really refers to mid-range densities. Another complete scan is made and a new set of eleven readings result. These are stored in memory as MDMAX (measured maximum density). With this information now in memory, the computer can now perform eleven separate linear regressions and generate eleven unique slopes and intercepts for each of the eleven emitter/detector pairs. In this way, as raw data is generated, it is operated on by a first order calibration routine. The final densities viewed later by the user, will have been corrected from wet, IR density to equivalent dry, diffuse density.

Figure 19 shows the relationship between the different variables involved. Calculations performed within the computer

are in a 32-bit floating point format. Please consult the Intel Floating Point Arithmetic Library User's Manual for details.




Figure 19: Software Calibration Plot.

MDMAX = Measured maximum density.

ADMAX = Actual maximum density.

MDO = Measured zero density.

Recall that eleven slope and intercept pairs must be generated for all eleven emitter/detector pairs.

From the plot:

$$MD = (AD) \times (SLOPE) + INTERCEPT$$

and the intercept is simply:

$$INTERCEPT = MDO$$

From another point on the plot:

$$MDMAX = (ADMAX) \times (SLOPE) + MDO$$

Solving for the slope:

$$\text{SLOPE} = (\text{MDMAX} - \text{MDO}) / \text{ADMAX}$$

Lastly, corrected densities can be calculated from:

$$\text{AD} = (\text{MD} - \text{MDO}) \times (\text{ADMAX} / (\text{MDMAX} - \text{MDO}))$$

Software Calibration Procedure:

- a). This calibration is to be performed after the emitter calibration. If the same type of film base is to be run for this processing session, and no emitter calibration was performed, turn the machine on and allow a 10 minute warm-up before proceeding.
- b) Open the chamber door. Make sure the chamber is absolutely clean and dry. Any moisture in the chamber will have an adverse effect on the calibration. Insert the calibration strip into the developer channel. The calibration sample is any material of uniform, known, IR density, close to the mid-point of the expected density range being investigated. Ideally, it should be cut to fit inside the development cavity, up against the IR emitters. It is important that it be placed in the channel so it will not affect the closing of the chamber door. Apply slight pressure to make sure that the strip is seated firmly in the channel.
- c). In total darkness, or under safelight conditions if the film will allow, unroll about eight inches of the film under study. Punch two registration holes into one end of the film strip using the special punch provided. Mount the film into the development chamber by aligning the two pins at the rear of the chamber into the holes

just punched. Be sure to mount the film emulsion side up, toward the developer channel. Hold the other end of film strip flat against the chamber bottom, close the chamber and seal it with the locking pin.

- d.) If the machine has just been turned on, the prompt "EMITTER CAL?" will be on the display. If performing an emitter calibration, complete the calibration and depress the "E" key until this same prompt appears. To get to the software calibration mode, enter zero and depress "E". The system should respond with: "ENTER D-MAX". The user should now enter the value of the calibration sample to be used. Remember that the decimal point is implied. That is, a density of 1.20 is entered as 120. After entering D-max, depress "E". The next prompt will be: "READ D-MAX". The next depression of the "E" key will read the density of the calibration sample plus the density of the dry film base. The value of this reading is assigned the calibration values just entered by the user.
- e). The system will now be prompting: "READ ZERO". Open the chamber and carefully remove the calibration sample. Try not to disturb the film sample while doing this. Reseal the chamber door and press the "E" key. The system has now just read the density of the film sample and assigned it the value of zero density.
- f). Open the chamber and remove the film base sample. The prompt: "ENTER MONTH" should now be on the display. At this point, the device is 100% calibrated and ready to go.

C.) Film Sample Preparation:

This device was originally designed to measured the density of a film sample that had been exposed to a standard 21-step sensitometric tablet. Alignment proved difficult, however, so a Kodak #5 step tablet is to be used for making all sensi exposures. This tablet has a step increment of about .30 density units, and has been mounted in a standard Kodak 101 sensitometer head. Note that any step tablet can be used, just as long as the step spacing is exactly 10.00 mm, otherwise poor registration of the step exposure to the detectors will give poor results.

Before making any exposures, make sure that this sensitometer head has been obtained! The procedure for making sensitometric exposures is as follow:

- a). In complete darkness, unroll about eight inches of the film under study. Use a safelight if the film will allow. Cut the film as close as possible to a right angle with the edge of the film. Punch two registration holes into the film with the punch provided.
- b). Align the holes in the film to the pins mounted on the end of the sensitometer head. Mount the film on the head with the emulsion facing the light source. Make the exposure. As with all sensitometer work, the best filter combination will have to be determined by calculation or trial.
- c). Mount the exposed sample in the development chamber. The emulsion side always faces up (toward the emitters). Close the chamber and seal it with the locking pin. The film sample is now ready for processing. Remember that the room lights have to remain off while the film is being loaded and processed.
- d.) The correct response to the "ENTER MONTH" prompt is to enter the two digit month code followed by an "E"

depression. This is repeated for the "ENTER DAY" and "ENTER YEAR" prompts.

- e.) The prompt: "MANUAL MODE?" should now be on the display. To select manual mode, enter any number from one to ten followed by an "E" depression. The printer will respond with the heading: "*** MANUAL MODE SELECTED ***". To select automatic, enter zero then a "E" depression. A more detailed description of manual mode and automatic is contained in Section E: Fluid Transport System.

D.) Data Entry:

Other than calibration, the user must also program the IR densitometer with certain process information. The device must be told for how long it is to process film samples and when to take density scans. This is done in the following manner:

- a). During the programming stage, the user will receive the prompt: "PROCESS 1=?". The device is asking the user how long to run the first process. There are three tanks and therefore the machine is capable of running three processes. This prompt is answered by entering any number in the range of 0 to 127 minutes. The smallest increment of time is one minute. Invalid entries are flagged as errors (see error message table).
- b). As this system has the ability to generate an incredible amount of data, a very selective methodology for the programming of scan acquisition was needed. This was accomplished by breaking each process into a number of smaller intervals. Then for each interval, a rate at which scan data will be acquired is defined. This rate is expressed as the number of scans to be taken per minute, with the range being from 1 scan/minute to 60 scans/minute. It is important to note that the system defines scan rates on a minute by minute basis. This means that for an entered scan rate of 32 scans/minute, the system will round this down to 30 scans/minute and take a scan every $60 / 30 = 2$ seconds. This was done because the shortest time for which a scan rate can be defined is one minute.
- c). After entering the process length, the system will prompt for the interval length. If an interval length of zero is entered, the device will assume that the user wishes to make no scans during this process and therefore jumps the user to the next process length entry. If

scans are to be made, the user should enter an interval length over which data is desired. At all times it should be remembered that the system can store only 187 scans and great care should be exercised not to exceed this value. Also, the entered interval cannot be greater than the process length. Once the interval length has been entered, a scan rate must be entered. The number of scans generated during any interval is simply the interval length times the scan rate. If invalid interval or rate information is entered, error messages will be issued.

- d). Interval/rate data can be entered until:
 - 1.) The number of generated scans equals 187.
 - 2.) Ten interval/rate pairs have been entered.
 - 3.) The sum of the interval lengths equals the process length.
 - 4.) An interval length of zero minutes is entered.
- e). This procedure is repeated for each of the three processes. If manual mode has been selected, the user is allowed to program only one process as the FTS is under total user control.

Every effort has been made to design the most efficient user interactive operating system as possible. There are, however, some ways a user can get into trouble. The most common mistake to be made is to use up all scans before the process is finished. This is done by entering an interval/rate combination that generates many scans. As an example, an interval length of 3 minutes and a rate of 60 scans per minute will produce 180 scans. This does not leave many scans for anything else. Remember that there are a total of 187 scans available and not 187 scans for each process. Once interval/rate information has been programmed, there is no way to clear it and start all over again. Either the device must be turned off and re-started or the programmed data has to be run through. The point is, be

careful when programming scan data!

At the completion of processing, the system will return control to the user with the prompt: "PROCESS #=?". Please refer to the prompt descriptions contained in Section F for directions.

E.) Fluid Transport System (FTS):

The Fluid Transport System, or FTS for short, consists of any component having to do with the transportation of chemistry to and from the development chamber. This system has two modes of operation, automatic and manual. There are certain aspects of this system that require careful attention, and will be repeated through-out this section. The first and foremost is; never run the device with the processing tanks empty! To do so will cause excessive wear of the pump gears and shorten pump life. Considerations of the two modes of operations are listed below.

Automatic:

During the programming of the device, the user will receive the prompt: "MANUAL MODE?". To by-pass manual mode and select automatic, the user must enter a zero followed by an "E" depression. Once in automatic mode, the computer will handle the complete operation of the inlet and outlet valves, and the pump. What is not under computer control, however, is the speed of the pump and the detection of fluids within the tanks. These two areas are always the user's responsibility.

To set the pump speed, locate the small knob mounted top dead center over the pump housing. The pump itself is located just below the chamber and right above the chemistry tanks. Refer to figure 20 for flowrate details. Position "0" is off. Position "8" is maximum speed. Best results are obtained at position "7", but the speed selected is up to the user. Tanks should be filled to just below the mounting screws, (about 2 liters) with care being taken to avoid spills onto the control valves. When programming for process lengths while in the automatic mode, the user can run each tank for 127 minutes each. Any combination of tanks can be run, but no tank can be run more than once, and




FIGURE 20: FLOWRATE AS A FUNCTION OF PUMP SETTING.

the order of usage is always tank 1, tank 2, tank 3. When the system begins to process, it has no way of knowing if the door to the development chamber is open or not and will always assume that it is closed. There is, however, an interlock switch on the chamber to interrupt the A.C. power to the pump if the door is left open, but the computer does not know this and will always try to start the pump. Even if the door is open, the inlet and outlet valves will open and the programmed timing sequence will begin. To see what exactly is turned on, refer to the FTS relay control box mounted under the cart and to the right. There are three pairs of LEDs and a single red LED at the bottom. The two leftmost red LED pairs are the control indicators for tank 1 inlet and outlet valves. The center two pairs are for tank 2 and the two on the right are for tank 3. The red LED by itself on the bottom is for the pump. Note that the layout of these LEDs is duplicated on the front panel of the control cabinet. This will be detailed below.

Manual:

To select manual mode, answer the "MANUAL MODE?" prompt by entering any number from 1 to 9 and depress the "E" key. The printer will print ** MANUAL MODE SELECTED ** and the manual FTS control buttons on the front panel of the control cabinet become active. When in manual mode, it is the user's responsibility to control the switching of tanks and the operation of the pump. This is to allow for the user to switch between any tank as often as desired. Tanks can be mixed, drained, new chemistry can be added, etc. The main purpose is to allow the user complete versatility in the processing of film samples.

Another reason for this manual control is for the draining and cleaning of tanks when finished with the device. To drain a tank, the procedure is as follows:

- a.) Locate the manual drain valve. It is the extra length of tubing that extends out from the top row of valves. Open the valve and place the free end of the tubing into a suitable container. See figure 21 for a schematic diagram of this system.
- b.) Depress one of the tank outlet buttons on the front panel. The valve should open with a loud snap. No more than one tank out button should be on at a time. If the manual drain valve is open, do not open any of the tank inlet valves, otherwise the tank will take forever to drain. Now that the bottom of a tank is open, depress the PUMP ON switch. Remember that the pump will only go on if the chamber door is sealed. Also, always keep a piece of film in the chamber when draining tanks, as this makes a better seal in the chamber and it is less likely to leak.
- c.) The tank should now be draining through the drain valve. Listen carefully for a sudden increase in pump speed. This is the best indication of when the tank is empty. Quickly turn off the pump, then the tank out valve. Repeat for any other tanks that need draining.

Important Note: NEVER NEVER!! try to operate the pump without having at least one outlet valve open and one inlet valve open. The drain valve can take the place of a inlet valve, however. The point is that the pump must always have a source from which to draw fluid and someplace to put the fluid. The pump can be damaged if made to draw from sealed tanks or empty tanks. Also, hoses may burst if the pump is not given an outlet destination. Therefore, ALWAYS make sure that the proper valves are open before turning on the pump. It is very easy to make this mistake so, BE CAREFUL!!!




Figure 21: Fluid Transport System Hydraulic Schematic.

Tank Cleaning:

The best method for cleaning the tanks is also the easiest. Place the device into manual mode. Fill all three tanks with clean, lukewarm water. Set the pump speed to about "5". Load a clean film sample in the development chamber. Open tank 1 inlet valve and outlet valve. Start the pump and run it for about 30 seconds. Now open the drain valve and close tank 1 inlet valve. The tank should now drain. Repeat this for the remaining two tanks. Remember that when the pump speed begins to increase, the tank is empty and the pump must be turned off immediately.

F.) Prompt Descriptions and Error Messages:

This section gives a brief description of every prompt the user will encounter when operating the In-process IR densitometer. The prompts are listed in their order of appearance to help guide the user through all steps of device usage.

<u>Prompt</u>	<u>User Action</u>
"EMITTER CAL?"	Device wishes to know if user would like to perform a hardware emitter output calibration. See calibration section for details. To perform the calibration, enter any number, then press "E". System will respond with "EMITTER 1,2,3..." for each "E" depression. Eleven depressions return the user to original prompt. To exit, or skip calibration, enter zero and press "E".
"ENTER D-MAX"	Enter the density value for the sample used as calibration standard. See calibration section for details. Density values are entered to the nearest hundredth with the decimal point implied. Example; 1.20 is entered as 120. Press "E" after entering correct value.
"READ D-MAX"	This is the prompt for the actual reading of the calibration standard. Place film stock sample and calibration strip into chamber and seal. Press "E" to take reading.
"READ ZERO"	Remove calibration strip from chamber, but leave film stock sample. Seal chamber door and press "E" to take zero reading. See calibration section for complete explanation of calibration routine.
"ENTER MONTH"	System is prompting for information so it can print a date header. Enter two digit month

code followed by "E" depression.

"ENTER DAY" Enter two digit day code, followed by "E".

"ENTER YEAR" Enter two digit year code, followed by "E".
system will print date header in the format:
xx/yy/zz.

"MANUAL MODE?" System is asking user if manual or automatic control of the FluidTransportSystem(FTS) is desired. Any number entered followed by the "E" key will enable front panel buttons and put FTS under manual (user) control. A zero followed by "E" will place FTS under automatic (computer) control. See FTS description for details. If manual mode is selected, system will print: ** MANUAL MODE SELECTED **.

"PROCESS 1=?" Enter 1st process length in minutes. Legal range is zero to 127 minutes. Process length can only be entered in increments of one minute. See data entry section for details. After depressing "E" key, system will print:
Process 1 = xyz minutes.

Note: An entered process length of zero minutes will jump user to "PROCESS 2=?" prompt. However, if manual mode has been selected, system will display "BEGIN" query instead and assume no scan data is to be taken.

"INTERVAL 1=?" Enter interval length in minutes; then press the "E" key. For the programming of scan acquisitions, the user must break up the entered process length into individual intervals. For each interval, the user must then enter a scan acquisition rate (see below). The sum of all the entered intervals cannot be greater than the total process length. If the sum is greater, an error prompt is issued.

"RATE 1=?"

See error prompts for details.

For the just-entered interval, a scan acquisition rate must now be entered. Legal rate entries are 1 scan/minute to 60 scans/minute (these translate to 60 seconds per scan to 1 second per scan). The system has the ability to store 187 complete scans. A running sum is kept internally as the user programs intervals and rates. The number of scans any combination generates is simply the rate times interval length.

If an entered rate value, in conjunction with the just entered interval, generates a scan sum greater than 187, the system prompts "RATE TOO BIG". If the scan sum exactly equals 187, the system prompts: "OUT OF SCANS" and prints out as follows:
Interval 01 = ab minutes

Rate 01 = xy scans/minute

"INTERVAL 2=?"

Interval/rate combinations may be entered until:

1.) The sum of interval lengths equals process length.

2.) Scan sum equals 187. Note that even though this is stored internally, it is not available for user display, and it is therefore the users responsibility to keep track of the scan sum independently of the machine.

3.) Ten interval/rate pairs have been entered.

Note: Entering an interval length of 00 minutes will jump the user to the next process. In the case of manual mode, only one process is allowed.

"PROCESS 2=?"

"PROCESS 3=?"

Same methodology applies for Process 2 and 3. As an example, if all scans are programmed in Process 1 and the "OUT OF SCANS" prompt is obtained, the remaining processes can still be run, however no scan acquisition data can be programmed.

"BEGIN"

At this point, system is ready to begin processing. User should now check to see that tanks are filled, chamber door is sealed, etc. Processing will begin when enter key is depressed. If manual mode has been selected, it is up to the user to operate the manual FTS controls. Scan acquisition will still be under computer control, however. If manual mode has not been selected, the FTS will be under complete computer control also.

Note: When processing is started, all displays will go blank (except the digital thermometer). The user can do nothing until processing is complete. If the processing must be interrupted, the main power switch is the only way to stop the device.

"PROCESS #=?"

When processing is complete, system returns with this prompt. User is now to enter process number from which it is desired to view scan data. Legal entries are 1, 2, or 3 only, followed by "E". If a process number is entered for which no scan acquisition data has been programmed, the error prompt "NO DATA" is displayed.

"TIME WANTED"

First set of double zeros are the number of minutes into selected process user wishes to view data. Legal entry range is 00 to 99 minutes, followed by the "E" key.

"TIME WANTED" Second set of zeros are the number of seconds to be added to the number of minutes already entered to form total time into process user wishes to view scan data. Legal entry range is 00s to 99 seconds, followed by "E".

Note: Time entries of 00 min 00 sec are not allowed and will return user to "PROCESS #=?" prompt.

"TIME ELAPSED" The times at which scans are made depends on how the scan acquisition data was programmed. Therefore it is possible for the user to enter a time at which no scan was made. The system will locate the closest scan made to the elapsed time wanted and display this actual scan time along with the scan data on the density display.

Data View Mode :

Once the first set of scan data is displayed, the special function keys become operative. Their functions are outlined below:

UP: Depressing this key will advance user to next acquired scan. Density display will show this new information and prompt display will be updated to show change in elapsed time. If already showing last scan in process, system will display: "NO MORE DATA".

DWN: This key will cause system to display scan acquired prior to present scan already on display. Note that time steps between scans are determined by scan rate information programmed earlier by the user. If very first scan in process is already being displayed, pressing "DWN" will cause system to display: "NO MORE DATA".

PRT: This key will cause the line printer to generate a copy of the scan data presently being displayed. The process number along with the elapsed time are also printed as the header. This key has absolutely no effect on the display contents.

NXT: Depressing this key will simply return the user to the prompt: "PROCESS #=?". In this fashion, desired scans can be arrived at directly, rather than inched up to with the UP and DWN keys.

"E"- "E" Two "enter" key-strokes in succession will return the user to "MANUAL MODE" query. This is to be done only when the user has completed the inspection of process results and wishes to make a new process run. Any number entered followed by an "E" will then enable the FTS manual control. This is to allow the user to change or alter the contents of the process tanks if so desired. A zero followed by the "E" key will jump the

user back to the "ENTER MONTH" prompt. An entire new process run can now be programmed. All calibration values obtained earlier are retained.

Note: Once the user has entered two "E" key strokes, there is no way to return to the data view mode except to run a complete new process. Also, the double "E" stroke will only work when the message: "TIME ELAPSED" is on the prompt display.

Error Messages :

The following error prompts are returned by the system when the user makes an illegal or invalid keyboard entry. Error prompts will appear after an "E" depression and remain on the display for approximately five seconds. After this time, the system will reprompt the user for the correct entry. For critical errors, the re-prompt will be flashing on and off to alert the user that some careful thought is required.

<u>Error prompt:</u>	<u>Meaning:</u>
"PROCESS FULL"	Not really an error, but just informing the user that the sum of entered intervals now equals entered process length, and no more interval entries will be allowed.
"RATE TOO BIG"	Entered rate in conjunction with just entered interval will generate more scans than system is capable of storing (187). Solution: re-enter smaller rate.
"OUT OF SCANS"	A reminder to the user that all scans have been used up. Further interval/rate entries will not be allowed.
"NO DATA"	Response to user request for data from a process where none has been acquired.
"NO MORE DATA"	Response generated as user views a block of scan data and comes to the top or bottom of that data block.

Note: Certain types of invalid entries will generate no error message. What will happen, however, is that the number field of the prompt display will be re-set to all zeros, and the prompt will begin to flash. In this way, the user cannot proceed in system programming until a legal entry is made.

G.) Printer Maintenance:

There are only two types of service operations the IR densitometer user may attempt on the 40-column line printer. They are to change the ribbon cassette and to install a new roll of paper. For other problems, it is recommended that the printer manufacturer be consulted (NCR corporation).

Important note: Always disconnect the A.C. power cord before attempting to service the printer. A.C. line voltage is present within the control cabinet and a severe electrical shock could result from the inadvertent touching of one of these connections.

To install a new ribbon cassette:

1. Remove the old ribbon cassette from the printer. This is done by gently squeezing the two plastic locking tabs at the center edges of the cassette together and lifting the cassette off the printer.
2. Before installing the new cassette, locate the small white ribbon advance knob on the underside of the cassette. Slowly turn the knob in the direction indicated by the arrow until any slack in the ribbon is removed. To replace the cassette on the printer, again squeeze the locking tabs together and set the cassette back onto its support tongs.

To install a new roll of paper:

1. Remove the ribbon cassette from the printer as described above.

2. Lift out the two vent plugs located just behind the printer on top of the control cabinet.
3. Slide the entire top panel out the back of the cabinet. Pull the panel out squarely as to prevent it from jamming. Lift slightly as the front end passes over the printer mechanism.
4. Step behind the device and view the paper feed spool located just behind the printer. Remove the wood spool from it's supports by lifting straight up on it at both ends.
5. Slide off the two plastic spacers and discard the cardboard core from the old roll. Insert the wool spool into a fresh roll of paper. The curl of the paper should be facing you. The long plastic spacer goes on the left and the short one on the right. Replace the entire assembly by guiding the slots in the spool into the U-shaped support holes.
6. Important!! Check now to see that the paper roll rotates freely on the spool and that it does not bind. If the roll is binding, it is because the wood spool was not replaced squarely in it's supports. Remove and replace if necessary.
7. While still standing behind the device, locate the feed roller release lever. This is a white plastic lever on the right hand side of the printer. Gently pull this lever towards you and note how the platen lifts away from the feed roller. While holding the platen open, feed the end of the paper under the metal roller opposite the rubber feed roller. Guide the paper out the top of the printer until the paper is feeding squarely

off the supply roll. Release the lever and tear off the excess paper you just pulled through.

8. Replace the top panel, again lifting the front edge slightly as it passes over the printer. Replace the two vent plugs. Before replacing the ribbon cassette, remove any slack by turning the white knob that is under the cassette in the direction indicated.
9. Reconnect the A.C. power.

References

1. L .J. Fortmiller and T. H. James, RPS Proceedings of the Centenary Conference, London, Sept. 1953.
2. J. Hughes, B.S. Thesis, Rochester Institute of Techonolgy, 1964.
3. J. Hisler and A. Casinelli, B.S. Thesis, Rochester Institute of Technology, 1970.
4. T.L. Beaupre and R.R. Jasper, B.S. Thesis, Rochester Institute of Technology, 1971.
5. D.A. Turbide and M.T. Williams, B.S. Thesis, Rochester Institute of Technology, 1972.
6. T.L. Beaupre, R.R. Jasper, D.A. Turbide, and M.T. Williams (presented by Dr. B.H. Carroll), Photogr. Sci. Eng., 18, 535 (1974).
7. A. Spitzak, J.SMPTE, 75, 103 (1966).
8. C.J. Keemink and G.J. Van der Wildt, J. of Applied Photogr. Eng., 2, No. 1, 49 (1976).
9. M. Piskacek, M.S. Thesis, Rochester Institute of Technology, June 1977.
- 10 C.R. Berry, Photogr. Sci. Eng., 16, No. 5, 349 (1972).

11. Processing Chemicals and Formulas, Data book J1, Eastman Kodak Co. (1973).
12. Personal Communication with Dr. B.H. Carroll, March, 1982.

Appendix 1 : Input/Output Port and Peripheral Device Addressing

1.) Line Printer Data and Control Lines:

The 40-column line printer is considered to be I/O mapped I/O. That is to say, all operations with the printer are made through three I/O ports. Port addresses and bit designations are listed below. Data is clocked into the printer by the printer clock input line and the clock signal itself is derived from the Timer Out pin located on the expansion 8155 (chip A17). The Printer Reset input is derived from the buffered 8085A Reset Out signal. See the NCR users manual and the SDK-85 users manual for more complete details.

PORT 23H : Define as Input.

Bit C5 : Not Used
Bit C4 : Not Used
Bit C3 : Not Used
Bit C2 : Motor Jam
Bit C1 : Printer Busy
Bit C0 : Not Used

PORT 22H : Define as Output.

Bit B7 : Not Used
Bit B6 : Not Used
Bit B5 : Not Used
Bit B4 : Not Used
Bit B3 : Not Used
Bit B2 : Not Used
Bit B1 : Not used
Bit B0 : Printer Write

PORT 21H : Define as Output.

Bit A7 : Data 7
Bit A6 : Data 6
Bit A5 : Data 5
Bit A4 : Data 4
Bit A3 : Data 3
Bit A2 : Data 2
Bit A1 : Data 1
Bit A0 : Data 0

2.) A/D Converter Control and Data Lines:

All the Input/Output ports listed below have been optically isolated from the peripherals they communicate with. The devices used are Hewlett-Packard 2731 opto-isolators. For more details, consult the system schematics in the appendix. The purpose of the isolation is prevent noise crossover between the digital and analog systems.

PORT 2BH : Define as Output.

Bit C5 : Not Used
Bit C4 : Not Used
Bit C3 : Not Used
Bit C2 : Not Used
Bit C1 : Not Used
Bit C0 : Start Conversion

PORT 2AH : Define as Input.

Bit B7 : Conversion Status
Bit B6 : Not Used
Bit B5 : Not Used
Bit B4 : Not Used
Bit B3 : Not Used
Bit B2 : Not Used
Bit B1 : Conversion Bit 1 (MSB)
Bit B0 : Conversion Bit 2

PORT 29H : Define as Input.

Bit A7 : Conversion Bit 3
Bit A6 : Conversion Bit 4
Bit A5 : Conversion Bit 5
Bit A4 : Conversion Bit 6
Bit A3 : Conversion Bit 7
Bit A2 : Conversion Bit 8
Bit A1 : Conversion Bit 9
Bit A0 : Conversion Bit 10 (LSB)

3.) FTS Control Output Bit Assignments:

The following bit assignments are all defined as outputs. The multiplexer address lines are optically isolated from the analog system to prevent noise cross-talk. All output bits are buffered by 74LS240s. These bits are wired in parallel with the buttons on the control panel (manual).

PORT 00H

Bit A7 : Not Used
 Bit A6 : Not Used
 Bit A5 : Not Used
 Bit A4 : Not Used
 Bit A3 : Multiplxer Address Line A3
 Bit A2 : Multiplxer Address line A2
 Bit A1 : Multiplxer Address Line A1
 Bit A0 : Multiplxer Address Line A0

PORT 08H

Bit A7 : Not Used
 Bit A6 : Not Used
 Bit A5 : Not Used
 Bit A4 : Not Used
 Bit A3 : Not Used
 Bit A2 : Not Used
 Bit A1 : FTS Automatic Enable (0 = "enable")
 Bit A0 : FTS Manual Enable (0 = "enable")

PORT 09H

Bit B7 : Not Used
 Bit B6 : Pump (1 = "on", 0 = "off")
 Bit B5 : Valve 3 Out
 Bit B4 : Valve 3 In
 Bit B3 : Valve 2 Out
 Bit B2 : Valve 2 In
 Bit B1 : Valve 1 Out
 Bit B0 : Valve 1 In

4.) 33-Digit Density Display/Keyboard System:

The 33-digit density display system is controlled by the expansion 8279 keyboard/display controller. This controller has an internal 16 byte/32 digit display RAM that is written to and read from by the 8085A CPU. For a more complete understanding on how the 8279 works, refer to the Intel MCS-85 User's Manual #9800366E.

Display Section:

The 8279 can be considered as memory mapped I/O. It is communicated with via the CPU as if it were a simple memory location. The two addresses of concern are:

Command Address : B801H

Data Address : B800H

Each byte in the display RAM is configured as : BBBBAAAA, where the four hi-order bits make up the B section and the four lo-order bits make up the A section. The hardware is set up to display each section as a BCD digit. Therefore, sixteen bytes of RAM with two nibbles each give a total of thirty-two digits. The physical display addresses are shown in the following list. The top left most digit is hard-wired to always read zero.

0. A0 A1

A2 A3 A4

A5 A6 A7

A8 A9 AA

AB AC AD

AE AF B0

B1 B2 B3

B4 B5 B6

B7 B8 B9

BA BB BC

BD BE BF

Digit addresses in relation to
their physical position.
(top view).

Keyboard Section:

The 8279 keyboard section is used in the encoded scan mode with 2-key lock out. When reading the FIFO, keystroke addresses are mapped as shown below:

CFH	CBH	C7H	C3H	
CEH	CAH	C6H	C2H	Keystroke addresses in
CDH	C9H	C5H	C1H	relation to actual key
CCH	C8H	C4H	C0H	pad location (top view).

5.) Twelve Character Alphanumeric display:

The twelve character 5x7 dot matrix alphanumeric display is controlled by the MTX-A1 display controller. The MTX-A1 is also set up as memory mapped I/O and therefore is treated as a single memory location. All commands and display data are written to what is known as the control address.

Control Address : BF00H

The MTX-A1 also contains a 32-character RAM, and is very similar in operation to the 8279. For complete details of operation, see the MTX-A1 User's Manual.

The relationship between character address and physical location is shown below:

B A 9 8 7 6 5 4 3 2 1 0

The display is written to from right to left and any excess characters roll off the left end and are lost.

6.) Five Digit Control Display:

These five display are simply wired in parallel with the displays mounted on the SDK-85. There is a small eight position

DIP switch located on the SDK-85 board that is used to switch between the two display systems. These displays are controlled by the Basic 8279 Display Controller. See Chapter Five of the SDK-85 System Design Kit User's Manual #9800451B for a very complete description on how these displays are used.

Data Address : 1800H

Command Address : 1900H

7.) I/O Port and Memory Map:I/O Map:

<u>Port Address</u>	<u>Port Name</u>	<u>Chip</u>
00H	Basic ROM Port A	A14
01H	Basic ROM Port B	A14
02H	Basic ROM DDRA	A14
03H	Basic ROM DDRB	A14
08H	Expansion ROM Port A	A15
09H	Expansion ROM Port B	A15
0AH	Expansion ROM DDRA	A15
0BH	Expansion ROM DDRB	A15
20H	Basic RAM CSR	A16
21H	Basic RAM Port A	A16
22H	Basic RAM Port B	A16
23H	Basic RAM Port C	A16
24H	Basic Lo-Order Timer Control	A16
25H	Basic Hi-Order Timer Control	A16
28H	Expansion RAM CSR	A17
29H	Expansion RAM Port A	A17
2AH	Expansion RAM Port B	A17
2BH	Expansion RAM Port C	A17
2CH	Expansion Lo-Order Timer Control	A17
2DH	Expansion Hi-Order Timer Control	A17

Memory Map:

<u>Active Address Range</u>	<u>Selected Device</u>
0000-07FF	8755 Basic UVEPROM
0800-0FFF	8755 Expansion UVEPROM
1000-17FF	Not Used
1800-1FFF	8279 Basic Controller
2000-27FF	8155 Basic RAM
2800-2FFF	8155 Expansion RAM
3000-7FFF	Not Used
8000-87FF	2142 RAMs (4)
8800-8FFF	2114 RAMs (4)
9000-97FF	Not Used
9800-9FFF	2716 UVEPROM
A000-AFFF	2732 UVEPROM
B000-B7FF	2716 UVEPROM
B800-BBFF	8279 Expansion Controller
BC00-BFFF	MTX-A1 Display Controller
C000-FFFF	Not Used

See SDK-85 User's Manual for more complete details on memory allocation.

Appendix 2 : Operating System Program Listings :

The most difficult portion of this thesis involved the development of the operating system control software. This appendix contains a short description of what each program does, followed by the complete listing of the assembled programs. The code is Intel 8080/8085 assembly language. All programs were developed on an Intel Inteltec MDS-220 development system. All software calculations are performed by routines contained in the Intel Floating-Point Arithmetic Library (FPAL). The operations provided are addition, subtraction, multiplication, division value comparison, negation, clearing to zero, absolute value, conversion between decimal and binary floating-point numbers and conversion between floating-point and 32-bit signed integer formats. All operations are 32-bit single-precision (positive number range approximates 1.2×10^{-38} to 3.4×10^{38}). The 32-bit single-precision format is fully described in the 8080/8085 Floating-Point Arithmetic Library user's manual, available from Intel Corporation, Santa Clara, California, order # 9800452-03.

Control of the IR densitometer is divided between seven different programs or modules. Each module is developed and debugged separately from one another. Once all programs are debugged, they are linked together to form one complete program. For this system, the linked programs require about 10K of program memory to run.

These program listings are provided to allow qualified persons to change or up-grade the IR densitometer operating system when needed. If more information is required, please consult the numerous Intel manuals concerning this subject.

1. Module CALIB:

This module performs all the initializations required to operate the IR densitometer. For example, certain memory locations and registers will contain unknown variables upon first power up. This program then loads these locations with the proper contents. This program also contains the sub-routines that; scan the 11 emitter/detector pairs, manage the output of data onto the 33-digit display, and perform the data calibration functions.

2. Module START:

Module START is a warm start routine in that after a process has been run, certain parameters have to be re-initialized. START also contains all the routines that control the 40-column line printer.

3. Module DATAIN:

This module interacts directly with the densitometer user to program certain process information. All timing and rate information is entered through this program and is also checked for validity. The majority of message tables are stored within this module, along with all delay routines.

4. Module SUBPAK:

SUBPAK contains all the routines needed to operate the 16 pad keyboard. When the user depresses a key, subpak determines which key has been pressed and loads its address and value (if it is a numeric key) into a series of memory locations called KSTOR. From there, it is up to the calling program to determine what to do next.

5. Module GODEV:

This program starts by displaying "BEGIN" on the prompt display. From there it simply controls the timing of the Fluid Transport System, and the timing of scan acquisition. GODEV also contains a memory allocation routine used to store density data as it is generated by the A/D converter. All timing functions are based on a 1 Hz clock input to RST 7.5.

6. Module DFIND:

This module is entered when processing is completed. The user enters a time at which density data was acquired, and the routine will calculate the memory address at which that data is stored. A routine to convert minutes to seconds is also contained here.

7. Module UPORDN:

Control of the special function keys is defined here. This program functions very similar to DFIND in that memory addresses are calculated depending on what density data is desired. Exit from this program will return the user to START where a new process run can be made.

```

LOC  OBJ      LINE      SOURCE STATEMENT
1 ;Date 7/26/81 revision 5, 8/11/81 revision 6, 8/29/81 revision 7.
2 ;
3 ; The In-Process IR Densitometer is controlled by seven program
4 ; modules. The module names are: CALIB, START, DATAIN, SUBPAK,
5 ; CODEV, DFIND, and UPORDN.
6 ;
7 ;     *** All programs by Steven p. Cox ***
8 ;
9 ;*****
10 ;
11 ;     Main start routine for In-process IR Densitometer.
12 ;     This program prompts the user through all steps of the
13 ;     start up calibration routine. After calibration, program
14 ;     control is given to START and then to DATAIN routine.
15 ;
16 ;*****
17 ;
18 ;     NAME      CALIB
19 ;
20 ;     PUBLIC  COMM1,DATA1,COMM2,DATA2,ALPHA,PA1,PB1,PC1,RAMPT
21 ;     PUBLIC  RPA2,RPB2,RPA1,POS2,POS1,RPOS1,CLEARA,CLEARB,SCAN,DSHOW
22 ;     PUBLIC  DRAM,MXDRAM,TMPNUM,TDSTOR
23 ;
24 ;     EXTRN   XDSPY2,XSTOR1,XSTOR2,FSET,FERHND,FCLR,FLOAD,FQFB2D
25 ;     EXTRN   XDSPY3,SUM3,FLTDS,FSTOR,FDIV,FMUL,FNEG,FSUB,FADD,START
26 ;     EXTRN   DELAY,DWRITE,SUM2,DEB
27 ;
28 ;
29 ;     Below are the variable names used in naming I/O port, display
30 ;     command and data addresses.
31 ;
32 ;
1900 33      COMM1  EQU      1900H    ; Basic 8279 command address.
1800 34      DATA1 EQU      1800H    ; Basic 8279 data address.
8801 35      COMM2  EQU      8801H    ; Expansion 8279 command address.
8800 36      DATA2  EQU      8800H    ; Expansion 8279 data address.
0020 37      CS1    EQU      20H      ; Basic RAM C/S REGISTER.
0021 38      PA1    EQU      21H      ; BASIC RAM PORT A.
0022 39      PB1    EQU      22H      ; BASIC RAM PORT B.
0023 40      PC1    EQU      23H      ; BASIC RAM PORT C.
0024 41      TILOB  EQU      24H      ; LO-ORDER BYTE OF BASIC TIMER.
0025 42      TIHOB  EQU      25H      ; HI-ORDER BYTE OF BASIC TIMER.
0026 43      CS2    EQU      26H      ; EXPANSION RAM C/S REGISTER.
0029 44      PA2    EQU      29H      ; EXPANSION PORT A
002A 45      PB2    EQU      2AH      ; EXPANSION PORT B.
002B 46      PC2    EQU      2BH      ; EXPANSION PORT C.
002C 47      TILOB  EQU      2CH      ; LO-ORDER BYTE OF EXPANSION TIMER.
002D 48      TIHOB  EQU      2DH      ; HI-ORDER BYTE OF EXPANSION TIMER.
8F00 49      ALPHA  EQU      8F00H    ; MTX-A1 CONTROL ADDRESS.
0000 50      RPA1   EQU      00H      ; BASIC ROM PORT A.
0001 51      RPB1   EQU      01H      ; BASIC ROM PORT B.
0002 52      RDDA1  EQU      02H      ; BASIC ROM DDR FOR PORT A.
0003 53      Rddb1  EQU      03H      ; BASIC ROM DDR FOR PORT B.

```

112

LOC	OBJ	LINE	SOURCE STATEMENT
0008		54	RPA2 EQU 00H ; EXPANSION ROM PORT A.
0009		55	RPR2 EQU 09H ; EXPANSION ROM PORT B.
000A		56	RDDA2 EQU 0AH ; EXPANSION ROM DDR FOR PORT A.
000B		57	RDDR2 EQU 0BH ; EXPANSION ROM DDR FOR PORT B.
0000		58	DRAM EQU 0000H ; START ADDRESS OF DENSITY RAM.
0001		59	MDMAXT EQU 0001H ; START OF MEASURED D-MAX DENSITY RAM.
0017		60	MDOT EQU 0017H ; START OF MD0 RAM.
		61 ;	
		62	CSEG
		63 ;	
0000	31FF20	64	CALIB: LXI SP, 20FFH ; STACK AT END OF EXPANSION RAM.
0003	210019	65	LXI H, COMM1 ; LOAD DISPLAY COMMAND ADDRESS.
0006	3600	66	MVI M, 00H ; 8 DIGIT, LEFT ENTRY, 2 KEY LOCK-OUT.
0008	CD0000	E 67	CALL DELAY
000B	36DC	68	MVI M, 0DCH ; CLEAR CONTROL DISPLAY RAM TO ZEROS.
000D	2101R0	69	LXI H, COMM2
0010	36A3	70	MVI M, 0A3H ; TURN DENSITY DISPLAY OFF.
0012	CD0000	E 71	CALL DELAY
0015	3E40	72	MVI A, 40H
0017	D32D	73	OUT T2HOR ; PROGRAM HI-ORDER BYTE OF EXPANSION TIMER.
0019	3E06	74	MVI A, 06H
001B	D32C	75	OUT T2LOR ; PROGRAM LO-ORDER BYTE OF EXPANSION TIMER.
001D	3ECC	76	MVI A, 0CCH ; PC2=OUT PR2=IN PA2=IN
001F	D320	77	OUT CS2 ; SET-UP EXPANSION I/O & START TIMER.
0021	3E00	78	MVI A, 00H ; SET CONVERT LINE TO LO.
0023	D32R	79	OUT PC2
0025	3E03	80	MVI A, 03H ; SET UP I/O FOR BASIC RAM.
0027	D320	81	OUT CS1 ; PRINTER CONTROL LINES.
0029	3EFF	82	MVI A, 0FFH
002B	D30A	83	OUT RDDA2 ; SET ALL BIT ON EXPAN. ROM PORT A TO OUTS.
002D	D30R	84	OUT RDDR2 ; SAME FOR EXPAN. ROM PORT B.
002F	3EFD	85	MVI A, 0FDR
0031	D308	86	OUT RPA2 ; ENABLE FTS CONTROL TO AUTO MODE.
0033	AF	87	XRA A
0034	D309	88	OUT RPR2 ; TURN OFF ALL VALVES AND PUMP.
0036	3E0F	89	MVI A, 0FH ; SET BASIC ROM PORT A TO OUTPUT.
0038	D302	90	OUT RDDA1 ; THESE ARE MUX ADDRESS LINES.
003A	3E00	91	MVI A, 00H ; PRESET MUX ADDRESS LINES TO 0000R.
003C	D300	92	OUT RPA1
003E	2100RF	93	LXI H, ALPHA ; LOAD 12 CHARACTER CONTROL ADDRESS.
0041	36E1	94	MVI M, 0E1H ; LOAD ALL BLANKS TO DISPLAY RAM.
0043	CD0000	E 95	CALL DELAY
0046	36C3	96	MVI M, 0C3H
0048	3E93	97	MVI A, 93H ; INITIALIZE PARAMETERS USED IN KDSPY2.
004A	326500	D 98	STA POS2
004D	3E04	99	MVI A, 04H
004F	326400	D 100	STA POS1
0052	3E64	101	MVI A, 64H
0054	326600	D 102	STA RPOS1
0057	3EDF	103	MVI A, 0DFH
0059	320000	D 104	STA CLEARA ; INITIALIZE ENTRY CLEAR TO KDSPY2.
005C	320100	D 105	STA CLEARB ; INITIALIZE EXIT CLEAR FROM KDSPY2.
005F	CD0000	E 106	CALL DELAY
0062	219903	C 107	ECAL: LXI H, EMITM ; PROMPT " EMITTER CAL? ".

LOC	ORJ	LINE	SOURCE STATEMENT
0065	040D	108	MVI B, 0DB
0067	CD0000	E 109	CALL DWRITE
006A	CD0000	E 110	CALL KDSPY2 ; ENTER WILL BYPASS ROUTINE.
006D	3A0000	E 111	LDA SUM2
0070	FE00	112	CPI 00H
0072	CAAA00	C 113	JZ INIT
0075	3E01	114	MVI A, 01H ; ANYTHING ELSE WILL TURN ON EMITTER 1.
0077	4F	115	MOV C, A ; C REG IS MUX ADDRESS COUNTER.
0078	320000	E 116	SETIT: STA DER ; DER CONTAINS DISPLAY NUMBER.
007H	79	117	MOV A, C
007C	D300	118	OUT RPA1 ; TURN ON EMITTER ADDRESSED BY C.
007E	C5	119	PUSH B
007F	219003	C 120	LXI H, EMIT ; DISPLAY " EMITTER (DER) ".
0082	060A	121	MVI B, 0AH
0084	CD0000	E 122	CALL DWRITE
0087	3A0000	E 123	LDA DER
008A	C430	124	ADI 30H ; CONVERT DER TO ASCII.
008C	3200BF	125	STA ALPHA ; SEND TO ALPHA DISPLAY.
008F	CD0000	E 126	CALL KDSPY2 ; WAIT.
0092	C1	127	POP B
0093	3A0000	E 128	LDA DER
0096	3C	129	INR A ; INCREMENT NUMBER COUNTER
0097	0C	130	INR C ; INCREMENT ADDRESS COUNTER.
0098	FED3	131	CPI 0D3H ; COMPARE TO " C ".
009A	CA4200	C 132	JZ ECAL
009D	FE0A	133	CPI 0AR
009F	CAA500	C 134	JZ FIX
00A2	C37800	C 135	JMP SETIT
00A5	3ED1	136	FIX: MVI A, 0D1H
00A7	C37800	C 137	JMP SETIT
		138 ;	
		139 ;	INITIALIZE FPAL:
		140 ;	
00AA	3E00	141	INIT: MVI A, 00H ; RESET MUX LINES TO 0000R.
00AC	D300	142	OUT RPA1
00AE	010800	D 143	LXI B, FPR
00B1	C5	144	PUSR B
00B2	010000	145	LXI B, 00H
00B5	CD0000	E 146	CALL FSET
		147 ;	
00B8	21A503	C 148	GDMAI: LXI H, DMAI ; PROMPT " ENTER D-MAX "
00BB	060E	149	MVI B, 0EH
00BD	CD0000	E 150	CALL DWRITE ; WRITE MESSAGE.
00C0	CD0000	E 151	CALL KDSPY3 ; ENTER D-MAX CALIBRATION VALUE (ADMAI).
00C3	3A0000	E 152	LDA SUM3 ; SUM3 CONTAINS ENTRY AS XYZ, IMPLIED X.YZ
00C6	321A00	D 153	STA ADMAI ; STORE ADMAI AS 32-BIT INTEGER.
00C9	AF	154	XRA A
00CA	321B00	D 155	STA ADMAI+1
00CD	321C00	D 156	STA ADMAI+2
00D0	321D00	D 157	STA ADMAI+3
00D3	010800	D 158	LXI B, FPR ; NOW CONVERT TO 32-BIT FLOATER.
00D6	111A00	D 159	LXI D, ADMAI
00D9	CD0000	E 160	CALL FLTDS
00DC	CD0000	E 161	CALL FSTOR ; STORE ADMAI AS FLOATER.

LOC	OBJ	LINE	SOURCE STATEMENT
00DF	21BA03	C 162	RDMAX: LXI H, RDDEM ; PROMPT " READ D-MAX ".
00E2	060C	163	MVI B, 0CR
00E4	CD0000	E 164	CALL DWRITE
00E7	CD0000	E 165	CALL KDSPLY2 ; CONTINUE ON RECEIPT OF ENTER KEY.
00EA	210080	166	LXI H, DRAM ; INITIALIZE MIDRAM.
00ED	220200	D 167	SHLD MIDRAM
00F0	212600	D 168	LXI H, SLOPE ; INITIALIZE SLOPE POINTER.
00F3	225200	D 169	SHLD SLOPEP
00F6	CD0D01	C 170	CALL SCAN ; 1ST 22 BYTES OF MIDRAM NOW CONTAIN MDMAX.
00F9	21B003	C 171	CZERO: LXI H, DZERO ; PROMPT " READ ZERO ".
00FC	060B	172	MVI B, 0BH
00FE	CD0000	E 173	CALL DWRITE
0101	CD0000	E 174	CALL KDSPLY2 ; CONTINUE ON RECEIPT ON ENTER KEY.
0104	CD0D01	C 175	CALL SCAN ; 2ND 22 BYTES OF MIDRAM CONTAIN MD0.
0107	211780	176	CALSLP: LXI H, MD0T ; INITIALIZE START ADDRESS OF 22 BYTE MD0 RAM.
010A	222200	D 177	SHLD MD0P
010D	210180	178	LXI H, MDMAIT ; INITIALIZE START ADDRESS OF MDRAM.
0110	222400	D 179	SHLD MDMAIP
0113	7D	180	MXSLOP: MOV A, L ; AT THIS POINT L CONTAINS LOB OF MDMAIP.
0114	FE17	181	CPI 17H
0116	CA0A01	C 182	JZ FINCAL ; IF EQUAL, THEN HAVE CALCULATED 11 SLOPES.
0119	2A2200	D 183	LHLD MD0P
011C	7E	184	MOV A, M ; PLACE LOB OF MD0 INTO A REG.
011D	321E00	D 185	STA TFLOAT
0120	23	186	INX H
0121	7E	187	MOV A, M ; PLACE HOB (2 MSBITS) INTO A.
0122	321F00	D 188	STA TFLOAT+1
0125	AF	189	XRA A
0126	322000	D 190	STA TFLOAT+2
0129	322100	D 191	STA TFLOAT+3
012C	010000	D 192	LXI B, FPR
012F	111E00	D 193	COVRTS: LXI D, TFLOAT ; CONVERT MD0 TO FLOATER, STORE IN TFLOAT.
0132	CD0000	E 194	CALL FLTDS
0135	CD0000	E 195	CALL FSTOR
0138	2A2400	D 196	LHLD MDMAIP ; LOAD LOB OF MDMAX INTO A REG.
013B	7E	197	MOV A, M
013C	320400	D 198	STA TMPNUM
013F	23	199	INX H
0140	7E	200	MOV A, M
0141	320500	D 201	STA TMPNUM+1
0144	AF	202	XRA A
0145	320600	D 203	STA TMPNUM+2
0148	320700	D 204	STA TMPNUM+3
014B	110400	D 205	LXI D, TMPNUM ; B SHOULD STILL CONTAIN FPR.
014E	CD0000	E 206	CALL FLTDS ; CONVERT MDMAX TO FLOATER, LEAVE IN FAC.
0151	111E00	D 207	LXI D, TFLOAT
0154	CD0000	E 208	CALL FSUB ; FAC = MDMAX - MD0 .
0157	CD0000	E 209	CALL FSTOR ; STORE RESULT IN TFLOAT.
015A	111A00	D 210	LXI D, ADMAX
015D	CD0000	E 211	CALL FLOAD ; LOAD FAC WITH ENTERED ADMAX.
0160	111E00	D 212	LXI D, TFLOAT
0163	CD0000	E 213	CALL FDIV ; FAC = 1/M = ADMAX/(MDMAX - MD0).
0166	2A5200	D 214	LHLD SLOPEP ; LOAD HL WITH CONTENTS OF SLOPEP.
0169	EB	215	XCHG ; EXCHANGE WITH CONTENTS OF DE.

115

LOC	OBJ	LINE	SOURCE STATEMENT
016A	CD0000	E 216	CALL FSTOR ; STORE 4-BYTE SLOPE AT LOCATION POINTED
		217	; TO BY SLOPE POINTER (SLOPEP).
016D	2A5200	D 218	INKS: LHL SLOPEP ; INCREMENT SLOPE COUNTER BY 4.
0170	23	219	INX H
0171	23	220	INX H
0172	23	221	INX H
0173	23	222	INX H
0174	225200	D 223	SHLD SLOPEP
0177	2A2200	D 224	LHL MDOP ; INCREMENT MD0 POINTER BY 2.
017A	23	225	INX H
017B	23	226	INX H
017C	222200	D 227	SHLD MDOP
017F	2A2400	D 228	LHL MDMAXP ; INCREMENT MDMAX POINTER BY 2.
0182	23	229	INX H
0183	23	230	INX H
0184	222400	D 231	SHLD MDMAXP
0187	C31301	C 232	JMP MXSLOP ; GO SEE IF MORE SLOPES TO BE CALCULATED.
		233 ;	
		234 ;	AT THIS POINT ALL SLOPES ARE CALCULATED, BEADY TO MAKE DENSITY
		235 ;	
018A	C30000	E 236	FINCAL: JMP STABT ; GO TO ENTER DATE ROUTINE.
		237 ;	
		238 ;	SUBROUTINES FOLLOW:
		239 ;	
		240 ;	*****
		241 ;	
		242 ;	*** SCAN ***
		243 ;	
		244 ;	SCAN PERFORMS ONE COMPLETE SCAN OF THE 18 DENSITOMETER BY
		245 ;	READING EACH OF THE 11 IR DENSITOMETERS ONE AFTER ANOTHER.
		246 ;	THE ROUTINE ASSUMES THAT MXDRAM, A 16-BIT ADDRESS WHERE DENSITY
		247 ;	VALUES ARE TO BE STORED, HAS AREADY BEEN INITIALIZED TO SOME VALUE.
		248 ;	
		249 ;	*****
		250 ;	
018D	2101B8	251	SCAN: LXI H, COMM2 ; BLANK DENSITY DISPLAY.
0190	36DF	252	MVI H, 0DFH
0192	36A3	253	MVI H, 0A3H
0194	CD0000	E 254	CALL DELAY
0197	CD0000	E 255	CALL DELAY
019A	3E01	256	MVI A, 01H ; INITIALIZE SCAN COUNTER,
019C	325500	D 257	STA SCNCNT
019F	FE0C	258	CKGCNT: CPI 0CH ; HAVE 11 PASSES BEEN MADE?
01A1	CAD701	C 259	JZ BESET
01A4	D300	260	OUT WPA1 ; PLACE SCNCNT CONTENTS ON MUX ADDRESS LINES.
01A6	CD0000	E 261	CALL DELAY ; LEAVE THAT EMITTER/DETECTOR PAIR ON 2MSEC
01A9	3E01	262	MVI A, 01H ; BRING CONVERT LINE HI
01AB	D32B	263	OUT PC2
01AD	160F	264	MVI D, 0FH ; EXTEND PULSE JUST A BIT.
01AF	15	265	LP2: DCR B
01B0	C2AF01	C 266	JNZ LP2
01B3	3E00	267	MVI A, 00H ; BRING CONVERT LINE LOW, START CONVERSION.
01B5	D32B	268	OUT PC2
01B7	163F	269	BUSY: MVI D, 3FH ; WAIT JUST OVER 25 MICROSECONDS.

LOC	DBJ	LINE	SOURCE STATEMENT
01B9	15	270	LP1: DCR B
01BA	C2B901	C 271	JNZ LP1
01BD	2A0200	D 272	LHLD NIDRAM ; INCREMENT DRAM POINTER BY ONE.
01C0	23	273	INX H
01C1	DB29	274	IN PA2 ; BRING IN 1ST 8 BITS OF CONVERSION.
01C3	77	275	MOV M, A
01C4	23	276	INX H ; INCREMENT DRAM POINTER BY ONE MORE.
01C5	DB2A	277	IN PB2 ; BRING IN REMAINING 2 MSBITS.
01C7	E603	278	ANI 03H ; MASK OUT HI-ORDER 6 BITS.
01C9	77	279	MOV M, A ; STORE WHERE POINTED BY NIDRAM.
01CA	220200	D 280	SHLD NIDRAM ; UPDATE NIDRAM.
01CD	3A5500	D 281	LDA SCNCNT
01D0	3C	282	INR A ; INCREMENT SCAN COUNTER BY ONE.
01D1	325500	D 283	STA SCNCNT
01D4	C39F01	C 284	JMP CKGCNT
01D7	3E00	285	RESET: MVI A, 00H ; RESET MUX ADDRESS LINES TO ZERO.
01D9	D300	286	OUT RPA1
01DB	C9	287	RET
		288 ;	
		289 ;	*****
		290 ;	
		291 ;	*** DSHOW ***
		292 ;	
		293 ;	DSHOW TAKES ONE 22-BYTE BLOCK OF DENSITY MEASUREMENTS AND OPERATES ON
		294 ;	IT WITH A LINEAR REGRESSION CALIBRATION ROUTINE. THE RESULTS ARE THEN
		295 ;	DISPLAYED ON THE 33-DIGIT DISPLAY SYSTEM. DSHOW MUST BE PASSED RAMPT,
		296 ;	(RAM POINTER) TO TELL IT WHICH BLOCK OF DATA IS WANTED.
		297 ;	THIS START ADDRESS IS FOUND FROM: DSTART = (RAMPT X 22) + 1 + DRAM
		298 ;	
		299 ;	SOME CONVENTIONS USED IN THIS SYSTEM ARE:
		300 ;	
		301 ;	1.) DENSITY READINGS BEGIN AT THE TOP OF THE CHAMBER AND DEPEND.
		302 ;	
		303 ;	2.) THE LO-DENSITY PRODUCING END OF THE SENSI STRIP WILL BE
		304 ;	PLACED AT THE TOP OF THE CHAMBER.
		305 ;	
		306 ;	3.) VALUES WILL DISPLAYED AND PRINTED WITH THE LOWEST DENSITY
		307 ;	READINGS ON TOP AND THE HIGHEST AT THE BOTTOM.
		308 ;	
		309 ;	4.) THE LOWEST DENSITY READING OF ANY DATA BLOCK WILL MAP TO
		310 ;	THE LOWEST OR START ADDRESS OF THAT DATA BLOCK.
		311 ;	
		312 ;	*****
		313 ;	
01DC	210000	314	DSHOW: LXI H, 0000H ; ZERO HL REGISTER PAIR.
01DF	3A5400	D 315	LDA RAMPT
01E2	FE00	316	MX: CPI 00H
01E4	CAEF01	C 317	JZ OVER ; CALCULATE (RAMPT X 22)
01E7	111600	318	LXI D, 0016H ; LOAD DE WITH 22.
01EA	19	319	DAD D ; HL = HL + DE .
01EB	3D	320	DCR A
01EC	C3E201	C 321	JMP NI
01EF	EB	322	DVER: XCHG ; STORE (RAMPT X 22) IN DE.
01F0	210000	323	LXI H, DRAM ; LOAD HL WITH DENSITY START ADDRESS.

117

LOC	OBJ	LINE	SOURCE STATEMENT
01F3	23	324	IMX M ; HL = DRAM + 1.
01F4	19	325	DAD D ; HL = DRAM + 1 + (RAMPT X 22) .
01F5	225400	D 326	SHLD DSTART
01F8	211780	327	LXI M, MDOT ; RE-INITIALIZE MD0 POINTER FOR USE IN CORECT.
01FB	222200	D 328	SHLD MD0P
01FE	212400	D 329	LXI M, SLOPE ; RE-INITIALIZE SLOPE POINTER FOR USE IN SAME.
0201	225200	D 330	SHLD SLOPEP
0204	2101B8	331	LXI M, COMM2 ; LOAD 8279 COMMAND ADDRESS.
0207	3490	332	MVI M, 90H ; SET DISPLAY RAM TO 0000, A1.
0209	CD0000	E 333	INHIBB: CALL DELAY
020C	34A5	334	MVI M, 0A5H ; BLANK AND WRITE INHIBIT B SEC OF DISPLAY.
020E	CD0000	E 335	CALL DELAY
0211	CD4C02	C 336	CALL CORECT ; FILL LO-DEM STEP, ACCOUNT FOR FIXED ZERO.
0214	2100B8	337	LXI M, DATA2
0217	3A5900	D 338	LDA TDSTOR+1
021A	77	339	MOV M, A
021B	CD0000	E 340	CALL DELAY
021E	3A5A00	D 341	LDA TDSTOR+2
0221	77	342	MOV M, A
0222	CD0000	E 343	CALL DELAY
0225	CD7C03	C 344	CALL DSEND
0228	CD7C03	C 345	CALL DSENB
022B	CD7C03	C 346	CALL DSEND
022E	CD7C03	C 347	CALL DSENB
0231	CD4C02	C 348	CALL CORECT
0234	2100B8	349	LXI M, DATA2 ; ACCOUNT FOR SWITCH BETWEEN SEC A AND SEC B.
0237	3A5800	D 350	LDA TDSTOR
023A	77	351	MOV M, A
023B	CD0000	E 352	CALL DELAY
023E	3A5900	D 353	LDA TDSTOR+1
0241	77	354	MOV M, A
0242	CD0000	E 355	CALL DELAY
0245	2101B8	356	LXI M, COMM2
0248	3490	357	MVI M, 90H ; SET RAM TO 0000.
024A	CD0000	E 358	INHIBA: CALL DELAY
024D	34A8	359	MVI M, 0A0H ; WRITE INHIBIT A SECTION OF RAM.
024F	CD0000	E 360	CALL DELAY
0252	2100B8	361	LXI M, DATA2
0255	3A5A00	D 362	LDA TDSTOR+2
0258	77	363	MOV M, A
0259	CD0000	E 364	CALL DELAY
025C	CD7C03	C 365	CALL DSEND
025F	CD7C03	C 366	CALL DSENB
0262	CD7C03	C 367	CALL DSEND
0265	CD7C03	C 368	CALL DSENB
0268	CD7C03	C 369	CALL DSEND
026B	C9	370	RET
		371 ;	
		372 ;	*****
		373 ;	
		374 ;	*** CORECT ***
		375 ;	
		376 ;	THIS ROUTINE RECEIVES AN ADDRESS PASSED TO IT CALLED DSTART WHICH IS
		377 ;	THE START ADDRESS OF A 22 BYTE BLOCK OF DATA CONTAINING 11 DENSITY

LOC	OBJ	LINE	SOURCE STATEMENT
		378 ;	READINGS. CORECT THEM TAKES A 2-BYTE, 10-BIT CONVERSION AND PERFORMS
		379 ;	A LINEAR REGRESSION CORRECTION ON THE VALUE. THE RESULT IS CONVERTED
		380 ;	TO SIMPLE BCD FOR DISPLAY, AND RETURNED IN TDSTOR AS :
		381 ;	
		382 ;	TDSTOR+2 LSD X
		383 ;	TDSTOR+1 Y
		384 ;	TDSTOR MSD X " X.YZ " FORMAT
		385 ;	
		386 ;	THE EXPONENT RETURNED FROM FQFB2D (FPAL ROUTINE) IS TESTED FOR
		387 ;	RANGE AND IF IN ERROR, ERROR CHARACTERS ARE RETURNED TO CALLING PROGRAM.
		388 ;	
		389 ;	*****
		390 ;	
026C	2A2200	D 391	CORECT: LHLD MD0P
026F	7E	392	MOV A, M ; CONVERT MD0 TO 32-BIT FLOATER.
0270	321E00	D 393	STA TFLOAT
0273	23	394	INX H
0274	7E	395	MOV A, M
0275	321F00	D 396	STA TFLOAT+1
0278	AF	397	XRA A
0279	322000	D 398	STA TFLOAT+2
027C	322100	D 399	STA TFLOAT+3
027F	010800	D 400	LXI B, FPR ; CONVERT
0282	111E00	D 401	LXI D, TFLOAT
0285	CD0000	E 402	CALL FLTDS
0288	CD0000	E 403	CALL FSTOR ; TFLOAT = FLOATER MD0.
028B	2A5600	D 404	FLOATS: LHLD BSTART
028E	7E	405	MOV A, M
028F	320400	D 406	STA TMPNUM ; CONVERT MD TO 32-BIT FLOATER.
0292	23	407	INX H
0293	7E	408	MOV A, M
0294	320500	D 409	STA TMPNUM+1
0297	AF	410	XRA A
0298	320600	D 411	STA TMPNUM+2
029B	320700	D 412	STA TMPNUM+3
029E	110400	D 413	LXI D, TMPNUM
02A1	CD0000	E 414	CALL FLTDS ; CONVERTED MD RESIDES IN FAC.
02A4	111E00	D 415	DOCOR: LXI D, TFLOAT
02A7	CD0000	E 416	CALL FSUB ; FAC = (MD - MD0)
02AA	2A5200	D 417	LHLD SLOPEP
02AD	EB	418	XCHG
02AE	CD0000	E 419	CALL FMUL ; FAC = (MD - MD0) X 1/SLOPE.
02B1	3E03	420	RESULT: MVI A, 03H ; CONVERT RESULT TO DECIMAL FORMAT.
02B3	325E00	D 421	STA DLNGTH
02B6	216100	D 422	LXI H, DECML
02B9	225F00	D 423	SHLD DADDB
02BC	010800	D 424	LXI H, FPR
02BF	115B00	D 425	LXI D, DSIGN
02C2	CD0000	E 426	CALL FQFB2D
02C5	8A5B00	D 427	SCLTST: LDA DSIGN ; TEST FOR NEGATIVE RESULT.
02C8	FE2B	428	CPI 2BH ; "2BH" IS ASCII FOR " + ".
02CA	C2DF02	C 429	JNZ ITNEC
02CD	3A5C00	D 430	LDA DSCALE ; GET DECIMAL EXPONENT.
02D0	FEFE	431	CPI 0FEH

LOC	OBJ	LINE	SOURCE STATEMENT
02D2	CAED02	C 432	JZ SHFT4 ; SHIFT POINT 4 PLACES.
02D5	FEFF	433	CPI OFFH
02D7	CAFF02	C 434	JZ SHFT3 ; SHIFT POINT 3 PLACES.
02DA	FE00	435	CPI 00H
02DC	CA1603	C 436	JZ SHFT2 ; SHIFT POINT 2 PLACES.
02DF	3E00	437	ITNEG: MVI A, 00H ; LOAD TDSTORS WITH ERROR CODE FOR DISPLAY.
02E1	325800	D 438	STA TDSTOR ; WHICH IS ALL ZEROS.
02E4	325900	D 439	STA TDSTOR+1
02E7	325A00	D 440	STA TDSTOR+2
02EA	C32E03	C 441	JMP INCREM ; GO TO INCREMENT ROUTINE.
02ED	AF	442	SHFT4: XRA A ; LOAD 1ST & 2ND PLACES W/ZERO.
02EE	325800	D 443	STA TDSTOR
02F1	325900	D 444	STA TDSTOR+1
02F4	3A6100	D 445	LDA DECML ; GET LSD.
02F7	D630	446	SUI 30H ; CONVERT FROM ASCII TO BCD.
02F9	325A00	D 447	STA TDSTOR+2
02FC	C32E03	C 448	JMP INCREM
02FF	AF	449	SHFT3: XRA A ; LOAD 1ST PLACE W/ ZERO ONLY.
0300	325800	D 450	STA TDSTOR
0303	3A6100	D 451	LDA DECML
0306	D630	452	SUI 30H
0308	325900	D 453	STA TDSTOR+1
030B	3A6200	D 454	LDA DECML+1
030E	D630	455	SUI 30H
0310	325A00	D 456	STA TDSTOR+2
0313	C32E03	C 457	JMP INCREM
0316	3A6100	D 458	SHFT2: LDA DECML
0319	D630	459	SUI 30H
031B	325800	D 460	STA TDSTOR
031E	3A6200	D 461	LDA DECML+1
0321	D630	462	SUI 30H
0323	325900	D 463	STA TDSTOR+1
0326	3A6300	D 464	LDA DECML+2
0329	D630	465	SUI 30H
032B	325A00	D 466	STA TDSTOR+2
		467 ;	
032E	2A5600	D 468	INCREM: LHLD DSTART ; INCREMENT ALL POINTERS USED IN THIS CALL.
0331	23	469	INX H
0332	23	470	INX H
0333	225600	D 471	SHLD DSTART
0336	2A2200	D 472	LHLD MD0F
0339	23	473	INX H
033A	23	474	INX H
033B	222200	D 475	SRLD MD0F
033E	2A3200	D 476	LHLD SLOPEF
0341	23	477	INX B
0342	23	478	INX H
0343	23	479	INX H
0344	23	480	INX H
0345	225200	D 481	SHLD SLOPEF
0348	3A5800	D 482	NIBBLE: LDA TDSTOR ; SET TDSTOR NIBBLE A = NIBBLE B.
034B	17	483	RAL
034C	17	484	RAL
034D	17	485	RAL

120

LOC	OBJ	LINE	SOURCE STATEMENT
034E	17	486	RAL
034F	E4F0	487	ANI 0F0H
0351	47	488	MOV B, A
0352	3A5000	D 489	LDA TDSTOR
0355	B0	490	ORA B
0356	325000	D 491	STA TDSTOR
0359	3A5900	D 492	LDA TDSTOR+1
035C	17	493	RAL
035D	17	494	RAL
035E	17	495	RAL
035F	17	496	RAL
0360	E4F0	497	ANI 0F0B
0362	47	498	MOV B, A
0363	3A5900	D 499	LDA TDSTOR+1
0366	B0	500	ORA B
0367	325900	D 501	STA TDSTOR+1
036A	3A5A00	D 502	LDA TDSTOR+2
036D	17	503	RAL
036E	17	504	RAL
036F	17	505	RAL
0370	17	506	RAL
0371	E4F0	507	ANI 0F0H
0373	47	508	MOV B, A
0374	3A5A00	D 509	LDA TDSTOR+2
0377	B0	510	ORA B
0378	325A00	D 511	STA TDSTOR+2
037B	C9	512	RET
		513 ;	
		514 ;	*****
		515 ;	
		516 ;	** DSEND **
		517 ;	
		518 ;	DSEND SIMPLY SENDS THE 3 BYTE OF BCD DATA CONTAINED IN TDSTOR TO
		519 ;	THE 33-DIGIT DISPLAY. IT IS UP TO THE CALLING PROGRAM TO SET THE
		520 ;	DISPLAY RAM ADDRESS TO WHICH THE DATA IS TO BE WRITTEN. AUTO
		521 ;	INCREMENT IS ALSO NOT ASSUMED.
		522 ;	
		523 ;	*****
		524 ;	
037C	CD6C02	C 525	DSEND: CALL CORECT
037F	2100B8	526	LXI H, DATA2
0382	3A5000	D 527	LDA TDSTOR
0385	77	528	MOV M, A
0386	CD0000	E 529	CALL DELAY
0389	3A5900	D 530	LDA TDSTOR+1
038C	77	531	MOV M, A
038D	CD0000	E 532	CALL DELAY
0390	3A5A00	D 533	LDA TDSTOR+2
0393	77	534	MOV M, A
0394	CD0000	E 535	CALL DELAY
0397	C9	536	RET
		537 ;	
		538 ;	DISPLAY MESSAGES FOLLOW:
		539 ;	

121

LOC	OBJ	LINE	SOURCE STATEMENT
0398	20	540	EMIT: DB 20H ; SPACE
0399	05	541	EMITM: DB 05H ; E
039A	0D	542	DB 0DH ; M
039B	09	543	DB 09H ; I
039C	14	544	DB 14H ; T
039D	14	545	DB 14H ; T
039E	05	546	DB 05H ; E
039F	12	547	DB 12H ; R
03A0	20	548	DB 20H ; SPACE
03A1	03	549	DB 03H ; C
03A2	01	550	DB 01H ; A
03A3	0C	551	DB 0CH ; L
03A4	3F	552	DB 3FH ; ?
03A5	05	553	DMAX: DB 05H ; E
03A6	0E	554	DB 0EH ; M
03A7	14	555	DB 14H ; T
03A8	05	556	DB 05H ; E
03A9	12	557	DB 12H ; R
03AA	20	558	DB 20H ; BLANK
03AB	04	559	DB 04H ; D
03AC	20	560	DB 20H ; BLANK
03AD	0D	561	DB 0DH ; M
03AE	01	562	DB 01H ; A
03AF	18	563	DB 18H ; X
03B0	20	564	DZERO: DB 20H ; BLANK
03B1	12	565	DB 12H ; R
03B2	05	566	DB 05H ; E
03B3	01	567	DB 01H ; A
03B4	04	568	DB 04H ; D
03B5	20	569	DB 20H ; BLANK
03B6	1A	570	DB 1AH ; Z
03B7	05	571	DB 05H ; E
03B8	12	572	DB 12H ; R
03B9	0F	573	DB 0FH ; O
03BA	12	574	RDDEN: DB 12H ; R
03BB	05	575	DB 05H ; E
03BC	01	576	DB 01H ; A
03BD	04	577	DB 04H ; D
03BE	20	578	DB 20H ; BLANK
03BF	04	579	DB 04H ; D
03C0	20	580	DB 20H ; BLANK
03C1	0D	581	DB 0DH ; M
03C2	01	582	DB 01H ; A
03C3	18	583	DB 18H ; X
03C4	20	584	TEST: DB 20H ; BLANK
03C5	04	585	DB 04H ; D
03C6	05	586	DB 05H ; E
03C7	0E	587	DB 0EH ; M
03C8	20	588	DB 20H ; BLANK
03C9	12	589	DB 12H ; R
03CA	05	590	DB 05H ; E
03CB	01	591	DB 01H ; A
03CC	04	592	DB 04H ; D
		593	;

LOC OBJ LINE SOURCE STATEMENT

122

```

594      DSEC
595 ;
0000    596 CLEARA: DS      01H      ; CLEAR CODE FOR ENTRY TO KDSPY2.
0001    597 CLEARB: DS      01H      ; CLEAR CODE FOR EXIT FROM KDSPY2.
0002    598 NIDRAM: DS      02H      ; STORAGE LOCATION FOR NEXT DEN RAM POINTER.
0004    599 THPNUM: DS      04H      ; TEMPORARY 32-BIT INTEGER STORAGE.
0006    600 FPR: DS         12H      ; 16 BYTE ALLOCATION FOR FLOATING POINT RECORD.
001A    601 ADMAX: DS      04H      ; FLOATER ADMAX.
001E    602 TFLOAT: DS     04H      ; FLOATER TEMP. STORAGE.
0022    603 MDOP: DS        02H      ; MEASURED ZERO DENSITY POINTER.
0024    604 MDMAXP: DS     02H      ; MEASURED D-MAX ADDRESS POINTER.
0026    605 SLOPE: DS        2CH      ; 44 BYTE BLOCK OF RAM FOR 11 SLOPE VALUES.
0032    606 SLOPEP: DS      02H      ; SLOPE ADDRESS POINTER STORAGE.
0034    607 RAMPT: DS       01H      ; RAM POINTER USED IN FINDING DATA BLOCKS.
0035    608 SCNCNT: DS      01H      ; SCAN COUNTER USED IN SCAN SUB-ROUTINE.
0036    609 DSTART: DS      02H      ; ADDRESS STORAGE FOR USE BY DSHOW AND CORECT.
0038    610 TDSTOR: DS      03H      ; TEMPORARY STORAGE FOR DISPLAY RESULTS.
003B    611 DSIGN: DS       01H      ; CONTROL BLOCK USED BY FQFB2D.
003C    612 DSCALE: DS     02H
003E    613 DLNGTH: DS     01H
003F    614 DADDR: DS      02H
0041    615 DECML: DS       03H
0044    616 POS1: DS        01H      ; PARAMETERS USED IN KDSPY2.
0045    617 POS2: DS        01H
0046    618 RPOS1: DS       01H
619 ;
620      END

```

PUBLIC SYMBOLS

ALPHA	A R000	CLEARA	D 0000	CLEARB	D 0001	COMM1	A 1900	COMM2	A R001	DATA1	A 1000	DATA2	A R000
DRAM	A 0000	DSHOW	C 01DC	NIDRAM	D 0002	PA1	A 0021	PR1	A 0022	PC1	A 0023	POS1	D 0044
POS2	D 0045	RAMPT	D 0054	RPA1	A 0000	RPA2	A 0000	RPR2	A 0009	RPOS1	D 0046	SCAN	C 010D
TDSTOR	D 0050	THPNUM	D 0004										

EXTERNAL SYMBOLS

0ER	E 0000	DELAY	E 0000	DWRITE	E 0000	FADD	E 0000	FCLR	E 0000	FDIV	E 0000	FERHND	E 0000
FLOAD	E 0000	FLTDS	E 0000	FMUL	E 0000	FNEG	E 0000	FQFB2D	E 0000	FSET	E 0000	FSTOR	E 0000
FSUB	E 0000	KDSPY2	E 0000	KDSPY3	E 0000	KSTOR1	E 0000	KSTOR2	E 0000	START	E 0000	SUM2	E 0000
SUM3	E 0000												

USER SYMBOLS

ADMAX	D 001A	ALPHA	A R000	BUSY	C 01B7	CALIB	C 0000	CALSLP	C 0107	CKGCNT	C 019F	CLEARA	D 0000
CLEARB	D 0001	COMM1	A 1900	COMM2	A R001	CORECT	C 024C	COVRTS	C 012F	CS1	A 0020	CS2	A 0028
DADDR	D 005F	DATA1	A 1000	DATA2	A R000	DER	E 0000	DECML	D 0041	DELAY	E 0000	DLNGTH	D 005E
DMAX	C 03A5	DOCOR	C 02A4	DRAM	A 0000	DSCALE	D 005C	DSEND	C 037C	DSHOW	C 01DC	DSIGN	D 005B
DSTART	D 0056	DWRITE	E 0000	DZERO	C 03B0	ECAL	C 0042	EMIT	C 0390	EMITH	C 0399	FADD	E 0000
FCLR	E 0000	FDIV	E 0000	FERHND	E 0000	FINCAL	C 010A	FIX	C 00A5	FLOAD	E 0000	FLOATS	C 020B
FLTDS	E 0000	FMUL	E 0000	FNEG	E 0000	FPR	D 0000	FQFB2D	E 0000	FSET	E 0000	FSTOR	E 0000
FSUB	E 0000	GDMAX	C 00B0	GZERO	C 00F9	INCREM	C 032E	INHIRA	C 024A	INHIRE	C 0209	INIT	C 00AA
INKS	C 016D	ITNEG	C 02DF	KDSPY2	E 0000	KDSPY3	E 0000	KSTOR1	E 0000	KSTOR2	E 0000	LP1	C 01B9
LP2	C 01AF	MDOP	D 0022	MDOT	A 0017	MDMAXP	D 0024	MDMAXT	A 0001	NIRBLE	C 0340	NX	C 01E2
NIDRAM	D 0002	MXSLOP	C 0113	OVER	C 01EF	PA1	A 0021	PA2	A 0029	PR1	A 0022	PR2	A 002A
PC1	A 0023	PC2	A 002B	POS1	D 0044	POS2	D 0045	RAMPT	D 0054	RDDA1	A 0002	RDDA2	A 000A

123

RDDB1	A 0003	RDDB2	A 000B	RDDEN	C 03BA	RDMAZ	C 00DF	RESET	C 01D7	RESULT	C 02B1	RPA1	A 0000
RPA2	A 0008	RPE1	A 0001	RPE2	A 0009	RPOS1	D 0066	SCAN	C 018D	SCLTST	C 02C5	SCMCNT	D 0055
SETIT	C 0078	SHFT2	C 0316	SHFT3	C 02FF	SHFT4	C 02ED	SLOPE	D 0026	SLOPEP	D 0052	START	E 0000
SUM2	E 0000	SUM3	E 0000	T1HOB	A 0025	T1LOB	A 0024	T2HOB	A 002D	T2LOB	A 002C	TDSTOR	D 0058
TEST	C 03C4	TFLOAT	D 001E	TMPNUM	D 0004								

ASSEMBLY COMPLETE, NO ERRORS

LOC	OBJ	LINE	SOURCE STATEMENT
		1	;*****
		2	;
		3	Warn start routine for in-process IR Densitometer.
		4	Come here after calibration routine is complete.
		5	The first part of this routine prompts for the process
		6	date, then prints the date information on the 48 column
		7	printer.
		8	;
		9	;*****
		10	;
		11	NAME START
		12	;
		13	PUBLIC DWRITE,PSEND,DELAY,START
		14	;
		15	EXTRN COMM2,KDSPY2,KSTOR2,KSTOR1,ALPHA,COMM1,DATAIN,PA1,PB1,PC1
		16	;
		17	CSEG
		18	;
0000	210000	E 19	START: LXI H, COMM2 ; CLEAR DENSITY DISPLAY.
0003	36DF	20	MVI M, 0DFH
0005	36A3	21	MVI M, 0A3H
0007	21E900	C 22	LXI H, MONTH ; SEND "ENTER MONTH:" TO DISPLAY.
000A	060D	23	MVI B, 0DH ; LOAD # OF CHARACTERS.
000C	CDA400	C 24	CALL DWRITE ; WRITE MESSAGE.
000F	CD0000	E 25	CALL KDSPY2 ; GET KEYBOARD ENTRIES.
0012	3E10	26	MVI A, 10H ; CLEAR PRINTER CONTROLS.
0014	CDC300	C 27	CALL PSEND ; SEND COMMAND TO PRINTER BUFFER.
0017	210C01	C 28	LXI H, DATE ; DATE POINTS TO "PROCESS DATE:"
001A	0E12	29	MVI C, 12H ; LOAD CHARACTER COUNTER.
001C	0D	30	PMESS1: DCR C ; DECREMENT CHARACTER COUNTER.
001D	3E00	31	MVI A, 00H ; PREPARE TO TEST COUNTER.
001F	B9	32	CMF C ; IF C=0, ZERO FLAG WILL SET.
0020	CA2800	C 33	JZ DONE ; IF NO MORE TO SEND, JUMP TO DONE.
0023	7E	34	MOV A, M ; PUT CHARACTER POINTED TO BY HL INTO A.
0024	CDC300	C 35	CALL PSEND ; SEND IT TO PRINTER.
0027	23	36	INX H ; INCREMENT TEXT POINTER.
0028	C31C00	C 37	JMP PMESS1
002B	3A0000	E 38	DONE: LDA KSTOR2 ; SEND MONTH DATA TO PRINTER.
002E	C630	39	ADI 30H ; ADD 30 HEX TO GET ASCII CHARACTER.
0030	CDC300	C 40	CALL PSEND
0033	3A0000	E 41	LDA KSTOR1
0036	C630	42	ADI 30H
0038	CDC300	C 43	CALL PSEND
003B	3E2F	44	MVI A, 2FH ; SEND " / " TO PRINTER.
003D	CDC300	C 45	CALL PSEND
0040	210000	E 46	LXI H, ALPHA
0043	36E1	47	MVI M, 0E1H ; BLANK ALPHA DISPLAY.
0045	CDDC00	C 48	CALL DELAY
0048	21F500	C 49	LXI H, DAY ; PROMPT "ENTER DAY:"
004B	060C	50	MVI B, 0CH ; LOAD CHARACTER COUNTER.
004D	CDA400	C 51	CALL DWRITE ; WRITE TO DISPLAY.
0050	CD0000	E 52	CALL KDSPY2 ; GET 2 DIGIT KEY ENTRY.
0053	3A0000	E 53	LDA KSTOR2

125

LOC	OBJ	LINE	SOURCE STATEMENT
0054	C430	54	ADI 30H
0058	CDC300	C 55	CALL PSEND ; SEND DAY INFO TO PRINTER.
005B	3A0000	E 56	LDA KSTOR1
005E	C430	57	ADI 30H
0060	CDC300	C 58	CALL PSEND
0063	3E2F	59	MVI A, 2FH ; SEND " / " TO PRINTER.
0065	CDC300	C 60	CALL PSEND
0068	210000	E 61	LXI H, ALPHA
006B	36E1	62	MVI M, 0E1H ; BLANK ALPHA DISPLAY.
006D	210001	C 63	LXI H, YEAR ; PROMPT "ENTER YEAR:".
0070	060D	64	MVI B, 0DH ; LOAD CHARACTER COUNTER.
0072	CDA400	C 65	CALL DWRITE ; WRITE TO DISPLAY.
0075	CD0000	E 66	CALL KDSPLY2 ; GET 2 DIGIT YEAR ENTRY FROM KEYBOARD.
0078	3A0000	E 67	LDA KSTOR2 ; SEND YEAR DATA TO PRINTER.
007B	C430	68	ADI 30H
007D	CDC300	C 69	CALL PSEND
0080	3A0000	E 70	LDA KSTOR1
0083	C430	71	ADI 30H
0085	CDC300	C 72	CALL PSEND
0088	210000	E 73	FINISH: LXI H, COMM1 ; LOAD CONTROL DISPLAY COMMAND ADDRESS.
008B	36DC	74	MVI M, 0DCH ; BLANK THE CONTROL DISPLAY.
008D	3E17	75	MVI A, 17H ; " START TO PRINT " COMMAND.
008F	CDC300	C 76	CALL PSEND ; PRINT THE DATE READER.
0092	3E17	77	MVI A, 17H ; SPACE NOW 3 BLANK LINES ON PRINTER.
0094	CDC300	C 78	CALL PSEND
0097	3E17	79	MVI A, 17H
0099	CDC300	C 80	CALL PSEND
009C	3E17	81	MVI A, 17H
009E	CDC300	C 82	CALL PSEND
00A1	C30000	E 83	JMP DATAIN ; GO TO DATAIN WHICH CONTAINS DATA ENTRY.
		84 ;	
		85 ;	SUBROUTINES FOLLOW:
		86 ;	
		87	DWRITE: ; HL CONTAINS STARTING ADDRESS OF TEXT STRING TO BE SENT TO DISPLAY
		88	; B REGISTER CONTAINS THE STRING LENGTH. ALPHA DISPLAY WILL ALWAYS
		89	; BE SET TO BLANKS BEFORE MESSAGE IS WRITTEN.
		90	;
00A4	E5	91	PUSH H ; STORE COPY OF HL ON STACK.
00A5	210000	E 92	LXI R, ALPHA
00A8	36E1	93	MVI M, 0E1H
00AA	CDDC00	C 94	CALL DELAY
00AD	3680	95	MVI M, 00H ; SET CURSOR POSITION TO 0.
00AF	CDDC00	C 96	CALL DELAY
00B2	E1	97	POP H ; RESTORE COPY OF HL REGISTERS.
00B3	05	98	AGAIN: DCR B ; DECREMENT CHARACTER COUNTER.
00B4	3E00	99	MVI A, 00H ; PREPARE A FOR CHARACTER COUNT TEST.
00B6	B8	100	CMP B ; IF B=0, ZERO FLAG WILL SET.
00B7	C8	101	RZ ; IF B=0, THEN ALL DONE AND RETURN.
00B8	7E	102	MOV A, M ; GET CHARACTER POINTED TO BY HL.
00B9	320000	E 103	STA ALPHA ; SEND IT TO ALPHA DISPLAY.
00BC	23	104	INX H ; INCREMENT TEXT POINTER.
00BD	CDDC00	C 105	CALL DELAY
00C0	C3B300	C 106	JMP AGAIN ; GO DO IT AGAIN ASSHOLE!
		107 ;	

126

LOC	OBJ	LINE	SOURCE STATEMENT
		108	;
		109	PSEND: ; THIS ROUTINE IS USED TO SEND EITHER ONE BYTE OF DATA OR ONE
		110	; COMMAND TO THE 40 COLUMN LINE PRINTER. THE PRINTER BUSY STATUS
		111	; LINE IS POLLED AND A RETURN IS NOT MADE UNTIL A NOT BUSY STATUS
		112	; IS READ. THE BYTE TO BE SENT MUST BE PLACED IN THE A REGISTER
		113	; BEFORE CALLING PSEND. (PRINTSEND).
		114	;
00C3	0300	E 115	OUT LOW PA1 ; PUT A ONTO PRINTER DATA LINES.
00C5	3E01	116	MVI A, 01H ; BRING WR LINE HI.
00C7	D300	E 117	OUT LOW PB1
00C9	3E00	118	MVI A, 00H ; BRING WR LINE LOW.
00CB	D300	E 119	OUT LOW PB1
00CD	3E01	120	MVI A, 01H ; BRING WR LINE HI AGAIN.
00CF	0300	E 121	OUT LOW PB1
00D1	CDDC00	C 122	CALL DELAY
00D4	DB00	E 123	PLOOK: IN LOW PC1 ; BRING IN PRINTER STATUS LINE.
00D6	0F	124	RRC
00D7	0F	125	RRC ; ROTATE BUSY BIT INTO CARRY POSITION.
00D8	D0	126	RNC ; IF NOT BUSY, NO CARRY AND RETURN.
00D9	C3D400	C 127	JMP PLOOK ; OTHERWISE, KEEP LOOKING!
		128	;
		129	;
		130	DELAY: ; DELAY IS CALLED WHENEVER AN EXTERNAL DEVICE SUCH AS THE MTX-A1,
		131	; OR THE PRINTER REQUIRES EXCESS TIME TO SET ITSELF UP.
		132	;
00DC	1402	133	MVI D, 02H
00DE	1EFF	134	LOOP2: MVI E, 0FFH
00E0	1D	135	LOOP1: DCR E
00E1	C2E000	C 136	JNZ LOOP1
00E4	15	137	DCR D
00E5	C2DE00	C 138	JNZ LOOP2
00E8	C9	139	RET
		140	;
		141	;
		142	;
		143	MONTH: DB 05H ; E
00EA	0E	144	DB 0EH ; M
00EB	14	145	DB 14H ; T
00EC	05	146	DB 05H ; E
00ED	12	147	DB 12H ; R
00EE	20	148	DB 20H ; BLANK
00EF	0D	149	DB 0DH ; M
00F0	0F	150	DB 0FH ; O
00F1	0E	151	DB 0EH ; M
00F2	14	152	DB 14H ; T
00F3	08	153	DB 08H ; H
00F4	3A	154	DB 3AH ; :
00F5	20	155	DAY: DB 20H ; BLANK
00F6	05	156	DB 05H ; E
00F7	0E	157	DB 0EH ; M
00F8	14	158	DB 14H ; T
00F9	05	159	DB 05H ; E
00FA	12	160	DB 12H ; R
00FB	20	161	DB 20H ; BLANK

127

LOC	OBJ	LINE	SOURCE STATEMENT
00FC	04	162	DB 04H ; D
00FD	01	163	DB 01H ; A
00FE	19	164	DB 19H ; Y
00FF	3A	165	DB 3AH ; :
0100	20	166	YEAR: DB 20H ; BLANK
0101	05	167	DB 05H ; E
0102	0E	168	DB 0EH ; M
0103	14	169	DB 14H ; T
0104	05	170	DB 05H ; E
0105	12	171	DB 12H ; R
0106	20	172	DB 20H ; BLANK
0107	19	173	DB 19H ; Y
0108	05	174	DB 05H ; E
0109	01	175	DB 01H ; A
010A	12	176	DB 12H ; R
010B	3A	177	DB 3AH ; :
010C	00	178	DATE: DB 00H ; BLANK
010D	00	179	DB 00H ; BLANK
010E	00	180	DB 00H ; BLANK
010F	50	181	DB 50H ; P
0110	72	182	DB 72H ; R
0111	6F	183	DB 6FH ; O
0112	63	184	DB 63H ; C
0113	65	185	DB 65H ; E
0114	73	186	DB 73H ; S
0115	73	187	DB 73H ; S
0116	00	188	DB 00H ; BLANK
0117	44	189	DB 44H ; D
0118	61	190	DB 61H ; A
0119	74	191	DB 74H ; T
011A	65	192	DB 65H ; E
011B	3A	193	DB 3AH ; :
011C	00	194	DB 00H ; BLANK
		195 ;	
		196	END

PUBLIC SYMBOLS

DELAY C 00DC DWRITE C 00A4 PSEND C 00C3 START C 0000

EXTERNAL SYMBOLS

ALPHA E 0000	COMM1 E 0000	COMM2 E 0000	DATAIN E 0000	KDSPY2 E 0000	KSTOR1 E 0000	KSTOR2 E 0000
PA1 E 0000	PB1 E 0000	PC1 E 0000				

USER SYMBOLS

AGAIN C 00B3	ALPHA E 0000	COMM1 E 0000	COMM2 E 0000	DATAIN E 0000	DATE C 010C	DAY C 00F5
DELAY C 00DC	DONE C 002B	DWRITE C 00A4	FINISH C 0088	KDSPY2 E 0000	KSTOR1 E 0000	KSTOR2 E 0000
LOOP1 C 00E0	LOOP2 C 00DE	MONTH C 00E9	PA1 E 0000	PB1 E 0000	PC1 E 0000	PLOOK C 00D4
PMESS1 C 001C	PSEND C 00C3	START C 0000	YEAR C 0100			

ASSEMBLY COMPLETE, NO ERRORS

LOC	OBJ	LINE	SOURCE STATEMENT
		1	*****
		2	;
		3	THIS PROGRAM INTERACTS WITH THE IR DENSITOMETER USER TO
		4	ALLOW FOR THE ENTRY OF THE THREE PROCESS TIME LENGTHS.
		5	ALSO TO BE ENTERED WITH THIS PROGRAM ARE THE INTERVAL-
		6	RATE COMBINATIONS USED IN THE DENSITY SCAN MEASUREMENTS.
		7	USER INFORMATION WILL BE PRINTED ON THE LINE PRINTER AS
		8	IT IS ENTERED.
		9	;
		10	*****
		11	;
		12	NAME DATAIN
		13	;
		14	PUBLIC DATAIN, LDELAY, BINASC, DEB, IRCNT, PLBYTE
		15	PUBLIC IRBRAM, STABCT, MANFLC, MANMOD
		16	;
		17	EXTRN ALPHA, DELAY, PSEND, PWRITE, DWRITE, KDSPY2, KDSPY3
		18	EXTRN SUM2, SUM3, DEVNUM, PNTNUM, RPA2, CODEV
		19	;
		20	CSEC
		21	;
0000	3EFD	22	DATAIN: MVI A, 0FDH ; RESET FTS TO AUTO MODE.
0002	D300	23	OUT LOW RPA2
0004	AF	24	XRA A ; ZERO A REG.
0005	320000	25	STA DEB ; ZERO DEB (DATA ENTRY BYTE).
0008	320200	26	STA STAB ; ZERO SCAN TABULATION BYTE.
000B	210700	27	LXI H, IRCNT ; ZERO I/R COUNT.
000E	23	28	INX H
000F	77	29	MOV M, A
0010	23	30	INX H
0011	77	31	MOV M, A
0012	23	32	INX H
0013	77	33	MOV M, A
0014	210300	34	LXI H, PLBYTE ; ZERO PROCESS LENGTH BYTES.
0017	23	35	INX R
0018	77	36	MOV M, A
0019	23	37	INX H
001A	77	38	MOV M, A
001B	23	39	INX B
001C	77	40	MOV M, A
001D	322A00	41	STA MANFLC ; MANFLAG = 00H MEANS AUTO MODE.
0020	210703	42	LXI H, MANMOD ; PROMPT "MANUAL MODE?"
0023	040D	43	MVI B, 0DH ; LOAD TEXT COUNTER.
0025	CD0000	44	CALL DWRITE ; TO ALPHA DISPLAY.
0028	CD0000	45	CALL KDSPY2 ; SUM2 = 0 MEANS AUTO, ANYTHING ELSE = MANUAL.
002B	3A0000	46	LDA SUM2
002E	FE00	47	CPI 00H
0030	CA5300	48	JX FRONT ; IF NOT MANUAL, GO WITH NORMAL ENTRY ROUTINE.
0033	3EFE	49	MVI A, 0FEH ; ENABLE FTS CONTROL FOR MANUAL SELECT.
0035	D300	50	OUT LOW RPA2
0037	219303	51	LXI H, PNTMAN ; PRINT "MANUAL MODE SELECTED."
003A	0E22	52	MVI C, 22H ; LOAD TEXT COUNTER.
003C	CD0000	53	CALL PWRITE ; WRITE TO PRINTER.

LOC	OBJ	LINE	SOURCE STATEMENT
003F	3E17	54	MVI A, 17H ; PRINT BUFFER CONTENTS.
0041	CD0000	E 55	CALL PSEND
0044	3E17	56	MVI A, 17H ; SPACE TWO LINES.
0046	CD0000	E 57	CALL PSEND
0049	3E17	58	MVI A, 17H
004B	CD0000	E 59	CALL PSEND
004E	3EFF	60	MVI A, 0FFH
0050	322A00	D 61	STA MANFLC ; SET MANFLC = FFH TO MEAN MANUAL MODE SELECTED
0053	AF	62 FRONT:	XRA A ; ZERO A REG AGAIN.
0054	320100	D 63	STA ITAB ; ITAB MUST BE ZEROED EVERY PASS THRU I/R IN.
0057	3A0000	D 64	LDA DEB
005A	E60F	65	ANI 0FH ; I/R COUNTER MUST BE RESET TO 0.
005C	47	66	MOV B, A ; STORE COPY OF A REG.
005D	212B00	D 67	LXI H, STABCT ; FORM SUCESSIVE SCAN TABS, DEPENDING ON
0060	05	68	ADD L ; NUMBER OF PROCESSES DESIRED.
0061	6F	69	MOV L, A
0062	3A0200	D 70	LDA STAB
0065	77	71	MOV M, A
0066	78	72	MOV A, B ; RETREIVE A.
0067	C601	73	ADI 01H ; INCREMENT LO-ORDER NIBBLE OF DEB.
0069	320000	D 74	STA DEB ; UP-DATE DEB.
006C	FE04	75	CPI 04H ; TEST TO SEE IF PROCESS ENTRY IS DONE.
006E	CA0000	E 76	JZ CODEV ; IF ZERO SETS, JUMP TO START PROCESS ROUTINE.
0071	CD0000	E 77	CALL DEVNUM ; PROMPT "PROCESS [LO-DEB] = ?".
0074	CD0000	E 78 OMT:	CALL KDSPY3 ; GET THREE DIGIT PROCESS ENTRY.
0077	3A0000	E 79	LDA SUM3 ; TEST FOR ALLOWABLE PROCESS LENGTH.
007A	FE7F	80	CPI 7FH ; 127 MINUTES IS MAX. LENGTH FOR PROCESS.
007C	DA8A00	C 81	JC PASS ; IF CY SETS, THEN SUM3 < 7FH.
007F	3ECB	82	MVI A, 0CBH ; COME HERE IF GREATER, AND SET ERROR MODE.
0081	320000	E 83	STA ALPHA ; START ALPHA DISPLAY BLINKING.
0084	CD0000	E 84	CALL DELAY
0087	C37400	C 85	JMP OMT ; GO AND GET ANOTHER PROCESS ENTRY.
008A	210300	D 86 PASS:	LXI H, PLBYTE ; LOAD HL WITH STARTING PLB ADDRESS.
008D	3A0000	D 87	LDA DEB ; GET LO-ORDEB PORTION OF DATA ENTRY BYTE.
0090	E60F	88	ANI 0FH ; MASK OUT HI-ORDER NIBBLE.
0092	05	89	ADD L ; ADD LOB OF PLBYTE TO PROCESS ENTRY COUNT.
0093	6F	90	MOV L, A ; PLACE RESULT INTO L REG.
0094	3A0000	E 91	LDA SUM3 ; GET ENTERED PROCESS LENGTH FROM STORAGE.
0097	77	92	MOV M, A ; STORE VALUE IN PLBYTE(1,2,3).
0098	CD0000	E 93	CALL PNTNUM ; PRINT "PROCESS [LO-DEB] = [KSTOR1,2,3] MIN"
009B	3A0000	E 94	LDA SUM3
009E	FE00	95	CPI 00H ; TEST TO SEE IF ENTERED PROCESS LENGTH = 0.
00A0	CAE001	C 96	JZ MEET ; IF SUM3=0, NO NEED TO GO TO INTERVAL ENTRY.
		97	; SO RETURN FOR NEXT PROCESS ENTRY.
00A3	3A0200	D 98 RATEIN:	LDA STAB
00A6	FEBA	99	CPI 0BAH ; COMPARE CURRENT STAB TO MAX. SCAN 0.
00AB	CAA201	C 100	JZ SCNOUT ; IF STAB = MAX. 0, THEN PROMPT "OUT OF SCANS"
00AB	3A0000	D 101	LDA DEB ; GET DATA ENTRY BYTE(DEB)
00AE	C610	102	ADI 10H ; INCREMENT HI-ORDER NIBBLE(INTERVAL/RATE).
00B0	320000	D 103	STA DEB ; STORE UPDATED VERSION OF DEB.
00B3	E6F0	104	ANI 0F0H ; MASK OUT LO-ORDER NIBBLE.
00B5	FEA0	105	CPI 0A0H ; ALLOW ONLY 10 INTERVAL/RATE ENTERIES.
00B7	CAE001	C 106	JZ NEXT ; IF 10 I/R REACHED, GO TO NEXT PROCESS ENTRY.
00BA	CDE702	C 107	CALL IRVNUM ; PROMPT "INTERVAL [HI-DEB] = ?".

LOC	OBJ	LINE	SOURCE STATEMENT
00BD	CD0000	E 108	OMTA: CALL KDSFY2 ; GET TWO DIGIT INTERVAL ENTRY.
00C0	3A0000	E 109	LDA SUM2 ; TEST VALUE OF SUM2.
00C3	FE00	110	CPI 00H ; IS IT = TO 0?
00C5	CAE001	C 111	JZ NEXT ; IF ZERO, TEST FOR MANUAL MODE.
00C8	CDD301	C 112	CALL RAMCNT ; GET 2(SUM OF IRCNTS)
00CB	3A0000	D 113	LDA DEB ; IRBRAM = INTERVAL RATE BYTE RAM.
00CE	E6F0	114	ANI 0F0H ; MASK TO GET HI-ORDER NIBBLE ONLY.
00D0	0F	115	RRC ; ROTATE 3 TIMES TO GET IT INTO LOW POSITION.
00D1	0F	116	RRC
00D2	0F	117	RRC
00D3	00	118	ADD B ; A NOW = 2(HI-DEB) + 2(SUM OF IRCNTS).
00D4	210B00	D 119	LXI H, IRBRAM
00D7	05	120	ADD L
00D8	4F	121	MOV L, A ; HL = 2(HI-DEB) + 2(SUM OF IRCNTS) + IRBRAM
00D9	3A0000	E 122	LDA SUM2 ; SUM2 = LAST INTERVAL ENTRY.
00DC	77	123	MOV M, A ; STORE IT AT HL POINTER.
00DD	3A0100	D 124	LDA ITAB ; LOAD A WITH WITH RUNNING INTERVAL SUM, ITAB.
00E0	06	125	ADD M ; A = SUM2 + ITAB
00E1	320100	D 126	STA ITAB ; UPDATE ITAB.
00E4	47	127	MOV B, A ; STORE COPY OF ITAB IN B REG.
00E5	3A0000	E 128	LDA SUM3 ; SUM3 SHOULD STILL CONTAIN PLB(LO-DEB).
00E8	B8	129	CMP B ; A-SUM3 COMPARED TO B=ITAB.
00E9	DA0501	C 130	JC EDBLNK ; IF CY=1, A(B, SUM3<ITAB, ERROR.
00EC	C21B01	C 131	JNZ LESS ; IF NO CY AND NO ZERO, ITAB<PROCESS LENGTH.
00EF	3EE1	132	MVI A, 0E1H ; IF COME HERE, THEN ITAB=SUM3.
00F1	320000	E 133	STA ALPHA ; CLEAR ALPHA DISPLAY.
00F4	CD0000	E 134	CALL DELAY
00F7	21FE01	C 135	LXI H, PFULL ; HL POINTS TO PROCESS FULL MESSAGE.
00FA	060D	136	MVI B, 0DH ; LOAD CHARACTER COUNTER.
00FC	CD0000	E 137	CALL DWRITE ; WRITE TO ALPHA DISPLAY.
00FF	CD7403	C 138	CALL LDELAY ; PAUSE FOR 5 SEC.
0102	C31B01	C 139	JMP LESS ; GO TO RATE ENTRY ROUTINE.
0105	3A0000	E 140	EDBLNK: LDA SUM2 ; LAST INTERVAL ENTRY TO BE SUBTRACTED
0108	47	141	MOV B, A ; FROM ITAB, THEN ITAB UPDATED.
0109	3A0100	D 142	LDA ITAB
010C	90	143	SUB B
010D	320100	D 144	STA ITAB
0110	3ECB	145	MVI A, 0CBH ; CAUSE ALPHA DISPLAY TO BLINK.
0112	320000	E 146	STA ALPHA
0115	CD0000	E 147	CALL DELAY
0118	C3BD00	C 148	JMP OMTA ; GO AND GET ANOTHER INTERVAL ENTRY.
011B	CD3B03	C 149	LESS: CALL RATNUM ; PROMPT "RATE (HI-DEB) = ?"
011E	CD0000	E 150	OMTYA: CALL KDSFY2 ; GET 2-DIGIT RATE ENTRY.
0121	3A0000	E 151	LDA SUM2 ; PUT RATE ENTRY INTO A REG.
0124	FE00	152	CPI 00H ; RATE ENTRIES OF 00 ARE NOT ALLOWED.
0126	CA3001	C 153	JZ INVALID
0129	47	154	MOV B, A ; STORE COPY IN B REG.
012A	3E3C	155	MVI A, 3CH ; LOAD A WITH MAX. RATE (60 SCANS/MIN.).
012C	B0	156	CMP B ; A=MAX. RATE CMP TO B=ENTERED RATE.
012D	D23B01	C 157	JNC PASSA ; IF CY SETS, RATE>MAX. SO BLINK DISPLAY.
0130	3ECB	158	INVALID: MVI A, 0CBH ; LOAD BLINK CODE.
0132	320000	E 159	STA ALPHA
0135	CD0000	E 160	CALL DELAY
0138	C31E01	C 161	JMP OMTYA ; GO AND GET ANOTHER RATE ENTRY.

131

LOC	OBJ	LINE	SOURCE STATEMENT
013B	CDD301	C 162	PASSA: CALL RAMCNT
013E	3A0000	D 163	LDA DEB
0141	E6F0	164	ANI 0F0H ; MASK OUT LO-ORDER NIBBLE.
0143	0F	165	RRC ; ROTATE 3 TIMES TO GET INTO LOW POSITION.
0144	0F	166	RRC
0145	0F	167	RRC
0146	00	168	ADD B
0147	D601	169	SUI 01H ; SUBTRACT ONE TO GET ONE BELOW ILB.
0149	210B00	D 170	LXI H, IRBRAM
014C	05	171	ADD L ; ADD 2(SUM OF IRCNTS) - 1 + IRBRAM +2[HI-DEB]
014D	6F	172	MOV L, A ; HL = IRBRAM - 1 + 2[HI-DEB] +2(SUM OF IRCNTS)
014E	3A0000	E 173	LDA SUM2 ; SUM2 IS SPB[HI-DEB].
0151	77	174	MOV M, A ; STORE SUM2 (RATE) AT THIS HL LOCATION.
		175 ;	
		176 ;	NOW CALCULATE PRODUCT OF ILB[HI-DEB] X SPB[HI-DEB]. THAT IS,
		177 ;	INTERVAL LENGTH TIMES SCAN RATE EQUALS TOTAL SCANS GENERATED.
		178 ;	
0152	AF	179	XRA A ; ZERO A REG.
0153	57	180	MOV D, A
0154	5F	181	MOV E, A ; ZERO DE REG. PAIR.
0155	7E	182	MOV A, M ; A = SPB (LAST RATE ENTRY).
0156	23	183	INX H ; HL NOW POINTS TO ILB (LAST INTERVAL ENTRY)
0157	FE00	184	CPI 00H ; IF SCAN RATE ENTRY = 0, JMP TO SMORE.
0159	CA6401	C 185	JZ PASS2
015C	46	186	HERE: MOV B, M ; B = ILB.
015D	13	187	MORE: INX D ; DE IS MULTI. COUNTER.
015E	05	188	DCR B
015F	C25D01	C 189	JNZ MORE ; DE GETS INCREMENTED BY ILB.
0162	3D	190	DCR A
0163	C25C01	C 191	JNZ HERE ; DE GETS INCREMENTED BY ILB, SPB TIMES.
0166	3A0200	D 192	PASS2: LDA STAB
0169	FE00	193	CPI 00H ; IF SCAN TABULATION = 00H, JUMP OVER SMORE.
016B	CA7301	C 194	JZ PASS3
016E	13	195	SMORE: INX D
016F	3D	196	DCR A
0170	C26E01	C 197	JNZ SMORE ; DE GETS INCREMENTED STAB TIMES.
		198 ;	
		199 ;	DE NOW EQUALS THE # OF SCANS THAT WOULD BE GENERATED W/RATE ENTERED
		200 ;	
0173	BA	201	PASS3: CMP D ; A=0. IF D>0, THEN RATE ENTRY TOO BIG.
0174	C2B201	C 202	JNZ RATBIG
0177	43	203	MOV B, E ; E CONTAINS 1 BYTE SCAN TABULATION.
0178	3EBA	204	MVI A, 0BAH ; MAX. SCANS AVAILABLE = 186.
017A	B0	205	CMP B
017B	CA9301	C 206	JZ NOMORE ; IF =, PROMPT "OUT OF SCANS", UPDATE STAB.
017E	DAB201	C 207	JC RATBIG ; IF CY SETS, RATE TOO BIG.
0181	70	208	MOV A, B ; A = SCAN TABULATION.
0182	320200	D 209	STA STAB
0185	CD2202	C 210	PRINT: CALL PNTIR ; COME HERE IF LESS THAN MAXIMUM.
		211	; AND PRINT INTERVAL/RATE INFO ON PRINTER.
0188	3A0200	D 212	LDA STAB ; TEST TO SEE IF WE HAVE SCANS REMAINING.
018B	FEBA	213	CPI 0BAH ; COMPARE STAB TO MAX. NUMBER.
018D	BAA300	C 214	JC RATEIN ; CY WILL SET IF SCANS REMAINING.
0190	C3E001	C 215	JMP NEXT ; IF NONE LEFT, GO TO NEXT PROCESS ENTRY.

LOC	OBJ	LINE	SOURCE STATEMENT
0193	7B	216	NOMORE: MOV A, E ; PUT SCAN TABULATION INTO A REG.
0194	320200	D 217	STA STAB
0197	CD2202	C 218	CALL PMTIR ; PRINT I/R INFO.
019A	3A0000	D 219	LDA DEB
019D	C610	220	ADI 10H ; INCREMENT [HI-DEB]
019F	320000	D 221	STA DEB ; UPDATE DEB
01A2	E5	222	SCNOUT: PUSH H
01A3	210A02	C 223	LXI H, OUTINT ; START ADDRESS OF "OUT OF SCANS"
01A6	040D	224	MVI B, 0DH ; LOAD TEXT COUNTER.
01A8	CD0000	E 225	CALL DWRITE ; WRITE TO DISPLAY.
01AB	E1	226	POP H
01AC	CD7403	C 227	CALL LDELAY ; PAUSE FOR 5 SECONDS.
01AF	C3E001	C 228	JMP NEXT ; GO STORE IRCNTS THEN GO TO NEXT PROCESS.
01B2	211602	C 229	RATBIG: LXI H, TOOBIG ; START ADDRESS OF "RATE TOO BIG".
01B5	040D	230	MVI B, 0DH ; LOAD TEXT COUNTER.
01B7	CD0000	E 231	CALL DWRITE
01BA	CD7403	C 232	CALL LDELAY ; PAUSE FOR "5 SEC.
01BD	3ECB	233	MVI A, 0CBH ; LOAD BLINK CODE.
01BF	320000	E 234	STA ALPHA
01C2	CD0000	E 235	CALL DELAY
01C5	C31E01	C 236	JMP OMTYA ; GO AND GET ANOTHER RATE ENTRY.
01C8	3A2A00	D 237	MANTEST: LDA MANFLG ; TEST TO SEE IF IN MANUAL MODE.
01CB	FEFF	238	CPI 0FFH
01CD	CA0000	E 239	JZ CODEV ; RECALL MANFLG=FFH MEANS MANUAL MODE.
01D0	C35300	C 240	JMP FRONT ; SKIP MANUAL MODE TEST FOR NOW.
		241 ;	
		242 ;	*****
		243 ;	
		244 ;	THIS SUBROUTINE TAKES IRBRAM ADDRESS AND DOES THIS:
		245 ;	IRBRAM POINTER = 2(SUM OF IRCNTS) RETURNED IN B REG.
		246 ;	
		247 ;	*****
		248 ;	
01D3	AF	249	RAMCNT: XRA A ; ZERO A REG.
01D4	210700	D 250	LXI H, IRCNT
01D7	23	251	INX H
01D8	06	252	ADD M ; ADD PROCESS 1 I/R COUNT TO A.
01D9	23	253	INX H
01DA	06	254	ADD M ; ADD PROCESS 2 I/R COUNT TO A.
01DB	23	255	INX H
01DC	06	256	ADD M ; ADD PROCESS 3 I/R COUNT TO A.
01DD	07	257	RLC ; MULTIPLY SUM OF I/R COUNTS BY 2.
01DE	47	258	MOV H, A
01DF	C9	259	RET
		260 ;	
01E0	3A0000	D 261	NEXT: LDA DE0
01E3	E60F	262	ANI 0FH ; GET PROCESS COUNT.
01E5	210700	D 263	LXI H, IRCNT ; START ADDRESS OF I/R COUNT STORAGE.
01E8	05	264	ADD L
01E9	4F	265	MOV L, A ; HL NOW POINTS TO IRCNT [LO-DEB].
01EA	3A0000	D 266	LDA DE0
01ED	E6F0	267	ANI 0F0H ; GET I/R COUNT.
01EF	0F	268	RRC
01F0	0F	269	RRC

LOC	OBJ	LINE	SOURCE STATEMENT
01F1	0F	270	RRC
01F2	0F	271	RRC
01F3	FE00	272	CPI 00B
01F5	CAFA01	C 273	JZ STORE
01F8	D601	274	SUI 01B
01FA	77	275	STORE: MOV M, A ; STORE I/R COUNT WHEN HL POINTS.
01FB	C3C801	C 276	JMP NANTST
		277 ;	
		278 ;	MESSAGE TABLES FOLLOW.
		279 ;	
01FE	10	280	PFULL: DB 10H ; P
01FF	12	281	DB 12H ; R
0200	0F	282	DB 0FH ; O
0201	03	283	DB 03H ; C
0202	05	284	DB 05H ; E
0203	13	285	DB 13H ; S
0204	13	286	DB 13H ; S
0205	20	287	DB 20H ; SPACE
0206	06	288	DB 06H ; F
0207	15	289	DB 15H ; U
0208	0C	290	DB 0CH ; L
0209	0C	291	DB 0CH ; L
020A	0F	292	OUTINT: DB 0FH ; O
020B	15	293	DB 15H ; U
020C	14	294	DB 14H ; T
020D	20	295	DB 20H ; SPACE
020E	0F	296	DB 0FH ; O
020F	06	297	DB 06H ; F
0210	20	298	DB 20H ; SPACE
0211	13	299	DB 13H ; S
0212	03	300	DB 03H ; C
0213	01	301	DB 01H ; A
0214	0E	302	DB 0EH ; N
0215	13	303	DB 13H ; S
0216	12	304	TOOBIG: DB 12H ; R
0217	01	305	DB 01H ; A
0218	14	306	DB 14H ; T
0219	05	307	DB 05H ; E
021A	20	308	DB 20H ; SPACE
021B	14	309	DB 14H ; T
021C	0F	310	DB 0FH ; O
021D	0F	311	DB 0FH ; O
021E	20	312	DB 20H ; SPACE
021F	02	313	DB 02H ; B
0220	09	314	DB 09H ; I
0221	07	315	DB 07H ; G
		316 ;	
		317 ;	SUBROUTINES FOLLOW:
		318 ;	
		319 ;	*****
		320 ;	
		321 ;	THIS ROUTINE TAKES THE CURRENT CONTENTS OF THE HI-ORDER NIBBLE
		322 ;	OF DEB (INTERVAL/RATE COUNTER) AND PRINTS IT ON THE LINE PRINTER
		323 ;	IN THE FOLLOWING FORMAT:

LOC	OBJ	LINE	SOURCE STATEMENT
		324 ;	
		325 ;	" INTERVAL (HI-DEB) = 1LB (HI-DEB) MINUTES "
		326 ;	
		327 ;	" RATE (HI-DEB) = SPB (HI-DEB) SCANS/MINUTE "
		328 ;	
		329 ;	WHEN CALLED, IT IS ASSUMED THAT HL IS POINTING TO 1LB (HI-DEB).
		330 ;	
		331 ;	*****
		332 ;	
0222	ES	333	PMTIR: PUSH H ; SAVE COPY OF HL AS DESTROYED IN CALLS.
0223	21BA02 C	334	LXI H, D1TAB ; START ADDRESS OF "INTERVAL".
0224	0E0E	335	MVI C, 0EH ; LOAD TEXT COUNTER.
0228	CD0000 E	336	CALL PWRITE
022B	3A0000 D	337	LDA DER
022E	E6F0	338	ANI 0F0H ; MASK OUT LOW-ORDER NIBBLE.
0230	0F	339	RRC ; ROTATE HIGH TO LOW.
0231	0F	340	RRC
0232	0F	341	RRC
0233	0F	342	RRC
0234	CD2203 C	343	CALL BINASC ; CONVERT I/R COUNTER TO 2 ASCII CHARACTERS.
0237	7A	344	MOV A, D ; PUT 10'S PLACE INTO A REG.
0238	CD0000 E	345	CALL PSEND
023B	79	346	MOV A, C ; PUT 1'S PLACE INTO A REG.
023C	CD0000 E	347	CALL PSEND
023F	3E20	348	MVI A, 20H ; SEND A SPACE.
0241	CD0000 E	349	CALL PSEND
0244	3E3D	350	MVI A, 3DH ; SEND "=" TO PRINTER.
0246	CD0000 E	351	CALL PSEND
0249	3E20	352	MVI A, 20H ; SEND ANOTHER SPACE.
024B	CD0000 E	353	CALL PSEND
024E	E1	354	POP H ; RETRIEVE COPY OF RL.
024F	7E	355	MOV A, H ; HL STILL POINTS TO 1LB(HI-DEB).
0250	ES	356	PUSH H ; SAVE COPY AGAIN.
0251	CD2203 C	357	CALL BINASC ; CONVERT INTERVAL LENGTH TO ASCII.
0254	7A	358	MOV A, D ; GET 10'S.
0255	CD0000 E	359	CALL PSEND
0258	79	360	MOV A, C ; GET 1'S.
0259	CD0000 E	361	CALL PSEND
025C	21C602 C	362	LXI H, MINUTE ; LOAD START ADDRESS OF "MINUTES".
025F	0E0A	363	MVI C, 0AH ; LOAD TEXT COUNTER.
0261	CD0000 E	364	CALL PWRITE
0264	3E17	365	MVI A, 17H ; NOW PRINT PRINTER BUFFER CONTENTS.
0266	CD0000 E	366	CALL PSEND
		367 ;	
		368 ;	NOW PRINT RATE INFORMATION
		369 ;	
0269	21CE02 C	370	LXI H, D2TAB ; START ADDRESS OF "RATE".
026C	0E0E	371	MVI C, 0EH ; LOAD TEXT COUNTER.
026E	CD0000 E	372	CALL PWRITE
0271	3A0000 D	373	LDA DER
0274	E6F0	374	ANI 0F0H ; MASK OUT LO-ORDER NIBBLE.
0276	0F	375	RRC
0277	0F	376	RRC
0278	0F	377	RRC

LOC	OBJ	LINE	SOURCE STATEMENT
0279	0F	378	RRC
027A	CD2203	C 379	CALL BINASC ; CONVERT I/R COUNT TO ASCII.
027D	7A	380	MOV A, D ; GET 10'S.
027E	CD0000	E 381	CALL PSEND
0281	79	382	MOV A, C ; GET 1'S.
0282	CD0000	E 383	CALL PSEND
0285	3E20	384	MVI A, 20H ; SEND SPACE.
0287	CD0000	E 385	CALL PSEND
028A	3E3D	386	MVI A, 3DH ; SEND "=" TO PRINTER.
028C	CD0000	E 387	CALL PSEND
028F	3E20	388	MVI A, 20H ; SEND ANOTHER SPACE.
0291	CD0000	E 389	CALL PSEND
0294	E1	390	POP H
0295	2B	391	DCX H ; INCREMENT HJ TO SPBHI-DEBI.
0296	7E	392	MOV A, H ; GET SPB.
0297	CD2203	C 393	CALL BINASC ; CONVERT TO ASCII.
029A	7A	394	MOV A, B
029B	CD0000	E 395	CALL PSEND
029E	79	396	MOV A, C
029F	CD0000	E 397	CALL PSEND
02A2	21DA02	C 398	LXI H, SCMMIN ; START ADDRESS OF "SCAN/MINUTES".
02A5	0E0E	399	MVI C, 0EH ; LOAD TEXT COUNTER.
02A7	CD0000	E 400	CALL PWRITE
02AA	3E17	401	MVI A, 17H ; PRINT PRINTER BUFFER CONTENTS.
02AC	CD0000	E 402	CALL PSEND
02AF	3E17	403	MVI A, 17H ; SPACE ONE BLANK LINE.
02B1	CD0000	E 404	CALL PSEND
02B4	3E17	405	MVI A, 17H ; SPACE ONE BLANK LINE.
02B6	CD0000	E 406	CALL PSEND
02B9	C9	407	RET ; RETURN TO CALLING PROGRAM.
		408 ;	
		409 ;	MESSAGE TABLES FOLLOW:
		410 ;	
02BA	20	411	D1TAB: DB 20H ; SPACE
02BB	20	412	DB 20H ; SPACE
02BC	20	413	DB 20H ; SPACE
02BD	49	414	DB 49H ; I
02BE	4E	415	DB 4EH ; M
02BF	74	416	DB 74H ; T
02C0	65	417	DB 65H ; E
02C1	72	418	DB 72H ; R
02C2	76	419	DB 76H ; V
02C3	61	420	DB 61H ; A
02C4	6C	421	DB 6CH ; L
02C5	20	422	DB 20H ; SPACE
02C6	20	423	MINUTE: DB 20H ; SPACE
02C7	6D	424	DB 6DH ; H
02C8	69	425	DB 69H ; I
02C9	6E	426	DB 6EH ; M
02CA	75	427	DB 75H ; U
02CB	74	428	DB 74H ; T
02CC	65	429	DB 65H ; E
02CD	73	430	DB 73H ; S
02CE	20	431	D2TAB: DB 20H ; SPACE

136

LOC	OBJ	LINE	SOURCE STATEMENT
02CF	20	432	DB 20H ; SPACE
02D0	20	433	DB 20H ; SPACE
02D1	20	434	DB 20H ; SPACE
02D2	20	435	DB 20H ; SPACE
02D3	20	436	DB 20H ; SPACE
02D4	20	437	DB 20H ; SPACE
02D5	52	438	DB 52H ; R
02D6	61	439	DB 61H ; A
02D7	74	440	DB 74H ; T
02D8	65	441	DB 65H ; E
02D9	20	442	DB 20H ; SPACE
02DA	20	443	SCNMIN: DB 20H ; SPACE
02DB	73	444	DB 73H ; S
02DC	63	445	DB 63H ; C
02DD	61	446	DB 61H ; A
02DE	6E	447	DB 6EH ; M
02DF	73	448	DB 73H ; S
02E0	2F	449	DB 2FH ; /
02E1	6D	450	DB 6DH ; M
02E2	69	451	DB 69H ; I
02E3	6E	452	DB 6EH ; M
02E4	75	453	DB 75H ; U
02E5	74	454	DB 74H ; T
02E6	65	455	DB 65H ; E
		456 ;	
		457 ;	*****
		458 ;	
		459 ;	THIS ROUTINE TAKES THE VALUE OF THE HI-ORDER NIBBLE STORED AT
		460 ;	DEB AND DISPLAYS IT ON THE ALPHA DISPLAY AS:
		461 ;	
		462 ;	" INTERVAL [HI-DEB] = ? "
		463 ;	
		464 ;	CALLS DWRITE AND DELAY.
		465 ;	
		466 ;	*****
		467 ;	
02E7	3EE1	468	IBVNUM: MVI A, 0E1H ; BLANK ALPHA DISPLAY.
02E9	320000	E 469	STA ALPHA
02EC	CD0000	E 470	CALL DELAY
02EF	211903	C 471	LXI H, IBVTAB ; LOAD START ADDRESS OF MESSAGE.
02F2	060A	472	MVI B, 0AH ; LOAD TEXT COUNTER.
02F4	CD0000	E 473	CALL DWRITE ; WRITE TO DISPLAY.
02F7	3A0000	D 474	LDA DEB
02FA	E6F0	475	ANI 0F0H ; MASK OUT LO-ORDER NIBBLE.
02FC	0F	476	RRC ; ROTATE H0N INTO L0N.
02FD	0F	477	RRC
02FE	0F	478	RRC
02FF	0F	479	RRC
0300	C630	480	ADI 30H ; CONVERT TO ASCII.
0302	320000	E 481	STA ALPHA
0305	CD0000	E 482	CALL DELAY
0308	3E3D	483	MVI A, 3DH ; SEND "=" TO ALPHA DISPLAY.
030A	320000	E 484	STA ALPHA
030D	CD0000	E 485	CALL DELAY

LOC	OBJ	LINE	SOURCE STATEMENT
0310	3E3F	486	MVI A, 3FH ; SEND "?".
0312	320000	E 487	STA ALPHA
0315	CD0000	E 488	CALL DELAY
0318	C9	489	RET
		490 ;	
		491 ;	MESSAGE TABLE FOLLOWS:
		492 ;	
0319	09	493	IRVTAB: DB 09H ; I
031A	0E	494	DB 0EH ; M
031B	14	495	DB 14H ; T
031C	05	496	DB 05H ; E
031D	12	497	DB 12H ; R
031E	16	498	DB 16H ; V
031F	01	499	DB 01H ; A
0320	0C	500	DB 0CH ; L
0321	20	501	DB 20H ; SPACE
		502 ;	
		503 ;	*****
		504 ;	
		505 ;	THIS ROUTINE CONVERTS A SINGLE BYTE BINARY NUMBER INTO TWO
		506 ;	4-BIT BCD VALUES.
		507 ;	
		508 ;	A REG: CONTAINS BINARY NUMBER TO BE CONVERTED.
		509 ;	D REG: RETURNS WITH BCD TENS PLACE IN ASCII.
		510 ;	C REG: RETURNS WITH BCD ONES PLACE IN ASCII.
		511 ;	
		512 ;	*****
		513 ;	
0322	1630	514	BINASC: MVI D, 30H ; LOAD D WITH ASCII ZERO.
0324	C601	515	ADI 01H ; OFF-SET NUMBER TO BE CONVERTED BY ONE
0326	060A	516	BEGIN: MVI B, 0AH ; LOAD B WITH LOOP COUNTER.
0328	3D	517	TENCMNT: DCR A ; SUBTRACT ONE FROM A REG.
0329	CA3403	C 518	JZ ONECNT ; IF A=0 JMP TO ONE'S PLACE.
032C	05	519	DCR B ; DECREMENT LOOP COUNTER.
032D	C22803	C 520	JNZ TENCMNT ; DO NOT LEAVE LOOP UNTIL CYCLED 10 TIMES.
0330	14	521	INR D ; IF THRU TEN TIMES, INCREMENT 10'S COUNTER.
0331	C32403	C 522	JMP BEGIN ; GO AND CYCLE AGAIN.
0334	3E0A	523	ONECNT: MVI A, 0AH ; DETERMINE HOW MANY ONE'S REMAIN.
0336	90	524	SUB B
0337	C630	525	ADI 30H ; CONVERT ONE'S TO ASCII.
0339	4F	526	MOV C, A ; STORE RESULT IN C REG.
033A	C9	527	RET
		528 ;	
		529 ;	*****
		530 ;	
		531 ;	THIS ROUTINE TAKES THE CURRENT VALUE OF HI-DEB AND DISPLAYS
		532 ;	IT ON THE ALPHA DISPLAY IN THE FORMAT:
		533 ;	
		534 ;	" RATE [HI-DEB] = ? "
		535 ;	
		536 ;	*****
		537 ;	
033B	3EE1	538	RATNUM: MVI A, 0E1H ; BLANK ALPHA DISPLAY.
033D	320000	E 539	STA ALPHA

LOC	OBJ	LINE	SOURCE STATEMENT
0340	CD0000	E 540	CALL DELAY
0343	216D03	C 541	LXI H, RATTAB ; START ADDRESS OF "RATE"
0346	0600	542	MVI B, 08H ; LOAD TEXT COUNTER.
0348	CD0000	E 543	CALL DWRITE ; WRITE TO DISPLAY.
034B	3A0000	D 544	LDA DEB ; GET I/R COUNT.
034E	E6F0	545	ANI 0F0H ; MASK OUT PROCESS COUNT (LOW).
0350	0F	546	RRC ; ROTATE NON INTO LOW POSITION.
0351	0F	547	RRC
0352	0F	548	RRC
0353	0F	549	RRC
0354	C430	550	ADI 30H ; CONVERT TO ASCII.
0356	320000	E 551	STA ALPHA
0359	CD0000	E 552	CALL DELAY
035C	3E3D	553	MVI A, 3DH ; SEND "=" TO ALPHA DISPLAY.
035E	320000	E 554	STA ALPHA
0361	CD0000	E 555	CALL DELAY
0364	3E3F	556	MVI A, 3FH ; SEND "!"
0366	320000	E 557	STA ALPHA
0369	CD0000	E 558	CALL DELAY
036C	C9	559	RET
		560 ;	
		561 ;	MESSAGE TABLE FOLLOWS:
		562 ;	
036D	20	563	RATTAB: DB 20H ; SPACE
036E	20	564	DB 20H ; SPACE
036F	12	565	DB 12H ; R
0370	01	566	DB 01H ; A
0371	14	567	DB 14H ; T
0372	05	568	DB 05H ; E
0373	20	569	DB 20H ; SPACE
		570 ;	
		571 ;	*****
		572 ;	
		573 ;	THIS DELAY ROUTINE GIVES "3 SECONDD PAUSE WHEN CALLED.
		574 ;	
		575 ;	*****
		576 ;	
0374	3E10	577	LDELAY: MVI A, 10H
0376	16FF	578	LOOOP3: MVI D, 0FFH
0378	1EFF	579	LOOOP2: MVI E, 0FFH
037A	1D	580	LOOOP1: DCR E
037B	C27A03	C 581	JNZ LOOOP1
037E	15	582	DCR D
037F	C27003	C 583	JNZ LOOOP2
0382	3D	584	DCR A
0383	C27403	C 585	JNZ LOOOP3
0386	C9	586	RET
		587 ;	
		588 ;	MESSAGE TABLES FOLLOW:
		589 ;	
0387	0D	590	MANMOD: DB 0DH ; M
0388	01	591	DB 01H ; A
0389	0E	592	DB 0EH ; H
038A	15	593	DB 15H ; U

LOC	OBJ	LINE	SOURCE STATEMENT
030B	01	594	DB 01H ; A
030C	0C	595	DB 0CH ; L
030D	20	596	DB 20H ; SPACE
030E	0D	597	DB 0DH ; M
030F	0F	598	DB 0FH ; O
0390	04	599	DB 04H ; D
0391	05	600	DB 05H ; E
0392	3F	601	DB 3FH ; ?
0393	20	602	PNTMAN: DB 20H ; SPACES
0394	20	603	DB 20H
0395	20	604	DB 20H
0396	20	605	DB 20H
0397	20	606	DB 20H
0398	20	607	DB 20H
0399	20	608	DB 20H
039A	2A	609	DB 2AH ; *
039B	2A	610	DB 2AH ; *
039C	20	611	DB 20H ; SPACE
039D	4D	612	DB 4DH ; M
039E	61	613	DB 61H ; A
039F	6E	614	DB 6EH ; N
03A0	75	615	DB 75H ; W
03A1	61	616	DB 61H ; A
03A2	6C	617	DB 6CH ; L
03A3	20	618	DB 20H ; SPACE
03A4	4D	619	DB 4DH ; M
03A5	6F	620	DB 6FH ; O
03A6	64	621	DB 64H ; D
03A7	65	622	DB 65H ; E
03A8	20	623	DB 20H ; SPACE
03A9	53	624	DB 53H ; S
03AA	65	625	DB 65H ; E
03AB	6C	626	DB 6CH ; L
03AC	65	627	DB 65H ; E
03AD	63	628	DB 63H ; C
03AE	74	629	DB 74H ; T
03AF	65	630	DB 65H ; E
03B0	64	631	DB 64H ; D
03B1	20	632	DB 20H ; SPACE
03B2	2A	633	DB 2AH ; *
03B3	2A	634	DB 2AH ; *
		635 ;	
		636	DSEC
		637 ;	
0000		638	DEB: DS 1H ; 1 BYTE OF RAM FOR DATA ENTRY BYTE.
0001		639	ITAB: DS 1H ; 1 BYTE OF RAM FOR INTERVAL TABULATION.
0002		640	STAB: DS 1H ; 1 BYTE OF RAM FOR SCAN TABULATION.
0003		641	PLBYTE: DS 4H ; 4 BYTE OF RAM FOR PROCESS LENGTH BYTES.
0007		642	IRCNT: DS 4H ; 4 BYTE OF RAM FOR # OF I/R ENTRIES PER PROCESS.
000B		643	IRBRAM: DS 1FH ; 30 BYTES OF RAM FOR I/R BYTES
002A		644	MANFLG: DS 1H ; MANUAL FLAG BYTE.
002B		645	STABCT: DS 04H ; 0=0,1=01 SCN TOTAL,2=#1+#2 SCN TOT,3=TOTAL SCANS.
		646 ;	
		647	END

LOC OBJ LINE SOURCE STATEMENT

140

PUBLIC SYMBOLS

BIMASC C 0322	DATAIN C 0080	DEB D 0000	IRBRAM D 000B	IRCNT D 0007	LDELAY C 0374	MANFLG D 002A
MANMOD C 0387	PLBYTE D 0003	STABCT D 002B				

EXTERNAL SYMBOLS

ALPHA E 0000	DELAY E 0000	DEVNUM E 0000	DWRITE E 0000	CODEV E 0000	KDSPY2 E 0000	KDSPY3 E 0000
PNTNUM E 0000	PSEND E 0000	PWRITE E 0000	RPA2 E 0000	SUM2 E 0000	SUM3 E 0000	

USER SYMBOLS

ALPHA E 0000	BEGIN C 0326	BIMASC C 0322	D1TAB C 02BA	D2TAB C 02CE	DATAIN C 8000	DEB D 0000
DELAY E 0000	DEVNUM E 0000	DWRITE E 0000	EDBLNK C 0105	FRONT C 0053	CODEV E 0000	HERE C 015C
INVALID C 0130	IRBRAM D 000B	IRCNT D 0007	IRVNUM C 02E7	IRVTAB C 0319	ITAB D 0001	KDSPY2 E 0000
KDSPY3 E 0000	LDELAY C 0374	LESS C 011B	LOOOP1 C 037A	LOOOP2 C 0378	LOOOP3 C 0376	MANFLG D 002A
MANMOD C 0387	MANTST C 01C8	MINUTE C 02C4	MOBE C 015D	NEXT C 01E0	NOMORE C 0193	OMT C 0074
OMTA C 00BD	OMTYA C 011E	ONECNT C 0334	OUTINT C 020A	PASS C 008A	PASS2 C 0166	PASS3 C 0173
PASSA C 013B	PFULL C 01FE	PLBYTE D 0003	PNTIR C 0222	PNTMAN C 0393	PNTNUM E 0000	PRINT C 0185
PSEND E 0000	PWRITE E 0000	RMCNT C 01D3	RATBIG C 01B2	RATEIN C 00A3	RATNUM C 033B	RATTAB C 036D
RPA2 E 0000	SCNMIN C 02DA	SCNOUT C 01A2	SMORE C 016E	STAB D 0002	STABCT D 002B	STORE C 01FA
SUM2 E 0000	SUM3 E 0000	TENCNT C 0328	TOOBIG C 0216			

ASSEMBLY COMPLETE, NO ERRORS

LOC	OBJ	LINE	SOURCE STATEMENT
		1	;*****
		2	;
		3	THIS PROGRAM MODULE CONTAINS MOST OF THE IMPORTANT
		4	SUBROUTINES USED BY THE IN-PROCESS IR DENSITOMETER.
		5	EACH SUBROUTINE IS PROCEEDED BY A BRIEF DESCRIPTION
		6	OF WHAT ITS FUNCTION IS AND THE PARAMETERS THAT MUST
		7	BE PASSED TO IT.
		8	;
		9	;*****
		10	;
		11	;
		12	NAME SUBPAK
		13	;
		14	PUBLIC KDSPY2,KDSPY3,KSTOR1,KSTOB2,KSTOR3,FMTNUM,DEVNUM
		15	PUBLIC PWRITE,SUM2,SUM3,LOOKUP,FIFOCH
		16	;
		17	EXTRN COMM1,DATA1,COMM2,DATA2,PSEND,DELAY,ALPHA,DWRITE
		18	EXTRN DEB,POS1,POS2,CLEARA,CLEARB,RPOS1
		19	;
		20	CSEC
		21	;
		22	;*****
		23	;
		24	THIS SUBROUTINE WILL ALLOW FOR THE ENTRY OF TWO NUMERIC DIGITS
		25	FROM THE CONTROL PANEL KEYBOARD. MISTAKES MADE ON THE KEYBOARD
		26	ARE CLEARED USING THE "CE" OR CLEAR ENTRY KEY. WHEN THE CORRECT
		27	DIGITS SELECTED BY THE USER ARE ON THE DISPLAY, THE "E" OR ENTER
		28	KEY WILL INPUT THE SELECTIONS TO THE COMPUTER. KEY ENTRIES ARE
		29	RIGHT ENTRY ONLY AND ARE RETURNED TO THE COMPUTER IN MEMORY
		30	LOCATIONS KSTOR1 AND KSTOR2.
		31	;
		32	;*****
		33	;
0000	CD0202	C 34	KDSPY2: CALL FIFOCH ; EMPTY THE FIFO.
0003	0400	35	MVI B, 00H ; SET DIGIT COUNTER TO 0.
0005	78	36	MOV A, B ; ZERO A REGISTER.
0006	320000	D 37	STA KSTOR1
0009	320100	D 38	STA KSTOR2 ; ZERO TEMP. KEY ENTRY REGISTERS.
000C	210000	E 39	LXI H, COMM1 ; WRITE ZEROS TO POS. 3 & 4 ON CONTROL DSP.
000F	3A0000	E 40	LDA CLEARA ; CLEAR DISPLAY RAM ACCORDING TO CLEAR CODE.
0012	77	41	MOV M, A
0013	CD0000	E 42	CALL DELAY
0016	3A0000	E 43	LDA POS2
0019	77	44	MOV M, A
001A	CD0000	E 45	CALL DELAY
001D	210000	E 46	LXI H, DATA1 ; 0279 DATA ADDRESS.
0020	3E0C	47	MVI A, 0CH ; LOAD DIGIT CODE FOR ZERO (0).
0022	77	48	MOV M, A
0023	77	49	MOV M, A ; WRITE 2 ZEROS TO DISPLAY.
0024	210000	E 50	KLOOK: LXI H, COMM2 ; STATUS WORD ADDRESS OF EXP. 0279.
0027	7E	51	MOV A, M ; PUT IT INTO A REG.
0028	1403	52	ANI 03H ; MASK OUT 5 HI-ORDER BITS.
002A	FE01	53	CPI 01H ; IF SOMETHING IN FIFO, ZERO WILL SET.

LOC	OBJ	LINE	SOURCE STATEMENT
002C	C22400	C 54	JNZ KLOOK ; IF NOT, KEEP LOOKING.
002F	3E40	55	MVI A, 40H ; SET FOR FIFO READ.
0031	77	56	MOV M, A ; SELECT FIFO AS READ SOURCE.
0032	210000	E 57	LXI H, DATA2 ; DATA ADDRESS OF EXPANSION 0279.
0035	7E	58	MOV A, M ; PUT FIFO CONTENTS INTO A REG.
0036	E60F	59	ANI 0FH ; MASK OUT HI-ORDER NIBBLE.
0038	FE00	60	CPI 00H ; IF CE KEY, ZERO FLAG WILL SET.
003A	CA0000	C 61	JZ KDSPLY2 ; START OVER IF CE KEY PUSHED.
003D	FE0F	62	CPI 0FH ; IF ENTER KEY, ZERO WILL SET.
003F	CAA200	C 63	JZ NOBLNK ; BEFORE RETURN, STOP ALPHA FROM BLINKING.
0042	21F301	C 64	LXI H, NUMTAB ; GET DECIMAL KEY VALUE.
0045	05	65	ADD L ; ADD KEY VALUE TO L REGISTER.
0046	4F	66	MOV L, A ; HL NOW POINTS TO ABSOLUTE KEY VALUE.
0047	7E	67	MOV A, M ; GET ACTUAL DECIMAL KEY VALUE.
0048	FEFF	68	CPI 0FFH ; IF ERROR KEY, ZERO FLAG WILL SET.
004A	CA2400	C 69	JZ KLOOK ; IF ERROR, GO AND READ KEYBOARD AGAIN.
004D	04	70	INR B ; INCREMENT DIGIT COUNTER.
004E	4F	71	MOV C, A ; STORE A IN C REG.
004F	3E01	72	MVI A, 01H ; PREPARE TO TEST DIGIT COUNTER.
0051	B8	73	CMP B ; IF DIGIT COUNTER=01, ZERO FLAG WILL SET.
0052	CA5E00	C 74	JZ STOR1 ; IF 01, GO TO KSTOR1 ROUTINE.
0055	3E02	75	MVI A, 02H ; TEST TO SEE IF = TO 02.
0057	B8	76	CMP B
0058	CA7300	C 77	JZ STOR2 ; IF ZERO SETS, GO TO KSTOR2 ROUTINE.
005B	C32400	C 78	JMP KLOOK ; MUST BE 3RD KEY ENTRY. NOT ALLOWED!!!!
005E	79	79	STOR1: MOV A, C ; RETRIVE A (DECIMAL KEY VALUE).
005F	320000	D 80	STA KSTOR1
0062	CDE201	C 81	CALL LOOKUP ; NOW GET VALUE FOR DISPLAY, RET IN C REG.
0065	210000	E 82	LXI H, COMM1
0068	3A0000	E 83	LDA POS1 ; DISPLAY POS. 1 NON-AI.
006B	77	84	MOV M, A
006C	210000	E 85	LXI H, DATA1
006F	71	86	MOV M, C ; SEND CHARACTER TO DISPLAY POS. 1.
0070	C32400	C 87	JMP KLOOK ; GO LOOK FOR NEXT KEY ENTRY.
0073	79	88	STOR2: MOV A, C ; RETRIVE A.
0074	320100	D 89	STA KSTOR2 ; STORE DECIMAL KEY VALUE IN KSTOR2.
0077	CDE201	C 90	CALL LOOKUP ; GET DISPLAY CHARACTER.
007A	210000	E 91	LXI H, COMM1
007D	3A0000	E 92	LDA RPOS1 ; READ DISPLAY RAM POS. 1
0080	77	93	MOV M, A
0081	210000	E 94	LXI H, DATA1
0084	56	95	MOV D, M ; READ VALUE INTO D REG.
0085	210000	E 96	LXI H, COMM1
0088	3A0000	E 97	LDA POS2 ; PREPARE TO WRITE TO POS. 2 (AI).
008B	77	98	MOV M, A
008C	210000	E 99	LXI H, DATA1
008F	72	100	MOV M, D ; WRITE POS. 1 TO POS. 2.
0090	71	101	MOV M, C ; NOW WRITE 2ND KEY ENTRY TO POS. 1
0091	3A0100	D 102	LDA KSTOR1
0094	57	103	MOV D, A
0095	3A0000	D 104	LDA KSTOR1
0098	320100	D 105	STA KSTOR1 ; REVERSE CONTENTS OF KSTOR 1 & 2.
009B	7A	106	MOV A, B
009C	320000	D 107	STA KSTOR1

143

LOC	OBJ	LINE	SOURCE STATEMENT
009F	C32400	C 100	JMP KLOOK ; GO LOOK FOR E OR CE KEY ENTRY.
00A2	210000	E 109	NOBLNK: LXI H, COMM1
00A5	3A0000	E 110	LDA CLEARB ; BLANK DISPLAY PER CODE IN CLRCOD.
00A8	77	111	MOV M, A
00A9	3EC3	112	MVI A, 0C3H ; CODE FOR NO-BLINK
00AB	320000	E 113	STA ALPHA ; CONTROL WORD.
00AE	CD0000	E 114	CALL DELAY
00B1	3A0000	D 115	LDA KSTOR1 ; WANT TO GET A SINGLE BYTE 0 FROM KSTOR1,2.
00B4	2E00	116	MVI L, 00H ; ZERO L AS L WILL CONTAIN RUNNING SUM.
00B6	0E01	117	MVI C, 01H ; C BECOMES 1'S FACTOR.
00B8	09	118	ONES2: DAD B ; ADD BC TO HL.
00B9	3D	119	DCR A ; DECREMENT ONES'S COUNTER.
00BA	FE00	120	CPI 00H ; IF COUNTER = 0, ZERO WILL SET.
00BC	C2B000	C 121	JNZ ONES2 ; IF NO ZERO, KEEP ADDING ONES.
00BF	3A0100	D 122	LDA KSTOR2 ; NOW DO IT FOR 10'S PLACE.
00C2	0E0A	123	MVI C, 0AH ; C IS NOW 10'S FACTOR.
00C4	09	124	TENS2: DAD B ; ADD BC TO RL.
00C5	3D	125	DCR A ; DECREMENT 10'S COUNTER.
00C6	FE00	126	CPI 00H
00C8	C2C400	C 127	JNZ TENS2
00CB	7D	128	MOV A, L ; PLACE RUNNING SUM INTO A REG.
00CC	320300	D 129	STA SUM2 ; STORE RUNNING SUM AT SUM2.
00CF	C9	130	RET ; RETURN TO MAIN PROGRAM.
		131 ;	
		132 ;*****	
		133 ;	
		134 ;	THIS IS KDSPY3. IT IS LIKE KDSPY2 EXCEPT THAT 3 KEY ENTRIES
		135 ;	ARE ALLOWED RATHER THAN ONLY TWO.
		136 ;	
		137 ;*****	
		138 ;	
00D0	CD0202	C 139	KDSPY3: CALL FIFOCR ; EMPTY THE FIFO.
00D3	0600	140	MVI B, 00H ; SET DIGIT COUNTER TO 0.
00D5	70	141	MOV A, B ; ZERO A REGISTER.
00D6	320000	D 142	STA KSTOR1
00D9	320100	D 143	STA KSTOR2
00DC	320200	D 144	STA KSTOR3 ; ZERO TEMP. KEY ENTRY REGISTERS.
00DF	210000	E 145	LXI H, COMM1 ; WRITE ZEROS TO POS. 2 & 3 & 4 ON CONTROL DSP.
00E2	36DF	146	MVI M, 0DFH
00E4	CD0000	E 147	CALL DELAY
00E7	3692	148	MVI M, 92H
00E9	210000	E 149	LXI H, DATA1 ; 0279 DATA ADDRESS.
00EC	3E0C	150	MVI A, 0CH ; LOAD DIGIT CODE FOR ZERO (0).
00EE	77	151	MOV M, A
00EF	77	152	MOV M, A
00F0	77	153	MOV M, A ; WRITE 3 ZEROS TO DISPLAY.
00F1	210000	E 154	KLOOK3: LXI H, COMM2 ; STATUS WORD ADDRESS OF EXP. 0279.
00F4	7E	155	MOV A, M ; PUT IT INTO A REG.
00F5	E603	156	ANI 03H ; MASK OUT 5 HI-ORDER BITS.
00F7	FE01	157	CPI 01H ; IF SOMETHING IN FIFO, ZERO WILL SET.
00F9	C2F100	C 158	JNZ KLOOK3 ; IF NOT, KEEP LOOKING.
00FC	3E40	159	MVI A, 40H ; SET FOR FIFO READ.
00FE	77	160	MOV M, A ; SELECT FIFO AS READ SOURCE.
00FF	210000	E 161	LXI H, DATA2 ; DATA ADDRESS OF EXPANSION 0279.

144

LOC	OBJ	LINE	SOURCE STATEMENT
0102	7E	162	MOV A, M ; PUT FIFO CONTENTS INTO A REG.
0103	E60F	163	ANI 0FH ; MASK OUT HI-ORDER NIBBLE.
0105	FE00	164	CPI 00H ; IF CE KEY, ZERO FLAG WILL SET.
0107	CAD000	C 165	JZ KDSPT3 ; START OVER IF CE KEY PUSHED.
010A	FE0F	166	CPI 0FH ; IF ENTER KEY, ZERO WILL SET.
010C	CAAA01	C 167	JZ NOBLN3 ; BEFORE RETURN, STOP ALPHA FROM BLINKING.
010F	21F301	C 168	LXI H, NUMTAB ; GET DECIMAL KEY VALUE.
0112	05	169	ADD L ; ADD KEY VALUE TO L REGISTER.
0113	4F	170	MOV L, A ; HL NOW POINTS TO ABSOLUTE KEY VALUE.
0114	7E	171	MOV A, M ; GET ACTUAL DECIMAL KEY VALUE.
0115	FEFF	172	CPI 0FFH ; IF ERROR KEY, ZERO FLAG WILL SET.
0117	CAF100	C 173	JZ KLOOK3 ; IF ERROR, GO AND READ KEYBOARD AGAIN.
011A	04	174	INR B ; INCREMENT DIGIT COUNTER.
011B	4F	175	MOV C, A ; STORE A IN C REG.
011C	3E01	176	MVI A, 01H ; PREPARE TO TEST DIGIT COUNTER.
011E	08	177	CMP B ; IF DIGIT COUNTER=01, ZERO FLAG WILL SET.
011F	CA3101	C 178	JZ STOR13 ; IF 01, GO TO KSTOR1 ROUTINE.
0122	3E02	179	MVI A, 02H ; TEST TO SEE IF = TO 02.
0124	B8	180	CMP B
0125	CA4401	C 181	JZ STOR23 ; IF ZERO SETS, GO TO KSTOR2 ROUTINE.
0128	3E03	182	MVI A, 03H ; IS IT 3RD KEY?
012A	08	183	CMP B
012B	CA6F01	C 184	JZ STOR33 ; IF ZERO SETS, GO TO KSTOR3 ROUTINE.
012E	C3F100	C 185	JMP KLOOK3 ; MUST BE 4TH KEY ENTRY. NOT ALLOWED!!!!
0131	79	186	STOR13: MOV A, C ; RETRIVE A (DECIMAL KEY VALUE).
0132	320000	D 187	STA KSTOR1
0135	CDE201	C 188	CALL LOOKUP ; NOW GET VALUE FOR DISPLAY, RET IN C REG.
0138	210000	E 189	LXI H, COMM1
013B	3604	190	MVI M, 84H ; DISPLAY POS. 1 NON-A1.
013D	210000	E 191	LXI H, DATA1
0140	71	192	MOV M, C ; SEND CHARACTER TO DISPLAY POS. 1.
0141	C3F100	C 193	JMP KLOOK3 ; GO LOOK FOR NEXT KEY ENTRY.
0144	79	194	STOR23: MOV A, C ; RETRIVE A.
0145	320100	D 195	STA KSTOR2 ; STORE DECIMAL KEY VALUE IN KSTOR2.
0148	CDE201	C 196	CALL LOOKUP ; GET DISPLAY CHARACTER.
014B	210000	E 197	LXI H, COMM1
014E	3644	198	MVI M, 64H ; READ DISPLAY RAM POS. 1
0150	210000	E 199	LXI H, DATA1
0153	7E	200	MOV A, M ; READ VALUE INTO A REG.
0154	210000	E 201	LXI H, COMM1
0157	3693	202	MVI M, 93H ; PREPARE TO WRITE TO POS. 2 (A1).
0159	210000	E 203	LXI H, DATA1
015C	77	204	MOV M, A ; WRITE POS. 1 TO POS. 2.
015D	71	205	MOV M, C ; NOW WRITE 2ND KEY ENTRY TO POS. 1
015E	3A0100	D 206	LDA KSTOR2 ; REVERSE KSTOR1 & KSTOR2.
0161	57	207	MOV D, A
0162	3A0000	D 208	LDA KSTOR1
0165	320100	D 209	STA KSTOR2
0168	7A	210	MOV A, D
0169	320000	D 211	STA KSTOR1
016C	C3F100	C 212	JMP KLOOK3 ; GO LOOK FOR 3RD KEY ENTRY.
016F	79	213	STOR33: MOV A, C ; RETRIVE A.
0170	320200	D 214	STA KSTOR3 ; STORE DECIMAL KEY VALUE IN KSTOR3.
0173	CDE201	C 215	CALL LOOKUP ; GET DISPLAY CHARACTER.

LOC	OBJ	LINE	SOURCE STATEMENT
0176	210000	E 216	LXI M, COMM1
0179	3663	217	MVI M, 63H ; READ DISPLAY RAM POS. 2 (ADDRESS=03H).
017E	210000	E 218	LXI M, DATA1
017F	7E	219	MOV A, M ; READ POS.2 VALUE INTO A REG.
0182	210000	E 220	LXI M, COMM1
0184	3664	221	MVI M, 64H ; READ DISPLAY RAM POS. 1 (ADDRESS=04H).
0187	210000	E 222	LXI M, DATA1
0188	56	223	MOV D, M ; READ POS.1 VALUE INTO D REG.
018B	210000	E 224	LXI M, COMM1
018C	3692	225	MVI M, 92H ; PREPARE TO WRITE TO POS.3 (A1).
018D	210000	E 226	LXI M, DATA1
0190	77	227	MOV M, A ; WRITE 1ST KEY ENTRY TO POS.3 (ADDRESS=02H).
0191	72	228	MOV M, D ; WRITE 2ND KEY ENTRY TO POS.2 (ADDRESS=03H).
0192	71	229	MOV M, C ; WRITE 3RD KEY ENTRY TO POS.3 (ADDRESS=04H).
0193	3A0200	D 230	LDA KSTOR3
0196	57	231	MOV D, A
0197	3A0100	D 232	LDA KSTOR2
019A	320200	D 233	STA KSTOR3
019D	3A0000	D 234	LDA KSTOR1
01A0	320100	D 235	STA KSTOR2
01A3	7A	236	MOV A, D
01A4	320000	D 237	STA KSTOR1 ; SCHUFFLE KSTORs AROUND ARIT.
01A7	C3F100	C 238	JMP KLOOK3 ; GO LOOK FOR E OR CE KEY ENTRY.
01AA	3EC3	239 NOBLN3:	MVI A, 0C3H ; CODE FOR NO-BLINE
01AC	320000	E 240	STA ALPHA ; CONTROL WORD.
01AF	CD0000	E 241	CALL DELAY
01B2	3A0000	D 242	LDA KSTOR1 ; WANT TO GET A SINGLE BYTE # FROM KSTOR1,2,3.
01B5	2E00	243	MVI L, 00H ; ZERO L AS L WILL CONTAIN RUNNING SUM.
01B7	0E01	244	MVI C, 01H ; C BECOMES 1'S FACTOR.
01B9	09	245 ONES3:	DAD B ; ADD RC TO HL.
01BA	3D	246	DCR A ; DECREMENT ONES'S COUNTER.
01BB	FE00	247	CPI 00H ; IF COUNTER = 0, ZERO WILL SET.
01BD	C2B901	C 248	JNZ ONES3 ; IF NO ZERO, KEEP ADDING ONES.
01C0	3A0100	D 249	LDA KSTOR2 ; NOW DO IT FOR 10'S PLACE.
01C3	0E0A	250	MVI C, 0AH ; C IS NOW 10'S FACTOR.
01C5	09	251 TENS3:	DAD B ; ADD RC TO HL.
01C6	3D	252	DCR A ; DECREMENT 10'S COUNTER.
01C7	FE00	253	CPI 00H
01C9	C2C501	C 254	JNZ TENS3
01CC	3A0200	D 255	LDA KSTOR3 ; NOW DO IT FOR HUNDREDS PLACE.
01CF	0E64	256	MVI C, 64H ; C IS NOW 100'S FACTOR.
01D1	09	257 HUNS3:	DAD B
01D2	3D	258	DCR A
01D3	FE00	259	CPI 00H
01D5	C2D101	C 260	JNZ HUNS3
01D8	7D	261	MOV A, L ; PLACE RUNNING SUM INTO A REG.
01D9	320400	D 262	STA SUM3 ; STORE RUNNING SUM AT SUM2.
01DC	210000	E 263	LXI M, COMM1
01DF	36DF	264	MVI M, 0DFH
01E1	C9	265	RET ; RETURN TO CALLING PROGRAM.
		266 ;	
		267 ;	SUBROUTINES FOLLOW:
		268 ;	
01E2	21E901	C 269	LOOKUP: LXI M, DTABLE ; STARTING ADDRESS OF DISPLAY TABLE.

146

LOC	OBJ	LINE	SOURCE STATEMENT
01E5	05	270	ADD L ; ADD KEY VALUE (0-9) TO LOB OF DTABLE.
01E6	6F	271	MOV L, A ; MOVE RESULT TO L REGISTER.
01E7	4E	272	MOV C, M ; PUT DISPLAY RESULT INTO C REGISTER.
01E8	C9	273	RET
		274 ;	
		275 ;	LOOK-UP TABLES USED IN FINDING KEY VALUES FOLLOW:
		276 ;	
01E9	0C	277	DTABLE: DB 0CH ; 0 FOR CONTROL DISPLAY.
01EA	9F	278	DB 9FH ; 1
01EB	4A	279	DB 4AH ; 2
01EC	0B	280	DB 0BH ; 3
01ED	99	281	DB 99H ; 4
01EE	29	282	DB 29H ; 5
01EF	20	283	DB 20H ; 6
01F0	0F	284	DB 0FH ; 7
01F1	00	285	DB 00H ; 8
01F2	09	286	DB 09H ; 9
		287 ;	
01F3	FF	288	NUMTAB: DB 0FFH ; FF IS ERROR CODE FOR INVALID KEY ENTRY.
01F4	03	289	DB 03H
01F5	06	290	DB 06H
01F6	09	291	DB 09H
01F7	00	292	DB 00H
01F8	02	293	DB 02H
01F9	05	294	DB 05H
01FA	00	295	DB 00H
01FB	FF	296	DB 0FFH ; ERROR
01FC	01	297	DB 01H
01FD	04	298	DB 04H
01FE	07	299	DB 07H
01FF	FF	300	DB 0FFH ; ERROR
0200	FF	301	DB 0FFH ; ERROR
0201	FF	302	DB 0FFH ; ERROR
		303 ;	
		304 ;	*****
		305 ;	THIS SUBROUTINE EMPTIES THE KEYBOARD FIFO.
		306 ;	*****
		307 ;	
0202	210000	E 308	FIFOCR: LXI H, COMM2 ; SET FOR FIFO READ AUTO-INCREMENT.
0205	3441	309	MVI M, 41H
0207	210000	E 310	EMPTY: LXI H, COMM2
020A	7E	311	MOV A, M ; GET FIFO STATUS WORD.
020B	E403	312	ANI 03H ; MASK OUT 5 HI-ORDER BITS.
020D	FE00	313	CPI 00H
020F	C0	314	RZ ; IF ZERO SETS, THEN FIFO EMPTY SO RETURN.
0210	210000	E 315	LXI H, DATA2
0213	7E	316	MOV A, M
0214	C30702	C 317	JMP EMPTY
		318 ;	
		319	DSEG
		320 ;	
0000		321	KSTOR1: DS 1
0001		322	KSTOR2: DS 1
0002		323	KSTOR3: DS 1

147

LOC	OBJ	LINE	SOURCE STATEMENT
0003		324	SUM2: DS 1
0004		325	SUM3: DS 1
		326	;
		327	CSEG
		328	;
		329	*****
		330	;
		331	THIS SUBROUTINE WRITES A MESSAGE TO THE LINE PRINTER.
		332	PWRITE (PRINTER WRITE) REQUIRES THE FOLLOWING INFORMATION.
		333	C REG: CONTAINS THE LENGTH OF THE TEXT STRING TO
		334	BE WRITTEN TO THE PRINTER.
		335	HL REG: CONTAINS THE STARTING ADDRESS OF THE TEXT
		336	STRING.
		337	THIS ROUTINE CALLS PSEND.
		338	;
		339	*****
		340	;
0217	0D	341	PWRITE: DCR C ; DECREMENT CHARACTER COUNTER.
0218	3E00	342	MVI A, 00H ; PREPARE TO TEST COUNTER.
021A	89	343	CMP C ; IF C=0 ZERO FLAG WILL SET.
021B	C8	344	RZ ; IF NO MORE TO SEND, THEN RETURN.
021C	7E	345	MOV A, H ; PUT CHARACTER POINTED TO BY HL INTO A.
021D	CD0000	E 346	CALL PSEND ; SEND IT TO PRINTER.
0220	23	347	INX H ; INCREMENT TEXT POINTER.
0221	C31702	C 348	JMP PWRITE
		349	;
		350	*****
		351	;
		352	THIS SUBROUTINE TAKES THE LO-ORDER NIBBLE STORED AT DEB
		353	(FOR DATA ENTRY BYTE) AND DISPLAYS IT IN THE FORMAT OF
		354	" PROCESS (LO-DEB) = ? " ON THE ALPHA DISPLAY.
		355	;
		356	*****
		357	;
0224	3EE1	358	DEVNUM: MVI A, 0E1H ; LOAD BLANK CODE.
0226	320000	E 359	STA ALPHA
0229	CD0000	E 360	CALL DELAY
022C	215202	C 361	LXI H, DEVTAB ; LOAD STARTING ADDRESS OF MESSAGE.
022F	0609	362	MVI B, 09H ; LOAD TEXT COUNTER.
0231	CD0000	E 363	CALL DWRITE ; WRITE MESSAGE ON ALPHA DISPLAY.
0234	3A0000	E 364	LDA DEB ; LOAD DATA ENTRY BYTE (DEB).
0237	E60F	365	ANI 0FH ; MASK TO GET PROCESS ENTRY COUNTER.
0239	C630	366	ADI 30H ; CONVERT TO ASCII.
023B	320000	E 367	STA ALPHA ; SEND TO DISPLAY.
023E	CD0000	E 368	CALL DELAY
0241	9E3D	369	MVI A, 3DH ; SEND "=" TO ALPHA DISPLAY.
0243	320000	E 370	STA ALPHA
0246	CD0000	E 371	CALL DELAY
0249	3E3F	372	MVI A, 3FH ; SEND "?" TO DISPLAY.
024B	320000	E 373	STA ALPHA
024E	CD0000	E 374	CALL DELAY
0251	C9	375	RET ; RETURN TO CALLING PROGRAM.
		376	;
		377	MESSAGE TABLE FOLLOWS:

LOC	OBJ	LINE	SOURCE STATEMENT
		378 ;	
0252	10	379	DEV TAB: DB 10H ; P
0253	12	380	DB 12H ; R
0254	0F	381	DB 0FH ; O
0255	03	382	DB 03H ; C
0256	05	383	DB 05H ; E
0257	13	384	DB 13H ; S
0258	13	385	DB 13H ; S
0259	20	386	DB 20H ; SPACE
		387 ;	
		388 ;	*****
		389 ;	
		390 ;	THIS SUBROUTINE TAKES THE VALUE OF THE LO-ORDER NIBBLE OF DEB
		391 ;	AND THE CURRENT CONTENTS OF KSTOR3,2,1 AND PRINTS THEM ON THE
		392 ;	LINE PRINTER. THE MESSAGE FORMAT IS:
		393 ;	" PROCESS [LO-DEB] = [KSTOR3,2,1] MINUTES "
		394 ;	CALLS PWRITE AND PSEND.
		395 ;	
		396 ;	*****
		397 ;	
025A	CDBE02	C 398	PNTNUM: CALL STAR ; PRINT A FULL ROW OF " ***** "
025D	21CF02	C 399	LXI H, PENTRY ; STARTING ADDRESS OF MESSAGE.
0260	0E11	400	MVI C, 11H ; LOAD MESSAGE LENGTH.
0262	CD1702	C 401	CALL PWRITE ; WRITE MESSAGE TO PRINTER.
0265	3A0000	E 402	LDA DEB
0268	E60F	403	ANI 0FH ; MASK TO GET PROCESS ENTRY COUNTER.
026A	C630	404	ADI 30H ; CONVERT IT TO ASCII.
026C	CD0000	E 405	CALL PSEND ; SEND IT TO PRINTER.
026F	3E20	406	MVI A, 20H ; SEND BLANK.
0271	CD0000	E 407	CALL PSEND
0274	3E3D	408	MVI A, 3DH ; SEND "=" TO PRINTER.
0276	CD0000	E 409	CALL PSEND
0279	3E20	410	MVI A, 20H ; SEND ANOTHER BLANK.
027B	CD0000	E 411	CALL PSEND
027E	3A0200	D 412	LDA KSTOR3 ; NOW SENDS ACTUAL VALUES TO PRINTER.
0281	FE00	413	CPI 00H ; USE LEADING ZERO BLANKING.
0283	C28E02	C 414	JNZ NOBLK1
0286	3E20	415	MVI A, 20H
0288	CD0000	E 416	CALL PSEND
028B	C39302	C 417	JMP DIC2
028E	C630	418	NOBLK1: ADI 30H ; HUNDREDS NOT EQUAL TO 0. CONVERT ASCII.
0290	CD0000	E 419	CALL PSEND
0293	3A0100	D 420	DIC2: LDA KSTOR2 ; GET 10'S PLACE.
0296	C630	421	ADI 30H ; CONVERT TO ASCII.
0298	CD0000	E 422	CALL PSEND
029B	3A0000	D 423	LDA KSTOR1 ; GET 1'S PLACE.
029E	C630	424	ADI 30H ; CONVERT
02A0	CD0000	E 425	CALL PSEND
02A3	21DF02	C 426	LXI H, PMIN ; END BY SENDING "MINUTES" TO PRINTER.
02A6	0E09	427	MVI C, 09H ; LOAD CHARACTER COUNT.
02AB	CD1702	C 428	CALL PWRITE
02AB	3E17	429	MVI A, 17H ; NOW PRINT PRINTER BUFFER CONTENTS.
02AD	CD0000	E 430	CALL PSEND
02B0	CDBE02	C 431	CALL STAR ; UNDERLINE MESSAGE WITH " ***** "

149

LOC	OBJ	LIME	SOURCE STATEMENT
02B3	3E17	432	MVI A, 17H ; NOW SPACE 2 BLANK LINES.
02B5	CD0000 E	433	CALL PSEND
02B0	3E17	434	MVI A, 17H
02BA	CD0000 E	435	CALL PSEND
02BD	C9	436	RET ; RETURN TO CALLING PROGRAM.
		437 ;	
		438 ;	ROUTINE TO GENERATE ONE FULL ROW OF " ***** "
		439 ;	
02BE	0E27	440	STAR: MVI C, 27H ; LOAD COLUMN COUNT WITH 39 COUNTS.
02C0	3E2A	441	UNMAS: MVI A, 2AH ; LOAD A WITH ASCII FOR " * ".
02C2	CD0000 E	442	CALL PSEND ; SEND IT TO PRINTER.
02C5	0D	443	DCR C
02C6	C2C002 C	444	JNZ UNMAS
02C9	3E17	445	MVI A, 17H ; START TO PRINT COMMAND.
02CB	CD0000 E	446	CALL PSEND
02CE	C9	447	RET
		448 ;	
		449 ;	MESSAGE TABLES FOLLOW:
		450 ;	
02CF	20	451	PENTRY: DB 20H ; SPACE
02D0	20	452	DB 20H ; SPACE
02D1	20	453	DB 20H ; SPACE
02D2	20	454	DB 20H ; SPACE
02D3	20	455	DB 20H ; SPACE
02D4	20	456	DB 20H ; SPACE
02D5	20	457	DB 20H ; BLANK
02D6	20	458	DB 20H ; BLANK
02D7	50	459	DB 50H ; P
02D8	72	460	DB 72H ; R
02D9	4F	461	DB 4FH ; O
02DA	43	462	DB 43H ; C
02DB	45	463	DB 45H ; E
02DC	73	464	DB 73H ; S
02DD	73	465	DB 73H ; S
02DE	20	466	DB 20H ; BLANK
02DF	20	467	PMIN: DB 20H ; BLANK
02E0	4D	468	DB 4DH ; M
02E1	49	469	DB 49H ; I
02E2	4E	470	DB 4EH ; N
02E3	75	471	DB 75H ; U
02E4	74	472	DB 74H ; T
02E5	45	473	DB 45H ; E
02E6	73	474	DB 73H ; S
		475 ;	
		476	END

PUBLIC SYMBOLS

BEVNUM C 0224	FIFO CR C 0202	KDSPY2 C 0000	KDSPY3 C 00D0	KSTOR1 D 0000	KSTOR2 D 0001	KSTOR3 D 0002
LOOKUP C 01E2	PMTNUM C 025A	PWRITE C 0217	SUM2 D 0003	SUM3 D 0004		

EXTERNAL SYMBOLS

ALPHA E 0000	CLEARA E 0000	CLEARB E 0000	COMM1 E 0000	COMM2 E 0000	DATA1 E 0000	DATA2 E 0000
DEB E 0000	DELAY E 0000	DWRITE E 0000	POS1 E 0000	POS2 E 0000	PSEND E 0000	RPOS1 E 0000

150

USER SYMBOLS

ALPHA E 0000	CLEARA E 0000	CLEARB E 0000	COMM1 E 0000	COMM2 E 0000	DATA1 E 0000	DATA2 E 0000
DEB E 0000	DELAY E 0000	DEVNUM C 0224	DEVTAB C 0252	DIG2 C 0293	DTABLE C 01E9	DWRITE E 0000
EMPTY C 0307	FIFOCH C 0202	HUMS3 C 01D1	KDSPY2 C 0000	KDSPY3 C 00D0	KLOOK C 0024	KLOOK3 C 00F1
KSTOR1 D 0000	KSTOR2 D 0001	KSTOR3 D 0002	LOOKUP C 01E2	MOBLK1 C 028E	MOBLN3 C 01AA	MOBLNK C 00A2
KUMTAB C 01F3	ONES2 C 00B8	ONES3 C 01B9	PENTRY C 02CF	PMIN C 02DF	PMTNUM C 025A	POS1 E 0000
POS2 E 0000	PSEND E 0000	PWRITE C 0217	RPOS1 E 0000	STAR C 02BE	STOR1 C 005E	STOR13 C 0131
STOR2 C 0073	STOR23 C 0144	STOR33 C 016F	SUM2 D 0003	SUM3 D 0004	TENS2 C 00C4	TENS3 C 01C5
UNMAS C 02C0						

ASSEMBLY COMPLETE, NO ERRORS

LOC	OBJ	LINE	SOURCE STATEMENT
		1	*****
		2	;
		3	THIS PROGRAM PERFORMS ALL THE TIMING FUNCTIONS
		4	USED IN THE PROCESSING OF FILM WITH THE IR DEN-
		5	SITOMETER. SCAN TIMING IS ALSO COMPUTED HERE, BUT
		6	THE ACTUAL SCANS ARE TO BE MADE BY ANOTHER PROGRAM
		7	CALL SCAN.
		8	THIS PROGRAM IS EXITED WHEN ALL PROCESSING IS DONE.
		9	;
		10	*****
		11	;
		12	NAME CODEV
		13	;
		14	PUBLIC CODEV,INVERT,MINCNT,MCNT
		15	;
		16	EXTRN IRCNT,IRBRAM,PLBYTE,LDELAY,DELAY,ALPHA,DWRITE,KDSPY2
		17	EXTRN RPB2,RPA2,DFIND,SCAN
		18	;
		19	CSEG
		20	;
0000	F3	21	CODEV: DI ; DISABLE INTERRUPT SYSTEM.
0001	210000	E 22	LXI H, ALPHA
0004	36E1	23	MVI M, 0E1H ; BLANKL ALPHA DISPLAY.
0006	CD0000	E 24	CALL DELAY ; GIVE ALPHA TIME TO CLEAR.
0009	219001	C 25	LXI H, PROCES ; START OF "BEGIN?" MESSAGE.
000C	060A	26	MVI B, 0AH ; LOAD TEXT COUNTER.
000E	CD0000	E 27	CALL DWRITE ; WRITE "BEGIN?" ON ALPHA DISPLAY.
0011	CD0000	E 28	CALL KDSPY2 ; "ENTER" KEY WILL START PROCESSING.
0014	06FF	29	MVI B, 0FFH ; INITIALIZE SECOND AND SCAN TIME COUNTERS.
0016	48	30	MOV C, 8
0017	AF	31	XRA A ; ZERO A.
0018	320600	D 32	STA SCNNUM ; SET SCAN COUNTER TO ZERO.
001B	320000	D 33	STA MINCNT ; SET MINUTE COUNTER TO ZERO.
001E	3E01	34	MVI A, 01H ; SET UP TASK COUNTERS.
0020	320100	D 35	STA MCNT ; M IS 1/R COMBO POINTER.
0023	320200	D 36	STA MCNT ; M IS PROCESS 0 POINTER.
0026	3E00	37	TIMING: MVI A, 00H ; SET FINTAG (FINISHED TAG) OFF.
0028	320700	D 38	STA FINTAG
002B	210000	E 39	LXI H, IRCNT ; HL CONTAINS START ADDRESS OF IRCNTS.
002E	3A0200	D 40	LDA MCNT
0031	85	41	ADD L
0032	4F	42	MOV L, A ; HL = IRCNT + MCNT.
0033	7E	43	MOV A, M ; A = IRCNTn.
0034	320300	D 44	STA TIRCNT ; TIRCNT NOW HAS COPY OF IRCNTn.
0037	210000	E 45	TLOAD: LXI H, IRBRAM ; HL CONTAINS START ADDRESS OF IRBRAM.
003A	3A0100	D 46	LDA MCNT
003D	85	47	ADD L
003E	4F	48	MOV L, A ; HL POINTS TO RATEm.
003F	E5	49	PUSH H
0040	CD1D01	C 50	CALL INVERT ; CONVERT Y SCANS/MIN TO Z SECONDS/SCAN.
0043	E1	51	POP H
0044	320400	D 52	STA TRATE ; COPY OF RATEm NOW IN TRATE.
0047	23	53	INX H ; HL NOW POINTS TO INTERVALm.

LOC	OBJ	LINE	SOURCE STATEMENT
0040	7E	54	MOV A, M
0049	320500	D 55	STA TINTR ; COPY OF INTERVALn NOW IN TINTR.
004C	3A0700	D 56	RESET: LDA FINTAG ; CHEC STATUS OF PROCESS FINISHED TAG.
004F	FEFF	57	CPI 0FFH ; WHEN FINTAG=FFH, THEN LAST SCAN.
0051	CA7700	C 58	JZ RELOAD
0054	3ED0	59	MVI A, 0D8H ; RESET AND UN-MASK RST 7.5, & SET SOD=1
0056	30	60	SIM
0057	20	61	PEND: RIM
0058	E640	62	ANI 40H ; MASK TO GET RST 7.5 PENDING BIT.
005A	FE00	63	CPI 00H
005C	CA5700	C 64	JZ PEND ; WHEN ZERO DOESNT SET, THEN INTERRUPT.
005F	210000	E 65	LXI H, ALPHA ; BLANK ALPHA DISPLAY WHEN PROCESS BEGINS.
0062	36E1	66	MVI M, 0E1H
0064	0C	67	INR C ; INTERRUPT RECEIVED.
0065	04	68	INR B ; SECCNT AND SCNSEC INCREMENTED BY 1.
0066	70	69	MOV A, B
0067	D670	70	SUI 70H ; IF B=60, ZEWRO FLAG WILL SET.
0069	C20400	C 71	JNZ NOT60
006C	3A0300	D 72	LDA TIRCNT ; COME HERE IF SECOND COUNT = 60.
006F	FE00	73	CPI 00H ; IF ZERO SETS, THEN TIRCNT = 0.
0071	CA7B00	C 74	JZ INC
0074	3A0500	D 75	LDA TINTR ; IF NOT ZERO, THEN DECREMENT TINTR.
0077	3D	76	DCR A
0070	320500	D 77	STA TINTR ; UP-DATE TINTR.
007B	3A0000	D 78	INC: LDA MINCNT ; INCREMENT MINUTE COUNTER.
007E	3C	79	INR A
007F	320000	D 80	STA MINCNT
0082	0600	81	MVI B, 00H ; RESET SECCNT TO 0.
0084	210000	E 82	NOT60: LXI H, PLBYTE ; LOAD HL WITH START ADDRESS OG PLBYTE.
0087	3A0200	D 83	LDA NCNT
008A	05	84	ADD L
008B	6F	85	MOV L, A ; HL NOW POINTS TO PLBn.
008C	3A0000	D 86	LDA MINCNT ; TEST FOR END OF PROCESS.
008F	BE	87	CMP M
0090	C2A000	C 88	JNZ NOTFIN ; IF ZERO SETS, PLBn = MINCNT, PROCESS DONE.
0093	3EFF	89	MVI A, 0FFH ; PROCESS DONE, SET FINTAG ON, DO 1 MORE SCAN.
0095	320700	D 90	STA FINTAG
0098	3A0000	D 91	LDA MINCNT ; TEST FOR ZERO PROCESS LENGTHS.
009B	FE00	92	CPI 00H
009D	CAEA00	C 93	JZ SKIP
00A0	3A0200	D 94	NOTFIN: LDA NCNT ; TURN ON TANKn AND PUMP.
00A3	FE01	95	CPI 01H
00A5	C2AF00	C 96	JNZ TANK2
00A8	3E43	97	MVI A, 43H ; TURN ON TANK 1 AND PUMP.
00AA	D300	E 98	OUT LOW RPB2 ; EXPANSION ROM PORT B.
00AC	C3BF00	C 99	JMP NOWON
00AF	FE02	100	TANK2: CPI 02H
00B1	C2BB00	C 101	JNZ TANK3
00B4	3E4C	102	MVI A, 4CH ; TURN ON TANK 2 AND PUMP.
00B6	D300	E 103	OUT LOW RPB2
00B8	C3BF00	C 104	JMP NOWON
00BB	3E70	105	TANK3: MVI A, 70H ; TURN ON TANK 3 VALVES AND PUMP.
00BD	D300	E 106	OUT LOW RPB2
00BF	3A0300	D 107	NOWON: LDA TIRCNT ; TEST TO SEE IF I/R COMBOS REMAINING.

LOC	OBJ	LINE	SOURCE STATEMENT
00C2	FE00	100	CPI 00H ; IF TIRCNT = 0, ZERO WILL SET.
00C4	CA4C00	C 109	JZ RESET ; GO WAIT FOR NEXT INTERRUPT.
00C7	3A0400	D 110	LDA TRATE ; DO WE MAKE A SCAM?
00CA	FEFF	111	CPI 0FFH ; TEST FOR ZERO RATE FLAG.
00CC	CAD300	C 112	JZ ZRATE ; IF ZERO RATE, BYPASS SCAM TEST.
00CF	B9	113	CMP C
00D0	CC7A01	C 114	CZ TSCAN ; YES. GO DO SCAM.
00D3	3A0500	D 115	ZRATE: LDA TINTR ; TEST FOR END OF INTERVAL.
00D6	FE00	116	CPI 00H
00D8	C24C00	C 117	JNZ RESET ; IF NOT END, GO AWAIT NEXT INTERRUPT.
00DB	3A0100	D 118	LDA MCNT
00DE	C602	119	ADI 02H ; MCNT = MCNT + 2.
00E0	320100	D 120	STA MCNT
00E3	3A0300	D 121	LDA TIRCNT ; ONE I/R SET COMPLETE.
00E6	3D	122	DCR A ; DECREMENT TIRCNT.
00E7	320300	D 123	STA TIRCNT
00EA	0E00	124	SKIP: MVI C, 00H ; RESET SCNSEC (C REG) TO ZERO.
00EC	3A0700	D 125	LDA FINTAG ; CHECK TO SEE IF PROCESS DONE.
00EF	FEFF	126	CPI 0FFH ; RECALL FFH = PROCESS DONE.
00F1	CAF700	C 127	JZ RELOAD ; GO AND TEST FOR REMAINING PROCESSES.
00F4	C33700	C 128	JMP TLOAD ; GET NEW TEMPORARY I/R VALUES.
00F7	0B00	E 129	RELOAD: IN LOW RPB2
00F9	E63F	130	ANI 3FH ; TURN OFF ONLY THE PUMP.
00FB	0300	E 131	OUT LOW RPB2
00FD	D5	132	PUSH D ; SAVE CONTENTS OF DE REGISTER.
00FE	CD6D01	C 133	CALL MDELAY
0101	CD6D01	C 134	CALL MDELAY ; GIVE PUMP TIME TO STOP BEFORE CLOSING VALVES.
0104	D1	135	POP D ; RETRIVE CONTENTS OF DE REGISTER.
0105	AF	136	XRA A ; ZERO A.
0106	D300	E 137	OUT LOW RPB2 ; NOW TURN OFF VALVES.
0108	320800	D 138	STA MINCNT ; RESET MINCNT TO ZERO.
010B	06FF	139	MVI B, 0FFH ; RESET SECOND AMS SCAM COUNTERS.
010D	40	140	MOV C, B
010E	3A0200	D 141	LDA MCNT
0111	3C	142	INR A ; INCREMENT PROCESS COUNTER.
0112	320200	D 143	STA MCNT
0115	FE04	144	CPI 04H ; HAVE THREE PROCESSES BEEN MADE?
0117	CA0000	E 145	JZ DFIND ; IF YES, GO TO DATA OUT ROUTINE.
011A	C32600	C 146	JMP TIMING ; OTHERWISE, GO GET NEXT SET OF PROCESS INFO.
011D	7E	147	INVERT: MOV A, M ; TEST TO SEE IF A ZERO RATE IS PASSED.
011E	FE00	148	CPI 00H
0120	C22601	C 149	JNZ OK ; IF > 0, PASS BY.
0123	3EFF	150	MVI A, 0FFH ; OTHERWISE TAG THE RATE.
0125	C9	151	RET
0126	212F01	C 152	OK: LXI H, TABLE ; LOAD HL WITH START ADDRESS OF TABLE.
0129	5F	153	MOV E, A
012A	1600	154	MVI D, 00H ; DE = RATE SCANS/MINUTE.
012C	19	155	DAD D ; HL = HL + DE.
012D	7E	156	MOV A, M ; HL POINTS TO RATE SEC/SCAN.
012E	07	157	RLC ; MULTIPLY RESULT IN A BY 2.
012F	C9	158	TABLE: RET
0130	3C	159	DB 3CH ; FOR RATE ENTRY OF 1 SCAN/MINUTE.
0131	1E	160	DB 1EH ; FOR ENTRY OF 2 SCAN/MINUTE.
0132	14	161	DB 14H ; FOR 3 SCANS/MINUTE.

LOC	OBJ	LINE	SOURCE STATEMENT
0133	0F	162	DB 0FH ; 4 SCANS/MINUTE.
0134	0C	163	DB 0CH ; 5 SCANS/MINUTE.
0135	0A	164	DB 0AH ; 6 S/M.
0136	09	165	DB 09H ; 7
0137	08	166	DB 08H ; 8
0138	07	167	DB 07H ; 9
0139	06	168	DB 06H ; 10
013A	05	169	DB 05H ; 11
013B	05	170	DB 05H ; 12
013C	05	171	DB 05H ; 13
013D	04	172	DB 04H ; 14
013E	04	173	DB 04H ; 15
013F	04	174	DB 04H ; 16
0140	04	175	DB 04H ; 17
0141	03	176	DB 03H ; 18
0142	03	177	DB 03H ; 19
0143	03	178	DB 03H ; 20
0144	03	179	DB 03H ; 21
0145	03	180	DB 03H ; 22
0146	03	181	DB 03H ; 23
0147	02	182	DB 02H ; 24
0148	02	183	DB 02H ; 25
0149	02	184	DB 02H ; 26
014A	02	185	DB 02H ; 27
014B	02	186	DB 02H ; 28
014C	02	187	DB 02H ; 29
014D	02	188	DB 02H ; 30
014E	02	189	DB 02H ; 31
014F	02	190	DB 02H ; 32
0150	02	191	DB 02H ; 33
0151	02	192	DB 02H ; 34
0152	02	193	DB 02H ; 35
0153	02	194	DB 02H ; 36
0154	02	195	DB 02H ; 37
0155	02	196	DB 02H ; 38
0156	02	197	DB 02H ; 39
0157	02	198	DB 02H ; 40
0158	01	199	DB 01H ; 41
0159	01	200	DB 01H ; 42
015A	01	201	DB 01H ; 43
015B	01	202	DB 01H ; 44
015C	01	203	DB 01H ; 45
015D	01	204	DB 01H ; 46
015E	01	205	DB 01H ; 47
015F	01	206	DB 01H ; 48
0160	01	207	DB 01H ; 49
0161	01	208	DB 01H ; 50
0162	01	209	DB 01H ; 51
0163	01	210	DB 01H ; 52
0164	01	211	DB 01H ; 53
0165	01	212	DB 01H ; 54
0166	01	213	DB 01H ; 55
0167	01	214	DB 01H ; 56
0168	01	215	DB 01H ; 57

LOC	OBJ	LINE	SOURCE STATEMENT
0169	01	216	DB 01H ; 57
016A	01	217	DB 01H ; 58
016B	01	218	DB 01H ; 59
016C	01	219	DB 01H ; 60
016D	16FF	220	MDELAY: MVI D, 0FFH ; MEDIUM LENGTH DELAY.
016F	1EFF	221	LOOP2: MVI E, 0FFH
0171	1D	222	LOOP1: DCR E
0172	C27101	C 223	JNZ LOOP1
0175	15	224	DCR D
0176	C26F01	C 225	JNZ LOOP2
0179	C9	226	RET
017A	9A0600	D 227	TSCAN: LDA SCNNUM ; SIMPLY INCREMENT SCAN COUNTER
017D	3C	228	INR A ; UNTIL SCAN ROUTINE IS WORKED OUT.
017E	320600	D 229	STA SCNNUM
0181	0E00	230	MVI C, 00H ; RESET SCNSEC TO ZERO.
0183	3E40	231	MVI A, 40H ; SET SOD LINE ON.
0185	30	232	SIM
0186	CD0000	E 233	CALL DELAY ; WAIT A BIT.
0189	3EC0	234	MVI A, 0C0H ; SET SOD LINE OFF.
018B	30	235	SIM
018C	CD0000	E 236	CALL SCAN
018F	C9	237	RET
		238 ;	
0190	20	239	PROCES: DB 20H ; BLANK
0191	20	240	DB 20H ; BLANK
0192	20	241	DB 20H ; BLANK
0193	02	242	DB 02H ; B
0194	05	243	DB 05H ; E
0195	07	244	DB 07H ; G
0196	09	245	DB 09H ; I
0197	0E	246	DB 0EH ; M
0198	3F	247	DB 3FH ; ?
		248 ;	
		249	DSEG
		250 ;	
0000		251	MINCNT: DS 1 ; MINUTE COUNTER.
0001		252	MCNT: DS 1 ; I/R COMBO COUNTER.
0002		253	NCNT: DS 1 ; PROCESS # COUNT.
0003		254	TIRCNT: DS 1 ; TEMPORARY COPY OF I/R COUNT.
0004		255	TRATE: DS 1 ; TEMPORARY COPY OF CURRENT RATE VALUE.
0005		256	TINTR: DS 1 ; TEMPORARY COPY OF CURRENT INTERVAL VALUE.
0006		257	SCNNUM: DS 1 ; RUNNING SUM OF SCANS MADE.
0007		258	FINTAG: DS 1 ; USED IN MAKING LAST SCAN OF ANY PROCESS.
		259 ;	
		260	END

PUBLIC SYMBOLS

CODEV C 0000 INVERT C 011D MINCNT D 0000 NCNT D 0002

EXTERNAL SYMBOLS

ALPHA E 0000 DELAY E 0000 DFIND E 0000 DWRITE E 0000 IRBRAM E 0000 IRCNT E 0000 KDSPY2 E 0000
 LDELAY E 0000 PLBYTE E 0000 RPA2 E 0000 RPB2 E 0000 SCAN E 0000

USER SYMBOLS

ALPHA E 0000	DELAY E 0000	DFIND E 0000	DWRITE E 0000	FINTAG D 0007	CODEV C 0000	INC C 007B
INVERT C 011D	IRBRAM E 0000	IRCNT E 0000	KDSPY2 E 0000	LDELAY E 0000	LOOP1 C 0171	LOOP2 C 016F
NCMT D 0001	MDELAY C 016D	MINCNT D 0000	NCMT D 0002	MOT60 C 0084	MOTFIN C 00A0	NOWON C 00BF
OK C 0126	PEND C 0057	PLBYTE E 0000	PROCES C 0190	RELOAD C 00F7	RESET C 004C	RPA2 E 0000
RPE2 E 0000	SCAN E 0000	SCNNUM D 0004	SKIP C 00EA	TABLE C 012F	TANK2 C 00AF	TANK3 C 00BB
TIMING C 0026	TINTR D 0005	TIRCNT D 0003	TLOAD C 0037	TRATE D 0004	TSCAN C 017A	ZRATE C 00D3

ASSEMBLY COMPLETE, NO ERRORS

LOC	OBJ	LINE	SOURCE STATEMENT
		1	*****
		2	;
		3	THIS MODULE PROMPTS THE USER TO ENTER THE PROCESS NUMBER OF A
		4	PARTICULAR SCAN SERIES IN WHICH HE IS INTERESTED IN VIEWING.
		5	;
		6	ONCE GIVEN A PROCESS NUMBER, THE ROUTINE CHECKS TO SEE IF ANY DATA
		7	WAS INDEED TAKEN DURING THAT PROCESS. IF SO, THE USER IS PROMPTED
		8	TO ENTER THE TIME AT WHICH HE THINKS A SCAN WAS MADE. THE PROGRAM
		9	SEARCHES UNTIL IT FINDS THE CLOSEST SCAN TO THE ENTERED TIME. THIS
		10	DENSITY BLOCK IS DISPLAYED ALONG WITH THE ACTUAL TIME IT WAS TAKEN.
		11	;
		12	*****
		13	;
		14	NAME DFIND
		15	;
		16	PUBLIC DFIND,SHOW,MINCON,ITEMP,TIMSUM,OFFSET,MIN
		17	;
		18	EXTRN DWRITE,KDSPY2,SUM2,DEB,IRCNT,STABCT,RAMPT,POS1,POS2,RPOS1
		19	EXTRN CLEARA,CLEARB,ALPHA,NCNT,IRBRAM,INVERT,DSHOW,MINCNT,COMM1
		20	EXTRN BINASC,LOOKUP,LDELAY,DATA1,UPORDN,PLBYTE,THPNUM
		21	;
		22	CSEC
		23	;
0000	3E8D	24	DFIND: MVI A, 0DH ; UN-MASK RST 6.5 FOR USE IN GETTING TO MONITOR
0002	30	25	SIM ; SET INTERRUPT MASK.
0003	FR	26	EI ; ENABLE INTERRUPTS. (RST 6.5 ONLY)
0004	214302	27	LXI H, TWANTM ; PROMPT "PROCESS #=?"
0007	060C	28	MVI B, 0CH ; LOAD TEXT COUNT.
0009	CD0000	29	CALL DWRITE
000C	CD0000	30	CALL KDSPY2 ; WAIT FOR PROCESS # ENTRY.
000F	3A0000	31	LDA SUM2 ; VALID ENTRIES ARE 1 2 3 ONLY.
0012	FE04	32	CPI 04H
0014	D20000	33	JNC DFIND ; INVALID ENTRIES RETURN YOU TO DFIND.
0017	FE00	34	CPI 00H
0019	CA0000	35	JZ DFIND
001C	320000	36	STA DEB ; STORE PROCESS ENTRY IN DEB.
001F	210000	37	LXI H, IRCNT ; LOAD HL WITH START ADDRESS OF IRCNTS.
0022	85	38	ADD L
0023	6F	39	MOV L, A ; HL POINTS TO IRCNT(DEB).
0024	7E	40	MOV A, H
0025	FE00	41	CPI 00H ; IF SELECTED PROCESS HAS 0 IRCNT, "NO DATA"
0027	CC1B02	42	CZ MODATA
		43	;
		44	DETERMINE NOW HOW MANY SCANS MADE BEFORE ARRIVAL AT PROCESS WANTED.
		45	;
002A	0E00	46	MVI C, 00H
002C	3A0000	47	LDA DEB
002F	FE01	48	CPI 01H ; IF DEB=1, START AT BOTTOM OF DRAM.
0031	CA4300	49	JZ ADDED
0034	47	50	MOV B, A
0035	3A0100	51	LDA STABCT+1
0038	4F	52	MOV C, A
0039	78	53	MOV A, B

LOC	OBJ	LINE	SOURCE STATEMENT
003A	FE02	54	CPI 02H ; IF DEB=2, LOAD C WITH PROCESS 1 SCAN CNT.
003C	CA4300	55	JZ ADDED
003F	3A0200	56	LDA STABCT+2
0042	4F	57	MOV C, A
0043	3E01	58	ADDED: MVI A, 01H ; ADD 1 TO OFFSET MDMAX AND MD0 BLOCKS.
0045	81	59	ADD C
0046	320000	60	STA RAMPT ; RAMPT = SUM OF SCANS + 2.
0049	214E02	61	GETTW: LXI H, TWMESS ; PROMPT "TIME WANTED?"
004C	060C	62	MVI B, 0CH
004E	CD0000	63	CALL DWRITE
0051	3E91	64	MVI A, 91H ; SET KDSPY2 FOR MINUTE ENTRY.
0053	320000	65	STA POS2
0056	3E82	66	MVI A, 82H
0058	320000	67	STA POS1
005B	3E62	68	MVI A, 62H
005D	320000	69	STA RPOS1
0060	3EDF	70	CLEAR: MVI A, 0DFH ; SET CLEAR CODE TO CLEAR ON ENTRY.
0062	320000	71	STA CLEARA
0065	3ECC	72	MVI A, 0CCH ; SET FOR NO CLEAR ON EXIT.
0067	320000	73	STA CLEARB
006A	CD0000	74	CALL KDSPY2 ; GET MINUTE ENTRY.
006D	3A0000	75	LDA SUM2
0070	CD2C02	76	CALL MINCON ; CONVERT MINUTES TO SECONDS.
0073	220000	77	SHLD TWANT ; TWANT= MINUTE ENTRY IN SECONDS.
0076	3E94	78	MVI A, 94H ; SET KDSY2 FOR SECOND ENTRY.
0078	320000	79	STA POS2
007B	3E85	80	MVI A, 85H
007D	320000	81	STA POS1
0080	3E65	82	MVI A, 65H
0082	320000	83	STA RPOS1
0085	3ECC	84	CLEAR: MVI A, 0CCH ; SET FOR NO CLEAR ON ENTRY TO KDSPY2.
0087	320000	85	STA CLEARA
008A	3EDF	86	MVI A, 0DFH ; SET CLEAR CODE TO CLEAR ON EXIT.
008C	320000	87	STA CLEARB
008F	CD0000	88	CALL KDSPY2 ; GET SECONDS ENTRY.
0092	2A0000	89	LHLD TWANT
0095	3A0000	90	LDA SUM2 ; ADD SECONDS TO TWANT.
0098	4F	91	MOV C, A
0099	0600	92	MVI B, 00H
009B	89	93	DAD B
009C	220000	94	SHLD TWANT
009F	7C	95	MOV A, H ; TEST TO SEE IF TWANT = 0.
00A0	FE00	96	CPI 00H
00A2	C2AB00	97	JNZ TWNOT0
00A5	7D	98	MOV A, L
00A6	FE00	99	CPI 00H
00A0	CA4900	100	JZ GETTW ; IF IT IS, GET ANOTHER TWANT ENTRY.
00AB	210000	101	TWNOT0: LXI R, PLBYTE ; PREPARE TO TEST TWANT AGAINST LENGTH
00AE	3A0000	102	LDA DEB ; OF SELECTED PROCESS.
00B1	85	103	ADD L
00B2	6F	104	MOV L, A ; HL POINTS TO SELECTED PLBYTE(DEB)
00B3	7E	105	MOV A, H
00B4	CD2C02	106	CALL MINCON ; CONVERT PLBYTE TO SECONDS.
00B7	3A0100	107	LDA TWANT+1 ; COMPARE ROB OF TWANT TO ROB OF PLBYTE.

LOC	OBJ	LINE	SOURCE STATEMENT
00BA	BC	108	CMP H
00BB	DAD400	C 109	JC OVER ; A CY MEANS NO TWANT (NO PLBYTE.
00BE	C2C800	C 110	JNZ BLINK ; NO ZERO MEANS NO TWANT) NO PLBYTE.
00C1	3A0000	D 111	LDA TWANT ; NOW TEST LOBS.
00C4	BD	112	CMP L
00C5	DAD400	C 113	JC OVER ; A CY MEANS LO TWANT (LO PLBYTE.
00C8	CAD400	C 114	JZ OVER ; A ZERO MEANS = AND THATS O.K. TOO.
00CB	3ECB	115 BLINK:	MVI A, 0CBH ; BLINK ALPHA AS TWANT) PLB(DEB).
00CD	320000	E 116	STA ALPHA
00D0	CD0000	E 117	CALL LDELAY
00D3	C30000	C 118	JMP DFIND ; GO START OVER AGAIN.
00D6	AF	119 OVER:	XRA A ; INITIALIZE N (I/R COUNTER).
00D7	320000	E 120	STA MCNT
00DA	47	121	MOV H, A
00DB	4F	122	MOV L, A ; ZERO RUNNING SUM : TIMSUM.
00DC	220200	D 123	SHLD TIMSUM
00DF	4F	124	MOV C, A
00E0	210100	E 125	LXI H, IRCNT+1 ; SUM IRCNTS TO GET START ADDRESS IN
00E3	3A0000	E 126	LDA DEB ; IRRAM SPACE
00E6	FE01	127 SUMIB:	CPI 01H ; NOTE THAT IF PROCESS 1 IS SELECTED, YOU
00E8	CAF500	C 128	JZ NOSUM ; WISH TO START AT THE VERY BOTTOM OF IRRAM.
00EB	47	129	MOV B, A
00EC	79	130	MOV A, C
00ED	06	131	ADD M
00EE	4F	132	MOV C, A
00EF	78	133	MOV A, B
00F0	3D	134	DCB A
00F1	23	135	INX H
00F2	C3E600	C 136	JMP SUMIR
00F5	79	137 NOSUM:	MOV A, C ; RETRIEVE IRCNT SUM.
00F6	87	138	ADD A ; DOUBLE IT.
00F7	C601	139	ADI 01H ; ADD 1 TO GET OFF IRCNT.
00F9	320400	D 140	STA OFFSET ; OFFSET MUST BE ADDED TO IRRAM TO PUT
00FC	210000	E 141	LXI H, IRRAM ; YOU IN THE PROCESS IR BLOCK YOU'VE
00FF	5F	142	MOV E, A ; SELECTED.
0100	1600	143	MVI D, 00H ; DE = OFFSET.
0102	19	144	DAD D ; HL = IRRAM + OFFSET.
0103	23	145	INX H ; INCREMENT TO MOVE FROM R TO 1.
0104	7E	146	MOV A, M ; MOVE INTERVAL INTO A REG.
0105	CD2C02	C 147	CALL MINCON ; CONVERT TO SECONDS.
0108	220400	D 148	SHLD ITEMP
010B	2A0000	D 149 WRDTST:	LHLD TWANT ; IS TIMSUM >= TWANT ?
010E	3A0300	D 150	LDA TIMSUM+1 ; TEST BOBS FIRST.
0111	BC	151	CMP H
0112	DA1F01	C 152	JC TSLSTV
0115	C29A01	C 153	JNZ SHOW
0118	3A0200	D 154	LDA TIMSUM ; TEST LOBS LAST.
011B	BD	155	CMP L
011C	D29A01	C 156	JNC SHOW
011F	2A0400	D 157 TSLSTV:	LHLD ITEMP ; IS TIMSUM = ITEMP ?
0122	3A0300	D 158	LDA TIMSUM+1 ; TEST BOBS FIRST.
0125	BC	159	CMP H
0126	C23001	C 160	JNZ ADDWRD
0129	3A0200	D 161	LDA TIMSUM ; TEST LOBS LAST.

LOC	OBJ	LINE	SOURCE STATEMENT
012C	BD	162	CMP I
012D	CA5001	C 163	JZ INCM
0130	3A0000	E 164	ADDWRD: LDA RAMPT ; COME HERE IF TIMSUM NOT EQUAL TO ITEM.
0133	3C	165	INR A
0134	320000	E 166	STA RAMPT ; RAMPT = RAMPT + 1
		167 ;	
		168 ;	CALCULATE R(N) POINTER BY: IRBRAM + OFFSET + 2(NCNT)
		169 ;	
0137	3A0000	E 170	LDA NCNT
013A	07	171	ADD A ; DOUBLE NCNT
013B	47	172	MOV B, A
013C	3A0400	D 173	LDA OFFSET
013F	00	174	ADD B ; A = OFFSET + 2(NCNT).
0140	0400	175	MVI B, 00H
0142	4F	176	MOV C, A ; BC = OFFSET + 2(NCNT).
0143	210000	E 177	LXI H, IRBRAM
0146	09	178	DAD B ; HL = IRBRAM + OFFSET + 2(NCNT).
0147	CD0000	E 179	CALL INVERT ; CONVERT RATE FROM SCAN/MIN TO SEC/SCAN.
014A	0F	180	RRC ; DIVIDE RESULT BY 2.
014B	4F	181	MOV C, A ; C = RATE N IN SEC/SCAN.
014C	0400	182	MVI B, 00H
014E	2A0200	D 183	LHLD TIMSUM
0151	09	184	DAD B ; HL = TIMSUM + R(N)
0152	220200	D 185	SHLD TIMSUM
0155	C30B01	C 186	JMP WRDTST
0158	3A0000	E 187	INCM: LDA NCNT ; N = N + 1
015B	3C	188	INR A
015C	320000	E 189	STA NCNT
015F	210000	E 190	LXI H, IRCNT ; IRCNT(DEB) = N ?
0162	3A0000	E 191	LDA DEB
0165	05	192	ADD I
0166	4F	193	MOV L, A ; HL POINTS TO IRCNT(DEB).
0167	3A0000	E 194	LDA NCNT
016A	BE	195	CMP M
016B	CA8B01	C 196	JZ LASSCN
016E	07	197	ADD A ; DOUBLE NCNT.
016F	47	198	MOV B, A ; STORE COPY OF NCNT IN B REG.
0170	3A0400	D 199	LDA OFFSET
0173	00	200	ADD B ; A = 2(NCNT) + OFFSET.
0174	3C	201	INR A ; A = 2(NCNT) + OFFSET + 1
0175	4F	202	MOV C, A
0176	0400	203	MVI B, 00H ; BC = 2(NCNT) + OFFSET + 1
0178	210000	E 204	LXI H, IRBRAM
017B	09	205	DAD B ; HL = IRBRAM + 2(NCNT) + OFFSET + 1
017C	7E	206	MOV A, M
017D	CD2C02	C 207	CALL NINCON ; CONVERT I(N) TO SECONDS.
0180	EB	208	XCHG ; STORE THID IN DE.
0181	2A0400	D 209	LHLD ITEMP
0184	19	210	DAD D ; HL = ITEMP + I(N).
0185	220400	D 211	SHLD ITEMP ; ITEMP = ITEMP + I(N).
0188	C33001	C 212	JMP ADDWRD
018B	210000	E 213	LASSCN: LXI H, NCNT ; N = N-1.
018E	35	214	DCR M
018F	216502	C 215	LXI H, LSCNM ; PROMPT "LAST SCAN"

LOC	OBJ	LINE	SOURCE STATEMENT
0192	060B	216	MVI B, 0BH
0194	CD0000	E 217	CALL DWRITE
0197	CD0000	E 218	CALL LDELAY
019A	CD0000	E 219	SHOW: CALL DSHOW ; DISPLAY DENSITY BLOCK POINTED TO BY RAMPT
		220 ;	
		221 ;	NOW DISPLAY CONTENTS OF TIMSUM AS: MIN:SEC .
		222 ;	
019D	AF	223	XRA A ; ZERO A AND SET CY TO 0.
019E	320900	D 224	STA MIN ; ZERO MINUTE COUNTER.
01A1	2A0200	D 225	LHLD TIMSUM ; STORE COPY OF TIMSUM IN MINUTE.
01A4	220000	E 226	SHLD TMPNUM ; STORE COPY FOR LATER RETRIVAL.
01A7	220700	D 227	SHLD MINUTE
01AA	2A0700	D 228	CHANGE: LHLD MINUTE ; STORE COPY OF MINUTE IN TIMSUM.
01AD	220200	D 229	SHLD TIMSUM
01B0	063C	230	MVI B, 3CH ; B=60 (60 SECONDS IN A MINUTE ASSHOLE!).
01B2	3A0700	D 231	LDA MINUTE ; LOB FIRST: MINUTE = MINUTE - 60.
01B5	98	232	SBB B ; A = MINUTE - 60.
01B6	320700	D 233	STA MINUTE
01B9	D2C801	C 234	JNC NOBROW ; NO CY IMPLIES NO NEED TO BORROW.
01BC	0600	235	MVI B, 00H ; IF CY, BORROW 1 FROM HOB OF MINUTE.
01BE	3A0000	D 236	LDA MINUTE+1
01C1	98	237	SBB B
01C2	320000	D 238	STA MINUTE+1
01C5	FAD201	C 239	JM SHOWTS ; IF HOB NEGATIVE, THEN EXIT.
01C8	3A0900	D 240	NOBROW: LDA MIN ; INCREMENT MINUTE COUNTER.
01CB	3C	241	INR A
01CC	320900	D 242	STA MIN
01CF	C3AA01	C 243	JMP CHANGE
01D2	215902	C 244	SHOWTS: LXI H, TPAST ; PROMPT "TIME ELAPSED"
01D5	060D	245	MVI B, 0DH
01D7	CD0000	E 246	CALL DWRITE
01DA	3A0900	D 247	LDA MIN ; CONVERT MIN TO ASCII.
01DD	CD0000	E 248	CALL BINASC
01E0	7A	249	MOV A, D ; D CONTAINS 10'S PLACE ASCII.
01E1	D630	250	SUI 30H ; CONVERT TO SIMPLE BCD.
01E3	41	251	MOV B, C ; HIDE C FOR A SEC.
01E4	CD0000	E 252	CALL LOOKUP ; CONVERT BCD TO DISPLAY CODE.
01E7	51	253	MOV D, C ; D NOW CONTAINS 10'S PLACE MINUTE CODE.
01E8	78	254	MOV A, B ; DO NOW FOR 1'S PLACE MINUTE COUNT.
01E9	D630	255	SUI 30B
01EB	CD0000	E 256	CALL LOOKUP
01EE	210000	E 257	LXI H, COMM1
01F1	3691	258	MVI M, 91H ; SET CONTROL DISPLAY FOR POS. 1 (A1).
01F3	210000	E 259	LXI H, DATA1
01F6	72	260	MOV M, D ; SEND 10'S PLACE MINUTE CODE TO DISPLAY.
01F7	71	261	MOV M, C ; SEND 1'S PLACE MINUTE CODE TO DISPLAY.
01F8	36FF	262	MVI M, 0FFH ; SEND BLANK TO SEPERATE MIN FROM SEC.
01FA	E5	263	PUSH H ; SAVE DISPLAY ADDRESS.
01FB	3A0200	D 264	LDA TIMSUM ; DO NOW FOR SECOND COUNT.
01FE	CD0000	E 265	CALL BINASC
0201	7A	266	MOV A, D ; D CONTAINS 10'S PLACE SECOND COUNT.
0202	D630	267	SUI 30B
0204	41	268	MOV B, C
0205	CD0000	E 269	CALL LOOKUP ; CONVERT 10'S BCD TO DISPLAY CODE.

162

LOC	OBJ	LINE	SOURCE STATEMENT
0200	51	270	MOV D, C
0209	70	271	MOV A, B
020A	B438	272	SUI 30H
020C	CD0000	E 273	CALL LOOKUP
020F	E1	274	POP H ; RETREIVE DISPLAY ADDRESS.
0210	72	275	MOV M, D ; SEND 10'S PLACE SECOND COUNT.
0211	71	276	MOV M, C ; SEND 1'S PLACE SECOND COUNT TO DISPLAY.
		277 ;	
0212	1A0000	E 278	LHLD TMPNUM ; RESTORE TIMSUM TO ITS ORIGINAL VALUE.
0215	220200	D 279	SHLD TIMSUM
		180 ;	
		181 ;	AT THIS POINT DENSITY VALUES ALONG WITH TIME ELASPED SHOULD BE
		182 ;	ON THE TWO DISPLAYS.
		183 ;	
0218	C30000	E 284	JMP UPORDM ; GO TO SELECT DATA UP, DATA DOWN OR RETURN.
		285 ;	
		186 ;	*****
		287 ;	
		188 ;	SUBROUTINES FOLLOW
		289 ;	
		190 ;	*****
		291 ;	
		192 ;	
		193 ;	*****
		294 ;	
		195 ;	** NODATA **
		296 ;	
		197 ;	THIS SUBROUTINE DISPLAYS " NO DATA " ON THE ALPHA DISPLAY.
		298 ;	IT ALSO POPS THE STACK AND PLACES A NEW RETURN ADDRESS OF
		199 ;	DFIND ONTO THE STACK.
		300 ;	
		301 ;	*****
		302 ;	
021B	113A01	C 303	NODATA: LXI H, NODATM ; START ADDRESS OF "NO DATA"
021E	060A	304	MVI B, 0AH
0220	CD0000	E 305	CALL DWRITE
0223	CD0000	E 306	CALL LDELAY
0226	E1	307	POP H
0227	210000	C 308	LXI H, DFIND
022A	E5	309	PUSH H
022B	C9	310	RET
		311 ;	
		312 ;	*****
		313 ;	
		314 ;	** MINCON **
		315 ;	
		316 ;	THIS ROUTINE CONVERTS MINUTES TO SECONDS.
		317 ;	
		318 ;	PASS: A CONTAINS THE NUMBER OF MINUTES.
		319 ;	
		320 ;	RETURN: HL CONTAINS THE RESULT OF A X 60 IN SECONDS.
		321 ;	
		322 ;	*****
		323 ;	

LOC	OBJ	LINE	SOURCE STATEMENT
022C	210000	324	MINCON: LXI H, 0000H
022F	013C00	325	LXI B, 003CH ; LOAD B WITH 40.
0232	FE00	326	MINADD: CPI 00H
0234	C0	327	RZ
0235	09	328	DAD B
0236	3D	329	DCR A
0237	C33202	C 330	JMP MINADD
		331 ;	
		332 ;	DISPLAY MESSAGES FOLLOW:
		333 ;	
023A	20	334	NODATM: DB 20H ; BLANK
023B	20	335	DB 20H ; BLANK
023C	0E	336	DB 0EH ; M
023D	0F	337	DB 0FH ; O
023E	20	338	DB 20H ; BLANK
023F	04	339	DB 04H ; D
0240	01	340	DB 01H ; A
0241	14	341	DB 14H ; T
0242	01	342	DB 01H ; A
0243	10	343	TWANTM: DB 10H ; P
0244	12	344	DB 12H ; R
0245	0F	345	DB 0FH ; O
0246	03	346	DB 03H ; C
0247	05	347	DB 05H ; E
0248	13	348	DB 13H ; S
0249	13	349	DB 13H ; S
024A	20	350	DB 20H ; BLANK
024B	23	351	DB 23H ; 0
024C	3D	352	DB 3DH ; =
024D	3F	353	DB 3FH ; ?
024E	14	354	TWMESS: DB 14H ; T
024F	09	355	DB 09H ; I
0250	0D	356	DB 0DH ; M
0251	05	357	DB 05H ; E
0252	20	358	DB 20H ; BLANK
0253	17	359	DB 17H ; W
0254	01	360	DB 01H ; A
0255	0E	361	DB 0EH ; M
0256	14	362	DB 14H ; T
0257	05	363	DB 05H ; E
0258	04	364	DB 04H ; D
0259	14	365	TPAST: DB 14H ; T
025A	09	366	DB 09H ; I
025B	0D	367	DB 0DH ; M
025C	05	368	DB 05H ; E
025D	20	369	DB 20H ; BLANK
025E	05	370	DB 05H ; E
025F	0C	371	DB 0CH ; L
0260	01	372	DB 01H ; A
0261	13	373	DB 13H ; S
0262	10	374	DB 10H ; P
0263	05	375	DB 05H ; E
0264	04	376	DB 04H ; D
0265	20	377	LSCNM: DB 20H ; BLANK

164

LOC	OBJ	LINE	SOURCE STATEMENT
0266	0C	378	DB 0CH ; L
0267	01	379	DB 01H ; A
0268	13	380	DB 13H ; S
0269	14	381	DB 14H ; T
026A	20	382	DB 20H ; BLANK
026B	13	383	DB 13H ; S
026C	03	384	DB 03H ; C
026D	01	385	DB 01H ; A
026E	0E	386	DB 0EH ; M
		387 ;	
		388 DSEC	
		389 ;	
0000		390 TWANT: DS 02H ; TIME WANTED IN SECONDS.	
0002		391 TIMSUM: DS 02H ; 2 BYTE RUNNING SUM OF RATES (ALSO IN SECONDS).	
0004		392 ITEMP: DS 02H ; LENGTH OF I(N) IN SECONDS.	
0006		393 OFFSET: DS 01H ; USED IN JUMPING I/R BLOCKS.	
0007		394 MINUTE: DS 02H ; PARTNER TO TIMSUM.	
0009		395 MIN: DS 01H ; MINUTE COUNTER.	
		396 ;	
		397 END	

PUBLIC SYMBOLS

DFIND	C 0000	ITEMP	D 0004	MIN	D 0009	MINCON	C 022C	OFFSET	D 0006	SHOW	C 019A	TIMSUM	D 0002
-------	--------	-------	--------	-----	--------	--------	--------	--------	--------	------	--------	--------	--------

EXTERNAL SYMBOLS

ALPHA	E 0000	BINASC	E 0000	CLEARA	E 0000	CLEARB	E 0000	COMM1	E 0000	DATA1	E 0000	DER	E 0000
DSHOW	E 0000	DWRITE	E 0000	INVERT	E 0000	IRBRAM	E 0000	IRCNT	E 0000	KDSPY2	E 0000	LDELAY	E 0000
LOOKUP	E 0000	MINCNT	E 0000	NCNT	E 0000	PLYTE	E 0000	POS1	E 0000	POS2	E 0000	RAMPT	E 0000
RPOS1	E 0000	STABCT	E 0000	SUM2	E 0000	TMPNUM	E 0000	UPORDN	E 0000				

USER SYMBOLS

ADDED	C 0043	ADDWRD	C 0130	ALPHA	E 0000	BINASC	E 0000	BLINK	C 00CB	CHANGE	C 01AA	CLEAR	C 0060
CLEARA	E 0000	CLEARB	E 0000	CLEARC	C 0085	COMM1	E 0000	DATA1	E 0000	DER	E 0000	DFIND	C 0000
DSHOW	E 0000	DWRITE	E 0000	GETTV	C 0049	INCN	C 0158	INVERT	E 0000	IRBRAM	E 0000	IRCNT	E 0000
ITEMP	D 0004	KDSPY2	E 0000	LASSCN	C 0108	LDELAY	E 0000	LOOKUP	E 0000	LSCNM	C 0265	MIN	D 0009
MINADD	C 0232	MINCNT	E 0000	MINCON	C 022C	MINUTE	D 0007	NCNT	E 0000	NORROW	C 01C8	NODATA	C 021R
MODATH	C 023A	MOSUM	C 00F5	OFFSET	D 0006	OVER	C 00D6	PLYTE	E 0000	POS1	E 0000	POS2	E 0000
RAMPT	E 0000	RPOS1	E 0080	SHOW	C 019A	SHOWTS	C 01D2	STABCT	E 0000	SUM2	E 0000	SUMIR	C 00E6
TIMSUM	D 0002	TMPNUM	E 0000	TPAST	C 0259	TSLSW	C 011F	TWANT	D 0000	TWANTH	C 0243	TWMESS	C 024E
TWNOT0	C 00AB	UPORDN	E 0000	WRDTST	C 0108								

ASSEMBLY COMPLETE, NO ERRORS

LOC OBJ LINE SOURCE STATEMENT

```

1 ;*****
2 ;
3 ;     UPORDN  (UP OR DOWN) IS A KEYBOARD READ/COMMAND ROUTINE THAT
4 ;     PERFORMS 5 DIFFERENT FUNCTIONS DEPENDING WHAT KEY IS DEPRESSED.
5 ;
6 ;     THESE FUNCTIONS ARE:
7 ;
8 ;     KEY      FUNCTION
9 ;
10 ;    UP        THIS MOVES USER UP TO NEXT SEQUENTIAL DENSITY SCAN.
11 ;
12 ;    DOWN      THIS MOVES USER DOWN TO PRECEEDING DENSITY SCAN.
13 ;
14 ;    PRT       THIS WILL PRINT ALL ELEVEN DENSITY VALUES ON THE
15 ;              40-COLUMN PRINTER THAT ARE CURRENTLY DISPLAYED.
16 ;
17 ;    NXT       RETURNS USER TO ENTER PROCESS PROMPT AS TO ALLOW
18 ;              FOR VIEWING DATA IN A DIFFERENT NUMBER PROCESS.
19 ;
20 ;    ENTER-ENTER  TWO ENTERS IN A ROW WILL RETURN USEB TO PROCESS
21 ;              INFORMATION ENTRY ROUTINE TO ALLOW FOR ANOTHER PROCESS.
22 ;
23 ;*****
24 ;
25 ;     NAME      UPORDN
26 ;
27 ;     PUBLIC    UPORDN
28 ;
29 ;     EXTRN     ITEMP,TIMSUM,RAMPT,NCNT,OFFSET,IRBRAM,DEB,FIFOCR,COMM2
30 ;     EXTRN     DATA2,MANFLC,START,DFIND,INVERT,SHOW,IRCNT,MINCON,DWRITE
31 ;     EXTRN     DRAM,MIDRAM,RPOS1,POS1,POS2,CLEARA,CLEARB,MIN
32 ;     EXTRN     TDSTOR,PWRITE,PSEND,COMM1,DATA1,MANMOD,SUM2,KDSPY2,RPA2
33 ;
34 ;     CSEC
35 ;
0000 CDF902  C 36 UPORDN: CALL  RESET          ; RESET RESETS PARAMETERS USED IN KDSPY2.
0003 CD0000  E 37          CALL  FIFOCR        ; CLEAR FIFO.
0006 210000  E 38 KLOOK: LXI   H,          COMM2  ; #279 STATUS WORD ADDRESS.
0009 7E      39          MOV    A,          H
000A E403    40          ANI    03H          ; MASK OUT 5 HI-ORDER BITS.
000C FE01    41          CPI    01H          ; IF SOMETHING IN FIFO, ZERO WILL SET.
000E C20400  C 42          JNZ   KLOOK        ; IF NOT, KEEP LOOKING.
0011 3440    43          MVI    H,          40H    ; SET FIFO AS READ SOURCE.
0013 210000  E 44          LXI   H,          DATA2 ; DATA ADDRESS OF EXPANSION #279.
0016 7E      45          MOV    A,          H    ; PUT FIFO CONTENTS INTO A REG.
0017 E40F    46          ANI    0FH          ; MASK OUT HI-ORDER NIBBLE.
0019 FE0F    47          CPI    0FH          ; IS IT ENTER KEY?
001B C22000  C 48          JNZ   NEKEY
001E 3A0000  E 49          LDA    MANFLC      ; CHECK FOR TWO CONCURRENT ENTERS.
0021 3C      50          INR    A
0022 320000  E 51          STA    MANFLC
0025 FE02    52          CPI    02H
0027 CA5903  C 53          JZ     CLEAN        ; YES, 2 ENTERS; GO TO MANUAL MODE TEST.

```

166

LOC	ORJ	LINE	SOURCE STATEMENT
002A	C30400	C 54	JMP KLOOK ; GO LOOK FOR ANOTHER ENTRY.
002D	47	55	NEKEY: MOV B, A
002E	AF	56	TRA A
002F	320000	E 57	STA MANFLG ; SECOND DEPRESSION NOT ENTER, RESET MANFLG.
0032	70	58	MOV A, B
0033	FE0E	59	CPI 00H
0035	CA4A00	C 60	JZ MEXITUP ; THIS IS INCREMENT UP ROUTINE.
0030	FE0D	61	CPI 00H
003A	CAC000	C 62	JZ MEXITDN ; THIS DECREMENT DOWN ROUTINE.
003D	FE0C	63	CPI 00H
003F	CAS201	C 64	JZ PRINT ; THIS PRINTS CURRENT DISPLAY CONTENTS.
0042	FE00	65	CPI 00H
0044	CA0000	E 66	JZ DFIND ; GO TO SELECT NEW PROCESS #.
0047	C30000	C 67	JMP UPORDN ; ALL OTHER ENTRIES INVALID, GO LOOK AGAIN.
004A	2A0000	E 68	NEXTUP: LRLD ITEM ; IS TIMSUM = ITEM?
004D	3A0100	E 69	LDA TIMSUM+1
0050	BC	70	CMP R
0051	C25R00	C 71	JNZ NOTEQ ; IF ROR'S NOT EQUAL, THEN LEAVE.
0054	3A0000	E 72	LDA TIMSUM
0057	BD	73	CMP L
0050	CA0300	C 74	JZ IMXM
005R	3A0000	E 75	NOTEQ: LDA RAMPT ; COME HERE IF TIMSUM NOT EQ TO ITEM.
005E	3C	76	INR A
005F	320000	E 77	STA RAMPT ; RAMPT = RAMPT + 1.
		78 ;	
		79 ;	CALCULATE R(N) POINTER BY: IRRHAM + OFFSET + 2(NCNT)
		80 ;	
0062	3A0000	E 81	LDA NCNT
0065	07	82	ADD A ; DOUBLE NCNT
0066	47	83	MOV R, A
0067	3A0000	E 84	LDA OFFSET
006A	00	85	ADD R ; A = OFFSET + 2(NCNT).
006B	0600	86	MVI R, 00H
006D	4F	87	MOV C, A ; RC = OFFSET + 2(NCNT).
006E	210000	E 88	LXI H, IRRHAM
0071	09	89	DAD R ; RL = IRRHAM + OFFSET + 2(NCNT).
0072	CD0000	E 90	CALL INVERT ; CONVERT RATE FROM SCAN/MIN TO SEC/SCAN.
0075	0F	91	RRC ; DIVIDE RESULT BY 2.
0076	4F	92	MOV C, A
0077	0600	93	MVI B, 00H
0079	2A0000	E 94	LRLD TIMSUM
007C	09	95	DAD B ; RL = TIMSUM + R(N).
007D	220000	E 96	SRLD TIMSUM ; UPDATE TIMSUM.
0080	C30000	E 97	JMP SHOW ; DISPLAY DEMSITIES POINTED TO BY RAMPT.
0083	210000	E 98	INXM: LXI R, IRCNT ; IRCNT(DEB) = N ?
0086	3A0000	E 99	LDA DER
0089	05	100	ADD L
008A	4F	101	MOV L, A ; RL POINTS TO IRCNT(DEB).
008B	56	102	MOV D, H
008C	15	103	DCR D
008D	3A0000	E 104	LDA NCNT
0090	BA	105	CMP D
0091	CAR500	C 106	JZ NOMDATA ; IF EQUAL, GO DISPLAY "NO MORE DATA"
0094	3C	107	INR A ; IF NOT EQUAL, ITEM = ITEM + 1(N).

LOC	OBJ	LINE	SOURCE STATEMENT
0095	320000	E 108	STA NCNT
0098	87	109	ADD A ; DOUBLE NCNT.
0099	47	110	MOV B, A ; GENERATE ITEMP = ITEMP + I(N).
009A	3A0000	E 111	LDA OFFSET
009D	00	112	ADD B ; A = OFFSET + 2(NCNT).
009E	3C	113	INR A ; A = OFFSET + 2(NCNT) + 1.
009F	4F	114	MOV C, A
00A0	0600	115	MVI B, 00H ; BC = OFFSET + 2(NCNT) + 1.
00A2	210000	E 116	LXI H, IRBRAM
00A5	09	117	DAD B ; HL = OFFSET + 2(NCNT) + 1 + IRBRAM.
00A6	7E	118	MOV A, M
00A7	CD0000	E 119	CALL MINCON ; CONVERT I(N) TO SECONDS.
00AA	EB	120	CHG ; STORE THID IN DE REGISTER.
00AB	2A0000	E 121	LHLD ITEMP
00AE	19	122	DAD D
00AF	220000	E 123	SHLD ITEMP ; ITEMP = ITEMP + I(N).
00B2	C35B00	C 124	JMP NOTEQ
00B5	217303	C 125	NMDATA: LXI H, NMDATH ; DISPLAY "NO MORE DATA"
00B8	060D	126	MVI B, 0DH ; LOAD TEXT COUNTER.
00BA	CD0000	E 127	CALL DWRITE
00BD	C30000	C 128	JMP UPORDN ; LOOK FOR MORE KEY COMMANDS.
00C0	AF	129	NEXTDN: XRA A ; CONVERT ITEMP BACK TO MINUTES.
00C1	320000	E 130	STA MIN
00C4	2A0000	E 131	LHLD ITEMP
00C7	EB	132	CHG ; PUT COPY OF ITEMP INTO DE.
00C8	7B	133	OMT: MOV A, E ; LOB OF ITEMP INTO A.
00C9	063C	134	MVI B, 3CH ; B=60.
00CB	98	135	SBB B
00CC	5F	136	MOV E, A ; UPDATE E REG.
00CD	D2D800	C 137	JNC NOBARO
00D0	7A	138	MOV A, D ; NOW DO HOB'S.
00D1	0600	139	MVI B, 00H
00D3	98	140	SBB B
00D4	57	141	MOV D, A ; UPDATE D REG.
00D5	FADF00	C 142	JM DONE ; WHEN HIB GOES NEGATIVE, THEN DONE.
00D8	210000	E 143	NOBARO: LXI H, MIN ; INCREMENT MINUTE COUNTER.
00DB	34	144	INR M
00DC	C3C800	C 145	JMP OMT
00DF	3A0000	E 146	DONE: LDA NCNT ; PERFORM LITEMP = ITEMP = I(N) IN MINUTES.
00E2	87	147	ADD A ; GET I(N) FIRST.
00E3	47	148	MOV B, A
00E4	3A0000	E 149	LDA OFFSET
00E7	00	150	ADD B ; A=2(NCNT) + OFFSET.
00E8	3C	151	INR A ; A = 2(NCNT) + OFFSET + 1.
00E9	4F	152	MOV C, A
00EA	0600	153	MVI B, 00H
00EC	210000	E 154	LXI H, IRBRAM
00EF	09	155	DAD B ; HL = IRBRAM + 2(NCNT) + OFFSET + 1.
00F0	3A0000	E 156	LDA MIN ; PLACE MINUTE VERSION OF ITEMP IN A.
00F3	96	157	SUB M ; A = ITEMP - I(N) ; ALL IN MINUTES.
00F4	CD0000	E 158	CALL MINCON ; CONVERT RESULT TO SECONDS.
00F7	220000	D 159	SHLD LITEMP ; LITEMP = ITEMP - I(N), NOW IN SECONDS.
00FA	EB	160	CHG ; PLACE COPY IN DE.
00FB	2A0000	E 161	LHLD TIMSUM ; TIMSUM = LITEMP ???

LOC	OBJ	LINE	SOURCE STATEMENT
00FE 7C		162	MOV A, H
00FF BA		163	CMP D ; TEST MOB'S FIRST.
0100 C21A01	C	164	JNZ DECKS ; GO TO DECKS IF MOB'S NOT EQUAL.
0103 7D		165	MOV A, L
0104 BB		166	CMP E ; TEST LOB'S.
0105 C21A01	C	167	JNZ DECKS ; GO TO DECKS IF LOB'S NOT EQUAL.
0106 3A0000	E	168	NE00: LDA NCNT ; N = 0 ?????
010B FE00		169	CPI 00H
010D CAB500	C	170	JZ NMDATA ; IF NOT EQUAL TO 0, THEN N = N-1.
0110 3D		171	DCR A
0111 320000	E	172	STA NCNT
0114 2A0000	D	173	LHLD LITEMP ; SET ITEMP = LITEMP.
0117 220000	E	174	SHLD ITEMP
011A 210000	E	175	DECKS: LXI H, RAMPT ; RAMPT = RAMPT - 1.
011D 35		176	DCR H
011E 3A0000	E	177	LDA NCNT
0121 07		178	ADD A ; NOW DO TIMSUM = TIMSUM - R(N)
0122 47		179	MOV B, A
0123 3A0000	E	180	LDA OFFSET
0126 00		181	ADD B ; A = 2(NCNT) + OFFSET
0127 0600		182	MVI B, 00H
0129 4F		183	MOV C, A
012A 210000	E	184	LXI H, 1RBRAM
012D 09		185	DAD B ; HL = 1RBRAM + OFFSET + 2(NCNT)
012E CD0000	E	186	CALL INVERT
0131 0F		187	RRC
0132 47		188	SUBBS: MOV B, A ; NOW DO ACTUAL SUBTRACTION.
0133 2A0000	E	189	LHLD TIMSUM
0136 7D		190	MOV A, L
0137 98		191	SBB B
0138 4F		192	MOV L, A ; UPDATE L WITH RESULT.
0139 D24101	C	193	JNC OUTIT
013C 7C		194	MOV A, H ; NOW DO MOB'S.
013D 0600		195	MVI B, 00H
013F 98		196	SBB B
0140 47		197	MOV H, A ; UPDATE H WITH RESULT.
0141 220000	E	198	OUTIT: SHLD TIMSUM ; TIMSUM = 0 ?????
0144 FE00		199	CPI 00H ; RECALL A = MOB OF TIMSUM.
0146 C20000	E	200	JNZ SHOW
0149 7D		201	MOV A, L
014A FE00		202	CPI 00H
014C CAB500	C	203	JZ NMDATA
014F C30000	E	204	JMP SHOW
0152 3E16		205	PRINT: MVI A, 16H ; SELECT 4 EXTRA DOT ROWS OF SPACE.
0154 CD0000	E	206	CALL PSEND
0157 3E04		207	MVI A, 04H
0159 CD0000	E	208	CALL PSEND
015C 218903	C	209	LXI H, PRONUM ; PRINT " PROCESSI(DEV).
015F 0E0A		210	MVI C, 0AH
0161 CD0000	E	211	CALL PWRITE
0164 3A0000	E	212	LDA DEB
0167 C630		213	ADI 30H ; CONVERT PROCESS NUMBER TO ASCII.
0169 CD0000	E	214	CALL PSEND
016C 219103	C	215	LXI H, TPAST ; LOAD NOW " XXITIMEXELASPEDIX " "

LOC	OBJ	LINE	SOURCE STATEMENT
016F	0E11	216	MVI C, 11H
0171	CD0000	E 217	CALL PWRITE
0174	210000	E 218	LXI H, COMM1 ; LOAD CONTROL DISPLAY COMMAND ADDRESS.
0177	3671	219	MVI M, 71H ; GET MINUTES ELAPSED.
0179	210000	E 220	LXI H, DATA1
017C	7E	221	MOV A, M
017D	ES	222	PUSH H ; SAVE DATA ADDRESS.
017E	CD2603	C 223	CALL REVERT ; CHANGE CONTROL DISPLAY CODE TO BCD.
0181	79	224	MOV A, C
0182	C630	225	ADI 30H ; CONVERT TO ASCII.
0184	CD0000	E 226	CALL PSEND ; SEND TO PRINTER BUFFER.
0187	E1	227	POP H
0188	7E	228	MOV A, M ; GET 1'S PLACE MINUTE COUNT.
0189	ES	229	PUSH H
018A	CD2603	C 230	CALL REVERT
018D	79	231	MOV A, C
018E	C630	232	ADI 30H
0190	CD0000	E 233	CALL PSEND
0193	3E20	234	MVI A, 20H ; SPACE
0195	CD0000	E 235	CALL PSEND
0198	3E4D	236	MVI A, 4DH ; H
019A	CD0000	E 237	CALL PSEND
019D	3E49	238	MVI A, 49H ; I
019F	CD0000	E 239	CALL PSEND
01A2	3E4E	240	MVI A, 4EH ; M
01A4	CD0000	E 241	CALL PSEND
01A7	E1	242	POP H
01A8	7E	243	MOV A, M ; GET BLANK BETWEEN MINUTES AND SECONDS.
01A9	3E20	244	MVI A, 20H
01AB	CD0000	E 245	CALL PSEND ; SEND BLANK TO PRINT BUFFER.
01AE	7E	246	MOV A, M ; GET 10'S PLACE SECOND COUNT.
01AF	ES	247	PUSH H
01B0	CD2603	C 248	CALL REVERT
01B3	79	249	MOV A, C
01B4	C630	250	ADI 30H
01B6	CD0000	E 251	CALL PSEND
01B9	E1	252	POP H
01BA	7E	253	MOV A, M ; GET 1'S PLACE SECOND COUNT.
01BB	CD2603	C 254	CALL REVERT
01BE	79	255	MOV A, C
01BF	C630	256	ADI 30H
01C1	CD0000	E 257	CALL PSEND
01C4	3E20	258	MVI A, 20H
01C6	CD0000	E 259	CALL PSEND ; SEND BLANK.
01C9	3E73	260	MVI A, 73H ; S
01CB	CD0000	E 261	CALL PSEND
01CE	3E45	262	MVI A, 45H ; E
01D0	CD0000	E 263	CALL PSEND
01D3	3E43	264	MVI A, 43H ; C
01D5	CD0000	E 265	CALL PSEND
01D8	3E17	266	MVI A, 17H ; PRINT LINE BUFFER CONTENTS.
01DA	CD0000	E 267	CALL PSEND
01DD	3E01	268	MVI A, 01H ; INITIALIZE STEP COUNTER.
01DF	320200	D 269	STA STEP

LOC	OBJ	LINE	SOURCE STATEMENT
01E2	AF	270	XRA A
01E3	320000	E 271	STA TDSTOR ; SET MSD OF STEP 1 = 0
01E6	210000	E 272	LXI H, COMM2
01E9	3670	273	MVI M, 70H ; SET FOR RAM READ AT 000 (AI) .
01EB	210000	E 274	LXI H, DATA2
01EE	7E	275	MOV A, M ; READ POSITION 000.
01EF	1F	276	RAR
01F0	1F	277	RAR
01F1	1F	278	RAR
01F2	1F	279	RAR
01F3	E60F	280	ANI 0FB
01F5	320100	E 281	STA TDSTOR+1
01F8	7E	282	MOV A, M
01F9	1F	283	RAR
01FA	1F	284	RAR
01FB	1F	285	RAR
01FC	1F	286	RAR
01FD	E60F	287	ANI 0FH
01FF	320200	E 288	STA TDSTOR+2
0202	CD8F02	C 289	FIRST: CALL DATPRT ; PRINT 1ST STEP ON PRINTER.
0205	210000	E 290	LXI H, DATA2
0208	CD6B02	C 291	CALL LOADA ; 3RD.
020B	CD6B02	C 292	CALL LOADA ; 5TH.
020E	CD6B02	C 293	CALL LOADA ; 7TH.
0211	CD6B02	C 294	CALL LOADA ; 9TH.
0214	7E	295	SECB: MOV A, M ; GET 11TH STEP AND CHANGE OVER TO SECTION B.
0215	1F	296	RAR
0216	1F	297	RAR
0217	1F	298	RAR
0218	1F	299	RAR
0219	E60F	300	ANI 0FR
021B	320000	E 301	STA TDSTOR
021E	7E	302	MOV A, M
021F	1F	303	RAR
0220	1F	304	RAR
0221	1F	305	RAR
0222	1F	306	RAR
0223	E60F	307	ANI 0FH
0225	320100	E 308	STA TDSTOR+1
0228	210000	E 309	LXI H, COMM2
022B	3670	310	MVI M, 70R
022D	210000	E 311	LXI H, DATA2
0230	7E	312	MOV A, M ; SECTION B, RAM ADDRESS 000
0231	E60F	313	ANI 0FH
0233	320200	E 314	STA TDSTOR+2
0236	CD8F02	C 315	CALL DATPRT ; PRINT 11TH STEP.
0239	210000	E 316	REST: LXI H, DATA2
023C	CD5302	C 317	CALL LOADB ; 13TH.
023F	CD5302	C 318	CALL LOADB ; 15TH.
0242	CD5302	C 319	CALL LOADB ; 17TH.
0245	CD5302	C 320	CALL LOADB ; 19TH.
0248	CD5302	C 321	CALL LOADB ; 21ST.
024B	3E17	322	MVI A, 17H ; PRINT ONE BLANK LINE TO SEPERATE DATA BLOCKS.
024D	CD0000	E 323	CALL PSEND

171

LOC	OBJ	LINE	SOURCE STATEMENT
0250	C30000	C 324	JMP UPORDM
		325 ;	
		326 ;	*****
		327 ;	
		328 ;	*** LOADB ***
		329 ;	
		330 ;	READS 3 DISPLAY LOCATIONS FROM DENSITY DISPLAY, STORES THEM
		331 ;	IN TDSTOR,+1,+2 AND THEN PRINTS IT BY CALLING DATPRTR.
		332 ;	
		333 ;	*****
		334 ;	
0253	7E	335	LOADB: MOV A, M ; ASSUME HL IS LOADED WITH DATA2 ADDRESS.
0254	E60F	336	ANI 0FH
0256	320000	E 337	STA TDSTOR
0259	7E	338	MOV A, M
025A	E60F	339	ANI 0FH
025C	320100	E 340	STA TDSTOR+1
025F	7E	341	MOV A, M
0260	E60F	342	ANI 0FH
0262	320200	E 343	STA TDSTOR+2
0265	E5	344	PUSH H ; SAVE CONTENTS OF HL.
0266	CD8F02	C 345	CALL DATPRT
0269	E1	346	POP H
026A	C9	347	RET
		348 ;	
		349 ;	*****
		350 ;	
		351 ;	*** LOADA ***
		352 ;	
		353 ;	READS 3 DISPLAY LOCATIONS, BUT THIS TIME FROM SECTION B AND MUST
		354 ;	ROTATE RESULTS 4 TIMES TO GET H0M IN LOW POSITION.
		355 ;	
		356 ;	*****
		357 ;	
026B	7E	358	LOADA: MOV A, M ; AGAIN ASSUMES HL CONTAINS DATA2 ADDRESS.
026C	1F	359	RAR
026D	1F	360	RAR
026E	1F	361	RAR
026F	1F	362	RAR
0270	E60F	363	ANI 0FH ; MASK OUT HI-ORDER NIBBLE.
0272	320000	E 364	STA TDSTOR
0275	7E	365	MOV A, M
0276	1F	366	RAR
0277	1F	367	RAR
0278	1F	368	RAR
0279	1F	369	RAR
027A	E60F	370	ANI 0FH
027C	320100	E 371	STA TDSTOR+1
027F	7E	372	MOV A, M
0280	1F	373	RAR
0281	1F	374	RAR
0282	1F	375	RAR
0283	1F	376	RAR
0284	E60F	377	ANI 0FH

LOC	OBJ	LINE	SOURCE STATEMENT
0206	320200	E 378	STA TDSTOR+2
0209	E5	379	PUSH H
020A	CD8F02	C 380	CALL DATPRT
020D	E1	381	POP H
020E	C9	382	RET
		383 ;	
		384 ;	*****
		385 ;	
		386 ;	*** DATPRT ***
		387 ;	
		388 ;	PASSED: STEP, TDSTOR,+1,+2 ALL ARE IN BCD.
		389 ;	STEP IS INCREMENTED BY 2 WHEN EXITED.
		390 ;	
		391 ;	PRINTS: " STEP (STEP) = (TDSTOR).(+1)(+2) "
		392 ;	
		393 ;	*****
		394 ;	
020F	217F03	C 395	DATPRT: LXI H, STEPM ; SEND XXXXSTEPXX TO PRINTER.
0292	0E0B	396	MVI C, 0BH ; LOAD TEXT COUNTER.
0294	CD0000	E 397	CALL PWRITE
0297	3A0200	D 398	LDA STEP ; STEP IS IN FORMAT: MON = 10'S, LON = 1'S.
029A	47	399	MOV B, A
029B	E60F	400	ANI 0FH ; MASK TO GET ONE'S PLACE.
029D	FE0B	401	CPI 0BH
029F	C2A902	C 402	JNZ DONOT ; IF NO ZERO, THEN 1'S = 9 OR LESS.
02A2	3E10	403	MVI A, 10H ; SET A = 00010000H
02A4	00	404	ADD B ; INCREMENT 10'S PLACE
02A5	E6F0	405	ANI 0F0H ; RESET ONE'S PLACE TO ZERO.
02A7	3C	406	INR A ; ADD 1 TO GET FROM 10D TO 11D.
02A8	47	407	MOV B, A
02A9	78	408	DONOT: MOV A, B
02AA	320200	D 409	STA STEP ; UPDATE STEP.
02AD	1F	410	RAR
02AE	1F	411	RAR
02AF	1F	412	RAR
02B0	1F	413	RAR ; ROTATE TO GET 10'S PLACE INTO LON.
02B1	E60F	414	ANI 0FH ; SET RON=0.
02B3	C430	415	ADI 30H ; CONVERT TO ASCII.
02B5	CD0000	E 416	CALL PSEND ; SEND TO PRINTER.
02B8	70	417	MOV A, B
02B9	E60F	418	ANI 0FH
02BB	C430	419	ADI 30H
02BD	CD0000	E 420	CALL PSEND
02C0	3E20	421	MVI A, 20H ; SEND BLANK TO PRINTER.
02C2	CD0000	E 422	CALL PSEND
02C5	3E3D	423	MVI A, 3DH ; SEND EQUAL SIGN (=).
02C7	CD0000	E 424	CALL PSEND
02CA	3E20	425	MVI A, 20H ; SEND ONE MORE SPACE.
02CC	CD0000	E 426	CALL PSEND
02CF	3A0000	E 427	LDA TDSTOR ; GET MOST SIGNIFICANT DENSITY DIGIT.
02D2	C430	428	ADI 30H ; CONVERT TO ASCII.
02D4	CD0000	E 429	CALL PSEND
02D7	3E2E	430	MVI A, 2EH ; SEND " " TO PRINTER BUFFER.
02D9	CD0000	E 431	CALL PSEND

LOC	OBJ	LINE	SOURCE STATEMENT
02DC	3A0100	E 432	LDA TDSTOR+1
02DF	C630	433	ADI 30H
02E1	CD0000	E 434	CALL PSEND
02E4	3A0200	E 435	LDA TDSTOR+2
02E7	C630	436	ADI 30H
02E9	CD0000	E 437	CALL PSEND
02EC	3E17	438	MVI A, 17H
02EE	CD0000	E 439	CALL PSEND ; NOW PRINT BUFFER CONTENTS.
02F1	210200	D 440	LXI H, STEP
02F4	3E02	441	MVI A, 02H
02F6	06	442	ADD M
02F7	77	443	MOV M, A
02F8	C9	444	RET
		445 ;	
		446 ;	A SIMPLE SUBROUTINE FOR THOSE OF YOU WHO ARE RUBBER PEOPLE.
		447 ;	
02F9	3E00	448	RESET: MVI A, 00H ; RESET MANFLC TO OFF.
02FR	320000	E 449	STA MANFLC
02FE	212C00	450	LXI H, 002CH ; RESET DRAM POINTER.
0301	220000	E 451	SHLD MIXDRAM
0304	3E93	452	MVI A, 93H
0306	320000	E 453	STA POS1
0309	3E04	454	MVI A, 04H
030R	320000	E 455	STA POS1
030E	3E64	456	MVI A, 64H
0310	320000	E 457	STA RPOS1
0313	3EDF	458	MVI A, 0DFH
0315	320000	E 459	STA CLEARA
0318	320000	E 460	STA CLEARR
031R	3E16	461	MVI A, 16H ; RESET PRINTER TO 0 EXTRA DOT FEED.
031D	CD0000	E 462	CALL PSEND
0320	3E00	463	MVI A, 00H
0322	CD0000	E 464	CALL PSEND
0325	C9	465	RET
		466 ;	
		467 ;	*****
		468 ;	
		469 ;	*** REVERT ***
		470 ;	
		471 ;	REVERT RECEIVES CONTROL DISPLAY CODE IN A REGISTER AND CONVERTS
		472 ;	IT TO THE EQUIVALENT RCD CODE AND RETURNS IT IN C REGISTER.
		473 ;	
		474 ;	*****
		475 ;	
0326	FE0C	476	REVERT: CPI 0CH ; DISPLAY CODE FOR ZERO.
0328	0E00	477	MVI C, 00H ; 00H
032A	C0	478	RZ
032R	FE9F	479	CPI 9FH ; 01
032D	0E01	480	MVI C, 01H ; 01
032F	C0	481	RZ
0330	FE4A	482	CPI 4AH ; 02
0332	0E02	483	MVI C, 02H ; 02
0334	C0	484	RZ
0335	FE0R	485	CPI 0BH ; 03

174

LOC	OBJ	LINE	SOURCE STATEMENT
0337	0E03	486	MVI C, 03H ; 03
0339	C0	487	R2
033A	FE99	488	CPI 99H ; 04
033C	0E04	489	MVI C, 04H ; 04
033E	C0	490	R2
033F	FE29	491	CPI 29H ; 05
0341	0E05	492	MVI C, 05H ; 05
0343	C0	493	R2
0344	FE28	494	CPI 28H ; 06
0346	0E06	495	MVI C, 06H ; 06
0348	C0	496	R2
0349	FE0F	497	CPI 0FH ; 07
034B	0E07	498	MVI C, 07H ; 07
034D	C0	499	R2
034E	FE08	500	CPI 08H ; 08
0350	0E08	501	MVI C, 08H ; 08
0352	C0	502	R2
0353	FE09	503	CPI 09H ; 09
0355	0E09	504	MVI C, 09H ; 09
0357	C0	505	R2
0358	C9	506	RET ; NO MATCH WILL RETURN ANY WAY.
		507	*****
		508	;
		509	** CLEAN **
		510	;
		511	CLEAN ALLOWS MANUAL USE OF FTS SYSTEM SO USER MAY CHANGE SOLUTIONS
		512	WITHOUT HAVING TO RE-CALIBRATE DEVICE. DEVICE WILL REMAIN IN MANUAL
		513	MODE UNTIL ENTER KEY IS PRESSED.
		514	;
		515	*****
		516	;
0359	210000	E 517	CLEAN: LXI R, MANMOD ; PROMPT "MANUAL MODE?"
035C	060D	518	MVI B, 0DB
035E	CD0000	E 519	CALL DWRITE
0361	CD0000	E 520	CALL KDSPLY2 ; ENTER 0 AND RETURN TO WARM START.
0364	3A0000	E 521	LDA SUM2 ; ANYTHING ELSE PUTS YOU INTO MANUAL.
0367	FE00	522	CPI 00R
0369	CA0000	E 523	J2 START
036C	3EFE	524	MVI A, 0FEH ; SET FOR MANUAL CONTROL.
036E	0300	E 525	OUT LOW RPA2
0370	C35903	C 526	JMP CLEAN
		527	;
		528	DISPLAY MESSAGES.
		529	;
0373	0E	530	NMDATM: DB 0EH ; N
0374	0F	531	DB 0FH ; O
0375	20	532	DB 20H ; BLANK
0376	0D	533	DB 0DH ; M
0377	0F	534	DB 0FH ; O
0378	12	535	DB 12H ; R
0379	05	536	DB 05H ; E
037A	20	537	DB 20H ; BLANK
037B	04	538	DB 04H ; D
037C	01	539	DB 01H ; A

175

LOC	OBJ	LINE	SOURCE STATEMENT
037D	14	540	DB 14H ; T
037E	01	541	DB 01H ; A
037F	20	542	STEPM: DB 20H ; SPACE
0380	20	543	DB 20H ; SPACE
0381	20	544	DB 20H ; SPACE
0382	20	545	DB 20H ; SPACE
0383	20	546	DB 20H ; SPACE
0384	53	547	DB 53H ; S
0385	74	548	DB 74H ; T
0386	65	549	DB 65H ; E
0387	70	550	DB 70H ; P
0388	20	551	DB 20H ; SPACE
0389	50	552	PRONUM: DB 50H ; P
038A	52	553	DB 52H ; R
038B	4F	554	DB 4FH ; O
038C	43	555	DB 43H ; C
038D	45	556	DB 45H ; E
038E	53	557	DB 53H ; S
038F	53	558	DB 53H ; S
0390	20	559	DB 20H ; SPACE
0391	20	560	TPAST: DB 20H ; SPACE
0392	20	561	DB 20H ; SPACE
0393	20	562	DB 20H ; SPACE
0394	54	563	DB 54H ; T
0395	69	564	DB 69H ; I
0396	6D	565	DB 6DH ; M
0397	65	566	DB 65H ; E
0398	20	567	DB 20H ; SPACE
0399	45	568	DB 45H ; E
039A	6C	569	DB 6CH ; L
039B	61	570	DB 61H ; A
039C	73	571	DB 73H ; S
039D	70	572	DB 70H ; P
039E	65	573	DB 65H ; E
039F	64	574	DB 64H ; D
03A0	20	575	DB 20H ; SPACE
03A1	3A	576	DB 3AH ; :
03A2	20	577	DB 20H ; SPACE
		578 ;	
		579 DSEC	
		580 ;	
0000		581	LITEMP: DS 02H ; 2 BYTES OF RAM FOR LOWER INTERVAL STORAGE IN NEXTDM.
0002		582	STEP: DS 01H ; 1 BYTE USED AS STEP COUNTER IN PRINT ROUTINE.
		583 ;	
		584	END

PUBLIC SYMBOLS

UPORDN C 0000

EXTERNAL SYMBOLS

CLEARA E 0000	CLEARB E 0000	COMM1 E 0000	COMM2 E 0000	DATA1 E 0000	DATA2 E 0000	DEB E 0000
DFIND E 0000	DRAM E 0000	DWRITE E 0000	FIFOGR E 0000	INVERT E 0000	IRBRAM E 0000	IRCNT E 0000
ITEMP E 0000	KDSPY2 E 0000	MANFLG E 0000	MANMOD E 0000	MIN E 0000	MINCOM E 0000	NCNT E 0000

NXDRAM E 0000	OFFSET E 0000	POS1 E 0000	POS2 E 0000	PSEND E 0000	PWRITE E 0000	RAMPT E 0000
RPA2 E 0000	RPOS1 E 0000	SHOW E 0000	START E 0000	SUM2 E 0000	TDSTOR E 0000	TIMSUM E 0000

USER SYMBOLS

CLEAN C 0359	CLEARA E 0000	CLEARB E 0000	COMM1 E 0000	COMM2 E 0000	DATA1 E 0000	DATA2 E 0000
DATPRT C 028F	DEB E 0000	DECKS C 011A	DFIND E 0000	DONE C 00DF	DOMOT C 02A9	DRAM E 0000
PWRITE E 0000	FIFOGR E 0000	FIRST C 0202	INKN C 0003	INVERT E 0000	IRBRAM E 0000	IRCNT E 0000
ITEMP E 0000	KDSPY2 E 0000	KLOOK C 0006	LITEMP D 0000	LOADA C 026B	LOADB C 0253	MANFLC E 0000
HAMMOD E 0000	MIN E 0000	MINCON E 0000	MCNT E 0000	MEKEY C 002D	MEQ0 C 0108	NEITDN C 00C0
NEXTUP C 004A	NMDATA C 00B5	NMDATM C 0373	NOBARO C 00D8	NOTEQ C 005B	NXDRAM E 0000	OFFSET E 0000
MT C 00C8	OUTIT C 0141	POS1 E 0000	POS2 E 0000	PRINT C 0152	PROMUM C 0389	PSEND E 0000
PWRITE E 0000	RAMPT E 0000	RESET C 02F9	REST C 0239	REVERT C 0326	RPA2 E 0000	RPOS1 E 0000
SECB C 0214	SHOW E 0000	START E 0000	STEP D 0002	STEPM C 037F	SUBBS C 0132	SUM2 E 0000
TDSTOR E 0000	TIMSUM E 0000	TPAST C 0391	UPORDN C 0000			

ASSEMBLY COMPLETE, NO ERRORS

Appendix 3: Raw Data

The first section of data that appears here is for EK Commercial Film, type 6127 developed in D-76. The raw data is reproduced directly from the printer tapes that are output by the IR densitometer. After this first set was collected, the electronics that control the printer broke down and the unit had to be returned to the manufacturer for repair. The remaining data from the Fine Grain Release Positive experiments was reproduced by hand from the 33-digit LED density display. The IR densitometer will still work without the printer, but it certainly saves a lot of work for the user.

At the time of this writing, the control electronics still had not yet been returned from the NCR Corporation.

Process 1 = 05 minutes

178

Run 1 12/10/81

Interval 01 = 01 minutes
Rate 01 = 60 scans/minute

Interval 02 = 04 minutes
Rate 02 = 10 scans/minute

Process 2 = 01 minutes

Process 3 = 00 minutes

PROCESS 1 Time Elapsed 00 min 50 sec

Step 01 = 0.19

Step 03 = 0.19

Step 05 = 0.19

Step 07 = 0.22

Step 09 = 0.23

Step 11 = 0.28

Step 13 = 0.33

Step 15 = 0.41

Step 17 = 0.51

Step 19 = 0.64

Step 21 = 0.75

PROCESS 1 Time Elapsed 00 min 19 sec

Step 01 = 0.21

Step 03 = 0.21

Step 05 = 0.20

Step 07 = 0.22

Step 09 = 0.22

Step 11 = 0.22

Step 13 = 0.22

Step 15 = 0.22

Step 17 = 0.21

Step 19 = 0.23

Step 21 = 0.22

PROCESS 1 Time Elapsed 01 min 30 sec

Step 01 = 0.19

Step 03 = 0.19

Step 05 = 0.22

Step 07 = 0.30

Step 09 = 0.40

Step 11 = 0.53

Step 13 = 0.68

Step 15 = 0.86

Step 17 = 1.06

Step 19 = 1.31

Step 21 = 1.51

PROCESS 1 Time Elapsed 00 min 23 sec

Step 01 = 0.20

Step 03 = 0.20

Step 05 = 0.20

Step 07 = 0.21

Step 09 = 0.21

Step 11 = 0.22

Step 13 = 0.21

Step 15 = 0.22

Step 17 = 0.21

Step 19 = 0.23

Step 21 = 0.23

PROCESS 1 Time Elapsed 02 min 00 sec

Step 01 = 0.19

Step 03 = 0.21

Step 05 = 0.26

Step 07 = 0.39

Step 09 = 0.53

Step 11 = 0.71

Step 13 = 0.90

Step 15 = 1.13

Step 17 = 1.38

Step 19 = 1.68

Step 21 = 1.96

PROCESS 1 Time Elapsed 02 min 30 sec

Step 01 = 0.19
Step 03 = 0.22
Step 05 = 0.30
Step 07 = 0.46
Step 09 = 0.64
Step 11 = 0.86
Step 13 = 1.08
Step 15 = 1.33
Step 17 = 1.62
Step 19 = 1.98
Step 21 = 2.38

PROCESS 1 Time Elapsed 04 min 00 sec

Step 01 = 0.19
Step 03 = 0.24
Step 05 = 0.24
Step 07 = 0.64
Step 09 = 0.90
Step 11 = 1.16
Step 13 = 1.44
Step 15 = 1.72
Step 17 = 2.09
Step 19 = 2.66
Step 21 = 3.76

PROCESS 1 Time Elapsed 03 min 00 sec

Step 01 = 0.19
Step 03 = 0.23
Step 05 = 0.33
Step 07 = 0.53
Step 09 = 0.74
Step 11 = 0.98
Step 13 = 1.22
Step 15 = 1.50
Step 17 = 1.81
Step 19 = 2.22
Step 21 = 3.00

PROCESS 1 Time Elapsed 05 min 00 sec

Step 01 = 0.21
Step 03 = 0.26
Step 05 = 0.41
Step 07 = 0.69
Step 09 = 0.99
Step 11 = 1.28
Step 13 = 1.58
Step 15 = 1.88
Step 17 = 2.31
Step 19 = 3.32
Step 21 = 3.76

Process 1 = 05 minutes

180

Run 2 12/10/81

Interval 01 = 01 minutes
Rate 01 = 60 scans/minute

Interval 02 = 04 minutes
Rate 02 = 10 scans/minute

Process 2 = 01 minutes

Process 3 = 00 minutes

PROCESS 1 Time Elapsed 00 min 50 sec

Step 01 = 0.18

Step 03 = 0.18

Step 05 = 0.19

Step 07 = 0.20

Step 09 = 0.22

Step 11 = 0.26

Step 13 = 0.33

Step 15 = 0.41

Step 17 = 0.52

Step 19 = 0.67

Step 21 = 0.77

PROCESS 1 Time Elapsed 00 min 19 sec

Step 01 = 0.19

Step 03 = 0.19

Step 05 = 0.20

Step 07 = 0.21

Step 09 = 0.20

Step 11 = 0.20

Step 13 = 0.20

Step 15 = 0.20

Step 17 = 0.19

Step 19 = 0.21

Step 21 = 0.21

PROCESS 1 Time Elapsed 01 min 30 sec

Step 01 = 0.18

Step 03 = 0.19

Step 05 = 0.22

Step 07 = 0.29

Step 09 = 0.39

Step 11 = 0.53

Step 13 = 0.68

Step 15 = 0.86

Step 17 = 1.07

Step 19 = 1.33

Step 21 = 1.53

PROCESS 1 Time Elapsed 00 min 23 sec

Step 01 = 0.19

Step 03 = 0.19

Step 05 = 0.19

Step 07 = 0.20

Step 09 = 0.20

Step 11 = 0.20

Step 13 = 0.20

Step 15 = 0.20

Step 17 = 0.19

Step 19 = 0.23

Step 21 = 0.23

PROCESS 1 Time Elapsed 02 min 00 sec

Step 01 = 0.18

Step 03 = 0.19

Step 05 = 0.26

Step 07 = 0.37

Step 09 = 0.52

Step 11 = 0.70

Step 13 = 0.89

Step 15 = 1.13

Step 17 = 1.38

Step 19 = 1.69

Step 21 = 1.97

PROCESS 1 Time Elapsed 02 min 30 sec

Step 01 = 0.18

Step 03 = 0.21

Step 05 = 0.30

Step 07 = 0.45

Step 09 = 0.63

Step 11 = 0.85

Step 13 = 1.08

Step 15 = 1.33

Step 17 = 1.62

Step 19 = 1.99

Step 21 = 2.39

PROCESS 1 Time Elapsed 04 min 00 sec

Step 01 = 0.19

Step 03 = 0.24

Step 05 = 0.39

Step 07 = 0.62

Step 09 = 0.88

Step 11 = 1.16

Step 13 = 1.43

Step 15 = 1.72

Step 17 = 2.09

Step 19 = 2.67

Step 21 = 3.76

PROCESS 1 Time Elapsed 03 min 00 sec

Step 01 = 0.18

Step 03 = 0.22

Step 05 = 0.33

Step 07 = 0.51

Step 09 = 0.73

Step 11 = 0.97

Step 13 = 1.21

Step 15 = 1.49

Step 17 = 1.81

Step 19 = 2.24

Step 21 = 2.98

PROCESS 1 Time Elapsed 05 min 00 sec

Step 01 = 0.20

Step 03 = 0.26

Step 05 = 0.41

Step 07 = 0.68

Step 09 = 0.98

Step 11 = 1.27

Step 13 = 1.55

Step 15 = 1.87

Step 17 = 2.30

Step 19 = 3.75

Step 21 = 3.76

Process 1 = 05 minutes

182

Run 3 12/10/81

Interval 01 = 01 minutes

Rate 01 = 60 scans/minute

Interval 02 = 04 minutes

Rate 02 = 10 scans/minute

Process 2 = 01 minutes

Process 3 = 00 minutes

PROCESS 1 Time Elapsed 00 min 50 sec

Step 01 = 0.19

Step 03 = 0.19

Step 05 = 0.19

Step 07 = 0.21

Step 09 = 0.24

Step 11 = 0.29

Step 13 = 0.35

Step 15 = 0.44

Step 17 = 0.55

Step 19 = 0.69

Step 21 = 0.81

PROCESS 1 Time Elapsed 00 min 19 sec

Step 01 = 0.21

Step 03 = 0.21

Step 05 = 0.21

Step 07 = 0.21

Step 09 = 0.22

Step 11 = 0.22

Step 13 = 0.22

Step 15 = 0.22

Step 17 = 0.22

Step 19 = 0.23

Step 21 = 0.23

PROCESS 1 Time Elapsed 01 min 30 sec

Step 01 = 0.19

Step 03 = 0.19

Step 05 = 0.22

Step 07 = 0.30

Step 09 = 0.41

Step 11 = 0.55

Step 13 = 0.71

Step 15 = 0.90

Step 17 = 1.11

Step 19 = 1.35

Step 21 = 1.59

PROCESS 1 Time Elapsed 00 min 23 sec

Step 01 = 0.20

Step 03 = 0.21

Step 05 = 0.20

Step 07 = 0.21

Step 09 = 0.21

Step 11 = 0.22

Step 13 = 0.21

Step 15 = 0.22

Step 17 = 0.22

Step 19 = 0.24

Step 21 = 0.26

PROCESS 1 Time Elapsed 02 min 00 sec

Step 01 = 0.19

Step 03 = 0.21

Step 05 = 0.26

Step 07 = 0.38

Step 09 = 0.54

Step 11 = 0.73

Step 13 = 0.93

Step 15 = 1.16

Step 17 = 1.42

Step 19 = 1.73

Step 21 = 2.03

PROCESS 1 Time Elapsed 02 min 30 sec

Step 01 = 0.19
Step 03 = 0.22
Step 05 = 0.31
Step 07 = 0.46
Step 09 = 0.65
Step 11 = 0.87
Step 13 = 1.11
Step 15 = 1.37
Step 17 = 1.65
Step 19 = 2.02
Step 21 = 2.46

PROCESS 1 Time Elapsed 04 min 00 sec

Step 01 = 0.19
Step 03 = 0.24
Step 05 = 0.39
Step 07 = 0.64
Step 09 = 0.90
Step 11 = 1.18
Step 13 = 1.45
Step 15 = 1.74
Step 17 = 2.12
Step 19 = 2.69
Step 21 = 3.76

PROCESS 1 Time Elapsed 03 min 00 sec

Step 01 = 0.19
Step 03 = 0.23
Step 05 = 0.33
Step 07 = 0.53
Step 09 = 0.75
Step 11 = 1.00
Step 13 = 1.24
Step 15 = 1.52
Step 17 = 1.84
Step 19 = 2.27
Step 21 = 3.23

PROCESS 1 Time Elapsed 05 min 00 sec

Step 01 = 0.21
Step 03 = 0.26
Step 05 = 0.41
Step 07 = 0.69
Step 09 = 0.99
Step 11 = 1.29
Step 13 = 1.58
Step 15 = 1.89
Step 17 = 2.33
Step 19 = 3.49
Step 21 = 3.76

Fine Grain Release Positive in D-19.

1/29/82

Temp: 70°F

Filter: 1.20 N.O.

Sens.ometer current = .82 amp

Calibration Standard: .70

Time:

	<u>7 sec</u>	<u>14</u>	<u>30</u>	<u>45</u>	<u>60</u>	<u>1:30</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>
①										
	.14	.16	.17	.17	.17	.17	.19	.23	.27	.32
	.14	.17	.17	.17	.18	.18	.20	.24	.28	.33
	.14	.17	.18	.18	.18	.18	.20	.24	.28	.33
	.16	.19	.19	.19	.19	.20	.22	.26	.31	.36
	.16	.19	.20	.20	.20	.22	.24	.28	.34	.40
	.16	.19	.20	.22	.24	.26	.29	.37	.44	.51
	.14	.18	.24	.28	.32	.39	.46	.58	.72	.84
	.16	.22	.38	.48	.57	.74	.87	1.12	1.32	1.49
	.19	.32	.64	.85	1.00	1.25	1.47	1.81	2.08	2.27
	.18	.41	.91	1.27	1.50	1.87	2.15	2.52	2.74	2.94
	.20	.53	1.17	1.65	1.93	2.35	2.60	2.94	3.14	3.08
②										
	.13	.16	.17	.17	.16	.17	.18	.21	.25	.30
	.13	.16	.16	.16	.16	.16	.17	.21	.25	.30
	.13	.17	.18	.18	.17	.18	.19	.23	.27	.31
	.14	.19	.19	.19	.19	.19	.21	.25	.29	.34
	.14	.17	.18	.19	.19	.20	.22	.26	.31	.36
	.17	.22	.23	.24	.25	.28	.31	.37	.44	.51
	.14	.19	.25	.28	.32	.39	.45	.58	.70	.80
	.17	.26	.39	.49	.58	.72	.86	1.07	1.26	1.43
	.18	.33	.59	.80	.95	1.21	1.42	1.76	2.03	2.24
	.18	.44	.86	1.22	1.44	1.80	2.09	2.50	2.75	2.90
	.19	.53	1.10	1.58	1.87	2.28	2.56	2.85	2.98	3.06

F.G.R.P. in D-19 (Cont)

1/29

Time:

③

<u>7sec</u>	<u>14</u>	<u>30</u>	<u>45</u>	<u>60</u>	<u>1:30</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>
.13	.16	.17	.17	.17	.17	.19	.22	.27	.33
.13	.16	.18	.17	.17	.17	.19	.24	.28	.33
.13	.16	.17	.17	.17	.17	.19	.23	.28	.33
.14	.18	.19	.19	.19	.20	.22	.26	.31	.36
.14	.18	.19	.20	.20	.21	.23	.28	.34	.40
.15	.19	.21	.22	.24	.26	.30	.37	.44	.51
.16	.22	.28	.32	.36	.43	.50	.63	.76	.87
.17	.25	.41	.51	.60	.75	.89	1.13	1.32	1.50
.14	.29	.60	.82	.98	1.26	1.48	1.83	2.09	2.31
.16	.41	.89	1.26	1.49	1.88	2.17	2.55	2.82	2.93
.19	.53	1.17	1.64	1.93	2.36	2.61	2.97	3.15	3.08

④ Developer and rinse water not replaced for this last run.

.16	.19	.20	.19	.20	.20	.22	.26	.31	.36
.14	.17	.17	.18	.18	.19	.20	.25	.30	.36
.16	.18	.19	.20	.20	.20	.22	.26	.32	.38
.16	.19	.20	.20	.20	.21	.23	.28	.33	.38
.16	.18	.19	.20	.20	.22	.23	.30	.36	.43
.17	.19	.22	.24	.25	.28	.32	.40	.48	.56
.18	.22	.29	.34	.38	.46	.54	.68	.82	.94
.19	.27	.44	.56	.66	.82	.96	1.23	1.43	1.61
.17	.33	.68	.91	1.08	1.35	1.59	1.94	2.20	2.40
.20	.48	1.05	1.39	1.64	2.02	2.29	2.69	2.92	3.07
.24	.62	1.37	1.81	2.08	2.50	2.75	3.03	3.09	3.12

Fine Grain Release Positive in D-76

1/29/82

Temp: 70°F

Filter: .78 N.O.

Sensitometer Current = .82 amp

Calibration Standard = .70

Time:

①

<u>20sec</u>	<u>30</u>	<u>45</u>	<u>1:00</u>	<u>1:30</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>
.19	.18	.18	.17	.16	.16	.13	.12	.10
.20	.19	.19	.18	.17	.16	.14	.13	.10
.16	.16	.16	.16	.15	.14	.13	.12	.10
.18	.18	.17	.18	.18	.18	.18	.17	.16
.17	.17	.20	.22	.25	.28	.31	.33	.34
.19	.21	.26	.32	.39	.45	.54	.60	.65
.21	.25	.35	.45	.59	.70	.85	.96	1.05
.20	.28	.44	.60	.81	.96	1.18	1.35	1.48
.25	.37	.60	.79	1.10	1.31	1.59	1.80	1.97
.29	.47	.74	.98	1.37	1.63	2.00	2.25	2.42
.29	.50	.82	1.11	1.54	1.88	2.30	2.58	2.72

②

.17	.17	.16	.16	.15	.13	.12	.10	.08
.21	.20	.19	.18	.17	.16	.14	.12	.10
.19	.19	.17	.17	.16	.16	.14	.12	.10
.21	.20	.20	.20	.20	.19	.19	.18	.16
.17	.17	.19	.20	.23	.26	.30	.32	.33
.18	.19	.25	.31	.38	.44	.53	.59	.64
.19	.23	.35	.44	.59	.69	.85	.97	1.07
.24	.32	.48	.63	.83	.98	1.20	1.37	1.50
.27	.39	.62	.81	1.12	1.32	1.60	1.83	1.99
.28	.44	.73	.97	1.37	1.64	2.01	2.28	2.47
.33	.54	.85	1.14	1.57	1.90	2.33	2.59	2.74

1/29

Time:

③

<u>20sec</u>	<u>30</u>	<u>45</u>	<u>1:00</u>	<u>1:30</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>
.17	.17	.16	.15	.14	.13	.11	.10	.08
.16	.16	.15	.14	.13	.13	.12	.10	.08
.16	.16	.15	.14	.13	.13	.12	.10	.09
.18	.17	.16	.16	.16	.17	.16	.16	.15
.16	.16	.19	.21	.24	.26	.30	.31	.33
.21	.22	.27	.32	.40	.46	.54	.60	.65
.21	.25	.36	.46	.60	.70	.86	.97	1.07
.22	.29	.46	.60	.82	.97	1.19	1.37	1.50
.25	.38	.59	.79	1.10	1.31	1.59	1.82	1.98
.29	.47	.75	.98	1.37	1.65	2.02	2.27	2.46
.30	.51	.83	1.12	1.56	1.90	2.33	2.60	2.80

④

Developer for this and rinse water not replaced last run.

.19	.18	.17	.16	.15	.14	.12	.10	.09
.19	.19	.17	.17	.16	.15	.13	.12	.10
.17	.16	.16	.15	.14	.13	.13	.11	.10
.20	.20	.19	.19	.19	.19	.19	.18	.17
.19	.19	.20	.23	.26	.29	.32	.34	.36
.21	.23	.29	.35	.42	.48	.57	.63	.68
.23	.28	.40	.50	.64	.73	.89	1.02	1.11
.25	.35	.52	.67	.88	1.03	1.25	1.42	1.54
.26	.41	.65	.86	1.16	1.36	1.64	1.87	2.02
.30	.50	.79	1.04	1.43	1.70	2.07	2.31	2.50
.33	.58	.91	1.19	1.64	1.97	2.39	2.63	2.85

12/11/81

Dry, fixed out white light density readings
Taken on TD-504.

Eastman Kodak Commercial Film developed in
D-76 @ 70°F. all readings after 5 minutes
of development.

Step \ Run	1	2	3
1	.08	.09	.09
2	.12	.12	.12
3	.23	.23	.23
4	.43	.43	.44
5	.71	.71	.71
6	1.02	1.02	1.03
7	1.34	1.33	1.35
8	1.68	1.67	1.68
9	2.02	2.02	2.03
10	2.40	2.42	2.42
11	2.72	2.73	2.73

1/30/82

Dry, fixed out, white lite density
Readings taken on TD-504.

Fine Grain Release Positive in D-19:

Step/Run	1	2	3	4
1	.29	.26	.30	.29
2	.29	.27	.30	.30
3	.30	.28	.31	.31
4	.32	.29	.33	.33
5	.35	.33	.36	.37
6	.48	.44	.49	.52
7	.88	.80	.88	.94
8	1.61	1.51	1.61	1.72
9	2.54	2.41	2.55	2.66
10	3.36	3.24	3.36	3.46
11	3.94	3.80	3.91	4.01

Fine Grain Release Positive in D-76:

Step/Run	1	2	3	4
1	.04	.04	.05	.05
2	.04	.05	.05	.06
3	.06	.06	.07	.07
4	.12	.12	.12	.13
5	.29	.29	.29	.32
6	.61	.60	.62	.66
7	1.06	1.05	1.07	1.12
8	1.55	1.54	1.56	1.62
9	2.04	2.04	2.05	2.13
10	2.63	2.61	2.62	2.70
11	3.05	3.03	3.05	3.15

Appendix 4: Infrared Emitting Diode LED55C* Specifications:

Absolute Maximum Ratings (@ 25°C):

Reverse Voltage	V_R	3	volts
Forward Current	I_F	100	mA
Power Dissipation	P_T	170	mW

Optical Characteristics (@ 25°C):

Total Power Output ($I_F=100\text{mA}$)	P_O	5.4	mW
Peak Emmision Wavelength ($I_F=100\text{mA}$)		940	nm
Spectral Bandwidth 50%		60	nm
Rise Time 0-90% of output		300	nsec

* From Optoelectronics Manual, General Electric Company,
Electronics Park, Syracuse, N.Y. 13201.




Figure 22: Power Output vs. Input Current for LED55C




Figure 23: Typical Radiation Pattern for LED55C

Appendix 5: UDT-450 Specifications* :

<u>Parameter</u>	<u>Symbol</u>	<u>Typical Value</u>	<u>Units</u>
Responsivity	R	.4 @850 nm	amp/watt
Active Area	A	.05	cm ²
Active Dia.	D	.10	inch
Output Resistance	Z _{out}	100	ohms
slew Rate		1	volt/microsec
Unit Gain Bandwidth		1	MHz
Supply Voltage		+/-15	volts
Supply Current		3	mA
Offset Voltage Drift w/ Temperature		+/-25	microvolt/°C
Light Range	L	10 ⁻² -5x10 ⁻¹²	watts/cm ²
Feedback resistor			
Range	R _f	.5K-50M	ohms
Frequency response	f	dc-10 ⁶	Hz
Spectral Range		350-1100	nm

* UDT-450 Data Sheet, United Detector Technology, 2644 30th st.
Santa monica, CA. 90405.




Figure 24: Spectral Response of UDT-450




Figure 25: Output voltage of UDT-450 as a function of incident energy

Appendix 6: Model 757N Logarithmic Ratio Amplifier* :

Transfer Function:

voltage mode
$$e_0 = -K \log_{10}(e_1/e_2 \times R_2/R_1)$$

Accuracy:

log conformity $\pm 1\%$, max.
1nA to 10mA Relative to Input

Input Specs:

maximum current ± 10 mA
worst case offset ± 85 microvolts/ $^{\circ}$ C max.
voltage

Rise Time:

increasing input 250 microsecond max.
decreasing input 600 microsecond max.

Power Supply:

rated performance ± 15 V dc
operating $\pm (12 \text{ to } 18)$ V dc
current, quiescent ± 8 mA

Mechanical:

case size 1.5" x 1.5" x 0.4"
weight 21 grams

* Data sheet, Analog Devices, P.O. Box 280, Norwood, Mass. 02062.

Appendix 7: SDK-85 Specifications* :

Central Processor

CPU : 8085A.

Instruction Cycle : 1.3 microsecond.

T_{cy} : 330 ns.

Memory

ROM : 4K bytes using 2 8755s.

RAM : 512 bytes with 2 8155s.

Addressing : Expandable to 64K bytes by use of additional buffers and decoders.

Input/Output

Parallel : 38 lines (expandable to 76).

Serial : Through SID/SOD ports of 8085. Software generated baud rate.

Baud rate : 110

Interfaces

Bus : All signals TTL compatible.

Parallel I/O : All signals TTL compatible.

Serial : 20 mA current loop.

Interrupts

Three levels : (RST 7.5) - Time base input.

(RST 6.5) - TTL input.

INTR - TL input.

* SDK Users Manual, order number 9800451B, Intel Corporation, 3065 Bowers Ave., Santa Clara, CA. 95051.

DMA

Hold request : Jumper selectable. TTL compatible input.

Physical Characteristics

Width : 12.0 in.

Height : 10.0 in.

Depth : 0.50 in.


Weight : Approx. 12 oz.

Electrical Characteristics (DC Power Required)


V_{CC} : +5V +/- 5% @ 1.3 amp.

Environmental

Operating temperature : 0-55 °C



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ON XX DEC .02 XXX DEC .005 ANGLES 10°	MATERIAL	LOG AMPLIFIER
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE FOR MFG TOLERANCES INCLUDING ROUND PINCHES AND DRILLED HOLES. SEE SPECS	FINISH	
DRAWN Steve Cox	CHECKER	SIZE B
APPROVED SC 11/13/81	GROUP CODE	SCALE None
		SHEET 1
		OF 1



ITEM	DESCRIPTION
S1-S6	0m121CBA3AD 2 Pole Push Button Rcty
S7	0m121FBA3AD 2 Pole Push Button
L1-L8	0m121LA36 6.3 Volt 20/5

REVISIONS			
ZONE	REV	ECO NO.	DATE

FTS MANUAL OVERRIDE		SIZE C	DRAWING NO.	SCALE 1/16" = 1"	SHEET 1 OF 1
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES		XX DEC TOLERANCES UNLESS OTHERWISE SPECIFIED			
XX DEC TOLERANCES UNLESS OTHERWISE SPECIFIED		ANGLES 1°			
COMMERCIAL TOLERANCES GOVT FIN		TOLERANCES INCLUDING PUNCHING AND DRILLED HOLES			
SEE NOTES		DRAWN: S.C.O.X			
CHECKER		GROUP CDDI			
APPROVED: 1a/11/80		2			

DATE 10-10-80

4


3

2

1

REVISIONS			
ZONE	REV	ECO NO.	DATE
B4	1	-	9/15/81
C3	2	-	2/26/82

APPD
S. Cox



TO J4 ON
CPU BOARD

U1-AD571JD
U2,3-74S240

TO
TERMINAL BLOCKS

200


UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ON XX DEC XX DEC ANGLES ±.02 ±.005 ±.1°		MATERIAL
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE FOR MFG. TOLERANCES, INCLUDING ROUND PUNCHED AND DRILLED HOLES. SEE SPI-ES		FINISH
NEXT ASSY.	USED ON	CHECKER
DRAWN S. COX		GROUP CODE
APPROVED 12/10/80		SCALE A 100 X
SHEET 1 OF 1		DRAWING NO. B
A/D CONVERTER BOARD ASSEMBLY		SIZE B

4

3

2

1



201

U1-ANALOG DEVICES
AD571JD

U2,U3 SN74S240
LINE DRIVERS

CARD EDGE
NUMBERS

V Start Conversion

TO ANALOG
GND

Gnd

ANALOG
IN

BNC
CONNECTOR

A/D BOARD

SIZE B DRAWING NO


SCALE 1/000

SHEET 1 OF 1

REVISIONS *			
ZONE	REV	ECO NO	DATE
B3	1		3/24/81
			S Cox

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ON ANGLES XX DEC XX DEC ANGLES ± .02 ± .005 ± .1°		MATERIAL	
COMMERCIAL STOCK SIZE FOR MFG. TOLERANCES, INCLUDING ROUND PUNCHED AND DRILLED HOLES. SEE SPECS		FINISH	
NEXT ASSY USED ON		DRAWN Steve Cox	
		APPROVED 9/14/81	
		CHECKER	
		GROUP CODE	


REVISIONS				
ZONE	REV	ECO NO.	DATE	APPD
Part Description				
D13-D17		5080-7740		H.P.
D1-D12		TIL-305		T.I.
U ₁ -U ₁₂		74LS164N		
Q ₁ -Q ₂		DH3467C N		
J1		26-Pin Header		
All resistors are 1/4 watt				



202

D 13

TOP VIEW		12 CHARACTER DISPLAY ASSEMBLY	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.		MATERIAL	
XX DEC TOLERANCES ON ANGLES ± .02 ± .005 ± .1°		FINISH	
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE FOR MFG TOLERANCES INCLUDING ROUND PUNCHED AND DRILLED HOLES. SEE SPECS		CHECKER	
DRAWN S. COX		GROUP CODE	
APPROVED 12/9/80		SCALE APPROX	
NEXT ASSY.		SIZE B	
USED ON		DRAWING NO.	
		SHEET OF	




203

PART	DESCRIPTION
D1-D12	5x7 ALPHA DISPLAY
L1-L12	506-770 11-P 7 SEGMENT DISPLAY
U1-U12	74LS08 2-INPUT 4-OUTPUT SHIFT REGISTER
Q1-Q2	D 3467 QUAD PNP TRANSISTOR

12 CHARACTER DISPLAY	
MATERIAL	FINISH
CONTROL PANEL	UNLESS OTHERWISE SPECIFIED
UNIT ASSEMBLY	UNLESS OTHERWISE SPECIFIED
APPROVED	12/8/80
CHECKER	GROUP CODE
DRAWN	SCALE
SIZE	1
DATE	1
ZONE	1
REV	1
ED NO	1
DATE	1
APPRO	1

REVISIONS				
ZONE	REV	ECO NO.	DATE	APPD

U_1 to U_{12} are Hewlett Packard
HCPL-2730 dual Channel
Optocouplers, mounted 2
per 16 pin DIP socket



SDK-95

204

[illegible]

REVISIONS			
ZONE	REV	ECO NO.	DATE

Parts:


- R1 to R22 - 2.7 KΩ 1/4 W
- R23 to R44 - 3.3 KΩ 1/4 W
- U1 to U12 - H/P HCPL-2730
- J4 - 26 pin header

SDK-85 Expansion RAM
8155

Port:	Address:	Status:
PA0-7	29 H	Input
PB0-7	2A H	Input
PC0-5	2B H	Output
C/S	28 H	D.N.A.


205

+5VDC } From SDK-85
GND



UNLESS OTHERWISE SPECIFIED		MATERIAL	
DIMENSIONS ARE IN INCHES		FINISH	
TOLERANCES UNLESS OTHERWISE SPECIFIED		DRAWING NO	
FRACTIONS .XX DEC .XXX DEC .10		SCALE (None)	
COMMERCIAL TOLERANCES GOVERN		CHECKER	
PUNCHED AND DRILLED HOLES		GROUP CODE	
SEE NOTES		APPROVED 12/5/80	
DRAWN S C O X		1	
NEXT ASBY.		2	
USED ON		3	
A/D CONDUCT		4	
A/D CONDUCT		5	
A/D CONDUCT		6	
A/D CONDUCT		7	
A/D CONDUCT		8	
A/D CONDUCT		9	
A/D CONDUCT		10	
A/D CONDUCT		11	
A/D CONDUCT		12	
A/D CONDUCT		13	
A/D CONDUCT		14	
A/D CONDUCT		15	
A/D CONDUCT		16	
A/D CONDUCT		17	
A/D CONDUCT		18	
A/D CONDUCT		19	
A/D CONDUCT		20	
A/D CONDUCT		21	
A/D CONDUCT		22	
A/D CONDUCT		23	
A/D CONDUCT		24	
A/D CONDUCT		25	
A/D CONDUCT		26	
A/D CONDUCT		27	
A/D CONDUCT		28	
A/D CONDUCT		29	
A/D CONDUCT		30	
A/D CONDUCT		31	
A/D CONDUCT		32	
A/D CONDUCT		33	
A/D CONDUCT		34	
A/D CONDUCT		35	
A/D CONDUCT		36	
A/D CONDUCT		37	
A/D CONDUCT		38	
A/D CONDUCT		39	
A/D CONDUCT		40	
A/D CONDUCT		41	
A/D CONDUCT		42	
A/D CONDUCT		43	
A/D CONDUCT		44	
A/D CONDUCT		45	
A/D CONDUCT		46	
A/D CONDUCT		47	
A/D CONDUCT		48	
A/D CONDUCT		49	
A/D CONDUCT		50	

OPTO-ISOLATED I/O




ITEM	TYPE	DESCRIPTION	QTY
Q1-Q8	DM3467CN	QUAD PNP Transistor	8
U1 and U2	SN7448	B.C.D to 7 Segment Decoder	2
U3 and U4	SN7448	Hex Inverting Drivers	2
U5	SN7445	1 of 16 Decoder	1
D1-D33	E083-7730	HP Seven Segment Display	33
UA	SN7442	BCD to Decimal Decoder	1
J1		2 Pin Molex	1
J2		26 Pin Right Angle Header	1

NOTES:

1. DATA FROM THE 8279 IS SENT AS 8 BITS MADE UP OF 2 NYBBLES, A & B.
2. U1 DECODES NYBBLE A AND IN TURN DRIVES U2 AND Q1-Q4. DIGITS DI-DIG ARE THEREFORE WRITTEN TO THRU NYBBLE A.
3. SIMILARLY, D17-D32 ARE WRITTEN TO THRU NYBBLE B.
4. D33 IS HARD-WIRED TO DISPLAY ZERO.
5. RESISTORS ARE IN OHMS AND RATED $\frac{1}{4}$ WATT.

206


33 DIGIT DISPLAY ASSEMBLY	
DRAWING NO. C	
SCALE	SHEET OF 1
MATERIAL	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES XX DEC XXX DEC 1.00 1.000 ANGLE 1.0	FINISH
COMMERCIAL TOLERANCES GOVERN DIMENSIONS UNLESS OTHERWISE SPECIFIED. INCLUDE ROUND PUNCHED AND DRILLED HOLES. ELECTRIC	
DRAWN S. Cox	CHECKER
APPROVED R. G. ISO	GROUP CODE



33 DIGIT DENSITY DISPLAY SYSTEM	
SIZE	DRAWING NO.
C	C
SCALE	1 OF 2
APPROVED	
DRAWN	
CHECKER	
GROUP CODE	
MATERIAL	
XX DEC. XXX DEC. ANGLES	
TOLERANCES ON	
DIMENSIONS ARE IN INCHES UNLESS OTHERWISE SPECIFIED	
TOLERANCES ON MATERIAL STOCK SIZE FOR MFG	
PUNCHED HOLES SHALL BE 0.015 IN.	
HARDWIRED ZERO	
D.33	

U1-74LS154
U2,3-7414
U4,5-9664


TO 11 IREDS IN.
CHAMBER TOP

11 x 100 1/4 W
Resistors

208

TOP VIEW

[illegible]



209


U1-SN74154 1 OF 16 DECODER
U2,U3-SN7414 INVERTER
U4,U5-SN7414 TI LED DRIVERS
DI-11-LED55C GE IRED

GREY CABLE TO
CHARTER BOTTOM

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES		MATERIAL	
XX DEC	XX DEC	XX DEC	XX DEC
XX DEC	XX DEC	XX DEC	XX DEC
COMMERCIAL TOLERANCES GOVERN		FINISH	
MATERIAL STOCK SIZE FOR MFG		CHECKER	
TOLERANCES IN PARENTESIS		GROUP CODE	
PUNCHED AND DULLED		DRAWN	
SEE DETAIL		APPROVED	
9/11/81		Cox	
1		2	
1		3	
1		4	
1		5	
1		6	
1		7	
1		8	
1		9	
1		10	
1		11	
1		12	
1		13	
1		14	
1		15	
1		16	
1		17	
1		18	
1		19	
1		20	
1		21	
1		22	
1		23	
1		24	
1		25	
1		26	
1		27	
1		28	
1		29	
1		30	
1		31	
1		32	
1		33	
1		34	
1		35	
1		36	
1		37	
1		38	
1		39	
1		40	
1		41	
1		42	
1		43	
1		44	
1		45	
1		46	
1		47	
1		48	
1		49	
1		50	
1		51	
1		52	
1		53	
1		54	
1		55	
1		56	
1		57	
1		58	
1		59	
1		60	
1		61	
1		62	
1		63	
1		64	
1		65	
1		66	
1		67	
1		68	
1		69	
1		70	
1		71	
1		72	
1		73	
1		74	
1		75	
1		76	
1		77	
1		78	
1		79	
1		80	
1		81	
1		82	
1		83	
1		84	
1		85	
1		86	
1		87	
1		88	
1		89	
1		90	
1		91	
1		92	
1		93	
1		94	
1		95	
1		96	
1		97	
1		98	
1		99	
1		100	
1		101	
1		102	
1		103	
1		104	
1		105	
1		106	
1		107	
1		108	
1		109	
1		110	
1		111	
1		112	
1		113	
1		114	
1		115	
1		116	
1		117	
1		118	
1		119	
1		120	
1		121	
1		122	
1		123	
1		124	
1		125	
1		126	
1		127	
1		128	
1		129	
1		130	
1		131	
1		132	
1		133	
1		134	
1		135	
1		136	
1		137	
1		138	
1		139	
1		140	
1		141	
1		142	
1		143	
1		144	
1		145	
1		146	
1		147	
1		148	
1		149	
1		150	
1		151	
1		152	
1		153	
1		154	
1		155	
1		156	
1		157	
1		158	
1		159	
1		160	
1		161	
1		162	
1		163	
1		164	
1		165	
1		166	
1		167	
1		168	
1		169	
1		170	
1		171	
1		172	
1		173	
1		174	
1		175	
1		176	
1		177	
1		178	
1		179	
1		180	
1		181	
1		182	
1		183	
1		184	
1		185	
1		186	
1		187	
1		188	
1		189	
1		190	
1		191	
1		192	
1		193	
1		194	
1		195	
1		196	
1		197	
1		198	
1		199	
1		200	
1		201	
1		202	
1		203	
1		204	
1		205	
1		206	
1		207	
1		208	

ITEM	TYPE	DESCRIPTION	QTY
U ₁	HCPL-2730	HP Opto-Isolator	2
U ₂ -U ₄	7415240	Bus Driver, Inverting	3
P ₁		26 Pin Header	1
P ₂ -P ₃		16 Pin DIP Socket	2
P ₄		2 Pin molex	1


Resistors are Rated 1/4 watt



TOP VIEW

210

[illegible]




211

REVISIONS			
ZONE	REV	DATE	APPD


UNLESS OTHERWISE SPECIFIED		MATERIAL
DIMENSIONS ARE IN INCHES		
TOLERANCES ON		ANGLES
FRACTIONS		1°
DECIMALS		0.005
COMMERCIAL TOLERANCES GOVERN		FINISH
MATERIALS AND FINISHES		
TOLERANCES INCLUDING HOLES		
PUNCHED AND DRILLED HOLES		
SEE NOTES		
DRAWN	S. COX	CHECKER
APPROVED	12/6/89	GROUP CODE

DRIVER BOARD	
SIZE	C
DRAWING NO.	
SCALE	1/8" = 1"
SHEET	1

2-pin header



Section A-A



2 Required

214

		UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ON ANGLES XX DEC. XXX DEC ANGLES ± .02 ± .005 ± 1°		MATERIAL 316 Stainless Steel	FINISH LOCK	SIZE B	DRAWING NO. 4R	SCALE 2:1	SHEET OF
		COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE. FOR MFG TOLERANCES, INCLUDING ROUND PUNCHED AND DRILLED HOLES. SEE SPIES							
		DRAWN S. Cox		CHECKER		4/8/80	GROUP CODE		
		APPROVED							
NEXT ASSY.		Chamber	Front (1R)						
		USED ON							

1

2

3

4


1

2

3

4

M.H.S. 10-10-1980



215
2 Required

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. TOLERANCES ON XX DEC. $\pm .02$ XX DEC. $\pm .005$ ANGLES $\pm 1^\circ$		MATERIAL 316 Stainless Steel	
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE FOR MFG. TOLERANCES, INCLUDING ROUND PUNCHED AND DRILLED HOLES. SEE SPIES		FINISH	
DRAWN S. Cox		CHECKER	
APPROVED 4/10/80		GROUP CODE	
Chamber Front USED ON		HINGE	
NEXT ASSY.		SIZE B	
		DRAWING NO. 3R	
		SCALE 2:1	
		SHEET OF	


4

3

2

1

--#27(.144) Drill
C6K For 6-32 Socket Hd



216

REVISIONS			
ZONE	REV	ECO NO.	DATE

MATERIAL		316 Stainless Steel	
FINISH		CHAMBER MOUNT	
DRAWN		S. Cox	
APPROVED		4/17/80	
CHECKER		9R	
GROUP CODE		1	
SCALE		Full	
SHEET		OF	
SIZE		B	


4

3

2

1


#27(.144) Drill
csk For 6-32 Socket Hd.





217

NEXT ASSY.		Chamber Front		USED ON	
DRAWN		S. COX		APPROVED	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. XXDEC. TOLERANCES ON ANGLES +.02 XXXDEC. TOLERANCES ON ANGLES +.005 +.10		MATERIAL 316 Stainless Steel		Cover Plate	
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE. FOR MFG. TOLERANCES, INCLUDING ROUND PUNCHED AND DRILLED HOLES, SEE SPI.E5		FINISH		DRAWING NO. 7R	
CHECKER		GROUP CODE		SCALE Full	
SHEET		OF		A	

Section A-A




15 Required




UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES.		MATERIAL	
XXDEC. +.02	XXDEC. +.006	316 Stainless Steel	
ANGLES ± 10		FINISH	
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE. FOR MFG. TOLERANCES INCLUDING ROUND, PUNCHED AND DRILLED HOLES, SEE SPI-E6		CHECKER	
DRAWN S. Cox		GROUP CODE	
APPROVED 4/17/80		SIZE A	
NEXT ASSY.		DRAWING NO. 8R	
Chamber		SCALE 6%1	
Back		SHEET OF	
USED ON		lens Tube	

REV ECO NO. DATE APPD




1/4 x 28 Hex Head Bolt 3/8 inch long

REV		ECO NO.		DATE		APPD	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. XXDEC. TOLERANCES ON ANGLES + .02 XXDEC. + .005 + 10				MATERIAL 316 Stainless Steel			
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE FOR MFG. TOLERANCES, INCLUDING ROUND PUNCHED AND DRILLED HOLES, SEE SPI-E6				FINISH			
DRAWN S, Co X				CHECKER			
APPROVED				GROUP CODE			
NEXT ASSY.		USED ON		SIZE A		DRAWING NO. 6R	
				SCALE 2:1		SHEET OF	
Cable Support							



1/2 x 20 Hex Head Bolt 1/2 inch Long

REV		ECO NO.		DATE		APPD	
UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. XXDEC. TOLERANCES ON ANGLES + .02 - .005 ± .10				MATERIAL 316 Stainless Steel			
COMMERCIAL TOLERANCES GOVERN MATERIAL STOCK SIZE. FOR MFG. TOLERANCES, INCLUDING ROUND PUNCHED AND DRILLED HOLES, SEE SP1-E5				FINISH			
DRAWN S. Cox				CHECKER			
APPROVED				GROUP CODE			
NEXT ASSY.		USED ON		Cable Support			
				SIZE A		DRAWING NO. 5R	
				SCALE 2:1		SHEET OF	




LOCK PIN

SCALE: FULL

DR. BY: M. PISKACEK

NOV. 18, 1976




INPUT PIPE

SCALE: 2/1

DR. BY: M. PISKACEK

OCT. 28, 1976




OUTPUT PIPE


SCALE: 2/1

DR BY: M PISKACEK


OCT. 28, 1973



SECTION A-A



$\frac{1}{32} \times 45^\circ$ CHAM-TYP



SECTION A-A

VITA

Steven Paul Cox was born a U.S. citizen in Los Angeles, California, in 1955. Raised in Burbank, California, he attended John Burroughs High School and graduated in 1973. The next five years were spent as a freelance photographer in the Los Angeles area, camera salesman, and a photography instructor for the Burbank Adult School Program. During these five years, he also obtained his B.S.E.E. from the California State University at Northridge.

During his last summer in California before coming to Rochester in 1978, he worked for Rockwell International on the Clinch River Fast Breeder Reactor project. Finding no joy in government contracts, he left Los Angeles to pursue his Masters of Science in Photographic Science and Instrumentation at the Rochester Institute of Technology. While attending RIT, Mr. Cox has worked part time for the Itek Corporation in the design of microprocessor controlled graphic arts cameras. In the spring of 1981, this thesis was presented at the annual SPSE conference in New York City.

Mr. Cox hopes to some day return to his native California.