



Western Washington University
Western CEDAR

Salish Sea Ecosystem Conference

2014 Salish Sea Ecosystem Conference
(Seattle, Wash.)

May 2nd, 10:30 AM - 12:00 PM

Use of Flowcam Technology for Phytoplankton Monitoring in Central Puget Sound

Gabriela Hannach

King County (Wash.). Department of Natural Resources and Parks, gabriela.hannach@kingcounty.gov

Lyndsey Sandwick

King County Environmental Lab

Follow this and additional works at: <https://cedar.wwu.edu/ssec>



Part of the [Terrestrial and Aquatic Ecology Commons](#)

Hannach, Gabriela and Sandwick, Lyndsey, "Use of Flowcam Technology for Phytoplankton Monitoring in Central Puget Sound" (2014). *Salish Sea Ecosystem Conference*. 52.
<https://cedar.wwu.edu/ssec/2014ssec/Day3/52>

This Event is brought to you for free and open access by the Conferences and Events at Western CEDAR. It has been accepted for inclusion in Salish Sea Ecosystem Conference by an authorized administrator of Western CEDAR. For more information, please contact westerncedar@wwu.edu.

Use of FlowCAM Technology for Phytoplankton Monitoring in Central Puget Sound

Gabriela Hannach
Lyndsey Sandwick

King County Environmental Lab
King County Department of Natural Resources and Parks
Seattle, Washington




King County


What is FlowCAM?

Imaging particle analyzer

- Flow cell + microscope + camera
 - each particle in field of view is imaged, analyzed, saved
 - single file contains all images from individual particles
- Pattern recognition software



FlowCAM in action



Why FlowCAM?

- Puget Sound phytoplankton: significant data gap
- Need ability to quantitatively characterize Puget Sound phytoplankton
- King County long term marine monitoring program
- Qualitative phytoplankton since 2008
- FlowCAM:
 - Population changes over temporal and spatial scales
 - Interaction with multiple environmental variables (complex data)
 - Biodiversity/Biomass/Productivity
 - Trophic structure
 - Rapid identification and counts of HAB species




Sample Preparation

- Fresh live sample is best
- Minimal preparation
 - Filter to separate size fractions
100 μm , 300 μm
 - Large size fraction using microscopy or FlowCAM at lower magnification
 - Dilution or concentration for certain samples



FlowCAM Hardware




Magnification	Flow cell depth
4x	300 μm
10x	100 μm
20x	50 μm

Run time: 10-60 minutes

Yield: 1000-5000 images

Data Acquisition Modes

- Autoimage Mode: Captures particle images at a user-defined rate; images everything in field of view.
- Trigger Mode: Captures particle images only if a fluorescence signal is detected; images only chlorophyll containing cells.



Particle or not?



Edge Trace:

The perceived edge of the particle is shown in blue.

Binarization:

Pixels that are considered to be part of the particle are shown in purple.








Image Data Processing

- Sort images using simple data filters
 - e.g. size, shape, color
- Sort images using complex pattern recognition
 - Goal: Auto-sorting into taxa by image recognition
 - Build image libraries (taxa-specific)
 - 'training sets' for pattern recognition
 - Used in advanced filter development
 - Based on combinations of large number of measured particle properties
 - Continuous development, limited capability, too many taxa...

Particle Properties

measured and stored for each captured image


- ID #
- Area (ABD)
- Area (Filled)
- **Aspect Ratio**
- Average blue
- Average green
- Average red
- Calibration factor
- Calibration image
- Camera
- Capture X
- Capture Y
- Ch1 area
- Ch1 peak
- Ch1 width
- Ch2 area
- Ch2 peak
- Ch2 width
- Ch2/Ch1 ratio
- Circle fit
- **Circularity**
- Circularity (Hu)
- Compactness
- Convex perimeter
- Convexity
- D[4,3](ABD)
- D[4,3](ESD)
- Date
- **Diameter (ABD)**
- Diameter (ESD)
- Edge Gradient
- Elongation
- Feret angle max
- Feret angle min
- Fiber curl
- Fiber straightness
- Filter score
- Geodesic length
- Geodesic aspect ratio
- Geodesic thickness
- Image file
- Image height
- Image width
- Image X
- Image Y
- Intensity
- **Length**
- Particles per chain
- Perimeter
- Ratio blue/green
- Ratio red/blue
- **Ratio red/green**
- Roughness
- Scatter area
- Scatter peak
- Scatter width
- Sigma intensity
- Source image
- Sum intensity
- Symmetry
- Time
- Timestamp
- Transparency
- **Volume (ABD)**
- Volume (ESD)
- Width



2014 Sampling


Pt. Jefferson
East Passage
Dockton Pk.

Point Wells
Westpoint TP Outfall
Elliott Bay
South TP Outfall
Fauntleroy




2014 FlowCAM Data


Central Basin: North



Central Basin: South







Quartermaster Harbor



Size Fraction (μm)


April 10, 2014

Taxa classification using FlowCAM ($<100 \mu\text{m}$ fraction)



Analysis of Harmful Algal Blooms

July 2013: *Heterosigma* bloom implicated in
Quartermaster Harbor fish kill




12-17 μm ABD (area based diameter) size fraction

Technical Challenges

- Flow cell size limit- currently using 100 μm and 300 μm
 - Exchanging objectives takes time - Not practical for each sample
 - Magnification-specific libraries are needed for image recognition
- Image resolution
- Long run times for sparse samples (>60min)
- Image libraries have significant limitations (natural heterogeneity)
- Taxa sorting is limited by capability of the software
 - Lots of manual sifting and sorting
- Data management
 - Data volume, processing speed of large files, primitive output



20 μ m



20 μ m

Acknowledgments

Special thanks to the

King County Environmental Lab Field Unit



20 μ m



20 μ m

