### Marshall University Marshall Digital Scholar

**Chemistry Faculty Research** 

Chemistry

Fall 9-2011

### An Integrated Systems Approach to Deconstructing Glycosylation

L.K.Mahal

John F. Rakus Marshall University, rakus@marshall.edu

Kanoelani Pilobello

P. Agrawal

Follow this and additional works at: http://mds.marshall.edu/chemistry\_faculty Part of the <u>Organic Chemistry Commons</u>

### **Recommended** Citation


Mahal, L. K., Rakus, J., Pilobello, K., Agrawal, P. (2011, September). An integrated systems approach to deconstructing glycosylation. Poster session presented at the NIH Director's Pioneer Award Symposium, Bethesda, MD.

This Article is brought to you for free and open access by the Chemistry at Marshall Digital Scholar. It has been accepted for inclusion in Chemistry Faculty Research by an authorized administrator of Marshall Digital Scholar. For more information, please contact zhangj@marshall.edu.



# **AN INTEGRATED SYSTEMS APPROACH TO DECONSTRUCTING GLYCOSYLATION** Lara K. Mahal, Ph.D. Rakus, J., Pilobello, K., Agrawal, P.


## **Lectin Microarray Results**




**CONCLUSION:** The lectin microarray platform is suitable to identify distinct, carbohydrate-based cell surface features which differentiate certain cancer tissue types.

**GOAL:** To deconstruct the regulation of the complex, non-template driven synthesis of carbohydrates using high-throughput microarray methods

**STRATEGY:** Isolate matched glycomic and transcriptomic samples from the NCI-60 panel, or metabolically perturbed cell lines and subject to a "systems biology" approach of lectin and genomic microarray followed by the decomposition and integration of the data sets



New York University, New York, NY



-We have used a "systems" approach to identify regulatory effects on the glycome

-We are now expanding our analysis to include glycomic analysis of the entire NCI-60 panel and genomic analysis with our custom designed arrays -We are perturbing the metabolic state of select cell lines in order to determine the effect on glycome expression

## **Decomposition Analysis**

	COLON	LEUKEMIA	MELANOMA
ec	eigen-celltype 2	APA, Con A, DSA, GNA, LcH, LPA, NPA, PAA, SNA (vector), UDA, HHL, MAL-II, CVN, SVN, GRFT, ACA, PSL, TJA-I, AMA, Calsepa, IRA, RS-IIL	High mannose, α2,6 sialic acid ( <i>N</i> -linked)
	eigen-celltype <b>3</b>	VRA, PTL-II, HAA	Terminal GalNAc ( <i>O</i> -linked)
c)	eigen-celltype <b>4</b>	TKA, PapGII	Glycolipid

**CONCLUSION:** Decomposing the lectin and gene expression data reveals 3 unique glycosylation patterns. N-linked expression is correlated, O-linked inversely positively. Each eigen-celltype has process specific gene/glycosylation patterns (not shown).

### **Future Directions**

