

1981

Interpretation of GC-MS data by reverse search and relative retention indices

Christopher Hein

College of William and Mary - Virginia Institute of Marine Science

Follow this and additional works at: <https://scholarworks.wm.edu/etd>



Part of the [Analytical Chemistry Commons](#)

Recommended Citation

Hein, Christopher, "Interpretation of GC-MS data by reverse search and relative retention indices" (1981). *Dissertations, Theses, and Masters Projects*. Paper 1539617520.
<https://dx.doi.org/doi:10.25773/v5-hk63-yh63>

This Thesis is brought to you for free and open access by the Theses, Dissertations, & Master Projects at W&M ScholarWorks. It has been accepted for inclusion in Dissertations, Theses, and Masters Projects by an authorized administrator of W&M ScholarWorks. For more information, please contact scholarworks@wm.edu.

INTERPRETATION OF GC-MS DATA BY REVERSE SEARCH
AND RELATIVE RETENTION INDICES

A Thesis

Presented to

The Faculty of the School of Marine Science
The College of William and Mary in Virginia

In Partial Fulfillment

Of the Requirements for the Degree of
Master of Arts

by

Christopher Hein

1981


ProQuest Number: 10626350

All rights reserved

INFORMATION TO ALL USERS

The quality of this reproduction is dependent upon the quality of the copy submitted.

In the unlikely event that the author did not send a complete manuscript and there are missing pages, these will be noted. Also, if material had to be removed, a note will indicate the deletion.



ProQuest 10626350

Published by ProQuest LLC (2017). Copyright of the Dissertation is held by the Author.

All rights reserved.

This work is protected against unauthorized copying under Title 17, United States Code
Microform Edition © ProQuest LLC.

ProQuest LLC.
789 East Eisenhower Parkway
P.O. Box 1346
Ann Arbor, MI 48106 - 1346

APPROVAL SHEET


This Thesis is submitted in partial fulfillment of
the requirements for the degree of


Master of Arts


Christopher Hein

Approved, December 1981


Rudolph Bieri, Ph.D.


Chi Wu Su, Ph.D.


Robert Lukens


Craig Smith, Ph.D.



David Evans, Ph.D.

TABLE OF CONTENTS

	Page
ACKNOWLEDGEMENTS	iv
LIST OF FIGURES	v
ABSTRACT	vii
INTRODUCTION	2
EXPERIMENTAL	13
RESULTS	14
Reference Mass Spectra	14
Relative Retention Index	15
Library Structure	15
Reference Spectra Acquisition	17
Data Control	19
Reverse Search Program	19
Analysis Parameters	20
GC-MS File Parameters	25
Reference Record Access	26
Retention Prediction	27
Reference Evaluation in Unknown MS	29
Interpretation of Observations	33
Analysis Output	42
DISCUSSION	50
CONCLUSION	56
REFERENCES	57
VITA	59
APPENDIX	
Program ONE	A1
Subroutine OPEN1	A7
Function CHNGE	A8
Subroutine GET	A9
Program TWO	A11
Subroutine GRAPH	A18
Program THREE	A19
Subroutine OPEN3	A27

ACKNOWLEDGEMENTS

I wish to thank Professor R. H. Bieri for his kind guidance and constant support. I am also indebted to Dr. Kent Cueman for providing an introduction into the access of the DuPont GCMS data System, to Buster Blystone and Robert Lukens for generous assistance in programming, to P. C. Gillette of Case Western University for helpful mathematical discussions, and to Molly Lutcavage for assistance in the production of this thesis.

Drafts and final copy of this report were prepared by the VIMS Report Center.

LIST OF FIGURES

Figure	Page
1. Illustration of Linear Combination of Reference Mass Spectra to Explain an Observed Mass Spectrum	9
2. Evaluation of the Occurrence of a Reference Mass Spectrum in a Series of Unknown Mass Spectra by the Method of Pesyna et al. (14)	11
3. Logical Structure of Program ONE	21
4. Decoded Reference Record	28
5. Search for 1-Methyl Naphthalene in MS#73 of Sample STATION J G3.2	30
6. Search for 1-Methyl Naphthalene in MS#84 of Sample STATION J G3.2	31
7. Attributable Intensity within the Search Window for 1-Methyl Naphthalene	32
8. Reconstructed Gas Chromatograph for Sample STATION J G3.2	34
9. Search for Fluorene in MS#216 of Sample STATION J G3.2	36
10. Search for Fluorene in MS#218 of Sample STATION J G3.2	37
11. Attributable Intensity within the Search Window for Fluorene	39
12. Reconstructed Gas Chromatogram of STATION J G3.2 with Illustrations of Intensity Attributable to Fluorene and 1-Methyl Phenanthrene	40
13. Attributable Intensity within the Search Window for 1-Methyl Phenanthrene	41
14. Reference Data for Two Structural Isomers of Molecular Weight 192	43
15. Analysis Printout of Program ONE, for Evaluation of STATION J G3.2 vs Library of Aromatic Compounds	44
16. Analysis Printout of Program ONE, for Evaluation of STATION J G3.2 vs Library of Pesticides and PCBs	47

Figure	Page
17. Illustration of the Relationship Between Intensity, Mass and Time in a GC-MS Series of Measurements During Elution of a Peak	53

ABSTRACT

A technique was developed to interpret GC-MS data from complex environmental samples. Two methods were combined into this technique. The first is a method to evaluate the occurrence of a given reference mass spectrum in an unknown GC-MS data set. The second is a method to optimize the use of gas chromatographic information via relative retention indices. The combination of these two methods provides a thorough, rapid technique to evaluate the potential of occurrence for any number of reference mass spectra and their associated retention indices.

INTERPRETATION OF GC-MS DATA BY REVERSE SEARCH
AND RELATIVE RETENTION INDICES

INTRODUCTION

Gas chromatography-mass spectrometry (GC-MS) is a technique used to determine the identity and concentration of components present in complex mixtures of organic compounds. The gas chromatography employed in such analyses may produce several hundred separable peaks. Subsequent mass spectrometry provides characterization of the compounds present in the gas chromatographic effluent. GC-MS analysis produces data on the retention time, and mass spectral (mass, intensity) characteristics of each compound present.

GC-MS data must be interpreted to provide meaningful results about the various compounds in a sample mixture. Gas chromatographic and mass spectral characteristics of a given compound are reproducible within limits. Compound identification requires that both GC retention and mass spectral parameters of the unknown compound match those of a standard compound. Identification by GC retention may be limited by the simultaneous or unresolvable elution of two or more compounds. Mass spectral identification is often limited by the similarity of fragmentation patterns in structural isomers. Ideally, GC-MS data analysis is only limited by the simultaneous elution of compounds with similar mass spectra. When this situation occurs, other methods of analysis must be employed.

Factors other than this theoretical limitation in GC-MS data interpretation provide far greater barriers to the complete interpretation of GC-MS data, especially from complex mixtures. As

the number of compounds increases in a sample, the probability for overlapping elution from the GC column also increases. Simultaneous elution produces superimposed mass spectra which in turn increases the difficulty of mass spectral interpretation. Simultaneous elution of various compounds remains a significant problem despite the fact that modern GC-MS analysis incorporates the best available high resolution chromatography. The data produced from an hour long GC-MS analysis may require several days of manual interpretation. Thus, manual interpretation of GC-MS data is an inefficient and time consuming process for routine analyses.

Computer programs have been developed to assist in the interpretation of GC-MS data. The goals of such programs include reducing the time of analysis and increasing the success of identifications. Interpretative programs have in most cases been used on mass spectral information alone. With increasingly complicated applications of GC-MS analysis, efforts have been made to pair retention information with mass spectral data to increase the probability of correctly identifying a compound.

The first computer programs for GC-MS analysis were based on reduced spectra and simulated manual interpretative procedures. Significant features of an unknown spectrum are selected with various criteria. They may then be compared with a library of features selected from reference mass spectra. This strategy of comparing an

unknown spectrum against known characteristics has in historical perspective become known as a 'forward search'.

Work in the field of computerized GC-MS analysis began with two papers co-authored by Crawford and Morrison (1,2). The concept presented in the first paper was to reduce unknown spectra to a small set of most significant mass intensities for comparison against a library of known spectra that had been reduced in similar fashion. The second paper demonstrated that a mass spectrum may be classified according to chemical type via condensed mass spectra.

In their first approach Crawford and Morrison took the six largest peaks of an unknown spectrum and compared this subset with the six largest peaks from 3200 known spectra. The ten library compounds with the greatest similarity to the unknown spectra were reported as possible identifications for the unknown compound. (1) A major problem with this concept is the frequent occurrence of an unknown mass spectrum matching the mass spectra of more than one compound. In this case unpredictable results occur including misidentifications and elimination of the correct spectrum.

By summing fragment intensities of a mass spectrum at intervals of 14 mass units, Crawford and Morrison realized that the resulting condensed mass spectrum was often indicative of a compound's functional group (2). This second strategy supposes that the compounds under investigation as a whole or in part contain saturated aliphatic chains, where 14 mass units signify the loss of CH_2

fragments. Difficulties with this method again relate to the simultaneous presence of an unknown number of spectra, leading to incorrect interpretation.

A number of scientists subsequently developed these concepts and combinations of these concepts for the purpose of identifying unknown spectra. Efficient use of computer time, number of elements versus success in identification, and determining the profile of the average condensed spectra for all functional groups were investigated. (3-7) Two later developments attempted to reduce the frequent problem of superimposed mass spectra. Dromey et al. (8) analyzed mass chromatogram peak shape, and Blaisdell et al. (9) developed a method attempting to separate unknown mixtures of mass spectra into contributing sets.

In 1975, Abramson first pointed out that by reversing the traditional order of library and unknown comparison, the problem of mixed or superimposed spectra could be avoided (10). This strategy had been utilized in part by Pereira et al. (11), and Sweeley et al. (12). In 1976, factor analysis was shown to be applicable to GC-MS data interpretation by Ritter et al. (13). This promising technique is still under development for use in applied GC-MS data analysis.

If a mass spectrum is generated by the simultaneous presence of more than one compound, it represents the simple algebraic sum of all contributing spectral components at their respective concentrations (10). Identification of a known compound in an unknown spectrum would

therefore require that all ions of the known spectrum be present in the unknown mass spectrum.

Pesyna et al. (14) presented a simple two step procedure to evaluate the contribution of a known mass spectrum in an unknown mass spectrum. Intensities of the unknown mass spectrum are divided by the corresponding relative intensities for all masses in the reference mass spectrum, yielding a series of normalized ratios expressed by the variable P below:

$$P_m = \frac{I_m}{RI_m}$$

where I is unknown spectrum intensity, RI is known spectrum relative intensity, and the subscript m indicates the mass, specific for both unknown and reference mass spectrum. The smallest observed P value determines the relative number of times the reference spectrum could be removed from the unknown mass spectrum before a fragment required for the presence of a reference mass spectrum would be exhausted. The intensity of the reference spectrum in the unknown mass spectrum is then obtained by multiplication of the relative number of reference spectra found to fit in the unknown P_{min} , by the sum of RI in each basic reference spectrum. This method assumes that at least one mass in the unknown spectrum is solely from the known compounds, and requires that all elements of the known be present in the unknown mass spectrum. The authors describe a number of additional factors to compensate for distortions in mass spectra from various instruments. The purpose of their study was to test large numbers of known mass

spectra for their possible occurrence in GC-MS data with utmost efficiency.

Gas chromatography (GC) is by itself a powerful analytical method. Separation of a complex mixture by GC leads to the elution of the various components at various points in time. Retention time is subject to variations in chromatographic parameters, thus retention data are usually modified from absolute retention times (injection to elution) and are referred to as relative retention indices. Kovats (15) found that by measuring the locations of eluting peaks in a gas chromatogram relative to a series of n-alkane marker compounds, the effects of variations in absolute retention are minimized. Such retention markers are either coinjected or they may be ubiquitous and clearly identifiable. Any point on the chromatogram that falls within the range of the markers may be expressed via a relative retention index (RRI), adequately defined for linear temperature programming by

$$RRI = \left[\frac{(t_x - t_n) \times 100}{(t_{n+1} - t_n)} \right] + 100n$$

n = number of carbon atoms in n-alkane marker, $(t_x - t_n)$ = retention time difference between unknown and nearest marker with $t_n < t_x$, $(t_{n+1} - t_n)$ = retention time difference of the nearest marker compounds to the unknown.

A relative retention index system will be most accurate if both retention markers and unknown have very similar chromatographic

behavior. Bieri (16) and Lee (17) found it necessary to define an aromatic retention index (ARI) system for work with aromatic compounds. Within one type of compound, the choice of markers is not critical for a retention index system, although a standardized system would facilitate interlaboratory comparison of results. The ARI system described by Bieri used 7 markers (naphthalene, biphenyl, phenanthrene, chrysene, pyrene, perylene and benzo(ghi)perylene) that are given retention indices of 0, 100, 200,...600, respectively.

Recognizing the importance of gas chromatography, Blaisdell and Sweeley (18) developed a GC-MS data interpretive program that incorporated RRI as an identifying parameter. This program observed the characteristic retention of the reference compounds under analysis, and only considered such compounds in their possible retention zones. Each unknown mass spectrum was best accounted for with linear combinations of the reference mass spectra under consideration. The result of this analysis was the printing of concentrations for up to ten library compounds for each unknown spectrum.

Expression of an observed mass spectrum by linear combinations of reference mass spectra takes the form of:

$$MS(\text{observed}) = [MS(N) * C(N)] + [MS(M) * C(M)] + \dots [MS(\text{un}) * C(\text{un})]$$

$$\text{where } MS(X) = [M(i)*RA(i)] + [M(j)*RA(j)] + [M(k) * RA(k)] \dots$$

(for X = reference mass spectrum N, M,...) and

M(x) is any mass occurrence in the spectrum of compound (X)

$RA(x)$ is a relative abundance associated with $M(x)$,

(for $x = \text{mass } i, j, \dots$)

$C(X)$ is an absolute intensity factor associated with $MS(X)$

$MS(\text{un})$ is unexplained mass spectrum

$C(\text{un})$ is intensity of unexplained mass spectrum

Blaisdell and Sweeley's approach was to adjust $C(N)$, $C(M)$, $C(L)$... in order to minimize $C(\text{un})$. This concept is illustrated below in Figure 1.




Figure 1. Diagram Illustrating the Strategy of Linear Combination of Reference Mass Spectra to Explain an Observed MS.

Success of this strategy requires prior knowledge of the compounds contributing to the observed mass spectrum. Approximately 30,000 reference mass spectra have been established. Many essentially similar mass spectra occur in similar retention zones. Identification of a compound may be beyond the scope of a given GC-MS analysis under

these circumstances. The linear combination of mass spectra is a method to optimize the interpretation of an observed mass spectrum with available data. The technique does have a drawback in that it can mask the limitations of a GC-MS analysis.

A compound undergoing gas chromatographic analysis will elute in one concentration vs time peak. In GC-MS analysis, a compound may be located by total ion current, however such measurement is no better than gas chromatographic analysis. Mass spectral analysis in conjunction with gas chromatography allows a greater differentiation between compounds in a complex gas chromatographic separation. Drawing of a compound's fragment intensity vs time allows the designation of time locations where the compound may exist. A common technique to search for the occurrence of a reference spectrum is to select a unique mass for such analysis. A thorough technique to measure the occurrence of a reference mass spectrum was presented by Pesyna et al as previously discussed. This technique determines the area (in scan time *ion intensity) in which a reference mass spectrum may exist in unknown GC-MS data.

The approach of Pesyna et al to GC-MS data interpretation is illustrated in Figure 2. In this manner, a two dimensional plot of the occurrence of a reference mass spectrum may be extracted from three dimensional data.

If more than one peak forms in the allowable intensity for a reference mass spectrum, additional analysis parameters are required.




Figure 2. Evaluation of the Occurrence of a Known Mass Spectrum in in a Series of Unknown GC-MS Data via the Method of Pesyna et al. (14).

An additional analysis parameter is available using gas chromatographic retention. A compound may potentially occur within only a limited retention range. If more than one peak may occur in a mass spectrum's intensity vs time plot, within the valid retention region, all GC-MS analysis parameters are satisfied. More definite analysis requires further knowledge of the sample under GC-MS analysis, or another analytical technique.

Pesyna and co-workers' technique would be optimized by integration with a relative retention index system. Such a combination would represent maximized utilization of available data,

and illustrate the limitations of such data. The intent of this study is to implement this combination of reverse mass spectral search with retention indices.

EXPERIMENTAL

Raw GC-MS data were generated with a DuPont 21-094B mass spectrometer and data system. The mass spectrometer scanning speed was modified to one scan cycle per 2.3 seconds, representing the fastest possible scanning speed without calibration difficulties. A Varian 2700 gas chromatograph, modified for capillary operation, was interfaced to the mass spectrometer. The DuPont data system incorporated a Hewlett Packard 2100A computer with a 16K 16 bit word memory.

A series of programs were developed in Fortran IV to provide the DuPont data system with reverse search capabilities. The programs were developed in a DOS III environment, using disk files of library data. A HP 7900A dual 1.5M word disk drive, with removable top disk, was employed with the HP 2100A. Reverse search programs and library files were retained permanently on removable disks, but were copied to the lower disk for analyses. In this manner, various removable disks of GC-MS data could be interfaced with the analysis program. Printouts of analytical results were generated during an analysis on a Tektronix 4012 terminal and 4610 hard copy unit.

RESULTS

An automated interpretative system was developed to supplant manual GC-MS data analysis. A reverse search strategy was determined to be preferable to forward search strategy because of the highly complex mixtures of mass spectra that are encountered in the environmental samples under investigation. The procedure allows an evaluation of the presence of a reference spectrum in unknown mass spectra, in conjunction with a relative retention index system for optimal use of qualitative GC data.

Six individuals systems had to be created before the reverse search could be implemented. These included: 1. a reference mass spectra system, 2. a retention index system, 3. a library for storage and retrieval of reference data, 4. a program to purify reference spectra from GC-MS data (Program THREE) for the collection of reference data, 5. a program to control data in the reference file (Program TWO), and finally, 6. the construction of the desired program itself (called Program ONE). Description of these systems follows.

Reference Mass Spectra

The various fragments of a given compound occur in constant intensity ratios to one another under reproduced mass spectral parameters. It is common practice to express the contributions of different fragments in relative abundance (RA) units. The most abundant fragment (the base peak) was assigned a relative abundance of 100. Other fragments were then expressed in percentage points

relative to the base peak. The reference mass spectrum is thus composed of pairs of information on fragment mass and RA. This data may be used to describe a compound's mass spectrum and thereby identify the occurrence of the compound in GC-MS data. Within the linear operating range of the instrument, the absolute intensity of a mass spectrum is proportional to the amount of compound analysed, thus quantification is possible.

Relative Retention Index

Relative retention index (RRI) systems have been routinely used by this laboratory for manual interpretation. The ARI system described by Bieri (16) and the KRI system described by Kovats were incorporated into the reverse search program to facilitate comparison of data. The version of the reverse search program used with the KRI system required modifying the format of the RRI, dropping the least significant figure and adding a digit in the 10^4 place. The systems are otherwise entirely comparable. In the following, emphasis will be on the ARI system.

Library Structure

In order to perform a reverse search for a series of aromatic compounds, records describing the corresponding reference compounds must be available for use by the computer. This was accomplished by creating a disk file that is read by the computer during analysis. This file was designed to be read in records of 128 integer words, one record for each reference compound. The reference file is ordered by

increasing ARI, thus the first record is of naphthalene (ARI=0) the last record is for benzo(ghi)perylene (ARI=600).

Each record in the reference file has 100 integer words reserved for mass and relative abundance data. Each of the 100 words contains two pieces of data: a fragment mass and the corresponding relative abundance. The fragment mass is limited to a range of 1 to 640 amu, and the relative abundance data must be expressed by numbers between 1 and 100. These data are compacted into a single word to reduce both reference file length, and to reduce the length of time necessary to read the disk. The fragment mass and relative abundance data are compacted into a 16 bit integer word (with range of expression from -32768 to 32767) in the following manner: The relative abundance data are reduced by 1 resulting in a value between 0 and 99. Fragment mass is reduced by 320, yielding a possible value of -319 to 320. This reduced fragment mass is then multiplied by 100 and added to the reduced relative abundance. These efforts produce a compacted number describing mass spectral fragment in the range of -31900 (mass of 1 with an RA of 1) to 32099 (mass of 640, RA of 100).

Decoding of the coded fragment data is accomplished by reversal of the compacting procedure while taking the absolute value of the RA remnant. Word 101 of the reference compound record stores the number of fragments in the reference mass spectrum. Word 102 records the ARI of the reference compound multiplied by 10 to allow one digit past the decimal point in integer number. (In use with the KRI system, the multiplication by 10 must be abandoned in order to allow the

additional descriptor in the 10^4 position of these values). Words 104-123 are used to record up to 40 alphanumeric characters of the compound name. Words 103, and 124-128 are available for parameters possibly useful in the future.

The library disk file was created with the Hewlett Packard Extended File Management Package (EFMP). The file was created on Pack Number (PN) 501, using 1800 sectors of disk area. One file, named ARILI, was set up on PN501 using 1200 of the available 1800 sectors. Thus, the limitation of storage space is 1200 reference records per disk. Access to the reference library may be gained with a Fortran program via HP EFMP EXEC commands, as will be discussed later.

Reference Spectra Acquisition

A computer program (Program THREE) was designed to facilitate compilation of reference data. Although it is possible to enter a reference record manually via Program TWO, it is more convenient to gather reference data from GC-MS records. Program THREE recognizes the phenomenon of mass chromatographic (MC) peaks that form due to the occurrence of a compound in a GC separation, and thus is capable of recognizing which masses may be relevant to the mass spectrum of the compound. Dromey et al. (8) described this technique for spectrum purification. The method is defeated when mass chromatographic peaks are distorted due to common fragment masses of neighboring compounds. Distorted mass chromatographic peak profiles make MC peak time to

compound assignment ambiguous, thereby increasing the difficulty of spectral purification.

Program THREE requests the location of the MS nearest the center in which all ions of a reference compound peak. The computer will then select the masses that form intensity peaks during the 5 mass spectra centered around the designated mass spectrum. For this select group of masses, exact peak center location is stored along with the ion counts (averaged for the three spectra nearest the peak center). The operator is then presented with a table of peak center time versus summed abundance data, and asked to define the time interval for inclusion into the reference spectrum. The mass spectrum that may be included into the designated time location is presented to the operator. If the operator rejects the reference mass spectrum, a new time interval may be designated, or a new set of GC-MS data may be selected for an improved mass spectrum. If the operator accepts the mass spectrum, the computer will request the compound's name and ARI. Finally the new reference spectrum is inserted by the computer into the appropriate position in the reference record file.

The majority of reference spectra in this study were extracted from identified unknowns in archived GC-MS data. A smaller number of reference spectra (approximately 100) were taken from the analyses of standards.

Data Control

A second program (Program TWO) was designed to edit the data on the reference disk file. Program TWO allows the user to display, alter, input, or delete reference data. The display function will plot the reference mass spectrum, name and ARI of a specified compound. This allows checking the record and referencing the library data. The alter function of Program TWO allows changing the ARI or compound name in a reference record. It also allows the addition and deletion of fragments in a reference mass spectrum. The input function of program TWO provides for a manual input of a reference record, thus supplementing Program THREE. The input function may be involved in cases where reference compounds are unavailable, but reference data are published. The delete function allows removal of an unwanted record, and subsequent repacking of the reference file over the vacancy.

Reverse Search Program

The program named ONE was created to perform the desired goal of this project: reverse search of unknown GC-MS data with regard to retention. The program is in the most general terms an interaction point between reference records and unknown GC-MS data with periodic status reports of the interaction sent to the user. This program may be divided into three sections. In the first section, certain variables must be defined before the analysis may be performed. In the second section, a repetitive cycle is carried out once for the

interpretation of each reference compound. The third and final section of Program ONE is the closure of the files accessed during analysis. The three divisions, and their respective steps are outlined in Figure 3.

The central interpretative part of Program ONE includes the following basic steps: In the first step, a reference record of a compound is read from disk to computer memory. The second step involves calculating where the compound could occur in the unknown data and thus define a window (a series of unknown mass spectra) for further investigation. The third step involves evaluation of the mass spectra in the analysis window for intensity attributable to the reference mass spectrum. The fourth step is to interpret these intensities for a peak profile (a requirement of GC analysis). Finally, the reference compound name is printed along with any observed peak data. The computer is then free to return back to the first step for the next compound under investigation, until an end of file mark is reached for either the reference or the GC-MS files.

Reverse Search Analysis Parameters

Early development work on the reverse search program demonstrated that some uncertainty existed in both relative index and reference mass spectra. Slight changes in ARI and reference MS are inevitable and must be accommodated. Parameters were designed into Program ONE to deal with these uncertainties in both reference and unknown GC-MS data.




Figure 3. Logical Structure of Program ONE.

Three analysis parameters must be defined upon initiation of Program ONE: IWND0, MINCD and N2REJ. Analysis parameters IWND0 and MINCD are factors that allow for uncertainty in ARI and MS data, respectively. The third parameter, N2REJ, was necessary to allow for the presence of minor irrelevant fragments in the reference mass spectra. Values required for this third parameter would decrease to zero if the irrelevant fragments could be recognized and eliminated from the reference data. These three parameters are designed to permit successful analysis using reference and GC-MS data of various qualities. Increasing the value of these parameters allows greater differences to exist in the matching of reference and unknown data.

The variable IWND0 (mnemonic for integer half window) defines the number of spectra to be analyzed on either side of the unknown spectrum of predicted occurrence for the reference compound. Acceptable IWND0 values are from 1 to 12 spectra. Two factors are of importance to the selection of an optimum IWND0 value. First, a wide enough range must be adopted to insure analysis of the compound over the region where it may be significantly represented. Second, the range for analysis must be limited to the region where the compound may possibly occur.

The maximum standard deviation of ARI observed in a 27 compound standard mixture analysed 14 times on one column was 0.63 ARI units (20). Working at a 95% confidence level would necessitate observation of ± 3 standard deviations from the average ARI, or ± 2 ARI units. Two ARI units may translate to a maximum scan time equivalent of 3

spectra, if evaluated within the ARI range of 100 to 200 where each ARI unit corresponds to about 1.4 scans.

The allowable range of IWND0 values is in excess of the statistically required range. This is necessary to allow full base-to-base peak analysis, and to allow flexibility with other systems that may be less accurate than the one used. Experimentation with values of IWND0 has led to common use of 4 spectra for this parameter, thus 9 unknown spectra are evaluated for each reference compound.

Ideally, the fragments produced by the mass spectrometry of a compound will occur in constant intensity relationships with one another. In practice, relative fragment abundances may vary by several percent. In the case of the system used in this experiment, the accuracy with which fragment ions can be measured varies with the square root of the total number of counts (21). Quantification via the technique of Pesyna et al. is based on mass spectral intensities. Error in this quantification is thus proportional to the square root of the absolute intensity of the mass spectral fragments. The variable MINCD was designed into Program ONE to allow a threshold minimum count for the determination of intensity attributable to a reference spectrum contained in an unknown spectrum. Limitation of the quantification to ions with counts of 100 ions or more will provide a relative error of less than 10% with each of the allowed fragments in the unknown data. It is desirable to use as many fragments as possible to define the occurrence of a known spectrum in

unknown GC-MS data, thus insuring a greater resistance to interference from neighboring compounds with masses in common with the reference compound. An optimum MINCD value therefore depends on two factors; first, the error associated with the ion counts used in quantification, and second, the number of intensities that are usable for the quantification. The number of intensities usable for quantification is a value that is determined by the number of fragments in the reference MS, and the number of fragments in the unknown GC-MS data that fall above the MINCD threshold intensity. Experimentally, a MINCD value of 10 has been found to compromise between the error of the single measurement and the accuracy afforded from multi-measurement compound definition.

The variable N2REJ (for the "number at which to reject") was designed to permit successful operation of Program ONE with varying amounts of irrelevant fragments in the reference MS. Irrelevant fragments occur in reference mass spectra due to the possible difficulty in separating relevant from irrelevant fragments. Experience has demonstrated that the majority of irrelevant fragments are of low RA. Fragments of low RA are often invaluable parameters for defining a compound, hence they are included in this analysis. The problem of impurities in principle is simple to solve: a comparison of reference mass spectra from several standard analyses, or from very pure standards, should provide superior reference spectra. This problem was not treated because of the time that would be involved in such multiple standard analyses. The value of N2REJ

defines the sum of missing ions in the reference spectrum (in percent of the total reference mass spectrum) that defines the absence of a reference mass spectrum in an unknown mass spectrum. This places greater importance on the more abundant fragments, i.e. if N2REJ is set at 8%, the absence of four 2% fragments or one 8% (or any other combination totalling N2REJ% or more) terminates the analysis for the mass spectrum of the reference compound in the unknown mass spectrum. N2REJ may be described as controlling qualitative analysis by mass spectrometry, with greater values for this parameter allowing greater dissimilarity between the reference and unknown MS with a positive identification.

GC-MS File Parameters

Some information about the GC-MS analysis data set is required before a reverse search can be performed by Program ONE. This includes defining the GC-MS data set that is to be analyzed and establishing the retention times for the marker compounds. The GC-MS data files of the DuPont 21-094B data system are accessible via Hewlett Packard Executive Extended File Management Package commands. The GC-MS data files are all recorded on Pack Number 111, with a five character (two letters, three numbers) file name assigned at the initiation of the GC-MS analysis. This file name must be given to Program ONE, along with specific information necessary to locate the marker compounds in the GC-MS data set.

Reference Record Access

Reference data may be read from the reference file, as may GC-MS data, by a Fortran program via the use of Hewlett Packard Executive Extended File Management Package commands, or HP EFMP EXEC. Three such commands are the minimum necessary sequence for such access, and are presented below:

```
CALL EXEC (24,1,IOPNTB,N,ITRBUF,NOTRB,2,IERROR)
```

```
CALL EXEC (24,4,IFNA,IPKN,M,ISC,K,IERROR)
```

```
CALL EXEC (24,6,IFNA,IRCN,IBUF,IERROR)
```

The first command initiates EFMP abilities by setting up opened file table buffers, track buffers, etc. The second command opens the desired file (IFNA) on the desired pack (IPKN) to allow access to the file by record number in EFMP. The third command creates a copy of the desired record number (IRCN) into a given buffer (IBUF). Complete explanation of EFMP EXEC details are available in a HP DOSIII manual. In practice, checks are made for error messages during EFMP commands. The first two commands listed above need to be executed only once until the file is closed. The effort to gain EFMP access to the files (GCMS and Library) in the programs of this project has been reduced by combining the first two commands into a subroutine called OPEN1 (See Appendix) that only requires the desired file name and pack number be specified.

Once the desired reference record is read into a buffer, it must be decoded for use. This includes unpacking of the mass-relative

abundance data and restoration of the true ARI value, as described in the Library Structure section. A typical reference file record is presented in decoded form in Figure 4.

Retention Prediction

The next step in the reverse search analysis is to predict where in the GC-MS data the reference compound could occur. This is performed by converting the characteristic RRI of the reference compound into the equivalent retention time of the GC-MS analysis. A Fortran function CHNGE (See Appendix) was designed to interconvert ARI and retention time (RT). Use of CHNGE requires that the RT and ARI of the retention standards are placed into buffers during initiation of Program ONE for the interconversion calculations. The buffers must be held in common between the main program and CHNGE. The CHNGE function is invoked by the general form:

$$Y=CHNGE (X,R2S)$$

whereby, if R2S is positive, X is converted from an RRI value into an RT value and result is stored in Y. If R2S is negative, X is converted from an RT value into a corresponding ARI value and the result is stored in Y. If an X is sent for conversion that falls out of range, a negative CHNGE value is returned, where a -1 specifies before, and -2 specifies after, possible conversion range.

Uncertainty in compound retention and compound peak width requires that a number of mass spectra be investigated before and after the predicted retention for the reference compound. The number

of spectra that must be analyzed before and after the predicted retention time is designated by the analysis parameter IWND0. The window for reverse search analysis is thus defined by the mass spectra between and including (CHNGE (ARI,1)-IWND0) to (CHNGE(ARI,1)+IWND0).

Reference Evaluation in Unknown MS

The mass spectra that are included within the window for reverse search analysis must be analyzed one at a time. This is accomplished by a method adapted from the work of Pesyna et al. (14). This method is illustrated with a computer printout of the analysis procedure (Figures 5 to 7) available by using a version of Program ONE called Program HOW, altered to illustrate the intermediate numeric values of the reverse search. Figure 5 shows the results of an analysis of the first spectrum (#78) in the search window for the compound 1-methyl naphthalene.

A negative result was determined in the analysis presented in Figure 5 due to the absence of the reference mass spectrum in the unknown mass spectrum within the specified value of N2REJ. A positive result was later determined in the search for this compound as illustrated in Figure 6. The intermediate step between analysis and user printout is presented in Figure 7.

The search for 1-methyl naphthalene was performed in a relatively simple region of the GC-MS data from analysis of sample STATION J. The reconstructed gas chromatogram of this area of the GC-MS

1-METHYL NAPHTHALENE

REFERENCE SPECTRUM	ION COUNT	P RATIO
MASS	FROM SAMPLE	
	M.S.# 78	
143	22.	1.7
142	80.	.8
141	67.	.8
140	.	IONCOUNT MISSING, SUM RA
139	10.	1.0
126	5.	N.A., IONCOUNT < MINCD
116	8.	N.A., IONCOUNT < MINCD
115	20.	.8
113	.	IONCOUNT MISSING, SUM RA
91	32.	18.0
89	.	IONCOUNT MISSING, SUM RA
87	.	IONCOUNT MISSING, SUM RA
1	.	IONCOUNT MISSING, SUM RA

ALLOWABLE PERCENTAGE OF TOTAL REFERENCE MASS SPECTRUM (3X) MISSING
REFERENCE SPECTRUM IS THUS DEFINED NOT PRESENT IN THIS MASS SPECTRUM

Figure 5. Search for 1-Methyl Naphthalene in MS #73 of Sample STATION J G3.2

1-METHYL NAPHTHALENE

REFERENCE SPECTRUM MSS	R.A.	ION COUNT FROM SAMPLE M.S.#	84	P RATIO
143	13	289.		22.2
142	99	2228.		22.5
141	81	1878.		23.2
140	5	117.		23.4
139	10	239.		23.9
126	1	22.		22.0
116	2	49.		24.6
115	24	429.		17.9
113	1	20.		20.0
91	2	35.		17.5
89	2	32.		18.0
87	1	15.		15.0
78	1	.		IONCOUNT MISSING, SUM RA
77	1	28.		26.0
76	1	15.		15.0
75	1	30.		30.0
74	1	20.		20.0
71	6	143.		23.8
70	5	75.		15.0
69	3	54.		18.0
65	1	19.		19.0
63	4	68.		17.0
62	1	23.		23.0
57	4	37.		9.2
51	3	44.		14.7

SUM OF IONS IN SAMPLE
SPECTRUM ATTRIBUTABLE - SUM OF R.A.s x MINIMUM P
TO REFERENCE COMPOUND 2525. 273. 9.2

Figure 6. Search for 1-Methyl Naphthalene in MS #84 of Sample STATION J G3.2

1-METHYL NAPHTHALENE

SUM OF IONS IN SAMPLE
SPECTRUM ATTRIBUTABLE
TO REFERENCE COMPOUND

SPECTRUM
77
78
79
80
81
82
83
84
85
86
87
88
89

710.
.
.
84.
.
.
1911.
2525.
262.
.
.
.
.

PEAK ANALYSIS
MAXIMUM AREA
80.0 84.
> 83.7 4698.

NOTE: > ASSIGNED TO PEAK NEAREST REFERENCE R.I.

Figure 7. Attributable Intensity within the Search Window for 1-Methyl Naphthalene in Sample STATION J G3.2

analysis is presented in Figure 8-A. A bracket is drawn over the search window of scans 77 to 89.

A second example will be presented from a region of GC-MS data that is more complex. This example will center on the search for fluorene, which is known to occur at an ARI of about 152. This reference retention index is changed into a scan location that the search window is centered on. The search window is then defined as the predicted retention time plus and minus the value of the IWND0 parameter established in program initiation. This window is drawn over scans 211-223 in Figure 8-B.

Each unknown mass spectrum in the search window is analysed for the occurrence of the reference mass spectrum. The results of this mass spectral evaluation for two of the unknown mass spectra are presented in Figures 9 and 10. The fluorene search is carried out in a region of high background, as apparent in Figure 8-B. Two major components, and several minor components are obviously present. A number of compounds contribute to the background.

Interpretation of Observations

As the spectra in the analysis window are analyzed, the computer program records the number of ions attributable to the reference spectrum. When all spectra have been analyzed, the sum of ions attributable to the reference compound versus the spectrum number is analyzed (for the occurrence of peaks). If one or more peaks are found, their exact centers are found along with equivalent ARI, peak

SEQUEN 21 PAGE 8
GCID EP 59 STATION J G3.2
#SCANS 150 HRDCPY NO
%SCALE 70 REZERO NO
BASE 2251x2xx 2




Figure 8A. Reconstructed Gas Chromatogram for Sample STATION J G3.2

SEQUEN 21 PAGE 9
GCID EP 59 STATION J G3.2




Figure 8B. Reconstructed Gas Chromatogram for Sample STATION J G3.2, Con't.

FLUORENE

REFERENCE SPECTRUM MASS	R.A.	ION COUNT FROM SAMPLE M.S.# 216	P RATIO
168	1	48.	48.0
167	16	265.	16.6
166	100	1246.	12.5
165	88	1154.	13.1
164	10	113.	11.3
163	13	191.	14.7
162	1	23.	23.0
139	6	102.	17.0
115	3	130.	43.3
113	1	40.	40.0
89	1	29.	29.0
87	1	21.	21.0
86	1	18.	18.0
83	8	582.	72.7
82	15	272.	18.1
81	4	304.	76.0
75	1	28.	28.0
74	1	25.	25.0
69	4	428.	107.0
63	3	63.	21.0
62	1	30.	30.0
51	1	49.	49.0

SUM OF IONS IN SAMPLE SPECTRUM ATTRIBUTABLE - TO REFERENCE COMPOUND 3164.	SUM OF R.A.s x MINIMUM P
	280.
	11.3

Figure 9. Search for Fluorene in MS #216 of Sample STATION J G3.2

FLUORENE

REFERENCE SPECTRUM MASS	R.A.	ION COUNT FROM SAMPLE M.S.# 218	P RATIO
168	1	12.	12.0
167	16	85.	5.3
166	100	664.	6.6
165	88	547.	6.2
164	10	72.	7.2
163	13	113.	8.7
162	1	30.	30.0
139	6	63.	10.5
115	3	93.	31.0
113	1	20.	20.0
89	1	38.	38.0
87	1	9.	N.A., IONCOUNT < MINCD
86	1	10.	10.0
83	8	248.	31.0
82	15	104.	6.9
81	4	291.	72.7
75	1	4.	N.A., IONCOUNT < MINCD
74	1	4.	N.A., IONCOUNT < MINCD
69	4	338.	84.5
63	3	29.	9.7
62	1	7.	N.A., IONCOUNT < MINCD
51	1	30.	30.0

SUM OF IONS IN SAMPLE SPECTRUM ATTRIBUTABLE TO REFERENCE COMPOUND 1487. SUM OF R.A.s x MINIMUM P 280. 5.3

Figure 10. Search for Fluorene in MS #218 of Sample STATION J G3.2

boundaries, peak area (in scan time *ion count), and peak with closest ARI to the reference compound. The intermediate step between analysis and user printout is presented in Figure 11.

The intensity attributable to the spectrum of fluorene within its search window is summarized in Figure 11. A graphical output of this data is presented on the reconstructed gas chromatogram in Figure 12. Small peaks were observed at either end of the major fluorene peak. The fluctuation in intensity has three possible causes. First, mass spectral distortions over the elution interval of the compound may have led or contributed to apparent fluctuations in the ideal mass spectrum. Second, the presence of a complex mixture of unresolved compounds (or 'noise') may have led or contributed to the observed results. Third, a compound exhibiting an essentially similar mass spectrum may have eluted in the vicinity of fluorene. This possibility often arises, as illustrated in the next sample for the analysis of 1-methyl phenanthrene.

The reverse search for 1-methyl phenanthrene occurred in a very complex chromatographic region of the GC-MS analysis of sample STATION J. This region is characterized by severe peak overlap, and a large background to peak ratio (approximately 5:1). Manual GC-MS interpretation of regions similar to this require a great deal of skill. The intermediate reverse search results (from Program HOW) are summarized in Figure 13. This data is presented in graphical form in Figure 12.

FLUORENE

SUM OF IONS IN SAMPLE
SPECTRUM ATTRIBUTABLE
TO REFERENCE COMPOUND

SPECTRUM
211
212
213
214
215
216
217
218
219
220
221
222
223

31.
59.
56.
165.
176.
3164.
6459.
1487.
283.
140.
118.
182.
182.

PEAK ANALYSIS
MAXIMUM AREA
212.4 113.
> 216.9 11961.
222.5 423.

NOTE: > ASSIGNED TO PEAK NEAREST REFERENCE R.I.

30) FLUORENE	152.1	149.6	212.4	212	213	221	224	113
		> 151.9	216.9	213	221	224		11961
		154.8	222.5	221	224			423

Figure 11. Attributable Intensity within the Search Window for
Fluorene in Sample STATION J G3.2

SEQUEN 21 PAGE 5
GCID EP 59 STATION J G3.2




Figure 12. Reconstructed Gas Chromatogram of STATION J G3.2 with Illustrations of Intensity Attributable to Fluorene and 1-Methyl Phenanthrene. Full Scale Represents an Intensity of 62,000 Ions. The Area found Attributable to Fluorene is Represented as a Black Peak Centered around Scan 217. The Area found Attributable to 1-Methyl Phenanthrene is Represented as a Series of Black Peaks Centered around Scan 369.

1-METHYL PHENANTHRENE

SUM OF IONS IN SAMPLE
SPECTRUM ATTRIBUTABLE
TO REFERENCE COMPOUND

SPECTRUM
363
364
365
366
367
368
369
370
371
372
373
374
375

2172.
1227.
1825.
2095.
989.
931.
3705.
4602.
3988.
3653.
1541.
843.
778.

PEAK ANALYSIS

MAXIMUM	AREA
365.7	5988.
> 370.1	19576.

NOTE: > ASSIGNED TO PEAK NEAREST REFERENCE R.I.


Figure 13. Attributable Intensity within the Search Window for 1-Methyl Phenanthrene in Sample STATION J G3.2

Two major peaks were observed in the search window for the reference spectrum of 1-methyl phenanthrene. This is due to the occurrence of structural isomers of methyl phenanthrene in the analysis window. The structural isomers of a compound often exhibit similar mass spectra and retention. The reference mass spectrum of 1-methyl phenanthrene, and one other compound referred to as me-178 is presented in Figure 14. The exact identity of this me-178 (an abbreviation for a methyl derivative of a compound of M.W. 178) was not established, although it is known to be a methyl phenanthrene or methyl anthracene from a previous GC-MS analysis of a mixture of these compounds. The mass spectra of the two compounds, in addition to the retention index, is very similar. Five isomers of me-178 compounds are presently in the reverse search library. Three more may potentially be defined for methyl phenanthracenes/anthracenes.

Analysis Output

The final step for the reverse search analysis for any compound is to provide a concise output of the computer interpretation. An example of the data output for analysis of aromatic compounds is presented in Figures 15A-C. Following the report on each compound, the computer proceeds to the next reference compound until an end of file mark is reached to either the GC-MS data or the reference library. A second example is presented in Figure 16 A-B for the analysis of pesticides and polychlorinated biphenyls in the same GC-MS data set.

>>>>>> LIBRARY NAME:ARILI FILE# 69
 RRI 240.8 <X>ME-178



>>>>>> LIBRARY NAME:ARILI FILE# 71
 RRI 244.1 <X>1-METHYL PHENANTHRENE




Figure 14. Reference Data for Two Isomers of Structure CH₃-X, where X Represents an Aromatic Compound of Molecular Weight 178.

PROGRAM ONE - REVERSE GC-MS SEARCH INCORPORATING RRI. REVISION 10/81/
81 (x)

----- ANALYSIS * PARAMETERS -----
SEARCH WINDOW +- 6 SCANS FROM PREDICTED R.I.
REJECT SCAN WHEN 3% REFERENCE SPECTRUM MISSING
DETN CONTRIB VIA 10 OR MORE INTENSE FRAGMENTS

----- GC-MS FILE PARAMETERS -----
GCID:EP059 STATION J G3.2
<X> RRI'S> 0 100 200 300 400 500
<X> 0SCAN> 27.0 117.7 308.8 447.5 573.0 705.7 805.8

REC#	COMPOUND NAME	ARI	ARI	SCAN	START#	END#	AREA
1)	NAPHTHALENE (R.S.)	.0	-1.0	27.0	25	29	17522
			3.3	30.0	29	31	199
2)	BENZOTHIOPHENE	4.0	3.3	30.0	29	31	265
3)	UNKN (SIM 3MECYHEXENON	15.2	8.8	35.0	35	36	81
			18.7	44.0	43	45	21
4)	BENZOTHIAZOLE	19.8	-	-	-	-	-
5)	4-PHENYL-3-BUTEN-2-ONE	38.9	-	-	-	-	-
6)	2-METHYL NAPHTHALENE	53.8	53.8	75.8	74	77	3937
			58.4	80.0	79	81	79
7)	1-METHYL NAPHTHALENE	62.6	58.4	80.0	79	81	84
			62.5	83.7	82	86	4698
8)	4-PHENYL-3-BUTEN-2-ONE	89.8	89.3	108.0	107	109	62

Figure 15A. Analysis Printout of Program ONE, for Evaluation of STATION J G3.2 with a Reference Library of Aromatic Compounds.

26) 4-PHENYLDECANE	147.5	145.2	204.2	203	207	33413
	>	147.3	208.1	207	211	2016
		150.2	213.7	211	215	3212
27) C3-NAPHTHALENE	147.9	147.8	209.0	205	216	5702
28) ME-ACENAPHTHYLENE	150.8	149.6	212.4	211	213	149
	>	151.8	216.8	213	221	9515
29) 3-PHENYLDECANE	151.9	149.9	213.1	211	217	17603
	>	152.5	218.0	217	220	1770
		154.4	221.7	220	223	2281
30) FLUORENE	152.1	149.6	212.4	212	213	113
	>	151.9	216.9	213	221	11961
		154.8	222.5	221	224	423
31) C3-NAPHTHALENE	153.9	152.0	217.0	215	218	623
	>	153.7	220.3	218	223	2158
		155.6	224.0	223	225	138
		156.6	225.8	225	227	97
32) C2-BIPHENYL	158.5	155.6	224.0	224	225	559
		156.7	226.0	225	227	547
	>	158.2	228.9	227	231	3233
		160.0	232.3	231	236	1697
33) 2-PHENYLDECANE	161.3	159.3	231.0	229	233	5508
	>	162.2	236.6	233	239	3677
		163.9	239.8	239	241	873
34) BIPHENYL-AMINE	163.6	162.0	236.1	235	237	1125
	>	163.2	238.4	237	240	2007
		164.7	241.3	240	244	2599

Figure 15B. Analysis Printout of Program ONE, for Evaluation of STATION J G3.2 with a Reference Library of Aromatic Compounds, Con't.

63)	ME-DIBENZOTHIOPHENE	222.9	> 220.6	337.4	335	343	10795
					335	338	794
					338	344	5752
					344	347	1022
64)	1-PHENYL NAPHTHALENE	227.5	> 226.7	345.8	343	351	3480
					350	353	361
65)	ME-DIBENZOTHIOPHENE	229.7	> 229.6	349.8	344	352	5892
					352	356	1957
66)	ME-DIBENZOTHIOPHENE	232.0	> 232.0	349.4	348	352	5606
					352	357	3984
					357	360	2239
67)	ME-178	236.8	> 234.3	356.3	355	357	370
					357	364	19759
					364	366	2883
68)	ME-178	238.3	> 234.3	356.3	356	357	430
					357	360	10130
					360	364	19228
					364	368	7543
69)	ME-178	240.8	> 238.1	361.6	360	364	18801
					364	368	7159
					368	372	16012
70)	4-H-CYCLOPENTA(d,e,f)P	242.4	> 238.1	361.7	362	364	1469
					364	371	19483
					371	374	3493
71)	1-METHYL PHENANTHRENE	244.1	> 241.0	365.7	364	368	5988
					368	376	19576

Figure 15C. Analysis Printout of Program ONE, for Evaluation of STATION J G3.2 with a Reference Library of Aromatic Compounds, Con't.

PROGRAM ONE - REVERSE GC-MS SEARCH INCORPORATING RRI. REVISION 10/21/81

ANALYSIS * PARAMETERS
 SEARCH WINDOW +- 6 SCANS FROM PREDICTED R.I.
 REJECT SCAN WHEN 3% REFERENCE SPECTRUM MISSING
 METH CONTRIB VIA 10 OR MORE INTENSE FRAGMENTS

GC-MS FILE PARAMETERS
 GOLD:EP059 STATION J G3.2
 <X> RRI'S> 200 300 -1 0 0
 <X> ESCAN> 308.8 447.5 -1.0 .0 .0 .0

LIBRARY DATA		ARI		OBSERVED PEAK DATA		AREA
REC#	COMPOUND NAME	ARI	ARI	SCAN	START#-END#	
24)	CL3-BIPHENYL	203.9	201.2	310.5	308 312	338
			> 203.7	313.9	312 316	1023
			206.2	317.4	316 319	397
			207.9	319.8	319 321	235
25)	DELTA BHC	213.3	210.0	322.7	322 324	567
			211.8	325.2	324 326	550
			> 213.2	327.2	326 329	853
			215.3	330.0	329 331	552
			216.9	332.2	331 334	615
26)	CL3-BIPHENYL	215.3	211.9	325.3	325 328	466
			214.6	329.1	328 330	491
			> 216.0	331.0	330 335	841
			219.4	335.8	335 337	241
27)	ALPHA ENDOSULFAN	219.5	-	-	-	-
28)	TRICHLOROBIPHENYL	227.0	227.4	346.7	343 351	3741

Figure 16A. Analysis Printout of Program ONE, for Evaluation of STATION J G3.2 with a Reference Library of Pesticides and Polychlorinated Biphenyls.

29) CL3-BIPHENYL	227.1 >	226.8	346.0	343	347	807
		228.2	347.9	347	350	676
		231.0	351.8	350	353	863
30) CL3-BIPHENYL	232.3 >	228.0	347.7	348	350	386
		231.8	352.9	350	355	1827
		234.7	356.9	355	358	615
		235.9	358.6	358	360	265
31) CL3-BIPHENYL	238.5 >	234.7	356.9	357	358	353
		237.7	361.1	358	364	1381
		240.4	364.8	364	366	357
		242.4	367.7	366	369	426
32) HEPTACHLOR	240.9	-	-			
33) CL3-BIPHENYL	243.0 >	240.4	364.8	364	366	363
		242.7	368.0	366	369	499
		244.5	370.5	369	371	362
		245.9	372.5	369	373	743
		246.7	373.5	369	375	1038
34) CL4-BIPHENYL	245.9	-	-			
35) TETRACHLOROBIPHENYL	252.7 >	253.6	383.1	380	385	1801
		256.2	386.8	385	388	738
36) CL4-BIPHENYL	254.4 >	253.7	383.2	382	385	521
		256.4	387.1	385	389	381
37) CL4-BIPHENYL	256.7 >	253.0	383.2	382	385	564
		256.4	387.0	385	389	404
38) CL4-BIPHENYL	258.4 >	256.4	387.1	385	388	328
		257.7	388.8	388	390	226

Figure 16B. Analysis Printout of Program ONE, for Evaluation of STATION J G3.2 with a Reference Library of Pesticides and Polychlorinated Biphenyls, Con't.

DISCUSSION

The presence of an unresolved complex mixture (UCM) in some circumstances defeats the ability of Program ONE to establish the absence of a given reference compound. This is usually due to intensity occurring at all masses during the window of potential occurrence for a given compound. It presently represents a severe limitation with the interpretation of the data. An analogy is the ability to filter a known signal from a noise signal. It is possible to define the maximum potential contribution of the signal to the noise. However, it is impossible to definitively say if the given known signal, or some combination of other signals, contributed to the noise. An unresolved complex mixture is composed of the mass spectra of a finite (but unknown) number of components at various but relatively small concentrations. It is presently impossible to make definitive statements about the composition of the UCM from an unknown sample. Two approaches may be taken in this case. First, chemically separate the components into less complex mixtures, and perform analyses on the fractions. A second approach is to determine the potential intensity attributable to a reference spectrum. If alarming results are produced (i.e. for the concentration of a highly toxic compound), only then would further efforts be expended on a more definitive chemical separation.

A second problem encountered in this reverse search program is the acquisition of high quality reference spectra. Reference spectra quality is a factor that controls the quality of the reverse search

analysis. There are three factors in the quality of a reference spectrum: 1) the inclusion of irrelevant fragments, 2) the exclusion of relevant fragments, and 3) improper assignment of relative abundance relationships for ions of different m/e . Factors 1 and 2 affect the qualitative analysis, and all three factors affect the quantitative analysis performed by Program ONE.

The inclusion of superfluous fragments in a reference mass spectrum will cause the consistent rejection of this compound in the search, except in cases where the same impurity fragments occur along with the reference compound in the unknown data. Such impurity fragments may be due to a compound that co-eluted with the reference compound used for a reference spectrum (in the cases where a pure reference compound was unavailable). Such an impurity may occur intermittently with the reference compound, hence complicating recognition and removal of excess fragments. The analysis variable N2REJ was designed to allow analysis to continue in the event of the above problem. Analysts are free to set the value as low as the quality of the reference mass spectra will permit. If an excessive amount of impurities is present in the reference mass spectrum, the reverse search program will erroneously report the corresponding reference name as not present. A second possible result from a mixed reference spectrum is the reporting of the area where the mixture occurs, with the corresponding (but incorrect) compound name.

Exclusion of relevant fragments from a reference mass spectrum causes a variable degradation in the qualitative, and constant

degradation in the quantitative analysis of the reference compound in the unknown data. Qualitative analysis is affected by such missing fragments when the fragment is a definitive factor in determining the occurrence or absence of the reference compound. This is most likely to be the case if the fragment is in the high mass range, as such fragments tend to be more definitive (22). The evaluation of reference against unknown mass spectra is thus carried out from high to low mass in order to increase analysis efficiency. Quantitative analysis is affected in that the missing fragment abundances will not be included in the area found attributable to the reference compound. This quantitative problem will manifest itself in a consistent percent error in the area attributed to the reference compound.

The improper assignment of relative ion abundance relationships in a reference mass spectrum will also affect the quantification of a reference compound. This factor represents a significant difficulty in producing good quantitative results. In a particular instrument and at a constant partial pressure, the relative ion abundances in a mass spectrum remain constant. In the GC-MS separation, the concentration of a sample is changing as a peak is eluted and this change can distort the observed mass spectrum (23). This problem is illustrated in Figure 17. A reference spectrum may account for only a fraction of the intensity exhibited by the mass spectrometry of the same compound during a GC-MS analysis. Assignment of the true average relationships in a reference mass spectrum will provide the most complete utilization of spectra that are distorted in favor of the




Figure 17. Illustration of the Relationship Between Intensity, Mass and Time in a Gas Chromatographic-Mass Spectrometric Series of Measurements Over the Elution of a Gas Chromatographic Peak. Mass Range on the System used in this Project was 517 to 51 amu. Scan Cycle Time was 2.3 Seconds. Data was Gathered in Each Scan Cycle During a 1.0 Second Interval. The Actual Measurement Period is Indicated by Thick Black Lines.

lower weight fragments during initial buildup of a peak's concentration, and distorted by magnification of the intensities of the higher weight fragments during subsidence of the peak's concentration (assuming higher weight fragments are sampled a finite amount of time before lower fragment masses during collection of the mass spectrum). Improperly assigned relative abundance relationships in a reference mass spectrum may produce unpredictable results, including odd peak shapes of intensity attributable to observed GC-MS data.

It may be possible mathematically to correct the distortions of sampling on an observed unknown mass spectrum. If such a solution is possible, it will be a function of GC peak width, sampling frequency, mass-time relation within a scan, fragment intensity change, and potentially other factors including knowledge of components in the GC peak. Solution of this problem was not attempted, hence the results of Program ONE are described as semiquantitative. Optimization of these semiquantitative results may be possible by creating a response factor for each reference spectrum, plus a certainty factor for the number of observed mass spectra used in quantification and hence express how closely the observed is likely to conform to the averaged reference mass spectrum.

It is possible to optimize the relative abundance relationships in a reference mass spectrum. This was observed by collection of several reference mass spectra from several GC-MS analyses of the same compound. Different quantitative results were provided by these

reference spectra when used to analyze test GC-MS data. The reference mass spectrum that consistently produced the largest quantitative value of the test data was selected as the best representation of the reference mass spectrum. In many cases, two reference spectra produced similar (differing by only a few percent) quantitative results and often neither was consistently more efficient in use of test data. Efforts have been made to develop a routine for optimization of reference mass spectra (24). Such efforts reduce the amount of work required in selecting optimum reference data.

CONCLUSION

Program ONE has demonstrated itself as a useful tool to quickly evaluate the potential occurrence of a series of reference spectra in unknown GC-MS data. The qualitative data thus provided are limited by the quality of the reference data, and in some cases, the complexity of the GC-MS data. The results of Program ONE are frequently too complex to be arrived at manually without extended tedious calculations. The quantitative results are expressed as scan time *ion count and are admittedly only semi-quantitative due to the approximate nature of reference mass spectra.

REFERENCES

1. Crawford, L. R., Morrison, J. D. *Anal. Chem.*, 140, 1464 (1968).
2. Crawford, L. R., and Morrison, J.D. *Anal. Chem.*, 40, 1469 (1968).
3. Knock, B., Smith, I., Wright, D. and Ridley, R. *Anal. Chem.*, 42,1516. (1970).
4. Hertz, H. S., Hites, R. A. and Biemann, K. *Anal. Chem.*, 43,681 (1971).
5. Smith, D. H., Gray N. A. B., Pillinger, C. T., Kimble, B. J. and Eglinton, G. *Adv. Org. Geochem.*, 249 (1971).
6. Heller, S. R. *Anal. Chem.*, 44, 1951 (1972).
7. Smith, D. H. *Anal. Chem.*, 44, 536 (1972).
8. Dromey, R. G. Stefik, M. J., Rindfleisch, T. C. and Duffield, A. M. *Anal. Chem.*, 48,1368 (1976).
9. Blaisdell, B. E. and Sweeley, C. C. *Anal. Chim. Acta*, 117, 1 (1980).
10. Abramson, F. P. *Anal. Chem.*, 47, 45 (1975).
11. Pereira, W. E., Hoyano, Y., Reynolds, W. E., Summons, R. E. and Duffield, A. M. *Anal. Biochem.*, 55, 236 (1973).
12. Sweeley, C. C., Young, N. C., Holland, J. F. and Gates, S. C. *J. Chromatography*, 99, 507 (1974).
13. Ritter, G. L., Lowry, S. R., Isenhour, T. L. and Wilkins, C. L. *Anal. Chem.*, 48 591 (1976).
14. Pesyna, G. M., Venkataraghavan, R., Dayringer, H. E. and McLafferty, F. W. *Anal. Chem.*, 48, 1362 (1976).
15. Kovats, *Helv. Chim. Acta*, 41, 1915 (1958).
16. Bieri, R. Abstracts of 30th Southeastern Regional American Chemical Society Meeting, November 8-10 (1978).
17. Lee, M. L., Vassilaros, D. L., White, C. M., and Novotny, M. *Anal. Chem.*, 51, 768 (1979).
18. Blaisdell, B. E. and Sweeley, C. C. *Anal. Chim. Acta*, 117, 17 (1980).
19. Tunnicliff, D. D. and Wadsworth, P. A. *Anal Chem.*, 37, 1082 (1965).

20. Personal Communication with C. L. Smith.
21. Faul, H., Ed. 1954. Nuclear Geology, a Symposium on Nuclear Phenomena in the Earth Sciences, John Wiley & Sons, Inc., New York. "Fundamental Considerations, Instruments and Techniques of Detection and Measurement" (Ch. 1).
22. Pesyna, G. M., McLafferty, F. W. Venkataraghavan, R. and Dayringer, H. E. Anal. Chem., 47, 1161 (1975).
23. Junk, G. A. Int. J. Mass Spectrom. Ion. Phys., 8 1 (1972).
24. Blaisdell, B. E. Anal Chem., 49, 180 (1977).

VITA

Christopher Steven Hein

Born in New Haven, Connecticut, June 4, 1957. Graduated from Amity High School in Woodbridge, Connecticut, 1975, B.A. in Chemistry College of Wooster, 1979. M.A. Candidate, College of William and Mary, School of Marine Science, 1979-1981, with a concentration in Chemical Oceanography.

In August 1979, the author entered the Virginia Institute of Marine Science as a research assistant in the Department of Ecology and Pollution.

```

0001 FTM4
0002
0003 PROGRAM ONE
0004 CCCCC PROGRAM ONE REVERSE SEARCHES GCMS DATA FOR COMPOUNDS
0005 CCCCC STORED IN LIBRARY CALLED ARLI, IN EFMP FILE PMS01.
0006 CCCCC PMS01 READ IN 128 INTEGER WD. RECORDS, WHERE WD 1-100
0007 CCCCC ARE COMPACTED MASS-REL.ABUND. PAIRS (M-1,650)(R.A.-1,100)
0008 CCCCC WD 101-NUMBER OF MASS-R.A. PAIRS USED,
0009 CCCCC WD 102: RRI*10, WDS 104-123: COMPOUND NAME(ALPHA),
0010 CCCCC WD 103, WDS 124-128 ARE OPEN FOR FUTURE USES.
0011 COMMON CAT(10),MAT(10),IBUF(128),MAIN(100),CAIN(100)
0012 DIMENSION LIRA(128),TALY(25),IGCID(3),LIRRI(3),NOTE(3)
0013 EXTERNAL EXEC,OPEN,ERASE,CHNGE,GET
0014 CCCCC RELOCATABLE SUBROUTINES FOR THIS VERSION SAVED AS,
0015 CCCCC OPEN1,PLOT,CHNG5,GET5
0016 DATA IOUT/1/,IN/6/,KSD/5/,KSCAN/1/,LUN/2/,ITRK/190/
0017 &,ISECT/0/,LIRRI(1)/2HAR/,LIRRI(2)/2HIL/,LIRRI(3)/2HI /
0018 &,NOTE(1)/2H >/,NOTE(2)/2H /
0019 CALL ERASE(1)
0020 1 WRITE(IOUT,111)
0021 111 FORMAT('***** PROGRAM ONE: REVERSE GC-MS DATA SEARCH',
0022 &' USING RETENTION INDICES. *****',/,10X,'DEFINE:',/,
0023 &' IUNDO - # SCANS TO SEARCH FROM PREDICTED OCCURRENCE',/,
0024 &' N2REJ - # REF M.S. ABSENCE FOR COMPOUND'S ABSENCE',/,
0025 &' MINCD - MIN. IONCOUNT TO BE USED IN CONTRIBUTION DETN.',/,
0026 &,' ** LIMITATIONS **',/, ' IUNDO MIN=1 MAX=12',/,
0027 &' N2REJ MIN=1 MAX=100',/, ' MINCD MIN=1 NO MAX',/,
0028 &' NOW INPUT IUNDO,N2REJ,MINCD')
0029 READ(IN,X)IUNDO,N2REJ,MINCD
0030 IF((IUNDO.GT.12).OR.(IUNDO.LT.1).OR.(N2REJ.LT.1).OR.
0031 (N2REJ.GT.100).OR.(MINCD.LT.1)) GO TO 1
0032 2 WRITE(IOUT,7)
0033 7 FORMAT(' ENTER EXISTING GCID AS FIVE CHARACTERS ')
0034 CCCCC MAKE SURE DESIRED GCMS ANALYSIS IS ON DISK.
0035 READ(IN,8) IGCID
0036 8 FORMAT(3A2)

```

```

0036 CALL OPEN(111,IGCID,IERROR)
0037 IF(IERROR)2,10,4
0038
0039 4 PAUSE 1
0040 10 WRITE(IOUT,5)
0041 CCCCC INPUT RRI DATA FOR LOCATING LIBRARY COMPOUNDS.
0042 5 FORMAT('INPUT 1ST STANDARD'S RRI, IE 0, OR 100 ...')
0043 READ(IN,X)MAT(1)
0044 WRITE(IOUT,9)
0045 9 FORMAT(' INPUT SCAN TIMES FOR THE RETENTION STANDARDS,
0046 &,/, " SMALLEST TO LARGEST. INPUT -1 AFTER LAST SCAN TIME.')
0047 READ(IN,X) CAT(1)
0048 DO 6 I=2,10
0049 READ(IN,X) CAT(I)
0050 IF(CAT(I).EQ.-1) GO TO 20
0051 IF(CAT(I).LE.CAT(I-1)) GO TO 10
0052 6 MAT(I)=MAT(I-1)+100
0053 20 MAT(I)=-1
0054 DO 21 L=(I+1),10
0055 MAT(L)=0
0056 21 CAT(L)=0.0
0057 CALL EXEC(24,6,IGCID,1,IBUF,IERROR)
0058 CALL EXEC(24,5,IGCID,0,IERROR)
0059 IF(IERROR.NE.0)STOP 2
0060 CALL ERASE(1)
0061 WRITE(IOUT,100)IUNDO,NREJ,MINCD,IGCID,(IBUF(I),I=3,29)
0062 &,(MAT(I),I=1,7),(CAT(I),I=1,7)
0063 100 FORMAT('<X> PROGRAM ONE - REVERSE GC-MS SEARCH',
0064 &' INCORPORATING RRI. REVISION 10/21/81 <X>',/,
0065 &,'----- ANALYSIS X PARAMETERS -----',/,
0066 &'SEARCH WINDOW +--,13," SCANS FROM PREDICTED R.I.',/,
0067 &'REJECT SCAN WHEN",13,"X REFERENCE SPECTRUM MISSING',/,
0068 &'DETN CONTRIB VIA",15," OR MORE INTENSE FRAGMENTS',/,
0069 &,'----- GC-MS FILE PARAMETERS -----',/,
0070 &,'GCID:",3A2.8X,27A2,
&'<X> RRI'S",7I8,/, <X> #SCAN",7F8.1,2/,

```

```

0071      &* /-----LIBRARY DATA-----\ /-----OBSERVED*
0072      &* PEAK DATA-----\ /, 'REC# COMPOUND NAME', 13X, 'ARI',
0073      &* ARI SCAN START#-END# AREA', /, )
0074      CCC >>>>> 999 IS MAIN LOOP, DO ANALYSIS FOR LIBR-TH COMP. <<<<<<
0075      DO 999 LIBR=1,1200
0076      CALL OPEN(501,LIRRI,IERROR)
0077      CCCCCC STOP 10 MEANS LIBRARY NOT ON DISK, OR OPEN DAMAGED.
0078      IF(IERROR.NE.0) STOP 10
0079      CCCCCC COPY LIBR-TH LIBRARY RECORD INTO LIRA(.).
0080      CALL EXEC(24,6,LIRRI,LIBR,LIRA,IERROR)
0081      CCCCCC STOP 77 INDICATES END OF LIBRARY REACHED.
0082      IF(LIRA(102).EQ.-10) STOP 77
0083      IF(IERROR.EQ.0) GO TO 24
0084      WRITE(IOUT,1000)IERROR
0085      PAUSE 3
0086      24 CALL EXEC(24,5,LIRRI,0,IERROR)
0087      ARI=LIRA(102)/10.
0088      CCCCCC CHANGE LIBRARY ARI TO RELEVANT SCAN#.
0089      ICENT=CHNGE(ARI,1)
0090      IF(ICENT+1)30,999,32
0091      CCCCCC ICENT RETURNS -1, BEFORE 1ST RRI, GO TO NEXT COMPOUND.
0092      CCCCCC ICENT RETURNS -2, AFTER LAST RRI, END ANALYSIS (STOP 7)
0093      30 STOP 7
0094      CCCCCC SETUP SCAN RANGE TO SEEK COMPOUND IN.
0095      32 IBSCN=ICENT-IUNDO
0096      IF(IBSCN.LE.0)IBSCN=1
0097      IESCN=ICENT+IUNDO
0098      LASTM=LIRA(101)
0099      CCCCCC STOP 11 MEANS # M/E-R.A. POINTS IN LIBRARY IS>100(BAD).
0100      IF(LASTM.GT.100) STOP 11
0101      NRM=LIRA(101)+1
0102      CCCCCC NRM-NUMBER REF. MASSES +1 FOR BACKWARDS SEARCH LOOP (42)
0103      SUMRA=0.
0104      CCCCCC NOW DECODE LIBRARY MASS-REL.AB PAIRS.
0105      DO 26 I=1,LASTM

```

```

0106 MAIN(I)=320-(LIRA(I)/100)
0107 CCCCCC PAUSE 4 MEANS LIBRARY MASSES NOT IN SEQUENTIAL ORDER.
0108 IF((I.NE.1).AND.(MAIN(I).LE.MAIN(I-1))) PAUSE 4
0109 LIRA(I)=1+IABS(LIRA(I))-IABS((LIRA(I)/100)*100)
0110 26 SUMRA=SUMRA+LIRA(I)
0111 IREJ=(N2REJ*SUMRA)/100
0112 CCCCCC OPEN UP THE GCMS FILE TO READ DATA.
0113 CALL OPEN(11,IGCID,IEROR)
0114 IF(IEROR.NE.0)PAUSE 5
0115 CCCCCC CLEAN OUT COMPOUND CONCENTRATION SCOREBOARD
0116 DO 33 I=1,25
0117 33 TALLY(I)=0.
0118 CCCCCC NOW LOOK AT EACH SCAN YOU WANT TO.
0119 DO 36 IWANT=IBSCN,IESCN
0120 CCCCCC ESTABLISH RELATIVE SCAN FOR 'SCOREBOARDS'.
0121 IRSCN=1+IWANT-IBSCN
0122 CCCCCC NZIPS IS $M/E MISSING, IF NZIPS>N2REJ: FORGET SCAN.
0123 NZIPS=0
0124 CALL GET(IGCID,KSD,KSCAN,KFND,IWANT,LASTM,IEROR)
0125 IF(IEROR.EQ.0) GO TO 38
0126 IF(IEROR.EQ.21) GO TO 45
0127 WRITE(IOUT,1000)IEROR
0128 PAUSE 6
0129 38 DIGMIN=10000.
0130 DO 42 J=1,LASTM
0131 I=NRM-J
0132 CCCCCC LOOK AT SPECTRUM STARTING WITH HIGH MASS END, MORE EFFIC?
0133 IF(CAIN(I).EQ.0.) GO TO 34
0134 IF(CAIN(I).LT.MINCD) GO TO 42
0135 DSIG=CAIN(I)/LIRA(I)
0136 CCCCCC DSIG=POTENTIAL QUANTIFYING PARAMETER FOR MASSES>MINCD.
0137 IF(DSIG.LT.DIGMIN)DIGMIN=DSIG
0138 CCCCCC DIGMIN=LIMITING DSIG, USED TO QUANTIFY KNOWN IN SPECTRUM.
0139 GO TO 42
0140 34 NZIPS=NZIPS+LIRA(I)

```

```

0141 CCCCC IF TOLERANCE OF MISSING IONS IS EXCEEDED, FORGET SCAN.
0142 IF(NZIPS.GT.IREJ) GO TO 41
0143 42 CONTINUE
0144 CCCCC 5T CONC OF LIBR COMPOUND IN TALLY().
0145 TALLY(IRSCL)=SUMRAXDIGMIN
0146 CCCCC IF FIRST SCAN OF AN ANALYSIS, SAVE ITS ADDRESS.
0147 41 IF(IWANT.NE.IBSCN) GO TO 36
0148 KSD=KFND
0149 KSCAN=IWANT
0150 36 CONTINUE
0151 CALL EXEC(24,5,IGCID,0,IEROR)
0152 PART TWO, INTERPRET THE DATA
0153 CCCCC LOOK FOR VALLIES AND PEAKS IN CONCENTRATION (TALY)
0154 1000 FORMAT(" EFMP ERR#",I4)
0155 45 LFENET=1+IWANT-IBSCN
0156 NPKS=0
0157 CCCCC 1ST DET'N WHERE PKS AREA LOCATED. STORE IN CAIN(NPK).
0158 DO 44 J=2,(LFENET-1)
0159 IF((TALY(J).LE.TALY(J-1)).OR.(TALY(J).LT.TALY(J+1)))GOTO44
0160 NPKS=NPKS+1
0161 CAIN(NPKS)=J
0162 44 CONTINUE
0163 IF(NPKS.EQ.0) GO TO 998
0164 SARI=1000.
0165 CCCCC >>>> WORK ON EACH PEAK INDIVIDUALLY <<<<< CCCCC
0166 DO 51 I=1,NPKS
0167 CCCCC FIND BEGINNING OF PK. DEFAULT TO BEG. OF WINDOW.
0168 DO 46 J=1,(CAIN(I)-1)
0169 K=1+CAIN(I)-J
0170 IF((TALY(K-1).GT.TALY(K)).OR.(TALY(K).EQ.0.)) GO TO 47
0171 46 CONTINUE
0172 CCCCC FIND END OF PK. DEFAULT TO END OF WINDOW.
0173 47 DO 48 L=CAIN(I),LFENET
0174 J=L
0175 IF((TALY(J+1).GT.TALY(J)).OR.(TALY(J).EQ.0.)) GO TO 49

```



```

0176 48 CONTINUE
0177 CCCCC STORE ABSOLUTE LOC OF PK ST IN CAIN(15+): END IN (30+
0178 49 CAIN(I+15)=K+IBSCN-1
0179 CAIN(I+30)=J+IBSCN-1
0180 CCCCC INOW INTEGRATE FROM VALLEY TO VALLEY VIA TRAPEZOIDAL METH.
0181 CAIN(I+45)=(TALY(K)+TALY(J))/2.
0182 NN=J-K-1
0183 DO 50 M=1,NN
0184 50 CAIN(I+45)=CAIN(I+45)+TALY(K+M)
0185 CCCCC DET'N EXACT PEAK LOCATION, STORE IN CAIN(80+NPK)
0186 X=CAIN(I)
0187 X0=(X-1.)*X(X-1.)
0188 X1=X*X
0189 X2=(X+1.)*X(X+1.)
0190 B=((X1-X2)*TALY(X-1.))+((X2-X0)*TALY(X))+((X0-X1)*TALY(X+1)
0191 &))/2.
0192 A=(TALY(X+1.)+TALY(X-1.)-(2.*TALY(X)))/2.
0193 SCAN=-B/(2.*A)+IBSCN-1
0194 CAIN(60+I)=SCAN
0195 CCCCC DET'N OBSERVED ARI. ALSO FIND MIN DELTA ARI.
0196 CAIN(75+I)=CHNGE(SCAN,-1)
0197 IBUF(I)=2
0198 DELRI=ABS(CAIN(75+I)-ARI)
0199 IF(DELRI.GE.SARI) GO TO 51
0200 CCCCC CLOSEST PK TO ARI BASE REL LOCS IN IBPK & IEPK.
0201 IBPK=K
0202 IEPK=J
0203 SARI=DELRI
0204 CCCCC MAKE PREPARATIONS FOR PUTTING " >" @ CLOSEST PEAK'S ARI.
0205 IBUF(MOLD)=2
0206 MOLD=I
0207 IBUF(I)=1
0208 51 CONTINUE
0209 CCCCC REPORT WHAT YOU KNOW.
0210 WRITE(IOUT,101)LIBR,(LIRA(I),I=104,114),ARI,(NOTE(IBUF(I)),

```

```

0211      &CAIN(75+I),CAIN(60+I),CAIN(15+I),CAIN(30+I),CAIN(45+I),
0212      &I-1,NPKS)
0213      101 FORMAT(14," " ,11A2,F6.1,20(1A2,F6.1,F7.1,I6,I6,I10,/,34X))
0214      GO TO 999
0215      998 WRITE(IOUT,102)LIBR,(LIRA(I),I=104,114),ARI
0216      102 FORMAT(14," " ,11A2,F6.1,7X,"-",6X,"-",/,)
0217      999 CONTINUE
0218      END
0219      ENDS
XXXX LIST END ***
6
:LI,S,1,SOPN1

```

```

0001 FTN4,L SUBROUTINE OPEN(IPKN,IFNA,IERROR)
0002 C SUBROUTINE OPEN DEFINES AND OPENS EXTENDED FILE MANAGEMENT
0003 C PACKAGE (EFMP) PACK NUMBER AND FILE NAME. IF SERIOUS
0004 C ERROR OCCURS, OPEN ABORTS AND REPORTS ERROR IN IERROR.
0005 C IF FILE NOT FOUND, THIS IS REPORTED AND IERROR=-3 RETND
0006 C FOR GCMS, IPKN=111,IFNA=IGCID FOR LIBR, IPKN=501,IFNA=ARILI
0007 DIMENSION ITRBUF(256),IOPNTB(128),NOTRB(2),IFNA(3)
0008 NOTRB(1)=1
0009 NOTRB(2)=1
0010 CALL EXEC(24,1,IOPNTB,128,ITRBUF,NOTRB,2,IERROR)
0011 IF((IERROR.NE.0).AND.(IERROR.NE.27)) GO TO 2
0012 CALL EXEC(24,4,IFNA,IPKN,1,0,1,IERROR)
0013 2 IF(IERROR.EQ.3) IERROR=-3
0014 IF(IERROR.EQ.25) IERROR=0
0015 IF(IERROR.GT.0)WRITE(1,1)IERROR
0016 1 FORMAT(' OPEN EFMP ERROR #',I3)
0017 RETURN
0018 END
0019 ENDS
0020 XXXX LIST END ***

```

:LI,S,1, SCHN5

```
0001 FTN4,L
0002 FUNCTION CHNGE(TI,IR2S)
0003 CHNGE CONVERTS SCAN TO RRI, OR RRI TO SCAN, DEPENDING ON R2S
0004 C IF R2S NEG.,SCN>>RRI IF R2S POS.,RRI>>SCN
0005 C IF TI>LAST STD (OR R2S OMITTED) -2 RETND.
0006 C IF TI<1ST STD: -1 RETND. CHNGE REQUIRES COMMON CAT(10) &
0007 C MAT(10). SCN TIMES IN CAT, RRI IN MAT. -1 TERMINATES CAT LIST.
0008 C CAT(N)'S SCN TIME MUST CORRESPOND TO MAT(N)'S RRI.
0009 COMMON CAT(10),MAT(10),II(128),MM(100),CCC(100)
0010 CHNGE=-1
0011 IF(IR2S)5,8,6
0012 5 IF(TI.LT.CAT(1)) RETURN
0013 GO TO 3
0014 6 IF(TI.LT.MAT(1)) RETURN
0015 3 DO 1 I=2,8
0016 IF(CAT(I).EQ.-1) GO TO 8
0017 IF(IR2S)7,9
0018 7 IF(TI.LE.CAT(I)) GO TO 2
0019 GO TO 1
0020 9 IF(TI.LE.MAT(I)) GO TO 4
0021 1 CONTINUE
0022 8 CHNGE=-2
0023 RETURN
0024 2 CHNGE=((TI-CAT(I-1))/(CAT(I)-CAT(I-1)))*100.)+MAT(I-1)
0025 RETURN
0026 4 CHNGE=((TI-MAT(I-1))/100.)*(CAT(I)-CAT(I-1))+CAT(I-1)
0027 RETURN
0028 END
0029 ENDS
0030 XXXX LIST END XXXX
●
```

```
:LI,S,1,SGET5
```

```

0001 FTN4,L
0002 SUBROUTINE GET(IFNA,KSD,KSCAN,KFND,IWANT,LASTM,IEERROR)
0003 COMMON XXX(10),III(10),IBUF(128),MAIN(100),CAIN(100)
0004 DIMENSION IFNA(3)
0005 ITMA=1
0006 IEERROR=-1
0007 IF(KSCAN.GT.IWANT) RETURN
0008 ISCAN=KSCAN-1
0009 K=KSD-1
0010 DO 10 I=1, LASTM
0011 10 CAIN(I)=0.
0012 CCCCCCCCCC GET GCMS SUBROUTINE CCCCCCCCCCCCCCCCCCCCCC
0013 C IONCOUNTS FOR MASSES IN MAIN(I) ARE STORED IN CAIN(I).
0014 C IFNA THE FILE NAME OF GCMS RUN
0015 C KSD A KNOWN FILE SECTOR, WITH...
0016 C ISCAN A KNOWN SCAN NUMBER BEGINNING WITH SECTOR KSD
0017 C IWANT THE SCAN NUMBER DESIRED FOR IONCOUNT LOADING
0018 C IWANT MUST BE GREATER OR EQUAL TO ISCAN
0019 C LASTM NUMBER OF MASSES TO FIND IONCOUNTS FOR.
0020 CCCCCCCCCCCCCC THESE VARS ARE RETURNED W/ANSWERS CCCCCCCCCCCCCC
0021 C KFND THE ITH RECORD THAT SCAN=IWANT BEGINS ON.
0022 C IEERROR IF ZERO, EVERYTHINGS O.K.
0023 C IF NEGATIVE, PROBLEM WITH GET SUBPROGRAM
0024 C -1) IWANT<ISCAN
0025 C -2) KSD IS NOT THE FIRST SECTOR OF A SCAN,
0026 C AND THEREFOR MAY NOT BE A VALID ADDRESS
0027 C IF POSITIVE, EFMP ERROR (IF 21, END OF GCMS DATA)
0028 C MASSES IN MAIN(I) MUST BE ORDERED F/LOW TO HIGHEST.
0029 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

```

```

0030 17 K=K+1
0031 KA=-1
0032 CALL EXEC(24,6,IFNA,K,IBUF,IERROR)
0033 IF(IERROR.NE.0)RETURN
0034 IF((IBUF(4).NE.0).AND.(ISCAN.LT.IWANT)) GO TO 17
0035 IF(((KSCAN-1).EQ.ISCAN).AND.(K.NE.KSD)) GO TO 13
0036 C THIS CHECKS THAT KSD IS FIRST SECTOR OF KSCAN.
0037 IF(IBUF(4).NE.0) GO TO 19
0038 ISCAN=ISCAN+1
0039 IF(ISCAN.LT.IWANT) GO TO 17
0040 IF(ISCAN.GT.IWANT) RETURN
0041 KFND=K
0042 AFCTR=2*IBUF(3)
0043 KA=5
0044 18 KA=KA-2
0045 19 KA=KA+2
0046 IF(KA.GT.128) GO TO 17
0047 IF(IBUF(KA).EQ.-1) RETURN
0048 IF(MAIN(ITMA)-IBUF(KA))14,11,19
0049 11 CAIN(ITMA)=IBUF(KA+1)
0050 IF(CAIN(ITMA).LT.0)CAIN(ITMA)=CAIN(ITMA)+65535.
0051 CAIN(ITMA)=CAIN(ITMA)*AFCTR
0052 14 ITMA=ITMA+1
0053 IF(ITMA.LE.LASTM) GO TO 18
0054 RETURN
0055 13 IERROR=-2
0056 RETURN
0057 C COMPILED VERSION SAVED AS GET5, APPX 650 WDS LONG.
0058 END
0059 ENDS
0060 **** LIST END ****

```

```

0001      FTN4
0002
0003      PROGRAM TWO
0004      DIMENSION IOPNTB(128),ITRBUF(128),RA(128),MON(128),LIRRI(3)
0005      EXTERNAL OPEN,ERASE,TPLLOT,GRAPH
0006      CC COMPILED VERS. OF SUBROUTINES CALLED OPEN1,PLOT,GRAPH.
0007      DATA LIRRI(1)/2HAR/,LIRRI(2)/2HIL/,LIRRI(3)/2HI /,
0008      &IOUT/1/,IN/6/
0009      GO TO 11
0010      2 CALL OPEN(501,LIRRI,IERROR)
0011      IF(IERROR.NE.0)STOP 13
0012      CALL ERASE(1)
0013      WRITE(IOUT,133)
0014      133 FORMAT(" WHAT RANGE OF FILE NUMBERS DO YOU WISH LISTED",
0015      &" IN THE DIRECTORY",/, "INPUT 1ST, 2ND. IF 1ST=-1, ALL.")
0016      READ(IN,X)L1,L2
0017      IF((L1.GT.0).AND.(L2.GT.L1)) GO TO 183
0018      L1=-1
0019      L2=1200
0020      183 CONTINUE
0021      WRITE(IOUT,14)
0022      14 FORMAT("FILE# ARI $M-RA COMPOUND NAME")
0023      DO 10 I=L1,L2
0024      CALL EXEC(24,6,LIRRI,I,MON,IERROR)
0025      IF(IERROR.EQ.21)GO TO 10
0026      IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0027      ARI=MON(102)/10.
0028      WRITE(IOUT,12)I,ARI,MON(101),(MON(L),L=104,123)
0029      12 FORMAT(I4," ",F6.1,2X,15,2X,20A2)
0030      IF(MON(102).EQ.-10) GO TO 11
0031      10 CONTINUE
0032      11 WRITE(IOUT,15)
0033      CALL OPEN(501,LIRRI,IERROR)
0034      IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0035      15 FORMAT(" INPUT FUNCTION DESIRED",/, "1) LOOK AT A SPECTRUM",
0036      &/, "2) DELETE A SPECTRUM",/, "3) MODIFY A LIBRARY SPECTRUM",

```

```

0036      &./,"4) MANUALLY INPUT A SPECTRUM",/,"5) LIST DIRECTORY",/
0037      &"6) STOP & QUIT?"
0038      READ(IN,X) ICHOI
0039      IF(ICHOI-2)16,17,303
0040      CCCCCCCCCCCC MANUAL INPUT SECTION CCCCCCCCCCCCCZZZZ
0041      303 IF(ICHOI-4)18,304,305
0042      305 IF(ICHOI-6)2,594
0043      594 STOP 67
0044      304 WRITE(IOUT,306)
0045      306 FORMAT("INPUT COMPOUND NAME")
0046      READ(IN,31)(MON(L),L=104,123)
0047      WRITE(IOUT,310)
0048      310 FORMAT("INPUT ARI")
0049      READ(IN,X)ARI
0050      MON(102)=ARI*10
0051      WRITE(IOUT,311)
0052      311 FORMAT("HOW MANY MASS-REL.ABUN POINTS DO YOU HAVE?")
0053      READ(IN,X)MON(101)
0054      DO 312 I=1,MON(101)
0055      GO TO 308
0056      307 WRITE(IOUT,309)
0057      309 FORMAT("ILLEGAL ENTRY, CHECK ORDER AND NUMBER ")
0058      308 WRITE(IOUT,313)I
0059      313 FORMAT("INPUT MASS,RA FOR",I3,"TH PAIR")
0060      READ(IN,X)ITRBUF(I),RA(I)
0061      IF((ITRBUF(I).LT.51).OR.(ITRBUF(I).GT.517).OR.(RA.LT.1).
0062      &OR.(RA.GT.100)) GO TO 307
0063      IF((I.NE.1).AND.(ITRBUF(I).LE.ITRBUF(I-1))) GO TO 307
0064      INTEN=RA(I)
0065      MON(I)=((320-ITRBUF(I))*100)+(INTEN-1)
0066      IF(ITRBUF(I).GT.320)MON(I)=MON(I)-(2*(INTEN-1))
0067      312 CONTINUE
0068      CALL OPEN(501,LIRRI,IERROR)
0069      IF(IERROR.NE.0)PAUSE
0070      CC LOOK FOR APPROPRIATE SPOT TO PUT THE NEW ENTRY.

```



```

0106
0107
0108
0109
0110
0111
0112
0113
0114
0115
0116
0117
0118
0119
0120
0121
0122
0123
0124
0125
0126
0127
0128
0129
0130
0131
0132
0133
0134
0135
0136
0137
0138
0139
0140

I=1
IF(ITRBUF(NUOB).GT.299)I=2
NUOB=NURO
IF(I.EQ.2)WRITE(IOUT,3207)
3207 FORMAT(6/,14X,"350",11X,"400",11X,"450",11X,"500")
WRITE(IOUT,2222)
WRITE(IOUT,2223)(ITRBUF(L),RA(L),L=1,NUOB)
2222 FORMAT(/,7("MASS RA "))
2223 FORMAT(7(I4,14,2X))
70 CALL GRAPH(I)
DO 71 I=1,NUOB
MF=50
IYLO=600
IF(ITRBUF(I).LT.299)GO TO 72
MF=300
IYLO=450
72 IXLO=((ITRBUF(I)-MF)*4)+12
IYHT=RA(I)+IYLO
CALL TPLOT(1,0,IXLO,IYLO)
71 CALL TPLOT(1,1,IXLO,IYHT)
READ(IN,X)DUMMY
GO TO 11
CCCCCCCCCCCC DELETE SECTION CCCCCCCCCCCCCC
17 WRITE(IOUT,21)
21 FORMAT(" WHAT FILES IS TO BE DELETED? INPUT -1 TO STOP.")
READ(IN,X) IDST
IF(IDST.EQ.-1) STOP 47
NDST=IDST+1
DO 22 INSTR=NDST,1200
INSTR=INSTR-1
C NOW PICK UP SECTOR INSTR AND PUT IT DOWN ON INSTR.
CALL EXEC(24,6,LIRI,INSTR,MON,IERROR)
IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
CALL EXEC(24,8,LIRI,INSTR,MON,IERROR)
IF((MON(102).EQ.-10).OR.(IERROR.EQ.21)) GO TO 23

```

```

0141 22 CONTINUE
0142 23 CALL EXEC(24,5,LIRRI,0,IERROR)
0143 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0144 GO TO 11
0145 666 STOP 45
0146 CCCCCCCCCC MODIFY OR STOP SECTION CCCCCCCCCCCCCC
0147 18 WRITE(IOUT,24)
0148 24 FORMAT(' WHAT FILE# DO YOU WISH TO MODIFY?, INPUT -1',
0149 &' TO STOP')
0150 READ(IN,X) IFILN
0151 IF(IFILN.LT.1) GO TO 2
0152 CALL OPEN(501,LIRRI,IERROR)
0153 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0154 CALL ERASE(1)
0155 CALL EXEC(24,6,LIRRI,IFILN,MON,IERROR)
0156 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0157 IF(IERROR.NE.0) GO TO 2
0158 WRITE(IOUT,25)
0159 25 FORMAT('DO YOU WISH TO,/,,'1) CHANGE ARI',/,
0160 &'2) CHANGE COMPOUND'S NAME',/,,'3) CHANGE M/R.A ?')
0161 READ(IN,X)ICHOI
0162 IF(ICHOI-2)26,27,28
0163 26 ARI=MON(102)/10.
0164 WRITE(IOUT,29)ARI
0165 29 FORMAT(' RRI=',F5.1,' INPUT NEW RRI')
0166 READ(IN,X)ARI
0167 MON(102)=ARI*10.
0168 GO TO 312
0169 27 WRITE(IOUT,30)(MON(L),L=104,123)
0170 30 FORMAT(' NAME=',20A2,' INPUT NEW NAME')
0171 READ(IN,31)(MON(L),L=104,123)
0172 31 FORMAT(20A2)
0173 CC MAKE SURE NEW ENTRY'S ARI IS PROPERLY PLACED.
0174 WRITE(IOUT,494)
0175 494 FORMAT(' DESTROY OLD NAME'S FILE? 1=YES, 2=NO')

```

```

0176 READ(IN,X)IDES
0177 IF(IDES.NE.1) GO TO 123
0178 CALL EXEC(24,8,LIRRI,IFILN,MON,IERROR)
0179 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0180 CALL EXEC(24,5,LIRRI,0,IERROR)
0181 GO TO 11
0182 CCCCCCCCCC CHAGE MASS-REL.ABUN PAIR SECTION CCCCCCCCCCCCCCCC
0183 28 CONTINUE
0184 CALL OPEN(501,LIRRI,IERROR)
0185 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0186 CALL ERASE(1)
0187 CALL EXEC(24,6,LIRRI,IFILN,MON,IERROR)
0188 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0189 DO 169 I=1,MON(101)
0190 ITRBUF(I)=320-(MON(I)/100)
0191 RA(I)=1+IABS(MON(I))-IABS((MON(I)/100)*100)
0192 WRITE(IOUT,222)
0193 WRITE(IOUT,223)(ITRBUF(L),RA(L),L=1,MON(101))
0194 WRITE(IOUT,34)
0195 34 FORMAT(' DO YOU WANT TO 1) ADD',/,16X,'2) DELETE A M-RA'
0196 ' ,. PAIR',/,16X,'3) STOP EDITING?')
0197 READ(IN,X)ICHOI
0198 IF(ICHOI-2) 35,36,37
0199 37 GO TO 11
0200 35 WRITE(IOUT,38)
0201 38 FORMAT('INPUT MASS, REL-ABUN TO ADD INTO LIBRARY')
0202 READ(IN,X)MAS,REL
0203 IF((MAS.LT.51).OR.(MAS.GT.517).OR.(REL.LT.1).OR.(REL.GT.100)
0204 &) GO TO 34
0205 DO 40 I=1,MON(101)
0206 IF(MAS.EQ.ITRBUF(I)) GO TO 28
0207 IF(ITRBUF(I).GT.MAS) GO TO 41
0208 40 CONTINUE
0209 41 CONTINUE
0210 MON(101)=MON(101)+1

```

```

0211 DO 42 J=I,MON(101)
0212 KMA=ITRBUF(J)
0213 KRA=RA(J)
0214 ITRBUF(J)=MAS
0215 RA(J)=REL
0216 MAS=KMA
0217
0218 42 REL=KRA
0219 DO 168 I=1,MON(101)
0220 MON(I)=$((320-ITRBUF(I))*100)+(RA(I)-1)
0221 IF(ITRBUF(I).GT.320) MON(I)=MON(I)-(2*(RA(I)-1))
0222 168 CONTINUE
0223 GO TO 676
0224 36 WRITE(IOUT,101)
0225 101 FORMAT('WHAT IS THE MASS OF THE PAIR YOU WISH TO DELETE?')
0226 READ(IN,X)MAS
0227 DO 102 I=1,MON(101)
0228 IF(ITRBUF(I).EQ.MAS) GO TO 103
0229 102 CONTINUE
0230 GO TO 28
0231 103 DO 104 J=I,MON(101)
0232 104 MON(J)=MON(J+1)
0233 676 CONTINUE
0234 CALL EXEC(24,8,LIRRI,IFILN,MON,IERROR)
0235 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0236 CALL EXEC(24,5,LIRRI,0,IERROR)
0237 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0238 GO TO 28
0239 33 CALL EXEC(24,8,LIRRI,IFILN,MON,IERROR)
0240 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0241 CALL EXEC(24,5,LIRRI,0,IERROR)
0242 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0243 GO TO 18
0244 1000 FORMAT(' EFMP ERROR #',I4)
0245 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC

```

```

0246 CCCCCC PROGRAM TWO REQUIRES OPEN1, PLOT AND GRAPH
0247 CCCCCC SUBROUTINES TO BE ATTACHED WHEN THIS
0248 CCCCCC PROGRAM IS LOADED.
0249 CCCCCC
0250 END
0251 ENDS
*** LIST END ***
@
:LI,S,1,SGRAPH

```

```

0001 FTN4,L SUBROUTINE GRAPH(NPLOTS)
0002 IFIF=4
0003 ITEN=1
0004 IFIU=4
0005 DO 400 MA=12,1004,4
0006 IFIU=IFIU+1
0007 MAN=0
0008 IF(IFIU.NE.5) GO TO 401
0009 MAN=MAN+5
0010 IFIU=0
0011 ITEN=ITEN+1
0012 IF(ITEN.NE.2) GO TO 401
0013 ITEN=0
0014 MAN=MAN+5
0015 IFIF=IFIF+1
0016 IF(IFIF.NE.5) GO TO 401
0017 IFIF=0
0018 MAN=MAN+5
0019
0020 401 CONTINUE
0021 CALL TPLOT(1,0,MA,601)
0022 CALL TPLOT(1,1,MA,(600+MAN))
0023 IF(NPLOTS.NE.2) GO TO 400
0024 CALL TPLOT(1,0,MA,449)

```

```

0025      CALL TPLOT(1,1,MA,(450+MAN))
0026      CONTINUE
0027      CALL TPLOT(1,0,7,700)
0028      CALL TPLOT(1,1,7,600)
0029      CALL TPLOT(1,1,1015,600)
0030      CALL TPLOT(1,1,1015,700)
0031      IF(NPLOTS.NE.2) RETURN
0032      CALL TPLOT(1,0,7,550)
0033      CALL TPLOT(1,1,7,450)
0034      CALL TPLOT(1,1,1015,450)
0035      CALL TPLOT(1,1,1015,550)
0036      RETURN
0037      END
0038      ENDS
0039      XXXX LIST END XXXX
0040      *
0041      :LI,S,1,STHR

0001      PROGRAM THREE
0002      COMMON ITRBUF(256),IOPNTB(128)
0003      DIMENSION RA(256),CAIN(50,5),MON(256),ABUN(85),IBUF(128)
0004      1, LAST(128),AMASS(40),MAIN(50),IGCID(3),LIRRI(3)
0005      EQUIVALENCE (LAST,AMASS)
0006      EQUIVALENCE (RA,CAIN)
0007      EQUIVALENCE (MON,ABUN), (MON(171),IBUF)
0008      DATA IOUT/1/,IN/6/,KSD/5/,KSCAN/1/,LUN/2/,ITRK/190/
0009      1, ISECT/0/,LIRRI(1)/2HAR/,LIRRI(2)/2HIL/,LIRRI(3)/2HI /
0010      CC LUN-2 MEANS LOWER DISK IS USED AS SCRATCH AREA & GCMS IS TOP.
0011      2 WRITE(IOUT,1001)
0012      READ(IN,1002)IGCID
0013      CALL OPEN(111,IGCID,IERROR)
0014      IF(IERROR)2,1,3
0015      3 PAUSE 1
0016
0001      FTN4

```

```

0017 1 WRITE(IOUT,700)
0018 700 FORMAT(" INPUT SCAN NUMBER CLOSEST TO COMPOUND MAX.")
0019 READ(IN,X)ICENT
0020 IF((ICENT.LT.2).OR.(ICENT.GT.1570)) GO TO 1
0021 IBSCN=ICENT-2
0022 IESCN=ICENT+2
0023 4 WRITE(IOUT,702)
0024 702 FORMAT("INPUT DESIRED MASS RANGE: LO,HI (LIMITS:51,517)")
0025 READ(IN,X)M1,M2
0026 IF((M2.LE.M1).OR.(M1.LT.51).OR.(M2.GT.517))GO TO 4
0027 CC DO 100 MASSES AT A TIME FOR EACH OF THE SCANS OBSERVED.
0028 DO 300 MA=M1,M2,51
0029 NUOB=0
0030 MB=MA+50
0031 IF(MB.GT.M2)MB=M2
0032 CC LOAD DESIRED MASSES INTO MAIN(), COUNT HOW MANY U/LASTM.
0033 LASTM=0
0034 DO 12 I=MA,MB
0035 LASTM=LASTM+1
0036 12 MAIN(LASTM)=I
0037 CC CLEAR OUT ION INTENSITY BUFFER CAIN VIA RA.
0038 DO 13 I=1,256
0039 13 RA(I)=0.
0040 IRS=0
0041 CC NOW GO LOOKING FOR INTENSITIES FOR SCANS NEAR CENTER.
0042 DO 10 IWANT=IBSCN,IESCN
0043 IRS=IRS+1
0044 CC THIS IS GET SUBROUTINE
0045 ITMA=1
0046 K=KSD-1
0047 ISCAN=KSCAN-1
0048 17 K=K+1
0049 KA=-1
0050 CALL EXEC(24,6,IGCID,K,IBUF,IERROR)
0051 IF(IERROR.EQ.0) GO TO 16

```

```

0052 WRITE(IOUT,1000)IERROR
0053 PAUSE 2
0054 16 IF((IBUF(4).NE.0).AND.(ISCAN.LT.IWANT)) GO TO 17
0055 IF((KSCAN-1).EQ.ISCAN).AND.(K.NE.KSD)) PAUSE 3
0056 CC PAUSE 3 INDICATES KSTOR IS NOT A BEGINNING SCAN ADDRESS.
0057 IF(IBUF(4).NE.0) GO TO 19
0058 ISCAN=ISCAN+1
0059 IF(ISCAN.LT.IWANT) GO TO 17
0060 IF(ISCAN.GT.IWANT) PAUSE 4
0061 CC PAUSE 4 INDICATES ERROR IN GET ROUTINE
0062 KFND=K
0063 AFCTR=2*XI BUF(3)
0064 KA=5
0065 18 KA=KA-2
0066 19 KA=KA+2
0067 IF(KA.GT.128) GO TO 17
0068 IF(IBUF(KA).EQ.-1) GO TO 5
0069 IF(MAIN(ITMA)-IBUF(KA))14,11,19
0070 11 CAIN(ITMA,IRS)=IBUF(KA+1)
0071 IF(IBUF(KA+1).LT.0)CAIN(ITMA,IRS)=CAIN(ITMA,IRS)+65535.
0072 CAIN(ITMA,IRS)=CAIN(ITMA,IRS)*AFCTR
0073 14 ITMA=ITMA+1
0074 IF(ITMA.LE.LASTM) GO TO 18
0075 CCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCCC
0076 5 IF(IRS.NE.1) GO TO 10
0077 KSD=KFND
0078 KSCAN=IWANT
0079 10 CONTINUE
0080 DO 15 KK=1,85
0081 15 IBUF(KK)=0
0082 DO 202 I=1, LASTM
0083 DO 100 J=2,4
0084 Y0=CAIN(I,(J-1))
0085 Y1=CAIN(I,J)
0086 Y2=CAIN(I,(J+1))

```



```

0037 IF((Y1.LE.Y0).OR.(Y1.LT.Y2)) GO TO 100
0038 IF((Y0.EQ.0).AND.(Y2.EQ.0)) GO TO 100
0039 CC IF DOWN HERE, MAX AT J.
0090 X0=J-1
0091 X1=J
0092 X2=J+1
0093 B=((X1XX1)-(X2XX2))XV0)+((X2XX2)-(X0XX0))XV1)
0094 B+((X0XX0)-(X1XX1))XV2))/2.
0095 A=(Y2+Y0-(2.XV1))/2.
0096 XMP=B/(2.XA)
0097 YMP=((XMP-X1)*(XMP-X2))/2.)XV0)+((XMP-X0)*(XMP-X2))
0098 B+((-Y1))+((XMP-X0)*(XMP-X1))/2.)XV2)
0099 IXMP=(XMPX10)-14
0100 IF((IXMP.GT.0).AND.(IXMP.LT.37)) GO TO 201
0101 WRITE(IOUT,666)(MAIN(N),(CAIN(LN,M),M=1,5),N=1,LASTN)
0102 666 FORMAT(' M/E=',I4,4X,'RSN(1-5):',5F9.0)
0103 WRITE(IOUT,666)MAIN(I)
0104 PAUSE 6
0105 CC PAUSE 6 MEANS PK. OBSERVED OUT OF RANGE.
0106 GO TO 100
0107 CC IXMP REPS. FRACTION BETWEEN IBSCN-ICENT-IESCN.
0108 201 NUOB=NUOB+1
0109 IF(NUOB.GT.85) PAUSE 7
0110 IBUF(NUOB)=MAIN(I)+(IXMPX1000)
0111 AMASS(IXMP)=AMASS(IXMP)+YMP
0112 ABUN(NUOB)=CAIN(I,(J-1))+CAIN(I,J)+CAIN(I,(J+1))
0113 100 CONTINUE
0114 202 CONTINUE
0115 ISECT=ISECT+2
0116 300 CALL BURIT(MON,256,LUN,ITRK,ISECT,0)
0117 CC NOW PRESENT AMASS DATA, ASK OPERATOR WHAT MC PK TIMES TO INCL.
0118 888 NUOB=0
0119 CALL ERASE(1)
0120 CC GENERATE REAL T.U.S
0121 DO 50 I=1,36

```

```

0122 50 ABUN(I)=(I/10.)+IBSCN
0123 WRITE(IOUT,52)ABUN(IBEG),ABUN(IEND)
0124 52 FORMAT("FORMAR MC PK RANGE",2F7.1,/,4("TU SCAN IONSUN X "))
0125 54 WRITE(IOUT,54)(L,ABUN(L),AMASS(L),L=1,36)
0126 54 FORMAT(4(I2,F6.1,I7," X "))
0127 WRITE(IOUT,56)
0128 56 FORMAT(" INPUT DESIRED INCLUSION RANGE FOR MC PKs",/,
0129 &" INPUT TU#1,TU#2. IF TU#1=-1, PROG HALTS.")
0130 READ(IN,X)IBEG,IEND
0131 IF( (IBEG.EQ.-1) ) STOP 7
0132 IF( (IEND.LT.IBEG).OR.(IEND.GT.36).OR.(IBEG.LT.1) ) GO TO 888
0133 WRITE(IOUT,58)
0134 58 FORMAT("WHEN SPECTRA IS PRESENTED",/, "RETN 1 FOR M/E",
0135 &,/, & R.A. TABLE NEXT TIME",/, " 0 FOR NO M/E TABLE",
0136 &,/, " 7 TRANSFERS SPECTRUM DATA TO DISK FOR LIBRARY USE"
0137 &,/, " NOW INPUT LARGEST DESIRED MASS FOR SPECTRUM")
0138 READ(IN,X)LAMA
0139 IF(LAMA.GT.M2)LAMA=M2
0140 WRITE(IOUT,59)MIN
0141 59 FORMAT(" MINIMUM R.A.",I4," INPUT NEW FILTER.")
0142 READ(IN,X)MIN
0143 IF( (MIN.LT.1).OR.(MIN.GT.100) ) MIN=1
0144 CALL ERASE(1)
0145 WRITE(IOUT,57) IGCID,ABUN(IBEG),ABUN(IEND),MIN,M1,LAMA
0146 57 FORMAT(5X,"<X> PROGRAM PARE: MS CLEANUP VIA MC PROFILE",
0147 &" ANALYSES <X>",/, "QC ID:",3A2,20X,"SCAN RANGE:",2F8.1,
0148 &, "MIN RA.",I3,"X",30X,"MASS RANGE:",2I8,7/,14X,"100",11X,
0149 &"150",11X,"200",11X,"250")
0150 DO 61 I=1,256
0151 ITRBUF(I)=0
0152 61 RA(I)=0.
0153 I100=IBEG*1000
0154 I200=(IEND*1000)+518
0155 BIG=0.
0156 DO 60 K=2,ISECT,2

```

```

0157 CALL BREAD(MON,256,LUN,ITRK,K,0)
0158 CC LOAD AAIN( )&MAIN( ) VIA TIM WITH TIMES,INTENS&MASSES.
0159 DO 64 I=1,85
0160 IF( (IBUF(I).LT.I100).OR.(IBUF(I).GT.I200) )GO TO 64
0161 CC IF IT COMES DOWN THIS WAY, MASS FOUND IN TIME RANGE.
0162 IF(NUOB.LT.257) GO TO 62
0163 PAUSE 11
0164 GO TO 888
0165 CC PAUSE 11 MEANS TOO MANY M/E PEAKS (>256).
0166 62 NUOB=NUOB+1
0167 ITRBUF(NUOB)=IBUF(I)-((IBUF(I)/1000)*1000)
0168 IF(ITRBUF(NUOB).EQ.ITRBUF(NUOB-1))PAUSE 10
0169 CC PAUSE 10 MEANS WIDE ENOUGH INCLUSION RANGE TO SEE MASS PX 2X.
0170 IF(ITRBUF(NUOB).GT.LAMA) GO TO 65
0171 RA(NUOB)=ABUN(I)
0172 IF(RA(NUOB).GT.BIG)BIG=RA(NUOB)
0173 64 CONTINUE
0174 60 CONTINUE
0175 65 IF(NUOB.EQ.0) GO TO 888
0176 NURO=0
0177 X1=100./BIG
0178 CC CODE MASS(ITRBUF( )) & RA(RA( )) INTO 1 INTERGER WORD (MON( ))
0179 DO 68 I=1,NUOB
0180 INTEN=RA(I)*X1
0181 IF(INTEN.LT.MIN) GO TO 68
0182 NURO=NURO+1
0183 MON(NURO)=((320-ITRBUF(I))*100)+(INTEN-1)
0184 IF(ITRBUF(I).GT.320)MON(NURO)=MON(NURO)-(2*(INTEN-1))
0185 68 CONTINUE
0186 CC DECODE MON( ) FOR PLOT, RESTORE SELECTED M-RA INTO ITR&RA.
0187 DO 69 I=1,NURO
0188 ITRBUF(I)=320-(MON(I)/100)
0189 69 RA(I)=1+IABS(MON(I))-IABS((MON(I)/100)*100)
0190 I=1
0191 IF(ITRBUF(NUOB).GT.299)I=2

```

```

0192 NUOB=MURO
0193 IF(I.EQ.2)WRITE(IOUT,3207)
0194 FORMAT(6/,14X,"350",11X,"400",11X,"450",11X,"500")
0195 3207 IF(NY.NE.1) GO TO 70
0196 WRITE(IOUT,2222)
0197 WRITE(IOUT,2223)(ITRBUF(L),RA(L),L=1,NUOB)
0198 2222 FORMAT(/,7("MASS RA "))
0199 2223 FORMAT(7(I4,I4,2X))
0200 70 CALL GRAPH(I)
0201 DO 71 I=1,NUOB
0202 MF=50
0203 IYLO=600
0204 IF(ITRBUF(I).LT.299)GO TO 72
0205 MF=300
0206 IYLO=450
0207 72 IXLO=((ITRBUF(I)-MF)*4)+12
0208 IYHT=RA(I)+IYLO
0209 CALL TPLOT(1.0,IXLO,IYLO)
0210 71 CALL TPLOT(1.1,IXLO,IYHT)
0211 READ(IN,X) NY
0212 IF((NY.NE.7).OR.(NUOB.GT.100))GO TO 888
0213 CALL ERASE(1)
0214 MON(101)=NUOB
0215 80 WRITE(IOUT,82)
0216 82 FORMAT("INPUT RRI, THEN COMPOUND NAME:")
0217 READ(IN,X)X1
0218 MON(102)=X1*10
0219 READ(IN,84)(MON(L),L=104,123)
0220 84 FORMAT(20A2)
0221 CC PREPARE TO ENTER NEW COMPOUND IN LIBRARY.
0222 DO 88 I=1,128
0223 88 LAST(I)=MON(I)
0224 CALL OPEN(501,LIRRI,IERROR)
0225 IF(IERROR.NE.0)PAUSE
0226 CC LOOK FOR APPROPRIATE SPOT TO PUT THE NEW ENTRY.

```

```

0227 DO 500 J=1,2000
0228 CALL EXEC(24,6,LIRRI,J,IBUF,IERROR)
0229 CC READ JTH, IF EOF OR RRR>NEW COMPOUNDS FALL THROUGH
0230 IF((IBUF(102).NE.-10).AND.(IBUF(102).LE.LAST(102)))GOTO500
0231 CALL EXEC(24,8,LIRRI,J,LAST,IERROR)
0232 IF(IERROR.NE.0)WRITE(IOUT,1000)IERROR
0233 CC PUT LAST WHERE IBUF WAS, THEN PUT IBUF IN LAST.
0234 IF(LAST(102).EQ.-10)GO TO 501
0235 DO 503 I=1,128
0236 503 LAST(I)=IBUF(I)
0237 500 CONTINUE
0238 501 CALL EXEC(24,5,LIRRI,0,IERROR)
0239 STOP 2
0240 CC PROGRAM THR REQUIRES OPEN3, PLOT & GRAPH IN LOADING.
0241 1000 FORMAT("EFMP ERROR NO.",I3)
0242 1001 FORMAT("ENTER AN EXISTING GCID AS FIVE CHARACTERS")
0243 1009 FORMAT(8F8.1)
0244 1002 FORMAT(3A2)
0245 END
0246 ENDS
0247 *** LIST END ***
●

```

!LI,5,1,SOPN3

```
0001 FTN4,L SUBROUTINE OPEN(IPKN,IFNA,IERROR)
0002 C SUBROUTINE OPEN DEFINES AND OPENS EXTENDED FILE MANAGEMENT
0003 C PACKAGE (EFMP) PACK NUMBER AND FILE NAME. IF SERIOUS
0004 C ERROR OCCURS, OPEN ABORTS AND REPORTS ERROR IN IERROR.
0005 C IF FILE NOT FOUND, THIS IS REPORTED AND IERROR=-3 RETND
0006 C FOR GCMS, IPKN=111,IFNA=IGCID FOR LIBR, IPKN=501,IFNA=ARILI
0007 COMMON ITRBUF(256),IOPNTB(128)
0008 DIMENSION IFNA(3),NOTRB(2)
0009 NOTRB(1)=1
0010 NOTRB(2)=1
0011 CALL EXEC(24,1,IOPNTB,128,ITRBUF,NOTRB,2,IERROR)
0012 IF((IERROR.NE.0).AND.(IERROR.NE.27)) GO TO 2
0013 CALL EXEC(24,4,IFNA,IPKN,1,0,1,IERROR)
0014 2 IF(IERROR.EQ.3) IERROR=-3
0015 IF(IERROR.EQ.25) IERROR=0
0016 IF(IERROR.GT.0)WRITE(1,1)IERROR
0017 1 FORMAT(' OPEN EFMP ERROR #',I3)
0018 RETURN
0019 END
0020 ENDS
0021 *** LIST END ***
```