Fungal Genetics Reports

Volume 25

Article 7

Storage of slime strains

C. P. Selitrennikoff *Merck and B., Ins.*

Follow this and additional works at: https://newprairiepress.org/fgr



This work is licensed under a Creative Commons Attribution-Share Alike 4.0 License.

Recommended Citation

Selitrennikoff, C. P. (1978) "Storage of slime strains," *Fungal Genetics Reports*: Vol. 25, Article 7. https://doi.org/10.4148/1941-4765.1727

This Technical Note is brought to you for free and open access by New Prairie Press. It has been accepted for inclusion in Fungal Genetics Reports by an authorized administrator of New Prairie Press. For more information, please contact cads@k-state.edu.

Storage of slime strains

Abstract

Storage of slime strain

Selitrennikoff, C. P.

Storage of slime strains.

The slime variant of N. crassa (FGSC #326; fz;sg;os-1, arg-1, cr-1, aur) can be maintained by repeated passage on-lititid or aggr-solidified medium and con be stored frozen in 10% dimethyl sulfoxide (Creighton and Trevithick (1973) Neurospora News], 20: 32) or as a component of a heterokaryon (Nelson et al. (1975) Neurospora News], 22: 15-16). However, I have found that petri dish and slant cultures of slime strains can be frozen in situ, stored at -70°C and subsequently thawed and revived. Simply, petri dishes and/or slants containing Nelson's medium B 17.5% Sorbitol (w/v), 1.5 Sucrose (w/v), 1X Vagel's Salts) solidified with 1.5% agar (and appropriate supplements) are inoculated and incubated for 5-10 days at 28°C. Petri dishes are wrapped en masse with aluminum foil (slants are sealed with parafilm) and placed in a -70° C freezer. To revive stored strains, dishes and slants ore allowed to thaw completely at room temperature and cell masses transferred to fresh dishes (or slants) with the gid of a rubber policeperson. Alternatively dishes or slants con be flooded with medium B and the liquid used as an inoculum for fresh agar-solidified medium. Thus for slime and two derivatives strains (slime strains containing cys-1) or inl) have been stored for four months and all cultures subsequently revived, Longer storage periods are currently being tested. - - Department of Boric Microbiology, Merck and b., Ins., Rahway, New Jersey 07065.