

The Utilization of High Frequency Percussive Ventilation to Reduce V-V Extracorporeal Oxygenation Membrane Support (Poster).

Kenneth Miller MEd, RRT-NPS

Lehigh Valley Health Network, Kenneth.Miller@lvhn.org

Kimberly D. Smith RRT-ACCS

Lehigh Valley Health Network, Kimberly_D.Smith@lvhn.org

James K. Wu MD

Lehigh Valley Health Network, james.wu@lvhn.org

Dave Marth RRT

Lehigh Valley Health Network, David.Marth@lvhn.org

Samantha K. Summers BSN

Lehigh Valley Health Network, samantha_k.summers@lvhn.org

See next page for additional authors

Follow this and additional works at: <https://scholarlyworks.lvhn.org/patient-care-services-nursing>



Part of the [Nursing Commons](#)

Published In/Presented At

Miller, K. Smith, K. Lindauer, L. Wu, J. Marth, D. Summers, S. Bennett, M. (2016, Oct). *The Utilization of High Frequency Percussive Ventilation to Reduce V-V Extracorporeal Oxygenation Membrane Support*. Poster Presented at: American Association of Respiratory care, San Antoni, Texas.

This Poster is brought to you for free and open access by LVHN Scholarly Works. It has been accepted for inclusion in LVHN Scholarly Works by an authorized administrator. For more information, please contact LibraryServices@lvhn.org.

Authors

Kenneth Miller MEd, RRT-NPS; Kimberly D. Smith RRT-ACCS; James K. Wu MD; Dave Marth RRT;
Samantha K. Summers BSN; Margaret Bennett BSN; and Lisa Lindauer MD

The Utilization of High Frequency Percussive Ventilation to Reduce V-V Extracorporeal Oxygenation Membrane Support

Kimberly Smith, BA, RRT-ACCS, Kenneth Miller MEd, MSRT, RRT-ACCS, Lisa Lindauer, MD, James K. Wu, MD, David Marth, RRT, Samantha Summers, BSN, Margret Bennett, BSN
Lehigh Valley Health Network, Allentown, Pennsylvania

Goals

Goal of Venous-Venous ECMO:


- Stabilization of gas exchange
- Minimize the risk of ventilator induced injury

Goal of Mechanical Ventilation During ECMO:


- Maintain lung recruitment
- Provide lung protection
- Augmentation of gas exchange if needed



VDR-4



VDR waveform



V-V ECMO



Maquet Rotoflow

Materials/Methods and Results

High Frequency Percussive Ventilation:

- Technologically classified as:
 - pneumatically driven
 - time cycled
 - pressure limited
 - bi-phasic percussive delivery
 - high frequency venturi flow interrupter
 - **exhalation is passive**


Role of HFPV (VDR) in ECMO Patients:

- Maximumized ECMO settings
- Retained secretions
- Inadequate lung inflation
- Constant poor pulmonary mechanics CLT<10cmH₂O
- Provides lung protection in unilateral lung disease process


Results						
Patient	Pre VDR ECMO FIO ₂ %	Pre VDR ECMO LPM	Post VDR ECMO FIO ₂ %	Post VDR ECMO LPM	VDR FIO ₂ %	VDR* PIP/PEEP cm/h20
1	100%	6lpm	70%	4lpm	40%	40/16
2	100%	7lpm	60%	5lpm	50%	38/20
3	100%	8lpm	80%	6lpm	50%	40/16
4	100%	7lpm	60%	5lpm	50%	38/18
5	100%	7lpm	100%	6lpm	90%	40/22
6	9lpm	9lpm	60%	6lpm	40%	40/18
7	100%	6lpm	70%	4lpm	50%	34/16
8	100%	6lpm	90%	4lpm	50%	38/18
9	100%	8lpm	60%	7lpm	40%	36/16
10	100%	7lpm	80%	5lpm	50%	40/18
11	100%	8lpm	90%	4lpm	50%	34/16
12	100%	7lpm	80%	5lpm	50%	40/16
13	100%	8lpm	60%	5lpm	40%	40/16
14	100%	8lpm	60%	6lpm	40%	38/14
15	100%	9lpm	100%	8lpm	90%	40/22

*Prior to placing on HFPV a pressure/ volume tool measurement was performed to determine starting airway pressure and PEEP parameters to set on the VDR.

Pre/POST HFPV



Pre VDR



Post VDR 24 hours

Discussion and Conclusions

DISCUSSION

- HFPV provides both an endobronchial wedge via the percussive rate and an oscillatory plateau via the convective rate.
- Provides an internal mucokinesis and maintains a patent airway
- With this ventilator strategy lower pressures and oxygen delivery can be employed and ECMO parameters can be often reduced.
- Benefit of lower ECMO FIO₂ and sweep increases the ability to match patient's hemodynamic demand by the delivery of a higher blood flow via the ECMO.
 - Less ECMO chatter
 - Less supplemental fluid replacement

CONCLUSIONS

- HFPV can help ECMO maintain gas exchange for patients at a lower FIO₂ and sweep settings.
- More research needs to be conducted to determine this ventilator strategy effect on morbidity and mortality.

© 2016 Lehigh Valley Health Network

610-402-CARE LVHN.org