

FITEXC – DIFFRACTION PROFILE FITTING PROGRAM RUN IN MS EXCEL

PETR VERTAT^{a,b,*}, JAN DRAHOKOUPIL^{a,c}

^a Department of Solid State Engineering, FNSPE CTU in Prague, Trojanova 13, 120 00, Prague 2, Czech Republic

^b Department of Magnetic Measurements and Materials, Institute of Physics AS CR, Na Slovance 2, 182 21, Prague 8, Czech Republic

^c Department of Material Analysis, Institute of Physics AS CR, Na Slovance 2, 182 21, Prague 8, Czech Republic

* corresponding author: vertat@fzu.cz

ABSTRACT. We present our custom made diffraction fitting program *FitExc* for the basic treatment of the 1D diffraction profiles, in particular, the *Intensity*(2θ) of X-ray, synchrotron or neutron diffraction measurements. In the current state, the program allows for simultaneous fitting of up to three profiles in one data set and is intended to be used mainly for precise lattice parameters determination or fitting of overlapped peaks. The entire program is built in the common *MS Excel* environment, it does not require any installation and is free to use.

KEYWORDS: X-ray diffraction, diffraction profile fitting, MS Excel.

1. INTRODUCTION

X-ray crystallography is a powerful tool with many applications: it is often used to determine the lattice parameters, the complete structure of solids, to examine the phase composition, residual stresses or texture. For multiple applications the precise examination of measured 1D *Intensity*(2θ) diffraction profiles represents a crucial part of data processing.

There are various programs for detailed diffraction profile fitting available. Many of them offer complex solutions implementing the Rietveld refinement process [1] (such as *X'Pert Highscore Plus* [2], *TOPAS* [3] or *FullProf* [4]). On the other hand, there are many applications, where quick and precise fitting of few diffraction peaks and their closest neighbourhood is sufficient. Even in this field, a few programs are used, such as *Fityk* [5]. Unfortunately, the majority of the programs for X-ray diffraction (XRD) data analysis does not, to our knowledge and experience, allow the user to fully control or customize all the aspects of the fitting in an easily accessible way, nor are they user-friendly or quick enough for easier applications and batch processing. Although the simple fitting of a few diffraction maxima might also be done manually (for example using the *Origin* software [6]), the approach is often clumsy and require wider knowledge of such environment.

We developed a simple user-friendly diffraction profile fitting program *FitExc* that does not suffer from mentioned disadvantages, yet provides precise results for basic applications. The program, first introduced in [7], is based on the *Microsoft Excel* environment, which is immediately available on the majority of computers today. Excel's interface together with *Visual Basic for Applications*-based macros create

a lightweight but powerful tool for the diffraction peak fitting. Another advantage is the immediate possibility to implement user's own scripts for various kinds of analyses such as macroscopic residual stress determination or the analysis of the crystallite size and strain using the Williamson-Hall plot.

2. PROGRAM INTERFACE AND FITTING FUNCTIONS

In the attempt to make *FitExc* as user-friendly as possible, the clickable interface with buttons, drop-down lists and check-boxes was created over the hidden worksheets used for the calculations. Actions are macro-powered and the overall approach is straightforward.

First, user is asked to select the file to be loaded. PANalytical *.xrdml* files are supported as well as standard *.xy*, *.xye* or *.txt* files (multiple columns in the ordinary textfile). Next, the fitting function and the radiation spectrum can be specified. Standard Cauchy, Gauss, Pearson VII and pseudo-Voigt functions are implemented, as well as their split variants for dealing with possible asymmetry. As for the radiation, user may choose from different X-ray tube presets or define his own spectrum containing up to 7 components with desired wavelength, relative intensity and width (Figure 1) [8]. Therefore, the program might be used not only for XRD, but also for neutron or synchrotron diffraction data fitting.

On the main fitting page, where the measured data and fitted function are plotted (Figure 2), the entire data processing is controlled. In the current state, up to three profiles might be fitted for one loaded dataset, in the future, the ability to use more than three profiles is to be added. Each of the profiles is defined by

Settings

loaded file: *H:\sample_400+040.xrxml* type: *xrxml*

file does contain multiple data sets
 number of data set: of 73 done !

fitting function: Cauchy pseudo-Voigt split pseudo-Voigt
 Gauss Pearson VII split Pearson VII

Radiation spectrum:

	include	rel. int.	λ (Å)	width
K α 1	<input checked="" type="checkbox"/>	1	1.789	1
K α 2	<input checked="" type="checkbox"/>	0.5	1.79289	1
component 3	<input type="checkbox"/>			
component 4	<input type="checkbox"/>			
component 5	<input type="checkbox"/>			
component 6	<input type="checkbox"/>			
component 7	<input type="checkbox"/>			

preset: ▼

initial values: try to determine initial values

use ";" instead of "," for solver commands
 fit positions first (more stable fit, sometimes)

FIGURE 1. Program interface: Initial settings worksheet shown after the data file is chosen. Here the basic settings can be specified, such as the spectrum of the radiation, function for fitting or initialisation of the new peak parameters.

function	formula ($f(2\theta)$)
Cauchy	$A_1/(1 + A_3(2\theta - A_2)^2)$
Gauss	$A_1 \exp(-A_3(2\theta - A_2)^2)$
Pearson VII	$A_1/(1 + A_3(2\theta - A_2)^2)^{A_4}$
pseudo-Voigt	$A_4 \text{ Gauss} + (1 - A_4) \text{ Cauchy}$

TABLE 1. Implemented fitting functions with following parameters: intensity at maximum A_1 , angular position of the maximum A_2 , the relative width A_3 and the shape parameter A_4 [9].

four parameters A_{1-4} , where A_1 is the intensity at maximum, A_2 the 2θ position of the maximum, A_3 is the relative width parameter and A_4 is the shape parameter. The simple formulas for such defined functions are denoted in Table 1 [9].

In the case of asymmetric profiles, instead of one A_3 and A_4 parameter, two sets of parameters (A_{31} , A_{41} and A_{32} , A_{42}) are used – one pair for the left and one pair for the right side of the peak. For background fitting, the constant, linear and quadratic forms were implemented so far. The program also calculates the integral intensity of the fitted reflections using analytical formulas where possible and numerical approach for the Pearson VII based functions. The loaded 2θ might be restricted if needed.

As one of the advantages, the user is immediately able to fix some parameters that will not be changed

during fitting. This makes the program very effective for fitting of complex patterns and multiple overlapped peaks.

Fitting algorithm is based on general reduced gradient non-linear least squares method. It is implemented using the *Solver (solver.xlam)* add-in and the standard worksheet environment of *MS Excel* run in the background.

To evaluate the goodness of the fit the standard discrepancy factors R_P , R_{WP} , R_{exp} and χ^2 are calculated and their values might be compared with those before fit.

3. ADDITIONAL FUNCTIONALITY

For processing of multiple data, our program allows loading of multiple files. These are then fitted automatically one after another. In case that *.xrxml* file containing multiple scans is loaded, the batch feature is enabled automatically providing the tool to fit the datasets sequentially. This feature has already been used for determination of the temperature development of selected lattice parameters of magnetic shape memory (MSM) alloys [10]. As our program enables determination of very small changes of the lattice parameters, their temperature evolution is clearly observable even in narrow temperature intervals (Figure 3).


Fitting algorithm is stable enough even for fitting of multiple peaks with great intensity differences

Fitting


profile parameters:

	profile 1	profile 2	profile 3	FIX
A1	6741.1363	7181.989	0	<input type="checkbox"/>
A2	73.593341	74.14202	81.07095	<input type="checkbox"/>
A3	15.831932	23	300	<input checked="" type="checkbox"/>
A4	0.4714673	0.726366	0.5	<input type="checkbox"/>

backgr.: type: linear
 parameters: $a + b(x-A2)$
 a: 104.678772
 b: 27.0485703



data loaded from 2θ: 72.00000 to: 76.00000
 min: 55.02636 max: 84.97295



factors included: Lorentz - Polarisation factor

file: H:\sample_400+040.xrdml
 scan: 20 of 73

fit 1 fit for the best χ^2

	R_p	R_{wp}	R_{exp}	χ^2
current:	1.170	1.646	0.978	2.830
previous:	1.371	1.906	0.982	3.762

copy results back to loading
 recover param. before fit

profile	int (num, rough)	%
profile 1	22500	56%
profile 2	17376	44%
profile 3	0	0%

profile	int (analytical)	%
profile 1	6343	57%
profile 2	4823	43%
profile 3	0	0%

FIGURE 2. Program interface: The main worksheet for data fitting. User is able to start the fitting, preset or check the fitted parameters, fix them if needed and to see the measured and fitted data including the difference curve.




FIGURE 3. Application example: precisely determined lattice parameters a and b of Ni-Mn-Ga MSM alloy and their evolution while heating. The sudden change denotes the structural transformation [10].

(so far even with two orders difference the results were reliable). To review a situation, the user may immediately choose between the linear and logarithmic scale of the main graph by clicking a single button.

Once the good fit is obtained, all the results including the basic settings might be stored in the easily accessible results worksheet. The procedure for easy recovery of the parameters of the fit to the fitting worksheet is also implemented.

4. CONCLUSIONS

We developed a new simple diffraction profile fitting program *FitExc*. This lightweight and user-friendly program can be run on any PC with *MS Excel* (2003 or later) installed. It is intended mainly for the precise fitting of the single reflections (especially their 2θ positions) based on selected or user-defined spectra of the radiation source and for fitting of overlapped peaks. This allows the precise determination of the lattice parameters, observation of their evolution or study of the residual stresses.

The program is free to use. To get your own copy, please visit <http://people.fjfi.cvut.cz/vertapet/>.

ACKNOWLEDGEMENTS

This work was supported by the Grant Agency of the Czech Technical University in Prague, grant No. SGS16/245/OHK4/3T/14, and Czech Science Foundation, grants No. 16-00043S and 17-00062S. Authors would like to thank O. Heczko for useful remarks and motivation - the tremendous amount of tasks requiring the development of our program :-).

REFERENCES

- [1] H.M. Rietveld: A profile refinement method for nuclear and magnetic structures. *J. Appl. Cryst.* **2** 65–71, 1969. doi: 10.1107/S0021889869006558
- [2] T. Degen, M. Sadki, E. Bron, U. Koenig, G. Nenert: The HighScore suite. *Powder Diffraction* **29** (Supplement S2):S13–S18, 2014. doi: 10.1017/S0885715614000840
- [3] DIFFRAC.SUITE TOPAS (Bruker, Karlsruhe, Germany), <https://www.bruker.com/products/x-ray-diffraction-and-elemental-analysis/x-ray-diffraction/xrd-software/topas.html>

- [4] J. Rodríguez-Carjaval: Recent advances in magnetic structure determination by neutron powder diffraction. *Physica B* **192** (1–2):55–69, 1993.
doi: 10.1016/0921-4526(93)90108-I
- [5] M. Wojdyr: Fityk: a general-purpose peak fitting program. *J. Appl. Cryst.* **43** 1126–1128, 2010.
doi: 10.1107/S0021889810030499
- [6] Origin (OriginLab, Northampton, MA).
- [7] P. Vertat: Diffraction profile fitting in MS Excel. *5. studentska vedecka konference fyziky pevných látek, Nove Hradky, Ceske vysoke uceni technicke v Praze*, 2015.
- [8] V. K. Pecharsky, P. Y. Zavalij: Fundamentals of Powder Diffraction and Structural Characterization of Materials, *Kluwer Academic Publishers: Boston/Dordrecht/London*, 2nd ed., p. 115, 2009.
- [9] R. Kuzel: Informace v praskovem difraktogramu a jeho zpracovani. *Materials Structure* **10** (1a) 18–19, 2003.
- [10] L. Straka, J. Drahokoupil, P. Vertat, J. Kopecek, M. Zeleny, H. Seiner, O. Heczko: Orthorhombic intermediate phase originating from {110} nanotwinning in Ni 50.0 Mn 28.7 Ga 21.3 modulated martensite. *Acta Mater.* **132** 335–344, 2017.
doi: 10.1016/j.actamat.2017.04.048