

# **Brain Computer Interfaces**




**Paul Gnanayutham**

# Rom Houben



# Structure of Brain




# Brain

1. Cerebrum
2. Cerebellum
3. Brainstem

# Cerebrum

- The cerebrum is the section where thoughts are created and memory is stored.
- Injury to the cerebrum can leave a person fully aware of their surroundings but unable to react to any events happening in the surroundings




# Lobes of Cerebrum

- The cerebrum also has five lobes:
  - **Frontal lobe:** motor cortex, which creates alpha brain waves
  - **Occipital lobe:** visual cortex that effects the visual perception
  - **Temporal lobe:** contains the cranial nerve and auditory cortex, damage can result in deafness
  - **Parietal lobe:** primary somatosensory cortex. Damage to this area of the brain affects the ability to use bio-potentials to manipulate a Brain-Body Interface
  - **Insular lobe:** affects emotion and damage to this region may affect a person's ability to relax

# Cerebellum


- Controls facial muscle co-ordination
- Damage to this area affects the ability to control eye movements and other facial muscle movements
- Will affect signals needed by Brain-Body Interfaces





# Brainstem

- Controls basic functions such as eating, respiration, heart rate
- The cranial nerves that carry the signals to control facial movements originate in the brainstem






# Bio-Potentials

- Electroencephalography (EEG) - amplitude can vary between 10 - 100  $\mu\text{V}$  when measured on the scalp or forehead, covers a frequency spectrum of 1 - 30 Hz and is divided into five classes.
- Electromyography (EMG) - 18 Hz, signals have an amplitude range of 0.2 - 2000  $\mu\text{V}$ .
- Electrooculargraphy (EOG) - 1.1 - 6.25 Hz, signals have an amplitude range of 1 - 4 mV.
- Slow Cortical Potentials (SCP) - frequency range 1 - 2 Hz, can be positive or negative, signals can be 5 - 8  $\mu\text{V}$
- Steady-State Visual Evoked Potential (SSVEP)/ Steady State Visual Evoked Responses (SSVER) - uses the 4 to 35 Hz frequency range, transfers data at high data transfer rates (68 bits/s) and occur at 100 - 1000 ms after the stimuli.


# Bio-Potentials

- P300 - displays a brain wave with positive amplitude, peaking at around 300 ms after task-relevant stimuli. This signal occurs in the delta (0.5 - 4 Hz) and theta (4 - 7 Hz) frequency range.
- N400 - displays a brain wave with negative amplitude, peaking at around 400 ms triggered by unexpected linguistic stimuli.
- Electrocochleography (ECoG) - - 300 - 1000  $\mu\text{V}$  amplitude and has a frequency of 40 Hz
- Low Frequency Asynchronous Switch Design (LF-ASD) - signals in the 1 - 4 Hz frequency range, with an amplitude of 10 - 100  $\mu\text{V}$ .
- Local Field Potential (LFP) - synchronous oscillations in the 15 - 30 Hz frequency range and have an amplitude of 6  $\mu\text{V}$ .
- Neuroprosthetic
- Motor Function


# Electroencephalography (EEG)




# Electromyography (EMG)



# Electrooculargraphy (EOG)



# P and N signals



# Brain Body Interfaces – Examples (these need some control movement above the neck)

- HeadMouse™
- Tonguepoint™
- Eye-tracking
- Thought Translation Device

Many brain-injured  
are so impaired that  
they cannot use any  
devices in this  
category



A	B	C	D	E	F
G	H	I	J	K	L
M	N	O	P	Q	R
S	T	U	V	W	X
Y	Z	1	2	3	4
5	6	7	8	9	space




# **Types of Brain - Body Interfaces for Traumatic Brain-Injured**


Two types:

- Invasive
- Non-invasive

# Invasive BBI



# Non-Invasive BBI




# Cyberlink™

Cyberlink™ used for our research, is a non invasive brain-body actuated control technology that combines eye-movement, facial muscle and brain wave bio-potentials detected at the user's forehead. This type of interface is used for non-verbal, quadriplegic, who have no way of communicating apart from using the brain waves.


# Cyberlink™ - Chosen as the best device for this research



# B1 thru B3 - Eye Movement




# B4 thru B6 - Alpha Brain Waves






# B7 thru B10 - Beta Brain Waves



# B11 - Muscle




# **Traumatic Brain Injury and Brain-body Interfaces (BBI)**

- Cyberlink can pick up various unwanted bio-potentials
- Such uncontrollable, erratic movements cause users frustration and fatigue
- Bringing the cursor back under control takes considerable effort, and may be impossible


# Research with BBI for the Traumatic Brain-Injured

- Issues with Brain-Injured Participants
  - Cognitive abilities are often not assessed
  - Unable to respond or really comatose
  - Difficulty of tracking participants over long periods
  - Etc.,


# Apparatus



# Personalised Tiling – Target Test



# Personalisation



A screenshot of a 'Change Profile' dialog box. The dialog has a title bar with the text 'Change Profile' and a close button. It contains two main sections: 'Settings' and 'Profile Appearance'. The 'Settings' section has six rows, each with a label, a text input field, and a unit description. The 'Profile Appearance' section has four rows, each with a label and a button with an ellipsis. At the bottom are 'OK' and 'Cancel' buttons.

Settings		
Wait Time on a Tile	<input type="text" value="250"/>	(In milli Seconds)
Blinking Speed	<input type="text" value="5"/>	(Times per Second)
Wait Time on a Target	<input type="text" value="1"/>	(In Seconds)
Wait Time at Start Area	<input type="text" value="15"/>	(In Seconds)
Application Launch Time	<input type="text" value="60"/>	(In Seconds)
Target Reach Time	<input type="text" value="15"/>	(In Seconds)

Profile Appearance	
Change Target	<input data-bbox="1062 878 1124 928" type="button" value="..."/>
Change Starting Area	<input data-bbox="1062 949 1124 999" type="button" value="..."/>
Change Background	<input data-bbox="1062 1021 1124 1071" type="button" value="..."/>
Change Target Positions	<input data-bbox="1062 1092 1124 1142" type="button" value="..."/>

OK Cancel



# Personalised Discrete Acceleration interface



- Personalised interface
- Discrete Acceleration
- Programmable targets for application launch or switching devices

# MIND OVER MATTER IS THE KEY TO POWERING THE LATEST MUST-HAVE GADGETS ABOUT TO HIT THE SHELVES

## Thought you were in control?

**A** WHOLE new generation of games is about to hit the shops, with some already claiming to be the "must have" toy of the year. So, what's so special about these new gadgets? Well, they can read your mind...

BY ROB PINNIGER


types of brain waves but the toys generally work by detecting two or three of them — alpha, beta or theta waves. The idea is to change your brain waves by either concentrating hard or attempting to relax and keep calm.

The toys can detect just how fiercely you are concentrating or how completely you are relaxing and use that information to drive the mechanism in the toy — a motor or fan for example.


As a result, higher concentration or deeper relaxation means the motor spins faster and your performance in the game is better. But be careful, as soon as you stop to think "yes, I'm winning", you lose your focus and your performance collapses.

Here we round up the pioneering games that could change the way we play...

### Unlocking the power of the mind



Can't get you out of my head: EEG gaming devices can read your brain waves



### TRICKS OF THE MIND

THE Mindflex requires you to concentrate to levitate the ball; the harder you concentrate the harder the motor which blows the ball in the air works and the higher the ball flies. The game also allows you to build up obstacle courses through which you have to manoeuvre your ball by a combination of levitating it and rotating the game board. This means you can build fiendish courses for others to attempt, try for quickest times around the course and so on, which should mean that the game keeps you amused. It's fun, easy to pick up and will hold your interest.

£49, amazon.com (pre-order for October, UK release unconfirmed)

9/10



### NEURAL PC

THE NIA, or Neural Impulse Actuator, is a headset designed to connect to your PC to allow you to control games without the use of a mouse or keyboard. That sounds pretty revolutionary but don't expect to sit at your desk thinking "move left", "swap weapon" and so on. You will need to assign certain muscular motions and mental states to specific game actions. For example, raising your eyebrows might be "fire" and concentrating might move your character forward. It requires extensive calibration and training but for serious gamers, the reward of fast reaction times can be worth it.

£79.99, amazon.co.uk

### REFLEX GAMING

SIMILAR to the NIA, the Emotiv EPOC is a control system for PC games. It uses 14 electrodes, compared with NIA's three. The extra hardware means that the EPOC is likely to offer finer levels of control but could also mean a longer period of training before you can use the device. Emotiv has decided not to sell the headsets direct to the public but instead to license them to game developers to be integrated into future releases.

£304 for a programmer's development kit, emotiv.com

8/10




### JEDI MIND CONTROL

THE idea with the Star Wars Force Trainer is to achieve higher and higher levels of concentration in order to move the Jedi Training Remote up and down the Training Tower. It's easy to get it going but if your concentration weakens, the ball drops back to the bottom. There are 15 levels of difficulty, with all interaction voiced by Yoda. It's fun, until the novelty has worn off.

£99.99, funfilledtoys.co.uk (Pre-order for delivery next month)

6/10



### BRAIN WAVE RACE

MINDBALL is a serious bit of kit and the price (£11,365) reflects this. It's a table, along which runs a track. You sit at the end of the table and attempt to move a ball along the track, forcing it closer to your opponent by relaxing your mind. The more you relax, the further it goes. The device comes with a monitor that shows measured brain waves and claims it has the health benefit of teaching you how to overcome anxiety.

Intended more for exhibitions. £11,365, or £595 to hire, vivifyeye.com/mindball/index.html

7/10




# Cyberlink

Control a computer with only your brain



<http://www.youtube.com/watch?v=ZKA2-zvinQc&feature=related>

# OCZ NIA Neural Impulse Actuator



<http://www.youtube.com/watch?v=P0NvIk-tTVk>

# Neurosky



[http://www.youtube.com/watch?v=wNr3yGcI\\_V8&feature=related](http://www.youtube.com/watch?v=wNr3yGcI_V8&feature=related)

# Other Brain Computer Games Consoles

Mind Balls

Mind Flex

**Jedi Mind Control**

**Any Ideas for  
collaborative Research?**



# Possible Future Research Opportunities

- Using Bio-Potentials as a diagnostic tool
- Rehabilitation games for brain injured
- Designing Interfaces to existing applications using BCIs
- Accessing web applications using bio-potentials
- Using a BCI to control a robot

# Possible Future Research Opportunities

- Designing an interface that can detect a range of bio-potentials, so that the appropriate/available one for that particular brain injury, can be chosen for the BCI
- A BCI that had probes that can fit in any part of a body to obtain bio-potentials not just from the head

**Any Questions?**